

GC Consumables Product Guide





GC Consumables Product Guide



Contents

	GC CONSUMABLES	P. 4	SAMPLE INTRODUCTION		SEPARATION	
	SAMPLE PREPARATION		Syringes	P. 8	Connections	P.14
	MEPSTM	P. 5	GC Inlet Liners	P.11		
			Septa	P.13	TROUBLESHOOTING	P.16
			Accessories	P.13		

GC CONSUMABLES

GC Consumables are designed to complement your Shimadzu GC-2010 and GC-2014 System. Shimadzu understands your chromatography analysis does not end with the selection of the GC column - the combination of components selected for your instrument make an important contribution to successful chromatography. Shimadzu GC Consumables are designed to provide you with the best separation possible.

SAMPLE INTRODUCTION – Syringes

Shimadzu's syringes for both manual and autosampler use incorporate Diamond Syringe Technology offering significantly improved levels of durability, clarity and accuracy. Features of the Diamond Syringe Technology include:

- Longest working life in the industry
- Improved solvent resistance
- Greater operational temperature range
- Reduced dead volume and carryover

SAMPLE INTRODUCTION – Septa

Many chromatography problems are caused as a result of the wrong septa material or inappropriate handling of the septa. Shimadzu septa are selected to contribute low bleed and optimum sealing for many injections.

SAMPLE INTRODUCTION – GC Inlet Liners

The GC inlet liner is where the sample is introduced and vaporized into the gaseous phase. The geometry of each of Shimadzu's inlet liners is important and the correct choice of inlet liner can significantly improve the performance of the chromatography. Inlet liner deactivation, quartz wool quantity and position are essential to ensure reproducible and accurate sample introduction for each sample type.

SEPARATION – GC Connections (Ferrules and Unions)

Shimadzu provides an array of connection types for use in the GC, each selected to ensure the best connection solution is achieved. The correct selection of the connection type will eliminate dead volumes, leaks during temperature cycling and problems with mismatched tubing sizes.




SEPARATION – GC Capillary Columns

Shimadzu offers a range of stationary phases including non-polar and polar functionalities.

SAMPLE PREPARATION – MEPS™

MEPS™ is Micro Extraction by Packed Sorbent and is a development for sample preparation and handling. MEPS is the miniaturization of conventional SPE packed bed devices from milliliter bed volumes to microliter volumes. The MEPS approach to sample preparation is suitable for reversed phases, normal phases, mixed mode or ion exchange chemistries. MEPS is available in a variety of common SPE phases.

MEPS BARREL INSERT AND NEEDLE



The MEPS Barrel Insert and Needle (BIN) assembly contains the stationary phase, and is built into the syringe needle.

SAMPLE SIZE AND SENSITIVITY

Sample volumes may be as little as 10 µL, or by taking multiple aliquots of 100 µL or 250 µL, samples of 1 mL or larger may be concentrated.


AUTOMATION

Extract samples and make injections on-line using a single device, reducing sample processing times and the need for operator intervention.

SORBENT LIFE

BIN life is dependent on the specific matrix being analyzed. For example, C18 analysis of whole plasma samples is conservatively 25-100 samples before the BIN needs to be changed. BIN life of cleaner samples is significantly longer.

Chromatogram of 22 pesticides standard solution (40 ng/mL)




CARRYOVER

The small quantity of phase in the MEPS BIN is easily and effectively washed between samples to reduce the possibility of carryover. This washing process is not practical with off-line SPE devices. With automation of MEPS, washing occurs while the previous sample is running.

FLEXIBLE AND EASY TO USE

The dimensions of the sorbent bed ensure performance remains identical to conventional SPE devices when used for extraction of similar samples. The AOC-MEPS system was developed to incorporate MEPS into a process automation workflow, combining sample preparation and analysis in a single platform

AOC-MEPS System



When automated by AOC-MEPS system, the injection volume is much larger compared to the commonly used GC injection volumes of 1 - 2 μL . With AOC-MEPS, the typical injections are 50 to 200 μL of elution solvent into the analytical instrument.

A large volume injection volume technique should be employed that removes the solvent volume from inside the injector unit while condensing the target compounds.

PTV Inlet Liner for AOC-MEPS System (AOC-20i + GC-2010 + PTV configuration)



OD (mm)	ID (mm)	Length (mm)	P/N
3.5	2.5	95	221-74830-09

AOC-20i

MEPS™ Syringe for AOC-20i

Syringe Volume	PTFE Tipped Plunger	Needle Length (mm)	Needle Gauge	Needle OD (mm)	Needle ID (mm)	Needle Tip	P/N	Spare Parts	P/N
100 µL	✓						221-74830-01	Plunger Pk 1	221-74830-10

Phase	Qty per pack	P/N
C18	5	221-74830-03
Silica	5	221-74830-04
C8 + SCX	5	221-74830-05
C2	5	221-74830-06
C8	5	221-74830-07
PDVB	5	221-74830-02
SDVB	5	221-74830-11
HDVB	5	221-74830-12
Development Kit - contains 1 each of C18, C8, Silica, C8+SCX, C2	1 of each phase	221-74830-08

AOC-5000


MEPS™ Syringes for AOC-5000

Syringe Volume	PTFE Tipped Plunger	Needle Length (mm)	Needle Gauge	Needle OD (mm)	Needle ID (mm)	Needle Tip	P/N	Spare Parts	P/N
100 µL	✓						221-75202-01	Plunger Pk 1	221-75202-03
250 µL	✓						221-75202-02	Plunger Pk 1	221-75202-04

Phase	Qty per pack	MEPS for 100 µL Syringe P/N	MEPS for 250 µL Syringe P/N
C18	5	221-75198-01	221-75199-01
Silica	5	221-75198-02	221-75199-02
C8 + SCX	5	221-75198-03	221-75199-03
C2	5	221-75198-04	221-75199-04
C8	5	221-75198-06	221-75199-06
Development Kit - contains 1 each of C18, C8, Silica, C8+SCX, C2	1 of each phase	221-75198-05	221-75199-05

SAMPLE INTRODUCTION – Syringes

Shimadzu Diamond syringes are the result of technological advancements in materials, design, and engineering. Designed to meet the ever increasing levels of sensitivity required by today's analyses, Shimadzu Diamond syringes give you a new level of accuracy and precision.



LONGER LIFE

Shimadzu Diamond syringes have a longer life. The improved solvent resistance and maximized operational temperature range along with the smoothest available internal glass surface ensure you receive the longest lifetime from your Shimadzu syringe.

SUPERIOR PERFORMANCE AND ROBUSTNESS

Shimadzu Diamond syringes have superior performance and robustness with unsurpassed levels of operational strength and durability. Potential for contamination is significantly reduced by the near-zero syringe dead volume and minimized adhesive in the flow path.

REDUCED CARRYOVER

Engineering enhancements have eliminated areas where fluid can become trapped and potentially cause carryover, improving accuracy, precision and analysis results.

Manual Autosampler or Instrument Syringes

If a syringe is being used by hand, a manual syringe should be selected. If a syringe is installed in an AOC autosampler then choose the appropriate syringe and volume to suit your instrument and application.

Shimadzu autosampler syringes are specifically designed to meet instrument dimensional specifications, have an accuracy of better than $\pm 1\%$ and are designed for precise, worry-free overnight sampling.

Needle Tip Styles

Cone: GC Autosampler



The cone shaped needle tip is specially developed to withstand multi injection demands and improve septum lifetime when used with the AOC autosampler. The cone design effectively "parts" the septum during piercing instead of cutting it, as would a bevel needle.

Bevel: Manual GC



The standard general purpose needle tip style supplied with many Shimadzu syringes is a 20° bevel tip. It is the preferred option for manual injection where piercing the septum in exactly the same place is difficult. The bevel tip is designed for optimum septum penetration and prevention of septum coring.

Dome:



This style needle is recommended for use with predrilled septa. The tip is rounded and polished to help septum penetration.

LC: HPLC



These needles are used for LC and HPLC valve injection and have a 90° square tip with rounded and polished edges. This eliminates damage to the valve's rotor seal and stator face. This needle tip style is a good choice for general liquid dispensing.

Side Hole Dome:



Samples are filled and dispensed through the side hole eliminating septum plugging of the needle. Ideal for large volume gas injection. The solid domed tip minimizes septum damage.

Valves



OPEN CLOSE

The push-button valve attaches directly to any luer lock 1 mL – 100 mL Shimadzu syringe.



OPEN CLOSE





The push-button valve attaches to any luer lock 5 mL – 100 mL Shimadzu syringe.

Syringe Volume	PTFE Tipped Plunger	Needle Length (mm)	Needle Gauge	Needle OD (mm)	Needle ID (mm)	Needle Tip	P/N	Spare Parts	P/N
Syringes for AOC-20i									
5 µL		42	23	0.63	0.11	Cone	221-75173		
10 µL	✓	42	23	0.63	0.11	Cone	221-74469	Plunger Pk 2	221-75173-01
10 µL		42	23	0.63	0.11	Cone	221-34618		
10 µL		42	23 – 26	0.63 – 0.47	0.11	Dome	221-37282-02		
10 µL	✓	42	23	0.63	0.11	Cone	221-75174	Plunger Pk 2	221-75174-02
								Needle Pk 2	221-75174-01
50 µL		42	23	0.63	0.24	Cone	221-45243		
250 µL	✓	42	23	0.63	0.24	Cone	221-45244	Plunger Pk 2	221-45244-01
MEPSTM Syringe for AOC-20i									
100 µL	✓						221-74830-01	Plunger Pk 1	221-74830-10
Syringes for AOC-5000									
10 µL		50	26	0.47	0.11	Cone	221-75175		
10 µL		50	23	0.63	0.11	Cone	221-75175-01		
10 µL		80	26	0.47	0.11	Cone	221-75176		
10 µL		80	23	0.63	0.11	Cone	221-75176-01		
10 µL	✓	50	26	0.47	0.11	Cone	221-75175-02		
10 µL	✓	50	23	0.63	0.11	Cone	221-75175-03		
25 µL	✓	50	23	0.63	0.24	Cone	221-75177		
1 mL	✓	50	26	0.47	0.15	Side-Hole Dome	221-75178		
1 mL	✓	50	23	0.63	0.15	Side-Hole Dome	221-75178-01		
2.5 mL	✓	50	26	0.47	0.15	Side-Hole Dome	221-75179		
2.5 mL	✓	50	23	0.63	0.15	Side-Hole Dome	221-75179-01		
MEPSTM Syringes for AOC-5000									
100 µL	✓						221-75202-01	Plunger Pk 1	221-75202-03
250 µL	✓						221-75202-02	Plunger Pk 1	221-75202-04
Manual Syringes									
5 µL		50	26	0.47	0.11	Bevel	221-75170		
10 µL		50	26	0.47	0.11	Bevel	670-12552-01		
10 µL	✓	50	26	0.47	0.11	Bevel	221-75170-01	Plunger Pk 2	221-75170-02
10 µL		51	22	0.028"	0.17	LC	670-12554-01		
25 µL		50	25	0.5	0.2	Bevel	670-12510-31		
25 µL	✓	50	25	0.5	0.2	Bevel	221-75171	Plunger Pk 1	221-75171-01
25 µL		51	22	0.028"	0.37	LC	670-12554-02		
50 µL		50	25	0.5	0.2	Bevel	670-12510-36		
50 µL	✓	50	25	0.5	0.2	Bevel	221-75172	Plunger Pk 1	221-75172-03
50 µL		51	22	0.028"	0.37	LC	670-12554-03		
100 µL		50	25	0.5	0.2	Bevel	670-12510-18		
100 µL	✓	50	25	0.5	0.2	Bevel	221-75172-01	Plunger Pk 1	221-75172-04
100 µL		51	22	0.028"	0.37	LC	670-12554-04		
250 µL		50	25	0.5	0.2	Bevel	670-12510-19		
250 µL	✓	50	25	0.5	0.2	Bevel	221-75172-02	Plunger Pk 1	221-75172-05
250 µL		51	22	0.028"	0.37	LC	670-12554-05		
500 µL		50	25	0.5	0.2	Bevel	670-12510-20		
500 µL		51	22	0.028"	0.37	LC	670-12554-06		
1 mL	✓	Luer Lock					221-54778-01	Plunger Pk 1	221-54778-11
5 mL	✓	Luer Lock					221-54778-02	Plunger Pk 1	221-54778-12
10 mL	✓	Luer Lock					221-54778-03	Plunger Pk 1	221-54778-13
25 mL	✓	Luer Lock					221-54778-04	Plunger Pk 1	221-54778-14
50 mL	✓	Luer Lock					221-54778-05	Plunger Pk 1	221-54778-15
100 mL	✓	Luer Lock					221-54778-06	Plunger Pk 1	221-54778-16
Luer Lock Needles									
		50	23	0.63	0.32	Bevel	221-54778-51		
		50	19	1.07	0.65	Bevel	221-54778-52		
		50	14	2.1	1.6	Bevel	221-54778-54		
Syringe Valves									
5 mL to 100 mL		Luer Lock					221-54778-49		
1 mL to 100 mL		Luer Lock					221-54778-50		

SAMPLE INTRODUCTION – GC Inlet Liners





The GC inlet functions as the interface between the syringe and the GC capillary column, where the sample is introduced, vaporized, mixed with carrier gas and transferred to the column. Shimadzu instruments offer several types of inlets - split, splitless, programmable temperature vaporization (PTV) and on-column.

The inlet liner prevents the sample contacting the metal walls of the injector block. Inlet liner geometry and packing materials enable the inlet liner to achieve greater heated surface area; this additional surface area can often improve sample vaporization. Conversely, choosing the wrong inlet liner geometry can significantly decrease the reproducibility and quality of analysis.

Injection Technique	Sample Types	Inlet Liner Geometry	Function
Splitless	Trace Level Analyses/ Active Compounds	Taper 	A bottom taper focuses sample onto the head of the column and minimizes sample contact with metal parts of the inlet. The addition of quartz wool to your inlet liner promotes mixing of analytes, aids the vaporization of liquid samples, and works as a trap to collect non-volatile residue in the sample (i.e. protects capillary column from 'dirty' samples).
Split	General Purpose/ Concentrated Samples/ Dirty Samples	FocusLiner™ 	Ensures quartz wool remains in the correct position in the liner. Excellent reproducibility results from the wiping of the sample from the syringe needle and the prevention of droplet formation. Minimizes high molecular weight discrimination.
Splitless	Trace Level Analyses/ Dirty Samples/ Wide Boiling Point Range	Taper FocusLiner™ 	Bottom taper focuses sample onto the head of the column and minimizes contact with metal parts of the inlet. Excellent reproducibility results from the wiping of the sample from the syringe needle and the prevention of droplet formation. Minimizes high molecular weight discrimination.
Direct	Trace Level Analyses/ Active Compounds	Direct Taper 	Direct inlet liners facilitate maximum transfer of sample by connecting directly to the GC column and inhibiting sample degradation due to hot metal components inside inlet.
Split/Splitless	General Purpose/ Concentrated Samples/ Dirty Samples (only if quartz wool is present)/ Gaseous Samples (also purge & trap, headspace)	Straight 	Straight inlet liners facilitate higher split flows. Narrow bore straight inlet liners facilitate fast GC work. Small injection volumes of less than 0.5 µL are best used with a narrow bore. Narrow bore straight inlet liners improve focussing of gaseous samples (purge, trap & headspace).

Inlet Liner Deactivation

Every batch of inlet liners are tested for inertness using the EPA 8081B method. This standard method ensures that each batch of inlet liners has less than 3 % Endrin breakdown from a 1 ppm injection.

Description & Geometry Sketch	OD (mm)	ID (mm)	Length (mm)	Pack Size	P/N
GC-2010					
LINERGC20103.4MMIDGNWWOOLPKT5 	5.0	3.4	95	5	221-75193
LNRGC20103.4MMIDGNPKT5 	5.0	3.4	95	5	221-75194
GC-2014					
LNRGC20143.4MMIDGNPKT5 	5.0	3.4	95	5	221-41444-05
LNRGC20142.6MMIDPKT5 	5.0	2.6	95	5	221-41544-05
LNRGC20142.6MMIDTAPPKT5 	5.0	2.6	95	5	221-41599-05
LINERGC20143.4MMIDWWOOLGNPKT5 	5.0	3.4	95	5	221-75195
GC-2010/GC-2014					
LNRGC2010/20143.4MMIDTAPFOC(ontowool)PKT5 	5.0	3.4	95	5	221-75187
LNRGC2010/20143.4MMIDFOC(ontowool)PKT5 	5.0	3.4	95	5	221-75188
LNRGC2010/20143.4MMIDFOC(intowool)PKT5 	5.0	3.4	95	5	221-75189
LNRGC2010/20143.4MMIDPKT5(StraightThrough) 	5.0	3.4	95	5	221-75190
LNRGC2010/20143.4MMIDTAPFOC(intowool)PKT5 	5.0	3.4	95	5	221-75191
LNRGC2010/2014SINGLETAPPKT5 	5.0	3.4	95	5	221-75192
LNRGC2010/20140.75MMIDSPMEPKT5 	5.0	0.75	95	5	221-75196
GLASSINSERTSPLITLESS/WBIWWOOLPKT5 	5.0	3.4	95	5	221-75197

SAMPLE INTRODUCTION – Septa



- High temperature silicone
- Excellent durability, resealing and solvent and tear resistant
- Injection temperature up to 350 °C

Septum Type	Description	Pack Size	P/N
GC-2010/GC-2014			
Enduro Blue	High thermal stable material in blue color	50	221-75180

GC CONSUMABLES

SAMPLE INTRODUCTION – Accessories

Additional Accessories

Description	P/N
Capillary Ceramic Tube Cutter (3 pc)	221-75181
Stainless Steel Nut/ SSNE-16-012S (Pk 5)	670-11009

SAMPLE PREPARATION

SAMPLE INTRODUCTION

SEPARATION

TROUBLESHOOTING

SEPARATION - Connections

Ferrules are available in a variety of different materials, shapes and sizes depending on their use, the instrument and the size of the capillary column being used. Probably the most important but difficult aspect of choosing a ferrule is the selection of the material type. The table below will help you choose the appropriate ferrule material for your application.

When choosing ferrules ensure you consider the following:

- 1) The material that best suits your application.
- 2) The connection type you want.

The following selection table will assist with your decision.

Ferrule Material Type	Graphite	Graphite Vespel®	SilTite™ Metal
			
Features	<ul style="list-style-type: none">• Easy to use.• Forms a stable seal.• Soft material.• Porous to oxygen.• Can be reused.• Forms a soft grip with capillary column.• Low emissions.	<ul style="list-style-type: none">• A composite of graphite and Vespel®.• Mechanically robust.• Hard material, long lifetime.• Forms a strong grip with capillary column.• Cannot be reused with another capillary column.• Requires re-tightening.	<ul style="list-style-type: none">• Specifically developed to overcome the problems associated with the use of 100% graphite and composite ferrules.• Strong seal on capillary columns.• Leak free - The ferrule and nut expand and contract at the same rate eliminating any chance of leaks with temperature cycling.• Nut does not need re-tightening after initial temperature cycles.
Suitable Uses	<ul style="list-style-type: none">• Column to injector connection.• Non-mass spectrometer detectors (FID, ECD, TCD and NPD).	MS interfaces, although even with a good seal will leak air compared to SilTite™ ferrules.	Ideal for MS interfaces due to leak free seal.
Not Suitable For	Connecting columns to mass spectrometers, as porous to oxygen.	High temperature applications.	
Risks	<ul style="list-style-type: none">• Can leave residue inside your column.• Can extrude into the injector or detector if it is over-tightened.	If not re-tightened after installation and temperature cycles of the GC, air may enter the column or detector decreasing sensitivity of the analysis and possibly degrading the column as well as components of the system.	Over-tightening of the seal can introduce leaks into the system. Follow the recommended installation instructions to avoid this problem.
Operating Temperature	Upper limit of 450 °C	Upper limit of 325 °C	No temperature limit in GC use.

Ferrule Material Type	Ferrule Code/Description	Pack Size	P/N
Graphite	GFF-505-050	10	221-75182
Graphite Vespel®	GVF16-005	10	670-15003-04
	GVF16-008	10	670-15003-07
	GVF16-004	10	670-15003-03
Metal	SilTite™ 0.25 mm ID Column	10	221-72563-04
	SilTite™ 0.32 mm ID Column	10	221-72563-05
	SilTite™ 0.53 mm ID Column	10	221-72563-08
	SilTite™ 1/32" ID Column	10	221-75200-04
	SilTite™ Kit 10/32" 0.25 mm ID Column	1	221-75200
	SilTite™ Kit 10/32" 0.32 mm ID Column	1	221-75200-01
	SilTite™ Kit 10/32" 0.53 mm ID Column	1	221-75200-02
	SilTite™ Kit 1/32" ID Column	1	221-75200-03
	SilTite™ Nut 10/32" 0.8 mm ID Column	5	221-75186

TROUBLESHOOTING

The purpose of this information is to help you troubleshoot the performance of your chromatography - your system manual is an excellent guide to help you troubleshoot the performance of the system.

The separation of structurally diverse analytes is often complicated by chance coelutions with other analytes or with matrix related compounds. Often the column is blamed, but while such coelutions make analysis difficult, they do not necessarily indicate a faulty column, poor chromatography or method design. No single column or method will resolve all compounds that can be chromatographed, so selecting a column that matches the needs of the application is an important first step. Rather than attempting to modify a method to resolve coeluting peaks, selecting a column with subtly different selectivity can achieve this aim without significant changes to established methods, elution orders or run times.

The following are good starting points to not only develop the ideal method for your chromatography, but also a guide on where you can troubleshoot to improve your desired separation.

Phase Selection


- Select the least polar phase that will perform the separation you require. Non-polar stationary phases separate analytes predominantly by order of boiling point. Increase the amount of phenyl and/or cyanopropyl content in the phase, and the separation is then influenced more by differences in dipole moments or charge distributions.
- To separate compounds that differ more in their hydrogen bonding capacities (for example aldehydes and alcohols), polyethylene glycol type phases are best suited, such as CBP-20.

Column Diameter

- The smaller the diameter, the greater the efficiency, and better resolution. Fast columns (0.1 mm ID) are used for faster analysis because the same resolution can be achieved in a shorter time.

Film Thickness

- For samples with a variation in solute concentration, a thicker film column is recommended. This will reduce the possibility of broad overloaded peaks coeluting with other compounds of interest. If the separation of two solutes is sufficient and co-elution is still unlikely, even with large differences in concentration, then a thinner film can be used.
- The greater the film thickness the greater the retention of solutes, therefore the higher the elution temperature. As a rule, doubling the film thickness results in an increase in elution temperature of approximately 15-20 °C under isothermal conditions. Using a temperature program, the increase in elution temperature is slightly less.



- Maintain phase ratio among different ID columns to yield similar chromatography.
- Columns should be conditioned to the maximum continuous temperature unless specified. When conditioning columns with a film thickness > 1 µm at the maximum operation temperature, it is recommended to do the initial conditioning without a connection to the detector to minimize contamination from the siloxane bleed.

$$\beta = \frac{id}{4d_f}$$

where

β = phase ratio


id = column internal diameter (µm)

d_f = film thickness (µm)

Formula to calculate Phase Ratio.

Column Length

- Always try to select the shortest column length that will provide the required resolution for the application. If the maximum column length available is being used and resolution of the sample mixture is still inadequate then try changing the stationary phase or internal diameter.
- Resolution is proportional to the square root of the column efficiency; therefore, doubling the column length will only increase the resolving power of the column by approximately 40%.



Effect of Length.

Injection

The function of the GC sample inlet is to introduce a representative portion of the sample as a narrow band onto the chromatographic column – failure to achieve this objective will significantly reduce the separation capability of the GC column. Because most samples are liquids, an essential feature of the common GC inlet type is that the sample and solvent are vaporized prior to reaching the column. These GC inlets are known as ‘Vaporizing’ injectors. Below are several tips on the injector but you should also consider the function of the injection process and how it influences the chromatography. If the injection is not performed correctly, you will end up with poor accuracy and poor precision. If you have an autosampler fitted, you should follow the recommendations in your user manual to determine whether the autosampler is performing to specification.


The Syringe



- Syringe life is influenced by the quality of the sample and the frequency of use – a blunt needle will potentially deposit septa material into the inlet liner and severely impact your chromatography.
- For the best syringe life, ensure your syringe is rinsed five to ten times with clean solvent after use. Syringe washing also helps to eliminate carryover but remember to discard the initial washes.
- For optimal reproducibility and accuracy, the smallest volume injected from any syringe should be no less than 10 % of its total capacity, for example: the smallest recommended injection volume from a 10 µL syringe would be 1 µL.

The Injector

- If you are doing trace analysis work or working with high injection port temperatures, replace the septum regularly.
- If the injection port temperature is not specified in the method, 250 °C is usually the recommended temperature.
- Choose the appropriate inlet liner for your application – review the inlet liner geometries in this Product Guide (page 11) and suggested applications.
- If the inlet liner is not specified, it is recommended to use a FocusLiner™ with quartz wool and bottom taper.
 - The wool provides additional surface area for complete volatilization of the sample to minimize sample discrimination.
 - It traps non-volatile components and septum particles from reaching the column.
 - It wipes any sample from the syringe needle, thereby increasing reproducibility and preventing sample residue build-up at the septum.
- When choosing the inlet liner ensure the inlet liner volume is larger than the volume of vaporized sample. If the vapor volume exceeds the liner's inner volume, there will be an overflow of sample vapor from the liner, resulting in contamination of the inlet system which in turn leads to carryover (ghost peaks) and poor run-to-run reproducibility.
- The inlet liner is a consumable and its life is influenced by the quality of the sample and the frequency of use – once the inlet liner is contaminated with septa particles, it needs to be replaced. Recommended replacement is every 100 - 200 injections.
- When peak shape deteriorates, replace the inlet liner immediately and remove ~ 30 cm from the front of the capillary column.
- When replacing the inlet liner it is the ideal time to replace the septa and the inlet liner sealing O-ring.



Liquid Volume (µL)	Inlet Temperature (°C)	Inlet Pressure (psi)	Gas Volume (µL)
1	250	10	399
2	250	10	798
2	250	34.7	399
5	200	45	746

Dichloromethane as solvent

GC Connections

Ferrules are the smallest, inexpensive and probably the most easily forgotten components that are used in every gas chromatograph. Yet, without ferrules, the leak-free sealing that is required at the detector and injector of a GC system would be impossible to achieve.

- Always cut the column **after** passing through a graphite ferrule – this eliminates the risk of graphite particles entering into the column and impacting your chromatography.
- Always re-tighten Graphite / Vespel® ferrules by ¼ to ½ turn after the first 2-3 oven temperature cycles – Graphite / Vespel® ferrules tend to “creep” during temperature cycling.
- For GC-MS connections always use SilTite™ metal ferrules to ensure a permanent leak-free connection.

High Baseline Level	Cause	Poor carrier gas quality
	Solution	Use high quality gas, install gas traps. Replace septa and insert a new inlet liner
	Cause	Exceeded the maximum temperature of the capillary column phase
	Solution	Lower the programmed maximum temperature
Loss of Sensitivity	Cause	System has become active
	Solution	Replace the inlet liner with a new deactivated liner, cut 50 cm off the front of the column and re-install
	Cause	Splitless conditions have changed
	Solution	Re-check solvent and method temperatures. Consider introducing an internal standard.
Loss of Peak Resolution	Cause	Manual injection technique or operator has changed
	Solution	Ensure technique is consistent
	Cause	Dirty column - if samples are dirty, non-volatile materials have deposited onto the column causing a change in polarity
	Solution	Cut 50 cm off the front of the column and re-install
Split Peaks	Cause	Poor injection technique
	Solution	Increase the manual plunger depression speed
	Cause	The column has been inserted too far into the injector
	Solution	Reposition the column according to the manufacturers recommendation
Fronting Peaks	Cause	Too much sample injected onto column
	Solution	Dilute the sample or increase the split ratio
Tailing Peaks	Cause	Glass inlet liner has active surface
	Solution	Replace with a fully deactivated inlet liner
	Cause	Graphite ferrule contamination in the start of the column
	Solution	Remove 5 cm off the front of the column and re-install
Broad Peaks	Cause	Split gas flow is too low
	Solution	Increase split flow or use the 'solvent effect' to focus peaks
	Cause	The column is contaminated
	Solution	Cut 50 cm off the front of the column and re-install
	Cause	Mass spectrometer sampling rate is too low
	Solution	Increase sampling rate or reduce the number of ions detected in SIM mode
	Cause	Co-elution of peaks
	Solution	Change the temperature program or polarity of the column
Shifting Retention Times	Cause	Leaking septum
	Solution	Tighten septum cap or replace septa
	Cause	Dirty column - if samples are dirty, non-volatile materials have deposited onto the column causing a change in polarity
	Solution	Cut 50 cm off the front of the column and re-install
No Peaks	Cause	Syringe is blocked or leaking around plunger
	Solution	Inspect syringe and change if damaged
	Cause	Column is blocked
	Solution	Cut 5 cm off the front and back ends of the column and re-install
Ghost peaks are appearing	Cause	Syringe has become contaminated from previous run
	Solution	Ensure syringe has been thoroughly washed with solvent between injections.
	Cause	Septum bleed
	Solution	Will appear as discrete peaks in a thermal gradient and disappear in an isothermal run. Replace septum
	Cause	A backflush event has occurred
	Solution	Inject twice the amount of pure solvent (repeat if necessary)



Shimadzu Corporation

www.shimadzu.com/an/

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

For Research Use Only. Not for use in diagnostic procedures.
The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

© Shimadzu Corporation, 2013

Printed in Japan 3655-05331-10AIT