

User's Guide and Reference Manual

Version 6.10

Jack T. Holladay, M.D., M.S.E.E., F.A.C.S.

HICSOAP Professional and Deluxe User's Manual

Contents

1.	Introduction.....	4
2.	System Requirements	5
3.	Getting Started	6
4.	Program and Surgeon Setup	7
4.1	Program Setup.....	7
4.2	Surgeon Lens Setup	8
4.3	Surgeon Piggy-Back Lens And Toric Calc Setup.....	10
4.4	Advanced Setup.....	13
5.	Main Menu	14
6.	Patient Screens.....	16
6.1	Patient Locator	16
6.2	PreOp Data	18
4.2.1	'Other Calculations'	20
4.2.2	Correct AL.....	21
4.2.3	Alternate K.....	21
6.3	IOL Calculations	26
6.4	Toric IOL Planner	27
5.4.1	Toric Lens Incision Optimizer	29
6.5	Surgery Data Entry.....	30
6.6	PostOp Data.....	31
6.7	Toric PostOp BackCalc	32
6.8	Individual Surgically Induced Refractive Change (SIRC)	32
6.9	Individual Back Caclulations	34
7.	Surgical Outcomes Analysis Program (SOAP).....	35
7.1	SOAP Prediction Error	37
5.1.1	Frequency Distribution (%) of Eyes	37
5.1.2	Cumulative Distribution (%) of Eyes.....	37
5.1.3	Summary Statistics	38
7.2	SOAP-SIRC Changes Report	38
5.2.1	Surgically Induced Spherical Equivalent (SIRC) – Upper left graph.....	39
5.2.2	Actual PreOp Cylinder and Axis – Lower left graph.....	39
5.2.3	Surgically Induced Cylinder and Axis (SIRC) – Upper right graph.....	39
5.2.4	Actual PostOp Cylinder and Axis – Lower right graph.....	39
5.2.5	Buttons under Graphs	39

HICSOAP Professional and Deluxe User's Manual

7.3	SOAP Potential Outcomes	40
5.3.1	Frequency Distribution (%) of Eyes – Upper Left Graph	40
5.3.2	Cumulative Distribution (%) of Eyes – Upper Right Graph.....	40
5.3.3	Frequency Distribution (%) of VA – Lower Left Graph	40
5.3.4	Cumulative Distribution (%) of VA – Lower Right Graph.....	40
5.3.5	Buttons under graphs.....	41
8.	Utilities	42
8.1	Database.....	42
8.2	Documentation.....	43
8.3	Remove Registration	44
A.	Sample Printouts	45
A.1	IOL Calc Report.....	45
A.2	Toric IOL Report.....	46

HICSOAP Professional and Deluxe User's Manual

1. Introduction

This document describes both the professional and the deluxe editions. They are identical except that the professional edition has additional functionality not found in the deluxe edition:

1. Networking: Network setup is simple with all of the security protections now present no longer causing problems for multiple users on the same station.
2. Multiple Users accessing the same database simultaneously..
3. Network downloads from biometers such as the IOL MASTER (Zeiss) or LenStar (HAAG STREIT)
4. Run under CompuLink and other practice management systems.

2. System Requirements

The program runs on any IBM compatible computer that is running Windows XP or newer.

For best results, the Display should be set to a Desktop area of no less than **1024 X 768** pixels (or larger) with the Font Size set to Small Fonts. If a smaller desktop area (640 X 480) or Large Fonts are chosen, then the entire screen cannot be seen at one time and scroll bars are necessary.

3. Getting Started

This section has moved to a new (improved and expanded) document: **HICSOAP Getting Started**

4. Program and Surgeon Setup

4.1 Program Setup

Holladay IOL Consultant Setup

Program/Clinic Setup | Surgeon Lens Setup | Surgeon Piggy-back Lens Setup | Advanced Settings

This software is registered for use at the following Clinic:

Test Clinic

Technicians

Visual Acuity Standard:

Verify your Computer's Clock

Printout alignment

Keratometer Settings:

Other Clinic Settings

OK Help

- **Clinic Name:** The 'clinic name' will be highlighted in blue and can be changed at any time.
- **Technicians:** A technician's name may be added and edited. PreOp records allow a technician to be entered.
- **Visual Acuity Notation:** Select from the four choices for visual acuity notation: **feet, metric, decimal or LogMAR.**
- **Verifying Correct Computer Clock Time:** The clock is used to calculate age, time between visits, etc. **Please make sure the computer clock time is correct at all times.**
- **Printer Alignment Adjustment:** The printed reports can be shifted on the page by using this feature. This feature is valuable if the left margin or top margin needs to be larger because the printout is placed in a chart with holes punched on the left side or top of the page. Positive values shift the printout down or to the right and negative values shift the printout up or to the left.
- **Keratometer Settings:** The refractive index used in the clinic keratometers must be selected.
NOTE: Both IOLMaster and Lenstar export the index of the refraction which is used by the program.

HICSOAP Professional and Deluxe User's Manual

- **Default Surgeon:** A default surgeon for the PreOp exam may be entered here. When a default Surgeon (other than None) is selected, then this Surgeon is entered automatically for new PreOp exams only.
- **Default Technician:** A default technician for the PreOp exam may be entered here. When a default Technician (other than None) is selected, then this Technician is entered automatically for new PreOp exams only.
- **Vertex Distance (VTX) Default:** The default value for the vertex distance of all refractions may be entered here. The average value for spectacles and trial frames is **12 mm**. The average vertex distance for the phoropter is **14 mm**. If the refraction exceeds ± 4 diopters, the vertex should be measured to avoid computation errors.

4.2 Surgeon Lens Setup

The screenshot shows the 'Holladay IOL Consultant Setup' dialog box with the 'Surgeon Lens Setup' tab selected. The interface includes a 'Surgeon' dropdown menu with 'Surgeon1, Example' selected, and 'Edit...' and 'New...' buttons. Below this is a spinner for 'Total Number of Lenses to appear in your IOL Power Calc report:' set to 4. A 'Search...' button is next to the text 'Select Lens for Active Lens List, or'. On the left, there are dropdowns for 'Lens Manufacturer:' (ACRIMED), 'Lens Model:' (11 C-11-BC), and 'Procedure:' (Std Phaco), along with 'Add >>' and '<< Remove' buttons. A text area for 'Manufacturer's Comments:' is at the bottom left. The 'Active Lens List' is a table with columns for ID, Manufacturer, Model, and Formula. At the bottom of the list are 'Move Up', 'Move Down', and 'Enter Constant' buttons. A large 'Add Lens not in data base' button is at the bottom center. The dialog has 'OK', 'Cancel', and 'Help' buttons at the bottom right.

ID	Manufacturer	Model	Formula
1	AMO	Veriseye VRSM5US	Phakic
2	B&L	P506UV	Std Phaco
3	Alcon	SN60T3-4-5	Std Phaco
4	AMO	Veriseye VRSM6US	Phakic
	AMO	SI-30NB	Std Phaco
	Alcon	MA30BA	Std Phaco
	B&L	Baikoff3	Phakic ACL
	B&L	L122UV	AC Backup
	Staar	Fyodorov1	Phakic ICL

- **New/Edit Surgeon:** Click the **New** button to add a **New surgeon** or the **edit** button to edit an existing surgeon. The **Last Name**, **First Name** and **Primary Formula** must be entered (**Holladay 2 Formula is recommended**). If the **Active** box is checked this surgeon will appear in the active group of surgeons, alphabetized at the beginning of the surgeon pull down list. **Inactive** surgeons appear alphabetically in a list after the **active** surgeons.
- **Number of lenses to be printed on IOL Calculation report from active lens list:** The program can print up to **4 IOL models** on each page on the IOL calculation screen and on a single page report.

HICSOAP Professional and Deluxe User's Manual

Any number of lenses may be chosen and the default is 4. **If more than 4 lenses are on the active lens list then multiple pages will be printed with 4 lenses per page.** There is no limit to the number of pages that may be printed.

- **Adding lenses to surgeon's Active Lens List:** To add a lens to the surgeon's **Active Lens List**, choose the manufacturer and the lens model from the pull down screens. The **procedure** defaults to **Std Phaco**, but may be changed to any description. For example, one may choose **AC Backup** for an anterior chamber backup lens or **ECCE** for a lens that is used for a large incision extracapsular cataract extraction. The procedure code is particularly valuable when the same lens model is used for different procedures, requiring a different lens constant. *If the lens constant is the same, there is no need to create a separate procedure group.* Once the Manufacturer, Lens Model, and Procedure have been chosen, click the **Add>>** button to place the lens in the **Active Lens List**. There is no limit to the number of active lenses that can be added. All Active Lenses will be displayed and printed with four lenses per page.
- **Removing lenses from the surgeon's Active Lens list: A Lens Model with no data may be removed from the surgeon's Active Lens List. If the Lens Model has data, it cannot be removed.** To remove a lens model with no data, highlight the lens with a single click, then click the remove button.
- **Changing the lens order number of printed lenses:** The order of lenses selected for printing may be rearranged by using the **Move Up/Move Down** buttons below the **Active Lens List**. Simply highlight a lens move it up or down to the desired position in the list.
- **Over-riding the Manufacturer's lens constant with a Surgeon-Entered lens constant:**

Surgeon-Entered Lens Constant

Surgeon: Surgeon1, Example
Lens Manuf.: Alcon/Cilco
Lens Model: SN60T3-4-5
Procedure: Std Phaco

Please enter your constant:

	Constant for Ultrasound AL	Constant for Optical Biometer AL
Holladay 2 or R	<input type="text"/> ACD	<input type="text"/> ACD
Holladay 1	<input type="text"/> SF	<input type="text"/> SF
SRK-T, I, & II	<input type="text"/> A-Constant	<input type="text"/> A-Constant
Hoffer-Q	<input type="text"/> ACD	<input type="text"/> ACD
Binkhorst I & II	<input type="text"/> ACD	<input type="text"/> ACD

If the A-constant and ACD do not agree, the A-constant is usually the better value and is for in-the-bag placement.
When they do not agree the manufacturer's ACD is usually ~0.30 mm less than the value converted from the A-constant and is actually for in-the-sulcus.

OK Cancel Help

The program automatically retrieves the manufacturer's lens constant from the Data base for IOL

HICSOAP Professional and Deluxe User's Manual

calculations. If the surgeon already has a personalized lens constant, he may override the manufacturer's lens constant by clicking the **'Enter Constant'**. Select the formula that the surgeon has been using and enter the personalized constant. This value overrides the manufacturer's constant and will be displayed on the IOL calculation screen and printout as **"Entered Cst." in place of "Manuf. Cst."**. The program automatically converts the **Surgeon-Entered** personalized constant for any of the formulas to the appropriate constant for the formula that has been chosen on the IOL Calculation Screen. **If the A-constant and ACD from the manufacturer do not agree, the A-constant is usually the better value and is for in-the-bag placement. When they do not agree the manufacturer's ACD is usually ~0.30 mm less than the value converted from the A-constant and is actually for in-the-sulcus.** A **"Personalized Lens Constant" (PLC)** is calculated after each postoperative data entry and overrides either of the previous constants when statistically significant (~ 30 to 40 cases). **PLC** will appear on the printout and the **PLC constant** may never be overridden.

- **View/Edit Lens PROPERTIES in Data base:** All of the properties of existing IOL's in the Data base may be viewed by clicking on the View/Edit Lens PROPERTIES in Data base button. Select the manufacturer and lens model from the pull down lists and the properties of the lens are displayed. This feature is particularly useful for determining the available powers or characteristics of a specific lens model.
- **Add lenses not in Data base:** Although every lens available at the time of printing is currently in the Data base, it is possible to Add Lens not in Data base by clicking this button. All of the lens properties must be completed, except for a lens constant. When one lens constant is entered, the other two fill in automatically. This feature is especially helpful for lens constant conversions. If the A-constant and ACD do not agree, the A-constant is usually the better value and is for in-the-bag placement. When they do not agree the manufacturer's ACD is usually ~0.30 mm less than the value converted from the A-constant and is actually for in-the-sulcus. The other properties are used in various calculations, such as optic diameter and vault, so all properties must be entered before a lens can be added. Lenses "Added" to the Data base by the user can be "edited" at any time on this Screen. Editing lens properties in the program's "ORIGINAL" DATA BASE can never be performed, ONLY LENSES THAT HAVE BEEN ADDED BY THE USER. (Exception: The Range of available dioptric powers may be changed by the user on any lens in the Database - original or added). Note: The format for the range of power of an IOL is min power, max power and dioptric steps; [e.g. (10.0,25.0,0.5;25.0,34.0,1.0) indicates that the lens is available from 10.0 to 25.0 in 0.5 D steps and from 25.0 to 34.0 in 1.0 D steps]. The min and max power are separated by commas and the step size is followed by a semicolon.
- **Search Lens Database:** The Lens Search button will display the lens search dialog, when the lens database may be searched based on various criteria.

4.3 Surgeon Piggy-Back Lens And Toric Calc Setup

Piggy-Back lenses are used when the required implant exceeds the power available for a given lens style. The choices below direct the program on how to handle this situation. This choice can be changed at any time on this screen by choosing the Lens Setup button on any screen. If both lenses are to be placed in-the-bag, then this should be chosen on the preoperative screen. If the posterior piggyback lens goes in-the-bag and the anterior piggy-back lens goes in-the-sulcus, then in-the-sulcus

HICSOAP Professional and Deluxe User's Manual


should be selected on the PreOp screen. **Note: Piggy-Back lenses are calculated using the Holladay 2 Formula ONLY.**

The screenshot shows the 'Holladay IOL Consultant Setup' dialog box with the 'Surgeon Piggy And Toric Lens Setup' tab selected. The 'Surgeon' dropdown menu is set to 'Surgeon1, Example'. The 'Piggy preferences' section contains a note: '*** NOTE: Piggy-Back lenses are calculated using the HOLLADAY 2 FORMULA ONLY ! ***'. Below the note are six radio button options: 'No piggyback lens', 'Fix Posterior lens power to max available power' (which is selected), 'Fix Posterior lens power to' (with an empty text box and 'D' label), 'Use two equal IOL's when power exceeds max available power', 'Use two equal IOL's when power exceeds' (with an empty text box and 'D' label), and 'Fix Front IOL to' (with an empty text box and 'D' label). The 'Toric Calc preferences' section includes 'Residual Refraction Mode' with 'Plus' selected and 'Minus' unselected, 'Print Inverted text on the Toric report:' with a checked checkbox, and 'Surgery Induced Astigmatism Temporal:' set to '0.2 D' and 'Superior:' set to '0.3 D'. At the bottom right are 'OK', 'Cancel', and 'Help' buttons.

- **No Piggy-Back Lens:** The surgeon does not want piggy-back calculations performed.
- **Fix Posterior lens power to max available power:** The surgeon wants the program to use the maximum available power from the manufacturer for a given lens model to be used as the posterior piggy-back lens. This is the **standard choice**, since most experts recommend placing the highest available power for the posterior lens and the remaining power for the anterior lens. Since the maximum power for each lens model may be different, the posterior piggy-back lens power on the screen and printout may be different for each lens model.
- **Fix Posterior Lens Power to:** Same as number 2, but the user has chosen to fix the power of the posterior IOL to a value other than the maximum available power. This selection is useful if the user wants to override the maximum power of a lens model in the Data base. For example, if the manufacturer makes a 34 diopter lens, but the maximum power available in the current inventory is 32, the value can be set to 32. **This value will be used for all lens models in the IOL calculation.**
- **Use Two Equal IOL's when power exceeds max available:** The surgeon wants piggy-back calculations performed using two lenses that are as close to equal power as possible and would like this to occur when the necessary lens power exceeds the maximum available power from the manufacturer.


HICSOAP Professional and Deluxe User's Manual

- **Use Two Equal IOL's when power exceeds:** Same as 4 above, but the maximum value is set by the user.
- **Fix Front IOL to:** The surgeon has chosen to fix the power of the anterior IOL. This is often helpful when the front lens required is very low, i.e., 1 to 5 diopters, and the lenses are not readily available. The front lens can be set to the lowest available power, such as 4 or 6 diopters.
- **Toric Calculation Residual Refraction Mode:** If Plus is selected, then all residual refractions in the toric reports will be using the plus cylinder mode. If Minus is selected, then they will be displayed/printed using the minus cylinder mode.
- **Toric Calculations Reversed Text:** If the checkbox is checked, then the Toric calculation report will print the IOL and Incision Angles using both normal and upside/down text:



- **Surgical Induced Astigmatism Temporal/Superior:** These values are used in the preop toric optimizer. Please refer to section 6.4 below for more details.


4.4 Advanced Setup



- **Database File Configuration:** The database file used by the program is specified. Normally the program will pick a default database file, and for most users that may be sufficient. Advanced users may move this file to a different location for networking or enterprise backup.
- **Automatic Database Backup:** The program maintains automatic backups over several days/weeks. The backup should not be turned off unless an alternate backup scheme is implanted in the clinic.
- **Soap Graph Printing Mode:** Some printers may not be capable of printing some of the SOAP graphs. If this is the case, the alternate or the 'screen dump' modes may be tried in that order.
- **Other Options:** Additional options are specified in this section and the default options would be appropriate for most clinics. It is strongly recommended to enable the first option (check for newer version), so that you always run the latest version of the software.

HICSOAP Professional and Deluxe User's Manual

5. Main Menu




- **Enter New Patient Data:** Displays the Patient Locator screen for data input of a new patient's demographic data.
- **Locate Patient Record:** Displays the Patient Locator screen to search for an existing patient.
- **Program & Surgeon Setup:** Displays the Program Setup, Surgeon Lens Setup and Surgeon Piggy-Back Lens Setup screens for initially setting up the program.
- **SOAP Formula and Lens Constant Analysis:** Displays the Formula and Lens Constant Analysis Table at the bottom of the screen.
- **SOAP Prediction Error Report:** Displays the Prediction Error Report Screen.
- **SOAP SIRC Changes Report:** Displays the SIRC Report screen.
- **SOAP Potential Outcomes:** Displays the Potential Outcomes screen.
- **Utilities Database:** Displays the Data base Operations screen for Backup and Restore Data base operations. If the shift key is depressed when the Data Operations Button is clicked, the Program and Data base Versions are displayed.
- **Utilities Documentation:** Displays the documentation list page

HICSOAP Professional and Deluxe User's Manual

- **Utilities Send Database to Tech Support:** Sends a copy of the database to tech support. Before sending the database all patient last names are removed.
- **Utilities Remove Registration:** The program registration may be removed from the computer.
- **Quit:** Displays the Exiting IOL screen for Backup & Exit or Exit when leaving the program. Daily backups are recommended or immediately after entering large amounts of data.
- **Help:** Displays the Help Screen for the Main Menu.

6. Patient Screens

6.1 Patient Locator



- **Show EXAMPLE Patients:** The Program has 3 EXAMPLE Surgeons and several EXAMPLE patients to illustrate many features of the program. The EXAMPLE patients may be “hidden” from the patient list by un-checking this box.
- **Add New Patient:** This button will popup the new patient dialog where the demographic data may be entered. Selecting the **PreOp Data** button displays the **PreOp Data** Screen for data entry.
- **Import from IOL Master or Lenstar:** Data that has been exported from the **IOL Master** to a 1.44 Meg floppy disk may be directly imported into the program.
- **Existing Patient List:** All of the existing patients can be found on the Patient Locator Screen. A patient may be located by Last Name, First Name or Chart Number. Once a patient has been selected, all of their PreOp, Surgical and PostOp visits are displayed in the window below. Any visit may be highlighted and when “double clicked” will display the visit.

HICSOAP Professional and Deluxe User's Manual

- **Edit Patient:** The Edit Patient button is activated by highlighting an item in the Patient List. If the Edit Patient button is clicked, the data entry screen for that item will be displayed and editing of the data allowed.
- **Delete Patient/Exam:** The delete buttons are activated by highlighting a Patient or Visit in the Patient Locator window. If the Delete button is clicked, the highlighted item will be deleted. The most recent entry for an eye must be deleted before any previous older data can be deleted.
- The active **Patient Locator, PreOp, IOL Calc, Toric IOL Calc, Surgery data, PostOp Exam, Back Calculation and SIRC Calculation Tabs** at the top of the screen display the respective data entry screens for the patient. If subsequent data has already been entered, the screen will be "read only". For example, if surgery data has already been entered for a patient's eye, the **PreOp Data** screen will not allow data entry ("read only"). **Data cannot be edited unless all subsequent data has been deleted.**
- **Export/Print Patient List:** Pressing one of these buttons displays the **Selection Criteria Screen**. After the selection criteria have been set the program generates a list of patients that satisfy all of the criteria selected which is exported to a .csv file or printed (**Example Report R – page 65**).

HICSOAP Professional and Deluxe User's Manual

6.2 PreOp Data

The screenshot displays the 'Holladay IOL Consultant - a' software window. The patient information at the top includes Name: Patient, Normal; ID: 11AA; DOB: 13-Oct-1924; Gender: F. The interface is divided into two main columns for the right eye (OD) and left eye (OS). Each column contains several sections: 'Calculation Type and Target Refraction' with radio buttons for 'Standard IOL (Cataract Removal)' and 'Other Calculations ->', and checkboxes for 'Bifocal IOL Power Calc.' and 'Toric IOL Power Calc.'; 'Pre existing condition' with checkboxes for 'Keratoconus', 'Silicone in Vitreous Cavity', 'Scleral Buckle', and 'Other Pathology'; 'Refraction and K's and Horizontal WTW' with input fields for Refractive Error (Ref), BCVA, UCVA, and Horizontal White-to-White distance; 'Axial Length (mm)' with input fields for AL(US) and AL(Optical), and radio buttons for 'UltraSound (US)' and 'Optical (IOLM, Lenstar)'; and 'Anterior Segment (mm)' with input fields for Phakic ACD and Phakic Lens Thickness. A 'Copy' button is located between the two eye columns. At the bottom, there are buttons for 'Main Menu', 'Help', 'Quit', 'Setup...', 'Print', '<- Pat. Locator', and 'IOL Calc. ->'. The status bar at the bottom left shows 'Ready' and the bottom right shows the date '10/26/20'.

This screen is used to enter all the PreOp information. Incorrect or missing items are highlighted in red, while optional values are highlighted in yellow.

The program allow for either standard IOL calculations (cataract removal) or advanced calculations grouped under the 'Other Calculations' button – see section below for more information.

It is necessary to enter all the pre-existing conditions for best accuracy. The following conditions influence the IOL Calculations and Personalization of the Lens Constant.

- **Keratoconus:** Necessary for the program to differentiate a normal steep K-reading in a small anterior segment from a patient with keratoconus.
- **Previous RK, PRK, ALK, LASIK:** When this box is checked the **Alternate K** button appears. Select the **Alternate K** button and the **Alternate Average K-value Source** screen appears – see below for details.
- **Silicone in Vitreous Cavity:** Changes the IOL calculations to compensate for **Silicone in the Vitreous Cavity**. Although **convexoplano IOLs (flat posterior surface) are recommended since they minimize this effect**, the program calculates for any lens selected including biconvex. The increase in lens power usually ranges from 3 to 6 diopters for convexoplano lenses and from 4

HICSOAP Professional and Deluxe User's Manual

to 10 diopters for biconvex lenses. **Only PMMA lenses are recommended.** **Note: Selecting this option does not correct the axial length for Silicone in the vitreous. Axial lengths measured with Silicone in the vitreous cavity are unreliable, even when corrected for ultrasound speed in silicone.** It is much more accurate to measure the eye without Silicone and use this value for the calculation. If both eyes have Silicone and no pre-Silicone axial lengths were measured, a useful approximation is to use 23.5 mm as a nominal axial length and add or subtract 1 mm for every 3 diopters of myopia or hyperopia, respectively.

- **Scleral Buckle:** A scleral buckle artificially lengthens the eye and makes the measured axial length inappropriate for use in calculating the size of the anterior segment.
- **Other Pathology:** Used for flagging a patient that is not to be included in the personalization of a lens constant. When this checkbox is pressed, a popup dialog is displayed that is used to enter more detailed information:

Other Pathology or Complications

No Yes (If yes, check all that apply)

<input type="checkbox"/> Wound leak	<input type="checkbox"/> Raised IOP requiring treatment
<input type="checkbox"/> Cystoid macular edema	<input type="checkbox"/> Pupillary block
<input type="checkbox"/> Macular degeneration	<input type="checkbox"/> Fibrin in pupil
<input type="checkbox"/> Optic atrophy	<input type="checkbox"/> Distorted pupil
<input type="checkbox"/> Hyphema	<input type="checkbox"/> Cortical remnants
<input type="checkbox"/> Flat anterior chamber	<input type="checkbox"/> Nuclear remnants
<input type="checkbox"/> Vitreous in anterior chamber	<input type="checkbox"/> IOL dislocation
<input type="checkbox"/> Vitreous to wound	<input type="checkbox"/> IOL optic decentration if > 1.0 mm
<input type="checkbox"/> Secondary glaucoma	If yes,
<input type="checkbox"/> Vitritis	<input type="checkbox"/> Optic atrophy
<input type="checkbox"/> Cyclitic membrane	<input type="checkbox"/> Retinal detachment
<input type="checkbox"/> Anterior synechiae	<input type="checkbox"/> Diabetic retinopathy
<input type="checkbox"/> Posterior synechiae	<input type="checkbox"/> Other
<input type="checkbox"/> Significant pigment deposits on IOL	

NOTE: Eyes with pathology will not use personalized constants

OK Cancel

NOTE: The Vertex Distance (VTX) must be specified for any refraction entered into the Data base. The Vertex Distance (VTX) Default may be set on the Program Setup Screen. The average value for

HICSOAP Professional and Deluxe User's Manual

spectacles and trial frames is 12 mm. The average vertex distance for the phoropter is 14 mm. If the refraction exceeds ± 4 diopters, the vertex should be measured to avoid computation errors.

4.2.1 'Other Calculations'

Selecting 'Other Calculations' enables the 'Click Here' button which, when clicked, pops up the following dialog when additional advance setting may be entered:

OS(Left) Pre Op Special Data

Phakic Eye
AC, PC or ICL Lens Calc. from Ref. (no Lens Removal)

Aphakic Eye (AC IOL always in Angle)
For PC lens, calc. for: in Bag in Sulcus

PseudoPhakic Eye (AC IOL only in Angle)
New PC lens: in Bag in Sulcus

Standard IOL or Exchange IOL (lens calc from AL)
 Secondary Piggy-Back IOL (lens calc from Ref.)

IOL Material in eye: Acrylic IOL Pwr in eye:

NOTE: IOL MATERIAL AND POWER of the lens in the eye must be entered, to activate "Correct AL" button.
Use a nominal value of 21D for the lens power if it is unknown.

For Aphakic and Pseudophakic eyes results will be improved if the Preop Phakic Refraction is entered. If the Preop Phakic refraction is unknown, please leave the following fields blank.

Preop Phakic Refraction: X VTX:

OK Cancel


- **Phakic Eye:** AC, PC & ILC Phakic IOL's.
- **Aphakic Eye:** AC or PC IOL's. Select whether secondary IOL is to be placed in Bag or in Sulcus.
- **PseudoPhakic Eye**
 - Select whether **New PC lens** is to be placed **in bag** or **in sulcus**.
 - Select **Exchange IOL Lens Calc** or **Secondary Piggy-Back IOL Lens Calc**.

HICSOAP Professional and Deluxe User's Manual

- Select the intended location of the **New PC lens, in Bag** or **in Sulcus**.
- Enter **IOL Material in Eye** and **IOL Power in eye**. If the **material and power** of the IOL in the eye are entered, the **Correct AL** button appears on the **PreOp Data Screen** that allows the surgeon to correct the pseudophakic axial length measurement.
- **PreOp Phakic Refraction**: If this refraction is available, the program will use it to improve the calculations

4.2.2 Correct AL

When the 'Correct AL' button is pressed the following screen is displayed:



Choose the velocity setting that corresponds to the velocity used to take the pseudophakic axial length measurement.

4.2.3 Alternate K

When the Previous RK, PRK, ALK, LASIK box is checked the Alternate K button appears. An alternate K may be calculated using different methods.

HICSOAP Professional and Deluxe User's Manual

Average K from Refractive data

Alternate Average K-value Source - Post-Keratorefractive Surgery Cornea

Select K-value to be used for IOL power Calculation:

Calculated average K-value from Refractive data 32.13

Calculated average K-value from Rigid Trial CL NA

Measured average K from preOp data screen 43.75

Surgeon-Entered K-Value NA

Average K from Refraction

a. Refraction and K-values prior to keratorefractive surgery (e.g., Rk, PRK, ALK, LASIK):

SPH	CYL	VTX	K1	K2
-12	-2	12	43.25	43.50

b. Stable Ref. after keratorefractive surgery, preferably before cataract formation:

SPH	CYL	VTX
-0.25	0.50	12

NOTE:

Lowest K value tends towards myopic error.

Highest K value tends towards hyperopic error.

OK Cancel Help

Allows the surgeon to calculate an average K-value from the data taken before and after the refractive surgery. The K-values before the refractive surgery and the refraction before and after the refractive surgery are needed.

HICSOAP Professional and Deluxe User's Manual

Average K from Rigid Trial Contact Lens

Alternate Average K-value Source - Post-Keratorefractive Surgery Cornea

Select K-value to be used for IOL power Calculation:

Calculated average K-value from Refractive data 32.13

Calculated average K-value from Rigid Trial CL 33.51

Measured average K from preOp data screen 43.75

Surgeon-Entered K-Value NA

Average K from Rigid Trial CL

a. Refraction without Rigid Contact Lens:

SPH	CYL	VTX
-0.50	-1.00	12

b. OverRefraction w/Rigid Contact Lens:

SPH	CYL	VTX
-1.50	0	12

c. Contact Lens BC in mm or D: 35 (D/mm)

d. Contact Lens Power: -1 D

NOTE:
Lowest K value tends towards myopic error.
Highest K value tends towards hyperopic error.

OK Cancel Help

Allows the surgeon to calculate the K-value from the base curve and power of a rigid trial contact lens and the refraction with and without the contact lens

HICSOAP Professional and Deluxe User's Manual

Measured Average K from PreOp Screen

Alternate Average K-value Source - Post-Keratorefractive Surgery Cornea

Select K-value to be used for IOL power Calculation:

- Calculated average K-value from Refractive data 32.13
- Calculated average K-value from Rigid Trial CL 33.51
- Measured average K from preOp data screen 43.75
- Surgeon-Entered K-Value NA

Measured average K: 43.75

NOTE:
Lowest K value tends towards myopic error.
Highest K value tends towards hyperopic error.

OK Cancel Help

Allows the surgeon to use the actual measured values from the PreOp Data Screen

HICSOAP Professional and Deluxe User's Manual

Surgeon Entered K

Alternate Average K-value Source - Post-Keratorefractive Surgery Cornea

Select K-value to be used for IOL power Calculation:

- Calculated average K-value from Refractive data 32.13
- Calculated average K-value from Rigid Trial CL 33.51
- Measured average K from preOp data screen 43.75
- Surgeon-Entered K-Value 33.50

Surgeon-Entered K

33.50 D

Because the K-readings prior to kerato-refractive surgery are not always available, & Contact Lens Trial is not always possible, the surgeon may choose to enter the central corneal power obtained from topography (Eff. RP).

NOTE:

Lowest K value tends towards myopic error.
Highest K value tends towards hyperopic error.

OK Cancel Help

Allows the surgeon to put in an average K-value directly. This value may be an average of the previous values or one from another source such as a corneal topography instrument.

6.3 IOL Calculations

Patient Name: **Nanophthalmos, Young** ID: **11AA** DOB: **13-Oct-1986** Gender: **M**

Patient Locator | PreOp Exam | **IOL Calc** | Toric IOL Planner | Surgery Data | PostOp Exam | SIRC Calculation | Back Calculation

OD Surgeon: **Surgeon1, Example** PreOp Exam date: 10/13/1996
 Ref: +12.50 -1.25 X 180 VTX: 12.00 Target Rx: -0.50
 BCVA: 20/80 UCVA: CF Target Add: 0.00
 K1: 54.62 @ 90 HWTW: 9.50 mm Avg. K: 53.44 n: 1.3375
 K2: 52.25 @ 0 ACD: 2.90 mm Adj. K: 53.44
 AL(US): 15.09 mm Lens Thick: 4.60 mm Adj. AL: 15.09

Formula: **Holladay II**

Lens #1: B&L P506UV Procedure: Std Phaco MFG ACD(US): 5.26	Lens #2: Alcon SN60T3-4-5 Procedure: Std Phaco MFG ACD(US): 5.20
IOL Ref.	IOL Ref.
26.50a+30.00p 0.20	25.50a+30.00p -0.05
27.00a+30.00p -0.15	26.00a+30.00p -0.40
27.50a+30.00p -0.50	26.13a+30.00p -0.50
27.50a+30.00p -0.50	26.50a+30.00p -0.76
28.00a+30.00p -0.86	27.00a+30.00p -1.13

Lens #3: AMO SI-30NB
Procedure: Std Phaco
MFG ACD(US): 4.62

Lens #3: AMO SI-30NB Procedure: Std Phaco MFG ACD(US): 4.62	Lens #4: Alcon MA30BA Procedure: Std Phaco MFG ACD(US): 5.49
IOL Ref.	IOL Ref.
25.50a+28.00p -0.11	26.50a+30.00p 0.14
26.00a+28.00p -0.47	27.00a+30.00p -0.21
26.00a+28.03p -0.50	27.41a+30.00p -0.50
26.50a+28.00p -0.85	27.50a+30.00p -0.56
27.00a+28.00p -1.22	28.00a+30.00p -0.92

Notes:
 The probability of the HWTW occurring in the population is < 0.1%.
 The probability of the AXIAL LENGTH occurring in the population is < 0.1%.
 The probability of the K-READING occurring in the population is < 0.1%.
 The probability of the ACD occurring in the population is < 2.3%.
 The probability of the REFRACTION occurring in the population is < 0.1%.

OS Surgeon: **Surgeon1, Example** PreOp Exam date: 10/13/1996
 Ref: +12.50 -1.00 X 180 VTX: 12.00 Target Rx: -0.50
 BCVA: 20/100 UCVA: CF Target Add: 0.00
 K1: 54.50 @ 90 HWTW: 9.50 mm Avg. K: 53.75 n: 1.3375
 K2: 53.00 @ 0 ACD: 2.95 mm Adj. K: 53.75
 AL(US): 15.31 mm Lens Thick: 4.72 mm Adj. AL: 15.31

Formula: **Holladay II**

Lens #1: B&L P506UV Procedure: Std Phaco MFG ACD(US): 5.26	Lens #2: Alcon SN60T3-4-5 Procedure: Std Phaco MFG ACD(US): 5.20
IOL Ref.	IOL Ref.
24.50a+30.00p 0.15	23.50a+30.00p -0.03
25.00a+30.00p -0.19	24.00a+30.00p -0.38
25.44a+30.00p -0.50	24.16a+30.00p -0.50
25.50a+30.00p -0.54	24.50a+30.00p -0.74
26.00a+30.00p -0.90	25.00a+30.00p -1.10

Lens #3: AMO SI-30NB
Procedure: Std Phaco
MFG ACD(US): 4.62

Lens #3: AMO SI-30NB Procedure: Std Phaco MFG ACD(US): 4.62	Lens #4: Alcon MA30BA Procedure: Std Phaco MFG ACD(US): 5.49
IOL Ref.	IOL Ref.
25.00a+26.00p 0.18	24.50a+30.00p 0.16
25.50a+26.00p -0.18	25.00a+30.00p -0.19
26.00a+25.93p -0.50	25.44a+30.00p -0.50
26.00a+26.00p -0.55	25.50a+30.00p -0.54
26.50a+26.00p -0.92	26.00a+30.00p -0.89

Notes:
 The probability of the HWTW occurring in the population is < 0.1%.
 The probability of the AXIAL LENGTH occurring in the population is < 0.1%.
 The probability of the K-READING occurring in the population is < 0.1%.
 The probability of the ACD occurring in the population is < 3.1%.
 The probability of the REFRACTION occurring in the population is < 0.1%.

Main Menu | Help | Quit | Lens Setup... | Print | <- PreOp. | Toric Calc -> | Surgery ->>

Ready 10/26/20

The patient's demographic and preoperative data are shown at the top of the screen with the surgeon's primary default formula displayed in the pull down box. If a different formula is desired, click on the **Formula** pull down menu and select the desired formula.

The calculated values for the selected formula change instantly. In the example IOL Calculation Screen above, Lens #1 for the right eye requires a **Primary Piggy Back IOL Calculation** because the IOL power necessary for this patient's desired postoperative refractive target exceeds the maximum power for that lens style (e.g. the P506UV maximum power is 30 D and the power needed is greater). The program is showing that you need a 27.5 D lens as the **FRONT (anterior)** Piggy Back IOL and a 30.0 D lens as the **BACK (posterior)** IOL. If the in-the-bag box is checked then the calculation is for both of the lenses in-the-bag (only recommended for PMMA lenses or rigid lenses) or if in-the-sulcus is checked, it is for the FRONT lens in-the-sulcus and the BACK lens in-the-bag.

The bold value for each IOL in the box is the IOL power nearest to the Target Refraction on the PreOp Data Screen.

HICSOAP Professional and Deluxe User's Manual

Values for **Bifocal IOL Add Power** are shown under notes at the bottom of the screen and printout.

NOTE: Previous versions of the program also displayed ideal toricity in the comments section (if selected). This version, however, has a new advanced Toric Calculation tab that should be used instead.

The EXACT bifocal IOL Add Power is calculated from the desired “add” entered on the preoperative screen at the vertex specified. The power shown on the IOL Calculation under screen notes is the “add” necessary on the IOL to achieve the specified add at the spectacle plane and vertex.

Click the **Print** button to print out the results. Click the **PreOp** button to return to PreOp Data Screen. Clicking on the **Surgery** button displays the Surgery data entry screen. Clicking on the **Lens Setup** button displays the Program Setup, Surgeon Lens Setup and Surgeon Piggy-Back Lens Setup Screens.

6.4 Toric IOL Planner

RIGHT

Surgeon: **Surgeon3, Example** PreOp Date: **11/24/2011**

Ref: **-3.00 -2.00 X 180** VTX: **12.00** Target Rx: **0.63**

BCVA: UCVA: Adj. AL:

Flat K: **39.00 @ 180** Steep K: **41.00 @ 90** n: **1.3375**

Avg. K: **40.00** Adj. K: **40.00**

Toric Lens from Active Lens List: **Tecnis Toric ZCTxxx (Std Phaco)** Optimize...

Surgically Induced Astigmatism: **0.50** D

Incision Location: **110** °

Formula: **Holladay II**

AMO Tecnis Toric ZCTxxx
Procedure: Std Phaco
SRG Entrd ACD(Opt): 5.60

IOL SEQ	SEQ Ref.
23.50	1.08
24.00	0.73
24.14	0.63
24.50	0.37
25.00	0.01

Lens Res. Astigm.

ZCT000	+1.65 D x 84°
ZCT150	+0.58 D x 84°
ZCT225	+0.05 D x 84°
ZCT300	+0.48 D x 174°

Incision @110° IOL @84°

IOL Placement Axis: **84°**

IOL Ideal Toricity: **2.25D @IOL**

Residual Ref.: **+0.706 +0.050D**

Post Incision K: **40.824 @84° 39.176 @174°**

LEFT

Surgeon: **Surgeon3, Example** PreOp Date: **11/24/2011**

Ref: **+0.00 +0.00 X 0** VTX: **12.00** Target Rx: **-0.50**

BCVA: **20/20** UCVA: **20/50** Adj. AL:

Flat K: **43.50 @ 45** Steep K: **45.00 @ 135** n: **1.3375**

Avg. K: **44.25** Adj. K: **44.25**

Toric Lens from Active Lens List: **Tecnis Toric ZCTxxx (Std Phaco)** Optimize...

Surgically Induced Astigmatism: **0.50** D

Incision Location: **328** °

Formula: **Holladay II**

AMO Tecnis Toric ZCTxxx
Procedure: Std Phaco
SRG Entrd ACD(Opt): 5.60

IOL SEQ	SEQ Ref.
29.00	0.07
29.50	-0.29
29.79	-0.50
30.00	-0.65
30.50	-1.01

Lens Res. Astigm.

ZCT000	+1.07 D x 129°
ZCT150	+0.02 D x 39°
ZCT225	+0.56 D x 39°
ZCT300	+1.11 D x 39°

Incision @328° IOL @129°

IOL Placement Axis: **129°**

IOL Ideal Toricity: **1.48D @IOL**

Residual Ref.: **-0.657 +0.016D x39°**

Post Incision K: **44.787 @129° 43.713 @39°**


Main Menu Help Quit Lens Setup... Print <- IOL Calc. Surgery ->

Ready 12/17/20

The *Toric IOL Planner* allows the User to determine the Ideal Toricity and Axis of Placement for a Toric IOL using the K's and Surgically Induced Astigmatism from the cataract incision. Confirm that the K-readings and the axis for the flat and steep K are correct before starting. The magnitude and axis of the manifest refraction are irrelevant, since the crystalline lens will be removed. Values for a patient may not be saved if Surgical or Post Op data have already been entered.

HICSOAP Professional and Deluxe User's Manual

Enter the magnitude of the Surgically Induced Astigmatism (SIA). For current small incision cataract surgery this value is usually between zero and 0.50 D. For a given surgeon, if Pre and Post Op K's have been entered, the mean value (centroid) may be found on the Aggregate SIRC Changes from the Main Menu as shown below. The SIRC from K's button and Minus Cylinder must be checked to obtain the correct value. In this example the magnitude of the SIA is -0.20 D. The USER would enter 0.20 D, ignoring the sign. The meridional location of the cataract incision must be entered. As soon as the magnitude and axis of the cataract incision have been entered, a BLACK ARCUATE LINE is shown at the Limbus on the drawing.




A Drop down list of available Toric IOLs is shown. The available Toric IOLs must be in the Active IOL List and on the IOL Calc Screen to be displayed in the drop down list. Select the Toric IOL to be used and the IOL Formula to be used (Holladay 2 recommended). Once this has been done, the resulting Steep Axis of the Cornea and correct axis for the IOL are shown. Below the picture is the 1) Axis of Placement, 2) Ideal Toricity of the IOL, and 3) the remaining Residual Refraction with the nearest available Toric IOL. All available Toricities are shown with the nearest to the Ideal Toricity highlighted in Bold Print and the Residual Refraction for each available toricity.

HICSOAP Professional and Deluxe User's Manual

5.4.1 Toric Lens Incision Optimizer

Clicking the 'Optimize' button in the preop toric screen will display the toric incision optimizer dialog, which graphs the residual astigmatism as a function of the incision angle:



The SIA values at 90 and 180 are pre-populated with the values entered in the Surgeon Toric setup screen (see section 4.3 above)

The radio buttons on the right specify the incision locations that minimize the residual astigmatism. Select one of these points and press OK to copy the incision angle into the Preop toric screen incision angle.

HICSOAP Professional and Deluxe User's Manual

6.5 Surgery Data Entry

The screenshot displays the 'Holladay IOL Consultant - a' software window. The patient information bar shows: Name: Nanophthalmos, Young; ID: 11AA; DOB: 13-Oct-1986; Gender: M. The 'Surgery Data' tab is active, with other tabs including Patient Locator, PreOp Exam, IOL Calc, Toric IOL Planner, PostOp Exam, SIRC Calculation, and Back Calculation. The interface is split into two columns for the right eye (OD) and left eye (OS). Each column contains the following fields: Surgery Date (11/13/1996), Surgeon (Surgeon1, Example), Lens Model Used & Procedure (4 Alcon MA30BA Std Phaco), and No. of Lenses Implanted this procedure. For the OD eye, '1 Lens' is selected with a Lens Power of 0 D, and '2 Lenses (Piggy-back)' is also selected with a Posterior lens of 30 D and an Anterior lens of 27 D. For the OS eye, '1 Lens' is selected with a Lens Power of 0 D, and '2 Lenses (Piggy-back)' is also selected with a Posterior lens of 30 D and an Anterior lens of 25 D. At the bottom of each column are checkboxes for 'Surgery Complication' and 'Corneal Astigmatism Procedure (LRI, AK, etc)'. The bottom of the window features buttons for 'Main Menu', 'Help', 'Quit', 'Setup...', 'Print', '<- IOL Calc.', and 'PostOp. ->'. The status bar at the bottom left shows 'Ready' and the bottom right shows the date '10/26/20'.

All the items in the screen must be entered. **Surgery Complication** should be selected if this patient is to be ignored in the personalized constant calculations. Corneal Astigmatism Procedures are also ignored for personalization and SIRC calculations.

Clicking the Print button prints the Surgery Data Screen. Clicking IOL Calc button returns to the IOL Calculation Screen. Clicking PostOp Data displays the PostOp Data Screen. Clicking on the Setup button displays the Program Setup, Surgeon Lens Setup and Surgeon Piggy-Back Lens Setup Screens.

HICSOAP Professional and Deluxe User's Manual

6.6 PostOp Data

The screenshot displays the 'Holladay IOL Consultant - a' software window. At the top, patient demographics are shown: Name: Nanophthalmos, Young; ID: 11AA; DOB: 13-Oct-1986; Gender: M. Below this, a navigation bar includes tabs for Patient Locator, PreOp Exam, IOL Calc, Toric IOL Planner, Surgery Data, PostOp Exam (selected), SIRC Calculation, and Back Calculation. The main area is split into two columns for OD (Ocular Dominant) and OS (Ocular Subordinate). Each column contains surgery information (Surgery Date: 11/13/1996, Surgeon: Surgeon1, Example, Lens: MA30BA, Procedure: Std Phaco, Implant Power: 57.0 (27.0a+30.0p) for OD and 55.0 (25.0a+30.0p) for OS) and a 'PostOp #1' section with fields for Date, Sphere, Cylinder, Vertex, BCVA, UCVA, K1, K2, and Complications. Navigation arrows are present above the PostOp #1 section in each column. At the bottom, there are buttons for Main Menu, Help, Quit, Setup..., Print, <- Surgery, SIRC ->, and Back Calc ->.

Patient demographics and surgery information is displayed at the top of the screen.

Check PostOp Complication if the patient's eye is to be ignored for personalized constant calculations.

Up to 3 postoperative exams for OD and OS are shown on the postoperative exams screen. If there are more than 3 postoperative exams, clicking the bold arrows allows the user to move forward and backward chronologically through the exams

Clicking the Setup button displays the Setup Screens. Clicking Surgery button returns to the Surgery Data Entry Screen. Clicking Print button prints a log of the postoperative data. Clicking BackCalcs button displays the Back Calculations screen. Clicking SIRC Tab at the top displays the SIRC (Surgically Induced Refractive Change) from keratometry and refraction.

6.7 Toric PostOp Back Calculator

The screenshot displays the 'Toric PostOp BackCalc' window for a patient named '_Toric, Calc'. The interface is split into two main columns for the 'RIGHT' and 'LEFT' eyes. Each column contains a summary of patient and surgical data, followed by calculated values for IOL Toricity & Axis, Ideal Placement Axis, and Predicted Residual Refraction. Below the text are four diagrams, each showing a circular IOL with a red 'Current' axis and a green 'Ideal' axis. The diagrams for the RIGHT eye show a 30-degree clockwise rotation from 120° to 90°. The diagrams for the LEFT eye show no rotation necessary, with the current axis at 96° and the ideal axis at 34°.

The **Toric PostOp Back Calculator** allows the surgeon to determine the exact amount that the Toric IOL needs to be rotated to achieve the smallest residual astigmatism. There are two methods available.

Method #1 uses 1) **PostOp K's** and 2) **Postop Refraction** to determine the a) the **toricity of the IOL** and its **current placement axis**, b) **ideal placement axis**, c) the **amount of clockwise rotation necessary to achieve the ideal placement axis**, and d) **predicted residual refraction at the ideal placement axis**. The blue line in the diagram illustrates the ideal axis of placement for the IOL.

Method #2 uses 1) **observed IOL meridian** at the slit lamp and 2) **PostOp Refraction** to determine the a) the **toricity of the IOL** and its **current axis**, b) **ideal placement axis**, c) the **amount of clockwise rotation necessary to achieve the ideal axis**, and d) **predicted residual refraction at the ideal placement**. The blue line in the diagram illustrates the ideal axis of placement for the IOL.

When the **required PostOp data** above has **NOT** been entered, an "Enter Missing Data" Button is present for Method #1 and/or #2 in place of the illustration. When this button is clicked, the PostOp EXAM screen appears and the missing data can be entered.

HICSOAP Professional and Deluxe User's Manual

It is prudent to measure both the **PostOp K's** and **observed IOL meridian** along with the **PostOp Refraction** as a double check of the amount to rotate the IOL to the **ideal axis of placement**. The two methods should agree to within a few degrees. When they do not match, something is wrong and measurements should be repeated until they agree.

As with the **Toric IOL Planner (forward calculation)**, the equations for the **back calculation are exact** and utilize the **effective lens position (ELP)** and **IOL Power and Toricity**. Approximation calculators that use a constant ratio to convert the cylinder at the IOL plane to the cylinder at the corneal plane will not get the same answer, especially for low and high power IOLs. If a CALCULATOR does not require the SEQ power of the IOL as input data, the calculator uses the approximation method.

6.8 Individual Surgically Induced Refractive Change (SIRC)

The screenshot displays the 'SIRC Calculation' window in the Holladay IOL Consultant software. The patient information at the top includes Name: Nanophthalmos, Young, ID: 11AA, DOB: 13-Oct-1986, and Gender: M. The interface is split into two columns for the right eye (OD) and left eye (OS).

OD (Right Eye) Data:

- Surgery Date: 11/13/1996
- Surgeon: Surgeon1, Example
- Lens: MA30BA Procedure: Std Phaco
- Implant Power: 57.0 (27.0a+30.0p)
- Post Op Date: 12/13/1996
- PreOp Rx: SPH 12.50, CYL -1.25, AXIS X 180
- PostOp Rx: SPH -0.43, CYL 0.00, AXIS X 90
- SIRC from Ref. (plus (+) Cyl form): 11.68, 1.25 @ 90
- SIRC from Ref. (minus (-) Cyl form): 12.93, -1.25 @ 0
- SIRC from Ref. ("cross" Cyl form): 11.68 @ 0 & 12.93 @ 90
- PreOp K's: K1: 54.62 @ 90, K2: 52.25 @ 0
- PostOp K's: K1: 54.25 @ 90, K2: 52.25 @ 0
- SIRC from K's (plus (+) Cyl form): -0.37, 0.37 @ 0
- SIRC from K's (minus (-) Cyl form): 0.00, -0.37 @ 90
- SIRC from K's ("cross" Cyl form): -0.37 @ 90 & 0.00 @ 0

OS (Left Eye) Data:

- Surgery Date: 11/13/1996
- Surgeon: Surgeon1, Example
- Lens: MA30BA Procedure: Std Phaco
- Implant Power: 55.0 (25.0a+30.0p)
- Post Op Date: 12/13/1996
- PreOp Rx: SPH 12.50, CYL -1.00, AXIS X 180
- PostOp Rx: SPH -0.40, CYL 0.00, AXIS X 0
- SIRC from Ref. (plus (+) Cyl form): 11.90, 1.00 @ 90
- SIRC from Ref. (minus (-) Cyl form): 12.90, -1.00 @ 0
- SIRC from Ref. ("cross" Cyl form): 11.90 @ 0 & 12.90 @ 90
- PreOp K's: K1: 54.50 @ 90, K2: 53.00 @ 0
- PostOp K's: K1: 54.25 @ 0, K2: 52.25 @ 90
- SIRC from K's (plus (+) Cyl form): -2.25, 3.50 @ 0
- SIRC from K's (minus (-) Cyl form): 1.25, -3.50 @ 90
- SIRC from K's ("cross" Cyl form): -2.25 @ 90 & 1.25 @ 0

The bottom of the window contains buttons for 'Main Menu', 'Help', 'Quit', 'Setup...', 'Print', '<- Surgery', and 'Back Calc ->'. The status bar at the bottom left shows 'Ready' and the bottom right shows the date '10/26/20'.

If postoperative K-values and refractions have been entered the program calculates the Surgically Induced Refractive Change (SIRC). The SIRC values from the refraction are shown at the top of the screen and from the K-values at the bottom of the screen.

The SIRC is displayed in "+" cylinder form, "-" cylinder form, and "cross" cylinder form. The "+" cylinder form is most helpful for compressing or steepening procedures, the "-" cylinder form is most helpful for

HICSOAP Professional and Deluxe User's Manual

relaxing or flattening procedures, and the “cross” cylinder form is most helpful for seeing the actual change in each meridian.

6.9 Individual Back Calculations

The screenshot displays the 'Holladay IOL Consultant - a' software window. The patient information at the top includes Name: Nanophthalmos, Young, ID: 11AA, DOB: 13-Oct-1986, and Gender: M. The 'Back Calculation' tab is selected, showing data for both the OD (Right Eye) and OS (Left Eye).

Parameter	OD (Right Eye)	OS (Left Eye)
Surgery Date	11/13/1996	11/13/1996
Surgeon	Surgeon1, Example	Surgeon1, Example
Lens	MA30BA	MA30BA
Procedure	Std Phaco	Std Phaco
Implant Power	57.0 (27.0a+30.0p)	55.0 (25.0a+30.0p)
Post Op Date	12/13/1996	12/13/1996
Formula	Holladay II	Holladay II
Actual Values:		
Implanted Lens Pwr	57.00	55.00
PostOp Ref. Sph.Eqv	-0.428	-0.404
PreOp AL	15.09	15.31
PreOp Avg K	53.44	53.75
Lens Constant(LC)	5.491 MFG	5.491 MFG
BackCalculated Values (Calculated from 4 measured values)		
From:		
IOL Power	58.658 D. (K, AL, LC, Ref)	56.551 D. (K, AL, LC, Ref)
Ref	0.708 D. (K, AL, LC, IOL)	0.647 D. (K, AL, LC, IOL)
Axial Length	15.300 mm (K, LC, IOL, Ref)	15.511 mm (K, LC, IOL, Ref)
Avg. K-value	54.590 D. (AL, LC, IOL, Ref)	54.817 D. (AL, LC, IOL, Ref)
LC	5.185 (K, AL, IOL, Ref)	5.202 (K, AL, IOL, Ref)
Prediction Error	1.136 D. (Pred Ref - PO Ref)	1.051 D. (Pred Ref - PO Ref)

At the bottom of the window, there are buttons for 'Main Menu', 'Help', 'Quit', 'Setup...', 'Print', '<- Sirc', and 'Pat. Locator ->'. The status bar shows 'Ready' and the date '10/26/20'.


The Back Calculation Screen very helpful for working up a refractive surprise following IOL implantation. The Actual Values are shown on the top of the screen and the Back Calculated Values are shown at the bottom of the screen. Each Back Calculated value shown is calculated from the four other Actual Values at the top of the screen. Since the Avg. K-value, Axial Length, and Refraction can be checked postoperatively, the ACD and IOL Power become the only unknown parameters. If the ACD appears to be correct by comparing the other formulas or with direct measurement, then the power of the mislabeled lens in the eye is the back calculated IOL Power. The back calculated values can be obtained using any of the formulas.

HICSOAP Professional and Deluxe User's Manual

7. Surgical Outcomes Analysis Program (SOAP)

SOAP provides the surgeon outcome analysis tools. The following sections describe in detail the different analysis tools available.

All the SOAP screens share a common toolbar that control several parameters used in the analysis:



- **Surgeon:** The pull down menu on the left of the tool bar allows the surgeon to be changed at any time.
- **Lens:** **All lenses** or a **specific lens model** may be chosen from the pull down menu. If "**All lenses**" is chosen then three columns in the statistics at the bottom of the report are blank: Personalized ACD, Eqv. Form. Constant and Current Form. Constant.
- **Formula:** The IOL formula used may be changed at any time from the pull down menu.
- **AL:** Either Ultrasonic (US) or Optical (IOL Master or Lenstar) may be selected.
- **Step Size:** The step size for the graphs may be chosen to be 0.125 D, 0.250 D or 0.500 D. This allows the user to see the results as details as 1/8 diopter steps.

In addition, all the SOAP screens share a common list of buttons near the bottom of the screen:



HICSOAP Professional and Deluxe User's Manual

- **Selection Criteria:** Displays the selection criteria pop up screen.

Soap Selection Criteria

Use enough cases to satisfy standard statistical criteria, or ...

Use the most recent cases, up to

Use cases from: to:

That satisfy the following conditions:

Include Only PostOp Results from to days Post-Op.

Include Only PostOp Results with BCVA of or better.

Axial Length - Lower Limit: mm - Upper Limit: mm.

IOL Power - Lower Limit: D - Upper Limit: D.

Avg. K-value- Lower Limit: D - Upper Limit: D.

PreOp Ref. - Lower Limit: D - Upper Limit: D.

PostOp Ref. - Lower Limit: D - Upper Limit: D.

Pred. Error - Lower Limit: D - Upper Limit: D.

Previous Rk, Lasik: No Yes

Include cases within limits only

Include cases outside limits only

Exclude PreOp Pathology

Exclude Surgery Complication

Exclude PostOp Complication

Technician:

OK Cancel


This screen allows the user to choose the number of cases, periods of time and various other conditions that may be of specific interest. The prediction error can look at the cases within the interval or outside the interval chosen. This is especially useful for looking at the outliers with large prediction errors.

- **Export Selected Patients:** The patients found by the “Selection Criteria” will have their data exported as a comma separated variable (*.CSV) format file which can be imported by EXCEL, etc.
- **Print Selected Patients:** A two page summary of the preoperative, surgical and postoperative data is printed for the Surgeon, Lens(es) and the Selection Criteria chosen.

HICSOAP Professional and Deluxe User's Manual

- **Print Current Page:**– The entire page is printed, including graphs and the summary statistics. See “Example A” in the Appendix.

7.1 SOAP Prediction Error



5.1.1 Frequency Distribution (%) of Eyes

This graph demonstrates the **Spherical Equivalent Prediction Error (D)** for the surgeon, lens(es) and formula selected for each interval of prediction error as a percentage of the total number. In the example above the surgeon has 37% of his cases with zero prediction error.

5.1.2 Cumulative Distribution (%) of Eyes

This graph demonstrates the **Spherical Equivalent Prediction Error (D)** for the surgeon, lens(es) and formula selected as a cumulative percentage. The surgeon can determine his cumulative absolute error for any level of prediction error. In the example above the surgeon has 85% of his cases less than or equal to 0.50 D.

HICSOAP Professional and Deluxe User's Manual


5.1.3 Summary Statistics

The **summary statistics window** near the bottom of the screen yields the numerical results for the Surgeon, Lens(es) and Formula selected. These statistics are similar to the values given on the Formula and Lens Constant Analysis Screen in the Holladay IOL Consultant.

Excluded Cases: Patients who are outliers and have been excluded from the Formula and Lens Constant Analysis Table (Summary Statistics Window) may be listed. The patient's name, chart number, eye, and reason for exclusion are listed. The possible reasons for exclusion are listed below:

- PO period less than 14 days.
- BCVA (best corrected visual acuity) below criteria selected (usually selected as less than 20/50 – “outliers”).
- Back calculated lens constant for this case is more than 2 standard deviations from the mean.
- Preoperative pathology, surgical or postoperative complication were selected on the respective screen.

7.2 SOAP-SIRC Changes Report



HICSOAP Professional and Deluxe User's Manual

5.2.1 Surgically Induced Spherical Equivalent (SIRC) – Upper left graph

This graph displays the Frequency Distribution (%) of the SIRC as a function SIRC spheroequivalent. This is the actual spheroequivalent change from preoperative to postoperative.

5.2.2 Actual PreOp Cylinder and Axis – Lower left graph

This graph displays the actual preoperative astigmatism on a double angle plot. The centroid (mean value), standard deviation and ring values (step size) are shown on the plot. **WTR** is the abbreviation for “With-the-Rule” astigmatism and **ATR** is the abbreviation for “against-the-rule” astigmatism. The ellipse (oval) or circle represents one standard deviation (67% of the data for a normal distribution).

5.2.3 Surgically Induced Cylinder and Axis (SIRC) – Upper right graph

This graph displays the actual Surgically Induced Cylinder and Axis on the double angle plot. The centroid, ring values and standard deviation are all displayed.

5.2.4 Actual PostOp Cylinder and Axis – Lower right graph

This graph displays the actual PostOp astigmatism cylinder and axis on the double angle plot. The centroid, ring values and standard deviation are all displayed.

5.2.5 Buttons under Graphs

5.2.5.1 SIRC from K

The SIRC from K calculates the surgically induced refraction change based on the preoperative and postoperative K-readings. This is the normal choice and reflects what changes the surgeon has made in the cornea.

5.2.5.2 SIRC From Ref

The SIRC from refraction calculates the surgically induced refraction change based on the preoperative and postoperative refractions. This value would reflect any changes in the pre- and postoperative refraction. Since this value is affected by the power and cylinder of the crystalline lens and the power of the IOL, it often has limited value in cataract surgery.

5.2.5.3 Plus (+) or Minus (-) Cylinder Notation

If the Plus(+) radio button is selected then all values are displayed in plus cylinder notation. If the Minus(-) radio button is selected then all values are displayed in minus cylinder notation. Usually, steepening procedures are best displayed in Plus(+) cylinder notation and flattening procedures are best displayed in Minus (-) cylinder notation.

7.3 SOAP Potential Outcomes



5.3.1 Frequency Distribution (%) of Eyes – Upper Left Graph

This graph displays the Spherical Equivalent Prediction Error (D) with the “Optimized Surgeon Constant” (Opt SC) in solid bars and with the “manufacturers recommended constant” in hatched bars.

5.3.2 Cumulative Distribution (%) of Eyes – Upper Right Graph

This graph displays the cumulative distribution of the Spherical Equivalent Prediction Error (D). The line with boxes represents the distribution with the “Optimized Surgeon Constant” and the line with triangles with the “manufacturers recommended constant”.

5.3.3 Frequency Distribution (%) of VA – Lower Left Graph

This graph displays the visual acuity’s at various acuity levels. The solid bars are the actual acuity’s and the hatched bars are the theoretical acuity’s that would be predicted from the patients’ refractions.

5.3.4 Cumulative Distribution (%) of VA – Lower Right Graph

This graph displays the cumulative distribution of the visual acuity’s at various acuity levels. The line with boxes is the actual acuity’s and the line with triangles is the theoretical acuity’s that would be predicted from the patients’ refractions.

HICSOAP Professional and Deluxe User's Manual

5.3.5 Buttons under graphs

5.3.5.1 Distance BSCVA

When this radio button is chosen all acuity values are displayed using the actual Best Spectacle Corrected Visual Acuity (BSCVA).

5.3.5.2 Distance UCVA

When this radio button is chosen all acuity values are displayed using the actual or theoretical Uncorrected Visual Acuity (UCVA) values.

5.3.5.3 Optimize Surgeon Constant

When this box is checked the theoretical values for visual acuity for patients' refractive error are displayed with the "optimized surgeon constant". When this is not checked the theoretical values using the "manufacturer's recommended constant" are displayed.

5.3.5.4 Zero Induced Astigmatism


When this box is checked, the theoretical visual acuity's reflect the results if the surgeon had not induced any astigmatism on any patient.

5.3.5.5 Zero PostOp Astigmatism

When this box is checked, the theoretical visual acuity's reflect the results if the surgeon had reduced every patient's postoperative astigmatism to zero.

8. Utilities

8.1 Database



Patient Database File: Indicates the database file currently open.

Backup: Backs up the database by first displaying the SAVE AS dialog box to specify the directory and drive to be used for the backup operation.

Restore: Restores from a list of automatic backups made by this program or restores from another backup. If more than one backup file is created on the same day, the older file will be replaced by the new file. The most recent backup file should be used.

Move Database: Moves the database to a different location (drive and folder). The program will continue to open the database at the new location when the program is started.


Repair and Compact Database: Repairs the database. This is useful when the program is not working right and should be done at least once a year.

HICSOAP Professional and Deluxe User's Manual

Open Other database: This feature allows the user to Open a different HICSOAP database. When this button is clicked the File Open Window below appears with a “Select a file from the following list” for existing Databases. If the Database has been opened before, it will appear in the list to select. If the Database has never been opened before, but is a *.hdf file, the “select a new file to open” radio button should be clicked and an Open File dialog box will appear from which you may choose the file.


Open Database Folder: Clicking this simply opens the folder containing the current database using Windows Explorer.

8.2 Documentation



Several PDF files documenting how to use HICSOAP are included with the program. Simply click the buttons to open the associated PDF file.

8.3 Remove Registration



We allow 3 active installations for each purchased license. We keep track of the number of active installations via the registration code. So if you no longer want to use the software on a computer then you need to remove the registration.

Remove Registration: Click this to remove registration. Doing so generates a registration removal number that you should email to support@hicsoap.com so that we can record your registration removal in our database.

System Information: The System Information button starts the Microsoft System Information program.

HICSOAP Professional and Deluxe User's Manual

A. Sample Printouts

A.1 IOL Calc Report

HicSoapPro - Iol Calc Report		Report Date: 10/26/2010	
Patient: Patient, Normal		ID: 11AA	Date of Birth: 13-Oct-1924 Sex: Female
Pre-Op. Data			
Surgeon: Surgeon1, Example 10/13/1996		Surgeon: Surgeon1, Example 10/13/1996	
Refraction: -0.50 +0.50 X 180		Refraction: -1.00 +1.00 X 180	
Vertex: 12.00	AL(US): 23.45	Vertex: 12.00	AL(US): 23.45
BCVA: 20/70	Adj.AL: 23.45	BCVA: 20/100	Adj.AL: 23.45
UCVA: 20/80	HorW-t-W: 11.70	UCVA: 20/200	HorW-t-W: 11.70
K1: 43.50 @90	Phakic Lens Th.: 4.70	K1: 43.75 @90	Phakic Lens Th.: 4.70
K2: 44.00 @0	TargetRef: -0.25	K2: 43.75 @0	TargetRef: -0.25
AverageK: 43.75	n: 1.3375	AverageK: 43.75	n: 1.3375
AdjustedK: 43.75	TgtAdd: -3.00	AdjustedK: 33.50	TgtAdd: 0.00
Additional Data			
Eye Status: Phakic New PC Lens: in bag		Eye Status: Phakic New PC Lens: in bag	
PreOp Pathology: No Prev. Rx...: No Keratoplasty: No Scleral Buckle: No Silicone In Vitreous Cavity: No		PreOp Pathology: No Prev. Rx...: No Keratoplasty: No Scleral Buckle: No Silicone In Vitreous Cavity: No	
Formula: Holladay II		Formula: Holladay II	
Lens #1: B&L P506UV Procedure: Std Phaco MFG ACD(US): 5.26		Lens #2: Alcon SN60T3-4-5 Procedure: Std Phaco MFG ACD(US): 5.20	
<u>IOL</u> <u>Ref.</u>	<u>IOL</u> <u>Ref.</u>	<u>IOL</u> <u>Pred. Ref.</u>	<u>IOL</u> <u>Pred. Ref.</u>
20.50 0.20	20.50 0.12	9.50a+26.50p 0.12	5.50a+30.00p 0.23
21.00 -0.14	21.00 -0.22	10.00a+26.50p -0.25	6.00a+30.00p -0.14
21.16 -0.25	21.04 -0.25	10.00a+26.51p -0.25	6.15a+30.00p -0.25
21.50 -0.48	21.50 -0.57	10.50a+26.50p -0.62	6.50a+30.00p -0.52
22.00 -0.83	22.00 -0.92	11.00a+26.50p -1.00	7.00a+30.00p -0.89
Lens #3: AMO SI-30NB Procedure: Std Phaco MFG ACD(US): 4.62		Lens #4: Alcon MA30BA Procedure: Std Phaco MFG ACD(US): 5.49	
<u>IOL</u> <u>Ref.</u>	<u>IOL</u> <u>Ref.</u>	<u>IOL</u> <u>Pred. Ref.</u>	<u>IOL</u> <u>Pred. Ref.</u>
19.00 0.41	21.00 0.18	25.00a+8.00p 0.48	9.00a+27.00p 0.46
19.50 0.05	21.50 -0.15	25.50a+8.00p 0.10	9.50a+27.00p 0.10
19.92 -0.25	21.64 -0.25	26.00a+7.95p -0.25	10.00a+26.98p -0.25
20.00 -0.31	22.00 -0.49	26.00a+8.00p -0.29	10.00a+27.00p -0.27
20.50 -0.67	22.50 -0.84	26.50a+8.00p -0.68	10.50a+27.00p -0.63
IOL Consultant Notes			
The probability of this IOL is < 0.1% Anterior Lens was set to minimum available			
Surgeon/Technician Notes			
HicSoapPro Ver 2010.1010		Test Eye Clinic	
		Page 1 of 1	

HICSOAP Professional and Deluxe User's Manual

A.2 Toric Iol Report

HicSoapPro - Toric Iol Report		Report Date: 10/26/2010		
Patient: Patient, Normal		ID: 11AA	Date of Birth: 13-Oct-1924 Sex: Female	
Pre-Op. Data				
Surgeon: Surgeon1, Example 10/13/1996 Refraction: -0.50 +0.50 X 180 AL(US): 23.45 Vertex: 12.00 Adj.AL: 23.45 BCVA: 20/70 HorW-t-W: 11.70 UCVA: 20/80 PhakioACD: 3.10 K1: 43.50 @90 Phakio Lens Th: 4.70 K2: 44.00 @0 TargetRef: -0.25 AverageK: 43.75 n: 1.3375 TgtAdd: -3.00 AdjustedK: 43.75		OD OS	Surgeon: Surgeon1, Example 10/13/1996 Refraction: -1.00 +1.00 X 180 AL(US): 23.45 Vertex: 12.00 Adj.AL: 23.45 BCVA: 20/100 HorW-t-W: 11.70 UCVA: 20/200 PhakioACD: 3.10 K1: 43.75 @90 Phakio Lens Th: 4.70 K2: 43.75 @0 TargetRef: -0.25 AverageK: 43.75 n: 1.3375 TgtAdd: 0.00 AdjustedK: 33.50	
Additional Data				
Eye Status: Phakic New PC Lens: in bag PreOp Pathology: No Prev. Rx...: No Keratocoids: No Scleral Buckle: No Silicone In Vitreous Cavity: No		Eye Status: Phakic New PC Lens: in bag PreOp Pathology: No Prev. Rx...: No Keratocoids: No Scleral Buckle: No Silicone In Vitreous Cavity: No		
<u>Lens</u> SA60T3 (1.50D @IOL) SA60T4 (2.25D @IOL) SA60T5 (3.00D @IOL)		<u>Res. Refraction</u> 0.17 -0.79 D x 0° 0.43 -1.31 D x 0° 0.69 -1.82 D x 0°		
HicSoapPro Ver 2010.1010		Test Eye Clinic		