

INFUSCALE

Version 7

Instructions for Use

TABLE OF CONTENTS

1. INTRODUCTION	5
1.1 ENHANCEMENTS TO VERSION 6.5.....	5
1.2 GENERAL INFORMATION ON THE INFUSION STANDARD.....	6
1.3 THE GRAVIMETRIC MEASUREMENT PROCEDURE	7
2. INSTALLATION.....	8
2.1 EQUIPMENT REQUIREMENTS	8
2.2 SOFTWARE INSTALLATION.....	8
2.3 CONNECTING THE BALANCE TO THE PC.....	9
2.4 BALANCE SETTINGS	10
2.4.1 <i>Mettler Analysis-Balance AE 100</i>	11
2.4.2 <i>Mettler Precision Balance PM4000</i>	12
2.4.3 <i>Mettler-Balances with SICS Serial Interface</i>	12
2.4.4 <i>Sartorius Balances</i>	14
3. MAIN MENU AND SETUP	15
3.1 TERMINAL.....	16
3.2 PRINTER SETUP.....	17
3.3 TIMER-CALIBRATION	17
3.4 OPTIONS.....	18
4. INFUSION PUMPS ACCORDING TO EN 60601-2-24	19
4.1 VOLUMETRIC INFUSION CONTROLLERS, VOLUMETRIC PUMPS AND SYRINGE PUMPS (50.102)	19
4.2 DRIP RATE INFUSION CONTROLLERS AND DRIP RATE INFUSION PUMPS (50.103)	19
4.3 INFUSION PUMPS FOR AMBULATORY USE TYPE 1 (50.104)	20
4.4 INFUSION PUMPS FOR AMBULATORY USE TYPE 2 (50.105)	20
4.5 INFUSION PUMPS FOR AMBULATORY USE TYPE 3 (50.106)	20
4.6 INFUSION PUMPS FOR AMBULATORY USE TYPE 4 (50.107)	20
4.7 INFUSION PUMPS FOR AMBULATORY USE TYPE 5 (50.108)	21
4.8 OTHER MEDICAL PUMPS	21
4.9 OVERVIEW OF INFUSION PUMPS ACCORDING TO EN 60601-2-24.....	21
4.10 CLASSIFICATION OF PUMPS	21
5. PERFORMING THE MEASURING	23
5.1 THE MEASURING SETUP.....	23
5.2 DATA ACQUISITION.....	24
5.2.1 <i>Entering data</i>	24
5.2.2 <i>Control of Data Transfer</i>	25
5.2.3 <i>Terminal</i>	26
5.2.4 <i>Start of measurement</i>	27
5.2.5 <i>Event marker</i>	30
5.2.6 <i>Storage of Graphics on File and to Clipboard</i>	30
5.3 SIMULATION.....	30
6. EVALUATION OF THE MEASURED DATA	31
6.1 THE TRUMPET CURVE	33
6.1.1 <i>The Meaning of the Trumpet Curve according to EN 60601-2-24</i>	33
6.1.2 <i>Calculation of the Trumpet Curve according to EN 60601-2-24, 50.102</i>	33
6.1.3 <i>Calculation of the Trumpet Curve according to EN 60601-2-24, 50.104</i>	35
6.1.4 <i>Calculation of the Trumpet Curve according to EN 60601-2-24, 50.105</i>	36
6.1.5 <i>Calculation of an optional Trumpet Curve</i>	37
6.1.6 <i>Modified Trumpet Curve and Constancy Index</i>	38
6.1.7 <i>Other Commands</i>	39
6.1.8 <i>Data Output</i>	40

6.2.	BOLUS EVALUATION	44
6.2.1.	<i>Evaluation according to EN 60601-2-24, 50.106</i>	44
6.2.2.	<i>Bolus Evaluation with adjustable Time Range</i>	45
6.2.3.	<i>Other Commands</i>	46
6.2.4.	<i>Data Output</i>	48
7.	MEASURING ACCURACY	52
7.1.	EVAPORATION.....	52
7.2.	INTERROGATION ERROR	54
7.2.1.	<i>Mettler Balance AE100 (Continuous Operation)</i>	54
7.2.2.	<i>Sartorius Balance MC1 (automatic data transmission)</i>	55
7.2.3.	<i>Sartorius Balance MC1 (request mode)</i>	56
8.	SYMBOLS AND FORMULA	57

Use of INFUSCALE

It is not permitted to make copies or to sell Infuscale. It is permissible, however to make copies of the program disc to secure the data.

Liability

We endeavour to supply a perfect product but we cannot guarantee a perfect running of the software under all conditions. Despite of an exact program validation we cannot guarantee the reaching of a special purpose of use. Liability for any damage during use of this program is excluded as far as lawful.

1. Introduction

The Infuscale program serves for the exact evaluation of the delivery accuracy of infusion pumps and controllers by the gravimetric method according to the international standard IEC601-2-24 (EN60601-2-24). It operates with electronic balances and runs on IBM compatible computers on Windows 95, 98, ME, NT, 2000, XP and 7 operating systems.

INFUSCALE provides powerful features including:

- Acquisition of balance data and storage on file
- Storage of event markers
- Operation with Mettler-, Sartorius-, Ohaus-, Chyo-, Precisa-, and Kern balances
- Graphical and numerical display of data and event marks
- Evaluation of data according to the specifications of the harmonized standard EN60601-2-24 / IEC601-2-24 or to the user's own requirements
- Evaluation according to MDA Device Bulletin DB2003(02)
- Calculation of Trumpet Curves
- Calculation of the Constancy Index
- Graphical display of weight and flow curves
- Storage of graphics as BMP or JPG files and to the Clipboard
- Display of values as tables and test reports
- Screen printout, table printout, printout of test report
- Storage of test report as Microsoft Word and Excel file, and in the data file formats: TXT, WMF, HTML and CSV.
- Storage of table as Excel data file

INFUSCALE has been used since many years in testing a large number of different pumps with different attributes demonstrating that it is a reliable tool. Data files stored with prior versions of INFUSCALE can be read and evaluated with the present version.

1.1 Enhancements to version 6.5

- Calculation of the modified trumpet curve and the constancy index
- Improved export to Microsoft Excel
- Improved entry of comments
- Improved setting of table print out parameters
- Storage of the name of the technician within the data file
- Optional operation of two balances on one PC.

1.2. General Information on the Infusion Standard


In recent years work has been done by a working group of the VDE and IEC (International Electrotechnical Committee) on a new international standard for pumps and controllers. This worldwide recognised commission has completed the Norm IEC 601-2-24 „Particular requirements for safety of infusion pumps and controllers“, First edition 1998-02. In the same year the standard has been harmonised according to the European Directive on Medical Devices 93/42/EEC as EN 60601-2-24.

As the graphics described in the above mentioned standards have to be part of any user manual for infusion pumps, a measuring system is necessary that fulfils those requirements.

The main task of INFUSCALE is to perform - together with the required hardware - and document accuracy measurements of infusion pumps according to sections 50.102, 50.104 - 50.108 of the standard. In addition INFUSCALE is useful during the development and the periodic safety checks of infusion pumps.

1.3. The Gravimetric Measurement Procedure

The gravimetric method is used in determining the accuracy of delivery of infusion pumps. During the measuring procedure the quantity of liquid delivered by the pump under test is weighed in equal time intervals and transmitted to a computer, that stores it together with a timestring on a data file. In a second process the accuracy data are calculated and documented.



Picture 1 – Measuring setup

The measuring setup (picture 1) is made according to IEC standard requirements (see pictures 104a and 104b of the standard). It consists of the pump with disposable accessories, an electronic balance and the connected digital computer. The measuring is controlled and documented by INFUSCALE. Depending on the type of pump (See chapter 50.102 - 50.108 of the standard) the pump and the measuring devices are set and started. The pump delivers liquid (water for medical use according to ISO class III) into the container located on the balance plate. The weight of the container increases in proportion to the delivered volume and is in regular intervals transmitted to the computer and stored on file. The menu option for performing this procedure is started by clicking on “Data Acquisition” / “Measurement”.

Note: Nutrition pumps must be tested with the specified nutrition solution instead of water. Before starting measuring, the density of the nutrition solution must be determined and entered into the respective field.

After the measuring is finished, the data can be analysed by calling up the menu item “Evaluation”.

With pumps of group 50.102 the evaluation according to EN 60601-2-24 is very easy. You just have to load the data file and click on key **102**. Prerequisite is a time interval of 30 seconds between two data points and a duration of measurement of > 2 hours which is equivalent to at least 241 measuring points.

With other pumps - i.e. pumps with bolus delivery or pumps with non continuous delivery characteristics - an automatic analysis is not recommended. In those cases the user's analysis is supported by different options in the menu “Evaluation”.

2. Installation

2.1. Equipment Requirements

- Electronic balance with RS232-interface (uni- or bidirectional).
- Picture 2 shows a list of balances supported by INFUSCALE.
- PC with Windows 95, 98, NT, 2000, XP or 7 operating system
- Colour Graphic Display of at least 800/600 pixel resolution.
- RS232 computer interface installed as COM1 .. COM4
- Printer installed on LPT1 or LPT2. All printers which are accepted by Windows can be used.
- Connection line between computer and balances
- Accessories (container, cannula, oil)

2.2. Software installation

Insert the program-CD into the CD-drive.

Start the software installation with "Drive:\Setup-E.EXE" in the "run" mode and press the OK key. The installation routine transfers all program- and data files from CD to the computer and creates an icon symbol on the desktop.

2.3. Connecting the Balance to the PC

The balance can be connected to any of the serial ports COM1 to COM4 of your computer. The specific port is selected in the INFUSCALE setup procedure (see chapter 3). As COM1 is normally used for the mouse, we recommend the selection of COM2 for connecting the balance.

We recommend to start without using handshake, and the balance set to AUTO (continuous) mode.

The following lines must at least be connected between the balance and the specified COM-port.

1. GND
2. Data line: Balance \Rightarrow PC
3. Data line: PC \Rightarrow Balance

It could happen, that the connection of the data lines is initially incorrect and that no data appear on the computer screen. In this case most often it is sufficient to cross the lines 2 and 3. Please check the settings of baudrate, parity, data bits and stopbits on balance and computer (see chapter 3). Set the handshake option of the electronic balance to OFF and the transmission mode to continuous (CONT or AUTO).

2.4. Balance Settings

The next table gives an overview of balances supported by INFUSCALE.

Manufacturer	Series	Auto	Abfr	Default parameters
Mettler	AE, AM, PM	X	X	2400 Bd, 7E1
Mettler SICS	AB-S, AG, AT, AX, MT, MX, PB-S, PG-S, PR, SB, SG, SR, UMT, UMX	X	X	2400 Bd, 7E1
Sartorius	AC, BA, BL, BP, LA, LC, LP, MC, RC, SC, CP	X	X	9600 Bd, 7E1
Ohaus GA	GA110, GA160		X	4800 Bd, 7E1
Ohaus GT	GT-Series	X	X	9600 Bd, 7E1
Chyo	MJ-Series	X		4800 Bd, 8N1
Precisa	Series 300S	X		9600 Bd, 7E1

Picture 2 - Balances and setting parameters

In column 3 and 4 the operation modes of the balances are marked.

„auto“ stands for continuous data transmission from the electronic balance to the computer.

„abfr“ stands for data transmission on request by the computer.

On some balances both transmission modes can be selected. INFUSCALE must be configured according to the mode selected at the balance (see chapter 3)

Column 5 of picture 2 contains the default parameters preselected by the balance manufacturer e.g. baudrate 9600 bd, word length 7 bit, 1 stop bit and even parity for Sartorius balances. In the Setup screen these default parameters can be automatically taken over from an internal table by pressing

Preset values

However any other value can be selected if necessary for proper communication. It is important that the settings of computer and balance correspond.

The next chapters give additional information for specific balance setups.

2.4.1. Mettler Analysis-Balance AE 100

integration time: step 2 = 3 sec (normal setting)
stability detector: step 2 (normal setting)

Data interface option 011:

In accordance with the installation instructions for option 011 for the Mettler AE-Balance the following settings have to be made:

Data transmission mode: send continuously
Baudrate: 2400 Baud
Data bits: 7
Stop bit: 1
Parity: parity enable, even parity
No protocol

Mark "Mettler" and "auto" in the Setup screen, and check the correspondence of the fixed parameters between computer and balance interface (see chapter 3).

With the correct setting of the COM-parameters the balance is continuously sending data to the PC. See next sequence for example:

```
S 0.00 g <--- weight 0.0 gramms
S 0.00 g
S 0.00 g
SD 0.5 g <- 7.2 gramms put on balance, result not stable
SD 2.6 g
SD 4.8 g
SD 6.1 g
SD 6.7 g
SD 7.0 g
SD 7.1 g
SD 7.1 g
SD 7.2 g
SD 7.20 g
SD 7.20 g
SD 7.20 g
S 7.20 g <---- result stable
S 7.20 g
S 7.20 g
S 7.20 g
S 7.20 g
S 7.20 g
```

Picture 3 - Data from a Mettler AE balance

During measuring INFUSCALE only accepts stable results.

2.4.2. Mettler Precision Balance PM4000

Please make the following settings in the sector "I-FACE" of the balance configuration register.

data transmission mode: "Cont" all values are transmitted automatically
 Baudrate: 2400 (default)
 Parity: E = even parity (default)
 Data bits: 7
 Stop bits: 1
 Interval (pause): 0 sec
 Calibration mark: AU = on (or off)

Mark "Mettler" and "auto" in the Setup screen, and check the correspondence of the fixed parameters between computer and balance interface (see chapter 3).

2.4.3. Mettler-Balances with SICS Serial Interface

This section covers all balances of model types:

AB-S, AG, AT, AX, MT, MX, PB-S, PG-S, PR, SB, SG, SR, UMT, UMX

The following table gives an overview of port connections to the PC.

Balance type	Port	PC-Connection
AB-S	RS-232C	Cable
AG	LocalCAN	Interface-Box with cable
AT	RS-232C	Cable
AX	RS-232C	Cable
MT	RS-232C	Cable
MX	RS-232C	Cable
PB-S	RS-232C	Cable
PG-S	RS-232C	Cable
PR	LocalCAN	Interface-Box with cable
SB	RS-232C	Cable
SG	LocalCAN	Interface-Box with cable
SR	LocalCAN	Interface-Box with cable
UMT	RS-232C	Cable
UMX	RS-232C	Cable

Picture 4 - Mettler balances with SICS Interface connection

Balances with RS-232C port are connected according to the following scheme:

Balance-Pin 13 (white)	⇒	Pin 5 of 9 pol. PC-Plug	GND
Balance-Pin 12 (brown)	⇒	Pin 2 of 9 pol. PC-Plug	data from Balance
Balance-Pin 2 (green)	⇒	Pin 3 of 9 pol. PC-Plug	data to Balance
Connect Pin 1,7 and 8 of 9 pol. PC-plug			

Picture 5 - Connection scheme

These balances can be operated in two different modes:

a) continuous data transmission to the computer:

The settings in sector "I-FACE" of the balance configuration register of the balance are:

Data output condition: CONT
Baudrate: 2400 (default)
Parity: E = even parity (default)
Data bits: 7
Stop bits: 1
Handshake: OFF
End of line mode: CRLF

b) data transmission on request by the computer:

Data output condition: ALL
Baudrate: 2400 (default)
Data bits: 7
Stop bits: 1
Parity: E = even parity (default)
Handshake: OFF
End of line mode: CRLF

The selection in the SETUP screen must be either "Mettler S", "auto" or "Mettler S", "abfr".

2.4.4. Sartorius Balances

All Sartorius balances of types **AC, BA, BL, BP, LA, LC, LP, MC, RC, SC, CP** have a RS-232C-port. They can be operated in two different modes:

a) continuous data transmission to the computer:

data output condition:	auto without standstill
Baudrate:	9600 (default)
Data bits:	7
Stop bits:	1
Parity:	E = even parity (default)
stopbits:	1
handshake:	hardware handshake with 2 figures after CTS

b) data transmission on request by the computer:


data output condition:	without standstill
Baudrate:	9600 (default)
Data bits:	7
Stop bits:	1
Parity:	E = even parity (default)
stop bits:	1
handshake:	hardware handshake with 2 figures after CTS

In accordance with the balance configuration the setup of INFUSCALE must be either "Sartorius", "auto" or "Sartorius", "abfr".

3. Main Menu and Setup


After starting INFUSCALE the main menu appears on the screen. The options “Data Acquisition”, “Evaluation”, “Setup” and “Help” can be selected from here.

Please click on Setup now:



Picture 6 - Main menu

At the very first call of the Setup routine, INFUSCALE tries to find an available serial port. If a suitable port has been found, the Setup screen appears:



Picture 7 – Setup Screen

If INFUSCALE has been installed on Windows 98 the checkbox **Win 98** must be marked. In case of Windows NT, 2000, XP or 7, this field remains empty. At next, the connected electronic balance is set by clicking on its name in the list "Balances".

A click on key



selects the preset values "Baudrate", "Data bits", " Stop bits" and "Parity" of the highlighted balance according to the table in picture 2. In the frame "COM-Port" the computer port to which the electronic balance is connected can be selected.

3.1. Terminal

After clicking on




a VT100 compatible terminal screen opens.

This terminal observes the data coming in from the balance to the PC in "auto" mode. These raw data can be recorded on file (extension "CAP") and played back afterwards.

The following functions are implemented:


Start:	Start of data display
Stop:	Stop of display
Capture:	Entering the name of the capture file and start recording
Play:	Open and play a capture file
Exit:	Terminal switch off.

On activating the capture function, a yellow light appears on the right side . Clicking on  switches recording on and off.


The seven green lights on the right side reflect the state of the port lines. Red means an active state.

3.2. Printer Setup

The printer setup can be made in frame "Print"




Picture 8 - Printer setup window

by clicking on . The field next to it shows the set printer resolution, which is in this example 300 dpi.

3.3. Timer-Calibration

The internal timer controls data acquisition very accurately. Its setting normally does not need to be changed. However, for calibration purposes, there is the possibility to accelerate or delay this timer by 5% compared to the internal PC-Clock.



Picture 9 - Timer calibration window

The timer can be adjusted by moving the track bar to the right or to the left side. A key-click on "Set" sets the new value. Clicking on sets the timer to 100%.

Alternatively, an automatic timer calibration can be performed. The duration of this process is determined by setting the time in the field next to "Duration (min)". The test starts after clicking on "Start". In the green field the time difference between timer time and internal PC-time is continuously displayed in seconds. The deviation in percents is shown in the field above. When the test has ended or been interrupted with "Stop", the optimal timer value is displayed next to "Set". Clicking on the Set-key sets the new value.

3.4. Options

In frame "Options" the following selections can be made:



Picture 10

- Log:** During data acquisition an additional log file is created. This file documents the data acquisition process in detail. Its extension is ".LOG".
- OWarning:** If this field is marked, an overwrite warning is displayed after starting the measuring in DATA ACQUISITION/MEASUREMENT if the set filename already exists.
- Check:** If this field is marked, a horizontal scroll bar appears in table "Balances". By moving to the right, all other fields of this table can be viewed. They contain the setting parameters of each balance.
- Par:** The 4 fields are related to the data format of the set balance.
- Field 1: Number of balance type
 Field 2: balance index
 Field 3: First read digit of a received data string starting from left
 Field 4: Number of read digits.
- If balance Custom auto is marked, the user can enter values into field 3 and 4 (see above). This allows the adaptation of already not supported balances to INFUSCALE.
- Win98** If INFUSCALE has been installed on Windows 98 this checkbox must be marked. In case of Windows NT, 2000, XP or 7, this field remains empty.

Key  leads back to the main menu.

4. Infusion Pumps according to EN 60601-2-24

The following chapters refer to the respective chapters of the standard EN 60601-2-24.

4.1. Volumetric Infusion Controllers, Volumetric Pumps and Syringe Pumps (50.102)

The following measuring conditions are defined in the standard:

Measuring interval:	30 sec
Set delivery rate:	<p>The measuring must be made at the following set rates:</p> <p>a) Volumetric pumps and controllers:</p> <ul style="list-style-type: none"> - minimum adjustable rate, but at least 1 mL/h - intermediate rate: 25 mL/h <p>b) Syringe Pumps:</p> <ul style="list-style-type: none"> - minimum adjustable rate, but at least 1 mL/h - intermediate rate: 5 mL/h
Measuring time:	More than 2 hours

For further details please refer to the standard or to picture 11. The evaluation with INFUSCALE is described in chapter 6.1.2.

4.2. Drip Rate Infusion Controllers and Drip Rate Infusion Pumps (50.103)

The standard does not recommend the gravimetric method, but a drip counting method.

4.3. Infusion Pumps for Ambulatory Use Type 1 (50.104)

Ambulatory pumps with continuous delivery characteristics are concerned.

Measuring interval:	900 s = 15 min
Set delivery rate:	<p>the measuring must be made at the following set rates :</p> <ul style="list-style-type: none"> - minimum adjustable rate - intermediate rate specified as typical by the manufacturer

For further details please refer to the standard or to picture 11. The evaluation with INFUSCALE is described in chapter 6.1.3.

4.4. Infusion Pumps for Ambulatory Use Type 2 (50.105)

Ambulatory pumps with non continuous delivery are concerned, i.e. repetitive delivery of equal boli during the whole delivery period. The measuring must be made with 20 and subsequent 100 shot cycles at intermediate rate as specified by the manufacturer. Further details are described in the standard and in picture 11. The evaluation with INFUSCALE is described in chapter 6.1.4.

4.5. Infusion Pumps for Ambulatory Use Type 3 (50.106)

Ambulatory Infusion and Syringe Pumps with Bolus Equipment are concerned (pumps that are able to deliver single, discrete boli)

The measuring must be made with 25 consecutive boli at the minimum setting either demanded manually or by program. Further details are described in the standard and in picture 11. The evaluation with INFUSCALE is described in chapter 6.2..

4.6. Infusion Pumps for Ambulatory Use Type 4 (50.107)

Combination of non continuous delivery and bolus delivery. Testing shall be done according to subclauses 50.104, 50.105 and 50.106 of the standard as appropriate. Please refer to the standard or to picture 11. The evaluation with INFUSCALE is described in chapters 6.1.3., 6.1.4. and 6.2..

4.7. Infusion Pumps for Ambulatory Use Type 5 (50.108)

Only for parenteral nutrition or patient controlled infusion. Use subclause 50.102 to 50.106 of the standard as appropriate. For further details please refer to the standard or to picture 11. The evaluation with INFUSCALE is described in chapters 6.1.2., 6.1.3., 6.1.4. and 6.2..

Note: Nutrition pumps must be tested with the specified nutrition solution instead of water. Before starting measuring, the density of the nutrition solution must be determined and entered into the respective field.

4.8. Other medical pumps

EN 60601-2-24 does not refer to implanted pumps or other pump systems in medical use. However INFUSCALE can be used for these pumps as well.

4.9. Overview of Infusion Pumps according to EN 60601-2-24

Sect.	Type	Min. rate ≥1ml/h	Intermed. rate	Intermed. rate +100mmHg	Intermed. Rate -100mmHg	Intermed. Rate + container 0.5 m below	Bolus min.	Bolus max.	Measur. Interval
50.102	Vol. Contr.	x	25ml/h		x				30 s
50.102	Vol. Pump	x	25ml/h	x	x	x	x	x	30 s
50.102	Syr. Pumpe	x	5ml/h	x	x	x	x	x	30 s
50.103	Drip Contr.	x	20/min						1 min
50.103	Drip Pump	x	20/min						1 min
50.104	Amb. Type1	x	x						15 min
50.105	Amb. Type2		x						
50.106	Amb. Type3						x	x	
50.107	Amb. Type4	x	x				x	x	
50.108	Amb. Type5	x	x				x	x	

Picture 11 - Set delivery rates according to IEC601-2-24, sections 50.102 - 50.108

4.10 Classification of Pumps

The British Health Ministry has published recommendations for users of medical infusion pumps. These publications are based on technical evaluations, user assessments, incident reporting and on own evaluations. Device bulletin 2003(02) addresses problems in the use and selection of infusion systems and the training of users. It raises awareness of the nature of infusion systems, their advantages and risks.

Regarding clinical application, pumps have been divided into three major categories according to the potential infusion risks. These categories A, B, and C are shown in the next table with a list of the performance parameters critical to each. The categories have been selected on the principle that, in general, the greater the potential risks associated with therapies, the higher the performance needed and the more important are the safety features.

Therapy Category	Therapy description	Patient group	Performance data
A	Drugs with narrow therapeutic margin	Any	Good long-term accuracy
	Drugs with short half-life (< 5 min)	Any	Good short-term accuracy
	Any infusion given to neonates	Neonates	Rapid alarm after occlusion Small occlusion bolus Able to detect very small air embolus Small flow rate increments Good bolus accuracy Rapid start-up time
B	Drugs, other than those with a short half-life (> 5 min)	Any except neonates	Good long-term accuracy Alarm after occlusion
	TPN Fluid maintenance Transfusions	Volume sensitive except neonates	Small occlusion bolus Able to detect small air embolus
	Application of Diamorphin (long half life of metabolite)	Any except neonates	Small flow rate increments Bolus accuracy
C	TPN Fluid maintenance Transfusions	Any except volume sensitive or neonates	Long-term accuracy Alarm after occlusion Small occlusion bolus Able to detect air embolus Incremental flow rates


Therapy categories and performance parameters

Minute-to-minute variability of flow becomes very important where drugs with short half-lives are being administered. The importance derives from the need to prevent undesirable fluctuation of effect-site concentration of the drug. There is documented evidence that minute-to-minute variability of flow can cause variation of physiological parameters and consequent difficulties in management, in both adults and neonates, where the half-life of the drug is short.

Short-term accuracy is expressed by the concept of **constancy index**. This is the shortest period during the pump's steady-state operation over which measurement of output consistently falls within $\pm 10\%$ of the mean rate. The data are derived from the flow tests performed over 24 hours at 1ml/h. Flow is recorded at 30 second intervals over the final 18 hour period. And the average rate compared with flow over each short period.

The constancy index of the pump should be less than or equal to the half-life of the drug used. Syringe pumps have a shorter constancy index than volumetric pumps i.e. a low constancy index indicates good short-term flow rate accuracy.

5. Performing the Measuring



Picture 12 - Measuring setup

5.1. The Measuring Setup

The measuring setup consists of the following components:


- Balance with RS232 interface
- Windows computer with RS232 interface COM1 .. COM4
- connection line between computer and electronic balance
- liquid container
- silicon oil for the oil film
- cannula (G18 or G21)

The balance must be installed perfectly plain, in a vibration free environment and must be calibrated. The container of the infusion pump that is to be tested must be filled with water (corresponding to ISO class III). The connected infusion line with a long cannula must be primed. In the liquid container placed on the measuring plate, a little amount of water is filled in with an oil film on the top (about 0.5 cm). This film reduces the evaporation of the water during the long measuring period.


The tip of the cannula is placed directly under the oil film, so that a continuous water pillar between the infusion container, the infusion line and the liquid container on the balance is maintained. Thus an influence of the measuring accuracy by the drop formation on the end of the cannula is prevented.

5.2. Data Acquisition

Data acquisition is started by clicking on :



The following screen appears:



Picture 13 - Measurement screen

5.2.1 Entering data

In the upper part of the screen the following parameters can be entered:

- Pump: Name of the pump under test
- Technician: Name of the technician performing the measuring
- Serial no.: Serial number of the pump under test
- Disposable: Used disposable
- Set flow rate: Set flow rate of the tested pump

Meas. Interval: Time interval between two measuring points in seconds

Evaporation: Evaporation rate in mL/h
 Density Density of the test fluid in g/mL
 Total time: Measuring time in minutes (0 or - = no time limitation)

When entering the parameters please regard the conditions, described in the IEC601-1-24 standard, and referenced in chapter 4 of this manual. By pressing TAB or ENTER the next field is highlighted for input.

The shortest time interval is 1 second. Please note, that such a short time period will increase the measuring error (see chapter 7).

Right to the "Pump" field the set balance is displayed.

Next to it is a green light that changes to red, if a new data line is transmitted to the PC.

The "Target" field shows name and path of the INFUSCALE file to be stored. By clicking on this field a new file name can be entered.

5.2.2. Control of Data Transfer

The area displayed on the lower left part of the screen serves for measuring control.

Start time	act.time	Meas. time	Time diff.	Event no.	Event time
16:01:51	16:01:55	16:01:55	00:00:04	1	00:00:02
Index	Time string	Weight	mean rate	mean error	Volume
4	00:00:04,00	0.00000	0.0000	-	0.0000
[32][32][32][32][32][32][32][48][46][48][32][32][32][32]			0.0	COM1 9600 7E1	

Picture 14 - Displayed values during data acquisition


On the bottom line the data received from the balance are shown in ASCII code. For instance [32] represents the Space-letter, [48] stands for zero. Right to it the same data appear in normal readable form. In "auto" mode these two fields allow checking whether correct data are transmitted from the balance to the PC.

5.2.3. Terminal



By clicking on



the screen of a VT100 compatible terminal appears. It allows the observe the incoming balance data and the RS232-control signals. In mode "abfr" also the signals that are transmitted from the PC to the balance are displayed.



Picture 15 - Terminal display during data acquisition

These raw data can be recorded on file (extension "CAP") and played back afterwards. On activating the capture function with , a yellow light appears on the right side. Clicking on  switches recording on and off.

The seven green lights on the right side reflect the state of the port lines. Red means an active state.

5.2.4. Start of measurement

When all parameters are entered, the measurement can be started by clicking on



If the Tare option Tare was been marked, the computer sends out a Tare-command which sets the balance to zero.

Clicking on



stops the measurement.

Below the 7 file header lines, the created file contains in each line the time and the measured weight (see next picture). Additionally, events that have been marked by pressing the blank key on the PC keyboard, are stored in a separate line. The name of the technician is automatically stored in the last line.

Example:

```
TIME INTERVAL (SEC): 30
SET FLOW RATE(ML/H): 20.0000
EVAPORATION (ML/H): 0.0000
DENSITY (G/ML): 0.9982
NAME OF PUMP : Example 5
SERIAL NUMBER : -
DISPOSABLE : -
13:26:33,25 0.00000
13:27:03,30 0.02041
13:27:33,29 0.12840
13:28:03,28 0.28762
13:28:33,27 0.45998
13:29:03,26 0.64032
13:29:33,25 0.80385
13:30:03,29 0.96463
13:30:33,28 1.12811
13:31:03,27 1.29109
13:31:33,26 1.44685
13:32:03,25 1.59950
13:32:33,29 1.74942
13:33:03,28 1.90455
13:33:33,27 2.06292
13:34:03,26 2.23061
usw.
```

Picture 16 – Infuscale data file

During the program run the following values are shown on the screen:

Upper line:


1. Starting time
2. Actual PC-time
3. Measuring time
4. Time difference between measuring time and starting time
5. Event number
6. Event time

Middle line:

1. Index
2. Stored time string
3. Stored weight
4. Calculated mean flow rate since start of measurement (mean rate)
5. Percentage deviation of the mean flow rate in relation to the set delivery rate
6. Actual delivered volume in mL

In the middle windows weight and flow rate are plotted in two graphs against time. By clicking on the keys below, the number of data points shown can be set. E.g. if key 100 is pressed, 100 measuring points are displayed. When exceeding this number, the curve is scrolling to the left. By clicking on "all" the maximum number of 10000 points can be seen.

When the key Time is pressed down, the time format of the x-axis changes from minutes to hh:mm:ss. 1 minute is displayed as 00:01:00.



Picture 17 – Data acquisition

5.2.5. Event marker

During measurement an event can be marked by pressing the blank key of the PC-keyboard. Such an event is displayed in the control area (see 5.2.2.), and is stored in the INFUSCALE file. Events are displayed on evaluation screens and tables.

5.2.6. Storage of Graphics on File and to Clipboard.


When the cursor is pointing on a graphic, a click on the right mouse key shows the following pop-up menu:



Clicking the left mouse key on one of the menu items stores the graphic on data file in either BMP or JPG format, or copies it to the clipboard from where it can be easily transferred to other programs e.g. Microsoft Word.

5.3. Simulation

The option "Simulation" does not perform a real data acquisition. It simulates the measuring by using an already existing INFUSCALE data file as data source instead of a balance. It can be selected by clicking on "Data acquisition" and "Simulation":




The source file can be opened by clicking on the input area on the top left side. The further operation is identical to chapter 5.2.


By marking the field x 10 the speed of simulation can be increased 10 fold.

6. Evaluation of the measured Data


To evaluate the measured data please return to the main menu and click on "Evaluation".




The following screen appears:



Picture 18 – Screen “Evaluation”

At this moment, the keys on the right side can only be seen in outlines without a function. However if an INFUSCALE file is loaded by clicking on , their view changes to normal. In the frame "EN 60601-2-24, 50" only those keys are activated, which lead to a reasonable evaluation according to the IEC-standard. This will be explained later in this chapter.




For opening another file, please click on the field right to "filename". A selection window appears immediately in which the file can be marked and opened.




Picture 19 – File has been loaded

The list on the left side contains the raw file data and the line index. In the center part of the screen the file parameters, which have been entered before starting the measurement, are displayed.

The name of the technician can be entered if necessary.

The keys   and  can be used anyway. Clicking on these keys lead the user directly into the program modules "Trumpet Curve", "Modified Trumpet Curve and Constancy Index" or "Bolus Evaluation". The calculations are made without restrictions regarding the minimum duration of measurement (> 2 hours).

 leads back to the main menu.

6.1. The Trumpet Curve

6.1.1. The Meaning of the Trumpet Curve according to EN 60601-2-24

In the trumpet curve the minimum and maximum deviation of the measured flow rate, which is averaged in 5 respectively 6 different observation windows, is plotted against the set delivery rate of the pump. The curve shows the mean flow rate deviation per time interval.


The trumpet curve is a useful information for the physician who has to decide whether a certain infusion pump can be used together with an intended drug.

Special drugs with short elimination half life demand high standards concerning delivery continuity. If a remedy with an elimination half life of e.g. one minute is used, flow variations of 15%/minute will cause corresponding changes of the plasma level. Thus a predictable, constant plasma level and a respective drug efficacy is hard to get.

Example:

Intravenous given insulin has an elimination half life of 15 minutes. Variations of the delivery rate of $\pm 15\%$ in the average, during a period of 40 minutes will cause similar deviations in the plasma level and, as a result, unreliable efficacy.

6.1.2. Calculation of the Trumpet Curve according to EN 60601-2-24, 50.102


Click on key , if it is activated. This is the case, if volumetric infusion controllers or syringe pumps are to be evaluated (see chapter 4.1). The time interval has to be exactly 30 s and the duration of the measurement must be more than 2 hours.

After clicking, INFUSCALE evaluates the data according to section 50.102 of the IEC-standard. For each observation period of 2, 5, 11, 19 and 31 minutes the program calculates the:

- minimum flow rate Q_{pmin}
- maximum flow rate Q_{pmax}
- maximum negative deviation of mean flow rate E_{pmin}
- maximum positive deviation of mean flow rate E_{pmax}
- percentage error A resp. B related to the set delivery rate.


These values are plotted as „Trumpet Curve“.

The next picture shows the Trumpet Curve of the second hour of example file BSP2.PMP.




Picture 20 - Trumpet Curve of the second hour

In the window right to the flow rate graph the calculated values are displayed as numbers, e.g. the maximum and minimum flow rates (38.60 resp. 10.94). The total error A (resp. B) is shown in the window right to the Trumpet Curve (A=0.07 %).

Note the row of keys on the bottom left: 

These keys allow the setting of the deviation range of the Trumpet Curve. The standard setting is $\pm 15\%$ deviation.


Key  on the upper left side of the flow rate graph changes the scaling of the ordinate. When this key is not pressed, the Y-range of the curve is identical to the range of Y-values of the data file. When this key is pressed, the maximum Y of the curve is rounded to the next even value above the maximum Y-value of the data file.

Clicking on  calculates and displays the Trumpet Curve for the last hour.

6.1.3. Calculation of the Trumpet Curve according to EN 60601-2-24, 50.104

Key **104** is activated when infusion pumps for ambulant use of type 1 are to be evaluated. (see chapter 4.3). The measuring interval in this case is 15 minutes. The specifications of section 50.104 of the IEC-standard apply.

Clicking on this key opens the input area




The startup period T1 can be entered here.

Clicking on



terminates the time input. INFUSCALE evaluates the file data now and draws the respective graphs. Each point of the flow rate curve over T1 represents the average of two data points as required in the IEC standard. The observation periods of the trumpet curve calculated over T2 are 15, 60, 150, 330, 570 and 930 minutes.



Picture 21 - Trumpet curve according to section 104

6.1.4. Calculation of the Trumpet Curve according to EN 60601-2-24, 50.105

Key **105** is activated when infusion pumps for ambulant use of type 2 are to be evaluated. (see chapter 4.4). The specifications of section 50.105 of the IEC-standard apply.

Clicking on this key opens the input area


T1 from minute to S = min

The startup period T1 can be entered here as well as the measuring interval S.

Clicking on




terminates the input. INFUSCALE evaluates the file data now, calculates n (n = 30 minutes divided by S), and draws the respective graphs. Each point of the flow rate curve over T1 represents the average of n data points as required in the IEC standard. The observation periods of the trumpet curve calculated over T2 are 1S, 2S, 5S, 11S, 19S and 31S depending on the time interval S.



Picture 22 - Trumpet curve according to section 105

6.1.5. Calculation of an optional Trumpet Curve

By clicking on  INFUSCALE can calculate and draw a Trumpet Curve for any time range, which is in the limits of the acquired data. The observation periods are 2, 5, 11, 19 and 31 minutes. The calculations are made according to section 102 of the standard, however, without restrictions regarding the minimum duration of measurement (> 2 hours).

The following input area opens immediately after clicking, in which the user can enter the desired time range for calculation.



The image shows a horizontal input field with a grey background. On the left, it is labeled "T1 from Minute" followed by a small blue icon of a trumpet. To the right of this is a white text input box containing the number "1". Further right is the word "to" followed by another white text input box containing the number "186". On the far right of the input area is a small grey button with a white right-pointing arrow.

Clicking on




terminates the input and starts the evaluation

6.1.6 Modified Trumpet Curve and Constancy Index

Short-term accuracy is expressed by the concept of constancy index. The principle is that the constancy index of the pump should be less than or equal to the half-life of the drug used.

The data from which the constancy index is calculated are derived from the flow tests performed over 24 hours at 1ml/h. Flow is recorded at 30 second intervals over the final 18 hour period. And the average rate is compared with flow over each observation window from 1,2,3,4,... 29,30 and 31 minutes. The constancy index is the shortest time period during the pump's steady-state operation over which measurement of output consistently falls within $\pm 10\%$ of the mean rate. The maximum positive and negative deviations of flow rate in each time window form the points of the modified trumpet curve. In contrast to the classical trumpet curve according to EN 60601-2-24, this curve is based on 31 observation intervals instead of 5 or 6.



Picture 23 - Modified trumpet curve and constancy index

In this example the constancy index is 4 as shown in the upper edge of the trumpet curve. Graphically the index is represented by a vertical broken line. In this example E_{pmax} is $+8.2116\%$, E_{pmin} is -8.2599% .

6.1.7. Other Commands

6.1.7.1. Set Time Range



This key opens an input area for the readjustment of the time baseline (X-axis) of the graphics.



Clicking on the key next to "Minute" sets the range back to a maximum.

6.1.7.2. Set Flow Rate Range



This key opens an input area for the readjustment of the range of flow rate (Y-axis) of the graphics.

Clicking on the key next to "Rate" sets the range back to a maximum.

6.1.7.3. Zoom



Clicking on this key displays an enlarged area of the flow rate graph. To determine its position and size please pull a rectangle with pressed left mouse key from the upper left corner to the lower right corner of the enlarged graphic area you want to display.

6.1.7.4. Unzoom



After clicking on this key, the original graphic before zooming is displayed.

6.1.7.5. Entering Comment lines



After clicking on this key the following window opens:

1.	
2.	
3.	
4.	

It contains 4 lines for entering comments, which are printed out in the graphic- as well as in the test report printout.

Clicking on **X** hides all comment lines. The OK-key below stores all lines into the comment memory.

6.1.8. Data Output

6.1.8.1. Printout of the Trumpet Curve



After clicking this key a printout of the flow rate graph and the Trumpet Curve is made.

6.1.8.2. Screen Print




After clicking on this key, a screen print is made.

6.1.8.3. Table Output



After clicking on this key, the file header, the comment lines, the single measured weight values and the calculated flow rate values are displayed in three tables. The column "Error[%]" in table 3 shows the percentage error of flow rate in relation to the set delivery rate for each data point. In columns Q1 - Q62 the mean flow rates are contained. At the end of table 3 the calculations according to the standard are shown.



Picture 24 - Table display Trumpet Curve


You can scroll through the table by pressing the left mouse key and move the mouse pointer up and down the scroll bar. For printing table 3 please click on . A screen print can be made by clicking on . The table data can be stored in an Excel file by clicking on .

The keys are used for loading and storing the comment lines.



6.1.8.4. Test Report



After clicking on this key, the INFUSCALE test report is displayed.




Picture 25 - Test report of trumpet curve evaluation

Clicking on  prints out the test report. After clicking on  the test report is stored on file in TXT, WMF, HTML, CSV, RTF, or XLS format. A data transfer to Microsoft Word and Excel is easily possible.

6.1.8.5. Storage of Graphics on File and to Clipboard.

When the cursor is pointing on a graphic, a click on the right mouse key shows the following pop-up menu:



Clicking the left mouse key on one of the menu items stores the graphic on data file in either BMP or JPG format, or copies it to the clipboard from where it can be easily transferred to other programs e.g. Microsoft Word.

Note: Auxiliary lines are not stored.

6.2. Bolus Evaluation

6.2.1. Evaluation according to EN 60601-2-24, 50.106

After clicking on key **106** the screen EVALUATION is called up, which is suitable for the assessment of pumps with bolus characteristics (section 50.106 of the standard). For this purpose weight and flow rate are simultaneously calculated and displayed.

In the following example the data file BSP6.PMQ has been loaded in screen EVALUATION before clicking on key **106**. The following graphic appears:




Picture 26 - Bolus evaluation


In the upper part the weight of the delivered liquid is displayed as a function of time with the weight at the beginning of the measurement on the lower left side and in the upper right corner the final value. On the right margin of the weight curve the minimum and the maximum values are shown as numbers (0,000 bzw. 44,73).

The steps, which represent the delivery of limited liquid volumes (boli) with a high flow rate, can be clearly identified. In this example a total of 9 liquid boli have been delivered.

The lower part of the screen shows the flow rate curve. The step increase of the flow rate during a weight step and its decrease in the time between two boli can be seen very clearly.


On the right margin of the flow rate curve the minimum and the maximum rates are shown as numbers (-0,2 resp. 109,2).

Key  on the upper left side of the graphs changes the scaling of the ordinate. When this key is not pressed, the Y-range of the curve is identical to the range of Y-values of the data file. When this key is pressed, the maximum Y of the curve is rounded to the next even value above the maximum Y-value of the data file.

If key  is pressed, two dotted lines are drawn in the flow rate curve marking the minimum and maximum flow rates.

The lists on the right side of the screen display the calculated values as numbers.

6.2.2. Bolus Evaluation with adjustable Time Range

By clicking on key  in screen "Evaluation" an INFUSCALE-file with adjustable time range can be calculated and displayed. The following input area opens immediately after clicking, in which the user can enter the desired time range for calculation.



T1 from Minute to 

Clicking on



terminates the input and start the evaluation.

6.2.3. Other Commands

6.2.3.1. Set Time Range



This key opens an input area for the readjustment of the time baseline (X-axis) of the graphics.



Clicking on the key next to "Minute" sets the range back to maximum.

6.2.3.2. Set Weight Range



This key opens an input area for the readjustment of the weight range (Y-axis) of the graphics.

Clicking on the key next to "Weight" sets the range back to maximum.

6.2.3.3. Set Flow Rate



This key opens an input area for the readjustment of the range of flow rate (Y-axis) of the graphics.

Clicking on the key next to "Rate" sets the range back to a maximum.

6.2.3.4. Zoom



After clicking on this key an enlarged area of the flow rate graph can be displayed. To determine its position and size please pull a rectangle with pressed left mouse key from the upper left corner to the lower right corner of the enlarged graphic area you want to display.

6.2.3.5. Unzoom



After clicking this on key, the original graph before zooming is displayed.


6.2.3.6. Curve Measuring



After clicking on this key, the graphic curves can be measured. Holding the left mouse key pressed, at first the left vertical line that indicates the lower measuring limit is drawn to the point of curve where the calculation shall start. Then, in the same way, move the upper measuring limit line.

INFUSCALE calculates all relevant values between these two borders and displays it in the two listing areas at the right side of the screen.

Example:



Picture 27 - Bolus measuring

In this example at first an enlarged area of the graph based on file BSP6.PMQ between minute 60 and 80 was drawn after clicking on "Time". Then a measuring between minute 64,01 and 70,00 has been made. The upper list shows, that the mean flow rate between the two lines is 49,4039^omL/h. The bolus between the lines has a volume of 4.9404 mL. The maximum flow rate between the vertical lines is 103,4081 mL/h, and 105,9 mL/h in the total displayed area.

6.2.3.7. Entering Comment lines



After clicking on this key the following window opens:

1.	
2.	
3.	
4.	

X
OK

It contains 4 lines for entering comments, which are printed out in the graphic- as well as in the test report printout.

Clicking on **X** hides all comment lines. The OK-key below stores all lines into the comment memory.

6.2.4. Data Output

6.2.4.1. Printout of the Graphs



After clicking on this key a printout of the two graphs (weight and flow rate) is made.

6.2.4.2. Screen Printout



After clicking on this key, a screen printout is made.

6.2.4.3. Table Output



After clicking on this key, the file header, the comment lines, the single measured weight values and the calculated flow rate values are displayed in three tables. At the end of table 3 the evaluation results are shown.

The screenshot shows a software window titled "Table display" with a file path "C:\Infu701e\data\Bsp2.pmp" and "Points" set to "372". It contains three tables:

Index	String	Value
1	MESSINTERVALL (SEC):	30
2	FOERDERRATE (ML/H):	30.000
3	VERDUNSTUNG (ML/H):	0.1200
4	DICHTE (G/ML):	0.9980
5	NAME DER PUMPE :	Rollen
6	SERIENNUMMER :	12678
7	EINMALARTIKEL :	-
8	TECHNIKER :	-

Index	Time (min)	Time (h:m:s)	vWeight [g]	Rate [ml/h]	Error [%]	vWeight	Rate
0	0.0	00:00:00	0.00000	0.00000	-100.0	min	
1	0.5	00:00:30	0.09000	10.94164	-63.5		min
2	1.0	00:01:00	0.37000	33.78733	12.6		
3	1.5	00:01:30	0.54000	20.56088	-31.5		
4	2.0	00:02:00	0.83000	34.98974	16.6		
5	2.5	00:02:30	1.02000	22.96569	-23.4		
6	3.0	00:03:00	1.33000	37.39455	24.6		
7	3.5	00:03:30	1.53000	24.16810	-19.4		
8	4.0	00:04:00	1.83000	36.19214	20.6		
9	4.5	00:04:30	2.08000	30.18012	0.6		
10	5.0	00:05:00	2.34000	31.38253	4.6		
11	5.5	00:05:30	2.63000	34.98974	16.6		
12	6.0	00:06:00	2.84000	25.37050	-15.4		

Picture 28 - Table display, Bolus evaluation


You can scroll through the table by pressing the left mouse key and moving the mouse pointer up and down the scroll bar. For printing out table 2 please click on . A screen print can be made by clicking on . The table data can be stored on an Excel file by clicking on .

The keys are used for loading and storing the comment lines.



6.2.4.4. Test Report



After clicking on this key, the INFUSCALE test report is displayed




Picture 29 - Test report of bolus evaluation

Clicking on  prints out the test report. After clicking on  the test report is stored on file in TXT, WMF, HTML, CSV, RTF, or XLS format. A data transfer to Microsoft Word and Excel is easily possible.

6.2.4.5. Storage of Graphics on File and to Clipboard.

When the cursor is pointing on a graphic, a click on the right mouse key shows the following pop-up menu:



Clicking the left mouse key on one of the menu items stores the graphic on data file in either BMP or JPG format, or copies it to the clipboard from where it can be easily transferred to other programs e.g. Microsoft Word.


Note: Auxiliary lines are not stored.

7. Measuring Accuracy

In this chapter we want to give you an overview about the influence of evaporation and the interrogation error, which results from the communication and data processing between electronic balance and computer. Other failures, like the cannula effect, that is based on the lift of the cannula, the influence of temperature and geographically latitude to the exactness of measuring are not as important and will be neglected here.


7.1. Evaporation

The next two graphics show the effect of the evaporation:



Picture 30 - Evaporation without oil film

In this picture the evaporation container was a bucket, containing 4 litres of water without an oil film. Within 67 hours 122.75 mL water evaporated, which can be observed in the upper graphic as falling weight curve. Respectively the evaporation curve below moves between -2.284 and -1.443 . This is equivalent to a mean evaporation rate of 1.83 mL/h.



Picture 31 - Evaporation with oil film

The volume of the container used for the next picture was about 100 mL with an oil film on the liquid surface. Within 870 minutes only 0.02 mL of water evaporated. The equivalent mean evaporation rate is 0.0014 mL/h.

The comparison of these pictures demonstrate, that evaporation can lead to large measuring errors. Therefore, it is recommended to work always with an oil film, when the measuring has to be very exact or during the measuring of pumps with low flow rates. Before this, the evaporation rate has to be measured in a first prospective run. This rate has to be entered into the program before starting the measuring (see chapter 5.2.1.).

7.2. Interrogation Error

Provided, the operation system has been set, as described in section 3.1, the timer error of the computer is one to two orders of magnitude below the interrogation error and can, therefore, be neglected.

7.2.1. Mettler Balance AE100 (Continuous Operation)

During continuous operation a data string is sent periodically to the computer by the Mettler balance AE100. The time distance between two data transmissions is maximally 440 ms. The interrogation by the computer is done whenever the set time interval has run out, and is not synchronised, so that the moment of interrogation may be delayed for a maximum of 440 ms. This means that the computer takes a weight value that was applied maximally 440ms before the correct moment of interrogation. With an average measuring interval of 30 sec the error is maximal 1.5%.

During the determination of the trumpet curve 4 to 62 measuring values are averaged. So the error will be reduced to about 0.38%.

The next table (picture 23) shows the max. interrogation error and the average errors for the smallest observation window $p = 4$ of the trumpet curve dependent from the time interval. The shown maximum errors are extreme values.

Time interval	max. interrogation error	Mean error (p=2 min)
1 s	44%	11%
3 s	14.6%	3.7%
6 s	7.3%	1.8%
10 s	4.4%	1.1%
30 s	1.5%	0.38%
60 s	0.73%	0.18%
100 s	0.44%	0.11%
300 s	0.15%	0.04%
600 s	0.07%	0.02%
900 s	0.05%	0.01%

Picture 32 - - Maximum error dependent of the time interval, Mettler auto

7.2.2. Sartorius Balance MC1 (automatic data transmission)

The Sartorius balance MC1 transmits periodically a data string to the computer (Sartorius auto). The time interval between two consecutive data transmissions is max. 220 ms. The interrogation by the computer is done whenever the set time interval has run out and is not synchronised, so that the moment of interrogation may be delayed by maximally 220 ms. This means that the computer takes a weight value that was applied maximally 220 ms before the correct moment of interrogation.

The next table (picture 24) shows the maximum interrogation error and the average errors for the smallest observation window $p = 4$ of the trumpet curve dependent from the time interval.

Time interval	Max. interrogation error	mean error (p=2 min)
1 s	22%	5.5%
3 s	7.3%	1.8%
6 s	3.6%	0.9%
10 s	2.2%	0.6%
30 s	0.73%	0.18%
60 s	0.36%	0.09%
100 s	0.22%	0.06%
300 s	0.07%	0.02%
600 s	0.04%	0.01%
900 s	0.02%	0.005%

Picture 33 - Maximum error dependent of the time interval, Sartorius auto

7.2.3. Sartorius Balance MC1 (request mode)

In the request mode (Sartorius abfr) the Sartorius balance MC1 is requested by the computer periodically after the time interval has run out. As response a data string is sent back to the computer. The maximum measured error is 60 ms in a test situation. In normal operation a maximal error of 30 ms can be expected.

The next table (picture 25) shows that this mode is more accurate than the automatic data transmission mode at least for measurements with short time intervals.

Time interval	max. interrogation error	mean error (p=2min)
1 s	6%	1.5%
3 s	2%	0.5%
6 s	1%	0.25%
10 s	0.6%	0.15%
30 s	0.2%	0.05%
60 s	0.1%	0.03%
100 s	0.06%	0.02%
300 s	0.02%	0.005%
600 s	0.01%	0.002%
900 s	0.005%	0.001%

Picture 34 - Maximum error dependent of the time interval, Sartorius abfr

8. Symbols and Formula

T = Total measuring time in minutes

$$T = t_k - t_1 \text{ [min]}$$

k = Number of measuring points

r = Set delivery rate [mL/h]

i = Measuring index $i = 1$ to k

Δt = Time interval [min] between two consecutive measuring points

t_i = Measuring time [min] at point i

p = Observation window index with $p = 4, 10, 22, 38, 62$

$t_{p,i}$ = Observation window with

$$t_{p,i} = t_{i+p} - t_i = p \cdot \Delta t$$

W_i = Weight [g] at the measuring point i corrected by the evaporation loss

d = Density of water = 0.998 [g/mL] at 20 °C

Q_i = Flowrate [mL/h] between two consecutive measuring points

$$Q_i = \frac{60 \cdot (W_i - W_{i-1})}{d \cdot \Delta t}$$

$W_{p,i}$ = Weight [g] delivered during the observation window $t_{p,i}$

$$W_{p,i} = W_{i+p} - W_i$$

W = Total weight [g] delivered during the total time T corrected by the evaporation loss

$$W = W_k - W_0$$

Q = Average flowrate [mL/h] during the total time T

$$Q = \frac{60 \cdot W}{d \cdot T}$$

$Q_{p,i}$ = Average flowrate [mL/h] during the observation window $t_{p,i}$

$$Q_{p,i} = \frac{60 \cdot W_{p,i}}{d \cdot t_{p,i}}$$

Q_{pmax} = Maximum $Q_{p,i}$ [mL/h] for each observation window $t_{p,i}$ and each value of p

Q_{pmin} = Minimum $Q_{p,i}$ [mL/h]] for each observation window $t_{p,i}$ and each value of p

E_{pmax} = Highest positive deviation of Q_{pmax} for each p, expressed as percentage deviation from the set rate r.

$$E_{pmax} = \frac{Q_{pmax} - r}{r} \cdot 100 \%$$

E_{pmin} = Highest negative deviation of Q_{pmax} for each p, expressed as percentage deviation from the set rate r.

$$E_{pmin} = \frac{Q_{pmin} - r}{r} \cdot 100 \%$$

A = Percentage total error of the average flow rate Q, related to the set rate r during the total time T.

$$A = \frac{Q - r}{r} \cdot 100 \%$$