

Revision No.: ZJ0003

Issue Date: Jul 1st, 2015

Avian influenza virus N8 Real Time RT-PCR Kit User Manual

For Research Use Only

REF MBS598228 - Instrument III, IV

For use with ABI Prism®7000/7300/7500/7900/Step One Plus; iCycler iQ™4/iQ™5; Smart Cycler II; Bio-Rad CFX 96; Rotor Gene™ 6000; Mx3000P/3005P; MJ-Option2/Chromo4; LightCycler®480 Instrument



1. Intended Use

Avian influenza virus N8 real time RT-PCR kit is used for the detection of gene N9 of avian influenza in human nasal and pharyngeal secretions and bird feces by using real time PCR systems.

2. Principle of Real-Time PCR

The principle of the real-time detection is based on the fluorogenic 5' nuclease assay. During the PCR reaction, the DNA polymerase cleaves the probe at the 5' end and separates the reporter dye from the quencher dye only when the probe hybridizes to the target DNA. This cleavage results in the fluorescent signal generated by the cleaved reporter dye, which is monitored real-time by the PCR detection system. The PCR cycle at which an increase in the fluorescence signal is detected initially (Ct) is proportional to the amount of the specific PCR product. Monitoring the fluorescence intensities during Real Time allows the detection of the accumulating product without having to re-open the reaction tube after the amplification.

3. Product Description

Highly pathogenic avian influenza (HPAI) caused by certain subtypes of influenza A virus in animal populations, particularly chickens, poses a continuing global human public health risk. Avian influenza virus N8 real time RT-PCR kit contains a specific ready-to-use system for the detection of gene N8 of avian influenza virus by Reverse Transcription Polymerase Chain Reaction (RT-PCR) in the real-time PCR system. The master contains Super Mix for the specific amplification of the virus RNA. The reaction is done in one step real time RT-PCR. The first step is a reverse transcription (RT), during which the avian influenza virus RNA is transcribed into cDNA. Afterwards, a thermostable DNA polymerase is used to amplify the specific gene fragments by polymerase chain reaction. Fluorescence is emitted and measured by the real time systems' optical unit during the PCR. The detection of amplified avian influenza virus DNA fragment is performed in fluorimeter channel FAM with the fluorescent quencher BHQ1.

4. Kit Contents

| Ref. | Type of reagent | Presentation | 25rxns |
|------|-----------------------|---------------|--------|
| 1 | N9 Super Mix | 1 vial, 480µl | |
| 2 | RT-PCR Enzyme Mix | 1 vial, 28µl | |
| 3 | Molecular Grade Water | 1 vial, 400µl | |
| 4 | N9 Positive Control | 1 vial, 30µl | |

Analysis sensitivity: 1×10^3 copies/ml;

Note: Analysis sensitivity depends on the sample volume, elution volume, nucleic acid extraction methods and other factors. If you use the RNA extraction kits recommended, the analysis sensitivity is the same as it declares. However, when the sample volume is dozens or even hundreds of times greater than elution volume by some concentrating method, it can be much.

5. Storage

- All reagents should be stored at -20°C. Storage at +4°C is not recommended.
- All reagents can be used until the expiration date indicated on the kit label.
- Repeated thawing and freezing (> 3x) should be avoided, as this may reduce the sensitivity of the assay.
- Cool all reagents during the working steps.
- Super Mix should be stored in the dark.

6. Additionally Required Materials and Devices

- Biological cabinet
- Vortex mixer
- Cryo-container
- Sterile filter tips for micro pipets
- Disposable gloves, powderless
- Refrigerator and Freezer
- Desktop microcentrifuge for "eppendorf" type tubes (RCF max. 16,000 x g)
- Real time PCR system
- Real time PCR reaction tubes/plates
- Pipets (0.5µl – 1000µl)
- Sterile microtubes
- Biohazard waste container
- Tube racks

7. Warnings and Precaution

- Carefully read this instruction before starting the procedure.
- For in vitro diagnostic use only.
- This assay needs to be carried out by skilled personnel.
- Clinical samples should be regarded as potentially infectious materials and should be prepared in a laminar flow hood.
- This assay needs to be run according to Good Laboratory Practice.
- Do not use the kit after its expiration date.
- Avoid repeated thawing and freezing of the reagents, this may reduce the sensitivity of the test.
- Once the reagents have been thawed, vortex and centrifuge briefly the tubes before use.
- Prepare quickly the Reaction mix on ice or in the cooling block.
- Set up two separate working areas: 1) Isolation of the RNA/ DNA and 2) Amplification/ detection of amplification products.
- Pipets, vials and other working materials should not circulate among working units.
- Use always sterile pipette tips with filters.
- Wear separate coats and gloves in each area.
- Do not pipette by mouth. Do not eat, drink, smoke in laboratory.
- Avoid aerosols

8. Sample Collection, Storage and transport

- Collected samples in sterile tubes;
- Specimens can be extracted immediately or frozen at -20°C to -80°C.
- Transportation of clinical specimens must comply with local regulations for the transport of etiologic agents

9. Procedure


9.1 RNA-Extraction

Different brand RNA extraction kits are available. You may use your own extraction systems or the commercial kit based on the yield. For the RNA extraction, please comply with the manufacturer's instructions. The recommended Extraction kit is as follows:

| Nucleic Acid Isolation Kit | Cat. Number | Manufacturer |
|---|-----------------|--------------|
| RNA Isolation Kit | ME-0010/ME-0012 | ZJ Biotech |
| QIAamp Viral RNA Mini extraction Kit (50) | 52904 | QIAGEN |

9.2 RT-PCR Protocol

The Master Mix volume for each reaction should be pipetted as follows:



- The volumes of Super Mix and Enzyme Mix per reaction multiply with the number of samples, which includes the number of controls, standards, and sample prepared. Molecular Grade Water is used as the negative control. For reasons of unprecise pipetting, always add an extra virtual sample. Mix completely then spin down briefly in a centrifuge.
- Pipet 20µl Master Mix with micropipets of sterile filter tips to each of the Real time PCR reaction plate/tubes. Separately add 5µl RNA sample, positive and negative controls to different reaction plate/tubes. Immediately close the plate/tubes to avoid contamination.
- Spin down briefly in order to collect the Master Mix in the bottom of the reaction tubes.
- Perform the following protocol in the instrument:

| | |
|--|----------|
| 45°C for 10min | 1cycle |
| 95°C for 15min | 1cycle |
| 95°C for 15sec, 60°C for 1min (Fluorescence measured at 60°C) | 40cycles |

| Selection of fluorescence channels | |
|------------------------------------|---------------------|
| FAM | Target Nucleic Acid |

- If you use ABI Prism® system, please choose "none" as passive reference and quencher.

10. Threshold setting: just above the maximum level of molecular grade water.

11. Quality control: Negative control and positive control must be performed correctly, otherwise the sample results is invalid.

| Channel | | Ct value |
|-------------------------------------|--|----------|
| Control | | FAM |
| Molecular Grade Water | | UNDET |
| Positive Control(qualitative assay) | | ≤35 |

12. Data Analysis and Interpretation

The following sample results are possible:

| | Ct value | Result Analysis |
|----|----------|---|
| | FAM | |
| 1# | UNDET | Below the detection limit or negative |
| 2# | ≤38 | Positive |
| 3# | 38~40 | Re-test; If it is still 38~40, report as 1# |

For further questions or problems, please contact our technical support _____