

Assisted Reproductive Technologies





Founded in 1987, Cryo Bio System is a division of IMV Technologies, the world leader in animal reproduction biotechnologies, based in L'Aigle, Normandy – France.

Cryo Bio System is exclusively dedicated to human life sciences and biodiversity.

Combined with IMV Technologies Group's expertise in cryopreservation of biological samples, Cryo Bio System provides innovative solutions to the scientific community through a range of High Security storage products for biorepositories, biobanks and biological resources centres in various fields and applications:

- Cancer research, epidemiological studies
- Blood transfusion haemovigilance
- Regenerative medicine: cell and genetic therapy, stem cells
- Genomics, proteomics
- Medical and pharmaceutical research
- Medico-legal storage
- Military medicine
- Biodiversity: animal and vegetal germplasm conservation
- Assisted Reproductive Technologies (ART).

Our international presence through subsidiaries (Italy, Netherlands, United States, China and India), and a network of distributors in over 60 countries offers customers both high quality products and expert technical support.



This catalog demonstrates our commitment to the Assisted Reproductive Technologies (ART) field with a complete product line that meets, in particular, all requirements of cryopreservation.

The solutions we have designed for humans make it possible to deep-freeze sperm, embryos, ovarian cortex, as well as to vitrify oocytes and embryos at various stages. We also have developed the unique concept of long-term high-security cryopreservation, an essential guarantee of storage safety and quality for all samples.

C O N T E N T S

ICI / IUI / IVF ▼

Recisperme	4
Isolating transport container	4
Makler® counting chamber	5
Leja counting chamber	5
ICSI, holding and biopsy pipettes	6
Aspiglaire	7
Craft™ pump	7
Intracervical insemination catheter for CBS™ High Security straws	8
Intra-uterine insemination catheter with plunger	9

Vitrification ▼

HSV High Security Vitrification kit	10
Opening device for HSV kit	11

Cryopreservation ▼

CBS™ High Security cotton-plugged sperm straw	12
CBS™ High Security hydrophobic plugged embryo straw	14
CBS™ High Security unplugged tissue 2 mL straw	15
PETG cotton-plugged sperm straw	16
PETG cotton-plugged embryo straw	16

Accessories ▼

LabPal printer	17
Micro-aspirator	17
CBS™ visotubes and goblets	18
Visotubes and goblets	19
Canisters	20

Equipment ▼

SIDE	21
SYMS II / SYMS	22
PACE	24
MAPI	25
Controlled rate freezers from Planer	26



Recisperme

SPERM COLLECTION VESSEL

007426 Pack of 20 Recisperme



DESCRIPTION / USE

Spermogram, sperm culture, spermocytogram (insemination, in vitro fertilization, diagnosis).

Sterilized by irradiation.

Individually packaged in box of 20 units.

TECHNICAL SPECIFICATIONS

The innovative design of the Cryo Bio System Recisperme allows sperm to be collected directly into a round bottom tube, eliminating the need for sample transfer and enhancing sample quality by avoiding tricky handling.

Furthermore, thanks to a graduated tube inside the device the sperm volume can be accurately read (contents 10 mL).

Preassembled sterile unit comprised of:

- a funnel-shaped upper part;
- a graduated 10 mL tube with a round base and blue stopper;
- a support; and
- an adhesive cover.



INSTRUCTIONS FOR USE

- After ejaculation directly into the vessel, the graduated tube is held in place to facilitate flow;
- The adhesive cover can be put on to guarantee that the sample remains sterile;
- Then the funnel is removed to avoid potential contamination, and the tube is closed with the stopper and sent to the laboratory for analysis.

SHELF-LIFE

Three years from date of manufacture.

Isolating transport container

015937 Bag of 5 isolating transport containers



DESCRIPTION / USE

For transportation of prepared sperm prior to intra-uterine insemination.

High density Polyethylene container with Polyurethane foam insert precut to hold one tube.

All commercially available 5 to 15 mL round or conical bottom culture tubes fit the size of the precut holder.

TECHNICAL SPECIFICATIONS

Length: 120 mm

Diameter: 65 mm.

Two hours prior to use, the container can be placed at 37°C. This ensures that tubes filled with prepared sperm are at an optimal temperature for transportation from the laboratory to the clinic where the insemination will be performed, guaranteeing sperm mobility by preventing thermal shock.

Guarantees sperm mobility by avoiding thermal shock.

INSTRUCTIONS FOR USE

- Open the container by pulling up the lid;
- Place the tube into the opening situated in the middle of the foam insert;
- Close the container by putting the lid back on.

The container is now ready to be transported to the clinic for insemination.

If necessary, the container can be labeled for identification.

Maximum duration of transportation: 15 minutes.



Makler® counting chamber

DESCRIPTION / USE

Rapid sperm analysis.

Counting chamber and cover glass for counting and analysis of spermatozoa.

TECHNICAL SPECIFICATIONS

Reusable: easily cleaned with a non-bleach disinfectant solution. The depth of 10 microns eliminates blurring and allows sperm to move freely. The applied sample is observed in one focal plane.

A 1 mm² grid, subdivided into 100 squares, each one of 0.1 x 0.1 mm, is on the cover glass, eliminating the need for insertion and removal of a grid into the microscope eye.

INSTRUCTIONS FOR USE

Easy-to-use: sperm count performed on undiluted specimen.

The number of spermatozoa counted in any strip of 10 squares of the grid indicates their concentration in millions/mL.

No additional factors are necessary for calculation.

Observation of color fringes at the four contact points, provides a self-controlled test for accuracy. The cover glass can never be raised by the applied sample.

Repeated use is possible with complete accuracy and without further calibration.

- 016261** Makler® counting chamber (complete with cover glass and grid)
- 016262** Extra cover glass with grid




LEJA counting chamber

- 019319** Box of 25 disposable Leja slides – 2 chambers



- 019320** Box of 25 disposable Leja slides – 4 chambers



DESCRIPTION / USE

High quality disposable counting chambers for semen analysis. Due to their special characteristics, Leja slides will both save time and increase the level of accuracy and precision of semen analysis.

All slides are covered with a special coating to prevent air bubble formation and to prevent semen from sticking to the chamber surface. The resin and ink used are both non-toxic.

The Leja slides have an excellent low limit of quantification and can handle high sperm concentrations. This has a direct clinical impact: you will be better prepared to determine the thresholds for the various modalities of Assisted Reproduction

Technologies (waiting for natural conception, intrauterine insemination, in-vitro fertilization or intracytoplasmic sperm injection) and to decide which treatment is optimal for a couple.

The high level of accuracy and precision of the Leja chambers will allow your results to be compared with other laboratories' scientific data concerning the relationship between semen parameters (concentration, motility, viscosity) and fertility.



TECHNICAL SPECIFICATIONS

Slides with 2 or 4 chambers
Depth: 20 microns

ICSI, holding and biopsy pipettes

022166	ICSI pipette – angle 20°	10
022167	ICSI pipette – angle 30°	10

ICSI pipettes



DESCRIPTION / USE

ICSI injection pipettes are used to perform intracytoplasmic sperm injection (ICSI).

TECHNICAL SPECIFICATIONS

The unique shape and taper of our ICSI pipettes provide for smooth movement of its content. A few sperm can be loaded for injecting a few oocytes at one time. The short and sharp fine spike at the tip allows for easy, and non-traumatizing penetration of the oocyte.

The 50V series is our classic design with an internal diameter of 5 μm , a spiked tip bevelled at 30°. They are supplied with an angle of 20° or 30° with a tip to elbow length of 550 μm .

Individually packaged in box of 10 units.

SHELF-LIFE

Three years from date of manufacture.

022168	Holding pipette – outer diameter 100 μm – angle 20°	10
022169	Holding pipette – outer diameter 100 μm – angle 30°	10
022170	Holding pipette – outer diameter 120 μm – angle 20°	10
022171	Holding pipette – outer diameter 120 μm – angle 30°	10

Holding pipettes



DESCRIPTION / USE

Holding pipettes are used to hold an oocyte or an embryo in position during ICSI and other micro-manipulation procedures.

TECHNICAL SPECIFICATIONS

The classic series has an internal diameter of 20 μm and an outer diameter of 100 μm or 120 μm .

All our holding pipettes are supplied with a tip angle of 20° or 30° with a tip to elbow length of 650 μm .



Individually packaged in box of 10 units.

SHELF-LIFE

Three years from date of manufacture.

022172	Pipette for blastomere biopsy – angle 30°	10
--------	---	----

Blastomere biopsy pipettes



DESCRIPTION / USE

Blastomere biopsy pipettes are designed to remove blastomeres from embryos for preimplantation genetic diagnosis (PGD).

TECHNICAL SPECIFICATIONS

They are supplied with a tip angle of 30° and a tip to elbow length of 550 μm .

Individually packaged in box of 10 units.

SHELF-LIFE

Three years from date of manufacture.

Aspiglaire

ENDOCERVICAL ASPIRATOR FOR HÜHNER TEST

CE marked



DESCRIPTION / USE

Hühner test (post-coital sample collection).
Mucus viscosity evaluations.
Endo- or exocervical mucus sampling.
Sterilized by irradiation.
Individually packaged in quantities of 1 or 5 units
in box of 20 or 40 units.

TECHNICAL SPECIFICATIONS

The Aspiglaire catheter has two components:

- a clear, flexible sheath with a round tip; and
- an inner blue plunger.

No syringe needed.
Length: 280 mm
External diameter: 2.92 mm

INSTRUCTIONS FOR USE

- Slowly insert the catheter into the cervical canal until reaching the endocervical sampling site;
- Aspirate the cervical mucus by pulling the plunger which automatically stays in the desired position to allow the mucus to be drawn into the tube;
- Gently remove the catheter from the cervix;
- Apply the sample on a slide for examination.

SHELF-LIFE

Three years from date of manufacture.

- 007447** Box of 20 Aspiglaire (individually packaged)
007448 Box of 40 Aspiglaire (8 cells containing 5 units)



Rocket Craft™ Pump

CE marked

- 005553** Rocket Craft™ pump - Vacuum pump for oocyte recovery
015621 Pedal for Craft™ pump
006553 Craft™ pump tubing
018204 Craft™ pump overflow reservoir
018268 Specific cap for Craft™ pump reservoir
019270 Luer-Lock tubing adaptor for Craft™ pump tubing

DESCRIPTION / USE

Oocyte aspiration.
The Craft™ pump is specifically designed to provide a smooth low volume vacuum at a pre-determined negative pressure.
12-month warranty.

TECHNICAL SPECIFICATIONS

Vacuum is activated by a foot operated switch controlled by the surgeon performing oocyte collection.
Dimensions: 20 x 14 x 12 cm
Weight: 3.2 kg
Vacuum: 0 - 400 mm Hg



Intracervical insemination catheter

FOR CBS™ HIGH SECURITY SPERM STRAWS

CE marked

015821 Box of 20 sterile ICI catheters for CBS™ High Security sperm straws



DESCRIPTION / USE

Artificial intracervical insemination with thawed sperm.

For use with CB™ 0.3 and 0.5 mL straws.

Sterilized by irradiation.

Individually packaged in box of 20 units.

TECHNICAL SPECIFICATIONS

The intracervical insemination catheter consists of 3 parts:

- a flexible PVC outer sheath with safety sleeve.
Length: 24.8 cm
Internal diameter: 3.7 mm
External diameter: 4.1 mm
- a white PETG guide to position the straw in the outer sheath.
Length: 13 cm
Internal diameter: 2.5 mm
External diameter: 3.1 mm
- and a white PETG plunger to expel the contents of the straw.
Length: 28 cm
External diameter: 1.5 mm

INSTRUCTIONS FOR USE

- Thaw the straw and release the safety plug by pinching it between finger and thumb;
- Decontaminate the straw with a bactericidal agent;
- Cut off the seal at the plug side of the straw;
- Place the straw in the clear sheath, sample side first;
- Make sure the opening of the straw is correctly positioned in the sheath;
- Introduce the catheter into the cervical canal and insert the plunger inside the guide;
- Gently push away the plug of the straw and inject the sperm very slowly.

SHELF-LIFE

Three years from date of manufacture.

Intra-uterine insemination catheter with plunger

CE marked

007446 Box of 20 sterile IUI catheters with plunger

DESCRIPTION / USE

Intra-uterine insemination with processed sperm.
Sterilized by irradiation.
Individually packaged in box of 20 units.

TECHNICAL SPECIFICATIONS

The intra-uterine insemination catheter consists of 2 parts:

- a clear flexible PVC sheath equipped with
- a thin and flexible cannula with a soft tip for easy entry into the uterine cavity.

No syringe needed.

A plastic plug guarantees sterility when the processed sperm must be transported to the clinic prior to insemination.

Length: 28 cm

External diameter: 3 mm

External diameter of the tip: 1 mm

INSTRUCTIONS FOR USE

- Pull on the plunger to draw the required volume of sperm into the catheter;
- Insert a speculum;
- Introduce the catheter into the cervical canal, then through the internal cervical orifice;
- Slowly inject the sperm into the uterine cavity;
- Keep the patient lying down for a few minutes.



SHELF-LIFE

Three years from date of manufacture.



Vitrification

HSV High Security Vitrification kit



DESCRIPTION / USE

Vitrification of oocytes and embryos.

The HSV High Security Vitrification kit is comprised of 3 parts:

- a High Security ionomeric resin straw;
- a capillary tube with a pre-formed gutter attached to a colored handling rod; and
- a blue plastic insertion device.

Sterilized by irradiation.

Packaged in a peel off blister pack (4 kits).

IDENTIFICATION

The HSV kit can be labeled using LN2 resistant labels.

FILLING AND SEALING

- Prepare the identification label (liquid nitrogen resistant) for the HSV straw and apply it approximately 20 mm (0.8 inch) from the flared end of the straw;
- Connect the end of the blue plastic insertion device to the colored end of the handling rod;
- Prepare the sample for vitrification according to laboratory protocol;
- Using a micro-pipette, carefully deposit the sample into the gutter one millimeter from the end. The drop holding the sample must be under 0.5 μ L. Maximum of 2 oocytes or embryos. In this case, use a droplet by sample;
- Immediately place the pre-formed gutter with handling rod in the flared end of the straw. Push until the rectangular portion of the blue handler comes into contact with the flared end of the straw;
- Slightly pinch the straw between your thumb and finger and remove the insertion device;
- While still holding the straw in place, seal the open end using a SYMS sealer;
- Hold the straw using tweezers in the area of the handling rod;
- Quickly plunge the entire straw into liquid nitrogen vertically;
- Gently stir the straw in liquid nitrogen for a few seconds so as to avoid formation of an isolating air bubble layer around the straw.

Pack of 4 sterile HSV High Security Vitrification kits

019916	red	019917	white	019918	yellow
019919	purple	019920	blue	019921	green
022406	Opening device for HSV kit				

CE marked
FDA 510(k) premarket notification
clearance # K092398
for 4-8 cell and blastocyst stage embryos



STORAGE IN LIQUID NITROGEN

Cryo Bio System carries a line of accessories for storage and management of HSV kits in liquid nitrogen.



THAWING

Thawing of the specimen should be done in a controlled environment, such as a laminar flow hood.

- Prepare the thawing / dilution media;
- Identify the HSV straw;
- Transfer the straw from the storage container to a transport dewar filled with liquid nitrogen;
- Lift the straw enough to expose the colored handling rod. Make sure the end with the sample remains immersed in the liquid nitrogen;
- While holding the straw, use the opening device for HSV kit to section the straw;
- Immediately (within 2 seconds), plunge the gutter into the first dilution media.

SHELF-LIFE



Three years from date of manufacture.



Opening device for HSV kit

Cryopreservation

CBS™ High Security cotton-plugged sperm straws



DESCRIPTION / USE

Storage and preservation of sperm.

The special design and the mode of filling and sealing of CBS™ High Security straws for cryopreservation of human biological samples makes them fully hermetic, especially in liquid nitrogen.

Manufactured from biocompatible materials, they are used in medically assisted procreation techniques and particularly for sperm preservation. The CBS™ straw is a clear, flexible tube made of ionomeric resin.

Length: 133 mm

Internal diameter: 2.25 mm

Sterilized by irradiation.

Model CBS™ 0.3 mL

This straw is composed of two distinct parts separated by a sliding white cotton safety plug. The first part with a working volume of 0.3 mL holds the sperm sample, while the second permits introduction of a colored identification tube and serves as a handle in liquid nitrogen.

Model CBS™ 0.5 mL

Same technical specifications as the 0.3 mL straw, except that the working volume is 0.5 mL, and the outer part of the cotton safety plug is colored.

IDENTIFICATION

The CBS™ 0.3 mL and 0.5 mL straws carry tamper-proof barcode and/or alphanumeric identification on an outer jacket that can be removed after thawing, guaranteeing sample traceability.

FILLING AND SEALING

Straws are filled individually by means of a yellow sterile nozzle and a micro-aspirator.

The impermeability of the seal is only guaranteed if the CBS™ straws are sealed with the SYMS or SYMS II sealing systems.

For large series (sperm banks), the CBS™ straws are filled and sealed on an automatic machine (PACE) or on an automatic filling, sealing and printing machine (MAPI).

STORAGE IN LIQUID NITROGEN

Cryo Bio System carries a line of accessories for storage and management of CBS™ straws in liquid nitrogen.

RECOMMENDATIONS FOR THAWING CBS™ STRAWS

After verifying the location of the desired straw, its identification is confirmed by the color of the internal identification tube (CBS™ 0.3 mL), its plug (CBS™ 0.5 mL) and/or its outer jacket.

The straw is gently removed from the storage visotube, taking care to keep the part containing the sperm in liquid nitrogen. Once it is checked, the straw is reheated by placing it in warm water (maximum 37°C) for a few seconds. After drying, it is ready to use.

SHELF-LIFE

Three years from date of manufacture.



CE marked
FDA 510(k) premarket notification clearance # K002595

010287 CBS™ High Security sperm straw 0.3 mL with white cotton plug*			20
010288			100
CBS™ High Security sperm straw 0.3 mL with white jacket and colored fixed rod			75
017797 red	017798 white	017799 yellow	
017800 blue	017801 green	017802 purple	
014651 CBS™ High Security sperm straw 0.5 mL with white cotton plug*			20
014650			100
CBS™ High Security sperm straw 0.5 mL with colored cotton plug*			100
016584 red	016611 yellow	016612 green	
016613 blue	016614 orange	016615 grey	
008656 Bag of 10 sterile yellow filling nozzles			
Bag of 10 weighted identification rods 30 mm			
019021 yellow	019022 green	019023 blue	
019024 orange	019025 white	019026 red	
Identification jacket			50
010272 red	010273 white	010274 orange	
010275 green	010276 blue	010277 yellow	
014939 clear			
014175 Set of 2 indelible black markers			4
016730 Sterile connection tubing for syringe (individually packaged)			1
016731 Sterile connection tubing for syringe (packaged by 5)			5
010284 Sterile emptying device			20




13



Cryopreservation

CBS™ High Security hydrophobic plugged embryo straw

CE marked
FDA 510(k) premarket notification clearance # K002595



DESCRIPTION / USE

Storage and preservation of embryos.

The special design and the mode of filling and sealing of CBS™ High Security straws for cryo-preservation of human biological samples makes them fully hermetic, in ultra low temperatures.

Manufactured from biocompatible materials, they are used in medically assisted procreation techniques and particularly for embryo preservation. The CBS™ embryo straw is a clear, flexible tube made of ionomeric resin.

Length: 133 mm
Internal diameter: 2.25 mm

Its disposable, clear filling tip allows introduction of the embryo without contaminating the opening.

Two distinct parts separated by a hydrophobic plug guarantee embryo sterility during introduction and recovery after thawing.

The first part holds the embryo and the adjacent dilution solvent (working volume 0.3 mL), while the second part allows insertion of a colored identification rod which also facilitates easy identification in liquid nitrogen.

Sterilized by irradiation.

Individually packaged in blister of 20 units.

IDENTIFICATION

The CBS™ High Security embryo straw carries a tamper-proof colored identification rod to guarantee sample traceability.

FILLING AND SEALING

The CBS™ High Security embryo straw is filled by aspiration using a micro-aspirator according to a specific protocol.

The impermeability of the seal is only guaranteed if the CBS™ straws are sealed with the SYMS or SYMS II sealing systems.

STORAGE IN LIQUID NITROGEN

Cryo Bio System carries a line of accessories for storage and management of CBS™ straws in liquid nitrogen.


RECOMMENDATIONS FOR THAWING CBS™ STRAWS

After verifying the location of the desired straw, its identification is confirmed by the color of the identification rod.

SHELF-LIFE

Three years from date of manufacture.

Weighted identification rod 40 mm						10
010266	yellow	010267	green	010268	blue	
010269	orange	010270	white	010271	red	



CBS™ High Security unplugged tissue 2 mL straws

CE marked

018960 CBS™ High Security tissue 2 mL straw

20



DESCRIPTION / USE

Cryopreservation of tissue fragments.
Sterilized by irradiation.
Individually packaged in blister of 20 units.

IDENTIFICATION

The CBS™ High Security tissue straws can be identified by using LN2 resistant labels.

FILLING AND SEALING

Filling and sealing of the straw should be done in a controlled environment, such as a laminar flow hood.

- Prepare the identification label and apply it on the outside of the CBS™ High Security tissue 2 mL straw;
- Holding the straw near the jaws, seal one extremity of the straw with the sealer SYMS equipped with the special guide for CBS™ High Security tissue 2 mL straws or SYMS II;
- Using a tong, carefully deposit the tissue fragment(s) into the straw, according to laboratory protocol (fragments or wrapped pieces);
- Holding the straw near the jaws, seal the other extremity of the straw with the SYMS or SYMS II sealing systems.

STORAGE IN LIQUID NITROGEN

Cryo Bio System carries a line of accessories for storage and management of CBS™ straws in liquid nitrogen.

SHELF-LIFE

Three years from date of manufacture.

Cryopreservation

PETG cotton-plugged sperm straws

CE marked



010261 Sterile PETG 0.25 mL sperm straw
014550 Sterile PETG 0.5 mL sperm straw

DESCRIPTION / USE

Freezing and storage of human sperm.
PETG sperm straws are clear flexible straws made of a copolymer.

- 0.25 mL straw
Length: 133 mm
Internal diameter: 1.6 mm.
- 0.5 mL straw
Length: 133 mm
Internal diameter: 2.5 mm.

They are composed of two distinct compartments separated by a sliding safety plug.

- The larger compartment is designed to receive sperm and has an effective volume of 0.21 mL or 0.47 mL.
- The smaller compartment is designed as for aspiration at the time of filling.

Sterilized by irradiation.

Packaged in bags of 100 units.

IDENTIFICATION

The straws can be printed by ink jet printer or identified by special LN2 labels or indelible marker for storage in liquid nitrogen.

FILLING AND SEALING

For storage of sperm in preparation of intra-uterine insemination or in vitro fertilization, the PETG straws are filled with a micro-aspirator (or with a syringe and an adapter) and closed with short ID plugs of 13 mm in length.

STORAGE IN LIQUID NITROGEN

Cryo Bio System carries a line of accessories for storage and management of PETG straws in liquid nitrogen.

Shelf-life


Three years from date of manufacture.

007342 Sterile obturation and identification rods (13 mm) yellow

20

PETG cotton-plugged embryo straws

CE marked



DESCRIPTION / USE

Freezing and storage of human embryos.
PETG embryo straws are flexible or rigid clear tubes made of copolymer.

The PETG embryo straws are closed at one end by two woven plugs enclosing a dose of powder that polymerizes when in contact with liquid.

Their useful volume is:

- 0.15 mL for the 91 mm long straws and
- 0.25 mL for the 133 mm long straws.

The internal diameter is 1.6 mm for both models.

Sterilized by irradiation.

Packaged in blisters of 4 x 5 units.

IDENTIFICATION

The straws are identified by different colored plugs.

FILLING AND SEALING

For embryo cryopreservation, the straws are filled with the use of a micro-aspirator.

After the loading procedure, the open extremity of the straw is closed by an obturation and an identification colored plug which is 42 mm in length for 91 mm straws, and 13 mm for 133 mm straws.

Shelf-life

Three years from date of manufacture.

006433 Blister of 20 rigid embryo straws 91 mm

006578 Blister of 20 rigid embryo straws 133 mm

014102 Blister of 20 flexible embryo straws 91 mm

014103 Blister of 20 flexible embryo straws 133 mm

Blister of 20 obturation and identification rods (42 mm)

007437 yellow

007438 purple

007440 red

007442 orange

014484 pink

014485 blue

014486 grey

014487 green

014488 royal blue

014489 white

007342 Sterile obturation and identification rods (13 mm) yellow

20

LabPal

DESCRIPTION / USE

Identification of CBS™ High Security straws and HSV High Security Vitrification kits.

Brady's laboratory labels are designed to withstand extreme laboratory environments.

12-month warranty.

Warning

Using our LN2 resistant labels, it is strongly recommended to overlap the label on the straw.

TECHNICAL SPECIFICATIONS

Print method: Thermal transfer (durable, non-smearing)

Multiple line printing

Weight: 0.55 kg with batteries and tape cartridge included

L x W x H: 228.6 x 114.3 x 63.5 mm

Power source: 6 AA-size batteries or AC with optional adapter

Operating temperature: 4° - 49°C

FCC Class A approved (printer)

019808 LabPal printer



- 019810** AC adaptor 220V for LabPal printer
- 019811** Hard case for LabPal printer
- 019812** Polyester ribbon for LabPal printer LN2 resistant
- 019813** Nylon tissue ribbon for LabPal printer LN2 resistant
- 019859** Clear polyester ribbon for LabPal printer for identification rods for CBS™ High Security straws

Micro-aspirator

014498 Micro-aspirator (delivered with one extra unit of tubing)

014500 Extra nozzle tubing for micro-aspirator

DESCRIPTION / USE

The micro-aspirator is used to manually fill all straws manufactured by Cryo Bio System without any risk of contamination.

- Slowly turn the thumb wheel towards the adaptor until reaching the correct filling height.
- Remove the straw just underneath the adaptor.

No maintenance required.

INSTRUCTIONS FOR USE

- Hold the aspirator in a horizontal position.
- Carefully push the straw into the adaptor and aspirate the liquid. The tube length should be about 10 cm.



CBS™ visotubes and goblets






DESCRIPTION / USE
Freezing and storage of biological materials in straws.

CBS™ High Security straws are placed in goblets for storage in freezers. Goblets come in 8 different colors and their sub-compartments (called visotubes) in 12 different colors.

In the Daisy configuration, goblets hold 12 visotubes of different colors and can hold a maximum of 168 straws of 0.3 ml or 0.5 mL.

Diameter: 65 mm (67 mm with cover)
Goblet height: 117 mm
Total height of Daisy goblet (with cover): 135 mm
Visotubes height: 120 mm
Height of visotube diameter 17.5 mm: 133 mm



		Round visotube ø 17.5 mm	Triangular visotube	Goblet with pre-attached cover	Cover for goblet ø65 mm (height 40 mm) (height 70 mm)		Daisy goblet* (bag of 1) (bag of 5)	
C O L O R S	Clear	006924	-	021224	021366	006321	015144	015152
	Black	-	018153	-	-	-	-	-
	Brown	-	018147	-	-	-	-	-
	Red	-	018146	021228	-	-	015146	015154
	Green	-	018145	021226	-	-	015147	015155
	Blue	-	018148	021225	-	-	015151	015159
	Grey	-	018144	-	-	-	018141	018142
	Purple	-	018152	021229	-	-	015150	015158
	Yellow	-	018143	021227	-	-	015148	015156
	Pink	-	018149	-	-	-	-	-
	Orange	-	018150	-	-	-	-	-
	Pistachio	-	018151	-	-	-	015149	015157
PACKAGING		100	100	10	50	10	1	5
C A P A C I T Y	0.25 mL classical straws	55	45	820	-	-	540	540
	0.5 mL classical straws	20	20	365	-	-	240	240
	0.3 mL & 0.5 mL CBS™ straws	14	14	225	-	-	168	168

Visotubes and goblets




DESCRIPTION / USE

Freezing and storage of biological materials in straws.

Arrangement using color coding of straws prior to storage in liquid nitrogen.

- **Standard goblet Ø 35 mm**
Diameter: 35 mm
Height: 122 mm
- **Standard goblet Ø 65 mm**
Diameter: 65 mm
Height: 117 mm
- **Polygonal, hexagonal and round visotubes**
Height: 120 mm

All visotubes are also available in bulk units of 300 (singular color/diameter); please contact us for further information.



		Polygonal visotube	Hexagonal visotube	Round visotube ø 7.1 mm	Round visotube ø 9.2 mm	Round visotube ø 10 mm	Round visotube ø 12 mm	Round visotube ø 13 mm	Goblet ø 35 mm	Goblet ø 65 mm
C O L O R S	White	020470	020486	020492	020499	020505	020511	020517	020523	006418
	Blue	020472	020487	020493	020500	020506	020512	020518	-	006755
	Green	020482	020488	020494	020501	020507	020513	020519	-	006754
	Yellow	020483	020489	020495	020502	020508	020514	020520	-	006756
	Orange	020484	020490	020496	020503	020509	020515	020521	-	-
	Red	020485	020491	020497	020504	020510	020516	020522	-	006397
	Purple	-	022141	-	-	-	-	-	-	006757
	Pistachio	-	-	-	-	-	-	-	-	006753
PACKAGING		10	10	10	10	10	10	10	10	10
C A P A C I T Y	0.25 mL classical straws	30	90	7	10	14	22	27	200	820
	0.5 mL classical straws	15	36	2	5	5	8	10	94	365
	0.3 mL & 0.5 mL CBS™ straws	8	22	1	3	4	7	7	75	225
	in a goblet ø 35 mm	4	1	15	10	8	5	5	-	-
	in a goblet ø 65 mm	15	7	60	37	31	21	19	-	-

Canisters

Canister with handle and white cane for classical goblets ø 65 mm without cover					10
007147	2-shelf	007148	3-shelf	007149	4-shelf
007150	5-shelf	007151	6-shelf		
Canister with handle and white cane for Daisy goblets					1
008644	2-shelf	008645	3-shelf	008646	4-shelf
007198	5-shelf	007199	6-shelf		
Plastic cane (color flag)					20
007328	white	015693	yellow	015694	salmon
015695	pink	015696	dark green	015697	turquoise
015698	orange	015699	purple	015706	grey
014175	Set of 2 indelible black markers				4
005524	Short tongs (250 mm), stainless steel				1
005523	Long tongs (700 mm), stainless steel				1

DESCRIPTION / USE

Cryo Bio System canisters are the sleeves that hold the goblets. They are compatible with all types of liquid nitrogen and vapor freezers, and also low-temperature mechanical freezers.

Five lengths of canisters are available, with a capacity of 2 to 6 goblets stored on top of each other. A colored flag can be placed on the lifter to facilitate retrieval.

Canisters are placed in the nitrogen freezer before filling the container with liquid nitrogen. Goblets are introduced and retrieved one by one with the metal lifter; all manipulation of goblets is in the cold environment of the neck of the nitrogen freezer to ensure the stable low temperature of the samples.

Diameter: 70 mm.



SIDE

AUTOMATIC PRINTING MACHINE FOR CBS™ HIGH SECURITY STRAWS

CE marked

DESCRIPTION / USE

An automatic system for ink-jet printing of biological samples in CBS™ High Security straws.

SIDE is compatible with CBS™ 0.3 or 0.5 mL open or sealed straws and can automatically print bar codes (code 128) and/or alphanumerical identification on the jacket of the straws.

SIDE is a flexible identification system enabling the laboratory to:

- print the requested number of straws for a sample;
- print the same code on the straws as the sample to ensure the link and
- print straws after sealing to ensure sterility of the content.

SIDE guarantees a one-to-one identification relationship between the straw and the sample code.

The SIDE machine consists of 3 parts:

- a straw handler;
- a PC with pilot software;
- and an ink-jet printing system.

SIDE machine cycle steps:

- Hopper feeding of straws;
- Ink-jet printing of tamper proof bar code and corresponding text of up to 10 digits and control key (code 128);
- Manual transfer to their storage address.

The ink-jet printing system requires active evacuation of its solvent vapors.

Setup, maintenance and shutdown: 15 minutes per day.

SIDE pilot software

The SIDE software features bar code reading or keyboard entry of codes to be printed and export of production files in different formats.

TECHNICAL SPECIFICATIONS

- SIDE
L x W x H: 500 x 450 x 190 mm
Weight: 20 kg
SIDE is provided with a Linx or Domino printer.
Power supply: 110 or 220 V / 60 or 50 Hz

SIDE automatic printing machine for CBS™ High Security straws

- 019651 220 V / with Linx printer
- 019652 110 V / with Linx printer
- 018802 220 V / with Domino printer
- 018803 110 V / with Domino printer





TECHNICAL SPECIFICATIONS


Setup, maintenance and downtime time are approximately 5 minutes per day.

Potentially contaminated parts of the system can be cleaned with ethanol or non-corrosive decontamination fluids.

L x W x H: 363 x 276 x 152 mm

Weight: 6 kg

Power supply: 110 or 220 V / 60 or 50 Hz



23

SYMS

Sealing system for CBS™ High Security straws

CE marked

DESCRIPTION / USE

Sealing system: Cryo Bio System has developed a thermal pulse sealing device for impermeable sealing of CBS™ High Security straws in ionomeric resin.

The sealing device is electronically operated by a push-button or foot pedal as desired.

INSTRUCTIONS FOR USE

- After the filling procedure with microaspirator, gently remove the transparent or yellow filling nozzle;
- Once this free part of the CBS™ straw is filled and sealed, there is no longer any risk of contamination;

- Insert an inner identification rod, or outer colored identification jacket;
- Seal the other end of the straw;
- Freeze according to a method appropriate for the type of product.

TECHNICAL SPECIFICATIONS

L x W x H: 230 x 220 x 170 mm

Weight: 6 kg

Power supply: 110 or 220 V / 60 or 50 Hz



016399 SYMS sealer for CBS™ High Security straws

009848 Spare kit for SYMS sealer (Teflon strips, fuses...)

PACE


SEMI-AUTOMATIC FILLING AND SEALING SYSTEM FOR CBS™ HIGH SECURITY STRAWS

CE marked

PACE semi-automatic filling and sealing system for CBS™ High Security straws

018451 PACE 220 V / 50 Hz

018450 PACE 110 V / 60 Hz



DESCRIPTION / USE

A semi-automatic system for packaging of biological samples in CBS™ High Security straws.

PACE automatically fills and thermally seals both ends of CBS™ 0.3 or 0.5 mL straws from an open primary sample tube.

CBS™ PACE system is the ideal packaging system for sperm banks or medium to large-size collection protocols.

PACE uses CBS™ straws that are pre-identified with bar codes and/or alphanumerical codes or just before printing with SIDE.

Different blood or sperm products from the same individual may be processed in a pre-defined order with the same filling nozzle.

The PACE system consists of 2 parts:

- the straw distribution and filling unit and
- its control box.

PACE machine cycle steps:

- Hopper feeding sample straws individually;
- Manual positioning of primary sample tube and filling tubing;
- Straw filling by aspiration with built-in pumps;
- Controlled thermal autogenic sealing at both ends of straw;
- Manual transfer of straws to their storage-address (visotube and goblet).

The throughput of straws is fully dependent on the protocol (e.g. number of straws filled per primary tube and the number of different biological products).

Setup, maintenance and shutdown: 15 minutes per day.

All potentially contaminated parts can be cleaned with ethanol or non-corrosive decontamination fluids.

The PACE system can be placed on the laboratory bench or under laminar flow with sterile straws and filling tubing and only needs power supply.

TECHNICAL SPECIFICATIONS

• PACE

L x W x H: 380 x 300 x 370 mm
Weight: 25 kg

• Control box

L x W x H: 410 x 365 x 155 mm
Weight: 14 kg

Power supply: 110 or 220 V / 60 or 50 Hz

A PACE can be configured for CBS™ 0.3 or 0.5 mL straws.

007454	Sterile CBS™ injection tubing (by 5 x 10)	50
018620	Sterile long blue CBS™ injection tubing (by 5 x 5)	25
007451	Sterile CBS™ aspiration tubing (by 5 x 10)	50

MAPI

AUTOMATIC FILLING, SEALING AND PRINTING SYSTEM FOR CBS™ HIGH SECURITY STRAWS

CE marked

DESCRIPTION / USE

An automatic system for packaging and identification of biological samples.

From an open primary sample tube, MAPI automatically fills CBS™ 0.3 or 0.5 mL straws, thermally seals both ends, prints bar codes and/or alphanumerical identification on the jacket and in-line validates the process by immediate bar code reading.

The MAPI system consists of the filling robot and its control box, a special ink-jet type printer and a PC with the MAPI pilot software. The automate is placed on a special laboratory bench at ergonomic height.

MAPI machine cycle steps:

- Hopper feeding of straws;
- Bar code reading or key board entry of primary tube code;
- Manual positioning of primary sample tube and filling nozzle;
- Controlled straw filling by aspiration with built-in pumps;
- Controlled thermal autogenic sealing at both ends of straw;
- Ink-jet printing of tamper-proof bar code and corresponding text of up to 10 digits and control key (code 128);
- Reading of straw bar code for validation;
- Selection of validated straws;
- Manual transfer of straws to their storage address (visotube and goblet).

The throughput of straws is fully dependent on the protocol (e.g., number of straws filled per primary tube and number of different biological products).

Setup, maintenance and shutdown: 15 minutes per day.


All potentially contaminated parts can be cleaned with ethanol or non-corrosive decontamination fluids.

The ink-jet printing system requires active evacuation of its solvent vapors.

MAPI pilot software

The MAPI pilot software is installed on the PC provided with the system to pilot the machine cycle of the MAPI system: distribution, filling, sealing, printing and re-reading of straws. The code to be printed on the straw can be either a bar code read from the primary tube, keyboard entered or programmed fixed text.

007218 MAPI automatic filling, sealing and printing system for CBS™ High Security straws



007454	Sterile CBS™ injection tubing (by 5 x 10)	50
018620	Sterile long blue CBS™ injection tubing (by 5 x 5)	25
007451	Sterile CBS™ aspiration tubing (by 5 x 10)	50

TECHNICAL SPECIFICATIONS

• MAPI

L x W x H: 760 x 260 x 470 mm
Weight: 35 kg

• Support table

L x W x H: 1200 x 600 x 900 mm
Height with machine: 1260 mm
Weight with machine: 102 kg

• Control box

L x W x H: 430 x 380 x 190 mm
Weight: 14 kg

• Printer Domino control box

L x W x H: 650 x 450 x 375 mm
Weight: 20 kg

Suction unit for DOMINO connectable to an extractor duct, not included (diameter 80 mm)

Power supply: 110 or 220 V / 60 or 50 Hz



Controlled rate freezers FROM PLANER

016225 Planer Freezer Kryo 360-1.7



016130 Planer Freezer Kryo 360-3.3



018134 Planer Freezer Kryo 560-16



DESCRIPTION / USE

Kryo 360-1.7 / Kryo 360-3.3

Fully featured biological freezer specifically designed for cryopreservation of sperm, oocytes, ovarian tissue and human embryos, in straws or ampoules.

- Chamber volume: 1.7 or 3.3 liters
- CBS™ sperm straw capacity: 45 or 90
- CBS™ embryo straw capacity: 45
- Controlled range: +40°C to -180°C
- Cooling rates: -0.01 to -50°C / min
- Controlled heating rates: 0.01 to 10°C / min
- Horizontal or vertical freezing
- System controller: MRV in 15 positions, radially located
- System pump: LNP4-C
- System Dewar: LAB 20 or 30
- PC software: Delta T

12-month warranty.

INDICATIONS

- Controller displays demand, sample and chamber temperatures, programme stage and current temperature graphic
- Menu driven controller, simple to program and operate
- Horizontal or vertical operation
- Compact design
- Standard operating features:
 - Start above ambient
 - Controlled heating
 - Data Printing (integral printer)
 - Comms port for PC connection
 - Fast cooling rates
 - Multiple safety features

DESCRIPTION / USE

Kryo 560-16

Fully featured biological freezer specifically designed for cryopreservation of bone marrow, stem cells, skin, cord blood, sperm, oocytes, ovarian tissue and human embryos, in straws, ampoules or bags.

- Chamber volume: 16 liters
- CBS™ straw capacity: 456
- Capacity: 11 x 250 / 500 ml bags, vertically in chamber
- Controlled range: +30°C to -180°C
- Cooling rates: -0.01 to -50°C / min
- Controlled heating rates: 0.01 to 10°C / min
- System controller: MRV
- System pump: LNP4-C
- System Dewar: LAB 30
- PC software: Delta T

12-month warranty.

INDICATIONS

- Controller displays demand, sample and chamber temperatures, programme stage and current temperature graphic
- Menu driven controller, simple to program and operate
- Protocol stage trigger on sample or chamber temperature or time
- Unique forced laminar flow system ensures efficient, even cooling
- Standard operating features:
 - Start above ambient
 - Controlled heating
 - Data printing (integral printer)
 - Comms port for PC connection
 - Fast cooling rates
 - Multiple safety features

Coming soon!

CBS™ Co-XStraw™

The new generation of fully colored CBS™ High Security 0.3 & 0.5 mL straws



27

CBS™ High Security tube

A revolutionary High Security 1.2 mL sterile container for long term storage of sperm and ovarian tissue in liquid nitrogen. Safe sealed as CBS™ High Security straws.



**Cryo Bio System France**

Groupe IMV Technologies
ZI n°1 Est
61300 L'Aigle
France
Tel.: +33 (0)233 346 464
Fax: +33 (0)233 341 198

Cryo Bio System Netherlands

Amperestraat 4
8912 AD Leeuwarden
Netherlands
Tel.: +31 58 2123279
Fax: +31 58 2129813

Cryo Bio System Italy

Via Stefano Merli
29122 Piacenza
Italy
Tel.: +39 0523 36 00 10
Fax: +39 0523 38 07 23

Cryo Bio System U.S.A.

11725 95th Avenue North
Maple Grove, MN 55369
USA
Tel.: +1 763 488 1881
Fax : +1 763 488 1888
Toll free within the USA: 1 800 DIAL IMV
(1 800 342 5468)

Cryo Bio System India

Plot No 750, Phase V, Udyog Vihar
Gurgaon 122016 (Haryana)
India
Tel.: +91 124 400 44 30/31/32
Fax: +91 124 400 21 72

Cryo Bio System China

Room 2402, Beijing Anzhen Plaza
N°2 Andingmenwai Street
Dongcheng District, Beijing
China
Tel.: +86 10 64481509
Fax: +86 10 64282447

www.cryobiosystem.com

e-mail: contact@cryobiosystem-imv.com