

We send you attached video and will explain how it works.

Video Link: [http://v.youku.com/v\\_show/id\\_XNDk1MzI4MjA4.html](http://v.youku.com/v_show/id_XNDk1MzI4MjA4.html)

Frist, Turn on the machine and select Lauage option: English, Chinese and etc.

Second, enter in the main interface like this:



the M1 M2 M3 is to adjust the flash times

M1 is 2 times flash/minutes , M2 is 3 times/minutes M3 is 5 times/minutes

when we adjust the button in the "PULSE" then choose the 3 modes

when we adjust it to the continuous we do not need to adjust the M1 M2 M3

"ENERGY" when we click this, then choose the number then choose "+ or -"

so we can change the eneregy for every light.

The number 1 2 3 4 5 6 1 2 3 4. it controlled all the 12 lipo light.

And how to put all the laser pads on body?

Pls use the bandage on your body and insert the pads in the bandage.

Like this picture:



We suggest that 20-25 mins/Each time of treatment and 2-3 times/Week.

1.What is the pulse for?

In the picture named screen.jpg,it mentioned that the pulse means that the light will be flashed and will not work continuously.

And when you choose the pulse button,then you will have three options for mode.

And the M1 M2 M3 is to adjust the flash times.

M1 is 2 seconds flash, M2 is 4seconds flash, M3 is 6 seconds flash.

2.When do I use pulse?

The choose pulse or continous mode,it might be determined by your reference.

Some people need stronger treatment,then he will use continous mode.

And sometimes,the exchange between continous and pulse will have better results.

3.When do I turn up the energy?

This might be determined by different person.Some skin can stand higher heat.Some are not.

So just listen to your client about their confortableness.

When he say that he need higher intensity for better results,you can try to increase it slowly from low to high grandually.

# Lipo Laser Fat System

## User manual



### Warning

Dear user :

Thanks for using our company's products, in order to fully make use of the products, we sincerely suggests follow:

- ◆◆ Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆ Please follow the instruction to fit and operate the machine correctly.
- ◆◆ Please do not remove and change any accessories of the machine.
- ◆◆ It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆ It forbids others to change and maintain the machine except the professional people.
- ◆◆ Please contact us if there is any trouble during operation. This product does not belong to medical, please don't use for disease treatment !

### Work Therapy

LLLT(Low Level Laser Therapy) has been used for many years, low level laser emits low level laser energy lipo laser treatment is completely noninvasive, safe, have no side effect for fat loss

It penetrate into skin surface and stimulates the fat cell membranes, changing their permeability, the fat cells reduce their overall size and intracellular fat is released then the fatty trigly cerides flow out of the disrupted cell membranes and into the interstitial space, where they gradually pass through the body\'s natural metabolic functions with no harmful physiological effects, this process is resulting in inch loss of patients

## Features

- 1 Pain free, no surgery and have no side effect
- 2 Complete safely and efficiently
- 3 Non-invasive, fat reduction and body shaping
- 4 Cold low laser



## Treatment procedure

- 1 Cleanse the skin, Put the Lipolaser probe to body where you want lost fat. then use a belt keep the pad to body.



- 2 You relax in a comfortable position while a trained Lipo Laser technician places two multi-laser paddles and two smaller lasers to the target areas. Then lay back, relax, read, or whatever you like.

## Comparison Before and after



## FAQ

1. Is LipoLaser Safe?

-Yes, Lipolaser treatments are completely non-invasive and clients feel nothing. Lipo Laser, low level laser technology has been studied and used for over 30 years for many biomedical purposes. CE approval for pain and inflammation, doctors can use it for fat loss .

2. Do I have to do other things to lose the fat?

-Only for bigger results - 1 inches without changes, 5 inches with major lifestyle changes, like diet and exercise we have support if you choose to change your life while doing the sessions

3. What do you do to keep the fat off?

-Normal things you need to keep weight off diet and exercise, drinking water, etc.

4. Where Does the Fat Go?

-That fat can be burned as energy in the body, or normal metabolic processing just like when you exercise and release fat

5. How Quickly Will I See Results

-Many clients see measurable results after their first session.

6. What Kind of Results are Produced?

-1-5 inches average-the most lost in 9 sessions is 23 inches.People who make minor lifestyle changes, are on the higher end of the average, and those who dont, the lower end.

7.How many times do I need it ?

-3times per week, 9times per course of treatment.

8. Forbidden group

-1.People with heart disease or high blood pressure, or configured cardiac pacemaker.

2.Patients with acute inflammation, asthma, deep vein thrombosis, thyroncus, cancer.

3.People with hemorrhagic disease, trauma or who is bleeding.

4.Pregnant women

5.Medical Plastic parts, or parts with artificial in filling inside.

6.People with metal inside the body(include metallic tooth)

7.Patient with skin inflammation or with edema.

8.People with an abnormal immune system.

9.Numb or insensitive to heat.

## **Technical Parameters**

Laser wave length: Laser diode AlGalnp 635-650nm

Maxinum output power:52W

intensity:3.0W/CM2(Belly),3.6W/CM2(Thigh),2.4W/CM2(Arm),1.5W/CM2(Face)

N.W:8Kg

G.W:12Kg

Dimension of the machine:37x31x21cm

Dimension with packing:57x51x36cm

## **Packing Include**

6xBig paddle with 8 led diode laser//Each (Total 48 Diode laser)

4xSmall paddle with 1 led diode laser//Each (Total 4 Diode laser)

1xMetal holder for diode paddles

1xPower line(We offer 100-240V,AU/EU/UK/AU plug)