

# ADC Read Software User Manual


UTA-HEP/LC 0024

*Shashwat Udit*  
*University of Texas at Arlington*  
*August 25, 2008*

## 1. Logging on to HEPPC 24

1-1. Log on to HEPPC 24 using user name and password on side of the computer.


1-2. Alternatively, log on to HEPPC24 using remote desktop feature from your own computer



**Figure 1: Remote Desktop Logon**

- a. Open up remote desktop connection, as shown by Figure 1, on your computer, under programs> accessories> communications. This is for Windows.
- b. Enter in the IP Number 129.107.77.54
- c. Log on using user name and password of HEP PC 24

## 2. Open up LabVIEW



**Figure 2: Desktop with Icons**


2-1 To read only one ADC Click on the Desktop Icon titled One ADC Read.

2-2 To read data from all 4 ADCs go to the Desktop Icon entitled All ADC Read.

2-3 To make changes to the LabVIEW Code or. Write new programs, go to Start> Programs> National Instruments>LabVIEW 8.0s


## 3. Set User Inputs

3-1 Set user inputs for One ADC Read.



**Figure 3: Front Panel showing user Inputs**

- The data collected will be stored in a text file. Write the file path in the input titled File to Write at the Top.
- Select the ADC Number. Enter 1 for the ADC in slot 17, 2 for the one in slot 18, 3 for the one in Slot 19, or 4 for the ADC in Slot 20. If any other number is entered, or nothing at all, the ADC in Slot 17 will be read.
- The last input is Number of Events (x16) The number typed in, multiplied by 16, is the number of events that are stored in the array at a time before they are outputted. For example, if 3 is typed in, when 48 events are stored in the array, the program will write them out to the graphs and the text file. If one is typed in, the program will only wait for 16 events to be stored before it outputs them. To see why this is, look in the attached program overview. Outputting data will causes the program to be a bit slower.



**Figure 4: Front Panel with Run Buttons Circled**

- d. To run the program, there are a couple buttons at the top of the display, as shown in Figure 4. The four to the left are all that are normally needed. The straight, right facing arrow is the run button. The program will execute once through, that is, it will execute until it outputs data one time. The two arrows forming a loop to the right of that button are the run continuously button. The program will run until the user hits either the stop or the pause button, which are to the right of the run continuously button, and stop and pause the program respectively.

### 3-2 Set user inputs for All ADC Read.

- a. The data collected will be stored in a text file. Write the file path in the input titled File to Write at the Top.
- b. The last input is Number of Events (x16) The number typed in, multiplied by 16, is the number of events that are stored in the array at a time before they are outputted. For example, if 3 is typed in, when 48 events are stored in the array, the program will write them out to the graphs and the text file. If one is typed in, the program will only wait for 16 events to be stored before it outputs them. Outputting data will cause the program to be a bit slower.
- c. To run the program, there are a couple buttons at the top of the display, as shown in Figure 4. The four to the left are all that are normally needed. The straight, right facing arrow is the run button. The program will execute once through, that is, it will execute until it outputs data one time. The two arrows forming a loop to the right of that button are the run continuously button. The program will run until the user hits either the stop or the pause button, which are to the right of the run continuously button, and stop and pause the program respectively.

#### 4. Data Output

File	Edit	Format	View	Help					
2859	334	364	344	385	380	398	378		
2741	333	364	345	386	379	398	378		
2770	334	365	345	386	382	401	378		
2806	334	365	344	385	380	398	379		
2766	332	363	344	384	380	400	376		
2755	332	365	344	385	380	398	378		
2767	334	365	344	386	379	398	378		
2775	332	364	344	385	380	398	378		
2757	333	364	346	385	379	398	379		
2730	333	364	345	385	381	400	378		
2770	334	364	345	384	381	400	378		
2741	334	364	344	384	380	398	379		
2770	333	364	344	386	380	399	379		
2773	333	364	343	384	380	398	379		
2737	334	364	344	386	381	400	380		
2714	332	363	344	385	380	400	379		

**Figure 5: Text File**

##### 4-1. Graph Output

- a. Each channel outputs to a waveform chart, with event number being the X-AXIS value, and ADC Counts represented on the Y-Axis.

##### 4-2 Text File Output

- a. The values are written to a text file, with each event being a row, and each channel a column, as shown in Figure 5. As more and more events are written, more rows are added to the text file.
- b. Text can be appended to an existing text file. To create a text file, open up My Computer, select the folder in which you wish for it to be stored, and add new text file.