


# URS-BPPV6 Series

## *Vacuum Compatible Precision Rotation Stages*



**Newport®**  
Experience | Solutions

## USER'S MANUAL

*For Motion, Think Newport*

# Warranty

Newport Corporation warrants this product to be free from defects in material and workmanship for a period of 1 year from the date of shipment. If found to be defective during the warranty period, the product will either be repaired or replaced at Newport's discretion.

To exercise this warranty, write or call your local Newport representative, or contact Newport headquarters in Irvine, California. You will be given prompt assistance and return instructions. Send the instrument, transportation prepaid, to the indicated service facility. Repairs will be made and the instrument returned, transportation prepaid. Repaired products are warranted for the balance of the original warranty period, or at least 90 days.

## Limitation of Warranty

This warranty does not apply to defects resulting from modification or misuse of any product or part.

---

### CAUTION

**Warranty does not apply to damages resulting from:**

- **Incorrect usage:**
    - **Load on the rotation stage greater than maximum specified load.**
    - **Carriage speed higher than specified speed.**
    - **Improper grounding.**
      - **Connectors must be properly secured.**
      - **When the load on the rotation stage represents an electrical risk, it must be connected to ground.**
    - **Excessive or improper cantilever loads.**
  - **Modification of the rotation stage or any part thereof.**
- 



### CAUTION

**Please return equipment in the original (or equivalent) packing.**

**You will be responsible for damage incurred from inadequate packaging if the original packaging is not used.**

This warranty is in lieu of all other warranties, expressed or implied, including any implied warranty of merchantability or fitness for a particular use. Newport Corporation shall not be liable for any indirect, special, or consequential damages.

No part of this manual may be reproduced or copied without the prior written approval of Newport Corporation.

This manual has been provided for information only and product specifications are subject to change without notice. Any changes will be reflected in future printings.

# Table of Contents

Warranty .....	ii
EU Declaration of Conformity.....	v
Definitions and Symbols.....	vi
Warnings and Cautions .....	vi
Warnings .....	vii
Cautions .....	viii
<b>1.0 — Introduction .....</b>	<b>1</b>
<b>2.0 — Description .....</b>	<b>2</b>
2.1 Design Details .....	2
<b>3.0 — Characteristics .....</b>	<b>3</b>
3.1 Definitions .....	3
3.2 Mechanical Specifications.....	4
3.3 Load Specification Definitions .....	4
3.4 Load Characteristics and Stiffness .....	4
3.5 Stage Weights .....	5
<b>4.0 — Drive .....</b>	<b>5</b>
4.1 Specifications.....	5
<b>5.0 — Motor .....</b>	<b>5</b>
5.1 Motor Characteristics.....	5
5.2 Command Signals.....	5
5.3 URS50BPPV6 Sensor Position .....	6
5.4 Pinouts.....	6
<b>6.0 — Connection to Newport Controllers .....</b>	<b>8</b>
6.1 Warnings on Controllers .....	8
6.2 Connection.....	9
6.3 Cables .....	9
<b>7.0 — Connection to Non-Newport Controllers .....</b>	<b>10</b>
<b>8.0 — Disabling of Limit Switches.....</b>	<b>11</b>
<b>9.0 — Dimensions .....</b>	<b>11</b>
9.1 URS50BPPV6 Rotation Stage.....	11
9.2 URS75BPPV6 Rotation Stage.....	12

---

<b>10.0 — Maintenance .....</b>	<b>13</b>
10.1 Maintenance.....	13
10.2 Repairing .....	13
10.3 Calibration.....	13

---

Service Form .....	15
--------------------	----

# EU Declaration of Conformity

## URS-BPPV6 Series

### EC Declaration of Conformity

following Annex II-1A  
of Directive 2006/42/EC on machinery



**The manufacturer:**

MICRO-CONTROLE Spectra-Physics,  
1 rue Jules Guesde ZI. Bois de l'Épine - BP189  
F-91006 Evry FRANCE

**Hereby declares that the machinery:**

Description: " URS-BPPV6 "  
Function: Precision rotation stages  
Models: URS50BPPV6 / URS75BPPV6.

- the technical file of which was compiled by:  
Mr Dominique DEVIDAL, Quality Director,  
MICRO-CONTROLE Spectra-Physics, Zone Industrielle - B.P.29  
F-45340 Beaune La Rolande France
- complies with all the relevant provisions of the Directive 2006/42/EC on machinery.
- complies with all the relevant provisions of the Directive 2004/108/EC relating to electro-magnetic compatibility.
- was designed and built in accordance with the following harmonised standards:  
NF EN 61326-1:2006 « Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements »  
NF EN 55011:2007 Class A
- was designed and built in accordance with the following other standards:  
NF EN 61000-4-2  
NF EN 61000-4-3  
NF EN 61000-4-4  
NF EN 61000-4-5  
NF EN 61000-4-6

### ORIGINAL DECLARATION

Done in Beaune La Rolande on 09 June 2011  
Dominique DEVIDAL  
Quality Director

DC1-EN rev:A

# Definitions and Symbols

The following terms and symbols are used in this documentation and also appear on the product where safety-related issues occur.

**General Warning or Caution**




The exclamation symbol may appear in warning and caution tables in this document. This symbol designates an area where personal injury or damage to the equipment is possible.

**European Union CE Mark**



The presence of the CE Mark on Newport Corporation equipment means that it has been designed, tested and certified as complying with all applicable European Union (CE) regulations and recommendations.




ATTENTION

**This stage is a Class A device. In a residential environment, this device can cause radioelectric interferences. In this case, suitable measurements must be taken by the user of this device.**


**Warnings and Cautions**

The following are definitions of the Warnings, Cautions and Notes that may be used in this manual to call attention to important information regarding personal safety, safety and preservation of the equipment, or important tips.



WARNING

**Situation has the potential to cause bodily harm or death.**



CAUTION

**Situation has the potential to cause damage to property or equipment.**

NOTE

**Additional information the user or operator should consider.**

# Warnings



---

**WARNING**

The rotation of objects of all types carries potential risks for operators. Ensure the protection of operators by prohibiting access to the dangerous area and by informing the personnel of the potential risks involved.

---

---

**WARNING**

Do not use this stage when its motor is emitting smoke or is unusually hot to the touch or is emitting any unusual odor or noise or is in any other abnormal state.

Stop using the stage immediately, switch off the motor power and then disconnect the electronics power supply.

After checking that smoke is no longer being emitted contact your Newport service facility and request repairs. Never attempt to repair the stage yourself as this can be dangerous.

---

**WARNING**

Make sure that this stage is not exposed to moisture and that liquid does not get into the stage.

Nevertheless, if any liquid has entered the stage, switch off the motor power and then disconnect the electronics from power supply.

Contact your Newport service facility and request repairs.

---



---

**WARNING**

Do not insert or drop objects into this stage, this may cause an electric shock, or lock the drive.

Do not use this stage if any foreign objects have entered the stage. Switch off the motor power and then disconnect the electronics power supply.

Contact your Newport service facility for repairs.

---

**WARNING**

Do not place this stage in unstable locations such as on a wobbly table or sloping surface, where it may fall or tip over and cause injury.

If this stage has been dropped or the case has been damaged, switch off the motor power and then disconnect the electronics power supply.

Contact your Newport service facility and request repairs.

---

**WARNING**

Do not attempt to modify this stage; this may cause an electric shock or downgrade its performance.

---

**WARNING**

Do not exceed the usable depth indicated on the mounting holes (see section “Dimensions”). Longer screws can damage the mechanics or cause a short-circuit.

---



# Cautions

---

**CAUTION**

Do not place this stage in a hostile environment such as X-Rays, hard UV or in a vacuum environment less than  $10^{-6}$  hPa.

---

**CAUTION**

Do not place this stage in a location affected by dust, oil fumes, steam or high humidity. This may cause an electric shock.

---

**CAUTION**

Do not leave this stage in places subject to extremely high temperatures or low temperatures. This may cause an electric shock.

- Operating temperature: +10 to +35 °C.
  - Storage temperature: -10 to +40 °C (in its original packaging).
- 

**CAUTION**

Do not move this stage if its motor power is on.



Make sure that the cable to the electronics is disconnected before moving the stage. Failure to do so may damage the cable and cause an electrical shock.

---

**CAUTION**

Be careful that the stage is not bumped when it is being carried. This may cause it to malfunction.

---

**CAUTION**

When handling this stage, always unplug the equipment from the power source for safety.

---

**CAUTION**

When the carriage is in its end-of-run position, it is strongly recommended not to go beyond this point by using the manual knob as this may damage the stage mechanism.

---

**CAUTION**

Contact your Newport service facility to request cleaning and specification control every year.


---




# Vacuum Compatible Precision Rotation Stages URS-BPPV6 Series

## 1.0 Introduction

This manual provides operating instructions for the URS-BPPV6 rotation stage that you have purchased.



*URS50BPPV6 rotation stage.*



*URS75BPPV6 rotation stage.*

### RECOMMENDATION

We recommend you read carefully the chapter “Connection to electronics” before using the URS-BPPV6 rotation stage.

## 2.0 Description

The URS-BPPV6 rotation stages features an enhanced, lower profile package providing easier top down mounting that is directly compatible with our vacuum compatible stages.

When using with our XPS or ESP301 motion controllers with high micro-step capability, low noise operation and very small incremental motions are guaranteed. The URS-BPPV6 rotation stage does not use encoder feedback, but calculates position by the number of commanded steps and micro-steps. For this purpose, the stepper motor is directly attached to the worm screw with a proprietary bellow coupling with high torsional stiffness and no gear or belt drive in between. Furthermore, the high output torque of the stepper motor minimizes the risk of lost steps and provides optimum motion sensitivity with good linearity between commanded micro-step and the actual motion of the stage.

The URS-BPPV6 rotation stages features a proprietary 4-contact point ball bearing. The unique, 2-piece design of the URS-BPPV6 provides a low-profile compact stage with exceptional stiffness, high reliability and outstanding wobble and eccentricity. The tilted worm screw arrangement allows for 4 symmetric mounting holes as compared to other designs that only feature 2 or 3 mounting holes. As a result, the URS-BPPV6 stages provides a higher stiffness for better support of higher or off-centered loads.

Additionally, the flexible preloading system for the worm gear has been improved to guarantee a backlash-free operation with an MTBF of 20,000 hours.

### 2.1 Design Details

Base Material	Hardened steel with aluminum body
Bearings	Large diameter ball bearings
Drive Mechanism	Grinded worm gear with self-compensating preload
Worm Gear Ratio	URS50BPPV6: 1:80 URS75BPPV6: 1:90
Feedback	<b>URS50BPPV6: Motor Mounted Rotary Encoder, 48 pts/rev.</b> URS75BPPV6: None
Limit Switches	URS50BPPV6: None URS75BPPV6: Two independently adjustable optical limit switches
Origin	Optical, fixed at position 0°. Typical 0.04° repeatability
Manual Adjustment	Via Allen wrench (included) at the end of the worm screw
Motor	URS50BPPV6: UE16PPHSV6 URS75BPPV6: UE41PPV6
Cable Length	Sub-D15F/Sub-D25M cable, length 3 m Not Vacuum compatible
MTBF	20000 h at 25% load and with a 30% duty cycle



#### NOTE

**This product complies with the RoHS directive  
(Restriction of Hazardous Substances)**

## 3.0 Characteristics

### 3.1 Definitions

Specifications of our products are established in reference to ISO 230 standard part II “Determination of the position, precision and repeatability of the machine tools with CNC”.

This standard gives the definition of position uncertainty which depends on the 3 following quantities:

#### **(Absolute) Accuracy**

Difference between ideal position and real position.

#### **On-Axis Accuracy**

Difference between ideal position and real position after the compensation of linear error sources.

Linear errors include: cosine errors, inaccuracy of screw or linear scale pitch, angular deviation at the measuring point (Abbe error) and thermal expansion effects. All Newport motion electronics can compensate for linear errors.

The relation between absolute accuracy and on-axis accuracy is as follow:

$$\text{Absolute Accuracy} = \text{On-Axis Accuracy} + \text{Correction Factor} \times \text{Travel}$$

#### **Repeatability**

Ability of a system to achieve a commanded position over many attempts. Can be specified as uni-directional or bi-directional.

#### **Reversal Value (Hysteresis)**

Difference between actual position values obtained for a given target position when approached from opposite directions.

#### **Minimum Incremental Motion (Sensitivity)**

The smallest increment of motion a device is capable of delivering consistently and reliably.

#### **Resolution**

The smallest increment that a motion device can be commanded to move and/or detect.

#### **Eccentricity**

Displacement of the geometric center of a rotation stage from the rotation axis in the plane defined by bearings.

#### **Wobble**

Tilt of rotation axis during rotation of a stage.

The testing of on-axis accuracy, repeatability, and reversal error are made systematically with our test equipment in an air-conditioned room (20 °C ±1 °C).

The test is comprised of 4 cycles, each taking 21 data points in the forward and reverse directions, for a total of 164 data points.

#### **Guaranteed Specifications**

Guaranteed maximum performance values are tested per the Newport's A167 metrology test procedure. For more information, please consult the metrology tutorial section in the Newport catalog or on the Newport web-site at [www.newport.com](http://www.newport.com)

### 3.2 Mechanical Specifications



	URS50BPPV6	URS75BPPV6
Travel Range (°)	360 continuous	360 continuous <sup>(1)</sup>
Resolution (°)	0.0014 <sup>(2)</sup>	0.0002 <sup>(2)</sup>
Minimum Incremental Motion (°)	0.0015	0.0002
	(0.0005 in open loop mode)	
Uni-directional Repeatability, guaranteed (°)	0.004 or ±0.002	0.004 or ±0.002
Bi-directional Repeatability, guaranteed (°)	0.006 or ±0.003	0.014 or ±0.007
Absolute Accuracy (°)	0.04 or ±0.02	0.03 or ±0.01
Max. Speed (°/s)	4 <sup>(3)</sup>	20 <sup>(3)</sup>
Wobble, guaranteed (μrad)	50 or ±25	
Eccentricity, guaranteed (μm)	6 or ±3	3 or ±1.5

<sup>(1)</sup> With disabled limit switches.

<sup>(2)</sup> Equal to 1/2 of a full step for URS50BPPV6, and 1/100 of a full step for URS75BPPV6.

<sup>(3)</sup> URS50BPPV6: 1.6 °/s; URS75BPP: 8 °/s with an SMC100PP controller.



#### CAUTION

To reach specifications stated, stages must be fixed on a plane surface with a flatness of 5 μm.

### 3.3 Load Specification Definitions

#### Normal Load Capacity (Cz)

Maximum load a rotation stage can move while maintaining specifications.

This value is given with speed and acceleration specified for each rotation stage, and with a load perpendicular to bearings.


#### Off-Centered Load (Q)

Maximum cantilever-load a rotation stage can move:  $Q \leq Cz / (1 + D/a)$

D: Cantilever distance.

a: Construction parameter.

### 3.4 Load Characteristics and Stiffness



	URS50BPPV6	URS75BPPV6
Cz (N)	50	100
a (mm)	20	25
$k\alpha$ (μrad/Nm)	100	30
Mz (Nm)	±0.25	±0.5

$$Q, \text{ Off-center load} \quad Q \leq Cz / (1 + D/a)$$

Cz, Normal center load capacity on bearings

D, Cantilever distance in millimeters

a, Construction parameter

$k\alpha$ , Transversal stiffness

Mz, Nominal Torque

3.5 Stage Weights

Weights indicated into the below table are values for stages without any cable.

	Weight [lb (kg)]
URS50BPPV6	1.5 (0.7)
URS75BPPV6	3.7 (1.7)
3-meter Cable	0.66 (0.3)

4.0 Drive

4.1 Specifications

	Full-Step Resolution <sup>(1)</sup> (°)	Speed (°/s)	Motor
URS50BPPV6	0.0028	4	UEPPHSV6
URS75BPPV6	0.02	20 <sup>(2)</sup>	UE41PPV6

<sup>1)</sup> When used with Newport motion controllers, this motor is driven in a dynamic micro-stepping mode (software commutation).


<sup>2)</sup> 1.6 °/s when used with a SMC100PP controller for URS50BPPV6, 4 °/s for URS75BPPV6.

5.0 Motor

5.1 Motor Characteristics


Motor	Angle by Step (°)	Current (A)	Resistance (Ω)	Inductance (mH)	Newport Utilization
UE16PPHSV6	15	0.18	12.5	5.5	Mini-Step (bipolar)
UE41PPV6	1.8	0.85	5.7	11.3	Micro-Step (bipolar)

5.2 Command Signals



End-of-Run and Mechanical Zero are TTL type: 5 V ±5%, 16 mA max.

5.3 URS50BPPV6 Sensor Position



CAUTION


“End-of-Run” and “Mechanical Zero” are active signals and should not be connected to any other source. Use appropriate TTL type receivers.

5.4 Pinouts

The 15-pin Sub-D connection for URS-BPPV6 rotation stages is given in the following table:

URS50BPPV6 UE16PPHSV6		URS75BPPV6 UE41PPV6	
1	+ Phase 1	1	+ Phase 1
2	+ Phase 2	2	+ Phase 2
3	Mechanical Zero	3	Mechanical Zero
4	N.C.	4	– End-of-Run
5	Ground	5	Ground
6	Encoder Phase /A	6	N.C.
7	Encoder Phase /B	7	N.C.
8	N.C.	8	N.C.
9	– Phase 1	9	– Phase 1
10	– Phase 2	10	– Phase 2
11	N.C.	11	+ End-of-Run
12	+5 V	12	+5 V
13	Encoder Phase A	13	N.C.
14	Encoder Phase B	14	N.C.
15	N.C.	15	N.C.

The 25-pin Sub-D connection at the end of the cable is given in the following table:



URS50BPPV6 UE16PPHSV6		URS75BPPV6 UE41PPV6	
1	+ Phase 1	1	+ Phase 1
2	+ Phase 1	2	+ Phase 1
3	– Phase 1	3	– Phase 1
4	– Phase 1	4	– Phase 1
5	+ Phase 2	5	+ Phase 2
6	+ Phase 2	6	+ Phase 2
7	– Phase 2	7	– Phase 2
8	– Phase 2	8	– Phase 2
9	N.C.	9	N.C.
10	N.C.	10	N.C.
11	N.C.	11	N.C.
12	N.C.	12	N.C.
13	Mechanical Zero	13	Mechanical Zero
14	Ground	14	Ground
15	N.C.	15	N.C.
16	Ground	16	Ground
17	N.C.	17	+ End-of-Run
18	N.C.	18	– End-of-Run
19	Encoder Phase A	19	N.C.
20	Encoder Phase B	20	N.C.
21	+5 V	21	+5 V
22	Ground	22	Ground
23	Encoder Phase /A	23	N.C.
24	Encoder Phase /B	24	N.C.
25	N.C.	25	N.C.



## 6.0 Connection to Newport Controllers

---

### NOTE

Visit [www.newport.com](http://www.newport.com) for compatible Newport controllers.

---

#### 6.1 Warnings on Controllers

Controllers are intended for use by qualified personnel who recognize shock hazards and are familiar with safety precautions required to avoid possible injury. Read the controller user's manual carefully before operating the instrument and pay attention to all written warnings and cautions.

---

### WARNING

Disconnect the power plug under the following circumstances:

- If the power cord or any attached cables are frayed or damaged in any way.
  - If the power plug is damaged in any way.
  - If the unit is exposed to rain, excessive moisture, or liquids are spilled on the unit.
  - If the unit has been dropped or the case is damaged.
  - If you suspect service or repair is required.
  - Whenever you clean the electronics unit.
- 

### CAUTION

To protect the unit from damage, be sure to:

- Keep all air vents free of dirt and dust.
  - Keep all liquids away from the unit.
  - Do not expose the unit to excessive moisture (>85% humidity).
  - Read this manual before using the unit for the first time.
- 



### WARNING

All attachment plug receptacles in the vicinity of this unit are to be of the grounding type and properly polarized.

Contact your electrician to check your receptacles.

---

### WARNING

This product is equipped with a 3-wire grounding type plug.

Any interruption of the grounding connection can create an electric shock hazard.

If you are unable to insert the plug into your wall plug receptacle, contact your electrician to perform the necessary alterations to ensure that the green (green-yellow) wire is attached to earth ground.

---

### WARNING

This product operates with voltages that can be lethal.

Pushing objects of any kind into cabinet slots or holes, or spilling any liquid on the product, may touch hazardous voltage points or short out parts.

---

## 6.2 Connection

On each packing of stage is represented a label which indicates its name and its serial number.



---

### WARNING

**Always turn the controller's power OFF before connecting to a stage.**

---

Stages may be connected to the rear panel motor connectors any time prior to power-up with the supplied cable assemblies.

---

### NOTE



URS-BPPV6 are ESP compatible stages. Enhanced System Performance is Newport's exclusive technology that enables Newport ESP motion controllers to recognize the connected Newport ESP stage and upload the stage parameters. This ensures that the user can operate the motion system quickly and safely.

---

## 6.3 Cables

Our URS-BPPV6 rotation stages are delivered equipped with a **MSCABLE-3** standard Sub-D15F/Sub-D25M 3-meter cable. This cable is not vacuum compatible.

---

### WARNING

The MSCABLE-3 cable supplied is not designed for using in a vacuum environment. The customer has the responsibility to link the rotation stage to the bulkhead coupling with a vacuum compatible cable and connect the cable supplied between the controller and the bulkhead coupling (the Sub-D15 connector must be removed).



---

### WARNING

This cable is shielded correctly. For a correct operation, make sure to lock connectors (ground continuity provided by cables).

---

---

### WARNING

Keep the motor cables at a safe distance from other electrical cables in your environment to avoid potential cross talk.

---

7.0

Connection to Non-Newport Controllers



WARNING

Newport takes no responsibility for improper functioning or damage of an URS-BPPV6 rotation stages when used with any non-Newport electronics.

WARNING

Newport guarantees the “CE” compliance of the URS-BPPV6 rotation stages only if are used with Newport cables and electronics.

Nevertheless, the figure below indicates the recommended wiring when an URS-BPPV6 rotation stages are used with non-Newport controllers.



“Encoder” and “Index Pulse” are “differential pair” type output signals. Using these signals permits a high immunity to noise. Emission circuits generally used by Newport are 26LS31 or MC3487. Reception circuits to use are 26LS32 or MC3486.


8.0

Disabling of Limit Switches

NOTE

The URS50BPPV6 stage is not equipped with limit switches.

The URS75BPPV6 stage is equipped with  $\pm 160^\circ$  optical limit switches. This limit switch can be disabled for continuous  $360^\circ$  rotation.




NOTE

On request, URS75BPPV6 limit switches can be set in factory at the desired position. Please contact our sales engineers.


9.0

Dimensions

9.1 URS50BPPV6 Rotation Stage



9.2 URS75BPPV6 Rotation Stage



## 10.0 Maintenance

---

### RECOMMENDATION

It is recommended to contact Newport Service who can guide you with the appropriate maintenance for your application.

---

### 10.1 Maintenance

The URS-BPPV6 stages require no particular maintenance. Nevertheless, this is a precision mechanical device that must be kept and manipulated with precaution.

---

### PRECAUTIONS

The URS-BPPV6 stages must operate, and be stocked in a clean environment, without dust, humidity, solvents or other substances.

---

### RECOMMENDATION

It is recommended to return your stage to Newport Service after every 2000 hours of use for lubrication.

If your URS-BPPV6 stage is mounted on a workstation and cannot be easily dismantled, please contact Newport Service for further instructions.

---

### 10.2 Repairing



### CAUTION

Never attempt to disassemble a component of the stage that has not been specified in this manual.

Disassembly of a non-specified component can cause a malfunction of the stage.

---

If you observe a malfunction in your stage, please immediately contact us to make arrangements for a repair.

---



### CAUTION

All disassembly attempts or repair of stage without authorization will void your warranty.

---

### 10.3 Calibration



### CAUTION

It is recommended to return your stage to Newport once a year for a recalibration to its original specifications.

---





# Service Form

## Your Local Representative

Tel.: \_\_\_\_\_

Fax: \_\_\_\_\_

Name: \_\_\_\_\_

Return authorization #: \_\_\_\_\_

(Please obtain prior to return of item)

Company: \_\_\_\_\_

Address: \_\_\_\_\_

Date: \_\_\_\_\_

Country: \_\_\_\_\_

Phone Number: \_\_\_\_\_

P.O. Number:

P.O. Number: \_\_\_\_\_ Fax Number: \_\_\_\_\_

**Item(s) Being Returned:**

Model #: \_\_\_\_\_

Serial #: \_\_\_\_\_

Description: \_\_\_\_\_

Reasons of return of goods (please list any specific problems): \_\_\_\_\_



# Newport®

Experience | Solutions

**Visit Newport Online at:  
[www.newport.com](http://www.newport.com)**

## **North America & Asia**

Newport Corporation  
1791 Deere Ave.  
Irvine, CA 92606, USA

### **Sales**

Tel.: (800) 222-6440  
e-mail: [sales@newport.com](mailto:sales@newport.com)

### **Technical Support**

Tel.: (800) 222-6440  
e-mail: [tech@newport.com](mailto:tech@newport.com)

### **Service, RMAs & Returns**

Tel.: (800) 222-6440  
e-mail: [rma.service@newport.com](mailto:rma.service@newport.com)

## **Europe**

MICRO-CONTROLE Spectra-Physics S.A.S

1, rue Jules Guesde – Bât. B  
ZI Bois de l'Épine – BP189  
91006 Evry Cedex  
France

### **Sales & Technical Support**

Tel.: +33 (0)1.60.91.68.68  
e-mail: [france@newport-fr.com](mailto:france@newport-fr.com)

### **Service & Returns**

Tel.: +33 (0)2.38.40.51.55