

**CIRCA Scientific Temperature Monitoring System**

**Model CS-1000**



**CIRCA Scientific, LLC**

**14 Inverness Drive East, Suite H240**

**Englewood, CO 80112**

**(303) 951-8767**

## *Table of contents:*

---

- 1. General warnings and cautions ... Page 3*
- 2. Indications for use ... Page 4*
- 3. Description ... Page 4*
- 4. Alarm System Overview ... Page 7*
- 5. Setup instructions ... Page 8*
- 6. Check-out instructions ... Page 9*
- 7. Operating instructions ... Page 10*
- 8. Shut-Down ... Page 12*
- 9. Cleaning and Disinfection ... Page 14*
- 10. Maintenance ... Page 14*
- 11. Storage and Transportation ... Page 15*
- 12. Troubleshooting ... Page 16*
- 13. Accessories ... Page 18*
- 14. Technical Information ... Page 19*
- 15. Symbols Key ... Page 21*
- 16. Electromagnetic Compatibility ... Page 22*

# CIRCA Scientific Temperature Monitor User Manual

---

## *1. General warnings and cautions:*

---

- Rx Only: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- The monitor is designed for use with CIRCA Scientific interconnect cables and temperature probes only. Incompatible components can result in degraded performance and could lead to damage to the unit.
- No modification of this equipment is allowed.
- Never pour any liquid into an opening. This may cause fire or electrical shock.
- Do not cover the openings. Overheating may occur as the openings on the enclosure are for air convection.
- If the equipment is not in use, disconnect it from the power source to avoid damage by transient overvoltage.
- Do not remove cover. No user serviceable parts inside. There is a danger of shock if incorrectly serviced. There is a danger of explosion if battery is incorrectly replaced. Refer servicing to qualified personnel.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.
- To avoid product damage, do not repair or replace power supply and USB Cable. Refer servicing to qualified personnel.
- This equipment needs special precautions regarding EMC (Electromagnetic Compatibility) and needs to be put into service according to the EMC information provided in section 16 “Electromagnetic Compatibility”.
- Do not connect additional accessory equipment to the analog and digital interfaces for signal input or output. Personnel who connect additional equipment configure a medical system; degraded performance and damage to the unit could occur.
- Part of defibrillation proof protection is provided by the S-Cath™ temperature probe. Do not use with any other applied part.

**Caution** is a statement that alerts the user to the possibility of a problem with the device associated with its use or misuse. Such problems include device malfunction, device failure, damage to the device or damage to other property.

**Warning** is a statement that alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

# CIRCA Scientific Temperature Monitor User Manual

---

## 2. Indications for use:


---

Display continuous temperature measurement ( $^{\circ}\text{C}$ ) from 12-sensor temperature probe.

## 3. Description:

---


The CIRCA Scientific Temperature Monitoring System consists of a touch-screen monitor, bracket, and power supply. The monitor is to be used with CIRCA Scientific interconnect cable and temperature probe only. The interconnect cable and temperature probe are sold separately.



[Figure 1: System Equipment]

# CIRCA Scientific Temperature Monitor User Manual

The monitor displays 12 temperature probe sensor readings (°C), the maximum temperature of all sensors, and contains an alarm system with user-selected levels. The monitor features the following three screens:




[Figure 2: Display Screens]

# CIRCA Scientific Temperature Monitor User Manual

---

Note: numbered channels 1 through 12 correspond to each sensor located within the S-Curve of the temperature probe. Channel number begins with sensor #1 located on distal end.



[Figure 3: Temperature Probe Sensor Location]

# CIRCA Scientific Temperature Monitor User Manual

## 4. Alarm System Overview

The alarm system provides visual and audible feedback when a sensor temperature is equal to and above the user-selected value.

The alarm system is intended to provide operator feedback regarding temperature compared to user-selected levels only; it does not provide physiological alarm conditions.


The visual and audible signal is triggered when a sensor temperature is equal to and above the set value. Alarm levels are set by the operator (see Operating Instructions section below for setting instructions).

WARNING is a lower priority than ALARM.


WARNING is identified by YELLOW.

ALARM is identified by RED.

AUDIO OFF inactivates the audible signal, only until another event occurs. Another event includes maximum temperature increasing from warning to alarm level, or temperature drops below warning or alarm value, then increases and is equal and above set value.



[Figure 4: Warning Signal]




[Figure 5: Alarm Signal]

# CIRCA Scientific Temperature Monitor User Manual

## 5. Setup instructions:

The operator is responsible for checking the compatibility of the monitor, interconnect cable, and temperature probe before use. Ensure only CIRCA Scientific components and equipment is connected.

- S1) Mount monitor to standard IV pole or similar rigid structure. Secure by tightening knob located at the bottom of the bracket.
- S2) Adjust monitor to desired position (tilt position and angle) using lock levers located on bracket. (*Note: intended position of the Operator to observe the Alarm Signal is within an approximate distance of 4.5 feet (1.4 meters).*)
- S3) Connect interconnect cable to monitor (connection located on back of monitor) by aligning snap-fit connectors and pushing firmly.
- S4) Connect temperature probe to interconnect cable by aligning snap-fit connectors and pushing firmly.
- S5) Plug power cord into power outlet. **Warning:** To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.



[Figure 6: Setup Connections]



# CIRCA Scientific Temperature Monitor User Manual




---

## 6. Check-out Instructions

---

- C1) Turn monitor on by momentarily pushing Power Button located on back of monitor.
- C2) During start-up, monitor will deliver a test audible alarm signal (3 beeps). If no audible signal is heard, sound system is defective. Refer service of monitor to qualified personnel.
- C3) Verify temperatures are displayed on monitor. If no temperature displays, verify connections are fully seated and resolve any error messages displayed on monitor.

The following messages will display if cables are not connected or unit does not initialize:

Message	Resolution
	<p>Touch “OK” button.</p> <p>Ensure USB Cable (back of monitor) is firmly seated in ports. Reconnect cable.</p> <p>Proceed with Shut Down from Windows screen.</p> <p>Turn monitor back on by momentarily pushing Power Button on back of monitor.</p>
	<p>Ensure temperature probe is firmly seated to interconnect cable and interconnect cable is firmly seated to monitor.</p> <p>When properly connected, message will terminate and temperature readings display.</p>
	<p>Touch “OK” button.</p> <p>Ensure USB Cable (back of monitor) is firmly seated in ports. Reconnect cable.</p> <p>Proceed with Shut Down from Windows screen.</p> <p>Turn monitor back on by momentarily pushing Power Button on back of monitor.</p>


[Table 1: Message Overview]

- C4) Proceed with operating instructions below.

# CIRCA Scientific Temperature Monitor User Manual


## 7. Operating instructions

Screen navigation and set-up is by touch screen. Different screens may be selected by touching the tab near the top of the screen.



[Figure 7: Display Screens]

01) Select “Set Up” tab to set desired alarm system levels and graph y-axis minimum and maximum.




[Figure 8: Set Up Screen]

- a. Settings – touch “+” button to increase value, touch “-“ to decrease value. (Note: warning temperature cannot be set equal to or greater than alarm temperature. Graph maximum and minimum cannot be set below alarm and warning levels.)

# CIRCA Scientific Temperature Monitor User Manual

---

- b. Alarm and Warning Audible Levels – touch “ALARM” or “WARNING” button to hear and update audible tone and volume settings. (Note: warning volume cannot be set equal to or greater than alarm volume.) Touch bar graph level to desired setting and then “UPDATE” button to save and exit.



[Figure 9: Tone and Volume Setting]

Note: all user selected values are retained; the last settings used before shut down will be those recalled upon start-up.


- O2) Once set up is complete, select desired screen, either “Maximum Temperature” or “Channel Graphs” in which to view temperatures measured by the temperature probe.
- O3) Disconnect probe from interconnect cable by grasping connectors. Do not pull on cable or probe wire to disconnect. (Note: after disconnecting, the message “No valid temperature data. Please check probe.” will display.)
- O4) Prepare temperature probe for patient use per temperature probe’s instructions for use.
- O5) Once temperature probe is set in patient and connected to interconnect cable, temperatures will display.

## 8. Shut-Down

Turn monitor off by following 3 steps:


### Step 1:

Go to “Set Up” Screen and select “Stop” Button:



After Selecting “Stop”, Program will exit to Windows Screen.


### Step 2: Select “Start” – “Shutdown”:



# CIRCA Scientific Temperature Monitor User Manual

---

Step 3: Select “OK” and unit will shut down:



## *9. Cleaning and Disinfection*

---

Disconnect from power before cleaning.

Use a damp cloth. Do not use liquid or spray detergents for cleaning.

Wipe outer surfaces with damp cloth and let dry. Do not rinse, soak, wash or sterilize.

If disinfection is required, apply non-abrasive/non-corrosive disinfection fluid to disposable cloth, wipe outer surfaces, and then let dry.

## *10. Maintenance*

---

### Routine Visual Inspection

Visually inspect the monitor and all accessories at least once before each use. Inspect the power cord, cables, and monitor for damage, wear, and loose components. Particular attention should be made to the power cord and interconnect cable for insulation damage such as cuts, brittleness, cracking, and bare spots. Do not use if equipment appears damaged.

### Annual Safety Inspection and Accuracy Test

The following safety inspection and accuracy test should be performed by hospital's equipment service department at least every 12 months:

- Inspect the equipment and accessories for mechanical damage.
- Inspect all labels and markings for legibility.
- Check temperature accuracy with the following test:

Equipment: 10K Resistor Test Part (available from CIRCA Scientific)

1. Connect Test Part to Interconnect Cable.
2. Turn unit on.
3. Observe temperature readings (12) displayed on the monitor.
4. Calculate the measurement error for each individual Output Temperature measurement by subtracting reading from 25.0°C.
5. Ensure that the measurement error is not greater than 0.3°C.

Do not use equipment if the inspections or test reveal a defect. The equipment has no serviceable parts. **Warning:** to prevent shock, do not remove cover. Refer servicing to qualified personnel.

## *11. Storage & Transport*




---

Do not leave this equipment in an environment where the storage temperature may go below -20°C (-4°F) or above 60°C (140°F). This could damage the equipment.

If shipping equipment, pack in original carton and packing materials. If original packing material is not available, cover monitor, pack with foam, and ship in sturdy box to prevent damage during transport.



# CIRCA Scientific Temperature Monitor User Manual

## 12. Troubleshooting

Problem	Resolution
<p>Message</p> 	<p>Touch “OK” button.</p> <p>Ensure USB Cable (back of monitor) is firmly seated in ports. Reconnect cable.</p> <p>Proceed with Shut Down from Windows screen.</p> <p>Turn monitor back on by momentarily pushing Power Button on back of monitor.</p>
<p>Message</p> 	<p>Ensure temperature probe is firmly seated to interconnect cable and interconnect cable is firmly seated to monitor.</p> <p>Replace probe and/or interconnect cable if no temperature data can be displayed after verifying connections.</p>
<p>Message</p> 	<p>Touch “OK” button.</p> <p>Ensure USB Cable (back of monitor) is firmly seated in ports. Reconnect cable.</p> <p>Proceed with Shut Down from Windows screen.</p> <p>Turn monitor back on by momentarily pushing Power Button on back of monitor.</p>
<p>Monitor does not boot up and BIOS Configuration reset to default.</p>	<p>The computer inside the monitor is provided with a battery-powered real-time clock circuit. The battery has no power.</p> <p><b>Warning:</b> There is a danger of explosion if battery is incorrectly replaced. Refer servicing to qualified personnel.</p>



# CIRCA Scientific Temperature Monitor User Manual

Problem	Resolution
<p>“----“ is displayed for an individual sensor temperature reading and an individual sensor graph (under “Channel Graphs” display) is orange.</p> 	<p>No action required.</p> <p>The four dashes (‘----’) and orange-colored graph indicate an individual sensor wire has failed.</p>
<p>Frozen Screen: Monitoring indicator does not flash.</p> 	<p>Shut down unit and re-start.</p>

[Table 2: Troubleshooting Overview]

# CIRCA Scientific Temperature Monitor User Manual

## 13. Accessories


The monitor is designed for use with CIRCA Scientific interconnect cables and temperature probe only. The interconnect cables and temperature probes are sold separately.

Contact CIRCA Scientific or Authorized Distributor for ordering information.

The “Recorder” Port located on the back of the monitor is intended to connect with the 0 to 5V analog input of a recording device that is compliant with IEC 60601-1. Do not connect this port to any other device.

**Caution:** Use of accessories or replacement parts not supplied by CIRCA Scientific can result in degraded performance and could lead to damage to the unit.

**Warning:** Use of cables and accessories other than those supplied by CIRCA Scientific may result in increased emissions or decreased immunity of the equipment.






[Figure 10: Equipment Illustration]

# CIRCA Scientific Temperature Monitor User Manual

---

## 14. Technical Information

---

Classification	Class I Defibrillation-Proof Type BF Applied Part Continuous Operation
Software	Revision level 1.1
Electrical (Mains)	Mains Supply Voltage: 100-240V AC  Mains Supply Frequency: 47-63 Hz Mains Rated Input: 1.35A
Electrical (Power Supply)	Output Voltage: 12V DC  Max Current: 3.75A Max Use only Sinpro model MPU50-105 power supply.
Electrical (Monitor Power Input)	Input Voltage: 12-19 V DC  Current: 3.75-2.63A
Electrical (Safety and Electromagnetic Compatibility)	Meets IEC 60601-1:1988 +A1:1991 + A2:1995 Meets IEC 60601-1-2:2007
User Settings	Alarm and Warning Temperature Levels in 0.1°C increments Alarm and Warning Tone and Volume 1 – 10 in single digit increments Graph y-axis Minimum and Maximum in 1°C increments Setting range = 0°C to 90°C
Measurement Display	Update rate = 50 milliseconds Graph time span = 60 seconds Accuracy = ± 0.3°C Precision = 0.1°C
Alarm System	Intended position of the Operator to observe the Alarm Signal is within an approximate distance of 4.5 feet (1.4 meters). Alarm Signal Sound Pressure Range = 45 to 85 dB

## CIRCA Scientific Temperature Monitor User Manual

---
















---

Physical	Dimensions (monitor): 10.2" W x 7.8" H x 3.25" D 258 W x 199 H x 83 D (mm) Weight: 5.5 lbs. (2.6 kg)
Disposal	No special precautions are required. Dispose of equipment per hospital policy. EU Only: Products affected by the directive Waste of Electrical and Electronic Equipment (WEEE). These products are not to be discarded together with non-electrical or non-electronic products.

[Table 3: Technical Information]

# CIRCA Scientific Temperature Monitor User Manual

## 15. Symbols Key (Equipment Markings)

	Alternating current		Dangerous voltage
	Direct current		Caution
	“ON” / “OFF” (push-push) NOTE Each position, "ON" or "OFF", is a stable position.		EU Only: Products affected by the directive Waste of Electrical and Electronic Equipment (WEEE). These products are not to be discarded together with non-electrical or non-electronic products.
	Centre Positive. Indicates that the center (tip) of the output plug is Positive (+) and the barrel of the output plug is Negative (-).		For indoor use only
	Defibrillation-Proof Type BF Applied Part		Manufacturer
	Consult Instructions for Use		Date of Manufacture
Rx Only	Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician		Serial Number
	Caution: part of defibrillation proof protection is provided by the S-Cath™ temperature probe. Do not use with any other applied part.		Catalogue Number

[Table 4: Symbols Key]

# CIRCA Scientific Temperature Monitor User Manual

---

## 16. Electromagnetic Compatibility

---

- **Warning:** Use of cables and accessories other than those supplied and sold by CIRCA Scientific may result in increased emissions or decreased immunity of the equipment.
- **Warning:** The equipment should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the monitor should be observed to verify normal operation. Normal operation is considered as absence of unusual, erratic variations in temperature readings.
- The equipment may be affected by portable and mobile RF (Radio Frequency) communications equipment.

**Table 5 – Guidance and manufacturer’s declaration – Electromagnetic Emissions**

The CIRCA Scientific CS-1000 Temperature Monitoring System is intended for use in the electromagnetic environment specified below.

The customer or the user of the CIRCA Scientific CS-1000 Temperature Monitoring System should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The CIRCA Scientific CS-1000 Temperature Monitoring System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The CIRCA Scientific CS-1000 Temperature Monitoring System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

# CIRCA Scientific Temperature Monitor User Manual

**Table 6 – Guidance and manufacturer’s declaration – Electromagnetic Immunity**

The CIRCA Scientific CS-1000 Temperature Monitoring System is intended for use in the electromagnetic environment specified below.

The customer or the user of the CIRCA Scientific CS-1000 Temperature Monitoring System should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0,5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CIRCA Scientific CS-1000 Temperature Monitoring System requires continued operation during power mains interruptions, it is recommended that the CIRCA Scientific CS-1000 Temperature Monitoring System be powered from an uninterruptible power supply or a battery.

# CIRCA Scientific Temperature Monitor User Manual

---

---

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: $U_T$ is the a.c. mains voltage prior to application of the test level.			



# CIRCA Scientific Temperature Monitor User Manual

**Table 7 – Guidance and manufacturer’s declaration – Electromagnetic Immunity**


The CIRCA Scientific CS-1000 Temperature Monitoring System is intended for use in the electromagnetic environment specified below.

The customer or the user of the CIRCA Scientific CS-1000 Temperature Monitoring System should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the CIRCA Scientific CS-1000 Temperature Monitoring System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ <p>80 to 800 MHz</p> $d = 2.3 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF</p>

# CIRCA Scientific Temperature Monitor User Manual

---

		<p>transmitters, as determined by an electromagnetic site survey<sup>(a)</sup> should be less than the compliance level in each frequency range<sup>(b)</sup>.</p> <p>Interference may occur in the vicinity of equipment that include RF transmitters or that apply RF electromagnetic energy for diagnosis, i.e. equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
--	--	---

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>(a)</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CIRCA Scientific CS-1000 Temperature Monitoring System is used exceeds the applicable RF compliance level above, the CIRCA Scientific CS-1000 Temperature Monitoring System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CIRCA Scientific CS-1000 Temperature Monitoring System.

<sup>(b)</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Note: If high frequency surgical equipment is used at the same time and interferes with the operation of the CIRCA Scientific CS-1000 Temperature Monitoring System, additional measures may be necessary, such as re-orientation of cables, relocation, and/or connecting the hospital grade power cable into a different grounded receptacle or separate grounded power source.

## CIRCA Scientific Temperature Monitor User Manual

**Table 8 – Recommended separation distances between portable and mobile RF communications equipment and the CIRCA Scientific CS-1000 Temperature Monitoring System**

The CIRCA Scientific CS-1000 Temperature Monitoring System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CIRCA Scientific CS-1000 Temperature Monitoring System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CIRCA Scientific CS-1000 Temperature Monitoring System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter <b>W</b>	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12 meters (4.7 inches)	0.12 meters (4.7 inches)	0.23 meters (9.1 inches)
0.1	0.38 meters (15.0 inches)	0.38 meters (15.0 inches)	0.73 meters (28.7 inches)
1	1.2 meters (3.9 feet)	1.2 meters (3.9 feet)	2.3 meters (7.6 feet)
10	3.8 meters (12.5 feet)	3.8 meters (12.5 feet)	7.3 meters (24.0 feet)
100	12 meters (39.4 feet)	12 meters (39.4 feet)	23 meters (75.5 feet)

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.