PA001GECA00 MS9662 Pen Needle HumaPen LUXURA® Pen Parts (sold separately) MS9663 Cartridge Lilly 3 mL Insulin Dose Paper Holder Cartridge Knob Tab **INSULIN DELIVERY DEVICE** For use only with Lilly insulin (Humalog® or Humulin®) 3 mL Cartridges (100 units/mL). Your HumaPen LUXURA® may be a different colour from the one shown, but it will work in the same way. Rubber Cartridge Injection Outer www.lilly.ca

HumaPen LUXURA is designed to enable you to reliably deliver an accurate dose of your insulin in an easy to use manner. You can inject from 1 to 60 units of any Humalog® or Humulin® insulin from a 3 mL cartridge (100 units/mL). You can dial your dose one unit at a time. If you dial too many units, you can correct the dose without wasting any insulin.

Read this manual completely and follow the directions carefully before using the pen. Failure to follow the instructions may result in too much or too little insulin being delivered.

DO NOT SHARE YOUR PEN OR NEEDLES AS THIS MAY RISK TRANSMISSION OF INFECTIOUS AGENTS.

If any of the parts of your HumaPen LUXURA appear broken or damaged, DO NOT USE. Contact your healthcare professional or call Lilly at 1-888-545-5972 for a replacement pen. HumaPen LUXURA is not recommended for the blind or visually impaired without the assistance of a sighted individual trained to use it. Always carry a spare source of insulin that you are trained on how to inject and could immediately use in case your pen is lost or damaged.

INSERT THE INSULIN CARTRIDGE

Important Notes

- Your healthcare professional has prescribed the type of insulin best for you.
- Any changes in insulin should be made only under medical supervision. HumaPen LUXURA is for use only with Lilly insulin (Humalog or Humulin) 3 mL cartridges.
- Do not use other brands of insulin cartridges.
- Read and follow the instructions provided in your Humalog or Humulin 3 mL insulin cartridge Information for the Patient Leaflet.
- Before each injection, read the cartridge label and be sure the pen contains the correct Humalog or Humulin 3 mL insulin cartridge.
- The colour of the pen is not intended to indicate insulin type.
- For more information on HumaPen LUXURA and Humalog or Humulin insulin, please contact your healthcare professional or Lilly at 1-888-545-5972.

Frequently Asked Questions about Inserting the Insulin Cartridge

1. Why doesn't the screw move out when there is no cartridge in the pen?

The screw may not move out when you push the injection button unless there is a cartridge in the pen. This feature of the pen allows you to easily push the screw into the pen when replacing a cartridge. Once a cartridge is installed, the screw will move out when the injection button is pushed.

2. What should I do if I can't attach the cartridge holder to the pen body?

Check that the insulin cartridge is fully inserted into the cartridge holder. Then carefully line up the cartridge holder with the pen body and screw together until

INJECT THE DOSE

Important Notes

- You must PUSH the injection button straight down for the dose to be delivered. You will NOT receive your insulin by turning the dose knob.
- Do not attempt to change the dose while injecting.

Frequently Asked Questions about Injecting

- 1. Why is it difficult to push the injection button when I try to inject?
- Your needle may be clogged. Try attaching a new needle.
- Pushing the injection button down quickly may make the button harder to push. Pushing the button more slowly may make it easier.
- Using a larger diameter needle will make it easier to push the injection button during injection. See your healthcare professional to determine which needle is best for you.
- Your injection button may become harder to push if the inside of your pen gets dirty with insulin, food, drink or other materials.

If none of the above steps resolves the problem, your pen may need to be replaced.

2. Why doesn't the dose knob go to zero when I inject my dose?

This can happen if the insulin cartridge does not have enough insulin left in it for your entire dose. The number in the dose window is the amount you did NOT receive. To get the rest of your dose, remove the needle, change the insulin cartridge, attach a new needle and prime the pen. Then complete your dose by dialing the amount that you did **NOT receive** and inject only this amount.

3. Why is insulin leaking from the needle after I finish my dose?

You may have removed the needle from your skin too quickly; remember to push and hold the injection button and slowly count to 5 before removing the needle. Make sure you see a 0 in the dose window to confirm you received the complete dose.

CARE AND STORAGE

- Keep the pen away from moisture, dust, extreme hot or cold temperatures, and direct sunlight.
- Do not store the pen in a refrigerator. • Wipe the pen cap, pen body and case with a damp cloth to clean them. Do not use
- alcohol on the pen body or dose window. • Do not use alcohol, hydrogen peroxide, bleach, cover in liquid or apply lubrication such
- as oil, as this could damage the pen.
- Refer to the insulin cartridge Information for the Patient Leaflet for complete insulin storage instructions.

USE AS DIRECTED

HumaPen LUXURA should only be used as described in this manual. The manufacturer will not be responsible for any problems with the device if these instructions are not followed. Failure to follow the instructions may result in breakage of the pen and too much or too little insulin being delivered.

PRIME THE PEN (BEFORE EVERY INJECTION)

Plunger

Body

Important Notes

- HumaPen LUXURA is suitable for use with BD pen needles.
- The directions regarding needle handling are not intended to replace local, healthcare professional or institutional policies.
- Do not store your pen with the needle attached, use a new needle for each injection. This will help ensure sterility. It will also help prevent leakage of insulin, keep out air bubbles, and reduce needle clogs.

Frequently Asked Questions about Priming

1. Why is it important to prime before every injection?

If you do NOT prime, you may get too much or too little insulin. Priming helps to ensure that the pen and needle are working properly. Once the pen is properly primed, a stream of insulin will flow from the needle. You may need to prime several times before you see a stream of insulin.

2. Why can it take several attempts to prime when a new cartridge is inserted?

There may be a gap between the screw and the cartridge plunger. Repeating the priming steps will move the screw out to touch the cartridge plunger. Once the end of the screw pushes the cartridge plunger out, insulin will flow from the needle.



- Repeat the priming steps until a stream of insulin is seen.
- If you are still unable to see a stream of insulin flow from the needle refer to Question 3.

3. Why should I prime to a stream, and why isn't the pen ready for use if I see a drop of insulin on the needle tip?

Priming moves the screw into contact with the cartridge plunger and gets the air out of the cartridge. When you are priming the pen:

- You may see a drop of insulin on the tip of the needle when you first attach it. This only shows that the needle is attached and not clogged. You must still prime the
- You may also see no flow at all. This may be because the screw is moving forward to close a gap between the screw and the cartridge plunger.
- If there is air in the cartridge, the insulin may sputter or drip until all the air is
- removed. - Insulin will flow as a stream only when the pen is properly primed.
- If the injection button is hard to push, the needle may be clogged. Attach a new needle. Repeat the priming steps until a stream of insulin is seen.

If you are still unable to see a stream of insulin out of the needle, do NOT use the pen. Contact your healthcare professional or Lilly at 1-888-545-5972 for assistance or to obtain a replacement.

4. What should I do if I have an air bubble in the cartridge?


Priming your pen will remove air bubbles. Point the pen up, and tap the cartridge gently with your finger so any air bubbles can collect near the top. Repeat the priming steps until a stream of insulin is seen. A small air bubble may remain in the cartridge after completion of the priming steps. If you have properly primed the pen, this small air bubble will not affect your insulin dose.

REPLACEMENT

Your HumaPen LUXURA has been designed to be used for up to 6 years after first use. Record the date your pen was first used here: __/__/_. Go directly to a pharmacy to get a new HumaPen LUXURA when your pen has been used for 6 years.

If you have any questions or concerns with your HumaPen LUXURA, contact your healthcare professional or Lilly at 1-888-545-5972 for assistance.





Distributed by: Eli Lilly Canada Inc. 3650 Danforth Avenue Toronto, Ontario M1N 2E8

Manufactured by: Eli Lilly and Company Pharmaceutical Delivery Systems Lilly Corporate Center Indianapolis, IN 46285 USA

Authorized Representative in the EU: Eli Lilly and Company Ltd. Erl Wood Manor, Windlesham, Surrey

GU20 6PH UK

HumaPen LUXURA meets the current dose accuracy and functional requirements of ISO 11608-1:2000 with Lilly 3 mL cartridges. Humalog, Humulin, and HumaPen LUXURA are registered trademarks owned by Eli Lilly and Company; used under license. Copyright © 2011 Eli Lilly and Company. All rights reserved. Document Revision Date: 26 July 2011

 ϵ 0344

PA001GECA00

