

STSTM

SUPER TURBINE SUCTION Dual Stage Dry Vacuum System with the Cyclonic Action Separator (CASTM)

STS-15TM with 20 GALLON CASTM
P/N 54130

STS-30TM with 20 GALLON CASTM
P/N 54132



USER'S MANUAL



**AIR
TECHNIQUES**
INC.

TABLE OF CONTENTS

<u>Description</u>	<u>Page</u>
Congratulations	2
Warranty	3
On-Line Warranty Registration	3
Safety Instructions	3
Sizing Guide	4
Key Parts Identification	4
Specifications	5
Site Requirements	5
Product Specifications	5
Dimensions	5
STS System Configurations	6
Operating Information	7
Pipe Flow Schematic	7
Installation Information	8
Installation Procedures	9
Single STS (STS-15) PVC Hose Installation	9
Dual STS (STS-30) Hose Installation	10
Electrical Connections	12
Maintenance	13
Trouble Shooting	16
Installation Accessory Kits	17
Accessories/Options	18

CONGRATULATIONS

Congratulations on the purchase of your new **STSTTM** Dual Stage Dry Vacuum System (Super Turbine Suction). The **STSTTM** is a dry vacuum pump that produces high-volume air flow with multiple users online. The **STSTTM** is a medical dry vacuum pump which is designed for use in a dental facility. The **CASTM** (Cyclonic Action Separator) tank will ensure that no liquids or foams enter the vacuum pump. The relief valve is easily accessible on the **CASTM**. The patented vacuum relief valve maintains a constant uniform vacuum. The powerful 3 phase motor, with a highly reliable contactor and powerful transformer can be depended upon to start every time. Your **STSTTM** vacuum system and **CASTM** separator tank are easily installed and maintained.

This manual provides operation, installation, and maintenance instructions for the support of the **STSTTM** Dual Stage Dry Vacuum System. Review and follow the guidelines included in this User Manual to ensure that the system provides the highest level of service.

The **STS™** and **CAS™** are warranted to be free from defects in material and workmanship from the date of installation for a period of thirty-six (36) months.

Any item returned to our factory through an Air Techniques Authorized Dealer, will be repaired or replaced at our option at no charge provided that our inspection shall indicate it to have been defective. Dealer labor, shipping and handling charges are not covered by this warranty.

This warranty does not apply to damage due to shipping, misuse, careless handling or repairs by other than authorized service personnel. Warranty is void if equipment is installed or serviced by other than dealer service personnel authorized by Air Techniques. Air Techniques, Inc. is not liable for indirect or consequential damages or loss of any nature in connection with this equipment.

This warranty is in lieu of all other warranties expressed or implied. No representative or person is authorized to assume for us any liability in connection with the sale of our equipment.

ON-LINE WARRANTY REGISTRATION

Quickly and easily register your new **STS™** on-line. Just have your product model and serial numbers available. Then go to the Air Techniques website, **www.airtechniques.com**, click the **warranty link** and complete the registration form. This on-line registration ensures a record for the warranty period and helps Air Techniques keep you informed of product updates and other valuable information.

SAFETY INSTRUCTIONS

Use of the **STS™** not in conformance with the instructions specified in this manual may result in permanent failure of the unit.

WARNING: To prevent fire or electrical shock, do not expose this appliance to rain in or moisture.

All user serviceable items are described in the maintenance section.

Manufacturing date code on serial number label is in the format Month YYYY.

ATTENTION USERS:



Alerts users to important Operating and Maintenance instructions. Read carefully to avoid any problems.



Warns users that uninsulated voltage within the unit may be of sufficient magnitude to cause electric shock.



Indicates the ON and OFF position for the Equipment power switch.



Indicates protective Earth Ground for the Equipment power switch.



Indicates type B equipment in accordance with IEC 601-1



Warns users of hot surfaces. There is a danger of burns. Work near these surfaces only after they have cooled down.



MEDICAL ELECTRICAL EQUIPMENT

WITH RESPECT TO ELECTRICAL SHOCK, FIRE, MECHANICAL
AND OTHER SPECIFIED HAZARDS ONLY
IN ACCORDANCE WITH UL-60601-1, CAN/CSA C22.2 NO.601.1
66CA

SIZING GUIDE

Choosing the correct size **STS™** to meet practice depends on the number of air users and the anticipated air demand. To assure optimum operation, the demands should not exceed the number of air handpiece users shown below. Each chart lists the number of simultaneous High Volume Evacuators (HVEs) and Saliva Ejectors (SEs) that can be used in respective systems.

STS-15 SYSTEM (with 1 20-gallon CAS)					
HVE's	+	SE's	HVE's	+	SE's
15	+	0	7	+	16
14	+	2	6	+	18
13	+	4	5	+	20
12	+	6	4	+	22
11	+	8	3	+	24
10	+	10	2	+	26
9	+	12	1	+	28
8	+	14	0	+	30

STS-30 SYSTEM (2 STS-15's in Tandem with 1 20-gallon CAS)											
HVE's	+	SE's	HVE's	+	SE's	HVE's	+	SE's	HVE's	+	SE's
30	+	0	22	+	16	14	+	32	6	+	48
29	+	2	21	+	18	13	+	34	5	+	50
28	+	4	20	+	20	12	+	36	4	+	52
27	+	6	19	+	22	11	+	38	3	+	54
26	+	8	18	+	24	10	+	40	2	+	56
25	+	10	17	+	26	9	+	42	1	+	58
24	+	12	16	+	28	8	+	44	0	+	60
23	+	14	15	+	30	7	+	46			

KEY PARTS IDENTIFICATION




Figure 1. STS-15 and STS-30 with 20-Gallon CAS Main Parts Location

SPECIFICATIONS

SITE REQUIREMENTS

ELECTRICAL	SINGLE (STS-15)	DUAL (STS-30)
Minimum Circuit Breaker Rating	30A	30A (Qty 2)
Wire Size AWG Minimum Gauge	#10	#10AWG (Qty 2)
Receptacle	Hardwired	Hardwired
PLUMBING		
Exhaust Vent Pipe	3-Inch Metal Pipe	3-Inch Metal Pipe (Qty 2)
Suction Line		
Riser Diameter ID	1/2 Inch	1/2 Inch
Branch Line Diameter ID Minimum	1-1/2 to 2 Inches	2 to 3 Inches
Main Line Diameter ID	2 to 3 Inches	3 to 4 Inches
End Fitting at STS	1-1/2 Inch FNPT	1-1/2 Inch FNPT
Drain Line	1-1/2 Inch Schedule 40 Pipe	1-1/2 Inch Schedule 40 Pipe
Wash-Out Line	3/4-Inch Garden Hose	3/4-Inch Garden Hose
Ambient Temperature	40 to 104°F (5 to 40°C)	40 to 104°F (5 to 40°C)

PRODUCT SPECIFICATIONS


ELECTRICAL	SINGLE (STS-3)	DUAL (STS-6)
Voltage (Minimum/Maximum)	200/250	200/250
Full Load	15 Amps	30 Amps
Starting Load	100 Amps	200 Amps
Frequency	60 Hz	60 Hz
Maximum Vacuum	17 InHg	17 InHg
Preset Vacuum Level	10 InHg	10 InHg
CAS TM SPECIFICATIONS		
Working Liquid Capacity:	20 Gallons plus adequate capacity for foam	20 Gallons plus adequate capacity for foam
Tank Material	304 Stainless Steel	304 Stainless Steel
ENVIRONMENTAL CONDITIONS		
Operating Temperature	40 to 104°F (5 to 40°C) with PVC vent pipe	40 to 104°F (5 to 40°C) with PVC vent pipe
Storage Temperature	0 to 150°F (-18 to 66°C)	0 to 150°F (-18 to 66°C)
Relative Humidity	90% (no condensation)	90% (no condensation)
Exhaust Fan Requirements	800 CFM Minimum	1600 CFM Minimum
IEC601-1 CLASSIFICATION		
Class 1, Type B, Transportable, Continuous Operation		
Equipment not suitable for use in the presence of flammable anaesthetic mixture(s).		
Protection against ingress of liquids -Ordinary		

DIMENSIONS

	SINGLE (STS-15)				DUAL (Two STS-30 Pumps Stacked)			
	Weight	Width	Depth	Height	Weight	Width	Depth	Height
STS	235 lbs.	24 in.	29 in.	28 in.	470 lbs.	24 in.	29 in.	56 in.
20-Gallon CAS	70 lbs.	22 in.	16 in.	62 in.	70 lbs.	22 in.	16 in.	16 in.


STS SYSTEM CONFIGURATIONS

STS-3/STS-6 Configuration Dimensions. Figure 2 shows the space requirements for STS-15 and STS-30 model configuration installations. Please note that all CAS heights can be increased or decreased by 8 inches using the adjusting holes in the 4-rail support assembly.



20 GALLON CAS, P/N 54170
(Used for STS-15 and STS-30)


Detail A-
Typical wall bracket
installation.




The 20 gallon CAS **must always** be bolted to the wall with brackets provided in the Tank Accessory Kit. See page 17.

20 gallon CAS tanks **cannot** be stacked on top of a STS-15 pump.

See
Detail A



STS-15, P/N 54110 Configuration



2 STS-15 Pumps Stacked Together Used in
STS-30 P/N 54132 Configurations

Figure 2. STS-15 and STS-30 with 20 Gallon CAS Configuration Dimensions



Any time the power to the STS™ is turned OFF the CAS™ tank will automatically drain.

- The STS™ may be turned “ON/OFF” from a single, convenient location within the dental office using a Remote Control Panel. Remote wiring must be done by a licensed electrician in accordance with local codes.
- The vacuum level is factory preset at 10 in Hg (inches of mercury). This is the reading on the gauge when all HVE's and SE's are CLOSED. If this setting needs to be adjusted contact your dealer to readjust the setting.
- The unit is capable of running continuously. To conserve electricity, the system may be turned off when not in use.
- The CAS™ separation tank has been designed to collect the fluids evacuated during a normal operating day. If an excessive amount of fluids are collected in the CAS™, the protective mechanism in the CAS™ will interrupt the vacuum flow in order for the tank to automatically drain. This process takes approximately 60 seconds. To restore the vacuum to full operation turn OFF the power to the STS™ for a minimum of 10 seconds and then turn back ON.
- Turn the power OFF at the end of the day. This will drain collected liquids in the CAS™ separation tank.

PIPE FLOW SCHEMATIC

Figure 3 shows the overall functional flow of STS systems with a CAS tank. Inlet and outlet connections to the tank is made at the top of the tank. This arrangement provides flexibility in installation, makes maintenance easier and adds to overall tank efficiency.




Figure 3. STS with 20 Gallon CAS Functional Flow Diagram

INSTALLATION INFORMATION

TEMPERATURE CONTROL

The STS unit must be used in a controlled temperature environment. Maintain equipment room temperature between 40 and 105 degrees Fahrenheit. Adequate forced ventilation must be provided across the unit by placing an appropriate exhaust fan opposite an equivalent air intake vent. The fan should be higher than the associated intake vent. Recommended minimum exhaust fan requirements are 800 CFM for STS-15 and 1600 CFM for STS-30

ELECTRICAL

Before connecting the STS to the facilities plumbing complete the following steps:

1. Measure and record the VAC between each leg and ground. If all measurements are equal skip to step 3, if not continue with the next step.
2. Two of the three legs must be equal to use the STS-15. The two equal legs must supply power to the transformer (consult the STS ELECTRICAL SCHEMATIC, Page 12). The transformer is powered with the black and red 18 gauge wires.
3. Three phase motors can spin in the clockwise and counterclockwise direction. Momentarily switch the motor "ON" and ensure that air is flowing out of the motors exhaust fitting at the rear of the unit. If the unit is running properly continue with the PLUMBING section below, if not continue to the next step.
4. Turn the power switch "OFF". Switch any two of the three incoming power lines and repeat step 3.

PLUMBING

For new installations it is recommended to follow the following guidelines:

1. Suction line from the operatories to be a minimum of 2 inches in diameter and must be sloped (1/4 inch minimum for every 10 feet) toward the separation tank. The suction line should not have any sharp right angle bends.
2. The suction line should be connected to the CASTTM separation tank with a short run of 1-1/2 inch diameter flexible tubing.
3. The STSTM system can replace turbine type pumps with larger diameter piping in existing installations.
4. The drain on the base of the separation tank must be connected to an open floor drain capable of handling 20 gallons per minute. Drain pipe size 1-1/2 inch schedule 40 minimum.
5. The drain line should be a short run with a minimum slope of 1/4 inch for every 10 feet toward the drain (avoid any sharp right angle bends).
6. The vent line should be a 3-inch diameter pipe for a single STSTM and two 3-inch diameter pipes for a dual system. The vent should be sloped 1/4 inch for every 10 feet towards the pump. Vent line must be capable of handling vapors and liquids and needs to terminate in a 2-inch end fitting.
7. The outside vent must be protected from rain and animals.
8. A 2-inch flexible air exhaust hose is provided to connect to the 3-inch diameter exhaust (vent), which must terminate with a 2-inch end fitting at the pump. Hose clamps are provided to secure hose to exhaust (vent) pipe.
9. Wash-out port uses a 3/4-inch garden hose connection, affixed to the CASTTM.
10. Refer to pages 9 through 11 for mounting and securing instructions.

SINGLE STS (STS-15) PVC HOSE INSTALLATION.

Note: Each kit provides the required hoses, clamps adapters and wall brackets. No schedule 40 PVC pipe is included. Refer to page 17 for the contents of each kit. If more than 15 feet of hose is needed, order P/N 54118 (order by the foot).

Installation Accessory Packs. Figure 4 shows the pipe and hose connections required for all STS configuration installations using accessory packs as follows:

CAS™ Tank Accessory Pack, P/N 54135 - used to make the following connections:

Between the suction line and CAS tank. See item (1) on Figure 4.

Between the CAS tank and STS pump. See item (2) on Figure 4.

Between the CAS tank and drain line. See item (3) on Figure 4.

STS™ Pump Accessory Pack, P/N 54468 - supplied with each STS pump, used to connect the pump to the facility vent line. See item (4) on Figure 4.

Connection Procedure. Installation of a single STS model configurations consists of making connections between a single 20 gallon CAS tank, P/N 54170, and a STS pump, P/N 54110. Refer to Figure 4 for the connection diagram and perform the following procedure.

1. Refer to Figure 2 and determine the footprint dimensions and connection requirements for the installation. Place the CAS tank and STS pump in position.
2. Install right and left side wall brackets to the CAS and secure the CAS tank to the wall by bolting with supplied lag bolts.
3. Measure and record distance between each connection point.
4. Cut the supplied hose to the length required for each connection.
5. Using industry standard techniques, install pipe to hose adapters as shown in Pipe Connection Diagram Figure 4.
6. Install hoses and secure with associated clamps.

LEGEND

Included in supplied accessory packs.

* Not Included

REQUIRED CONNECTIONS

- (1) CAS tank suction input to facility suction line
- (2) CAS tank pump outlet to STS pump inlet
- (3) CAS tank drain to facility drain line
- (4) STS pump to the facility vent line

Notes:

The CAS accessory kit comes with 15 feet of 1-1/2" Black Flex PVC hose. Cut hose to the length required for the following connections.

1. Suction line to CAS
2. Pump to CAS
3. CAS to Drain.

If more than 15 feet of hose is needed, please order P/N 54118 (order per foot).




Figure 4. STS Pump and CAS Tank Pipe Connection Diagram

INSTALLATION PROCEDURES

DUAL STS (STS-30) PVC HOSE INSTALLATION.

Note: Each kit provides the required hoses, clamps and adapters. No schedule 40 PVC pipe is included. Refer to page 17 for the contents of each kit. If more than 15 feet of hose is needed, order P/N 54118 (order by the foot)

Installation Accessory Packs . Figure 5 shows the pipe and hose connections required for a dual STS configuration installation while Figure 4 shows the other necessary connections. The PVC hose installation is accomplished using accessory packs as follows:

CAS™ Tank Accessory Pack, P/N 54135 - used to make the following connections:

Between the suction line and CAS tank. See item (1) on Figure 4.

Between the CAS tank and drain line. See item (2) on Figure 4.

Between the CAS tank and combined STS pumps. See item (A) on Figure 5.

STST™ Pump Accessory Pack, P/N 54468 - used to connect each pump to the facility vent line. See item (4) on Figure 4.

STST™ Dual System Accessory Kit; P/N 54400 - used to connect two STS pumps together either side-by-side or stacked. See Figure 5.

Connection Procedure. Installation of the dual STS model configuration consists of making connections between a single 20 gallon CAS tank P/N 54170, and two STS pumps, P/N 54110. The associated pumps can be stacked or installed side-by-side depending on the size of the installation site. When stacking the two STS pumps, refer to Figure 6 showing the method to secure the pumps. In any case, the hose connections between the suction line and CAS tank, the CAS tank and drain line and pump to the facility vent line are identical as shown by Figure 5. Use the STST™ Dual System Accessory Kit; P/N 54400 to connect the two STS pumps together for connection to the CAS tank as shown by Figure 5. Refer to Figures 4, 5 and 6 and perform the following procedure.

1. Install check valves supplied in Dual System Accessory Kit; P/N 54400, on both pumps as shown by Figure 5.
2. Refer to Figure 2 and determine the installation footprint dimension and connection requirements. Place the CAS tank and STS pumps in position.
3. If installing side-by-side configuration, proceed to step 6. If stacking, perform step 4.

Note: Make sure to bolt or secure installation to the floor or wall when stacking STS pumps.

4. Stack the two STS pumps as shown by Figure 6 and secure the pumps together using supplied hardware.
5. Refer to Figure 2 and secure the CAS tank to the wall using supplied hardware.
6. Measure and record distance between each connection point as shown in Figure 5.
7. Cut the supplied hose to the length required for each connection.
8. Using industry standard techniques, install pipe to hose adapters as shown in Figures 4 and 5.
9. Install hoses and secure with associated clamps.

INSTALLATION PROCEDURES




Figure 5. Pump Combining Using Dual System Accessory Kit, P/N 54400




Figure 6. Stacking Two STS Pumps Hardware Detail

ELECTRICAL CONNECTIONS



Remove all power to the system prior to working within the electrical box. Contacting high voltage can cause serious injury or even death



All systems must be wired directly from an electrical box that complies with local electrical codes

Refer to Figure 7 when making STS-15 and STS-30 electrical connections. A Boost/Buck Transformer must be installed if the facility power is below the minimum 200V required prime input power.

Note: If the STS pump is spinning counter clockwise and producing pressure at the outlet, turn unit OFF and disconnect from power source. Switch any two of the three incoming power lines in the 4X4 Handy Box.

See INSTALLATION INFORMATION, page 8.




Figure 7. Electrical Schematic for STS-15 and STS-30

Preventive Maintenance

Whenever a service technician fulfills a repair call at the customer site routine checks should also be performed to detect general overall wear, and replacement of parts should be made if necessary before a failure causes a prolonged shut-down. This preventive maintenance program will aid in dependable equipment operation and help reduce breakdown.

Initial Maintenance

The **CASTM** uses patented technology to provide cyclonic action cleaning within the tank. After installation, clean the vacuum lines with Clean StreamTM Cleaner. This is especially necessary when a new system is being installed into existing dental system piping. Using Clean StreamTM Cleaner helps the **STSTM** system to remove any built up deposits in the piping system.

Daily Maintenance Recommendation

Further, to maintain the cleanliness of the **CASTM**, including all the vacuum lines and tubing in the dental system, Air Techniques recommends that the dental office personnel be instructed to use Clean StreamTM Cleaner daily as part of the overall preventive maintenance program.

Yearly Maintenance

Although the **CASTM** uses patented technology to provide cyclonic action cleaning within the tank during operation of the **STSTM** system, the tank does need to be inspected and cleaned yearly to insure proper operation. Procedures to inspect and clean the new and old style **CASTM** tanks are provided by the following pages. Perform the procedures corresponding to the installed tank style to remove any solids trapped in the tank.

Maintenance-Free STS Pumps

All STS pumps are designed for maintenance-free operation. The pump features a powerful permanent split capacitor motor, with a highly reliable contactor and powerful transformer. The motor is completely water and oil-free and provides a dependable operation requiring no scheduled maintenance.

MAINTENANCE

Yearly Tank Maintenance Procedure. Refer to Figure 8 and perform the procedures to remove any solids trapped in the tank.

1. Turn OFF the power to the **STS™**.
2. Remove the vacuum gauge and bushing assembly.
3. Attach the male garden hose adapter to the wash-out port in the **CAS™** (Be sure to include washers).
4. Attach the the water supply to the garden hose adapter (use female garden hose adapter if required).
5. Run water through the **CAS™** for approximately 5 minutes.
6. Remove the male garden hose adapter from the **CAS™** and re-install the gauge and bushing assembly using teflon tape.




Figure 8. Stainless CAS Washout Port Detail



Over or under adjustment of the Vacuum Relief Valve can degrade overall system operation. Adjustments should only be made to keep the suction level at the factory-set 10 InHg level, which is adequate to provide optimum operation.

Vacuum Relief Valve Adjustment

The **STS™** system vacuum level is factory set at 10 InHg (inches of Mercury) as shown by the Vacuum Gauge. This suction level is more than adequate to provide a properly sized system service for a multiple-user dental facility. Whenever the suction level varies above or below the factory set point, adjustments can be made as necessary via the Vacuum Relief Valve. Refer to Figure 9 and adjust the system suction level by performing the following:

1. Access Vacuum Relief Valve by removing Vacuum Relief Valve Cap from the Vacuum Connection Assembly.
2. Hold Adjusting Nut with a 7/16-inch open-end wrench.
3. Using a flat screwdriver, adjust suction level as follows:
 - a. Increase suction by turning adjusting screw clockwise in no more than 1/4 increments.
 - b. Decrease suction by turning adjusting screw counterclockwise in no more than 1/4 increments.

Vacuum Relief Valve Cleaning

A dirty or clogged Vacuum Relief Valve degrades the **STS™** system suction level. Clean the Vacuum Relief Valve by removing Vacuum Relief Valve Cap and carefully pulling the valve from the Vacuum Connection Assembly. Blow out accumulated solid deposits using clean low pressure compressed air.




Figure 9. Vacuum Relief Valve Adjustment and Cleaning Detail

TROUBLESHOOTING

PROBLEM	POSSIBLE CAUSE	POSSIBLE SOLUTIONS
1. No suction.	a. Pump not turned on. b. Pump not running. c. CAS separator tank full. d. CAS hooked up backwards. e. Drain check valve clogged. f. Clogged drain. g. Kinked or collapsed suction hose. h. Separator tank is full and will not drain.	a. Turn pump on.. b. Call your authorized Air Techniques dealer for repair service. c. Shut unit "OFF" for 10 seconds then turn back "ON". d. Connect hose from STS to swivel tee fitting on CAS. e. Call your authorized Air Techniques dealer for repair service. f. Call your local plumber. g. Check the suction line from the unit to the separation tank and the separation tank to the operatory line. If clogged, collapsed or kinked call your authorized Air Techniques dealer for repair service. h. Tank must be hooked up to an open drain. If hooked to an open drain and tank won't drain call your authorized Air Techniques dealer for repair service.
2. Poor or low suction	a. Restricted air exhaust. b. Restricted air suction.	a. Check air exhaust pipe to make sure it conforms to specifications. Check and clear possible restrictions in exhaust line. b. Check the suction line from the unit to the separation tank and the separation tank to the operatory line. If clogged , collapsed or kinked call your authorized Air Techniques dealer for repair service.
3. Excessive suction	a. Relief valve set high or stuck closed. b. Relief valve screen clogged.	a. Call your authorized Air Techniques dealer for repair service. b. Call your authorized Air Techniques dealer for repair service.
4. Pump does not run.	a. Site circuit breaker is "OFF". b. Pump circuit breaker is "OFF". c. Low voltage circuit breaker is open. d. Low voltage remote switch turned "OFF", or not connected properly. e. Electrical problem.	a. Turn "ON" the site circuit breaker. b. Turn "ON" the pump circuit breaker. c. If the white section of the circuit breaker is visible, it is tripped. Flip this section back in to reset breaker. d. Make sure remote switch is turned "ON" or if not using remote, switch yellow and orange wires twisted together. e. Call your authorized Air Techniques dealer for repair service.

INSTALLATION ACCESSORY KITS

STSTM Pump Accessory Pack; P/N 54468, Contents -

Part No.	Component Description	Quantity
57253	Barbed Adapter, 1-1/2-Inch MNPT by 1-1/2 Inch	1
56057	2 1/4-Inch Diameter Hose, 17 Inch Long, 1 Inch Wide Cuff	1
89324	Hose Clamp, 1-9/16- 2-1/2 Inch Maximum by 1/2 Inch	2
57234	1-1/2-Inch Street Elbow 1-1/2 Inches NPT	1
57169	1.31 to 2.25-Inch Diameter Hose Clamp	1
54105	STS User's Manual	1



CAS™ Tank Accessory Pack; P/N 54135, Contents -

Part No.	Component Description	Quantity
54118	1-1/2 Diameter; Black PVC Hose	15 FT
57169	1.31 to 2.25-Inch Diameter Hose Clamp	5
57253	Barbed Adapter, 1-1/2-Inch MNPT X 1-1/2 Inch	4
57234	1-1/2-Inch Street Elbow 1-1/2" NPT	1
31453	1/4-20 X 1-3/4-Inch HEX Head Bolt, Grade 5	4
30610	1/4-Inch Plated Flat Washer	12
30222	1/4-Inch Plated Split Lock Washer	4
30049	1/4-20, Plated HEX Nut	4
A1028	Right-Hand Side CAS Tank Bracket	1
A1029	Left-Hand Side CAS Tank Bracket	1
30124	1/4 X 2-Inch, HEX Head Lag Bolt	4



STSTM Dual System Accessory Kit; P/N 54400, Contents -

Part No.	Component Description	Quantity
54129	Check Valve; 1-1/4 Inch to 1-1/2 Inch, Modified. Includes: Connector Adapter, 54234 PVC Pipe Length, 54291 Barbed P Connector Adapter, 57253	2
54128	PVC Tee PCONN, 1-1/2 Inch FNPT All Sides: D. Includes: Barbed P Connector Adapter, 57253	1
55222	PVC P Connector Adapter, 1-1/2 Inch MNPT X 1-1/2 Inch SPG	2
54118	1-1/2 Diameter; Black PVC Hose	6 FT
57169	Hose Clamp, 1.31 - 2.25 Inch Diameter.	5
31454	1/4-20 X 2-3/4-Inch, HEX Head Bolt, Grade 5, Plated	4
30958	1/4-Inch Flat Washer, 18-8 ST	8
30920	1/4-Inch Split Lock Washer, 18-8 ST	4
30049	1/4-20, Plated HEX Nut	4
54069	Instruction, Dual System Assembly	1



ACCESSORIES/OPTIONS

Accessories/Equipment Options. The following lists the ordering number and description for accessory components available to maintain the **STS™** product family. Contact an authorized Air Techniques' dealer for information.

<u>Part No.</u>	<u>Description</u>
54400	Installation Kit, Dual System
54061	Washout Connector Washer
56200	Vacuum Equalizer
54360	Pre-Installation Guide, STS and CAS
54199	Kit, Remote Vacuum Relief Valve Assembly

Air Techniques and ALLPRO Imaging are leading manufacturers of fine dental, medical and veterinary equipment from air and vacuum systems and X-ray film processors, to an impressive line of new products incorporating the most recent technological advances. These new products, vital components of the innovative professional practice, include intraoral cameras, digital imaging systems, which utilizes phosphor plate technology and, most recently, an intraoral digital X-ray system using sensor technology.

Air Techniques and ALLPRO Imaging have been manufacturing quality products for the dental, medical and veterinary professional since 1962.

Air Techniques and ALLPRO Imaging products are distributed only through authorized dealers. Refer to www.airtechniques.com or www.allproimaging.com to find a dealer in your area.

- ☐ Accent™ Intraoral Digital X-ray Image System
- ☐ Acclaim® Intraoral Digital Video Camera System
- ☐ Acclaim® USB Only Intraoral Digital Video Camera System
- ☐ AirStar®
- ☐ A/T 2000® XR
- ☐ Guardian™ Amalgam Collector
- ☐ Peri-Pro®
- ☐ Provecta 70™
- ☐ Rinsendo Root Canal Disinfection System
- ☐ ScanX®
- ☐ STS™
- ☐ VacStar™



1-800-AIR-TECH
(1-800-247-8324)
www.airtechniques.com

- ☐ 100 Plus
- ☐ 2010
- ☐ Medscope
- ☐ Provecta V
- ☐ ScanX® 12
- ☐ ScanX® DVM
- ☐ ScanX® NDT
- ☐ ScanX® 12 EV
- ☐ ScanX® 14 Portable
- ☐ ScanX® NDT Portable
- ☐ ScanX® 14 In-Counter



1-800-AIR-TECH
(1-800-247-8324)
www.allproimaging.com