




## Data Input YDI 50 Z

Dedicated Keyboard Rated to IP 65 with LCD  
Installation and Operating Instructions



sartorius



- |  |   |
|--|---|
| <p><b>1</b> LCD<br/>LCD Anzeige<br/>Affichage LCD</p> <p><b>2</b> DIP switches and screw for adjusting the contrast of the LCD<br/>DIP-Schalter und Stellschraube für den Kontrast der LCD-Anzeige<br/>Commutateurs DIP et vis de réglage pour le contraste de l'affichage LCD</p> <p><b>3</b> COM B<br/>COM B<br/>COM B</p> <p><b>4</b> COM A<br/>COM A<br/>COM A</p> <p><b>5</b> Manufacturer's label<br/>Typenschild<br/>Plaque signalétique</p> <p><b>6</b> CF key (clear function)<br/>CF-Taste (Clear-Funktion)<br/>Touche CF (Clear-Function)</p> <p><b>7</b> CE key (clear entry)<br/>CE-Taste (Clear Entry)<br/>Touche CE (Clear Entry)</p> | <p><b>8</b> Numeric keypad<br/>Zahlenblock<br/>Bloc numérique</p> <p><b>9</b> Function keys<br/>Funktionstastenfeld<br/>Touches de fonctions</p> <p><b>10</b> PRINT key<br/>PRINT-Taste<br/>Touche d'impression PRINT</p> <p><b>11</b> ON/OFF switch<br/>EIN/AUS-Schalter<br/>Commutateur Marche/Arrêt (ON/OFF)</p> <p><b>12</b> Fuse<br/>Netzsicherung<br/>Fusible de secteur</p> <p><b>13</b> Interface for connecting the balance/scale<br/>Waagenschnittstelle<br/>Interface balance</p> <p><b>14</b> Voltage selector<br/>Spannungswähler<br/>Sélecteur de voltage</p> |
|--|---|

# Contents

	Page
<b>About the Product (Warranty)</b>	4
Storage and Shipping Conditions	4
<b>Equipment Supplied</b>	5
<b>How to Set Up the YDI 50 Z Data Input Keyboard</b>	6
How to Adjust the Voltage Selector to Your Local Line Voltage	6
How to Connect the Data Interface Cable and the Power Cord	6
Adjusting the Data Transmission Rate of the Balance/Scale Interface	7
Turning the YDI 50 Z Data Input Keyboard On and Off	7
<b>Overview of the DIP Switch Settings</b>	8
<b>How to Adjust the Contrast of the LCD</b>	8
<b>Keypad Description</b>	9
<b>Adapting the YDI 50 Z Data Input Keyboard to Individual Requirements</b>	9
<b>Configuration Options</b>	10
<b>Additional Literature</b>	12
<b>Specifications</b>	13
<b>CE Marking</b>	14

## About the Product (Warranty)

The YDI 50 Z Data Input dedicated keyboard is a programmable terminal that can be connected to all MP-8 and MC1 Balances and Scales that are equipped with an interface (MC1 Balances and Scales must have a YDO 01 interface).

**Read these “Installation and Operating Instructions” carefully before operating your Data Input keyboard.**

Do not miss out on the benefits of our full warranty. Please complete the warranty registration card, indicating the date of installation, and return the card to your Sartorius dealer.


### Storage and Shipping Conditions

Storage temperature:  $-40^{\circ}\text{C} \dots +70^{\circ}\text{C}$   
 $-40^{\circ}\text{F} \dots +158^{\circ}\text{F}$

After unpacking the Data Input keyboard, please check it immediately for any visible damage as a result of rough handling during shipment.

**Save all parts of the packaging and the box because you may need to ship your equipment. Before you pack your equipment to ship it, unplug all connected cables to prevent damage.**

## Equipment Supplied



- Power cord
- Interface cable for connection to an MP-8 balance/scale
- Fuse
- Keyboard overlays (for the function keys)


# How to Set Up the YDI 50 Z Data Input Keyboard

Set up the Data Input dedicated keyboard in a suitable area.

If you plan on connecting peripherals (such as a balance or scale with an interface, a printer, etc.), make sure that the ambient conditions are also suitable for these devices.


## How to Adjust the Voltage Selector to Your Local Line Voltage

The Data Input keyboard has been factory-set to 220 volts.



To adjust the voltage selector to the local line voltage rating that you use, unscrew the cover plate from the base of the YDI 50 Z keyboard. Then adjust the voltage selector to your local line voltage rating. Write the new voltage setting on the manufacturer's label.

Whenever you change the operating voltage to **100/120 V**, make sure to exchange the installed **T 80 mA** fuse for the **T 160 mA** fuse that comes with the Data Input keyboard.



The fuse **(12)** is located on the side panel of the Data Input keyboard.

## How to Connect the Data Interface Cable and the Power Cord

Connect the data interface cable to the balance/scale first (**with the unit unplugged**), then to the YDI 50 Z Data Input keyboard. Plug the power cord of the Data Input keyboard into a properly installed wall outlet.

**Important Note:**


If you connect your Data Input keyboard to a power supply which does not have a protective grounding conductor, be sure to ground the unit.

For information on startup of the balance/scale, please refer to the installation and operating instructions supplied with it.

**Adjusting the Data Transmission Rate of the Balance/Scale Interface**

To adapt the Data Input keyboard and the balance/scale, it is necessary to change the data transmission rate in the balance/scale operating program from 1,200 baud to 9,600 baud (see the installation and operating instructions of the balance or scale you are using). For this purpose:


**change code 224 to 227.**

**Turning the YDI 50 Z Data Input Keyboard On and Off**

You can turn the YDI 50 Z on and off using the ON/OFF switch (11).

"0": no power


"1": power is on

**Important Note:**

Always unplug the Data Input keyboard and the balance/scale from line power before you connect or disconnect devices to or from the interface ports.




## Overview of the DIP Switch Settings



The DIP switches **(2)** are located behind a cover plate on the left side of the Data Input keyboard.

- DIP switch S1.1 = OFF (V24 mode)
- DIP switch S1.2 = OFF (V24 mode)
- DIP switch S1.1 = ON (Sartonet)
- DIP switch S1.2 = ON (Sartonet)
- DIP switch S1.3 = not assigned
- DIP switch S1.4 = not assigned
- DIP switch S2.1 =  
(program-dependent switch setting\*)
- DIP switch S2.2 =  
(program-dependent switch setting\*)
- DIP switch S2.3 =  
(ON = data can be transmitted to the LCD)
- DIP switch S2.3 =  
(OFF = data cannot be transmitted to the LCD)
- DIP switch S2.4 =  
(program-dependent switch setting\*)

### How to Adjust the Contrast of the LCD



Located between the DIP switches **(2)** is a small screw which you can use to adjust the contrast of the LCD.


\* Please refer to the program description for information on the switch setting.



# Configuration Options


## COM A and COM B in the Sartonet Mode

In the SARTONET mode, up to 32 units can be connected in parallel to one data interface cable (SARTONET system for networking). A YDI 50 Z Data Input keyboard functions as the "MASTER," and all other connected devices are referred to as "SLAVES" (please see "Additional Literature" on page 12).




## COM A and COM B in the RS-232C Mode

If you operate COM A and COM B in the RS-232C mode, you can directly connect other devices, such as a computer, a printer, a floor scale, a recorder or a spectrometer (please see "Additional Literature" on page 12).



**COM A in the Sartonet Mode and COM B in the RS-232C Mode**

You can connect a printer to COM B and simultaneously work with COM A in the Sartonet mode.



**Additional Literature:**

SARTONET – Description and User’s Manual

Interface Description for Sartorius Balances and Scales (MP 8-4) 98647-000-03

# Specifications

Model	YDI 50 Z
Power requirements	110/120/220/240 V (selectable), approx. 30 VA, 50–60 Hz
Fuse	T 80 mA (160 mA for 110/120 V)
Dimensions of the housing	W 260 mm, D 160 mm (210 mm incl. plug), H (max.) 110 mm
Weight	3 kg
Temperature range	+10°C to +40°C
ROM for programs	64 KB
Data memory	32 KB
LCD	20 characters, 2 lines, backlit
Protection rating	IP 65

# CE-Marking

## CE

The **CE** marking affixed to the equipment indicates that the equipment meets the requirements of the following Directive(s) issued by the Council of the European Union:

### **89/336/EEC “Electromagnetic compatibility (EMC)”**

Applicable European Standards:

Limitation of emissions:	<b>EN 50081-1</b>	Residential, commercial and light industry
	<b>EN 50081-2</b>	Industrial environment
Defined immunity to interference:	<b>EN 50082-1</b>	Residential, commercial and light industry
	<b>EN 50082-2</b>	Industrial environment

#### **Important Note:**

The operator shall be responsible for any modifications to Sartorius equipment and for any connections of cables or equipment not supplied by Sartorius and must check and, if necessary, correct these modifications and connections. On request, Sartorius will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

### **73/23/EEC “Electrical equipment designed for use within certain voltage limits”**

Applicable European Standards:

<b>EN 60950</b>	Safety of information technology equipment including electrical business equipment
<b>EN 61010</b>	Safety requirements for electrical equipment for measurement, control and laboratory use Part 1: General requirements

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.

## Sartorius AG

✉ 37070 Goettingen, Germany

🏠 Weender Landstrasse 94-108, 37075 Goettingen, Germany

☎ (+49/551) 308-0, 📠 (+49/551) 308-32 89

Internet: <http://www.sartorius.com>

Copyright by Sartorius AG, Goettingen, Germany.  
All rights reserved. No part of this publication  
may be reprinted or translated in any form or by any means  
without the prior written permission of Sartorius AG.

The status of the information, specifications and  
illustrations in this manual is indicated by the date  
given below. Sartorius AG reserves the right to  
make changes to the technology, features,  
specifications and design of the equipment  
without notice.

Status: May 1997, Sartorius AG, Goettingen, Germany

Printed in Germany on paper that has been bleached without any use of chlorine - W1A000 - KT  
Publication No.: WYD6068+97054

---

The Sartorius logo consists of the word "sartorius" in a bold, lowercase, sans-serif font. A vertical line passes through the letter "o", which is highlighted with a yellow circle.