

Because we care.

Track Light with Monitor Installation Instructions

Applies to Models:


153829 (-000 thru -008) 153830 (-000 thru -008) 153906 (-010 thru -018) 153907 (-010 thru -018)

Language of origin: English



WARNING

Equipment is not suitable for use in the presence of a flammable anesthetic mixture.



61-3/4"

Equipment Alert

This installation requires at least two service technicians.

.61-3/4" Possible_ **Trolley Travel**

Note

The Light has a 46° Dead Zone at the foot end of the track. The Monitor also has a 24° Dead Zone as shown here.

The Dead Zones move along the track length with each trolley.





Equipment Alert

The Light and Monitor Flex Arms can come into contact with one another.

Note

In addition to a template on the shipping box, a full size paper template, kit 002-1260-00, is also available.






Recommended Ceiling Support Structure

Note

Midmark does not supply the ceiling support structure; these illustrations show suggestions only. Consult a qualified professional to install the appropriate support structure. The weights listed apply only to the track system, not the ceiling support. To prevent the flex arms from drifting, the track assembly must be level and not deflect more than 1/16" over a 12" span.




Equipment Alert

The A/V Cable and Facility Power line require two 2-1/2" minimum diameter access holes thru the ceiling. Visually confirm that you will have room to feed the cables thru the backer board access hole before mounting the track to the ceiling.

WARNING

Use only the Lag Bolts (3 1/2" x 5/16") supplied with the track assembly. Do not substitute a lower grade. All 6 Lag Bolts and Washers must be used to secure track to support structure and they each must be anchored in a supporting stud. Pre-drill pilot holes in structure with a 7/32" diameter drill bit to prevent cracking the support structure. Drill holes 2 3/8" deep.



Service Tip

Leave the plastic wrap on the Trolleys while mounting the track to prevent them from rolling up and down the track.



Equipment Alert

Track with Light and Monitor Weight (Max.)......157 Lbs. (71 kgs)

Recommended Ceiling Support Structure...continued

Note

Midmark does not supply the ceiling support structure; these illustrations show suggestions only. Consult a qualified professional to install the appropriate support structure. The weights listed apply only to the track system, not the ceiling support. To prevent the suspension arms from drifting, the track assembly must be level and not deflect more than 1/16" over a 12" span.



Equipment Alert

The A/V Cable and Facility Power line require two 2-1/2" minimum diameter access holes thru the ceiling. Visually confirm that you will have room to feed the cables thru the backer board access hole before mounting the track to the ceiling.



Equipment Alert

Track with Light and Monitor Weight (Max.)......157 Lbs. (71 kgs)


For Existing 8-5/8" On **Construction with** Center Joists running parallel Pilot Hole Power and A/V Cable access holes Existing Ceiling 1/4" Fender Washer 3-1/2" x 5/16" For installations Lag Screws on suspended ceilings Ceilina Support


Support Structure built to Local Codes



WARNING

Use only the Lag Bolts (3-1/2" x 5/16") supplied with the track assembly. Do not substitute a lower grade. All 6 Lag Bolts and Washers must be used to secure track to support structure and they each must be anchored in a supporting stud. To prevent cracking of support structure, pre-drill pilot holes in structure using a 7/32" diameter drill bit. Drill holes 2-3/8" deep.






Level

DA210601i

Video Source Access Hole



Equipment Alert

Disconnect power from the Track Light Monitor with the use of a wall mounted light switch may cause issues with display settings on the Personal Computer. While power is still being applied to the TLM, the PC has the ability to detect and communicate with the monitor. Communication between the PC and the monitor of the TLM can happen even if the monitor has been placed in Standby Mode using the power button on the monitor or the remote control. If a wall mounted light switch is used to cut power to the TLM, the PC may revert back to the default display settings of a single monitor. To avoid this, please do not install a wall mounted light switch to cut power to the TLM. If a wall monitor light switch must be installed for service reasons, please contain switch inside of a locked cover.

Note


Recommendation, to properly connect multiple Track Light Monitors together on one coax cable connection, please read the following:

• Do not use a splitter to connect multiple TLM's together without an amplifier. Adding splitters decreases the signal power, plus long runs of coax cable further increases the power


Input-

- Proper connection of splitters and amplifiers should occur as follows: SPLITTER AMPLIFIER SPLITTER AMPLIFIER. Do not connect a TLM directly from a splitter output. Best performance is achieved when the TLM's are directly connected to the output of an amplifier.
- Use a Bidirectional Cable-TV Amplifier to split and amplify the cable signal power up so that each output is as strong as or stronger than the original signal. Amplifying the signal helps the signal propagate down a long cable and still maintain enough power at the monitor or TV to display a nice picture.
- It is recommended to use an Bidirectional CableTV Amplifier that comes with an adjustable gain control knob. Gain knob to be used to properly adjust signal for TLM.

Examples of Bidrectional Cable -TV Amplifier Set-up




Monitor and Light Down Tubes



Light Down Tube Weldment

Monitor Down Tube



| Ceiling Height inches (cm) | Tape Color | Coiler Bowl Assembly | Monitor Suspension Tube | | Light Suspension Tube | |
|----------------------------|---------------|-------------------------|---------------------------------|-------------------|--|--------------------|
| | | 029-4871- <u>0X</u> | 052-0678- <u>0X</u> -216 | | 030- <u>184</u> 5- <u>0X</u> -216 | |
| | (see Note) | Dash No. | Dash No. | Tube Length | Dash No. | Tube Length |
| 96*/98 in. (244/249 cm) | Red | -00 | -00 | 5.69 in. (15 cm) | -00 | 11.71 in. (30 cm) |
| 99/102 in. (251/259 cm) | Orange | -01 | -01 | 8.69 in. (22 cm) | -01 | 15.71 in. (40 cm) |
| 103/106 in. (262/269 cm) | Yellow | -02 | -02 | 9.69 in. (25 cm) | -02 | 19.71 in. (50 cm) |
| 107/110 in. (272/279 cm) | Green | -03 | -03 | 10.69 in. (27 cm) | -03 | 23.71 in. (60 cm) |
| 111/114 in. (282/290 cm) | Blue | -04 | -04 | 14.69 in. (37 cm) | -04 | 27.71 in. (70 cm) |
| 115/118 in. (292/300 cm) | Violet | -05 | -05 | 18.69 in. (47 cm) | -05 | 31.71 in. (80 cm) |
| 119/122 in. (302/310 cm) | Gray | -06 | -06 | 22.69 in. (58 cm) | -06 | 35.71 in. (90 cm) |
| 123/126 in. (312/320 cm) | Brown | -07 | -07 | 26.69 in. (68 cm) | -07 | 39.71 in. (101 cm) |
| >126 in. (320 cm) | Black | -08 | -08 | 30.69 in. (78 cm) | -08 | 43.71 in. (111 cm) |

Note:

The "Tape Color" Column corresponds to the correct tape color on the VGA and Coax Harnesses for a specific suspension tube. That color will be exposed directly beneath the coiler bowl assembly. If you should ever move this unit to another location with a different ceiling height, you will need different down tubes and the harness will have to be retracted or extended in the coiler bowl until the correct tape color is exposed.

^{*} There are no Down Tubes available for ceilings lower than 96 " in height

Light Flex Arm Connections and Installation

Note


Have an assistant support the Flex Arm as you make the connections necessary for installation.




Equipment Alert

To ensure quiet operation, apply a thin coat of lubricant to the inside of the Light Down Tube and to the outer surface of the bearing shaft.

Step 1: Connect the light cable connector to the pivot arm. bend any excess cable as shown and push it up into the down tube.



Step 3: Rotate flex arm as needed to align nearest two mating slots and insert the key. Lower collar over key.




Step 2: Slip the flex arm into down tube. Align the flex arm with the track, light facing the head end. (Set screw hole will be 180° across from flex arm.)

Tension Adjustments Neck Cover Arm Pivot (Long Neck) **Adjustment Screws** To adjust tension at Light pivot joint... A) Remove screw cover as shown. B) Loosen set screw. C) Loosen or tighten light pivot adjustment screw. D) Tighten set screw. To adjust tension at Arm / Yoke pivot joints... A) Remove neck cover as shown (see Note). **Light Pivot** Adjustment Screw B) Loosen or tighten appropriate screws to adjust tension. Note: Some halogen models have a short neck on the lighthead. It is not Yoke Pivot Adjustment Screws necessary to remove the neck cover on these models. Set Screw **Short Neck** Arm Pivot Adjustment Screw Yoke Pivot Adjustment Screw To adjust Flex Arm tension... KA974400i A) Remove end cap and top cover as shown. B) Loosen or tighten adjustment nut as desired.


DA3633i


Adjustment Nut



Flex Arm Cable Connections

Step 1: Route the A/V cable into the flex arm channel as shown. The cable connections are staggered for a better fit along the length of the channel.







Note: The screw can be installed on either side of the channel.

Step 2: Install the bracket (two screws) over the cable to hold it in place.

Bracket

Step 4: Replace the bottom down tube cover, lower and screw the top cover onto it.





Side Rail Installation


Note

The side rails are notched at each end to prevent them from hitting the mounting brackets on the track.


End Cap


Step 2: Use a rubber mallet to lightly tap the side rails into the space between the track extrusion and the backer board.

Step 1: Remove the 2 screws and the end cap from the front end of the track.




Step 3: Secure the end cap on the front end of the track with 2 screws.





Monitor and Light Flex Arm Tension Adjustments

Monitor Flex Arm Tension Adjustments



Monitor Specification:

Refer to Manfaturer's User Manual for monitor specifications.

Note

Refer to: **www.exorvision.com** for additional information on service and warranty for the monitor and remote control.

The Model and Serial Number will be required.

Light Flex Arm Tension Adjustments

