

chenomx nmr
suite 7.0



Chenomx NMR Suite 7.0

Copyright © 2001-2010 Chenomx Inc.

Table of Contents

Notices and Trademarks	vii
1. Introduction	1
1.1. Installation	1
Windows (32-bit)	1
Windows (64-bit)	1
Mac OS X	1
Linux (32-bit)	1
Linux (64-bit)	2
1.2. Activation and Licensing	2
1.3. Adding the Chenomx Reference Compounds	3
2. Handling Samples and Spectra	4
2.1. Sample Preparation	4
Internal Standard Solution	4
2.2. Spectral Data Acquisition	5
Pulse Sequences	5
Required NMR Parameters	5
2.3. Spectral Data Processing	5
Data Processed with Other Software	6
3. Common Tasks	8
3.1. Processing a Spectrum	8
3.2. Profiling a Spectrum	8
3.3. Identifying Compounds	9
3.4. Profiling Overlapped Regions in a Spectrum	11
3.5. Determining Compound Concentrations	11
4. The Basics	13
4.1. Opening Files	13
4.2. Importing Spectra	13
4.3. Saving Files	14
4.4. Closing Files	15
4.5. Exiting Modules	15
4.6. Undo and Redo	16
4.7. Change Columns	16
4.8. Sidebar	16
Legend	17
Reference Panel	18
Processing History	18
Compound Sets	18
Simulation	19
4.9. Spectrum Overlays	19
4.10. Spectrum View Tools	20
Show Entire Spectrum	20
Set Zoom	20
Undo/Redo Zoom	20
Auto Zoom	21
Increase/Decrease Vertical Zoom	21
Increase/Decrease Horizontal Zoom	21
Spectrum View	21
Spectrum Thumbnail	22
Select Region	22
4.11. Export Spectrum Image	23
4.12. Copy Spectrum Image to Clipboard	24

4.13. Display Options	24
4.14. Preferences	25
5. Processor	28
5.1. Overview	28
5.2. Quick Reference Guide	29
5.3. Information Tools	30
CSI Editor	31
pH Editor	32
Spectrum Details	34
Processing History (Sidebar)	34
5.4. Processing Tools	35
Phasing	36
Baseline Correction	37
Line Broadening	38
Reference Deconvolution	39
Water Deletion	40
Reverse Spectrum	41
5.5. Importing and Exporting Data	41
Convert Spectra	41
Send to Profiler	42
6. Profiler	43
6.1. Overview	43
6.2. Quick Reference Guide	44
6.3. Profiling Tools	46
Compound Table	46
Custom Colors	48
Pinned Compounds	49
Starred Compounds	50
Quick Searches	50
Concentrations	52
Transforms	53
Fit Automatically	55
Enforce Transform Windows	55
Maximum Concentrations	56
Scale Concentrations	57
Cluster Navigator	58
6.4. Compound Set Tools	59
Compound Sets (Sidebar)	59
Open Library Manager	59
6.5. Importing and Exporting Data	60
Export Profiled Data	60
Export Profiled Compounds	61
Import Profile	62
Spectral Binning	62
Send to Processor	64
6.6. Information Tools	64
Legend (Sidebar)	64
Reference Panel (Sidebar)	65
Spectrum Details	65
Concentration Units	66
7. Spin Simulator	67
7.1. Overview	67
7.2. Quick Reference Guide	68

7.3. Simulating Tools	70
New Simulation	70
Spin Systems	70
Spin Definitions	71
Spin Navigator	73
J-Modifiers	73
7.4. Information Tools	75
Simulation (Sidebar)	75
Simulation Details	76
8. Compound Builder	77
8.1. Overview	77
8.2. Quick Reference Guide	78
8.3. Building Tools	80
New Compound	80
Import Simulation	80
Add Peak	81
Selecting Peaks and Clusters	82
Delete Selected Peaks	83
Adjusting Peaks	83
Grouping Peaks as a Cluster	84
Optimize Selected Peak Shapes	85
Generate Cluster for Region	86
Transform Windows	86
pH Sensitivities	88
Manage Cluster IDs	90
Cluster Navigator	90
8.4. Information Tools	91
Reference Panel (Sidebar)	91
Compound Details	91
Information Panel	92
9. Library Manager	94
9.1. Overview	94
9.2. Quick Reference Guide	95
9.3. Library Tools	96
Compound Table	96
Quick Searches	96
Add Compounds	97
Export Compounds	98
Send to Compound Builder	99
Remove Selected Compounds	99
9.4. Compound Set Tools	100
Compound Sets (Sidebar)	100
Compound Sets	101
Smart Compound Sets	102
Rename Compound Sets	103
Remove Compound Sets	103
10. Tutorials	105
10.1. File Location Note	105
10.2. Processor	105
Introduction	105
Processing a Spectrum	106
10.3. Profiler	110
Introduction	110

Basic Profiling Techniques	110
Advanced Profiling Techniques	112
10.4. Targeted Profiling	116
Identifying Compounds	116
Profiling Clusters	117
10.5. Spin Simulator	122
Introduction	122
Simple Simulation - Valine	122
Complex Simulation - Proline	125
10.6. Compound Builder	128
Introduction	128
Simple Compound - Valine	129
Complex Compound - Proline	131
10.7. Building pH Sensitive Compounds	133
Introduction	133
Building a pH Sensitive Compound	133
10.8. Library Manager	134
Introduction	134
Adding Compounds	134
10.9. Using Expressions	135
Chemistry-Related Expressions	138
Glossary	141
Index	156

Notices and Trademarks

Developed by Chenomx Inc.

© 2001-2010 Chenomx Inc. No copying or reproduction allowed without express written permission.

Technical Support

If you have any questions or concerns, please contact Chenomx Technical Support at:

Phone: +1 780 432 0033
 Fax: +1 780 432 3388
 E-mail: <support@chenomx.com>

When you contact Chenomx Technical Support, you may be asked to provide an application log. The application log contains data that can help us track down issues that you may have experienced while using Chenomx NMR Suite.

To save an application log file

1. Open the Save Application Log dialog.
2. Locate the folder in which you would like to save the application log file.
3. In the File Name field, type a name for the application log file. If you do not specify a different name, the file is saved as `application.log`.
4. Click **Save**.

Tips and Tricks

- In the unlikely event that Chenomx NMR Suite encounters an unrecoverable error and crashes, you will be asked to save an application log the next time you run the software. We strongly recommend that you do so, and send the resulting file to Chenomx Technical Support. The information contained in the log file will help us to improve future versions of Chenomx NMR Suite.

Chenomx is located at:

800 - 10050 112 Street NW
 Edmonton, Alberta
 T5K 2J1
 Canada

Disclaimer of Warranties and Liabilities

The software described in this document is furnished under a License Agreement. The software may be used or copied only in accordance with the terms of such agreement. No part of this publication may be reproduced, stored in a retrieval system, or transmitted in any form or any means electronic or mechanical including photocopying and recording, except as permitted in the License Agreement.

The information contained in this document is subject to change without notice. Chenomx Inc. assumes no responsibility for any errors, omissions, or inaccuracies whatsoever. Chenomx Inc. shall not be held liable for any damages, including special or consequential

Menu Location	Help > Save Application Log...
Icon(s)	
Hotkey(s)	Ctrl-Shift-A Cmd-Shift-A

damages, arising out of the use of the information contained in this document and the software described herein, even if Chenomx Inc. has been advised in advance of the possibility of such damages. Chenomx Inc. makes no warranties with regard to this material including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Trademarks

Chenomx NMR Suite is a trademark of Chenomx Inc.

Microsoft, Windows, Windows XP, Windows Vista and Windows 7 are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Mac and Mac OS are trademarks of Apple Computer, Inc., registered in the U.S. and other countries.

Java is a trademark of Sun Microsystems, Inc in the United States and other countries. Copyright © 1992-2010 Sun Microsystems, Inc.

Linux is a registered trademark of Linus Torvalds.

Other brand and product names are trademarks or registered trademarks of their representative holders and should be noted as such.

Chapter 1: Introduction

Welcome to Chenomx NMR Suite! Before you begin analyzing spectra with Chenomx NMR Suite, you need to install and activate the software.

1.1 Installation

Chenomx NMR Suite requires a Java Runtime Environment (JRE). Please download and install the latest version of Java from Oracle [<http://java.com/>]. The Windows installers for Chenomx NMR Suite will install a JRE automatically, unless you already have a more recent JRE installed.

For Mac OS X, you must have at least OS X 10.5 (Leopard) installed on a 64-bit Intel Macintosh and install the latest version of Java supplied by Apple [<http://www.apple.com>].

After installing the software, you need to activate it. Please see “Activation and Licensing” on page 2 for more details.

Windows (32-bit)

1. Run the installer, named `ChenomxNmrSuite-7.0-windows-setup-x86.exe`.
2. Follow the instructions displayed on screen.
3. You have now installed Chenomx NMR Suite!

Windows (64-bit)

1. Run the installer, named `ChenomxNmrSuite-7.0-windows-setup-x86_64.exe`.
2. Follow the instructions displayed on screen.
3. You have now installed Chenomx NMR Suite!

Mac OS X

1. Open the DMG file, named `ChenomxNmrSuite-7.0-osx.dmg`.
2. Drag and drop the “Chenomx NMR Suite” folder to your desired installation location.
3. You have now installed Chenomx NMR Suite!

Linux (32-bit)

1. Copy the installer script, named `ChenomxNmrSuite-7.0-linux-x86.sh`, to your desired installation folder.
2. Run `ChenomxNmrSuite-7.0-linux-x86.sh`.
3. The install script creates a new folder named `ChenomxNmrSuite`; you can run Chenomx NMR Suite by running the `ChenomxNmrSuite.sh` script located in this folder.
4. You have now installed Chenomx NMR Suite!

Linux (64-bit)

1. Copy the installer script, named `ChenomxNmrSuite-7.0-linux-x86_64.sh`, to your desired installation folder.
2. Run `ChenomxNmrSuite-7.0-linux-x86_64.sh`.
3. A new folder named `ChenomxNmrSuite` is created; you can run `Chenomx NMR Suite` by running the `ChenomxNmrSuite.sh` script located in this folder.
4. You have now installed `Chenomx NMR Suite`!

1.2 Activation and Licensing

Chenomx NMR Suite uses an activation system designed to verify that Chenomx software products are legitimately licensed. After you install the `Chenomx NMR Suite` software, you need to activate it to ensure that you can use all of the modules that you purchased, and that you can open, edit and save changes to `Chenomx` file types (`.cnx`, `.xcpd`, `.cnxlib`, `.xss`).


To activate Chenomx NMR Suite

1. Open any module of `Chenomx NMR Suite`. If you have not previously activated, or if your license has expired, the `Chenomx Software Activation` dialog appears. Otherwise, on the `Help` menu, click **Software Activation...**
2. Click **Request License File...** to create a new e-mail message in your system's default e-mail client, requesting activation of `Chenomx NMR Suite`.

Enter your name in the generated message, and send the message. If you do not have access to e-mail, call Chenomx at +1 780 432 0033.
3. You will receive a license file (`.cnxlic`) from a support staff member. Save the license file to a folder on your computer.
4. Open any module of `Chenomx NMR Suite`. The `Chenomx Software Activation` dialog appears again.
5. Click the **Install License File** button.
6. Browse to the folder in which you saved the license file, select the file, and click the **Open** button.
7. You have now activated `Chenomx NMR Suite`. Click the **Accept** button to run the module that you selected.

If you have not purchased `Chenomx NMR Suite`, or if you have not yet received your license file, you can still use the software in evaluation mode. When evaluating `Chenomx NMR Suite`, you have complete access to every module, with the following exceptions:

- Saving profiled `.cnx` files, using **Send to Processor**, exporting profiled data and performing spectral binning in `Profiler`. Note that you can still process and save `.cnx` files in `Processor`.
- Using **Send to Profiler** in `Processor`.
- Saving `.xcpd` files in `Compound Builder`



To evaluate Chenomx NMR Suite

1. Open any module of Chenomx NMR Suite; the Chenomx Software Activation dialog appears.
2. Click **Evaluate**.
3. You can now run Chenomx NMR Suite in Evaluation mode. Click the **Accept** button to run the module that you selected.

You can view your activation details at any time by clicking **Help > Software Activation** in any module.

Tips and Tricks

- If you do not have a default e-mail client configured, or if a message is not automatically generated when you click **Request License File...**, you can compose the message yourself. To ensure a timely response, include your name, the version and edition of Chenomx NMR Suite that you would like to activate, and the activation key.
- The activation key is unique to the computer used to generate it, and you cannot use it to activate the software on other computers. If you are installing Chenomx NMR Suite on multiple computers, you need to send activation keys for each to Chenomx.
- Chenomx NMR Suite is licensed per seat, meaning that one license can be used on one computer.
- Multiple users can use a one-seat license. To enable multiple users to share a license on the same computer, select the license file as described above, check the box marked **License extends to other users on this computer**, then click the **Accept** button.
- Installing a license for all users may require running Chenomx NMR Suite with administrator privileges. However, once the license is installed, all users (administrator or not) can run Chenomx NMR Suite.
- New major releases (for example, version 5.x to 6.x) require reactivation, while minor releases (for example, version 5.0 to 5.1) do not.

1.3 Adding the Chenomx Reference Compounds

To add the Chenomx Reference Compounds

1. Open the Library Manager module of Chenomx NMR Suite.
2. On the **Library** menu, click **Add Compounds...**
3. Select Chenomx Reference Compounds, and click **Next**.
4. Select the frequencies for which you would like to install Chenomx Reference Compounds, and click **Finish**.

Reference Compound Sets corresponding to the Chenomx Reference Compounds that you added are created automatically.

Chapter 2: Handling Samples and Spectra

Accurate analysis of your NMR spectra with Chenomx NMR Suite requires preparing good samples and acquiring high-quality spectra.

2.1 Sample Preparation

Preparing your samples properly helps you get the best possible results from your analyses with Chenomx NMR Suite.

Internal Standard Solution

Your samples must contain an internal standard for accurate quantification using Chenomx NMR Suite. The internal standard solution recommended for use with Chenomx NMR Suite contains the following components:

To prepare the Chenomx-recommended internal standard solution:

1. To the desired volume of 99.9% D₂O, add 5 mM 2,2-Dimethyl-2-silapentane-5-sulfonate (DSS) as a chemical shift indicator (CSI) and 0.2% w/v sodium azide (NaN₃) to inhibit bacterial growth.

Remember that DSS is hygroscopic in the solid form. For accurate concentrations from a Chenomx NMR Suite analysis, you must determine the final DSS concentration accurately.

2. For optional pH measurement in Processor, add one or more of 100 mM imidazole, 100 mM creatinine and 20 mM DFTMP.
3. Adjust the resulting solution to about pH 6.5.

Add the internal standard as a 1 in 10 (10% v/v) spike to each of your samples. For example, add 70 µL of the internal standard to 630 µL of a sample to be analyzed for 700 µL of spiked sample.

Tips and Tricks

- Instead of DSS, you may use 3-trimethylsilylpropionate (TSP or TMSP) or formate as a CSI. However, we do not recommend TSP or formate, as their chemical shifts are sensitive to pH.
- For accurate sample analysis with Chenomx NMR Suite, the pH of your samples should be between 4 and 9, and ideally close to 7 for manual analysis.
- Spectra that are poorly collected, or collected from samples that are poorly prepared, may be difficult or even impossible to convert and analyze accurately.
- If you choose not to add any pH indicators to your samples, you need to supply a measured pH value for proper analysis with Chenomx NMR Suite. More accurate pH values for your samples allow you to make better use of pH-sensitive signatures to profile spectra of the samples.

2.2 Spectral Data Acquisition

For the best results from your analysis with Chenomx NMR Suite, you must acquire high-quality spectra. Pay particular attention to proper shimming, phasing and baseline correction to produce spectra that you can accurately analyze and interpret using Chenomx NMR Suite.

Chenomx NMR Suite can automatically detect DSS or TSP signals in most spectra. However, automatic detection relies on the water peak being centered in the spectrum. For best results, ensure that the water peak is centered in your spectra.

Pulse Sequences

Using Chenomx NMR Suite, you can quantitatively analyze spectra produced by NOESY pulse sequences with t_{mix} of 100 ms. You can use other pulse sequences for identification purposes, but quantification based on alternate pulse sequences may be unreliable.

Required NMR Parameters

For the best results from Chenomx NMR Suite when analyzing your data, use the following parameters for NMR data collection:

- Temperature = 25 °C
- Acquisition time = 4 s
- Relaxation delay = 1 s
- Spectral width = 12 ppm

Tips and Tricks

- Chenomx NMR Suite is capable of automatically detecting DSS in your spectra. However, automatic detection relies on the water peak in your aqueous ^1H NMR spectrum being centered in the spectrum.

2.3 Spectral Data Processing

Chenomx NMR Suite can import raw data from a number of common sources. Supported raw file formats include:

- Varian (`.fid` files)
- Bruker (`.fid` files)
- JEOL (`.jdf` files)
- NMRPipe (`.fid` files)

To analyze a spectrum with Chenomx NMR Suite, you must first process your raw spectrum file using either the included Processor module or one of the following software packages.

- VNMR and VNMRj (Varian `phasefile` files)
- XWIN-NMR (Bruker `1r` files)
- Delta (JEOL `.jdf` files)
- NMRPipe (`.ft2` files)

- ACD/Labs software package. Export your spectra to the JCAMP-DX format (.jdx files).

Zero-Filling

- To ensure maximum compatibility with the Chenomx Reference Compounds, zero-fill your data.
- When zero-filling, use the nearest power of 2 that is at least twice as large as the number of points acquired in the original spectrum (e.g., for 16000 points in the original spectrum, use zero-filling of 32K points).

To do this in Processor, select **Automatic** in the import dialog under zero-filling.

Data Transformation

- Do not apply any transformations like sine bell or Gaussian functions to your spectra.

Basic Processing

- For the best results, you must properly phase and baseline correct all of your spectra.

Line Widths

- For the best results, your spectra should have CSI linewidths of no more than 1.5 Hz. Use an appropriate combination of reference deconvolution and line broadening to adjust the linewidths.

Data Processed with Other Software

If you would like to process your spectra using software packages other than the Processor module, you need to ensure that the Fourier-transformed (frequency domain) data from your processed spectra are available for Chenomx NMR Suite to import.

Spectra processed from most third-party processing software require no additional preparation, including XWIN-NMR from Bruker, Delta from JEOL, and NMRPipe. Spectra processed elsewhere and converted to JCAMP-DX format also require no additional preparation.

If you have spectra that you processed using VNMR or VNMRJ, and you would like to retain the VNMR processing when you import the spectra into Chenomx NMR Suite, you will need to copy the phasefile into the fid folder. Note that it is possible to automate this process using VNMR, if it is something you expect to use often. Consult your VNMR documentation for further information on automating tasks.

To manually copy the Varian phasefile into the fid folder

1. Open the spectrum in VNMR and process the spectrum.
2. Change experiments or quit VNMR, making note of the experiment number in which you have opened the spectrum (e.g., exp2).
3. In a Solaris terminal window, copy the `datdir` folder into the `mydata.fid` folder where `mydata.fid` is your Varian spectrum. `mydata.fid` is the folder that you need to access to open the spectrum using Chenomx NMR Suite.

Use the following commands:

```
cd ~/vnmr/sys/data/mydata.fid
```

```
cp -r ~/vnmrsys/expn/datdir .
```

where *expn* is the experiment number noted previously, and *mydata.fid* is the fid folder you want to access with Chenomx NMR Suite. In this example, *expn* would be exp2.

Chapter 3: Common Tasks

Chenomx NMR Suite can help you to perform many useful tasks. You can find overviews of some common tasks in this chapter. Detailed descriptions of each function mentioned here appear in later chapters.

3.1 Processing a Spectrum

Before you analyze a spectrum in Chenomx NMR Suite, whether in Profiler or in Compound Builder, you need to process the spectrum.

To process a spectrum in Processor

1. Open the spectrum (see “Opening Files” on page 13 and “Importing Spectra” on page 13). If necessary, reverse the spectrum plot.

Reverse Spectrum flips the spectrum in the frequency domain. Reversing is only necessary for some Bruker spectra.
2. Adjust the phasing (see “Phasing” on page 36).
3. Apply baseline correction (see “Baseline Correction” on page 37) as needed.
4. Define the CSI parameters (see “CSI Editor” on page 31).
5. (Optional) Apply reference deconvolution to correct lineshape distortions (see “Reference Deconvolution” on page 39), or apply line broadening to improve the signal-to-noise ratio (see “Line Broadening” on page 38).
6. If you applied reference deconvolution or line broadening, readjust the CSI parameters to take into account any newly-introduced changes in lineshape (see “CSI Editor” on page 31).
7. Determine the sample pH using imidazole, creatinine or DFTMP (see “pH Editor” on page 32).
8. (Optional) Delete the water peak from the spectrum (see “Water Deletion” on page 40).

Water deletion may not be necessary or desirable for all samples.
9. Save the spectrum (see “Saving Files” on page 14), or jump to Profiler (see Chapter 6, *Profiler* on page 43) to continue with your analysis.

3.2 Profiling a Spectrum

Profiling a spectrum using targeted profiling is a very flexible task. You can vary the number of compounds profiled in a particular spectrum to suit a variety of timeframes, research goals, or other requirements. One of the most important aspects of your approach to profiling a spectrum should be maintaining a consistent profiling strategy, as this gives you more consistent and more easily-interpreted results.

For more details on specific profiling techniques and how to apply them to your spectra, please see “Targeted Profiling” on page 116.

To profile a spectrum using targeted profiling

1. Process the spectrum (see “Processing a Spectrum” on page 8).
2. (Optional) Create a compound set containing the compounds that you would like to analyze (see “Compound Set Tools” on page 100).
3. Open the spectrum in Profiler (see “Opening Files” on page 13).
4. Select a compound set containing compounds appropriate for the sample. If you created a compound set above, select the compound set that you created.
5. Identify and profile prominent signals in the spectrum (see “Identifying Compounds” on page 9).

Some compounds in the spectrum are very easy to pick out, including CSI compounds like DSS or TSP and pH indicators like imidazole, creatinine or DFTMP. Depending on the type of sample that you are analyzing, there may be other compounds that also stand out.

6. Systematically review the remaining signals in the spectrum; identify and profile as needed.

You can use information and external links from the Reference Panel (see “Reference Panel (Sidebar)” on page 65) to help establish the presence of particular compounds in the samples that you are analyzing.

3.3 Identifying Compounds

Chenomx NMR Suite provides computer-assisted methods to profile spectra. To identify compounds in a sample, you need to match the compound signatures from your library with patterns in the sample spectrum. There are several ways that you can identify compounds in a sample.

The Selection Line must be enabled via the Legend to view individual compound signatures. It is enabled by default.

To identify specific compounds

1. Select at least one compound set in the Compound Sets view (see “Compound Sets (Sidebar)” on page 59).
2. Locate a compound in the Compound Table that you believe is present in the sample and select it (see “Compound Table” on page 46).
3. For each cluster, compare the shape of the preview line (light blue) and the spectrum line (black) and decide whether it is reasonable that the preview line could be contributing to the spectrum line in the displayed region.
4. If the shape appears reasonable, click the preview line, and adjust the resulting shape by determining appropriate heights and transforms for the clusters of the compound (see “Concentrations” on page 52 and “Transforms” on page 53).
5. If the shape does not appear reasonable, continue by selecting another compound in the Compound Table.

To identify possible compounds causing a specific peak

1. Select at least one compound set in the Compound Sets view (see “Compound Sets (Sidebar)” on page 59).
2. Right-click on the peak, and click **Filter for Compounds Near x.xx ppm**.

3. In the Compound Table, select one of the listed compounds to display the preview line and Cluster Navigator for the compound (see “Cluster Navigator” on page 58).
4. Compare the shape of the preview line (light blue) and the spectrum line (black) and decide whether or not it is reasonable that the preview line could be contributing to the spectrum line in the displayed region.
5. If the shape appears reasonable in the region near the cluster that you are investigating, use the Cluster Navigator to verify that the other clusters for the compound also appear to match reasonably well (see “Cluster Navigator” on page 58).
6. If all clusters appear reasonably well-matched, apply an automatic fit, then adjust the concentration and transforms for the compound as needed (see “Concentrations” on page 52, “Transforms” on page 53 and “Fit Automatically” on page 55).
7. If the shape does not appear reasonable, continue evaluating the rest of the compounds in the list.
8. When you have finished with the filtered list, clear the frequency filter to continue your analysis with the full compound list (see “Quick Searches” on page 50).
9. If none of the compounds in the list of potential candidates provides a reasonable shape near the peak that you are investigating, consider the possibility that the compound responsible for the peaks is not in your compound library.

To identify all compounds contributing to a particular region

1. Select at least one compound set in the Compound Set Panel (see “Compound Sets (Sidebar)” on page 59).
2. Select a region in the Spectrum View (see “Select Region” on page 22).
3. Filter the Compound Table for compounds in the selected region (see “Quick Searches” on page 50).
4. In the Compound Table, select one of the listed compounds.
5. Compare the shape of the preview line (light blue) and the spectrum line (black) and decide whether or not it is reasonable that the preview line could be contributing to the spectrum line in the displayed region.
6. If the shape appears reasonable in the region near the peak that you are investigating, use the Cluster Navigator to verify that the other clusters for the compound also appear to match reasonably well (see “Cluster Navigator” on page 58).
7. If all clusters appear reasonably well-matched, apply an automatic fit, then adjust the concentration and transforms for the compound as needed (see “Concentrations” on page 52, “Transforms” on page 53 and “Fit Automatically” on page 55).
8. If the shape does not appear reasonable, continue evaluating the remainder of the compounds in the list.
9. When you have finished with the filtered list, clear the frequency filter to continue your analysis with the full compound list (see “Quick Searches” on page 50).

3.4 Profiling Overlapped Regions in a Spectrum

For all but the simplest spectra, you need to profile compounds whose signatures overlap. Approaches to profiling these overlap regions are as varied as are the number of people profiling them, but the general guidelines below may help you to develop your own strategy. As with your more general approach to profiling spectra, a consistent strategy for profiling overlap regions gives you more consistent results.

For more details on this and other profiling strategies, see “Targeted Profiling” on page 116.

To profile compounds in an overlapped region of the spectrum

1. Identify likely compounds contributing to an overlapped region (see “Identifying Compounds” on page 9).
2. Evaluate each compound for the presence of additional, non-overlapped peaks elsewhere in the spectrum (see “Cluster Navigator” on page 58).
 - a. If they are present:
 - i. Profile the non-overlapped peaks to help establish a concentration range for the compound (see “Concentrations” on page 52).
 - ii. Look for clusters of peaks in the overlapped region with some clearly visible peaks. Such clusters occur mostly at the edges of the overlap region, but remember to look throughout the region, as some peaks stand out nicely between the peaks of another cluster.
 - iii. Use these partially-exposed peaks to establish an optimal frequency for the cluster, using transforms (see “Transforms” on page 53).
 - b. If they are not present:
 - i. Look for clusters of peaks in the overlapped region with some clearly visible peaks.
 - ii. Use these partially-exposed peaks to establish an optimal frequency for the cluster, using transforms (see “Transforms” on page 53).
3. Once you have approximately profiled most of the groups in the overlap region, fine-tune the heights and frequencies of each cluster in the region to obtain an optimal shape. You can use the subtraction line to guide your modifications, as a well-matched cluster will leave a subtraction line that resembles a normal spectrum with the effects of the cluster removed.
4. In overlap regions, profile multi-peak clusters before single-peak clusters. Also, if urea is present, you may want to profile it last, to minimize its interference with other compounds in the spectrum.

3.5 Determining Compound Concentrations

Once you have identified compounds present in a sample, you need to determine their concentrations based on the data in the spectrum. If you have automatically fit a compound, you may need to refine the suggested concentration.

For more details on techniques for determining optimal profiles for compounds in experimental spectra, see “Targeted Profiling” on page 116.

To determine the concentration of a compound

1. If the subtraction line (green) is not currently visible, turn it on (see “Legend” on page 17).
2. Set a concentration for the compound so that you can see the clusters clearly (see “Concentrations” on page 52).
3. Use the Cluster Navigator to locate a cluster that you can use to approximate the compound concentration (see “Cluster Navigator” on page 58).

When choosing a cluster with which to establish a compound's concentration, use a cluster that overlaps other clusters as little as possible.

4. Click and drag the height control for the compound and observe the subtraction line. Adjust the height of the compound until the subtraction line under the cluster approximates a normal spectrum. You may need to move (transform) the cluster side to side somewhat to properly adjust the subtraction line (see “Subtraction Line” on page 118, “Concentrations” on page 52 and “Transforms” on page 53).
5. If any of the clusters for the current compound overlap with compounds that you have already profiled, you may need to readjust those compounds.

Chapter 4: The Basics

Chenomx NMR Suite includes many advanced features that help you process and analyze spectra, create custom compound signatures, and manage compound signatures and compound sets in your library. Underlying all of these advanced features is a set of more basic ones that are fundamental to the proper use of Chenomx NMR Suite. The basics include file operations, navigation, display options and preferences that directly affect your day-to-day use of Chenomx NMR Suite.

4.1 Opening Files

You can open and overlay Chenomx NMR Suite spectra (.cnx) in Processor and Profiler, and overlay them in Compound Builder and Spin Simulator. You can open Chenomx compound files (.xcpd) in Compound Builder, and you can open Chenomx spin simulation files (.xss) in Spin Simulator and import them in Compound Builder.

Menu Location	File > Open
Icon(s)	
Hotkey(s)	Ctrl-O Cmd-O

How do I open a file?

1. Open a module.
2. Open the file selection dialog.
3. Locate the folder containing the spectrum.
4. Select the file and click **Open**.
5. If the file you are opening requires importing, select import options as necessary.

Tips and Tricks

- If you are opening a Varian spectrum named `Example.fid`, the raw data is stored in the file named `Example.fid\fid`, while the processed data is stored in the file named `Example.fid\datdir\phasefile`.
- If you are opening a Bruker spectrum named `Example\1`, the raw data is stored in the file named `Example\1\fid`, while the processed data is stored in the file named `Example\1\data\1\1r`.

4.2 Importing Spectra

When you open a spectrum in a supported non-Chenomx format in Processor or Profiler, you need to select import options for conversion into the Chenomx spectrum format (.cnx). You can import raw spectrum data in Varian (fid), Bruker (fid), JEOL (.jdf) or NMRPipe (.fid) file formats, and processed spectra in Varian (phasefile), Bruker (1r), JEOL (.jdf), NMRPipe (.ft2) or JCAMP-DX (.jdx) file formats.

When you import processed spectra, your import options are somewhat limited, since most of the processing that you may want to apply is assumed to have already taken place.

Table 4.1. Import NMR Data Options

CSI	Specify the Chemical Shape Indicator (CSI) used in the sample. You can choose DSS, TSP or formate.
Concentration (mM)	Enter the concentration of the specified CSI in millimoles per litre.
pH: Automatic	Attempt to automatically detect the pH of the sample, subject to the specified uncertainty, using signals from imidazole, creatinine and DFTMP. You must have at least one of these pH indicators present in the sample to attempt automatic pH detection. If the pH is not detected, it will be set to 7.00; you can change it manually after the spectrum is imported.
pH	Specify a pH value manually, subject to the specified uncertainty. Useful if your sample does not contain any of the pH indicators.
Advanced Details	Access advanced import options.
Defaults	Restore all import settings for the current file type to their default values.
Advanced Options	
Line Broadening (Hz)	Apply an exponential FID weighting (line broadening) function to your data. You can enter any value between 0 and 1000 Hz.
Zero Fill To	Enable zero-filling. Automatic zero-fills to the nearest power of two that is at least twice as large as the number of points acquired in your spectrum. You can also choose a specific number of data points, if you prefer.
Phasing (Varian only)	Automatically apply phase adjustments to your spectrum. Auto applies automatically-determined phase angles, while Import VNMR Phasing imports the phase angles found in the Varian procpa file. Only available when you import Varian raw data (fid).
Auto Phasing	Apply automatically-determined phase angles to your spectrum. Only available when you import non-Varian raw data
Auto Baseline Correction	Automatically apply baseline correction to your spectrum. Linear (Drift) applies a simple two-point linear baseline correction, while Spline applies a cubic spline baseline correction based on automatically-determined breakpoints.
Delete Water Region	Remove the water region from the spectrum during the import process.
Reverse Spectrum (Bruker only)	Reverse your spectrum horizontally. Only available when you import Bruker spectrum files (fid).

4.3 Saving Files

Saving a file ensures that any changes you have made are preserved for future reference. When you are working on a complex task, save often.

Menu Location	File > Save/Save As...
Icon(s)	
Hotkey(s)	Ctrl-S Cmd-S

How do I save a file?

1. Click **Save**.
2. If you have not previously specified a location and filename, do so now.

How do I save a file with a different file name?

1. Click **Save As**.
2. Locate the folder in which you would like to save the spectrum.
3. In the File Name field, type the new name of the spectrum file.
4. Click **Save**.

Tips and Tricks

- Saving compound files in Compound Builder or simulation files in Spin Simulator preserves only the mathematical peak descriptions created during the fitting process. Compound Builder and Spin Simulator do not retain spectrum overlays in the saved files.

4.4 Closing Files

When you finish working with a file, you can close it and open a new one without exiting the program.

Menu Location	File > Close
Icon(s)	
Hotkey(s)	Ctrl-W Cmd-W

How do I close a file?

1. Click **Close**.
2. If you have made changes to the file and they are not saved, choose whether you would like to save the changes.
3. If necessary, choose a location and filename for the file.

Tips and Tricks

- You cannot close a file that has unsaved changes without having the option of saving the changes. If there are no unsaved changes, the file is simply closed without further interaction.

4.5 Exiting Modules

When you are finished working with a module, you can exit the module.

Menu Location	File > Exit
Icon(s)	

How do I exit a module?

1. Click **Exit**.
2. If you have made changes to the currently open file and they are not saved, choose whether you would like to save the changes.
3. If necessary, choose a location and filename for the file.

Tips and Tricks

- You cannot exit a module when a file is open with unsaved changes, without having the option of saving the changes. If there are no unsaved changes, the module exits without any further interaction.

4.6 Undo and Redo

You can undo recent changes to the currently open file, and redo recently undone changes.

How do I undo a recent change?

1. Click **Undo**.

How do I redo a recently undone change?

1. Click **Redo**.

Tips and Tricks

- Only the last 100 changes are available to be undone. If the change that you want to undo is more than 100 steps back, you cannot undo the change.
- You can only redo changes that you have just undone. If you undo a change, then make other changes, you can no longer redo the change that was undone.

Menu Location	Edit > Undo/Redo
Icon(s)	
Hotkey(s)	Ctrl-Z / Ctrl-Y Cmd-Z / Cmd-Y

4.7 Change Columns

In modules that have Compound Tables (Profiler and Library Manager), you can specify which of the available columns are visible in the Compound Table.

How do I change the visible columns in the Compound Table?

1. Open the Change Columns dialog.
2. Select the check boxes for the columns you would like to display, and clear the check boxes for the columns you would like to hide.
3. Click OK.

Tips and Tricks

- You can quickly display or hide a particular column by right-clicking in the header row of the Compound Table and clicking the column that you would like to display or hide.
- The Concentration, Maximum, Autofit Clusters, Matched Clusters, Color, Pinned, Starred and Transforms columns are only available in Profiler.

Menu Location	View > Change Columns...
Icon(s)	

4.8 Sidebar

The Sidebar is a dedicated area of the screen used to display frequently-used information related to the Spectrum View. Several views are available in the Sidebar, depending on the module currently in use. Available views include Legend (all modules except Library Manager), Reference Panel (Profiler and Compound Builder), Compound Sets (Profiler), Simulation (Spin Simulator) and Processing History (Processor).

Related Topics

- “Information Tools” on page 30, in Chapter 5: *Processor*
- “Information Tools” on page 64, in Chapter 6: *Profiler*
- “Information Tools” on page 75, in Chapter 7: *Spin Simulator*
- “Information Tools” on page 91 in Chapter 8: *Compound Builder*

Menu Location	View > Show Sidebar
Icon(s)	
Hotkey(s)	Ctrl-Shift-S Cmd-Shift-S

Legend

The Legend is an interactive display of the various lines currently appearing in the Spectrum View. You can use the Legend to toggle the display of each line, change their colors, or even remove overlays or individual compound lines from the display completely.

How do I toggle the display of a line?

1. Switch to the Legend view in the Sidebar.
2. Select the check box next to a line to display it, or clear the check box to remove it from the Spectrum View.

How do I change a line color using the Legend?

1. Click the colored box beside an entry in the Legend.
2. Click one of the preconfigured colors in the color chooser, or click **More...** to choose from a larger selection of colors.

How do I remove an overlay using the Legend?

1. Click the red X beside an overlay entry in the Legend.

Tips and Tricks

- The topmost section of the Legend (Profiler Lines in Profiler, Compound Builder Lines in Compound Builder, and so on) contains the default lines that are always available in the current module.
- Changing the color of a compound line using the Legend in Profiler is equivalent to setting a custom color using the Compound Table. The custom color is saved with the spectrum, and will be used when the spectrum is reopened. All other color changes made via the Legend are temporary, and are reset to default colors the next time you run Chenomx NMR Suite.
- You can change the default colors of various lines in a module using the preferences for that module.
- Each of the default lines for a module can be toggled on or off temporarily using the Legend. To change the default state of these lines, use the Preferences for the module that you would like to change.

Related Topics

- “Spectrum Overlays” on page 19
- “Preferences” on page 25
- “Compound Table” on page 46
- “Custom Colors” on page 48

Menu Location	View > Sidebar > Legend
Icon(s)	
Hotkey(s)	Ctrl-Shift-L Cmd-Shift-L

Reference Panel

Menu Location	View > Sidebar > Reference Panel
Icon(s)	
Hotkey(s)	Ctrl-Shift-R Cmd-Shift-R

The Reference Panel Sidebar view displays information about specific compounds. The displayed information includes a structure image, molecular formula and weight, alternate names and links to external reference sites. You can use this information to verify the identity of the compound that you are fitting and better understand the signature and behaviour of the compound in your spectra.

The compound represented in the Reference Panel in Profiler is the compound that you currently have selected in the Compound Table. In Compound Builder, it is the compound matching the Chenomx Compound ID that you have specified in the Compound Details.

Tips and Tricks

- Some sections in the Reference Panel are collapsed by default, including Alternate Names and Alternate CAS Registry. To expand them, click the + symbol next to the section name. Click the - symbol to collapse the section again.
- Each of the External Database References are links to the indicated database. Click a link to open the associated database entry for the current compound in your system's web browser.
- Chenomx NMR Suite retrieves all Reference Panel information from a Chenomx server named reference.chenomx.com, to ensure that you always have the most accurate information. If you access the internet through an HTTP proxy server, or have no internet access, you may want to change the internet preferences for Chenomx NMR Suite.
- If Chenomx NMR Suite cannot access the internet, the content of the Reference Panel will instead come from a local cached copy.

Related Topics

- “Preferences” on page 25
- “Profiling Tools” on page 46
- “Building Tools” on page 80
- “Compound Details” on page 91

Processing History

The Processing History view in Processor tracks processing effects that you add to a spectrum, recording details about each effect and letting you reverse them as needed.

Related Topics

- “Processing History (Sidebar)” on page 34

Compound Sets

The Compound Sets view in Profiler lets you modify the contents of the Compound Table to include any combination of the compound sets that you have created in Library Manager. The panel contains a list of every compound set that is currently defined. Contents of the checked compound sets appear in the Compound Table.

Related Topics

- “Compound Sets (Sidebar)” on page 59

Simulation

The Simulation view lets you select, view and edit elements of the currently open spin simulation, including spin systems, spin definitions and J-modifiers.

Related Topics

- “Simulation (Sidebar)” on page 75

4.9 Spectrum Overlays

You can overlay Chenomx spectra (.cnx) in any module with a Spectrum View. Overlaying a spectrum in Compound Builder or Spin Simulator can aid the fitting process, while overlays in Profiler and Processor can provide useful visual references to other spectra in a dataset.

How do I overlay spectra?

1. Click **Add Spectrum Overlays (Set Spectrum Overlay** in Spin Simulator and Compound Builder).
2. Locate the folder containing the spectra that you would like to overlay, and select the spectra.
3. Click **Overlay**.

How do I remove a single spectrum overlay?

1. In the Legend view of the Sidebar, click the red X beside the overlay that you would like to remove.

How do I remove all spectrum overlays?

1. Click **Remove All Spectrum Overlays (Clear Spectrum Overlay** in Spin Simulator and Compound Builder)).

Tips and Tricks

- You can also remove all spectrum overlays by clicking on the red X beside the Overlay Lines section in the Legend.
- You can only overlay spectra in the Chenomx NMR Suite format (.cnx). If the spectrum you would like to use is in a different, supported format, you need to import and save it in Processor before using it as an overlay.
- In Compound Builder and Spin Simulator, you can only have one spectrum overlay at a time. In Profiler and Processor, you can have up to fifty spectrum overlays.
- In Compound Builder, if you overlay a spectrum whose frequency does not match the magnet frequency specified in the Compound Details, you are notified of the mismatch, and given the chance to update the Compound Details to match the spectrum overlay.

Related Topics

- “Sidebar” on page 16
- Chapter 5, *Processor* on page 28
- “Compound Details” on page 91

Menu Location	Tools > Add Spectrum Overlays Tools > Set Spectrum Overlay
Icon(s)	
Hotkey(s)	Ctrl-Shift-O Cmd-Shift-O

Menu Location	Tools > Remove All Spectrum Overlays Tools > Clear Spectrum Overlay
Icon(s)	

4.10 Spectrum View Tools

The Spectrum View contains the spectrum as a graph of intensity vs. frequency. In addition to the simple zoom and scroll that you may be familiar with in other software, Chenomx NMR Suite includes a number of specialized tools for manipulating the Spectrum View.

Show Entire Spectrum

Show Entire Spectrum lets you quickly expand the current view, setting both the vertical and horizontal scales to their largest possible values.

Menu Location	View > Zoom > Show Entire Spectrum
Icon(s)	
Hotkey(s)	Ctrl-Minus Cmd-Minus

Tips and Tricks

- Show Entire Spectrum is not available when the full spectrum is already visible.
- Show Entire Spectrum is applied by default when you first open any spectrum, unless you have set specific zoom settings in the Preferences.

Related Topics

- “Preferences” on page 25

Set Zoom

You can set the Spectrum View to display specific frequency (horizontal) and height (vertical) ranges. Displaying specific ranges can be useful for preparing illustrations or focusing on a particular area of the spectrum.

Menu Location	View > Zoom > Set Zoom
Icon(s)	
Hotkey(s)	Ctrl-M Cmd-M

How do I zoom to specific frequency and height ranges?

1. Open the Set Zoom dialog.
2. Enter the horizontal range (X-Axis) in ppm that you would like to display, or select Entire Axis to zoom out to the full spectrum width.
3. Enter the vertical range (Y-Axis) in standard intensity (si) that you would like to display, or select Entire Axis to zoom out to the full spectrum height.
4. Click OK.

Tips and Tricks

- When you open the Set Zoom dialog, the values that appear initially correspond to your current view. If you change the values in only one dimension, your current settings in the other dimension remain intact.

Undo/Redo Zoom

Just as Undo and Redo let you recall recent actions, Undo Zoom and Redo Zoom let you recall recent views of the spectrum.

Menu Location	View > Zoom > Undo/Redo Zoom
Icon(s)	
Hotkey(s)	[or]

Tips and Tricks

- You can access up to twenty recent views through Undo Zoom.
- Redo Zoom is not available until you have used Undo Zoom at least once.

Menu Location	View > Zoom > Auto Zoom
Icon(s)	
Hotkey(s)	Ctrl-F Cmd-F

Auto Zoom

Auto Zoom lets you quickly adjust the vertical scale of the current view so that the tallest peak in the currently displayed frequency range fits in the viewing pane, without changing the displayed frequency range.

Tips and Tricks

- You can apply Auto Zoom to the Spectrum View at any time by clicking anywhere on the spectrum with the middle mouse button.

Menu Location	View > Zoom > Increase/Decrease Vertical Zoom
Icon(s)	
Hotkey(s)	+ (plus) or - (minus)

Increase/Decrease Vertical Zoom

Increase/Decrease Vertical Zoom lets you incrementally adjust the vertical scale of the current view, without changing the displayed frequency range.

Tips and Tricks

- Decrease Vertical Zoom is not available when the full intensity range is visible.

Menu Location	View > Zoom > Increase/Decrease Horizontal Zoom
Icon(s)	
Hotkey(s)	* (asterisk) or / (slash)

Increase/Decrease Horizontal Zoom

Increase/Decrease Horizontal Zoom lets you incrementally adjust the horizontal scale of the current view, without changing the displayed intensity range.

Tips and Tricks

- Decrease Horizontal Zoom is not available when the full frequency range is visible.

Spectrum View

Processor offers a number of useful ways to change the current view by interacting directly with the Spectrum View.

How do I scroll in the Spectrum View?

- Use the horizontal and vertical scroll bars appearing at the bottom and right sides of the Spectrum View.
- Place the mouse cursor anywhere in the Spectrum View, then hold down the **Shift** key and use the scroll wheel on your mouse to scroll horizontally.

How do I zoom in the Spectrum View?

- Click and drag directly on the spectrum. A red zoom box appears to show you the region you are about to zoom in on.
- Select a region using the Select Region tool, and double-click in the region, or right-click in the region and select **Expand Region**. The region expands horizontally to fill the Spectrum View.
- Select a region using the Select Region tool, and middle-click in the region, or right-click in the region and select **Auto Zoom Region**. The region expands horizontally to fill the Spectrum View, and an Auto Zoom is applied.
- Place the mouse cursor anywhere in the Spectrum View, and use the mouse scroll wheel to adjust the Vertical Zoom.
- Place the mouse cursor anywhere in the Spectrum View, press and hold the **Control** key (**Command** for Mac OS X) and use the mouse scroll wheel to adjust the Horizontal Zoom.

- Press + (plus) or - (minus) to adjust the Vertical Zoom.
- Press * (asterisk) or / (slash) to adjust the Horizontal Zoom.
- Middle-click anywhere in the Spectrum View to apply an Auto Zoom.

Spectrum Thumbnail

The Spectrum Thumbnail appears in the bottom right corner of modules with a Spectrum View, and displays a miniature version of the Spectrum View. You can interact with the Spectrum Thumbnail to change the view appearing in the Spectrum View.

How do I toggle the display of the Spectrum Thumbnail?

1. Click **Show Spectrum Thumbnail**.

How do I use the Spectrum Thumbnail?

- Click and drag in the Spectrum Thumbnail. A red zoom box appears to show you the region you are about to zoom in on, and remains to indicate the region currently displayed in the Spectrum View.
- Click and drag inside the zoom box to move it. The Spectrum View updates to show the new region.
- Click and drag one of the solid borders of the zoom box to move that border, or click and drag one of the corners to move the size and shape of the displayed region.
- Click anywhere in the Spectrum Thumbnail outside the zoom box to show the entire spectrum in the Spectrum View. The results are the same as using the Show Entire Spectrum zoom tool.
- Double-click anywhere in the Spectrum Thumbnail to show the entire spectrum in the Spectrum View. The results are the same as using the Show Entire Spectrum tool.

Related Topics

- “Show Entire Spectrum” on page 20

Select Region

You can use the Select Region tool as an alternate zoom function, to specify areas with which to filter the Compound Table, and to obtain information about the selected region itself, including the region's width and the area under various lines within the region.

How do I select a region?

1. Activate Select Region.
2. Click and drag across a region of the Spectrum View to select the region.
3. Double-click in the selected region to expand the region horizontally to fill the Spectrum View.

Or

Middle-click in the selected region to expand the region horizontally to fill the Spectrum View and apply an Auto Zoom.

How do I quickly select a region?

1. Press and hold down the **Shift** key.

Menu Location	View > Show Spectrum Thumbnail
Hotkey(s)	Ctrl-Shift-T Cmd-Shift-T

Menu Location	View > Select Region View > Zoom > Select Region
Icon(s)	
Hotkey(s)	Ctrl-T Cmd-T

- Click in the Spectrum View where you would like to start the selection, then drag to the end point.

Tips and Tricks

- You can also activate Select Region by double-clicking anywhere in the Spectrum View.
- You must have a spectrum, simulation or compound currently displayed to select a region.
- Additional information about the selected region appears in the top right corner of the Spectrum View. The information displayed always includes the region width in Hz and ppm. Depending on the module, the display can also include the areas under the compound line, overlay line, simulation line, spectrum line or sum line within the selected region, as well as the percentage of the selected area fit by the sum line (Profiler only).
- In Profiler, you can use the selected region to filter the contents of the Compound Table via Quick Searches. Right-click on the selected region, then click **Filter for Compounds in Selected Region**.
- The selected region remains in the Spectrum View until you select a new region using the Select Region tool, or expand the region to fill the Spectrum View by double- or middle-clicking inside the region.

Related Topics

- “Auto Zoom” on page 21
- “Compound Table” on page 46
- “Quick Searches” on page 50

4.11 Export Spectrum Image

When you prepare publications or presentations based on your work in Chenomx NMR Suite, you may need to include images illustrating the fitting process or particular regions of the spectrum. You can export spectrum images from any module with a Spectrum View into a variety of bitmap and vector formats.

Menu Location	File > Export > Spectrum Image...
Icon(s)	

How do I export an image of the current spectrum?

- Start the Export Spectrum Image wizard.
- Choose a file format. You can specify bitmap (.gif, .png) or vector formats (.svg, .pdf, .eps).
- Choose a folder and filename for the exported image.

Tips and Tricks

- If you plan to distribute the image via e-mail or post it to a Web page, use Portable Network Graphics (.png), Graphics Interchange Format (.gif) or Scalable Vector Graphics (.svg).
- If you need an image that others can view or print on multiple computers, use Portable Document Format (.pdf).
- If you are preparing images for professional typesetting or publication, use Encapsulated PostScript (.eps).

- Choose **Do not shade area under lines** to ensure that none of the peak shading effects are retained in the exported image. Disabling shading can be useful when you export an image of an area containing many clusters, as shading can create a very heavy look.

Related Topics

- “Spectrum View Tools” on page 20
- “Display Options” on page 24

4.12 Copy Spectrum Image to Clipboard

You can use Copy Spectrum Image to Clipboard to quickly make an image of the current spectrum available to paste directly into another application, without creating an intermediary file.

Menu Location	Edit > Copy Spectrum Image to Clipboard
Icon(s)	

Tips and Tricks

- The resolution of the spectrum image copied to the clipboard varies based on the size of the window on your screen in pixels. If adjusting combinations of your window size and screen resolution do not let you produce images of a high enough quality, consider using Export Spectrum Image.

Related Topics

- “Spectrum View Tools” on page 20
- “Export Spectrum Image” on page 23
- “Display Options” on page 24

4.13 Display Options

There are a number of ways in which you can customize modules in Chenomx NMR Suite to display information as you prefer. You can toggle the display of graph elements in the Spectrum View, and you can change scales. You can find all of the following options in the **View** menu.

Table 4.2. View Menu Options

X Scale	Changes the units represented on the horizontal scale in the Spectrum View. Available units are ppm and Hz.
Y Scale	Changes the units represented on the vertical scale in the Spectrum View. Available units are si (standardized intensity) and ni (normalized intensity).
Sidebar > Compound Sets	Switches the Sidebar to the Compound Sets view. Available in Profiler.
Sidebar > Legend	Switches the Sidebar to the Legend view. Available in Processor, Profiler, Spin Simulator and Compound Builder.
Sidebar > Processing History	Switches the Sidebar to the Processing History view. Available in Processor.
Sidebar > Reference Panel	Switches the Sidebar to the Reference Panel view. Available in Profiler and Compound Builder.

Sidebar > Simulation	Switches the Sidebar to the Simulation view. Available in Spin Simulator.
Show Sidebar	Toggles the display of the Sidebar. Available in every module except Library Manager.
Show Spectrum Thumbnail	Toggles the display of the spectrum thumbnail in modules with a Spectrum View.

Tips and Tricks

- You can change the colors of the various graph elements noted above, as well as that of the background and zoom box, in the Preferences for the appropriate module.

Related Topics

- “Preferences” on page 25

4.14 Preferences

In Preferences, you can specify default folders, default settings for the Spectrum View, colors of various graph elements and some module-specific settings.

Menu Location	File > Preferences
Icon(s)	

How do I change preferences?

- Open the Preferences module.
- Select the module for which you would like to change preferences.
- Make your desired changes.
- Click **OK**.

How do I use a proxy server with Chenomx NMR Suite?

- Switch to the General module in Preferences.
- In the Internet Preferences section, select **Proxy connection to the Internet**.
- Enter the IP address of your proxy server, and the port that Chenomx NMR Suite should use to access it.

How do I turn off animated zooming?

- Switch to the General module in Preferences.
- Clear the check box labeled **Enable Zoom Animation**.

How do I change colors?

- Click the swatch box corresponding to the color you would like to change.
- Click a new color on the color palette.
- Click **OK** to accept the selected color.

How do I set the initial zoom?

- Select a module tab.
- Select the **Initial Zoom** check box.
- In the Initial Zoom fields, enter the default values (in ppm) for the range of the spectrum to display.

How do I specify the default file folder?

1. Select a module tab.
2. Click the **Always use this folder** button.
3. Click the default folder browse button (...).
4. Locate the desired folder.
5. Click **Select** to select the specified folder.

Tips and Tricks

- When you click on a color swatch to change a color, the simple color chooser that appears initially lets you choose from a selection of 15 common, well-distinguished colors. Colors in this selection that are already used in the current module's preferences are indicated with a faint white checkmark, while the currently selected color is indicated with a strong white checkmark.
- When you change a color in Preferences, you are setting a new default color, meaning that the element whose color you set will appear in that color every time you use Chenomx NMR Suite.
- Avoid light line colors with a light (or white) background, as they are difficult to see on screen or in a printout.
- The low and high frequencies for the initial zoom must be at least 1 ppm apart.

Related Topics

- "Spectrum View Tools" on page 20

Table 4.3. Preferences Overview


General Preferences	
Internet Preferences	Lets you configure Chenomx NMR Suite to access the Internet through a proxy server, if needed, or to disable connecting to the Internet entirely.
Zooming Preferences	Lets you change the color of the zoom box, and toggle the zooming animations used when you change zoom settings in the Spectrum View.
Graph Preferences	
Show Spectrum Line	Toggles the default display of the spectrum line in the Spectrum View in Processor and Profiler, and lets you set its default color.
Show Sum Line	Toggles the default display of the sum line in the Spectrum View in Profiler, and lets you set its default color.
Show Compound Line	Toggles the default display of the compound line in the Spectrum View in Compound Builder, and lets you set its default color.
Show Simulation Line	Toggles the default display of the simulation line in the Spectrum View in Spin Simulator, and lets you set its default color.

Show Subtraction Line	Toggles the default display of the subtraction line in the Spectrum View in every module except Library Manager, and lets you set its default color.
Selection Line	Changes the default color of the selection line in the Spectrum View in Profiler, Compound Builder and Spin Simulator.
Overlay Lines	Changes the default color of spectrum overlay lines in the Spectrum View in Profiler, Processor, Compound Builder and Spin Simulator.
CSI Line	Changes the default color of the CSI line in the CSI Editor's Spectrum View in Processor.
Processing Tools	Changes the default color of the processing tools in Processor, including the baseline correction breakpoints, reference deconvolution peak, adjustment triangles, etc.
Initial Zoom (ppm)	Sets the frequency range of the Spectrum View that initially appears when you open a file.
File Preferences	
Display Recent Files	Sets the number of files displayed when you click File , then Open Recent .
Default File Folder	Sets the default folder for saving and opening files. Options are to use the last used folder, or always use the specified folder.

Chapter 5: Processor

The Processor module lets you convert spectra into the Chenomx NMR Suite format (.cnx) and perform basic processing. Currently, Processor can import raw data in Varian (.fid), Bruker (.fid), JEOL (.jdf) and NMRPipe (.fid) formats and processed data in Varian (phasefile), Bruker (1r), JEOL (.jdf), NMRPipe (.ft2) and JCAMP-DX (.jdx) formats. Regardless of the source format, Processor saves spectra as spectrum files (.cnx).

5.1 Overview



- ❶ Sidebar View Buttons. Let you change the currently visible Sidebar view.
- ❷ Processing History. Displays a list of every processing layer that you have applied to the current spectrum, and lets you remove any processing layer that you have applied using the Processor module.
- ❸ Information Panel. Displays information relating to the current stage of spectrum processing. Also, contains the controls used during pH determination, CSI adjustments and other processing steps.
- ❹ Spectrum View Tools. Let you navigate the spectrum more easily.
- ❺ Spectrum View. Displays the spectrum as well as guidelines used during pH determination, CSI adjustments and other processing steps.
- ❻ Spectrum Thumbnail. Displays a thumbnail image of the current spectrum. You can also use the thumbnail to navigate the spectrum at various zoom settings.
- ❼ Status Bar. Displays the current cursor position in the units specified by the **View > X Scale** and **View > Y Scale** settings. Also displays information about recent user actions.

5.2 Quick Reference Guide

	Linux, Windows	Mac OS X
File Menu		
 Open...	Ctrl-O	Cmd-O
 Open Recent	---	---
 Close	Ctrl-W	Cmd-W
 Save	Ctrl-S	Cmd-S
 Save As...	---	---
 Export > Spectrum Image...	---	---
 Spectrum Details...	---	---
 Send to Profiler	Ctrl-J	Cmd-J
 Preferences...	---	---
 Exit	---	---
Edit Menu		
 Undo	Ctrl-Z	Cmd-Z
 Redo	Ctrl-Y	Cmd-Y
 Copy Spectrum Image to Clipboard	---	---
View Menu		
X Scale	---	---
Y Scale	---	---
 Zoom > Show Entire Spectrum	Ctrl-Minus	Cmd-Minus
 Zoom > Undo Zoom	[	[
 Zoom > Redo Zoom	]	]
 Zoom > Set Zoom...	Ctrl-M	Cmd-M
 Zoom > Select Region	Ctrl-T	Cmd-T
 Zoom > Auto Zoom	Ctrl-F	Cmd-F
 Zoom > Increase Vertical Zoom	+	+
 Zoom > Decrease Vertical Zoom	-	-
 Zoom > Increase Horizontal Zoom	*	*

	Linux, Windows	Mac OS X
 Zoom > Decrease Horizontal Zoom	/	/
Edit Mode > CSI	---	---
Edit Mode > pH	---	---
 Sidebar > Legend	Ctrl-Shift-L	Cmd-Shift-L
 Sidebar > Processing History	Ctrl-Shift-H	Cmd-Shift-H
Show Sidebar	Ctrl-Shift-S	Cmd-Shift-S
Show Spectrum Thumbnail	Ctrl-Shift-T	Cmd-Shift-T
Processing Menu		
 Phasing...	Ctrl-Alt-1	Cmd-Alt-1
 Baseline Correction	Ctrl-Alt-2	Cmd-Alt-2
 Line Broadening...	Ctrl-Alt-3	Cmd-Alt-3
 Advanced > Reference Deconvolution...	Ctrl-Alt-4	Cmd-Alt-4
 Advanced > Water Deletion	Ctrl-Alt-5	Cmd-Alt-5
 Advanced > Reverse Spectrum	Ctrl-Alt-6	Cmd-Alt-6
Tools Menu		
 Add Spectrum Overlays...	Ctrl-Shift-O	Cmd-Shift-O
 Remove All Spectrum Overlays	---	---
 Find CSI Automatically	---	---
 Find pH Automatically	---	---
 Convert Spectra...	---	---
Help Menu		
 NMR Suite Help...	F1	F1
 Save Application Log...	Ctrl-Shift-A	Cmd-Shift-A
 Software Activation...	---	---
 About Chenomx NMR Suite...	---	---

5.3 Information Tools

In addition to processing the spectrum itself, you can add extra information to a spectrum that lets you analyze the spectrum in Profiler or create compound signatures and simulations

in Compound Builder and Spin Simulator. You can also review processing layers that you have added, and remove them if necessary.

CSI Editor

The chemical shape indicator (CSI) settings form the basis for all fits performed in Profiler and Compound Builder. The accuracy of these settings also have a direct impact on the accuracy of profiled spectra (.cnx), exported data (.txt), compound files (.xcpd) and spin simulations (.xss). The CSI editor lets you view, define and change these settings for the current spectrum.

When you import a spectrum, Processor attempts to automatically locate the CSI, and determine its width and height. You can change these automatically determined parameters using the CSI editor if necessary.

Menu Location	View > Edit Mode > CSI
---------------	------------------------

Menu Location	Tools > Find CSI Automatically
Icon(s)	

How do I modify the CSI settings?

1. Switch to the CSI Editor.
2. Select the CSI present in your sample using the **Chemical Shape Indicator** drop-down box. Valid CSIs for use in Chenomx NMR Suite include DSS, TSP and formate.
3. Type the concentration of the CSI in your sample in the **Concentration (mM)** text entry box.
4. *Automatic:* Click the **Find Automatically** button.
5. *Manual:*
 - Click and drag the CSI peak to match it with the spectrum line (black).
 - Adjust the width by dragging the width control (red triangles to either side of the peak) or typing a number directly into the **Width (Hz)** text entry box.
 - Use the up and down arrow keys to adjust peak height, and the left and right arrow keys to adjust peak position.

Tips and Tricks

- If Processor cannot find the CSI automatically when you use DSS or TSP, check the phasing and baseline correction of your spectrum and try again. If Processor still cannot find your CSI, you need to use the manual tools to locate it instead.
- To adjust the CSI settings for your spectrum, Processor must be in the CSI Edit Mode.
- Processor cannot automatically find formate as a CSI, due to potential overlap with other compounds.
- If you have selected DSS as the CSI, then the indicated center position always reads 0 ppm. If you use TSP or formate instead, the center position varies based on the current pH value for the spectrum.
- If you are using formate as a CSI, finding the CSI automatically does not work.
- If the automatic method does not work, you need to define the CSI manually.
- If you would like to adjust the height or width of the CSI peak while keeping its total area constant, check the **Preserve CSI Area** box.
- Define CSI parameters in the CSI Editor only after you have finished applying all necessary processing layers. In particular, if you change the line broadening or add

reference deconvolution, you need to set new CSI parameters to ensure accurate quantification and generally good fits.

- Finding the DSS or TSP peak automatically relies on having the water peak centered in the spectrum. The formate peak cannot be located automatically.

Related Topics

- Chapter 2, *Handling Samples and Spectra* on page 4
- “Importing Spectra” on page 13
- “pH Editor” on page 32
- “Spectrum Details” on page 34
- “Processing Tools” on page 35

pH Editor

Signal position and shape for compounds in aqueous solution can be very sensitive to pH. The accuracy of the pH that you set for a spectrum may have an impact on the quality and accuracy of fits that you may perform in Profiler, Compound Builder or Spin Simulator based on the spectrum. The pH editor lets you calculate or set the pH for the current spectrum.

If your sample contains any combination of creatinine, imidazole or DFTMP, Processor may be able to calculate the spectrum pH. If your sample does not contain these compounds, or if you prefer to use an externally measured pH, you may enter a pH value manually.

While importing a spectrum, Processor attempts to calculate the pH of the sample automatically using the frequencies of the imidazole, creatinine and DFTMP peaks. You can refine or change this calculated pH in the pH Editor. If Processor cannot calculate the pH of your sample, you need to set it manually.

Whether you use automatic or manual methods to set the spectrum pH, you can also set an uncertainty for the pH value. The pH uncertainty helps to determine the size of the transform windows for pH-sensitive signatures in Profiler.

$$pH = pKa + \log_{10} \left(\frac{\delta_{min} - \delta}{\delta - \delta_{max}} \right)$$

Figure 5.1. Equation to calculate spectrum pH from the position of a single peak. pKa is the acid dissociation constant of the indicator, δ is the chemical shift of the peak, and δ_{max} and δ_{min} are the upper and lower bounds on the chemical shift.

Processor calculates the spectrum pH as a function of the peak positions of the indicator compounds. The formula used for DFTMP calculates spectrum pH based on the observed position of a single peak (see Figure 5.1 on page 32), while the formula used for both imidazole and creatinine uses the distance between the two observed peaks (see Figure 5.2 on page 32).

$$pH = pKa + ml + \log_{10} \left(\frac{\Delta_{min} - \Delta\delta}{\Delta\delta - \Delta_{max}} \right)$$

Figure 5.2. Equation to calculate spectrum pH from the positions of two peaks. pKa is the acid dissociation constant of the indicator, ml is a correction for ionic strength, $\Delta\delta$ is the chemical shift difference between the peaks, and Δ_{max} and Δ_{min} are the upper and lower bounds on the chemical shift difference.

The single peak equation for DFTMP derives from a similar equation published recently (see citation below). The two peak equation was developed internally by Chenomx.

MD Reily, LC Robosky, ML Manning, A Butler, JD Baker, and RT Winters. **2006.** *DFTMP, an NMR Reagent for Assessing the Near-Neutral pH of Biological Samples.* *J Am Chem Soc.* 128(38): 12360-12361

The predicted frequencies of the imidazole, creatinine and DFTMP peaks, based on the pH value that you currently have set, appear in the Spectrum View as a pair of blue lines at high frequency (imidazole), a pair of green lines at mid frequency (creatinine) and a single orange line at low frequency (DFTMP). You can use these lines during manual pH determination by aligning them with the corresponding peaks for creatinine, imidazole and DFTMP, when these compounds are present.

How do I set the spectrum pH?

1. Switch to the pH Editor.
2. Set the desired pH using any one of the following methods:
 - Click the **Find Automatically** button.
 - On the **Tools** menu, click **Find pH Automatically**.
 - Click and drag the pH slider bar to set the determined frequencies of imidazole, creatinine and DFTMP as close as possible to the actual frequencies of the peaks in the spectrum (if known). The resulting pH appears in the text box beside the slider bar.
 - Select the pH slider control, and use the arrow keys to make fine adjustments to the pH in increments of 0.01.
 - Type the desired numeric pH into the text entry box to the right of the slider bar.

Tips and Tricks

- Setting an accurate spectrum pH helps Profiler choose more accurate starting positions and transform windows for clusters.
- When aligning the displayed lines to peaks in your spectrum, it is not necessary to have the lines match the peak positions precisely. A common situation is to arrive at two pH values separated by 0.01, one of which places the indicator lines on one side of the peak, while the other places them on the other side of the peak. In such cases, you may use either value.
- Once you have made coarse adjustments to the pH, you can zoom in to particular indicator lines to make finer adjustments. Select the **Fine Tune** check box to let you make more detailed adjustments over a smaller pH range.
- If you have selected TSP or formate as the CSI, then changing the pH changes the position of the horizontal axis, since the frequencies of the TSP and formate peaks are pH-sensitive.
- Processor must be in the pH Edit Mode to adjust the pH for your spectrum.
- You must save your spectrum to retain the changes that you have made to the spectrum pH.

Related Topics

- Chapter 2, *Handling Samples and Spectra* on page 4
- “Importing Spectra” on page 13
- “CSI Editor” on page 31

Menu Location	View > Edit Mode > pH
---------------	-----------------------

Menu Location	Tools > Find pH Automatically
Icon(s)	

- “Spectrum Details” on page 34

Spectrum Details

Spectrum details describe the acquisition parameters for the spectrum. Your selected CSI and its concentration, along with the spectrum pH that you set, also appear in the Spectrum Details. Since most of these parameters come directly from the imported spectrum, you cannot modify most of the spectrum details, except for the **Frequency Domain Size** and **Comment** fields.

The **Frequency Domain Size** field displays the number of points in the Chenomx spectrum (.cnx) in the frequency domain, after spectrum data has been imported and zero-filling applied. You can contrast this with the **Time Domain Size**, which displays the number of point originally present in the imported spectrum, in the time domain.

The **Comment** field initially consists of the contents of any text descriptor file associated with the original spectrum, but you can replace or expand it as desired.

Menu Location	File > Spectrum Details
Icon(s)	

How do I view and edit details of the current spectrum?

1. Open Spectrum Details.
2. Click anywhere in the **Comment** field to place a text cursor in the field.
3. Make any desired changes.
4. Click the **Accept** button to apply your changes or click the **Cancel** button to discard them.

Tips and Tricks

- You must save the spectrum to retain changes that you have made to the spectrum comments.
- You can only change the CSI settings and spectrum pH using the CSI and pH editors.
- The value of the **Frequency Domain Size** field is set by the amount of zero-filling that you select when you convert a spectrum. If you selected "Automatic," the frequency domain size will be the nearest power of two that is at least twice the time domain size.

Related Topics

- Chapter 2, *Handling Samples and Spectra* on page 4
- “Importing Spectra” on page 13
- “CSI Editor” on page 31
- “pH Editor” on page 32

Processing History (Sidebar)

Processing spectra in Processor does not alter your original spectral data. Every processing effect that you add to a spectrum acts as a layer modifying the display of your original data. All of these effects are completely reversible, and are tracked in the Processing History Sidebar view.

While processing a spectrum, you may want to review the processing layers that you have already applied, or you may need to reverse the effects of one or more of the layers that you have already applied. The Processing History maintains a list of every processing layer

that you have applied to the current spectrum using Processor. You can review this list to find out the parameters used in each layer, like phase angles or reference deconvolution parameters, and step back through past layers to reverse their effects.

When you import a processed spectrum, Processor imports phasing, baseline correction and other processing that you may have applied in other processing software. The imported processing effects appear in the Processing History as a layer called Preprocessing. The Preprocessing details, appearing at the bottom of the panel, indicate the format from which you imported the spectrum. The Preprocessing layer is not reversible; you cannot remove it.

How do I review and remove processing layers?

1. Switch to the Processing History view in the Sidebar.
2. Click the name of a processing layer to review its associated parameters; they appear in the Details box below the list.
3. Click the **Clear Last** button to remove the most recent processing layer. Click the button again to remove the next most recent layer, and so on.

Or

Click the **Clear All** button to remove all of the processing layers listed in the Processing History.

Tips and Tricks

- The Processing History only tracks processing effects, including phasing, baseline correction, reference deconvolution, water deletion and line broadening. Changes to spectrum metadata, such as the spectrum details or CSI and pH settings, do not appear in this panel.
- If you have just accidentally removed a layer in the Processing History, you can use the **Undo** button to restore it.

Related Topics

- Chapter 2, *Handling Samples and Spectra* on page 4
- “Importing Spectra” on page 13
- “Legend” on page 17
- “Processing Tools” on page 35
- “Convert Spectra” on page 41

5.4 Processing Tools

Processor applies each processing step as a layer, altering the way that you see and interact with your spectrum, but leaving the original raw data untouched. The Processing History records each layer along with its associated parameters, like phase angles, baseline breakpoints and linewidth adjustments. You can review the layers, and even remove them if necessary, restoring your spectrum to earlier state(s).

When you import a processed spectrum, processing that you have performed in other software appear in the Processing History as a layer called Preprocessing; the details of the Preprocessing layer indicate the format from which the processed data was imported. The Preprocessing layer is not reversible, and cannot be removed. If you need to work with the

Menu Location	View > Sidebar > Processing History
Icon(s)	
Hotkey(s)	Ctrl-Shift-H Cmd-Shift-H

spectrum as it appeared before processing in other software, import the raw data (fid, .jdf, .fid) directly.

Phasing

Phase correction, or phasing, lets you correct phase shifts that may have occurred during data acquisition. You can recognize phase shifts in a spectrum as an asymmetric baseline on either side of peaks and peak clusters in the spectrum. In some cases, you may see whole peaks inverted, pointing down from the baseline instead of up, while extreme phase shifts may also add a periodic oscillation or "rolling" to the baseline.

Most spectra require some degree of phase correction. Adjust phasing before applying any other processing layers to your spectrum, like baseline correction or reference deconvolution.

While importing a spectrum, you may choose to have Processor automatically apply phasing adjustments. After importing a spectrum, you can apply automatic or manual phase correction. In either case, a Phasing layer appears in the Processing History; you can remove the layer if necessary.

How do I apply phase correction to my spectrum?

1. Switch to Phasing mode.
2. *Automatic:* Click the **Auto** button to have Processor automatically calculate phase angles. Review the spectrum after automatic phasing to ensure that the results are acceptable. If necessary, you can use the manual controls to adjust the phase angles before accepting them.
3. *Manual:* Click and drag the **Zero-order Phase** and **First-order Phase** sliders, or type angles in the boxes beside each slider. Use the **Normal**, **Fine** or **Very Fine** modes to select the scale of changes made with the sliders.
4. Click the **Accept** button to apply your changes, or click the **Cancel** button to discard them.

Tips and Tricks

- Only phase your spectrum once. If you do not like the results of phasing, undo it or remove it using the Processing History, then try again.
- Always phase your spectrum before applying baseline correction or any other processing layers.
- Select one of the phase angle controls, and use the arrow keys to make fine adjustments to the phase angle. The arrow keys adjust in increments of 1° in normal mode, 0.1° in fine mode, and 0.01° in very fine mode.
- Properly phase all of your spectra, as good phasing is essential to accurately analyzing your data using Chenomx NMR Suite.
- Your spectrum is properly phased when it has a baseline that describes a smooth curve across the entire width of the spectrum.
- In spectra of aqueous samples, the water peak may be distorted relative to the rest of the spectrum; do not use the water peak to determine phasing. Instead, try to obtain smooth baseline curves to either side of the water peak.
- You will usually need to adjust the zero-order phase angle more than the first-order phase angle, but remember that the requirements of a specific spectrum may vary.

Menu Location	Processing > Phasing
Icon(s)	
Hotkey(s)	Ctrl-Alt-1 Cmd-Alt-1

- During a phasing session, you can clear your current adjustments by clicking the **Reset** button. After you have accepted your changes, you can clear previous adjustments through the Processing History panel.

Related Topics

- “Importing Spectra” on page 13
- “Processing History (Sidebar)” on page 34
- “Baseline Correction” on page 37
- “Reference Deconvolution” on page 39
- “Convert Spectra” on page 41

Baseline Correction

Baseline correction lets you remove distortions in the baseline of the spectrum. You can recognize distorted baselines as a net curvature or slant of the regions of your spectrum that contain no signal (“noise” regions). Common forms of baseline distortions include “smiles” (outer regions of the spectrum turned up), “frowns” (outer regions of the spectrum turned down) and “wings” (left and right edges of the spectrum turned up or down, with the middle region flat). If you notice a periodic oscillation or “rolling” of the baseline, check the phasing of your spectrum before proceeding with baseline correction.

Most spectra require some degree of baseline correction, particularly if you need accurate quantification in Profiler. Apply baseline correction after phasing, but before any other processing layers, like reference deconvolution.

While importing a spectrum, you may choose to have Processor automatically apply either simple linear baseline correction (drift correction) or a cubic spline baseline correction based on automatically-determined breakpoints. After importing a spectrum, you can apply automatic or manual baseline correction. In either case, a Baseline Correction layer appears in the Processing History panel; you can remove this layer if necessary.

For more direct control over the baseline correction applied to your spectrum, you can specify endpoints (linear fit) or breakpoints (cubic spline fit) manually.

How do I apply baseline correction to my spectrum?

1. Switch to Baseline Correction mode.
2. *Automatic (Linear)*: Click the **Auto Linear** button to automatically determine a linear baseline adjustment.
Or
Automatic (Spline): Click the **Auto Spline** button to automatically determine a set of breakpoints for a cubic-spline-based baseline adjustment.
3. *Manual*: Click on the spectrum to set breakpoints for a cubic spline fit. Click and drag existing breakpoints (including the endpoints) to move them. Hold down **Control** or **Command** to place breakpoints anywhere in the spectrum, instead of directly on the spectrum line itself. Hold down **Shift** while clicking on an existing breakpoint to delete it.
4. Click the **Accept** button to apply your changes, or click the **Cancel** button to discard them.

Menu Location	Processing > Baseline Correction
Icon(s)	
Hotkey(s)	Ctrl-Alt-2 Cmd-Alt-2

Tips and Tricks

- Only apply baseline correction after you are satisfied with the phasing in your spectrum.
- Even if you intend to manually adjust baseline correction, you may find it helpful to use the **Auto Linear** function first. **Auto Linear** lets Processor first locate the "natural" endpoints of the spectrum, giving you a useful starting point for any more complex baseline fit.
- When setting breakpoints, be aware of the amount of curvature in the area you are trying to fit. You need to set more breakpoints to fit regions with very strong curvature, such as the areas near the water peak, than you would for a relatively flat region.

Related Topics

- "Importing Spectra" on page 13
- "Processing History (Sidebar)" on page 34
- "Phasing" on page 36
- "Reference Deconvolution" on page 39
- "Convert Spectra" on page 41

Line Broadening

Line broadening is a mathematical operation that multiplies your fid by an exponential function before the fid is Fourier transformed. Line broadening effectively increases the linewidth in your spectrum while averaging out instrument noise.

You may not need to apply line broadening to every spectrum. If you choose to use line broadening, apply it after phasing and baseline correction.

While importing a spectrum, you may choose to have Processor automatically apply line broadening of 0 to 1000 Hz. After importing a spectrum, you can apply line broadening manually. In either case, a Line Broadening layer appears in the Processing History; you can remove the layer if necessary.

How do I apply line broadening to my spectrum?

1. Switch to Line Broadening mode.
2. Enter a value between 0 and 1000 Hz.
3. Click the **Accept** button.

Tips and Tricks

- Only apply line broadening after you are satisfied with the phasing and baseline correction in your spectrum.
- If you intend to apply reference deconvolution to your spectrum, you can add line broadening at the same time. See "Reference Deconvolution" on page 39 for details.

Related Topics

- "Importing Spectra" on page 13
- "Processing History (Sidebar)" on page 34
- "Phasing" on page 36

Menu Location	Processing > Line Broadening...
Icon(s)	
Hotkey(s)	Ctrl-Alt-3 Cmd-Alt-3

- “Baseline Correction” on page 37
- “Reference Deconvolution” on page 39
- “Convert Spectra” on page 41

Reference Deconvolution

Reference deconvolution is a method of reconstructing an ideal spectrum by removing lineshape distortions based on the shape of a reference peak. To successfully apply reference deconvolution to a spectrum, you need to find a single, isolated peak, well-separated from the signals in the rest of the spectrum. If you have added DSS or TSP to your sample, the methyl signal from either of these is an ideal choice.

Reference deconvolution is a linear process, using direct and inverse Fourier transforms, so quantitative relationships among compounds in the spectrum remain intact. Applying reference deconvolution to a spectrum assumes that all lineshapes in the spectrum are affected by systematic distortions in the same way.

Many spectra can benefit from the use of reference deconvolution, but it is not required. It is particularly useful in reducing the effects of varying shimming techniques, or when you need accurate quantification from regions with strong overlap of signals.

You need to apply reference deconvolution manually after you have imported your spectrum. After you have applied reference deconvolution, a Reference Deconvolution layer appears in the Processing History panel; you can remove this layer if necessary.

How do I apply reference deconvolution to my spectrum?

1. In the CSI editor, fit the red CSI peak to your spectrum. See “CSI Editor” on page 31 for details.
2. Switch to Reference Deconvolution mode.
3. Zoom in to the highlighted peak in the Spectrum View.
4. If you are using DSS or TSP as a CSI, select the **Use DSS/TSP Satellites** check box. Otherwise, clear the check box.
5. Enter the desired target linewidth.
6. Click the **Accept** button to apply your changes, or click the **Cancel** button to discard them.
7. In CSI Edit mode, refit the red CSI peak to your spectrum to account for changes introduced by reference deconvolution. See “CSI Editor” on page 31 for details.

Tips and Tricks

- Only apply reference deconvolution after you are satisfied with the phasing and baseline correction in your spectrum.
- The target linewidth text box initially contains the native linewidth of your spectrum, based on your CSI definition.
- Choosing a target linewidth larger than the native linewidth in your spectrum is equivalent to applying line broadening. Choosing a value smaller than the native linewidth may provide some resolution enhancement, at the expense of a reduced signal-to-noise ratio and possible spectral artifacts.

Menu Location	Processing > Advanced > Reference Deconvolution
Icon(s)	
Hotkey(s)	Ctrl-Alt-4 Cmd-Alt-4

Related Topics

- “Importing Spectra” on page 13
- “Processing History (Sidebar)” on page 34
- “Phasing” on page 36
- “Baseline Correction” on page 37
- “Line Broadening” on page 38
- “Convert Spectra” on page 41

Water Deletion

Menu Location	Processing > Advanced > Water Deletion
Icon(s)	
Hotkey(s)	Ctrl-Alt-5 Cmd-Alt-5

In many spectra of biofluids or cell extracts, the water peak may be very large. In some cases, it is so large that seeing any other peaks while completely zoomed out is difficult. Water deletion is a simple processing layer that masks data points in the frequency domain across a range associated with water. It is designed to prevent the height of the water peak from dominating the vertical range of the spectrum display.

Using water deletion is optional. Ensure that no peaks of interest occur close to the water peak before you apply water deletion, as such peaks may be masked by the water deletion layer. Apply water deletion after any other processing layers, including phasing, baseline correction and reference deconvolution.

While importing a spectrum, you may choose to have Processor apply water deletion automatically. After importing a spectrum, you can apply water deletion manually. In either case, a Water Deletion layer appears in the Processing History panel; you can remove this layer if necessary.

Tips and Tricks

- Water deletion will not work if you have not defined the CSI for the spectrum, as the location of the water peak is determined relative to your CSI.
- If your spectrum needs to be reversed, but is not yet reversed, water deletion may not work as intended. If necessary, reverse the spectrum before applying water deletion.

Related Topics

- “Importing Spectra” on page 13
- “CSI Editor” on page 31
- “Processing History (Sidebar)” on page 34
- “Reverse Spectrum” on page 41
- “Convert Spectra” on page 41

Reverse Spectrum

Menu Location	Processing > Advanced > Reverse Spectrum
Icon(s)	
Hotkey(s)	Ctrl-Alt-6 Cmd-Alt-6

Reverse spectrum reverses the order of the data points in your spectrum in the frequency domain, effectively flipping the spectrum horizontally.

You only need to use reverse spectrum with certain types of spectra from Bruker spectrometers. If you have trouble locating the CSI peak (when you have a CSI in the sample), or even making any sense of your spectrum after importing it, you may need to reverse the spectrum.

While importing a spectrum, you may choose to have Processor reverse the spectrum automatically. After importing a spectrum, you can reverse the spectrum manually. In either case, a Reverse Spectrum layer appears in the Processing History panel; you can remove this layer if necessary.

Tips and Tricks

- If you notice that a spectrum from a particular spectrometer must be reversed, it is likely that all spectra from that instrument, especially those acquired with the same parameters, need to be reversed.
- If it is necessary, apply Reverse Spectrum first, as other processing layers may be difficult to apply correctly when the spectrum is backwards.

Related Topics

- “Importing Spectra” on page 13
- “CSI Editor” on page 31
- “Processing History (Sidebar)” on page 34
- “Convert Spectra” on page 41

5.5 Importing and Exporting Data

Processor can help you to efficiently import raw spectral data from various sources, convert them into Chenomx spectra (.cnx) and begin working with them.

Convert Spectra

For larger projects, you may want to convert groups of spectra into Chenomx spectra (.cnx), processing them with a consistent set of parameters to ensure consistent analysis results.

How do I convert a group of spectra?

1. Start the Convert Spectra wizard.
2. Choose the source spectra to convert. The source spectra must be in one of the supported file formats (“Spectral Data Processing” on page 5)
3. If the source spectra include multiple file types, select the type that you would like to convert.
4. Click **Conversion Settings** to choose processing options to apply to each spectrum.
5. Specify a target folder for the converted spectra.

Tips and Tricks

- If the source folder contains more than one type of spectrum file, you need to select one type to convert in the current session, since different conversion options may

Menu Location	Tools > Convert Spectra...
Icon(s)	

apply to each available file type. You can do additional sessions to convert other file types.

- The Conversion Settings available through the Convert Spectra wizard are the same as those available when you import individual spectra (“Importing Spectra” on page 13).
- Storing the converted spectra next to the source files means that each new spectrum file appears in the same folder as the spectrum file from which it was created. You may find it more useful to have Processor place the converted spectra in a different folder.
- If you specify a folder other than the source folder for converted spectra, you get a folder containing all of the converted spectra, with no subfolders created. Select **Mirror source folder structure** to create a folder structure in the target folder that is the same as that in the source folder, and store the resulting converted spectra within the new folder structure.

Related Topics

- “Importing Spectra” on page 13
- “Processing History (Sidebar)” on page 34
- “Processing Tools” on page 35

Send to Profiler

When you finish processing a spectrum in Processor, your next task is often to start profiling the spectrum in Profiler. Send to Profiler closes the current spectrum in Processor, and opens it in Profiler, retaining any changes you have made in Processor.

Menu Location	File > Send to Profiler
Icon(s)	
Hotkey(s)	Ctrl-J Cmd-J

Tips and Tricks

- Any changes that you have made to the spectrum in Processor are retained when you send the spectrum to Profiler, and are visible if you bring the spectrum back into Processor. However, you still need to save the spectrum, either in Processor or Profiler, to permanently save those changes to the spectrum file.


Related Topics

- Chapter 6, *Profiler* on page 43

Chapter 6: Profiler

The Profiler module helps you to identify and quantify compounds in a sample spectrum, as well as apply spectral analysis techniques like spectral binning. You can fit any available compound to the current spectrum, and you can modify the compounds displayed using compound sets and filters.

6.1 Overview



- 1 Sidebar View Buttons. Let you change the currently visible Sidebar view.
- 2 Legend. Sidebar view letting you customize the lines appearing in the Spectrum View.
- 3 Cluster Navigator. Lets you navigate among the clusters of the currently selected compound. You can enable or disable the Cluster Navigator using the **View** menu.
- 4 Compound Table. Displays compounds, their calculated concentrations, and their fitting status. You can access signatures of these compounds in the Spectrum View by selecting them. As well, you can view additional compound information in the Reference Panel by double-clicking on a compound name in this table. You can change the default concentration units displayed in the Compound Table using the **View** menu.
- 5 Spectrum View Tools. Let you navigate the spectrum more easily.
- 6 Spectrum View. Displays the spectrum as well as the sum line. You can modify the display by focusing in on a particular area or displaying and modifying compound signatures.
- 7 Quick Search. Lets you filter the Compound Table to focus on a particular compound or group of compounds.

- ⑧ Spectrum Thumbnail. Displays a thumbnail image of the current spectrum. You can also use the thumbnail to navigate the Spectrum View at various zoom settings.
- ⑨ Status Bar. Displays the current cursor position in the units specified by the **View > X Scale** and **View > Y Scale** settings. Also displays information about recent user actions.

6.2 Quick Reference Guide

	Linux, Windows	Mac OS X
File Menu		
 Open...	Ctrl-O	Cmd-O
 Open Recent	---	---
 Close	Ctrl-W	Cmd-W
 Save	Ctrl-S	Cmd-S
 Save As...	---	---
 Send to Processor	Ctrl-J	Cmd-J
 Open Library Manager	Ctrl-B	Cmd-B
 Import Profile...	---	---
 Export > Profiled Data (.txt)...	---	---
 Export > Profiled Compounds (.pack)...	---	---
 Export > Spectrum Image...	---	---
 Spectrum Details...	---	---
 Preferences...	---	---
 Exit	---	---
Edit Menu		
 Undo	Ctrl-Z	Cmd-Z
 Redo	Ctrl-Y	Cmd-Y
 Copy Spectrum Image to Clipboard	---	---
Enforce Transform Windows	Ctrl-Shift-X	Cmd-Shift-X
View Menu		
Concentration Units	---	---
X Scale	---	---
Y Scale	---	---

		Linux, Windows	Mac OS X
	Zoom > Show Entire Spectrum	Ctrl-Minus	Cmd-Minus
	Zoom > Undo Zoom	[	[
	Zoom > Redo Zoom	]	]
	Zoom > Set Zoom...	Ctrl-M	Cmd-M
	Zoom > Select Region	Ctrl-T	Cmd-T
	Zoom > Auto Zoom	Ctrl-F	Cmd-F
	Zoom > Increase Vertical Zoom	+	+
	Zoom > Decrease Vertical Zoom	-	-
	Zoom > Increase Horizontal Zoom	*	*
	Zoom > Decrease Horizontal Zoom	/	/
	Sidebar > Legend	Ctrl-Shift-L	Cmd-Shift-L
	Sidebar > Reference Panel	Ctrl-Shift-R	Cmd-Shift-R
	Sidebar > Compound Sets	Ctrl-Shift-K	Cmd-Shift-K
	Show Sidebar	Ctrl-Shift-S	Cmd-Shift-S
	Show Spectrum Thumbnail	Ctrl-Shift-T	Cmd-Shift-T
	Show Cluster Navigator	Ctrl-Shift-N	Cmd-Shift-N
	Change Columns...	---	---
Compound Menu			
Selected Compound Submenu			
	Clear > Concentration	---	---
	Clear > Transforms	---	---
	Clear > Custom Color	---	---
	Clear > Pin	---	---
	Clear > Star	---	---
	Clear > Everything	Delete	Delete
	Fit Automatically...	Space	Space
	Set to Maximum Concentration	---	---
	Previous Cluster	Q	Q

	Linux, Windows	Mac OS X
 Next Cluster	E	E
Starred Compounds and All Compounds Submenus		
 Clear > Concentrations	---	---
 Clear > Transforms	---	---
 Clear > Custom Colors	---	---
 Clear > Pins	---	---
 Clear > Stars	---	---
 Clear > Everything	---	---
 Set to Maximum Concentrations	---	---
 Scale Concentrations...	---	---
Tools Menu		
 Add Spectrum Overlays...	Ctrl-Shift-O	Cmd-Shift-O
 Remove All Spectrum Overlays	---	---
 Filter Compounds...	---	---
 Spectral Binning...	---	---
Help Menu		
 NMR Suite Help...	F1	F1
 Save Application Log...	Ctrl-Shift-A	Cmd-Shift-A
 Software Activation...	---	---
 About Chenomx NMR Suite...	---	---

6.3 Profiling Tools

Profiler provides you with a variety of tools to help you extract compound data from your spectra. The data extraction process, called profiling, includes selecting compounds from your Library, setting concentrations, determining transforms, and finding and using extended information about compounds that you would like to profile.

Compound Table

The Compound Table summarizes the results of your profiling work with the current spectrum. It contains entries for compounds that are available to profile in your spectrum, as well as those that you have already profiled, along with their maximum and profiled concentrations. You can modify the display of the Compound Table to show as much or as little information as you like, to better suit your preferred approach to profiling.

The currently selected compound in the Compound Table is the active compound in the Spectrum View, letting you see its preview line and adjust its concentration and transforms.

You can filter the current contents of the Compound Table using Quick Searches. Filtering can help you quickly find the exact compound that you would like to profile. You can filter the table based on text search, regions of the spectrum, or specific locations in the spectrum.

Using the Compound Table, you can set custom display colors for compounds, pin compounds so that they are visible even when you are not working on them and toggle stars on compounds to generally mark them for future reference.

There are many columns available for display in the Compound Table. You can choose to display any number of the available columns, with a minimum of one visible at any time, and you can drag the column headers to position the columns in any order that you like. You can sort the Compound Table on any column by clicking its column header. Click the header again to reverse the sort order.

Tips and Tricks

- You must select at least one compound set in the Compound Sets view to see any compounds in the Compound Table.
- When you have nothing selected in the Compound Table, pressing the up or down arrow keys will select the first compound in the Compound Table.
- After typing in the Quick Search box, press the down arrow key to navigate to the first compound in the Compound Table.
- When you select a compound without a concentration in the Compound Table, a preview line for the compound appears in the Spectrum View, scaled vertically to suit your current zoom level. You can click on the preview line to set the concentration for the selected compound to match the height of the preview line.
- When you select a compound with a status icon, the status area to the left of the Quick Search box shows a detailed status message. You can also see the detailed status message in the tooltip that appears when you hover the mouse cursor over the status icon.
- When you hover the mouse cursor over some parts of the Compound Table, tooltips appear containing more details. For example, hovering over the optional matched clusters column for a compound shows you which clusters are matched.
- You can set compound concentrations directly in the Compound Table, by double-clicking on a compound's concentration column and typing in a value.
- You can perform many common compound operations through the context menu, available by right-clicking on a compound's row in the table. Operations accessible through the context menu include setting the compound to its maximum concentration, automatically fitting it to the spectrum, and clearing concentrations, transforms, custom colors, stars and pins.
- You can copy the contents of a cell to the Clipboard by right-clicking on the cell and selecting **Copy [cell name] to Clipboard**. You can then paste the contents as text to any other program.
- Compounds that you have modified from their default state (including concentrations, transforms, custom colors, pins and stars) will appear in the Profiled Compounds set, and appear in bold in the Compound Table.
- Press **Delete** with a compound selected in the Compound Table to clear all concentrations, transforms, custom colors, pins and stars for the selected compound.

- To make more space for the Compound Table, you can hide the Spectrum Thumbnail by pressing **Ctrl-Shift-T**, or clicking **View > Show Spectrum Thumbnail**.

Related Topics

- “Change Columns” on page 16
- “Sidebar” on page 16
- “Preferences” on page 25
- “Custom Colors” on page 48
- “Pinned Compounds” on page 49
- “Starred Compounds” on page 50
- “Quick Searches” on page 50
- “Concentrations” on page 52
- “Maximum Concentrations” on page 56
- “Compound Set Tools” on page 59
- “Export Profiled Data” on page 60
- “Information Tools” on page 64
- “Concentration Units” on page 66

Custom Colors

When working with multiple compounds, it can become difficult to distinguish which compounds are contributing which parts to an overlapped region. You can set specific compounds to appear in colors of your choice whenever they are visible in the Spectrum View.

How do I set a custom color for a compound?

1. Select the compound in the Compound Table.
2. Click in the Custom Color column (🌈) to access the color chooser.
3. Click one of the preconfigured colors in the color chooser, or click **More...** to choose from a larger selection of colors.

How do I clear a custom color from a compound?

1. Select the compound in the Compound Table.
2. Click **Clear Custom Color**, or click in the Custom Color column (🌈) and click **Clear Custom Color**.

Menu Location	Compound > Selected Compound > Clear > Custom Color
Icon(s)	

Tips and Tricks

- You can also set or clear custom colors for compounds visible in the Legend, by clicking on the colored box beside the compound in the Legend.
- Setting a custom color for a compound adds it to the Profiled Compounds set, even if you have not set concentrations or transforms for the compound.
- Custom colors are stored with the spectrum when you save it as a .cnx file. The colors you select are still associated with the compounds when you next open the file, until you clear them.

- You can export custom colors to a file using Export Profiled Data.

Related Topics

- “Legend” on page 17
- “Preferences” on page 25
- “Compound Table” on page 46
- “Pinned Compounds” on page 49
- “Starred Compounds” on page 50
- “Legend (Sidebar)” on page 64

Pinned Compounds

It is often useful to display several different compounds at the same time while profiling, especially when you are working in heavily overlapped area. Pinning a compound ensures that it remains visible regardless of which compound you have selected in the Compound Table. In combination with custom colors and exporting spectrum graphics, pinning can also help you to prepare snapshots of your profile that are suitable for publication.

How do I pin a compound?

- Click the pin for the compound in the Pinned column (📌) of the Compound Table.
- Click the pin next to the selection line in the Legend (📌).

How do I clear a pinned compound?

1. Select the compound in the Compound Table.
2. In the Compound menu, click **Selected Compound > Clear > Pin**, clear the check box for the compound in the Pinned column, click the red pin next to the selection line in the Legend, or click the red X beside the compound's entry in the Legend.

Menu Location	Compound > Selected Compound > Clear > Pin
Icon(s)	

Tips and Tricks

- Pinning a compound adds it to the Profiled Compounds set, even if you have not set concentrations or transforms for the compound.
- Pin settings are stored with the spectrum when you save it as a .cnx file. The compound that you pin remains pinned when you next open the file, until you clear the pins.
- You can export the 'pinned' status of all compounds to a file using Export Profiled Data.

Related Topics

- “Legend” on page 17
- “Compound Table” on page 46
- “Custom Colors” on page 48
- “Starred Compounds” on page 50
- “Legend (Sidebar)” on page 64

Starred Compounds

Starring compounds lets you mark them for future reference, or to perform group operations on them, including setting them to maximum concentrations, scaling concentrations, and clearing concentrations, transforms, custom colors, pins or stars. Stars are saved in the spectrum file, and remain visible when you next open the file, or until you clear them.

How do I star a compound?

1. Click the grey star for the compound in the Starred column (★) of the Compound Table.

How do I clear a starred compound?

1. Select the compound in the Compound Table.
2. In the Compound menu, click **Selected Compound > Clear > Star**, or click the yellow star for the compound in the Starred column (★).

Menu Location	Compound > Selected Compound > Clear > Star
Icon(s)	

Tips and Tricks

- Starring a compound adds it to the Profiled Compounds set, even if you have not set concentrations or transforms for the compound.
- Stars are stored with the spectrum when you save it as a .cnx file. The compound that you star remains starred when you next open the file, until you clear the stars.
- You can export the 'starred' status of all compounds to a file using Export Profiled Data.

Related Topics

- “Compound Table” on page 46
- “Custom Colors” on page 48
- “Pinned Compounds” on page 49
- “Concentrations” on page 52
- “Transforms” on page 53
- “Maximum Concentrations” on page 56
- “Scale Concentrations” on page 57

Quick Searches

When you are profiling a spectrum, you may have several compound sets active. The number of compounds appearing in the Compound Table can make it difficult to locate a particular compound to profile. Quick Searches let you temporarily limit the number of compounds appearing in the Compound Table based on search criteria. The search criteria can include text searches based on compound name, location, compound IDs, compound comments or all of these. You can also search by frequency, listing all compounds appearing at a specific frequency, or within a range of frequencies.

How do I limit the compounds visible in the Compound Table?

1. Click the binoculars icon at the left edge of the search box to choose the type of search to use (Compound Name, Compound Comment, Compound ID, Location or All). If you do not specify a search type, the default is All.

Icon(s)	
Hotkey(s)	Ctrl-L Cmd-L

2. Type your search terms in the Quick Search box at the top left of the Compound Table. The contents of the Compound Table update to show all matches to your current entry.
3. Click the red X at the right edge of the search box to cancel the search and restore the full list of compounds. Clearing the search field also resets the search type to All.

How do I search for compounds near a single frequency using the Spectrum View?

1. Right-click on the spectrum, and click **Filter for Compounds Near x.xx ppm**.
2. Click the red X at the right edge of the search box to cancel the search and restore the full list of compounds.

How do I search for compounds in a frequency range using a selected region?

1. Use the Select Region tool to select a region of the spectrum.
2. Right-click on the selected region, and click **Filter for Compounds in Selected Region**.
3. Click the red X at the right edge of the search box to cancel the search and restore the full list of compounds.

Filter Compounds is an alternate method for setting up a Quick Search based on a frequency or range that you type, rather than based on clicked locations in the Spectrum View.

Menu Location	Tools > Filter Compounds...
Icon(s)	

How do I filter the compounds in the Compound Table based on a specific frequency or range?

1. Open the Filter Compounds dialog.
2. Enter the frequency or range that you want to use to filter the Compound Table.
3. Click OK.

Tips and Tricks

- Use Quick Searches to temporarily limit the number of compounds displayed. Use compound sets created in Library Manager to define a more permanent working set of compounds for profiling a particular type of spectra.
- You can click repeatedly anywhere under the sum line (red) in the Spectrum View to cycle through the compounds contributing to the sum line at the frequency that you clicked.
- When you right-click anywhere under the sum line (red) in the Spectrum View, the resulting context menu shows a list of compounds contributing to the sum line at the frequency that you clicked. Select a compound from the list to make it the active compound.
- You can perform frequency and frequency range searches directly from the Quick Search box. For example, type **@2.89** to find compounds appearing at 2.89 ppm, or type **@3.0 to 7.0** to find compounds appearing between 3.0 and 7.0 ppm.
- You can reset the Quick Search box and clear all filters at any time by pressing **Esc**.
- While editing the Quick Search box, press the down arrow to navigate to the Compound Table.

Related Topics

- “Spectrum View Tools” on page 20
- “Compound Table” on page 46
- “Compound Set Tools” on page 59

Concentrations

While profiling a spectrum, you measure compound concentrations by determining heights of compounds signatures that best match peak heights in the current spectrum. You can only adjust the concentration of the compound currently selected in the Compound Table. Profiler determines the concentration of a compound using the standardized intensities (si) and standardized widths (sw) of its component peaks. When you adjust the height of a compound, the heights of all peaks and clusters corresponding to the compound scale proportionately.

When you have set a compound concentration to something other than unknown (--), Profiler adds it to the Profiled Compounds set at the top of the Compound Sets view. If you set a profiled compound's concentration to unknown (--), Profiler removes the compound from the Profiled Compounds set.

How do I set compound concentrations?

1. Select a compound in the Compound Table. If you do not see any compounds in the Compound Table, select additional compound sets or use Quick Searches to find one.
2. Click and drag the blue triangle on the vertical axis of the graph upward until the compound signature reaches the desired height. If the compound's concentration is already set, you can click and drag any peak or cluster vertically to change the concentration.

How do I clear a concentration?

1. Select a compound with a concentration in the Compound Table. If you do not see any compounds with concentrations in the Compound Table, make sure that you have selected the Profiled Compounds set. If you still do not see any compounds with concentrations, then the current spectrum has no concentrations set.
2. In the Compound menu, click **Selected Compound > Clear Concentration**, or right-click on the compound in the Spectrum View or the Compound Table and click **Clear Concentration**. The compound's concentration in the Compound Table appears as unknown (--).

Tips and Tricks

- You can have Profiler suggest concentration and transform values for the currently selected compound by clicking **Compound > Selected Compound > Fit Automatically**, or by pressing the space bar.
- You can double-click on a compound's concentration in the Compound Table and type a value directly to set the concentration. Type a number from 0 to 5000 mM (or equivalent concentrations in μM or mg/dL), or type "--" to set the concentration to unknown.
- When you have selected a compound in the Compound Table, you can click its shape in the Spectrum View, then use the W and S keys or up and down arrow keys to adjust its concentration up and down. Hold down **Shift** to allow finer adjustments.

- Once you have set a concentration for a compound, it appears in the Profiled Compounds set. Selecting that set shows only compounds that you have profiled in the Compound Table.
- You can hover the mouse cursor over any cluster, or its transform adjustment triangle on the horizontal axis, or the compound's concentration adjustment triangle on the vertical axis, to see the concentration for the current compound, even when the Concentration column is not visible in the Compound Table. The concentration appears in the top right corner of the Spectrum View.
- To clear all of the concentrations for the current spectrum, select **Compound > All Compounds > Clear Concentrations**.
- To clear concentrations only for starred compounds, select **Compound > Starred Compounds > Clear Concentrations**.
- All compound concentrations start out as unknown (--).
- If a compound's concentration is unknown (--) and you drag it up and then back down to the baseline, the concentration is set to 0.00 mM. To return the concentration to unknown, you must clear it.
- When you Export Profiled Data, compounds with unknown concentrations (--) are not exported unless other spectra in the exported set have those compounds profiled, or you have chosen to include compounds with no concentrations.

Related Topics

- "Reference Panel" on page 18
- "Spectrum View Tools" on page 20
- "Compound Table" on page 46
- "Starred Compounds" on page 50
- "Quick Searches" on page 50
- "Transforms" on page 53
- "Maximum Concentrations" on page 56
- "Scale Concentrations" on page 57
- "Compound Set Tools" on page 59
- "Export Profiled Data" on page 60
- "Reference Panel (Sidebar)" on page 65
- "Concentration Units" on page 66

Transforms

The frequencies of clusters in a compound signature that you are using may be different from those in the actual spectrum, due to changes in solution pH, ion concentrations, and other matrix effects. Transforms let you move the clusters from their starting locations to compensate for these effects and better fit the spectrum.

Transforms are normally limited to regions called transform windows. Profiler displays the transform window for a particular cluster as triangles on the spectrum baseline when the cluster is selected. You can move clusters beyond these limits if necessary.

How do I set transforms?

1. Select a compound in the Compound Table. If you do not see any compounds in the Compound Table, select additional compound sets or use Quick Searches to find one.
2. Click and drag the blue triangle for a cluster along the horizontal axis in the Spectrum View. If the compound already has a concentration set, you can click and drag any cluster horizontally to change its transform.

How do I clear transforms for a compound?

1. Select a transformed compound in the Compound Table. If you do not see any transformed compounds in the Compound Table, make sure the Profiled Compounds set is selected. If you still do not see any transformed compounds, then the current file may not contain any transformed compounds.
2. In the Compound menu, click **Selected > Clear Transform**, or right-click on the compound in the Spectrum View or the Compound Table and click **Clear Transform**. The compound's clusters return to their default locations.

Tips and Tricks

- You can have Profiler suggest concentration and transform values for the currently selected compound by clicking **Compound > Selected Compound > Fit Automatically**, or by pressing the space bar.
- When you have selected a compound in the Compound Table, you can click one of its clusters in the Spectrum View, then use the left and right arrow keys to adjust the cluster's transform. Hold down **Shift** to allow finer adjustments.
- If the transform window does not allow you to move a cluster far enough to fit your spectrum, you can disable enforcement of the transform windows by clicking **Edit > Enforce Transform Windows**.
- You can hover the mouse cursor over any cluster, or its transform adjustment triangle on the horizontal axis to see the transform applied to the cluster. The information appears in the top right corner of the Spectrum View, and visually as a thick blue line connecting the starting position of the cluster and its current position. If the cluster is not transformed, the display reads "Cluster is not transformed", and no visual indicator appears.
- When you transform a cluster outside its transform window, its transform adjustment triangle turns red. Also, when you hover over the cluster or its triangle, the visual indicator for the transform is red, as is the text in the top right corner of the Spectrum View describing the extent of the transform.
- To clear all of the transforms for the current spectrum, select **Compound > All Compounds > Clear Transforms**.
- To clear transforms only for starred compounds, select **Compound > Starred Compounds > Clear Transforms**.
- You can see which compounds have been transformed by turning on the **Transform** column in the Compound Table. The **Transform** column will also indicate whether the compound has any clusters that have been transformed beyond their transform windows.
- You can export a file describing the transform status of all compounds using Export Profiled Data.

Related Topics

- “Spectrum View Tools” on page 20
- “Compound Table” on page 46
- “Starred Compounds” on page 50
- “Quick Searches” on page 50
- “Concentrations” on page 52
- “Compound Set Tools” on page 59
- “Transform Windows” on page 86

Fit Automatically

Profiling a spectrum manually involves a lot of effort in setting initial concentrations and transforms. You can have Profiler suggest a concentration and transforms for the currently selected compound, allowing you to focus your efforts on tasks like refining your overall analysis of the spectrum, or resolving overlapped compounds.

How do I fit a compound automatically?

Menu Location	Compound > Selected Compound > Fit Automatically
Icon(s)	

1. Select the compound in the Compound Table. If you do not see the compound in the Compound Table, select additional compound sets or use Quick Searches to find it.
2. On the Compound menu, click **Selected Compound > Fit Automatically**, or right-click on the compound and select **Fit Automatically** from the context menu.

Tips and Tricks

- Although Profiler will often arrive at acceptable values, always review the results of an automatic fit before moving on to other compounds.
- Fit Automatically yields more reliable results from properly processed spectra. This includes well-defined CSI parameters, and accurate spectrum pH values.

Related Topics

- “CSI Editor” on page 31
- “pH Editor” on page 32
- “Concentrations” on page 52
- “Transforms” on page 53

Enforce Transform Windows

Transform windows appear in compound signatures to guide profiling by restricting the amount that you can transform each cluster. All Chenomx Reference Compounds have had transform windows set to allow the flexibility to fit most common biofluids. You can choose to have the transform windows enforced, preventing you from moving any cluster outside its transform window, or not enforced, allowing you to move any cluster beyond its transform window.

Menu Location	Edit > Enforce Transform Windows
Hotkey(s)	Ctrl-Shift-X Cmd-Shift-X

How do I enforce transform windows?

1. On the Edit menu, check **Enforce Transform Windows**.

How do I allow clusters to move outside their transform windows?

1. On the Edit menu, clear **Enforce Transform Windows**.

Tips and Tricks

- Clusters that are outside their transform window appear in the Transforms column of the Compound Table as **Outside Windows**. Note that you may need to turn on the Transforms column in the Compound Table.
- When you transform a cluster outside its transform window, its transform adjustment triangle turns red. Also, when you hover over the cluster or its triangle, the visual indicator for the transform is red, as is the text in the top right corner of the Spectrum View describing the extent of the transform.
- Transform windows for pH sensitive compound signatures are determined in part based on the pH uncertainty for the current spectrum. If you find that you frequently need to override the transform windows for a spectrum, check the spectrum pH, and consider increasing the pH uncertainty value.
- If, while Enforce Transform Windows is on, you try to move a cluster that has been set outside its transform window, the cluster will first snap to its nearest transform window edge. After that, you can move it normally, within its transform window.

Related Topics

- “pH Editor” on page 32
- “Compound Table” on page 46
- “Transforms” on page 53
- “Cluster Navigator” on page 58
- “Transform Windows” on page 86

Maximum Concentrations

The maximum concentration of a compound in Profiler is an estimate of the upper limit of the compound's concentration. Profiler calculates this value as the maximum concentration for the compound that does not result in any peak in the compound's signature exceeding the height of the spectrum line, or any cluster moving outside its transform window. The value is made available as a guideline only.

How do I set a compound to its maximum concentration?

1. Select the compound in the Compound Table. If you do not see the compound in the Compound Table, select additional compound sets or use Quick Searches to find it.
2. On the Compound menu, click **Selected Compound > Set to Maximum Concentration**, or right-click on the compound and select **Set to Maximum Concentration** from the context menu.

Menu Location	Compound > Selected Compound > Set to Maximum Concentration
Icon(s)	

Menu Location	Compound > Starred Compounds > Set to Maximum Concentrations
Icon(s)	

Menu Location	Compound > All Compounds > Set to Maximum Concentrations
Icon(s)	

How do I set starred compounds to their maximum concentrations?

1. Star at least one compound.
2. On the Compound menu, click **Starred Compounds > Set to Maximum Concentrations**.

How do I set all compounds to their maximum concentrations?

1. On the Compound menu, click **All Compounds > Set to Maximum Concentrations**.

Tips and Tricks

- Maximum concentrations in Profiler refer to estimates that are subject to error, based on assumptions that may not coincide with your interpretation of the spectrum. Verify any results that you obtain based on maximum concentrations.
- The method used to determine maximum concentrations is simpler than that used for automatic fits. Concentrations calculated by the two methods may differ.

Related Topics

- “Compound Table” on page 46
- “Starred Compounds” on page 50
- “Quick Searches” on page 50
- “Concentrations” on page 52
- “Fit Automatically” on page 55
- “Compound Set Tools” on page 100

Scale Concentrations

Scale Concentrations lets you modify the concentrations of all profiled compounds in the current spectrum by a common multiplication factor. You can use this to compensate for dilution effects, errors in CSI concentration or other systematic influences on compound concentrations.

Menu Location	Compound > All Compounds > Scale Concentrations...
Icon(s)	

Menu Location	Compound > Starred Compounds > Scale Concentrations...
Icon(s)	

How do I scale the concentrations of all compounds by a common factor?

1. On the Compound menu, click **All Compounds > Scale Concentrations...**
2. Enter the multiplication factor you would like to use. To increase concentrations, use a factor greater than one. To decrease them, use a factor less than one.
3. Click OK.

How do I scale the concentrations of starred compounds by a common factor?

1. Star at least one compound.
2. On the Compound menu, click **Starred Compounds > Scale Concentrations...**
3. Enter the multiplication factor you would like to use. To increase concentrations, use a factor greater than one. To decrease them, use a factor less than one.
4. Click OK.

Tips and Tricks

- If you enter a multiplication factor that would result in a concentration higher than 5000 mM for any compounds, those compounds are set to 5000 mM. All other compounds scale with the exact factor entered.
- If you set the multiplication factor to zero, all of the profiled compounds are set to 0 mM.
- Scale Concentrations does not affect compounds with unknown concentrations (--).

Related Topics

- “Starred Compounds” on page 50
- “Concentrations” on page 52
- “Compound Set Tools” on page 59

Cluster Navigator

The Cluster Navigator lets you move quickly among the clusters of the currently displayed compound. It appears in the top left corner of the Spectrum View as an ordered list of the approximate positions of all the clusters of the compound in ppm. When you click on a cluster position, the view zooms to provide an optimal view of the cluster, including the larger of the full transform range of the cluster or the full width of the cluster, and the full vertical extent of any peaks within the horizontal range. You can also use the Cluster Navigator to move sequentially through the list of clusters in either direction.

Menu Location	View > Show Cluster Navigator
Hotkey(s)	Ctrl-Shift-N Cmd-Shift-N

How do I display the clusters of a compound in the Spectrum View?

1. Make sure that the Cluster Navigator is visible, and select a compound in the Compound Table.
2. Click a cluster position to zoom directly to that cluster. Click the left or right arrow to move sequentially through the list of clusters.

Tips and Tricks

- Use the Q and E keys to display the next or previous cluster.
- Clusters that appear green in the Cluster Navigator are matched, meaning they closely match the spectrum line as determined by comparing the sum line (red) and the spectrum line (black) near each cluster.
- You can click on a cluster and drag along the Cluster Navigator to display multiple clusters in the Spectrum View.

Related Topics

- “Spectrum View Tools” on page 20
- “Compound Table” on page 46
- “Concentrations” on page 52
- “Transforms” on page 53
- “Enforce Transform Windows” on page 55
- “Compound Set Tools” on page 59

6.4 Compound Set Tools

Compound sets are groups containing compounds found in your Library. You can use compound sets that you have created in Library Manager to make the task of profiling your spectra easier, by placing compounds that interest you the most at your fingertips.

Compound Sets (Sidebar)

Menu Location	View > Sidebar > Compound Sets
Icon(s)	
Hotkey(s)	Ctrl-Shift-K Cmd-Shift-K

The Compound Sets view lets you modify the contents of the Compound Table to include any combination of the Compound Sets that you have created in Library Manager. The panel contains a list of every compound set that is currently defined. Contents of the checked compound sets appear in the Compound Table.

Tips and Tricks

- Reference Compound Sets, containing any Chenomx Reference Compounds in your library, appear automatically in the Compound Sets view. Their icons are colored green.
- When you filter the Compound Table by frequency or range, only compounds from the currently active compound sets are considered. You can use this, for example, by creating a compound set containing only compounds that appear in urine samples. Then, when this urine compound set is the only one active, any quick searches that you perform will only return compounds that appear in urine samples as potential hits.
- A compound set called Profiled Compounds always appears at the top of the Compound Sets view in Profiler. The Profiled Compounds set includes all compounds in the current spectrum that have been modified, based on concentrations, transforms, custom colors, pins and stars.
- The Profiled Compounds set updates continuously, so if you clear the concentration of a compound, it immediately disappears from Profiled Compounds. In this case, you can still select the compound from its original compound set.
- If you change the contents of a compound set in Library Manager, your changes immediately appear in Profiler, including modified selection rules for Smart Compound Sets, manual changes to Compound Sets by adding or removing compounds, and created or deleted compound sets.

Related Topics

- “Profiling Tools” on page 46
- “Open Library Manager” on page 59

Open Library Manager

While profiling a spectrum, you may want to define a new compound set that contains a smaller list of compounds, or modify an existing compound set to focus more directly on your current project. To access and edit compound sets, you can open the Library Manager module from Profiler.

Menu Location	File > Open Library Manager
Icon(s)	
Hotkey(s)	Ctrl-B Cmd-B

Tips and Tricks

- Compound sets that you change in Library Manager update immediately in Profiler.
- Profiler does not close when you open Library Manager. Also, you can leave Library Manager open for future reference while you continue your work in Profiler.

Related Topics

- “Compound Set Tools” on page 59
- Chapter 9, *Library Manager* on page 94

6.5 Importing and Exporting Data

You can export data from your profiled spectra for further analysis, bin the spectra for analysis by more traditional techniques, or create images of the spectra for use in publications and presentations.

Export Profiled Data

When you have finished profiling a spectrum or set of spectra, you often need the resulting data in a spreadsheet or statistical analysis package. Export Profiled Data lets you create a text file containing all of the relevant details from profiled spectra, in a format that you can easily open in software like Microsoft Excel or Umetrics SIMCA-P.

How do I export data from one or more spectra?

1. Start the Export Profiled Data wizard.
2. Choose the source spectra to export, and specify a target folder and filename for the exported text file. You can specify the currently open spectrum, other individual spectrum files or one or more folders containing spectrum files as the source.
3. Select the type of data that you would like to export. You can export concentrations in various units, other profiling results like matched clusters and transform summaries, or annotations like stars, pins and custom colors.
4. Choose the compounds for which you would like to export data. You can use custom compound sets that you have created in Library Manager to filter the output to specific compounds, or you can export all profiled compounds in all of the spectra.
5. Choose additional data that you would like to include in the exported file. Additional data can include spectrum pH values, Chenomx compound IDs for the exported compounds, or InChI or SMILES strings for the exported compounds.
6. Choose the default layout (files in rows and compounds in columns) or transposed layout (compounds in rows, files in columns) for the output file, depending on how you intend to use the file.

Tips and Tricks

- If you do not explicitly specify a target file and folder, an appropriately-named file is created in the folder containing your source spectra. Possible file names include:
 - concentrations.txt
 - matched-clusters.txt
 - transform-summaries.txt
 - stars.txt
 - pins.txt
 - colors.txt

Menu Location	File > Export > Profiled Data (.txt)...
Icon(s)	

- To filter the exported compounds without modifying the source spectra, create a Compound Set or Smart Compound Set in Library Manager, and only export compounds in your custom Compound Set.
- To set spectrum pH values, open the spectrum in Processor, and use the pH editor. The value that you set in the pH editor appears in the exported file when you choose spectrum pH values as extra data to export.
- Use the default layout to export files that you intend to open with statistical analysis packages. Use the transposed layout to export files that you intend to open in a spreadsheet program, or if you prefer the transposed layout.
- If you intend to open the exported file in Microsoft Excel 2003 or earlier, limit the number of compounds to less than 253 for the Default layout, or the number of files to less than 253 for the Transposed layout, since these versions of Excel do not support large numbers of columns.

Related Topics

- “Profiling Tools” on page 46
- “Compound Set Tools” on page 59
- “Concentration Units” on page 66

Export Profiled Compounds

When you collaborate with other people on a project, you may have different compounds in your library than your collaborator has in theirs. Exporting profiled compounds creates a compound pack (.pack) containing the set of actual compound signatures that you used to fit spectra for the project. You can exchange compound packs with other users of Chenomx NMR Suite, and import new packs via Library Manager. You can then be certain that you and your collaborator are always profiling spectra for that project using the same collection of compound signatures.

Menu Location	File > Export > Profiled Compounds (.pack)...
Icon(s)	

How do I export a compound pack from the current spectrum?

1. Open the Export Profiled Compounds dialog.
2. Enter a name for the exported compound pack. The name you enter is used to label the associated Compound Set when the compound pack is imported into another library via Library Manager.
3. Click OK.
4. In the resulting file chooser, locate the folder in which you would like to save the compound pack.
5. In the File Name field, type the new name for the compound pack file. The specified name can be different from the name of the compound pack itself, but you may prefer to keep them the same.
6. Click **Save**.

Tips and Tricks

- You and your collaborator each need a valid licensed copy or evaluation edition of Chenomx NMR Suite, version 5.0 or later, installed to make use of compound packs.

- If you want to export all of the compounds in a compound set, regardless of whether they have been fit in a particular spectrum, use **Export as Compound Pack** in Library Manager.

Related Topics

- “Profiling Tools” on page 46
- “Compound Set Tools” on page 59
- “Add Compounds” on page 97
- “Export Compounds” on page 98
- “Compound Set Tools” on page 100

Import Profile

When you are profiling a series of spectra from the same study, the spectra are often quite similar to one another. While it is possible to profile each spectrum from scratch, it can be more efficient to copy a profile from an earlier spectrum in the series and adjust it to match the current spectrum. Importing profiles is an important technique for profiling large series of similar spectra efficiently.

How do I import a profile from another spectrum?

1. Open the Import Profile dialog.
2. Select a profiled spectrum (.cnx) from which you would like to import the profile.
3. Click Import.

Tips and Tricks

- When you import a profile, the imported data completely replaces any profile of the current spectrum that you may have done before the import.

Related Topics

- “Profiling Tools” on page 46

Spectral Binning

Statistical approaches to analyzing large numbers of spectra often involve data reduction methods like spectral binning to reduce the number of variables that must be considered. Spectral binning divides spectra into a series of regions, or bins; subsequent analysis involves the integrated area of these bins instead of the raw spectral data. Profiler lets you perform spectral binning on collections of spectra, exporting the results to a text file format that is easily read by spreadsheet programs and statistical analysis software packages.

How do I bin spectra?

1. Start the Spectral Binning wizard.
2. Choose the source spectra (.cnx) for binning, and specify a target folder for the resulting binned data. You can specify a series of spectrum files or one or more folders containing spectrum files as the source.

If you would like to repeat a binning session that you have previously saved, you can specify a binning configuration file (**BinningConfiguration.xml**) to import.

Menu Location	File > Import Profile...
Icon(s)	

Menu Location	Tools > Spectral Binning...
Icon(s)	

3. Choose a binning method to apply to the specified spectra. You can use a defined bin size in ppm across a range of the spectrum in ppm or a defined number of bins. You can also use bins based on the transform windows of compounds in a selected Compound Set ("targeted binning").

If you would like to reuse a predefined set of bins, you can specify a bin definition file (`BinDefinitions.txt`) to load.

4. Choose regions of the spectrum to exclude from binning, if desired. You can specify any regions of the spectrum using start and end values in ppm.
5. Choose a binning target. You can bin spectra based on the sum line (standard binning) or the subtraction line (residual binning).
6. Choose a normalization method. You can normalize the bins relative to the total area under the spectrum line or to units of standardized area (sa), where 1 sa is the area under a theoretical DSS methyl peak at 0.50 mM.

Tips and Tricks

- Binning can take a long time, especially on older computers or particularly large datasets. If you have never tried binning before, try a small dataset of five or ten spectra first, before moving on to larger datasets. The time needed for the smaller dataset will help you to more accurately estimate the time needed to complete a binning session.
- If you choose one or more folders as the source, binning only includes spectra both directly in the specified folders and in subfolders of the specified folders.
- If you do not explicitly specify a target folder, output files are created in the folder containing your source spectra.
- Binning based on number of bins divides the spectrum into the specified number of bins after excluding the regions that you have defined, if any.
- To do residual binning, each of the source spectra must have been analyzed in Profiler and must contain at least one compound with a defined concentration.
- Excluding regions from a binning dataset can help you to remove the influence of interfering or confounding signals.
- Normalization based on total area gives the area of a particular bin as a fraction of the total area of all bins, except area in any of the excluded regions that you may have defined.
- Spectral binning is subject to several numerical limits:
 - Each binning session can include a maximum of 5000 spectra.
 - The range that you define when using number of bins or size of bins as a binning method must be at least 1 ppm wide, and can be as large as 1998 ppm (-999 to 999 ppm).
 - You can specify a maximum of 1000000 (one million) bins or a minimum of 1.
 - You can specify a minimum bin width of 0.00001 ppm and a maximum bin width equal to the width of the range that you have chosen. For example, if you are binning across a range of 2 to 10 ppm, your maximum bin width is 8 ppm.

Related Topics

- "Profiling Tools" on page 46

- “Compound Set Tools” on page 59

Send to Processor

Menu Location	File > Send to Processor
Icon(s)	
Hotkey(s)	Ctrl-J Cmd-J

As you work with a spectrum in Profiler, you may notice that the phasing, baseline correction or other aspects of the spectrum's lineshapes are not properly adjusted. Similarly, you may import a spectrum directly into Profiler, only to realize that it requires phasing, baseline or other correction. Send to Processor closes the current spectrum in Profiler, and opens it in Processor, retaining any changes you have made in Profiler.

Tips and Tricks

- Any changes that you have made to the spectrum in Profiler are retained when you send the spectrum to Processor, and are visible when you bring the spectrum back into Profiler. However, you still need to save the spectrum, either in Processor or Profiler, to permanently save those changes to the spectrum file.

Related Topics

- Chapter 5, *Processor* on page 28

6.6 Information Tools

While profiling spectra, you can view acquisition and sample information connected to the current spectrum, and change the concentration units that appear in the Compound Table and any exported profiles. The Sidebar lets you change the lines displayed in the Spectrum View, select compound sets to use in fitting spectra, and view extended information about the compounds in your Library

Legend (Sidebar)

The Legend is an interactive display of the various lines currently appearing in the Spectrum View. In Profiler, you can use the Legend to pin or unpin compound lines to the spectrum, change their colors, or make a pinned compound active for fitting.

How do I pin a compound using the Legend?

1. Switch to the Legend view in the Sidebar.
2. Select a compound in the Compound Table.
3. In the Legend view, click the gray pin icon (📌) appearing to the right of the selection line.

How do I unpin a compound using the pinned lines in the Legend?

1. Switch to the Legend view in the Sidebar.
2. Locate a pinned compound in the Legend view.
3. Click the red X (✖) appearing beside the pinned compound's name.

How do I unpin a compound using the Legend?

1. Switch to the Legend view in the Sidebar.
2. Locate a pinned compound in the Legend view.
3. Click the red X (✖) appearing to the right of the pinned compound's name.

Menu Location	View > Sidebar > Legend
Icon(s)	
Hotkey(s)	Ctrl-Shift-L Cmd-Shift-L

How do I unpin a compound using the selection line in the Legend?

1. Switch to the Legend view in the Sidebar.
2. Select a compound in the Compound Table.
3. Click the red pin () appearing beside the selection line.

How do I make a pinned compound active using the Legend?

1. Switch to the Legend view in the Sidebar.
2. Click a pinned compound in the Legend view.

Tips and Tricks

- The compound currently selected in the Compound Table will always appear in the Legend as the selection line.
- You can adjust the display of compound lines through the Legend view in the Sidebar, including setting custom colors, pinning compounds and selecting a compound for use in profiling.

Related Topics

- “Legend” on page 17
- “Compound Table” on page 46
- “Custom Colors” on page 48
- “Pinned Compounds” on page 49

Reference Panel (Sidebar)

Menu Location	View > Sidebar > Reference Panel
Icon(s)	
Hotkey(s)	Ctrl-Shift-R Cmd-Shift-R

The Reference Panel Sidebar view displays information specific to the compound that you have selected in the Compound Table. This includes a structure image, molecular formula and weight, alternate names and links to external reference sites. You can use this information to verify the identity of the compound that you are fitting and better understand the signature and behaviour of the compound in your spectra.

Tips and Tricks

- You can open the Reference Panel for a specific compound by double-clicking the compound in the Compound Table. Double-click again to close the Reference Panel.

Related Topics

- “Reference Panel” on page 18
- “Preferences” on page 25
- “Profiling Tools” on page 46

Spectrum Details

Menu Location	File > Spectrum Details...
Icon(s)	

Spectrum Details lets you view acquisition parameters for the current spectrum, including sweep width, magnet frequency, pulse sequence and so on, as well as user-defined properties such as sample pH and comments that have been added in Processor.

Tips and Tricks

- The information in Spectrum Details is read-only in Profiler. To change sample pH or add or modify comments for the spectrum, send the spectrum to Processor.

Related Topics

- Chapter 5, *Processor* on page 28

Concentration Units

Menu Location	View > Concentration Units
---------------	----------------------------

You can change the units of concentration used in the Compound Table. Available units include millimoles per litre (mM), micromoles per litre (μ M) and milligrams per decilitre (mg/dL).

Tips and Tricks

- The maximum allowed concentration for any compound is 5000 mM (5 M). When you select mg/dL as the concentration unit, the maximum value varies per compound, since the units incorporate the molecular weight of the compound. However, the maximum value is always equivalent to 5000 mM.
- Concentrations in mg/dL are only available for compounds whose molecular weight has been specified, and the accuracy of such concentrations depends on the accuracy of the molecular weight provided. To change the molecular weight associated with a compound, open it in Compound Builder and edit the Compound Details.

Related Topics

- “Profiling Tools” on page 46
- “Open Library Manager” on page 59
- “Export Profiled Data” on page 60
- “Compound Details” on page 91
- “Send to Compound Builder” on page 99


Chapter 7: Spin Simulator

Spin Simulator is a simple yet powerful tool for creating simulations of NMR spectra, based on user-defined spin systems and reference spectra. You can use these simulations as starting points for creating your own compound signatures, especially if you need to create signatures at multiple magnet frequencies, or you want to create robust signatures with a strong theoretical backing.

To use Spin Simulator effectively, you must have a basic understanding of NMR theory, and know how to interpret NMR spectra, especially J-coupling. Simulating complex compounds requires at least this basic theoretical knowledge. Chenomx experts can provide compound fitting services and develop specialized compound libraries for you. For more information about customizing compound libraries and related services, please contact Chenomx Contract Services.

You can only use Spin Simulator if you have purchased the Professional edition of Chenomx NMR Suite, or if you are using the Evaluation edition of Chenomx NMR Suite.

7.1 Overview



- ❶ Sidebar View Buttons. Let you change the currently visible Sidebar view.
- ❷ Simulation. Sidebar view that displays Spin Definition Tabs for the currently selected spin system. Also lets you select spin systems (via a drop-down menu) and spin definitions (via the tabs).
- ❸ J-Value Editor. Lets you vary the J-value for the currently selected J-modifier.

- ④ J-Modifier List. Lets you add, select, edit or remove J-modifiers associated with the currently selected spin definition.
- ⑤ Simulation Details. Display information about the current simulation, including the simulation frequency and numbers of spin systems, transitions, and J-modifiers represented in the simulation.
- ⑥ Spin Navigator. Lets you navigate among the spin definitions of the currently open simulation. You can enable or disable the Spin Navigator using the **View** menu.
- ⑦ Spectrum View Tools. Let you navigate the spectrum more easily.
- ⑧ Spectrum View. Displays the spectrum overlay as well as the simulation line and the currently selected spin definitions.
- ⑨ Spectrum Thumbnail. Displays a thumbnail image of the current spectrum. You can also use the thumbnail to navigate the Spectrum View at various zoom settings.
- ⑩ Status Bar. Displays the current cursor position in the units specified by the **View > X Scale** and **View > Y Scale** settings. Also displays information about recent user actions.

7.2 Quick Reference Guide

	Linux, Windows	Mac OS X
File Menu		
 New Simulation...	Ctrl-N	Cmd-N
 Open...	Ctrl-O	Cmd-O
 Open Recent	---	---
 Close	Ctrl-W	Cmd-W
 Save	Ctrl-S	Cmd-S
 Save As...	---	---
 Export > Spectrum Image...	---	---
 Preferences...	---	---
 Exit	---	---
Edit Menu		
 Undo	Ctrl-Z	Cmd-Z
 Redo	Ctrl-Y	Cmd-Y
 Copy Spectrum Image to Clipboard	---	---
Select All	Ctrl-A	Cmd-A
View Menu		
X Scale	---	---
Y Scale	---	---
 Zoom > Show Entire Spectrum	Ctrl-Minus	Cmd-Minus

	Linux, Windows	Mac OS X
 Zoom > Undo Zoom	[	[
 Zoom > Redo Zoom	]	]
 Zoom > Set Zoom...	Ctrl-M	Cmd-M
 Zoom > Select Region	Ctrl-T	Cmd-T
 Zoom > Auto Zoom	Ctrl-F	Cmd-F
 Zoom > Increase Vertical Zoom	+	+
 Zoom > Decrease Vertical Zoom	-	-
 Zoom > Increase Horizontal Zoom	*	*
 Zoom > Decrease Horizontal Zoom	/	/
Sidebar > Legend	Ctrl-Shift-L	Cmd-Shift-L
Sidebar > Simulation	Ctrl-Shift-D	Cmd-Shift-D
Show Sidebar	Ctrl-Shift-S	Cmd-Shift-S
Show Spectrum Thumbnail	Ctrl-Shift-T	Cmd-Shift-T
Show Spin Navigator	Ctrl-Shift-N	Cmd-Shift-N
Show J-Value Editor	Ctrl-Shift-J	Cmd-Shift-J
Simulation Menu		
 Add New Spin System	---	---
 Delete Current Spin System	---	---
 Add New Spin Definition...	Ctrl-D	Cmd-D
 Delete Current Spin Definition	---	---
 Add New J-Modifier...	Ctrl-J	Cmd-J
 Edit Selected J-Modifier...	---	---
 Delete Selected J-Modifiers	---	---
 Uncouple Selected Split	---	---
 Group/Ungroup Selected Splits	---	---
 Previous Spin Definition	Q	Q
 Next Spin Definition	E	E
Tools Menu		

	Linux, Windows	Mac OS X
 Set Spectrum Overlay...	Ctrl-Shift-O	Cmd-Shift-O
 Clear Spectrum Overlay	---	---
NMR Suite Help...		
 NMR Suite Help...	F1	F1
 Save Application Log...	Ctrl-Shift-A	Cmd-Shift-A
 Software Activation...	---	---
 About Chenomx NMR Suite...	---	---

7.3 Simulating Tools

Spin Simulator provides you with a variety of tools to help you create accurate spin simulations, by adding, manipulating and removing spin definitions and coupling constants, grouping spin definitions into spin systems, selecting simulation algorithms, and viewing summary details about your current simulation.

Spin simulation is a computationally intensive process, and takes more time as the simulation you are creating becomes more complex. Generally, the amount of time required to recalculate a simulation increases more or less linearly with the number of spin definitions (protons) in the simulation, and exponentially with the number of distinct J-couplings. Spin Simulator monitors the overall complexity of your simulation as you work, and if you attempt to take an action that would result in unusable simulation times, the action is cancelled, and you see a warning dialog indicating that the resulting simulation would be too complex.

New Simulation

New Simulation creates a new simulation file (.xss). The new simulation file has the default simulation frequency of 800 MHz.

Menu Location	File > New Simulation...
Icon(s)	
Hotkey(s)	Ctrl-N Cmd-N

Tips and Tricks

- The new simulation file does not have a spectrum overlay. Add a spectrum overlay before you continue creating a simulation.
- You can overlay only one spectrum at a time in Spin Simulator.
- You need to save the new file to retain changes that you have made.

Related Topics

- “Spectrum Overlays” on page 19
- “Simulation Details” on page 76

Spin Systems

If you need to model many protons in a simulation, you may want to use spin systems to organize them. A spin system in Spin Simulator is a set of spin definitions that you have defined as only interacting with other spin definitions within the set. Chemically, this may be loosely interpreted as a set of protons that couple only to each other.

Menu Location	Simulation > Add New Spin System...
Icon(s)	

How do I add a new spin system?

1. Click **Add New Spin System**.
2. If this is the first spin system you have added, a new dropdown list appears in the Simulation view, above the Spin Definition Tabs, listing the available spin systems. Additional spin systems also appear in this dropdown list.

How do I select a different spin system?

1. In the Simulation view, select the spin system to which you would like to switch from the drop-down list.

Alternately, select a spin definition belonging to the other spin system in the Spectrum View.

Menu Location	Simulation > Delete Current Spin System...
Icon(s)	

How do I delete a spin system?

1. Select a spin system from the drop-down list in the Simulation view.
2. Click **Delete Current Spin System**.

Tips and Tricks

- Deleting a spin system, deletes all spin definitions and J-modifiers associated with the spin system as well.

Related Topics

- “Spin Definitions” on page 71
- “Spin Navigator” on page 73
- “Simulation (Sidebar)” on page 75
- “Simulation Details” on page 76

Spin Definitions

Spin Simulator creates computer models of the interactions among the members of a spin system. Within Spin Simulator, you may consider a member of a spin system as equivalent to a single proton, or a group of protons that are both chemical shift equivalent and magnetic-equivalent. Spin definitions described the members of spin systems, and are the basic building blocks of the models created by Spin Simulator.

A spin definition consists of three basic pieces of information: The number of protons represented by the spin definition, the center of the spin definition in ppm (equivalent to a chemical shift in an NMR spectrum), and a width in sw (analogous to a linewidth at half-height in an NMR spectrum).

You can organize spin definitions into spin systems, as described elsewhere in this chapter (see “Spin Systems” on page 70).

Menu Location	Simulation > Add New Spin Definition...
Icon(s)	
Hotkey(s)	Ctrl-D Cmd-D

How do I add a new spin definition?

1. Open the Add New Spin Definition dialog.
2. Enter a number of protons, starting center in ppm, and starting width in sw. You can change the entered values later.
3. Click OK.

How do I edit a spin definition?

1. Select the spin definition by clicking on it in the Spectrum View, or selecting its Spin Definition Tab in the Simulation view.
2. Do any of the following:
 - Click and drag the blue triangle on the horizontal axis to change the center of the spin definition.
 - Click and drag the blue triangle on the vertical axis to change the height of the simulation.
 - Click and drag the spin definition directly in the Spectrum View to change its center and the simulation height simultaneously.
 - Click and drag the blue triangles on either edge of the spin definition (at around half the height of the tallest peak) to change the width of the spin definition.
 - Use the W and S keys or the up and down arrow keys to adjust the height of the whole simulation up and down.
 - Use the A and D keys or the left and right arrow keys to adjust the position of the selected spin definitions left and right.
 - Use the Q and E keys to display the next or previous spin definition.

Menu Location	Simulation > Delete Current Spin Definition...
Icon(s)	

How do I delete a spin definition?

1. Select a spin definition using the Simulation view or the Spin Navigator.
2. Click **Delete Current Spin Definition**.

Tips and Tricks

- A spin definition can describe no more than ten protons. However, given that a spin definition typically describes chemically and magnetically equivalent protons, you may find that a value greater than three (a methyl group) is not often necessary.
- A spin system can contain any number of protons and spin definitions, as long as the total number of transitions is less than 2000. If you need to create a simulation requiring more than 2000 transitions, you will need to use multiple spin systems.

The total number of transitions in the current simulation appears in the simulation details (see “Simulation Details” on page 76).

- Height adjustments always affect the whole simulation, not just the currently selected spin definition. If you need to adjust the heights more than is possible in Spin Simulator, save the simulation and open it in Compound Builder, then make your adjustments there.
- Deleting a spin definition also deletes any J-modifiers that you have associated with that spin definition.

Related Topics

- “Spin Systems” on page 70
- “Spin Navigator” on page 73
- “J-Modifiers” on page 73
- “Simulation (Sidebar)” on page 75

- “Simulation Details” on page 76

Spin Navigator

The Spin Navigator lets you quickly move among the spin definitions that you have created in the current simulation. It appears in the top left corner of the Spectrum View as a ordered list of the approximate positions of all the spin definitions of the simulation in ppm. When you click on a spin definition's position, the view zooms to provide an optimal view of the spin definition. The optimal view includes the full width of the cluster and the full vertical extent of any peaks within the horizontal range. You can also use the Spin Navigator to move sequentially through the list of spin definitions in either direction.

How do I display the spin definitions of a simulation in the Spectrum View?

Menu Location	View > Show Spin Navigator
Hotkey(s)	Ctrl-Shift-N Cmd-Shift-N

1. Make sure that the Spin Navigator is visible.
2. Click a spin definition's position to zoom directly to that spin definition. Click the left or right arrow to move sequentially through the list of spin definitions.

Tips and Tricks

- Click on a spin definition and drag along the Spin Navigator to include adjacent spin definitions in the Spectrum View.

Related Topics

- “Spin Systems” on page 70
- “Spin Definitions” on page 71
- “Simulation (Sidebar)” on page 75
- “Simulation Details” on page 76

J-Modifiers

When you create a new spin definition, it appears as a single, unsplit peak. J-modifiers let you split a spin definition into more complex shapes, to better fit real spectra. J-modifiers can be uncoupled, affecting only one spin definition, or coupled, linking the spin definition to another one via its coupling constant. Coupled J-modifiers result in more accurate simulations, as Spin Simulator can use them to calculate higher-order coupling effects, accounting for tenting and other complex behaviour.

How do I add an uncoupled J-modifier to a spin definition?

Menu Location	Compound > Add New J-Modifier...
Icon(s)	
Hotkey(s)	Ctrl-J Cmd-J

1. Select the spin definition to which you would like to add a J-modifier.
2. Open the **Add New J-Modifier...** dialog.
3. Enter a J-value for the new J-modifier, in Hz.
4. Select Uncoupled Split, and specify a split count for the new J-modifier.
5. Click OK.

How do I add a coupled J-modifier to a spin definition?

1. Select one of the two spin definitions that you would like to link with a coupled J-modifier.
2. Open the **Add New J-Modifier...** dialog.
3. Enter a J-value for the new J-modifier, in Hz.

Menu Location	Compound > Edit Selected J-Modifier...
Icon(s)	

4. Select Coupled, and select the spin definition to which you would like to link.
5. Click OK.

How do I edit an existing J-modifier?

1. Select a spin definition.
2. Select a J-modifier that you would like to edit.
3. Open the **Edit Selected J-Modifier...** dialog.
4. Make your desired changes.
5. Click OK.

How do I couple two spin definitions using splits that I have already defined?

1. Select one of the two spin definitions that you would like to link, and select the uncoupled split that will link it to the other spin definition.
2. Open the **Edit Selected J-Modifier...** dialog.
3. Select Coupled, and select the spin definition to which you would like to link.
4. Select the uncoupled split from the second spin definition that you would like to link to.
5. Click OK.

Menu Location	Compound > Delete Selected J-Modifier
Icon(s)	

How do I delete a J-modifier?

1. Select a spin definition.
2. Select a J-modifier that you would like to delete.
3. Click **Delete Selected J-Modifier**.

Menu Location	Compound > Uncouple
Icon(s)	

How do I uncouple a J-modifier?

1. Select a spin definition.
2. Select a J-modifier that you would like to uncouple.
3. Click **Uncouple**.

Menu Location	Compound > Ungroup Selected Splits
Icon(s)	

How do I ungroup J-modifiers?

1. Select a spin definition.
2. Select an uncoupled split with a split count greater than 1.
3. Click **Ungroup Selected Splits..**

Menu Location	Compound > Group Selected Splits
Icon(s)	

How do I group several J-modifiers?

1. Select a spin definition.
2. Select several uncoupled splits that you would like to combine. Hold down **Ctrl** to multi-select J-modifiers.
3. Click **Group Selected Splits..**

Tips and Tricks

- J-Modifiers can only link two spin definitions that are in the same spin system.

- When you delete a coupled J-modifier, you will delete it from both spin definitions to which it applies. If you only want to delete the J-modifier from one spin definition, uncouple it first.
- Uncoupling a coupled J-modifier will result in an uncoupled split on each of the two previously coupled spin definitions, each split having the same J-value as the original coupled J-value.
- Most operations involving J-modifiers are available by right-clicking on a J-modifier in the Simulation view.

Related Topics

- “Spin Systems” on page 70
- “Spin Definitions” on page 71
- “Simulation (Sidebar)” on page 75

7.4 Information Tools

While creating spin simulations, you can view and update spin systems, spin definitions and J-modifiers associated with the current compound, and view overall details on the current simulation.

Simulation (Sidebar)

The Simulation Sidebar view contains Spin Definition Tabs, which let you view and edit the properties of spin definitions, and create, edit or remove J-modifiers.

How do I edit spin definitions and J-modifiers using the Simulation view?

1. Select a spin definition by clicking it in the Spectrum View, or by selecting its tab in the Simulation view.
2. Click in any of the text fields and enter a new value to edit the number of protons, center or width of the spin definition.
3. Edit J-modifiers by doing any of the following:
 - Click **Add...** to add a new J-modifier to the current spin definition.
 - Click **Edit...** to edit the selected J-modifier in the current spin definition.
 - Click **Delete** to delete the selected J-modifier in the current spin definition.

Tips and Tricks

- Double-click on a J-modifier to change its coupling constant and coupling status, including which spin definition the current one is coupled to, if any.
- Only the Spin Definition Tabs for the currently selected spin system are visible at any given time.
- The tabs for spin definitions that are coupled to the currently selected one appear highlighted in bold in the Simulation view.

Related Topics

- “Spin Systems” on page 70
- “Spin Definitions” on page 71

Menu Location	View > Sidebar > Simulation
Icon(s)	
Hotkey(s)	Ctrl-Shift-D Cmd-Shift-D

- “J-Modifiers” on page 73

Simulation Details

The bottom center pane in the Spin Simulator window displays simulation details, information about the current simulation. Simulation details include the simulation frequency, numbers of spin systems and spin definitions (spins), and numbers of transitions present in the current simulation.

Tips and Tricks

- When you do not have a spectrum overlay, you can change the simulation frequency in the simulation details to explore the effects of spectrometer frequency on your simulated compound.
- When you overlay a spectrum, the simulation frequency is set to match that of the spectrum overlay.

Related Topics


- “Spectrum Overlays” on page 19
- “Spin Systems” on page 70
- “Spin Definitions” on page 71

Chapter 8: Compound Builder

The Compound Builder module lets you create compound signatures so that you can identify and quantify user-defined compounds using Chenomx NMR Suite. You can add the compound signatures that you create with Compound Builder to user-defined compound sets using Library Manager, and use the resulting sets in Profiler to help you profile spectra.

You can only use Compound Builder if you have purchased the Professional edition of Chenomx NMR Suite, or if you are using the Evaluation edition of Chenomx NMR Suite.

8.1 Overview



- 1 Sidebar View Buttons. Let you change the currently visible Sidebar view.
- 2 Reference Panel. Displays additional information for the compound indicated by the Chenomx Compound ID in the Compound Details, including structure, molecular weight, CAS-RN and so on.
- 3 Cluster and Peak Lists. Let you edit, select or remove clusters or peaks in the currently open compound. Also lets you define new clusters based on the currently selected peaks or clusters.
- 4 Cluster Navigator. Lets you navigate among the clusters of the currently open compound. You can enable or disable the Cluster Navigator using the **View** menu.
- 5 Spectrum View Tools. Let you navigate the spectrum more easily.
- 6 Spectrum View. Displays the spectrum as well as the sum line and the currently selected cluster(s) or peak(s).
- 7 Information Panel. Displays information relating to the currently selected cluster(s) and peak(s).

- ⑧ Spectrum Thumbnail. Displays a thumbnail image of the Spectrum View. The thumbnail can also be used to navigate the spectrum at various zoom settings.
- ⑨ Status Bar. Displays the current cursor position in the units specified by the **View > X Scale** and **View > Y Scale** settings. Also displays information about recent user actions.

8.2 Quick Reference Guide

	Linux, Windows	Mac OS X
File Menu		
 New Compound...	Ctrl-N	Cmd-N
 Open...	Ctrl-O	Cmd-O
 Open Recent	---	---
 Close	Ctrl-W	Cmd-W
 Save	Ctrl-S	Cmd-S
 Save As...	---	---
 Import Simulation...	---	---
 Export > Spectrum Image...	---	---
 Compound Details...	---	---
 Preferences...	---	---
 Exit	---	---
Edit Menu		
 Undo	Ctrl-Z	Cmd-Z
 Redo	Ctrl-Y	Cmd-Y
 Copy Spectrum Image to Clipboard	---	---
Select All	Ctrl-A	Cmd-A
Select by Cluster	Ctrl-B	Cmd-B
View Menu		
X Scale	---	---
Y Scale	---	---
 Zoom > Show Entire Spectrum	Ctrl-Minus	Cmd-Minus
 Zoom > Undo Zoom	[	[
 Zoom > Redo Zoom	]	]

	Linux, Windows	Mac OS X
 Zoom > Set Zoom...	Ctrl-M	Cmd-M
 Zoom > Select Region	Ctrl-T	Cmd-T
 Zoom > Auto Zoom	Ctrl-F	Cmd-F
 Zoom > Increase Vertical Zoom	+	+
 Zoom > Decrease Vertical Zoom	-	-
 Zoom > Increase Horizontal Zoom	*	*
 Zoom > Decrease Horizontal Zoom	/	/
 Sidebar > Legend	Ctrl-Shift-L	Cmd-Shift-L
 Sidebar > Reference Panel	Ctrl-Shift-R	Cmd-Shift-R
Show Sidebar	Ctrl-Shift-S	Cmd-Shift-S
Show Spectrum Thumbnail	Ctrl-Shift-T	Cmd-Shift-T
Show Cluster Navigator	Ctrl-Shift-N	Cmd-Shift-N
Compound Menu		
 Add Peak	Ctrl-K	Cmd-K
 Edit Selected Peaks...	Ctrl-E	Cmd-E
 Delete Selected Peaks	Delete	Delete
 Group Selected Peaks as Cluster	Ctrl-G	Cmd-G
 Transform Window > Edit for Selected Cluster...	Ctrl-Shift-E	Cmd-Shift-E
 Transform Window > Set Default on Selected Clusters	---	---
 Transform Window > Clear for Selected Clusters	---	---
 Transform Window > Set Default on All Clusters	---	---
 Transform Window > Clear for All Clusters	---	---
 Previous Cluster	Q	Q
 Next Cluster	E	E
Tools Menu		
 Set Spectrum Overlay...	Ctrl-Shift-O	Cmd-Shift-O
 Clear Spectrum Overlay	---	---

	Linux, Windows	Mac OS X
 pH Sensitivities...	Ctrl-Shift-P	Cmd-Shift-P
 Manage Cluster IDs...	---	---
 Optimize Selected Peak Shapes	Ctrl-I	Cmd-I
 Generate Cluster for Region	Ctrl-Shift-I	Cmd-Shift-I
Help Menu		
 NMR Suite Help...	F1	F1
 Save Application Log...	Ctrl-Shift-A	Cmd-Shift-A
 Software Activation...	---	---
 About Chenomx NMR Suite...	---	---

8.3 Building Tools

Compound Builder includes a variety of tools to help you build compound signatures, by adding, manipulating and removing peaks, grouping peaks into clusters and defining ranges within which clusters can move. You can also generate arbitrary clusters to match the overall lineshape of a region, or optimize peaks to better match the spectrum.

New Compound

New Compound creates a new compound file (.xcpd). The new compound file has the default Compound Details, including a magnet frequency of 800 MHz and a sample concentration of 5 mM.

Menu Location	File > New Compound...
Icon(s)	
Hotkey(s)	Ctrl-N Cmd-N

Tips and Tricks

- New compound files do not have a spectrum overlay. Add a spectrum overlay before you begin build a compound signature.
- You can overlay one spectrum at a time in Compound Builder.
- You need to save the new file to retain changes that you have made.

Related Topics

- "Spectrum Overlays" on page 19
- "Compound Details" on page 91

Import Simulation

Importing lets you use an existing spin simulation (.xss) or legacy compound album (.cca) to guide building a new compound signature.

How do I import a spin simulation or compound album?

1. Open the Import Compound dialog.
2. Locate the folder containing the simulation or album file that you would like to import.

Menu Location	File > Import Simulation...
Icon(s)	

3. Select the simulation or album file and click **Open**.

4. Choose a calibration method.

If you are calibrating to a spectrum file (.cnx), click the Browse button (...) and locate the spectrum file you would like to use.

The spectrum that you choose is overlaid in the Spectrum View.

5. Click **OK**.

Tips and Tricks

- Calibrating to the Default Environment instead of a spectrum file will set a magnet frequency of 800 MHz, a pH of 7.00, and CSI settings corresponding to DSS at 0.500 mM (height of 1 si, width of 1 sw).
- The magnet frequency and spectrum pH for the current spectrum overlay appear in the status bar at the bottom of the Compound Builder window. If there is no spectrum overlay, the displayed values correspond to the default environment.
- You must calibrate a signature to a real spectrum for accurate results when using the signature in an analysis.

Related Topics

- “Spectrum Overlays” on page 19
- Chapter 7, *Spin Simulator* on page 67
- “Compound Details” on page 91

Add Peak

Add Peak lets you add a new peak to the current compound signature. The new peak appears in the peaks list, gets added to the currently selected cluster and becomes the currently selected peak in that cluster.

Menu Location	Compound > Add Peak
Icon(s)	 
Hotkey(s)	Ctrl-K Cmd-K

How do I add a peak to the current signature?

1. Switch to Add Peak mode. The mouse cursor changes to a large crosshair.
2. Click somewhere on the spectrum to add a new peak. The point that you click defines both the position and the height of the new peak.

Tips and Tricks

- To create a new cluster, you can deselect all clusters in the cluster list, then create a new peak. The new peak appears in its own new cluster.
- The new peak appears in the peak list in numerical order with the currently displayed peaks, not necessarily at the beginning or end of the list.
- Although the position that you click on the spectrum sets the position and height of the new peak, you are free to change both after you have added the peak.
- Right-click anywhere in the Spectrum View that does not already contain a peak and select **Add a New Peak to Cluster** or **Add a New Peak to New Cluster** to add a new peak to a specific cluster or to a new one, respectively.
- When you right-click in the Spectrum View to add a new peak, the specific cluster to which you can add a peak is the currently selected cluster if only one is selected, or the closest cluster to the point that you clicked if multiple clusters are selected.

Related Topics

- “Selecting Peaks and Clusters” on page 82
- “Delete Selected Peaks” on page 83
- “Adjusting Peaks” on page 83
- “Grouping Peaks as a Cluster” on page 84
- “Optimize Selected Peak Shapes” on page 85
- “Generate Cluster for Region” on page 86
- “Transform Windows” on page 86
- “Cluster Navigator” on page 90
- “Information Panel” on page 92

Selecting Peaks and Clusters

Many actions that you may perform in Compound Builder require you to select a peak or cluster to act on. There are several ways to select peaks, unassociated groups of peaks or clusters.

Clicking on the graph selects either peaks or clusters. You can toggle between selecting peaks or clusters via the graph itself, but the Information Panel always lets you select either.

How do I select a cluster?

1. Turn on selection by cluster.
2. Click on the cluster in the Spectrum View. Alternately, click on the cluster's position in the Information Panel, or double click on any peak in the cluster.

How do I select a peak?

1. Turn off selection by cluster.
2. Click on the peak in the Spectrum View. Alternately, click on the peak's position in the Information Panel.

How do I select multiple peaks or clusters?

1. Turn selection by cluster on or off, as needed.
2. Hold down the **Ctrl** or **Cmd** key, then click several peaks or clusters, either in the Spectrum View, the Information Panel, or both.

Tips and Tricks

- When you have a cluster currently selected, only the peaks associated with that cluster are available to select via the Information Panel.
- When you have no clusters currently selected, no peaks are visible or available to select via the Information Panel.
- If the peak that you would like to select is not visible in the Information Panel, you may need to select a different cluster.

Related Topics

- “Add Peak” on page 81

Menu Location	Select by Cluster
Hotkey(s)	Ctrl-B Cmd-B

- “Delete Selected Peaks” on page 83
- “Adjusting Peaks” on page 83
- “Grouping Peaks as a Cluster” on page 84
- “Optimize Selected Peak Shapes” on page 85
- “Generate Cluster for Region” on page 86
- “Transform Windows” on page 86
- “Cluster Navigator” on page 90
- “Information Panel” on page 92

Delete Selected Peaks

As you build a compound signature, you may need to delete one or more peaks from the signature to refine or simplify the signature. You can delete one or more selected peaks from the current signature, including peaks in different clusters, whole clusters or multiple clusters.

How do I delete peaks?

1. Select the peaks or clusters that you would like to delete.
2. Click Delete Selected Peaks.

Tips and Tricks

- You can delete multiple peaks or clusters by selecting them via the Spectrum View or the Information Panel.
- Right-click on a selection in the Spectrum View or Information Panel and select **Delete Peaks** to delete the selected peaks.

Related Topics

- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Information Panel” on page 92

Adjusting Peaks

In most cases, the initial location of a peak does not accurately match the spectrum overlay. You need to modify the height or frequency of one or more peaks to arrive at a signature that satisfies your requirements.

How do I set peak properties to specific values?

1. Select the peaks or clusters that you would like to adjust.
2. Open the Edit Selected Peaks dialog.
3. Enter the desired center (position in ppm), amplitude (height in si) and width (in sw) to apply to all selected peaks or clusters.
4. Click OK.

How do I adjust peak properties?

1. Select the peaks or clusters that you would like to adjust.

Menu Location	Compound > Delete Selected Peaks
Icon(s)	
Hotkey(s)	Delete Delete

Menu Location	Compound > Edit Selected Peaks...
Icon(s)	
Hotkey(s)	Ctrl-E Cmd-E

2. Do any of the following:

- Click and drag the selected peaks or clusters in the Spectrum View to adjust their positions and heights. Use the blue triangles on the horizontal and vertical axes for fine control over position and height, respectively.
- Click and drag the blue triangles on either side of the selected peaks and clusters, at approximately half-height, to adjust their widths.
- Use the W and S keys or the up and down arrow keys to adjust the height of the selection up and down.
- Use the A and D keys or the left and right arrow keys to adjust the position of the selection left and right.
- Use the Q and E keys to display the next or previous cluster.

Tips and Tricks

- When you enter specific amplitudes and widths in the Edit Selected Peaks dialog, all of the selected peaks are set to the specified values, regardless of which clusters they are associated with.
- When you enter a new center in the Edit Selected Peaks dialog, the geometric center of the selected peaks becomes the value that you entered
- Double-click on a peak or cluster in the Information Panel to open the Edit Selected Peaks dialog for the double-clicked item.
- Right-click on a selection in the Spectrum View or Information Panel and select **Edit Selected Peaks** to edit properties for the selection.
- When you adjust the properties for a single peak, you can use mathematical expressions instead of single values to set the center, height and width of the peak. To insert an expression instead of a single value, right-click in the text box that you would like to change, then click **Expression....**
- To view or edit an existing expression, click the expression icon () in the text box that you would like to change.

Related Topics

- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Grouping Peaks as a Cluster” on page 84
- “Optimize Selected Peak Shapes” on page 85
- “Generate Cluster for Region” on page 86
- “Transform Windows” on page 86
- “Cluster Navigator” on page 90
- “Information Panel” on page 92
- “Using Expressions” on page 135

Grouping Peaks as a Cluster

During the building process, and especially when you are building more complex compounds, you may need to group peaks from one or more clusters into a cluster of their own. Grouping

ensures that all of the selected peaks move together when you use the compound signature for profiling, and that they obey the same transform window settings.

How do I group peaks as a cluster?

1. Select the peaks and clusters that you would like to group as a new cluster.
2. Click Group Selected Peaks as Cluster.

Tips and Tricks

- Right-click on a selection in the Spectrum View or Information Panel and select **Group Selected Peaks as Cluster** to group the selected peaks directly.
- After grouping peaks into a new cluster, you may want to use the Edit Selected Peaks dialog to set them all to a common width.

Related Topics

- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Adjusting Peaks” on page 83
- “Transform Windows” on page 86
- “Cluster Navigator” on page 90
- “Information Panel” on page 92

Optimize Selected Peak Shapes

Adding peaks to a signature and adjusting them to approximately model a particular region of the spectrum is fairly simple. However, the task of fine-tuning those peaks to accurately match the lineshape in that region can be tedious. Once you have approximated the shape of a region of the spectrum, you can use Compound Builder's optimization algorithms to rapidly obtain a more refined match to that region.

How do I optimize the shape of peaks in my compound signature?

1. Select the peaks and clusters that you would like to optimize.
2. Click **Optimize Selected Peak Shapes**.

Tips and Tricks

- You must select fewer than 100 peaks to optimize.
- Optimization will generally give better results if your starting shape is reasonably close to the actual shape of the spectrum. Simulating the compound in Spin Simulator can often give you a very good starting shape.

Related Topics

- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Adjusting Peaks” on page 83
- “Generate Cluster for Region” on page 86

Menu Location	Compound > Group Selected Peaks as Cluster
Icon(s)	
Hotkey(s)	Ctrl-G Cmd-G

Menu Location	Tools > Optimize Selected Peak Shapes
Icon(s)	
Hotkey(s)	Ctrl-I Cmd-I

Generate Cluster for Region

You may need to build a compound signature that will let you verify that the compound is present in a spectrum, without necessarily providing the full range of flexibility that a chemically accurate signature affords. You can create these limited-use signatures using automatically generated clusters. Generated clusters lack the chemical accuracy of signatures that are backed by spin simulations, but provide a quick model of the lineshape of a compound under specific conditions.

How do I generate a cluster to match a region of the spectrum?

1. Select the region of the spectrum that you would like to match with a generated cluster.
2. Click **Generate Cluster for Region**.

Tips and Tricks

- Generated clusters have no formal theoretical basis. They simply add a number of Lorentzian peaks that approximate the overall shape of the spectrum in the selected region. The generated peaks will not let you accurately account for significant changes in the spectrum based on magnet frequency, pH or other systematic variations.
- Signatures based on generated clusters may only be expected to give accurate results when profiling spectra acquired under conditions very similar to the source spectrum.
- Generated clusters are already optimized. Further optimizations will not generally give any better results.
- When you generate a cluster, Compound Builder will only use up to 100 peaks to match the selected region. If the region you are trying to match cannot be accurately modeled using 100 peaks, try selecting a smaller region, adding more peaks manually, or simulating the compound in Spin Simulator and importing the simulation.

Related Topics

- “Select Region” on page 22
- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Adjusting Peaks” on page 83
- “Optimize Selected Peak Shapes” on page 85

Transform Windows

You can limit the amount that you can transform (or move) individual clusters in a compound signature in Profiler using transform windows. Establishing transform windows for your own compounds improves the quality of profiles that you perform with your custom compound signatures.

If you set transform windows manually, ensure that they are wide enough to account for the degree of variability in cluster position that you expect to see in your samples. For example, if you routinely analyze spectra with widely varying pH, you will need larger transform windows than if your sample spectra are always close to the same pH.

The most common influences on the position of a cluster in a compound are pH and ion concentrations. To manually set accurate transform windows for a compound, consider acquiring spectra of the pure compound at the outer limits of the conditions you expect to

Menu Location	Tools > Generate Cluster for Region
Icon(s)	
Hotkey(s)	Ctrl-Shift-I Cmd-Shift-I

encounter when profiling the compounds in experimental mixtures. The additional spectra let you map the maximum extent of these effects on the cluster positions in your compounds.

In the absence of specific information about the behavior of a cluster, you can set a default transform window for the cluster. The default transform window is 0.006 ppm wide, centered on the current cluster center. For a cluster appearing at 1.000 ppm, the default transform window sets lower and upper limits of 0.997 and 1.003 ppm, respectively.

Menu Location	Tools > pH Sensitivities...
Icon(s)	

How do I set pH sensitive transform windows?

1. On the Tools menu, click **pH Sensitivities....**
2. Click **Load Data Points....**
3. Select the signatures that you would like to use to describe pH sensitive centers and transform windows for the compound.
4. Click **Set Data Points.**
5. Click **Optimize Lines.**
6. Review the clusters in the compound. Click **Accept** to apply the automatically generated expressions for cluster centers and transform windows to your compound..

Menu Location	Compound > Transform Window > Edit for Selected Cluster...
Icon(s)	

How do I set the transform window for a cluster?

1. Select the cluster for which you would like to set a transform window.
2. Open the Edit Transform Window dialog.
3. Select **Limit transforms to**, and enter the new lower and upper limits for the selected cluster.
4. Click OK.

How do I change the transform window for a cluster?

1. Select the cluster whose transform window you would like to change.
2. Open the Edit Transform Window dialog.
3. Enter the new lower and upper limits for the selected cluster.
4. Click OK.

How do I set the default transform window for a cluster?

1. Select the cluster for which you would like to set the transform window.
2. On the Compound menu, select **Transform Window > Set Default on Selected Cluster.**

Menu Location	Compound > Transform Window > Set Default on Selected Clusters
Icon(s)	

How do I set the default transform window for all clusters?

1. On the Compound menu, select **Transform Window > Set Default on All Clusters.**

Menu Location	Compound > Transform Window > Set Default on All Clusters
Icon(s)	

Menu Location	Compound > Transform Window > Clear for Selected Clusters
Icon(s)	

Menu Location	Compound > Transform Window > Clear for All Clusters
Icon(s)	

How do I clear the transform window for a cluster?

1. Select the cluster for which you would like to clear the transform window.
2. On the Compound menu, select **Transform Window > Clear for Selected Clusters**.

How do I clear the transform windows for all clusters?

1. On the Compound menu, select **Transform Window > Clear for All Clusters**.

Tips and Tricks

- When you set the transform windows for a cluster, you can use mathematical expressions instead of constants to set the maximum and minimum values. A special icon in the text box for a transform window () indicates that you may use an expression to set the transform window.
- When you create pH sensitive compounds, expressions for the transform windows are calculated automatically, based on the set of spectra that you select.
- A transform window applies to the center of a cluster, not the positions of its outermost peaks. Consider only cluster centers when you manually set transform windows for your compounds.
- When you manually set a transform window for a cluster, the cluster center must lie within the selected limits.
- After you have set a transform window for a cluster, if you try to move the cluster beyond either limit of the transform window, the limit moves to accommodate the new position of the cluster. If the transform window was an expression, adjusting it in this manner will replace the expression with a constant.
- You may need to acquire multiple spectra to determine appropriate transform windows for a compound. Remember that some compounds, and some clusters within a compound, are more sensitive than others to pH changes and other matrix effects.
- Setting default transform windows and clearing transform windows, described above for single clusters, can also apply to multiple selected clusters.

Related Topics

- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Adjusting Peaks” on page 83
- “Grouping Peaks as a Cluster” on page 84
- “Cluster Navigator” on page 90
- “Using Expressions” on page 135

pH Sensitivities

You can make compound signatures that you create in Compound Builder sensitive to the pH of the spectrum that you are profiling. The clusters of pH sensitive compounds will adapt their centers and transform windows based on the spectrum pH, ensuring that they start closer to their actual position in the spectrum and simplifying compound identification in complex regions of the spectrum.

Menu Location	Tools > pH Sensitivities...
Icon(s)	
Hotkey(s)	Ctrl-Shift-P Cmd-Shift-P

How do I make a compound signature pH sensitive?

1. Make sure that you have a set of signatures available to be used as data points on your pH curves. They must satisfy the following criteria:
 - All of the signatures must have the same number of clusters, in the same order.
 - Each signature must have its valid pH set to the pH value from the signature's originating sample, via the Compound Details dialog. Use a single value instead of a range.
2. Select a signature to use as a 'master' signature. This signature will become the pH sensitive signature.
3. Open the master signature file (.xcpd) in Compound Builder, and overlay the associated spectrum file (.cnx).
4. On the Tools menu, click **pH Sensitivities....**
5. Click **Load Data Points....**
6. Click **Optimize Lines**.
7. Review the resulting pH curves for each cluster by clicking the various tabs at the top of the pH Sensitivities dialog.
8. Click **Accept** to apply the automatically generated cluster center and transform window expressions to the clusters in your master signature, and set the valid pH of the master signature to the range defined by the pH values of your individual data points.

Tips and Tricks

- Clusters that appear green in the Cluster Navigator are pH sensitive, meaning that one or more of their heights, positions, and widths change with spectrum pH.
- If the clusters in your set of signatures appear in different orders in some of the files, use Manage Cluster IDs to rearrange the cluster IDs in the signatures so that they match.
- If you manually adjust any cluster after applying pH sensitivities to the compound, some or all of the automatically generated expressions may be replaced with constant values.
- If the signatures you use as data points are at a different spectrometer frequency than your master signature, you may need to apply an offset to ensure accurate pH curves. Click **Lock to Current Clusters** in the pH Sensitivities dialog to calculate the offset before reviewing and accepting the pH curves.
- If the pH of your master signature is close to a pKa of the compound, using **Lock to Current Clusters** may introduce unintended errors. Review the suggested pH curves carefully.

Related Topics

- "Adjusting Peaks" on page 83
- "Transform Windows" on page 86
- "Manage Cluster IDs" on page 90
- "Compound Details" on page 91

Manage Cluster IDs

Menu Location	Tools > Manage Cluster IDs...
Icon(s)	

How do I change the order of cluster IDs in a signature?

1. On the Tools menu, click **Manage Cluster IDs...**
2. Click and drag the clusters in the list to change their order.

Or

Select a cluster in the list and use the green up and down arrows to change the position of the cluster in the list.

3. Click OK.

Tips and Tricks

- Cluster IDs are used by the pH sensitivities dialog to determine which clusters should be considered as corresponding to the same signal across all of your selected data points. Make sure that the clusters appear in the proper order in each of your data points before applying pH sensitivities.
- A cluster that appears upfield from another cluster at one pH may appear downfield from the same cluster at a different pH. Consider overlaying all of your data point spectra (.cnx) in Processor or Profiler to better determine the behavior of clusters across the full pH range.
- If you have used the signature to profile spectra, changing the cluster IDs may make your existing fits of the compound invalid.

Related Topics

- “pH Sensitivities” on page 88

Cluster Navigator

The Cluster Navigator lets you quickly move among the clusters that you have created in the current compound signature. It appears in the top left corner of the Spectrum View as an ordered list of the approximate positions of all the clusters of the compound in ppm. When you click on a cluster position, the view zooms to provide an optimal view of the cluster. The optimal view includes the larger of the full transform range of the cluster or the full width of the cluster, and the full vertical extent of any peaks within the horizontal range. You can also use the Cluster Navigator to move sequentially through the list of clusters in either direction.

How do I display the clusters of a compound in the Spectrum View?

1. Make sure that the Cluster Navigator is visible.
2. Click a cluster position to zoom directly to that cluster. Click the left or right arrow to move sequentially through the list of clusters.

Tips and Tricks

- Clusters that appear green in the Cluster Navigator are pH sensitive, meaning that one or more of their heights, positions, and widths change with spectrum pH.
- Click on a cluster and drag along the Cluster Navigator to include adjacent clusters in the Spectrum View.

Menu Location	View > Show Cluster Navigator
Hotkey(s)	Ctrl-Shift-N Cmd-Shift-N

Related Topics

- “Add Peak” on page 81
- “Selecting Peaks and Clusters” on page 82
- “Adjusting Peaks” on page 83
- “Grouping Peaks as a Cluster” on page 84
- “Transform Windows” on page 86

8.4 Information Tools

While building compound signatures, you can access extended information about common compounds, view and update details associated with the current compound, and more precisely manage peaks and clusters in complex compounds.

Reference Panel (Sidebar)

Menu Location	View > Sidebar > Reference Panel
Icon(s)	
Hotkey(s)	Ctrl-Shift-R Cmd-Shift-R

The Reference Panel Sidebar view displays information specific to the compound matching the Chenomx Compound ID that you have specified in the Compound Details. The displayed information includes a structure image, molecular formula and weight, alternate names and links to external reference sites. You can use this information to better understand the structure and behaviour of the compound as you build its signature.

Related Topics

- “Reference Panel” on page 18
- “Preferences” on page 25
- “Building Tools” on page 80
- “Compound Details” on page 91

Compound Details

Menu Location	File > Compound Details...
Icon(s)	

Using Compound Details, you can associate information with compound signatures that you build. The information that you enter lets you use the signature with appropriate spectra, display correct information in the Reference Panel, calculate accurate concentrations, and manage the contents of your compound library effectively.

Table 8.1. Compound Details

General	
Compound Name	Name of the compound, as you would like it to appear in the Compound Table in Profiler
Compound Author	Author of the compound signature, as you would like it to appear in the Compound Table in Profiler
Magnet Frequency (MHz)	Frequency of the acquired NMR spectrum, in megahertz (MHz)
Sample Concentration (mM)	Concentration of the compound in the sample spectrum, in millimoles per litre (mM)
Last Modified	Date and time at which the signature file was last modified (and saved).

UUID	Universally Unique Identifier used to uniquely identify a particular signature in your compound library.
Reference	
Chenomx Compound ID	Chenomx compound ID associated with the compound
InChI Code	IUPAC International Chemical Identifier (InChI) associated with the compound (http://www.iupac.org/inchi/)
SMILES Code	Simplified Molecular Input Line Entry Specification (SMILES) associated with the compound
Advanced	
Molecular Weight (u)	Molecular weight of the compound in atomic mass units, not including counter-ions
Pulse Sequence	Pulse sequence used to acquire the sample spectrum that you would like to use as a reference (overlay) for the compound signature
Valid pH	pH value or range at which you expect the compound signature to be valid
Comment	
Comment	Text comments associated with the compound signature; you can use this to record any additional information that you consider relevant

Tips and Tricks

- For accurate quantification results, consider verifying the sample concentration that you enter using multiple independent techniques.
- Entering the correct Chenomx Compound ID ensures that information appropriate to the compound appears in the Reference Panel.
- If you intend to use concentration units of milligrams per decilitre (mg/dL) in Profiler, enter a molecular weight.
- Chenomx calculates molecular weights based on IUPAC atomic weights [<http://www.chem.qmul.ac.uk/iupac/AtWt/>], and rounds the result to four decimal places.
- You need to save the compound file to retain new or changed Compound Details.

Related Topics

- “Reference Panel” on page 18
- “Building Tools” on page 80
- “Reference Panel (Sidebar)” on page 91
- Chapter 9, *Library Manager* on page 94

Information Panel

The Information Panel is always visible at the bottom of the Compound Builder window. It contains current information about the compound signature as you work on it, listing clusters and peaks that you have created, describing the elements that you have currently selected and the environment on which your compound is based, and letting you access the

edit dialog for any clusters or peaks you may have selected. You can also select peaks and clusters directly in the Information Panel, with more precision than selecting them from the Spectrum View.

Tips and Tricks

- You can select multiple peaks and clusters in the lists in the Information Panel. Hold down the **Ctrl** or **Cmd** keys while clicking on the desired peaks and clusters.
- If you do not have a spectrum overlay, the Information Panel indicates that your compound signature is in the default environment. When you overlay a spectrum, the environment changes to match the overlay.
- The Information Panel displays a warning if your spectrum overlay does not match the information that you provided in the Compound Details. You need to update the Compound Details to match the overlay or choose a new overlay spectrum to remove the warning.


Related Topics

- “Spectrum Overlays” on page 19
- “Building Tools” on page 80
- “Compound Details” on page 91

Chapter 9: Library Manager

The Library Manager module lets you create and modify compound sets based on compounds in your library, as well as compounds that you have created in Compound Builder. You can use compound sets in Profiler to focus your profiling efforts on specific compounds, based on dynamic filters or predefined lists of compounds.

9.1 Overview



- 1 Compound Sets. Sidebar view letting you control which compound set appears in the Compound Table.
- 2 Compound Table. Displays information about compounds and their status in the library. You can send the signature of a compound to Compound Builder for editing by clicking **File > Send to Compound Builder** or by double-clicking on the compound.
- 3 Quick Search. Lets you filter the Compound Table to focus on a particular compound or group of compounds.
- 4 Status Bar. Displays information about recent user actions.

9.2 Quick Reference Guide

	Linux, Windows	Mac OS X
File Menu		
 New Compound Set...	Ctrl-N	Cmd-N
 New Smart Compound Set	Ctrl-Shift-N	Cmd-Shift-N
 Edit Smart Compound Set...	Ctrl-E	Cmd-E
 Rename Compound Set	Ctrl-R	Cmd-R
 Remove Compound Set	Delete	Delete
 Send to Compound Builder	Ctrl-J	Cmd-J
 Preferences...	---	---
 Exit	---	---
Edit Menu		
 Undo	Ctrl-Z	Cmd-Z
 Redo	Ctrl-Y	Cmd-Y
Select All	Ctrl-A	Cmd-A
View Menu		
 Change Columns...	---	---
Library Menu		
 Add Compounds...	---	---
Selected Compounds Submenu		
 Add to [Compound Set]	---	---
 Remove	Delete	Delete
 Export to > Compound Pack (.pack)...	---	---
 Export to > Compound Files (.xcpd)...	---	---
Help Menu		
 NMR Suite Help...	F1	F1
 Save Application Log...	Ctrl-Shift-A	Cmd-Shift-A
 Software Activation...	---	---
 About Chenomx NMR Suite...	---	---

9.3 Library Tools

Library Manager lets you browse and manage the contents of your library. Your library contains all of the compounds that you can use to profile a spectrum, at all magnet frequencies, pH ranges, and pulse sequences.

Compound Table

The Compound Table in Library Manager serves a different purpose from the one in Profiler. The most obvious difference is that concentrations (maximum or otherwise) are not available in Library Manager, since you are not profiling spectra. More generally, the Compound Table displays compounds in your library and lets you modify them.

When you first open Library Manager, the Compound Table displays all of the compounds in your library. As you select Compound Sets, the Compound Table updates to reflect the contents of the currently-selected set. You can use the Compound Table to help you build Compound Sets, via drag-and-drop and context menus. You can also use text-based Quick Searches to filter any list currently visible in the Compound Table.

There are many columns available for display in the Compound Table. You can choose to display any number of the available columns, with a minimum of one visible at any time, and you can drag the column headers to position the columns in any order that you like. You can sort the Compound Table on any column by clicking its column header. Click the header again to reverse the sort order.

Tips and Tricks

- You can add more compounds to a Compound Set by dragging them from the Compound Table onto the name of the Compound Set, or by right-clicking on the compound and selecting **Add to New Compound Set** or **Add to**. You cannot add compounds to Smart Compound Sets directly.
- You can edit a compound in the Compound Table by right-clicking on its entry and selecting **Send to Compound Builder**. As soon as you save your changes in Compound Builder, they are available in all Compound Sets containing the compound.
- You can also double-click a compound in the Compound Table to send it to Compound Builder.
-

Related Topics

- “Change Columns” on page 16
- “Quick Searches” on page 96
- “Send to Compound Builder” on page 99
- “Remove Selected Compounds” on page 99
- “Compound Set Tools” on page 100

Quick Searches

The number of compounds appearing in the Compound Table can make it difficult to locate a particular compound. Quick Searches let you temporarily limit the number of compounds appearing in the Compound Table based on search criteria. The search criteria can include text searches based on name, compound IDs, comments or all of these.

Icon(s)	
Hotkey(s)	Ctrl-L Cmd-L

How do I limit the compounds visible in the Compound Table?

1. Click the binoculars icon at the left edge of the search box to choose the type of search to use (Compound Name, Compound Comment, Chenomx Compound ID or All). If you do not select something here, the default search type is All.
2. Type your search terms in the Quick Search box at the top left of the Compound Table. The contents of the Compound Table updates to show all compounds matching your current entry.
3. Click the clear button at the right edge of the search box to cancel the search and restore the full list of compounds. Clearing a search also resets the search type to All.

Tips and Tricks

- Use Quick Searches to temporarily limit the number of compounds displayed. Use compound sets to define a more permanent working set of compounds for profiling a particular spectrum.

Related Topics

- “Compound Table” on page 96
- “Remove Selected Compounds” on page 99
- “Compound Set Tools” on page 100

Add Compounds

You can add compounds to your library from a variety of sources, including compound files, compound packs, Chenomx spectra, old Chenomx library files and ZIP (compressed) archives. Add Compounds lets you add compounds from all of these sources to your library, making the compound signatures available in Profiler.

Menu Location	Compounds > Add...
Icon(s)	

How do I add Chenomx Reference Compounds to my library?

1. Open the Add Compounds wizard.
2. Select Chenomx Reference Compounds, then click Next.
3. Select the frequencies for which you would like to add Chenomx Reference Compounds.
4. Click Finish to add the selected compounds.

How do I add other compounds to my library?

1. Open the Add Compounds wizard.
2. Select Importable Files, then click the Browse button (...).
3. Select any combination of compound files (.xcpd), compound packs (.pack), spectrum files (.cnx), old Chenomx library files (.cnxlib), ZIP files (.zip) or folders containing any of these file types, then click Choose Files.
4. Click Next.
5. If you like, choose how you would like compound sets created.
6. Click Finish to add the selected compounds.

Tips and Tricks

- When you add compounds, any compound sets that were created or modified during the addition process 'throb' for a couple seconds to indicate which sets were affected by the import.
- The addition process will not create duplicate Compound Sets. If the addition options that you selected would normally create a Compound Set with the same name as an existing set, the existing set will be updated with the added compounds.

Related Topics

- "Compound Table" on page 96
- "Compound Set Tools" on page 100

Export Compounds

When you collaborate with other people on a project, you may have different compounds in your library than your collaborator has in theirs. Exporting compounds allows you to exchange compound signatures with other users of Chenomx NMR Suite. You can then be certain that you and your collaborator are always fitting spectra for that project using the same set of compound signatures.

You can export compounds as a single compound pack file (.pack), or as collections of individual signature files (.xcpd).

How do I export compounds to a compound pack?

1. Select the compounds or compound set that you would like to export.
2. Open the **Export Compound Pack** dialog.
3. Enter a name for the exported compound pack. The entered name is used to label the associated compound set when the compound pack is added to another library via Library Manager.
4. Click **OK**.
5. In the resulting file chooser, locate the folder in which you would like to save the compound pack.
6. In the File Name field, type the new name for the compound pack file. The specified name can be different from the name of the compound pack itself, but you may prefer to keep them the same.
7. Click **Save**.

Menu Location	Compounds > Export Selected Compounds to > Compound Pack (.pack)...
Icon(s)	

How do I export compounds to individual compound files?

1. Select the compounds or compound set that you would like to export.
2. On the Compounds menu, click **Export Selected Compounds to > Compound Files (.xcpd)...**
3. In the resulting file chooser, select the folder in which you would like to save the compound files.
4. Click **Select**.

Menu Location	Compounds > Export Selected Compounds to > Compound Files (.xcpd)...
Icon(s)	

Tips and Tricks

- You and your collaborator each need a valid licensed copy or evaluation edition of Chenomx NMR Suite, version 5.0 or later, installed to make use of compound packs and signature files.
- If you only want to export compounds that were actually profiled in a particular spectrum, load the spectrum in Profiler and use **Export Profiled Compounds**.

Related Topics

- “Add Compounds” on page 97
- “Compound Set Tools” on page 100

Send to Compound Builder

Library Manager can help you to manage building and maintaining of your custom compound signatures. You can send any compound visible in the Compound Table to Compound Builder for modifications. When you save your changes, you can add the modified compound to your library.

How do I send a compound to Compound Builder?

1. Select a compound in the Compound Table.
2. Click **Send to Compound Builder**.
3. Make any desired changes to the compound signature in Compound Builder.
4. Save the changes in Compound Builder.
5. Close Compound Builder, or close the file.

Menu Location	File > Send to Compound Builder
Icon(s)	
Hotkey(s)	Ctrl-J Cmd-J

Tips and Tricks

- Make sure that you have the Version column turned on in the Compound Table, so you can see which compounds you have modified recently.
- You can send multiple compounds to Compound Builder by selecting them in the Compound Table, right-clicking on your selection, and selecting **Send to Compound Builder**.
- You can only open five Compound Builder windows at the same time. If you try to send more than five compounds to Compound Builder, only five of them open, and an error message appears indicating that the remaining compounds could not be opened.

Related Topics

- Chapter 8, *Compound Builder* on page 77
- “Compound Table” on page 96
- “Quick Searches” on page 96
- “Compound Set Tools” on page 100

Remove Selected Compounds

If you are certain that you no longer want a compound accessible in Chenomx NMR Suite, you can remove it from your library.

Menu Location	Library > Selected Compounds > Remove
Icon(s)	
Hotkey(s)	Delete Delete

How do I remove compounds from my library?

1. Select the compound that you would like to remove in the Compound Table.
2. Click **Remove Selected Compounds**.

Tips and Tricks

- When you have removed compounds from the library, you cannot recover them through Undo. If you still have the original compound files (.xcpd) or compound packs (.pack) available, you can add them again.

Related Topics

- “Compound Table” on page 96
- “Quick Searches” on page 96
- “Compound Set Tools” on page 100

9.4 Compound Set Tools

Compound Sets are groups containing compounds from your library that you can use in Profiler. You can create them either by selecting specific compounds (Compound Sets) or by defining rules for selecting compounds (Smart Compound Sets). Library Manager includes a number of tools to help you browse and manage compound sets.

Compound Sets (Sidebar)

The Compound Sets view in the Sidebar contains a list of the currently defined compound sets and lets you interact with them by displaying their contents, or editing them. The Compound Sets view also contains the Compounds set, letting you review the complete contents of your library.

How do I display the contents of a compound set?

1. Select the set in the Compound Sets view. The contents of the set appear in the Compound Table.

Tips and Tricks

- You can drag compounds onto a Compound Set in the Compound Sets view to add them to the set. You cannot drag compounds onto Smart Compound Sets.
- You can rename compound sets directly in the Compound Sets view. Click on the name of a selected set, type the new name, and press **Enter** to accept the new name.
- Several operations involving compound sets are available by right-clicking on them in the Compound Sets view.

Related Topics

- “Library Tools” on page 96
- “Export Compounds” on page 98
- “Compound Sets” on page 101
- “Smart Compound Sets” on page 102
- “Rename Compound Sets” on page 103

- “Remove Compound Sets” on page 103

Compound Sets

You may need to create a compound set that does not change as you add new compounds to your library. For example, analysis of an extended dataset may require profiling a fixed set of compounds to help ensure a consistent analysis, or you may want to profile only specific compounds, and cannot easily define rules that select the compounds you want. Compound Sets let you include arbitrary compounds on a per-compound basis, hand-picking the members of the set to suit your needs.

Although Compound Sets do not automatically include new compounds that you add to your library, some changes to the library still affect Compound Sets. If you remove a compound from the library, or replace a compound with a new version, the removal or update affects both Compound Sets and Smart Compound Sets.

How do I create an new, empty Compound Set?

1. Click **New Compound Set**.
2. Enter a name for the new Compound Set.
3. Start adding compounds from your library or other compound sets.

Menu Location	File > New Compound Set...
Icon(s)	
Hotkey(s)	Ctrl-N Cmd-N

How do I create a new Compound Set containing specific compounds?

1. Select a compound or group of compounds in the Compound Table.
2. Click **New Compound Set**.
3. Enter a name for the new Compound Set.

How do I add compounds to an existing Compound Set?

1. Select a compound or group of compounds in the Compound Table.
2. Click (e.g.) **Add to My Compound Set** to add the compounds to an existing set called "My Compound Set."

Alternately, drag the selected compounds onto the name of the Compound Set in the Compound Sets view.

How do I remove compounds from a Compound Set?

1. Select the Compound Set from which you would like to remove compounds in the Compound Sets view.
2. Select the compounds that you would like to remove in the Compound Table.
3. Click **Remove Selected Compounds from Compound Set**.

Tips and Tricks

- When you change the contents of any compound set, its name and icon in the Sidebar 'throbs' for a couple seconds to indicate which set has changed. If you have made a change that affects multiple compound sets, all of the affected sets will throb.
- When you add new compounds, you can choose to create Compound Sets containing the newly added compounds. If the Compound Sets that would normally be created already exist, the new compounds are added to the existing sets.

Menu Location	Compounds > Add Selected Compounds to [Compound Set]
---------------	--

Menu Location	Compounds > Remove Selected Compounds from > Compound Set
Hotkey(s)	Delete

- Right-click on a selected compound or group of compounds in the Compound Table, and click **Add to New Compound Set** to create a new set containing the selected compounds.
- You can only remove compounds from Compound Sets, and not Smart Compound Sets.

Related Topics

- “Library Tools” on page 96
- “Export Compounds” on page 98
- “Compound Sets (Sidebar)” on page 100
- “Smart Compound Sets” on page 102
- “Rename Compound Sets” on page 103
- “Remove Compound Sets” on page 103

Smart Compound Sets

In Library Manager, you can create Smart Compound Sets based on sets of rules that select compounds from your library. Smart Compound Sets automatically update to reflect changes in your library. If you add new compounds to your library, they immediately appear in any Smart Compound Sets with rules that match the new compounds.

Simple rules might be *"all compounds that I can use to profile 800 MHz spectra"* or *"all compounds with 'acetate' in their names."* However, you can also use Smart Compound Sets to select *"all compounds useful for profiling 600 MHz spectra, with a valid pH range from 5.54 to 8.39, added to the library after January 1, 2010, with a molecular weight less than 300 g/mol and whose name starts with the letter 'a'."* You can be as specific or as general as you like. You can also quickly merge results from existing Compound Sets, selecting *"compounds that are in Set A, but not in Set D."*

Smart Compound Sets can help you focus on just the compounds that you are interested in, without having to worry about updating the sets as you add new compounds to your library or update existing compounds.

How do I create a new Smart Compound Set?

1. Open the New Smart Compound Set dialog.
2. Enter a name for the new Smart Compound Set.
3. Enter the desired magnet frequency for the Smart Compound Set.
4. Use the drop-down menus to define any additional search filters that you would like to apply to the Smart Compound Set. Use the + buttons to add new filters.
5. Click OK.

How do I edit a Smart Compound Set?

1. Select a Smart Compound Set in the Compound Sets view.
2. Click **Edit Smart Compound Set**.
3. Make any desired changes to the name of the set, magnet frequency restriction, or filters.
4. Click OK.

Menu Location	File > New Smart Compound Set...
Icon(s)	
Hotkey(s)	Ctrl-Shift-N Cmd-Shift-N

Menu Location	File > Edit Smart Compound Set...
Icon(s)	
Hotkey(s)	Ctrl-E Cmd-E

Tips and Tricks

- You cannot manually add compounds to a Smart Compound Set. Only compounds in your library that match the filters that you define in the Smart Compound Set appear in the set. Use **Edit Smart Compound Set** to modify the compounds appearing in a Smart Compound Set.
- When you change the contents of any compound set, its name and icon in the Sidebar 'throbs' for a couple seconds to indicate which set has changed. If you have made a change that affects multiple compound sets, all of the affected sets will throb.
- Smart Compound Sets update automatically as the contents of your library change. When you add a new compound to your library, it appears in any Smart Compound Set whose filters it matches, even if you created the Smart Compound Set before you added the new compound.
- Editing the filters of a Smart Compound Set is the only way that you can modify the contents of a Smart Compound Set. You cannot directly add compounds to or delete compounds from a Smart Compound Set.

Related Topics

- "Library Tools" on page 96
- "Compound Sets (Sidebar)" on page 100
- "Compound Sets" on page 101
- "Rename Compound Sets" on page 103
- "Remove Compound Sets" on page 103

Rename Compound Sets

You can rename Compound Sets or Smart Compound Sets after you have created them.

How do I rename a compound set?

1. Select a compound set in the Compound Sets view.
2. Click **Rename Compound Set**.
3. Enter a new name for the set.
4. Press **Enter**.

Menu Location	File > Rename Compound Set...
Icon(s)	
Hotkey(s)	Ctrl-R Cmd-R

Tips and Tricks

- Click on the name of a selected compound set, or right-click on the name of a compound set and select **Rename** to rename the set.
- You can rename a Smart Compound Set via the **Edit Smart Compound Set** dialog.

Related Topics

- "Compound Sets (Sidebar)" on page 100
- "Compound Sets" on page 101
- "Smart Compound Sets" on page 102

Remove Compound Sets

You can remove compound sets when you no longer need them.

Menu Location	File > Remove Compound Set...
Icon(s)	

How do I remove a compound set?

1. Select a compound set in the Compound Sets view.
2. Click **Remove Compound Set**.

Tips and Tricks

- Right-click on the name of a compound set and select **Remove** to remove the set.

Related Topics

- “Compound Sets (Sidebar)” on page 100
- “Compound Sets” on page 101
- “Smart Compound Sets” on page 102

Chapter 10: Tutorials

Whether you are a brand new user of Chenomx NMR Suite or you have just upgraded from an earlier version, the tutorials in this chapter can help you use any of the available modules more effectively.

10.1 File Location Note

You can find all files mentioned in these tutorials in the `Samples\` folder in your Chenomx NMR Suite install folder

10.2 Processor

Introduction

You can use the Processor module to process your spectra in preparation for using them in Profiler, Spin Simulator or Compound Builder. Before processing your spectra, add spectrum metadata. The most important metadata to add is a CSI definition. Processor adds some spectrum metadata based on the CSI type and concentration that you specify in the Import NMR Data dialog when you import the spectrum. However, Profiler, Spin Simulator and Compound Builder use the CSI to determine lineshapes for any fit or simulation that you do, so it is very important to define the CSI shape properly.

Once you have added spectrum metadata, you can begin processing spectra. There are a variety of processing effects available, and it is possible to apply them in any order that you prefer. However, for the best results in any module, we recommend the following processing order:

1. Phasing

Unless you need to reverse the spectrum, phasing should always be your first processing step.

2. Baseline Correction (2-point linear or cubic splines)

Even if your spectrum appears to have a fairly flat baseline, consider applying a simple baseline correction layer, to ensure that the spectrum baseline aligns with that of the signatures in your library.

3. Reference Deconvolution (if necessary)

You will not need to apply reference deconvolution to every spectrum, but when you do, make sure that you have properly adjusted the phasing and baseline correction first.

4. Line Broadening (0 to 1000 Hz; if necessary)

Line broadening is available, but many ^1H NMR spectra will not need it. If you have CSI linewidths of 1 to 1.5 Hz or higher, do not apply line broadening.

5. Water Deletion

Variation in the performance of solvent suppression pulse sequences can produce water peaks that are still much larger than all the other signals in the spectrum. When you do not expect any signals of interest very close to the water peak

(4.68-4.88 ppm), deleting the water peak can greatly reduce the amount of zooming in you will need to do while profiling the spectrum.

Processing a Spectrum

Getting Started

1. Open Processor.
2. Open the spectrum `Sample.fid\fid`.
3. In the Import Varian Raw FID Dialog:
 - a. Select DSS as the Chemical Shape Indicator, and set its concentration to 0.484 mM.

Setting an accurate CSI concentration ensures that quantification in Profiler is as accurate as possible.

When you import your own spectra, use a CSI concentration appropriate to your samples.

- b. Expand the **Advanced Details** section.
- c. Set **Zero Fill to** Automatic, if it is not set already.

Zero-filling maximizes the visible resolution of the data that you collected, allowing the most accurate fits possible with your data.

- d. Clear the **Line Broadening** check box.

Unless your CSI linewidths are significantly less than 1 Hz, you do not need to add line broadening.

- e. Clear the **Phasing** check box.

For this tutorial, we will do manual phasing. Since multiple phasing steps are detrimental to spectrum quality, turn off the automatic phasing. This also disables the option to import phasing angles set in VNMR, which should only be retained if you know they are accurate.

- f. Clear the **Auto Baseline Correction** check box.

For this tutorial, we will manually correct the baseline.

- g. Clear the **Delete Water Region** check box.

For this tutorial, we would like to see what the spectrum looks like before the water region is deleted.

- h. A warning appears stating that Processor could not determine the pH automatically. For this tutorial, we will manually measure the spectrum pH.

Click the **OK** button.

Adjusting the Spectrum Phasing

4. Switch to Phasing mode.
5. Click the **Auto** button.

A good first step to manual phasing is to apply the automatic phasing, which often provides a reasonable starting point. For the best results, visually inspect the automatic results and refine them as needed.

6. Click on **Very Fine** to make small adjustments to the phasing. Drag the sliders for zero-order and first-order phase angles until most peaks in the spectrum appear symmetrical.

The water peak often does not phase properly with the rest of the spectrum, so focus on other signals in the spectrum instead. Phase angles for the sample spectrum are near 54° (zero-order) and -2° (first-order).

7. Click the **Accept** button to apply the changes to the spectrum and exit Phasing mode.

Applying Baseline Correction

8. Switch to Baseline Correction mode.
9. Click the **Auto Linear** button to have Processor locate the spectrum endpoints automatically.

The automatic linear baseline correction provides a useful starting point for the more involved spline correction.

10. Click the **Auto Spline** button to have Processor place an initial set of breakpoints on the spectrum for you.
11. Click in the spectrum viewing pane to add more breakpoints for a cubic spline baseline fit. Use more points where the baseline is strongly curved, and fewer points where it is more linear.

Breakpoints normally appear on the spectrum line. Hold down the **Control** or **Command** key while clicking to place breakpoints at arbitrary positions.

12. Continue placing breakpoints until the blue line follows the spectrum baseline as closely as possible.
13. When you are satisfied with your fit of the spectrum baseline, click the **Accept** button to apply your changes to the spectrum.

Adjusting the CSI Settings

14. Switch to the CSI Editor.
15. Zoom in on the single DSS peak at the far right side of the spectrum. Use the **Automatic Vertical Zoom** button to ensure the best view of the single peak.
16. Click the **Find Automatically** button to have Processor attempt to match the peak automatically.

If your spectrum has poor phasing or baseline correction, finding the CSI automatically will not work; readjust the phasing and baseline correction before proceeding.

When you process spectra using format as a CSI, the automatic fit will not work; you will need to use the manual fit instead.

17. If necessary, use the Adjustment Handle, width controls or keyboard controls to adjust the shape of the red overlaid peak to match the black spectrum line as closely as possible.

Many of the CSI parameters interact to affect the quality of the analysis that is possible with your spectra. Perhaps most obviously, the CSI peak position calibrates the chemical shift scale. In addition, the CSI concentration and peak height both

affect quantification accuracy, while an accurate CSI peak width allows you to better fit regions of the spectrum in which multiple clusters overlap.

18. While you are zoomed in to the CSI peak, apply reference deconvolution.

Reference deconvolution corrects your spectrum for some types of systematic errors that may have occurred during acquisition. It may not be appropriate for all spectra.

Applying Reference Deconvolution

19. Switch to Reference Deconvolution mode.
20. Leave the Target Linewidth at its initial value of around 0.6 Hz.

If you also want to apply line broadening to the spectrum, you may add the desired amount of line broadening to the target linewidth here. For example, to add 0.5 Hz of line broadening, you would change the target linewidth here to 1.1 Hz.

21. Select the **Use DSS/TSP Satellites** check box.

This option includes the ^{29}Si satellite peaks on either side of the DSS methyl peak in the reference deconvolution model.

22. Click the **Accept** button to apply your changes to the spectrum.
23. Switch to the CSI Editor.

The CSI definition needs adjustment following reference deconvolution, as peak shapes have changed throughout the spectrum. It is good practice to double check the fit of the CSI peak any time you change the processing of a spectrum.

24. Click the **Find Automatically** button.

Applying Line Broadening

25. Switch to Line Broadening mode.
26. Enter a value of 0.4 Hz.

When you analyze spectra using the Chenomx library, the ideal CSI linewidth is somewhere around 1.0 to 1.5 Hz, to match the spectra originally used to create the compound signatures. Since the spectrum started with a linewidth of 0.6 Hz, adding 0.4 Hz here brings the spectrum into the ideal range.

You could also have added this value (0.4 Hz) to the target linewidth during reference deconvolution, instead of applying it in a separate step.

27. Click the **Accept** button to apply your changes to the spectrum.
28. Switch to the CSI Editor.

Once again, since we have modified the processing of the spectrum, it is a good practice to verify the fit of the CSI peak.

29. Click the **Find Automatically** button.

Determining the Spectrum pH Automatically

30. Switch to the pH Editor.

Setting the spectrum pH is very important! Most of the Chenomx Reference Compounds are pH sensitive, meaning their starting positions and transform windows are affected by the spectrum pH.

If you set a spectrum pH that differs a lot from the actual sample pH, clusters will start farther than necessary from their real positions, and the transform windows may not allow you to move each cluster far enough to compensate without overriding them.

31. Click the **Find Automatically** button to have Processor attempt to determine the spectrum pH automatically.

Now that you have corrected the phasing and baseline of the spectrum, automatic pH determination will work.

The calculated pH is 8.37.

Determining the Spectrum pH Manually

32. Click and drag the pH slider bar to adjust the pH, or manually enter a pH value in the pH text entry box.

The blue, green and yellow lines, indicating the positions of the imidazole, creatinine and DFTMP peaks respectively, move as you adjust the pH. You can determine a coarse pH simply by adjusting these lines as close as possible to the corresponding peaks in your spectrum.

33. Using the coarse pH, and the information box to the left of the pH slider bar, note which pH indicator is most active in the pH range of the spectrum. In this case, it is imidazole.

If the coarse pH is near the transition between creatinine and imidazole (around 5.5), then both indicators are useful in refining the pH value.

34. Zoom in on the two peaks for imidazole. Click the **Automatic Vertical Zoom** button to ensure the best view of the two peaks.

The range of pH available on the slider bar changes as you zoom in to yield a finer-grained control over your pH adjustments. To make the entire range of 0-14 available, clear the **Fine Tune** check box.

35. Click and drag the pH slider bar to adjust the pH, or manually enter a pH value in the pH text entry box.

As before, the indicator lines move as you adjust the pH. To obtain a refined pH value, ensure that each line is as close as possible to its associated peak, keeping in mind that every change in pH can move both peaks.

Once again, the measured pH should be around 8.37.

36. Click the **Show Entire Spectrum** button to restore a view of the whole spectrum.

Once you have finished adjusting your spectrum in Processor, you can:

- Move the spectrum directly to Profiler for further analysis (using the **Jump to Profiler** button)
- Save the spectrum for later analysis in Profiler, for creating a simulation in Spin Simulator, or for creating a compound signature in Compound Builder.

10.3 Profiler

Introduction

The Profiler module is the primary tool in Chenomx NMR Suite for profiling sample solutions. You can profile samples by matching peak frequencies and heights contained in a compound library with frequencies and heights in the sample spectrum. To help you become familiar with the profiling process, the first part of the following tutorial provides a step-by-step guide to fitting a sample spectrum containing only twelve compounds, none of which have overlapping signatures. The second part provides some techniques for fitting spectra in which some compound signatures overlap.

Basic Profiling Techniques

The simplest way to start analyzing a spectrum is to profile easily identifiable compounds first, as they often have additional peaks in the busier regions of the spectrum. Profiling them first can help rule out other compounds later in the profiling process. In this simple spectrum, all of the compounds are quite well separated from each other, and should be relatively simple to profile.

Single Peak Compound - Acetate

1. Open Profiler.
2. Open the sample file called `Samples\Profiler\Basic Start.cnx`.
3. Zoom in on the single peak near 1.90 ppm.
4. Right-click on the peak, and click **Filter for Compounds Near 1.91 ppm**.
5. In the Compound Table, select the compound with the highest Maximum Concentration (acetate), as this is the most likely candidate for fitting the peak.

You may want to sort the Compound Table by Maximum Concentration to make selecting the highest Maximum Concentration simpler. Click the header of the Maximum Concentration column to sort by that column. Additional clicks on the header reverse the sort order.

6. Click **Compound > Selected Compound > Fit Automatically**, or press the space bar.

Profiler's automatic fit suggests concentrations and transforms for the currently selected compound. Although it will often arrive at acceptable values, always review the results of an automatic fit before moving on to other compounds.

7. Adjust the acetate signature (blue) as needed to match it to the corresponding peak in the spectrum (black).

You can click and drag the controls on the horizontal and vertical axes to move the peak along only one axis at a time, or use the arrow keys to move the peak up, down, left and right.

8. Clear the quick search box above the Compound Table by clicking the red "X" or pressing **Esc**.
9. Click **Show Entire Spectrum**.

Single Peak Compound - Pyruvate

10. Zoom in on the single peak near 2.40 ppm.

11. Right-click on the peak, and click **Filter for Compounds Near 2.36 ppm**.
12. In the Compound Table, select pyruvate, since it has the highest Maximum Concentration.
13. Use **Fit Automatically** to suggest a concentration and transforms for pyruvate (press the space bar).
14. Adjust the pyruvate signature (blue) as needed to match it to the corresponding peak in the spectrum (black).
15. Clear the quick search box.
16. Click **Show Entire Spectrum**.

Water Suppression Effects - Profiling Alanine

17. Zoom in on the doublet near 1.50 ppm.
18. Right-click on the doublet, and click **Filter for Compounds Near 1.47 ppm**.
19. In the Compound Table, select alanine.
20. Use **Fit Automatically** to suggest a concentration and transforms for alanine (press the space bar).
21. Adjust the height of the alanine doublet (blue) as needed to match the corresponding doublet in the spectrum (black).

You should not normally increase the height of a signature above the black line, so use the right-most peak to determine the appropriate height.

22. Move to the second cluster of alanine (3.8 ppm) using the Cluster Navigator, and adjust its frequency as needed to match the spectrum.

Since the second cluster is much closer to water than the first one, it is more sensitive to side effects of the water suppression pulse used to acquire the spectrum. Thus, the cluster at 1.5 ppm is a good choice for determining the concentration of alanine, while the cluster at 3.8 ppm must be fit using only frequency adjustments.

For more details on handling water suppression effects, see “Water Suppression” on page 120.

23. Clear the quick search box.
24. Click **Show Entire Spectrum**.

Profiling the Remaining Compounds

25. Using the techniques discussed above, profile the remaining nine compounds in the spectrum. You can check your results for all of the compounds with the following list of concentrations. If you are unsure what the finished profile should look like, open the file `Basic_End.cnx` to view this spectrum with profiling completed.

Note that in one instance (3-indoxylsulfate), Profiler is unable to suggest a concentration or transforms via **Fit Automatically**, due to the atypical shape of the cluster at 7.3 ppm. For this compound, choose a concentration that provides the best match to the other clusters, then adjust frequencies of each cluster to best match the spectrum.

Some clusters are more reliable than others when determining concentrations of compounds that have multiple clusters. When you determine concentrations manually, use clusters

meeting as many of the following criteria as possible to set the concentration for such compounds:

- Minimal overlap. Ideally, the cluster should have no other clusters overlapping.
- More protons. If a cluster corresponds to three protons (such as a methyl group), the signal for that cluster is stronger than the signal for a cluster corresponding to only one proton.
- Unobscured. If a cluster appears near water (around 4.73 ppm), its signal may be obscured by the presence of the water peak. Baseline correction or solvent suppression during processing may further affect quantification. Use obscured clusters only when no alternative is available, and recognize that concentrations based solely on obscured clusters may be unreliable.

Table 10.1. Expected Concentrations for Basic Start.cnx

Compound	Conc (μM)	Compound	Conc (μM)
3-Indoxylsulfate	458.2	Formate	321.2
4-Hydroxy-3-methoxymandelate	49.1	Fumarate	155.6
Acetate	613.8	Niacinamide	564.4
Alanine	608.3	Pantothenate	94.7
Citrate	432.1	Pyruvate	598.9
DSS (Chemical Shift Indicator)	481.9	Urea	147094.0

Advanced Profiling Techniques

Real samples are rarely as simple as the previous example, and usually contain compounds whose signatures overlap. Although profiling becomes more challenging in these cases, the strategy remains the same; profile easily identifiable compounds first, and use the contributions of those compounds to refine your efforts at profiling the overlapped regions.

There is a certain artistry to efficiently profiling a complex spectrum. The best method for learning many of the techniques that you can use to profile such a spectrum is practice, profiling progressively more complex spectra until you develop a level of proficiency. More detailed discussions of profiling strategies are available elsewhere in this manual (see “Targeted Profiling” on page 116, “Profiling a Spectrum” on page 8 and “Profiling Overlapped Regions in a Spectrum” on page 11).

Getting Started

1. Open **Profiler**.
2. Open the file called `Samples\Profiler\Advanced Start.cnx`.

Although completely non-overlapping compounds are easy to profile, they cannot give any further insight into the shape of other compounds in the spectrum. Profile these compounds first, so that you can focus on the more challenging regions of the spectrum. Since we covered the relevant techniques in “Basic Profiling Techniques” on page 110, they will not be covered in detail here.

Profiling Compounds using Basic Profiling Techniques

3. Using the basic techniques discussed in “Basic Profiling Techniques” on page 110, profile the following twelve compounds.
 - 3-Indoxylsulfate

- 4-Hydroxy-3-methoxymandelate
- Acetate
- Citrate
- Dimethylamine
- DSS (Chemical Shape Indicator)
- Formate
- Fumarate
- Niacinamide
- Pyruvate
- Threonine
- Trimethylamine

Once you have profiled the non-overlapped compounds, a more systematic approach will help in profiling the remainder. For this tutorial, start at the right-most edge of the spectrum, and move to the left.

Fitting Compounds with Overlap

4. Zoom in on the two single peaks near 0.90 ppm.

These peaks correspond to methyl groups in the structure of pantothenate (the highest Maximum Concentration in this region of the spectrum). Methyl groups can be very useful in setting an appropriate concentration for a compound, as they are usually easy to locate.

5. Use **Fit Automatically** (press the space bar), then adjust the concentration of pantothenate to match the spectrum in this region.
6. Click the location of the next pantothenate cluster (2.4 ppm) in the Cluster Navigator.
7. Adjust the frequency of the cluster to match the signature (blue) to the spectrum (black). Although the left-most peak of this cluster overlaps with 2-oxoglutarate, the remaining two peaks provide enough information to establish an appropriate frequency for the cluster.

It is good practice to completely finish fitting an overlap region before moving on to the next one. Since the only other compound in this region is 2-oxoglutarate, it is the next logical compound to consider.

8. In the Compound Table, select 2-oxoglutarate, and use **Fit Automatically**.
9. Adjust the 2-oxoglutarate cluster near 2.43 ppm as needed to match the spectrum.

Notice that the sum line (red) at 2.417 ppm now matches the spectrum more closely than it did with only pantothenate adjusted.

Since 2-oxoglutarate is quite simple, you can finish profiling it before continuing with pantothenate.

10. Move to the remaining 2-oxoglutarate cluster (3.0 ppm) using the Cluster Navigator, and adjust the frequency of the cluster as needed to match the spectrum.
11. Select pantothenate again.

You can simply click on the sum line (red) to select the nearest compound. Click repeatedly to cycle through the contributing compounds in an overlap region.

12. Move to the next pantothenate cluster (3.4 ppm) using the Cluster Navigator, and adjust the frequency of the cluster at 3.39 ppm to fit the spectrum.
13. In the Compound Table, and use **Fit Automatically** on fucose and taurine, and then select taurine.
14. Move to the taurine cluster near 3.3 ppm, and adjust it as needed to match the spectrum.
15. Return to the taurine cluster near 3.4 ppm, and adjust its frequency as needed to match the spectrum.

In overlapped regions like this, the sum line (red) becomes a better indicator of the quality of your profile in the region, so check the match of the sum line to the spectrum before proceeding.

16. Select fucose, and use the Cluster Navigator to display the clusters near 1.2 ppm.
Click and drag across the series of clusters in the Cluster Navigator to view the entire group.
17. Adjust the concentration of fucose so that the clusters near 1.2 ppm match the spectrum. The cluster at 1.24 ppm provides the best reference to establish the concentration for fucose.
18. Return to the fucose cluster near 3.4 ppm, and adjust its frequency as needed to match the spectrum, using the sum line to evaluate your progress.
19. Select pantothenate again, and use the Cluster Navigator to view the cluster near 3.5 ppm.
20. Adjust the frequency of the pantothenate cluster at 3.51 ppm as needed to match the spectrum.
21. Move to the pantothenate cluster near 4.0 ppm, and adjust its frequency as needed to match the cluster to the nearest peak in the spectrum.

The profiled peak for pantothenate in this region does not completely match the area available, due to overlap with serine. Once you have profiled serine, this peak will match better.

22. In the Compound Table, use **Fit Automatically** on glycolate and serine, and then select serine.
23. Use the Cluster Navigator to display the serine cluster near 3.8 ppm, and adjust its concentration to better match the spectrum, using the sum line to evaluate your progress.
24. Use the Cluster Navigator to display the serine clusters at 4.0 and 3.9 ppm.
Remember that you can click and drag across several clusters to display all of them.
25. Adjust the frequencies of the remaining two serine clusters.

Notice that the sum line over the pantothenate cluster in this region now yields a much better match to the spectrum.

26. Select glycolate, and adjust the single glycolate peak as needed to match the spectrum.

27. Select fucose, display the two clusters nearest 3.7 ppm, and adjust their frequencies as needed to match the spectrum.
28. In the Compound Table, select alanine and use **Fit Automatically**.
29. Display the alanine cluster near 1.5 ppm, and adjust it as needed to match the spectrum.
Use the right-most peak to set the height for the cluster, much like you did for the basic spectrum.
30. Return to the alanine cluster near 3.8 ppm, and adjust its frequency as needed to match the spectrum, using the sum line to evaluate your progress.
31. Adjust the frequencies of each of the fucose clusters near 3.8 ppm. Try finding individual peaks in each cluster that match peaks in the spectrum, and lining them up.
32. Adjust the frequencies of each of the remaining fucose clusters as needed. If there is not enough of a pattern visible in the spectrum to which you can match a cluster, leave the cluster in its original position.
33. In the Compound Table, select urea and use **Fit Automatically**.
34. Adjust the urea cluster near 5.8 ppm as needed to match the spectrum.
35. Check the overall shape of the spectrum to verify that the sum line lies generally along the spectrum line.

You can check your results for all of the compounds with the following list of concentrations. If you are unsure what the finished profile should look like, open the file Advanced End.cnx to view this spectrum with profiling completed.

Table 10.2. Expected Concentrations for Advanced Start.cnx

Compound	Conc (μM)	Compound	Conc (μM)
2-Oxoglutarate	620.9	Fumarate	155.5
3-Indoxylsulfate	448.7	Glycolate	852.1
4-Hydroxy-3-methoxymandelate	50.2	Niacinamide	581.4
Acetate	616.3	Pantothenate	98.7
Alanine	614.2	Pyruvate	587.2
Citrate	452.7	Serine	1185.2
DSS (Chemical Shape Indicator)	481.9	Taurine	1196.3
Dimethylamine	655.7	Threonine	424.1
Formate	328.3	Trimethylamine	618.8
Fucose	618.2	Urea	332288.8

Bonus compounds: There are at least two additional compounds present in trace quantities in this spectrum. Identify and profile these compounds and check their concentrations here:

Suberate (60.6 μM)

Isopropanol (14.5 μM)

10.4 Targeted Profiling

Targeted profiling involves matching a series of compounds to an experimental spectrum. Practically, targeted profiling includes the two tasks of identifying a compound or compounds that you would like to profile, and matching the individual clusters of each compound to their respective regions of the spectrum. Identifying a compound to profile is simply a matter of displaying the compound and visually comparing its signature with visible patterns in the spectrum, remembering that you can adjust both height (concentration) and position (frequency) during the profiling process. There are several methods of filtering and sorting the compound list to include compounds that are most likely to appear in a given region of the spectrum.

Profiling individual clusters in a region of the spectrum requires some care and attention; the most important facet of any profiling technique is applying it consistently to all spectra in a dataset. Profiling can involve a variety of methods, depending on the intensity of the peaks, the presence of other peaks in the region, the degree of overlap with other peaks, and the presence of other clusters from the same compound elsewhere in the spectrum. Some of these methods are described below; you may discover other methods as you analyze more spectra.

These techniques refer specifically to manual profiling, as well as refining the results of automatic fits using **Fit Automatically**.

Identifying Compounds

There are a variety of methods that you can use to help you identify compounds that are present in any given spectrum. Some are implemented directly in the Profiler module, while others may involve some literature review or external reference.

Each of the following techniques is essentially a method of limiting the list of potential compounds for a given region of the spectrum. You must make the final decision as to whether or not a particular compound is represented in the spectrum by visually comparing the signature for the compound to the patterns visible in the spectrum. Remember that each cluster can be moved somewhat to either side of its starting position (inside its transform window), and the overall height of the compound can be adjusted.

Frequency Filters

You can use Filter Compounds to filter the Compound Table to include only compounds that can appear at a certain frequency, or in a certain range of frequencies (see "Quick Searches" on page 50). To quickly filter based on a single location, right-click on the spectrum, and click **Filter for Compounds Near x.xx ppm** to select a single frequency. Similarly, to filter based on a frequency range, select a region, then right-click inside the selected region and click **Filter for Compounds in Selected Region**.

The frequency filters display all compounds that have transform windows overlapping the selected frequency or frequency range. Since transform windows are a function of spectrum pH, the number of compounds displayed using the frequency filters varies based on your spectrum pH.

Text Filters

You can filter the Compound Table based on a simple text search, by typing the text to search for in the empty text box at the top right corner of the Compound Table. Filtering can help you look for a particular compound, or a group of compounds with similar names, like all compounds containing "acetate".

You can also use text filters to type in frequency searches (see “Frequency Filters” on page 116). To search for compounds near 1.90 ppm, type “@1.90” in the text box (without the quotes). To search for compounds from 2.75 to 4.30 ppm, type “@2.75 to 4.30” (without the quotes).

Compound Lists

A list of compounds known to occur in samples of the type that you are analyzing can be an invaluable tool to guide your profiling efforts. In many cases, the number of compounds that could be contributing to a specific region of the spectrum can be large, so it can be helpful to know some compounds that are more likely to be present. If you do not have such a list, and are expecting to analyze samples of a particular type on a regular basis, you may find it helpful to develop a list of your own based on representative samples.

For the Profiler sample files, the following compound list guided the profiling process:

2-Oxoglutarate	DSS (Chemical Shape Indicator)	Pyruvate
3-Indoxylsulfate	Formate	Serine
4-Hydroxy-3-methoxymandelate	Fucose	Taurine
Acetate	Fumarate	Threonine
Alanine	Glycolate	Trimethylamine
Citrate	Niacinamide	Urea
Dimethylamine	Pantothenate	

If you have a list of compounds common to a sample type, you may want to create a Compound Set in Library Manager and use it to perform initial profiles of the spectra in your dataset.

Test Mixtures

As you profile a series of spectra in a larger dataset, you will notice that the profiling some mixtures of compounds is difficult. You may want to try running controlled test mixtures of the compounds in question to establish their behaviour under your experimental conditions, and use information from the test mixtures to guide your profiles of the larger dataset. Test mixtures will help increase confidence and consistency in your profiling results.

Profiling Clusters

Once you have selected a compound to profile, you need to match the individual clusters of the compound in their respective regions of the spectrum. There are several common approaches to profiling clusters, and in a typical spectrum, more than one method will be necessary. Try to maintain a consistent approach to profiling spectra in a dataset, as this will help you obtain accurate and consistent results from your analysis.

Compound Line

The simplest way to match a cluster to a region of the spectrum is to adjust the cluster so that the compound line (blue by default, unless you have assigned custom colors) matches the spectrum line (black). The compound line is most useful for isolated, low-intensity clusters. If you profile all of the clusters for a compound using only the compound line, the measured concentration of the compound may be unreliable, particularly if other clusters are nearby; the identity can still be considered accurate if the cluster shapes are clearly

discernible. Using the compound line on tall, isolated clusters generally provides accurate identity and quantity.




Figure 10.1. Apart from these clusters in the aromatic region, none of the other peaks for tyrosine are clearly visible. Although it is clear that tyrosine is present, the quantity of 19 μM should be considered approximate.

Subtraction Line

For complex regions shapes, use the subtraction line to determine an optimal fit for the cluster, where possible. Adjust the cluster height and position to remove the influence of the adjusted cluster while approximating a normal spectrum with the subtraction line. Using the subtraction line is an advanced technique offering some of the most accurate shapes available, but requires more time than other, simpler techniques.




Figure 10.2. The overfit subtraction line (top center) shows a negative residual peak, and results in a concentration of 25 μM . The underfit subtraction line (bottom center) leaves a positive residual peak, and results in a concentration of 12 μM . The accurate profile (far right) results in a subtraction line that approximates a smooth baseline through the area, indicating that the signal has been effectively subtracted from the spectrum.




Figure 10.3. Inaccurate positions result in a combination of over- and underfits: A negative residual peak on the right suggests the cluster should be moved left, or downfield (top center). A negative residual peak on the left suggests the cluster should be moved right, or upfield (bottom center). Both identity and quantity can be considered reliable in the accurate profile (far right).

Sum Line

You will encounter situations in which the cluster you need to profile overlaps one or more clusters from other compounds; regions fitting this description are often called "busy regions". You often need to profile compounds that have overlapped clusters in groups, ensuring that the combined shape (sum line) matches the spectrum line.

First, identify likely compounds contributing to the overlapped region. Where available, use non-overlapped peaks, including those in other regions of the spectrum, to establish concentration ranges for the identified compounds. Use partially-exposed peaks to establish positions for each of the overlapped clusters. Starting with the approximate shape thus obtained, fine-tune the position and height of each of the overlapped clusters to optimize the overall shape in the region.

Generally, profile multiplets before singlets in overlap regions, as they are easier to position accurately. If all of the clusters for a compound are in overlap regions, the profile may be unreliable, especially if most or all of the clusters are low-intensity singlets.




Figure 10.4. Using frequency filters and sorting by maximum concentration lets you identify serine, glycolate and pantothenate as candidate compounds. Visual comparison of the signatures to the spectrum line confirms the identifications.




Figure 10.5. Profiling each compound individually using the basic method (simple match to spectrum line) yields reasonable concentrations, but the sum line doesn't quite match.




Figure 10.6. Considering all the compounds together, along with confirmation using other clusters of serine and pantothenate, gives a better match. Also, note the more than 20% change in pantothenate concentration!

Water Suppression

Spectra of aqueous samples are often acquired using a water suppressing pulse sequence, since the water signal would otherwise drown out many useful signals. The water suppression pulse can also partially suppress other resonances occurring near the frequency of the water signal. When you profile compounds with multiple clusters, use clusters well away from the water signal to determine compound quantities. This may result in a slight "overfit" (sum line higher than spectrum line) of the clusters closer to water, but gives you more consistent and accurate concentrations.




Figure 10.7. Since cluster 0 (right) is quite far from the water peak, it experiences little suppression. Profiling lactate based on cluster 0 results in a slight but acceptable overfit of cluster 1 (left).




Figure 10.8. Profiling lactate based on cluster 1 (left) incorrectly reduces the reported concentration by almost 11%.

Best Guess

When you profile a compound with multiple clusters, you may encounter a cluster that is difficult to profile accurately, due to a lack of prominent peaks in the signature or in the spectrum, or both. In this case, use minor adjustments to approximately match the sum line to the spectrum line; the goal here is mostly aesthetic, since very little information is available to produce an accurate profile. If all of the clusters for a compound must be profiled using this technique, the profiled compound (both identity and quantity) should be considered unreliable.




Figure 10.9. The large cluster near 3.2 ppm was sufficient to establish the identity and quantity of O-phosphocholine. Although there are no obvious peaks that can be used to align clusters 3 and 4 (top), it seems reasonable to move them over to fit more neatly under the nearby "hump" in the spectrum (bottom).

No Fit

When profiling some compounds, you may find that a cluster is not visible in the signature or the spectrum, or both, at concentrations consistent with the spectrum. Since no information is available to help you profile the cluster, the best course is to leave the cluster at its starting position. If all of the clusters for a compound fall into this category, the compound should be considered not detected in the sample.




Figure 10.10. The doublet near 1 ppm provides enough information to set a concentration for isoleucine, and the triplet near 0.92 ppm lets you confirm the identification. It is reasonable that isoleucine is present, at or near 9 μ M.




Figure 10.11. There is no clear information elsewhere in the spectrum to guide the placement of the remaining clusters, so they are left at their default positions. The initial assessment of isoleucine at 9 μ M is still consistent with the spectrum.

10.5 Spin Simulator

Introduction


You can use Spin Simulator to create spin simulations of compounds that you would like to add to your library. Simulations that you create in Spin Simulator can be imported into Compound Builder to create compound signatures that you can then add to your library.

You may find it is helpful to have the molecular structure of the compound you are simulating available for reference purposes.

You can only complete the Spin Simulator tutorial if you have purchased the Professional edition of Chenomx NMR Suite, or if you are using the Evaluation edition of Chenomx NMR Suite.

Simple Simulation - Valine

To become familiar with the various tools available in Spin Simulator, as well as with the basic process involved in creating a simulation, you will simulate a simple compound, valine.



Aside from exchangeable protons, like the carboxylic acid and amine protons, there are four types of protons expected to have visible signals for this compound, marked as 1 through 4 in the structure above. The two methyl groups (1 and 2) and one of the methine groups (3) will each couple to the remaining methine group (4). Thus, in the spectrum for valine, you would expect to see three doublets, each coupled to one complex cluster (rigorously, a quartet of quartets of doublets, or qqd).

Simulating Valine

1. Assign the spectrum, matching protons in the structure of the compound with signals in the spectrum. For valine, the assignments are as follows:
 - a. Methyl signals (1 and 2) appear at 0.978 and 1.031, each coupled to one of the the methine signals (4).
 - b. Methine signal (3) appears at 3.604 ppm, and is coupled to the other methine (4).
 - c. Methine signal (4) appears at 2.264 ppm, and is coupled to each of the other three signals.
2. Open **Spin Simulator**.
3. Overlay the valine spectrum, called `Samples\Spin Simulator\Valine.cnx`.
4. Add spin definitions for each of the signals from the assignment above, starting from the right-hand side of the spectrum (low numbers) and moving left (to higher numbers).

You now have SD0 and SD1 for the methyl groups (3 protons each, at 0.978 and 1.031 ppm), and SD2 and SD3 for the methine groups (1 proton each, at 2.264 and 3.604 ppm), with a tab for each spin definition in the Simulation view.
5. Zoom in to the two methyl signals near 1.0 ppm.

You can click and drag across the two corresponding entries in the Spin Navigator to provide a useful zoom level.
6. Measure the distance between the pairs of peaks in each doublet using the tape measure.

The right-hand doublet measures about 7.02 Hz, while the left-hand one measures about 7.07 Hz.
7. Select SD0 (0.978 ppm) by clicking on it in the graph, or clicking on its tab in the Simulation view.
8. Right-click on SD0 (blue) and click **Add J-Modifier...**
9. Enter a J-value of 7.02 Hz, select Uncoupled Split with a split count of 1, and click OK.

10. Adjust first spin definition (blue) to fit the spectrum (black). You can adjust freely by clicking and dragging the peak shape, or use either the controls on the horizontal and vertical axis or the up, down, left and right arrow keys to restrict your adjustment to specific dimensions and directions.
11. Adjust the width of the doublet using the controls to either side of the spin definition, so that it better matches the spectrum.

A width of about 1.2 to 1.3 Hz is sufficient.
12. Select SD1 (1.031 ppm).
13. Add a J-modifier to SD1.

Use a J-value of 7.07 Hz, and select Uncoupled Split with a split count of 1.
14. Adjust SD1 to fit the spectrum.

Once again, a width of 1.2 to 1.3 Hz is sufficient.
15. Zoom to SD3 (3.604 ppm) using the Spin Navigator or the Spin Definition pane.
16. Measure the distance between peaks in the doublet.

This doublet measures about 4.36 Hz.
17. Add a J-modifier to SD3.

Use a J-value of 4.36 Hz and select Uncoupled Split with a split count of 1.
18. Adjust SD3 to fit the spectrum.

This doublet has a width of around 1.1 Hz.
19. Zoom to SD2 (2.264 ppm) using the Spin Navigator or the Spin Definition pane.

The resulting view does not reveal much of the actual spectrum, so adjust it to get a better view.
20. Click **View > Zoom > Set Zoom**, and enter 2.23 to 2.30 ppm for the X Axis, and -0.05 to 0.87 si for the Y Axis.
21. Add a J-modifier to SD2. Select Coupled, to Spin Definition 0, using the uncoupled split of 7.02 Hz.

The spin definition is now linked to one of the methyl groups; notice that it has split to a quartet (n=3).
22. Add two more J-modifiers, one coupled to SD1 using the uncoupled split of 7.07 Hz, and the other coupled to Spin Definition 3 using the uncoupled split of 4.36 Hz.

Notice that with each added J-modifier, the spin definition becomes more complex, and matches more closely with the spectrum.
23. Adjust SD2 to better fit the spectrum.

This spin definition has a width of about 1.5 Hz.
24. Notice that changing the height of either spin definition changes the height of the whole simulation; you cannot obtain a perfect height match with all four spin definitions at the same time. The discrepancy is normal, and can be corrected in Compound Builder prior to adding the compound to your library (see "Simple Compound - Valine" on page 129 for details).


For now, simply choose one of the spin definitions with which to fit the simulation height.

25. Save the simulation file.

If you would like to review the sample simulation of valine, do not replace the file called **Samples\Spin Simulator\Valine.xss**.

Complex Simulation - Proline

Now that you have developed a basic familiarity with some of the tools available in Spin Simulator, you are ready to use them to simulate a more complex compound, proline.



The spectrum of proline is much more complex than the spectrum of valine. The protons at each position couple to all of the protons at adjacent positions. Also, the protons at each position are diastereotopic, meaning that they can couple with each other (geminal coupling), as well as each coupling to adjacent protons. While this means that splitting patterns will be based exclusively on doublets ($n=1$), it also means that you need to be aware of a lot more coupling interactions.

Simulating Proline

1. Assign the spectrum, matching protons in the structure of the compound with signals in the spectrum. Note that this is often one of the largest tasks in simulating a new compound. For proline, the assignments are as follows (protons labelled 'a' project above the plane of the page, while those labelled 'b' project below the plane):
 - a. The signal for proton 4b appears at 1.979 ppm, and couples to protons 3a, 3b, 4a, 5a and 5b (Spin Definition 0, or SD0).
 - b. The signal for proton 4a appears at 2.012 ppm, and couples to protons 3a, 3b, 4b, 5a and 5b (Spin Definition 1, or SD1).
 - c. The signal for proton 3b appears at 2.058 ppm, and couples to protons 2, 3a, 4a and 4b (Spin Definition 2, or SD2).
 - d. The signal for proton 3a appears at 2.339 ppm, and couples to protons 2, 3b, 4a and 4b (Spin Definition 3, or SD3).
 - e. The signal for proton 5a appears at 3.328 ppm, and couples to protons 4a, 4b and 5b (Spin Definition 4, or SD4).
 - f. The signal for proton 5b appears at 3.413 ppm, and couples to protons 4a, 4b and 5a (Spin Definition 5, or SD5).
 - g. The signal for proton 2 appears at 4.121 ppm, and couples to protons 3a and 3b (Spin Definition 6, or SD6).
2. Open **Spin Simulator**.
3. Overlay the proline spectrum, called **Samples\Spin Simulator\Proline.cnx**.

4. Add spin definitions for each of the protons listed in the assignment above, making sure that the number of protons for each is set to one. Although you do not need to add the spin definitions in the indicated order to complete the tutorial, doing so will make the tutorial easier to understand, since the Spin Definition (SD) notation appears throughout this tutorial for clarity.
5. Zoom in to SD6 (4.121 ppm).
6. Measure the two coupling constants represented in this region of the spectrum, using the spectrum line (black) as your reference. Specifically, find the distances from an outermost peak to the adjacent peak and from an outermost peak to the second peak in. The distances in this spectrum are 8.79 Hz and 6.52 Hz.
7. Add one J-modifier to SD6 for each of the measured distances. Use uncoupled splits, with split counts of 1 and J-values of 8.79 and 6.52 Hz.

You know that each of these couplings will link to single protons, so you can use a split count of 1 for each of them.

8. Adjust the frequency, height and width of SD6 to better fit the spectrum.
9. If you need to adjust the J-values for SD6, select the J-modifier in the Simulation view, and use the J-Value Editor appearing under the Spin Navigator to fine-tune the J-value.
10. Zoom in to SD4 and SD5.

Based on the assignment above, there are three coupling constants associated with each of these spin definitions. One can be easily measured as the distance between an outermost peak and the second peak in; this measurement yields around 11.58 Hz for both spin definitions, suggesting that this is the coupling constant linking these two spin definitions. In addition, the value itself (11.58 Hz) is consistent with the geminal coupling relationship between the protons responsible for these signals.

11. Add a J-modifier to SD5, and couple it to SD4 using a New Split, with a J-value of 11.58 Hz.

The remaining two coupling constants for SD5 are both near 7 Hz (the distance from an outermost peak to the adjacent one), but it is not clear what the precise values should be. You can set them to 7 Hz to approximate the correct value, and use the J-Value Editors to refine the values visually.

12. Add two J-modifiers to SD5 as uncoupled splits, each with a J-value of 7 Hz, and a split count of 1.
13. Use the J-Value Editor to fine-tune the J-modifiers that you just added. Aim for values that center each of the peaks in the spin definition (blue) on the corresponding peaks in the spectrum (black). The resulting values are 7.40 Hz and 6.67 Hz.
14. Adjust the frequency, height and width of SD5 to better fit the spectrum.

You can determine the remaining two coupling constants for SD4 using similar techniques. The approximate distances are 7.2 Hz.

15. Add two J-modifiers to SD4 as uncoupled splits, with J-values of 7.2 Hz each, and split counts of 1.
16. Use the J-Value Editor to fine-tune the J-modifiers. The resulting values are 7.44 and 6.88 Hz.

17. Adjust the frequency, height and width of SD4 to better fit the spectrum.

The remaining spin definitions are quite complicated. There is a lot of coupling among them, including some to spin definitions that you have already fit. Some of these links to solved spin definitions can help reduce the complexity of modeling the spin definitions that remain. Start with SD3, since it is the last one that is well-isolated from other signals.

Due to the complexity of this signal, it is not immediately clear which measurements you can use to extract coupling constants. However, one clue exists in the geminal relationship between the protons responsible for this signal and SD2. You observed a coupling constant of 11.58 Hz previously between geminal protons, so start by looking for a measurement of similar magnitude. Measuring from an outermost peak to a peak near 11.58 Hz gives about 8.7 or 13.7 Hz. The former value does not appear quite large enough for a geminal coupling, so use the latter to couple SD3 to SD2.

18. Add a J-modifier to SD3, and couple it to SD2 using a New Split with a J-value of 13.7 Hz.

SD3 should also be coupled to SD6, based on the original assignments. There are two J-modifiers available to create this link, with J-values of 8.79 and 6.52 Hz. However, based on the structure, you can estimate the dihedral angle between these two protons as larger than that between SD2 and SD6. A larger dihedral angle implies a larger coupling constant, according to the vicinal Karplus correlation, so use the larger of the two here.

19. Add a J-modifier to SD3, and couple it to SD6 using the larger of the two available uncoupled splits (8.79 Hz).

Finally, you know that this signal should be coupled to SD0 and SD1, but the magnitude of the coupling constants is not clear. You do have information from some of the other couplings in the molecule that suggest the coupling constants should be near 7 Hz, so use this as a starting point, and fine-tune using the J-Value Editor. Similarly, based on the dihedral angles among the protons, you can speculate that the coupling constant to SD1 will be larger than the coupling to SD0, and bias your adjustments accordingly.

20. Add two J-modifiers to SD3, one coupled to SD0, and the other to SD1, each with a J-value of 7 Hz.
21. Use the J-Value Editor to fine-tune the two new J-modifiers. Ensure that the J-value for the coupling to SD1 is larger than the J-value for the coupling to SD0.

Unfortunately, there is no entirely satisfactory combination of these two coupling constants that results in a good fit of the cluster. However, the coupling constant to SD2 is still speculative, so try adjusting that, as well.

22. Continue adjusting the J-modifiers, this time including the coupling to SD2 in your adjustments. As before, aim for values that center each of the peaks in the spin definition (blue) on the corresponding peaks in the spectrum (black). The resulting values are 12.99 Hz (to SD2), 7.38 Hz (to SD1) and 6.65 Hz (to SD0).
23. Adjust the frequency, height and width of SD3 to better fit the spectrum.

SD2 has a coupling pattern similar to SD3. You have already coupled it to SD3, and there is only one remaining J-modifier on SD6 to which you can couple. Also, the coupling constants to SD0 and SD1 should be of similar magnitude to those you just established for SD3, so again, start them at 7 Hz. Note that the Karplus correlation suggests that this time, the coupling constant to SD0 should be larger than the coupling to SD1.

24. Add a J-modifier to SD2, and couple it to SD6 using the only remaining uncoupled split (6.52 Hz).
25. Add two J-modifiers to SD2, one coupled to SD0, and the other to SD1, each with a J-value of 7 Hz.
26. Use the J-Value Editor to fine-tune the J-modifiers. Ensure that the J-value for the coupling to SD0 is larger than the J-value for the coupling to SD1. The resulting values are 7.57 Hz (to SD0) and 6.87 Hz (to SD1).
27. Adjust the frequency, height and width of SD2 to better fit the spectrum. Since the spin definition overlaps with others on the right-hand side, use the peaks on the left-hand side to guide your adjustments.

Most of the coupling constants relevant to the last two spin definitions (SD0 and SD1) have already been defined. There is yet another geminal relationship between the associated protons, so couple them using a J-value comparable to the other geminal couplings you have seen (about 13 Hz).

28. Add a J-modifier to SD1, and couple it to SD0 using a J-value of 13 Hz.

Since the signals for SD0 and SD1 in the spectrum both appear simpler than you might expect given that they are coupled to five other signals, it seems reasonable that the remaining two coupling constants for each spin definition will be similar. The uncoupled splits that are available are associated with SD4 (7.44 Hz and 6.88 Hz) and SD5 (7.40 Hz and 6.67 Hz). Thus, the two larger coupling constants must be associated with SD0 and the two smaller with SD1, or vice versa. A quick test of each alternative indicates that the former option yields a better approximation of the spectrum.

29. Add a J-modifier to SD1, and couple it to SD4 using the smaller of the available uncoupled splits (6.88 Hz).
30. Add a J-modifier to SD1, and couple it to SD5 using the smaller of the available uncoupled splits (6.67 Hz).
31. Add a J-modifier to SD0, and couple it to SD4 using the remaining uncoupled split (7.44 Hz).
32. Add a J-modifier to SD0, and couple it to SD5 using the remaining uncoupled split (7.40 Hz).
33. Use the J-Value Editor to fine-tune the J-modifier coupling SD0 to SD1. The resulting value is 13.23 Hz.
34. Adjust the frequency, height and width of SD1 and SD0 to better fit the spectrum.
35. Save the simulation file. If you would like to review the sample simulation of proline, do not replace the file called `Samples\Spin Simulator\Proline.xss`.

10.6 Compound Builder


Introduction

You can use the Compound Builder module to build compound signatures that you can add to your library. Although you can build compound signatures from scratch directly in Compound Builder, you can build complex compounds more accurately if you simulate them in Spin Simulator first. The simulation gives you a good starting point from which to build a high-quality compound signature.

You can only complete the Compound Builder tutorial if you have purchased the Professional edition of Chenomx NMR Suite, or if you are using the Evaluation edition of Chenomx NMR Suite.

Simple Compound - Valine

To become familiar with the various tools available in Compound Builder, you will build a compound signature for a simple compound, valine. You can use the simulation that you created during the Spin Simulator tutorial, or just use the Spin Simulator sample files provided. If you have not yet completed the Spin Simulator tutorial, you may find it helpful to do so before proceeding with Compound Builder.



The simulation of valine is a model representation of the expected spectrum for valine. Compound Builder helps you to turn this model spectrum into a more accurate representation of an experimental spectrum. For this purpose, the four main areas you need to consider are refining the model, simplifying the model and defining transform windows (and the associated valid pH range). For the best results when you add the compound to your library, you also need to add some descriptive metadata about the compound.

Creating a Compound Signature for Valine

1. Open **Compound Builder**.
2. Import the valine simulation.
If you have already completed the Spin Simulator tutorial, import the simulation that you created. Otherwise, import the sample simulation, called `Samples\Spin Simulator\Valine.xss`.
3. When prompted to calibrate the simulation environment, select Spectrum (.cnx), and click the Browse button (...). Select the valine spectrum, called `Samples\Compound Builder\Valine.cnx`.

The compound signature that you build will become part of your compound library. To ensure that the new compound is properly represented in your library, you need to add some basic information about the compound itself. Some information has been filled in automatically based on the calibration spectrum that you selected.

4. Open the Compound Details dialog.
5. In the General tab, enter your name as the Compound Author, and enter a Sample Concentration of 10 mM.
For better integration with compounds based on data from multiple spectrometers, including the Chenomx library, you may want to replace the specific magnet frequency obtained directly from the spectrum (799.81 MHz) with the more general 800 MHz.
6. In the Reference tab, enter the Chenomx Compound ID for valine, 215. If you have access to the InChI or SMILES codes for valine, you can enter them here as well.

7. In the Advanced tab, enter the molecular weight of valine to four decimal places, 117.1463. Also, enter a valid pH of 6.92. Enter 'noesy' as the pulse sequence.
8. Click **OK**.

You will need to set appropriate transform windows before adding the signature to your library. However, setting the default (0.006 ppm) transform windows for all clusters provides an easy way to ensure useful zoom levels when using the Cluster Navigator to move around the spectrum. Also, turning on the subtraction line can give you valuable insight into how well the compound fits the experimental spectrum.

9. Click **Compound > Transform Window > Set Default on All Clusters** to set a default transform window (0.006 ppm wide, centered on the cluster center) for every cluster in the compound.
10. Switch to the Legend sidebar, and turn on the subtraction line (green).

Several of the clusters in the imported simulation are quite simple. Focus on matching the simple clusters to the spectrum first, and then deal with the more complex cluster.

11. Move to the cluster near 3.6 ppm using the Cluster Navigator and select it.
12. Adjust the cluster (blue) to fit the spectrum (black). To further refine the shape, adjust the cluster so that the subtraction line is as flat as possible.

When you are adjusting simple clusters, you can get a very good approximate shape for the cluster automatically, using optimization.

13. Move to the clusters near 1.0 ppm using the Cluster Navigator.
14. Select the cluster near 0.98 ppm and click Optimize Selected Peak Shapes.
15. Review the optimized shape, and refine it as needed.

To obtain accurate quantification when you use this compound to analyze spectra in Profiler, it is helpful to recognize that peak heights have a larger effect on the quantification results than peak width or peak area. Thus, when you refine the shape of a cluster, it is important to ensure that the heights of the peaks match the experimental spectrum particularly well.

16. Select the cluster near 1.03 ppm, optimize it, and refine the shape as needed.
17. Move to the cluster near 2.3 ppm using the Cluster Navigator, select it, and optimize the shape.

In this case, optimization cannot provide a better solution than the one you started with. However, when you review the list of peaks in the cluster (click on each to view its contribution to the cluster), you may notice that there appears to be more peaks present than are strictly necessary. Simulation can produce peaks that do not contribute much to the compound at the current frequency. Part of building a good compound signature is recognizing these extraneous peaks and removing them.

18. Select the first peak in the peak list for the cluster. Using the up and down arrow keys, flip through the list and observe the contributions of each peak.

The first two peaks (2.234 and 2.240 ppm) appear normal, but the next two peaks (both very near 2.243 ppm) appear almost identical, but offset slightly. Near-identical pairs of peaks are excellent candidates for pruning, as the two peaks could be replaced with one without any significant change in the overall shape.

19. Find all of the near-identical pairs of peaks in the cluster, and remove one peak from each pair (6 peaks removed).

Further review reveals that some of the taller peaks in the spectrum are represented by groups of three peaks; one tall central peak flanked by two much shorter peaks. The flanking peaks can also be removed without reducing the shape quality.

20. Find all of the groups of three peaks as described above, and remove the two flanking peaks from each group (6 groups, 12 peaks removed).
21. Optimize the shape of the cluster again.


The results of the second optimization are much better, and require only minimal refinement. In addition, you have now reduced the complexity of the model (from 38 peaks to 20), while improving the overall quality of the fit.

Now that the compound matches the spectrum well, you can add proper transform windows to each cluster, build additional signatures for a pH sensitive compound, or make this signature pH sensitive using other signatures as data points.

22. Save the compound file. If you would like to review the sample signature for valine, do not replace the file called `Samples\Compound Builder\Valine.xcpd`.

Complex Compound - Proline

Now that you have developed a basic familiarity with the tools available in Compound Builder, you are ready to use them to build a compound signature for a more complex compound, proline.



Creating a Compound Signature for Proline

1. Open **Compound Builder**.
2. Import the proline simulation.
If you have already completed the Spin Simulator tutorial, import the simulation that you created. Otherwise, import the sample simulation, called `Samples\Spin Simulator\Proline.xss`.
3. Select Spectrum (`.cnx`), and click the Browse button (...). Select the proline spectrum, called `Samples\Compound Builder\Proline.cnx`.
4. Open the Compound Details dialog.
5. In the General tab, enter your name as the Compound Author, and enter a Sample Concentration of 10 mM.

For better integration with compounds based on data from multiple spectrometers, including the Chenomx library, you may want to replace the specific magnet frequency obtained directly from the spectrum (799.81 MHz) with the more general 800 MHz.

6. In the Reference tab, enter the Chenomx Compound ID for proline, 230. If you have access to the InChI or SMILES codes for proline, you can enter them here as well.
7. In the Advanced tab, enter the molecular weight of proline to four decimal places, 115.1305. Also, enter a valid pH of 6.93. Enter 'noesy' as the pulse sequence.
8. Click OK.
9. Set default transform windows on all clusters, and make sure that the subtraction line (green) is visible.

Several clusters are quite easy to adjust, and do not overlap with any other clusters, so refine these ones first.

10. Optimize each of the clusters near 4.1, 3.4, 3.3 and 2.3 ppm, respectively, and refine them as needed.

The cluster near 2.1 ppm overlaps slightly with the cluster near 2.0 ppm, but the overlap should not cause any particular difficulty. However, the cluster does include some extraneous peaks, so remove them before proceeding.

11. Move to the cluster near 2.1 ppm, select it, and optimize the shape.
12. Review the contributions of each of the peaks to the overall cluster shape. Delete peaks that appear particularly small relative to the peaks immediately adjacent (6 peaks removed).
13. Optimize the cluster, and refine it as needed.

The remaining two clusters are heavily overlapped, and optimizing them at this point yields unsatisfactory results. You will need to simplify both of these clusters before refining their shapes.

14. Review the contributions of each peak in the cluster at 1.979 ppm to the overall shape of the cluster.
 - a. For pairs of near-identical peaks, remove one peak from each pair (8 peaks removed).
 - b. For groups containing peaks much smaller than the immediately adjacent peaks, remove the small peaks (13 more peaks removed).
 - c. To evaluate the effectiveness of removing a particular peak from the cluster, remove it, and then select all the peaks contributing to the apparent peak in the spectrum. Optimize the subset of peaks, and review the results. If all of the remaining peaks contribute significantly to the optimized shape, then move on. Otherwise, undo the optimization, and remove another peak from the subset.
15. Review the contributions of each peak in the cluster at 2.011 ppm in the same way as described for the previous cluster (22 peaks removed).
16. Select both clusters, and optimize them together.
17. Refine the shapes of both clusters as needed.

Now that the compound matches the spectrum well, you can add proper transform windows to each cluster, build additional signatures for a pH sensitive compound, or make this signature pH sensitive using other signatures as data points.

18. Save the compound file. If you would like to review the sample signature for proline, do not replace the file called `Samples\Compound Builder\Proline.xcpd`.

10.7 Building pH Sensitive Compounds

Introduction

You can use Compound Builder to build pH sensitive compounds similar to the Chenomx Reference Compounds. pH sensitive compounds use expressions incorporating the spectrum pH in place of constant values for certain cluster properties, including cluster centers and transform windows. These expressions allow compounds to partially adapt to the spectrum before you profile it, reducing the time and effort needed to accurately identify compounds in the spectrum.

Before you can create a pH sensitive compound for your Library, you will need some data to work with.

Preparing to create a pH sensitive compound

1. Prepare one signature file to be made pH sensitive (the 'master' signature). You can use the valine signature file that you created in the Compound Builder tutorial (see "Simple Compound - Valine" on page 129) as a starting point.
2. Prepare a number of other signature files ('partial' signatures) spanning the pH range over which you would like your pH sensitive compound to be valid. For each of these:
 - Set the valid pH to the pH of the original sample.
 - Make sure each partial signature has the same number of clusters, in the same order, as the master signature. Rearrange the cluster IDs if necessary (see "Manage Cluster IDs" on page 90).
 - Make sure that the cluster centers in each partial signature are as accurate as possible, as these will be used along with the pH values to generate the necessary expressions.

Building a pH Sensitive Compound

Once you have built the necessary signatures, turning them into a pH sensitive compound signature is quite simple. Remember that the quality of the pH sensitive signature that you build is directly related to the quality of the source signatures, and in particular the master signature.

Creating a pH sensitive version of valine

1. Open **Compound Builder**.
2. Open the valine signature.

If you have already completed the Compound Builder tutorial, open the signature that you built. Otherwise, open the sample signature, called `Samples\Compound Builder\Valine.xcpd`.

3. Overlay the valine spectrum, called `Samples\Compound Builder\Valine.cnx`.

4. Open the pH Sensitivities dialog.
5. Click **Load Data Points....**
6. Select the folder called `Samples\Compound Builder\pH Data Points` and click **Set Data Points**.
7. Use the tabs at the top of the window to review each cluster. The red dots, representing the loaded data points, should appear in a reasonably shaped curve.
8. Click **Optimize Lines**.
9. Review each cluster again. Ensure that the blue line representing the cluster centers and the dotted gray lines representing the transform windows are reasonable, and pass through or close to the red data points.
10. When you are satisfied with the proposed fit, click **Accept**.
11. Save the compound file. If you would like to review the sample pH sensitive signature for valine, do not replace the file called `Samples\Compound Builder\Valine-pH.xcpd`.

10.8 Library Manager

Introduction

You can use the Library Manager module to add new compounds to your library, remove existing compounds, and create and edit compound sets using the compounds in your library. Sources of compounds that you can add to your library include the Chenomx Reference Compounds, compound files (.xcpd) created in Compound Builder, compound packs (.pack) from Chenomx or from other users of Chenomx NMR Suite, embedded libraries in legacy spectrum files (.cnx) and legacy Chenomx library files (.cnxlib).

Adding Compounds

To use compound signatures that you create using Compound Builder, you need to add them to your library. Adding compound files (.xcpd) makes the compounds available in your library, letting you use them in Profiler, add them to compound sets, and generally use them alongside the Chenomx library.

Adding Valine and Proline to Your Library

1. Open **Library Manager**.
2. Click **Library > Add Compounds...**
3. Select Importable Files, then click the Browse button (...) and navigate to the folder in which you saved the compound files that you created during the Compound Builder tutorial.

If you have not completed the Compound Builder tutorial, you can use the provided sample files, in the `Samples\Compound Builder` folder.

4. Select your compound files and click **Choose Files**.

If you are using the sample files, select `Valine-pH.xcpd` and `Proline.xcpd` and click **Choose Files**.

5. Click **Next**.

6. Click the **Create associated compound set:** check box, and enter a name for the new Compound Set that you would like the compound to appear in.

If you would not like to automatically add the compound to a new Compound Set, clear the check box.

7. Click **Finish**, and click **OK** in the resulting information dialog.

If you created a Compound Set, it appears in the Compound Set Table. Otherwise, you can find your newly added compounds in the **Compounds** set at the top of the Compound Set Table.

You can use the same method to add compound packs that you may receive from other Chenomx NMR Suite users to your library. Simply select one or more compound pack files (.pack) instead of selecting individual compound files. You can even select a folder containing a mixture of compound files (.xcpd), compound packs (.pack) and legacy files containing library information (.cnx and .cnxlib) to import from a variety of sources at the same time.

10.9 Using Expressions

Building advanced compound signatures in Compound Builder may involve using mathematical expressions to set peak centers, heights or widths, or cluster transform windows. Text fields in which mathematical expressions are allowed contain a button () indicating that expressions are supported for the field. You can click the button to open a simple expression editor, useful for editing more complicated, multiline expressions.

Numbers that you enter in expressions must be in US English style decimal notation (e.g., 6.4378 or -12.038203) or US computerized scientific notation (e.g., 2.0394E+4 or 1.0384E-2). You may add a "+" or "-" character before any number, indicating whether the following number is positive or negative; if you use neither, the number is assumed to be positive. Numeric results of calculations are stored and displayed in double precision (64-bit).

You can use supported text variables in an expression. Supported variables currently include "ph" (returns the pH associated with the current spectrum, as set in Processor), "true" (boolean, equivalent to "1") and "false" (boolean, equivalent to "0"). You cannot declare your own variables.

There are two types of mathematical operations available in expressions, "infix" operations and "prefix" operations. Infix operations involve an operator that occurs inside the two operands, like "5+3", where 5 and 3 are the operands and "+" is the operator. Prefix operations involve an operator that occurs before the operands, like "cos(3.12)" or "max(1,100)", where cos and max are operators and the contents of the parentheses are the operands.

Table 10.3. Infix Operations

Operation	Description	Examples	Result
Boolean Operations			
> or gt	Greater than	ph > 7.00 ph gt 7.00	False for spectrum pH of 6.00, true for 8.00; also false for 7.00
< or lt	Less than	ph < 7.00 ph lt 7.00	True for spectrum pH of 6.00, false for 8.00; also false for 7.00
>= or gte	Greater than or equal to	ph >= 7.00 ph gte 7.00	False for spectrum pH of 6.00, true for 8.00; also true for 7.00
<= or lte	Less than or equal to	ph <= 7.00 ph lte 7.00	True for spectrum pH of 6.00, false for 8.00; also true for 7.00
= or ==	Equal to; use with care when dealing with non-integer values	ph = 7.00 ph == 7.00	True for spectrum pH of 7.00, false for anything else
!=	Not equal to	ph != 7.00	False for spectrum pH of 7.00, true for anything else
Mathematical Operations			
+	Addition	3.202 + 5.47 1.204E-4 + 4.473E-4	8.672 5.677E-4
-	Subtraction	9.01 - 2.356 5.469E2 - 1.24E2	6.654 422.9
*	Multiplication	1.0438 * 3.389 5.469E2 * 1.248E2	3.5374382 68253.12
/	Division, returning a value with a decimal component	5 / 3 7.20 / 3.20	1.666666666666667 2.25
//	Integer division; rounds operands to the nearest integer, then returns the integer result of the division, throwing away the remainder	5 // 3 7.20 // 3.20	1 2
%	Integer modulo; rounds operands to the nearest integer, then returns the remainder after performing integer division	5 % 3 7.20 % 3.20	2 1
^	Exponentiation; the first operand is raised to the power of the second operand	5 ^ 3	125

Table 10.4. Prefix Operations

Operation	Description	Examples	Results
Logical Operations			
and(arg1, arg2)	Returns true if and only if both supplied values are true; returns false otherwise	and(ph != 7.00, ph lt 7.65)	True for spectrum pH of 6.00, false for 8.00
if(arg1, arg2, arg3)	Based on the criterion arg1, returns arg2 if arg1 is true; returns arg3 otherwise	if(ph lt 7.00, 2.103, 3.032)	2.103 for spectrum pH of 6.00, 3.032 for 8.00
not(arg1)	Returns true if the supplied value is false and false if it is true	not(ph = 8.00)	True for spectrum pH of 6.00, false for 8.00

Operation	Description	Examples	Results
<code>or(arg1, arg2)</code>	Returns false if and only if both supplied values are false; returns true otherwise	<code>or(ph lt 6, ph gt 8)</code>	True for spectrum pH of 4.00 or 9.00, false for 7.00
Mathematical Operations			
<code>abs(arg1)</code>	Returns the absolute value of the supplied argument	<code>abs(3.2084)</code> <code>abs(-1.391)</code>	3.2084 1.391
<code>ceil(arg1)</code>	Returns the smallest integer that is not less than the supplied value	<code>ceil(3.2084)</code> <code>ceil(-1.391)</code>	4 -1
<code>exp(arg1)</code>	Returns Euler's number (<i>e</i>) raised to the power of the supplied value	<code>exp(1)</code>	2.7182818284590455 (<i>e</i>)
<code>floor(arg1)</code>	Returns the largest integer that is not greater than the supplied value	<code>floor(3.2084)</code> <code>floor(-1.391)</code>	3 -2
<code>log(arg1)</code>	Returns the base-10 logarithm of the supplied value	<code>log(10)</code>	1
<code>logn(arg1)</code>	Returns the base- <i>e</i> (natural) logarithm of the supplied value	<code>logn(2.71828)</code>	0.999999327347282 (~1)
<code>max(arg1, arg2)</code>	Returns the greater of the two supplied values	<code>max(2.23, 4.30)</code> <code>max(2.23, -4.30)</code>	4.30 2.23
<code>min(arg1, arg2)</code>	Returns the smaller of the two supplied values	<code>min(2.23, 4.30)</code> <code>min(2.23, -4.30)</code>	2.23 -4.30
<code>pow(arg1, arg2)</code>	Returns the first supplied value raised to the power of the second supplied value; equivalent to the "^" infix operation	<code>pow(5, 3)</code>	125
<code>sqrt(arg1)</code>	Returns the square root of the supplied value	<code>sqrt(25)</code>	5
Trigonometric Operations			
<code>acos(arg1)</code>	For a value representing an angle in radians, returns its arc cosine in the range 0 through π (pi)	<code>acos(-1)</code> <code>acos(1)</code>	3.141592653589793 (π) 0
<code>asin(arg1)</code>	For a value representing an angle in radians, returns its arc sine in the range $-\pi/2$ through $+\pi/2$	<code>asin(1)</code> <code>asin(0)</code>	1.5707963267948966 ($\pi/2$) 0
<code>atan(arg1)</code>	For a value representing an angle in radians, returns its arc tangent in the range $-\pi/2$ through $+\pi/2$	<code>atan(1)</code> <code>atan(0)</code>	0.7853981633974483 ($\pi/4$) 0
<code>cos(arg1)</code>	For a value representing an angle in radians, returns its cosine	<code>cos(0)</code> <code>cos(3.14159)</code>	1 -0.9999999999964793 (~1)
<code>sin(arg1)</code>	For a value representing an angle in radians, returns its sine of the value	<code>sin(0)</code> <code>sin(3.14159)</code>	0 -2.65358979335273E-6 (~0)
<code>tan(arg1)</code>	For a value representing an angle in radians, returns its tangent of the value	<code>tan(0)</code> <code>tan(0.785398)</code>	0 0.9999996732051568 (~1)

Chemistry-Related Expressions

In addition to the basic operations described above, Chenomx NMR Suite supports a number of chemistry-related expressions, allowing a sort of shorthand notation for certain commonly encountered situations.

Cluster Center

The cluster center operation computes the center (chemical shift) of a cluster based on a given pH value and a sigmoid function with potentially more than one set of inflection points and "plateaus". For the purposes of this operation, inflection points or "pka" values are interpreted as acid dissociation constants for the compound being modeled, and the sigmoid "plateaus" describe the ppm locations (chemical shifts) of each pure species associated with the dissolved compound.

The general form of the cluster center operation is:

```
clusterCenter(ph, n, pka1, ..., pkaN, shift1, ..., shiftN+1)
```

where:

- `ph` is the pH of the solution the compound is dissolved in, represented by the spectrum pH for the currently overlaid spectrum
- `n` is the number of pka values for the modeled compound; must be at least 1.
- `pka1, ..., pkaN` are the compound's pka values (dissociation constants), in ascending order; there must be exactly `n` of these. These values are in the same units as pH.
- `shift1, ..., shiftN+1` are the pure species chemical shifts associated with the given pka values and the cluster being modelled. They are in the same order as the their associated pkas. Note that there must be exactly `n+1` of these. These values are expressed in ppm.

The following expression is an example of a valid use of the cluster center operation:

```
clusterCenter(ph, 3, 0.6902281, 6.851296, 9.289176, 4.3607736,
0.9786361, 0.9779636, 0.9253154)
```

The above example describes a cluster from a compound with 3 pKas of about 0.69, 6.85 and 9.29, and plateaus at about 4.361, 0.979, 0.978 and 0.925 ppm. The resulting cluster center is 0.9781524039197456 when the spectrum pH (`ph`) is 6.92

Transform Window

The transform window operation, performs the cluster center operation described above at several pH values and uses the results to determine error bounds for the computed center of the modeled cluster. In effect, it computes the cluster's transform window.

Transform windows are based on the underlying assumption that there is always an error in any pH value passed into the `clusterCenter` operation, termed `errorPh`. Since a cluster's center moves primarily with pH, `errorPh` can then be used to represent a boundary on the amount of error in our overall cluster center calculations.

More specifically, upper and lower error bounds are determined for the cluster center by recalculating the cluster center at `ph + pHError` and `ph - pHError`. If the cluster center function has any local extrema that fall within that range, the cluster center at each extremum is also calculated. The extrema must be provided as parameters to the transform window operation; see the example below.

The transform window operation returns either the highest or lowest of all the computed centers, depending on whether or not the "upper" parameter is set to true. You can also set a "minDeltaPPM" parameter, which ensures that the returned boundary is at least a minimum distance away from the unmodified clusterCenter value. This is useful for setting a "minimum" transform window on clusters that do not move much with pH.

The general form of the cluster center operation is:

```
transformWindow(ph, n, pka1, ..., pkaN, shift1, ..., shiftN+1, upper,
errorPh, minDeltaPpm, x, extrema1, ..., extremax)
```

where:

- The first four parts are the same as for the cluster center operation.
- `upper` is non-zero (true) if you want the operation to return the upper bound of the calculated transform window, and zero (false) if you want the lower.
- `errorPh` is the maximum amount of pH error to incorporate into the calculated transform window.
- `minDeltaPpm` is the minimum distance (in ppm) between the unmodified cluster center and the cluster center at $\text{pH} \pm \text{errorPh}$, or the minimum size of one half of the transform window.
- `x` is the number of extrema that the clusterCenter function (represented in the first 5 arguments) has.
- `extrema1...extremax` the extrema locations (in pH) of the clusterCenters extrema points. There must be "x" of these, but x may be zero.

The following expression is an example of a valid use of the transform window operation:

```
transformWindow(ph, 3, 0.6902281, 6.851296, 9.289176, 4.3607736,
0.9786361, 0.9779636, 0.9253154, 0, 0.5, 0.02, 0)
```

The above example refers to the same cluster as the cluster center example above. The expression calculates the lower transform bound of the cluster (since `upper = 0`), subject to a pH error of 0.5 and a minimum distance of 0.02 ppm, with no extrema. The resulting lower bound is 0.9581524039197455 when the spectrum pH is 6.92.

Bounded Cluster Center and Transform Window

The cluster center and transform window operation described above rely on pKas and chemical shifts. When the pKas and shifts are determined using observed data and an automatic fitting routine, the results of the pKa/shift method outside of a certain known pH range is not reliable, as they would require inappropriately extrapolating cluster centers from an arbitrary sigmoidal function.

The bounded forms of these operation are designed to handle such situations. The bounded forms are identical to the unbounded forms except for two extra arguments specifying the lower and upper pH boundaries within which the operation is considered valid.

The value returned by these bounded operations does not change at all outside the valid pH range. Results from pHs below the valid lower bound are the same as the result at the lower bound, and results from pHs above the valid upper bound are the same as the result at the upper bound.

The general form of the bounded cluster center operation is:

```
bClusterCenter(lowerBound, upperBound, ph, n, pka1, ..., pkaN, shift1,
..., shiftN+1)
```

where:

- lowerBound is the lowest pH at which the operation is considered valid.
- upperBound is the highest pH at which the operation is considered valid.
- The remaining arguments are the same as those found in the unbounded cluster center operation (see “Cluster Center” on page 138).

The following expression is an example of a valid use of the bounded cluster center operation:

```
bClusterCenter(4.07, 9.02, ph, 3, 0.6902281, 6.851296, 9.289176,
4.3607736, 0.9786361, 0.9779636, 0.9253154)
```

The above example describes the same cluster as in the cluster center example (see “Cluster Center” on page 138), now subject to a pH range of 4.07 to 9.02.

The general form of the bounded transform operation is:

```
bTransformWindow(lowerBound, upperBound, ph, n, pka1, ..., pkaN, shift1,
..., shiftN+1, upper, errorPh, minDeltaPpm, x, extrema1, ..., extremax)
```

where:

- lowerBound is the lowest pH at which the operation is considered valid.
- upperBound is the highest pH at which the operation is considered valid.
- The remaining arguments are the same as those found in the unbounded transform window operation (see “Transform Window” on page 138).

The following expression is an example of a valid use of the bounded cluster center operation:

```
bTransformWindow(4.07, 9.02, ph, 3, 0.6902281, 6.851296, 9.289176,
4.3607736, 0.9786361, 0.9779636, 0.9253154, 0, 0.5, 0.02, 0)
```

The above example describes the same transform bound as in the unbounded transform window example (see “Transform Window” on page 138), now subject to a pH range of 4.07 to 9.02.

Glossary

#

1r *See Bruker processed spectrum data..*

A

acquisition time	Refers to the period of active data collection in an NMR pulse sequence. The recommended acquisition time for spectra intended for use with Chenomx NMR Suite is 4 s.
activation	The process of licensing Chenomx NMR Suite for legitimate use, involving the installation of a Chenomx license file (.cnxlic).
adding compounds	Refers to expanding the compound library by adding signatures from one or more compound sources, including the Chenomx Reference Compounds, compound packs (.pack), compound files (.xcpd), and legacy Chenomx NMR Suite files (.cnx and .cnxlib).
animated zooming	Smooth transitions between zoom settings in the Spectrum View of Processor, Profiler, Compound Builder or Spin Simulator.
application log	Text file that Chenomx NMR Suite generates for troubleshooting purposes.
auto baseline correction	Option available on importing spectra into Processor or Profiler to automatically determine placement of breakpoints for spline-based baseline correction. Equivalent to selecting Auto Spline during baseline correction in Processor.
auto phasing	Option available on importing spectra into Processor or Profiler to automatically determine phase angles. Equivalent to selecting Auto during phasing in Processor.
auto zoom	Adjusts the vertical scale of the Spectrum View in Processor, Profiler, Compound Builder and Spin Simulator to include all visible data from available lines, including the spectrum, spectrum overlay, sum, compound, simulation, subtraction, CSI and pinned lines (where applicable).

B

baseline correction	Processing layer available in Processor. Also, the process of compensating for distortions in the non-signal (noise) regions of a spectrum.
best guess fit	Simple profiling technique involving small transforms of indistinct clusters to yield an aesthetic improvement in the profile, in the absence of concrete information on which to base a more accurate approach.

bin	In spectral binning, a region of the source spectrum that has been reduced to a single value representing the area under the curve in that region.
Bruker processed spectrum data	Processed (Fourier-transformed) data from spectra acquired using Bruker spectrometers and processed using the XWIN-NMR software package developed by Bruker.

C

chemical shift and shape indicator	Refers to a signal in the spectrum used to calibrate the chemical shifts, lineshapes and concentrations of compounds in your Library. Also refers to the compound responsible for the signal.
chemical shift equivalent	Describes a relationship within a molecule between nuclei that are interchangeable by symmetry operations or by a rapid process.
Chenomx compound ID	Arbitrary positive integer assigned to a compound by Chenomx staff, which uniquely identifies the compound in the Chenomx compound library.
Chenomx Reference Compounds	Compound signatures distributed with Chenomx NMR Suite, comprising the complete collection of validated signatures available from Chenomx at all supported frequencies (400, 500, 600, 700 and 800 MHz).
cluster	Group of peaks that you can move or otherwise adjust as a unit in Compound Builder and Profiler. Clusters can correspond to the multiplets familiar to NMR spectroscopists, or to arbitrary groupings created in Compound Builder.
Cluster Navigator	Graphical tool in Compound Builder and Profiler that lets you adjust the Spectrum View to include specific clusters present in the current compound.
.cnx	<i>See spectrum file..</i>
.cnxlic	<i>See license file..</i>
Compound Builder	Module in Chenomx NMR Suite that lets you prepare and modify compound signatures (.xcpd) for use in your Library.
compound details	Information about a compound signature created in Compound Builder. Some of the information appears directly in the reference panel, some provides access to other information for display in the reference panel, and some indexes the compound within a library.
compound file	Chenomx file format (.xcpd) designed to store the mathematical descriptions of compound signatures created with Compound Builder. The compound library of Chenomx NMR Suite consists of a pool of compound files stored in a folder that you can designate via Library Manager preferences.

compound library	A pool of compound files (.xcpd) stored in a common folder and indexed by Library Manager, used to profile spectra in Profiler.
compound line	In Compound Builder, the line in the Spectrum View indicating the overall shape of the compound signature, colored red by default.
compound line fit	Simple profiling technique involving matching the currently selected compound line (blue) to the spectrum line (black).
compound pack	Chenomx file format (.pack) designed to store a collection of compound files (.xcpd) for exchange among users of Chenomx NMR Suite.
compound set	List of compounds defined in Library Manager that appears in Profiler for use during profiling. Can refer to any set appearing in the Compound Sets view, or to fixed lists of compounds specifically as opposed to the automatically-updating lists known as smart compound sets.
Compound Sets view	Sidebar view appearing in Profiler and Library Manager that contains a list of available compound sets, based on the current compound library. Lets you modify the displayed contents of the Compound Table to aid in profiling spectra.
Compound Table	Table appearing in Profiler and Library Manager. Contains information on compounds in your Library, including names, compound IDs, valid magnet frequencies and pH ranges, InChI and SMILES strings, and so on. Also contains maximum and measured concentrations when available in Profiler.
coupled split	In Spin Simulator, a J-modifier that links two spin definitions via a common split count and J-value, resulting in more accurate calculation of complex coupling interactions. Equivalent to J-couplings in previous versions of Chenomx NMR Suite.
creatinine	pH indicator supported by Chenomx NMR Suite. Most accurate for pH 3.75 to 5.50.
CSI	<i>See chemical shift and shape indicator..</i>
CSI editor	Mode in Processor that lets you define the properties of the CSI for a spectrum, including the experimental concentration (in mM), peak center (in ppm), height (in si) and width (in Hz).
CSI line	In Processor, the line in the Spectrum View illustrating the current CSI settings as defined in the CSI editor.

D

default transform window	In Compound Builder, the limits applied to a cluster when such limits are not otherwise specified (e.g., through the menu item Compound > Transform Window > Set Default on Selected Clusters). Equivalent to a range of 0.006 ppm centered on the
---------------------------------	--

position of the cluster when the default is applied, or 0.003 ppm in either direction.

Delta

NMR software package developed by JEOL to run their spectrometers and process the resulting data.

DFTMP

Abbreviation of [difluoro(trimethylsilyl)methyl]phosphonate, a pH indicator supported by Chenomx NMR Suite. Most accurate for pH 5.00 to 7.50.

DSS

Abbreviation of 2,2-dimethyl-2-silapentane-5-sulfonate, the CSI recommended by Chenomx for most samples. All other CSIs in Chenomx NMR Suite are defined relative to DSS.

E

.eps

See encapsulated PostScript file..

encapsulated PostScript file

PostScript® file which contains image size and layout information according to Adobe's Document Structuring Convention. Used primarily for electronic and professional desktop publishing.

exporting spectrum graphic

In Profiler, refers to creating an image file showing the current Spectrum View. Used to generate images in various formats for use in presentations or publications.

exporting profiled compounds

In Profiler, refers to creating a compound pack file (.pack) or set of compound files (.xcpd) based on profiled compounds in the currently open spectrum.

exporting profiled data

In Profiler, refers to creating a tab-delimited text file describing the data stored in one or more profiled Chenomx spectrum files (.cnx). The resulting file is suitable for use in a variety of statistical analysis and spreadsheet software.

F

.fid

See NMRPipe spectrum data..

fid

See spectrum data..

fit automatically

In Profiler, function that provides suggested concentration and transforms for the currently selected compound.

fit percentages

In Profiler, the percentage of the area under the spectrum line in the selected region. Reported for the spectrum line and the sum line.

frequency domain size

The number of data points in a Chenomx spectrum, after importing spectrum data and applying zero-filling.

frequency filter

In Profiler, a quick search based on the position of peaks in the spectrum, in ppm. You can specify a single frequency or a range of frequencies, resulting in a list of all compounds having at least

one transform window that overlaps the specified frequency or range. A more limited variety of a text filter.

.ft2

See NMRPipe processed spectrum data..

G

generate cluster for region

In Compound Builder, an operation that automatically creates a series of peaks to match the selected region, by minimizing the average intensity of the subtraction line.

.gif

See graphics interchange format..

graphics interchange format

Raster graphic format with lossless compression, 8-bit palette color depth and simple transparency support. Suitable for image distribution on the internet.

grouping J-modifiers

In Spin Simulator, combines multiple uncoupled J-modifiers into a single one with the same total split count, and the average J-value for the group. For example, grouping three uncoupled splits, with split counts of 1, 2 and 3 and J-values of 4, 6 and 7 Hz, will give a single uncoupled split with a split count of 6 (1+2+3) and a J-value of 5.67 Hz (average of 4, 6 and 7 Hz).

group peaks as cluster

In Compound Builder, defines a collection of peaks as belonging to the same cluster, so that they will transform together, appear under the same entry in the Cluster Navigator and obey the same transform windows.

I

imidazole

pH indicator supported by Chenomx NMR Suite. Most accurate for pH 5.50 to 8.80.

InChI

See IUPAC international chemical identifier..

IUPAC international chemical identifier

A string of characters that uniquely represents a chemical compound based on its structure. Developed and approved by the International Union of Pure and Applied Chemists (IUPAC) and the National Institute of Standards and Technology (NIST).

J

Java runtime environment

Software based on a specification by Sun Microsystems that lets a computer run Java applications like Chenomx NMR Suite.

JCAMP-DX spectrum data

Spectrum data (unprocessed or processed) stored in the JCAMP-DX data interchange format.

.jdf

See JEOL Delta spectrum data..

.jdx

See JCAMP-DX spectrum data..

JEOL Delta spectrum data	Data from spectra acquired using JEOL spectrometers and processed (where applicable) using the Delta software package developed by JEOL.
J-modifier	In Spin Simulator, a parameter modifying a spin definition by splitting it to form more complex shapes. Composed of a J-value and a split count, and includes both coupled and uncoupled splits.
J-value	In Spin Simulator, defines the magnitude of the split defined by a J-modifier. Equivalent to a coupling constant for coupled splits.
J-value editor	In Spin Simulator, a graphical control that lets you vary the J-value for the currently selected J-modifier.

L

Legend	Sidebar view providing an interactive display of the various lines appearing in the Spectrum View
library	Pool of compound signatures (.xcpd) available for use in Chenomx NMR Suite. The Library replaces the Chenomx Compound Libraries available in previous versions of Chenomx NMR Suite.
Library Manager	Module in Chenomx NMR Suite that lets you manage your Library by adding or removing compounds, updating existing compounds with new versions, and creating, updating and deleting compound sets.
license file	Chenomx file format (.cnxlic) used to activate Chenomx NMR Suite.
line broadening	Processing layer available in Processor. A mathematical operation that multiplies a fid by an exponential function in the time domain, averaging out instrumental noise at the expense of increasing linewidths. The exponential function is represented by a value measured in Hz, approximately corresponding to the increase in linewidth that will result from the operation.

M

magnetic-equivalent	Describes a relationship within a molecule between nuclei that couple equally to any nucleus in the same spin system. Nuclei that are magnetic-equivalent are chemical shift equivalent by definition.
matched clusters	Calculated value visible in the Compound Table and exportable using Export Profiled Data, the matched clusters for a compound are clusters that closely matched the spectrum line, as determined by comparing the compound and spectrum line shapes near each cluster.
maximum concentration	Estimated upper limit of a compound's concentration, calculated automatically by Profiler. You can use maximum concentrations to guide analysis of a spectrum in Profiler.

N

ni	<i>See normalized intensity..</i>
NMRPipe	A free software package designed to process NMR spectra.
NMRPipe processed spectrum data	Data from spectra processed using the NMRPipe software package.
NMRPipe spectrum data	Data from spectra stored using the NMRPipe software package.
no fit	Simple profiling technique applied to clusters for which no information is available to guide profiling. The cluster is simply not modified.
normalization	In Profiler, refers to dividing the area of each bin in a spectral binning session by a common denominator, letting you compare binning results across multiple spectra. You can normalize bins relative to the total area under the spectrum line or to units of standardized area (sa).
normalized intensity	Refers to an intensity relative to that of the tallest peak in the spectrum, that is, the point with the largest intensity. Using this scale, if the tallest peak in the spectrum is the creatinine methyl signal, and the urea signal is half as tall as the creatinine signal, then the intensity of the urea signal is 0.5 ni.

O

optimize selected peak shapes	In Compound Builder, an operation that automatically adjusts the centers, heights and widths of the selected peaks to match the spectrum overlay, by minimizing the average intensity of the subtraction line.
overlay lines	Lines in the Spectrum View corresponding to spectra loaded for comparison to the currently loaded spectrum in Processor and Profiler, the current simulation in Simulator, or the current compound signature in Compound Builder.

P

.pack	<i>See compound pack..</i>
.pdf	<i>See portable document format..</i>
peak	Mathematical model of a signal appearing in a spectrum. Chenomx NMR Suite models all peaks as pure Lorentzian curves.
pH editor	Mode used in Processor to determine the sample pH for a spectrum, based on the approximate peak positions of one or more pH indicators.

pH indicator	Compound that you add a sample to let you detect the sample pH. Indicators supported by Chenomx NMR Suite, like creatinine, imidazole and DFTMP, have signals whose chemical shifts change with the pH of the sample. Sample pH can thus be determined by measuring the signals' chemical shifts.
phase correction	<i>See phasing..</i>
phasefile	Data from spectra processed using the VNMR or VNMRj software packages developed by Varian.
phasing	Processing layer available in Processor. Also, the process of choosing coefficients for the parts of a linear combination of real and imaginary spectra, derived from the source spectrum and varying linearly with frequency, to optimize the appearance of the source spectrum.
pinned compound	In Profiler, a compound with a check in the Pinned column of the Compound Table. The signatures of pinned compounds always appear in the Spectrum View, regardless of what is selected in the Compound Table.
pinned lines	In Profiler, the lines used to indicate pinned compounds. Colored to match the default selection line (blue by default), or the custom color for each compound, if one has been set.
.png	<i>See portable network graphics..</i>
portable document format	Open file format created by Adobe Systems Incorporated. Vector graphic format that stores document information in a device independent and resolution independent way. Suitable for multiplatform viewing and printing.
portable network graphics	Raster graphic format with lossless compression, 48-bit true-color depth, variable transparency and gamma correction. Suitable for image distribution on the internet.
preprocessing	In Processor, an entry that appears in the Processing History to indicate prior processing in third-party software packages. Appears whenever processed data is imported (1r, phasefile, .ft2, .jdf, .jdx).
preview line	In Profiler, the line in the Spectrum View displayed for a compound whose concentration is not set, used to estimate the likelihood of the compound being present in the spectrum.
Processing History	Sidebar view showing the processing layers that have been applied to the currently open spectrum. Also lets you review the details of each layer and remove layers as needed.
processing layer	In Processor, one of several modifications that you can make to imported spectrum data to refine the appearance of the spectrum and compensate for defects in acquisition. Includes phasing, baseline correction, line broadening, reference deconvolution, water deletion and reverse spectrum.

processing tools	In Processor, lines and shapes in the Spectrum View illustrating the settings of various processing layers, including phasing, baseline correction and reference deconvolution.
Processor	Module in Chenomx NMR Suite that lets you import NMR data from various sources, add various processing effects like phasing, baseline correction, line broadening and reference deconvolution, and save the results as a spectrum file (.cnx).
profile	Part of a spectrum file (.cnx) describing the results of analyzing the spectrum in Profiler. Includes the list of identified compounds, along with their determined concentrations and transforms. Represented in Profiler as the Compound Table.
profiled compounds	In Profiler, compounds that have been modified in some way from their original state, including concentrations, transforms, custom colors, pins and stars. Also, compounds that appear in the Profiled Compounds set in the Compound Sets view.
Profiler	Module in Chenomx NMR Suite that lets you analyze spectra by comparison with compound signatures from your Library, creating a profile of the spectrum.
proxy server	Used to mediate communication between a computer on one network and computers on other networks. Sometimes used to improve network operations by caching frequently-accessed data, or to improve network security by acting as a secured gateway through which other computers access other networks (like the Internet).
.ps	<i>See PostScript..</i>
pulse sequence	Specific series of radiofrequency bursts and pauses in an NMR experiment, designed to excite nuclei in a particular way before beginning to detect the resulting magnetization signal (during the acquisition time).

Q

quick search	Filter for the current view of the Compound Table based on entered text. Displayed results can match to any or all of compound name, compound comment, compound ID or location in the spectrum (in ppm).
---------------------	--

R

Reference Compound Sets	Compound sets containing all Chenomx Reference Compounds currently appearing in your library. The Reference Compound Sets are created and removed automatically by Chenomx NMR Suite as you add or remove Chenomx Reference Compounds from your library.
--------------------------------	--

reference deconvolution	Processing layer available in Processor. Lets you reconstruct an ideal spectrum by removing lineshape distortions caused by systematic variations in the magnetic field of the spectrometer.
Reference Panel	Sidebar view that displays additional data for the compound currently selected in the Compound Table in Profiler, or for the compound ID specified in the compound details in Compound Builder (if any).
relaxation delay	Refers to a pause preceding active data collection in an NMR pulse sequence, to ensure that residual signal from the radiofrequency pulse used to excite the sample is not visible in the spectrum. The recommended relaxation delay for spectra intended for use with Chenomx NMR Suite is 1 s.
residual binning	Variation of spectral binning in which the subtraction line is binned instead of the spectrum line.
reverse spectrum	Processing layer available in Processor. Inverts the order of the data points in a spectrum in the frequency domain, effectively flipping the spectrum horizontally. Only necessary for spectra from certain Bruker spectrometers.
RGB color model	Additive color model used to represent color as a combination of red, green and blue light for display on electronic devices (compare the subtractive CMYK color model used to reproduce color in print). Usually encoded for computer programs as a triplet of 8-bit unsigned integers, one for each color, represented in decimal (0 to 255) or hexadecimal (00 to FF) notation, where 0 or 00 represents no contribution from a color, while 255 or FF represents maximum contribution. In Profiler, colors are expressed as 6-digit hexadecimal numbers, two digits each for red, green and blue, in that order. Using this system, red is expressed as FF0000, green as 00FF00, and blue as 0000FF. More complex colors can also be expressed using this notation; 6495ED is cornflower blue, while FF69B4 is hot pink.
RGB Hex	<i>See RGB color model..</i>

S

sa	<i>See standardized area..</i>
scalable vector graphics	Open standard created by the World Wide Web Consortium. XML-based format for describing two-dimensional vector graphics. Suitable for image distribution on the internet.
select region	Tool that lets you temporarily mark a specific region of the spectrum in the Spectrum View. Allows more precise zooming, defining frequency ranges for quick searches, measuring areas under various lines, marking clusters, peaks or spin definitions for selection, marking a region in which to generate a new cluster, and so on.

selection line	Line in the Spectrum View in Compound Builder and Spin Simulator corresponding to the currently selected peak, cluster or spin definition, colored blue by default.
send to Compound Builder	In Library Manager, opens the compound file (.xcpd) of the selected compound in Compound Builder for modification or review.
send to Processor	In Profiler, opens the current spectrum in Processor to let you modify the processing layers, CSI settings, spectrum pH or spectrum details. The details of the current profile (concentrations, transforms, etc.) are retained, even if you then save and close the spectrum from Processor.
send to Profiler	In Processor, opens the current spectrum in Profiler to start (or continue) profiling. Any changes to the processing layers, CSI setting, spectrum pH, or spectrum details are retained, even if you then save and close the spectrum from Profiler.
si	<i>See standardized intensity..</i>
Sidebar	Dedicated area of the screen that displays frequently-used information related to the Spectrum View. Available views include Compound Sets (Profiler), Legend (all modules except Library Manager), Processing History (Processor), Reference Panel (Profiler and Compound Builder) and Simulation (Spin Simulator).
simplified molecular input line entry specification	A string of characters that uniquely represents a chemical compound based on its structure. Developed by David Weininger and extended by Daylight Chemical Information Systems, Inc.
Simulation	Sidebar view that lets you view and edit spin definition properties, and create, edit or remove J-modifiers
simulation details	In Spin Simulator, information about the current simulation, including simulation frequency in MHz, number of spin systems, number of transitions modeled and number of J-modifiers defined, including coupled and uncoupled splits.
simulation file	Chenomx file format (.xss) designed to store the parameters of simulations created with Spin Simulator.
simulation line	Line in the Spectrum View in Spin Simulator corresponding to the overall lineshape of the current spin simulation, colored red by default.
smart compound set	Compound set defined in Library Manager as a set of rules for dynamically selecting compounds from your Library.
SMILES	<i>See simplified molecular input line entry specification..</i>
spectral binning	Data reduction method widely used to prepare spectra for statistical analysis. Involves dividing spectra into pieces or bins, and considering only the integrated area in each bin for further analysis.

spectral width	The full range of frequencies represented in the current spectrum, in ppm.
spectrum details	Information about the currently open spectrum in Processor and Profiler. Includes acquisition parameters stored by the spectrometer, CSI concentration and spectrum pH determined in Processor, and text comments associated with the spectrum.
spectrum file	Chenomx file format (.cnx) designed to store raw spectrum data (fid), processing layers added using Processor and profiles created using Profiler. You can use spectrum files as overlays in any module with a Spectrum View.
spectrum line	Line in the Spectrum View in Processor and Profiler corresponding to the experimental data of the currently open spectrum, colored black by default.
spectrum overlay	Spectrum files (.cnx) appearing in the Spectrum View as supplements to the currently open file. In Spin Simulator and Compound Builder, overlays guide you in constructing accurate simulations and compound signatures. In Processor and Profiler, they let you compare the currently open spectrum with others in the same dataset or other datasets.
spectrum overlay line	Line in the Spectrum View in Compound Builder and Spin Simulator corresponding to the spectrum overlay, colored black by default.
spectrum thumbnail	Miniature version of the Spectrum View that you can use to quickly navigate or zoom to different regions of the spectrum, regardless of the currently viewed region or zoom level.
Spectrum View	Main display area in Processor, Profiler, Compound Builder and Spin Simulator. Contains the spectrum graph and several tools for manipulating the displayed information, including the Cluster Navigator or Spin Navigator, and controls for adjusting the heights and frequencies of displayed peaks, clusters or spin definitions.
spin definition	Basic element of a simulation created in Spin Simulator. Each spin definition consists of a number of protons, a center in ppm and a width in sw.
Spin Navigator	Adjusts the Spectrum View in Spin Simulator to include specific spin definitions in the current simulation. Analogous to the Cluster Navigator.
Spin Simulator	Module in Chenomx NMR Suite that lets you prepare and modify spin simulations (.xss).
spin system	Grouping of spin definitions in Spin Simulator that are uncoupled to any other spin definitions. Chemically, a spin system corresponds to a group of protons in a compound that do not interact with any other protons in the compound in a manner detectable by NMR.

spline	A curved line formed by two or more vertices or breakpoints and a mathematical formula that describes the curves between them.
split count	One of the two parameters that define a J-modifier. For coupled splits, the split count is predefined as the number of protons associated with the spin definition to which the current spin definition is coupled. For uncoupled splits, you may consider the split count as the number of virtual protons to which the spin definition would need to be coupled to show the desired splitting pattern.
standardized area	Refers to an area relative to that of a theoretical DSS methyl peak at a concentration of 0.50 millimoles per litre (mM).
standardized intensity	Refers to an intensity relative to that of a theoretical DSS methyl peak at a concentration of 0.50 millimoles per litre (mM). Using this scale, if you select DSS as your CSI, and indicate a DSS concentration of 1.00 mM, then the intensity of the DSS methyl signal in that spectrum is 2 si, since 1.00 mM is twice the standard concentration of 0.50 mM.
standardized width	Refers to a linewidth relative to that of a theoretical DSS methyl peak at a concentration of 0.50 millimoles per litre (mM).
starred compound	In Profiler, a compound for which you have toggled the Star column on (yellow).
subtraction line	Line in the Spectrum View in Profiler, Compound Builder and Spin Simulator showing the difference between the black line (spectrum or overlay line) and the red line (sum, compound or simulation line), colored green by default.
subtraction line fit	Advanced profiling technique giving very accurate results. For a subtraction line fit, you manipulate concentrations and transforms for a compound so that the subtraction line (green) approximates a normal spectrum, indicating that the influence of the compound has been efficiently subtracted from the original spectrum.
sum line	Line in the Spectrum View in Profiler representing the sum of all of the compound lines in the profile of the currently open spectrum, colored red by default.
sum line fit	Profiling technique needed for fitting overlapped regions. For a sum line fit, you manipulate clusters overlapping in a particular region so that the sum line (red) in that region best approximates the spectrum line (black). Other clusters for the same compounds that occur in other, potentially less overlapped regions will often contribute to this technique.
.svg	<i>See scalable vector graphics..</i>
sw	<i>See standardized width..</i>

T

targeted binning	Variation of spectral binning in which the bins are defined as the transform windows of each cluster of each compound in a compound set.
targeted profiling	Data reduction method embodied in the Profiler module, involving comparing a spectrum to compound signatures in your Library, thus reducing spectral data to a list of quantified compounds.
text filter	A quick search based on entered text. Filters the Compound Table by matching compound name, compound comment, compound ID or location in the spectrum (Profiler only, in ppm).
time domain size	The number of data points originally acquired for a spectrum, if it originated from a raw fid, or the number of points originally saved in a spectrum, if it was processed using a third-party software package.
TMSP	<i>See TSP.</i>
transform	Translation of a cluster in a compound signature along the frequency axis during profiling. Transforms are necessary to compensate for variation in cluster frequency due to pH, ionic strength and other matrix effects.
transform summaries	In Profiler, values for profiled compounds indicating whether clusters for each compound are within their respective transform windows, or outside them. May be exported using Export Profiled Data.
transform window	Limit placed on the extent of transforms possible on a particular cluster.
transitions	Lorentzian peaks corresponding to energy transfer among spins in a spin system. The result of a simulation in Spin Simulator is a series of transitions, the sum of which is the simulation line (red).
TSP	Abbreviation of 3-trimethylsilylpropionate, a CSI supported by Chenomx NMR Suite as an alternative to DSS. The chemical shift of TSP is sensitive to pH.

U

uncoupled split	In Spin Simulator, a J-modifier that splits a spin definition as though it were coupled to a number of protons equivalent to the spin count, by an amount equal to the specified J-value (in Hz).
universally unique identifier	Standard identifier commonly used in software to identify a piece of information (like a compound signature) in such a way that the same identifier is never used for another piece of information.

unprocessed spectrum data Spectrum data stored as it was recorded by the spectrometer, with no modifications.

UUID *See universally unique identifier..*

V

valid pH Optional parameter describing the pH at which a compound signature is expected to accurately reflect an experimental spectrum. Often correlates with the pH range used to determine transform windows.

vertical zoom Adjusts the scale of the Spectrum View in the vertical direction, without changing the displayed frequency range.

VNMR Older NMR software package developed by Varian to run their spectrometers and process the resulting data.

VNMRj NMR software package developed by Varian to run their spectrometers and process the resulting data.

W

water deletion Processing layer that masks data points in a spectrum in the frequency domain across a range associated with water. Do not confuse water deletion with water subtraction methods that operate in the time domain to remove the water signal itself.

water suppression NMR technique that reduces the contribution of the water signal to ^1H NMR spectra of aqueous solutions via a special pulse sequences.

X

.xcpd *See compound file..*

.xss *See simulation file..*

XWIN-NMR NMR software package developed by Bruker to run their spectrometers and process the resulting data.

Z

zero-filling Refers to adding zeros to a fid prior to Fourier-transforming it. Used by modern spectrometers to ensure that all spectra have a number of points that is a power of 2, enabling the use of fast Fourier-transform algorithms. Can also increase the digital resolution of a spectrum, to a certain extent.

Index

A

- ACD/Labs, 5
- acquisition time, 5
- activation
 - key, 2
 - license, 2
 - license file (.cnxlic), 2
 - multiple users, 2
 - software, 2
- add
 - Chenomx Reference Compounds, 3, 97
 - compound file (.xcpd), 97
 - compound pack (.pack), 97
 - J-modifier, 73, 75
 - old library file (.cnxlib), 97
 - peak, 81
 - selected compounds to [Compound Set], 101
 - spectrum file (.cnx), 97
 - spin definition, 71
 - spin system, 70
 - ZIP-compressed file (.zip), 97
- animated zoom, 25
- annotations
 - custom color, 46, 48
 - pin, 46, 49
 - star, 46, 50
- application log, vii
- auto zoom, 21
- automatic
 - baseline correction, 37
 - fit (single compound), 55
 - phasing, 36

B

- baseline correction
 - automatic, 37
 - breakpoint, 37
 - 'frowns', 37
 - 'rolling', 37
 - 'smiles', 37
 - 'wings', 37
- best guess fit, 121
- bin
 - label, 62
 - size, 62
 - width, 62
- bin definitions, 62
- binning configuration file, 62
- binning target, 62

- breakpoint, 37
- Bruker
 - processed data (1r), 13
 - spectrum file (fid), 13

C

- change columns, 16
- chemical shift, 4
- Chenomx file types
 - compound album file (.cca), 80
 - compound file (.xcpd), 13, 80
 - compound pack (.pack), 61
 - license file (.cnxlic), 2
 - simulation file (.xss), 13, 70, 80
 - spectrum file (.cnx), 13, 41
- Chenomx Reference Compounds
 - add, 3, 97
 - Reference Compound Sets, 59
 - transform windows, 55
- Chenomx Technical Support, vii
- clear
 - concentration, 52
 - transform, 53
 - transform window, 86
- cluster
 - generate for region, 86
 - group peaks as, 84
 - list, 92
 - navigator, 58, 90
 - select, 82
- cluster IDs
 - manage, 90
- comment, 34
- compound
 - add to [Compound Set], 101
 - album file (.cca), 80
 - author, 91
 - details, 91
 - file (.xcpd), 13, 80
 - filter in selected region, 50
 - filter near x.xx ppm, 50
 - ID, 91
 - magnet frequency, 91
 - name, 91
 - new, 80
 - pack (.pack), 61
 - pH sensitive, 88
 - remove, 99
 - remove from compound set, 101
 - sample concentration, 91
 - set, 59, 96, 100, 101, 102, 103
 - table, 46, 96
 - UUID, 91

- compound line fit, 117
- compound set
 - edit, 102
 - export compound files (.xcpd), 98
 - export compound pack (.pack), 98
 - new, 101, 102
 - remove, 103
 - remove compound from, 101
 - rename, 103
 - sidebar, 18, 59, 100
 - smart, 100, 102, 103
- concentration
 - clear, 52
 - export profiled data, 60
 - exported data (.txt), 31
 - maximum, 56
 - scale, 57
 - set to maximum, 56
 - units, 66
- convert
 - automatic, 41
 - multiple spectra, 41
 - settings, 41
 - spectra, 41
- copy spectrum image to clipboard, 24
- coupled split, 73
- creatinine, 4, 8, 32
- CSI
 - DSS, 4, 8, 31
 - editor, 31
 - find automatically, 31
 - formate, 4, 8, 31
 - internal standard, 4
 - preserve area, 31
 - TSP, 4, 8, 31
- custom color, 17, 46, 48

D

- data acquisition, 5
 - NMR parameters, 5
 - (see also NMR parameters)
 - pulse sequence, 5
- datdir, 6
- default
 - file folder, 25
 - transform window, 86
- default layout (see export)
- delete
 - J-modifier, 73, 75
 - peak, 83
 - spin definition, 71
 - spin system, 70
- Delta, 5

DFTMP, 4, 8, 32

display options

sidebar, 24

spectrum thumbnail, 24

subtraction line, 118

sum line, 119

X scale, 24

Y scale, 24

DSS (see CSI)

DSS/TSP satellites, 39

E

edit

J-modifier, 73, 75

peak, 83

smart compound set, 102

spin definition, 71

transform window, 86

encapsulated PostScript (.eps), 23

enforce transform windows, 53, 55

exclude region, 62

export

colors, 60

compound files (.xcpd), 98

compound pack (.pack), 98

compound set, 98

concentrations, 60

data, 31

default layout, 60

image options, 23

matched clusters, 60

pins, 60

profiled compounds, 61

profiled data, 60

spectrum image, 23

stars, 60

transform summaries, 60

transposed layout, 60

exported data (.txt), 31

expressions, 83, 86

F

file operations

close, 15

open, 13

save, 14

filter

compounds, 50

compounds in selected region, 50

compounds near x.xx ppm, 50

find CSI automatically, 31

find pH automatically, 32

first-order phase, 36

fit automatically, 55

folder

default, 25

formate (see CSI)

frequency domain size, 34

'frowns', 37

G

generate cluster for region, 86

graphics interchange format (.gif), 23

group

peaks as cluster, 84

selected splits, 73

H

horizontal zoom, 21

I

identify compound

based on a peak, 9

based on a region, 9

image file types

encapsulated PostScript (.eps), 23

graphics interchange format (.gif), 23

portable document format (.pdf), 23

portable network graphic (.png), 23

scalable vector graphic (.svg), 23

imidazole, 4, 8, 32

import

compound album file (.cca), 80

profile, 62

simulation file (.xss), 80

spectrum, 13

information panel, 82, 83

initial

zoom, 25

installation, 1

Linux, 1

Mac OS X, 1

Windows, 1

internal standard, 4

J

J-coupling, 75

J-modifier

add, 73, 75

coupled split, 73

delete, 73, 75

edit, 73, 75

group selected, 73

J-value, 73

J-value editor, 73

list, 73

split count, 73

uncoupled split, 73

ungroup selected, 73

J-value, 73

J-value editor, 73

Java, 1

JCAMP-DX spectrum file (.jdx), 13

JEOL Delta spectrum file (.jdf), 13

JRE (see Java)

K

keyboard shortcut

cluster adjustment, 84

cluster selection, 93

compound table, 47

concentration, 52

CSI peak, 31

display next cluster, 58, 84

display previous cluster, 58, 84

horizontal scroll, 21

peak adjustment, 84

peak selection, 93

phasing, 36

select by cluster, 82

select by peak, 82

select region, 22

spectrum pH, 33

spin definition, 72

transform, 54

zoom, 21

L

legend

sidebar, 17, 64

library manager

open, 59

license file (.cnxlic), 2

line broadening, 38

line color, 25

M

magnet frequency, 91

manage cluster IDs, 90

maximum concentration, 56

N

new

compound, 80

compound set, 101

simulation, 70

smart compound set, 102

next spin definition, 72

NMR parameters

- acquisition time, 5
- relaxation delay, 5
- spectral width, 5
- temperature, 5

NMRPipe

- processed data (.ft2), 13
- spectrum file (.fid), 13

no fit, 121

normalization, 62

O

optimize selected peak shapes, 85

output options, 62

overlay spectrum, 19

overview

- Compound Builder, 77
- Library Manager, 94
- Processor, 28
- Profiler, 43
- Spin Simulator, 67

P

peak

- add, 81
- delete, 83
- edit, 83
- group as cluster, 84
- list, 92
- optimize selected, 85
- select, 82

pH

- creatinine, 4, 8, 32
- DFTMP, 4, 8, 32
- editor, 32
- find automatically, 32
- imidazole, 4, 8, 32
- indicator, 4
- internal standard, 4
- sensitivities, 88
- spectrum, 32
- uncertainty, 8

pH sensitivities, 88

phasing

- automatic, 36
- first-order, 36
- zero-order, 36

pin, 46, 49

portable document format (.pdf), 23

portable network graphic (.png), 23

preprocessing, 34, 35

- ACD/Labs, 5

Delta, 5

Varian processed data (phasefile), 6

VNMR, 5

VNMRj, 5

XWIN-NMR, 5

preserve CSI area, 31

previous spin definition, 72

processed data

- Bruker (1r), 13
- JCAMP-DX (.jdx), 13
- JEOL Delta (.jdf), 13
- NMRPipe (.ft2), 13
- Varian (phasefile), 13

processing

- baseline correction, 37
- history, 34
- history (sidebar), 18, 34
- line broadening, 38
- overview, 8
- phasing, 36
- preprocessing, 34, 35
- reference deconvolution, 39
- reverse spectrum, 41
- supported spectra, 5

(see also supported file types)

water deletion, 40

profiled compounds, 46, 48, 49, 50, 52, 53, 57, 61

profiled data, 48, 49, 50, 52, 53, 60

profiling

- best guess fit, 121
- compound line fit, 117
- determining concentrations, 11
- identifying compounds, 9
- no fit, 121
- overlapped regions, 11
- overview, 8
- subtraction line fit, 118
- sum line fit, 119
- water suppression, 120

proline, 125, 131, 134

proxy server, 25

Q

quick reference guide

- Compound Builder, 78
- Library Manager, 95
- Processor, 29
- Profiler, 44
- Spin Simulator, 68

quick search, 50, 96

R

redo zoom, 20

Reference Compound Sets, 3, 59

reference deconvolution

- and line broadening, 39
- DSS/TSP satellites, 39
- target linewidth, 39

reference panel

- Compound Builder, 91
- Profiler, 65
- sidebar, 18, 65, 91

relaxation delay, 5

remove

- compound, 99
- compound from compound set, 101
- compound set, 103
- spectrum overlays, 19

rename compound set, 103

reverse spectrum, 41

'rolling', 37

S

sample concentration, 91

sample preparation, 4

scalable vector graphic (.svg), 23

scale concentrations, 57

select

- cluster, 82
- peak, 82
- region, 22

selected region, 22

send to

- Compound Builder, 99
- Processor, 64
- Profiler, 42

set

- to maximum concentration, 56
- zoom, 20

show entire spectrum, 20

sidebar, 16, 24

- Compound Builder, 91
- compound sets, 18, 59, 100
- legend, 17, 64
- processing history, 18, 34
- Profiler, 64
- reference panel, 18, 65, 91
- simulation, 19, 75
- Simulator, 75

simulation

- details, 76
- file (.xss), 13, 70, 80
- frequency, 76

- new, 70
- sidebar, 19, 75
- spin definition, 71
- spin system, 70, 71
- smart compound set, 100, 102, 103
- 'smiles', 37
- spectral binning
 - bin definitions, 62
 - bin label, 62
 - bin size, 62
 - bin width, 62
 - binning configuration file, 62
 - binning target, 62
 - exclude region, 62
 - normalization, 62
 - output options, 62
 - total area, 62
- spectral width, 5
- spectrum
 - Bruker (fid), 13
 - file (.cnx), 13, 41
 - JCAMP-DX (.jdx), 13
 - JEOL Delta (.jdf), 13
 - NMRPipe (.fid), 13
 - overlay, 19
 - pH, 32
 - remove overlays, 19
 - reverse, 41
 - thumbnail, 22
 - Varian (fid), 13
- spectrum image
 - copy to clipboard, 24
 - export, 23
- spectrum thumbnail, 22, 24
- Spectrum View, 21, 58, 73
- spin definition
 - add, 71
 - delete, 71
 - edit, 71
 - navigator, 73
 - panel, 70
- spin system
 - add, 70
 - delete, 70
- split count, 73
- star, 46, 50
- statistical analysis
 - export profiled data, 60
 - spectral binning, 62
- subtraction line fit, 118
- sum line fit, 119
- supported file types
 - Bruker processed data (1r), 13

- Bruker spectrum file (fid), 13
- JCAMP-DX spectrum file (.jdx), 13
- JEOL Delta spectrum file (.jdf), 13
- NMRPipe processed data (.ft2), 13
- NMRPipe spectrum file (.fid), 13
- Varian processed data (phasefile), 6, 13
- Varian spectrum file (fid), 13

T

- tape measure (see selected region)
- target linewidth, 39
- temperature, 5
- time domain size, 34
- total area, 31, 62
- transform window
 - clear, 86, 88
 - default, 86
 - edit, 86
 - enforce, 53, 55
- transposed layout (see export)
- TSP (see CSI)
- tutorial compounds
 - proline, 125, 131, 134
 - valine, 122, 129, 134

U

- uncoupled split, 73
- undo zoom, 20
- ungroup selected splits, 73
- UUID, 91

V

- valine, 122, 129, 134
- Varian
 - processed data (phasefile), 6, 13
 - spectrum file (fid), 13
- vertical zoom, 21
- VNMR, 5
- VNMRj, 5

W

- water
 - deletion, 40
 - suppression, 120
- 'wings', 37

X

- X scale, 24
- XWIN-NMR, 5

Y

- Y scale, 24

Z

- zero-order phase, 36
- zoom
 - animated, 25
 - auto, 21
 - box, 21, 22
 - horizontal, 21
 - initial, 25
 - redo, 20
 - set, 20
 - show entire spectrum, 20
 - undo, 20
 - vertical, 21