


WARNING

1. Do not use sling if a patient needs head or neck support.
2. Do not use sling if a patient suffers from back injury.
3. Do not use sling if a patient can not weight bear.
4. Please understand hoist operating instructions before using.
5. Check sling before use for signs of wear.

Standing Sling Lifting Procedures

KH400SFB 40kg—110kg
KH400SFLB 100kg—150kg

STEP 1: Undo waist strap (A) slide the body of the sling down the patients back until the base of the sling is level with the base of the spine. Fasten the velcro waist strap (A) snugly around the waist and clip buckle. Adjust strap on buckle (C) to ensure a snug fit.



STEP 2: Move the lifter into position directly in front of the patient. Ensure lifting arm is at lowest point. Legs of the lifter are gated to clear the legs of the chair and castor brakes are released (not in braking position). Patients feet are placed comfortably on the footplate and patients knees/shins are placed against the knee pad.



STEP 3: When using the KH400SFB or KH400SFLB use the coloured strap that fits comfortably over the attachment knobs. There are 2 recommended mounting points relevant to the patients height. The front knob is used for patients over 5'6" (165cm) and the back knobs are used for patients below 5'6" (165cm). Ensure both sides have the same coloured straps fitted

Front Knob

Rear Knob



CAUTION: Patients can use the hand rails to steady themselves. The operator/carer can also steady the patient by holding the side of the sling.

STEP 4: Start the lifting process using the slow speed. Always ensure the patient is comfortable, if at anytime the patient looks uncomfortable or unbalanced, stop immediately. Lower the patient and adjust the sling. Encourage the patient to lean slightly back and not forward, this will assist in the lifting process by not allowing the sling to ride up under the arms of the patient. Lift the patient to the desired position.



STEP 5: Ensure the patient is comfortable and then move to the desired position



STEP 6: Position the patient over the chair, bed or toilet and begin lowering. Always make sure the patient is comfortable and feeling safe.

CARE INSTRUCTIONS

1. Sling Care

The velcro strap should always be closed when not in use. Protect from direct sunlight. Do not puncture with any sharp objects. Repair any tear immediately. Do not lift a patient if there are any sling defects. Washing and drying temperature must not exceed 60 degrees.

2. Spot Cleaning

The sling is machine washable. Treat spills and stains as soon as possible. Gently scrape away any soils or mop away liquid from surface of fabric with damp cloth. Do not scrub with stiff brush. Dry in shade, away from direct heat. Dry thoroughly before use.




A Hills Company

4 Halley Rd
Balcatta WA 6021
Tel: 08 9240 2944
Fax: 08 9240 1469
Email: info@kerryequipment.com.au
Website: www.kerryequipment.com.au

PATIENT RECORDS FOR MOUNTING POINTS USING KH400SFB

KNOB 1: For patients over 5'6" (165cm)

KNOB 2: For patients below 5'6" (165cm)

[illegible]

USERS GUIDE

KH400SFB

STANDING SLING



FEATURES

Neoprene waist strap.
Padded for comfort.
Colour coded attachment straps.
Machine washable.



A  Hills Company