

# **GeneFishing™**

## **DEG Premix Kit**

### **User Manual**

- Version 5 (Japan Only)
- Published February 2005
- Catalog No. K1001, K1002  
K1021 ~ 1026

---

# Table of Contents

---

<b>1. Notices to Customers</b>	
1.1 Product Warranty and Liability-----	1
1.2 Safety Warning and Precautions-----	1
1.3 Notice to Customers-----	1
<b>2. Introduction</b>	
2.1 ACP™ Technology-----	2
2.2 GeneFishing™ DEG Premix kits	
2.2.1 Overview-----	4
2.2.2 Principle-----	5
2.2.3 Key Features-----	8
<b>3. Components</b>	
3.1 List of Components -----	10
3.2 Storage Conditions -----	12
3.3 Reagents and Equipments to be Supplied by User -----	12
<b>4. Protocol</b>	
4.1 Positive Control Experiments-----	13
4.2 Protocol for GeneFishing™ DEG Premix Kit	
4.2.1 Reverse Transcription-----	14
4.2.2 GeneFishing™ PCR-----	15
<b>5. Expected Results -----</b>	<b>17</b>
<b>6. Troubleshooting Guide -----</b>	<b>19</b>
<b>7. Ordering Information -----</b>	<b>21</b>

---

# 1 Notices to Customers

---

## 1.1 Product Warranty and Liability

Seegene guarantees the performance of all products as described when they are used according to the instructions in this manual. Any problem that occurs for reasons other than misuse should be reported to Seegene immediately. This warranty limits our liability for product replacement.

## 1.2 Safety Warning and Precautions

This product is limited to research use only. It is not recommended or intended for the diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

## 1.3 Notice to Customers

The PCR process is covered by patents owned by Hoffman-La Roche Inc. No license or immunity under any other patent is either expressed or implied by the sale of any Seegene product.

---


# 2 Introduction

---

## 2.1 ACP™ Technology

The specificity with which a primer anneals to its target sequences is the most critical factor for successful target-specific PCR amplification. **The principle of our ACP™ Technology (patent pending) focuses on an oligonucleotide primer that anneals with exquisite specificity to the intended template, therefore allowing only the target products to be amplified.** The primer is denoted **Annealing Control Primer (ACP™)**, and it has a unique tripartite structure. As shown in Figure 1, it has distinct 3'- and 5'-end portions separated by a regulator, which interact with the template in turn during a two-stage PCR.

ACP (Annealing Control Primer)



Designation	Function
a. Core sequence	Anneal at the 1 <sup>st</sup> PCR Stage (targeting)
b. Regulator	Universal sequence
c. Universal sequence	Anneal at the 2 <sup>nd</sup> PCR Stage

Fig. 1. ACP™ Structure

---

The 3'-end core, which is designated as 'a' in Figure 1, is the targeting portion that consists of a hybridizing sequence that is substantially complementary to a site on a nucleic acid template. The 5'-end portion (designated as 'b' in Figure 1) has a universal sequence. The regulator bridging the core and universal sequences of the ACP (designated as 'c' in Figure 1) plays a key role in controlling the annealing of each portion to the template.

The ACP system requires a two-stage PCR (ACP<sup>TM</sup>-based PCR) amplification to maximize the functions of each portion as follows:

**1<sup>st</sup> stage PCR to generate a specific PCR product**

During the 1<sup>st</sup> stage PCR, the conditions are set such that the 3'-end core portion of the ACP anneals to a specific site on the template. The annealing of the regulator to the template is not favored under these conditions. Consequently, the regulator interrupts the annealing of the 5'-end portion to the template, thereby restricting the primer-annealing portion to the 3'-end. The reaction equilibrium in the 1<sup>st</sup> stage PCR lies more favorably towards the specific annealing of the 3'-end portion sequence than the non-specific annealing under this annealing temperature. **Thus, the effect of the regulator on the 5'- and 3'-end portions of the ACP is one of the key features of the ACP<sup>TM</sup> Technology.**

**2<sup>nd</sup> stage PCR to amplify the 1<sup>st</sup> PCR product**

The product generated by the annealing and extension of the 3'-end portion sequence of the ACP in the 1<sup>st</sup> stage PCR is amplified in the 2<sup>nd</sup> stage PCR, which is conducted under high stringency conditions. This prevents further annealing of the 3'-end core portion of the ACP to the original template. **Instead, only the sequences at the 3'- and 5'-ends of the 1<sup>st</sup> PCR product act as donors of priming sequences for the amplification.** This results in the amplification of only the target product, with an amplification efficiency that is close to the theoretical optimum of a two-fold increase of product for each PCR cycle.

## 2.2 GeneFishing™ DEG Premix Kit

### 2.2.1 Overview




Fig. 2. Overview of the GeneFishing™ DEG Premix Kits

---

## 2.2.2 Principle

The GeneFishing™ DEG Premix Kits are designed to identify differentially expressed genes (DEGs) in two or more nucleic acid samples. The designing of the kit specifically focused on overcoming the disadvantages and limitations of current gene expression profiling-related methodologies, such as the microarray and differential display techniques. The GeneFishing™ DEG Premix Kits require the following three steps, which consist of reverse transcription (RT) and a two-stage PCR (GeneFishing™ PCR, see Figures 2 and 3, pages 4 and 7, respectively):


**Step 1:** RT is conducted using dT-ACP1 to **synthesize the first-strand cDNAs from the samples**. The 3'-end core portion of **dT-ACP1** consists of a hybridizing sequence that is complementary to the poly A region of mRNA transcripts. This results in the first-strand cDNA bearing the universal sequence of dT-ACP1 at its 5'-end.

**Step 2:** The first-strand cDNA resulting from the first step is diluted and placed in a PCR tube with an arbitrary ACP and dT-ACP2. The 3'-end core portion (10-mer) of the arbitrary ACP consists of a hybridizing sequence that is sufficiently complementary to a region on the first-strand cDNAs. A single PCR cycle (1<sup>st</sup> stage PCR) is conducted under conditions that permit only the arbitrary ACP to anneal with its 3'-end core to the first-strand PCR. These conditions do not allow the 3'-end core portion of dT-ACP2 to anneal to the first-strand cDNA. The result of the 1<sup>st</sup> stage PCR is the synthesis of second-strand cDNA that bears the complementary sequence of the universal sequence of dT-ACP1 on its 3'-end and the universal sequence of the arbitrary ACP on its 5'-end. **This step ensures that no artifacts and dT-ACP2/dT-ACP2 products are obtained.**

---

**Step 3:** The mixtures in the PCR tubes arising from step 2 are then subject to 2<sup>nd</sup> stage PCR under high stringency conditions that allow **both** dT-ACP2 and the arbitrary ACP to respectively anneal to the 3'- and 5'-ends of the arbitrary ACP-primed second-strand cDNA generated in step 2. These annealing events exclusively involve the universal sequences of the two primers since the 3'-end of the second-strand cDNA is complementary to the universal sequence of the arbitrary ACP while the 5'-end is recognized by the universal sequence of dT-ACP2. This results in the amplification of the targeted PCR products **ONLY**.






### First-strand cDNA synthesis

- ◆ First-strand cDNA is synthesized by RT using dT-ACP 1 (1-2).




## GeneFishing™ PCR



### First-stage PCR for second-strand cDNA synthesis

- ◆ Annealing is allowed at  $T_{\text{anneal}}=50^{\circ}\text{C}$ , which is high enough to prevent the core of dT-ACP2 from binding(3a).
- ◆ Only 10-mer core of arbitrary ACP is able to bind to the first-strand cDNA with 8-10 base pair matches (3b).
- ◆ Second-strand cDNA is synthesized in only one cycle (4).



### Second-stage PCR for amplifying the second-strand cDNA

- ◆ 35~40 cycles of PCR are performed at  $T_{\text{anneal}}=65^{\circ}\text{C}$ , which is stringent enough to prevent core sequences from annealing yet still permits the binding of universal sequences to each other.
- ◆ Arbitrary ACP and dT-ACP2 cannot anneal on the first-strand cDNA at  $65^{\circ}\text{C}$  (5a).
- ◆ The second-strand cDNA is amplified by priming the Arbitrary ACP and dT-ACP2 at 3'- and 5'-ends, respectively (5b).
- ◆ Only target PCR products are amplified (6).

Electrophoresis of PCR products on an agarose gel

  
NO false positives!

Fig 3. Flow chart of cDNA synthesis and the GeneFishing™ PCR

---

### **2.2.3 Key Features**

#### **(1) No false positives**

GeneFishing™ DEG Premix Kits identify Differentially Expressed Genes (DEGs) and never fail upon confirmation by Northern blot analysis (or RT-PCR)! The problem of false positives has been a major remaining bottleneck for current DEG-profiling methodologies, such as the microarray and differential display technologies. The lack of false positives allows the researcher to concentrate on authentic DEGs.

#### **(2) No PAGE required, as an agarose gel is sufficient**

The Annealing Control Primer (ACP) dramatically improves the specificity and sensitivity of PCR amplification and results in only specific PCR products. In addition, the GeneFishing™ Technology requires small amounts of starting material. The PCR products can be detected on standard ethidium bromide-stained agarose gels, and hence, obviating the use of PAGE gels with the attendant issues in handling and use. The bands shown on agarose gels by GeneFishing™ Technology are of adequate resolution to be detected by Northern blot analysis.

#### **(3) No expensive detection methods**

The use of Radioactive/fluorescent detection methods is limited in terms of its cost and risks. Such expensive detection methods are additional drawbacks of current gene expression profiling methods. Therefore, the detection of PCR products by EtBr staining of agarose gels is a key benefit of the GeneFishing™ DEG Premix Kit.

#### **(4) Guaranteed reproducibility**

The GeneFishing™ DEG Premix Kits provide highly reproducible PCR products time after time as all reagents are provided with the kit. This ensures that variations caused by the use of old, incompatible reagents are prevented.

---

**(5) No special skills required**

GeneFishing™ Technology is a robust PCR based technology and as such is simple to use.

**(6) Speedy and cost-effective**

The GeneFishing™ DEG Premix Kits enable researchers to identify authentic DEGs within 5 hrs and not to waste time on accessing false positives. In contrast, all other current DEG profiling methods require intensive downstream work and time to identify the authentic DEG candidates.

**(7) Wide range of PCR products**

Each GeneFishing™ PCR reaction generates a wide range of PCR products ranging from 150 bp to 2 kb, which not only increases the chances of identifying DEGs and but also provides more significant sequence information for the prediction of gene function.

---

# 3 Components

---

## 3.1 List of Components

- **dT- ACP1 (10  $\mu$ M):** For Reverse Transcription, 1 vial  
dT-ACP1: 5'-CTGTGAATGCTGCGACTACGATXXXXX(T)<sub>18</sub>-3'
  
- **dT- ACP2 (10  $\mu$ M):** For GeneFishing™ PCR, 1 vial  
dT-ACP2: 5'-CTGTGAATGCTGCGACTACGATXXXXX(T)<sub>15</sub>-3'
  
- **Arbitrary ACPs (5  $\mu$ M):** For GeneFishing™ PCR

<u>Cat. No.</u>	<u>Arbitrary ACPs</u>				<u>vial</u>	
K1001 :	Arbitrary ACP	1	-	Arbitrary ACP	10	10
K1002 :	Arbitrary ACP	11	-	Arbitrary ACP	20	10
K1021 :	Arbitrary ACP	1	-	Arbitrary ACP	20	20
K1022 :	Arbitrary ACP	21	-	Arbitrary ACP	40	20
K1023 :	Arbitrary ACP	41	-	Arbitrary ACP	60	20
K1024 :	Arbitrary ACP	61	-	Arbitrary ACP	80	20
K1025 :	Arbitrary ACP	81	-	Arbitrary ACP	100	20
K1026 :	Arbitrary ACP	101	-	Arbitrary ACP	120	20

**Note**

**Note:** The 3'-end core portion of each arbitrary ACP consists of a 10-mer of randomly selected nucleotides, while the 5'-end consists of the universal sequence. For example, the 20 arbitrary ACPs enclosed in the K1021 have the following sequences:

---

ACP1 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXGCCATCGACC-3'  
ACP2 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXAGGCGATGCC-3'  
ACP3 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXCCGGAGGATG-3'  
ACP4 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXGCTGCTCGCG-3'  
ACP5 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXAGTGCGCTCG-3'  
ACP6 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXGGCCACATCG-3'  
ACP7 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXCTGCGGATCG-3'  
ACP8 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXGGTCACGGAG-3'  
ACP9 : 5'-GTCTACCAGGCATTCGCTTCATXXXXXGATGCCGCTG-3'  
ACP10: 5'-GTCTACCAGGCATTCGCTTCATXXXXXTGGTCGTGCC-3'

ACP11: 5'-GTCTACCAGGCATTCGCTTCATXXXXXCTGCAGGACC-3'  
ACP12: 5'-GTCTACCAGGCATTCGCTTCATXXXXXACCGTGGACG-3'  
ACP13: 5'-GTCTACCAGGCATTCGCTTCATXXXXXGCTTCACCGC-3'  
ACP14: 5'-GTCTACCAGGCATTCGCTTCATXXXXXGCAAGTCGGC-3'  
ACP15: 5'-GTCTACCAGGCATTCGCTTCATXXXXXCCACCGTGTG-3'  
ACP16: 5'-GTCTACCAGGCATTCGCTTCATXXXXXGTCGACGGTG-3'  
ACP17: 5'-GTCTACCAGGCATTCGCTTCATXXXXXCAAGCCCACG-3'  
ACP18: 5'-GTCTACCAGGCATTCGCTTCATXXXXXCGGAGCATCC-3'  
ACP19: 5'-GTCTACCAGGCATTCGCTTCATXXXXXCTCTGCGAGC-3'  
ACP20: 5'-GTCTACCAGGCATTCGCTTCATXXXXXGACGTTGGCG-3'

- **Control cDNAs (Kidney, 10 ng/μl), 1 vial**  
These control cDNAs were synthesized from total RNAs isolated from adult mouse kidney tissues.
- **Control cDNAs (Liver, 10 ng/μl), 1 vial**  
These control cDNAs were synthesized from total RNAs isolated from adult mouse liver tissues.
- **Control ACP (5 μM):** For GeneFishing™ PCR, 1 vial

**Note**

- **2X SeeAmp™ ACP™ Master Mix (280 rxns), 4 vials**  
**Note:** K1001 and K1002 have 2 vials (140 rxns)

---

### 3.2 Storage Conditions

**Store all reagents below -20°C except the SeeAmp™ ACP™ Master Mix.**

In case of the SeeAmp™ ACP™ master mix, 4°C for short-term storage and -20°C for long-term storage are recommended. Avoid repetitive thawing as it may decrease the activity of the master mix. The master mix is proven to be stable for at least 3 months at 4°C.

### 3.3 Reagents and Equipments to be Supplied by the User

RNase-free H <sub>2</sub> O	Reverse transcriptase	2 mM dNTP
RNase inhibitor	Thermal cycler	Micro-centrifuge

---

# 4 Protocol

---

## 4.1 Positive Control Experiments

◆ Please consider the following in prior to your experiments.




If you are using GeneFishing™ DEG Premix kits for the first time, please establish the GeneFishing™ PCR conditions with the positive control experiments.

The positive control experiments should be conducted using the control cDNAs as templates and a primer set of Control ACP and dT-ACP2 (provided in the kit) according to the instructions in this protocol.

You can assume that the initial experiment conditions are set up when the pattern of the positive control displayed on the agarose gel is similar with the figure on the right.

*Refer the instructions on page 15.*

3 µl	Control cDNA-Kidney or Liver
2 µl	Control ACP
1 µl	dT-ACP2 (10 uM)
4 µl	Distilled water
10 µl	2X SeeAmp™ ACP™ Master Mix
<hr/>	
20 µl	Total volume



---

## 4.2 Protocol for GeneFishing™ DEG Premix Kits

### 4.2.1 Reverse Transcription

1. Add the following reagents to a RT tube on ice.

? $\mu$ l	Total RNA (3 $\mu$ g)
2 $\mu$ l	10 $\mu$ M dT-ACP1
? $\mu$ l	RNase-free water
<hr/>	
9.5 $\mu$ l	

**Note**

**Note:** Mix the reagents by tapping or pipetting.

**Note:** In order to identify differentially expressed bands between RNA samples, it is important to add equal amount of RNA to each tube.

2. Incubate the tube at 80°C for 3 min.
3. Chill the tube on ice for 2 min and spin the tube briefly.
4. Add the following reagents to the tube from step 3.

4 $\mu$ l	5X RT buffer
5 $\mu$ l	2 mM dNTP
0.5 $\mu$ l	RNase inhibitor (40 u/ $\mu$ l)
1 $\mu$ l	M-MLV reverse transcriptase (200 u/ $\mu$ l)
<hr/>	
20 $\mu$ l	Total volume

5. Incubate the tube at 42°C for 90 min.
6. Heat the tube at 94°C for 2 min.
7. Chill the tube on ice for 2 min and spin the tube briefly.
8. Dilute the first-strand cDNA by adding 80  $\mu$ l of DNase-free water.

**Note**

**Note:** Store all cDNA samples at -20°C until ready for use.



## 4.2.2 GeneFishing™ PCR

1. Add the following reagents to a PCR tube **on ice**.

3~7 $\mu$ l	Diluted first-strand cDNA (~50 ng)
2 $\mu$ l	5 $\mu$ M arbitrary ACP (one of the arbitrary ACPs)
1 $\mu$ l	10 $\mu$ M dT-ACP2
? $\mu$ l	Distilled water
10 $\mu$ l	2X SeeAmp™ ACP™ Master Mix
<hr/>	
20 $\mu$ l	Total volume

### Note

**Note:** Depending on the samples, different amounts of diluted first-strand cDNA can be used as templates for GeneFishing™ PCR. High amounts of the starting material (the first-strand cDNA) results in perfect reproducibility and amplification of rare mRNAs. It also permits the use of the ethidium bromide-stained agarose gel to detect differentially expressed products. We recommend using 3~5  $\mu$ l of the diluted first-strand cDNA as the template for the GeneFishing™ PCR.

### Note

**Note:** A primer set consisting of one of the ten arbitrary ACPs in the kit and dT- ACP2 should be used as the 5' and 3' primers in each PCR reaction.



**Important !!**

2. Place the tube in a preheated (94°C) thermal cycler.

**Note: It is important to preheat (94°C) the thermal cycler before placing the tube in the thermal cycler.**

3. Immediately commence the PCR reaction using the following program.

Segment	No. of cycles	Temperature	Duration
1	1	94°C	5 min
2	1	50°C	3 min
3	1	72°C	1 min
4	40	94°C	40 sec
		65°C	40 sec
		72°C	40 sec
5	1	72°C	5 min

### Note

**Note:** We recommend the GeneAmp PCR System 9700 of Applied Biosystems that has a heated lid.

---

4. Electrophorese 2~3  $\mu$ l of the PCR products on a 2% agarose gel containing EtBr.

**Note**

**Note:** If the band intensity of the DEGs of your sample is very weak, the intensity can be increased by raising the amount of starting material (diluted first-strand cDNA) or the ACP primer concentration in the GeneFishing™ PCR.

5. Extract the differentially expressed bands from the agarose gel.

**Note**

**Note:** The gel extraction kit (e.g. QIAquick Gel extraction kit (Qiagen, cat. No. 28704) or the GENECLEAN II kit (Q.BIOgene, cat. No. 1001-400)) is recommended to extract the differentially expressed bands from the agarose gel.

6. Clone the product into a TA cloning vector

**Note**

**Note:** We suggest that you avoid re-amplifying the extracted PCR product, because using the same primer set in re-amplification may generate smearing or non-specific products.” Even if the differentially expressed band of interest is faint, the extracted DNA should be directly used for cloning into a TA cloning vector. In this case, we strongly recommend to use an efficient TA cloning system (e.g., the Invitrogen, TOPO TA Cloning Kit).


# 5

## Expected Results


### Examples of GeneFishing™ DEG Premix Kits Applications

Various experiments to compare differentially expressed genes (DEGs) in two or more RNA samples from different organisms were conducted with the ACPs of the GeneFishing™ Premix DEG Kits in accordance with the instructions given in this User Manual. Below are examples of resulting data shown by agarose gel photographs. These agarose gels show the typical results that are generated by such experiments. Note that the number of PCR products can vary depending on the type of sample loaded.


#### A. Mouse Conceptus Tissues (E4.5 vs. E11.5 vs. E18.5)



#### B. Human Brain Tissues (Normal vs. Tumor)



### C. Human Stomach Tissues (Normal vs. Tumor)



**Fig. 4** Agarose gel photographs indicating various DEGs obtained by using GeneFishing™ Technology. **A.** An agarose gel photograph showing DEGs obtained from E4.5, E11.5 and E18.5 mouse embryos by using five arbitrary ACPs. **B.** An agarose gel photograph showing DEGs obtained from normal human brain and astrocytic tumor tissues by using four arbitrary ACPs. **C.** An agarose gel photograph showing DEGs obtained from normal human stomach and stomach tumor tissues by using three arbitrary ACPs. The DEGs are marked by white arrowheads. M represents the molecular weight marker obtained by using the Forever 100 bp Ladder Personalizer Kit (Seegene, Cat. No. M0100).

# 6 Troubleshooting Guide

Problems	Comments and suggestions
<b>No band</b>	<p><b>You may have a problem with the RT or PCR reaction.</b></p> <p><b>Please follow the instructions:</b></p> <ol style="list-style-type: none"> <li>Confirm the integrity of the PCR reaction by using our control cDNAs. If several bands arise from the control cDNAs but not from your experimental samples, you may have a problem with the RT reaction and/or the RNA quality. Make sure that Seegene's dT-ACP1 was added to the RT reaction and check the integrity of the RNA samples by formaldehyde-agarose gel electrophoresis.</li> <li>If there is no problem with the RT, check the concentration of the first-strand cDNA: too little cDNA may cause no or poorly intense bands. If this is the case, add more cDNA.</li> </ol>
<b>When comparing samples, you see the same band pattern</b>	<p>Your RNA samples or PCR reagents may be contaminated.</p> <ol style="list-style-type: none"> <li>If your RNA samples are contaminated with chromosomal DNA, treat your RNA with DNase I.</li> <li>If PCR products have cross-contaminated your reagents, replace the reagents.</li> </ol>
<b>Fewer than usual bands in the sample lanes (while the control cDNA experiment works)</b>	<ol style="list-style-type: none"> <li>Check the integrity of your total RNA.</li> <li>The number of PCR products may fluctuate depending on the type of samples. More bands may be obtained by decreasing the annealing temperature at the 1<sup>st</sup> stage PCR, although it may increase the chance of obtaining false positives as well.</li> </ol>
<b>Bad resolution on agarose gel</b>	<p>We recommend using 2% agarose gel and running the gel (12 X 14 cm) until the bromophenol blue dye has migrated 6 cm from the well.</p>

<b>Problems</b>	<b>Comments</b>
<b>Smearing only in experimental samples</b>	<p>You may have problems with the RNA quality or the PCR reaction itself.</p> <ol style="list-style-type: none"> <li>a. Clean up your RNA samples by phenol:chloroform extraction or use fresh RNA samples in the RT reaction.</li> <li>b. Check the integrity of the RNA by formaldehyde-agarose gel electrophoresis.</li> <li>c. Confirm the integrity of the PCR reaction by using the control cDNAs.</li> <li>d. It is important to place the PCR reaction on ice immediately before samples are placed in the thermal cycler.</li> </ol>
<b>Minimum RNA quantity required</b>	<p>When we performed RT using mouse conceptus total RNA to determine the lower limit of total RNA needed, we found that reproducible results could be obtained when at least 250 ng of total RNA was used.</p>
<b>Cloning of weak DEG band</b>	<p>If you failed to clone the extracted PCR product due to its poor band intensity, we recommend you try the following procedures.</p> <ol style="list-style-type: none"> <li>a. Repeat the GeneFishing™ PCR 2~3 times employing the same arbitrary ACP/dT-ACP2 combination used in the original PCR. Extract the DEG band in question from each gel and combine. Directly clone the DEG band.</li> <li>b. Repeat the GeneFishing™ PCR using the arbitrary ACP/dT-ACP2 combination at a 2~3 times higher concentration. This may increase the band intensity, allowing the cloning of the DEG bands after gel extraction.</li> </ol>

---

# 7 Ordering Information

---

## GeneFishing™ DEG Premix Kits

<b>Cat. No.</b>	<b>Products</b>	<b>Size</b>
K1001	GeneFishing™ DEG101 Premix Kits	10 X 12 rxns
K1002	GeneFishing™ DEG102 Premix Kits	10 X 12 rxns
K1021	GeneFishing™ DEG101 & 102 Premix Kits	20 X 12 rxns
K1022	GeneFishing™ DEG103 & 104 Premix Kits	20 X 12 rxns
K1023	GeneFishing™ DEG105 & 106 Premix Kits	20 X 12 rxns
K1024	GeneFishing™ DEG107 & 108 Premix Kits	20 X 12 rxns
K1025	GeneFishing™ DEG109 & 110 Premix Kits	20 X 12 rxns
K1026	GeneFishing™ DEG111 & 112 Premix Kits	20 X 12 rxns
K1040	GeneFishing™ DEG101 & 112 Premix Kits	Full-package

<b>Cat. No.</b>	<b>Products</b>	<b>Description</b>	<b>Size</b>
P1131	dT-ACP1	Reverse Transcription	20 rxns
P1132	dT-ACP2	GeneFishing™ PCR	40 rxns
P1001-1120	Arbitrary ACPs	GeneFishing™ PCR	40 rxns
E1010	2X SeeAmp™ ACP™ Master Mix	2X Master Mix	5 ml





