




User's Manual For 1000-DR Digital Conversion for the PC-1000 Imaging Software

Software version 6.0



Panoramic Corporation

4321 Goshen Road • Fort Wayne, IN 46818-1242
800-654-2027 • www.pancorp.com

Document Information

Document Version 1.7

Copyright (c) 2008 Oy Ajat Ltd.

Oy Ajat Ltd. may at its sole discretion from time to time provide updates to this Manual or other supporting documentation or withdraw them at any time without prior notice.

This document is the user's guide for the software provided with the SCAN300FPx digital panoramic imaging system. This manual provides a documentation of the software's features and main concepts. Specific equipment usage and safety precautions are considered to be out of the scope of this document and such topics should be considered prior to usage. Any and all reproduction, transfer, distribution or storage of the contents of this document in part or in whole and in any form without the prior written permission of Oy Ajat Ltd. is prohibited.

Specifications of SCAN300FP

SCAN300FP"x" SERIES COMMON FEATURES:

Single crystal CdTe(CdZnTe) -CMOS linear sensor

Pixel size: 100um (hundred micro meters)

Six tiles of approximately 25mmx6.4mm each.

Approximately maximum active area: 150mmx6.4mm

Approximately Useful active area: 146mmx5.6mm

Physical gap between adjacent tiles: 2 lines (=200um)

Maximum interpolated gap between last used pixels on each tile: 4 lines

Responsible Author: Matt Powell, Oy Ajat Ltd.

Approved by: K. Spartiotis, Oy Ajat Ltd.



Oy AJAT Ltd. declares that SCAN300FPC,FPBX,FPA,FPSC, FPI, FPS, FPG and its components are in conformity with the essential requirements of the Medical Device Directive 93/42/EEC, Annex II for Class IIa medical device. The product is marked with the CE marking and the Notified Body number 0434.



DISCLAIMER

THE CONTENTS OF THIS DOCUMENT ARE PROVIDED "AS IS". EXCEPT WHERE REQUIRED BY APPLICABLE LAW, NO WARRANTIES OF ANY KIND, EITHER EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE ARE MADE IN RELATION TO THIS DOCUMENT.

Oy Ajat Ltd. shall under no circumstance be responsible for any loss of data, income or any incidental, consequential, special or indirect damages however caused or occurred.

The software or related imaging devices may contain technology subject to export laws and regulations from the US and/or other countries. Any export activities contrary to law are strictly prohibited.

Table of Contents

1	INSTALLATION	7
1.1	COMPUTER PREREQUISITES	7
2	OPERATING GUIDELINES	8
2.1	3RD PARTY SOFTWARE AND UPDATES.....	8
3	SOFTWARE INSTALLATION	9
3.1	INSTALLATION	9
3.2	FIRMWARE.....	12
3.3	PANORAMIC EQUIPMENT SELECTION.....	13
4	SENSOR ALIGNMENT AND CALIBRATION.....	14
4.1	SENSOR ALIGNMENT.....	14
4.2	CALIBRATION	15
5	IMAGE MANAGEMENT.....	17
5.1	STORING IMAGES IN WINDOWS.....	17
5.2	IMPORTING INTO IMAGE MANAGEMENT SOFTWARE.....	18
5.3	TWAIN INTEGRATION.....	18
6	SOFTWARE OVERVIEW	19
6.1	MENUS OVERVIEW	19
6.1.1	<i>File</i>	<i>19</i>
6.1.2	<i>View.....</i>	<i>20</i>
6.1.3	<i>Patient</i>	<i>20</i>
6.1.4	<i>Image.....</i>	<i>20</i>
6.1.5	<i>Tools.....</i>	<i>20</i>
6.1.6	<i>Window.....</i>	<i>21</i>
6.1.7	<i>Help.....</i>	<i>21</i>
7	ADVANCED SETTINGS & CONFIGURATION.....	22
7.1	GENERAL	22
7.2	PANORAMIC UNIT.....	23
7.3	TWAIN	24
7.4	IMAGE.....	25
7.5	FILES AND SAVING	26
7.6	DICOM	27
7.7	ADVANCED	27
8	INTERFACE OVERVIEW.....	28
8.1	STATUS INDICATOR (TRAFFIC LIGHT).....	28
8.2	HISTOGRAM.....	28
8.3	TEMPERATURE.....	29
8.4	TOOLS.....	30
8.4.1	<i>Magnifier.....</i>	<i>30</i>
8.4.2	<i>Move Tool.....</i>	<i>30</i>
8.4.3	<i>Region Selection Tool.....</i>	<i>30</i>
8.4.4	<i>Highlight Tool</i>	<i>31</i>
8.4.5	<i>Enhance.....</i>	<i>32</i>
8.4.6	<i>Focus</i>	<i>32</i>
8.4.7	<i>Colors.....</i>	<i>33</i>
8.4.8	<i>User.....</i>	<i>33</i>
8.4.9	<i>DICOM.....</i>	<i>33</i>
8.4.10	<i>Undo.....</i>	<i>33</i>
8.5	PANORAMIC UNIT PROGRAM	34
8.6	SAVE TWAIN IMG.....	34
9	POWER UP AND SYSTEM CHECK.....	35

10	SETTING THE EXPOSURE LEVELS	36
11	TAKING A PANORAMIC IMAGE.....	37
12	TAKING A TMJ IMAGE	38
13	USING THE IMAGE WORKSPACE.....	39
13.1	HIDING THE TOOLBAR	39
13.2	WORKING WITH MAGNIFIED IMAGES	39
14	HISTOGRAM ADJUSTMENT & DIAGNOSTIC TOOLS.....	40
14.1	HISTOGRAM TOOL.....	40
14.2	IMAGE	42
14.3	USER IMAGE PROCESS	44
15	FOCUS IMAGING	45
15.1	INTRODUCTION	45
15.2	AUTO FOCUS (LOCAL).....	46
15.3	MANUAL FOCUSING.....	48
16	IMAGE REPOSITORY	49
16.1	SETTING UP THE IMAGE REPOSITORY	49
16.2	VIEWING A SAVED REPOSITORY IMAGE	49
17	PROBLEM SOLVING	50
18	APPENDIX - DENTAL SOFTWARE INTEGRATION.....	51
18.1	APTERYX XRAYVISION 3.5-3.9 & PROFESSOR SUNI.....	52
18.2	PROIMAGE 6	54
18.3	APTERYX XVLITE 2.0.....	55
18.4	MOGO 11.01.....	56
18.5	DEXIS 8 WITH DEXSCAN.....	57
18.6	KDIS 6.4 KODAK DENTAL IMAGING SOFTWARE	58
18.7	DENTRIX 11 & DENTRIX IMAGE 4.5.....	63
18.8	DENTALEYE 3.1	67
18.9	SCHICK CDR DICOM FOR WINDOWS 3.5.....	69
18.10	CLINIVIEW 6.1.2.....	70

Table of Figures

FIGURE 1-1 HARDWARE PREREQUISITES	7
FIGURE 3-1 SOFTWARE INSTALLATION	9
FIGURE 3-2 FULL INSTALLATION OR SOFTWARE ONLY OPTION	10
FIGURE 3-3 FRAME GRABBER DRIVER INSTALLATION	10
FIGURE 3-4 FINISH INSTALLATION	11
FIGURE 3-5 FIRMWARE UPDATE TOOL	12
FIGURE 3-6 FIRMWARE UPDATE	12
FIGURE 4-1 MECHANICAL ALIGNMENT TOOL	14
FIGURE 4-2 SUCCESSFUL CALIBRATION	16
FIGURE 5-1 INFRANVIEW BY INFRAN SKILJAN	18
FIGURE 7-1 TOOLS > SETTINGS > GENERAL	22
FIGURE 7-2 TOOLS > SETTINGS > PANORAMIC UNIT	23
FIGURE 7-3 TOOLS > SETTINGS > TWAIN	24
FIGURE 7-4 TOOLS > SETTINGS > IMAGE	25
FIGURE 7-5 TOOLS > SETTINGS > FILES AND SAVING	26
FIGURE 7-6 TOOLS > SETTINGS > DICOM	27
FIGURE 8-1 HISTOGRAM	28
FIGURE 8-2 TEMPERATURE PLOT	29
FIGURE 9-1 SYSTEM STATUS	35
FIGURE 10-1 EXPOSURE ON AN OBJECT	36
FIGURE 12-1 TMJ IMAGE	38
FIGURE 13-1 HIDING & SHOWING THE TOOLBAR	39
FIGURE 13-2 HIDDEN SLIDERS IN MAGNIFIED VIEW	39
FIGURE 14-1 HISTOGRAM ADJUSTMENT	40
FIGURE 14-2 IMAGE CONTROLS	41
FIGURE 14-3 NUMERIC VALUES FOR IMAGE CONTROLS	41
FIGURE 14-4 IMAGE PROCESS OPTIONS	42
FIGURE 14-5 USER ASSIGNED IMAGE PROCESSING	44
FIGURE 15-1 FOCUS IMAGING	45
FIGURE 15-2 AUTO FOCUS (LOCAL) EXAMPLE	46
FIGURE 15-3 FOCUSING AN IMAGE	47
FIGURE 15-4 MANUAL FOCUSING WINDOW	48

1 Installation

1.1 Computer Prerequisites

The computer running AJAT's software must fulfill the following requirements.

Item	Requirement
CPU	SSE2 compatible processor with performance equal to or more than in Intel Pentium 4 @ 3 GHz or Core 2
Memory	2 GB DDR-333 or higher
Chipset	NVidia NForce or Intel (VIA, SIS or ATI chipset are not supported)
Hard Drive	1 GB of free disk space for AJAT's software.
Operating System	Windows XP Professional SP2, Vista Business
PCI	One PCI 2.2 compatible 3v or 5v full height, bus master slot
Display Adapter	AGP or PCI Express connection or integrated non-PCI connected. PCI display adapters are not supported unless approved by AJAT separately (high definition medical displays)

Figure 1-1 Hardware Prerequisites

2 Operating guidelines

It is strongly recommended that a dedicated PC is used for the panoramic imaging application! All hardware devices not required for Panoramic Imaging use should be removed or disabled in BIOS and/or Windows' Device Manager.

The PC used for Panoramic Imaging application should be protected by an Uninterruptible Power Supply (UPS) or at least by a surge protector. Using a medical-grade PC with a high-quality power supply is highly recommended! Power supply failure in the Panoramic Imaging PC may irreversibly damage both SCAN300FP system and the panoramic unit and is not covered by the warranty!

The PC used for Panoramic Imaging application and the panoramic unit must be connected to the same power grid with common ground.

2.1 3rd Party Software and Updates

AJAT requires that upgrades (Service Packs or Hotfixes) to the Windows operating system are not installed after the installation of the sensor before AJAT has verified the compatibility.

AJAT requires that the computer used for capturing images from the SCAN300FP sensor doesn't have real-time virus protection or software firewalls active. Also, it should be verified that no background application (for example automatic update service for virus protection software) can activate during data collection.

3 Software Installation

Note: If this software will be used with an AJAT sensor, be sure to physically install the PCI frame grabber prior to beginning the software installation. The necessary drivers will be automatically installed during the setup process.

3.1 Installation

Double click the AJAT Panoramic installation package to begin the installation.

Select Next on all screens to perform a typical installation.




Figure 3-1 Software installation

The select components window provides the options for a full installation for systems with the AJAT sensor attached or software only for standalone workstations.




Figure 3-2 Full installation or software only option

Select Continue Anyway to install the frame grabber and drivers.




Figure 3-3 Frame grabber driver installation

On the last screen, select Yes > Finish.




Figure 3-4 Finish installation

After closing any open software, restart the computer.

3.2 Firmware

After successfully installing the AJAT Panoramic Imaging software and restarting the computer, please run the frame grabber firmware update tool.

Start the firmware update tool from Start>All Programs>Panoramic Imaging Software >Debug>PC2-CamLink firmware upgrade.

Select Yes to begin the update.




Figure 3-5 Firmware update tool

Select Update > StartPC2-CamLink update. Restart the computer when the update is finished.




Figure 3-6 Firmware update

3.3 Panoramic Equipment Selection

Start the AJAT Panoramic Imaging

Select Tools > Settings

Select the Manufacturer and model of the panoramic equipment to be used.

Click Ok.

Restart the Computer after making any selections.

4 Sensor Alignment and Calibration

Note: Please refer to the documentation provided with the sensor for any specific sensor installation, alignment or calibration procedures.

4.1 Sensor Alignment

Select Tools > Maintenance > Mechanical Alignment

The alignment window provides a graphical representation of the X-RAY beam and instructions to achieve optimal alignment.




Figure 4-1 Mechanical alignment tool

Sometimes the software is not able to fully detect the beam and may only display arrows without any numbers. In this case, the sensor should be moved 2-3 millimeters in the direction of the arrows.

The Instructions window provides information specific to the panoramic unit.

While the panoramic unit is producing radiation, you will see an illuminated area representing the position of the actual beam detected by the sensor.

The numbers located on the top and bottom of the window indicate the direction the sensor should be moved for proper alignment. (seen from behind the sensor)

Once the green "Alignment ok" is displayed, the alignment is complete.

Be sure to check the alignment again after replacing any covers.


4.2 Calibration

The software calibration tool is used after any alignment or hardware adjustments performed to the panoramic machine or sensor. Calibration is to be performed after the sensor has been properly aligned, securely mounted and the entire panoramic unit has been fully assembled with all covers. Any alignment parts that may obstruct the XRAY beam must be removed. This includes bite pieces, chin rests, and forehead alignments.

Calibration allows the imaging software to get to know a particular panoramic machine and makes software adjustments for optimum performance. For this process, you will need the Calibration Block provided with your sensor. On screen instructions are provided as the specific process varies depending on your panoramic equipment.

To perform a calibration click Tools> Maintenance> Calibration. Attach the supplied calibration block over the tube head covering the focal spot. Click Calibrate when ready.

The calibration process consists of several steps which include making three exposures. Follow the on screen directions for your particular panoramic unit.



The screenshot shows a software window titled "Panoramic calibration" with a blue header. Inside, there is a "Status" tab. Below the tab, a list of steps is displayed: Step 1: Exposure length measurement NOT NEEDED, Step 2: First X-ray exposure, Step 3: 2nd X-ray exposure, Step 4: 3rd X-ray exposure, Step 5: Calculate parameters, Step 6: Save results, and Step 7: Load calibration. Below the list, there are two progress bars: "Current step" and "Total progress", both of which are empty. At the bottom of the window, there is a "Calibrate" button.

Once the calibration is complete, the results will be displayed on the right. If all the categories are green, the calibration was successful. Yellow graphs indicate that maximum tolerances are close to being exceeded and should be checked. If the results are green or yellow, verify that the sensor installation was successful by taking a test image. Red graphs indicate that a maximum tolerance has been exceeded. Please refer to the Troubleshooting section. Never accept an installation if any of the fields is red.




Figure 4-2 Successful calibration

5 Image Management

The AJAT Panoramic Imaging software has been designed to adapt to any storage environment. Images may simply be saved to a Windows folder or integrated into an existing dental imaging application using the TWAIN interface.

5.1 Storing Images in Windows

Microsoft Windows XP contains all the necessary tools to save, view and print your panoramic images, if you decide not to use dental imaging software. The Windows Picture and Fax Viewer allow viewing and printing to an installed printer.

Create a folder for saving images.

Locate a place of your computer to store your images and right click > New > Folder. Name the folder Panoramic Images.

From the Panoramic Software, select File > Save Image as...

Locate the Panoramic Images folder.

Name the image. I.e. Doe, John

Available Image Types:

- Tagged Image File (TIF) 8 bit
- Tagged Image File (TIF) 16 bit
- Portable Network Graphics (PNG) 8 bit
- Portable Network Graphics (PNG) 16 bit
- DICOM (DCM) 8 bit
- DICOM (DCM) 16 bit
- Binary 32 bit
- Text (TXT)

Viewing and Printing

Locate the Panoramic Images folder

Open the desired image using the Windows Picture and Fax Viewer.

Select the printer icon at the bottom the window or type (CTRL+P).

5.2 Importing into Image Management Software

Most dental image management software has a feature that allows importing images that were acquired outside of its environment. Please consult with your software provider for further instructions.

5.3 TWAIN Integration

TWAIN is a common language that most software uses to communicate with imaging devices. If your software has support for a scanner, you will be able to use this interface for a seamless integration. Please consult with your software provider for further instructions. The TWAIN interface is versatile enough to allow even non dental software to be used as well, as seen in the illustration below.

Non-dental imaging software




Figure 5-1 InfranView by Infran Skiljan

6 Software Overview

6.1 Menus Overview

6.1.1 File

Open image file	Opens images saved as PNG,TIFF or DAT.
Open dataset	Opens a saved dataset.
Reopen	Lists the ten recent images for quick opening
Data repository	Allows reopening of automatically saved images or dataset.
Save Image As	Allows the saving of an acquired image. The formats are Tagged Image (TIF) 8 bit, Tagged Image (TIF) 16 bit, Portable Network, Graphics (PNG) 8 bit, Portable Network Graphics (PNG) 16 bit , DICOM (DCM) 8 bit, DICOM (DCM) 16 bit, Binary 32 bit and Text (TXT).
Save Dataset As	Allows you to save the full acquired data. The full dataset will be several hundred megabytes and can only be viewed in this software.
Save image as TWAIN	Returns an acquired image back to the software that called it. (TWAIN mode only)
Close image	Closes currently displayed image
Close dataset	Finalizes dataset
Exit	Closes the panoramic software

6.1.2 View

Show mean frame signal level	Displays the average signal strength of the last exposure for support purposes.
Show overview image	Enables or disables the small overview window used while viewing magnified images.
Show temperature plot	Enables or disables the sensor temperature status windows.
Show function bar	Enables or disables tools and status information displayed above the image area.

6.1.3 Patient

Patient information	Opens a dialog box for entering patient information.
---------------------	--

6.1.4 Image

Denoise Gaussian	Noise removal tool
Denoise Median	Smoothing tool
Sharpen unsharp mask	Sharpening tool
Edge enhancement horizontal	Sharpens horizontal edges
Edge enhancement vertical	Sharpens vertical edges
Apply auto-leveling	Automatically changes the grey levels to the best use of the display's dynamic range.

6.1.5 Tools

Settings	Opens the advance setting area
----------	--------------------------------

Dataset Tools	Advanced dataset handling tools
DICOM	Opens the DICOM saving and sending options
Maintenance	Sub menu that contains mechanical alignment, calibration, diagnostics tools and open log file tool.
Send error report	Can be configured to electronically send reports to support personnel.

6.1.6 Window

Tile horizontally	Displays multiple open images side by side.
Tile vertically	Displays multiple open images vertically stacked.
Tile square	Displays multiple images in equal square pattern.
Cascade	Displays multiple images as a cascade.
Dataset image	Displays the current dataset on top of all open images.

6.1.7 Help

About	Displays the software date and version.
-------	---

7 Advanced Settings & Configuration

The software's advanced settings and configurations are located under Tools>Settings. The items in this area are used to adjust the functionality of many features which generally will not need to be changed after the initial setup.

7.1 General

The general tab of the settings area is used for patient information handling. Additionally, there is an option for automatically maximizing the view when opening the software.




Figure 7-1 Tools > Settings > General

Patient information is generally used, when displaying a patient's name and ID number on an image is necessary or when not provided by existing dental image management software. In this window, a selection can be made to prompt for this information before or after an image is made.

The patient DB bridge allows other applications to pass patient information to be displayed on an image using a temporary UTF-8 XML file. The location of this file should be constant and indicated in the location shown above. This file is automatically deleted after each usage.

XML Example:


```
<?xml version="1.0"?>
<Examination>
  <Patient>
    <PatientID>12346</PatientID>
    <GivenNames>FIRST</GivenNames>
    <FamilyName>LAST</FamilyName>
  </Patient>
</Examination>
```

7.2 Panoramic unit

This tab is used to configure the panoramic software to work with your specific panoramic equipment. Expand the manufacturer name and select the correct supported model. The "Show advanced programs" option will enable additional control of panoramic and TMJ image capturing. The default option is not enabled and is the recommended setting.

Rotation direction should be selected to match the exposure direction of the panoramic unit.

Asymmetry reduction is used to remove some unnecessary artifacts from the edges of an image caused by the rotation start and stop. The effectiveness of this option depends on the panoramic unit and model.



Selected panoramic unit model

Manufacturer Panoramic Corporation

Model PC4000 US

Direction of rotation

☒ Left to right (Clockwise)

☐ Right to left (Counter clockwise)

☒ Show advanced programs

☐ Asymmetry reduction

Figure 7-2 Tools > Settings > Panoramic unit

7.3 TWAIN

The TWAIN selection screen allows the user to customize how the TWAIN functionality is implemented and offers several different GUI (Graphical User Interface) options. Full application mode and Automatic compatibility are selected by default.




Figure 7-3 Tools > Settings > TWAIN

7.4 Image

The image settings area provide selectable options that change the appearance of the acquired image to appear more like a traditional film. Orientation imprinting option will place an L & R marker in the corners of the image to indicate the direction of the exposure.

Patient information imprinting options control the display location of the patient information entered or passed by other application.

Image options are preset to the ideal configuration and generally should not need to be changed.




Figure 7-4 Tools > Settings > Image

7.5 Files and saving

The Files and saving tab provides options to customize the background data repository feature. This image protection feature is not intended to be a replacement of a solid daily method of backup, but more as a convenience. When enabled, every image or dataset acquired is automatically time stamped and saved to the repository

Options to remove the unnecessary white edges, only store unsaved images and disk space allowances may also be configured in this tab as well.




Figure 7-5 Tools > Settings > Files and saving

7.6 DICOM

The DICOM storage and printing servers may be configured in the section. This information should be provided and tested by an administrator.

DICOM file save folder selection allows a default file path to be set. When using the DICOM Save tool, this path will be used.




Figure 7-6 Tools > Settings > DICOM

7.7 Advanced

CPU Performance optimization is generally automatically performed after the initial installation and does not need to run again.

8 Interface Overview

8.1 Status Indicator (Traffic Light)

Green	The sensor and software are ready to acquire and image.
Yellow	The system is preparing to go to ready mode.
Red	The system is not ready. See Troubleshooting.

8.2 Histogram

A histogram or graphical representation of how many times a color occurs in an image. Use this tool to adjust the brightness and intensity for the overall image. Additional information such as gamma can be accessed by right clicking on the control window and selecting "show output control" or "show numeric fields".




Figure 8-1 Histogram

8.3 Temperature

The optional temperature view is a real time report of the current temperature of the sensor. This window can be toggled by selecting View>Show temperature plot. By default this control is not enabled.




Figure 8-2 Temperature plot

8.4 Tools

8.4.1 Magnifier

This is a spot magnifier that is controlled by the mouse while the left button is pressed. More magnification is available by right clicking on the tool.



8.4.2 Move Tool


This is used to move the image when magnified beyond the limits of the screen. When desired the entire image may be magnified by right clicking anywhere on the image and select a desired magnification. Use the move tool to navigate around the image.



Optionally, a small overview may be enabled for easier navigation while the image is magnified. This feature is enabled by selecting View>Show overview image.

8.4.3 Region Selection Tool

Used for selecting a region of interest during focusing. Select a starting position and click and hold the left mouse button. While still holding the left button, move the mouse to enlarge the box to the desired size.




8.4.4 Highlight Tool

The highlight tool is used by moving the mouse while pressing the left mouse button. This tool produces a high contrast filter, revealing information in difficult regions. Right click the mouse button to choose different sizes for the "High Light" square.




8.4.5 Enhance

Image enhancement tools are provided for more advance diagnostics. Denoiseing, Sharpening and Enhancing can be quickly applied to improve quality, overcome anatomical defects or enhance areas of interest. Additional information is provided in the Histogram Adjustment & Diagnostic Tools chapter.




8.4.6 Focus

The focus area contains several advanced focusing tools. The area selection tool is used to select an area of the image that is outside the default focal trough and then the auto focus can be applied. The manual focus or tilt tool may also be used for greater control. Detailed descriptions of these features are provided in the Focus Imaging Chapter.



8.4.7 Colors

The colors section offers four different color schemes as well as invert. These colors can be toggled on and off by left clicking the desired color.



8.4.8 User

This area is used to configure and apply user-defined image processing. Many times you will find yourself using the same image enhancement or other processing tools. To save time, the auto process button can be configured to automatically apply any combination of tools available in the image process area. Four additional configurable tools are available, quickly configurable to any preference.

8.4.9 DICOM

DICOM options include saving, printing, and sending to a server. Sending and printing require advanced configuration and must be configured by an administrator.

Save will automatically save an image in the DICOM (DCM) format to a preconfigured location and automatically name the image based on the patient information provided.

8.4.10 Undo


The Undo incrementally (once per click) removes image processing that has been applied by the user utilizing the "Enhance" area.

8.5 Panoramic unit program

This allows you to tell the software which direction your panoramic unit rotates during an exposure or select unit specific features such as TMJ.

8.6 Save TWAIN img

Returns an approved image to a program such as your practice management software. (TWAIN mode only)



9 Power Up and System Check

Note: Refer to the documentation provided with the sensor for any specific instructions for start up.

Note: Refer to the documentation provided by the manufacturer of the panoramic unit for proper power up procedures.

1. Open that AJAT Panoramic Software.
2. Power on the panoramic unit.
3. Verify that the Green Light is illuminated and that system is ready.



Figure 9-1 System status

10 Setting the Exposure Levels

Note: Exposure levels must be followed strictly according to the instructions of the panoramic equipment manufacturer.

Note: Taking several test images on inanimate objects is a good way to safely get comfortable with the panoramic units exposure levels.

Note: Only a person trained in dental radiography should be using the panoramic unit!




Figure 10-1 Exposure on an object

11 Taking a Panoramic Image

Select the correct panoramic rotation direction. The software must match the setting on the panoramic unit's control panel.


Position that patient as specified by the panoramic unit's manufacturer. Proper patient positioning is key to successful panoramic imaging.

Adjust the panoramic unit to the proper dosage according to the panoramic unit manufacturer's instructions.

Activate the panoramic machine.

Note: refer to the manufacturer's documentation for specific usage instructions.

During the exposure, if the traffic light turns to green or red, the exposure should be interrupted to prevent unnecessary exposure.



12 Taking a TMJ Image

The TMJ mode produces one image consisting usually of four short exposures. The TMJ image shows the patient's right closed, right open, left open and left closed temporomandibular joints.

Depending on the panoramic unit, it may be necessary to change the panoramic unit program in the software from panoramic to TMJ while others will remain in the automatic mode.

Manually changing to TMJ mode:

Locate the panoramic unit program box in the upper left corner of the screen. Click on the current mode button and select TMJ. Remember to set this back to panoramic before taking another exposure.

The TMJ sequence will automatically advance to the next position and will be complete when all four images are taken.




Figure 12-1 TMJ image

13 Using the image workspace

The image workspace consists of an area where an X-RAY is displayed, tools for diagnostic viewing, and information about the panoramic digital sensor. This section will provide an overview of the built in controls you may find useful while working with panoramic images.


13.1 Hiding the toolbar

The upper tool bar can be easily hidden to allow a panoramic image to be viewed at maximum size. Click the small arrow located on the bottom right area of the toolbar to collapse the bar. Click the arrow again or press ESC to expand.



Figure 13-1 Hiding & showing the toolbar

13.2 Working with magnified images



Method 1 - Images may be magnified by right clicking anywhere on the image and selecting a desired magnification. By default, view to fit is automatically selected.

While magnified, hidden right and bottom scroll bars can be used to navigate around the image. Move the mouse to the far right or bottom to reveal the hidden scroll bars. The blue scroll position is moved by left clicking and moving.

Method 2 - The optional overview image, enabled by selecting View>Show overview image, may be used to relocate your position on the image. To move, left click the small red rectangle and slide it to a new location. Entire image magnification is adjustable using the slider located at the bottom of the overview window.

Figure 13-2 Hidden sliders in magnified view

14 Histogram Adjustment & Diagnostic Tools

14.1 Histogram tool

The image adjustment tool provides a graphic histogram of an image, representing how many times a color occurs in an image. Errors in dosage levels can be easily noticed and corrected using the image adjustment tool.

The histogram is adjusted using three sliding arrows. The arrow to left represents the white end of the dynamic range and the arrow to the right, the black end. The arrow in middle is used to shift the histogram toward the black or white (gamma correction). The result of making a shift will be seen as brightening or darkening of the image.

To change the contrast of an image, black or white may be clipped by moving the far left and right arrows towards its opposite side. Normally shifting the histogram is sufficient and the contrast can be left alone. Below are some examples of shifting and clipping.




Figure 14-1 Histogram adjustment

To move any of the arrows, left click on the arrow and move to a new location. As the pointer is moved the image will be adjusted in real time.

Additional adjustment information is available by right clicking on the tool and selecting any desired information.




Figure 14-2 Image controls

Additional controls may be enabled, as shown below, to view or change the histogram numerically. Right click on the Image Controls windows and select the desired views.




Figure 14-3 Numeric values for image controls

14.2 Image

The image processing tab contains several enhancement tools for smoothing, sharpening and edge enhancing. These tools may be applied to an image in any order or combination and can be easily undone using the undo button. Below is a description of each tool.




Figure 14-4 Image process options

1) Denoise Lowpass

Applies lowpass filtering with a Gaussian kernel to the image. This tool helps to remove normal (white) noise from the image to make an image appear smoother or less noisy. Several levels of this filter are available by right clicking on the tool and selecting a different size.

2) Denoise Median

Applies median filtering to the image. Median filtering smoothes the image by removing the salt and pepper type of noise that appears as dots in the image. Several levels of this filter are available by right clicking on the tool and selecting a different size.

3) Sharpen

Sharpens an image by applying an unsharp mask method. Unsharp masking enhances boundaries in the image so that the image appears clearer. Additional degrees of sharpening may be selected by right clicking the tool and choosing the desired filter size.

4) Equal Vertical

The equal vertical tool adjusts the local brightness of an image so that brighter areas appear less obvious. This tool is useful for evening the appearance of a panoramic image and diminishing the bright artifacts caused by the spine (spine compensation).

5) Vertical Edge Enhancement

Vertical edge enhancement is similar to sharpening but only vertical edges are enhanced.

6) Horizontal Edge Enhancement

Vertical edge enhancement is similar to sharpening but only horizontal edges are enhanced.

14.3 User Image Process

The user image process area provides several user configurable buttons that can be used to apply frequently used image processing tools as desired. One additional tool can be configured to apply image processing automatically when images are acquired.

Use the following steps to configure user button:

- 1) Select and apply any smoothing, sharpening or edge enhancement desired.
- 2) Select the User Image Process tab
- 3) Right click the desired configurable button and select "Assign current image processing to desired button".
- 4) Type a name for the button and click ok.




Figure 14-5 User assigned image processing

15 Focus Imaging

15.1 Introduction

The SCAN300 panoramic imaging system utilizes cutting edge technologies to acquire ultra high resolution images. In addition, the imaging system allows changing locally the focal trough location within approximately a 3 centimeter range depending on the location.

The Ajat Panoramic software typically collects a dataset consisting of 100-400 mega-bytes of frames per exposure. The large dataset of information gives the software the ability to move through the image and find locally the most focused layers to build the panoramic image.

The SCAN300FP expands the focal trough and collects more diagnostic images.





Figure 15-1 Focus Imaging

15.2 Auto Focus (local)

Auto Focus is a tool for guiding the software to make additional focal adjustments that may help to further enhance the diagnostic quality of an image.

Auto Focus consists of two steps:

First the selection tool  is used. Select the area to be focused. It's recommended to select an area slightly larger in both vertical and horizontal direction to improve the results. This provides the software the necessary information to perform the refocus.

Once an area is selected, use the Auto Focus  button to make the change.





Figure 15-2 Auto focus (local) example

In the image below, the anterior region is noticeably out of focus. The selection tool is used to select an area, within the local area to be focused, that is considered to be out-of-focus. The auto focus tool is then used and the difference can be seen clearly in the second image. It should be remembered that the focusing alters the focal trough for the selected part and thus changes also the region below and above the selected area.



Unfocused image



Focused image

Figure 15-3 Focusing an image

At any point the Default Focus button can be used to undo all changes. If the desired focal result is not achieved after several attempts, it may mean that there is not enough information in the selected area and the profile should be reset to the "Default focus".

15.3 Manual Focusing

Manual image focusing works in the same way as the auto focus tool described above, but with the added ability to control in real time the focal point position as well as tilt.

Begin by using the region selection tool to select an area you would like to focus then select manual focus. The focus control window will now appear with two vertical sliding controls.




Figure 15-4 Manual focusing window

Move the vertical sliders until the desired focus is achieved. The graphical representation at the top of the slider will illustrate the current focal spot once the Apply button is selected.

The close button will exit this window and prompt to save the changes.

Any focusing can easily be undone by selecting the Default Focus button in the focus area.

16 Image Repository

The image repository is a safety feature that, when enabled, automatically saves acquired panoramic images or datasets to allow recovery in the event an image is lost or damaged. The repository is intended for convenience and as a last resort to image recovery. The repository is not a replacement for a full system daily backup.

16.1 Setting up the Image Repository

Select File Tools>Settings> Files and saving

To protect and image only until it's successfully save, select Remove backup image after save. If this options is not selected, the repository will automatically purge based on the space allocated and number of max images.

Make sure that there is sufficient free space on the hard drive to accommodate the maximum size repository and max number of images selected.

When finished, click Ok.

16.2 Viewing a saved repository image

To view or recover an image or dataset, select File > Data repository

Select a saved image. Be careful, once an image is deleted from the repository, it's permanently removed!

Images will be automatically deleted based on the settings selected when enabled. The delete and delete all should not need to be used.

17 Problem Solving

Installation

Q. Status indicator flashes Init Failed and sensor will not function.

A. This is normal condition that may occur periodically and is corrected by resetting the sensor.

Q. Status indicator displays Not Connected.

A. Verify that all cables are connected and that all equipment is powered on.

Image Quality

Q. Image is a blurry and the spine is very clear.

A. Verify that the software's rotation direction matches the panoramic machine.

Q. The image has poor tone quality or excessive background noise.

A. The image is under exposed. Increase the kVp or mA.

Q. Excessive dark areas in image or poor tone.

A. The image is over exposed. Reduce the kVp or mA.

Q. Image is solid black.

A. The image is excessively over exposed.

Q. Anterior teeth are out of focus, blurry and narrow.

A. Patient is positioned forward of the focal trough.

Q. Anterior teeth are out of focus, blurry and wide.

A. Patient is positioned too far back from the focal trough.

Sensor

Q. Solid red light on Frame Grabber.

A. Sensor is not powered on or sensor link cable is disconnected.

Q. Slow blinking light on Frame Grabber.

A. Make sure that the panoramic service is started. Restart computer and re-check.

Q. Sudden reduced image quality or lines occur.

A. Sensor may be misaligned. Check mechanical alignment.

18 Appendix - Dental Software Integration.

The following section provides specific support integration unto the following dental imaging programs*.


Apteryx XrayVision 3.5-3.9 & Professor Suni 3.5-3.9 - Done
ProImage 6
XVLite 2.1 & Lightyear Technology 2.
MOGO 11.01
Dexis 8 with DexScan
KDIS 6.4
Dentrix 11 & Denterix Image 4.5
DentalEye 3.1
Schick CDR DICOM for Windows 3.5
CliniView 6.1.2

* Any mention of dental imaging products, with the exception to Ajat branded products, is for personal reference only and no form of relation is implied.


18.1 Apteryx XrayVision 3.5-3.9 & Professor SUNI

Method One

Open the desired patient and select "Capture Single".




Make sure the "Scan300FP" is selected in the Scanned Image and click "Ok".



Take an exposure and select "Save TWAIN Image".


Select the center area of the tooth association window and then select "Ok".




Method Two

Open the desired patient.

Select File > Acquire > Capture Single Image.




Make sure the "Scan300FP" is selected in the Scanned Image and click "Ok".



Take an exposure and select "Save TWAIN Image".



Select the center area of the tooth association window and then select "Ok".



18.2 ProImage 6

Configuration

Select File > Select Scanner (Twain)

Select SCAN300 as the source.

Usage


Select File> Scan Image (Twain Input)

Take an exposure and select "Save TWAIN Image".

18.3 Apteryx XVLite 2.0

Select a patient.


Select Acquire Image.



Select Take and image for all teeth.



Make sure the "Scan300FP" is selected in the Scanned Image and click "Ok".



Take an exposure and select "Save TWAIN Image".

18.4 MOGO 11.01

Setup

Open MOGO Windows Image Editor.

Select File > Select TWAIN Device.



Select the SCAN300 and click Select.



Usage

Open MOGO Windows Image Editor

Select File > TWAIN Acquire.

Take an exposure and select "Save TWAIN Image".

18.5 Dexis 8 with DexScan

Note: DexScan requires an additional license to enable Twain acquisitions.

Note: The dexis.ini files must be edited and the value fits=yes added to the [SCAN300FP] section. The ini file is located in the Dexis root folder and can be opened with Windows Notepad. The section should look like the following.

```
[SCAN300FP]
fits=yes
frame=0
camframe=1
```

Select the SCAN300FP in the Settings menu.

Choose a patient.

Select the scanner icon.


Take an exposure and select "Save TWAIN Image".

18.6 KDIS 6.4 Kodak Dental Imaging Software


To start the program, click the KODAK icon on the desk top
After you have highlighted the patient you wish to open, click the OPEN file/folder

Once the patient file is OPEN, you will see the patient info displayed as it is below

Once the patient information is displayed, click on the IMAGING icon to proceed to the next screen.



Once the imaging section is OPEN, your screen will look similar to this




To acquire an image using the AJAT sensor follow these steps:

Click FILE on the tool bar to get a drop down menu


In the drop down menu, first click the TWAIN option

Once TWAIN is selected, you will see two options, Acquire and Select Source


First select the option "select source"



The next screen you will see should look like this
 If multiple TWAIN compliant devices are connected to the computer
 Look in the list and highlight and select the option "SCAN300FP"



Next click on File from the tool bar again to get the drop down menu
 This time select the TWAIN option, and the Acquire tab
 This will automatically launch the AJAT software




Once the AJAT program is launched, you will see the acquisition software over the top of the KDIS software
 Verify the stoplight is GREEN, and the PC-1000 is turned ON
 Position your patient and take your exposure
 After the image is captured your screen should look similar to this
 From this screen you can manipulate the image with the tools from the AJAT software



This is an example of the image after all modifications have been done and we are ready to save the image into KDIS software.

To SAVE the image, click the Save TWAIN Img button on the right side of the screen


This shows the image transferred directly into the KDIS software
Once in the KDIS software the operator would have the ability to use all of the tools within KDIS to manipulate and modify the image.



To save the image into KDIS, simply click the disc icon on the tool bar.



Once you have clicked the icon to save the image, you will see a box asking for comments to be attached to the image, and what classification for the image (Not classified, Extraoral or Intraoral) then click OK




Once the patient is selected and the name appears in the blue title bar, you will then have the ability to capture an image.

To capture an image, first you must get into the Imaging side of DENTRIX


To do this, please click on the icon which resembles an intraoral image.



The imaging module will then launch and automatically ask if you would like to select a template.
Please select Cancel for if you are going to take a Pan image.




After you have hit Cancel, you should have a blank screen.
In this screen, please click the Scanner Icon on the tool bar. This is the fifth button (5th) from the Left.



When you hit the Scanner button, a window will pop up to ask you to Select a Source.

Make sure the SCAN300FP is selected and press the Select button




This will OPEN the AJAT imaging software.

With the program open, please verify the Green light is lit, and verify the temperature setting of the camera in the black screen.

With the Green light you are ready to proceed to patient positioning, and capture you Panoramic image.

With your Pan image captured, you will see the image on the screen and have the ability to manipulate the image as needed in the AJAT program. When finished, click the Save TWAIN img to send the image to the patient file within DENTRIX



With the image in the DENTRIX program you will have the ability to manipulate the images with all of the tools available from Dentrix.

18.8 DentalEye 3.1

The DentalEye imaging software may be easily used the Ajat Panoramic Imaging software by using the TWAIN interface.


Note: The TWAIN plug-in must be installed in order to utilize this feature.

Configuration:


Open DentalEye and select the configuration button on the TWAIN plug-in.



Select DX Panoramic for the image type and select TWAIN source.




Choose the SCAN300FP 5.0 from the sources and click Select.




Open the Ajat Panoramic Imaging software. Click Tools>Settings>TWAIN.

Select Adobe Photoshop / Elements under application compatibility.



Usage:

Click the acquire TWAIN button in DentalEYE. The Ajat Panoramic Imaging software will open.




Once a panoramic image is acquired and adjusted as desired, Select the Save TWAIN Img.




The panoramic image will be returned to DentalEye for editing and saving as usual.

18.9 Schick CDR DICOM for Windows 3.5

Open the Schick software and select Start A New Exam.



Select a patient then select PANO 1x1 from the Select Series window.



Change the acquisition mode to scanner using the Mode button.



When the Scanner is selected, a scanner will be displayed in the acquisition window.

Click the scanner icon to launch the Panoramic Imaging software.


After the panoramic image acquired, select the Save TWAIN Image window. The image will be returned to the Schick software and edited as usual.



18.10 CliniView 6.1.2


Configuration:

Select File>TWAIN>Select Source. Select the SCAN300FP, then click Select.



Usage:

With a patient selected, select File>TWAIN>Capture>Panoramic.



After the panoramic image acquired, select the Save TWAIN Image window. The image will be returned to the Schick software and edited as usual.

