• COLLECTION SET CODE 712 CARDIO, CODE 712 • CIRCUITO DI RACCOLTA COD. 712 CARDIO, COD. 712 •
CIRCUIT DE RECOLTE CODE 712 CARDIO, CODE 712 • SAMMELSYSTEM CODE 712 CARDIO, CODE 712 •
CIRCUITOS DE RECOGIDA CÓDIGO 712 CARDIO, CÓDIGO 712 • SET DE COLHEITA CÓDIGO 712 CARDIO,
CÓDIGO 712 • ΣΕΤ ΣΥΛΛΟΓΗΣ ΚΩΔΙΚΟΣ 712 CARDIO, ΚΩΔΙΚΟΣ 712 • VERZAMELSET CODE 712 CARDIO, CODE
712 • UPPSAMLINGSSET KOD 712 KARDIO, KOD 712 • COLLECTION SET CODE 712 CARDIO, CODE 712 •
KERÄYSSETIT KOODI 712 CARDIO, KOODI 712 • OPPSAMLINGSSETT KODE 712 KARDIO, KODE 712 • TOPLAMA
SETİ (712 CARDIO KODLU VE 712 KODLU)

011670/00 - 06/2009

GB - ENGLISH - INSTRUCTIONS FOR USE	3
IT - ITALIANO - ISTRUZIONI PER L'USO	7
FR - FRANCAIS - INSTRUCTIONS D'UTILISATION	11
DE - DEUTSCH - GEBRAUCHSANWEISUNG	15
ES - ESPAÑOL - INSTRUCCIONES DE USO	19
PT - PORTUGÊS - INSTRUÇÕES DE UTILIZAÇÃO	23
GR - ΕΛΛΗΝΙΚΑ - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	27
NL - NEDERLANDS - GEBRUIKSAANWIJZING	31
SV - SVENSKA - BRUKSANVISNING	35
DK - DANSK - BRUGSANVISNING	39
FI - SUOMI - KÄYTTÖOHJEET	43
NO - NORSK - BRUKSANVISNING	47
TR - TÜRKÇE - KULLANMA TALİMATI	51


fig.1

GB	IT	FR	DE	ES	PT	GR	NL	SV	DK	FI	NO	TR
Fig.	Fig.	Fig.	Abb.	Fig.	Fig.	ΕΙΚ.	Fig.	Fig.	Fig.	Kuva	Fig.	Şek.

GB - ENGLISH - INSTRUCTIONS FOR USE

CATALOGUE N.	PRODUCT DESIGNATION	PRODUCT DESCRIPTION
04051	CODE 712 CARDIO	CARDIO COLLECTION SET WITH BT 894 AUTOTRANSFUSION RESERVOIR
04057	CODE 712	COLLECTION SET WITH BT 894 AUTOTRANSFUSION RESERVOIR

DESCRIPTION

The Collection Set Code 712 Cardio and Code 712 are composed of the following components, respectively:

COLLECTION SET CODE 712 CARDIO	COLLECTION SET CODE 712	DESCRIPTION
BT 894 Autotransfusion Reservoir	BT 894 Autotransfusion Reservoir	Container fitted with filter for removal of aggregates larger than 40 µm diameter from recovered blood, with blood outlet connector on the reservoir cover
Aspiration and Anticoagulation line	Aspiration and Anticoagulation line	Transfer line for anticoagulated blood aspirated in the operating field.
Vacuum Connection Line	Vacuum Connection Line	Reservoir connection line to vacuum outlet
Kit Cardio (containing): - Three-way adapter - Adapter Tube (1/4"-1/4") - Extension Line - Luer-Lock Adapter	Not present	Adapter for connection to the oxygenator

All these devices are single use, non toxic, non pyrogenic, supplied **STERILE** and individually packaged. Sterilized by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

INTENDED USE

The intended use for the Collection Set Code 712 Cardio and Code 712, are respectively:

Collection Set Code 712 Cardio

The Collection Set Code 712 Cardio is intended for sterile recovery and filtration of the blood shed **during cardio-surgical** operations. It is to be used in combination with the aspiration and anticoagulation line and with the vacuum line that are included in the set.

The Collection Set Code 712 Cardio is provided with the Cardio Kit, which includes a set of adapters for connection of the extra-corporeal circulation (ECC) circuit to the wash set of an autotransfusion equipment or to the BT 894 Autotransfusion Reservoir. The adapters contained in the Kit Cardio allow operators to extend the use of the autotransfusion techniques to the following procedures:

- Concentration of priming fluid during ECC
- Recovery of priming fluid residue from the ECC circuit
- Simultaneous use of two BT 894 Autotransfusion Reservoirs.

Collection Set Code 712

The Collection Set Code 712 is intended for sterile recovery and filtration of the blood shed **during non-cardio surgical** operations. It is to be used in combination with the aspiration and anticoagulation line and with the vacuum line that are included in the set.

The products listed above should be used in combination with medical devices listed in section "Medical devices for use with the Collection Sets Code 712 Cardio and Code 712".

SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious adverse reactions and potential safety hazards for practitioner and/or patient that can occur in the proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such cases.

CAUTION

Indicates any special care to be exercised by a practitioner for the safe and effective use of the device.

EXPLANATION OF THE SYMBOLS USED ON THE LABELS

For single use only (Do not reuse)

Batch code (number)
(reference for product traceability)

Use by (Expiry date)

Date of manufacture

Manufactured by

Sterile - Ethylene oxide sterilised


Non Pyrogenic

Contains PHTHALATE

Latex free

CAUTION: Federal law (usa) restrict this device to sale by or on the order of a physician

Warning: Do not resterilize.

Contents sterile only if package is not opened, damaged or broken

Catalogue (code) number

Attention, see instruction for use

This way up

Fragile; handle with care

Keep away from heat

Keep dry

ea Quantity

The following is general safety information with the aim of advising the operator in preparing to use the device. Also, specific safety information is given in the instructions for use at locations in the text where that information is relevant for correct operation.

WARNING

- The device must be used in accordance with the instructions for use provided in this manual.
- The device is intended to be used by professionally trained personnel.

- Do not use the device if it has been cracked, dropped or otherwise physically damaged.
- SORIN GROUP ITALIA is not responsible for problems arising from inexperience or improper use.
- FRAGILE, handle with care.
- Avoid any condition which may cause the blood temperature to exceed 37°C (98.6°F).
- Keep dry. Store at room temperature.
- Always apply and maintain a correct dose and accurate monitoring of the anticoagulant.
- For single use and for single patient use only: during use the device is in contact with human blood, body fluids, liquids or gases for the purpose of eventual infusion, administration or introduction into the body and due to its specific design it cannot be fully cleaned and disinfected at the end of use. Therefore the reuse on other patients might cause cross-contamination, infection and sepsis. In addition, the reuse increases the probability of product failure (integrity, functionality and clinical effectiveness).
- The device must not undergo any further processing.
- Do not resterilize.
- After use, dispose of the device in accordance with applicable regulations in force in the country of use.
- The device must only be used if STERILE.
- The device contains phthalates. Considering the nature of body contact, the limited contact duration and the number of treatments per each patient; the amount of phthalates which might be released from the device do not raise specific concerns about residual risks. Further information is available on request from Sorin Group Italia.
- For further information and/or in case of complaint contact SORIN GROUP ITALIA or the authorized local representative.

⚠ CAUTION

- For a detailed description of the circuits, refer to the user's manual of the SORIN GROUP ITALIA ELECTA cell separator
- Refer to the ELECTA user's manual for complete instructions on use of the device
- Federal law (U.S.A.) restricts this device to sale by or on the order of a physician

SET UP

HOLDERS SET UP

BT894 Holder (Catalogue N. 04108)

- Fit the BT894 holder to an I.V. pole, place the reservoir holder ring at working height and tighten the knob on the side.

Integrated reservoir holder on SORIN GROUP ITALIA ELECTA machine

- Turn the reservoir holder ring on the SORIN GROUP ITALIA ELECTA machine by 180° with respect to the closing position and place it at working height. (fig.1.1 and 1.2)

⚠ CAUTION

- For accurate measuring of the blood volume inside the autotransfusion reservoir, the reservoir arm should be raised at least a few cm.

DISPOSABLE SET UP

⚠ WARNING

- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.
- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- Do not use solvents such as alcohol, ether, acetone, etc. to prevent damage to the device.
- Make the connections using a suitable aseptic technique.
- Due to the possibility of operator exposure to bloodborne pathogens (such as HIV, hepatitis viruses, bacteria, Cytomegalovirus, etc.) when handling extracorporeal blood circuits, adequate precautions should be taken at all times to prevent the exposure to and transmission of such agents.
- Carefully observe this set for leaks before and during use. Leakage may result in loss of sterility or loss of blood and/or fluid. If leakage is observed

before or during use, replace leaking component or retighten leaking connection, as appropriate.

Set up of the BT 894 Autotransfusion Reservoir

1. Extract the product from the protective wrapping. Remove the protection lid.
2. Place the external ring of the reservoir cover into the guide on the inside of the holder so that they fit together correctly (fig.1.2)
3. From this initial position push the reservoir towards the holder until the locking clip is released returning to the rest position (fig.1.3).
4. Check the seal of the caps on the connectors, tightening them in their housings.

⚠ WARNING

- Before use, check that the reservoir is completely inserted in its housing on the holder. In the event of incorrect installation, any accidental knocks against the reservoir might cause it to detach from the holder.

Set up of the vacuum connection line

1. Using a sterile technique, remove the external wrapping and give the content to the surgical area.
2. Connect one end of the vacuum connection line to the vacuum pump filter connector on the rear of the unit.
3. Connect the opposite end of the vacuum connection line to the connector marked VACUUM on the reservoir. Position the clamp on the line in the proximities of the VACUUM connector.
4. If a vacuum pump is not available, connect the opposite end of the vacuum connection line to the vacuum source available in the operating field, by means of a vacuum regulator to control the aspiration level.

Set up of the aspiration and anticoagulation line

1. Using a sterile technique, remove the external wrapping and give the content to the surgical area.
2. Open the internal wrapping and take out the line. Connect the surgical aspirator cannula to the aspiration connector of the aspiration and anticoagulation line and give the opposite end to the operator outside the sterile area.
3. Connect the aspiration line coming from the operative field to one of the horizontal connectors marked SUCTION on the reservoir.
4. Close the flow regulator on the anticoagulation line and insert the spike in the anticoagulant solution container. Keep the line closed until the procedure starts.

Set up of the three-way adapter (for Cod. 712 Cardio only)

1. Using a sterile technique, remove the external wrapping and give the content to the surgical area.
2. If the reservoir is used during cardiac surgery or connected to a second reservoir, insert the three-way adapter, provided in the Kit Cardio, in the reservoir outlet connector marked WASH SET
3. Close the clamps on the three-way adapter branches waiting for subsequent connections. The branch ending with a male port is used to connect the reservoir to the WASH SET, the second branch to connect the extension line going to a circuit or oxygenator for ECC.

Set up of the extension line to a circuit or oxygenator for ECC (for Cod. 712 Cardio only)

1. Before adding the priming liquid to the ECC circuit, connect the 1/4" tube segment contained in the Kit Cardio to a 1/4" connector of the circuit.
2. Close the tubing segment with a clamp, keeping the protective cap in place to maintain sterility. During priming circulation remove the air from the circuit.
3. Take the extension line and connect the two ends to the tubing segment on the ECC circuit and to the free connector on the three-way adapter of the Reservoir or connect it directly to the autotransfusion wash set device.

INSTRUCTIONS FOR USE

Blood collection and anticoagulation

⚠ WARNING

- SORIN GROUP ITALIA recommends not to exceed a negative pressure of -100 mmHg (13 Kpa/0,13 Bar/1,9 psi) during recovery procedure as an excess vacuum increases hemolysis of the recovered blood
1. Adjust the aspiration pressure according to the vacuum source available:
 - 4.1. If the Electa vacuum pump is going to be used, adjust the aspiration pressure to -100 mmHg (13 KPa/ 0.13 bar/1.9 psi).
 - 4.2. If an external vacuum source is going to be used, close all lines (the aspiration line and the tree-way adapter line, if present s), apart from the vacuum line, with a clamp and adjust the aspiration pressure to -100 mmHg (13 KPa/ 0.13 bar/1.9 psi). Reopen the lines when the adjustment has been made, keeping the vacuum on.
 2. Open the aspiration and anticoagulation line and transfer at least 200 ml of solution to the reservoir before starting recovery. The anticoagulant solution containing heparin, ACD-A or CPDA-1, must be prepared beforehand in a quantity

sufficient for the duration of the procedure. If heparin is used, 30.000-60.000 I.U. of heparin per litre of sterile isotonic saline are sufficient to guarantee good anticoagulation of recovered blood.

- After the transfer of the anticoagulant to the reservoir, adjust the flow rate of the anticoagulant solution to between 60 and 100 drops a minute. The lower flow rate is suitable for ACD-A or CPDA-1 and the higher rate for heparin. In case of heavy bleeding, appropriately increase the rate.

NOTE: The flow of anticoagulant into the reservoir is manually controlled by the roller clamp on the anticoagulant line. This flow must be adjusted to the rate of blood collection from the surgical field. If the rate of collection varies without adjusting the flow of anticoagulant, the ratio of anticoagulant to blood can be too low or too high. The blood in the reservoir may clot if there is too little anticoagulant.

⚠ WARNING

- In the event of excess heparin in collection reservoir due to inappropriate ratios, the salvaged blood may contain residual heparin.
- In the event of decreased patient antithrombin III levels if using heparin anticoagulation, consult the physician in order to provide alternate anticoagulation.

⚠ CAUTION

- Failure to maintain adequate anticoagulation during blood collection can cause excessive clotting and possible clogging of the collection reservoir or centrifuge bowl.
 - Prior to collecting any blood, prime the reservoir with approximately 200 mL of the anticoagulant solution.
 - When an anticoagulant solution using ACD-A or CPDA-1 is used during blood recovery, do not use ringer lactate solutions in subsequent wash.
- The reservoir is equipped with a device for interrupting aspiration and preventing the blood collected from entering the vacuum line. Should this level regulating valve activate, follow the instructions below to reset aspiration to the operative field:
 - Temporarily interrupt aspiration and transfer the contents of the reservoir to the bowl until the level of the blood is below the safety valve.
 - Pull the yellow tab located on the connector marked VACUUM and check that the valve has dropped down in the guide.
 - Release the yellow tab and reset aspiration.

⚠ WARNING

- A pressure relief valve is provided to prevent implosion of the device if exposed to a sudden, deep vacuum
 - The negative pressure inside the BT 894 Autotransfusion Reservoir must not exceed -300 mmHg (40 KPa/0.4 bar/5.7 psi).
 - Check the cap seals of the connectors not in use, pushing them fully on.
 - Do not occlude the safety valve release present on the BT 894 reservoir with foreign bodies to prevent the potential possibility of implosion

Using the BT 894 Autotransfusion Reservoir connected to a wash set

- Make the connection between the reservoir and the washing circuit:
 - connect the washing circuit inlet connector to the connector marked WASH SET on the reservoir cover or to the three-way adapter as already described in the specific paragraph. In this case open the connection clamp to the washing circuit on the three-way adapter.

Using the BT 894 Autotransfusion Reservoir connected to a wash set and to a circuit or oxygenator for ECC (for cod. 712 cardio only)

- Connect the wash set connector to the corresponding one on the three-way adapter and open the clamp on the relevant branch and the clamp on the adapter-reservoir connection line.
- Keeping the three-way adapter branch for connection to the ECC circuit closed; proceed with recovery of blood lost during operation following the instructions above.

Priming fluid concentration during ECC (for cod. 712 cardio only)

- Close the 1/4" tube of the three-way adapter connected to the reservoir outlet connector with a clamp and open the clamps located on the three-way adapter branches. Remove the clamp used on the 1/4" tubing segment connected to the ECC circuit and open any others that may be present on the extension line.
- Whenever the extension line is directly connected to the wash set, ensure that all the clamps are open.
- Transfer the priming liquid to the Autotransfusion device following the instructions for use of the device, taking care not to create any negative pressure within the ECC circuit.
- During the return of the concentrated priming fluid, take particular care to avoid any air being routed to the ECC circuit.

Recovery of priming residue from the ECC circuit (for cod. 712 cardio only)

- If the BT 894 Autotransfusion Reservoir is used, it is advisable to remove the priming residue from the ECC circuit to the reservoir to keep the time the pump console remains in the operating theatre to a minimum.
- Close the aspiration and anticoagulation line originating in the operative field and keep both the clamps on the three-way adapter branches open. Open the clamp on the 1/4" tube segment connected to the ECC circuit.
- Transfer the priming residue to the reservoir from the ECC circuit until it is completely empty, appropriately adjusting the vacuum.
- After transferring the priming residue to the reservoir, reopen the aspiration and anticoagulation line, close the clamp on the three-way adapter branch connected to the ECC circuit and remove the connection to the circuit. The pump console is now free of connections and the washing process can start as normal.

⚠ CAUTION

- Wash the recovered blood with at least the volume of saline solution recommended in the instructions for use of SORIN GROUP ITALIA cell separators for autotransfusion, to guarantee removal of excess anticoagulant and other undesired components collected.

REINFUSING SALVAGED BLOOD

⚠ WARNING

- It is recommended that intraoperative and postoperative autologous blood be processed. The washing operation is necessary whenever the recovery is carried out in the presence of heavy haemolysis caused by surgery, or in cases of high coagulopathy or hepatic or renal failure.
- The decision to re-infuse the blood components obtained by this concentration-washing procedure is the responsibility of the attending physician.
- Do not reinfuse the patient's blood from the primary reinfusion bag when it is connected to the cell separator's wash circuit. Reinfusion from the primary reinfusion bag connected to the circuit could lead to air embolism to the patient.
- Do not reinfuse under pressure (i.e., do not use a blood pressure cuff on the secondary reinfusion bag). Reinfusion under pressure could lead to air embolism.
- To minimize the complications of particulate matter infusion, use of an in-line microaggregate filter on the patient reinfusion line is highly recommended. Failure to use an in-line filter may result in reinfusing particles potentially harmful to the patient.
- To reduce risk of air embolism, remove all air from the primary reinfusion bag before handing the bag over for reinfusion.
- The American Association of Blood Banks recommends the following guidelines for expiration of salvaged blood¹.
 - If not transfused immediately, units collected and processed under sterile conditions with a device for intraoperative blood collection that washes with 0.9% saline, USP, shall be stored under one of the following conditions prior to completion of transfusion:
 - At room temperature for up to 4 hours after completion of processing
 - at 1 °C to 6 °C for up to 24 hours, provided that storage at 1°C to 6 °C is begun within 4 hours after completion of processing
 - The transfusion of shed blood collected under postoperative or posttraumatic conditions shall be completed within 6 hours of initiating the collection.

¹ American Association of Blood Banks. *Standards for Perioperative Autologous Blood Collection and Administration*. 3rd ed. Bethesda, MD: American Association of Blood Banks, 2007.

USING THE BT894 ATS RESERVOIR QUICK REINFUSION PORT

The BT894 Autotransfusion Reservoir is provided with a pierceable port on the bottom, to gravity drain blood out of the reservoir. This port can be accessed with a spike connector and may be used to empty the reservoir in case needed, or to quickly reinfuse the unprocessed salvaged blood to the patient, in emergency situations.

Observe the following instructions to access the Quick Reinfusion Port:

- Shut off the vacuum source, pull the yellow tab on the port of BT 894 cover marked "PULL TO RESET", then raise the reservoir as high as needed
- Remove the cap on the Quick Reinfusion Port on the bottom of the reservoir, and insert a spike through the integral membrane
- In case of quick reinfusion, completely prime the line before connection to the patient

⚠ WARNING

- Limit the reinfusion of unprocessed blood to the surgical situations of real necessity. The decision and responsibility of reinfusing unprocessed blood belongs solely to the physician in charge. When reinfusing unprocessed blood, the use of a microaggregate filter on the patient reinfusion line is mandatory.

MEDICAL DEVICES TO BE USED WITH THE PRODUCTS

- The BT894 Holder (catalogue 04108) stand alone must be used to set up the reservoir separated from the ELECTA equipment.
- The Collection Set Code 712 Cardio and Code 712 must be connected to the SORIN GROUP ITALIA Wash Sets for removal of undesired substances collected from the surgical field.
- The connections to the product must be of a diameter which is compatible with the dimensions of the connectors on the device (3/8" or 1/4").
- The aspiration must be performed with SORIN GROUP ITALIA surgical aspirator or other aspirators having compatible technical characteristics.

⚠ CAUTION

- **The User should observe the warnings and cautions and follow Instruction for Use accompanying the separate device.**

RETURN OF USED PRODUCTS

Should the user be dissatisfied with anything related to the quality of the product, the product distributor or the authorized local SORIN GROUP ITALIA representative should be notified.

All parameters considered critical by the user must be reported with particular care and urgency. The following is the minimum information that should be provided:

- Detailed description of the event and, if pertinent, the conditions of the patient;
- Identification of the product involved;
- Lot number of the product involved;
- Availability of the product involved;
- All the indications the user considers useful in order to understand the origin of the elements of dissatisfaction.

SORIN GROUP ITALIA reserves the right to authorize, if necessary, recall of the product involved in the notification for assessment. If the product to be returned is contaminated, it must be treated, packed and handled in conformity with the provisions of the legislation in force in the country where the product was used.

⚠ CAUTION

- **It is the responsibility of the health care institution to adequately prepare and identify the product for return shipment. Do not return products that have been exposed to blood borne infectious diseases.**

ONLY for US customers

If for any reason the product must be returned to the manufacturer, a returned good authorisation (RGA) number is required from Sorin Group USA, Inc. prior to shipping. If the product has been in contact with blood or blood fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton to prevent damage during shipment, and it should be properly labelled with an RGA number and an indication of the biohazardous nature of the content in the shipment.

The shipping address for returned goods in the US is:

SORIN GROUP USA, Inc.

Returned CV Products
14401 West 65th Way
Arvada, CO 80004-3599
FAX (800) 323 4031

LIMITED WARRANTY

This Limited Warranty is in addition to any statutory rights of the Purchaser pursuant to applicable law.

SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected. SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or

any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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All these devices are single use, non toxic, non pyrogenic, supplied **STERILE** and packaged singly. Sterilized by ethylene oxide. The level of ethylene oxide residuals in the device is within the limits established by national regulations in the country of use.

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The Collection Set Code 712 Cardio is provided with the Cardio Kit, which includes a set of adapters for connection of the extra-corporeal circulation (ECC) circuit to the wash set of an autotransfusion equipment or to the BT 894 Autotransfusion Reservoir. The adapters contained in the Kit Cardio allow to extend use of the autotransfusion techniques to the following procedures:

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⚠ CAUTION

- For accurate measuring of the blood volume inside the autotransfusion reservoir, the reservoir arm raised a few cm.

DISPOSABLE SET UP

⚠ WARNING

- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the device if sterility cannot be guaranteed.
- Check the expiry date on the label attached. Do not use the device after the date shown.
- The device must be used immediately after opening the sterile packaging.
- The device must be handled aseptically.
- Carry out a visual inspection and carefully check the device before use. Transport and/or storage conditions other than those prescribed may have caused damage to the device.
- To prevent damage to the device, do not use solvents such as alcohol, ether, acetone, etc.
- Make the connections using a suitable aseptic technique.
- Due to the possibility of operator exposure to bloodborne pathogens (such as HIV, hepatitis viruses, bacteria, Cytomegalovirus, etc.) when handling extracorporeal blood circuits, adequate precautions should be taken at all times to prevent the exposure to and transmission of such agents.
- Carefully observe this set for leaks before and during use. Leakage may result in loss of sterility or loss of blood and/or fluid. If leakage is observed before or during use, replace leaking component or retighten leaking connection, as appropriate.

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⚠ WARNING

- SORIN GROUP ITALIA recommends not to exceed a negative pressure of -100 mmHg (13 Kpa/0,13 Bar/1,9 psi) during recovery procedure as an excess vacuum increases hemolysis of the recovered blood
1. Adjust the aspiration pressure according to the vacuum source available:
 - 4.1. If the Electa vacuum pump is going to be used, adjust the aspiration pressure to -100 mmHg (13 KPa/ 0.13 bar/1.9 psi).
 - 4.2. If an external vacuum source is going to be used, close all lines (the aspiration line and the tree-way adapter line, if present s), apart from the vacuum line, with a clamp and adjust the aspiration pressure to -100 mmHg (13 KPa/ 0.13 bar/1.9 psi). Reopen the lines when the adjustment has been made, keeping the vacuum on.
 2. Open the aspiration and anticoagulation line and transfer at least 200 ml of solution to the reservoir before starting recovery. The anticoagulant solution containing heparin, ACD-A or CPDA-1, must be prepared beforehand in a quantity sufficient for the duration of the procedure. If heparin is used, 30.000-60.000 I.U. of heparin per litre of sterile isotonic saline are sufficient to guarantee good anticoagulation of recovered blood.

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- If the BT 894 Autotransfusion Reservoir is used, it is advisable to remove the priming residue from the ECC circuit to the reservoir to keep the time the pump console remains in the operating theatre to a minimum.

- Close the aspiration and anticoagulation line originating in the operative field and keep both the clamps on the three-way adapter branches open. Open the clamp on the 1/4" tube segment connected to the ECC circuit.
- Transfer the priming residue to the reservoir from the ECC circuit until it is completely empty, appropriately adjusting the vacuum.
- After transferring the priming residue to the reservoir, reopen the aspiration and anticoagulation line, close the clamp on the three-way adapter branch connected to the ECC circuit and remove the connection to the circuit. The pump console is now free of connections and the washing process can start as normal.

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REINFUSING SALVAGED BLOOD

WARNING

- It is recommended that intraoperative and postoperative autologous blood be processed. The washing operation is necessary whenever the recovery is carried out in the presence of heavy haemolysis caused by surgery, or in cases of high coagulopathy or hepatic or renal failure.
- The decision to re-infuse the blood components obtained by this concentration-washing procedure is the responsibility of the attending physician.
- Do not reinfuse the patient's blood from the primary reinfusion bag when it is connected to the cell separator's wash circuit. Reinfusion from the primary reinfusion bag connected to the circuit could lead to air embolism to the patient.
- Do not reinfuse under pressure (i.e., do not use a blood pressure cuff on the secondary reinfusion bag). Reinfusion under pressure could lead to air embolism.
- To minimize the complications of particulate matter infusion, use of an in-line microaggregate filter on the patient reinfusion line is highly recommended. Failure to use an in-line filter may result in reinfusing particles potentially harmful to the patient.
- To reduce risk of air embolism, remove all air from the primary reinfusion bag before handing the bag over for reinfusion.
- The American Association of Blood Banks recommends the following guidelines for expiration of salvaged blood¹.
 - If not transfused immediately, units collected and processed under sterile conditions with a device for intraoperative blood collection that washes with 0.9% saline, USP, shall be stored under one of the following conditions prior to completion of transfusion:
 - At room temperature for up to 4 hours after completion of processing
 - at 1 °C to 6 °C for up to 24 hours, provided that storage at 1 °C to 6 °C is begun within 4 hours after completion of processing
 - The transfusion of shed blood collected under postoperative or posttraumatic conditions shall be completed within 6 hours of initiating the collection.

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SORIN GROUP ITALIA warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. SORIN GROUP ITALIA warrants that the medical device is capable of functioning as indicated in the current instructions for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, SORIN GROUP ITALIA cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performance and effectiveness of the device with damaging consequences for the patient, even though the specified instructions for use have been respected. SORIN GROUP ITALIA, whilst emphasizing the need to adhere strictly to the instructions for use and to adopt all the precautions necessary for the correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. SORIN GROUP ITALIA undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by SORIN GROUP ITALIA up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of SORIN GROUP ITALIA or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device except as expressly stated herein. SORIN GROUP ITALIA disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this Limited Warranty and in particular agrees, in the event of a dispute or litigation with SORIN GROUP ITALIA, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty is given as well as every dispute related to it or in any way connected to it as well as anything related to it or any dispute concerning this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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