<u>AGScan</u> USERS MANUAL

Version 0.1 / February 2006 Edition

Rémi Cathelin



Copyright © 2006 SIGENAE / INRA

Permission is granted to copy, distribute and/or modify this document under the terms of the GNU Free Documentation License, Version 1.2 or any later version published by the Free Software Foundation; with the Invariant Sections being just "Background", no Front-Cover Texts, and no Back-Cover Texts. A copy of the license is included in the section entitled "GNU Free Documentation License" at the end of this manual.



Revision table :

Version	Date	Major modifications	Author(s)
0.1	February 15th 2006	First release	R. Cathelin

Acknowledgment :

I wrote this guide in English having the aim to give access to AGScan to the largest possible scientific community. But English is not my mother language, so any remark about the orthography or the grammar is welcome.



INDEX

I. INTRODUCTION	4
I.1. A bit of history	
I.2. Main features of AGScan	
I.3. Vocabulary	4
II. GENERAL OVERVIEW	6
II 1 The process	6
II 2 FUE FORMATS	7
	7
II 4 INSTALLING AGSCAN ON YOUR COMPLITER	8
II 4 a Download	8
II 4 h packages contents	10
II.4.c. Launching AGScan	
III. PROCESSING STEPS	
	11
III.1. USER INTERFACE	
III.2. LOADING AN IMAGE FILE	
III.5. CREATING A GRID.	
III.4. ALIGNING A GRID ON AN IMAGE	
III.4 h. alianing manually a grid	
III.4.c. Global alignment algorithm	
III.4 d Local alianment algorithm	23
III 4 e Spot detection	25 26
III 5 THE TABLE	20 26
III 6 Oliantieving an Alignment	20
III 6 a quantifications	27
III 6 h analyze your results	29
III 6 c Export results	31
III.7. PROCESSING A BATCH.	32
III.7.a. Batch tab	32
III.7.b. Preparing the batch	
III.7.c.Save/Open a batch	
IV. PLUGINS	
IV 1 Presentation	37
IV.1. I RESENTATION. IV.2. Example: the "Memory Plugin"	
V. FREQUENTLY ASKED QUESTIONS	
VI. INDEX	40
VI 1 Word Index	۸۵
VI.2. Illustration index	
VII. GNU GENERAL PUBLIC LICENSE	44



I. Introduction

This manual is aimed to provide the necessary information to biologists using AGScan image processing software. The document is split in four parts :

- Introduction
- General overview
- Processing steps
- Plugins presentation
- Frequently asked questions

I.1. A bit of history

AGScan is a micro-array image processing software developed by the SIGENAE team. http://www.sigenae.org/

It is an evolution of BZScan : BZScan is a radioactive DNA micro-array image processing software based on the JAI library . It was developed by the TAGC (<u>http://tagc.univ-mrs.fr/</u>) from 2002 to 2005.

In 2005, the SIGENAE team decided that it was worthwhile to go on with the development of BZScan, and mainly by going through the new Sigenae version named AGScan.

This software is freely available using the following address: http://mulcyber.toulouse.inra.fr/projects/agscan/

You can also download the source code, report bugs and propose enhancements at the same address.

I.2. Main features of AGScan

The main features of AGScan are :

- 1. Grid conception
- 2. Able to use single or multiple (up to 3) channel images as input
- 3. Image manipulation and processing based on ImageJ
- 4. Manual and automatic global alignment
- 5. Automatic local alignment
- 6. Automatic quantification
- 7. Batch processing
- 8. Data export
- 9. Plugin support, to add new file formats, alignment or quantification methods...
- 10. Multi-language support

These different features will be presented in the rest of this document.

I.3. Vocabulary

Image : is a file produced by a micro-array scanner and sometimes modified (cut,...).

Grid : is first a file generated by AGScan describing the way the array has been spotted. From this it can be aligned on an image and be a part of an alignment. You will see later the importance of the grid structure.

Alignment : is composed of an image and a grid. It can have several states like 'not aligned', 'globally aligned' and 'locally aligned'.



Quantification : once the alignment is finished the quantification is the process of extracting one or more figures for each spot giving an intensity value for the spot.

Batch : is a group of alignments processed by AGScan from global alignment to quantification through local alignment.



II. General Overview

This chapter presents an overview of the whole process, some information about image file types and quality and gives the software install procedure.

$\tt II.1.$ The process

To start the process you need at least :

- a micro-array image (different formats can be used),
- a grid file (generated in AGScan see next chapter section),
- AGScan installed (see install procedure below)

Notice that AGScan works with 16bits-unsigned grayscale images.

By default, AGScan input image format is the TIFF format but ISAC format (INF/IMG) is also recognized by a plugin automatically provided into all AGScan versions.

The grid file can be generated with AGScan or provided by someone having generated it. To share a grid you have to be sure that the micro-array and the definition used to scan the image are the same.

The process can be split in five steps :

- experiment image(s) loading
- grid loading

- global alignment : This step can be done by AGScan or by the user. It corresponds to the general alignment of the designed grid on the image(s).

- **local alignment :** This step is done by AGScan. Starting from the global alignment AGScan tries to figure out what the small local changes are to fit the sub-grids as well as possible to the image.

- quantification : This step is done by AGScan. Starting from the local alignment coordinates AGScan calculates for each spot several quantification values. These values will be presented later on.



As presented before, AGScan is able to process an image and a grid without any human help. This is a very important and user asked AGScan feature.



II.2. File formats

AGScan can load and save several file formats depending on the element. See table hereunder.

Element	file extension	comments
image ¹	.tif	Standard format of compression of image without loss of information,
	.inf	This file is the information file of the corresponding ".img" file. Both files have to be in the same directory to be used by AGScan. This file type is named ISAC and is for example produced by the Fuji BAS5000 scanner.
grid	.grd ²	This file type is produced and used by AGScan. The format is kind of XML compliant.
experiment: Alignment / quantification	.zaf	This file type is produced and used by AGScan. It's an archive file that contains a ".ali" file (experiment data file), a manifest file and may contains the image(s), The ".ali" format is kind of XML compliant. ₃
batch	.bzb ⁴	This file type is produced and used by AGScan. The format is kind of XML compliant.

Notice that AGScan provides also non-openable export formats like a ".txt" output, easily usable with your favorite spreadsheet.

II.3. Image quality

Even if AGScan proposes methods for automatically align grids on image, the result depends of a lot of things like the quality of the image, the grid structure and parameters used for the alignment.

¹⁻ Be aware that the image files have to be a 16 bits format to be processed by AGScan.

²⁻ Be careful: BZScan and AGScan .grd grid files and .bzb batch files are incompatible!

³⁻ Never try to open a .zaf file: .ali files are not alone support

⁴⁻ Be careful: BZScan and AGScan .grd grid files and .bzb batch files are incompatible!



II.4. Installing AGScan on your Computer

II.4.a. Download

AGScan is freely available at:

http://mulcyber.toulouse.inra.fr/projects/agscan

•••	🛛 🏠 🖸	a http://mulcyber	r.toulouse.inra.fi	/projects/agscan/					v 0	ок 💽	
Démarrage 🔂 De	ernières nouvelles										
	INRA	B	ite Toulo	use P	iule	ybe	Search th	ne entire project	*		Sea
Home		My Page		Project Tree		Code Snippets	Pro	oject Openings		AGSca	an M
Summary	Admin	Forums	Tracker	Lists	Tasks	Docs	Surveys	News	SCM	5	Files
AGScan is a 1 purposes fluc develop func • Develop • Intende • License • Natural • Operat • Prograr • Topic: Registered: 2 Activity Perc View project	micro-array in micro-array tionalities). ad Audience: :: GNU Gener Language: E ng System: (noming Langua Bio-Informati 005-11-16 0 entile: 100% activity stati	mage processi images suppo s: 4 - Beta Developers, E al Public Licen English, French Windows, Linu age: Java cs 9:21 istics.	ng software rt, multi-lan ind Users/De ise (GPL) 1, German ix	. It is an evol guage menus esktop, Other	lution of BZ It is also Audience	Scan2. It is ba	used on Imag	eJ and I (user can Ré De Ch [V [R	Devel oject Adm mi CATHE ivelopers: ristophe k iew Memb equest to	oper In iins: LIN (lopp iers] join]	fo [~]

Just click on the "Files" subsection:

Backago	Roloaco & Notoc	Filonamo		Date	
Package	Release & Notes	Filename	Size	D/L Arch	Туре
AGScan - Full 🖾					
v270106			2006	-01-27 14:35	
AGScanFull270106.zip			19.5 MB	0 i386	.zip
v081205			2005	-12-08 00:00	
AGScanFull081205.zip			19.48 MB	0 i386	.zip
v061205			2005	-12-06 00:00	
AGScanFull061205.zip			19.48 MB	0 i386	.zip
v211105			2005	-11-21 00:00	
AGScanFull211105.zip			19.49 MB	0 i386	.zip
AGScan - Light 🖾					
v270106			2006	-01-27 14:38	
AGScanLight270106.zip			3.3 MB	0 i386	.zip
v081205			2005	-12-08 00:00	
AGScanLight081205.zip			3.28 MB	1 i386	.zip
v061205			2005	-12-06 00:00	
AGScanLight061205.zip			3.28 MB	0 i386	.zip
v211105			2005	-11-21 00:00	12000
AGScanLight211105.zip			3.29 MB	0 i386	.zip
Data examples 🖾					
Radio grids			2005	-11-25 16:15	
GRID_AGScan_example_8x4.qr	d 🖌		6 KB	25 (386)	Other
GRID_AGScan_generic_16x4.gr	d		6 KB	23 i386	Other

Mulcyber AGScan page



All versions of the software are available. They are named using the date of the release.

The software releases are available in two different packages: light and full versions. With the light version, Java Runtime Environment needs to be installed previously on your operating system. The full version does not need any other software, just to click on the icon!

Downloading AGScan is very easy.

First choose a release of the package. For each release, notes are available by clicking on the release name.

Next you can download the corresponding "AGScanxxxxx.zip" file.



Download AGScan



II.4.b. packages contents

The "AGScanFullxxxxx.zip" archive file contains the following directory:



AGScanFull270106

The "AGScanLightxxxxx.zip" archive file contains the following directory:

Package contents are similar except 3 elements added to the Full version:

- the install directory: contains a Java Runtime Environment
 - the "AGScan.cfg" file: the launcher configuration file
 - the "AGScan.exe" file : the icon to launch AGScan (on Windows OS)

Note: the Full version launcher is only usable for Windows OS.

The file structure of files is composed of:

- the "AGScanxxxxx.jar" file is the java executable archive file
- the languages directory containing language files available in AGScan
- the "lib" directory containing AGScan libraries and dependences used by the software (ij.jar for example). Notice that the subdirectory called "plugins" of this directory contains imageJ plugins that can be used by AGScan
- the "plugins" directory containing AGScan ".class" plugins files
- the "ChangeLog.txt" file containing the history of release notes
- the "parameters.properties" file containing the current user preferences of the software

II.4.c. Launching AGScan

On Linux OS, use the command line: java -Xss10m -Xmx420m -jar AGScanxxxxx.jar

On Windows OS:

If you choose a full version: just click on the AGScan.exe icon! You can also run the command line (or write it in a .bat file): "java.exe" -Xss10m -Xmx420m -jar "AGScanxxxxx.jar"

The ".jar" file name is the release name of AGScan.

Note: the "-Xss" and the "-Xmx" parameters correspond to the memory allocation given to the Java Virtual Machine(JVM) in order to support the application. If you encounter "java OutOfMemory errors" then increase the "-Xmx" value.

By default, AGScan full version is configured with 420Mo of memory allocated for the JVM.You can modify this value in the launcher configuration file ("AGScan.cfg"file).



III. Processing steps

This guided tour will walk us through the different parts listed hereunder :

- User interface,
- Loading an image file,
- Creating a grid,
- Aligning a grid on an image,
- The table,
- Quantifying an alignment,
- Analyzing quantification results,
- Processing a batch.

First it will start by presenting the user interface layout:

III.1. User interface



The user interface of the software wants to be simple, readable and as ergonomic as possible. It is composed of five elements:

- The "title window": presenting the used AGScan version, in the selected language.
- The "menu bar": containing the different menus of the application.
- The "toolbar": icons wich are shortcuts to some sub-menu actions.
- The "display window": the main window containing experiment view. Several elements may be displayed in the same time (images, grids, alignments). They are open in different tabs.
- The "status bar": displays current information like status and progress bar.

As we want the software to be simple to use we tried to place the menus in a logical order from left to right. We hope that menus are explicit enough to be easily memorized. Shortcuts to common actions have been added below the menus in order to open, align and quantify an image in three clicks! Shortcut icons are also duplicated next to each menu title to ease memorization. Notice that if you point your mouse on a shortcut, you will obtain a little description of the associated action.

File	Edition	Disp	lay	Image	Grid	A
0	Dpen	Þ		Image		
C	Ireate	•		Grid		-
	ave			Alignmen	t	
	ave as			Batch		



Example : the "Open Image" action menu (and related toolbar item)

The first menu you can use is the "Edition>Language" one. It allows to change the language of the application menus.

🖢 Langi	iage		×
Choice -			
	E	nglish	
	Fi	rench	
	G	erman	
Validate			
	ОК	Cancel	

Notice that you need to reload the application to apply a new language.

Note: in the current version, three languages are available: English, French and German.

There can be spelling mistakes. However it is easy to modify these faults. Indeed, all the words are stored in text files (named "messages_en.properties", "messages_fr.properties" and "messages_de.properties"). We can even easily generate new ones in other languages. For that, it is necessary to add the new "messages_xx.properties" file into the subdirectory of AGScan named "languages". It will be automatically recognized by the software.

Now Let us present the different steps needed to quantify an image.

III.2. Loading an image file

<u>Preliminary remark:</u> an experiment can be composed from 1 to 3 images, but in the following chapters we will refer to an experiment including only one image.

To open an image, use the "File>Open>Image" menu or the following short icon of the toolbar:









The image is loaded into the "Display window".

In the case of a multi-channels experiment, the RGB image is displayed :



Once opened, you can now modify it using the different widgets disposed in the "Image" menu.



They allow to invert the LUT, to modify the contrast and the brightness, to rotate the image, to flip it horizontally or vertically. Some of these options are very interesting when you have scanned your array the wrong way. You can also select a rectangular region and crop it to create a new sub-image. Sometimes due to low contrast or quality, the spots can be quiet difficult to locate. All these image tools can help you to visualize them. However, you have to know that the modification of the contrast and the brightness of the image doesn't affect pixel values of the image and isn't save.

Other general functionalities can be used in the "Display" menu.

An example: the "zoom":

Working with images, we sometimes need a general view and sometimes a very precise one. This can be obtained by using the zoom . Clicking on the left of the progress bar will zoom out of the image, clicking in the right will zoom in the image.

When the zoom is active, the icon background becomes gray. By clicking a second time on the icon, the zoom is disable.



Example: image with zoom action enabled.

Display Image Grid Align	The "display" tools are :
• Tabs	- "tabs" or "windows": the display mode. Each item (image, grid
Windows	alignment) is on a tab or on a distinct window, I abs are the defau
Q Zoom	display mode.
	- "zoom": disable/enable the zoom,
	- "zoom window". Allows to display a little window containing the
Display the ruler	entire image, the zoom factor applied and that show the part of the
✓ Display the state bar	image zoomed
Properties	- "Display the ruler": disable/enable a ruler. Disabled by default.
	- "Display the State bar" disable/enable the status bar Enable
	by default
	by default.
	- Properties : allows to display the properties (name, size) of the
	current element displayed into the current tab or window.

By default, the display mode is "tabs", the zoom isn't active, the ruler is disabled and the status bar is displayed. (see the image above).



On this example, the window mode is active, the status bar is hidden, the ruler is displayed and the zoom window is open (the current zoom factor is 3),

The next part will show you how to create a grid.

III.3. Creating a grid

The grid is a very important element for the image analyze. After opening an image, you have to load a grid on it and therefore you must have already define it.

Note that you can use the same grid for all your experiments that use the same spotting scheme.

So this part will show you how to create, open and modify a grid.

At first, choose the creation grid menu found in the "File>Create>New Grid" menu or in the "Grid>New Grid" menu.

A new tab and the "structure editor" will appear. The main window is divided into two parts:

- the **grid**: the graphical representation of the grid
- the table: a table with all the characteristics of the grid spots (number, name position). It will be describe later.

|--|

	GScan -	January	27, 20	005	Vers	ion		×
File	Edition	Display	Image	G	irid	Alignme	nt Quantification Batch Tools About	
				9	1	Q		
GRJ	D_113897	4828336.g	rd \					
	0000	0000		2			Structure editor	•
	Name	X	Y	R	С	Block	Number of levels : 1	
1	A01	187.5	187.5	1	1			
2	B01	562.5	187.5	1	2			2000
3	C01	937.5	187.5	1	з		Number of columns : 12	
4	D01	1312.5	187.5	i	4			
5	E01	1687.5	187.5	i	5		Height of lines (µm) : 375	
6	F01	2062.5	187.5	1	6		spots diameter(µm) : 250	
7	G01	2437.5	187.5	1	7		Starting position - Top Left	
8	H01	2812.5	187.5	1	8		Soot pames :	
9	I01	3187.5	187.5	1	9			
10	J01	3562.5	187.5	1	10		Close	
11	K01	3937.5	187.5	1	11	1	1 11 250.0	
Read	ly .						1x	0

By default, the editor appears on a "one level" grid, So we have to introduce you the notion of **level** in the grid. A grid can have two levels representing the block and the grid.

- level1 = the block level. In fact, the spotting is usually made using a block pattern. This is usual but
 not always the case. The level defines a block of spots. Information needed is the dimension of a
 block (number and size of columns and lines) and the diameter of a spot.
- level2 = the grid level. A grid can often be compared to a set of blocks. So if you have this second level, you just need to inform how many blocks you have vertically and horizontally. So you can also notify inter-blocks space size.

To sum up, a grid is seen by AGScan as a set of elements. The smallest element is the spot. When you have two levels, a set of spots forms a block and a set of block a grid .

In the next table the first line gives you number of levels used, the second gives you the structure of the first level and the third line gives you the structure of the second level. In the last line, you can see the graphical view of the type of spotting pattern you have used.

1 level	2 levels	2 levels
20 x 12 spots	12 x 12 spots	1 x 1 spots
-	4 x 2 blocks	10 x 10 blocks
	Examples of grids	

The grid, the table and the editor are linked. So, be aware that you need to move on an other box or hit <ENTER> (carriage return) after having entered a value in each text box of the "grid structure editor" in order to synchronize the view of the grid .

You have now define your grid and you can save it. The default proposed name of the saved file has no sense. It is just created randomly in order to be unique.

The resulting file has a ".grd" extension and is in XML format.

Notice that you cannot create a grid directly on a loaded image. You need before to load the image and after the grid above.

III.4. Aligning a grid on an image

III.4.a. Loading a grid

This is where all really begins. After having opened an image, you can load a grid on it. To do this, you can use the "Grid>Load" menu or the following icon :

The selected grid will be loaded on the image, placed on the top left and the tab is divided into two parts, in order to add the table of spots characteristics.

00016	States and the second	A DEAL PROPERTY	and the second se		-	Second International Internati					
24			2.62	12.83	100	-	8	22-32 TO 12-58	2.2.2.2.2.8 B.2.2.2		
											-
											-
	Name	x	Y	R	C	Block	Row	Column	Diameter	Computed Diameter	-
1	Name A01-A01	X 173.0	Y 173.0	R 1	С 1	Block	Row	Column	Diameter 300.0	Computed Diameter -1.0	
1	Name A01-A01 A01-B01	X 173.0 519.0	Y 173.0 173.0	R 1	C 1 2	Block	Row 1	Column 1 2	Diameter 300.0 300.0	Computed Diameter -1.0 -1.0	-
1 2 3	Name A01-A01 A01-B01 A01-C01	X 173.0 519.0 865.0	Y 173.0 173.0 173.0	R 1 1	C 1 2 3	Block 1 1 1	Row 1 1	Column 1 2 3	Diameter 300.0 300.0 300.0	Computed Diameter -1.0 -1.0 -1.0	
1 2 3 4	Name A01-A01 A01-B01 A01-C01 A01-D01	X 173.0 519.0 865.0 1211.0	Y 173.0 173.0 173.0 173.0 173.0	R 1 1 1	C 1 2 3 4	Block 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	Row 1 1 1	Column 1 2 3 4	Diameter 300.0 300.0 300.0 300.0	Computed Diameter -1.0 -1.0 -1.0 -1.0 -1.0	

The next step consists in correctly aligning the grid on the image.

You have several possibilities to do it. This step can be done manually or automatically by launching alignments algorithms.

All the alignment functionalities are available in the "Alignment" menu:



III.4.b. aligning manually a grid

How to manually place the grid?

You can either use the mouse and the "Move/Rotate/Resize" menu.

Let us present to you the "color code" of a grid: a loaded grid on an image appears in blue if it is not selected and orange if you select it.

To select in the grid, just click on it.

Note that you can easily only move a sub-grid by changing the selection mode:



These three icons of the toolbar allow to change the grid selection mode. By default, the first is enabled. It corresponds to the entire grid. If you select the second icon, you turn into the block mode



that allows you to only select a block. The third icon corresponds to the spot mode that allows to only select one spot.

A selected item (spot, block or the grid) can be moved using the mouse. For this, move it by maintaining it selected (left mouse button down).

Notice that the table view corresponding selected spots lines are yellow.



You can also deform the grid by placing your mouse at the corners and drag them.

A specific and useful tool allow to easily manually align a grid :the "Move/Rotate/Resize" menu. Click on it: a window appears, offering you the possibility to rotate, move and resize the grid of given values:



Example: the grid has been moved and rotated.



You can choose the unit of your rotation (degrees, radian), the unit of your translation or resizing (pixels, µm or spots).

Thanks to this tool, you can manually correctly place the grid on the image. In order to verify the position of each spot, you can apply a "spot detection" (presented below).

Note that it is possible to perform an automatic local alignment on a manually aligned grid.

It is interesting to know how manually place a grid but in the most cases, you can perform automatic global and local alignments.

In the next pages we will see how to perform these automatic alignments.

Let us first introduce the notions of "global" and "local" alignments. The global alignment places the grid (as good as possible) on the entire image. The local alignment searches to correctly align each block separately.

These two algorithms have been initially developed for the BZScan software. So they have a lot of parameters that can be modified and memorized.

III.4.c. Global alignment algorithm



The "Global alignment" sub-menu contains two different algorithms, the "sniffer" and the "edges detection" global algorithms. The default algorithm is the "sniffer" one (in bold in the menu). You can also call it thanks to the following toolbar icon :

If you click on "Default algorithm", it is launched without proposing you to change the parameters. To modify them, you need to choose the "Sniffer" menu. So the parameters windows appears:



	E Farin	1		
	t : [0.1] c:: [0.7	2 <u></u>	
	spot detection Alg	orithm : Detection	by 4 profiles 👻	
	Research are	a of the center (%)	: 40	
	Computed are	a of the center (%)	: 10	
		Quality threshold	: 30	
	1			
ocal alignment a	lgorithm paramete	rs		
Bootstrap :	sensitiv	vity (%) : 20	Iterations number :	5
1st run :	sensitiv	vity (%) : 30	Iterations number :	50
2nd run :	sensitiv	vity (%) : 20	Iterations number :	10
3rd run :	sensitiv	vity (%) : 10	Iterations number :	10
Thre	shold of change of se	nsitivity : 1.0E-7		
Vi2 Vsrifier la co	bhiż 1⁄3rence de la grille	a		
		1	0	
rementation for th	ne research of the an	igle (degrees) : 0	.2 📜	
Contras	t for the research of I	the angle (%) : 5	50 1	

Notice that all these parameters can be changed and saved. It is also always possible to restore the default parameters. So don't hesitate to change them. For more information about their meaning, please contact us.

"Sniffer" algorithm quick presentation:

At first, the angle of the image is calculated. Then (about 10% of the process), the algorithm try to parse all the image with a sniffer block to find all spots. To do this, it already searches to place the "sniffer block" (a defined block of spots) at 40% in width and in height of the image, the condition to place it is that the center spot must be a detected one. When the algorithm detects this spot, it considers this position as the first position of its route (pict1). The sniffer block moves after in all the direction of the image and flags every spot it finds. Each time it will find something considered a spot, it will put a little colored cross on its center (pict2). After having recovered the entire image, the algorithm computes how to place the grid in order to include the most of detected spots.

At the end of the process, you will have a yellow (selected) grid aligned on the image(pict3).



Some visual steps of the sniffer global alignment.



Our grid is now globally aligned.

The next step of the process is the local alignment.

III.4.d. Local alignment algorithm

A local alignment is an alignment that you must launch on a "already aligned" grid. It allows to improve the alignment quality.

In the actual version, two different local alignment algorithms are available.

Move/R	otate/Resize		0	♠ 83	88
Spots d Global a	etection alignment	•			
Local ali	ignment	•	Defaul	t algorithm	
			Alignm	ent of the e	ntire grid

The "Local alignment" sub-menu contains two different algorithms, the "entire grid" and the "block by block" local alignment algorithms. The default algorithm is the "block by block" one (in bold in the menu). You can also call it thanks to the following toolbar icon :

88

If you click on "Default algorithm", it is launched without proposing you to change the parameters. To modify them, you need to choose the "Alignment block by block" menu.

"Alignment block by block" algorithm quick presentation:

Each block is processed individually. The position of each spot is verified and adjusted. The block may be moved from the initial global alignment place.



Progression of the local algorithm

In this process, AGScan takes each block and tries to fit it as precisely as possible to the spots it finds underneath. Each block is processed separately and can therefore move from the previously calculated global position to a new one. It is quiet amazing to see how a block is progressing towards its end position. Once the block position is optimized, AGScan goes on with the next block.

The red circles indicate that the program identify a spot ; on the opposite, the green ones indicate that the spot is deduced due to the grid structure and the red spots locations.

Note that the coordinate values of the spots have changed during the local alignment process. You can see values in the alignment table.

The progress of the alignment is given by a progress bar. At the end of the process your display windows will look like this :







III.4.e. Spot detection

The spot detection is the method used to define if a spot is found (red) or estimate (green). A spot is recognize as one if its form fit the model.

You can call it on a selected a grid. Indeed, a selected grid becomes orange, so you don't see the spot nature.

About
0 83 8
Algorithm
s
fil

The "Spots detection" sub-menu contains two different algorithms, the "2 profiles" and the "4 profiles" algorithms. The default algorithm is the "4 profiles" one (in bold in the menu). You can also call it thanks to the toolbar icon :

88

If you click on "Default algorithm", it is launch without propose you to change parameters. To modify them, you need to choose the "4 profiles" menu.

The "2 profiles" and the "4 profiles" spot detection algorithms are similar. The difference is that profiles of spots used to search the center of the spot since the center of the research zone are diagonal lines (for the "2 profiles" detection) and diagonals + vertical + horizontal lines (for the "4 profiles" detection).



2 profiles

4 profiles

These two methods progress rapidly and give you the percentage of detected spot in the status bar:



$\tt III.5.$ The table

When a grid is open, the main window is split into two parts, giving you access to a table view of the spot informations.

This table is synchronized with the alignment. So coordinates have been updated by the move of the grid. However, no new columns are added.



· · · · ·		B	8899	ğč	ğ	5886	36	<u>soo</u>	2000
			8000	AC	ŚČ	ŎŎŎ	3 8	NOO Ç	<u>AZZZZ</u>
	Name	x	Y	R	С	Block	Row	Column	Diameter
86	A01-B08	3723.348	8476.989	8	2	1	8	2	300.0
87	A01-C08	4062.763	8455.912	8	3	4	8	3	300.0
88	A01-D08	4400.071	8391.421	8	4	1	8	4	300.0
		A selec	ted spot	and	its	table d	escrip	tion	

Base columns are:

- the first column is the number of the spot, determined by the spotting plan.
- Name: ID of each spot defined by AGScan: the first part is the block ID, the second part is the spot ID. The letter corresponds to the line position (A,B....) and the number is the column position (1,2,3,,,).
- X : X coordinate of the center of the spot in the image(in micrometer).
- Y : Y coordinate of the center of the spot in the image(in micrometer).
- R: row number in the entire grid.
- C: column number in the entire grid.
- Block: number of the current block.
- Row: row number in the current block.
- Column: column number in the current block.
- Diameter: constant diameter of the spot defined in the grid structure.

Note: the "Computed diameter" column is present but not used (values are "-1") until the quantification is processed.

The table and the alignment are interactive. You can select a spot in the table and see its location at the same time on the alignment: the spot (or the set of spots) are highlighted into the alignment (see the figure above).

So it's also possible to sort columns values, copy it, color some values...



The columns options can be active by clicking on the right click of the mouse. The interest of these options will be describe to you after the quantification part.

III.6. Quantifying an alignment

Now you have an aligned grid. So the next step is the quantification.

III.6.a. quantifications



Several computations are available. To launch a quantification, use the "Quantification>Compute" menu and select your action.

You can choose the type of quantification you want to call or select several of them with the "Selection" sub-menu or the following toolbar icon:



Compute	•	🔆 Selection
Export	×.	Image/Fixed diameter Quantification
Display computed diamete	ers	Image/Computed diameter Quantification
FG Total & Ratio		Fit/Fixed diameter Quantification
	100	Fit/Computed diameter Quantification
1 28 8000		Diameter
		Correction of overshining
1 P 25 CO		Correction by the fit
11 3000		QM

The selection menu opens the next window:

	Make your selection;
🖌 Image / C	onstant diameter Quantification
🖌 Image / V	ariable diameter Quantification
Fit / Cons	tant diameter Quantification
🖌 Fit / Varia	ble diametre Quantification
🖌 QM	
Fit Correc	tion
🖌 Overshini	ng correction
Select all	/ Deselect all
	OK Capcel

After selecting the actions you want to process, click on "OK".

Each calculation is done one after the other. You can follow advance thanks to the progress bar. Results are available into the table.

	Computed Diameter	Qtf. Fit/Const.	Qtf. Fit/Calc.	Qtf. Image/Const.	Qtf. Image/Calc.	QM	Fit Correction	Overshining Correction
Ē.	344.91243	18	22	33	42	0.6	0.0	0.0
	344.91243	21	27	32	46	0.6	0.0	0.0

Example of columns added, results of quantifications

Note: the "Computed diameter" column is filled during quantifications that need this value.

By selecting the data column title, you can change its position in the table.

Available Quantifications:

Four quantification are available in AGScan:

- Image/Constant diameter Quantification: for each spot, the quantification value is computed thanks to the image pixels intensities and the constant diameter given by the grid structure.
- Image/Variable diameter Quantification: for each spot, the quantification value is computed thanks to the image pixels intensities and a variable computed diameter.
- Fit/Constant diameter Quantification: the quantification value is computed thanks to the estimated intensity of the spot (with a fit function) and the constant diameter given by the grid structure.
- Fit/Variable diameter Quantification: the quantification is computed thanks to the estimated intensity of the spot (with a fit function) and its variable computed diameter.

The QM value is an indication of quality of the spots. Each spot has a value of QM ranging between 0 and 1 (the closer the value is to 1, the more it is reliable).

The "fit correction" and the "overshining correction" are two corrections values that can be used with quantification. Be careful, to apply them, you need a specific formula.

III.6.b. analyze your results

Thanks to the table columns options, you can sort the result values,

Example of a sorted column: click right on the column header and choose "Sort> Descending Order".

бшш	tage/Co *
	13759
	12058
	11290
	10656
	9888
	9840
	9275

An useful option is also the "color" tool. You can choose to give a color to some spots of the alignment in order to locate them.

<u>Example</u>: we want to color in orange spots that have an "image/constant" quantification value higher than 10000. To do that:

1/select the "Modify" menu of the "Qtf Image/Const" column.2/Check the "superior" case, fill the value (10000) and choose a color(orange).

Name .	Qtf. Image/Const.
Type :	Entier 👻
Default value :	-1
Values :	
if superior to 10000 If between	and

3/Click on "Apply" and "Close",

4/select the "Color spots" menu of the "Qtf Image/Const" column.

The text header of the column becomes red. Spots concerned by this condition will appear in orange on the aligned image.

In our example, only 4 spots are colored (as the sorted column shows us before).



Spot coloration example

You can also color spots between two values or lower than an other one.



III.6.c. Export results

You can save the experiment in every process step. An alignment and a quantification are considered as the same thing. In fact, there are two "aligned grid with a table". Only tables columns differ. So use the "save" menu to save it. The extension of an experiment file is ".zaf".

You can also export only the table values in order to use them into a spreadsheet. To do that, select the menu:

Quantification	Batch	Tools	Abo	ut
Compute	đ.		×.	
Export			Ð	📅 Text File

"Quantification>Export>Text File" menu

The file saved content looks likes:

Block	Row	Column	Diameter	Computed Diameter	Qtf. Fit/Const.	Qtf. Fit/Calc.
1	1	1	300.0	344.91243	18	22
1	1	2	300.0	344.91243	21	27
1	1	3	300.0	344.91243	48	60
1	1	4	300.0	344.91243	11	14
1	1	5	300.0	344.91243	67	82
1	1	6	300.0	346.0	4475	5307
1	1	7	300.0	271.35695	4654	4397
1	1	8	300.0	344.91243	120	140
1	1	9	300.0	344.91243	46	58
1	1	10	300.0	346.0	3617	4098
1	1	11	300.0	344.91243	37	44
1	1	12	300.0	344.91243	27	32
1	2	1	300.0	344.91243	24	29
1	2	2	300.0	344.91243	11	14
1	2	3	300.0	286.43438	342	330
1	2	4	300.0	344.91243	176	212
1	2	5	300.0	121.14123	935	191
1	2	6	300.0	295.18784	1720	1720
1	2	7	300.0	344.91243	125	146
1	2	8	300.0	344.91243	18	23
100		1 12	000.0	01101010		

Text file output example

Now, you have done a complete analyze of your image.

Do you know that this process can also be repeat automatically, using the batch mode?



III.7. Processing a batch

III.7.a. Batch tab

A batch allows to analyze automatically several experiments!

Important note: all these experiments are analyzed with the same grid.

To create a batch, you can use the "File>Create>New Batch" menu or the "Batch>New Batch". The associated icon of the toolbar is:

٥

The batch tab contains two parts:

💦 AGScan - January 27, 2005 Version		
File Edition Display Image Grid Alignmen	Quantification Batch Tools Ab	out
Image: Second		Right part
BATCH_1139324241609.bzb	ge 3 Status	Bight part Grid : Output directory : Save images with alignments Channels Number: 1 Add images Select al Job Images process Rotation : 0 +90° Pip : horizontal Vertical Vertical V Global Alignment V M Edges Detection Parameters V Incal Alignment Export results Entire Grid Parameters
		Launch Display the current alignment Erase
Ready	''E	atch" tab

The tab is split into two parts, the left part that displays experiments status and images open and the right part that contains the configuration of the batch.

Note that the left part is by default hidden. You have to increase it with the mouse to see it better.



III.7.b. Preparing the batch

Grid :	GRID_AGScan_example_8x4.grd	
Output directory :	:\Documents and Settings\rcathelin\test	2
	Save images with alignments	
	Channels Number: 1 1 Add images	Select all

Batch files configuration

In the top of the right part of the tab:

1/ Select the grid of the batch.

2/ Choose the output directory where results will be saved.

3/ Check or not the fact to save image into the archive created for each image. Notice that if the image isn't in the archive, you need to give its localization when you will open the experiment.

4/ Select the number of channels (1 to 3).

5/ Click on "Add Images" to select images.

Note: you can add images one by one or several in the same time in a multi-selection (only in "1 channel mode").

After choosing several images: the left part is filled:

	Image 1	Image 2	Image 3	Status
1	microarray_example.tif			WAITING
2	microarray_example2.tif			WAITING

Batch experiment information

You can now choose all the actions you want to process for each image thanks to the job right part.

Rotation : +90° ·90° ·180° Flip : horizontal vertical	 Image / Constant Diameter Quantification Image / Computed Diameter Quantification Fit / Constant Diameter Quantification Fit / Computed Diameter Quantification
Edges Detection Parameters Sniffer	QM Fit Correction Overshining Correction
Local Alignment Entire Grid Block per block	Export results Text

You just need to check or un-check process you want.

You can choose to display or not the current alignment processing by checking the last checkbox. To launch the batch, just click on "Launch".



😽 AGScan - January 27, 2	005 Version	Ř			
File Edition Display Imag	e Grid Alig	nment Quan	tification Batch T	ools About	
	1			•	88 ★ 🔕
BATCH_1139324241609.bzb \					
		1000	5		l 1980 - Star 🛓
and the second					+
		+ +	•		* * * *
	- -	+	• • •	S + +	+
		+ +	* + +		++
	+ +		* * *		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
Sela Station	++++		* +		•
Image 1	Image 2	Image 3	Status		cua 🦛
microarray_example_b			17%		Griu :
microarray_example2.tif			WAITING		Output directory :
	-	505055			
		99998		00/	
	Follo	wing a "Ba	tch" in progres	0 70 S	

When you choose to display the process, you can follow it into the same tab: indeed a new part appears in the top of the current tab. As for the left window, you have to stretch it with the mouse. The status of each image is updated.

In the statusbar of AGScan, the total percentage of the batch is given.

When the process is successfully finished, the image status is "DONE".

	Image 1	Image 2	Image 3	Status
1	microarray_example			DONE
2	microarray_example2			DONE

Batch progress finished

When the batch is finished, you can see results in the output directory.

For each experiment, you will have:

- the archive .zaf file
- the snapshot .tif file (rapid snapshot of the alignment)
- the export files (optional)

In our example, we have analyze two images and we asked for the text result file. So we have three files by image:



III.7.c.Save/Open a batch

You can save a batch configuration in order to launch it several times or to re-open and modify it. The extension of a batch file is a ".bzb". This file can only be used with AGScan. The format is kind of XML compliant.

Note:In the current version, it is not possible to record a batch in multichannel mode.

To open a batch, use the "File>Open>Batch" menu:



IV. Plugins

 ${\tt IV.1.}$ Presentation

One of the most interesting feature of AGScan is the plugins system. It allows to easily add "functionalities" to AGScan. For example, to add new file formats, alignment or quantification methods... The only thing to do is to add plugins files in the agscan directory named "plugins".

Some plugins are already available on the mulcyber following page:

Plugins 🖾			
SPlotMatrixPlugin	08/	12/2005 14:28	3
SPlotMatrixPlugin.zip	18 KB	14 i386	.zip
SPlotPlugin	08/	12/2005 14:23	1
SPlotPlugin.zip	6 KB	12 i386	,zip
SAlignmentSnapshotPlugin	06/12/2005 15:12		2
SAlignmentSnapshotPlugin.zip	З КВ	19 i386	.zip
SSurfacePlot_3DPlugin	28/11/2005 16:04		
SSurfacePlot_3DPlugin.zip	49 KB	17 i386	.zip
SSurfacePlot_3DPlugin_v2.zip	49 KB	0 i386	.zip
SQuantif1Plugin	28/11/2005 15:13		3
SQuantif1Plugin.zip	27 KB	17 i386	.zip
SFilesMergePlugin	24/11/2005 14:45		
SFilesMergePlugin.zip	4 KB	27 i386	.zip
SMemoryPlugin	24/11/2005 09:15		
SMemoryPlugin.zip	6 KB	22 i386	.zip

Plugins download page

On the website, you can click on the name of the plugin to watch its presentation and on the ."zip" file to download it.

IV.2. Example: the "Memory Plugin"

For example, we will present you how to download/install and run the SmemoryPlugin. This plugin allows to add to AGScan a little window that give you the current memory consumption.

At first, click on the SMemoryPlugin.zip file to download it.

The .zip file contains sources, classes and compile notes (version of javac and AGScan used) of the Memory Controler plugin. Note that there are one source file SMemoryPlugin.java and 3 .class files (caused by embedded classes).

To use this plugin, just put .class files into the plugins directory of AGScan.

A "Memory Controler" submenu is added into the tool menu.





"Memory plugin" window

Notice that different categories of plugins exist. So, each plugin menu appears into its menu (alignment, tool, quantification...). There are also other plugins like "format plugins" that allow to read/save other images formats and that doesn't create menu into the application. You know that they are installed by showing the different type of files you can open. All these inforamtion are available int o the note provided with the plugin.





V. Frequently asked questions

(to be completed)



VI. Index

VI.1. Word Index

-Xmx	10
.bzb	7, 36
.grd	7, 18
.jar	10
.txt	7
.zaf	7, 31, 35
A	
AGScan	1, 2, 4-13, 17, 24, 27, 29, 35-37
Alignment	4-7, 11, 15, 19, 21-24, 26, 27, 29, 31, 34, 35, 37, 38
В	
Batch	4, 5, 7, 11, 31-36
Block	17, 19-24, 27
BZScan	4, 21, 37
С	
ChangeLog.txt	10
D	
Display window	11, 14
Download	4, 8, 9, 37
E	
Export	4, 7, 31, 35
F	
File formats	4, 7, 37
Fit correction	29
Fuji BAS5000 scanner	7
Full	9, 10, 45
G	- , - , - , -
Global alignment	4-6. 21-23
Grid	4, 6, 7, 11, 15-27, 29, 31, 33
Grid structure	4
Image	4, 6, 7, 11-16, 18, 19, 21, 22, 27, 29-31, 33, 35, 42
ImageJ	4, 10
INF/IMG	6
Intensity	5, 29
Interface	11, 45
ISAC	6,7
J	
Java Runtime Environment	9, 10
L	-, -
Languages	10, 12
Launcher configuration file	10
Level1	17
Level2	17
Light	9
Linux OS	10
Local alignment	4-6. 21. 23. 24
М	
Menu bar	11
Mulcyber	4.8.37
0	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,
Overshining correction	2.9
P	

Packages contents	10
Plugin	4, 6, 37, 38
Plugins	4
Progress bar	11, 14, 24, 28
Properties	10, 12, 13, 15
Q	
QM	29
Quantification	4-7, 11, 27-29, 31, 37, 38
S	
Shortcuts	11, 12
SIGENAE	1, 4
Sniffer	21, 22
Spot	5, 6, 17, 20-24, 26, 27, 29, 30
Status bar	11, 15, 16, 26
Structure editor	16, 18
Т	
Tab	15, 16, 18, 32, 33, 35
Table	2, 7, 11, 16, 18, 20, 24, 26-29, 31
TIFF	6
Toolbar	11, 12, 19, 21, 23, 26, 28, 32
V	
Versions	6, 9, 46
W	
Windows OS	10
Z	
Zoom	14-16

VI.2. Illustration index

AGScan logo	1
Mulcyber AGScan page	8
Download AGScan	9
"Open" menu	.12
"Open" shortcut	.12
"Languages" window	.12
"Open image" shortcut	.12
"Open images" window	.13
"image parameters" window	.13
Radioactive experiment open	.13
Cy3/Cy5 experiment open	.14
"Image" menu	.14
"Zoom" shortcut icon	.14
A zoomed image	.15
"Display" menu	.15
"Windows" display mode & "Zoom window"	.16
"Grid Structure Editor" window	.17
Grid structures	.18
"Load a grid" shortcut icon	.18
"Alignment" menu	.19
"Grid elements selection" shortcut icons	.19
Sub-arid selection	.20
"Move/Rotate/Resize" arid	.20
"Global alignment" menu	.21
"Global alignment" shortcut icon	.21
"Global alignment parameters" window	.22
"Sniffer" local alignment algorithm	.22
"Local alignment" menu	.23
"Lobal alignment" shortcut icon	.23
Progression of the local algorithm	.24
An aligned grid	.25
"Spot detectiont" menu	.26
"Spot detection" shortcut icon	.26
Percentage of detected spots	.26
A selected spot and its table description	.27
Columns options	.27
"Quantification>Compute>Selection" shortcut icon	.28
"Quantification>Compute>Selection" menu	.28
"Quantifications selection" window	.28
Example of columns added, results of quantifications	.29
A sorted column	.29
"Modify a column" window	.30
Spot coloration example	.30
"Quantification>Export>Text File" menu	.31
Example of output text file contents	.31
"Batch" shortcut icon	.32
"Batch" tab	.32
"Batch" files configuration	.33
Batch experiment information	.33
"Batch" parameters configuration	.34
Following a "Batch" in progress	.35
Batch progress finished	.35

Elements of the output directory	36
"File>Open>Batch" menu	36
Plugins download page	37
"Memory plugin" window	
······································	

VII. GNU GENERAL PUBLIC LICENSE

Version 2, June 1991

Copyright (C) 1989, 1991 Free Software Foundation, Inc. 51 Franklin St, Fifth Floor, Boston, MA 02110-1301, USA

Everyone is permitted to copy and distribute verbatim copies of this license document, but changing it is not allowed.

Preamble

The licenses for most software are designed to take away your freedom to share and change it. By contrast, the GNU General Public License is intended to guarantee your freedom to share and change free software--to make sure the software is free for all its users. This General Public License applies to most of the Free Software Foundation's software and to any other program whose authors commit to using it. (Some other Free Software Foundation software is covered by the GNU Library General Public License instead.) You can apply it to your programs, too.

When we speak of free software, we are referring to freedom, not price. Our General Public Licenses are designed to make sure that you have the freedom to distribute copies of free software (and charge for this service if you wish), that you receive source code or can get it if you want it, that you can change the software or use pieces of it in new free programs; and that you know you can do these things.

To protect your rights, we need to make restrictions that forbid anyone to deny you these rights or to ask you to surrender the rights. These restrictions translate to certain responsibilities for you if you distribute copies of the software, or if you modify it.

For example, if you distribute copies of such a program, whether gratis or for a fee, you must give the recipients all the rights that you have. You must make sure that they, too, receive or can get the source code. And you must show them these terms so they know their rights.

We protect your rights with two steps: (1) copyright the software, and (2) offer you this license which gives you legal permission to copy, distribute and/or modify the software.

Also, for each author's protection and ours, we want to make certain that everyone understands that there is no warranty for this free software. If the software is modified by someone else and passed on, we want its recipients to know that what they have is not the original, so that any problems introduced by others will not reflect on the original authors' reputations.

Finally, any free program is threatened constantly by software patents. We wish to avoid the danger that redistributors of a free program will individually obtain patent licenses, in effect making the program proprietary. To prevent this, we have made it clear that any patent must be licensed for everyone's free use or not licensed at all.

The precise terms and conditions for copying, distribution and modification follow.

TERMS AND CONDITIONS FOR COPYING, DISTRIBUTION AND MODIFICATION

0. This License applies to any program or other work which contains a notice placed by the copyright holder saying it may be distributed under the terms of this General Public License. The "Program", below, refers to any such program or work, and a "work based on the Program" means either the Program or any derivative work under copyright law: that is to say, a work containing the Program or a portion of it, either verbatim or with modifications and/or translated into another language. (Hereinafter, translation is included without limitation in the term "modification".) Each licensee is addressed as "you".

Activities other than copying, distribution and modification are not covered by this License; they are outside its scope. The act of running the Program is not restricted, and the output from the Program is covered only if its contents constitute a work based on the Program (independent of having been made by running the Program). Whether that is true depends on what the Program does.

1. You may copy and distribute verbatim copies of the Program's source code as you receive it, in any medium, provided that you conspicuously and appropriately publish on each copy an appropriate copyright notice and disclaimer of warranty; keep intact all the notices that refer to this License and to the absence of any warranty; and give any other recipients of the Program a copy of this License along with the Program.

You may charge a fee for the physical act of transferring a copy, and you may at your option offer warranty protection in exchange for a fee.

2. You may modify your copy or copies of the Program or any portion of it, thus forming a work based on the Program, and copy and distribute such modifications or work under the terms of Section 1 above, provided that you also meet all of these conditions:

a) You must cause the modified files to carry prominent notices stating that you changed the files and the date of any change.

b) You must cause any work that you distribute or publish, that in whole or in part contains or is derived from the Program or any part thereof, to be licensed as a whole at no charge to all third parties under the terms of this License.

c) If the modified program normally reads commands interactively when run, you must cause it, when started running for such interactive use in the most ordinary way, to print or display an announcement including an appropriate copyright notice and a notice that there is no warranty (or else, saying that you provide a warranty) and that users may redistribute the program under these conditions, and telling the user how to view a copy of this License. (Exception: if the Program itself is interactive but does not normally print such an announcement, your work based on the Program is not required to print an announcement.)

These requirements apply to the modified work as a whole. If identifiable sections of that work are not derived from the Program, and can be reasonably considered independent and separate works in themselves, then this License, and its terms, do not apply to those sections when you distribute them as separate works. But when you distribute the same sections as part of a whole which is a work based on the Program, the distribution of the whole must be on the terms of this License, whose permissions for other licensees extend to the entire whole, and thus to each and every part regardless of who wrote it.

Thus, it is not the intent of this section to claim rights or contest your rights to work written entirely by you; rather, the intent is to exercise the right to control the distribution of derivative or collective works based on the Program. In addition, mere aggregation of another work not based on the Program with the Program (or with a work based on the Program) on a volume of a storage or distribution medium does not bring the other work under the scope of this License.

3. You may copy and distribute the Program (or a work based on it, under Section 2) in object code or executable form under the terms of Sections 1 and 2 above provided that you also do one of the following:

a) Accompany it with the complete corresponding machine-readable source code, which must be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,

b) Accompany it with a written offer, valid for at least three years, to give any third party, for a charge no more than your cost of physically performing source distribution, a complete machine-readable copy of the corresponding source code, to be distributed under the terms of Sections 1 and 2 above on a medium customarily used for software interchange; or,

c) Accompany it with the information you received as to the offer to distribute corresponding source code. (This alternative is allowed only for noncommercial distribution and only if you received the program in object code or executable form with such an offer, in accord with Subsection b above.)

The source code for a work means the preferred form of the work for making modifications to it. For an executable work, complete source code means all the source code for all modules it contains, plus any associated interface definition files, plus the scripts used to control compilation and installation of the executable. However, as a special exception, the source code distributed need not include anything that is normally distributed (in either source or binary form) with the major components (compiler, kernel, and so on) of the operating system on which the executable runs, unless that component itself accompanies the executable. If distribution of executable or object code is made by offering access to copy from a designated place, then offering equivalent access to copy the source code from the same place counts as distribution of the source code, even though third parties are not compelled to copy the source along with the object code. 4. You may not copy, modify, sublicense, or distribute the Program except as expressly provided under this License. Any attempt otherwise to copy, modify, sublicense or distribute the Program is void, and will automatically terminate your rights under this License. However, parties who have received copies, or rights,

from you under this License will not have their licenses terminated so long as such parties remain in full compliance.

5. You are not required to accept this License, since you have not signed it. However, nothing else grants you permission to modify or distribute the Program or its derivative works. These actions are prohibited by law if you do not accept this License. Therefore, by modifying or distributing the Program (or any work based on the Program), you indicate your acceptance of this License to do so, and all its terms and conditions for copying, distributing or modifying the Program or works based on it.

6. Each time you redistribute the Program (or any work based on the Program), the recipient automatically receives a license from the original licensor to copy, distribute or modify the Program subject to these terms and



conditions. You may not impose any further restrictions on the recipients' exercise of the rights granted herein. You are not responsible for enforcing compliance by third parties to this License.

7. If, as a consequence of a court judgment or allegation of patent infringement or for any other reason (not limited to patent issues), conditions are imposed on you (whether by court order, agreement or otherwise) that contradict the conditions of this License, they do not excuse you from the conditions of this License. If you cannot distribute so as to satisfy simultaneously your obligations under this License and any other pertinent obligations, then as a consequence you may not distribute the Program at all. For example, if a patent license would not permit royalty-free redistribution of the Program by all those who receive copies directly or indirectly through you, then the only way you could satisfy both it and this License would be to refrain entirely from distribution of the Program.

If any portion of this section is held invalid or unenforceable under any particular circumstance, the balance of the section is intended to apply and the section as a whole is intended to apply in other circumstances. It is not the purpose of this section to induce you to infringe any patents or other property right claims or to contest validity of any such claims; this section has the sole purpose of protecting the integrity of the free software distribution system, which is implemented by public license practices. Many people have made generous contributions to the wide range of software distributed through that system in reliance on consistent application of that system; it is up to the author/donor to decide if he or she is willing to distribute software through any other system and a licensee cannot impose that choice.

This section is intended to make thoroughly clear what is believed to be a consequence of the rest of this License.

8. If the distribution and/or use of the Program is restricted in certain countries either by patents or by copyrighted interfaces, the original copyright holder who places the Program under this License may add an explicit geographical distribution limitation excluding those countries, so that distribution is permitted only in or among countries not thus excluded. In such case, this License incorporates the limitation as if written in the body of this License.

9. The Free Software Foundation may publish revised and/or new versions of the General Public License from time to time. Such new versions will be similar in spirit to the present version, but may differ in detail to address new problems or concerns.

Each version is given a distinguishing version number. If the Program specifies a version number of this License which applies to it and "any later version", you have the option of following the terms and conditions either of that version or of any later version published by the Free Software Foundation. If the Program does not specify a version number of this License, you may choose any version ever published by the Free Software Foundation. If you wish to incorporate parts of the Program into other free programs whose distribution conditions are different, write to the author to ask for permission. For software which is copyrighted by the Free Software Foundation, write to the Free Software Foundation; we sometimes make exceptions for this. Our decision will be guided by the two goals of preserving the free status of all derivatives of our free software and of promoting the sharing and reuse of software generally.

NO WĂRRANTY

11. BECAUSE THE PROGRAM IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE PROGRAM, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE PROGRAM "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE PROGRAM IS WITH YOU. SHOULD THE PROGRAM PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

12. IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE PROGRAM AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE PROGRAM (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE PROGRAM TO OPERATE WITH ANY OTHER PROGRAMS), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.