OCS Small Animal Imaging Core Facility

Quick Users Manual

This is a self service facility
It only works if everyone follows the
procedures

How to Start

1. Signup user logbook

Login: IVISuser

Pwd: user


2. Double click on Living Image Software icon



3. Select your 2-letters lab user ID



4. Initialize the unit and WAIT for the GREEN light



This takes about 10 min

You can prepare your animals during this time

URL link to online manual: http://www.caliperls.com/assets/022/8234.pdf


Prepare your animals

Only 3 mice at a time can fit together in the IVIS chamber

- 1. Inject animals with chromiluminescent substrate
- 2. Place a pad of paper towels inside the anesthesia box before introducing your mice
- 3. Anesthetize your mice




4. Place your animal(s) in the IVIS chamber



- Nose in the plastic nose cone adaptor
- Make sure to plug the unused stations with rubber corks (nose cones and plugs are in the table drawer)
- 5. Take <u>ALL</u> the mice out of the IVIS chamber back into a cage asap

Bioluminescence Images

- 1. Place your <u>anesthetized</u> mice in the chamber
- 2. Put a check mark next to Luminescent



3. Auto exposure – At first try:


Exposure time: 1 min

Binning: Large or Medium

• F/Stop: 1

Adjust time if needed


4. Acquire data




5. Label you data file and save it in the "Data temporary folder" (desktop)

Add an X-Ray

- 1. Place your anesthetized mice in the chamber
- 2. Follow "Bioluminescence Images" procedure
- 3. Add a check mark next to X-Ray




4. Set up X-Ray mode



- Lock the door handle completely
- Emergency switch in the ON (out) position
- Key is turned ON.
- Push orange switch (light must be ON)

The X-ray source cannot be energized until these conditions have been fulfilled.

5. Acquire data




- 6. Label you data file and save it in the "Data temporary folder" (desktop)
- 7. Return the key to the OFF position

Fluorescence Images

- 1. Place your <u>anesthetized</u> mice in the chamber
- 2. Turn on the Fluorescence light source (lower shelf)



3. Put a check mark next to Fluorescent



- 3. Select excitation and emission filters
- 4. Acquire data
- 5. Label you data file and save it in the "Data temporary folder" (desktop)

When you are done

- 1. Take ALL the mice out of the IVIS chamber
- 2. Gently wipe clean all surfaces
 Use cleaning solution provided by DLAR
 Do not spray or pour cleaning solution into the chamber
- 3. Clean the anesthesia box
 Remove paper towels
 Wipe with cleaning solution

What to do with the images?

1. Close Living Image Software:

x in the upper right corner of the screen

2. DO NOT LEAVE YOUR IMAGES ON THIS COMPUTER

Option 1: Export in your lab SAN space

• Click on the "Connect Research SAN" icon (desktop)



- Follow the prompts
 - enter your SAN lab name (XXXIab)
 - login with Kerberos ID and password
 - drag your data onto your SAN space
 - **❖ LOGOUT: Click the "Disconnect Research SAN" icon**



Option 2: Copy your data on a USB flash drive

• USB port on the left side of the monitor

Option 3: Burn a CD

- 1. Delete your files from the computer
- 3. Logoff IVISuser





- 4. DO NOT Shutdown the computer
- 5. DO NOT turn off the main ON/OFF switch in the back of the instrument

Data Analysis Workstations

Analysis software is available at 4 workstations:

- Ehrman Medical Library (Graphic area)
- Smilow 308
- MSB 530
- HJD (location TBD)

2-letters lab user IDs

IA Aifantis lab

AE Erlebacher lab

• JH Hay lab

EH Hernando Lab

KS Kirsch lab

DL Levy lab

RN Novick lab

OO O'Connor lab

HB Samuels lab

JT Torres lab

HB Samuels lab

YW Radiology

If your lab does not have a user ID contact Sheenah Mische (Skirball 3). Please do not create your own lab ID.

PLEASE READ

- If you don't know ask!
 - ► <u>Luigi Scotto</u> (O'Connor's lab) ext: 3-9123 Smilow 1306
 - Maria Guijarro (Hernando's lab) ext: 3-9054 Smilow 3
- DO NOT TURN OFF the main ON/OFF switch in the back of the instrument
- DO NOT shutdown the computer
- All data files left on the computer will quickly be erased
- Leave the equipment clean
- Leave the room in order