PhIDO Manual

Phase identification & indexing from ED patterns





©Calidris, Manhemsvägen 4, S-191 45 Sollentuna, Sweden Tel & Fax: +46 8 6250041, http://www.calidris-em.com

1. Introduction

What can PhIDO be used for?

- *Phase identification:* PhIDO tells you which compound your sample corresponds to, by comparison of the observed d-spacings and angles of an electron diffraction pattern with a database of known substances.
- *Indexing*: PhIDO indexes automatically the electron diffraction pattern when the sample is identified and gives the zone axis and the indices of two reflections with the shortest reciprocal lattice vectors.

2. Installation

The PhIDO program is incorporated, together with ELD in the software package CRISP supplied on a CD. PhIDO will be installed automatically when CRISP/ELD is installed (see more details in the CRISP installation).

PhIDO is active when an ED pattern is loaded and a lattice refinement has been performed by ELD.

3. Preparation of database

1. Database format

- The information needed for PhIDO is: lattice type, crystal system and the cell parameters. This information is stored in a database file called PHIDO.TBL under the same directory as CRISP.
- The database of known substances is stored in ASCII form, allowing the user to modify and extend it. It can be edited with, for example, Notepad, or any other text editor. If you use a word processor such as Microsoft Word, you must save the file in plain text mode, not in the format which is used by the word processor.

An initial database containing 77 common phases found in steels is provided in PHIDO.TBL. The format of the database is as follows:

📋 Phido.tbl -	No	otep	ad					_ 🗆	×
<u>F</u> ile <u>E</u> dit <u>S</u> e	arch	ו <u>ד</u>	<u>H</u> elp						
L-NiM	1	2	3.520	3.520	3.630	90.000	90.000	90.000	
CR2N P31m	1	3	4.811	4.811	4.484	90.000	90.000	120.000	
Cr2N P31m	1	3	2.750	2.750	4.439	90.000	90.000	120.000	
FCC-Fe Fm3m	3	1	3.585	3.585	3.585	90.000	90.000	90.000	
BCC-Fe Im3m	2	1	2.870	2.870	2.870	90.000	90.000	90.000	
M23C6Fm3m	3	1	10.621	10.621	10.621	90.000	90.000	90.000	
M6C Fd3m	3	1	11.000	11.000	11.000	90.000	90.000	90.000	
FeCr P42/mnm	1	2	8.790	8.790	4.544	90.000	90.000	90.000	
FeMo P42/mnm	1	2	9.188	9.188	4.812	90.000	90.000	90.000	
LavesP63/mmc	1	3	4.740	4.740	7.720	90.000	90.000	120.000	
Chi I43m	2	1	8.900	8.900	8.900	90.000	90.000	90.000	
WC P6/mmn	1	3	2.904	2.904	2.835	90.000	90.000	120.000	
NbC Fm3m	3	1	4.431	4.431	4.431	90.000	90.000	90.000	
TiC Fm3m	3	1	4.328	4.328	4.328	90.000	90.000	90.000	
MoC Fm3m	3	1	4.274	4.274	4.274	90.000	90.000	90.000	-
4		-							
	Γ	Ţ							
Substance Latti	ce	Crys	tal a	b	С	ά	β	γ	
information type	e	syste	em						

• Substance Information

The first 12 characters are names which identify the substances.

• Lattice type

The numbers in this column correspond to the Bravais lattice types. They are defined on the right.

If the lattice is A or B centred, the lattice vectors have to be redefined (exchange the vectors a, b and c) to be converted into a C-centred one.

Number	Lattice Type
1	P, primitive
2	I, body centred
3	F, face centred
4	C, C-base centred
5	R, rhombohedral

• Crystal system

The numbers in this column correspond to the crystal systems. They are defined on the right.

• Lattice parameters

The following 6 columns (Columns 4 to 9) then give the lattice parameters, *a*, *b*, *c*, in Ångströms and α , β , γ in degrees.

Number	System
1	cubic
2	tetragonal
3	hexagonal
4	orthorhombic
5	monoclinic
6	triclinic

2. Modify the database

You should add your candidate substances in the database PHIDO.TBL. For example, here two phases (marked in the table below), Li2NaTa7O19 with a = 15.23, b = 23.57 and c = 3.84 Å and K-Nb-O with a = 23.73, b = 11.67 and c = 7.84 Å are added in the database PHIDO.TBL. Both have primitive lattices and are orthorhombic, so α , β and γ are all 90 degrees.

📋 Phido.tbl -	No	ter	bad					_ 🗆	×
<u>F</u> ile <u>E</u> dit <u>S</u> e	arch		<u>H</u> elp						
Li2NaTa7019	1	4	15.23	23.57	3.84	90.000	90.000	90.000	
K-ND-0	1	4	23.73	11.67	7.84	90.000	90.000	90.000	
L-NiM	1	2	3.520	3.520	3.630	90.000	90.000	90.000	
CR2N P31m	1	3	4.811	4.811	4.484	90.000	90.000	120.000	
Cr2N P31m	1	3	2.750	2.750	4.439	90.000	90.000	120.000	
FCC-Fe Fm3m	3	1	3.585	3.585	3.585	90.000	90.000	90.000	
BCC-Fe Im3m	2	1	2.870	2.870	2.870	90.000	90.000	90.000	
M23C6Fm3m	3	1	10.621	10.621	10.621	90.000	90.000	90.000	
M6C Fd3m	3	1	11.000	11.000	11.000	90.000	90.000	90.000	
FeCr P42/mnm	1	2	8.790	8.790	4.544	90.000	90.000	90.000	
FeMo P42/mnm	1	2	9.188	9.188	4.812	90.000	90.000	90.000	
LavesP63/mmc	1	3	4.740	4.740	7.720	90.000	90.000	120.000	
Chi I43m	2	1	8.900	8.900	8.900	90.000	90.000	90.000	
WC P6/mmn	1	3	2.904	2.904	2.835	90.000	90.000	120.000	
NbC Fm3m	3	1	4.431	4.431	4.431	90.000	90.000	90.000	
TiC Fm3m	3	1	4.328	4.328	4.328	90.000	90.000	90.000	-
4								<u> </u>	· //

4. Preparing input for PhIDO

1. How to digitize the ED patterns?

The ED patterns can be digitized for example with CRISP or other equipment. In order to know the cell parameter, one has to know the digitization constant, i.e. the number of pixels one millimiter in the ED negative corresponds to. This can be calculated from for example a millimeter paper or an ED pattern with known d-values digitized under the same magnification as the ED negative films (see ELD user manual 5-3 to 5-5).



ED patterns (ta-388.pcx and ta-453.pcx) from two different zones axes taken from the same crystal, with 6.5° angle between them (included in the supplied CD).

2. Prepare input data by ELD

• Information: before processing an ED pattern, information about the conditions under which the ED pattern was taken and digitized should be given in the Information dialogue box. Activate the box next to Diffraction pattern. Then edit the accelerating voltage, camera length, digitization and if necessary the aspect ratio X/Y. The ED constant will be calculated from these data. Save the now calibrated diffraction pattern back to the disk. (If you don't save, the calibration information will not be preserved.)

Information on <ta-388.pcx></ta-388.pcx>	×
<u>Microscope parameters</u>	Camera length
JEOL 2010HT	1000 mm
Acc Voltage 200 (kV)	-Digitization
Sph. Aberr. 1 (mm)	Along X 121613 (pix/mm)
Defocus (u) 578 (Â)	
Defocus (⊻) -578 (Â)	ED constant (d*R) C
Azimuth (x.u) 0 (*)	Along⊻ 305 (Â*pixel)
D <u>e</u> f. Spread 0 (Â)	Aspect Ratio
Convergence 0 (*)	Aspect X/Y 1 (x)
Specimen parameters	Correct
Beam tilt 0 (*)	
Specimen <u>t</u> ilt 0 (*)	I Diffraction pattern
T <u>h</u> ickness 1 (Â)	Update Close

- The parameter ASPECT X/Y present in the Information Panel gives the difference in calibration between the horizontal (x) and vertical (y) directions in the electron diffraction pattern. Some TEMs, scanners and CCD cameras may not have the same magnification along the x and y. If there is a significant difference in magnification along the x and y directions, this should be corrected by editing the Aspect Ratio and then pressing the Correct button.
- ED processing:
 - a. Indexing the ED pattern manually, using a primitive cell. You should select any three reflections which do not lie on the same line (will be marked by red (1), blue (2) and green (3) crosses, respectively) and give their indices in 2D under h and k.

IMPORTANT! The shortest reciprocal vector in the ED pattern should be indexed as 1 0 (U) and the second shortest as 0 1 (V).



b. Activate "**Unbend**", then click on <u>Refine</u>. ELD will find the position of each diffraction spot from the three indexed reflections and then refine the lattice parameters using all the diffraction spots.

5. Phase identification

1. Phase identification from one ED pattern

• When the lattice refinement has been performed by ELD with the ELD window still active (blue on the top of the ELD dialogue box), the PhIDO program can be started by clicking on ELD in the main menu and selecting Phase Identification:

ELD	O <u>p</u> tions	<u>W</u> indov	w <u>H</u> elp
ED	processin	g	Ctrl+E
ED	calibration	ì	
<u>P</u> ha	ase Identifi	action	
ED	I <u>R</u> ad. Distr	ibution	Ctrl+R

A Dialogue box **Phase Identification** will be open.

The PhIDO Dialogue contains two halves, Lattice 1 and Lattice 2, one for each ED pattern. PhIDO will search the database for substances which contain zones which would give the lattice or lattices observed. For each lattice, the lengths of the lattice vectors U and V (in $Å^{-1}$) and the angle between them (U^V in degrees) are specified, together with a tolerance (in percent and degrees, respectively) for defining a "hit" with a possible lattice in the database. The lattice vectors U and V and the angle between them (U^V) determined by ELD will be given automatically in the dialogue box in Lattice 1. The default tolerances are 2%, 2% and 2° for U, V and the U^V, respectively.

🖁 Phase Identification from d-spacing (ta 🗙						
Lattice 1 Lattice 2 /						
U 0.0424 1/Â ± 2.0 % U 0.0424 1/Â ± 2.0 %						
V 0.4199 1/Â ± 2.0 % V 0.4199 1/Â ± 2.0 %						
U^V 90.469 * ± 2.0 * U^V 90.469 * ± 2.0 *						
8 matches (2 phases) 1^2 90.000 * ± 2.0 *						
Find Save Clear list Close						
1 0 0 0 2 3 [0 3 -2] 0.469 0.000 0.001 💌						
; . Eile : CryCyleouie IMOCENTS-200 cou. Phase Identification						
; ; U1 V1 [Zone axes] U^Verr Verr Verr						
Li2NaTa7019 : Lattice=Orthorhombic , Simple (P) 0 -1 0 5 0 -1 [-1 0 -5] 0.469 0.000 0.001 0 -1 0 5 0 1 [1 0 -5] 0.469 0.000 0.001 0 1 0 5 0 1 [1 0 5] 0.469 0.000 0.001 0 1 0 5 0 1 [-1 0 5] 0.469 0.000 0.001 						
K-Nb-O : Lattice=Orthorhombic , Simple (P)						
1 0 0 0 -2 -3 [0 -3 2] 0.469 0.000 0.001 1 0 0 0 -2 3 [0 3 2] 0.469 0.000 0.001						
1 0 0 0 2 -3 [0 -3 -2] 0.469 0.000 0.001						
1 0 0 0 2 3 0 3 -21 0.469 0.000 0.001						

• Click on Find to search the database for hits. When the search is finished, the number of phases and the number of matches (zone axes) found will be given (here two phases Li2NaTa7O19 and K-Nb-O are found). The corresponding indices for the lattice vectors U and V and the zone axis for each match are given under the column U1, V1 and [Zone axes], respectively. The respective errors for the angle between the two vectors U and V are given under U^Verr (in degrees) and those for the two vectors U and V are given under Uerr and Verr (in Å⁻¹), respectively.

This list can be saved in a file with <u>Save...</u>. The list can be cleared with <u>Clear list</u>. Close the PhIDO dialogue box by clicking on <u>Close</u>.

PhIdo list				? ×
🔁 Sample Images	•	£	d	
Li-Ta-O		_		<u>S</u> ave
List files		•		Cancel
	Phido list Sample Images Li-Ta-Ol List files	Phildo list Childo list Li-Ta-Ol List files	Phido list Sample Images E	Phido list Sample Images Li-Ta-Ol List files

2. Phase identification from two ED patterns

• If you have two ED patterns taken from the same crystal and the angle between these two ED patterns are known, you can make the phase identification more accurately. Process the second ED pattern in the same way as the first one and run PhIDO:

Ta-453.pcx - 512x512x8 (1:1)	
	🛔 ELD on Ta-453.pcx 🔀
	Lattice Display
	$\begin{array}{c c c c c c c c c c c c c c c c c c c $
0.0	C 1 6 0 γ=89.602
	C 20 1 0K= 15
and the second	C 3 -6 0 V Unbend ZoomY
# ^{16,00}	Refine Auto Cent. Marks
	Estimation
	C Integration C Amp. Do all
	⊙ Shape fitting ⊙ Int.
	Estimate Protocol Synth
and the second	Remove 'red' reflections Save Clear list Close 3D>>
· · · · · · · · · · · · · · · · · · ·	

ED processing performed on ta-453.pcx. Note that the diffraction rows with weak reflections are from high order Laue zones and should be ignored when indexing, so the Reflection 2 (blue) is indexed as 0 1.

• Similar to the ED pattern ta-388.pcx, PhiDO found two phases from the database, Li2NaTa7O19 and K-Nb-O, which may give similar ED patterns as ta-453.pcx.

Cross the box next to Lattice 2. Fill in the parameters U, V and U^V for the second lattice, obtained from the other ED pattern (here the parameters are from ta-388.pcx). Give the angle between the two ED patterns in the box 1² (here 6.5°). Change the allowed tolerances if needed. Then click Find to search again. on Now phase only one Li2NaTa7O19 is found.

🛔 Phase Identification	n from d-spacing (ta 🗙
Lattice 1	Lattice 2 T
U 0.0423 1/A ± 2.0 %	U 0.0423 1/Å ± 2.0 %
V 0.3704 1/Â ± 2.0 %	V 0.3704 1/Â ± 2.0 %
U^V 90.398 * ± 2.0 *	U~V 90.398 • ± 2.0 •
8 matches (2 phases)	1^2 90.000 * ± 2.0 *
Find Save	Clear list Close
100041[01-	-4] 0.398 0.000 0.005
; File : C:\Crispwin\IMAGE\ta	-453.pcx, Phase Identification
U1 V1 [Zone axe	s] U^Verr Uerr Verr
Li2NaTa7019 : Lattice=Orth 0 -1 0 4 0 -1 [-1 0 -	orhombic , Simple (P) 4] 0.398 0.000 0.001
0 -1 0 4 0 1 [1 0 - 0 1 0 4 0 -1 [1 0	4] 0.398 0.000 0.001 4] 0.398 0.000 0.001
0 1 0 4 0 1 L -1 0 ;	4] 0.398 0.000 0.001
1 0 0 0 -4 -1 [0 -1 1 0 0 0 -4 1 [0 1	4] 0.398 0.000 0.005 4] 0.398 0.000 0.005
1 0 0 0 4 -1 [0 -1 - 1 0 0 0 4 1 [0 1 -	4] 0.398 0.000 0.005 4] 0.398 0.000 0.005
A Phase Identification	from d-spacing (ta ×
Lattice 1	Lattice 2 🔽
U 0.0423 1/Â ± 2.0 %	U 0.0424 1/Â ± 2.0 %
V 0.3704 1/Â ± 2.0 %	V 0.4199 1/Â ± 2.0 %
U^V 90.398 * ± 2.0 *	U^√ 90.469 * ± 2.0 *
4 matches (1 phases)	1^2 6.5000 * ± 2.0 *
Find Save	Clear list Close
; ; ; ; ; ; ; ; ; ; ; ; ; ; ; ; ; ; ; ;	
; U1 V1 [Zone ax.	1] U2 V2 [Zone ax
; Li2NaTa7019 : Orthorhombic	>, Simple (P)
0-1040-1[-10- 0-10401[10-	4] 0-1 0 5 0-1 [-1 0 4] 0-1 0 5 0 1 [1 0
01040-110	4] 0 1 0 5 0 -1 [1 0 4] 0 1 0 5 0 1 [-1 0

- The corresponding 3D indices for the reflections 1 0 (U1) and 0 1 (V1) in the first ED lattice are listed under U1 and V1. The corresponding zone axis is under [Zone ax. 1]. The corresponding 3D indices for the reflections 1 0 (U2) and 0 1 (V2) in the second ED lattice are listed under U2 and V2. The corresponding zone axis is under [Zone ax. 2].
- In general, different phases may give ED patterns with similar lattice vectors U and V, so the number of phases which may match an ED pattern will be more than one, especially if the ED patterns are not from the basic zone axes. Combining two ED patterns normally will give you a unique answer which is the correct.