

NeuroPORTER[®] Transfection Reagent

Instruction Manual

Catalog Number

T400150

T400750



Genlantis

A Division of Gene Therapy Systems, Inc.

10190 Telesis Court

San Diego, CA 92121 USA

Telephone: 858-457-1919

US Toll-free: 888-428-0558

Fax: 858-623-9494

E-mail: tech1@genetherapysystems.com

Web Site: <http://www.genetherapysystems.com>

PAGE INTENTIONALLY LEFT BLANK

Purchaser Notification

Limited License

The purchase price paid for the NeuroPORTER[®] Transfection Reagent kit by end users grants them a non-transferable, non-exclusive license to use the kit and/or its separate and included components (as listed in the Kit Contents section). This kit is intended **for internal research only** by the purchaser. Such use is limited to the transfection of nucleic acid into neuronal cells as described in the product manual. Furthermore, research only use means that this kit and all of its contents are excluded, without limitation, from resale, repackaging, or use for the making or selling of any commercial product or service without the written approval of Genlantis, a division of Gene Therapy Systems, Inc (“Genlantis”).

Separate licenses are available from Genlantis for the express purpose of non-research use or applications of the NeuroPORTER[®] Transfection Reagent. To inquire about such licenses, or to obtain permission to transfer or use the enclosed material, contact the Director of Licensing at Genlantis.

Purchasers may terminate this License at any time by returning all NeuroPORTER[®] Transfection Reagent material and documentation to Genlantis, or by destroying all NeuroPORTER[®] Transfection Reagent kit components. Purchasers are advised to contact Genlantis with the notification that a NeuroPORTER[®] Transfection Reagent kit is being returned in order to obtain a refund and/or to expressly terminate a research only license granted through the purchase of the kit(s).

This document covers in full the terms of the NeuroPORTER[®] Transfection Reagent research only license, and does not grant any other express or implied license. The laws of the State of California shall govern the interpretation and enforcement of the terms of this License.

Product Use Limitations

The NeuroPORTER[®] Transfection Reagent and all of its components are developed, designed, intended, and sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use. All care and attention should be exercised in the handling of the kit components by following appropriate research lab practices.

For more information, or for any comments on the terms and conditions of this License, please contact:

Director of Licensing
Genlantis, a Division of Gene Therapy Systems, Inc.
10190 Telesis Court, San Diego, CA 92121.
Telephone: 858-457-1919
Fax: 858-623-9494.
Email: licensing@genlantis.com

TABLE OF CONTENTS

| | Page |
|--|-------------|
| OVERVIEW | |
| Purchaser Notification | 3 |
| Kit Contents | 5 |
| Shipping and Storage | 5 |
| Related Products | 5 |
| Product Support | 5 |
| Introduction | 6 |
| | |
| METHODS AND PROCEDURES | |
| Transfection of Primary Rat Hippocampal Neurons | 6 |
| Transfection of Other Primary Neurons | 8 |
| Transfection of Neuronal Cell Lines | 10 |
| Transfection of Differentiated Post-Mitotic Neurons and Glial Cell Lines | 12 |
| | |
| APPENDIX | |
| Quality Control | 14 |
| Examples of Optimization of Transfection Conditions | 14 |

OVERVIEW

Kit Contents

Each NeuroPORTER[®] Transfection Reagent kit (Cat. No. T400150) contains sufficient material for 75 to 300 transfection reactions depending on the cell type. Each NeuroPORTER[®] Transfection Reagent kit (Cat. No. T400750) contains sufficient material for 375 to 1500 transfection reactions depending on the cell type. Each reaction is for transfecting 2 µg of DNA.

| Item | Description | Quantity |
|---|---|--|
| NeuroPORTER [®] Transfection Reagent | Dried NeuroPORTER [®] lipid film | 1 vial (Cat # T400150) 5 vials (Cat # T400750) |
| Hydration Buffer | Transfection grade hydration buffer used to hydrate NeuroPORTER [®] dried lipid film before transfection | 1.5 ml (Cat # T400150) 5 x 1.5 ml (Cat # T400750) |
| DNA Diluent | Solution for diluting DNA for optimal transfection efficiency in neuronal cell lines | 7.5 ml (Cat # T400150) 5 x 7.5 ml (Cat # T400750) |

Shipping and Storage

The NeuroPORTER[®] Transfection Reagent kit is shipped at room temperature. For maximum stability store all reagents at 4°C upon receipt. All components are stable for at least one year if stored properly.

Related Products

phCMV Mammalian Expression Vectors

| Product Name | Cat. No. | Quantity |
|-------------------|----------|----------|
| phCMV1 Vector Kit | P003100 | 25 µg |
| phCMV2 Vector Kit | P003200 | 25 µg |
| phCMV3 Vector Kit | P003300 | 25 µg |

For efficient transfection and high-level expression in non-neuronal cells

| Product Name | Cat. No. | Quantity |
|--|----------|----------------------------|
| GenePORTER [®] 2 Transfection Reagent | T202007 | 75 reactions (0.75 ml) |
| GenePORTER [®] 2 Transfection Reagent | T202015 | 150 reactions (1.5 ml) |
| GenePORTER [®] 2 Transfection Reagent | T202075 | 750 reactions (5 x 1.5 ml) |

Product Support

| | |
|--|--|
| Telephone: 858-457-1919 OR 888-428-0558 (US toll free) | Fax: 858-623-9494 |
| E-mail: tech1@genlantis.com | Web: http://www.genlantis.com |

For a complete list of international distributors, visit our web site at www.genlantis.com.

Introduction

NeuroPORTER[®] Transfection Reagent is a novel cationic lipid specially formulated for optimal transfection in neuronal cells, including primary neurons, differentiated post-mitotic neurons, neuronal cell lines, and glial cells. NeuroPORTER[®] Transfection Reagent is much easier to use than the traditional viral delivery method for transfecting DNA into neuronal cells. NeuroPORTER[®] Transfection Reagent is compatible with serum eliminates the need to change media following transfection. An included DNA Diluent is designed to facilitate DNA/lipid complex (lipoplex) formation and enhance the transformation efficiency in certain neuronal cells such as NT2 (not recommended for primary and differentiated neurons). Compared to other commercially-available transfection reagents, NeuroPORTER[®] provides superior transfection efficiency and minimized cytotoxicity. Cell type specific protocols are developed for NeuroPORTER[®] Transfection Reagents to ensure optimal transfection results.

METHODS AND PROCEDURES

1. Transfection of Primary Rat Hippocampal Neurons

- 1.1. Seed primary rat hippocampal cells in poly-D-lysine coated plates (Becton Dickinson Labware) in the numbers listed in Table 1 below using the following *Plating Medium*: *Neurobasal medium* (Invitrogen Cat. No. 21103-049) *supplemented with B27, 0.5 mM L-glutamine and 25 μM glutamate*. Incubate the cells at 37°C in 5% CO₂ for 72 hours.

Table 1: Suggested Cell Plating Numbers

| Tissue Culture Plate | Cell Number | Plating Medium Volume |
|----------------------|--------------------|-----------------------|
| 96-well | 15,000 cells/well | 0.125 ml |
| 24-well | 100,000 cells/well | 0.5 ml |
| 12-well | 200,000 cells/well | 1.0 ml |
| 6-well | 500,000 cells/well | 2.0 ml |

- 1.2. After 72 hours of incubation, remove ½ volume of the Plating Medium and replace with the following *Culture Medium*: *Neurobasal medium supplemented B27 and 0.5 mM L-glutamine (no 25 μM glutamate)*. Continue incubation for an additional 24 hours.
- 1.3. Hydrate the NeuroPORTER lipid vial at room temperature with 1.5 ml of the hydration buffer. Vortex for 30-60 seconds at top speed. Store the hydrated reagent at 4°C and vortex briefly before use.
- 1.4. Dilute the DNA and hydrated NeuroPORTER reagent with serum-free medium. (**do not use the DNA Diluent for primary neurons**) Refer to Tables 2 and 3 for recommended DNA, NeuroPORTER, and serum-free medium volumes for different tissue culture plates.

Table 2: Volumes of Transfection Reagents

| DNA (µg) | Serum Free Medium for DNA (µl) | NeuroPORTER (µl) | Serum Free Medium For NeuroPORTER (µl) |
|-----------------|---------------------------------------|-------------------------|---|
| 0.1 – 0.5 | 12.5 | 2.5 | 10 |
| 1 – 3 | 25 | 5 | 20 |
| 2 – 4 | 37.5 | 7.5 | 30 |
| 4 – 6 | 62.5 | 12.5 | 50 |

- 1.5. Add the DNA solution to the diluted NeuroPORTER Transfection Reagent. Mix by pipetting up and down several times. Incubate at room temperature for 10 minutes to allow the NeuroPORTER /DNA complexes to form. *Do not incubate for longer than 30 minutes.*

Table 3: Medium Volumes and DNA Amount for Various Culture Dishes

| Tissue Culture Dish | DNA (µg) | Serum Free Medium Volume (ml) | Total Transfection Volume (ml) |
|----------------------------|-----------------|--------------------------------------|---------------------------------------|
| 96-well | 0.1-0.5 | 0.1 | 0.125 |
| 24-well | 1-3 | 0.45 | 0.5 |
| 12-well | 2-4 | 0.925 | 1.0 |
| 6-well | 4-6 | 1.375 | 1.5 |

- 1.6. Remove the Plating Medium from the cells, and add the volume of serum-free medium indicated in Table 3 to each well.
- 1.7. Apply the DNA/NeuroPORTER complexes from step 1.5 to each well. The total transfection volume at this step is indicated in Table 3.
- 1.8. Gently mix the DNA/NeuroPORTER/serum-free medium by swirling, and place the cells in a 37°C incubator with 5% CO₂.
- 1.9. After two hours of incubation, add one additional volume of fresh Culture Medium containing 2X concentration of B27 onto the cells.
- 1.10. Perform assay for gene expression after 24-48 hours.

2. Transfection of Other Primary Neurons

- 2.1. Hydrate the NeuroPORTER lipid vial at room temperature with 1.5 ml of the hydration buffer. Vortex for 30-60 seconds at top speed. Store the hydrated reagent at 4°C and vortex briefly before use.
- 2.2. Dilute the hydrated NeuroPORTER reagent with serum-free medium. Refer to Table 4 for the appropriate volume of serum-free medium.

Table 4: Volumes of Transfection Reagents

| DNA (µg) | Serum Free Medium for DNA (µl) | NeuroPORTER (µl) | Serum Free Medium for NeuroPORTER (µl) |
|----------|--------------------------------|------------------|--|
| 0.5 | 12.5 | 2.5 | 10 |
| 1 | 20 | 5 | 15 |
| 2 | 40 | 10 | 30 |
| 4 | 55 | 20 | 35 |
| 6 | 70 | 30 | 40 |
| 8 | 110 | 40 | 70 |

Although NeuroPORTER has been optimized for specific cell culture conditions, optimization may be needed to achieve maximum transfection efficiency. The two critical variables are the ratio of NeuroPORTER reagent to DNA and the quantity of DNA used. For optimization of the ratio of NeuroPORTER reagent to DNA start by using 2.5 to 15 µl of reagent for each 1 µg of DNA. Use a fixed amount of DNA or vary the amount as suggested in the Appendix to optimize this ratio.

- 2.3. Dilute the DNA with the serum free medium (**do not use the DNA Diluent for primary neurons**). Refer to Table 4 for the appropriate volume of serum-free medium.

To obtain maximum efficiency in particular cells, some optimization may be needed. The two critical variables are the ratio of NeuroPORTER reagent to DNA and the quantity of DNA used. For optimization of the DNA quantity used, maintain a fixed ratio of NeuroPORTER reagent to DNA, and then vary the DNA quantity over a suggested range (see Table 5). See the Appendix for examples.

- 2.4. Add the DNA solution to the diluted NeuroPORTER Transfection Reagent. Incubate at room temperature for 5 to 10 minutes to allow the NeuroPORTER /DNA complexes to form.

Do not incubate the DNA solution with the NeuroPORTER Transfection Reagent for longer than 30 minutes

- 2.5. Add your complexes directly to the cells growing in serum-containing culture medium. Refer to Table 5 for suggested medium volumes.

Table 5: Medium Volumes and DNA Amount for Various Culture Dishes

| Tissue Culture Dish | DNA (μg) | Medium Volume (ml) |
|----------------------------|--------------------------------|---------------------------|
| 96-well | 0.1-0.5 | 0.2 |
| 24-well | 0.5-3 | 0.5 |
| 12-well | 1-4 | 1 |
| 6-well | 2-6 | 1.5 |
| 60 mm | 6-8 | 2.5 |
| 100 mm | 8-12 | 5 |

- 2.6. Add fresh growth media as needed 24 hours post transfection,. Depending on the cell type and promoter activity, the assay for the reporter gene can be performed 24 to 72 hours following transfection.

For some cell types, the old media can be replaced with fresh media at this step.

3. Transfection of Neuronal Cell Lines

- 3.1. Hydrate NeuroPORTER lipid film at room temperature with 1.5 ml of the hydration buffer. Vortex for 30-60 seconds at top speed. Store the hydrated reagent at 4°C and vortex briefly before use.
- 3.2. Dilute the hydrated NeuroPORTER reagent with serum-free medium. Refer to Table 6 for the appropriate volume of serum-free medium.

Table 6: Volumes of Transfection Reagents

| DNA (μg) | DNA Diluent (μl) | NeuroPORTER (μl) | Serum Free Medium for NeuroPORTER (μl) |
|--------------------------|----------------------------------|----------------------------------|---|
| 0.5 | 6.25 | 1.25 | 5 |
| 1 | 12.5 | 2.5 | 10 |
| 2 | 25 | 5 | 20 |
| 4 | 50 | 10 | 40 |
| 6 | 75 | 15 | 60 |
| 8 | 100 | 20 | 80 |

- 3.3. Dilute the DNA with the DNA Diluent and incubate 1 to 5 minutes at room temperature. Refer to Table 6 for the appropriate volume of DNA Diluent.

Do not incubate DNA with the DNA Diluent for longer than 5 minutes. Avoid vortexing the DNA diluent.

Although NeuroPORTER consistently delivers high transfection efficiencies, in order to obtain maximum efficiency in particular cell types, some optimization may be needed. The two critical variables are the ratio of NeuroPORTER reagent to DNA and the quantity of DNA used. For optimization, first maintain a fixed ratio of NeuroPORTER reagent to DNA, and then vary the DNA quantity over the suggested range. If necessary, optimize the ratio of NeuroPORTER reagent to DNA by using 1.25 to 12.5 μl of reagent for each 1 μg of DNA. Use a low DNA quantity to optimize this ratio. Following this process, cell number can also be optimized. See the Appendix for examples.

- 3.4. Add the DNA solution to the diluted NeuroPORTER Transfection Reagent. Incubate at room temperature for 5 to 10 minutes to allow the NeuroPORTER /DNA complexes to form.

Do not incubate the DNA solution with the NeuroPORTER Transfection Reagent for longer than 30 minutes

- 3.5. Add your complexes directly to the cells growing in serum-containing culture medium. Refer to Table 7 for suggested cell numbers for specific tissue culture dishes. Refer to Table 8 for appropriate medium volumes.

Cells plated the day before transfection should be 50% to 70% confluent on the day of transfection.

Table 7: Suggested Cell Culture Conditions for Transfection of Neuronal Cell Lines

| Tissue Culture Dish | Number of Cells / Well |
|---------------------|---------------------------|
| 96-well | 25-30 x 10 ³ |
| 24-well | 125-150 x 10 ³ |
| 12-well | 250-300 x 10 ³ |
| 6-well | 500-600 x 10 ³ |
| 60 mm | 1-1.5 x 10 ⁶ |
| 100 mm | 2.5-3 x 10 ⁶ |

Table 8: Medium Volumes and DNA Amount for Various Culture Dishes

| Tissue Culture Dish | DNA (µg) | Medium Volume (ml) |
|---------------------|----------|--------------------|
| 96-well | 0.1-0.5 | 0.2 |
| 24-well | 0.5-3 | 0.5 |
| 12-well | 1-4 | 1 |
| 6-well | 2-6 | 1.5 |
| 60 mm | 6-8 | 2.5 |
| 100 mm | 8-12 | 5 |

- 3.6. Add fresh growth media as needed 24 hours post transfection,. Depending on the cell type and promoter activity, the assay for the reporter gene can be performed 24 to 72 hours following transfection.

For some cell types, the old media can be replaced with fresh media at this step.

The same protocol can be used to produce stably transfected cells: 48 to 72 hours post transfection, put the cells in fresh medium containing the appropriate selection antibiotic. It is important to wait at least 48 hours before exposing the transfected cells to the selection media. For some cell types it may be necessary to wait as long as 4 to 5 days before applying the selection condition.

4. Transfection of Differentiated Post-Mitotic Neurons and Glial Cell Lines

- 4.1. Hydrate NeuroPORTER lipid film at room temperature with 1.5 ml of the hydration buffer. Vortex for 30-60 seconds at top speed. Store the hydrated reagent at 4°C and vortex briefly before use.
- 4.2. Dilute the hydrated NeuroPORTER reagent with serum-free medium. Refer to Table 9 for the appropriate volume of serum-free medium.

Table 9: Volumes of Transfection Reagents

| DNA (µg) | Serum Free Medium for DNA (µl) | NeuroPORTER (µl) | Serum Free Medium for NeuroPORTER (µl) |
|----------|--------------------------------|------------------|--|
| 0.5 | 15 | 5 | 10 |
| 1 | 25 | 10 | 15 |
| 2 | 50 | 20 | 30 |
| 4 | 75 | 40 | 35 |
| 6 | 100 | 60 | 40 |
| 8 | 150 | 80 | 70 |

- 4.3. Dilute the DNA with the serum free medium. Refer to Table 9 for the appropriate volume of serum-free medium.

Although NeuroPORTER consistently delivers high transfection efficiencies, in order to obtain maximum efficiency in particular cell types, some optimization may be needed. The two critical variables are the ratio of NeuroPORTER reagent to DNA and the quantity of DNA used. For optimization, first maintain a fixed ratio of NeuroPORTER reagent to DNA, and then vary the DNA quantity over the suggested range. If necessary, optimize the ratio of NeuroPORTER reagent to DNA by using 5 to 20 µl of reagent for each 1 µg of DNA. Use a low DNA quantity to optimize this ratio. Following this process, cell numbers can also be optimized. See the Appendix for examples.

- 4.4. Add the DNA solution to the diluted NeuroPORTER Transfection Reagent. Incubate at room temperature for 5 to 10 minutes to allow the NeuroPORTER /DNA complexes to form.

Do not incubate the DNA solution with the NeuroPORTER Transfection Reagent for longer than 30 minutes

- 4.5. Add your complexes directly to the cells growing in serum-containing culture medium. Refer to Table 10 for suggested cell number according to culture dishes size and cell types. Refer to Table 11 for appropriate medium volumes.

Cells plated the day before transfection should be 50% to 70% confluent on the day of transfection.

Table 10: Suggested Cell Culture Conditions for Transfection of Differentiated Neurons and Glial Cells

| Tissue Culture Dish | Cells / Well Diff. Neurons | Cells / Well Glial Cells |
|---------------------|-------------------------------|-----------------------------|
| 96-well | 35 x 10 ³ | 50 x 10 ³ |
| 24-well | 150 x 10 ³ | 200 x 10 ³ |
| 12-well | 300 x 10 ³ | 400 x 10 ³ |
| 6-well | 600 x 10 ³ | 800 x 10 ³ |
| 60 mm | 1.5 x 10 ⁶ | 2 x 10 ⁶ |
| 100 mm | 3 x 10 ⁶ | 4 x 10 ⁶ |

Table 11: Medium Volumes and DNA Amount for Various Culture Dishes

| Tissue Culture Dish | DNA (µg) | Medium Volume (ml) |
|---------------------|----------|--------------------|
| 96-well | 0.1-0.5 | 0.2 |
| 24-well | 0.5-3 | 0.5 |
| 12-well | 1-4 | 1 |
| 6-well | 2-6 | 1.5 |
| 60 mm | 6-8 | 2.5 |
| 100 mm | 8-12 | 5 |

- 4.6. 24 hours post transfection, add fresh growth media as needed. Depending on the cell type and promoter activity, the assay for the reporter gene can be performed 24 to 72 hours following transfection.

For some cell types, the old media can be replaced with fresh media at this step. Also, the same protocol can be used to produce stably transfected cells: 48 to 72 hours post transfection, put the cells in fresh medium containing the appropriate selection antibiotic. It is important to wait at least 48 hours before exposing the transfected cells to the selection media. For some cell types it may be necessary to wait as long as 4 to 5 days before applying the selection condition.

APPENDIX

Quality Control

To assure the performance of each lot of the NeuroPORTER reagent, we pre-qualify the chemical synthesis of NeuroPORTER lipid by mass spectrometry and thin layer chromatography. The final product is further tested by in vitro β -galactosidase transfection assay in NT2 neuronal precursor cell. Each lot shall have an acceptance specification of >70% of the activity of the Reference lot.


Examples of Optimization of Transfection Conditions

1. Optimization conditions for primary neuron transfection in 24-well plates

Follow the general protocol to prepare the DNA/NeuroPORTER complexes. We **do not recommend** using the DNA Diluent for primary neurons.

| Condition | DNA dilutions in serum free medium | NeuroPORTER dilutions in serum free medium | Total Volume | Final DNA Concentration |
|-----------|---------------------------------------|--|-------------------|----------------------------|
| 1 | 10 μg in 250 μl | 25 μl in 225 μl ($V_t = 250 \mu\text{l}$) | 500 μl | 20 $\mu\text{g}/\text{ml}$ |
| 2 | " | 50 μl in 200 μl ($V_t = 250 \mu\text{l}$) | " | " |
| 3 | " | 75 μl in 175 μl ($V_t = 250 \mu\text{l}$) | " | " |
| 4 | " | 100 μl in 150 μl ($V_t = 250 \mu\text{l}$) | " | " |
| 5 | " | 125 μl in 125 μl ($V_t = 250 \mu\text{l}$) | " | " |
| 6 | " | 150 μl in 100 μl ($V_t = 250 \mu\text{l}$) | " | " |

Add the appropriate volume of complexes solution directly to your cells as illustrated below.




2. Optimization conditions for neuronal cell line transfection in 24-well plates

Follow the general protocol to prepare the DNA/NeuroPORTER complexes. We **recommend** using the DNA Diluent for neuronal cell lines such as NT2.

| Condition | DNA Diluent | NeuroPORTER dilutions in serum free medium | Total Volume | Final DNA Concentration |
|-----------|-----------------|--|--------------|-------------------------|
| 1 | 10 µg in 125 µl | 12.5 µl in 112.5 µl (Vt = 125 µl) | 250 µl | 40 µg/ml |
| 2 | " | 25 µl in 100 µl (Vt = 125 µl) | " | " |
| 3 | " | 50 µl in 75 µl (Vt = 125 µl) | " | " |
| 4 | " | 75 µl in 50 µl (Vt = 125 µl) | " | " |
| 5 | " | 100 µl in 25 µl (Vt = 125 µl) | " | " |
| 6 | " | 125 µl NeuroPORTER | " | " |

Add the appropriate volume of complexes solution directly to your cells as illustrated below.



3. Optimization conditions for differentiated post-mitotic neurons and glial cell line transfection in 24-well plates

Follow the general protocol to prepare the DNA/NeuroPORTER complexes. We **do not recommend** using the DNA Diluent for differentiated post-mitotic neurons and glial cells.

| Condition | DNA dilutions in serum free medium | NeuroPORTER dilutions in serum free medium | Total Volume | Final DNA Concentration |
|-----------|------------------------------------|--|--------------|-------------------------|
| 1 | 10 µg in 250 µl | 50 µl in 200 µl (Vt = 250 µl) | 500 µl | 20 µg/ml |
| 2 | " | 75 µl in 175 µl (Vt = 250 µl) | " | " |
| 3 | " | 100 µl in 150 µl (Vt = 250 µl) | " | " |
| 4 | " | 125 µl in 125 µl (Vt = 250 µl) | " | " |
| 5 | " | 150 µl in 100 µl (Vt = 250 µl) | " | " |
| 6 | " | 200 µl in 50 µl (Vt = 250 µl) | " | " |

Add the appropriate volume of complexes solution directly to your cells as illustrated below.

