

TF Interaction Arrays

User Manual
Cat #: MA5000-MA5063

Panomics, Inc.

p300-TF Interaction Array User Manual

Copyright

© Copyright 2008, Panomics, Inc. All rights reserved.

Trademarks

Hyperfilm and ECL are trademarks of GE Healthcare, Inc., FluorChem is a registered trademark of AlphaInnotech, Inc. All other trademarks are of their respective owners.

Citing Protein/DNA Arrays in Publications

When describing a procedure for publication using this product, we would appreciate it if you would refer to it as the p300-TF Interaction Array Kit from Panomics

If a paper cites our product and is published in a research journal, the lead author(s) may receive a travel stipend for use at a technology conference or tradeshow by sending a copy of the paper to our technical support group at techsupport@panomics.com or via fax at (510) 818-2610.

Disclaimer

Panomics, Inc. reserves the right to change its products and services at any time to incorporate technological developments. This manual is subject to change without notice.

Although this manual has been prepared with every precaution to ensure accuracy, Panomics, Inc. assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.

Introduction

Introduction The human genome encodes approximately 35,000 genes. While some of these, called housekeeping genes, are constitutively expressed, most are not. Non-housekeeping genes are expressed in a particular type of cell, under specific environmental conditions, or at a certain stage of development. The expression of such genes must be tightly regulated.

Gene expression is regulated by a group of proteins known as transcription factors (TFs). An estimated 2,000 TFs are at work in any given human cell. By binding to specific DNA *cis*-elements, each TF contributes a different functional role in gene expression.

Just as TFs regulate gene expression, the TFs themselves are regulated in a number of ways. Known as combinatorial regulation, these different levels of regulation allow gene expression to be tightly controlled. TFs may be regulated by protein modification, such as phosphorylation; protein translocation; and/or protein interaction. In the latter category, TFs may interact with upregulators or downregulators in the basal transcription machinery, or with other transcription factors.

TF-TF **A novel method for detecting TF-TF interactions**

Traditionally, TF-TF interactions have been studied by co-immunoprecipitation and super-gel shift. But because these methods are notoriously time consuming and inefficient, they are not conducive to mapping the network of TF-TF interactions.

TF-TF Interaction Arrays enable you to determine how a particular TF interacts with multiple other TFs—all in one experiment. This technology is based on the non-covalent binding of DNA binding sequences to TF proteins. After precipitation with specific antibodies against a target TF, the co-precipitated TFs are determined by array analysis of attached DNA sequences. The advantage of this technology is its capability of high-throughput analyzing interactions of native TF proteins in a given sample. Figure 1 illustrates how this simple procedure works. First, incubate your nuclear extract with the provided set of biotin-labeled double-stranded oligonucleotide probes, which represent a known library of *cis*-elements. During the incubation step, these probes bind to their corresponding TFs in your nuclear extract. The next step is immunoprecipitation: using a specific antibody against your TF of interest, pull out that TF and associated TFs, along with corresponding *cis*-elements. Next, wash away free *cis*—an extracellular stimulus binds to a cell surface receptor. Binding triggers a signaling cascade, which migrates through the cytoplasm and into the nucleus. Once transmitted to the nucleus, the signal regulates transcription factors, which turn on and off genes that elements and nonspecific binding proteins. Finally, elute the *cis*-elements and hybridize them to the Protein/DNA Array membrane.

We currently offer three versions of the TF-TF Interaction Array. Version I is spotted with 54 different consensus sequences; Version II is spotted with 96 different consensus sequences; and Version III is spotted with 94 different consensus sequences. (Maps of Versions I, II and III are provided in Appendices B, C and D, respectively.) By

combining all of these versions, you can profile the interactions of 244 unique TFs. However, the TF-TF Interaction Array system can be easily expanded to include up to 2,000 consensus sequences.

p300 To initiate transcription, activated TFs must either bind to their *cis*-element sequences or interact with other TFs. To connect the TFs to one another and to the basal transcriptional machinery, secondary proteins called cofactors are usually required. The interactions between TFs, cofactors, and the components of the transcriptional machinery create a complex, multidimensional transcription network.

One group of proteins involved in this network is the p300/CREB-binding protein (CBP) family. Originally discovered through their binding to the adenoviral E1A and CREB proteins, p300 and CBP are now known to be potent transcriptional co-activators. More than 30 transcription factors-p53, ER, HIF1, Stat1, E2F-1, and Jun to name a few-have been identified that interact with p300/CBP (1-3).

To regulate TFs, p300/CBP integrates into the transcriptional apparatus, becoming a physical link between TFs, cofactors, and the basal transcription machinery. This interaction occurs through three distinct mechanisms: p300/CBP acts as a bridge to connect TFs directly to the basal transcription machinery (3-5); p300/CBP forms a scaffold for the assembly of multi-protein cofactor complexes (5-6); and p300/CBP increases access to the DNA by modifying chromatin structure with histone acetyltransferase (HAT) activity (5,7-8).

The distribution of p300/CBP among TFs is believed to be signal-dependent (5, 9). By responding to activated signal pathways, p300/CBP mediates a wide array of cellular processes, including proliferation, differentiation, and apoptosis (1,3-5). Mutations in p300/CBP genes have even been found in some human tumors, suggesting a tumor suppressor functionality (10--10). To fully understand the many physiological roles of p300/CBP and the relationships between signaling pathways and TFs, the intricate network of p300/CBP interactions must be dissected.

We currently offer two versions of the p300-TF Protein Interaction Array. Version I is spotted with 54 different consensus sequences; Version II is spotted with 96 different consensus sequences. (Schematic Diagrams of Versions I and II are provided in Appendices B and C, respectively.) By combining both versions, you can profile the interactions of p300 with 150 unique TFs.

SRC Hormones, vitamins, and other metabolites exert their cellular effects through binding to a specialized class of transcription factors known as nuclear receptors. By interacting directly with DNA response elements as either homo- or heterodimers, nuclear receptors transduce the effects of their ligands into the regulation of gene expression. In this manner, nuclear receptors play a significant role in critical physiological processes, such as development, reproduction, and metabolism.

Like most transcription factors, nuclear receptors depend on cofactors, or coregulator proteins, to regulate gene expression. Cofactors can act as a bridge to bring

nuclear receptors into close contact with chromatin or the basal transcription machinery to facilitate transcription. Alternatively, cofactors can form multi-protein complexes with nuclear receptors to alter their ability to bind DNA. Cofactors with negative effects on transcription are called corepressors, whereas cofactors with positive effects on transcription are called coactivators.

Steroid and nuclear hormone receptors are coactivated by interaction with a cofactor dubbed the steroid receptor coactivator (SRC). There are three splicing variants of this particular coactivator. While the functionality of each isoform differs, the SRC family shares histone acetyltransferase activity and characteristic structural elements, which include a basic helix-loop-helix motif, a PAS domain, and a nuclear localization signal.

The first SRC to be cloned and identified was steroid receptor coactivator 1 (SRC1). This 160 kDa cofactor enhances the transcriptional activity of the nuclear receptors for estrogen, progesterone, glucocorticoid, thyroid hormone, and retinoid X hormone. Research has shown that SRC1 can also bind to transcription factors that are not nuclear receptors, such as TBP and TFIIB. These findings suggest that SRC1 is a bridging coactivator that connects nuclear receptors to other transcription factors and an integral part of the regulation of hormone-induced cellular response. Traditionally, protein interactions have been studied by co-immunoprecipitation and super-gel shift. But because these methods are notoriously time consuming and inefficient, they are not conducive to mapping the network of interactions between the TFs and the SRC1 cofactor. The SRC1-TF Interaction Arrays enable you to determine how SRC1 interacts with multiple TFs—all in one experiment. This array-based system for detecting protein interactions is based on Panomics proprietary TF-TF Interaction Array technology.

We currently offer two versions of the SRC1-TF Protein Interaction Array. Version I is spotted with 54 different consensus sequences; Version II is spotted with 96 different consensus sequences; and Version III is spotted with 94 different consensus sequences. (Schematic Diagrams of Versions I, II, and III are provided in Appendices B, C, and D, respectively.) By using all three versions, you can profile the interactions of SRC1 with 244 unique TFs.

PPARalpha Peroxisome Proliferator-Activated Receptors, or PPARs, are nuclear proteins that compose a subset of the steroid receptor super-family. These proteins were first identified by their ability to bind peroxisome proliferators, fatty acid-like molecules that modulate hepatic lipid metabolism. It is now well established that PPARs are ligand-dependent TFs and are activated by PPs, endogenous fatty acids, and other metabolites.

When bound to their cognate ligands, PPARs form a heterodimer with another steroid receptor, retinoid receptor X (RXR). This protein complex then induces the expression of metabolism-related genes through binding to specific DNA response elements. PPARs are also known to interact with other TFs, including Sp1 and Smad4, to regulate transcription. Through these mechanisms, the PPAR family of nuclear receptors

is responsible for the cellular response to excess or shortage of dietary factors, an essential process for maintaining homeostasis.

Three isoforms of PPAR exist— α , β , and γ —each playing a separate role in systemic metabolism. The localization of each protein is illustrative of this diversity of function: PPAR α is found predominantly in the liver, heart, kidney, as well as other tissues where metabolism is elevated. PPAR γ , on the other hand, is found only in adipose tissue and tissues related to the immune system. Of the three, PPAR β is the most ubiquitously expressed.

Despite its abundance, the physiological role of PPAR β is largely unknown. PPAR α and PPAR γ , on the other hand, have been identified as regulators of critical biological processes. PPAR α is involved in the transcriptional regulation of enzymes that break-down fatty acids; PPAR γ promotes the creation of adipocytes, the fat cells that control lipid homeostasis and glucose metabolism. Disruption of PPAR function has been implicated in the development of several chronic diseases, including diabetes and cancer

NF κ B p50 More than 15 years after its discovery, nuclear factor- κ B (NF κ B) still commands the interest of researchers in diverse fields. This is due to a variety of factors, including NF κ B's unusual and rapid regulation, the wide range of genes that it controls, its central role in immunological processes, and its involvement in various diseases.

Regulation of NF κ B

NF κ B is a family of structurally and functionally related proteins that regulate the transcription of a wide variety of genes. The NF κ B proteins, which include p65 (Rel A), c-Rel, p105/50 (NF κ B1), p100/52 (NF κ B2), and p68 (Rel B), interact to form homodimers or heterodimers that bind DNA with different specificities.

Uninduced NF κ B is sequestered in the cytoplasm and is bound by I κ Bs, proteins that specifically inhibit NF κ B family members. NF κ Bs can be activated by a number of agents, including lipopolysaccharide, phorbol esters, IL-1, TNF α , protein synthesis inhibitors, and reactive oxygen intermediates. Treatment with these agents leads to the phosphorylation and degradation of I κ B, and subsequent translocation of induced NF κ B to the nucleus.

NF κ B target genes

NF κ B regulates the transcription of various target genes. Just a few examples of NF κ B target genes include interferon β , various cytokines, genes encoding cell adhesion molecules and p53.

The expression of NF κ B target genes depends not only on the activity of NF κ B itself, but also on a number of other factors, including NF κ B's interaction with other tran-

scription factors. Thus, NFkB target genes might not be induced even when NFkB is active. Furthermore, the pattern of NFkB target genes will be different under different physiological conditions. By profiling the expression pattern of NFkB target genes, we can begin to elucidate the NFkB-linked functions of those genes.

Both NFkB p50-TF and NFkB p65-TF Interaction Arrays come in three versions. Version I is spotted with 54 different consensus sequences; Version II is spotted with 96 different consensus sequences; and Version III is spotted with 94 different consensus sequences. (Schematic Diagrams of Versions I, II, and III are provided in Appendices B, C, and D, respectively.) By using all three versions, you can profile the interactions of NFkB with 244 unique TFs.

GRalpha Glucocorticoid receptor (GR) is a 94-kDa intracellular protein, which belongs to the family of steroid hormone receptors. GR mediates the effects of glucocorticoids and functions as a ligand-activated transcription factor playing an important role in the development, metabolism, and neurobiology of programmed cell death and anti-inflammation. In the absence of glucocorticoid, GR is bound by HSP 90 and localizes in the cytoplasm. Upon glucocorticoid binding, GR becomes hyperphosphorylated, dissociates from Hsp90, and translocates into the nucleus. In the nucleus, GR binds to a specific response element in the promoter region of a target gene and regulates gene expression. These responsive genes make up around 20 % of the human genome. As corresponding ligands, glucocorticoids are among the most widely prescribed class of drugs in the world.

GR-TF Protein Interaction Arrays come in three versions. Version I is spotted with 54 different consensus sequences; Version II is spotted with 96 different consensus sequences; and Version III is spotted with 94 different consensus sequences. (Schematic Diagrams of Versions I, II, and III are provided in Appendices B, C, and D, respectively.) By using all three versions, you can profile the interactions of GR with 244 unique TFs.




Figure 1: Schematic flow chart of the TF Interaction Array procedure

Materials

Materials Provided

Box 1: Array Membranes & Hybridization Reagents	
Interaction Array (either TF-TF, p300-TF, SRC1-TF, PPARalpha-TF, NFkB p50-TF or GRalpha-TF)	3 each
Hybridization Buffer	15 ml
Substrate Solution I	3 ml
Substrate Solution II	3 ml
Substrate Solution III	5 ml
1000x Streptavidin HRP Conjugate	60 µl
20X SSC	32 ml
20% SDS	15 ml
4X Wash Buffer	45 ml
10X Detection Buffer	25 ml
Distilled H ₂ O	RNase & DNase free; 300 µl

Box 2: Reaction Kit	
TF-TF Probe Mix	60 µl
Poly d(I-C)	50 µl
Control Nuclear Extract	10 µl
2X Blocking Buffer	30 ml
5X Binding Buffer	100 µl
5X IP Wash Buffer	2 x 1.5 ml
1X IP Dilution Buffer	1.5 ml
1X IP Elution Buffer	0.5 ml

* These reagents are provided as stock solutions. See Appendix A for instructions on how to dilute these solutions/buffers.

Depending on which kit you purchased, you will also receive:

p53 Goat Polyclonal, IgG (TF-TF)	10 µl; 200 µg/ml
p300 Antibody, rabbit polyclonal (p300-TF)	45 µl; 200 µg/ml
SRC1 Antibody, rabbit polyclonal (SRC1-TF)	45 µl; 200 µg/ml
PPARα Antibody: either Polyclonal OR Monoclonal (PPARalpha-TF)	45 µl; 200 µg/ml
NFkB p50 Antibody: either Polyclonal OR Monoclonal (NFkB p50-TF)	45 µl; 200 µg/ml
GR Antibody Polyclonal, Rabbit IgG (GRalpha-TF)	45 µl; 200 µg/ml

Additional Materials Required

Reagents & Solutions
Nuclear Extraction kit (e.g., Panomics Nuclear Extraction Kit, Cat. # AY2002)
Hybridization Wash I (2X SSC/0.5% SDS)*
Hybridization Wash II (0.1X SSC/0.5% SDS)*
Deionized H ₂ O
Antibody directed against your TF of interest.
[Optional] Normal IgG Antibody for your negative control during immunoprecipitation. (Choose a normal IgG antibody that is specific to your antibody of choice.)

See Appendix A for recipes.

Materials and Equipment
Dynabeads with either protein A (DynaL Biotech, Cat# 100.01) or protein G (DynaL Biotech, Cat# 100.03) depending on antibody used for immunoprecipitation. <i>Please do not use agarose or sepharose beads as these bind non-specifically to the DNA probes, providing false positive results.</i>
Magnetic separation stand (Promega, Cat. # Z5332)
0.5-ml and 1.5-ml tubes
Pipetman and tips
Microcentrifuge
Hybridization oven and bottles † (Stratagene, Cat.# 401030) † Note: skirted centrifuge tubes with screw caps may be used in place of the hybridization bottles (VWR, Cat. # 21008-480). Hybridization bottle dimensions: 150 mm x 35-mm diameter tubes.
Plastic containers (~4.5" x 3.5"; equivalent to the size of a container for 200-ml pipet tips)
Heating block or water bath
Shaker
Plastic sheet protectors or overhead transparencies
Hyperfilm ECL (Amersham, Cat.# RPN3114K) OR Chemiluminescence imaging system (e.g., FluorChem from Alpha Innotech Corp.)

Preparing Nuclear Extract from Cells or Tissues

Nuclear extracts can be prepared using the method described by Dignam *et al.* (1). Alternatively, you can use a commercially available kit, such as Panomics' Nuclear Extraction Kit (Cat. # AY2002).

Preparing TF-Bound DNA

- For each nuclear extract sample, combine the following components into a sterile 1.5-ml microcentrifuge tube. Prepare an additional sample for your negative control.

Nuclear Extract*	10 µl
Poly d(I-C)	10 µl
TF-TF Probe Mix	10 µl
5X Binding Buffer	15 µl
dH ₂ O (RNase, DNase free)**	30 µl
Total Volume	75 µl

*Use Control Nuclear Extract for positive control. We recommend using 30-100 µg of nuclear extract per reaction.

**Adjust the volume of dH₂O based on the concentration of your nuclear extract.

- Mix well by pipeting.
- Incubate samples at 15°C for 30 min. Afterwards, incubate on ice for an additional 30 min.

Immuno-precipitation & Elution

BEFORE YOU START: Chill all solutions to 4°C prior to use. Dilute solutions accordingly (see Appendix A). If you plan to proceed directly to overnight hybridization, we recommend that you first prehybridize your membranes (see Section 7.1, 7.2) before beginning with Section 6. **Estimated time: 2–3 hours.**

Using a specific antibody against your TF of interest, you will pull out that TF and any associated TFs, along with corresponding cis-elements. Wash away free cis-elements and nonspecific binding proteins; then, elute the cis-elements for use as probe.

NOTE: Antibodies bind with different affinities to protein A or protein G. We recommend that you use the chart below to choose the best affinity Dynabeads for your specific antibody. **Note:** Do not use agarose or sepharose beads as these bind non-specifically to DNA, providing false positive results.

Ig origin	Dynabeads
mouse IgG1	protein G
mouse Ig2a, Ig2b	protein A*
mouse IgG3	protein G
rabbit IgG1	protein A* or protein G
goat IgG	protein G

*Add 0.1% BSA to 1X IP Dilution Buffer and 1X IP Wash Buffer for protein stability.

1. To the tube prepared in step 5, add:

IP Dilution Buffer	200 μ l
Antibody for your TF*	(~200 ug/ml), 10 μ l

Note: The amount of antibody for immunoprecipitation is variable; adjust the amount accordingly.

*For your negative control, replace antibody with normal IgG. For your positive control, replace antibody with p53 Goat Polyclonal, IgG.

- 2 Mix well by gently tapping the tube. Place sample on rocking platform and rock gently for 90 min at 4°C.
- 3 . Wash Dynabeads
 - 3.1 Resuspend the Dynabeads thoroughly in vial.
 - 3.2 Transfer 75 μ l of the Dynabead solution to 1.5 ml microcentrifuge tube.
 - 3.3 Place tubes in magnetic stand for one minute (beads will precipitate out). Pipet off fluid and discard.
 - 3.4 Equilibrate the beads by resuspending them in 200 μ l of 1X IP Dilution Buffer.
 - 3.5 Place tube in magnetic stand for one minute. Carefully, remove and discard the supernatant without disturbing the beads attached to the wall of the tube.
4. Transfer the DNA/protein/antibody complex (prepared in Step 6.2) to the tube containing the washed beads. Rock gently on the rocking platform at 4°C for 1 hr.
5. Collect beads:
 - 5.1 Briefly spin tube in a microcentrifuge at 2,000 rpm for 2 sec.
 - 5.2 Place tube in magnetic stand for 2 min to collect beads. Keeping tubes on magnetic stand, pipet off fluid and discard (do not disturb pellet).
6. Wash beads with 1X IP Wash Buffer:
 - 6.1 Add 400 μ l of pre-chilled 1X IP Wash Buffer to tube containing beads. Completely resuspend the beads by gently tapping the tube. **IMPORTANT:** Gently **invert** the tube four times to clean the sides of the tube.
 - 6.2 Spin the tube in a microcentrifuge at 2,000 rpm for 2 sec.
 - 6.3 Place the tube in magnetic stand for 2 min to collect beads. Pipet off fluid

and discard. **Important: Be sure to pipet off as much fluid as possible in this step.**

6.4 Repeat Step 6.6.1–3 three times.

7. Elute probe from beads:

7.1 Preheat heat block to 100°C. Heat the Elution Buffer to dissolve precipitates before use.

7.2 Add 60 µl of 1X IP Elution Buffer to tube containing beads. Gently mix.

7.3 Incubate at 100°C for 5 min.

7.4 Transfer to ice immediately. Keep on ice for 2 min.

7.5 Spin tube in a microcentrifuge at 2,000 rpm for 5 sec.

7.6 Place tube in magnetic stand for 2 min to collect beads on ice (i.e., place magnetic stand on ice with tubes).

7.7 **Important:** Transfer the fluid to a clean tube and place tube on ice. This is your eluted probe. The beads can be discarded. If you do not plan to use the probe immediately, store the probe at –20°C. Prior to use, thaw probe, and repeat steps 6.7.3 and 6.7.4, to ensure complete denaturation of double-stranded DNA.

Hybridization In this Section, you will hybridize the eluted probe (prepared in Section 6) to the array membrane.

1. Place each array membrane into a hybridization bottle. To wet the membrane, fill bottle with deionized H₂O. Then, carefully decant water.
2. To each hybridization bottle that contains an array membrane, add 3–5 ml of prewarmed Hybridization Buffer (provided). Place each bottle in the hybridization oven at 42°C for 2 hr.
3. Add the eluted probe from step 6.7.7 to each hybridization bottle and hybridize at 42°C overnight.
4. Decant the hybridization mixture from each hybridization bottle, and wash each membrane:
 - 4.1 Add 50 ml of prewarmed Hybridization Wash I*, incubate at 48°C for 20 min in a rotating hybridization oven. Decant liquid and repeat wash.
 - 4.2 Add 50 ml of prewarmed Hybridization Wash II*, incubate at 48°C for 20 min in a rotating hybridization oven. Decant liquid and repeat wash.

* See Appendix A for recipes.

Detection **IMPORTANT:** Do not allow the membrane to dry during the detection. Dilute all solutions before use (see Appendix A). Pre-warm the 1X Blocking Solution and 1X Wash Buffer at 37°–50°C until all particulates are dissolved.

1. Using forceps, carefully remove each membrane from its hybridization bottle and transfer to a new container containing 20 ml of 1X Blocking Buffer; each membrane needs its own container. (We use a container that is equivalent to the size of a 200-ml pipet-tip container, ~4.5" x 3.5".)
2. Incubate at room temperature for 30 min with gentle shaking.
3. Transfer 1 ml of Blocking Buffer from each membrane container to a fresh 1.5-ml tube. To each, add 20 µl of 1000x Streptavidin-HRP Conjugate and mix well. Return each mixture to the appropriate container and incubate at room temperature for 20 min.
4. Decant the Blocking Solution containing the antibody. Wash each membrane three times at room temperature with 20 ml of 1X Wash Buffer, each 5 min.
5. Add 20 ml of 1X Detection Buffer to each membrane and incubate for 5 min at room temperature.
6. To prepare the substrate solution, mix 1 ml Substrate Solution I and 1 ml Substrate Solution II. Add 1 ml Substrate Solution III. Mix well.
7. Using a plastic sheet protector, overhead transparency film, or plastic wrap (whichever is readily available), place each membrane on a plastic sheet. Then, pipet 3 ml of substrate solution onto each membrane and overlay each with a second plastic sheet. Ensure that substrate is evenly distributed over the membrane with no air bubbles. Incubate at room temperature for 5 min.
8. Remove excess substrate by gently applying pressure over the top sheet. Using a paper towel, remove excess substrate that might be remaining on the surface of the sheets. Expose the membranes using either Hyperfilm ECL or a chemiluminescence imaging system, such as the FluorChem imager from Alpha Innotech Corp. In either case, we recommend that you try several different exposures of varying lengths of time (e.g., 2–5 min). Overexposure of your blot may result in excessive background

.Results & Analysis The main advantage of the TF-TF Interaction Array is that you can detect the interaction of a protein with multiple transcription factors simultaneously. Follow these guidelines to analyze your results:

1. Acquire the images using either x-ray film or a chemiluminescent imaging system.
2. Adjust the exposure time such that the biotin spots (along the right and bottom sides of the membrane) have equal signal intensity.
3. Compare the images between negative control (normal IgG) and experimental. All the spots that appear in your experimental should be positive; in contrast, the negative control should not have any spots. If the image of the negative control has some weak spots, then you'll need to check whether the same spots also appear in the experimental. If so, the spots in the experimental may represent false positives. To confirm your results, perform co-immunoprecipitation or super-gel shift assays.

Troubleshooting Guide

Problem	Cause	Recommendation
Uneven background	Substrate is not evenly distributed on the membrane.	Do not use thin cling wrap materials during the overlay procedure
		Re-detect with substrate that is evenly covering the membrane surface.
High Background	Incubation with substrate is too long.	Expose longer—10 to 15 min
		Incubation should not exceed 5 min.
Signal is too weak.	The yield of recovered DNA probe is low.	Confirm using control nuclear extract (provided).
		Check the concentration of nuclear extract.
No Signal	The antibody may be not suitable for IP	Test the antibody before applying to the array.

References TF-TF, SRC-TF, PPARalpha-TF, NFkB p50-TF, GRalpha-TF

1. Dignam, J.D., Lebovitz, R.M., and Roeder, R.G. (1983) Accurate transcription initiation by RNA polymerase II in a soluble extract from isolated mammalian nuclei. *Nucleic Acids Research* 11:1475–1489.

p300-TF

1. Janknecht, R. and Hunter, T. (1996) Transcription control: versatile molecular glue. *Curr Biol* 6:22-23.
2. Janknecht, R. and Hunter, T. (1996) MAP kinase-dependent transcriptional coactivation by Elk-1 and its cofactor CBP. *Biochem Biophys Res Commun* 228:831-837.
3. Shikama, N., Lyon, J., and La Thangue, N.B. (1997) The p300/CBP family: integrating signals with transcription factors and chromatin. *Trends Cell Biol* 7:230-236.
4. Goodman, R.H. and Smolik, S. (2000) CBP/p300 in cell growth, transformation,

- and development. *Genes Dev* **14**:1553-1577.
5. Chan, H.M. and La Thangue, N.B. (2001) p300/CBP proteins: HAT s for transcriptional bridges and scaffolds. *J Cell Science* **114**:2363-2373.
 6. Xu, L., Glass, C.K., and Rosenfeld, M.G. (1999) Co-activator and corepressor complexes in nuclear receptor function. *Curr Opin Genet Dev* **9**:140-147.
 7. Chan, H.M. and La Thangue, N.B. (2001) p300/CBP proteins: HAT s for transcriptional bridges and scaffolds. *J Cell Science* **114**:2363-2373.
 8. Ogryzko, V.V., Schiltz, R.L., Russanova, V. *et al.* (1996) The transcriptional co-activators p300 and CBP are histone acetyltransferases. *Cell* **87**:1107-1112.
 9. Korzus, E., Torchia, J., Rose, D.W., *et al.* (1998) Transcription factor-specific recruitment of co-activators and their acetyltransferase functions. *Science* **279**:703-707.
 10. Moran, E. (1993) DNA tumour virus transforming proteins and the cell cycle. *Curr Opin Genet Dev* **3**:63-70.
 11. Giles, R.H., Peters, D.J.M. and Breuning, M.H. (1998) Conjunction dysfunction: CBP/p300 in human disease. *Trends Genet* **14**: 178-183.

APPENDIX A:**Recipes &
instructions for
diluting stock
solutions****SECTION 6**

- **15 ml of 1X IP Wash Buffer**

To 12 ml of deionized H₂O, add 3 ml of 5X IP Wash Buffer (provided).

SECTION 7

- **300 ml of 2X SSC/0.5% SDS (Hybridization Wash I)**

To 262.5 ml of deionized H₂O, add 30 ml of 20X SSC (provided) and 7.5 ml of 20% SDS (provided). Mix well.

- **300 ml of 0.1X SSC/0.5% SDS (Hybridization Wash II)**

To 291 ml of deionized H₂O, add 1.5 ml of 20X SSC (provided) and 7.5 ml of 20% SDS (provided). Mix well.

SECTION 8

- **60 ml of 1X Blocking Buffer**

To 30 ml of deionized H₂O, add 30 ml of 2X Blocking Buffer (provided). Mix well and store at 4°C.

- **200 ml of 1X Wash Buffer**

To 150 ml of deionized H₂O, add 50 ml of 4X Wash Buffer (provided). Mix well.

- **100 ml of 1X Detection Buffer**

To 90 ml of deionized H₂O, add 10 ml of 10X Detection Buffer (provided). Mix well.

APPENDIX B:

Schematic diagram of the Interaction Array I

LEGEND

	1	2
A	AP1	AP1
B	AP1	AP1

NORMAL
DILUTED 1:10

A	AP1 ₍₁₎	AP1 ₍₁₎	AP2 ₍₁₎	AP2 ₍₁₎	ARE	ARE	Bmi3	Bmi3	C/EBP	C/EBP	CBF	CBF	CDP	CDP	c-Myb	c-Myb	AP1 ₍₁₎	AP1 ₍₁₎	A
B	AP1 ₍₁₎	AP1 ₍₁₎	AP2 ₍₁₎	AP2 ₍₁₎	ARE	ARE	Bmi3	Bmi3	C/EBP	C/EBP	CBF	CBF	CDP	CDP	c-Myb	c-Myb	AP1 ₍₁₎	AP1 ₍₁₎	B
C	CREB ₍₁₎	CREB ₍₁₎	E2F1	E2F1	EGR ₍₁₎	EGR ₍₁₎	ERE	ERE	Ets	Ets	Ets/PEA3	Ets/PEA3	FAST1	FAST1	GAS/SRE	GAS/SRE	AP2 ₍₁₎	AP2 ₍₁₎	C
D	CREB ₍₁₎	CREB ₍₁₎	E2F1	E2F1	EGR ₍₁₎	EGR ₍₁₎	ERE	ERE	Ets	Ets	Ets/PEA3	Ets/PEA3	FAST1	FAST1	GAS/SRE	GAS/SRE	AP2 ₍₁₎	AP2 ₍₁₎	D
E	GATA	GATA	GRE	GRE	HNF4	HNF4	IRF-1	IRF-1	MEF-1	MEF-1	MEF-2	MEF-2	MyelMax	MyelMax	NF-1	NF-1			E
F	GATA	GATA	GRE	GRE	HNF4	HNF4	IRF-1	IRF-1	MEF-1	MEF-1	MEF-2	MEF-2	MyelMax	MyelMax	NF-1	NF-1			F
G	NFATc	NFATc	NF-E1	NF-E1	NF-E2	NF-E2	NFκB	NFκB	Oct-1	Oct-1	p53	p53	Pax-5	Pax-5	Pbx1	Pbx1			G
H	NFATc	NFATc	NF-E1	NF-E1	NF-E2	NF-E2	NFκB	NFκB	Oct-1	Oct-1	p53	p53	Pax-5	Pax-5	Pbx1	Pbx1			H
I	Pit1	Pit1	PPAR	PPAR	PRE	PRE	RAR(DR5)	RAR(DR5)	RXR(DR1)	RXR(DR1)	SIE	SIE	Smad SBE	Smad SBE	Smad 3/4	Smad 3/4			I
J	Pit1	Pit1	PPAR	PPAR	PRE	PRE	RAR(DR5)	RAR(DR5)	RXR(DR1)	RXR(DR1)	SIE	SIE	Smad SBE	Smad SBE	Smad 3/4	Smad 3/4			J
K	Sp1	Sp1	SRE	SRE	Stat1	Stat1	Stat3	Stat3	Stat4	Stat4	Stat5	Stat5	Stat6	Stat6	TFIID	TFIID			K
L	Sp1	Sp1	SRE	SRE	Stat1	Stat1	Stat3	Stat3	Stat4	Stat4	Stat5	Stat5	Stat6	Stat6	TFIID	TFIID			L
M	TR	TR	TR(DR-4)	TR(DR-4)	USF-1	USF-1	VDR(DR3)	VDR(DR3)	HSE	HSE	MRE	MRE							M
N	TR	TR	TR(DR-4)	TR(DR-4)	USF-1	USF-1	VDR(DR3)	VDR(DR3)	HSE	HSE	MRE	MRE							N
0																			0
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	

Figure 2. Schematic diagram of the Interaction Array I. The genes on the array are spotted in duplicate: the first row is DNA spotted normally, the second row is DNA diluted 1:10. The dark grey columns along the right and bottom sides of the array indicate where biotinylated DNA has been spotted. These spots are intended for alignment. (Note that the notch is at the top, right-hand corner.)

Schematic
diagram of the
Interaction Array II

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
A	ADR1	ADR1	Ahr/ Arnt	Ahr/ Arnt	ANG-IRE	ANG-IRE	Antiox. RE	Antiox. RE	AP3	AP3	AP4 (RIP60)	AP4 (RIP60)	β -resp. element	β -resp. element	Cdx4/ NKX2	Cdx4/ NKX2	
B	CEFI	CEFI	CEF2	CEF2	CEF/ CRE	CEF/ CRE	COUP- TF	COUP- TF	CREB- BP1	CREB- BP1	c-Rel	c-Rel	DR-0	DR-0	DR-1	DR-1	
C	DR-2	DR-2	DR-3	DR-3	DR-4	DR-4	DR-5	DR-5	E47	E47	E4BP4	E4BP4	E4F ATF	E4F ATF	Elk1	Elk1	
D	EVI-1	EVI-1	FKHR (human)	FKHR (human)	FKhr (mouse)	FKhr (mouse)	Afxh (Foxo4)	Afxh (Foxo4)	Freac-2 (1)	Freac-2 (1)	Freac-2 (2)	Freac-2 (2)	Freac-4	Freac-4	Freac-7	Freac-7	
E	GAG	GAG	GAS	GAS	GATA-1	GATA-1	GATA- 1/2	GATA- 1/2	GATA-2	GATA-2	GATA-3	GATA-3	GATA-4	GATA-4	GATA-6	GATA-6	
F	Gfi-1	Gfi-1	HAS	HAS	HBS	HBS	HBS + HAS	HBS + HAS	HFH-1	HFH-1	HFH-2	HFH-2	HFH-3	HFH-3	HFH-8	HFH-8	
G	HIF-1	HIF-1	HNF-3	HNF-3	HNF-3 β	HNF-3 β	HNF-4	HNF-4	HNF-4/ COUP-TF	HNF-4/ COUP-TF	Ikaros	Ikaros	IRF-1/2	IRF-1/2	ISRE(1)	ISRE(1)	
H	ISRE(2)	ISRE(2)	ISRE (Transf)	ISRE (Transf)	L-III BP	L-III BP	KKLF	KKLF	MEF-2(1)	MEF-2(1)	MEF-2(2)	MEF-2(2)	MEF-3	MEF-3	MSPT	MSPT	
I	MT-Box	MT-Box	MUSFT	MUSFT	MZF1	MZF1	NF-Y	NF-Y	NKx-2.5	NKx-2.5	NPAS2	NPAS2	HBS/ xbp-1	HBS/ xbp-1	NZF-3	NZF-3	
J	ORE	ORE	p300	p300	PARP	PARP	Pax2	Pax2	Pax3	Pax3	Pax4	Pax4	Pax6	Pax6	Pax8	Pax8	
K	PPAR α	PPAR α	PPAR γ	PPAR γ	RREB(1)	RREB(1)	RREB(2)	RREB(2)	RSRF4	RSRF4	SAA	SAA	Skn	Skn	SRY	SRY	
L	Tax/ CREB	Tax/ CREB	TCF/LEF	TCF/LEF	XBP-1	XBP-1	XRE	XRE	ZIA	ZIA	ZIB	ZIB	ZIC	ZIC	ZID	ZID	
M																	
	1	2	3	4	5	6	7	8	8	10	11	12	13	14	15	16	17

Figure 3. Schematic diagram of the Interaction Array II. The dark grey columns along the right and bottom sides of the array indicate where biotinylated DNA has been spotted. These spots are intended for alignment. (Note that the notch is at the top, right-hand corner.)


Schematic
diagram of the
Interaction Array
III

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
A	ADD1	ADD1	AF1/ ARP1	AF1/ ARP1	AFP1	AFP1	AIC/ CBF	AIC/ CBF	AML1	AML1	AP2/ YY1	AP2/ YY1	AREB6	AREB6	ARP	ARP	A
B	ATF/ CRE	ATF/ CRE	ATF2	ATF2	CACC	CACC	CCAAT	CCAAT	CCAC	CCAC	CD28RC	CD28RC	CD28RC/ NF-IL2B	CD28RC/ NF-IL2B	Cdx2	Cdx2	B
C	CEA	CEA	CP1/CTF/ CBF	CP1/CTF/ CBF	CPE	CPE	CREB2	CREB2	CSBP	CSBP	CTCF	CTCF	E12	E12	EBP40/ 45	EBP40/ 45	C
D	EKLF(1)	EKLF(1)	EKLF(2)	EKLF(2)	Elf	Elf	H4TF	H4TF	HEN1	HEN1	HFH-8/ HNF3 LUN	HFH-8/ HNF3 LUN	HNF	HNF	HLF	HLF	D
E	HMG	HMG	HNF-1A	HNF-1A	HOX4C	HOX4C	HOXD8/9 /10	HOXD8/9 /10	HOXD8(1)	HOXD8(1)	HOXD8(2)	HOXD8(2)	HOXD9/ 10	HOXD9/ 10	ICSBP	ICSBP	E
F	Isl-1	Isl-1	KPFI	KPFI	KTP1	KTP1	LCR-F1	LCR-F1	LF-A1	LF-A1	LR1	LR1	LyF(1)	LyF(1)	LyF-1(2)	LyF-1(2)	F
G	LyF-1(2)	LyF-1(2)	MAZ	MAZ	MTB-Zf	MTB-Zf	MTF(1)	MTF(1)	MyoD	MyoD	MyoG	MyoG	MyT1	MyT1	MZF1(1)	MZF1(1)	G
H	MZF1(2)	MZF1(2)	NF-4FA	NF-4FA	NF-Gma (1)	NF-Gma (1)	NF-Gma (2)	NF-Gma (2)	NFIL-2	NFIL-2	NF-Y	NF-Y	NRF-1	NRF-1	p53(2)	p53(2)	H
I	p55	p55	PEA3 HIP	PEA3 HIP	PEBP	PEBP	PEBP2	PEBP2	PPUR(1)	PPUR(1)	PPUR(2)	PPUR(2)	Pur-1	Pur-1	PYR	PYR	I
J	RB	RB	REF1/2 /3(1)	REF1/2 /3(1)	RIPE3a1	RIPE3a1	SIF1	SIF1	SIF2	SIF2	SIF3	SIF3	SP1 ASP	SP1 ASP	SRF SAP	SRF SAP	J
K	SYR	SYR	TCE	TCE	TFE3	TFE3	Tf-LF	Tf-LF	TGT3	TGT3	TIF1	TIF1	TRE1/2	TRE1/2	v-Maf	v-Maf	K
L	WT1(1)	WT1(1)	WT1(2)	WT1(2)	WT1(3)	WT1(3)	X2BP	X2BP	XBP1/ X2BP	XBP1/ X2BP	ZNF174	ZNF174					L
M																	M
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17

Figure 4. Schematic diagram of the Interaction Array III. The dark grey columns along the right and bottom sides of the array indicate where biotinylated DNA has been spotted. These spots are intended for alignment. (Note that the notch is at the top, right-hand corner.)

APPENDIX E:**TF-TF Protein
Interaction Array
Grid**

Photocopy this page onto separate piece of paper or transparency film. (Note that the notch is closest to A24.)



APPENDIX F:

Stripping Procedure for Arrays

Note: We do not encourage stripping the array membranes more than two times.

Procedure

1. Wash membranes in 0.4M NaOH at 45°C for 30 min .
2. Wash membranes in 0.2M Tris-HCL, pH 7.6; 0.1X SSC, and 0.1% SDS at 45°C for 15 min .
3. To ensure that stripping was successful, run it through the standard chemiluminescence detection procedure as described in this user manual.
4. After detection, wash the membrane in 1X Washing Buffer at 42°C for 30 min.
5. Membranes are ready for hybridization or dry the membrane in an 80°C incubator for later use.

Contacting Panomics

For ordering information or technical support, contact the appropriate resource provided below according to your geographical location.

Location	U.S. Corporate Headquarters
Address	Panomics, Inc. 6519 Dumbarton Circle Fremont CA 94555 USA
Telephone	1.510.818.2600
FAX	1.510.818.2610
Email	info@panomics.com
Technical Support	techsupport@panomics.com or 1.877.726.6642 option 3
Ordering Information	orders@panomics.com

Location	European Headquarters
Address	Panomics Srl Via Sardegna 1 20060 Vignate Milano Italy
Telephone	+39.02.95.360.250
FAX	+39.02.95.360.992
Email	info_europe@panomics.com
Technical Support	techsupport_europe@panomics.com
Ordering Information	order_europe@panomics.com

Location	Asia Pacific Headquarters
Address	Panomics, Inc. 16F Gemdale Plaza Tower A, No. 91 Jianguo Road Beijing 100022 P.R. China
Telephone	+86.10.59208157
FAX	+86.10.59208111
Email	info_asia@panomics.com
Technical Support	techsupport_asia@panomics.com
Ordering Information	order_asia@panomics.com

For an updated list of FAQs and product support literature, visit our website at www.panomics.com.
