

Nanospray Flex Ion Source

Version 1.0

Getting Connected Guide

60053-97126 Revision A November 2010

DOCUMENTATION
SURVEY

© 2010 Thermo Fisher Scientific Inc. All rights reserved.

Nanospray Flex, DirectJunction, EASY-nLC, and EASY-Column are trademarks and Xcalibur is a registered trademark of Thermo Fisher Scientific Inc. in the United States.

The following are registered trademarks in the United States: Nanospray is a registered trademark of AB Sciex Pte. Ltd. Lemo is a registered trademark of Interlemon Holding S.A. Dino-Lite is a registered trademark of Anmo Electronics Corporation. Swagelok is a registered trademark of the Crawford Fitting Company.

Thermo Fisher Scientific San Jose has a non-exclusive license with the University of Texas for US Patent 5,572,023. Caprioli: Electro spray methods and apparatus for trace analysis.

All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Fisher Scientific Inc. provides this document to its customers with a product purchase to use in the product operation. This document is copyright protected and any reproduction of the whole or any part of this document is strictly prohibited, except with the written authorization of Thermo Fisher Scientific Inc.

The contents of this document are subject to change without notice. All technical information in this document is for reference purposes only. System configurations and specifications in this document supersede all previous information received by the purchaser.

Thermo Fisher Scientific Inc. makes no representations that this document is complete, accurate or error-free and assumes no responsibility and will not be liable for any errors, omissions, damage or loss that might result from any use of this document, even if the information in the document is followed properly.

This document is not part of any sales contract between Thermo Fisher Scientific Inc. and a purchaser. This document shall in no way govern or modify any Terms and Conditions of Sale, which Terms and Conditions of Sale shall govern all conflicting information between the two documents.

Release history: Revision A, November 2010

Software version: EASY-nLC 2.5 and later; Xcalibur 1.3 and later

For Research Use Only. Not for use in diagnostic procedures.

Regulatory Compliance

Thermo Fisher Scientific performs complete testing and evaluation of its products to ensure full compliance with applicable domestic and international regulations. When the system is delivered to you, it meets all pertinent electromagnetic compatibility (EMC) and safety standards as described in the next section or sections by product name.

Changes that you make to your system may void compliance with one or more of these EMC and safety standards. Changes to your system include replacing a part or adding components, options, or peripherals not specifically authorized and qualified by Thermo Fisher Scientific. To ensure continued compliance with EMC and safety standards, replacement parts and additional components, options, and peripherals must be ordered from Thermo Fisher Scientific or one of its authorized representatives.

FCC Compliance Statement

THIS DEVICE COMPLIES WITH PART 15 OF THE FCC RULES. OPERATION IS SUBJECT TO THE FOLLOWING TWO CONDITIONS: (1) THIS DEVICE MAY NOT CAUSE HARMFUL INTERFERENCE, AND (2) THIS DEVICE MUST ACCEPT ANY INTERFERENCE RECEIVED, INCLUDING INTERFERENCE THAT MAY CAUSE UNDESIRE OPERATION.



CAUTION Read and understand the various precautionary notes, signs, and symbols contained inside this manual pertaining to the safe use and operation of this product before using the device.

Notice on the Proper Use of Thermo Scientific Instruments

In compliance with international regulations: Use of this instrument in a manner not specified by Thermo Fisher Scientific could impair any protection provided by the instrument.

Notice on the Susceptibility to Electromagnetic Transmissions

Your instrument is designed to work in a controlled electromagnetic environment. Do not use radio frequency transmitters, such as mobile phones, in close proximity to the instrument.

For manufacturing location, see the label on the instrument.

WEEE Compliance

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC. It is marked with the following symbol:



Thermo Fisher Scientific has contracted with one or more recycling or disposal companies in each European Union (EU) Member State, and these companies should dispose of or recycle this product. See www.thermo.com/WEEERoHS for further information on Thermo Fisher Scientific's compliance with these Directives and the recyclers in your country.

WEEE Konformität

Dieses Produkt muss die EU Waste Electrical & Electronic Equipment (WEEE) Richtlinie 2002/96/EC erfüllen. Das Produkt ist durch folgendes Symbol gekennzeichnet:



Thermo Fisher Scientific hat Vereinbarungen mit Verwertungs-/Entsorgungsfirmen in allen EU-Mitgliedsstaaten getroffen, damit dieses Produkt durch diese Firmen wiederverwertet oder entsorgt werden kann. Mehr Information über die Einhaltung dieser Anweisungen durch Thermo Fisher Scientific, über die Verwerter, und weitere Hinweise, die nützlich sind, um die Produkte zu identifizieren, die unter diese RoHS Anweisung fallen, finden sie unter www.thermo.com/WEEERoHS.

Conformité DEEE

Ce produit doit être conforme à la directive européenne (2002/96/EC) des Déchets d'Equipements Electriques et Electroniques (DEEE). Il est marqué par le symbole suivant:



Thermo Fisher Scientific s'est associé avec une ou plusieurs compagnies de recyclage dans chaque état membre de l'union européenne et ce produit devrait être collecté ou recyclé par celles-ci. Davantage d'informations sur la conformité de Thermo Fisher Scientific à ces directives, les recycleurs dans votre pays et les informations sur les produits Thermo Fisher Scientific qui peuvent aider la détection des substances sujettes à la directive RoHS sont disponibles sur www.thermo.com/WEEERoHS.

Contents

	Preface	ix
	Related Documentation	ix
	Safety and Special Notices	ix
	Contacting Us	xi
	Disclaimer	xi
	Declaration of Conformity	xii
Chapter 1	Introduction	1
	Content of Nanospray Flex Ion Source Kit	2
	Unpacking the Ion Source and Preinstallation	3
	Installing the Ion Source on the Mass Spectrometer	5
Chapter 2	Electrical Connections	9
	Connecting the Monitor and Cameras	9
	Adjusting the Monitor and Cameras	11
Chapter 3	DirectJunction Adaptor	13
	Choosing a One-Column or Two-Column Configuration	14
	Choosing Stainless Steel or Glass Emitters	15
	Standard Configuration	16
	Assembling the Columns for Standard Configuration	16
	Connecting the Columns to a Stainless Steel Emitter	17
	Mounting the Column and Emitter Setup onto the DirectJunction Adaptor	18
	Connecting the LC System to the Column Setup	19
	Connecting the DirectJunction to HV	19
Chapter 4	Alternative Configurations for DirectJunction	21
	Stainless Steel Emitter and One-Column Setup	21
	Glass Emitter with Liquid Junction and One-Column Setup	22
	Glass Emitter with Liquid Junction and Two-Column Setup	22
	Packed Glass Emitter with Liquid Junction in One-Column Setup	23
	Packed Glass Emitter with Liquid Junction in Two-Column Setup	24
	Setup with a Long Analytical Column	25

Chapter 5	Nano Flex Ion Source Configuration	27
Chapter 6	Operation	33
	Positioning the Emitter Tip	33
	Maintaining Stable Spray	35
	Using the Optional Gas Connection	36
Appendix A	Accessories and Spare Parts	37
	Nano ES Spray Kit for Offline Analysis	37
	Accessories List	38
	Limited Warranty	38
	Additional Support	38
	Index	39

Preface

This guide provides information on the Thermo™ Nanospray Flex™ ion source, including setup, installation, and operation.

The Nanospray Flex ion source works in conjunction with these Thermo Scientific mass spectrometers:

- LTQ™ Series: LTQ, Orbitrap™, FT™, Velos™
- Exactive™
- TSQ™ Series: Quantum™, Quantum Access™, Quantum Discovery MAX™, Quantum Ultra™, Vantage™
- LCQ™ Deca XP MAX

Contents

- [Related Documentation](#)
- [Safety and Special Notices](#)
- [Contacting Us](#)
- [Disclaimer](#)
- [Declaration of Conformity](#)

Related Documentation

In addition to this guide, Thermo Fisher Scientific provides the following documents for the Nanospray Flex ion source available as PDF files:

- Thermo Xcalibur™ manual suite
- *Ion Max Source and Ion Max-S API Source Hardware Manual*

The software also provides Help.

Safety and Special Notices

Make sure you follow the precautionary statements presented in this guide. The safety and other special notices appear in boxes.

Safety and special notices include the following:



CAUTION Highlights hazards to humans, property, or the environment. Each CAUTION notice is accompanied by an appropriate CAUTION symbol.

IMPORTANT Highlights information necessary to prevent damage to software, loss of data, or invalid test results; or might contain information that is critical for optimal performance of the system.

Note Highlights information of general interest.

Tip Highlights helpful information that can make a task easier.



CAUTION

- Emitter tips are extremely sharp. Never touch the tips as this may cause injury.
- Always depressurize the LC-system with transfer line before removing emitters from the source head. Failure to do so may result in the emitter being ejected at high speed.



CAUTION

- Although the ion source is shielded, in certain cases there might be access to the emitter which is on high voltage (HV). Always turn off the HV before touching the source head.
- Never connect a power supply able to deliver more than 8kV and 100 μ A to this Nanospray source.
- Do not leave the source unattended while the spray voltage is on.



CAUTION The interface region, especially the cone of the mass spectrometer can be hot. Do not touch this area.

Contacting Us

There are several ways to contact Thermo Fisher Scientific for the information you need.

❖ To contact Technical Support

Phone	800-532-4752
Fax	561-688-8736
E-mail	us.techsupport.analyze@thermofisher.com
Knowledge base	www.thermokb.com

Find software updates and utilities to download at mssupport.thermo.com.

❖ To contact Customer Service for ordering information

Phone	800-532-4752
Fax	561-688-8731
E-mail	us.customer-support.analyze@thermofisher.com
Web site	www.thermo.com/ms

❖ To get local contact information for sales or service

Go to www.thermoscientific.com/wps/portal/ts/contactus.

❖ To copy manuals from the Internet

Go to mssupport.thermo.com, agree to the Terms and Conditions, and then click **Customer Manuals** in the left margin of the window.

❖ To suggest changes to documentation or to Help

- Fill out a reader survey online at www.surveymonkey.com/s/PQM6P62.
- Send an e-mail message to the Technical Publications Editor at techpubs-lcms@thermofisher.com.

Disclaimer

- This equipment, including spray capillaries, is designed specifically for creating ions in a mass spectrometer. (No other use is recommended.)
- This equipment is for laboratory use only and is not a medical device.
- This equipment must not be connected to any power supply delivering more than 8 kV and 100 μ A.
- Do not allow unauthorized or untrained operators to use this equipment.

- Any misuse will be the sole responsibility of the user/owner. Thermo Fisher Scientific assumes no implied or inferred liability for direct or consequential damages or injuries from this instrumentation if it is operated or used in any way other than that for which it is designed.
- Any legal disagreements must be settled before a Danish court of law.
- Only use this equipment if you agree with the conditions above.

Declaration of Conformity

European Safety Standards

Declaration of Conformity

Standard(s) to which conformity is declared:	IEC 61010-1/ EN61010-1:2001, Second edition - "Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory Use."
Manufacturer's Name:	Proxeon Biosystems A/S
Manufacturer's Address	Edisonsvej 4, DK-5000 Odense, Denmark
Type of Equipment:	Laboratory Instrumentation
Model Name:	Nanospray Flex ion source
Model Numbers:	ES071
Serial Number:	ES-001000 and later
Year of Manufacture:	2010–

I, the undersigned hereby declare that the equipment specified above conforms to the above Directive(s) and Standard(s).



Ole Vorm, Site Manager

October 11, 2010

Introduction

Nanospray Flex ion sources minimize the liquid flow rate (nL/min) and maintain excellent spray stability to ensure evaporation and ionization of liquid samples —the key to achieving the highest sensitivity.

Key benefits of the Nanospray Flex ion source include

- User-friendly design
- Single setup for all online nanoflow applications
- Ability to interface with online nanoscale LC separation techniques


Contents

- [Content of Nanospray Flex Ion Source Kit](#)
- [Unpacking the Ion Source and Preinstallation](#)
- [Installing the Ion Source on the Mass Spectrometer](#)

1 Introduction

Content of Nanospray Flex Ion Source Kit

Figure 1. Nanospray Flex ion source



Content of Nanospray Flex Ion Source Kit

The Nanospray Flex Ion Source Kit (P/N ES071) includes the following:


- Flexible ion source housing with precision XYZ manipulator and fittings
- DirectJunction™ adaptor ready for use with a variety of different column/emitter configurations for online analysis
- LiquidJunction MicroCross for one column setup with long packed emitters
- Dual camera setup including one LCD monitor
- Nanobore stainless steel emitters fitting the DirectJunction adaptor

- All necessary cables and connectors for the ion source
- *Nanospray Flex Ion Source Getting Connected Guide*

Unpacking the Ion Source and Preinstallation

The Nanospray source kit is packed in a single box. To fully protect the components during transport, the source mechanics, DirectJunction adaptor, and cameras are not assembled (see [Figure 2](#)).

Figure 2. Nanospray ion source components



❖ To unpack the ion source


1. After removing the foam from the top of the box, access the accessories for the Nanospray source.
2. Carefully lift up the next foam bar including the accessories to access the body of the Nanospray source.
3. Carefully take out the body of the Nanospray source and place it on a table to complete the assembly of the Nanospray source.

1 Introduction

Unpacking the Ion Source and Preinstallation

4. Pull back the drawer of the Nanospray source to access the arm on the XYZ manipulator (see [Figure 3](#)). The DirectJunction adaptor has to be mounted to this arm using the 5 mm. screw in the holder of the DirectJunction.

Figure 3. Assembly of the DirectJunction to the XYZ manipulator arm




5. Use the supplied screwdriver to mount the DirectJunction adaptor on the XYZ manipulator arm. Before tightening the 5 mm screw, position the DirectJunction with the back of the DirectJunction elevated about 20° above horizontal. Tighten the screw firmly.

The two white boxes in the box contain Dino-Lite™ cameras, which are located in the opening on the top and left side of the source body.

6. Loosen the small screw which is used to secure the cameras in place with the supplied screwdriver (see [Figure 4](#)).

Figure 4.



Small screw to fix
camera in place



Dino-Lite camera placed on
top of the source body

7. Remove the cameras and install the first one in the top of the source body. Note that the camera can only be placed one way, ensuring optimal alignment. Slowly tighten the screw until it touches the camera housing. Tighten by another half turn.
8. Repeat the procedure for the second camera, which is located on the left side of the source body.

You are now ready to install the Nanospray source.


Installing the Ion Source on the Mass Spectrometer

After assembling the body and accessories, you are ready to mount the completed source on the mass spectrometer.

1 Introduction

Installing the Ion Source on the Mass Spectrometer


Figure 5. Thermo Nanospray Flex ion source with DirectJunction



❖ To install the ion source on the mass spectrometer

1. Before removing the previous ion source, disconnect any columns and flow lines from the HPLC system.
2. Refer to the mass spectrometer user manual for guidance on how to remove the existing ion source.
3. Remove the existing ion source as follows:
 - a. Turn the two black handles on top of the source in the direction indicated in [Figure 6](#). The handles might give some resistance.
 - b. Carefully pull out the source.

Figure 6. Turning the black handles on top of the ion source




4. Take the Nanospray source and align the guiding rods of the source with the guiding holes of the mass spectrometer (see [Figure 7](#)).
5. Push the source gently towards the mass spectrometer. The top part with the connectors might give some resistance. If so, push a little harder in the top area of the source.
6. Close the two handles by turning them so that they point toward each other (you might feel some resistance).

The mass spectrometer is now in Standby mode and the source is mounted.

1 Introduction

Installing the Ion Source on the Mass Spectrometer

Figure 7. Guiding holes on the Nanospray source for mounting on the guiding rods of the mass spectrometer



Electrical Connections

For correct operation, use these procedures to connect the cameras and monitor electrically and to adjust the focus for sharp pictures.

Contents

- [Connecting the Monitor and Cameras](#)
- [Adjusting the Monitor and Cameras](#)



CAUTION When you push the XYZ manipulator onto the mounting frame of the source, the HV for the spray emitter switches on. For your protection, always retract the drawer of the source before working with the source head.

You must use the supplied HV cables and connectors in the Nanospray Flex ion source.


Connecting the Monitor and Cameras

After installing the Thermo Nanospray Flex ion source, you must set up the LCD monitor and Dino-Lite cameras.

❖ To connect and adjust the LCD monitor and Dino-Lite cameras

1. In the ion source kit, remove the monitor from the brown box and place it on top of the mass spectrometer by using the support arm on the back of the monitor (see [Figure 8](#).)


Figure 8. Back of LCD monitor



2. Connect the 12 V power supply connector to the monitor.
3. Locate the two BNC adaptors in the ion source kit and, using a quarter-turn, insert them into the two Video-IN ports on the back of the monitor (Figure 8).
4. Connect the two yellow video connectors from the Dino-Lite cameras onto the two BNC adaptors.
5. Choose and insert the appropriate connector to the power supply:
 - a. Locate the power supply for the cameras in the camera box.

The power supply comes equipped with an EU standard net connector. The kit box contains one U.K. and one U.S. connector.
 - b. Connect the low-voltage connectors of the two power supplies to the two Dino-Lite camera power connectors (see Figure 9).

Figure 9. Connection of the low-voltage connectors to the camera power connectors



6. Connect the power adaptors for the two Dino-Lite cameras and the power adaptor for the monitor to your net power sockets.


The cameras and the monitor are now ready for use.

Adjusting the Monitor and Cameras

❖ To adjust the LCD monitor and Dino-Lite cameras

1. Turn on the LCD monitor by pressing the POWER button on the front of the monitor (see [Figure 10](#)).

Figure 10. LCD Monitor setup and control panel




After a few seconds the picture from one of the two cameras appears. The picture is likely to be very blurred.

2 Electrical Connections


Adjusting the Monitor and Cameras

2. Adjust the focus with the focusing wheel (Figure 4).

Figure 11. Camera light source switch and focusing wheel



The adjustments produce a sharp image on the monitor:



3. Switch to the second camera by pressing the SOURCE button (see Figure 10 on page 11), and adjust the focus with the focusing wheel as above.
4. If the picture on the LCD monitor is too dark, use the camera's internal light to make adjustments.
 - a. To turn on the light, gently press the small button on the back of either camera (Figure 4).
 - b. Adjust the picture to achieve the sharpest image by using either camera light or both lights together.

DirectJunction Adaptor

DirectJunction is a configurable adaptor for the Thermo Nanospray Flex ion source (see [Figure 12](#)). It enables the flexible interchange between most online column and emitter configurations, and uses stainless steel or glass emitters and a one-column or two-column setup, with or without a liquid junction.

Contents

- [Choosing a One-Column or Two-Column Configuration](#)
- [Choosing Stainless Steel or Glass Emitters](#)
- [Standard Configuration](#)
- [Connecting the DirectJunction to HV](#)


3 DirectJunction Adaptor

Choosing a One-Column or Two-Column Configuration

Figure 12. DirectJunction adaptor

HV clamp for stainless steel or coated glass emitters

ZDV union for connecting the column and emitter



HV protection shield

Mounting bracket for the Nanospray ion source

Multipurpose adaptor for ZDV unions, liquid junctions, and connecting tees

Choosing a One-Column or Two-Column Configuration

The DirectJunction adaptor supports both a one- and two-column configuration.

A one-column configuration has the following advantages:

- It has fewer connections, which minimize any potential peak broadening that results from dead volumes.
- It enables MS analysis of compounds that elute during sample loading (as they might not bind to the column material) because the fluid path leads directly to the MS.

A two-column setup (typically a short pre/trap-column with a large inner diameter and a longer analytical column with a narrower inner diameter) has these advantages:

- A pre/trap-column provides an increased loading capacity and an increased loading rate when compared to loading directly onto the analytical column.
- It directs the loading solvent to waste so that sample “contaminants” do not enter the MS, reducing the need to clean the MS inlet.
- A pre/trap-column acts as a guard column by protecting the analytical column from particulate matter.

Choosing Stainless Steel or Glass Emitters

Stainless steel emitters are more robust than glass emitters, and help to maintain a stable and consistent spray for longer periods of time.

The DirectJunction adaptor enables high voltage (HV) to be applied directly to the stainless steel emitter through the small HV clamp at the front of the adaptor. By applying voltage across the entire emitter, you can maintain a constant optimal electrical contact with the liquid. The result is superior spray stability.

A cross section of scientists consider glass emitters to be more bio-inert than stainless steel emitters. This state means reducing the risk of non-specific adsorption of biomolecules and leads to slightly improved sensitivity. However, the degree of this adsorption depends on the chemical characteristics of the sample.

You can pull glass emitters to produce very small ID emitter openings. However, the very small opening at the tip often results in stability problems associated with blockage of the tip, and means that glass emitters rarely last as long as their steel counterparts.

Glass emitters usually have HV applied through a liquid junction. This means that only a small percentage of the bulk liquid flow is exposed to the electrode carrying the HV.


Tip When using glass emitters, do the following:

- Place the liquid junction at the high-pressure side of the column. Electrochemical processes occurring at the electrode can otherwise create gas that leads to spray instability.
- Use glass emitters with small ID emitter tips (<20 µm) to create back pressure in the emitter and avoid outgassing and consequent spray instability.

Standard Configuration

Figure 13 shows the recommended configuration for the EASY-nLC™ nano-flow liquid chromatography system and EASY-Columns™ with a stainless steel emitter installed on the DirectJunction prior to mounting onto the Nanospray Flex ion source.

Figure 13. Configuration for the EASY-nLC and EASY-Columns on DirectJunction using a stainless steel emitter



Assembling the Columns for Standard Configuration

❖ **To assemble the pre-column and analytical column**

1. Remove both red nuts from the ZDV union (SC600).
2. Screw the white blind plug into the union.


Figure 14. ZDV union with white blind plug (SC600)



3. Insert the pre-column (SC001) through a 2 cm (0.79 in.) 1/32 in. OD sleeve (F-385) with arrows on the label pointing away from the sleeve.

4. Insert the sleeved pre-column through a red nut.
5. Screw in the red nut while ensuring that the pre-column and sleeve are firmly pushed against the white blind plug.
6. Insert the other end of the pre-column through a 1/32 in. OD sleeve (F-385), ferrule, and nut from the venting tee (SC601).

Figure 15. ZDV union (SC600) and pre-column (SC001) assembly



7. Tighten the nut while ensuring that the pre-column and sleeve are firmly pushed against the center of the venting tee.
8. Install the analytical column on the other side of the venting tee using a 1/32 in. OD sleeve (F385). Check the label for the flow, and firmly push the sleeve and column against the center of the venting tee.

Figure 16. ZDV union (SC600), pre-column (SC001), venting tee (SC601), and analytical column (SC200) assembly



Connecting the Columns to a Stainless Steel Emitter

❖ To connect the columns to a stainless steel emitter

1. Remove both red nuts from a **second** ZDV union (SC600).
2. Insert the white blind plug into the union.
3. Insert the outlet end of the analytical column through a 1/32 in. OD sleeve (F385) and a red nut.
4. Screw the outlet end into the union, ensuring that the column and sleeve are firmly pushed against the white blind plug. Then, remove the white blind plug.
5. Insert the sleeved emitter (ES542) through a red nut.

6. Screw the sleeved emitter into the other end of the union, ensuring that the emitter and sleeve are firmly pressed against the column outlet and sleeve.

Figure 17. Stainless steel emitter with 1/32 inch sleeve (ES542)



Figure 18. ZDV union assembly between analytical column and stainless steel emitter




Mounting the Column and Emitter Setup onto the DirectJunction Adaptor

❖ **To mount the column and emitter setup onto the DirectJunction adaptor**

1. Mount the DirectJunction adaptor onto the Nanospray ion source.
2. Connect the red HV cable to the banana plug on the DirectJunction adaptor.
3. Open the front HV clamp on the DirectJunction adaptor.
4. Mount the ZDV union connecting the emitter and analytical column in the black holder onto the DirectJunction adaptor, and close the clamp over the emitter.
5. Adjust the position of the multipurpose adaptor on the steel rod, and mount the venting tee in the adaptor.

Figure 19. DirectJunction adaptor with two-column setup and steel emitter mounted on a Nanospray Flex ion source




Connecting the LC System to the Column Setup

❖ To connect the LC system to the column setup

1. Connect the “Column out” line from the EASY-nLC system to the pre-column through the ZDV union, ensuring the line is firmly pushed against the pre-column and sleeve.
2. Connect the “Waste in” line from the EASY-nLC system to the third port on the venting tee.


Figure 20. DirectJunction adaptor (as in [Figure 19](#)) but with “Column out” and “Waste in” flow lines connected from the EASY-nLC system



Connecting the DirectJunction to HV

Make sure that the high voltage connector of the DirectJunction adaptor is inserted into the Lemo™ HV socket on the underside of the Nanospray ion source drawer.

Figure 21.



Alternative Configurations for Direct Junction


Depending on the specific experimental requirements on columns, emitters, and connections, the following configurations might work well in different situations.

Contents

- [Stainless Steel Emitter and One-Column Setup](#)
- [Glass Emitter with Liquid Junction and One-Column Setup](#)
- [Glass Emitter with Liquid Junction and Two-Column Setup](#)
- [Packed Glass Emitter with Liquid Junction in One-Column Setup](#)
- [Packed Glass Emitter with Liquid Junction in Two-Column Setup](#)
- [Setup with a Long Analytical Column](#)

Stainless Steel Emitter and One-Column Setup

Figure 22. Stainless steel emitter with one column




Note EASY-nLC software version 2.5 and later require use of the “Waste in” line connected to the vent (SC601) upstream of the analytical column in one-column configurations.

4 Alternative Configurations for DirectJunction Glass Emitter with Liquid Junction and One-Column Setup

For earlier software versions or with other LC systems, you can replace the venting tee with a ZDV union (SC600).

Glass Emitter with Liquid Junction and One-Column Setup


Figure 23. Glass emitter with a liquid junction and one column



Note EASY-nLC software version 2.5 and later require use of the “Waste in” line connected to the vent (ES257) upstream of the analytical column in one-column configurations. For earlier software versions or with other LC systems, you can replace the Liquid Junction Cross (ES257) with a Liquid Junction Tee (ES258). Use small ID emitter tips (<math><20\ \mu\text{m}</math>) to prevent out gassing in the tip and subsequent spray instability.

Glass Emitter with Liquid Junction and Two-Column Setup


Figure 24. Glass emitter with a liquid junction and two columns



Note Use small ID emitter tips (<20 μm) to prevent out gassing in the tip and subsequent spray instability.

Packed Glass Emitter with Liquid Junction in One-Column Setup

Figure 25. Packed glass emitter with a liquid junction and one column



Note The packed glass emitter contains the stationary phase and has dual roles: as the analytical column and the emitter. The packed glass emitter is inserted through the ZDV union, which only serves as positional support.


Note EASY-nLC software version 2.5 and later require use of the “Waste in” line connected to the vent (ES257) upstream of the analytical column in one-column configurations. In earlier software versions or with other LC systems, you can replace the Liquid Junction Cross (ES257) with a Liquid Junction Tee (ES258).

4 Alternative Configurations for DirectJunction

Packed Glass Emitter with Liquid Junction in Two-Column Setup

Packed Glass Emitter with Liquid Junction in Two-Column Setup


Figure 26. Packed glass emitter with a liquid junction and two columns



Note The packed glass emitter contains the stationary phase and has dual roles: as the analytical column and the emitter. The packed glass emitter is inserted through the ZDV union, which only serves as positional support.

Setup with a Long Analytical Column

Figure 27. Long analytical column



Note You can easily fit long columns to the DirectJunction adaptor by coiling them between the ZDV union and the vent (ES257). [Figure 27](#) shows a glass emitter and a liquid junction, but combinations with a steel emitter, a two-column configuration, or both of these are also possible.

Note EASY-nLC software version 2.5 and later require use of the “Waste in” line connected to the vent (ES257) upstream of the analytical column in one-column configurations. In earlier software versions or with other LC systems, you can replace the Liquid Junction Cross (ES257) with a Liquid Junction Tee (ES258). Use small ID emitter tips (<20 μm) to prevent out gassing in the tip and subsequent spray instability.


Nano Flex Ion Source Configuration

When analyzing with nanospray sources, you only need to change a few mass spectrometer settings in the Xcalibur data system. The MS software will automatically recognize the ion source type.

❖ To specify configuration options for the ion source

1. Open the Instrument Configuration window from the computer desktop as follows:
 - For Xcalibur 2.0.7 or earlier versions, choose **Start > All Programs > Xcalibur > Instrument Configuration**. Or, double-click the **Instrument Configuration** icon.
 - For Xcalibur 2.1.0 or later versions, choose **Start > All Programs > Thermo Foundation 1.0 > Instrument Configuration**.

The [Thermo Foundation] Instrument Configuration window appears with a list of the installed available devices (see [Figure 28](#)).

Figure 28. Thermo Foundation Instrument Configuration window

2. Double-click the icon that represents your MS system (for example, LXQ MS).


A copy of the icon appears in the Configured Devices pane.

3. Complete the configuration of the mass spectrometer as follows:

- a. In the Configured Devices pane, double-click the **LXQ MS** icon (or the icon for your system).

The LXQ Configuration dialog box appears (see [Figure 29](#)).


Figure 29. LXQ Configuration dialog box



- b. Select **Ion Source** in the left pane if it is not already selected.
 - c. Ensure that the default source is **Nanospray** (the instrument recognizes the source automatically).
If the source is not correct, follow the instructions in the Xcalibur manual.
 - d. Click **OK** to close the window.
4. From your desktop, double-click the **LXQ Tune** icon (or the correct system name of your MS).

The Thermo Tune Plus window opens where you can analyze and choose specific settings (see [Figure 30](#)).


Figure 30. Thermo Tune Plus window



- From the main menu, choose **Setup** > **NSI Source**.

The NSI source dialog box opens.

Figure 31. NSI Source dialog box



6. Set the ion spray voltage for nanospray analysis. Pay attention to these nanospray-specific settings:

- For online analysis:

I Spray Voltage (kV): 1.40 to 2.40
(Start with **1.80** kV.)

- For offline analysis:

I Spray Voltage (kV): 0.70 to 1.20
(Start with **0.90** kV.)

7. Click **OK**.

Note Before you begin nanospray analysis, check that the drawer of the Nanospray Flex ion source is pushed completely toward the frame; otherwise, a message that the source is open appears and the MS will not scan.

The Nanospray Flex ion source is now be ready for analysis.

Operation

Both for initial setup and regular maintenance, use these procedures to monitor the position of the emitter tip and the spray voltage/current for optimal performance.

Contents

- [Positioning the Emitter Tip](#)
- [Maintaining Stable Spray](#)
- [Using the Optional Gas Connection](#)

Positioning the Emitter Tip

The correct position of the emitter tip is an important parameter in nano ES analysis.


❖ To optimize the position of the emitter tip

1. Position the spray capillary tip almost on-axis with the orifice by using the four knobs on the XYZ manipulator (see [Figure 32](#)).
2. As you make adjustments, observe the magnified inlet/orifice on the monitor.

6 Operation

Positioning the Emitter Tip

Figure 32. Adjustment knobs on the XYZ manipulator



CAUTION Emitter tips are extremely sharp. Never touch the tips as this may cause injury.




CAUTION

- Although the ion source is shielded, there is easy access to the source head.
- Always turn off the high voltage before touching the source head. Do not leave the source unattended while the spray voltage is on.
- Always mount the high voltage insulator correctly before turning on the high voltage.
- The interface region of the mass spectrometer can be hot. Do not touch this area.

Figure 33 illustrates some important settings whose values can vary depending on the mass spectrometer in use. Also, the optimal distance between the emitter tip and orifice depends on the flow rate, as does the optimal potential on the emitter.

Figure 33. Schematic setup of the ion source and the mass spectrometer interface



Values for stainless steel or fused silica emitters:
A = 3 to 5 mm and the potential on the emitter, 1400 to 2400 V

Maintaining Stable Spray

For best results, use stainless steel emitters for nano electrospray experiments. These emitters ensure straightforward initiation of the spray and help to maintain spray stability over extended periods (up to 1000 hours).

Unstable spray is detrimental to most analyses. [Table 1](#) lists the most common reasons for unstable spray and ways to correct it.

Table 1. Correcting unstable electrospray

Cause	Possible correction
The HV contact is interrupted.	Check/clean the gold contact in the HV clamp where the emitter is held in the DirectJunction adaptor. Check the red cable and connectors. Check if the mass spectrometer is supplying the correct high voltage.
Air bubbles in the emitter (if glass) might cause the emitter to “spit.”	Check if the solvents are degassed properly.
The emitter is blocked because of particulates in the sample, other small particles from the flow lines/valves, and so on.	(Stainless steel emitter only) If sonicating the emitter does not remove the blockage, replace the emitter.
There is a leak at some point in the liquid path.	Check if the correct flow is being delivered.


6 Operation

Using the Optional Gas Connection

Using the Optional Gas Connection

The drawer of the Nanospray Flex ion source is equipped with a gas connector, a Swagelok™ fitting for tubing with 1/8" (OD).

Figure 34.



Swagelok 1/8" gas connector, option for drying- or curtain gas.

In research analyses, use this optional gas connection for drying gas when you want to work in a controlled environment around the spray area.

Accessories and Spare Parts

This appendix provides ordering information for the Nano ES Spray Kit for offline analysis, and spare parts and accessories for the Thermo Nanospray Flex ion source.

For detailed descriptions and ordering information on the full range of spare parts and accessories for the Nanospray Flex ion source, visit www.proxeon.com.

Contents

- [Nano ES Spray Kit for Offline Analysis](#)
- [Accessories List](#)
- [Limited Warranty](#)
- [Additional Support](#)


Nano ES Spray Kit for Offline Analysis

Through offline analysis, you can extract the maximum amount of information from very limited amounts of sample. You can also average data extensively to improve the S/N-ratio and conditions optimized for MS/MS experiments. Using the Nano ES Spray Kit (see [Figure 35](#)), you can do the following:

- Work at low flow rates of 10–40 nL/min.
- Utilize nearly 100 percent of sample.
- Work effectively with sample volumes down to 300 nL.
- Avoid cross contamination by using disposable emitters.
- Spray from purely aqueous as well as purely organic solvents.

The Nano ES Spray Kit includes all necessary items for offline analysis. The kit is not included with the Nanospray Flex ion source but can be ordered separately.

Figure 35. ES259 offline Nano ES Spray Kit



Accessories List

This section lists the part numbers for the following Nanospray Flex ion source accessories:

DirectJunction	ES256
LiquidJunction with Cross 1/32 in.....	ES257
LiquidJunction with Tee 1/32 in.	ES258
Off-line source head kit.	ES259
One Dino-Lite video camera	ES216
One LCD Monitor 8 in.	ES217

Limited Warranty

- Thermo Fisher Scientific limits the warranty on repair of supplied Nano ES equipment and replacement, provided that instrument operation has been in accordance with the instructions outlined in the instruction manual.
- Abuse or misuse of the Nano ES equipment, or unauthorized repairs void this warranty.
- Limited warranty work is performed only at the factory, and the user bears the cost of shipment to and from the factory.
- The limited warranty is as stated above and no implied or inferred liability for direct or consequential damage is intended.

Additional Support

For further assistance, go to www.proxeon.com and click **Support and Service** in the left margin of the window.

Index

A

accessories 38
 adaptor
 See DirectJunction
 analysis, offline 37

B

BNC adaptor 10

C

camera
 adjusting 11
 Dino-Lite 9
 column and emitter setup 18
 column setup, long analytical 25
 compliance
 FCC iii
 regulatory iii
 WEEE v
 components, ion source 3
 configuring the ion source 27

D

declaration of conformity xii
 DirectJunction
 alternative configurations 21
 column and emitter setup 18
 description 13
 HV connection 19
 mounting 4
 Disclaimer xi

E

electromagnetic compatibility iii
 EMC compliance iii

emitter
 glass
 liquid junction, one-column setup 22
 liquid junction, two-column setup 22
 packed glass
 liquid junction, one-column setup 23
 liquid junction, two-column setup 24
 stainless steel and one-column setup 21
 emitter tip (caution) 34

F

FCC compliance iii

G

gas connector 36

H

HV cables and connectors 9

I

installation, ion source on MS 5
 instrument configuration window (Xcalibur) 27
 ion source
 removing existing 6
 ion source - MS schematic 35
 ion source kit, contents 2

M

mass spectrometers, ion source compatible ix
 monitor (LCD)
 adjusting 11
 connecting to cameras 9

O

offline analysis 37
 online/offline analysis 31

Index: P

P

positioning, emitter tip [33](#)
power supply and connectors [10](#)
preinstallation, ion source [3](#)

R

regulatory compliance [iii](#)

S

safety standards [iii](#)
spare parts [37](#)
spray capillary [33](#)
spray stability, maintaining [35](#)

U

unpacking the ion source
See preinstallation, ion source

W

warranty information [38](#)
WEEE compliance [v](#)

X

Xcalibur, settings for ion source [27](#)
XYZ manipulator [2](#), [33](#)
XYZ manipulator (caution) [9](#)