



• OwnIndex™ to reliably predict your VO₂max



• OwnCal™ the most accurate calorie consumption calculation




• OwnCode™ locks out possible interference from other heart rate monitors



• ExerciseProfiles™ to guide you through multi-phase exercise sessions



• Infra red data communication with your computer and Uplink™ feature



• Relaxation rate, beat to beat measurement

179291 GBR A

S810™
HEART RATE MONITOR



Quick Guide

Signal/Light

Stop

- exit the displayed mode
- return to the Time of day display



▲ Up

- move onward
- increase the selected value

▼ Down

- move backward
- reduce the selected value

OK

- start measuring the heart rate
- lock in your selection

How to start measuring your heart rate (BasicUse)

1. Place the transmitter around your chest and wet its electrodes. Wear the wrist receiver on your wrist.
2. Press and hold the OK button to start the measurement. The stopwatch and the exercise recording start. Your heart rate will appear in max. 15 seconds.

How to stop measuring your heart rate

1. Press the stop button. The stopwatch and other calculations stop. Heart rate measurement continues but exercise data is not recorded into a file anymore.
2. Press the stop button again. The receiver returns to the Time of day display.



The Polar S810 and Polar Precision Performance™ software offer you various possibilities on customising exercise sessions to fit your personal needs.

For example, you can build up exercise profiles to guide you through various training sessions. Choose the functions that suit you best: OwnCal calorie calculation, timers or recovery calculation, etc. Monitor your fitness level, and any long-term changes in it, by performing the quick and easy Polar Fitness Test™.

After training, you can review your exercise file and transfer it to your computer for analysis with the software. Or, you can transfer your heart rate to the computer in real time through an optional Polar IR Interface™. Even more thorough information about your body's recovery state you get by performing an overtraining test included in the software.

Please read this manual carefully to familiarise yourself with the Polar S810. Here is where you can find the information you need to use and maintain your equipment. The function road map on the front cover is a quick guide to the versatile features of the Polar S810; please tear it off.

Please refer to the glossary for special terms, display symbols and text explanations. The alphabetical index will help you quickly find answers to any questions that you might encounter while using your Polar S810.

Thank you for choosing Polar!



CONTENTS

A. Getting going	5	Setting phases on/off	25	Setting alarm	36	Exercise phase	47
Heart rate monitor parts	5	Setting interval calculation	26	Setting time of day	36	Interval period	47
Buttons and their functions	6	Setting interval type	26	Setting date	37	Recovery period	48
Easy start	8	Setting number of intervals		Turning reminders on/off	37	Interrupting an exercise profile ...	49
How to prepare the wrist receiver	8	and recoveries	28	Setting tips	38	Exercising with BasicSet	49
How to put the transmitter on	9	Turning heart rate limits on/off	28	Swap time 1 or time 2	38	Starting the timers	49
How to start measuring		Setting recovery calculation	28	Shortcut to watch settings	38	Swapping the heart rate limits	49
your heart rate	9	Setting heart rate limits	28	Swap measuring units	38	Starting the recovery calculation	50
How to stop measuring		Naming profiles	30	Swap date, user ID,		Stopping the exercise	50
your heart rate	10	Setting recording rate	31	user name and logo	38	Tips during exercise	50
After exercising	10	Turning functions on/off	32	C. Exercise	39	Swapping the exercise type	50
Modes and their functions	11	OwnCal calorie counter on/off	32	Heart rate measure mode	39	Swapping the recording rate	
B. How to enter settings	14	Fitness test on/off	32	Exercise recording mode	41	of the exercise	50
User information settings	16	Predicted maximum		Turning heart rate		Resetting the interval or	
Exercise settings	20	heart rate on/off	33	zone alarm on/off	41	recovery calculation of	
Selecting exercise type	22	Online recording to		illuminating the display	42	an exercise profile	50
Setting BasicSet exercise	22	computer on/off	33	Storing lap and split time	42	Interrupting an interval or	
Setting timers	23	Monitor settings	34	Pausing the exercise	43	recovery calculation of an	
Setting heart rate limits	23	Turning activity/button		Swapping the displayed		exercise profile	51
Turning heart rate limits on/off ...	23	sound on/off	34	exercise information	43	Interrupting a phase of	
Setting recovery calculation	24	Selecting measuring units	34	Resetting the stopwatch to zero	45	an exercise profile	51
Setting profile exercise	25	Turning Help on/off	35	Starting an exercise profile	46	Heart touch function	51
		Watch settings	36				



A. GETTING GOING


D. How to recall training information	52	F. Polar Fitness Test	65
Exercise time	55	OwnIndex	65
Recovery information of BasicSet exercise	55	Predicted maximum heart rate	67
Heart rate limits	55	Fitness test settings	68
Exercise time within, above and below the heart rate zones	55	Carrying out the test	68
Calorie expenditure	56	Update OwnIndex and HR _{max} -p	70
Cumulative calorie count	56	Fitness classes	71
Cumulative exercise time count	56	G. Care and maintenance	73
Profile exercise information	57	H. Precautions	75
Lap information	59	I. Frequently asked questions	79
Recorded heart rate samples	60	J. Technical specifications	83
Records file	60	K. Limited Polar Warranty	86
Resetting the cumulative counts to zero	60	L. Disclaimer	87
Deleting a file	61	M. Polar glossary	88
Deleting all files	61	Display symbols	90
E. Connection to computer	62	Display texts	91
Downloading training data to computer	62	Index	95
Uploading settings from computer	63		

Heart rate monitor parts and their functions

The elastic strap holds the transmitter comfortably, yet firmly, in the correct position.

The Polar T61™ transmitter detects your heart rate and transmits it to the wrist receiver.

The Polar S810 wrist receiver displays your heart rate and other data during exercise. Enter your personal settings into the receiver and analyse the exercise files after you're done.





Buttons and their functions

- Signal/ Light** Turns the beep on or off.
Turns the backlight on.
- Reset** Resets the heart rate monitor.
- Stop Return** - Stops measuring the heart rate.
Exits the displayed mode and returns to the previous mode level.
Returns to the Time of day display from any mode by pressing and holding the stop button.
- OK Accept** - Starts measuring the heart rate (start).
Enters the displayed mode shown on the lower level (start). Locks in your selection (ok). Records lap information (lap). The Help text (start, ok or lap) that appears on the display above the button indicates its use.



SIGNAL/
LIGHT

RESET



STOP

OK



Up

Moves to the following mode. Increases the selected value.



Down

Returns to the previous mode. Decreases the selected value.



Useful tips

- Briefly pressing a button (approximately 1 second) will do different things than pressing and holding the button for a longer period of time (approximately 2 to 5 seconds). If you press the button for a longer period of time, you can use shortcuts.
- The buttons are slightly stiffer than those of an ordinary watch to prevent them from being pressed accidentally.
- You can always return to the Time of day display from the Options setting or File recalling loops by pressing and holding the stop button.
- See the glossary for display symbols and text interpretation.

A 6

A 7





Easy start

How to prepare the wrist receiver

1. Since the display is blank when the receiver leaves the factory, you have to activate the wrist receiver by pressing the OK button twice, after which the Time of day display appears. This is a one-off procedure; once you have turned on your receiver you can not turn it off anymore.
2. You may start measuring your heart rate right away using the default settings. However, the following settings are ones that it would be useful to change as soon as possible: Time of day (see section Watch settings) and personal user settings (see section User information settings).
3. Polar S810 has various functions that you can use according to your personal training needs. For more information, see section Modes and their functions. You can either manually enter the settings by pressing the receiver buttons or you can define your settings using Polar Precision Performance software, after which the data should be transferred from your computer to the wrist receiver.
4. Wear the wrist receiver as you would wear an ordinary watch. If you go biking, we recommend attaching the wrist receiver to a Polar Bike Mount™.

A 8

How to put the transmitter on

1. Attach the transmitter to the elastic strap.
2. Adjust the strap length to fit snugly and comfortably. Secure the strap around your chest, just below the chest muscles, and buckle it.
3. Lift the transmitter off your chest and **moisten the grooved electrode areas** on the back.
4. Check that the wet electrode areas are firmly against your skin and the Polar logo is in a central, upright position.

How to start measuring your heart rate

1. Begin with the display that shows the time of day.
2. Keep the wrist receiver within 1 metre/ 3 feet of your transmitter.
3. Check that you are not near other people that have heart rate monitors, high voltage power lines, televisions, mobile phones or other sources of electromagnetic disturbance.

A 9





4. Press the OK button to start measuring your heart rate. A heart symbol will start flashing and your heart rate (beats per minute) will appear in a maximum of 15 seconds.
5. Press OK again. The stopwatch starts running and you can start exercising. The exercise data will only be stored in a file if you have turned the stopwatch on.

How to stop measuring your heart rate

1. Press the stop button. The stopwatch and other calculations stop. The heart rate measuring continues, although exercise data is no longer recorded.
2. Press the stop button again. The heart rate measuring stops. The receiver shows the Time of day display again.

After exercising

1. Carefully wash the transmitter with a mild soap and water solution.
2. Rinse it with pure water.
3. Dry the transmitter carefully with a soft towel.
4. Store the transmitter in a clean and dry place. Dirt impairs the elasticity and functioning of the transmitter. Sweat and moisture can keep the electrodes wet and the transmitter activated, which shortens battery life.

A 10



Modes and their functions

Your Polar S810 has six main modes: Time of day, Exercise, File, Options, Fitness test and Connection. The S810 modes offer you the following features:

1. TIME OF DAY

- 12 h or 24 h time, alarm, date and weekday
- two time options
- backlight

In addition, the following functions can be uploaded from a computer and utilised:

- 7 reminders to remind you with alarm sounds
- your personal logo, user number or user name

2. EXERCISE

There are three optional exercise modes:

- simple BasicUse function without settings or exercise guiding functions
- freely programmable Exercise Profile function with six exercise phases
- programmable BasicSet function with three adjustable target heart rate limits, timers and a recovery period

The following are in use during all three different exercise modes:

- three adjustable display options
- current exercise time (stopwatch)
- time of day
- lap and split time as well as the lap number

A 11



- current heart rate, average heart rate or a reading of your maximum heart rate in percentages
- calorie count for an exercise session
- graphic and numeric display options of relaxation rate

3. FILE

- up to 99 exercise files and one resettable records file, which keeps track of cumulative calories and exercise time
- date and starting time of the exercise
- total exercise time
- average, maximum and minimum heart rate readings for the exercise session
- recovery information of BasicSet exercise
- heart rate limits used during exercise
- time spent in, above and below your target zones during exercise
- OwnCal calorie expenditure
- cumulative calorie expenditure
- cumulative exercise time
- profile information on the exercise phases
- lap information: best lap time; up to 99 lap and split times; current, average and maximum heart rate for each lap
- heart rate samples recorded in the preset recording rate

4. OPTIONS

- BasicUse with no settings, BasicSet with three heart rate limits, timers and a recovery period or seven programmable exercise profiles with

exercise phases. Each phase may include heart rate limits, intervals and recovery periods.

- heart rate for 1 minute, 15 seconds, 5 seconds or each heart beat is recorded.
- OwnCal calorie counter, fitness test, predicted maximum heart rate and online measurement*
- user-specific information: weight, height, date of birth, gender, activity level, maximum heart rate, maximal oxygen uptake and baseline for relaxation rate
- activity/button sound, measurement units and Help function
- alarm, time of day, date and reminders

5. FITNESS TEST

In this mode you can:

- perform the Polar Fitness Test
- recall your latest test results; OwnIndex and HR_{max} -p values
- update your user information to reflect your most recent OwnIndex and HR_{max} -p values

6. CONNECTION

In this mode you can:

- download recorded exercise information to your computer for detailed analysis with Polar Precision Performance software version 3.0 or later
- upload your receiver settings from your computer to receiver with Polar Precision Performance software version 3.0 or later

* An optional Polar IR Interface is required.

A 12

A 13



B. HOW TO ENTER SETTINGS



This section describes how to manually enter settings by pressing the receiver buttons. However, Polar S810 offers you another and faster way as well: prepare your settings with Polar Precision Performance software and upload them to your receiver from your computer. See chapter **Uploading settings from computer**.

B 14




1. In the Time of day display, scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Scroll up or down until the desired setting mode is displayed.
4. Press OK to start the settings.

- Move forward and backward by using the up or down buttons.
- Press OK to start adjusting the value.
- Scroll up or down to adjust the value.
- Press OK to lock in your selection.
- Move backward to the previous setting level by pressing the stop button.
- To return to the Time of day display, wherever you are in the Options settings: Press and hold the stop button.
- The digits run faster if you press and hold the up or down button as you adjust the value. Release the button when approaching the desired value.
- On the display each function is presented by an abbreviated text. A flashing block on the display indicates how far you are in the setting mode.

B 15



User information settings



! In the User Set loop you can move forward only by pressing OK.

Measuring units depend on the settings in Options/ Monitor Set/
Unit 1 or Unit 2.

	Units 1	Units 2
Weight	kg	lbs
Height	cm	ft/ inch
Date of Birth setting order	day - month - year	month - day - year

1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Scroll up or down until USER SET is displayed.
4. Press OK to start the user information settings.

Weight (kg or lbs)

5. Scroll up or down to adjust your weight. Press OK.

Height (cm or ft/ inch)

6. Scroll up or down to adjust your height. Press OK.

B 16

Date of birth (Birthday)

7. Scroll up or down to adjust your day of birth. Press OK.
8. Month appears. Scroll up or down to select your month of birth. Press OK.
9. Year appears. Scroll up or down to adjust the year of birth. Press OK.

Units 2: setting order is month - day - year.

Sex

10. Scroll up or down to select your gender. Press OK.

Activity level (Activity)

11. Scroll up or down to select your activity level for Polar Fitness Test. Press OK.

Assess your long-term physical activity level. Do not change your activity level description if your regular exercise habits have changed during the last six months.

- Low You do not participate regularly in programmed recreational sport or heavy physical activity. E.g. you walk only for pleasure or occasionally exercise sufficiently to cause heavy breathing or perspiration.
- Middle You participate regularly in recreational sports. E.g. you run 5-10 km per week or spend 1/2-2 hours per week in comparable physical activity or, your work requires modest physical activity.
- High You participate regularly, at least 3 times a week, in heavy physical exercise. E.g. you run 10-20 km per week or spend 2-3 hours per week in comparable physical activity.

B 17





Top You participate regularly in heavy physical exercise at least 5 times a week.
E.g. you exercise to improve performance for competitive purposes.

Maximum heart rate value (HR_{max})

- Your age-predicted maximum heart rate value is displayed as a default setting when you set this value for the first time. **If you know your exact clinically tested current maximum heart rate value**, set the value by scrolling up or down.
- Press OK.

If you do not know this value, update it in the Fitness Test mode after having carried out the Polar Fitness Test.

Maximal oxygen uptake value (VO_{2max})

- The wrist receiver suggests value 35 for females and value 45 for males as a default setting when you set this value for the first time. **If you know your exact clinically tested current maximal oxygen uptake value**, set the value by scrolling up or down.
- Press OK.

If you do not know this value, update it in the Fitness Test mode after having carried out the Polar Fitness Test.

B 18

Baseline

16. Scroll up or down to select your relaxation rate baseline. Press OK.
USER SET appears.

The degree of relaxation rate is highly individual, therefore compare the value to your own baseline. To find out your personal baseline relaxation rate, you need to measure your relaxation rate in a recovered state. Follow the instructions below or refer to Polar Precision Performance software to easily determine your relaxation rate baseline.

	Time	Relaxation rate Example	Your value
1	0:15	32	_____
2	0:30	33	_____
3	0:45	35	_____
4	1:00	34	_____
5	1:15	33	_____
6	1:30	35	_____
7	1:45	36	_____
8	2:00	33	_____
		271	_____
		/ 8	/ 8
Baseline		34	_____

- Before the measurement, relax for 2-3 minutes.
- Start measuring your heart rate by pressing and holding the OK button.
- Press the up button until relaxation rate value (Rlx) in milliseconds is displayed in the upper row.
Or, if INT is displayed in the upper row press first the down button once and then up button until relaxation rate value (Rlx) is displayed.
- Write down your relaxation rate each 15 seconds for 2 minutes.
- Sum up the eight relaxation values and divide them by eight to obtain the average value. This is your baseline relaxation rate.

B 19





For getting the most of your relaxation rate measure your baseline relaxation rate:

- in that situation which best represents the time you will be using the relaxation rate.
- always in the same way (e.g., laying down), place and time.

Update the baseline value once a month or whenever you want to refer to your value in different situations.

To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.

Exercise settings

When you start setting exercise preferences, you first select whether you want to pre-programm an exercise profile or a BasicSet or exercise with the BasicUse option.

BasicUse (E0)

BasicUse is a simple mode for exercising. You do not have to do any settings before exercising with this Exercise mode. There are no heart rate limits in use in this mode. While exercising, you can see the following optional values on the display:

- current heart rate or average heart rate of the exercise
- stopwatch
- time of day

B 20



- split and lap times
- relaxation rate

Additionally, you can follow calories of your exercise session and percentage of your maximum heart rate, if you have set all the user information and turned the OwnCal calculator on.

Exercise profile (E1-E7)

Exercise profile will guide you through the exercise with up to six pre-programmable phases. You can name the exercise profiles, e.g., by the exercise type. Profiles can be used for versatile exercises where different phases are successive, e.g., in interval training.

During exercise, you can see the same optional values as in the BasicUse mode. Heart rate limits and up to 30 interval and recovery periods can be set for each phase.

BasicSet (E8)

BasicSet is a pre-programmable exercise set, but simpler than an exercise profile. There are no phases available although three timers and a recovery period help you split up the exercise session. The three adjustable heart rate limits allow you to monitor your workouts, securing the correct intensity. During exercise, you can see the same optional values as in the BasicUse mode.


B 21





Selecting exercise type

1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Press OK to start the exercise settings. BasicUse E0, Profil E1-E7 or BasicSet E8 is displayed.
4. Scroll up or down until the desired exercise type is displayed. Press OK.




If you chose BasicUse (E0): to return to the Time of day display, press and hold the stop button and skip the rest of the exercise settings.

If you chose Profil (E1-E7): continue to the step 5 of the chapter "Setting profile exercise/ Setting phases on/off".

If you chose BasicSet (E8): continue setting the functions for the BasicSet.

Setting BasicSet exercise



B 22



Setting timers

If you start from the Time of day display, repeat steps 1-5.

5. Scroll up or down until desired Timer is displayed. Press OK to start setting the timer.
 6. Scroll up or down to turn the timer on or off. Press OK. OFF: skip steps 7-8.
 7. Scroll up or down to adjust minutes. Press OK.
 8. Scroll up or down to adjust seconds. Press OK.
- Repeat until you have set all timers you want to use.

! You should select the timers on sequentially. E.g. you cannot set timer 3 on if timer 1 is off.

Setting heart rate limits for the BasicSet

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until desired Limits is displayed.
 6. Press OK to start setting the heart rate limits.
 7. Lim High appears. Scroll up or down to adjust your upper limit. Press OK.
 8. Lim Low appears. Scroll up or down to adjust your lower limit. Press OK.
- Repeat until you have set all heart rate limits you want to use.

Turning heart rate limits for the BasicSet on/off

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until desired Limits is displayed.
6. Press and hold the signal/light button to turn limits on or off. Release.

B 23





Setting recovery calculation for BasicSet

- A.** Recovery timer Ends your recovery as you reach the preset time.
- B.** Recovery heart rate Ends your recovery as you reach the preset heart rate.



If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until Recovery is displayed.
6. Press OK to start setting the recovery on or off. On/OFF starts to flash.
7. Scroll up or down to set the recovery on or off. Press OK. OFF: Skip steps 8-10. On: TIMER or Hr starts to flash.
8. Scroll up or down to select the recovery timer or heart rate. Press OK.

A. RecoTime is displayed.

9. Scroll up or down to adjust minutes. Press OK.
 10. Scroll up or down to adjust seconds. Press OK.
- Recovery appears.

Or

B. RecoHr is displayed.


9. Scroll up or down to adjust the heart rate. Press OK.
- Recovery appears.

To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.

B 24



Setting profile exercise



1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Press OK to start the exercise settings. BasicUse E0, Profil E1-E7 or BasicSet E8 is displayed.
4. Scroll up or down until the desired Profil E1-E7 is displayed. Press OK.


See an example of an exercise profile from chapter "Starting an exercise profile".

B 25





Setting Phases on/off




If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until the desired phase number is displayed.
6. Press OK to start setting the phase. On/OFF starts to flash.
7. Scroll up or down to set the phase on or off. Press OK. On: Int. and the phase number (P1) are displayed. OFF: continue setting the other phases from step 5 or to return to the Time of day display, press and hold the stop button and skip the rest of the exercise settings.

! You can not set the phase on, if you have not set the interval or recovery calculation on.

Setting interval calculation (Int.)



If you start from the Time of day display, repeat steps 1-7.

8. Start with the display showing Int. (= interval).
9. Press OK to start setting the interval. On/OFF starts to flash.
10. Scroll up or down to set the interval on or off. Press OK. On: Int.Type is displayed. OFF: continue setting the recovery or heart rate limits for the phase or to return to the Time of day display, press and hold the stop button and skip the rest of the exercise settings.

B 26



Setting interval type (Int. Type)

Select the interval type that will be used to end the interval calculation from three optional types:

- | | |
|-------------------------------------|--|
| A. Time-based interval | Ends your interval when you have used the selected time. |
| B. Heart-rate-based interval | Ends your interval as you reach the selected heart rate.
! When selecting heart-rate-based interval, recovery calculation must be set on. |
| C. Manual interval | End the interval by pressing OK. |

If you start from the Time of day display, repeat steps 1-10.

11. Start with the display showing Int. Type.
12. Press OK to start setting the interval type. TIMER, Hr or MANUAL starts to flash.
13. Scroll up or down to select the interval type that will be used to end the interval. Press OK.

If you chose manual interval type (MANUAL) skip steps 14 - 15.

- | | |
|-------------------------------------|---|
| A. Time-based interval | 14. Scroll up or down to adjust minutes. Press OK.
15. Scroll up or down to adjust seconds. Press OK. TIMER appears. |
| Or | |
| B. Heart-rate-based interval | 14. Scroll up or down to adjust the heart rate that will end your interval. Press OK. Hr appears. |

B 27






Setting the number of intervals and recoveries (Repeat)

If you start from the Time of day display, repeat steps 1-10.

You can set up to 30 intervals and recoveries (if set on) for each phase.

11. Scroll up or down until Repeat is displayed.
12. Press OK to start setting the number.
13. Scroll up or down to adjust the number. If you do not know in advance the precise number of the intervals select continuous interval by pressing and holding the up or down button until Cont  is displayed. If you choose Cont, the intervals and recoveries will succeed until you manually stop the exercise or until 30 intervals and recoveries have been made. Press OK.

Continue setting the recovery or heart rate limits by pressing the stop button. Int is displayed. To return to the Time of day display, press and hold the stop button and skip the rest of the exercise settings.

Turning heart rate limits for profile on/off

If you start from the Time of day display, repeat steps 1-7.

8. Scroll up or down until desired Limits is displayed.
9. Press and hold the signal/light button to turn limits on or off. Release.

Setting heart rate limits for profile

If you start from the Time of day display, repeat steps 1-7.


8. Scroll up or down until desired Limits is displayed.
 9. Press OK to start setting the limits.
 10. Lim High appears. Scroll up or down to adjust your upper limit. Press OK.
 11. Lim Low appears. Scroll up or down to adjust your lower limit. Press OK.
- Heart rate limits are used in the interval and recovery calculations.

B 28



Setting recovery calculation for profile (Reco)

- | | |
|-------------------------------|--|
| A. Recovery timer | Ends your recovery as you reach the preset time. |
| B. Recovery heart rate | Ends your recovery as you reach the preset heart rate. |



If you start from the Time of day display, repeat steps 1-7.

8. Scroll up or down until Reco and the number of the phase (P1) are displayed.
9. Press OK to start setting recovery on or off. On/OFF starts to flash.
10. Scroll up or down to turn the recovery on or off. Press OK.
On: TIMER/ Hr is flashing.
OFF: continue setting the interval calculation or heart rate limits for the phase or to return to the Time of day display, press and hold the stop button and skip the rest of the exercise settings.
11. Scroll up or down to select recovery timer or heart rate. Press OK.

- A.** RecoTime is displayed.
12. Scroll up or down to adjust minutes. Press OK.
 13. Scroll up or down to adjust seconds. Press OK.
- Reco appears.

Or

B 29





B. RecoHr is displayed.

12. Scroll up or down to adjust the value. Press OK.

Reco appears.

! If you set both heart-rate-based interval and recovery, adjust your recovery heart rate to be smaller than the interval heart rate.

To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.


Naming profiles

You can name the profiles, with seven letters, numbers or punctual marks.

The available digits are: 0-9, space, A-Z, a-z, - % / () * + : ?.


1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Press OK. E0 - E8 is displayed.
4. Scroll up or down until desired profiles (E1 - E7) is displayed.
5. Press and hold signal/light button. Adjustable letter is flashing.
6. Scroll up or down to select the desired mark. Press OK.
7. Repeat previous step until you have chosen all 7 letters.

To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.



Setting recording rate

1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Scroll up or down until MEMORY SET is displayed.
4. Press OK to start the memory settings. Rec. Rate is displayed.
5. Scroll up or down to select the desired recording rate. As you change the recording rate you can see the memory left in the lower row. Press OK.



The recording time depends on the selected recording rate. The use of lap recording and profile function reduces the maximum exercise times. When the memory becomes full, you can continue exercising but the exercise is no longer recorded.

Recording rate	Maximum exercise time
RR (Each heart beat is recorded)	100 % is approx. 30 000 heart beats. The memory left is displayed in percentages
5 s	44 h 13 min
15 s	99 h 59 min and total 120 h is approx. 5 days
60 s	99 h 59 min in one file and total 500 h is approx. 10 days


To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.





Turning functions on/off

1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Scroll up or down until FUNCTION SET is displayed.
4. Press OK to start the function settings. OwnCal is displayed.



! To set the following functions on you have to set your user information. If kg/lbs are displayed as you start setting the function on, the wrist receiver points out that you have not set your all user information.

OwnCal calorie counter on/off

If you start from the Time of day display, repeat steps 1-4.

5. Press OK to start setting the OwnCal function. On/OFF starts to flash.
6. Scroll up or down to set the OwnCal on or off. Press OK.

For getting the most accurate OwnCal results, update the VO_{2max} setting and the HR_{max} setting with clinically measured values (if available). If you do not know your clinically measured values perform Polar Fitness Test and update the OwnIndex and HR_{max} -p settings in Fitness Test mode.

Fitness test on/off

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until Fit. Test is displayed.
6. Press OK to start setting the fitness test. On/OFF starts to flash.
7. Scroll up or down to set the fitness test on or off. Press OK.

B 32

B 33

Predicted maximum heart rate on/off (HR_{max-p})

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until HR_{max-p} is displayed.
6. Press OK to start setting the HR_{max-p} . On/OFF starts to flash.
7. Scroll up or down to set the HR_{max-p} on or off. Press OK.

HR_{max-p} can be turned on only if Polar Fitness Test is on.

Online recording to computer on/off

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until Online is displayed.
6. Press OK to start setting the online recording. On/OFF starts to flash.
7. Scroll up or down to set the online recording on or off. Press OK.

To use this function you need an optional Polar IR Interface.


To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.





Monitor settings

1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Scroll up or down until MONITOR SET is displayed.
4. Press OK to start the monitor settings. Sound is displayed.



Turning the activity/button sound on/off

If you start from the Time of day display, repeat steps 1-4.

5. Press OK to start setting the sound. On/OFF starts to flash.
6. Scroll up or down to set sound on or off. Press OK.

Activity/button sound consists of the button sound and in the Measure and Exercise modes the following activity sounds: start, stop and lap recording, phase start, interval phase end and recovery calculation end sounds. Does not consist of heart rate zone alarms, watch alarm or fitness test alarm.

Selecting measuring units

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until Units is displayed.
6. Press OK to start setting the units. 1 or 2 starts to flash.
7. Scroll up or down to select the desired unit. Press OK.

The units settings affect the units in the user information settings and the order of your date of birth settings.

Units 1: kg, cm, day-month-year **Units 2:** lbs, feet, month-day-year

Turning Help on/off

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until Help is displayed.
6. Press OK to start setting the Help function. On/OFF starts to flash.
7. Scroll up or down to set the Help on or off. Press OK.

If you turn Help on, the flashing arrows guide you to use the correct buttons in the Options and File modes. During the exercise, when you are changing the middle row information, you can see the name of the function for few seconds.

To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.





Watch settings

1. In the Time of day display scroll up or down until OPTIONS is displayed.
2. Press OK to enter the Options mode. EXERCISE SET is displayed.
3. Scroll up or down until WATCH SET is displayed.
4. Press OK to start the watch settings. ALARM is displayed.



Setting alarm

If you start from the Time of day display, repeat steps 1-4.

5. Press OK to start alarm settings.
 -) On/OFF starts to flash.
6. Scroll up or down to set alarm on or off. Press OK.
 - 24 h: skip step 7.
7. 12 h: AM/ PM starts to flash. Scroll up or down to select AM or PM. Press OK.
8. The hours start to flash. Scroll up or down to adjust the hours. Press OK.
9. The minutes start to flash. Scroll up or down to adjust the minutes. Press OK.

OFF: skip steps 6- 9.

When the alarm comes on, it will sound for one minute or you can turn it off by pressing any of the five buttons.

Setting time of day

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until TIME1/ TIME 2 is displayed.
6. Press OK to start setting the Time of day.

B 36

7. Time is flashing. Scroll up or down to select time 1 or time 2. Press OK.
8. 12 h/ 24 h starts to flash.
 - Scroll up or down to select 12 h or 24 h time mode. Press OK.
 - 24 h: skip step 9.
9. AM/PM starts to flash. Scroll up or down to select AM or PM. Press OK.
10. The hours start to flash. Scroll up or down to adjust the hours. Press OK.
11. The minutes start to flash. Scroll up or down to adjust the minutes. Press OK.

Setting date

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until DATE is displayed.
6. Press OK to start adjusting the date.
7. Day appears. Scroll up or down to select the desired value. Press OK.
8. Month appears. Scroll up or down to select the desired value. Press OK.
9. Year digits start to flash. Scroll up or down to select the year. Press OK.

12 h time mode: month - day - year.

Turning reminders on/off

If you start from the Time of day display, repeat steps 1-4.

5. Scroll up or down until REMIND is displayed.
6. Press OK to start setting the reminders.
7. Scroll up or down until desired reminder is displayed. The number of the reminder is displayed in the corner of the display.
8. Press OK. On/OFF starts to flash.
9. Scroll up or down to set the reminder on or off. Press OK.

B 37





You can upload 7 reminders, each including time and name, to your wrist receiver from Polar Precision Performance software.

The reminder sounds an alarm only when time of day is displayed. When the reminder alarm comes on, it will sound for 15 seconds or you can turn it off by pressing any of the buttons except the signal/light button. Erase the reminder text by pressing any of the buttons except the signal/light button.

To continue your settings, press the stop button or to return to the Time of day display, press and hold the stop button.

Setting tips

Swap time 1 or time 2

In the Time of day display press and hold the down button. Time 1 or Time 2 is displayed in the uppermost row for a couple of seconds. Displayed time will now be used in all watch and alarm functions. If you chose time 2, "2" is displayed in the Time of day display.

Shortcut to watch setting loop

Press and hold the signal/light button for a couple of seconds to enter to watch setting loop. If you want to skip some parts and go directly to set a particular watch function, scroll up or down button until the desired watch function is displayed. Set the function as in the "Watch settings" section.

Swap measuring units

As you set your user information on weight and height press and hold the signal/light button until the units are swapped. Release the button.

Swap date, user ID, user name and logo

You are able to set a user ID number, user name and logo by Polar Precision Performance software. These settings are shown on the upper row of the Time of day display. You can swap between the user ID, user name, logo and date by pressing and holding the up button.

B 38

C. EXERCISE

Time of day



Heart rate measure mode



Exercise recording mode

You can use two modes in your exercise: Measure or Exercise mode. In the Measure mode you can see your heart rate but your exercise is not recorded. In the Exercise mode your exercise is recorded and stopwatch and other calculations are started.

Heart rate measure mode

1. Wear the transmitter and the wrist receiver as described in "Easy start" section.
2. To make sure that the code search will be successful:
 - Keep the wrist receiver within 1 metre/ 3 feet of your transmitter.
 - Check that you are not near to other people with heart rate monitors, high voltage power lines, televisions, mobile phones, cars, motor driven exercise equipment or other electromagnetic disturbance.



C 39






- Start from the Time of day display by pressing the OK to start measuring your heart rate. A heart symbol will start flashing and your heart rate in beats per minute will appear in max. 15 seconds. The frame around the heart symbol indicates that the heart rate reception is coded. The receiver returns automatically to the Time of day display within 5 minutes if there is no heart rate reception.

The following displays start to alternate as you enter the Measure mode:



*To use this function you need an optional Polar IR Interface.



Exercise recording mode

If you are continuing from the Measure mode, press OK.
Or, if you are starting from the Time of day display, press and hold OK.

! If BasicUse exercise starts, even if you have set an exercise profile on please check that you have set at least one of the exercise profile functions on.


Exercise information will be stored in an exercise file only when the stopwatch is running. The ongoing recording is indicated with a graphic bar that runs continuously on the display.

If there is 00 / - - reading on the display, there is no heart rate reception. Bring the wrist receiver up to your chest near the transmitter's Polar logo. The receiver starts looking for the heart rate signal again.

💡 In the beginning of your exercise, when it has lasted less than 60 minutes, the exercise time is displayed in minutes and seconds. When you have exercised more than 60 minutes, the exercise time is displayed in hours and minutes.

Turning the heart rate zone alarm on/off

Press and hold the signal/light button. When «)» appears on the display the zone alarm is on. You will hear a beeping sound with each heart beat if your heart beat is outside the heart rate limit. If you do not want to use the zone alarm the flashing heart rate reading indicates if your heart rate is outside the limits. This function can be used if you have set heart rate limits on in profile and BasicSet exercises.





In addition to the alarm you can see:

- ▲ if you are exercising above heart rate zone.
- ▼ if you are exercising below heart rate zone.

You have three optional Exercise modes depending on the selected exercise:

- BasicUse (E0)
- BasicSet (E8)
- Profiles (E1-E7)

The following functions can be used during all three optional Exercise modes.

Illuminating the display

Press the signal/light button.

After you use backlight once in the Measure mode, Polar S810 automatically illuminates the display when you start or stop the exercise session, when you store lap information and at the end of an interval and a recovery calculation.

Storing lap and split time

Press OK to store lap and split time.

First

Lap time
Average heart rate of the lap
Number of the lap



Then

Split time
18:00.0
158 01

C 42



Polar S810 automatically stores a lap time when you stop recording an exercise. Polar S810 can store 99 laps if you have not used the profile function. After storing 99 laps FULL appears for a few seconds every time you take a lap. You can continue taking laps but they will not be stored. The wrist receiver continues recording all other information of the exercise, except lap and profile information.

Pausing the exercise

Press the stop button.

The exercise recording, stopwatch and other calculations are paused. To continue them, press OK. Or, to return to the Time of day display, press and hold the stop button.



The receiver automatically returns to the Time of day display within 5 minutes if you forget to quit the heart rate measuring after having stopped the stopwatch and taken the transmitter off your chest.

Swapping the displayed exercise information

There are three display options, which you can swap between by pressing the down button. Once you have configured the display options, Polar S810 saves the settings. If the Help function is on, the name of the selected function appears for a few seconds.

Stopwatch — 5w21:38.0
Time of day — 8:50 10
Current heart rate — 155 02

Time of day — Time 8:50
Stopwatch — 21:38.0
Lap time — 20:10.0


Stopwatch — 5w21:38.0
Lap time — 20:10.0
Current heart rate — 155 02

C 43





! When you have set a profile exercise on, you will have one extra display. You can change the middle and lower rows for this display



See the alternative phase displays from "Starting an exercise profile/ Exercise phases" section. Configure the display option and choose the displayed information for each row.

You can configure the middle and upper rows during heart rate measurement.

Select the middle row information

1. Press the down button to select the desired display alternative. When you swap the middle row information, the upper and lower rows change also. Continue selecting the upper and lower rows.

Select the upper row information

2. Press the up button to select from the following: Time of day (Time), stopwatch (Sw), lap time (Lp), relaxation rate (RLX) or calorie expenditure (Cal), if OwnCal is set on.

Additionally you will have a graphical relaxation rate display in the upper row. Relaxation rate helps you recognise stressful situations. The degree of relaxation rate is highly individual, therefore compare the value to your own baseline.



The bar indicates the baseline of your relaxation rate

 Bar of dots indicates the current relaxation rate in relation to your baseline value. If you are in a more relaxed state than the baseline value, the bar of dots is on the right side of the baseline bar, and vice versa.

You cannot select the same displayed function to the upper and middle row i.e. the option on the middle row is no longer available for the upper row.


Select the lower heart rate information row

3. Press and hold the down button to select from the following: Current heart rate, average heart rate (avg) or the reading in percentages of your maximum heart rate (% max) if you have set your user information.

You can only change the displayed function of the lower row if the stopwatch is on.

Resetting the stopwatch to zero


1. Pause your exercise by pressing the stop button.
2. Press and hold the down button until the stopwatch is reset. "Reset" stops flashing as the stopwatch is reset. Release. If you are resetting a profile exercise repeat step 2.
3. Press OK to restart the stopwatch and recording or to return to the Time of day display, press and hold the stop button. As a result of resetting your stopwatch the recorded file is deleted.






Starting an exercise profile (E1-E7)

1. Start from the Time of day display. Press OK to enter the Measure mode.
2. If you want to swap the profile (E1-E7), press and hold the up button until the exercise is swapped. Release. Repeat until desired profile is displayed.
3. To start exercising and recording the exercise press OK.

Profile function automatically guides you through the exercise. In order to hear the activity sounds when the phase starts and stops make sure the activity/button sound is on. When the stopwatch starts running  indicates that profile function is on. As exercise ends the Polar S810 automatically starts a BasicUse exercise, which is recorded in the same file as the exercise profile.

Example of a exercise profile structure:



C 46



Exercise phase (P1-P6)

Each exercise profile may include up to six preprogrammable exercise phases. Each phase may include up to 30 interval and recovery periods, which alternate continuously through the phase.

Interval period

1. First
Preset heart rate limits (if set on).
2. Then

One of the following is displayed depending on your settings:
A. manual, **B.** time-based or **C.** heart-rate-based interval.

A. Manual interval

Number of set intervals. (C means you have set continuous intervals)

Current interval number.

Stopwatch

Phase number

If you selected manual interval, your interval will end as you press OK.

Or

B. Time-based interval

Countdown timer.

If you selected time-based interval, your interval will end as you reach the preset time.

Or

C. Heart-rate-based interval

The difference between the preset interval heart rate and current heart rate.

If you selected heart-rate-based interval, your interval will end as you reach the preset heart rate reading.



C 47





- 3. At the end of interval
- Interval duration.
- Average heart rate of the interval.

When your intervals end, the wrist receiver automatically continues following the preset function, recovery calculation or the following exercise phase.

Recovery period

One of the following is displayed depending on your settings:

A. time-based recovery or **B.** heart-rate-based recovery.

A. Time-based recovery



- Recovery countdown timer.
- Current heart rate.

Or

B. Heart-rate-based recovery



- The difference between the preset recovery heart rate and current heart rate.
- Current heart rate.

You will hear two beeps when the recovery calculation is done.



- At the end of recovery
- Recovery duration.
- Drop in your heart rate.

C 48



After the last phase has ended the receiver automatically continues exercise recording with BasicUse mode.

Interrupting an exercise profile

1. Press the stop button to pause the exercise.
2. Press and hold the stop button to return to the Time of day display.

Or

Press the stop button. BasicUse appears. You can continue measuring your heart rate in this mode (Measure mode) or to continue exercising with BasicUse and recording your exercise information press OK. Repeat steps 1 and 2 to return to the Time of day display.

Exercising with BasicSet (E8)

In this mode you can also use the functions in “Exercise recording mode” section.

Starting the timers

If timers are set on, they automatically start running once you start the stopwatch. Timers alternate in the following order: timer 1, 2, 3, 1, 2, 3, etc. At the end of timer 1 you will hear 1 beep, at the end of timer 2 you will hear 2 beeps and at the end of timer 3 you will hear 3 beeps. Timers stop running when you stop the stopwatch.

Swapping the heart rate limits

Press and hold the up button. Repeat until the desired limits are displayed.

C 49





Starting recovery calculation

To start the recovery calculation, press and hold the OK button. Depending on your setting RecoHr or RecoTime appears in the upper row of the display. See "Starting an exercise profile"/ "Recovery calculation" for the displayed information during the recovery calculation. After the recovery calculation is done the exercise is paused. You can continue the exercise by pressing the OK button after the recovery calculation. If you do another recovery calculation, the recovery information will be deleted in consequence.

Stopping the exercise

1. Press the stop button to pause the exercise.
2. Press and hold the stop button until time of day is displayed. Your lap and split time are automatically recorded, when you end the exercise.

Tips during the exercise

Swapping the exercise type (E0-E8)

In the Measure mode press and hold the up button. Release. Repeat until desired exercise is displayed.

Swapping the recording rate of the exercise (rr, 5 s, 15 s or 60 s)

In the Measure mode press and hold the down button. Release. Repeat until desired recording rate is displayed.

Resetting the interval or recovery calculation of an exercise profile

Press the stop button to pause the exercise.

C 50



Reset the calculation by pressing and holding the down button. Reset Int? (=interval) or Reset Re? (=recovery) is displayed until the calculation is reset. To restart the same calculation press the OK button or to reset the exercise, press and hold the down button.

Interrupting an interval or recovery calculation of an exercise profile

Press and hold the OK button during your profile exercise to interrupt the calculation. The following section automatically starts. E.g. if you interrupt an interval your recovery calculation (if set on) automatically starts.

Interrupting a phase of an exercise profile

Press the stop button to pause the exercise. Press and hold the OK button. The following phase automatically starts. If you have done all the preset phases the BasicUse exercise starts.

Heart touch function

The heart touch function is activated as you bring the wrist receiver near the transmitter's Polar logo.

You have three alternatives for the heart touch function. You need to select it by Polar Precision Performance software:

- Checking the current heart rate limits
- Swapping the three or four display options
- Storing lap time

During the heart touch, the backlight is on for a couple of seconds, if you have already used the backlight during the exercise.

C 51






D. HOW TO RECALL TRAINING INFORMATION



The Polar S810 starts saving exercise information to a file as you start the stopwatch in the Exercise mode. You can recall the information in the File mode. In the File mode there are up to 99 exercise files and a record file, which keeps track of the cumulative values.

D 52



Your first file is F1, then F2 and so on, until you have recorded maximum 99 exercise files. The bigger the file number the more recent it is.




Some information alternates on the display automatically. You can make the alternation faster by pressing the OK button , e.g. , the target heart rate zone information.

Exercise file

1. In the Time of day display scroll up or down until FILE is displayed.
2. Press OK to enter the File mode. The main information of your file is displayed.

Starting date and time alternate on the display.




3. Scroll up or down until the desired file is displayed.
4. Press OK to start recalling an exercise file. Exe. Time is displayed.

D 53





- Move forward and backward within file loops with up or down buttons.
- Press OK to enter a file or to enter a deeper level in the file loop.
- Press the stop button to move backward to the previous file level.
- Press and hold the stop button to return the Time of day display, wherever you are in the File recalling cycle.



D 54



Exercise time (Exe. Time)

The exercise time is the amount of time that you have exercised with the stopwatch running. Average, maximum and minimum heart rates of your exercise alternate on the display. To continue recalling the file, scroll up or down button.



Recovery information of BasicSet exercise (Reco)

- Recovery time.
- Recovery duration.
- Drop in your heart rate.

Or



- Recovery heart rate.
- Time it took reach the preset heart rate reading.
- Drop in your heart rate.

Heart rate limits (Limits)

Settings of the heart rate limits alternate on the display. Limits 1 are used for heart rate zone 1, limits 2 for zone 2, etc.

Exercise time within, above and below the heart rate zones

Indicates the exercise time, which you have spent within (InZone), above and below each heart rate limits during the entire exercise.

D 55





Calorie expenditure (KCal)

You will see how many kilocalories you have expended during the exercise.

Cumulative calorie count (Tot. KCal)

You will have a count for cumulative calories expended during several exercise sessions starting from the previous resetting of the value in the records file. This function allows you to follow the calorie count as a gauge of effective exercise e.g. during one week.

Cumulative exercise time count (Tot. Time)

This function counts cumulative exercise time expended during several exercise sessions starting from the previous resetting of the value in the records file. It allows you to follow the exercise time count as a gauge of amount of exercise e.g. during one week.




Your cumulative total time is displayed in hours and minutes until 99 hours 59 minutes is reached. After this your cumulative exercise time is displayed in hours (100 h) until 9999 hours is reached.

D 56



Profile exercise information (EXE. SET)



Start with the display showing EXE. SET.

Profile exercise time.

Profile number.

1. Press OK to start recalling the profile exercise information.
2. Scroll up or down to see the information on exercise phases.



Phase number.

Phase duration.

Average heart rate of intervals and recoveries.

3. Press the stop button to return to the previous file level and skip steps 4 and 5, or continue recalling more detailed information of the phase by pressing OK.

D 57






4. Scroll up or down between intervals and recoveries.


Interval information

Interval number


Split time of the interval



Recovery information



Or



5. To continue recalling information on the other phases, press the stop button or to continue recalling, e.g., lap information press the stop button twice.



D 58



Lap information (LAPS)

Start with the display showing LAPS and the number of recorded laps.

1. Press OK to start recalling the lap information.
2. Scroll up or down to see the following information:




Best Lap

Shortest lap time.

The number of the best lap.

The best lap information is displayed if you have stored at least 3 laps. The best lap can not be the last lap.



Laps

Split time.

Lap time.

Lap number.

The heart rate at the end of the lap,

average and

maximum heart rates of the lap.

Press the stop button to exit the LAPS display.

D 59






Recorded heart rate samples

Check from the sample file detailed heart rate information the exercise. Start with the display showing SAMPLES.

1. Press OK to start recalling the recorded samples of your exercise information. The following is displayed:
The RR reading in milliseconds (ms) or the sample number (S1), if RR recording rate was set off, the time when the heart rate sample was recorded and heart rate.
2. Scroll up or down the samples.
Press and hold the stop button to return to the Time of day display.



Records file

This file is updated each time your exercise is recorded to file.

1. In the Time of day display scroll up or down until FILE is displayed.
2. Press OK to enter the File mode.
3. Scroll up or down until Records is displayed.
4. Press OK to enter the records file. Tot. KCal is displayed.
You will have a count for cumulative calories expended during several exercise sessions starting from the previous resetting.
5. Scroll up or down to see the cumulative exercise time (Tot. Time) starting from the previous resetting.
You can now return to the Time of day display by pressing and holding the stop button or continue resetting your cumulative counts.

Resetting your cumulative counts to zero

Once your have reset the cumulative count you cannot retrieve it.

Start with the display showing Tot. KCal or Tot. Time in the records file.

1. Press OK to start resetting the desired total count. Reset appears and the value starts to flash.

D 60

2. Press and hold the down button. Reset starts to flash. You can still cancel the resetting by releasing the down button before the value is reset. If you are sure that you want to reset the value keep holding the down button until value is reset. Release.
3. Press the stop button to exit the reset display. Press and hold the stop button to return to the Time of day display.

Deleting a file

1. In the Time of day display scroll up or down until FILE is displayed.
2. Press OK to enter the File mode.
3. Scroll up or down until desired exercise file is displayed. Press and hold the signal/light button. DELETE FILE and file number are displayed.
4. Press OK. The receiver confirms if you want to delete a file (Are you sure).
5. Press OK to delete the selected file.
6. Press stop to exit the file. Wait appears. The receiver reorganises your file, which takes few seconds.
7. Press stop to return to the Time of day display.

Deleting all files

1. In the Time of day display scroll up or down until FILE is displayed.
2. Press and hold the signal/light button. DELETE FILES is displayed.
3. Press OK. The receiver confirms if you want to delete all files (Are you sure).
4. Press OK to delete all exercise files.
5. Press stop to return to the Time of day display.

E 61



E. CONNECTION TO COMPUTER

There is two-way communication between Polar S810 and a computer. Polar Precision Performance software offers an easy way to analyse training data. A training log makes it possible to plan future training. Wrist receiver settings are easy to upload from the software instead of entering them manually.



After preparing your settings you can upload options lock function. Your options can no longer be changed until you deactivate the options lock from the "user features" of your software.

Time of Day ↔ File ↔ Options ↔ Fit. Test ↔ Connect

Downloading training data to a computer

1. Start the software.
2. Connect the Polar IR Interface to the serial port (RS-232) of your computer.
3. In the Time of day display scroll up or down until CONNECT is displayed.
4. Point the infrared window of the receiver towards the infrared window of the interface. The maximum infrared transfer distance is 15 cm/ 6" and the infrared transfer angle ± 15 degrees.

E 62



5. Activate the heart rate monitor connection function of the software. Your receiver should now communicate with the software. The receiver automatically returns to the Time of day display after data transfer.

Interrupt the data transfer to the computer by pressing the stop button during data transfer. CONNECT is displayed.

! To be able to transfer the recorded training information via IR communication and to analyse it with Polar Precision Performance software version 3.0 or later, you need the Polar IR Interface, unless there is an IR data port and Windows® 95 operating system on your computer.

Uploading settings from computer

Your Polar S810 offers you the option of preparing your settings with Polar Precision Performance software. You are also able to set following functions with the software:

- a user ID, user name or logo for easier receiver recognition in multi-user situations
- the heart touch function to show the current heart rate limits, store a lap time or swap the display options during your exercise
- names and alarm times for 7 reminders
- options lock function on or off

E 63



You may upload settings from a computer to a wrist receiver either via IR communication or via audible Polar UpLink. For IR communication please see the previous chapter and for the Polar UpLink communication proceed as follows:

1. Start the software.
2. If the receiver is in the Time of day display scroll up or down until CONNECT is displayed.
3. Place the receiver less than 10 cm/ 4" from the loudspeaker.
4. Activate the upload function of the software. Do not move the receiver during the data transfer.
5. The receiver automatically starts to upload the information as UpLink and COM appear. You will hear a chirping sound while information is transferred to the receiver. The receiver automatically returns to the Time of day display after data transfer.

If the data transfer was unsuccessful bring the receiver closer to the loudspeakers or turn the volume level of the loudspeakers up.

Interrupt the data transfer to the computer by pressing the stop button during data transfer. CONNECT is displayed.

! To be able to upload the settings from Polar Precision Performance software via Polar UpLink, you need a personal computer with a sound card (compatible with SoundBlaster™) and dynamic loudspeakers or earphones.



F. POLAR FITNESS TEST



Polar Fitness Test is an easy, safe and quick way to estimate individual maximal aerobic power and to get predicted maximum heart rate value. Polar Fitness Test is targeted to healthy adults.

OwnIndex

OwnIndex is a value which is comparable to maximal oxygen uptake (VO_{2max}), a commonly used descriptor of aerobic fitness. Cardiovascular (aerobic) fitness relates to how well your cardiovascular system works to transport and utilize oxygen in your body. The stronger and more efficient your heart is, the better the cardiovascular fitness is. VO_{2max} is a good indicator of performance capability in endurance sports.





If you want to improve your cardiovascular fitness, it takes a minimum of 6 weeks on an average to see a noticeable change in OwnIndex. Less fit individuals see progress even more rapidly and for more fit individuals more time is needed. The better the cardiovascular fitness of an individual, the smaller the improvements in OwnIndex.

Cardiovascular fitness is best improved by exercise types which employ large muscle groups. Such activities include cycling, running, swimming, rowing, skating, cross-country skiing and walking.

To monitor your fitness progress start with measuring your OwnIndex a couple of times during the first two weeks, to get a baseline value. Thereafter, repeat the test approximately once a month. OwnIndex is based on resting heart rate, heart rate variability at rest, age, gender, height, body weight and self-assessed physical activity.



Predicted maximum heart rate (HR_{max-p})

The HR_{max-p} definition is carried out simultaneously with the Polar Fitness Test.

The HR_{max-p} score predicts your individual maximum heart rate value more accurately than the age-based formula (220-age). The age-based method provides a rough estimation and is thus not very accurate, especially for people who have been fit for many years or for older people. The most accurate way of determining your individual maximum heart rate is to have it clinically measured (in maximal treadmill or bicycle stress test) by a cardiologist or an exercise physiologist.

Your maximum heart rate changes to some extent in relation to your fitness. Regular endurance exercise tends to decrease HR_{max} . There may also be some variation according to the sport you participate in. For example, running $HR_{max} >$ cycling $HR_{max} >$ swimming HR_{max} .

HR_{max-p} gives the possibility of defining training intensities as percentages of maximum heart rate and of following maximum heart rate changes due to training, without an exhaustive maximal stress test. HR_{max-p} is based on resting heart rate, heart rate variability at rest, age, gender, height, body weight and maximal oxygen uptake, VO_{2max} (measured or predicted). The most accurate HR_{max-p} is obtained by entering your clinically measured VO_{2max} into the wrist receiver.





Fitness test settings

To be able to carry out the Polar Fitness Test you need to do the following settings in the Options mode:

- Set your personal user information and long-term physical activity level.
- Set the Polar Fitness Test on.
- Set the HR_{max}-p on if you want to get a prediction of your maximum heart rate.

Carrying out the test

To get reliable test results, the following basic requirements apply:

- You should be relaxed and calm.
- The test can take place anywhere – home, office, health club – as long as the testing environment is peaceful. No disturbing noises (e.g. television, radio or telephone), no other people talking to you.
- Keep the testing place, time of day and environment the same every time the test is repeated.
- Avoid eating a heavy meal or smoking 2-3 hours prior to testing.
- Avoid heavy physical effort, alcoholic beverages or pharmacological stimulants on the test day and the day before.

F 68



Time of Day ↔ File ↔ Options ↔ Fit. Test ↔ Connect


1. In the Time of day display scroll up or down until FIT. TEST is displayed. If you cannot find Fitness Test mode check that you have turned the function on in the Options mode.
2. Press OK to enter the Fitness Test mode. Your latest OwnIndex and the testing date are displayed.
3. Scroll up or down to see your latest HR_{max}-p and the testing date.
4. **Lie down and relax** yourself for 1-3 minutes before the test.

Start the test

It takes 3-5 minutes to carry out the test.

5. Press OK to start the test. The wrist receiver starts searching for your heart rate. TEST On and your heart rate are displayed. The test begins.
6. Lie relaxed. Keep your hands beside your body and avoid body movements. Don't communicate with other people.
7. Your current OwnIndex and the testing date are displayed.
8. If HR_{max}-p is on: Scroll up or down to see the HR_{max}-p result and the testing date.

Turn the fitness test sound signal on/off

The sound signal will alarm you at the end of the test. If the sound signal is on  appears on the display.

- To turn the sound signal on or off either before or during the test, press and hold the signal/light button.

F 69





Interrupt the test

- You can stop the test any time during the test by pressing the stop button. Failed TEST is displayed for a few seconds. The latest OwnIndex and HR_{max-p} are not replaced.

If the wrist receiver does not succeed in receiving your heart rate at the beginning or during the test, the test will fail. Check that the transmitter electrodes are wet and the elastic strap is snug enough.

Update OwnIndex and HR_{max-p}

You can update your OwnIndex and HR_{max-p} values to the Options mode/ User information. When you update these values regularly, the reading in percentages of your maximum heart rate and calorie calculation give more accurate individual information.

- In the Time of day display scroll up or down until FIT. TEST is displayed.
- Press OK to enter the Fitness Test mode.
- Scroll up or down until Update USER SET is displayed. This display appears if you have not updated your latest test results.
- Press OK to update the OwnIndex and HR_{max-p} values.
- To return to the Time of day display, press and hold the stop button.

! Your maximum heart rate value is updated if the HR_{max-p} is on.

F 70



Fitness classes

The Polar Fitness Test result, your OwnIndex, is most meaningful when comparing your individual values and changes in them. OwnIndex can also be interpreted according to your gender and age. Locate your OwnIndex on the table to find out your current cardiovascular fitness classification, compared with those who are of the same age and gender.

	Age (years)	1 (very poor)	2 (poor)	3 (fair)	4 (average)	5 (good)	6 (very good)	7 (excellent)
MEN	20-24	< 32	32-37	38-43	44-50	51-56	57-62	>62
	25-29	< 31	31-35	36-42	43-48	49-53	54-59	>59
	30-34	< 29	29-34	35-40	41-45	46-51	52-56	>56
	35-39	< 28	28-32	33-38	39-43	44-48	49-54	>54
	40-44	< 26	26-31	32-35	36-41	42-46	47-51	>51
	45-49	< 25	25-29	30-34	35-39	40-43	44-48	>48
	50-54	< 24	24-27	28-32	33-36	37-41	42-46	>46
55-59	< 22	22-26	27-30	31-34	35-39	40-43	>43	
60-65	< 21	21-24	25-28	29-32	33-36	37-40	>40	
WOMEN	20-24	< 27	27-31	32-36	37-41	42-46	47-51	>51
	25-29	< 26	26-30	31-35	36-40	41-44	45-49	>49
	30-34	< 25	25-29	30-33	34-37	38-42	43-46	>46
	35-39	< 24	24-27	28-31	32-35	36-40	41-44	>44
	40-44	< 22	22-25	26-29	30-33	34-37	38-41	>41
	45-49	< 21	21-23	24-27	28-31	32-35	36-38	>38
	50-54	< 19	19-22	23-25	26-29	30-32	33-36	>36
55-59	< 18	18-20	21-23	24-27	28-30	31-33	>33	
60-65	< 16	16-18	19-21	22-24	25-27	28-30	>30	


The classification is based on literature review (Shvartz & Reibold 1990) of 62 studies where VO_{2max} was measured directly in healthy adult subjects in the USA, Canada and 7 European countries. Reference: Shvartz E, Reibold RC: Aerobic fitness norms for males and females aged 6 to 75 years: a review. Aviat Space Environ Med; 61:3-11, 1990.

F 71





For various fitness classes we recommend the following:



Top athletes typically score OwnIndex values above 70 (men) and 60 (women). As high values as 95 can be reached by Olympic level athletes. OwnIndex is highest in sports that involve large muscle groups such as cross-country skiing and cycling.

F 72

G. CARE AND MAINTENANCE

Your Polar heart rate monitor is a high-tech instrument of superior design and workmanship and should be treated with care. The suggestions below will help you fulfill the warranty obligations and enjoy this product for many years to come.

Taking care of your Polar heart rate monitor

- Wash the transmitter regularly after use with a mild soap and water solution. Dry it carefully after washing with soft towel.
- Never store the transmitter wet. Sweat and moisture can keep its electrodes wet and the transmitter activated, which shortens the battery life.
- Store your Polar heart rate monitor in a cool and dry place. Do not store it in any kind of non-breathing material, such as a plastic bag or a sports bag if it is wet.
- Do not bend or stretch the transmitter. This may damage the electrodes.
- Do not dry the transmitter in any other way than with a towel. Hard-handed handling may damage the electrodes.
- Keep your Polar heart rate monitor out of extreme cold and heat. The operating temperature is -10 °C to 50 °C/ 14 °F to 122 °F.
- Do not expose the Polar heart rate monitor to direct sunlight for extended periods, such as by leaving it in a car.

G 73





Batteries

Polar T61 transmitter

The estimated average battery life of the Polar T61 transmitter is 2500 hours of use. Contact your distributor/ dealer for a transmitter replacement. Polar recycles used transmitters. See your Customer Care Charter for detailed instructions.

Polar S810 wrist receiver

The estimated average battery life of the Polar S810 wrist receiver is 2 years in normal use (2h/ day, 7 days a week). Note also that excessive use of the backlight and the alarm signals will use the battery more rapidly. Do not open the Polar wrist receiver yourself. To ensure the water resistance properties and the use of qualified components, the wrist receiver battery should be replaced only by an authorised Polar Service Centre. At the same time a full periodic check of the Polar heart rate monitor will be done.

Service

Should your Polar heart rate monitor need service, see Customer Care Charter to contact an authorised Polar Service Centre.

! Water resistance cannot be guaranteed after unauthorised service.



H. PRECAUTIONS

Using the Polar S810 in a water environment

Your Polar S810 is water resistant to 50 metres/ 160 feet.

Heart rate measurement in a water environment is technically demanding for the following reasons:

- Pool water with a high chlorine content and seawater may be very conductive and the electrodes of Polar transmitter may get short circuited and ECG signals cannot be detected by the transmitter unit.
- Jumping into the water or strenuous muscle movement during competitive swimming may cause water resistance that shifts the transmitter on the body to a location where it is not possible to pick up ECG signal.
- The ECG signal strength varies depending on the individual's tissue composition and the percentage of people who have problems in heart rate measuring is considerably higher in a water environment than in other use.





Polar heart rate monitor and interference

Electromagnetic interference

Disturbances may occur near high voltage power lines, traffic lights, overhead lines of electric railways, electric bus lines or trams, televisions, car motors, bike computers, some motor driven exercise equipment, cellular phones or when you walk through electric security gates.

Crosstalk

The Polar wrist receiver in non-coded mode of operation picks up transmitter signals within 1 metre/ 3 feet. Non-coded signals from more than one transmitter picked up simultaneously can cause an incorrect readout.

Exercise equipment

Several pieces of exercise equipment with electronic or electrical components such as LED displays, motors and electrical brakes may cause interfering stray signals. To try to tackle these problems, relocate the Polar wrist receiver as follows:

1. Remove the transmitter from your chest and use the exercise equipment as you would normally.
2. Move the wrist receiver around until you find an area in which it displays no stray reading or flashing of the heart symbol. Interference is often worst right in front of the display panel of the equipment, while the left or right side of the display is relatively free of disturbance.
3. Put the transmitter back on the chest and keep the wrist receiver in this interference-free area as far as it is possible.

4. If the Polar heart rate monitor still does not work with the exercise equipment, this piece of equipment may be electrically too noisy for wireless heart rate measurement.

Minimising possible risks in exercising with a heart rate monitor

Exercise may include some risk, especially for those who have been sedentary. To minimise possible risks consider the following advice:

Before beginning a regular exercise program consult your doctor

- if you are over 40 years of age and have not been following a regular exercise program during the past 5 years
- if you smoke
- if you have high blood pressure
- if you have high blood cholesterol
- if you have any signs or symptoms of any disease
- if you are recovering from a serious illness or a medical treatment e.g. surgery
- if you use a pacemaker or another implanted electronic device.

Note also that in addition to exercise intensity the heart rate can be affected by, amongst others, heart, circulation (particularly high blood pressure) and asthma and other breathing medications, as well as some energy drinks, alcohol and nicotine.

It is important to be sensitive to your body's reactions during the exercise. If you feel unexpected pain or tiredness regarding your exercise intensity level, it is recommended to stop the exercise or continue at a lighter intensity.

H 76

H 77





If you are allergic to any substance that comes into contact with the skin, check the Polar heart rate monitor materials from the “Technical specifications” chapter. It is recommended that you wear the transmitter against your bare skin to ensure flawless operation. However, if you wish to wear the transmitter over a shirt, moisten the shirt well under the electrodes.

Notice to persons with pacemakers, defibrillator or other implanted electronic device. Persons who have a pacemaker use the Polar heart rate monitor at their own risk. Before starting use, we always recommend an exercise test under doctor’s supervision. The test is to ensure the safety and reliability of the simultaneous use of the pacemaker and the heart rate monitor.



I. FREQUENTLY ASKED QUESTIONS

What should I do if...

...I don't know where I am in the Options or File cycle?

Press and hold the stop button until time of day is displayed.

...I cannot find the previous exercise's file?

You may have used the Measure mode during your exercise, which means that your exercise was not stored. To record the exercise you have to turn the stopwatch on by pressing the OK button in the Measure mode.

...I cannot find the expended calories from my exercise file?

Check that you have set your personal user information and turned the OwnCal function on. Check that your heart rate has reached over 90 bpm or above 60 % of your maximum heart rate during your exercise.

...I cannot turn my HR_{max}-p on?

Check that you have set your personal user information and turned the fitness test on before turning the HR_{max}-p on.





...there is no heart rate reading (- -)?

1. Check that the electrodes of the transmitter are moistened and that you are wearing it as instructed.
2. Check that you have kept the transmitter clean.
3. Check that you are not near other people that have heart rate monitors, high voltage power lines, televisions, mobile phones or other sources of electromagnetic disturbance.
4. Have you had a cardiac event which may have altered your ECG waveform? In this case consult your physician.

...the code search was not successful?

The frame around the heart symbol will disappear. Bring the receiver up to your chest near the transmitter's Polar logo. The receiver starts looking for a heart rate signal again. If you still do not succeed in code searching, but your heart rate is displayed and the heart symbol without a frame is flashing, you can start exercising but your wrist receiver may become disturbed by other heart rate monitors.

...heart symbol flashes irregularly?

1. Check that the wrist receiver is not further than 1 metre/ 3 feet from the transmitter.
2. Check that the elastic strap has not become loose during exercise.
3. Make sure that the electrodes of the transmitter are moistened.
4. Make sure that there is no other heart rate transmitter within 1 metre/ 3 feet.
5. Cardiac arrhythmia may cause irregular readings. In this case consult your physician.

...other person with heart rate monitor is causing interference?

Your training partner may have exactly the same code as you. In this case proceed as follows: Keep a distance from your training partner and continue your training session normally.

Or

1. Take the transmitter off your chest for 30 seconds. Keep a distance from your training partner.
2. Put the transmitter back and bring the receiver up to your chest near the transmitter's Polar logo. The receiver starts looking for a heart rate signal again. Continue your training session normally.

...heart rate reading becomes erratic or extremely high?



You may have come within range of strong electromagnetic signals which cause erratic readings. Check your surroundings and move further away from the source of disturbance.

...the display is blank or fading?

If the display is blank, activate the wrist receiver by pressing the OK button twice, after which the Time of day display appears. The first sign of a discharged battery is the battery symbol on the display. Also the digits fade when the backlight is used. Have the batteries checked.

...there are no reactions to any buttons?

Reset the Polar heart rate monitor. Resetting clears watch settings which will return to the default settings. User information and measuring units settings will be saved.





1. Press the RESET button with a pen tip. Display becomes full of digits.
If you don't press any button after reset within one minute the wrist receiver goes to the Time of day display.
2. Press any of the buttons once. The time of day is displayed.

...the battery of the wrist receiver must be replaced?

We recommend having all service done by an authorised Polar Service Centre. This is especially necessary to ensure the warranty remains valid and is not affected by incorrect repair procedures done by an unauthorised agent. Polar Service will test your wrist receiver for water resistance after battery replacement and make a full periodic check of the complete Polar heart rate monitor.

...the data transfer from receiver to computer was unsuccessful

1. Check the correct position of the wristwatch receiver with the interface.
2. Check that the surroundings are not too bright. This problem often occurs near windows.
3. Check that there is no obstruction between the infrared windows of the interface and the wristwatch receiver.



J. TECHNICAL SPECIFICATIONS

Polar heart rate monitor is designed to indicate the level of physiological strain and intensity in sports and exercise. No other use is intended or implied. Heart rate is displayed as number of heart beats per minute (bpm).

Polar transmitter

Battery type:	Built-in Lithium Cell
Battery life:	Average 2 500 hours of use
Operating temperature:	-10 °C to +50 °C/ 14 °F to 122 °F
Material:	Polyurethane

Elastic strap

Buckle material:	Polyurethane
Fabric material:	Nylon, polyester and natural rubber including a small amount of latex





Wrist receiver

Battery type: CR 2430
 Battery life: Average 2 years (2h/ day, 7 days/ week)
 Operating temperature: -10 °C to +50 °C/ 14 °F to 122 °F
 Water resistance: to 50 metres/ 160 feet
 Wrist strap material: Polyurethane

Back cover and wrist strap buckle material: Stainless steel complying with EU Directive 94/27/EU and its amendment 1999/C 205/05 on the release of nickel from products intended to come into direct and prolonged contact with the skin.

Watch accuracy better than $\pm 0,5$ seconds/ day at 25 °C/ 77 °F temperature.

Polar S810 is a Class 1 Laser Product.

Accuracy of heart rate measurement: $\pm 1\%$ or ± 1 bpm, whichever larger, definition applies to steady state conditions.

Default Settings

Time of day	24h	HR _{max}	220-age
Alarm	OFF	VO _{2max}	
Date of Birth	0	male	45
(adjustable range for year is 1921 - 2020)		female	35
Sex	Male	Baseline	0
Weight	0 (kg)	Activity sound	On
Height	0 (cm)	Units	1
Activity	Low	Help	On
		OwnCal	OFF

Fitness test	On
HR _{max-p}	OFF
Online	OFF
Exercise settings	BasicUse On
BasicSet function	OFF
Profile function	OFF
Heart touch function	Shows heart rate limits
Recording rate	5 s

Limit Values

Exercise time	99 h 59 min 59 s
Heart rate limits	30 - 240
Heart rate zone time	99 h 59 min 59 s
Recovery time	99 min 59 s
Calories	99 999 kcal
Cumulative calories	999 999 kcal
Cumulative exercise time	9 999 h
Laps recorded to memory	99
Files recorded to memory	99

J 84

J 85





K. LIMITED POLAR WARRANTY

- Polar Electro Oy warrants to the original consumer/ purchaser of this heart rate monitor that the product will be free from defects in material or workmanship for one year from the date of purchase. **Please keep the receipt which is your proof of purchase!**
- Warranty does not cover receiver battery, damage due to misuse, abuse, accidents or non-compliance with the precautions; improper maintenance, commercial use, cracked or broken cases.
- Warranty does not cover damage or consequential damage caused by service not authorised by Polar Electro Oy.
- During the warranty period the product will be either repaired or replaced at an authorised service centre free of charge.
- This warranty does not affect the consumer's statutory rights under applicable national or state laws in force, or the consumer's rights against the dealer arising from their sales/ purchase contract.



This CE marking shows compliance of this product with Directive 93/42/EEC.

Copyright © 2001 Polar Electro Oy, FIN-90440 KEMPELE, Finland. All rights reserved. No part of this manual may be used or reproduced in any form or by any means without prior written permission of Polar Electro Oy. Polar logotype is a registered trademark and Polar S810, Polar IR Interface and Polar Precision Performance are trademarks of Polar Electro Oy.

Sound Blaster is a trademark of Creative Technology, Ltd. Windows is a registered trademark of Microsoft Corporation.

K 86



L. DISCLAIMER

- The material in this manual is for informational purposes only. The products it describes are subject to change without prior notice, due to the manufacturer's continuous development program.
- Polar Electro Oy makes no representations or warranties with respect to this manual or with respect to the products described herein. Polar Electro Oy shall not be liable for any damages, direct or incidental, consequential or special, arising out of, or related to the use of this material or the products described herein.

This product is protected by one or several of the following patents:

FI 88223, DE 4215549, FR 92.06120, GB 2257523, HK 113/1996, SG 9591671-4, US 5491474, FI 88972, FR 92.09150, GB 2258587, HK 306/1996, SG 9592117-7, US 5486818, FI 96380, US 5611346, EP 665947, DE 69414362, FI 68734, DE 3439238, GB 2149514, HK 812/1989, US 4625733, FI 100924, FI 100452, US 5840039, FI 4150, FI 4069, DE 29910633, FI 104463.

Other patents pending.

L 87





M. POLAR GLOSSARY

Activity level: Your long-term physical activity level you need to assess in order to do the Polar Fitness Test.

Coded heart rate transmission: Polar T61 transmitter automatically locks in a code to transmit your heart rate to the wrist receiver. In a coded transmission the receiver accepts heart rate data only from the Polar T61 transmitter. Coding significantly reduces crosstalk caused by other users of heart rate monitors. However, it does not necessarily reduce all environmental interference.

Electrodes: The grooved areas of the transmitter, against the skin. Pick up the heart rate.

Exercise mode: The stopwatch is on and your exercise is recorded. You can select either the BasicUse, Profiles and BasicSet exercise modes.

M 88



Heart touch function: In the Exercise mode bring your receiver close to your transmitters Polar logo and the selected heart touch function is activated for a few seconds.

Maximal oxygen uptake: (maximal aerobic power, VO_{2max}) Is the maximal rate at which oxygen can be used by the body during maximal work. VO_{2max} is a good index of aerobic fitness.

Maximum heart rate: (HR_{max}) Is the highest number of heart beats per minute of a person.

Measure mode: Measuring your heart rate without recording the exercise.

Target heart rate zone: The area between the upper and lower target heart rate limits. Target zone selection is based on personal fitness goals.

M 89







Display symbols

- ♥ Indicates an ongoing heart rate measurement and flashes to the beat of your heart. A heart without frames indicates a non-coded heart rate transmission.
- ♥ Frame around the heart indicates a coded transmission of heart rate.
- ∞∞ Indicates there has not been heart rate reception for less than 5 minutes.
- - Indicates there has not been heart rate reception for at least 5 minutes.
-)) In the Time of day display the alarm is on.
Measure mode: beeps if you are exercising outside your target heart rate zone. Fitness Test mode: sounds an alarm at the end of the fitness test.
- ▮▮▮ Indicates profile exercise.
- ▮▮▮▮▮ Indicates continuous intervals and recoveries. See Cont in "display texts" chapter.
- ▲ Exercise mode: indicates the heart rate above the target zone.
File mode: indicates the time spent above the target zone.
- ▼ Exercise mode: indicates the heart rate below the target zone.
File mode: indicates the time spent below the target zone.
- ▮▮▮▮▮ The graphic bar consists of five blocks. A flashing block indicates in which level you are in the Options settings or File recalling loops.
- ▮▮▮▮▮ Continuously running bars indicate that the exercise is recorded.
- ⇒ Indicates low battery.



Display texts

- ALARM:** Indicates the alarm in watch settings.
- AM or PM:** Are used in 12h time mode. E.g. in 24h time mode 13:00 means 1:00 PM.
- AVG:** Together with a reading indicates average heart rate.
- Baseline:** Set your personal relaxation rate baseline.
- BasicUse:** Exercising without settings.
- BasicSet:** Preprogrammable exercise with three alternating timers and heart rate limits. Additionally you can set a time- or heart-rate-based recovery function that is turned on manually in the exercise situation.
- BestLap:** Indicates the shortest lap time information.
- CONNECT:** Communication mode between wrist receiver and computer.
- Cont:** Up to 30 intervals and recoveries succeed one another unless the phase is manually stopped.
- E0:** Indicates BasicUse exercise.
- E1-E7:** Indicates the profile exercise.
- E8:** Indicates the BasicSet exercise.
- EXE.SET:** Indicates the recorded exercise profile information in the File mode.
- Exe. Time:** Exercise duration
- FILE:** Indicates the file mode.
- FIT. TEST:** Indicates fitness test mode.
- FULL:** After having stored 99 laps (the use of Profiles reduces the amount) FULL appears in the display. You can continue taking lap times or Profiles but they won't be recorded into the file.
- HR_{max}-p:** Predicted maximum heart rate value.
- Int.:** Indicates the work period of a profile exercise.
- InZone/ Above/ Below:** Indicate the time spent within, above and below the target zone.
- KCal:** Indicates the energy expenditure of your exercise.





LAPS: Indicates how many lap times have been stored.

Lim High: The higher heart rate limit of the target zone.

Lim Low: The lower heart rate limit of the target zone.

Limits: Heart rate limits for target zones.

MAX: Indicates the highest heart rate.

MAX %: Indicates what percentage your current heart rate is of your maximum heart rate.

Mem full: Appears on the display when the memory is full. You can continue your exercise and taking lap times but they won't be recorded into the file.

MIN: Indicates minimum heart rate of exercise in the File mode.

ms: Millisecond

Online*: Online recording allows you to transfer your heart rate to computer in real time. Place the wrist receiver within the reception range, 15 cm/ 6", of the interface

and within range of the transmitter, 1 metre/ 3 feet.

OPTIONS: Indicates the Options mode.

OwnCal: Calculates energy expenditure during exercise in kilocalories. 1 kilocalorie (kcal) = 1000 calories (cal). You can follow how much energy you lose by exercising in one exercise session or in a day, week or even a year. The OwnCal can be used as a dose measure of exercise e.g. in nutrition planning. OwnCal calorie calculation starts when your heart rate reaches 90 bpm or is above 60% of your maximum heart rate, whichever is smaller. These limits are set to count only exercise calories. The higher the heart rate, the faster the energy expenditure. OwnCal is individually calibrated by your preset weight, maximal oxygen uptake (VO_{2max}) and maximum heart rate (HR_{max}). The most accurate OwnCal is received by entering your clinically measured (in

maximal treadmill or bicycle stress test) VO_{2max} and HR_{max} into the wrist receiver. Energy expenditure measurement is most accurate in continuous activities such as cycling and running.

OwnINDEX: Result of Polar Fitness Test and comparable score with maximal oxygen uptake (VO_{2max}).

P1-P6, exercise phases: Indicates pre-programmable exercise phases of the profile exercise. Each phase includes adjustable number of interval periods and a recovery period.

Profil: Indicates a preprogrammable exercise profile with six optional phases.

RecoHr, heart-rate-based recovery: You set the heart rate, which will end your recovery calculation. In the File mode you will see the drop in your heart rate and the recovery duration.

RecoTime, time-based recovery: You set the recovery time in the Options mode. In the File mode you



will see the drop in your heart rate and the recovery duration.

RLX, Relaxation rate: Relaxation rate is most useful when measuring your heart rate at rest or during light daily activities, such as office work. After exercise, e.g., in the following morning, the value can also indicate the state of your physical recovery. A resting heart rate significantly higher or a relaxation rate smaller than normal may be an indication of an incomplete recovery. Mental stress typically increases heart rate and decreases relaxation rate. Relaxation rate measurement is based on the heart rate variability, which indicates your autonomic nervous system's state. The degree of heart rate variability is highly individual and tends to decrease over time.

RR: R-R recording rate measures heart beat intervals, i.e. the time between successive heart beats in milliseconds. This information is also shown as instantaneous heart rate

*An optional IR Interface is required.





in beats per minute in recorded samples. Recording every interval accurately is more difficult than obtaining the averaged heart rate only, because R-R recording rate is less shielded from measurement disturbances. We recommend using contact gel (ECG gel) to ensure contact between your skin and the transmitter. It is useful to minimise the use of the heart touch function and backlight, because the heart rate receiver does not detect new heart beat intervals when the backlight is on. Also the use of activity/button sound may occasionally cause a loss of one R-R interval. Readings that are interpreted as incorrect in the heart rate data can be adjusted and corrected with Polar Precision Performance software.

Split Time: The elapsed time from the beginning of the exercise until storing the split time.

Sw: Indicates stopwatch.



Index

12h/ 24h time mode	36	Lap and split time storing	42, 51, 59
Activity/button sound On/OFF	34	Measuring units settings	34, 38
Alarm settings	36	Modes and their functions	11
Baseline for relaxation rate	19	Name profiles	30
BasicSet exercise	20, 49	Online* measurement to	
Buttons	6	computer	40, 92
Care and maintenance	73	Online* On/OFF	33
Code search	39, 80, 88	Options lock function	62
Date settings	37	OwnCal calorie expenditure	44, 92
Deleting a file	61	OwnCal On/OFF	32
Display configuration	43	OwnIndex	65, 93
Download training data		Pause an exercise	43
to a computer	62	Precautions	73
Easy start	8	Predicted maximum heart rate	67, 91
File	52	Profile exercise	46, 57, 93
Fitness test	65	Recall training information	52
Fitness test On/OFF	32	Record rate settings	31
Frequently asked questions	79	Records file	60
Functions during the exercise	41	Recovery during BasicSet	50, 55
Heart rate limits	51, 55	Recovery during exercise profile ...	48, 58
Heart rate limits settings		Recovery settings for BasicSet	24
for BasicSet	23	Recovery settings for	
Heart rate limits settings for Profile ...	28	exercise profile	29
Heart rate zone alarm On/OFF	41	Relaxation rate	19, 44, 93
Heart touch function	51	Reminder On/OFF	37
Help On/OFF	35	Reset cumulative counts	60
How to enter settings	14	Reset the stopwatch	45
HR max -p On/OFF	30	Reset the heart rate monitor	81
Interrupt an exercise phase	51	Return to the Time of day display	7
Interrupt an interval or		Selecting exercise type	20
recovery calculation	51	Setting interval calculation	26

* An optional Polar IR Interface is required.





Setting tips 38

Start a profile exercise 46

Start measuring your heart rate 9, 39

Stop a profile exercise 49

Stop measuring your heart rate 10, 50

Swap the exercise type 50

Swap the heart rate limits 49

Swap the time zones 1/ 2 38

Technical specifications 83

Time of day settings 36

Timer settings for BasicSet 23


Transmitter 5, 10

Update HR max -p and OwnIndex
to user information 70

Upload settings from computer 63


User information settings 16

Warranty 86





Road Map



Dear Customer

Congratulations on purchasing a Polar Heart Rate Monitor! We are proud to assist you in achieving personal fitness and performance goals. If you want to get even more out of your Polar S810 during training, check www.polar.fi for:

- the best product and training tips in Polar Personal Trainer
- the software updates and other value added service for your S-series heart rate monitor in Polar Product Support



Manufactured by:
 Polar Electro Oy
 Professorintie 5
 FIN-90440 KEMPELE

Tel. +358-8-520 2100
 Fax +358-8-520 2300
www.polar.fi