



# XcelGen | Plasmid Mini Kit

## User Guide

Cat No: XG1211-01

XcelGen

Ver.: 0.2/13

Revised Protocol

## Table of Contents

Introduction .....	02
Storage and Stability .....	02
Important Notes .....	02
Kit Contents .....	04
Before Starting .....	04
Safety Information .....	04
Plasmid Mini Kit Spin Protocol .....	05
Plasmid Mini Kit Spin/ Vacuum Protocol .....	07
Purification of Low-Copy-Number Plasmid/Cosmid .....	07
Purification of plasmid >12kb .....	07
Troubleshooting Guide .....	09
Related Products .....	10
Limited Use and Warranty .....	10

## Introduction

Key to the kit is our proprietary DNA binding systems that allow the high efficient reversible binding of DNA to the mini column while proteins and other impurities are removed by wash buffer. Nucleic acids are then eluted with sterile water or elution buffer.

This kit is designed for fast and efficient purification of plasmid DNA from 1 to 4ml of *E. coli* culture. The mini column has a plasmid DNA binding capacity of 50µg. The yield from 1 ml culture is typically around **8 to 12µg**.

The purified DNA is ready for downstream applications such as cloning/subcloning, RFLP, sequencing and transfection of robust cells such as HEK293 cells.

## Storage and Stability

Buffer A1 should be stored at 4°C once RNase A is added. All other materials can be stored at room temperature (22-25°C). All kit components are stable up to 12 months.

## Important Notes

**Plasmid Copy Numbers:** The yield of plasmid DNA depends on the origin of the replication and the size of the plasmid. **The protocols are optimized for high copy number plasmid purification. For low copy number plasmids, both the culture volume and the buffer volume need to be scaled up 2 times.** Please contact our customer service for further information and reference Table 1 for the commonly used plasmids.

**Table 1: Commonly used plasmids and expected yield.**

Plasmid	Origin	Copy Numbers	Expected Yield (µg/ml)
PSC101	PSC101	5	0.1-0.2
pACYC	P15A	10-12	0.4-0.6
pSuperCos	pMB1	10-20	0.4-1
pBR322	pMB1	15-20	0.6-1
pGEMR	Muted pMB1	300-400	6-7
pBluescriptR	ColE1	300-500	6-8
pUC	Muted pMB1	500-700	8-12

**Host Strains:** The strains used for propagating plasmid have significant influence on yield. Host strains such as Top 10 and DH5 $\alpha$  yield high-quality plasmid DNA. *endA*+ strains such as JM101, JM110, HB101, TG1 and their derivatives, normally have low plasmid yield due to either endogenous endonucleases or high carbohydrates released during lysis. We recommend transform plasmid to an *endA*- strain if the yield is not satisfactory. Please reference Table 2 for the *endA* information.

**Table 2: *endA* strains of *E. Coli*.**

<i>EndA</i> - Strains of <i>E. Coli</i>							
DH 5 $\alpha$	DH1	DH 21	JM106	JM 109	SK2267	SRB	XLO
TOP10	DH108	JM103	JM107	SK1590	MM294	Stbl2™	XL1-Blue
BJ5182	DH20	JM105	JM108	SK1592	Select96™	Stbl4™	XL10-Gold

<i>EndA</i> + Strains of <i>E. Coli</i>							
C600	JM110	RR1	ABLE® C	CJ236	KW251	P2392	BL21(DE3)
HB101	TG1	TB1	ABLE® K	DH12S™	LE392	PR700	BL21(DE3) pLysS
JM101	JM83	TKB1	HMS174	ES1301	M1061	Q358	BMH 71-18
All NM strains				All Y strains			

**Optimal Cell Mass (OD<sub>600</sub> x ml of Culture):** This procedure is designed for isolating plasmid grown in standard LB medium (Luria Bertani) for 12-16 hours to a density of OD<sub>600</sub> 2.0 to 3.0. If rich medium such as TB or 2xYT are used, make sure the cell density doesn't exceed 3.0 (OD<sub>600</sub>). A high ratio of biomass over lysis buffers result in low DNA yield and purity. The mini column has an optimal biomass of 10-15. For example, if the OD<sub>600</sub> is 3.0, the optimal culture volume should be 1-5 ml. For over amount of cell numbers, either reduce the biomass or scale up the volumes of Buffer A1, B1 and N1.

**Culture Volume:** Use a flask or tube 4 times bigger in volume than the culture medium to secure optimal condition for bacteria growth. Don't exceed the maximum culture volume suggested in the protocol. Incomplete lysis due to over amount of bacterial culture results in lower yield and less purity.

## Kit Contents

Product	XG1211-00	XG1211-01
Preps	4	50
Columns	4	50
Buffer A1	1.2 ml	15 ml
Buffer B1	1.2 ml	15 ml
Buffer N1	1.6 ml	20 ml
Buffer BL	1.8 ml	22 ml
Buffer KB	3 ml	30 ml
DNA Wash Buffer*	2 ml	15 ml
Elution Buffer	600 µl	10 ml
RNase A (20 mg/ml)	0.2 mg (10 µl)	1.5 mg (75 µl)
User Manual	1	1

## Before Starting

Prepare all components and get all necessary materials ready by examining this instruction booklet and become familiar with each steps and pay special attention to the followings.

## Important

- **RNase A: 20mg/ml.** It is stable for more than half a year when stored at room temperature. Spin down RNase A vial briefly. Add the RNase A solution to Buffer A1 and mix well before use. Store at 4°C.
- Add 8ml (XG1211-00) or 60 ml (XG1211-01) or 96-100% ethanol to each **DNA Wash Buffer** bottle before use.
- **Buffer B1** and **Buffer BL** precipitates below room temperature. It is critical to warm up the Buffer at 50°C to dissolve the precipitates before use.
- Keep the cap tightly closed for **Buffer B1** and **Buffer BL** after use.
- Carry out all centrifugation at room temperature.

## Safety Information

- Buffer N1 contains acidic acid, wear gloves and protective eyewear when handling.
- Buffer N1 and KB contains chaotropic salts, which may form reactive compounds when combines with bleach. Do not add bleach or acidic solutions directly to the preparation waste.

## Plasmid Mini Kit Spin Protocol

1. Inoculate **1-4ml** LB containing appropriate antibiotic with a fresh colony from a freshly streaked selective plate. Incubate at 37°C for 14-16 hours with vigorous shaking.

**Note:** Prolonged incubation (> 16 hours) is not recommended since the *E.coli* starts to lyse and the plasmid yields may be reduced.

**Note:** Do not grow the culture directly from the glycerol stock.

**Note:** This protocol is optimized for *E. coli* strain cultured in LB medium. When using TB or 2xYT medium, special care needs to be taken to ensure the cell density doesn't exceed 3.0 (OD<sub>600</sub>). Buffers need to be scaled up proportionally if over amount of cultures are being processed.

2. Harvest the bacterial culture by centrifugation for 1 min at 10,000 rpm. Pour off the supernatant and blot the inverted tube on a paper towel to remove residue medium. Remove the residue medium completely.

**Note:** Residue medium will cause, Poor cell lysis and thus lower DNA yield.

3. Add **250µl Buffer A1** (Add RNase A to **Buffer A1** before use) and completely resuspend bacterial pellet by vortexing or pipetting.

**Note:** Complete resuspension is critical for bacterial lysis and lysate neutralization.

4. Add **250µl Buffer B1**, mix gently by inverting the tube 10 times (do not vortex), and incubate at room temperature for 5 minutes.

**Note:** Do not incubate for more than 5 minutes.

**Note:** Buffer B1 precipitates (cloudy look) below room temperature. Warm up Buffer B1 at 50°C to dissolve precipitation before use.

5. Add **350µl Buffer N1**, mix completely by inverting/shaking the vial for 5 times and sharp hand shaking for 2 times.

**Note:** Incubating the lysate in ice for 1 min will improve the yield.

**Note:** It is critical to mix the solution well. If the mixture still appears conglobated, brownish or viscous, more mixing is required to completely neutralize the solution.

6. Centrifuge the lysate at 13,000 rpm for 10 minutes at room temperature.

**Note:** If the lysate doesn't appear clean, reverse the tube angle, centrifuge for 5 more minutes and then transfer the clear lysate to DNA column.

7. Add **400µl of Buffer BL** into the spin column (provided), incubate at room temperature for 2 minutes, centrifuge at 12000 rpm for 2 minutes and discard the flow through. The column is ready and will work well for binding DNA.
8. Carefully transfer the clear lysate into a DNA column with a collection tube, avoid the precipitations, spin at 13,000 rpm for 1 minute, discard the flow-through and put the column back to the collection tube.
9. Add **500µl Buffer KB** into the spin column, centrifuge at 13,000 rpm for 1 minute. Remove the spin column from the tube and discard the flow-through. Put the column back to the collection tube.

**Note:** This step is important to remove residual protein contaminations especially for endA+ strains and be highly recommended for high quality plasmid DNA.

10. Add **650µl DNA Wash Buffer** (Add ethanol to DNA wash buffer before use) into the spin column, centrifuge at 13,000 rpm for 1 minute at room temperature. Remove the spin column from the tube and discard the flow-through. Repeat step "10" to improve the recovery.
11. Reinsert the spin column, with the lid open, into the collection tube and centrifuge for 2 minutes at 13,000 rpm.

**Note:** Residual ethanol can be removed more efficiently with the column lid open. It is critical to remove residual ethanol completely.

12. Carefully transfer the spin column into a sterile 1.5ml microfuge tube and add **50-100µl** sterile ddH<sub>2</sub>O or **Elution buffer** into the center of the column and let it stand for 2 minutes. Elute the DNA by centrifugation at 13,000 rpm for 1 minute. Reload the eluate into the column and elute again.

**Note:** It is recommended to use elution buffer instead of ddH<sub>2</sub>O.

**Note:** If ddH<sub>2</sub>O is applied, please make sure the pH is no less than 7.0 (7.0-8.5 is preferred). NaOH could be used to adjust the pH of ddH<sub>2</sub>O.

**Note:** The DNA is ready for downstream applications such as cloning/subcloning, RFLP, library screening, in vitro transfection, sequencing, transfection of robust cells such as HEK293 cells.

**Note:** It is highly recommended to remove the endotoxin (XG1212-01) if the DNA is used for endotoxin-sensitive cell lines, primary cultured cells or microinjection.

13. The DNA concentration can be calculated as follows :-

$$\text{Concentration } (\mu\text{g/ml}) = \text{OD}_{260} \text{ nm} \times 50 \times \text{dilution factor}$$

## Plasmid Mini Kit Spin/Vacuum Protocol

1. Set up the vacuum manifold according to manufacture's instruction and connect the column to the manifold.
2. Carry out step 1-7 on Page 6 in previous protocol.
3. Carefully transfer the clear lysate to the DNA column and turn on the vacuum to allow the lysate pass through the column.
4. Add **500µl Buffer KB** into the spin column and allow the lysate pass through the column by vacuum.


**Note:** *This step is important to remove residual protein contaminations especially for endA+ strains and be recommended for high quality plasmid DNA.*

5. Add **650µl of DNA Wash Buffer** to the column and allow the vacuum to draw the liquid through the manifold. Turn o the vacuum. Repeat step "5" to improve the recovery.
6. Transfer the column, with the lid open, to a 1.5 ml collection tube and centrifuge at 13,000 rpm for 2 minutes.
7. Carefully transfer the spin column into a clean 1.5 ml microfuge tube and add **50-100µl (>50µl) Sterile ddH<sub>2</sub>O or Elution Buffer** into the column and let it stand for 2 minutes. Elute the DNA by centrifugation at 13,000 rpm for 1 minute. Reload the eluate into the column and elute again.

**Note:** *It is recommended to use Elution Buffer instead of ddH<sub>2</sub>O*

**Note:** *The DNA is ready for downstream applications such as cloning, RFLP, library screening, in vitro translation, sequencing and transfection of robust cells such as HEK293 cells.*





**Fig:** Agarose gel analysis of plasmid DNA purified with XcelGen Plasmid mini Kit

Lane 1: pbluescript II  
Lane 2: pUC 18  
Lane 3: pBR 322  
Lane 4: pGEM  
Lane 5: pET 43.1  
Lane L: 1 kb ladder

## Purification of Low-Copy-Number Plasmid/Cosmid

The yield of low copy number plasmid is normally around 0.1–1 µg /ml of overnight culture. For isolating low copy number or medium copy number plasmid DNA, use the following guideline:

- **Culture volume:** Use **2 x volumes** of the high copy number culture.
- Use **2 x volumes** of the **Buffer A1, Buffer B1 and Buffer N1**. Additional Buffers can be purchased from XcelGen.
- Use **same volume** of **DNA Wash Buffer** and **Elution Buffer**.

## Purification of plasmid > 12kb

For isolating plasmid DNA > 12 kb, use the following guideline:

- **Culture volume:** Use **2 x volumes** of the culture.
- Use **2 x volumes** of the **Buffer A1, Buffer B1 and Buffer N1**. Additional Buffers can be purchased from XcelGen.
- Use **same volume** of **DNA Wash Buffer** and **Elution Buffer**.
- Pre-warm the **Elution Buffer** at 65–70°C and let the column stand for 5 minutes after adding **Elution Buffer**.

## Troubleshooting Guide

Problems	Possible Reasons	Suggestions
Low Yield	Poor Cell lysis.	Resuspend pellet thoroughly by vortexing and pipetting prior to adding Buffer B1. Make fresh Buffer B1 if the cap had not been closed tightly. (Buffer B1: 0.2M NaOH and 1%SDS).
Low Yield	Bacterial culture overgrown or not fresh	Grow bacterial culture upto 12-16 hours. Spin down cultures and store the pellet at -20°C if the culture is not purified the same day. Do not store culture at 4°C over night.
Low Yield	Low copy-number plasmid	Increase culture volume and the volume of Buffer A1, B1, N1 as instructed on page 8.
No DNA	Plasmid lost in Host <i>E. coli</i>	Prepare fresh culture.
Genomic DNA contamination	Over-time incubation after adding Buffer B1.	Do not vortex or mix aggressively after adding Buffer B1. Do not incubate more than 5 minutes after adding Buffer
RNA contamination	RNase A not added to Buffer A1.	Add RNase A to Buffer A1.
Plasmid DNA floats out of wells while running in agarose gel, DNA doesn't freeze or smell of ethanol	Ethanol traces were not completely removed from column.	Make sure that no ethanol residue remains in the silicon membrane before elute the plasmid DNA. Re-centrifuge or vacuum again if necessary.
No phase partitioning after centrifugation	Temperature is lower than 23°C.	Make sure the temperature is greater than 23°C for centrifugation or incubate the sample at 60°C for 5 min and then perform centrifugation

## Related Products

1. EndoFree plasmid mini kit (XG1212-01)
2. 96 Well Plasmid Isolation Kit (XG1201-96)
3. DNA Gel/PCR Purification Miniprep kit (XG3511-01/XG3514)
4. PremixTaqV2.0 (XG334A)
5. Agarose (XGA-100)

## Limited Use and Warranty

This product is intended for *in vitro* research use only. Not for use in human. This product is warranted to perform as described in its labeling and in XcelGen's literature when used in accordance with instructions. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by XcelGen. XcelGen's sole obligation and purchaser's exclusive remedy for breach of this warranty shall be, at the option of XcelGen, to replace the products, XcelGen shall have no liability for any direct, indirect, consequential, or incidental damage arising out of the use, the results of use, or the inability to use it product.

For technology support or learn more product information, please visit our website at [www.xcelrisgenomics.com](http://www.xcelrisgenomics.com)

### **XcelGen** Quality Kits made by **Xperts**

- Plasmid DNA Isolation Kits • Genomic DNA Extraction Kits • RNA Extraction Kits • Polymerase • DNA Ladders • DNA Markers
- Premix Taq • dNTP's • RAPD kits • Agarose • Glycerol • Tms NA Stabilizers & RNA Protectant solutions

### **PrimeX** Oligo Synthesis & Purification Services

- 10 nmole • 25 nmole • 50 nmole • 100 nmole • 200 nmole • 1000 nmole

### **NXT** NGS Services

- Denovo Genome Sequencing • Whole Genome Resequencing • GBS/RAD Sequencing • Exome Sequencing • Amplicon Sequencing
- Whole Transcriptome Analysis/RNA-Sequencing • Small RNA Sequencing • Metagenomics • Metatranscriptomics
- ChIP Sequencing • Mitochondrial Sequencing • Next Generation Genomic Services on Illumina MiSeq
- Genotyping by Sequencing • Tilling/Scotilling using NGS • Genome Database development Services

### **NGS Bioinformatics**

- In silico Primer Design • Microarray Analysis • Metagenomics • Physical, Genetic and QTL mapping
- Assembly and annotation of prokaryotic and eukaryotic genome • Genome Mapping and SNP discovery
- Transcriptome discovery and analysis • sRNA analysis and discovery

### **XcelSeq** Sanger Sequencing Services

- Plasmid /PCR Sequencing Services • r-E. coli Culture Sequencing Services • Primer Walk Sequencing Services
- Microbial Identification Service • Multilocus Sequence Typing

### **Customised Services**

- SNP Genotyping by SNaPshot Assay • Microsatellite Genotyping • Golden Gate Assays and Arrays
- Gene Expression on Real Time PCR • Gene expression on Agilent / Microarray / Affymetix • Library construction