

X-Ray Fluorescence - Data Collection and Analysis

Matthias Zeller, Jerry Hunter, Sept. 19th 2008

A) Data Collection

Logging In:

The login and password to log into the program as a normal user are “ysu1” for both the username and password. If you need to use any feature only available with an administrator login ask one of the instrumentation scientists for assistance.

Loading Samples:

1. If using powders or liquids, the sample should fill the XRF sample cups about half way.
2. Place the 3.6 μm mylar film in between the two cups and press down fast and hard. *Try not to touch the mylar film with bare hands!*



3. If the sample is a solid, make sure the sample fits into the 40mm rings inside the machine. If it is too small place it in a plastic sample cup with mylar film such as used for powders or liquids.
4. Lift the chamber door and place the sample in **any** empty space.
5. Close chamber door.
6. For manual loading of samples, see Appendix 3.

Calibration:

1. It is important to calibrate the machine at the beginning of each day.
2. The screen on the XRF machine should be highlighted on **Loader**. If not highlight **Loader**.
3. To calibrate the machine select the space in which the copper disc is located (delete any experiment defined for this position).
4. Under **choose your application** press **recal**.
5. Highlight **copper disk calibration**. Press **Next**
6. On the loader screen, highlight the copper disk space.



7. Press **measure**.
8. The hourglass shows calibration in progress in the bottom right-hand corner. It will take about 8min for the calibration to complete.
9. While the calibration is running you can start setting up your sample(s). They will be queued and once the calibration disk is returned to its original spot the next samples will be measured.

Sample Analysis:

1. Highlight the spot on the touch screen that corresponds to where you will put your sample.
2. Under **Choose Application** select which application you want to run. (*the application you choose will depend on your sample and what results you are looking for*)

Description of Methods:




- *Elements Powder and Solid*- used to find out what elements are present in the sample.


- *Fluoride Pellet/Powder*- used to find what ratios of fluorides are contained in the sample. Use ‘pellet’ method if sample is in a pressed pellet, and powder if the sample is in powder form.
 - *Oxides Powders/Solid*- used to find what ratios of oxides are contained in the sample. Use “powder” if sample is in powder form and “solid” if sample is in solid form.
 - *Water*- used to find what elements are contained in a liquid sample.
 - More methods that are less often used are listed under **Other**
 - ***You are not limited to just these methods. Methods can be created and adjusted to accommodate different types of analysis.***
3. Enter the sample ID and press **next**.
 4. Under **Select the Preparation** choose **Mylar 3.6µm** for powders or, for solid samples, **Solid, SpectroBlend, or X-ray mix** (*this will depend on the type of solid. As a standard choose “Solid.”*)
 5. Press **complete and/or next**.
 6. If asked, enter the estimated weight of the sample. **Press Complete**. A picture of your sample will appear on the screen
 7. Highlight your sample on screen by touching it. **Press measure**. An hourglass appears to show your sample is being analyzed.
 - Starting and finishing time are displayed in the bottom right hand corner of the **loader** screen. A standard sample will take about 8-9min to run.
 - If more than one sample needs to be run, you can repeat the above procedure for any number of samples to analyze (up to 28).
 8. You can abort the procedure at any time by highlighting your sample on the loader screen and pressing **abort** on the screen.
 9. The spectrum of your sample can be viewed as it is running by pressing **graph**.

10. Once the sample is returned to its original spot and “ok” is displayed, the results can be viewed at the bottom of the screen. *It takes about one and half to two minutes for results to appear.*
11. Press **results** and press **last hour** to refresh the screen. Touch the screen to highlight your sample. The results appear at the bottom of the screen.

At this point you will have to “clean-up” your results because the default software has a tendency to miss-assign lines to wrong elements or to assign escape peak or sum peaks as real signals. The spectrum needs to be visually inspected and every line needs to be assigned to a certain element. Only these lines should then be used in the re-evaluation. Leaving out elements that exhibit X-ray lines that overlap with lines of actually present elements allows these to be calculated with higher accuracy. You can do a quick re-evaluation via the Touch Screen Interface (see chapter B) Re-evaluation of Data using the Touch Screen Interface). A thorough and reliable re-evaluation should be done via the program **Spectra EDX** on the flat screen computer to the left of the instrument (see section C) Re-evaluation of Data using SpectraEDX).

B) Re-evaluation of Data using the Touch Screen Interface

1. For a quick and non-quantitative analysis you can make a re-evaluation of the data using the touch screen interface. **The X-ray facility staff HIGHLY recommends to use Spectra EDX as described in section C) for sample evaluation!** If you decide to use the touch screen interface click **results** to make sure your sample is still highlighted or highlight it.
2. Click **graph**. You may need to refresh the screen by pressing 
3. You can identify the peaks on the graph by touching the center of the peak on the screen (you may have to zoom in on the peak first) and pressing 
4. Select and/or unselect X-ray lines by pressing , this will bring you to the periodic table to select or deselect X-ray lines of elements.
5. Change or add elements:



The periodic table is composed of buttons. A **red or orange button** indicates that the X-ray lines for this element are selected and a **grey button** indicates that they are deselected. The spectra series K, L and M can be switched on or off for all elements. The Energy Information area lists available lines and reports their energy in eV. Select the elements you are interested in.







6. By highlighting the elements in the table, you can choose which elements you want to keep or eliminate. This might be important as you might have additional information such as those certain elements are for sure absent. By leaving these out all elements that exhibit X-ray lines

that are close to the ones unselected will be calculated with higher accuracy. After changing elements, press the **re-eval** button to updated the spectrum. This may take a few minutes. Press **ok**. An updated graph appears with only the elements you specified.


7. To print your results press **print** in the bottom-right hand corner of the screen. .





WAIT UNTIL YOU HEAR THE CUTTING NOISE TO REMOVE THE PAPER! IF YOU TEAR TOO EARLY THE MACHINE WILL JAM!!!!!!

Summary of Graph Controls:

Auto-Scale		The Auto-Scale function scales the data so that the whole spectrum fits inside the Graph Area.
Live-Mode	 	When Live-Mode is on, the Graph Area will display either the current measurement (if a sample is being measured) or the last measurement (if no sample is being measured). If Live-Mode is off, the spectra from the selected sample in the Loader screen (automatic systems) or the selected sample from the Results screen is displayed.
Guess Line		This button automatically adds an element emission line. Simply position the Cross-Hairs near the peak of interest and press this button.
Add Lines		This button will display the periodic table page so that you can select specific element emission lines to be displayed in the Graph Area (see section "Adding Element Lines").
Zoom		This Thumbwheel is used to zoom in or out of the spectra. Position the Cross-hairs on the region of interest and do the following: <ul style="list-style-type: none"> • Move the Thumbwheel up to zoom out. • Move the Thumbwheel down to zoom in • Press and hold the thumbwheel to centre the region.

C) Re-evaluation of Data using SpectraEDX

1. Measure your sample using the touch screen interface
2. Switch over to the computer on the left (small flat screen). Open the program SpectraEDX. (If the path is lost, e.g. after a restart, you can define it again by highlighting the instrument and clicking connect)
3. Open the Evaluation option in SpectraEDX. Highlight the sample you want to analyze (you can get older/expired data via the "advanced search" option)
4. If there are more than one entry for your file (e.g. after a reevaluation via the touch screen) use the top entry (original data). (If you want to copy data to another drive use the windows explorer. Go to X:\SPECPlus\Data Files\TempSSD. Use copy and paste, do not drag or delete files)
5. Highlight the entry you want to analyze. Click "evaluation options" (right upper corner). Select a calibration (e.g. elements.evm for a chunk of metal) and a preparation (e.g. solid or mylar 3.6). Click "Apply & Close" to close the "evaluation options" window
6. Click the "Interactive Qualitative Mode" button (bottom of page). A new window opens. Select the "Elements Toolbar" option . A periodic system will pop up. Click on individual elements to see their XRF lines. Decide if an element is present or not by comparison of the spectrum with the XRF lines. Take into account that some lines might be "escape" or "sum peaks", and that all lines of an element (K, L and M) have to be found in the correct ratios for the element to be present. Continue until all lines are accounted for as either an element peak or an escape or sum peak. Make a note for elements that might be potentially hidden by e.g. an escape peak (For the later evaluation leave out palladium. The machine has a palladium source and Pd lines will always be present whether Pd is in the sample or not).

7. Write down all elements you wish to use for the analysis. When done you may want to print out the spectrum including the XRF lines for your records (hard copy and/or pdf) and label individual lines by element on the hard copy.
8. Click on the "alternative mode window" button . In the new window click on "initialize" , then "Show All Elements" . Highlight elements you do not want to have included in the analysis. Via a right mouse click select "set not present". After every change click "initialize" again. When only the elements of interest are left over click on "compute concentrations" .
9. You may want to write down the RWP and the R/R0 value (bottom bar of window) after the refinement is finished. Values for R/R0 should be below a value of 20 to be regarded as "reliable".
10. Print the results for your records to pdf.

D) Data analysis pitfalls

Artifacts - escape peaks and sum peaks, peak overlaps, and excessive dead time are all phenomena that an EDS user must be aware of. Modern analytical software used in processing energy dispersive X-ray spectra can generally take them into account -- but such software is not perfect. Also, many users will look at the raw spectra, where the software may or may not have labeled the artifacts.

Escape peaks - there is a statistical probability that some of the X-rays, generated in the sample and impacting the solid state detector (e.g., a SiLi device), will 'inadvertently' knock out Si K-shell electrons in the detector, reducing that X-ray's energy measured in the detector by the Si absorption edge energy (1.84 KeV). Say you're looking at something with lots of Fe (Ka of 6.40 KeV); the Si-escape peak of Fe Ka will appear at 4.56 KeV. You see this escape peak only for the major elements present.

Sum peaks - this phenomenon occurs where the count rates are moderate to high, when two X-rays impact the detector virtually instantaneously; the pulse created and measured is the sum of the two X-ray energies. Say you have a sample with lots of Si (Ka of 1.74 KeV) and Al (Ka of 1.487); a peak at 3.23 KeV is the sum peak, not to be assumed to be a K peak (Ka of 3.31 KeV).