INFINICYTTM, 1.1 VERSION

To improve the benefits that this software offers, InfinicytTM 1.1 has included the following:

1.	Th	e FL-based merger option has been deleted, now there is only a type of merger	1
2.	<u>Ne</u>	ew options and types of Graphs have been created to show information	2
3.	It i	s possible to create and save <u>Analysis Profiles</u> for successive uses	5
-	5.1.		~
		3.1.1. It is possible to use formulas to calculate results	2
,		3.1.2. Absolute counts from the population of reference selected by the user	9
-	5.2.	current compensation can be changed manually and saved if the origin FCS file contains this information internally	9
-	33	New option of creating consecutive gates to identify a population and save the	
•		sequence to apply it to similar files	10
2	3.4.	It is allowed to configure negative visibility of the parameters with Logical scale	11
	3.5.	New option to check all the parameters measured of a data file (Properties)	13
4.	Ne	w options to facilitate the Manual Analysis of the data file	
2	4.1.	It has been improved the way in which the user can save and find in his computer a	
		document already analyzed (.cyt)	14
2	4.2.	Option to choose the columns to visualize in the Analysis Window	14
2	4.3.	User can configure the mathematical transformation of the parameter SSC to	
		improve the separation of populations	15
2	4.4.	New options regarding printed report	15
5.	Ne	ew options regarding. Merger Module	
	51	New information and commands in the merger screen	17
	5 2	News regarding merger process controls	17
•			1/
6.	Ne	w options regarding File Export	19

1. <u>The FL-based merger option has been deleted.</u>

In the version 1.1 there is a unique merger option since it is possible to obtain the information that in the previous version was obtained with the FL-based merger option, from the new "Parameter Band Dot Plot" diagram.

1. Select the population of interest in the file resulting from the merger.



2. Open a new diagram type "Parameter Band Dot Plot" selecting as number of dimensions the number of parameters that you want to visualize from the merged file.

NEW DIAGRAM	\mathbf{X}							
Histogram 2D Dotplot DotPlot Parameter Band DotPlot 3D Graph APS Graph Density 3D DotPlot MultiDimensional								
Number of Dimensions:								
OK Cancel								

3. By right-clicking in the middle of the new Parameter Band Dot Plot, go to "Graph Configuration" menu and select "Axes Configuration". Through this command the user may select the parameters to visualize on each band of the new diagram:

	😢 PB-DP	
Axes Configuration	en e	10 PE ACTIVADO:FL2-H Log
Band 1: IL-2 PE ACTIVADD:FL2-H V LOGABITHMIC V	T THE MARK MARK TRADEWARK AND AND ALL ALL ADDRESS AND ADDRESS AND ADDRESS ADDRE	IFN PE ACTIVADO:FL2-H Log
Band 2: IL-4 PE ACTIVADO:FL2-H 💙 LOGARITHMIC 💙		
Band 3: IL-10 PE ACTIVADO:FL2-H V LOGARITHMIC V	· · · · · · · · · · · · · · · · · · ·	AFB FE ACTIVADO.FE2-A E0g
Band 5: TNFa PE ACTIVADO:FL2:H V LOGARITHMIC V		NFa PE ACTIVADO:FL2-H Log
Band 6: TNFb PE ACTIVADO:FL2:H V LOGARITHMIC V	a standard and a	-13 PE ACTIVADO:FL2-H Log
Band 7: GIFN PE ACTIVADU:FL2-H V LUGARITHMIC V Band 8: IP10 PE ACTIVAD0:FL2-H V LUGARITHMIC V		-10 PE ACTIVADO:FL2-H Log
FL4-Height:FL4-H IL-4 PE ACTIVADO:FL2-H IL-10 PE ACTIVADO:FL2-H ancel		-4 PE ACTIVADO:FL2-H Log
		-2 PE ACTIVADO:FL2-H Log

Parameter Band Dot Plot shows the legend of the analyzed parameters.

Each band represents the real expression of the parameter selected in the original file.

See detailed information about the different options related to the file merger process in chapter 4.2 of the User's Manual.

2. <u>New options and types of Graphs</u> have been created to show information.

See detailed information about types of graphic representations in chapter 5.3 of the User's Manual.

In the version 1.1 of InfinicytTM, <u>Parameter Band Histogram</u> and <u>Parameter Band Dot Plot</u> have been included as new type of graphics that allow to select as many different bands as parameters desired to be represented in the same graphic. The number of bands can be as many as the number of parameters in the data file. For example in a merged file, we could obtain a Parameter Band Histogram with "n" bands, one for each parameter.



The following Parameter Band Dot Plot shows on red the phenotypic expression of the plasma cells population previously selected.



Band 1:	FSC-H:FSC-Height	*	LINEAR	Y
Band 2:	SSC-H:SSC-Height	*	Exp-SSC Low	~
Band 3:	FL1-H:CD38 FITC	~	LOGARITHMIC	~
Band 4:	FL2-H:CD56 PE	~	LOGARITHMIC	~
Band 5:	FL3-H:CD19 PERCPCY5.5	*	LOGARITHMIC	~
Band 6:	FL4-H:CD45 APC	¥	LOGARITHMIC	Y

In some assays it will be very useful to visualize the data with <u>Density Representations</u> since they can help to better distinguish the different cellular populations. Version 1.1 of InfinicytTM allows visualizing the data as dot plots or density views in different types of graphs. To change to a density view, right-click on the diagram and select "Density"

Graph Configuration	Þ	1
Density	₽	✓ Contour
Assign Population	₽	Population
Undo/Redo	Þ	Configuration
Visibility	₽,	
Window Manipulation	۲	
Percentage of Events Displayed	۲	
Invert Axes		
Graph Info		

1. "Countour Density" considers all the events visualized to perform a density display.



2. "Population density" considers the population already classified to maintain the colours for the density representation. The density will be represented degrading the choosen colour.



3. "Configuration" allows the user to choose the scale of density to show in the graphic.



Using the new <u>Automatic Population Separator (APS) graph</u>, in the 1.1 version of InfinicytTM, the user can choose the parameters to use in order to obtain the best separation of populations. The software will take into account both the parameters selected and the rest of events independently. This option is available from the "Configure APS Parameter" window of the "Profile Configuration" menu and it is especially useful in merged files where we will select only the common parameters for the APS calculation.

Preferences	Configure APS Parameters		
Analysis Graphs Controls	Select all parameters		
• Statistics	FSC-A	LINEAR	~
 Stored Compensations Gates 	SSC-A	Exp-SSC Low	~
FCS file	FITC-A:CD22	LOGICAL	~
Current Compensation Negative Visibility Configure APS Parameters Properties	PE-A:CD23	LOGICAL	~
	PerCP-Cy5-5-A:CD19	LOGICAL	~
	PE-Cy7-A	LOGICAL	×
	APC-A:CD20	LOGICAL	~
	APC-Cy7-A	LOGICAL	¥
	Pacific Blue-A	LOGICAL	~
	AmCyan-A	LOGICAL	Ŷ

The aim of APS graph is for the user to be able to distinguish in these representations as many different populations as possible (especially useful when working with a high number of parameters) for a further characterization of these populations in the conventional 2D graphics. New in this version, whenever you gate and classify a population, the APS graph will calculate again the view taking into account selected populations. This way, it is possible to distinguish even more subpopulations among already-classified events.



3. It is possible to create and save <u>Analysis Profiles</u> for further uses.

In order to facilitate and standardize the analysis of similar samples, version 1.1 of InfinicytTM allows to configure, to save and to load for further uses different analysis profiles.



A profile is all the information about the software configuration. It is possible to create different profiles for different diseases, so the user can always work with the same graphs, population tree, controls, statistics, compensations, gates, report, etc with the same type of study.

It is recommended to see detailed information about Analysis Profiles in chapter 5.1.4 of the User's Manual.



3.1. Statistics

3.1.1. It is possible to use results calculation formula.

The result of the formula can be located in the populations tree, if the created formula relates different previously established populations (for example CD4/CD8 ratio), or in the columns relating to the statistics in the Analysis Window.

3.1.1.a) In order to place the result of the formula as a new population in the Population Tree, right-click on the first label of "Events" of the Population Tree and select "Add Function".



A new window will open where the user can write the name of the function desired and choose the formula.



Once the formula is finished, click OK and the new function will not be a new column with values for each population, but a unique value shown as a new population of the Population Tree.

🤨 ANA	LYSIS Normal peripheral blood sample						
<u>F</u> ile Sta	itistics <u>S</u> election Pr <u>o</u> file <u>W</u> indows <u>H</u> elp						
VIS.	Population	Events	T otal %	Partial %	Visibility %	Comments	
	Events	60000					Í
	Other Events	6529	10.88%	10.88%	0.00%		
		53471	89.12%	89.12%	100.00%		
	 Other Leucocytes 	0	0.00%	0.00%	0.00%		
	<mark>Eosinophils</mark>	1838	3.06%	3.44%	3.44%		
	Neutrophils	31582	52.64%	59.06%	59.06%		
	Monocytes	3074	5.12%	5.75%	5.75%		
	É⊷ <mark>_Lymphocytes</mark>	16977	28.30%	31.75%	31.75%		
	Other Lymphocytes	0	0.00%	0.00%	0.00%		
	T lymphocytes	12438	20.73%	73.26%	23.26%		
	f CD4/CD8 Ratio: 0.60	\triangleright					
	 Other T lymphocytes 	584	0.97%	4.70%	1.09%		
	CD4 T lymphocytes	4432	7.39%	35.63%	8.29%		
	CD8 T lymphocytes	7422	12.37%	59.67%	13.88%		
V	B Lymphocytes	2586	4.31%	15.23%	4.84%		
	🔲 🔷 Other B Lymphocytes	0	0.00%	0.00%	0.00%		
	kappa B lymphocytes	989	1.65%	38.24%	1.85%		
	Lambda B lymphocytes	1597	2.66%	61.76%	2.99%		
	NK cells	1953	3.26%	11.50%	3.65%		
			Gate:	0 events	Total %: 0.0	00 Visib	oility %: 0.00

3.1.1.b) In order to show the result of the formula in the columns relating to the statistics in the Analysis Window, go to Profile Configuration menu click Statistics and activate the "New Function" option. A new window will open, where the user can write the name of the function desired and choose the formula.

Preferences Analysis Graphs Controls	Show Statistics Mean V SD V CV V Median Mode New function V
Statistics Stored Compensations Gates FCS file Current Compensation	Name CD4/CD8 Rate
Negative Visibility Configure APS Parameters Properties	Value 0.00 V
	Number Data Population Data Subfunction Parameter Data Summation Value 0.00 \$
	Format Decimal separator Character separator in CSV file; Define decimal places 2
	Absolute Counts Events / µl of the reference population (beads 1,000



Once the formula is finished, click OK and the new function will be a new statistic column with values for each population and each parameter.

😢 AN	ANALYSIS Normal peripheral blood sample											
File SI	atistics Selection Profile Windows Help											
VIS.	Population	Events	Total %	Partial %		■ CE	FL4-H 03 APC	\bigcirc				
	Events	60000			Mean	SD	CV	CD4/CD8	1			
	 Other Events 	6529	10.88%	10.88%	659.50	1447.93	219.55	0.60	1			
	Eucocytes	53471	89.12%	89.12%	532.29	1010.19	189.78	0.60	10			
	 Other Leucocytes 	0	0.00%	0.00%	0.00	0.00	0.00	0.60				
	Eosinophils	1838	3.06%	3.44%	131.25	617.14	470.21	0.60				
	Neutrophils	31582	52.64%	59.06%	26.95	172.47	640.06	0.60	5			
	Monocytes	3074	5.12%	5.75%	272.82	831.14	304.64	0.60				
	E Lymphocytes	16977	28.30%	31.75%	1562.78	1194.11	76.4	0.60	3			
	Other Lymphocytes	0	0.00%	0.00%	0.00	0.00	0.00	0.60				
	🖨 T lymphocytes	12438	20.73%	73.26%	2127.68	865.95	40.70	0.60	i			
	 Other T lymphacytes 	584	0.97%	4.70%	3926.38	1336.96	34.05	0.60				
	CD4 T lymphocytes	4432	7.39%	35.63%	2122.35	769.34	36.25	0.60				
	CD8 T lymphocytes	7422	12.37%	59.67%	1989.34	699.17	35.15	0.60	1			
	B Lymphacytes	2586	4.31%	15.23%	16.31	113.44	695.37	0.60				
	 Other B Lymphocytes 	0	0.00%	0.00%	0.00	0.00	0.00	0.60				
	kappa B lymphocytes	989	1.65%	38.24%	14.70	138.40	941.79	0.60				
	Lambda B lymphocytes	1597	2.66%	61.76%	17.32	94.73	547.06	0.60	1			
V	NK cells	1953	3.26%	11.50%	12.82	28.95	225.89	0.60				
								\bigcirc				
<									>			
					Gate:	0 events	Total %: 0.00	Visibility %: (0.00			
					,				_			

3.1.2. Absolute counts from the population of reference selected by the user.

To calculate absolute counts there is no need to create the formula. The user can choose which population will be taken as reference for the calculation in the "Absolute Counts" option of the Statistics Window of the "Profile Configuration" menu.

Once the population which will be taken as reference for the calculations the number of event per and microlitre of this population is selected, the software will automatically show a new column in the Analysis Window with the absolute counts (events /µl) of each population.



3.2. <u>Current compensation can be changed manually and saved</u> if the origin FCS file contains this information internally.</u>

The new cytometric file acquisition programs can save as internal data the compensation of fluorescences used. In those cases, with the version 1.1 of InfinicytTM the user can modify the compensation values manually so that the changes will be applied to the file. In addition, the user can save the values of the new compensation and apply them to other files. These options are available in the Current Compensation Window of the Profile Configuration menu.

A B C D E F G B FITC-A:D: Statistics A - 10.84 3.49 0.32 0.06 0.00 D.00 B PE-A:CD23 C Per-CP-Cy5 B FITC-A:D: B PE-A:CD23 C Per-CP-Cy5 B PE-Cy7-A C C 0.32 0.78 - 16.73 2.64 5.97 0.00 D Per-Cy7-A CS file D 0.26 11.16 7.65 - 0.15 2.72 0.00 F APC-A:CD2 D 0.26 11.16 7.65 - 0.15 2.72 0.00 G Partine Blue	2
Controls A B C D E F G A FITC-A:CD3 Statistics Statistics A - 10.84 3.49 0.32 0.06 0.00 0.00 C Per-A:CD23 C Per-CP-Cy5 Statistics B 1.30 - 35.31 2.64 0.15 0.00 D D Per-Cy-Cy5 S Per-Cy7-Cy5 S F A F A F A A A A A A A A A B 0.32 0.06 D.00 D D Per-Cy7-Cy5 C Per-Cy7-Cy5 F A F A F A B A <t< td=""><td>2</td></t<>	2
A - 10.84 3.49 0.32 0.06 0.00 0.00 C PerALD2S Stored Compensations B 1.30 - 35.31 2.64 0.15 0.00 D.00 D<	
B 1.30	-5-A:CD19
C 0.32 0.78 16.73 2.64 5.97 0.00 E APC-A:CD2 C 0.76 D 0.26 11.16 7.65 0.15 2.72 0.00 F APC-CY7-A C Current Compensation E 0.00 0.00 1.34 0.14 10.01 0.00 G Pacific Blue	C
D 0.26 11.16 7.65 - 0.15 2.72 0.00 F APC-Cy7-A Monsteine Visibility E 0.00 0.00 1.34 0.14 - 10.01 0.00 G Pacific Blue	D
Merchie Wichildy E 0.00 0.00 1.34 0.14 - 10.01 0.00 G Pacific Blue	
	A
Configure 405 Parameters F 0.07 0.01 0.21 (2.12 ♦ 0.28 - 0.02 F AIRLyan-A	
G 0.12 0.00 0.38 0.14 0.09 6.74	
H 17.89 1.31 0.29 0.02 0.00 0.00 19.23	<u> </u>
Save As.	

Selecting "Stored Compensations" option of the Profile Configuration menu, the user can check, delete and change the current file compensation for one of the compensation already stored.

rrererences	Storeu	Sureu Compensations									
Graphs	Name	Name									
Controls	COMPEN	ISATION 1	L								
 Statistics 	COMPER	COMPENSATION 2									
Stored Compensations	COMPEN	ISATION 3	3								
Gates											
FCS file											
Current Compensation											
Negative Visibility											
Configure APS Parameters	a page 1										
Properties											
			1		1		1	1			
		A	В	C	D	E	F	G	A FITC-A		
	A	-	10.84	3.49	0.32	0.06	0.00	0.00	D PC-A C PerCP-CV5-5-A		
	В	В	1.30	-	35.31	2.64	0.15	0.00	0.00	D PE-Cy7-A	
	C	0.32	0.78	-	16.73	2.64	5.97	0.00	E APC-A		
	D	0.26	11.16	7.65	-	0.15	2.72	0.00	F APC-Cy7-A		
	E	0.00	0.00	1.34	0.14	-	10.01	0.00	G Pacific Blue-A		
	F	0.07	0.01	0.21	2.12	10.28	-	0.02	H AmCyan-A		
	G	0.12	0.00	0.38	0.14	0.09	6.74	-			
	н	17.89	1.31	0.29	0.02	0.00	0.00	19.23			
	<			1111				>			

When the user selects a compensation already stored and clicks OK, then the selected compensation will be applied to the file, but whenever a new file is opened, it will be always opened with its internal compensation.

3.3. New option for <u>creating consecutive gates to identify a population and save the sequence to apply it to similar files</u>

1) Draw consecutive gates to identify a population:



2) Go to "Selection" menu of the Analysis Window and select "Save Gate As..."

🤨 ANALYSIS	Increased CD4+ cells PB sample			
File Statistics	Selection Profile Windows Help			
VIS.	Undo Draw Region Ctrl+Z Redo Ctrl+Y	Events	Total %	Partial %
E 🖃	Delete All Selections (Ctrl+0	60000		
-	Select File	5116	8.53%	8.53%
Image: Construction of the second	Save Gates As	54884	91.47%	91.47%
	Load Gates	0	0.00%	0.00%
	Edit Gates	300	0.50%	0.55%
	Neutrophils	27721	46.20%	50.51%
	Monocytes	2655	4.43%	4.84%
	É <mark>_Lymphocytes</mark>	24208	40.35%	44.11%
	 Other Lymphocytes 	0	0.00%	0.00%
	🖅 – Tlymphocytes	21137	35.23%	87.31%
	<mark>⊕</mark> B Lymphocytes	955	1.59%	3.94%
	NK cells	2116	3.53%	8.74%

3) Write a name for the identified population.

Save Gate	
Gate Name	CD34+ CELLS
	OK Cancel

- 4) When you want to apply the same gating sequence to a new file, go to the "Selection" menu of the Analysis Window and click "Load Gates". Choose the appropriate gate from the list and then go to the Analysis Window to classify the population in the Population Tree.
- 5) The pre-defined gates are part of the current profile configuration, therefore clicking the option "Gates" of the "Profile Configuration" menu or clicking the option "Edit Gates" of the "Selection" menu, the software will open a window with options to delete them or change the color of the events and the contour.

3.4. It is allowed to <u>configure negative visibility of the parameters with Logical scale</u>.

Negative visibility represents the percentage of negative data range in the Logical scale which is shown in the graph. It corresponds to the part of negative axes relative to the positive axes. Thus, whenever the FCS file has been obtained from an instrument with

electronic compensation, negative values that are usually close to the axes may be moved to improve the visualization of the positive events.

Go to "Profile Configuration" menu and click "Negative Visibility" to change the percentage shown for each parameter. Depending on this percentage, the negative part of the scale of the parameter will be more or less represented.

Preferences Analysis Graphs Controls	Negative Visibility Percentage	3 (\$	
Statistics	Parameter	%	
 Stored Compensations Cater 	FSC-A	16.88%	
FCS file	SSC-A	0.01%	
Current Compensation Negative Visibility	FITC-A:CD22	0.06%	
	PE-A:CD23	0.05%	
Configure APS Parameters Departies	PerCP-Cy5-5-A:CD19	0.1%	
 Properties 	PE-Cy7-A	0.38%	
	APC-A:CD20	0.1%	
	APC-Cy7-A	0.22%	
	Pacific Blue-A	0.05%	
	AmCyan-A	0.06%	
	Apply to all:	0.1%	
		· · · · · · · · · · · · · · · · · · ·	

If we increased the negative visibility percentage applied for a certain parameter, the negative part of the scale of this parameter will be less represented although it continues being visible.



3.5. New option to <u>check all the parameters measured of a data file</u> (Properties)

It is possible to go to the information about the parameters measured on a given data file, by clicking "Properties" of the Analysis Window "File" menu, by clicking the "Properties" option of the "Profile Configuration" menu or with Crtl+I.

This is especially useful when working with files acquired with more than 3-4 fluorescences or with merged data files because in this window we will always have available the matrix with the information for all the parameters and tubes.

Preferences	Parameters Info	
- Analysis	Parameters\File	123
Graphs	A.722	
Controls	CD22EITC.4	
Statistics	CD23 PE-A	
 Stored Lompensations 	CD19PerCP-Cu5.5.6	
	PE-Cu7-A	
FLS he	CD20-APC-A	
Current Compensation	FSCA	
Casting APC Recently	FMC7-FITC-A	
Configure AFS Farameters	CD24 PF-A	
Flopendes	CD34-APC-A	R
	CD103FITC-A	
	CD25:PE-A	
	CD11CAPC-A	

Preferences	Parameters Info	
Analysis Graphs	Parameters	1
 Controls Statistics 	FSC-A SSC-A	R
Stored Compensations	FITC-A:CD22	R
Gates	PE-A:CD23	R
FCS file	PerCP-Cy5-5-A:CD19	R
Current Compensation	PE-Cy7-A	R
Negative Visibility	APC-A:CD20	R
Configure APS Parameters	APC-Cy7-A	R
Properties	Pacific Blue-A	R
	AmCyan-A	R
	< III	>

4. New options make the <u>Manual Analysis of the data file</u> easier

4.1. The way in which the <u>user can save and find in his computer a document already</u> <u>analyzed (.cyt)</u> has been improved

InfinicytTM software allows saving the files analyzed, containing the information of classification, statistics, graphics, report, profile, etc. and store them in the hard disk. Select "Save Analyzed File" option of the "File" menu when you want to store the analysis of a file. The software will ask you the directory where you want to save it and the file will be saved as **.cyt**.



The new feature in version 1.1 of $Infinicyt^{TM}$ is that the analyzed file will be saved in the folder of the hard disk that the user selects with the extension **.cyt**. The original data file will remain the same in the corresponding folder, because the **.cyt** file does not affect the original one.

Whenever the user wants to open an already analyzed file, then the **.cyt** file must be opened from the option "Open Analyzed File" of the "File" menu.

4.2. Option to choose the columns to visualize in the Analysis Window.

By default, all columns of the Analysis Window, except the statistics column, will be visible. Nevertheless, at any time, right-clicking on upper columns, this information can be hidden or shown.

🤨 AN/	ALYSIS Normal peripheral blood sample				
File Sta	atistics Selection Profile Windows Help				
VIS.	Population	Events	T otal %	Partial %	now Events Column
	Events	60000		• SI	now Total Column
	 Other Events 	3117	5.20%	5.20	now Partial Column
	🖃 Leucocytes	56883	94.81%	94.8	now visible Column
	 Other Leucocytes 	1654	2.76%	2.9	now Commence Column
	Eosinophils	1552	2.59%	2.73%	Iow Statistics
	Neutrophils	30075	50.12%	52.87%	
	Monocytes	5743	9.57%	10.10%	
	Lymphocytes	17859	29.76%	31.40%	
	🛶 Other Lymphocytes	1	0.00%	0.01%	
	🖳 T lymphocytes	13044	21.74%	73.04%	
		2632	4.39%	14.74%	
	NK cells	2182	3.64%	12.22%	
		Gate: 0) events	Total %: 0.00) Visibility %: 0.00

4.3.The user can <u>configure the mathematical transformation of the parameter SSC</u> to improve the separation of populations.

In Graph Configuration option, the possibility of selecting the values of the mathematical transformation of the SSC-parameter to apply in the data file, in order to improve the separation of populations, as been added.

The "Configurable EXP-SSC Low" option gives the possibility to expand the SSCparameter with the values desired by the user. When clicking on it, the SSC plot will appear with two scrolls in the lower part of the graph to regulate the part of the scale to amplify.





4.4. New options regarding the Printed Report

- **4.4.1.** It is allowed to copy, cut and paste on Printed Report components of the software and external images or texts.
- **4.4.2.** It is possible to choose the size of paper from the "Properties" option of the report "File" menu.
- **4.4.3.** Undo (Cntrl+Z) / Redo (Cntrl+Y) options have been added from the Report "Edit" menu.
- **4.4.4.** "Data" option has been added as component of the printed report.



4.4.5. Different options have been added to configure the populations, the parameters and the statistics to be shown on the population component of the report. In the version

1.1 of InfinicytTM it is possible to select more statistic data, choose different files and show comments if this option is selected in the "columns" menu.

٩	🖁 Populatio	n				×
Ł	Select only po	pulatio	ns with events		~	>
	VIS.		Population			Chaliatian
			Events		1	Statistics
		÷	Other Events		1	All Files
			Eosinophils			
	V		Neutrophils			
		±	Monocytes			
		.	Lymphocytes			
		.	Other Lymphocytes			Columns
			T Cells			SD 🔼
			Other T Cells			CV
		.	CD4+ T lymph.			Median
		±	CD8+ T lymph.			New function
		±	CD4-/CD8- T lymph.			Comments 💌
		±	CD4+/CD8+ T lymph.			Display Options
		÷	B cells			
		.	Other B cells			
		±	Карра			Show Background
		÷	Lambda			
		±	Precursors			Data font
		±	NK cells			Header font
		÷	CD34+ precursors		~	
			OK Ap	ply Cancel		

5. New options regarding the Merger Module:

5.1 New information and commands in the Merger Screen.

	1- C:1001.PATIENT X TUBE 01 2- C:1002.PATIENT X TUBE 02 3- C:1003.PATIENT X TUBE 03 4- C:1004.PATIENT X TUBE 04	
Configure	Info Parameter/File	
Select Files	SSC-A	
	CD22 FITC-A:FITC-A	
Source Allo	CD23 PE-A:PE-A	
Save As	CD19 PerCP-Cy5-5-A:PerCP-Cy5-5-A	
	CD20 APC-A: APC-A	5 V
Main Menu	CD103 FITC-A:FITC-A	
	OD25 PE-A:PE-A	
	OD11C APC-A:APC-A	
	FMC7 FITC-A:FITC-A	
	CD24 PE-A:PE-A	
	CD34 APC-A:APC-A	
	LAMBDA FITC-A.FITC-A	

Version 1.1 of InfinicytTM includes the following available commands in the Merger Screen:

- Configure: File merge options allow manipulating the controls that will warn you in case your common parameters are not similar enough to perform a merge process.
- Select files: this command opens the directory where the user can select the files to be merged. Once the data files are selected, in the upper part of the dialogue box appears the list of the files to merge and in the lower part of the dialogue box appears the matrix (parameters/file).
 - You can select and move the files on this list with the ▲ and ▼ cursors, in order to obtain the best sequence to represent the tubes in the graphics.
 - In the matrix, the user can add or not each parameter, in order to avoid parameters that are not interesting because there is no marker on it. With this, we will optimize the memory of the new data file reducing the number of parameters to be merged.
- Save as... when the selection of the data files is done and the parameters to merge are the correct ones, the user can save this information. The software will ask you for a new name and directory for the new FCS file containing all the parameters merged.
- To go back to the main screen of the program, click the Main Menu command.

5.2 New features regarding merger process controls.

See chapter 4.2.4 of the User's Manual to obtain detailed information about Coherence Control to determine if the user is performing the merge process with the optimal files, analyzing in different aspects the similarity of expression in the common parameters.

The version 1.1 of $Infinicyt^{TM}$ includes the following new features about the merger process controls:

a. The coherence information window about the differences relative to the global coherence control (Merge Control) has been simplified. When the program detects significant

differences between the distributions of the events of the different files in the common parameters space, it shows a warning message. The warning will indicate the percentage of difference (value 0 means that there are no differences between the files and value 100 indicates that the distribution of events relative to common parameters is completely different), and will also list, by significance, the parameters and files which make the differences.

Maximum difference	between files		
	26 % difference	between some files	
Files and parameter	s with significant differences	5	
	FILE 6:	FSC-A	
	FILE 5:	FSC-A	

b. In order to avoid errors in analysis, the option of correcting the parameters when a difference of no coherence between populations (Population Control) is detected has been deleted

6. New options regarding File Export

- **6.1. The available export file formats have been reduced to Standard FCS or TXT**. When FCS format is chosen, the data will be saved in the form of an FCS file, readable by any cytometric analysis software, but not by common worksheet or database programs. When it is chosen TXT format, the values will be saved in ASCII format as a data base, readable by worksheets or word-processing programs
- **6.2.** Some options on the advanced functions of the File export to avoid extra memory requirements resulting from unnecessary information have been added:
 - The user can export only the <u>parameters</u> which are necessary for the analysis. Initially, all the parameters are selected for export. If the user wants to eliminate some of these parameters, unclick the parameter to eliminate it from the export.
 - <u>Events</u>: if the user is trying to export a merged file, InfinicytTM will give the option to either export events proportional to each file maintaining the file number separation, or mixed in a file, without separation between original files.

	Basic Advant	3,2,1	F	estore Values	
	PnE format So	ource 🗸			
	Export	Parameter		PnB	PnB
	V	FSC-Height :FSC-H		1024	16 🗖
	V	SSC-Height:SSC-H		1024	16
	V	CD103 FITC:FL1-H		1024	16
Carriella		CD25 PE:FL2-H		1024	16
Save As		CD19 PERCPCY5.5:FL3-H		1024	16
		CD11C APC:FL4-H		1024	16
Main Menu	V	FMC7 FITC:FL1-H		1024	16
	V	CD24 PE:FL2-H		1024	16
	V .	CD34 APC:FL4-H		1024	16
	V .	CD43 FITC:FL1-H		1024	16
	V	CD79B PE:FL2-H		1024	16 🗸
	Events 2	50000 Pro	portional to each File 🔜		
		Pro	portional to each File		
		kdis.	ed in a file		

6.3.From the Main Screen the File Export option will consider the export of a whole file. If the user is interested on **exporting only one or more populations**, then, the "Export" option must be opened from the "File menu" of the Analysis Window. If the box relating to "Only visible populations" is enabled, only the visible populations in the Analysis Window will be considered for exporting. If it is not, all the events contained in the file will be exported.

	Basic Advanced			
	Only visible population:	s 🗹		
	Simulated Time			
Save As	Exp-SSC Low		SSC-A	*
Main Menu	1			
	Export File Format			
		Standard FC	S	
		Ο ΤΧΤ		