




MICROARRAY DESIGN USING CLONETRACKER DBTM

A Supplement to the CloneTracker DB™





Every effort has been made to ensure the accuracy of this document. However, BioDiscovery makes no warranties with respect to this documentation and disclaims any implied warranties or merchantability and fitness for a particular purpose. BioDiscovery shall not be liable for any errors or for any accidental or consequential damages in connection with the furnishing, performance or use of this document or the examples herein. The information in this document is subject to change.

Trademarks

CloneTracker DB™, ImaGene™, GeneSight™, and GeneDirector™ are trademarks of BioDiscovery, Inc. Other product names mentioned in this document may be trademarks or registered trademarks of their respective companies and the sole property of their respective manufacturers.

Copyright 2002 BioDiscovery, Inc. All Rights Reserved.


Written at BioDiscovery, Inc., Marina Del Rey, CA, 90292 Printed in the United States of America



4640 Admiralty Way, Suite 710, Marina del Rey, CA 92092 Phone 310.306.9310 Fax 310.306.9101 E-mail: support@biodiscovery.com

CONTENTS

CHAPTER 1 - Introduction to CloneTracker DB and its role in Array Design	
Microarray Experimental Design to Assist Normalization of Data	.6
CHAPTER 2 - CloneTracker DB terminology	.8
Experimental Terms	. 9
Array Designs	. 9
Arrays	10
Clones	10
Plates	10
Array Design Terms	10
Field	10
Layout Mode	11
Sector	11
Pin Configuration	11
Pin Spacing	11
Spot Staggering	12
CHAPTER 3 - Moving from CloneTracker DB v1.x to v2	13
Database Converter	14
User Interface Changes	15
Naming Convention Changes	16
CHAPTER 4 - Design examples	19
Example 1: Designing a cDNA macro Array— working with a nylon membrane	20
Example 2: Designing a a High density Array— working with glass slides	26
Example 3: Plate-to-plate-to-slide mapping	29
Example 4: Designing an array with replicated / Duplicated spots	32
Example 5: Exercise for practice	36
APPENDIX 1: FAQ - Frequently Asked Questions	38
APPENDIX 2: Solution to Exercise in Chapter 5	39
References	40
New Features	41



Introduction

Microarray projects involve acquisition and validation of large data sets. Each study typically requires iterations of series of processes, which start with experimental design and array fabrication, proceed with array scanning, image analysis, and finally gene expression data analysis.

The initial data resources required to design and fabricate gene expression microarrays that need to be documented include cDNA sequence data, cDNA clones information, experimental procedures, physical location of materials, and others. This information needs to be integrated within the data flow during array design, analysis, and interpretation of results (Figure 1-1).




Fig. 1-1

The software tool CloneTracker DB from BioDiscovery, Inc. (Marina Del Rey, CA; www.biodiscovery.com) has been designed to make management of array fabrication data easy and rewarding.

The CloneTracker DB software offers:

Laboratory information management. A management system that provides access interface to the underlying fabrication database. It helps the user monitor the overall fabrication process by allowing queries of fabrication data, updating data step by step according to the progress, or tracing a particular spot on a particular array back to the originating plate and well position.

Array design. An easy-to-use visual interface that allows the user to design or adjust the array layout based on the type of arrayer used. It can also export the settings for the arrayer or even generate complete programs to control the arrayer robot directly.

Utilities for complex data manipulation. The information flow in the fabrication process involves several transformations that cannot be handled using simple database operations. Examples: (i) transformation of clone data from a set of plates (plate data) into spot data on the array (array data); (ii) generation of fabrication data in the format recognized by the image analysis system ImaGene and the data-mining tool GeneSight; (iii) tracking clone transfer-processes from 96 well to 384 well plates to a DNA array, and back.

The primary goal of this document is to illustrate uses of CloneTracker DB with examples while explaining to existing users new features and terminology. Readers should review the CloneTracker DB User's Manual when using this software.


CloneTracker DB is useful not only as a virtual design tool to facilitate visualization before printing of the array and organizing data, but also to help the actual experimental design.

Microarray Experimental Design to Assist Normalization of Data

A designed experiment is a test or series of tests, in which purposeful changes are made to the input variables of a process or system so that we may observe and identify reasons for changes in the output response (Montgomery, 1991).

The basics principles of experimental design are replication, randomization, and blocking. Approaches for experimental design in microarray experiments are as follows:

- Use more than one cDNA spot for normalization purposes spot replication. CloneTracker DB offers an easy approach for creating duplicates, triplicates, and more spot replications. This information is submitted to the data analysis tool GeneSight (BioDiscovery) for data normalization procedures.
- Replicate the array on the same slide array replication. This is a blocking strategy (performing an experiment on a homogenous material, i.e. the same slide), which can be achieved by the click of the mouse with CloneTracker DB.
 - Use control spots. CloneTracker DB keeps track and can color-code control spots.
- Replicate the experiment on different slides experiment replication. CloneTracker allows you to use an array design that is automatically applied to many slides.
- Randomize the spot replicates on the array. You may use CloneTracker Replicates Editor to design your array.




Due to the flexibly associated with various parts of the arraying processes, CloneTracker DB has developed a unique terminology covering various design aspects. The terms listed below reoccur throughout the usage of CloneTracker DB. As a result, users are encouraged to have a working knowledge of their meaning. Most terms covered below are found throughout microarray design and clone management. For additional information, please consult the *CloneTracker DB v.2 User's Manual*.

EXPERIMENTAL TERMS

Array Designs

Array Designs are templates of spotted arrays, describing how a particular array is spotted (Figure 2-1) and from which plates the clones are taken (Figure 2-2). Array Designs are used to create multiple Arrays of the same template, or design. Some of the parameters used to describe an Array Design include the array dimensions (whether for a slide or a membrane), pin and spot configurations, pin and spot spacing, and sectors.



Example: If a particular set of genes (we denote any length of DNA whether an EST or cdna as a gene for simplicity) needs to be printed on hundreds of slides by a core facility then the set up applying to each slide is the array design. This blue print contains all the information about the dimension of the slide, how many pins are used for printing, what plates with DNA probes were utilized, etc.

Arrays

Arrays are particular instances of Array Designs corresponding to physical arrays. The plates used in the Array Design are the plates used in the Array. Additional information such as barcode information can be applied to each Array.

Clones

Clones are the biological materials contained in Plate wells. Therefore a clone may mean a PCR product, cdna, or other material to be printed. Clones may be imported into the database along with Plate Importation, or they may be added to existing plates that are not populated with Clone data.


Plates

Plates are small plastic containers consisting of 96, 384, and 1536 wells. These wells contain the material to be spotted onto the microarray slide. Within CloneTracker DB, the Plate Type indicates whether the plate is 96, 384, or 1536 well plate.

ARRAY DESIGN TERMS

Field

A field is the largest design element within a slide (Figure 2-3). A field typically consists of the arraying done by a single print head on the slide. For example, if the arrayer has a print head with 8 pins in a 2 x 4 configuration, the region of the slide containing the resulting printing is a field. For most experiments and arrayers, users will only have a single field. However, CloneTracker DB 2 fully supports the creation of multiple fields if necessary.



Common uses for multiple fields include:

Fig. 2-3

- Generating slides which contain arrays of different geometries: for example, some microarray core facilities procedures include a set of designated genes on each slide produced in the facility for calibration purposes. These genes are not within the field with the genes under investigation. Such sets of genes will be located in a separate field. CloneTracker DB 2 offers the flexibility and ease-to-use to create this field.
- Generating slides with duplicated arrays next to each other. This will be needed if you decide to replicate the whole array on the same slide as discussed in Lee et al. (2000).

Layout Mode

Regular - In Regular Mode (Figure 2-4) the printing is done in such a way that each pin prints the required number of spots from a single plate as to use all samples from all wells. See the examples in following chapters.



Fig. 2-4

Compact - In Compact Mode, users are free to design the rows and columns to be printed irregardless of the relationship to the plate. This allows for more complicated array design.

Sector

A sector is a region on the array containing spots printed by one pin. The source for a sector could be one or several plates depending on the Layout Mode. (Fig. 2-5)

Example: when printing in Regular Mode, a sector is the printing done by a single pin from a single plate. In this case, each sector is displayed as a different color within the software.

Please note that when the Layout Mode is switched to Compact, the notion of a sector increases in complexity. Under Compact Mode, a sector defines the rows and columns of printing; however, this may or may not be equal to the printing done from a single plate. Sector Size refers to the number of rows and columns of clones/spots placed on the slide by a single pin. In Regular Mode this is **FIXED** to a size which allows the print head to dip to every well of the plate. In Compact Mode, the sector size can be modified to virtually any values by the user. Due to the gap between pins on the pin head, researchers can spot several Sectors between each pin. CloneTracker DB allows this to be done by specifying the Number of Sector Rows and Number of Sector Columns.




Fig. 2-5

Pin Configuration

The arrangement and number of individual pins on the print head. CloneTracker allows you to use either hand-held devices with 96 pins or any arrayer with any combination of pins.

Pin Spacing

Pin Spacing determines whether the pins on the print head dip into wells next to one another or whether the pins skip a well in the middle. Typically, users should select Consecutive wells.

CloneTracker DB Terminology

Spot Staggering

Spot Staggering allows researchers to fit an increased number of spots within an array. Spot Staggering prints spots in a hexagonal structure minimizing the space between spots. Left staggered means the initial row is spotted to the left of the row below it. Right staggered means the opposite where the top row is position slightly to the right of the row below it. (Fig. 2-6)





Fig. 2-6



DATABASE CONVERTER

CloneTracker DB 2 includes a handy database conversion utility designed to transition data from **CloneTracker version 1.4.4 and higher** to CloneTracker DB version 2. User with older versions of CloneTracker DB will need to export plates from version 1 then recreate the slides within CloneTracker DB 2. Through use of the DB Converter, all information within the database is automatically transferred to the new MySQL database.

To convert the database, perform the following:

- 1. From within the Start Menu, click DB Converter within the CloneTracker DB 2 menu.
- 2. Read the initial message. click OK to continue. (Fig. 3-1)







Fig. 3-1

Fig. 3-2

- 3. The next window is the Database Transfer Window. For most users, the correct parameters will be provided. However, **the user must enter the source password**. By default, this password is: **sql**. There is no need to enter a Target password.
 - 4. Click **OK** to continue.
- 5. The information is transferred and Log window is then displayed. If the same data exists within the new database, the user is prompted what to do with the data. (Fig.3-2)
- 6. Upon completion, **"Logout... Succeeded"** is displayed in the Log window. The data is now available within CloneTracker DB 2.

USER INTERFACE CHANGES

CloneTracker DB 2 has included the best aspects of version 1 and applied several new options for the main user interface, which were required by a world-wide group of scientists employing CloneTracker DB in their research. The user interface in version 2, as it was in 1.x, is primarily geared to virtual array design. Within this virtual design, the key elements, slide, array geometry, clones and mapping information, are placed predominantly before the user. Please see diagram below for similarities.




CloneTracker DB 2 also adds several new features as outlined below. (Fig. 3-5)

When zoomed in, spotting orders become visible.

All Arraying Parameters conveniently displayed on main screen.

New Toolbar for quick access to important features.



Ability to edit individual "Replicate Patterns"


Fig. 3-5

Ability to control the mapping information that is displayed.

NAMING CONVENTION CHANGES

CloneTracker DB 2 uses a unique naming convention which allows for the description of highly complicated arraying design. Note, although there have been terminology changes between previous versions of CloneTracker DB and version 2, the same designs can still be created.

CloneTracker 1 used the notions of spot layout, subgrid and super grid to describe all the design elements of the slide. CloneTracker DB 2 also has a set of terms, sector, sub-grid and field, which perform the same task. The following information explains the relationship between naming conventions.






Fig. 3-6

Fig. 3-7

The smallest level of design in v1 is the spot layout. (Fig 3-6) In v2, this is the **Sector**. (Fig. 3-7) When arraying in Regular Mode, spot layout and sector represent the same thing.

Note: When performing more complicated arraying this is not true; however, such arraying is beyond the scope of this document.

In v1, all the spot layouts generated by a single pin are grouped into a subgrid. (Fig. 3-8) In v2, the sectors are grouped into a subgrid, (Fig. 3-9) but are loosely called simply Rows and Columns of Sectors. The subgrid term is not used much in v2 except to indicate to which direction printing should continue upon filling the allocated sector rows and columns.




Fig. 3-8


Fig. 3-9

The top level of design in v1 is the supergrid. The supergrid was simply all the subgrids printed within the slide. However, the problem of a supergrid is that it only allows for a single supergrid design per slide. For example, this prevents users from generating a slide where extra marker genes are placed outside of the supergrid, as is displayed below.

It also prevents different arraying designs, or shapes, to occur on a single slide. V2 addresses these issues by introducing the concept of a field.

Moving from CloneTracker DB version 1.x to version 2

Typically, a single field represents a the entire supergrid structure on the slide. For example, if the slide has a supergrid structure of 1x3, (Fig 3-10) then this design would be contained within a single field in v2. (Fig. 3-11) The two versions differ because v2 allows users to create multiple or different fields for additional designs within the same slide. (Fig. 3-12)








Fig. 3-11





Fig. 3-12



EXAMPLE 1

Designing a cDNA Macro Array: working with a nylon membrane

OBJECTIVE: Print a cDNA array on a nylon membrane with a hand held device, which has 96 pins.

Array Design

Specify the array design name and the name of the user under the Array Design tab. Let us call this array design "Membrane Mouse". The dimensions of this membrane are 20 cm x 20 cm. Therefore you may enter the corresponding dimensions in microns, i.e. 200,000 – (see frequently asked questions). You many choose margins depending on your experiment. For this example enter 10,000 for all margins. After all the information is specified, press the *Apply* button (Fig. 4.1.1)




Fig. 4.1.1

FIELD DESIGN

Click on the **Field Tab**. Specify the name of the field under the Field tab. For example, type in "Device 96 pins".

Step 1 - Plates

Initially, only the **Plates** tab is enabled (Figure 4.1.2). Depending on your microarray experiment, choose either 96 (8x12), 384 (16x24) or 1536 (32x48) well plates to be used for the final stage of printing by selecting from the **Plate Type** drop-down menu. For this experiment we need to choose 96 (8x12).




Fig. 4.1.2

Click on the Add Plate icon to bring up the Edit Plate List window (figure 4.1.3). Choose your choice of plate(s) from the Available Plates panel. For this membrane-printing project, select "8x12 (96-well plates)" and choose Plates 501, 502 and 503.

Click on the right arrow icon to bring your chosen plate(s) to the Selected Plates panel. Hit the OK button.




Fig. 4.1.3

If initially there are no plates listed or you would like to import plates then click on the Plates icon in the toolbar. Alternatively, select Plates under View menu (Figure 4.1.4).

Click on the "Import Plates From Text" icon on the right icon bar selection (Figure 4.1.5). A "Plates Importer" window shows up. Click on Browse and find the folder Samples (located in your CloneTracker folder). Choose 96 plates.txt file, which contains clone information for four 96 well plates called 501, 502, and 503, i.e. they are concatenated. A file can also contain information for one plate only. You may want to view or edit this file by using Excel. Once editing within Excel is complete, remember to save the file as a tab delimited text file for compatibility with CloneTracker.




Fig. 4.1.4




Fig. 4.1.5

In the **"Plates Importer"** click on the selection buttons (squares) for the columns you wish to import (Figure 4.1.6). In this case, select all columns. The selected columns will be colored in red. Click on the first row of plate 501 and select it. Scroll down to the last row of plate 503. Hold the **SHIFT** key on your keyboard to select all rows including the last one (row **h**, column **12** of plate 503). All the columns become highlighted. Click the Down Arrow button at the bottom of the "Plates Importer" window. The selected plates appear in the lower half of the Importer window (Figure 4.1.7). You may enter some remarks about each plate in the column termed "Info". Click OK. Now you may go to the Field Tab and click on the "Add Plates" icon. Select plate 501, 502, and 503. Click the right arrow. Now plate 501, 502, and 503 are added to the **Selected Plates** panel. And then click on the **0K** button.





Fig. 4.1.6






Fig. 4.1.7

Fig. 4.1.8

Step 2 - Pins

The Pins tab allows you to customize the desired pin configuration. You could modify the pin spacing by selecting either **Consecutive wells** or **Every other well**. Choose the 8x12 pin configuration for this experiment as shown if Figure 4.1.8.

Several well-dipping order patterns could be selected by clicking on the **Well-Dipping Order** drop-down menu. The most common dipping pattern used by most arraying machines is "Left to Right, Top-Down".

Step 3 - Spots

For simplicity purposes, we will not change the default values under the Spots tab. Researchers usually change the spot spacing and diameter to achieve special printing patterns. Depending on the type of arrayer machine, one is able to decrease the vertical and horizontal spot spacing to fit more spots into each sector area. Decreasing the spot diameter allows one to further decrease the spacing between each spot.

Step 4 - Sectors

A sector is most easily defined as the total number of spots printed by one pin from one plate in Regular mode (please see the following examples for further discussion).

The sector size specifies the layout of all those spots. In our case, we have 96 wells that are printed by 96 pins. There is only one spot in one sector. Therefore, the only sector size possible is 1x1 in Regular Mode. (96 wells/96 pins = 1). In fact, in this case one spot is equivalent to a sector.

The numbers of Sector Rows and Sector Columns define in what pattern the sector, as a geometrical entity, will be repeated. Enter the number of Sector Rows and Sector Columns and then change to another pair of numbers to see how the printing pattern alters. Let us choose 4x4.

Select a specific order for **Sector Ordering** from the pull down menu. The Sector Ordering refers to the order of clones that are spotted by one pin. If you choose Left-to-right, Top-down (and the well dipping order is Left to Right, Top-Down), the then the first spot (row one, column one) in the first sector will come from well a1 in plate 501, the second spot left to right (row one, column two) will come from well a1 in plate 502, and so on.

Let us consider the following design. Now your objective is to replicate the same clone (material from the same well) on the same array twice. Several research teams print the spots resulting from the same well adjacent to each other. For example, the material from well a1 plate 501 will be printed in row 1 column 1 and row 1, column 2. To achieve this we need to enter the plates in a different order. In the Plate importer enter the plate (in this case 501) twice as shown in Figure 4.1.9. Then do the same with plate 502, 503, and so on.

Additionally, some researchers like to leave blank spots for either background corrections or other reasons. Let us imagine that in row 4, column 4, for each pin you would like to leave a BLANK spot. The easiest way is to create a fake plate with all gene IDs = BLANK. Then import this plate in such an order after all the others so that the spots coming from this plate appear in row 4, column 4.




Fig. 4.1.9

Step 5 - General

For this project, we will stay with the Regular layout and keep the default settings under the General tab. Select the Regular layout mode ensures that the spots from one plate are spotted in a rectangular sector. The Compact layout mode will be explained in another project example.

Step 6 - Saving the Design

There are two ways to save the array design that you have created. By using either method, the Array Design can be easily recalled for future analysis or modification.

Saving to the CloneTracker Database

Click on File > Save on the menu bar (Figure 4.1.10). This will save the current Array Design to the CloneTracker database.

Saving to a Specific Location

Click on File > Export on the menu bar. The Save dialog box is displayed, where you can browse to the location of where you would like your array design exported to. Click the Save button and your slide set is saved in .xml format in the desired location.




Fig. 4.1.10

EXAMPLE 2

Designing a High-Density Array - working with glass slides

OBJECTIVE: Design a high-density array, which will be printed on a glass slide. This example will require 2x2 pins and 96 well plates.

Array Design

Specify the array design name and the name of the user under the **Array Design** tab. Let us call this array design "Slide Mouse". The dimensions of the slide by default are width $75,000\mu$ x height $25,000\mu$ (7.5 cm x 2.5 cm). The frosted part of the slide is $25,000\mu$. With these dimensions the slide width is horizontal to the viewer and the frosted side is to your right (Figure 4.2.1). If you wish to rotate the slide, then you can enter $75,000\mu$ for the height and $25,000\mu$ for the slide bottom margin. Let us leave the slide with the default settings.




Fig. 4.2.1

Field Design

Click on the Field Tab. Specify the name of the field under the Field tab. Enter 2x2 for Pin Configuration.

Step - 1 Selecting Plates

To select the plates you wish to use in CloneTracker, click on the Add Plates icon on the right, vertical toolbar in the Plates Panel. The Edit Plate List window is displayed. If there are no plates listed see Example 1 for details on importing plates.

Select plates 501, 502, and 503 and click on the right arrow to the right to import these plates. Select plate 501 again and import it so we can work with 4 plates in this experiment.

Step 2 - Pins

Select 2x2 pin configuration.

Step 3 - Sector

Click on the **Sector Size** tab. Your Sector Size pull down Menu gives you options, which are all combinations of 24, i.e. 1x24, 2x12, 3x8, 4x6, etc (Figure 4.2.2). The reason for this is that 24 is the number of times each pin will touch the slide from one plate during the printing. If this number is denoted by N then $N = [W/P] \times SR$, where W is the number of wells in the plate, P is the number of pins, and SR is the number of duplicate spots. In this case 96 wells / 4 pins (only one gene replicate) yields 24.

Let us choose sector size 3x8 (3 rows x 8 columns)

Enter the number of Sector rows = 3 and Sector columns = 2. Watch the array. Flip these numbers and see the design. In this way you will understand the architecture of the array and the flexibility of CloneTracker DB to create your designs.




Fig. 4.2.2

To understand the design of the finished array and to conceptualize what will be printed by one pin, we can compute a subgrid size. The subgrid is simply the total numbers of rows and columns of spots printed by a pin from several plates. The actual calculation is:

(Number of sector rows) x (number of rows in each sector) = total number of rows of each subgrid.

(Number of sector columns) x (number of columns in each sector) = total number of columns of each subgrid.

So in the above example, we have 3x3=9 rows and 2x8=16 columns in each subgrid.

Step 4 - Spot Information

To create a high-density array, we will need to change the information regarding the size and position of spots. Click on the Spots tab to access this information. First, let us change the spot staggering to **Staggered Right**. Upon selection, the array design to the right will be updated to reflect the new

design. To further increase the number of possible spots within the array, let us decrease the **Spot Space Horizontal** and **Spot Space Vertical**. The default values of 250 should be changed to 225. As a result, by manipulation the parameters under the **Spot** tab, users can design arrays for machines capable of high-density printing.

Exportation of a gene ID file

After generating an array design, you can export the gene ID file to your image analysis software. ImaGene from BioDiscovery automatically imports gene ID files generated by CloneTracker.

Select "Export Mapping" under the File menu in CloneTracker. Choose ImaGene format (Figure 4.2.3). You have the option to assign the name "blank" to all empty spots. Those are the spots where nothing is printed. This is particularly useful if you want to carry out your background corrections based on such spots for the data analysis. Alternatively, you can choose not to export any empty spots.

With the Browse button select the directory and folder where you would like to save the file. Assign a name to the file and click the save button.




Fig. 4.2.3

EXAMPLE 3

Plate-to-plate-to-slide mapping

OBJECTIVE: It is a common practice to perform a PCR reaction in 96 well plates and then transfer some volume of the PCR reaction from each well to a well in a 394 well plate. Subsequently, the printing of the microarray will be carried out by using the 384 well plate (Figure 4.3.1). This protocol allows for printing of high density arrays and storing some of the initial PCR material for future use. Therefore, we need to map the clone IDs from a 96 well plate to a 384 plate and then to the glass slide array.




Fig. 4.3.1

Alternatively, some protocols call for transfer of biological material from 384 well plates to 96 well plates, which are subsequently used as a source for printing the arrays. CloneTracker has a module, which assists in mapping clones in such protocols.

A. MERGING FOUR PLATES INTO ONE

Under the "View" pull down menu, select "Plates". Click the "Merge Four Plates into One" icon on the right, vertical toolbar of the Plates Panel to display the Plate Merger dialog box. Select the four plates that you want to merge from the "All Available Plates" list. Click the left arrow button to move these plates to the Four Plates list. For example, if you have the plates 501, 502, and 503 then transfer them to the left. Now select plate 501 again and transfer it to the left (Figure 4.3.2).




Fig. 4.3.2

Click the Merge button to merge the plate tables Figure 4.3.3). Specify a name for the merged plate under Plate ID (as an example, Mouse Scale Up) and click the OK button to save your changes. You should always enter notes about the created plate in the "Description" field. These notes prove always useful in future data analysis.

Return to the CloneTracker Plates Panel.




Fig. 4.3.3

Importing Merged Plate

In the Array Design Panel of CloneTracker, click on the Field tab. Under the Plates tab, select the type of plate you wish to import. In this project, we want to use the 384-well plate we had just created from the four 96-well plates, so we select 16x24 (384-well plates) as the Plate Type. Click on the Edit Plate List icon to bring up the Edit Plate List window. You will see that the newly created "Mouse Scale Up" plate is there. Choose your choice of plate(s) from the Available Plates list, click on the right arrow icon to bring your chosen plate(s) to the Selected Plates list. Hit the OK button.

B. SPLITTING ONE 384 WELL PLATE INTO FOUR 96 WELL PLATES

Some of your projects may require performing PC in a 394 well PCR plate and subsequently splitting the materials into 96 well plates (Figure 4.3.4). Under the "View" pull down menu select "Plates". Click the "Split One Plate into Four" icon on the right, vertical toolbar of Plates Panel to display the Split plates dialog box. Select the 384 button for plates to be split. Select the "Mouse Scale Up" plate. Click the right pointing arrow on the right. The 384 well plate is automatically split into four 96 well plates, which have the following names: "Mouse Scale Up-0", "Mouse Scale Up-1", "Mouse Scale Up-2", and "Mouse Scale Up-3".

If you decide to design an array with 96 well plates and repeat the plate import procedure as outlined in example 1, you will see that the newly created 96 well plates are available for mapping to the glass slide array.




Fig. 4.3.4

EXAMPLE 4

Designing an Array with Replicated / Duplicated Spots

As discussed in Chapter 1, replication is one of the basic principles of experimental design. Many researchers prefer to have replicates of their genes on the same array. There are a two ways to achieve this in CloneTracker.


The easiest way to duplicate spots is to copy the entire field to a different location on the slide. In the Array Design Panel, start off by loading some plate files. (See Example 1: Step 1 - Plates)

Copy Field

- 1. Select Copy Field from the Fields menu.
- 2. Click the Field Mode button on the toolbar.



3. Click on the new field and drag it into the desired position.






Fig. 4.4.1 – Duplicated Fields

Some researchers prefer to have their duplicate clones positioned differently in each field. Changing the Spot Ordering or Sector Ordering in the duplicated field easily does this.

Change Spot Ordering

- 1. Select the duplicated field tab (typically Field 2).
- 2. Select the Spot tab.
- 3. Click the drop-down triangle under Spot Ordering to change the ordering to something different than the original field.

Change Sector Ordering


- 1. Select the duplicated field tab (typically Field 2).
- 2. Select the Sector tab.
- 3. Click the drop-down triangle under Sector Ordering to change the ordering to something different than the original field.

The other way to duplicate spots is to make duplicates within the same sector as the original spots. This is called editing the Replicates Pattern. This tool is a little complicated, but it is very powerful.

Edit Replicates Pattern

- 1. Under the Field tab, click on the Spots sub-tab.
- 2. Click the Edit Replicates Pattern button to display the Replicates Editor dialog box.

In the left panel is an empty grid, and on the right panel is a grid full of spots. The left grid represents the layout of every sector in this field. A blank grid indicates a default layout.



Equals




Fig. 4.4.2 – A blank grid means a default layout

Indicate Spot Layout

- 1. On the right grid, click on the spot that you wish to replicate.
- 2. On the left grid, click on the location(s) that you wish to have the selected spot in.
- 3. Repeat steps 1 and 2 for all of your desired spot(s) to be replicated.

In the example below, we have indicated that we want the first clones to be duplicated three more times (marked in red). Any left grid boxes that you do not explicitly select are assumed to be empty.




Fig. 4.4.3 – Blank areas on the grid does NOT mean a default layout.

Fig 4.4.3 Blank areas on the grid does NOT mean a default layout.

Looking at the right grid above, you may wonder what happened to clones 1a2, 1a3, and 1a4. Since they are not on the grid, they will not be printed.

In order to have all the clones in the sector AND have duplicated spots, you must increase the size of the sector and add the missing clones back in (marked in green).

Change Sector Size

- 1. Enter the new number of rows to use in the Number of Rows field.
- 2. Enter the new number of columns to use in the Number of Columns field.
- 3. Press <ENTER> on your keyboard to apply the changes.




Fig. 4.4.4 – Expanded sector with missing clones added back in.

Now that we have our customized Replicates Pattern, we can apply the changes by clicking OK and exiting the Replicate Pattern Editor. Our new slide design now shows up on the Spotted Array window box.




Fig. 4.4.5 – New Slide Design after editing Replicates Pattern

After creating this customized Replicates Pattern, we may want to save our design and load it back in at a later time.

Export Replicates Pattern

- 1. Click the Export Replicates Pattern button to display the Save dialog box.
- 2. Select a destination folder for the replicates file.
- 3. Click the Save button to save the file to the designation location.

Import Replicates Pattern

- 1. Click the Import Replicates Pattern button to display the Open dialog box.
- 2. Locate and select a replicates file to import.
- 3. Click the Open button to display the selected file.

EXAMPLE 5

Exercise for practice

Using the examples outlined in the previous sections try to perform the following exercise by yourself. The proposed design configurations are given in an Appendix of this tutorial.

Let us try to design the following experiment. You will print your arrays with an arrayer equipped with 4x4 pins. You will use 384 well plates (the distance between spots is 240 um). The printing of the array sectors will be performed Right-to-Left, Top-to-Bottom. Genes will be replicated on this array by creating a duplicate meta-grid. The slide should be aligned vertically and there is no frosted area. You need to print the spots in sectors of 18x18 (Figure 4.5.1). Additionally, the objective is to have 12x4 such sectors.





Fig. 4.5.1

APPENDIX ES



APPENDIX I

FAQ - Frequently Asked Questions

- **Q:** Why does CloneTracker DB use microns instead of millimeters?
- A: Using microns allows CloneTracker DB to be highly accurate within it internal calculations. Also, as spot size and spacing is typically determined in microns, CloneTracker DB maintains unit consistency.
 - Q: Why does CloneTracker DB offer so many possible combination of arraying styles?
- **A:** CloneTracker DB 2 has been designed from the start as a comprehensive clone management system capable of performing virtually any task required of researchers. As a result, CloneTracker DB provides greater flexibility than any proprietary solution available.
 - Q: Is CloneTracker DB cross platform compatible?
- **A:** Currently, CloneTracker DB is only available for the windows platform. However, BioDiscovery is actively porting to Unix systems and the Mac OS X. Please contact BioDiscovery Support for the latest listing of supported OS.

APPENDIX II

Solution to Exercise in Chapter 5

Exercise 1 - Solution:

- 1. Open CloneTracker.
- 2. Click on the Array Design Tab.
- 3. Enter "Exercise 1" as the name for this design.
- 4. Set slide height to $75,000\mu$ and the width to $25,000\mu$ because we would like to have the slide positioned vertically.
 - 5. Assign 300µ to all margins. This will remove the frosted area.
 - 6. Click on the Field Tab.
 - 7. Select the Plates Tab.
 - 8. Select 384 as the plate type.
 - 9. Click the Add Plates button.
 - 10. Pick any plate from the Available Plate list.
 - 11. Click the Right Arrow button to add this plate to the Selected Plates list.
 - 12. Click OK. This will activate all other Field sub-tabs.
 - 13. Select the Pins Tab
 - 14. Pick 4x4 from the Pin Configuration list.
 - 15. Click on the sub-tab General.
- 16. Choose Compact layout mode. Because we need to print spots in a structure 18x18 = 324 we cannot use regular mode as explained earlier in one of the examples. If we use regular mode then we will have to choose a sector size equal to a combination of two numbers which have been multiplied and their result is 24 = 384 / (4x4).
 - 17. Enter 4,000 for Field Location X and 15,000 for Field location Y to position it.
 - 18. Click on the Sectors sub-tab.
 - 19. Type 18x18 in sector size (this size will NOT appear on the list).
 - 20. Select Left-to-Right, top-down sector ordering.
- 21. Type in 1 for sector rows and for sector columns. After we have created this array (field) with 4x4 sectors and each sector with 18 rows and 10 columns we need to print to more arrays on the same slide with the same configuration.
 - 22. Click on the Apply button to apply these changes to the current field.
- 23. Select Copy Field from the Fields menu at the top of the screen. CloneTracker DB will now identify the first Field as Field 1 and the new one as Field 2.
 - 24. Click on the Field 2 tab.
 - 25. Select the general sub-tab.
 - 26. Type 4,000 for Field location X and 53,000 for Field location Y.
 - 27. Click on the Apply button to apply these changes to the current field.

REFERENCES

Lee, M-L., Kuo, F.C., Whitemore, G.A., and Sklar, J. 2000.

Importance of replication in microarray gene expression studies: statistical methods and evidence from repetitive cDNA hybridizations. PNAS USA, 29:9834-9839.

Montgomery, D.C. 1991. **Design and analysis of experimental design.** John Wiley & Sons, New York, NY.

NEW FEATURES

Dynamically Structured Database - Unlike common databases, CloneTracker DB's has been designed to grow with new data and information. For example, virtually any type of additional clone data, whether assesion number or chromosomal info, can be added to the database during import.

Extremely Flexible File Formats supporting XML - Researchers are now free to import files in various formats with virtually any number of columns including, using the use of XML.

Enhanced Visual Manipulations - Zoom, move and align arrays with a mouse click.

Advanced Spot Duplication Tools - CloneTracker DB includes the replicates editor allowing for almost any form of spot replication pattern.

New Online Help System - Java based help system.

Enhanced Query Capability - New built in query tools allow researchers to retrieve information on a variety of microarray information.

Java Core - Java allows CloneTracker DB to be a robust, cross platform application with enhanced memory management.

Staggered Spotting - Users whose machine's support spot staggering are now supported within CloneTracker DB 2.

MySQL Database - CloneTracker DB includes a powerful, cross platform, relational database system.



