



i-limb™ ultra revolution User Manual

Part number: MA01141: Issue No. 1, April 2013



This document provides information for i-limb ultra revolution users and should be used in conjunction with advice from your clinical team.



This symbol signifies important information and is used throughout the manual.

A separate USB datadrive is included with your kit that contains all relevant product manuals. You may also refer to www.touchbionics.com to ensure the latest copy of this document.

Table of Contents

1	i-limb ultra revolution	1.1	Product Description
		1.2	Know Your Prosthesis
2	Socket Interface	2.1	Socket Interface and Control
		2.2	Battery Charging
3	biosim	3.1	biosim Overview
		3.2	biosim Connecting
		3.3	Navigating biosim
		3.3.1	Features
		3.3.2	Training
		3.3.3	Hand Health Check
		3.3.4	Exit
4	i-limb ultra revolution Coverings	4.1	Cover Options
		4.2	Donning the i-limb skin active cover
		4.3	Doffing the i-limb skin active cover
		4.4	Donning the i-limb skin natural cover
		4.5	Doffing the i-limb skin natural cover
		4.6	Wear and Care Guidelines for i-limb natural and custom covers
5	Support Information	5.1	Troubleshooting
		5.2	Frequently Asked Questions
		5.3	Warnings and Precautions
6	User Information	6.1	User Details
7	Appendix	7.1	Technical Information
		7.2	Component Compatibility
		7.3	Warranty

1.0 i-limb ultra revolution

1.1 Product Description

We are pleased that you and your clinical team have decided that the i-limb ultra revolution is the most appropriate prosthetic hand for your needs. You may have discussed your functional goals with your clinical team. This manual along with the training and support of your clinical team should help you understand how the i-limb ultra revolution will help you accomplish these goals.

Each digit of the hand has its own motor. This allows the digits to run until they meet the object being grasped, so the hand takes the shape of that object (compliant grip). Fig. 1 is an example of a hand holding a ball with compliant grip. The power for the hand comes from a battery which has been built into the socket of your prosthesis that your prosthetist has made especially for you. Initially you and your clinical team will develop control of opening and closing your hand. At this stage, you should be able to carry out a wide range of functional daily activities using your i-limb ultra revolution hand.

It takes time and practice to gain control and master the best way to do tasks with your prosthetic hand. The pace at which individuals gain this control varies, but we encourage you to be patient and seek training.

Your i-limb ultra revolution is covered under the Touch Care program. The duration of the program varies and you can contact your clinician regarding this program.

Figure 1.



1.2 Know Your Prosthesis

Each i-limb ultra revolution has a unique identifying number called the serial number as depicted on the image to the right (fig. 2). This is an important piece of information and we recommend that you make a note of it on page 34 of this manual.

The i-limb ultra revolution will not provide you with sensation and so you will not be able to feel things such as heat and moisture. The i-limb ultra revolution cannot absorb any shock as your own body is able to do.

If high degrees of force are applied, you may cause damage to your i-limb ultra revolution. If this does happen, then please contact your clinical team who will contact Touch Bionics customer care team to resolve the matter as quickly as possible.

To best protect the hand when carrying objects, the weight should be distributed evenly across the digits as close to the knuckles and palm of the hand as possible rather than out on the tips of the digits. The image to the right (fig. 3) is an example of how to correctly carry objects.

The i-limb ultra revolution is like any other electrical device so do not immerse in water and do not pour water over the device. When you have to clean the covering then use a damp cloth to remove any dirt.


If you have any specific questions about your hand then please contact your clinical team or visit the Touch Bionics web site www.touchbionics.com for the latest information and news.

Figure 2.



Figure 3.






2.0 Socket Interface

2.1 Socket Interface and Control

To fit the i-limb ultra revolution you have had a socket especially made for you by your clinical team. They will have made this to ensure it is comfortable and fits directly onto your residual limb. You may experience some redness but this should not last for very long after you remove the socket. If the socket is causing you discomfort or affecting you using your i-limb ultra revolution then we recommend that you contact your clinician.

As part of the fitting process, your prosthetist selected the best method for you to provide control of your device. This tells the processor (the brains of the hand) whether you want the hand to open or close. The most common control option is called an electrode (fig. 4). To the right you will see two examples of the type of electrode that may have been used to build your prosthesis. Electrodes pick up on the electrical activity that your body gives off when you contract a muscle. While gaining good control may take time and practice, your method of input should not cause pain. Talk with your prosthetist if you do not understand how your input control is working or if it is causing discomfort. There are also other potential control options that your prosthetist may have determined to be the best option for you. If so, they will describe to you how that option works.

Figure 4. Electrode Options



2.2 Battery Charging

As the i-limb ultra revolution is powered by a battery, you will need to charge it. How often you need to charge it will depend on how much you use the i-limb ultra revolution each day. You do not have to worry about over charging your battery.

First, remove the the socket from your arm and turn the hand off. The on/off switch is pictured on the image to the right. Hand is off when switch is in left position.

Only use the charger that was supplied with the your i-limb ultra revolution.

If you intend to travel outside of your home country please ensure you have a Touch Bionics charger that will work in the country to which you are traveling. Additional chargers are available from Touch Bionics. Please contact your clinician for details on how to obtain them.

For customers residing in parts of Europe and the United States, the charger pictured to the right (fig. 5) is used. The light display is either:

Red – rapid charge

Green – fully charged

Charging time is approximately:

1,300 mAh battery- 90 minutes

2,000 mAh battery- 180 minutes

Note: If the green light is on when you first plug in your device, ensure that the hand is off.

For customers residing in the UK, Australia and South Africa the charger illustrated to the right (fig. 6) is used. The light display is as follows:

Solid Amber – on standby

Slow flashing amber – pre-charge mode

Rapid flashing amber – Error

Slow flashing green – maintenance charge

Rapid flashing green – rapid charge

Solid green – fully charged

Charging time from full discharge is approximately:

1,300 mAh battery 180 minutes

2,000 mAh battery 180 minutes

To charge, insert the charger lead connector into the charge port. A “click” should be heard on connection.

Insert the charger into the power outlet.

To remove the charger lead connector from the i-limb ultra revolution, grip the connector and pull directly away from the port. Consult Warnings and Precautions for additional information.



On/Off Switch

Figure 5.



Figure 6.



As an alternative to charging directly from domestic power, a car charger (fig. 7) is provided with your i-limb ultra revolution.

Figure 7.



DO NOT PULL THE CABLE TO REMOVE THE LEAD.



To ensure the i-limb ultra revolution is continually functional, charge at the end of each day.



Switch the hand OFF to preserve battery power when not in use.



Replace the battery annually for optimal performance.
Contact your clinical team to order new batteries.

3.0 biosim

3.1 biosim Overview

biosim is a software program that wirelessly connects to i-limb ultra revolution using a Bluetooth® connection.

You can use biosim to select the grip patterns and hand features that you want to use, including the option to create your own custom gestures. Built-in training modes are included to help you achieve better results with your prosthesis. It also allows you to view the signals that are being used to control the hand, which can be useful if you are struggling with your control.

Your prosthetist will also have a version of biosim to make changes to how your i-limb ultra revolution functions.

To use the biosim mobile application with the i-limb ultra revolution you will need to download the biosim application (downloadable from the Apple® App Store) to the iPod® Touch supplied by Touch Bionics or compatible device: iPhone® 4s, iPhone® 5, iPad® with retina display or iPad® mini.

A PC or Mac loaded with the biosim software and used with the biosim Bluetooth handshake receiver may also be used.

PC requirements are as follows: Window XP, Windows Vista, Windows 7, or Windows 8; Microsoft.NET framework v3.5; USB port for Bluetooth® connector.

Mac requirements are as follows: Mac OS X Lion (10.7) and higher, USB port for Bluetooth connector.

For users of the biosim mobile application consult the iPod Quick Start Guide included with your iPod. You can also download the iPod Quick Start Guide from www.touchbionics.com.

Occasionally updates to the biosim software will be provided. An automated notification will advise you when this is needed. Follow the steps as advised in the notification.

Note: if you would like to utilize your prosthesis on devices such as smartphones or tablets, you may use a stylus or touchscreen gloves.

3.2 biosim Connecting

biosim can be downloaded at <http://www.touchbionics.com/biosimdownload>. After downloading and installing, open the biosim program by clicking on the biosim icon (fig. 8) which should be clearly visible on your computer's desktop. The biosim handshake USB (fig. 9) must also be inserted into a USB port to allow the Bluetooth® signal to be received from the i-limb ultra revolution prosthesis. The Handshake Bluetooth® receiver will pick up signals within a 10 meter range (32 feet).

The opening welcome screen will load and the first numerical icon "insert handshake" request will flash. On insertion of the Bluetooth receiver (biosim USB handshake) the icon will be constantly lit, the "handshake installed and ready" box will then self-tick and the second icon "Turn device off, then back on" will begin to flash. i-limb ultra revolution does not need to be turned off, then turned back on for this step. The third icon "Connecting to device" will be constantly lit. The tab marked "connect" must now be selected for the connection to be made, this may take up to 24 seconds.


Figure 8.



Figure 9.



If there is more than one i-limb device within range of the Bluetooth receiver then a box will appear listing all devices by serial number. The correct i-limb ultra revolution device can then be selected from the list.



3.3 Navigating biosim

Once you are connected the home screen will now be displayed, with five options as illustrated. This homepage can be accessed at any time from subsequent pages.


Features are where you will find the different grip patterns and gestures the hand can perform. You can also change what grip patterns are programmed from this page.

Training provides access to the training suite and a selection of games to help improve overall control.

Hand health check icon provides a quick and easy diagnostic check of the hand to ensure it is working properly.

Don or Doff icon positions the i-limb ultra revolution for putting on or taking off a cover.


Exit icon draws the session to a close, exiting the program.



3.3.1 Features



Click on the features icon to enter the features. The features are the hand positions and grips that can be programmed onto your i-limb ultra revolution and triggers are the muscle signals you give in order to enter the feature.







Precision Pinch Grip Options

Precision pinch grip options are best for picking up small items between and the thumb and index finger. There are 4 options available depending on how you want the other digits to perform while doing the pinch. The most popular is Thumb Precision Pinch Closed.

Feature	Picture Example	Description	Use	Task Examples
Standard Precision Pinch Opened		Middle, ring and little fingers remain fully opened and switch off. Both index finger and thumb will move to provide grip.	Allows for a wider opening than thumb precision. Aids with visualization or for pinching objects where the non-active digits may get in the way.	1. Returning cards or money to wallet 2. Picking up napkins 3. Folding laundry
Thumb Precision Pinch Opened		Middle, ring and little fingers remain fully opened and switch off. Thumb automatically moves to a partially closed position. Only index finger will move to provide grip against the fixed thumb.	Accuracy is improved when picking up an object by allowing you to place the thumb against the object to be pinched. Only the index finger moves to grasp the object. Ideal for repetitive tasks.	1. Pick up pencil or slim, long objects 2. Thread needle 3. Sort/Pick up medications
Standard Precision Pinch Closed		Middle, ring and little fingers automatically close and switch off. Both index finger and thumb will move to provide grip.	Will allow for better visualization in some tasks, especially when the working surface is not at eye level.	1. Slide small object from shelf over head 2. Pick up small object from floor
Thumb Precision Pinch Closed		Middle, ring and little fingers automatically close and switch off. Thumb automatically moves to a partially closed position. Only index finger will move to provide grip against the fixed thumb.	Can improve accuracy for picking an object by allowing you to place the thumb against the object to be pinched and only the index finger moves to grasp the object. Ideal for repetitive tasks.	1. Pick up and open sugar packet from a coffee stand 2. Pick up coins 3. Alternative way to tie shoes (also see "lateral grip")







Tripod Grip Options






Tripod Grip Options provide a grip where the thumb is meeting up with the index and middle digits. Typically the thumb is rotated between the index and middle to provide more stability when gripping slightly larger objects than those with the precision pinch grips. As a result of having a motor in each digit and the hand conforming to the shape of the object being held, the hand will automatically form a tripod grip when the object is small enough to only be gripped between the thumb, index, and middle digits. The advantage of using the tripod grip options is if you want the ring and small digits to stay either fully open or fully closed when gripping.






Feature	Picture Example	Description	Use	Task Examples
Standard 3 Jaw Chuck (Tripod) Opened		Ring and little fingers remain fully opened and switch off. Thumb, index and middle fingers will move to provide grip	Allows for larger objects to be grasped using the tripod pinch. Can also improve positioning. May also decrease compensation at shoulder for tasks where closed digits get in the way of the pinch.	1. Preparing food with ring and little finger clear, such as slicing cucumber or carrot 2. Holding cell phone (would manually rotate thumb laterally based on size of phone or setup custom grip)
Standard 3 Jaw Chuck (Tripod) Closed		Ring and little fingers automatically close and switch off. Thumb, index and middle fingers will move to provide grip.	Allows for the wider opening when using this feature with larger objects. Tripod grip allows greater stability for larger or round objects in comparison to precision pinch.	1. Open bag of chips 2. Pull up socks
Thumb 3 Jaw Chuck (Tripod) Opened		Ring and little fingers remain fully opened and switch off. Thumb automatically moves to a partially closed position. Only index and middle fingers will move to provide grip against the fixed thumb.	For smaller motion, but increased accuracy. Tripod grip allows greater stability for larger or round objects in comparison to precision pinch. Can also decrease elbow compensation when seated.	1. Grab spices from cupboard
Thumb 3 Jaw Chuck (Tripod) Closed		Ring and little fingers automatically close and switch off. Thumb automatically moves to a partially closed position. Only index and middle fingers will move to provide grip against the fixed thumb.	Can improve stability and control grasping larger objects. With thumb in a static position, your ability to align an object can improve particularly on small objects. Can also decrease elbow compensation when seated.	1. Repetitive factory type work moving pieces 2. Grasp pen for writing 3. Laundry 4. Pack bag

Additional Grip and Gesture Options

With these additional grip patterns and gestures, consider how often you perform different tasks and if having a program for that activity would make it easier/faster for you to perform. Since each digit has its own motor, you can put pressure against any digit to stop it from moving (called stalling the digit). Many of these additional grips can be performed by stalling one or more digits. If it is a task that is not performed often (maybe using index point to press the "Shift" key on a keyboard) then you may prefer to just stall the digit. However, if you frequently perform the task or want to go into the position without needing to put the hand against the table or your other hand, then programming the grip may be the better option.

Feature	Picture Example	Description	Use	Task Examples
Thumb Park Continuous		All four fingers remain fully open and switch off. Only the thumb will move.	For longer dressing period that will require more than the 1.5 seconds of thumb park quick or grasping light weight flat objects. Can also use stalling out digits to complete, putting pressure against index to little and closing thumb in to hand.	1. Put on jacket 2. Grasp flat objects, such as a book or a tablet computer or clipboard 3. Allows readjustment of the object opening letters or sorting paperwork.
Thumb Park Quick		Thumb automatically rotates into opposition and fully closes. Thumb remains active. All four fingers fully open and switch off. For 1.5 seconds only thumb will move, after which the hand automatically returns to normal function.	Dressing or putting on a coat. You can achieve this position by using the stall out capability of digits against a table if you do not want to program.	1. Put on jacket or shirt
Lateral Grip		Thumb automatically rotates into lateral position and partially closes. All four fingers fully close and switch off. Only thumb will move.	Holding onto plate, papers, CD or other flat objects. Improved stability with use of the side of the index finger.	1. Hold plate while serving food 2. Hold clipboard 3. Get card from ATM 5. Open ziplock bag 6. Tie Shoe Laces
Index Point		Thumb automatically rotates into lateral position and fully closes. Thumb switches off. Middle, ring and little fingers fully close and switch off. Index finger remains open and active	Used for pushing buttons, typing on computer, or simply pointing. Can be achieved with stall out.	1. Keyboarding (hitting shift key, CTRL, ALT, DEL, general typing) 2. Push Button for elevator
Customer Gesture		All digits automatically move to a user defined position with each digit option as open or closed.	Allows user to fully customize the hand. Many positions can be obtained using stall out if the user does not have enough available triggers to use this feature.	1. Boy scout salute 2. Sign Language gestures such as "I love you"
Custom Grip		Digits move to user defined position with options for degree of open or closed and option to keep selected digits "active".	Can be used to fully customize for use of specific tool or work or hobby. Multiple can be set up as "favorites" on iPod.	1. Customize grip to specific work or hobby type tool 2. Supporting and operating a camera

Feature	Picture Example	Description	Use	Task Examples
Rotate Thumb		Thumb and all four fingers fully open and switch off. Only thumb will rotate	No flexion or extension of the digits will occur, one signal will rotate thumb out to the side and the other will rotate into the palm	1. Used to position the thumb when you do not want to use the manual override. Can ONLY be used if the option to rotate the thumb when exiting the feature is turned off.
Open Palm		Thumb automatically rotates laterally and closes partially. All fingers remain open. All digits switch off.	For tasks where balancing an object on the open hand is useful.	1. Stabilize a tray of food 2. Balance a Bowl
Cylindrical		Thumb rotates automatically to a semi lateral position. All digits are active in this feature	Many handled type objects require this cylindrical grip. You can also use the manual override of the thumb to position manually rather than having to use this feature.	1. Broom Handle 2. Shovel or Rake 3. Hair Brush or Hair Dryer
Mouse		Hand forms shape appropriate for using a computer mouse.	Note that multiple steps are involved in this grip: Step 1: Enter the feature using either the trigger or "Quick grip" Step 2: Place the closed little finger against the side of the mouse. Step 3: Provide a close signal until the thumb rests against the other side of the mouse (the thumb must detect pressure, so you may also provide pressure with your other hand if the thumb will not stop) Step 4: The index digit will now close with the signal to click the left mouse button, relaxing signal will re-open.	1. To left click, use signal to close index digit 2. To right click use motion of the arm to put pressure on the mouse from the ring digit
Handshake		Thumb automatically rotates into semi lateral position and partially closes and switches off. All fingers remain fully open. When a close signal is given the fingers close, when this signal is relaxed the fingers open.	Note that with this feature the thumb moves to a slightly closed position and only the digits will move with the closed signal. As soon as you relax the signals the hand will automatically re-open.	1. Shaking hands

Feature	Picture Example	Description	Use	Task Examples
Grasp		Hand adopts partially closed position. When a close signal is given the hand closes at full speed. (no proportional control). When the close signal is relaxed the hand will fully open and then reposition into its partially closed starting position again.	For grasping objects where crushing the object is not a concern. Note that as soon as you stop giving a signal the hand will open automatically and reset to the slightly flexed position.	1. Sorting objects 2. Grasping large soft objects such as socks when packing a suitcase
One Finger Trigger		Thumb automatically rotates into opposition and partially closes. All fingers partially close. When a close signal is given the thumb, middle, ring and little fingers close. Once they have stalled around an object the next close signal will operate the index finger in squeeze mode (when a close signal given index will close, when signal relaxed index will open).	Note this feature requires multiple steps: Step 1: Enter the feature using either the trigger or "Quick grip". Step 2: Provide a close signal until the thumb, middle, ring, and small digits have fully closed around the object being held. Step 3: Providing a close signal will close the index digit onto the trigger and relaxing the signal will re-open the digit.	1. Pulling trigger on a small spray bottle
Trigger Two Finger		Thumb automatically rotates into opposition and partially closes. All fingers partially close. When a close signal is given the thumb, ring and little fingers close. Once they have stalled around an object the next close signal will operate the index and middle fingers in squeeze mode (when a close signal given index and middle will close, when signal relaxed index and middle will open).	Note this feature requires multiple steps: Step 1: Enter the feature using either the trigger or "Quick grip" Step 2: Provide a close signal until the thumb, ring, and small digits have fully closed around the object being held Step 3: Providing a close signal will close the index digit onto the trigger and relaxing the signal will re-open the index and middle digits	1. Pulling trigger on spray bottle
Thumb Trigger		Thumb automatically rotates to a semi lateral position and partially closes. All fingers partially close. When a close signal is given, fingers will close around the object, once fingers are stalled the thumb will operate in squeeze mode (when a close signal given thumb will close, when signal relaxed thumb will open).	Note this feature requires multiple steps: Step 1: Enter the feature using either the trigger or "Quick grip". Step 2: Provide a close signal until the index, middle, ring, and small digits have fully closed around the object being held. Step 3: Providing a close signal will close the thumb onto the trigger and relaxing the signal will re-open the thumb.	1. Pressing down on a small aerosol can button
Donning or Doffing a cover		Hand forms the proper shape for donning and doffing a cover.	When changing out covers for the hand to avoid too much pressure against the thumb.	1. Putting on covers 2. Taking off covers

Triggers

Your muscles create signals that you will use to control the i-limb ultra revolution. These signals are called triggers. Four potential triggers are available to activate features: Hold open, Co-contraction, Double impulse and Triple impulse.

- 1 Hold open: a prolonged open signal for the programmed amount of time, normally 2-3 seconds.
- 2 Co-contraction: Both the open and closed signals must be activated at the same time. Your clinical team may have a recommendation on how to do this based on your control, such as imagining making a tight fist, flicking your fingers, or forcefully straightening your arm.
- 3 Double impulse: two rapid open signals, separated by relaxation below the threshold. Impulses are very quick open signals. They do not need to be strong, just fast like blinking your eyes. Your prosthetist will have pre-set this at the time of fitting. To activate a double impulse the device needs to be fully opened prior to attempting activation.
- 4 Triple impulse: three rapid open signals, separated by relaxation below the threshold. Impulses are very quick open signals. They do not need to be strong, just fast like blinking your eyes. Your prosthetist will have pre-set this at the time of fitting. To activate a triple impulse the device needs to be fully opened prior to attempting activation.

Any of the above four potential triggers can be linked with any of the above features. Some users are able to perform activities of daily life (ADL's) without the use of any triggers. In many cases users initially set one or two triggers while familiarity and control is gained when learning how the i-limb ultra revolution works. If you have several features set up with different triggers and notice that you accidentally jump into a particular feature often, you may want to stop using that trigger until you can follow up with your prosthetist to determine why it is happening inadvertently. Your clinical team may be able to make adjustments or training recommendations to improve your control.

Note on Quick Grips:

In both the biosim software on your computer and on the app, you will find "Quick Grips". This is another way you can enter and exit features. By clicking the "Quick Grip" button for the selected grip you hand will enter that mode. To exit, click the "Quick Grip" again.


To Exit Features:

You may use 3 different options:

1. Do the same thing you did to get into the feature (either the same trigger or pressing the Quick Grip)
2. Choose to go from one feature to a different feature by either using the trigger linked to that feature OR select the Quick Grip of the feature you want to enter
3. Turn the hand off and back on - this option will always return the hand to the full hand normal operation

Linking Triggers with Features

In order to link a trigger with a feature simply click on to the desired feature (the feature will now be highlighted). Select and click the desired trigger that you would like to link with the highlighted feature. Both feature and trigger should now be highlighted and the descriptive label will now appear under the feature icon. The example illustration on the right indicates the Index Point Feature linked to the Hold Open Trigger.




By clicking on any of the features on the features page the available triggers and example sections will also appear. In order to preview a feature simply highlight the feature and click on example to the right of the screen. The section will then enlarge and provide a demonstration of the feature. The hand image can be rotated to improve visualization by holding the right mouse key and moving the cursor in the direction of the intended rotation. The left mouse key will zoom in and out while the right mouse key will allow you to drag the hand around on the screen to a different view point.




Additional Functions

Toward the bottom right of the screen the additional Global Options box can be used to give access to Vari-grip/pulsing, Rotate thumb on exit, and natural hand mode.




1 Vari-grip/pulsing: this mode provides additional grip force when the closed signal is sustained for the amount of time selected (default 0.5 seconds). This can also be customized between 250ms (0.25 seconds) to 3,000ms (3 seconds). Vari-grip is excellent for tasks such as tying shoes, opening bags (ziplock, chips, etc.), tearing paper or opening packages.



2 Rotate thumb on exit: you can choose to turn this off or on. If off the thumb will remain in the same position as it is located during the last feature used. If on the thumb will always move to the fully open (off to the side) position.

3 Natural Hand Mode: this mode allows you to select what position you would like the hand to move to when it is not in use. Once the hand is fully opened and you relax your signals for the selected amount of time, the hand will move to a more relaxed position. The timing for this position can be adjusted as well as how far opened or closed each digit moves.



3.3.2 Training








The training suite contains a variety of short exercises aimed at developing control of the i-limb ultra revolution device. The opening screen highlights a series of exercises which can be selected individually and in any order. Both open and closed signals can be practiced by selecting the appropriate box to the right. An indication of difficulty is provided by the 5 point scale on each module.

Real time signals: This graph allows you to view the signals you are providing to the hand. The red signals are your open signals and the blue are closed signals. Press the play button to start the graph and stop when finished. The dark lines with arrows to the left show how strong your signals must be before the hand starts to pay attention to them. To make the hand run at full speed, the signals should be about half way up the graph. Your signals may be too strong if they peak out at the very top of the graph. Ideally your open and close signals should operate independently with no interference between the two signals. If interference does occur then this can impact on your ability to control your hand. If your signals are either too low or too high or you are having difficulty controlling your hand, you should contact your prosthetist. The “done” button allows you to return to the training page.


Strength: This is to help you focus on generating a strong muscle signal one muscle at a time. At the right of the screen select either the open or close icon to work on that signal. Press the start button and then give the signal you have selected to try to get the ball above the line. See how many times you can successfully get the ball to rise. You can view the strength of your signals on the right by watching the red and blue rise. Red is open and blue is closed. Press “stop” and “done” when finished.

Speed Only: This is to help you generate quick muscle signals. Quick signals help with giving quick impulses. Typically the open signal is the one that needs to be quick. Work to make the ball rise to the top of the hill as quick as possible.


Strength and Separation: This game focuses on generating strong muscle signals while isolating the opposing muscle. Choose to work on either open or close then press “start”. Only let the balance scales move in the direction you have selected.




Co-contraction: Focuses on generating quick strong simultaneous contractions with both muscles. Press Start then try to get both balls to reach the top of the hill at the same time indicating both signals being activated. Mastering this will aid you to be successful with the co-contraction trigger.



Hold Open: Focuses on generating a strong open muscle signal and maintaining it at a high level for several seconds. After pressing "start" aim to hold the open signal above the bar for 2-3 seconds to be successful. Mastering this exercise will assist acquirement of the hold open trigger.




Speed and Separation: Focuses on generating quick, strong muscle signals while isolating the opposing muscle. This exercise lays the foundation for impulse control. Remember, impulses are the very quick signals that can be used to change from one grip pattern feature to another using a double or triple impulse trigger. The open signal should be focused on for this purpose.




3.3.3 Hand Health Check



The Hand Health Check screen provides a quick and easy check of the hand. Click on the "Run Hand Health Check" icon and the health check will begin. There must be an internet connection to run the Hand Health Check.



The i-limb ultra revolution hand will then go through a series of movements as each digit is checked, the process will run for approximately 8 seconds and provide basic feedback on each step. If there is a problem with your hand, please contact your clinician.



3.3.4 Exit



Exit

Use the Exit icon to fully exit biosim.

A prompt will be provided allowing you to save the settings for that session for future reference.

Do not power i-limb ultra revolution down until configurations are saved and biosim closes.

4.0 i-limb ultra revolution Coverings

4.1 Cover Options

The cover of the i-limb ultra revolution is very important. The i-limb ultra revolution hand should not be used without an approved cover that is well maintained.

A number of options are available to cover the hand. All covers require regular checks for wear and tear and some require routine maintenance. Covers will need replacement depending on usage.

The i-limb skin active cover is designed for the i-limb ultra revolution hand and will fully cover the hand. The palm and inner surface of the cover provide some frictional properties for slip resistance when gripping. The back of the glove is smooth to allow the hand to easily move through sleeves of clothing. You are eligible to receive replacement i-limb skin active covers if you are covered under the Touch Care program.




The i-limb skin natural cover is designed to be close to natural human anatomy. A color swatch is available to select the closest color match between your natural skin color and the color of your i-limb skin natural cover. There are ten color options for the i-limb skin natural cover, the glove may not be an exact match, but will be very close.

For individuals who desire the most realistic covering option possible, i-limb skin match is available, which is a fully customized silicone covering that is hand painted and exactly matches the skin tone and features.



4.2 Donning (putting on) the i-limb skin active Cover

- 1 There are two ways to position the i-limb ultra revolution to don the active skin cover. The first way is to manually position your i-limb ultra revolution. The 4 fingers should be fully open and the thumb closed with a 15mm space between thumb and index finger. The hand should then be switched off to ensure this position is maintained. The second way would be to select the “Don or Doff” button in biosim. By selecting this button, the i-limb ultra revolution will automatically position itself to the correct position.
- 2 Align the i-limb skin active cover with the digits and pull the cover down over the digits. Gently pull up and over the thumb, taking care not to cause any undue downward pressure on the thumb.
- 3 Pull the i-limb skin active cover over the hand to the wrist.
- 4 Individually maneuver the fingers of the i-limb skin active cover over the fingers until they are fully aligned. It is important not to pull the cover too tightly, the cover is designed to allow the hand to function optimally, if the cover is pulled too tight it will negatively affect hand function.



4.3 Doffing (removing) the i-limb skin active Cover

- 1 There are two ways to position the i-limb ultra revolution to doff the active skin cover. The first way is to manually position your i-limb ultra revolution. The 4 fingers should be fully open and the thumb closed with a 15mm space between thumb and index finger. The hand should then be switched off to ensure this position is maintained. The second way would be to select the "Don or Doff" button in biosim. By selecting this button, the i-limb ultra revolution will automatically position itself to the correct position.

Gently maneuver the i-limb skin active cover from the wrist up and over the base of the thumb taking care not to cause any undue downward pressure on the thumb.



- 2 Pull the fingertips of the cover upwards on each finger in turn to loosen the cover. Once loosened the cover can be gently pulled up and off the i-limb ultra revolution hand to remove.



4.4 Donning (putting on) the i-limb skin natural Cover



It is highly recommended that you visit your clinician to don and doff your i-limb skin natural cover. If your clinician is not available, it is highly recommended that you seek assistance.

- 1 There are two ways to position the i-limb ultra revolution to don the natural skin cover. The first way is to manually position your i-limb ultra revolution. The 4 fingers should be fully open and the thumb closed with a 15mm space between thumb and index finger. The hand should then be switched off to ensure this position is maintained. The second way would be to select the "Don or Doff" button in biosim. By selecting this button, the i-limb ultra revolution will automatically position itself to the correct position.
- 2 Spray the Silicone outside of the covering generously with the IPA (Isopropyl Alcohol) lubricant / cleaner supplied.
- 3 Roll the covering back on itself, from wrist to fingers until the opening for the fingers is showing.
- 4 Push the bat provided into each of the fingers and thumb in turn to straighten the lining. While doing this ensure the fingers of the cover align side by side and are not folding or overlapping each other. This will help the covering slide onto the extended i-limb ultra revolution fingers.




- 5 Spray into the opening at the wrist and on the outside of the inverted cover generously with the IPA lubricant / cleaner supplied. Ensure the IPA does not contact the fabric lining, it should only be sprayed onto the silicone.



- 6 Slide the covering onto the four fingers of the hand until only the aluminium knuckle is visible.



- 7 Place the thumb hole of the covering over the thumb tip. Hold the wrist section of the covering and pull on either side. In one motion pull over the thumb, ensuring it does not cause any undue pressure on the thumb.



- 8 Once the covering is over the digits and roughly half way down the thumb, take a firm hold of the cover at the wrist end and gently maneuver over the thumb and hand – take care not to cause undue pressure on the thumb.



- 9 Once glove is fully donned massage cover to ensure the lining has no folds or tightness which may affect the hand function. Position the thumb in opposition and open hand fully. Check if the thumb is able to fully open or if there is tightness in the web space restricting it. Close the hand and ensure a tip to tip grip can be achieved with the index finger. Then position the thumb in the lateral position and ensure it can open and close fully without restriction. Fully close the fingers. You will notice the little finger does not close as far as the others, it should close within 25mm of the palm. If the hand does not perform optimally in any of these positions massage the cover and lining further



Note: When using the i-limb natural cover, adjustments may be required in biosim to ensure that digit alignment during functional tasks is not sacrificed. Contact your clinician or Touch Bionics with questions regarding these adjustments.

4.5 Doffing (removing) the i-limb skin natural Cover

- 1 There are two ways to position the i-limb ultra revolution to doff the natural skin cover. The first way is to manually position your i-limb ultra revolution. The 4 fingers should be fully open and the thumb closed with a 15mm space between thumb and index finger. The hand should then be switched off to ensure this position is maintained. The second way would be to select the "Don or Doff" button in biosim. By selecting this button, the i-limb ultra revolution will automatically position itself to the correct position.

Start by gently pulling the i-limb skin natural cover from the i-limb ultra revolution. You should begin to feel the finger tips of the i-limb skin natural cover begin to loosen from the i-limb ultra revolution.

- 2 Slide the covering gently up and over the thumb and continue over the digits.

- 3 Remove the covering entirely by continuing to gently draw the fingers of the cover off the digits and off the i-limb ultra revolution.



4.6 Wear and Care Guidelines for i-limb skin natural and custom covers

Care of the covers is very important in order to maintain its use long-term. The following are some guidelines for wear and care.

Cleaning the Cover

Use warm water and plain soap regularly to clean the outside of the cover. Using the correct soap is very important - plain brand liquid dish washing soaps such as Dawn® or Ivory® are recommended. Clean the cover with medical grade rubbing alcohol once a week to help with disinfection. When cleaning you should not run the hand under the water, but instead use a damp cloth with the water and soap.



Some soaps contain additives such as oil, perfume, hand lotion, glycerin or aloe and should be avoided.

Cleaning Precautions

Oils and oily substances alone may not damage the cover. However, oily substances attract and create an adherence of dirt, grime and other chemicals which can lead to damage of the cover if it is not washed and cleaned regularly as described above. Many hair and body products contain oil and should only be used when the prosthesis is not being worn. Exposure to oils and oily substances may make the prosthesis look shiny. Do not use oil control cleansers to manage the shine, these are abrasive and will give the cover a tarnished appearance.

Body perspiration contains oils which can accumulate and damage the prosthesis if not removed by appropriate cleaning.

Do not attempt to repair the i-limb skin natural, or custom cover. If the cover is damaged, do not attempt to fix it. Contact your prosthetist. Damage can often be repaired, rendering the damage almost undetectable. Cutting, gluing or altering your cover may make it irreparable and will invalidate the warranty.

Staining

Although the cover is made of a durable material, it is possible for it to become stained. Some substances to avoid include blue Woolite®, permanent markers, some pens and inks (newspaper ink will not stain), unwashed blue jeans and fabrics, bleach, chronic exposure to cigar and cigarette smoke, carbon paper and 'carbonless' carbon paper, other products may also stain your cover. The list provides examples of products that can damage the cover, when using new products for the first time, use with caution. If the cover is stained, washing quickly with warm water and plain soap and then cleaning with alcohol is most likely to remove the stain.

Failure to keep the cover clean as instructed above will invalidate the warranty on the cover.

If cleaning does not remove the stain, contact your clinician or Touch Bionics for an evaluation and repair. Refer to the back page for contact information.

Accessing the Battery or Switch

If you need to access the charge port or switch, and it is covered by a glove or sleeve, spray the outside with medical rubbing alcohol to reduce friction, this will aid partial removal and will also help prevent damage caused by friction.

Nails

The nails of the i-limb skin natural cover are made of silicone and cannot be polished or painted as this will damage the cover. Attaching acrylic nails will invalidate the warranty.

If you have a custom high definition cover the nails may either be silicone or acrylic. Acrylic fingernails may be polished as usual. Use only non-acetone nail polish remover. Silicone nails should not be polished.

Jewelry and Gloves

Jewelry is not recommended to be worn with Touch Bionics covers. A silk or smooth fabric lined glove will be easily put on and taken off the hand. Gloves should not have bleeding dyes when damp or wet.

General Precautions

- The i-limb ultra revolution hand must be used with an approved Touch Bionics cover.
- Never put more than one cover on your i-limb ultra revolution hand.
- Always use the correct size cover and the correct orientation (left or right).
- Ensure the cover is fitted properly.
- Only use IPA provided as directed in the manual.

The i-limb ultra revolution warranty may be invalid if used outside of the recommendations.

Please consult the Cover Care Guide supplied for a complete overview of cover wear and care guidelines.

5.0 Support Information

5.1 Troubleshooting

Problem	Action
Does not operate	<p>Ensure the prosthesis is switched “on” at the switch on the hand and on any addition power switch you may have on your prosthesis</p> <p>Ensure the battery is charged</p> <p>Ensure the electrodes are making good contact with your skin</p> <p>Ensure the hand is engaged at the wrist</p>
Stops midway through an action	<p>Check that the i-limb ultra revolution device is engaged at the wrist</p>
Difficult to operate	<p>Ensure battery has good charge</p> <p>Review with your clinical practitioner</p>
Battery charge does not last a full day	<p>Fully charge the battery. This may take up to 2 hours</p> <p>Ensure that you are not holding excessive sustained muscle signals. Refer to section 3.3.2 to review proper signal control of the i-limb ultra revolution device</p>
Battery is not working	<p>Ensure the battery charger is working correctly and the appropriate lights are shown on your charger when charging - Refer to section 2.2 for proper charging</p> <p>Check that battery is charged</p>
Digits open when a closed signal is activated	<p>Contact your prosthetist</p>

5.2 Frequently Asked Questions

Will I be able to move each digit individually?	No, while there are a number of customized grip patterns that will allow you to just move certain digits for pinching an object or making a certain gesture, you will not have individual control of each digit.
How should I hold a fork or knife?	There are several options of how to hold utensils. The more contact points on the utensil, the more stability you will have with the utensil. Try rotating the thumb out to the side to grasp the handle of the utensil between the thumb and side of the index digit. To get different angles you can slide the utensil between different digits such as the index and middle or middle and ring.
How long will it take to learn to use the i-limb ultra revolution?	It varies from person to person. You should anticipate a learning curve. You have probably compensated to do activities without a prosthesis or learned to use a different device in the past. You must now train to best utilize the i-limb ultra revolution and fully understand all the functional benefits. This process can take several months before use of the device feels natural. If you are struggling with specific tasks, talk to your clinical team or go to www.touchbionics.com to see if there are any videos demonstrating the activity you are struggling to complete.
Can I use the i-limb ultra revolution for sporting tasks?	Depending on the activity you are performing the i-limb ultra revolution may or may not be the most appropriate tool. Talk with your prosthetist about your goals and what device may be best to help accomplish your goals.
My hand is beeping at me, what does that mean?	Three beeps in a row is the low battery warning signal. You should charge your hand soon after you hear this warning.
My battery does not seem to last all day, what can I do?	If the battery used to last all day and no longer does, you should talk to your prosthetist. Sometimes your signals have gotten stronger resulting in greater usage. If signal is present even when you are not using the hand that will drain the battery. You could try turning the hand off when you are not using it to save battery life. Your prosthetist may be able to discuss other options.
I had good control when I first got the hand, but now it seems like I drop items more easily, why is that?	The signals you are giving may have gotten stronger over time since you have been using them more. Now even a small muscle contraction may make the hand move. If you connect to biosim you may see your signals on the graph reach the top very easily as opposed to in the past. Contact your prosthetist.
What if the device gets wet in the rain?	You should dry it off as soon as you can and always ensure your cover is free from holes. You do not want water to get under the cover and sit on the electronics. Treat the hand like you would a cell phone and protect it from water.
How should I clean it/how often?	You should clean your prosthesis and the cover over your i-limb ultra revolution regularly. Discuss with your prosthetist the best method based on the material of your socket. Often times a damp cloth with a dish soap without any lotions is used to clean the cover. You can also use alcohol to clean any stains or prevent bacteria. After using the restroom, you can wash your hands using a damp paper towel with the soap rather than running the hand under the water.
What is the best way to tie my shoes?	Some people like to use the pinch mode where the index finger and thumb come together, but often times individuals feel they get a better grip and better ability to see the shoe lace if they pinch with the thumb against the side of the hand (called a lateral pinch).
How can I put a coat on if my device is already on?	To make the hand as slim line as possible, you can put pressure against the index, middle, ring, and small finger while giving a close signal or alternatively use the thumb park feature. Both of these will bring the thumb into the palm of the hand so it does not catch inside the coat sleeve. Some people put the prosthetic hand through the sleeve first and use their other hand to move up the coat sleeve to the hand to guide it down the sleeve without sticking. Others will put a grocery bag or other slick material over the hand and then just slide on the coat sleeve normally with the bag decreasing the friction making it easier.

5.3 Warnings and Precautions

i-limb ultra revolution

Do not use without an approved cover

Do not use under water

Do not use to operate heavy / industrial machinery

Do not use with machinery with moving parts that may cause personal injury or damage

Users must comply with local regulations on the operation of automobiles, aircraft, sailing vessels of any kind and any other motorized vehicle or device

Do not use for extreme activities that may cause injury to a natural hand e.g. rock climbing

Do not expose to excessive moisture, liquid, dust, vibration or shock

Do not expose to high temperatures

Do not expose to flames

Do not use or expose to explosive atmospheres

Do not disassemble componentry or modify in any way

Maintenance, repairs and upgrades may only be performed by qualified Touch Bionics technicians and technical partners

Do not use with a damaged cover

Damaged covers must be replaced or repaired by a qualified Touch Bionics technician or technical partner

Only approved Touch Bionics accessories and tooling may be used with the i-limb ultra revolution

If utilizing an i-limb ultra revolution multi-flex wrist, the hand should be locked if carrying an object that is hot

Batteries

Do not bend or exert excessive pressure on the battery

Do not pierce the battery

Do not disassemble the battery

Do not expose to high temperatures

Do not incinerate batteries

Do not alter battery terminal wires

Do not short circuit the battery

Do not store batteries inside a vehicle

Dispose of batteries in accordance with US, European or local regulations

Only use the appropriate Touch Bionics charger to charge Touch Bionics batteries

CAUTION: If the battery has visibly ballooned or swollen:

- discontinue the charging process immediately
- disconnect the battery
- remove to a safe area
- leave and observe for 15 minutes
- replace the battery
- do not re-use
- dispose of any leaking batteries in an appropriate manner

Failure to comply with the above guidelines will invalidate the warranty.



If you experience technical problems with the i-limb ultra revolution and your clinical team is not available, contact Touch Bionics as follows:

North American Customers
(Canada, Mexico & US)
Tel: +1 855 MYiLIMB (694 5462)
UK & Non-North American Customers
Tel: +44 1506 438 556

Driving of Motor Vehicles



The i-limb ultra revolution has the functional capability to assist a patient with driving a motor vehicle however due to factors including the differences in world-wide driving regulations and the variations in the level of ability between patients Touch Bionics is unable to provide definitive advice in respect of a patient with an i-limb device driving a motor vehicle.

Touch Bionics is aware that individuals have used the i-limb to drive a motor vehicle and our recommendations prior to a patient doing so would include the following:

- contacting the driving authority in your home location to obtain and understand the local regulations;
- working with the appropriate authorities to have your car modified to meet the local regulations for your respective disabilities as required;
- re-taking any mandatory driving test using your i-limb device to demonstrate your ability to operate a motor vehicle safely if required by local regulations;
- contacting your insurance provider and advise them that you will be using the i-limb device to drive a motor vehicle;
- ensuring that the device has a fully charged battery. Please note that the i-limb device will emit a low battery signal which will alert you if the battery requires to be charged;
- switching off the i-limb device. This is due to the possibility of involuntary muscle signals being generated; and
- moving the thumb into the lateral position to allow the i-limb to be removed from the steering wheel without opening the hand.

It is entirely the patient's responsibility to seek confirmation that they are physically and legally able to drive using the device and to the fullest extent permitted by law Touch Bionics shall under no circumstances whatsoever be liable to the patient or any other party as a result of or in connection with a patient with an i-limb device driving a motor vehicle.

6.0 User Information

6.1 User Details

Provision of the following basic information will enable easy identification if the prosthesis is returned to Touch Bionics customer service.

User Name:

Fitting Date:

Hand Purchase Date:

Hand Serial Number:

Prosthetist Name & Contact Information:

Therapist Name & Contact Information:

7.0 Appendix

7.1 Technical Information

i-limb ultra revolution	
Voltage	7.4 V (nominal)
Max. Current	5 A
Battery Capacity	Rechargeable lithium polymer 7.4 V (nominal); 2,000 mAh capacity; 1,300 mAh capacity
Max hand load limit (static limit)	90kg/198lb
Finger Carry Load (static limit)	32kg/71lbs
Time from open position to full power grip	1.2 seconds
Device Weight	Weight with QWD 507g (small) and 515g (medium) Weight with WD unit 443g (small) and 4g (medium)

7.2 Component Compatibility

1. General Safety

1.1 The i-limb ultra revolution device is an electrical device, which under certain circumstances could present an electrical shock hazard to the user. Please read the accompanying user manual thoroughly and follow directions stated in the manual to assure maximum safety during charging and operation.

1.2 EN 60601-1:2006

1.2.1 Protection against electrical shock – Class II

1.2.2 Degree of protection against electrical shock – Type BF provides additional protection against electric shock

1.2.3 Degree of protection against ingress of water (IEC 60529:2001) – IP40

1.2.4 Not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide

1.3 EMI/EMC

1.3.1 Compliance against standard EN 60601-1-2:2007

1.4 Radio Spectrum Matters (ERM)/Bluetooth

1.4.1 Compliance against standard EN 301 489-1 V1.8.1

1.5 EN 301 489-3 Clause 7.1

1.5.1 EN55022: 2006








1.6 Radiation emissions, Enclosure

1.6.1 EN 301-489-1 Clause 8.2 - Pass (30MHz to 6,000MHz)

1.7 Zones of Use

1.7.1 Not recommended in zones 0, 1, 20 and 21

N.B. See www.touchbionics.com for further information on EMC testing carried out on products within this manual.

	Refer to operating instructions
	Class II equipment – provides double Isolation to protect against electric shock
IP40	Degree of protection – IP40 Protection against penetration by solid particles with diameters larger than 1 mm. No special protection against penetration by water
	Batch/Lot/ID Number
	For i-limb ultra revolution devices: Each device has a guaranteed unique id number example: 0001:2012 The unique serial number for i-limb ultra revolution devices is a “R” with a 4 digit alpha / numeric number. The year of manufacture of the device is then added.
	WEEE Compliance
	Catalogue number
	Manufacturer
	Keep Dry

1. Customer Service/Contact Information:

Touch Bionics, Unit 3 Ashwood Court,
Oakbank Park Way, Livingston EH53 0TH, UK

Tel: Customer Service: +44 (0) 1506 445 415
Tel: General Enquiries: +44 (0) 1506 438 556

www.touchbionics.com

Touch Bionics, 35 Hampden Road
Mansfield MA 02048, USA

Tel: +1 855 MY iLIMB (694 5462)

www.touchbionics.com

7.3 Warranty

Limited Warranty for i-limb ultra revolution

Touch Bionics warrants that the i-limb ultra revolution will conform to its specifications and be free of defects in material and/or workmanship for twelve (12) to sixty (60) months (depending on package purchased) from the date of Touch Bionics invoice for the i-limb ultra revolution. This Limited Warranty applies only to an i-limb ultra revolution provided by Touch Bionics or an affiliate authorized by Touch Bionics to provide the i-limb ultra revolution. This Limited Warranty applies to all components including but not limited to fixtures, motors, bearings, and electronics. This Limited Warranty is governed by UK law and is not transferable.

Warranty:

Touch Bionics reserves the right to credit, repair or replace an “in-warranty” i-limb ultra revolution as its option. If required, replacements will be new products. The wearer shall report any defect claim to Touch Bionics directly or to the facility that provided the i-limb ultra revolution immediately upon discovering the defect, and, in any event, within the warranty period. The defective i-limb ultra revolution must be returned to Touch Bionics or any other Touch Bionics authorized representative. To find the nearest location, visit www.touchbionics.com or call +1-855-MY-iLIMB (US & Canada), or +44 (0) 1506 438 556 (International). The unit must be in assembled condition and include an approved covering when returned. The warranty is void if the i-limb ultra revolution is subjected to abuse, neglect, alteration, modification, improper repair and/or maintenance performed by anyone other than Touch Bionics or a Touch Bionics’ affiliate. Damage as the result of normal wear and tear including the result of fatigue is not covered during the warranty period. Damage resulting from installation of parts and accessories not compatible with the i-limb ultra revolution by anyone other than Touch Bionics or an affiliate is not covered, including use of non-Touch Bionics batteries. The warranty is void if damaged covers are not replaced or repaired in a timely manner by a Touch Bionics cosmesis facility or another facility authorized by Touch Bionics, or if an approved cosmetic covering is not worn at all times when the i-limb ultra revolution is being worn.

This is the exclusive remedy under this warranty, any and all other remedies that may otherwise be applicable are excluded, including, but not limited to, incidental or consequential damage or punitive damage to the maximum extent permitted by law. This is the only warranty made by Touch Bionics on the i-limb ultra revolution and components, and there are no warranties which extend beyond the description herein. Any warranties that may otherwise be implied by law including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose are extended.

This Limited Warranty gives the consumer specific legal rights. The consumer may also have other legal rights which vary from country to country, from state to state in the U.S, from province to province in Canada and from state to state in Mexico. Some countries and states may not allow the exclusion or limitation of incidental or consequential damages or warranties, so the above limitations or exclusions may not apply to you. If it is determined by a court of competent jurisdiction that a certain provision of this limited warranty does not apply, such determination shall not affect any other provision of this limited warranty and all other provisions shall remain in effect.



North American Customers
(Canada, Mexico & US)
Touch Bionics
35 Hampden Road
Mansfield MA 02048
USA
Tel: +1 855 MY iLIMB (694 5462)



International Customers
Touch Bionics
Unit 3, Ashwood Court
Oakbank Park Way
Livingston EH53 0TH
Scotland
Tel: +44 1506 438 556
Email: info@touchbionics.com

For address details and further information
please visit www.touchbionics.com

Third party products and brand names
may be trademarks or registered
trademarks of their respective owners

