



**KDS35 Small Animal
Ventilator**
User's Manual

WEEE/RoHS Compliance Statement


EU Directives WEEE and RoHS

To Our Valued Customers:

We are committed to being a good corporate citizen. As part of that commitment, we strive to maintain an environmentally conscious manufacturing operation. The European Union (EU) has enacted two Directives, the first on product recycling (Waste Electrical and Electronic Equipment, WEEE) and the second limiting the use of certain substances (Restriction on the use of Hazardous Substances, RoHS). Over time, these Directives will be implemented in the national laws of each EU Member State.

Once the final national regulations have been put into place, recycling will be offered for our products which are within the scope of the WEEE Directive. Products falling under the scope of the WEEE Directive available for sale after August 13, 2005 will be identified with a "wheelie bin" symbol.

Two Categories of products covered by the WEEE Directive are currently exempt from the RoHS Directive - Category 8, medical devices (with the exception of implanted or infected products) and Category 9, monitoring and control instruments. Most of our products fall into either Category 8 or 9 and are currently exempt from the RoHS Directive. We will continue to monitor the application of the RoHS Directive to its products and will comply with any changes as they apply.



- **Do Not Dispose Product with Municipal Waste**
 - **Special Collection/Disposal Required**

KDS35 Small Animal Ventilator

Contents

- Specifications.....2
- Parts and Accessories2
- General Information - Warranty and Repairs3
- General Safety Summary.....4
- Description and Features5
- Voltage Selection6
- Fuse Replacement6
- Pump Operation:
 - Connecting Animal7
 - Selecting Rate & Volume.....8
 - Setting Volume & Rate.....9
 - Use of Gas Mixture9
 - Creating PEEP.....9
- Cleaning Maintenance10
- Changing Piston & Cylinder.....11
- Repair Kit11

KDS35 Small Animal Ventilator

Specifications

Rate	Adjustable from 30 to 150 strokes/min while the ventilator is running
Phase Ratio	1:1 fixed
Dimensions, H x L x W	25 x 32.5 x 20 cm (10 x 13 x 8 in)
Weight	9.5 kg (21lbs)
Power	100-120 VAC/200-230 VAC, 0.7A/0.35A max., 50/60 Hz
Stroke Rate	Continuously variable from 30 to 150 strokes/min. while the ventilator is running, with LIGHT-EMITTING DIODE (LED) DISPLAY.
Volume	Adjustable from 0.2 to 5.0 cc/stroke or 1.5 to 30.0 cc/stroke while the ventilator is running
Environmental:	
Operating Temperature	0°C to 40°C (32°F to 104°F)
Operating Humidity	20% to 80% RH, non-condensing
Installation Category	II
Pollution Degree	I
Port:	
Part Number	5121-074
Size	
ID	.128" (3.25mm)
OD	.188" (4.78mm)
Catalog No.	78-0391

Parts and Accessories

78-0391	Small Animal Volume Controlled Ventilator
78-0395	Overhaul Kit for 5cc Small Animal Ventilator; contains O-rings, valve springs, lubricants, etc. to overhaul Ventilator
78-0396	Overhaul Kit for 30cc Small Animal Ventilator; contains O-rings, valve springs, lubricants, etc. to overhaul Ventilator
78-0397	Connection Kit

KDS35 Small Animal Ventilator

General Information - Warranty and Repairs

SERIAL NUMBER

All inquiries concerning our product should refer to the serial number of the unit. Serial numbers are located on the rear of the chassis.

CALIBRATION

All syringe pumps are designed and manufactured to meet their performance specifications at all rated voltages and frequencies. A calibration certificate is available upon request. Contact customer service for details and pricing.

WARRANTY

KD Scientific warrants this instrument for a period of one year from date of purchase. At its option, KD Scientific will repair or replace the unit if it is found to be defective as to workmanship or material.

This warranty does not extend to damage resulting from misuse, neglect or abuse, normal wear and tear, or accident. This warranty extends only to the original customer purchaser.

IN NO EVENT SHALL KD SCIENTIFIC BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES. Some states do not allow exclusion or limitation of incidental or consequential damages so the above limitation or exclusion may not apply to you.

THERE ARE NO IMPLIED WARRANTIES OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR USE, OR OF ANY OTHER NATURE. Some states do not allow this limitation on an implied warranty, so the above limitation may not apply to you.

If a defect arises within the one-year warranty period, promptly contact your local distributor or **KD Scientific, 84 October Hill Road Holliston, Massachusetts 01746-1388, 508-429-6809.** Goods will not be accepted for return unless an RMA (returned materials authorization) number has been issued by our customer service department. The customer is responsible for shipping charges. Please allow a reasonable period of time for completion of repairs, replacement and return. If the unit is replaced, the replacement unit is covered only for the remainder of the original warranty period dating from the purchase of the original device.

This warranty gives you specific rights, and you may also have other rights which vary from state to state.

REPAIR FACILITIES AND PARTS

KD Scientific stocks replacement and repair parts. When ordering, please describe parts as completely as possible, preferably using our part numbers. If practical, enclose a sample or drawing. We offer a complete reconditioning service.

CAUTION

This pump is not registered with the FDA and is not for clinical use on human or veterinary patients. It is intended for research use only.

KDS35 Small Animal Ventilator

General Safety Summary

Please read the following safety precautions to ensure proper use of your syringe pump. To avoid potential hazards and product damage, use this product only as instructed in this manual. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

TO PREVENT HAZARD OR INJURY:

Use Proper Line Cord

Use only the line cord shipped with the product and make sure line cord is certified for country of use.

Ground the Product

This product is grounded through the grounding conductor of the power cord. To avoid electric shock, the grounding conductor must be connected to earth ground. Before making any connections to the input or output terminals of the product, ensure that the product is properly grounded.

Make Proper Connections

Make sure all connections are made properly and securely.

Orient Equipment Properly

Do not position the equipment such that it is difficult to reach the disconnecting device.

Observe all Terminal Ratings

Review the operating manual to learn the ratings on all connections.

Use Proper Fuse

Use only specified fuses with product.

Avoid Exposed Circuitry

Do not touch any electronic circuitry inside of the product.

Do Not Operate with Suspected Failures

If damage is suspected on or to the product do not operate the product. Contact qualified service personnel to perform inspection.

Place Product in Proper Environment

Review the operating manual for guidelines for proper operating environments.

Observe all Warning Labels on Product

Read all labels on product to ensure proper usage.



CAUTION
Refer to Manual



*Protective Ground
Terminal*

KDS35 Small Animal Ventilator

Description & Features

The KDS35 is a volume controlled small animal/rodent respirator:

- Strong motor and motor speed control with LED (light emitting diode) read-out of rate (stroke rate or breaths/min).
- Vertical slide valve that has no "O"-Rings or passive check valves and requires no adjustments.
- Orientation allows for minimal system dead space and thus more precise tidal volume delivery to the animal

The KDS35 is a positive volume ventilator intended to meet the respiratory requirements of small laboratory animals weighing up to 5 kg (11 lb), covering rodents through rabbits.

This ventilator is furnished with two interchangeable piston/cylinder assemblies of 0.2 to 5.0 ml/breath (5 CC) and 1.5 to 30.0 ml/breath (30 CC). The smaller cylinder is factory installed, and the larger one is supplied separately. The mechanical design features a linkage that enables the piston to touch the end of the cylinder at each stroke regardless of the volume. This assures that all the air/gas mixture taken into the ventilator is expelled at each stroke.

A positive mechanical slide valve is actuated by a cam on the motor to synchronize inspiration and expiration with the motion of the piston. Any NON-EXPLOSIVE gas/air mixture can be used with the ventilator. Connection kit (78-0937) is suggested for this connection. A three-way slide valve is attached to the mechanism plate at the end of the cylinder and actuated by a cam on the motor.

This ventilator CANNOT BE USED IN AN EXPLOSIVE ATMOSPHERE OR WITH ANY EXPLOSIVE GASES SUCH AS ETHER OR CYCLOPROPANE, ETC. THESE VENTILATORS ARE NOT EXPLOSION-PROOF.

KDS35 Small Animal Ventilator

Voltage Selection


- 1) The 'KDS35' Ventilator is equipped with a 115/230V selector feature, which is incorporated into the power entry module located on the side of the unit (see figure below). The Ventilator is shipped from the factory with the voltage selector set at 115V. To change voltage from 115V to 230V, follow these steps:
 - a) Turn off main power, and disconnect power cord.
 - b) Using a small flathead screwdriver, pry open the fuse access door on the power entry module (see figure below).
 - c) Remove fuse holder by using the small flathead screwdriver.
 - d) Rotate the red fuse holder 180 deg. and install back into power entry module.
 - e) Close small access door firmly until fully closed. In power entry module window (see detail below) 230V should be seen.
 - f) Re-install power cord.
- 2) Turn on main power switch located on the rear panel. The display will now illuminate. The number displayed is the current strokes per minute rate.

Fuse Replacement

Key Information

- 1) Make sure power cord is disconnected from the main power supply before servicing the fuse.
- 2) Use only Type 3AG, 1/4" x 1-1/4", 0.75 amp, 250 volts; Type T (time delay) fuses (KD Scientific part # 5153-554 or equivalent).

Turn off power and remove power cord from power module. Use a straight blade screwdriver to pry open the access door. Remove the fuse holder and then remove the fuses from this holder as shown in the figure below. Replace fuses using only 0.75A amp, 250 volts, Type T fuses. Then replace the fuse holder.



KDS35 Small Animal Ventilator

Pump Operation:

CONNECTION TO THE ANIMAL

The vertical slide valve has four ports (see figure below):


- The top port (furthest from the chassis in new orientation) is the return channel from the animal, which leaves the Pump at the port on the side. It is at this top port that expired air can be collected or recycled.
- The middle port is the room air/gas channel to the animal.
- The bottom port (closest to the chassis in new orientation) is the room air or gas mixture intake.

NOTE: It is important that the tubes to and from the animal be as short as possible. The new orientation allows the animal to be as close as possible to the Ventilator for minimal dead space in the system.

The actual connection to the animal can be either:

- Tracheal cannula with y adapter
- Intubation cannula with y adapter

A selection of tracheal cannula and endotracheal tubes are listed in the KD Scientific Catalog and on our website www.kdscientific.com. It is important that the volume of tubing shared by inspired and expired air (dead space) be kept to a minimum by having the "Y" or "T" connector as close to the animal as possible.




KDS35 Small Animal Ventilator

Pump Operation: (continued)

SELECTION OF RATE AND VOLUME

Refer to the Ventilator Graph (below) to select the appropriate tidal volume and respiration rate as a function of animal body weight. As a rough guide an adult rat will weigh between 200 and 500 grams with tidal volumes between 0.6 to 1.25 cm³ with respiration rates between 66 to 210 breaths/min. The average rate is 100 breaths/min. For most rodents the 5 cc cylinder will be adequate. For larger animals it will be necessary to change the cylinder.



Ventilator Graph: Tidal Volume vs. Weight & Rate for Laboratory Mammals in Resting State.
(Apparatus dead space must be added)

KDS35 Small Animal Ventilator

Pump Operation: (continued)

SETTING VOLUME

The volume delivered is read by observing the excursion of the O-ring, on the piston, which is closest to the valve block. This piston O-ring excursion is relative to the volume scale on the outside of the cylinder wall. Volume change is accomplished while the Ventilator is running by turning "volume" knob in the desired direction.

SETTING RATE

Adjust the "Rate" control until the desired rate appears on the LED display. Avoid running the Ventilator at rates exceeding 150 breaths/min for periods in excess of 2 hours.

CAUTION - DO NOT OVER INFLATE THE ANIMAL. A GOOD RULE OF THUMB IS TO OBSERVE THE DEPTH OF RESPIRATION IN THE NORMAL CALM UNANESTHETIZED ANIMAL AND ADJUST THE VOLUME TO APPROXIMATE THIS RESPIRATION DEPTH WHEN THE ANIMAL IS ATTACHED TO THE VENTILATOR.

USE OF GAS MIXTURE

It is possible to introduce non-explosive air/gas mixtures into the lower room air port. However, the source of gas mixture is usually from a tank delivering gas at a constant rate while the Ventilator takes in air/gas intermittently. Some sort of buffer/expansion chamber must be provided between the tank and the Ventilator. KD Scientific offers a connection kit (73-0397) to depressurize the gas prior to entering the ventilator. **This Ventilator is not designed to take air/gas inputs under pressure.**

CREATING A POSITIVE END-EXPIRATORY PRESSURE (PEEP)

To create a PEEP pressure when using the ventilator, simply connect a length of tubing to the "Exhaust" port and submerge the other end in a column of water. There are two locations in which the clip and cylinder can be installed - depending on the chosen orientation. The depth of the end of the tubing (in cm) below the water surface will give the PEEP pressure in cmH₂O.

KDS35 Small Animal Ventilator

Cleaning & Maintenance

CLEANING


To clean the exterior surfaces, use a lint-free cloth to remove loose dust. Use care to avoid scratching the clear display window. For more efficient cleaning, use a soft cloth dampened with water or an aqueous solution of 75% isopropyl alcohol.

MAINTENANCE

No special maintenance is required other than periodic lubrication of the piston O-rings and the slide valve with the special non-toxic grease provided. To get to the slide valve for lubrication:

- 1) Remove the four screws holding the valve to the 1/4" aluminum black back plate.
- 2) Remove the screws in front that hold the cylinder and valve together.
- 3) Lift entire valve assembly up and out of Ventilator.
- 4) Pull stainless steel valve rod free of valve block.
- 5) Lubricate sparingly.

KDS35 RODENT VENTILATOR REPLACEMENT PARTS



LUBRICATION

This KDS35 Ventilator is shipped with the appropriate grease to lubricate the piston and O-rings in the cylinders and the inside of the slide valve.

KDS35 Small Animal Ventilator

Changing Piston and Cylinder


This Ventilator is furnished with a 5 cc piston and cylinder installed. A 30 cc cylinder and piston is also supplied for users that require larger tidal volumes.

NOTE: It is extremely important when changing piston and cylinders that they be aligned to be concentric and in line with the piston shaft.

The cylinder is fastened to the 4 port vertical slide valve by means of two screws located between the "room air inlet" port and the "to animal" port on the front of the slide valve.

For ease of mounting, there is an alignment bracket located below the valve body. When re-attaching a valve to the KDS35 Ventilator, the cut-out in the valve should sit square and parallel with the corner of the alignment bracket. This should ensure that the valve and piston shaft are in line with each other.

The vertical slide valve is fastened to the black 1/4" aluminum sub base by four screws located on the opposite side from the printed label. After mounting a valve; if coupling does not engage piston shaft, loosen cylinder mounting screws and adjust cylinder position until piston shaft enters coupling with ease.



KDS35 Small Animal Ventilator

Changing Piston and Cylinder (continued)

TO REMOVE THE CYLINDER AND PISTON

- 1) Set the Ventilator volume to maximum.
- 2) Run Ventilator at lowest speed and stop when the helical coupling is exposed.
- 3) Using a 0.050 Allen wrench (provided), disconnect the white piston from the coupling, leaving the coupling on the piston shaft.
- 4) Remove the two screws located in the front of the slide valve.
- 5) Remove piston and cylinder.
- 6) Be careful to save the O-ring that is recessed in the slide valve block between the screws. The O-ring provides the airtight seal between the plastic cylinder and the valve block.


TO INSTALL THE PISTON AND CYLINDER

- 1) Pick out the new piston and cylinder.
- 2) Using pliers, remove the piston and lubricate with the grease provided.
- 3) Insert the lubricated piston in the cylinder so that the rear of the piston is flush with the rear of the cylinder with the 1/8" diameter coupling stud protruding.
- 4) Fasten the cylinder to the valve block with the two screws, making sure the O-ring is in place. Do not tighten these two screws completely, but just enough to hold the cylinder in place.
- 5) The object of alignment is to have the stud at the end of the piston be exactly in line with its mating hole in the helical coupling. This insures that the piston moves freely within the cylinder without undue friction.
- 6) The two screws at the front of the valve block control left to right alignment while the four screws at the rear of the slide valve control up and down alignment. By adjustment of these 6 screws the proper alignment can be made. Remember that the screws interact with one another so the procedure may have to be repeated one or more times.
- 7) When proper alignment is achieved, the helical coupling should be connected. If alignment is correct, the coupling should slide onto the piston stud without bending of the piston rod.

KDS35 Small Animal Ventilator

Repair Kit

REPAIR INSTRUCTIONS FOR THE KDS35 RESPIRATOR





Declaration of Conformity

Manufacturer: *KD Scientific, Inc.
84 October Hill Rd.
Holliston, MA 01746-1388 U.S.A.
Phone: (508) 429-6809*

We herewith declare that the following product:

Product Name: *KDS35 Small Animal Ventilator*

Model Numbers: *78-0391*

To which this declaration relates, is in conformity with the applicable EC Directives, harmonized standards, and other normative requirements:

**Application of
Council Directive(s):** *73/23/EEC, 89/336/EEC*

**Standard(s) to which
Conformity is Declared:**

Safety: *EN 61010-1, 2nd Edition (2001)*

Emissions/Immunity: *EN 61326:1997 w/ A1:1998
EN 61000-4-2:1995
EN 61000-4-3:1996
EN 61000-4-4:1995
EN 61000-4-5:1995
EN 61000-4-6:1996
EN 61000-4-11:1994
EN 61000-3-2:2000 w/ A.14
EN 61000-3-3:1997*

Signature: 
Beth Bauman, V.P. of Operations

Date: 8/1/06