**ROTINA 420
ROTINA 420 R**

(DE)	Bedienungsanleitung	5
(EN)	Operating Instructions.....	31
(FR)	Mode d'emploi	56
(IT)	Istruzioni per l'uso	83


Fig. 1


Fig. 2 ROTINA 420


Fig. 3 ROTINA 420 R

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

Das bezeichnete Gerät, inklusive Zubehör entspricht den aufgeführten EG-Richtlinien und Normen.
The denoted device, including accessories corresponds to the listed EC guidelines and standards.
L'appareil désigné, y compris les accessoires, correspond aux directives CE et aux normes énumérées.
L'apparecchio designato, compresi gli accessori, è conforme alle direttive CE e alle norme citate.

Geräteart, Type of device, Type d'appareil, Tipo di apparecchio:

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Contents

1	Use according to specification	33
2	Residual risks	33
3	Technical specifications	33
4	Notes on safety	34
5	Symbol meanings	35
6	Delivery checklist	36
7	Unpacking the centrifuge	36
8	Initial operation	36
9	Interface (only for centrifuge with interface)	36
10	Opening and closing the lid	37
10.1	Opening the lid	37
10.2	Closing the lid	37
11	Installation and removal of the rotor	37
12	Loading the rotor	38
13	Aerosol tight sealing of angle rotors	38
14	Control and display elements	39
14.1	Control knob	39
14.2	Keys and setting options	39
15	Enter centrifugation parameters	40
15.1	Running time	40
15.2	Speed (RPM)	40
15.3	Relative centrifugal force (RCF) and centrifuging radius (RAD)	40
15.3.1	Centrifuging radius (RAD)	40
15.3.2	Relative centrifugal force (RCF)	41
15.4	Run-up and run-down parameters	41
15.4.1	Run-up step	41
15.4.2	Braking step	41
15.4.3	Brake switch-off speed	41
15.5	Temperature (only for centrifuge with refrigeration)	41
16	Programming	41
16.1	Inputting or changing programs	41
16.2	Calling up programs	42
16.3	Write protection for programs	42
16.4	Program linkage	42
16.4.1	Activating or deactivating program linkage	42
16.4.2	Linking programs or changing a program linkage	43
16.4.3	Calling up program linkage	43
16.5	Automatic temporary storage	43
17	Centrifugation	44
17.1	Centrifugation with preselected time	44
17.2	Continuous operation	44
17.3	Short-time centrifugation	45
18	Change the settings during the centrifugation run	45

19	Emergency Stop.....	45
20	Acoustic signal	45
21	Displayed centrifugation data after switch-on.....	46
22	Set temperature unit (only for centrifuge with refrigeration).....	46
23	Setting program locks.....	46
24	PIN (Personal Identification Number)	47
24.1	Setting or changing the PIN	47
24.2	Changing a forgotten PIN	47
25	Centrifuge address.....	47
26	RCF integral and B ramp functions	48
27	Query of operating hours and the number of centrifugation runs	48
28	Querying system information.....	48
29	Immediate display of centrifugation data after switch-on.....	48
30	Cooling (only in centrifuges with cooling)	48
30.1	Standby-cooling	48
30.2	Precooling of the rotor.....	49
31	Heating (only on centrifuges with heating/cooling option)	49
32	Relative centrifugal force (RCF)	49
33	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³	50
34	Rotor Identification	50
35	Emergency release	50
36	Maintenance and servicing.....	51
36.1	Centrifuge (housing, lid and centrifuging chamber).....	51
36.1.1	Surface cleaning and care	51
36.1.2	Surface disinfection	51
36.1.3	Removal of radioactive contaminants	51
36.2	Rotors and Attachments	52
36.2.1	Cleaning and care.....	52
36.2.2	Disinfection	52
36.2.3	Removal of radioactive contaminants	52
36.2.4	Trunnions.....	52
36.2.5	Rotors and accessories with limited service lives	53
36.3	Autoclaving	53
36.4	Centrifuge containers.....	53
37	Faults	54
38	Switch on automatic circuit breaker.....	55
39	Acceptance of the centrifuges for repair	55
40	Disposal	55
41	Anhang / Appendix.....	110
41.1	Rotoren und Zubehör / Rotors and accessories.....	110

1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm³. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen				
Model	ROTINA 420		ROTINA 420 R		
Type	4701	4701-01	4706, 4706-20, 4706-50	4706-01	
Mains voltage (± 10%)	200 – 240 V 1~	100 – 127 V 1~	200 – 240 V 1~	100–127V 1~	100 V 1~
Mains frequency	50 – 60 Hz	50 – 60 Hz	50 Hz	60 Hz	50 Hz
Connected load	750 VA	800 VA	1500 VA	1750 VA	
Current consumption	3.8 A	8.0 A	7.0 A	16.0 A	
Cooling medium	----		R 404A		
Max. capacity	4 x 600 ml				
Allowed density	1.2 kg/dm ³				
Speed (RPM)	15000				
Force (RCF)	24400				
Kinetic energy	24000 Nm				
Obligatory inspection (BGR 500)	yes				
Ambient conditions (EN / IEC 61010-1) – Set-up site – Altitude – Ambient temperature – Humidity – Excess-voltage category (IEC 60364-4-443) – Pollution degree	Indoors only Up to 2000 m above sea level 2°C to 35°C 5°C to 35°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C. II 2				
Device protection class	I				
Not suitable for use in explosion-endangered areas.					
EMC – Emitted interference, Interference immunity	EN / IEC 61326-1, Class B	FCC Class B	EN / IEC 61326-1, Class B	FCC Class B	
Noise level (dependent on rotor)	≤ 63 dB(A)		≤ 56 dB(A)		
Dimensions – Width – Depth – Height	506 mm 650 mm 423 mm		713 mm 654 mm 423 mm		
Weight	approx. 75 kg	approx. 84 kg	approx. 107.5 kg, 109 kg (4706-50)	approx. 117 kg	

4 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
- With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

- Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
 For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").
 In a biosafety system, droplets and aerosols are prevented from escaping by a bioseal (packing ring).
 If the hanger of a biosafety system is used without the lid, the packing ring must be removed from the hanger in order to prevent the packing ring from being damaged during the centrifugation run. Damaged packing rings must not be used to seal the biosafety system.
 Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 610101-2-020 standard.

For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.

- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
 - prescribed tests to BGV A1, BGR 500 are carried out by an expert.

5 Symbol meanings

Symbol on the machine:

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in this document:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the machine and in this document:

Beware of hot surface.

Nonobservance of this warning can lead to material damage and personal injury.

Symbol on the machine and in this document:

Plastic suspension can only be used in temperatures up to a maximum of 40°C / 104°F.

Nonobservance of this warning can lead to material damage and personal injury.

Symbol on the machine:

Equipotential: Connector (PE connector) for potential equalization (only for centrifuge with PE connector).

Symbol on the machine:

RS232 interface (only for centrifuges with RS232 interface).

Symbol on the machine:

Automatic circuit breaker (only for types 4701-01 and 4706-01).

Symbol in this document:

This symbol refers to important circumstances.

Symbol on the machine and in this document:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

6 Delivery checklist

- 1 Connecting cable
- 1 Hex. pin driver 2,5 mm
- 1 Hex. pin driver 5 mm
- 1 Lubricating grease for trunnions
- 1 Operating instructions
- 1 Notes on moving the equipment safely

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

7 Unpacking the centrifuge

- Lift the carton upward and remove the padding.

Do not lift by the front panel.

Observe the weight of the centrifuge, refer to chapter "Technical specifications".

Lift the centrifuge on both sides with an appropriate number of helpers and place it on the laboratory table.

8 Initial operation

- According to the laboratory instrument standards EN / IEC 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Do not place any object in front of the ventiduct.
Keep a ventilation area of 300 mm around the ventiduct.
- In the case of the centrifuge, type 4706-20, connect the nitrogen supply according to the enclosed instruction sheet AH4706-20XX.

It is mandatory that the centrifuge be connected according to the enclosed instruction sheet.
Make sure you observe the enclosed instruction sheet.

- Centrifuge with PE connector:
If required, connect the PE connector at the back of the device with an additional medical potential equalization system.
- Centrifuge with RS232 interface:
Connect the RS232 centrifuge interface with an RS232 connecting cable (not included in the scope of delivery) to the PC.
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Switch on the mains switch (switch position "I"). The LEDs in the keys are blinking.
The following displays appear one after the other:
 1. The centrifuge model
 2. The program version
 3. The rotorcode (rotor), the maximum rotor speed (Nmax) and the preset centrifuging radius (R) of the last rotor that was recognized by rotor recognition.
 4. **⚡ OPEN ⚡ OEFFNEN**
- Open the lid.
The centrifuge data of the last used program or program 1 will be displayed.
- Remove the transport safety device (see instruction sheet on "Moving the equipment safely").

9 Interface (only for centrifuge with interface)

As an option, the device can be equipped with an RS232 interface.

The RS232 interface is labeled with the  symbol.

The centrifuge can be controlled and data queried via this interface.

The LED in the  key lights up during data communication.

10 Opening and closing the lid

10.1 Opening the lid

The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Press the button **OPEN/STOP**. The lid unlocks via the motor and the left LED in the pushbutton **OPEN/STOP** extinguishes.

10.2 Closing the lid


Do not put your fingers between lid and housing.
Do not bang the lid shut.

If the left LED in the **OPEN/STOP** button flashes, press the **OPEN/STOP** button so that the motor-driven lid lock goes into the basic position (opened).

- Place the lid and lightly press down the front edge of the lid. The locking action is effected by motor. The left LED in the button **OPEN/STOP** lights up.

11 Installation and removal of the rotor


- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

12 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position. Hangers that are marked with a set number (e.g. S001/4) may only be used in the set.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded


Not permitted!

Rotor is not evenly loaded

- On certain hangers, the weight of the maximum load or the weight of the maximum load and the maximum weight of the completely loaded hanger is specified. This weight may not be exceeded. In case of exception, see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³". The weight specified for the maximum loading includes the total weight of adapter, frame, centrifuging container and content.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- Always fill the centrifuge containers outside of the centrifuge.
- No liquid should be allowed to enter the centrifugal chamber during filling and swinging out of the hangers.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.

13 Aerosol tight sealing of angle rotors

To ensure aerosol sealing the lid of an aerosol sealed angular rotor must be firmly closed.

This is performed with the help of the supplied spanner that is inserted through the borehole in the rotary handle. If the rotary handle does not have a borehole then the lid must be firmly closed by turning in clockwise direction.

For available aerosol tight angle rotors, please see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

14 Control and display elements

See figure on page 2.

Fig. 2, Fig. 3: Display and control panel

14.1 Control knob

For setting the individual parameters.
Turning anticlockwise reduces the value. Turning clockwise increases the value.

14.2 Keys and setting options

- Running time, parameters **t/hms**.
h: hours. Adjustable from 1 h to 99 h, in 1 hour increments.
m: minutes. Adjustable from 1 min to 59 min, in 1 minute increments.
s: seconds. Adjustable from 1 s to 59 s, in 1 second increments.
- Continual running " ∞ "

- Revolution, parameters **rpm**.
Adjustable from 50 rpm to a maximum rotor speed of the rotor (Nmax), in increments of 10. For maximum rotor speed, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

- Relative centrifugal force, parameter **RCF**.
The RCF is shown in brackets $\rangle \langle$. The LED in the button lights up.
The figure that is adjustable is the figure that produces a revolution speed between 50 rpm and the maximum rotor speed (Nmax). Adjustable in increments of 1.
- Centrifuging radius, parameters **RAD**.
Adjustable from 85 mm to 330 mm in 1 millimeter increments. For centrifuging radius, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

- Run-in and run-out parameters
- Run-in steps, parameters $\swarrow$.
Step 9 = shortest run-in time, ... Step 1 = longest run-in time.
- Braking step, parameters $\searrow$.
Step 9 = shortest run-out time, ... Step 1 = long run-out time, step 0 = non-braked run-out.
- Brake switch-off revolution, parameters **N Brake**.
Adjustable from 50 rpm to the maximal rotor speed (Nmax) in increments of 10.
Once the rotor speed has been reached, non-braked run-out begins.

- Temperature (only for centrifuge with refrigeration)
Adjustable in degrees Celsius (°C) or in degrees Fahrenheit (°F). For setting the temperature unit, see chapter "Temperatur-Einheit einstellen/Setting the temperature unit".
Parameters **T/°C** = degrees Celsius (°C).
Adjustable from -20°C to +40°C, in 1°C increments (for the heat/cool option, adjustable in the range from -20°C to +90°C).
Parameters **T/°F** = degrees Fahrenheit (°F).
Adjustable from 4°F to +104°F, in 1°F increments (for the heat/cool option, adjustable in the range from -4°F to +194°F).
The lowest temperature that can be reached depends on the rotor (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").

- Call up programs and program linkages; parameter **RCL**.
Programs: Program locations 1 to 99 and PREC. program linkages: Program locations A to Z.
- Save programs and program linkages; parameter **STO**.
99 programs can be saved (for centrifuges without refrigeration: program places 1 to 99, for centrifuges with refrigeration: program places 1 98 and PREC). The program place PREC (PRECOOLING) is reserved for the precooling program. Program place 0 is used for temporary storage for centrifuge data from the last centrifugation run. No programs can be saved at this program place.
It is possible to save 25 program linkages (program locations A to Z; there is no program location J). A program linkage can consist of 20 programs.
- Link programs; parameter **EDIT**.
- Start "Machine Menu" (keep the key held down for 8 seconds) and select the parameters from this menu.

- Start rotor precooling (only for centrifuges with refrigeration). The LED in the key is lit during the centrifugation run, as long as the rotor is revolving. Precooling of the rotor takes place automatically with the **PREC** (PRECOOLING) program.

- Start centrifugation run. The LED in the key is lit during the centrifugation run, as long as the rotor is revolving.
- Short-time centrifugation. The centrifugation run occurs while the key is kept pressed. The LED in the key is lit during the centrifugation run, as long as the rotor is revolving.
- Saving entries and changes.

- End centrifugation run. The rotor decelerates with the preset run-down parameters. The right-hand LED in the button lights up until the rotor is stationary. Once the rotor is stationary the left-hand LED flashes in the button. Pressing the button twice triggers the EMERGENCY STOP.
- Unlock the lid. The left-hand LED in the button goes out.
- Leave parameter input and the "Machine Menu".

15 Enter centrifugation parameters

If no key is pressed for 8 seconds long after the selection or during the input of parameters, the previous values will be shown in the display. The input of parameter then has to be executed again.

When you enter several parameters, the **START** key does not have to be pressed until you have made the settings for the last parameter.

If parameters are changed, the program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.

You can no discontinue entering parameters at any time by pressing the key **OPEN / STOP**. In this case the adjustments are not saved.

15.1 Running time

In order to set continual running, the minutes, seconds and hours have to be set to zero. Continual running is represented in the display by the following symbol, "∞".

- Press the **TIME** key. The parameters **t/hms** is displayed. The minutes (**m**) are shown in parentheses < >, and can be changed.
- Use the adjusting knob $\odot$ to set the value you want.
- Press the **TIME** key. The seconds (**s**) are shown in parentheses < > and can be changed.
- Use the adjusting knob $\odot$ to set the value you want.
- Press the **TIME** key. The hours (**h**) are shown in parentheses < > and can be changed.
- Use the adjusting knob $\odot$ to set the value you want.
- Press the **START** key, to have the settings applied to the display.

15.2 Speed (RPM)

- Press the **RPM** key. The parameters **RPM** is displayed.
- Use the adjusting knob $\odot$ to set the value you want.
- Press the **START** key, to have the settings applied to the display.

15.3 Relative centrifugal force (RCF) and centrifuging radius (RAD)

15.3.1 Centrifuging radius (RAD)

- Keep pressing the **RCF** key until the parameters **RAD** and **RCF** are displayed, e.g. **RAD = 146 -> RCF = 3695**. The LED is lit in the key.
- Use the adjusting knob $\odot$ to set the centrifuging radius you want. By changing the centrifuging radius, the value adjusts automatically to the RCF.
- Press the **START** key, to have the settings applied to the display.

15.3.2 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) depends on the radius of the centrifuge. Before setting the RCF, the radius of the centrifuge must be set. See chapter "Centrifuge radius (RAD)".

- Keep pressing the key until the parameters **RCF** and **RAD** are displayed, e.g. **RCF = 3695 RAD = 146**. The LED is lit in the key.
- Use the adjusting knob to set the RCF you want.
- Press the key, to have the settings applied to the display.

15.4 Run-up and run-down parameters

15.4.1 Run-up step

- Press the key as often as necessary until the parameter is shown.
- Use the adjusting knob to set the level you want.
- Press the key, to have the settings applied to the display.

15.4.2 Braking step

- Press the key as often as necessary until the parameter is shown.
- Use the adjusting knob to set the level you want.
- Press the key, to have the settings applied to the display.

15.4.3 Brake switch-off speed

- Press the key as often as necessary until the parameter **N Brake** is shown.
- Use the adjusting knob to set the value you want.
- Press the key, to have the settings applied to the display.

15.5 Temperature (only for centrifuge with refrigeration)

The temperature can be entered in degrees Celsius (°C) or in degrees Fahrenheit (°F). For setting the temperature unit, see chapter "Set temperature unit".

- Press the key. The following parameters are displayed: **T/°C** or **T/°F**.
- Use the adjusting knob to set the value you want.
- Press the key, to have the settings applied to the display.

16 Programming

If parameters are changed, the program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.

16.1 Inputting or changing programs

- Set the parameters you want (see chapter "Enter centrifugation parameters").
- Press the key as often as necessary until the **STO** parameter is shown.
- Use the adjusting knob to set the program place you want.

If a "+" is shown after the program location, then this data is write-protected. In this case, the write protection has to be removed first before saving (see the chapter on "Write protection for programs").

- Press the key to have the settings saved to the program place you want. Briefly, **Program store ...** is displayed as confirmation.

The previous program location data will be overwritten when the new data is saved.

If "**Protected !!**" is displayed, then the data at the program location is write-protected and it will not be saved.

16.2 Calling up programs

- Press the **◀PROG** key. The parameters **RCL** is displayed.
- Use the adjusting knob **○** to set the program place you want.

If a "+" is shown after the program location, then this data is write-protected.

- Press the **◀START** key. Briefly, **Program recall ...** is displayed as confirmation. The centrifugation data of the selected program place is shown.

16.3 Write protection for programs

The programs can be protected against unintentional modification.

The write protection can be activated/deactivated as follows when the rotor is at a standstill:

- Call up the desired program (see the chapter on "Calling up programs").
- Press the **◀PROG** key. The parameters **RCL** is displayed.
- Press and hold the **◀PROG** button for eight seconds.
After eight seconds, **Set Protection = 1-** (for example) appears in the display.
- Set **○ "+"** or **○ "-"** with the adjusting knob.
 - +** = program is write-protected,
 - = program is not write-protected.
- Press the **◀START** button to save the setting.

16.4 Program linkage

"Program linkage" can be used to link several centrifuge operations together.

A program linkage is only possible if this has been activated (parameter **Multi programs = on**; see the chapter on "Activating or deactivating program linkage").

16.4.1 Activating or deactivating program linkage

The program link can be activated/deactivated as follows when the rotor is at a standstill:

- Press and hold the **◀PROG** button for eight seconds.
After eight seconds, ***** Machine Menu ***** appears in the display.
- Press the **◀PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **◀START** key. The **SOUND / BELL = off/on** is displayed.
- Press the **◀PROG** key as often as necessary until the **Multi programs = off/on** is displayed.
- Set **off** or **on** with the adjusting knob **○**.
 - off** = program linkage deactivated,
 - on** = program linkage activated.
- Press the **◀START** key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **◀OPEN/STOP** key once to leave the "Settings" menu or press the **◀OPEN/STOP** key twice to leave the "Machine Menu".

16.4.2 Linking programs or changing a program linkage

It is possible to save 25 program linkages (program locations A to Z; there is no program location J). A program linkage can consist of no more than 20 programs.

The speed adjustment from one program to the next one is always done in a program linkage with the start-up parameter of the next program.

There cannot be any continuous operation programs in a program linkage.

No centrifuge parameters may be modified in a program linkage. A parameter modification is only possible in the individual programs.

The **TIME** button can be used during the centrifuge operation to call up the total run time of the program linkage

(e. g. **Σ=00:05:30**) and the run time of the currently running program (e.g. **t B.02=00:01:00**).

1. Press the **PROG** button repeatedly until the parameter **EDIT A...Z** is displayed.
2. Use the **◂** adjusting knob to set the desired program location where the program linkage is to be saved.
3. Press the **START** button. The program location of the program linkage and the first program of the program linkage will be shown (e.g. **EDIT B.01 = 01**).
4. Use the **◂** adjusting knob to set the first program of the program linkage.
5. Press the **PROG** button. The next program of the program linkage will be shown (e.g. **EDIT B.02 = END**).
6. Use the **◂** adjusting knob to set the next program of the program linkage.
7. Press the **PROG** button. The next program of the program linkage will be shown (e.g. **EDIT B.03 = END**).
8. Repeat steps 6 and 7 until all programs have been set.
9. Set **END** with the adjusting knob (turn adjusting knob anti-clockwise).

For program linkages which consist of 20 programs, **END** cannot be set after the 20th program.

10. Press the **START** button. **STO B**, for example, is displayed.
 11. Press the **START** button to save the program linkage.
- As confirmation, **Multi program store ..** will be shown briefly.

16.4.3 Calling up program linkage

- Press the **PROG** button repeatedly until the parameter **RCL A...Z** is displayed.
- Use the adjusting knob **◂** to set the program place you want.
Press the **START** key. Briefly, **Multi program recall ...** is displayed as confirmation.
- The centrifugation data of the first program of the program linkage will be shown.

16.5 Automatic temporary storage

The program place 0 serves as temporary storage for centrifugation data of the last centrifugation run that took place.

No programs can be saved at this program place.

After every start to a centrifugation run, the centrifugation data that is used for the run is automatically saved at program place "0" and can be accessed there.

17 Centrifugation

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

For centrifuges with the heating / cooling option, after a centrifugation run with a very high temperature (e.g. +90°C), it is necessary to wait until the cover of the centrifuge has cooled down to ambient temperature before carrying out a centrifugation run with cooling. Failure to comply with this rule may result in cracks in the cover.

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, and **IMBALANCE** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor (Nmax), it will not be possible to start a centrifugation run. The message **N > ROTOR MAX** will be displayed (see Chapter "Faults").

If **N > ROTOR MAX** in Prog : e.g. **5**, Runtime **00:00** in Prog : e.g. **5** or **Empty Program** is shown during program linkage, no centrifuge operation can be started (see the chapter on "Faults").

A centrifugation run can be stopped at any time by pushing the key **OPEN / STOP**.

All parameters can be selected and altered during the centrifugation run (see Chapter "Enter centrifugation parameters").

It is possible to switch over between RPM and RCF display at any time using the **RPM** and **RCF** keys. Switching is not possible if program linkages are being used. When working with the RCF display, it is necessary to enter the centrifuging radius.

If **OPEN** **OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If the rotor has been changed, no centrifugation run takes place and this display is shown, e.g. **Rotor 4 Nmax= 4500 R=184 mm** (see chapter "Rotor Identification").

Operation errors and malfunctions will be shown (see Chapter "Faults").

- Switch on the mains supply switch. Switch position **I**.
- Load the rotor and close the centrifuge lid.

17.1 Centrifugation with preselected time

- Set a run time, call up a program with preset time or a program linkage (see the chapters on "Enter centrifugation parameters", "Calling up programs" or "Program linkage").
- Press the **START** key. The LED in the **START** key blinks until the rotor has been imported, it is subsequently lit.
- Once the time has run up or if the centrifugation run is canceled by pressing the **OPEN / STOP** key, run-down takes place with the selected run-down parameter. The run-down parameter is displayed e.g. **9**. The right-hand LED in the **OPEN / STOP** key is lit. Once the rotor has come to a standstill, the LED in the **START** key and the right LED in the **OPEN / STOP** key switch off and the left LED in the **OPEN / STOP** key lights up.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the remaining time will be displayed.

17.2 Continuous operation

- Set the minutes, seconds and hours to "0" or start a continual run program (see chapter "Enter centrifugation parameter" or "Calling up programs").
- Press the **START** key. The LED in the **START** key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- Press the **OPEN / STOP** key to end the centrifugation run. Run-down is carried out with the selected run-down parameters. The run-down parameter is displayed e.g. **9**. The right-hand LED in the **OPEN / STOP** key is lit. Once the rotor has come to a standstill, the LED in the **START** key and the right LED in the **OPEN / STOP** key switch off and the left LED in the **OPEN / STOP** key lights up.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the time run will be displayed.

17.3 Short-time centrifugation

A short-time centrifugation is not possible if program linkages are being used.

- Keep the **START** key pressed. The LED in the **START** key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- Release the **START** key again to end the centrifugation run. Run-down is carried out with the selected run-down parameters. The run-down parameter is displayed e.g. 9. The right-hand LED in the **OPEN/STOP** key is lit. Once the rotor has come to a standstill, the LED in the **START** key and the right LED in the **OPEN/STOP** key switch off and the left LED in the **OPEN/STOP** key lights up.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the time run will be displayed.

18 Change the settings during the centrifugation run

It is not possible to modify settings during centrifuge operation if program linkages are being used.

The running time, the speed, the relative centrifugal force (RCF), the centrifuging radius, the run-up and run-down parameters as well as the temperature (for centrifuges with cooling only) can be changed while the centrifuge is running.

- Select the parameter you want and use the $\odot$ adjusting knob to change the value (see chapter "Enter centrifugation parameter").
- Press the **START** key. The changed setting is saved at program place "0" (see chapter "Automatic temporary storage"). The program place number is displayed in parentheses (0). The original program is not overwritten.

19 Emergency Stop

- Press the **OPEN/STOP** key twice.

During a not-stop, run-out is carried out with braking level 9 (shortest run-out). Braking level 9 is displayed. If braking level 0 has been preselected, run-out is carried out with braking level 9d. With braking level 9d, the run-out time is longer, technically speaking, than braking level 9.

20 Acoustic signal

The acoustic signal sounds:

- Upon the appearance of a disturbance in 2 second intervals.
- After completion of a centrifugation run and rotor standstill in 30 second intervals.

The acoustic signal is stopped by opening the lid or pressing any key.

The signal after completion of the centrifugation run can be activated or deactivated in the following manner, if the rotor is at standstill:

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Use the adjusting knob $\odot$ to set **off** or **on**.
- Press the **START** key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN/STOP** key once to leave the "Settings" menu or press the **OPEN/STOP** key twice to leave the "Machine Menu".

21 Displayed centrifugation data after switch-on.

After switch-on the centrifugation data from program 1, or from the last program that was used, is displayed. If the rotor is at standstill, this can be set as follows:

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Press the **PROG** key as often as necessary until **Start program = Last/First** is displayed.
- Use the adjusting knob  to set **Last** or **First**.
Last = the program used last, First = program 1.
- Press the **START** key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN / STOP** key once to leave the "Settings" menu or press the **OPEN / STOP** key twice to leave the "Machine Menu".

22 Set temperature unit (only for centrifuge with refrigeration)

The temperature can be entered in degrees Celsius (°C) or in degrees Fahrenheit (°F).

When the rotor is at standstill, the temperature unit must be set as follows:

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as necessary until the **->Settings** is displayed.
- Press the **START** key. The **SOUND / BELL = off/on** is displayed.
- Press the **PROG** key as often as necessary until **Temp Unit = Celsius/Fahrenheit** is displayed.
- Use the adjusting knob  to set **Celsius** or **Fahrenheit**.
- Press the **START** key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **OPEN / STOP** key once to leave the "Settings" menu or press the **OPEN / STOP** key twice to leave the "Machine Menu".

23 Setting program locks

The following program locks can be set:

- LOCK 1** **LOCK 1** is displayed in the "" field.
Programs can only be called up but not modified.
- LOCK 2** **LOCK 2** is displayed in the "" field.
No programs can be called up and modified.
The centrifuge can be controlled via the interface (only for centrifuge with interface).
- LOCK 3** No status display.
No program locking. Programs can be called up and modified.

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as necessary until the **->Change LOCK** is displayed.
- Press the **START** key. The lock status is displayed.
If no PIN has been entered, **LOCK = {3} confirm by START**, for example, will be displayed.
If a PIN has been entered, **LOCK = 3**, for example, will be displayed.
- Use the  adjusting knob to set the desired lock status.

If a PIN has been entered, **PIN = ---- confirm by START** will be shown. In this case, a valid PIN has to be entered with the  adjusting knob and then the **START** button has to be pressed before the lock status can be set.

- Press the **START** key to save the setting.
Briefly, **Store LOCK 2 ...** is displayed as confirmation, followed by **-> Change LOCK**.
- Press the **OPEN / STOP** key once to leave the "Change LOCK" menu or press the **OPEN / STOP** key twice to leave the "Machine Menu".

24 PIN (Personal Identification Number)

To prevent the program lock from being changed by unauthorized persons, a PIN can be set.

 No PIN is set in the factory.

24.1 Setting or changing the PIN

The PIN can be set as follows when the rotor is at a standstill:

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as necessary until the **->Change PIN** is displayed.
- Press the **START** key. **old PIN = ----** **<START>** is displayed.
- With the **0** adjusting knob, set the valid PIN.

 If the PIN is being set for the first time, then skip this step or set "0000".

Input assistance:

Press and hold the respective button.

1000 : Only the 1000s place of the PIN will be changed.

RCF : Only the 100s place of the PIN will be changed.

RPM : Only the 10s place of the PIN will be changed.

TIME : Only the 1s place of the PIN will be changed.

- Press the **START** key. **new PIN = ----** **<START>** is displayed.

 If an incorrect PIN is set, then **old PIN = ----** **<START>** will be displayed again. In this case, set the valid PIN with the **0** adjusting knob and then press the **START** button.

- With the **0** adjusting knob, set the new PIN.

 To delete the PIN, "0000" must be set.

- Press the **START** key to save the setting.
Briefly, **Store PIN ...** is displayed as confirmation, followed by **-> Change PIN**.
- Press the **OPEN/STOP** key once to leave the "Change PIN" menu or press the **OPEN/STOP** key twice to leave the "Machine Menu".

24.2 Changing a forgotten PIN

If the PIN is forgotten, a so-called help number can be called up. With the help of this number, the manufacturer can calculate a PIN which can be used to change the forgotten PIN. In case this is required, please contact your supplier.

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as necessary until the **->Change PIN** is displayed.
- Press the **START** key. **old PIN = ----** **<START>** is displayed.
- Press the **PROG** key. **Get HELP # no** is displayed.
- Set **yes** with the **0** adjusting knob.
- Press the **START** key. **Are you sure ? no** is displayed
- Set **yes** with the **0** adjusting knob.
- Press the **START** key. A help number is displayed, e.g. **HELP # = 5487**.
- Note this help number and use it to request the needed PIN.
- Set a new PIN with the PIN you receive (see the chapter on "Setting or changing the PIN").

25 Centrifuge address

 Ex works, the address is] = 29. Address is set.

26 RCF integral and B ramp functions

These functions are not activated as standard.

27 Query of operating hours and the number of centrifugation runs

The operating hours are divided up into internal and external operating hours.

Internal operating hours: Total time the device was switched on.

External operating hours: Total time of previous centrifugation runs.

The operating hours and the number of centrifugation runs can be queried as follows when the rotor is at standstill:

- Keep the **[PROG]** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **[PROG]** key as often as necessary until the **-> Operating Time** is displayed.
- Press the **[START]** key. The external operating hours are displayed, e.g. **OP Time ext = 0h25m**.
- Press the **[PROG]** key. The internal operating hours are displayed, e.g. **OP Time int = 1h36m**.
- Press the **[PROG]** key. The number of centrifugation runs is displayed, e.g. **Number of Starts = 10**.
- Press the **[OPEN/STOP]** key to leave the menu, "Operating Time" or press the key **[OPEN/STOP]** three times to leave the "Machine Menu".

28 Querying system information

Rotor information and the program version of the centrifuge as well as the frequency converter can be queried as follows when the rotor is at a standstill:

- Keep the **[PROG]** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **[PROG]** key as often as necessary until the **->Info** is displayed.
- Press the **[START]** key. The rotorcode (rotor), the maximum rotor speed (Nmax) and the preset centrifuging radius (R) of the last rotor that was recognized by the rotor recognition, e.g. **Rotor 4*: Nmax = 4500 R=184**.
The last recognized rotor is marked with a star (*). Use the adjusting knob  to have the information of the rotors permitted in the centrifuge displayed.
- Press the **[PROG]** key. The program version of the centrifuge is displayed, e.g. **SW version = 01.00**.
- Press the **[PROG]** key. The program version of the frequency converter is displayed e.g. **FC-SW-Version = 4**.
- Press the **[OPEN/STOP]** key twice to leave the "Info" menu or press the key **[OPEN/STOP]** three times to leave the "Machine Menu".

29 Immediate display of centrifugation data after switch-on.

- Switch on the power supply. (Switch setting I).
- The first time the display changes optically, (inverse display) press any key and keep it held down. The centrifugation data is displayed immediately.

30 Cooling (only in centrifuges with cooling)

The temperature set-point can be adjusted from -20°C to +40°C / -4°F to +104°F. On centrifuges with a heating/cooling option, the temperature set-point is settable from -20°C to +90°C / -4°F to +194°F. The lowest obtainable temperature is dependent on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").

30.1 Standby-cooling

When the rotor is at a standstill and when the lid is closed, the centrifugal chamber is cooled down to the preselected temperature if this is lower than 20°C / 68°F.

During the standby cooling, the preset temperature is displayed.

30.2 Precooling of the rotor

For quick precooling of the unloaded rotor and the accessories, it is a good idea to do a centrifugation run with the continuous operation setting and a speed of about 20% of the maximum rotor speed.

The precooling is done automatically with the **PREC** (PRECOOLING) program.

Precooling is not possible if program linkages are being used.

- Press the key. The LED in the key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- Press the key to end the centrifugation run. Run- down is carried out with the selected run- down parameters. The run-down parameter is displayed e.g. 9. The right-hand LED in the key is lit. Once the rotor has come to a standstill, the LED in the key switches off and **OPEN** **OEFFNEN** is displayed. The right-hand LED in the key also switches off, the left LED in the key starts blinking and continues blinking until the lid is opened.

During the centrifugation run, the speed of the rotor, or the resulting RCF value, the sample temperature and the remaining time or expired time are displayed.

31 Heating (only on centrifuges with heating/cooling option)

During operation of the centrifuge the centrifuge chamber can be heated as required to the preselected temperature. The heating is switched off when the rotor is stationary.

Risk of burning! The surface temperature of the heating element in the centrifugal chamber of the centrifuge can reach 500°C / 932°F. Do not touch the heating element.

Plastic suspension can only be used in temperatures up to a maximum of 40°C / 104°F.

32 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$RCF = \left(\frac{RPM}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad RPM = \sqrt{\frac{RCF}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.
For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

33 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

34 Rotor Identification

After every start of a centrifugation run the rotor utilised is identified.

After changing a rotor, the engine switches off and the rotorcode (rotor), the maximum rotor speed of the rotor (Nmax) and the preset centrifuging radius (R) of the rotor being used are displayed, e.g. **rotor 4 Nmax= 4500 R=184 mm**.

- Press the **OPEN/STOP** key to open the lid, or press the **START** key to start the centrifugation run. For centrifuges with refrigeration, precooling of the rotor can also be started by pressing the **PRE** key.

If, following a rotor change, the maximum speed of the rotor is less than the set speed, the speed is limited to the maximum speed of the rotor. In this case, the program place number is displayed in parentheses ().

35 Emergency release

During a power failure the lid cannot be unlocked by motor. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.
Open the lid only during rotor standstill.

See figure on page 2.

- Switch off the mains switch (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the hexagonal wrench key into the bore hole (Fig. 1, A) and carefully rotate by half a turn in clockwise direction until the lid can be opened.
- Pull the hexagon socket head wrench out of the drilling again.
- If the left LED in the **OPEN/STOP** button flashes after the centrifuge is switched on again, press the **OPEN/STOP** button so that the motor-driven lid lock goes into the basic position (opened) again.

36 Maintenance and servicing

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

36.1 Centrifuge (housing, lid and centrifuging chamber)

36.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

36.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants: ethanol, n-propanol, isopropyl alcohol, glutardialdehyde, quaternary ammonium compounds.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

36.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

36.2 Rotors and Attachments

36.2.1 Cleaning and care

- In order to prevent corrosion and material changes, rotors and accessories must be cleaned regularly with soap or a mild detergent and a damp cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents:
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- In the case of biosafety systems (for further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"), the packing rings must be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear. In order to prevent the packing ring from twisting when opening and closing the cover, the packing ring must be lightly rubbed with talcum powder or a rubber care product.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and the accessory parts have to be checked on a monthly basis for wear and corrosion. For swing-out rotors, it is important to check the area of the lifting lugs, for hangers, the grooves and the base should be checked for cracks.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

36.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants:
glutaraldehyde, propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

36.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

36.2.4 Trunnions

With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

36.2.5 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

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For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

The number of centrifuge operations can be queried; see the chapter on "Querying the operating hours and the number of centrifuge operations".

36.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium, suspension made of metal, lids with bio-containment as well as stands and reductions can be autoclaved at 121° C / 250°F (20 mins).

Otherwise you must ask the manufacturer.

No statement can be made about the degree of sterility.

The lids of the rotors and containers must be removed prior to autoclaving.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

We recommend that the packing rings of the bio-safety system be replaced after autoclaving.

36.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splitters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

37 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

Perform a MAINS RESET:

- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "I").


Message / fault		Cause	Remedy
No display	---	No voltage. Overvoltage protection tripped out. Triggering the automatic circuit breaker (only for types 5650-01, 5660-01 and 5660-51).	– Check supply voltage. – Switch automatic circuit breaker ON again; refer to Chapter "Switch on automatic circuit breaker" (only for types 5650-01, 5660-01 and 5660-51). – Mains switch ON.
TACHO - ERROR	1, 2, 96	Faulty speedometer. Motor, electronics defective.	– Open the cover. – Switch off the mains switch (switch position "0"). – Wait at least 10 seconds. – Turn the rotor vigorously by hand. – Switch on the mains switch again (switch position "I"). The rotor must turn during switch-on.
IMBALANCE	---	The rotor is unevenly loaded.	– Open lid. – Check the loading of the rotor, see chapter "Loading the rotor". – Repeat the centrifugation run.
CONTROL - ERROR	4.1 - 4.5, 6	Error in lid locking or lid closure.	– Perform a MAINS RESET.
N > MAX	5	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10.1, 10.3	Incorrect rotor coding	
MAINS INTERRUPT	---	Power failure during the centrifugation run. (The centrifugation run was not finished.)	– Open lid. – Push START button. – Repeat the centrifugation run if necessary.
VERSION-ERROR	12	Mismatch between electronic components Error / defect electronics	– Perform a MAINS RESET.
CONTROL-ERROR	22, 25.1-25.4	Error / defect electronics	
SER I/O - ERROR	31, 34, 36	Error / defect electronics	
° C * - ERROR	51, 53 - 55 97,98	Error / defect electronics	
° C * - ERROR	52	Temperature is too high in overspeed room. Error / defect electronics	
FU / CCI - ERROR	60, 61.2-61.20, 61.128-61,131, 62	Error / defect electronics / motor	
FU / CCI - ERROR	61.1	Supply voltage is too low Error / defect electronics / motor	– Check supply voltage. – Perform a MAINS RESET.

Message / fault		Cause	Remedy
SENSOR-ERROR	90	Error / defect electronics	– Perform a MAINS RESET.
SENSOR-ERROR	91 - 93	Error / defect unbalance sensor	
NO ROTOR OR ROTORCODE ERROR	---	No rotor installed. Faulty speedometer.	– Open lid. – Install rotor.
N > ROTOR MAX	---	Speed in the selected program greater than the maximum speed of the rotor.	– Check the set speed. Reduce the set speed
N > ROTOR MAX in Prog : e.g. 3		There is a program in the displayed program location whose speed is higher than the maximum speed of the rotor.	
Runtime 00:00 in Prog : e.g. 3		There is a continuous operation program in the displayed program location.	– Replace the continuous operation program in the program linkage with a program with preset time.
Empty Program		There is no program linkage saved on the displayed program location.	– Call up a program linkage.
FC INIT ERROR		Error / defect electronics	– Perform a MAINS RESET.
FC VERSION ERROR		Error / defect electronics	
FATAL EEPROM ERROR	1 - 5	Error / defect electronics	

38 Switch on automatic circuit breaker (only for types 4701-01 and 4706-01)

Switch off the mains switch and separate the centrifuge from the mains!

The automatic circuit breaker is located in the lower right side of the front panel.

- Push the plastic pin (A) of the automatic circuit breaker upwards.
- Reconnect the centrifuge to the mains supply.

39 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

40 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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








- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
- 3) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 4) Load only the inner hole circle. Not closeable with lid 4751
- 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.



4723	4750 + 4751								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ↙ 90°	 mit Bioabdichtung / with bio-containment 10)								
	4763				4764				
									
									
	Kapazität / capacity	ml	7,5 - 8,2	9 - 10	4 – 7	8,5 – 10	15	25	30
Maße / dimensions Ø x L	mm	15 x 92	16 x 92	16 x 75	16 x 100	17 x 100	24 x 100	26 x 95	25 x 90
Anzahl p. Rotor / number p. rotor		72	72	72	72	72	28	28	28
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	³⁾	4030	4030	4030	4030	4030	4030	4030	4030
Radius / radius	mm	178	178	178	178	178	178	178	178
 9 (97%)	sec	40							
 9	sec	40							
Temperatur / temperature	°C ¹⁾	2							
Probenerwärmung/Sample temp. rise	K ²⁾	10							

4723	4750 + 4751								
Ausschwingrotor 4-fach / Swing out rotor 4-times  90°	 mit Bioabdichtung / with bio-containment 10)								
	4764		4765		4766				
									
									
			0521	0548	0523	0547	0538	0549	
Kapazität / capacity	ml	30	30	50	75	100	85	94	85
Maße / dimensions Ø x L	mm	25 x 110	25 x 110	34 x 100	35 x 105	40 x 115	38 x 106	38 x 106	38 x 106
Anzahl p. Rotor / number p. rotor		28	28	16	16	12	12	12	12
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	³⁾	4030	4030	4007	4007	4007	4007	4007	4007
Radius / radius	mm	178	178	177	177	177	177	177	177
 9 (97%)	sec	40							
 9	sec	40							
Temperatur / temperature	°C ¹⁾	2							
Probenerwärmung/Sample temp. rise	K ²⁾	10							

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


- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
- 3) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4723	4750 + 4751								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ↙ 90°	 mit Bioabdichtung / with bio-containment ¹⁰⁾								
	4767	4768	4769	4770	4770 5)	4770 5)	4771	4772	
									
	0526	0530	0509	0513		0546	5127	---	
							 ¹²⁾	 ¹²⁾	
Kapazität / capacity	ml	100	250	15	50	50	50	250	400
Maße / dimensions Ø x L	mm	44 x 100	65 x 115	17 x 120	29 x 115	29 x 115	29 x 107	62 x 122	81 x 136
Anzahl p. Rotor / number p. rotor		8	4	52	20	20	20	4	4
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	³⁾	4007	3985	4166	4143	4143	4143	4166	4166
Radius / radius	mm	177	176	184	183	183	183	184	184
 9 (97%)	sec	40							
 9	sec	40							
Temperatur / temperature	°C ¹⁾	2							
Probenerwärmung/Sample temp. rise	K ²⁾	10							

4723	4750 + 4751												
Ausschwingrotor 4-fach / Swing out rotor 4-times  90°	 mit Bioabdichtung / with bio-containment ¹⁰⁾												
	4772	4773	4774	4775									
													
	---	2078	0536								0553	0501	
	 12)												
Kapazität / capacity	ml	400	1,5	2,0	12	2,7 - 3	2,6 – 2,9	4,9	1,6 - 5	4 - 7	4	5	6
Maße / dimensions Ø x L	mm	84 x 134	11 x 38		17 x 100	11 x 66	13 x 65	13 x 90	13 x 75	13 x 100	12 x 60	12 x 75	12 x 82
Anzahl p. Rotor / number p. rotor		4	96	60		84		84	84		84	84	
Drehzahl / speed	RPM	4500	4500	4500	4500	4500		4500	4500	4500	4500	4500	
RZB / RCF	³⁾	4166	4030	4143	4030	4030		4030	4030	4030	4030	4030	
Radius / radius	mm	184	178	183	178	178		178	178	178	178	178	
 9 (97%)	sec	40											
 9	sec	40											
Temperatur / temperature	°C ¹⁾	2											
Probenerwärmung/Sample temp. rise	K ²⁾	10											

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
- 3) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 5) Take the inserts out of the frame / adapter
- 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.

4723	4750 + 4751							
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°	 mit Bioabdichtung / with bio-containment ¹⁰⁾							
	4775		4776		4777			
								
	0578				Nalgene®		Nunc®	
								
Kapazität / capacity	ml	7	5	225	175	175	200	
Maße / dimensions Ø x L	mm	12 x 100	13 x 75	61 x 137	61 x 118	62 x 144	60 x 130	
Anzahl p. Rotor / number p. rotor		84	84	4	4	4	4	
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	
RZB / RCF	³⁾	4030	4030	4166	4166	4166	4166	
Radius / radius	mm	178	178	184	184	184	184	
 9 (97%)	sec	40						
 9	sec	40						
Temperatur / temperature	°C ¹⁾	2						
Probenerwärmung/Sample temp. rise	K ²⁾	10						

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
 3) Centrifuge containers of glass will not stand RCF values exceeding 4000
 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4723		4752							
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°									
									
		MTP	MS	CP	DWP	MTP	Microtest-platten / plate Terasaki	96-PCR-Platte / plate	PCR-Strips
									
Kapazität / capacity	ml	---	---	---	---	---	---	---	0,2
Maße / dimensions TxBxH / DxWxH	mm	86x128x15	86x128x46	86x128x22	86x128x44,5	86x128x17,5	59x84x11	82x124x20	---
Anzahl p. Rotor / number p. rotor		16	4	12	4	16	8	4	48 x 8
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF		3600	3600	3600	3600	3600	3600	3600	3600
Radius / radius	mm	159	159	159	159	159	159	159	159
 9 (97%)	sec	40							
 9	sec	40							
Temperatur / temperature	°C ¹⁾	2							
Probenerwärmung/Sample temp. Rise	K ²⁾	10							

4753		5051 + 5053							
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°		 							
		5227				5242	5243 ⁴⁾	5243	5247
							 2 x 6316		
			0501			0519		0521	0578
									
Kapazität / capacity	ml	5	6	2,7 – 3	4,5 – 5	25	50	50	7
Maße / dimensions Ø x L	mm	12 x 75	12 x 82	11 x 66	11 x 92	24 x 100	29 x 115	34 x 100	12 x 100
Anzahl p. Rotor / number p. rotor		80	80	80	80	20	8	8	80
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		2773	2773	2773	2773	2755	2755	2755	2755
Radius / radius	mm	155	155	155	155	154	154	154	154
 9 (97%)	sec	18							
 9	sec	16							
Temperatur / temperature	°C ¹⁾	-1							
Probenerwärmung/Sample temp. Rise	K ²⁾	7							

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



CP Kulturplatte /
Culture plate

DWP Deep Well Platte /
Deep well plate


MS Micronic System /
Micronic system


4753	5051 + 5053									
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°	 									
	5247-91	5248				5248-91		5249	5257	
	 6)					 6)				
	0578	0507	----	0518	0507	0518	0523	---	---	
										
Kapazität / capacity	ml	7	15	8,5 - 10	15	15	15	100	1,5	2,0
Maße / dimensions Ø x L	mm	12 x 100	17 x 100	16 x 100	17 x 100	17 x 100	17 x 100	40 x 115	11 x 38	
Anzahl p. Rotor / number p. rotor		80	48	48	48	48	48	4	160	
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF		2755	2755	2755	2755	2755	2755	2755	1950/2826	
Radius / radius	mm	154	154	154	154	154	154	154	109/158	
 9 (97%)	sec	18								
 9	sec	16								
Temperatur / temperature	°C ¹⁾	-1								
Probenerwärmung/Sample temp. rise	K ²⁾	7								

4753	5051 + 5053									
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°	 									
	5281	5258	5258	5259	5262	5264				
		 4)		 4)						
	2078	0536		0513	0526	0500				
										
Kapazität / capacity	ml	1,5	2,0	10	9 - 10	50	100	9	4 - 5,5	7,5 - 8,2
Maße / dimensions Ø x L	mm	11 x 38	15 x 102	16 x 92	29 x 115	44 x 100	14 x 100	15 x 75	15 x 92	
Anzahl p. Rotor / number p. rotor		64	44	44	8	4	48	48	48	
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF		2826	2755	2755	2826	2755	2773	2773	2773	
Radius / radius	mm	158	154	154	158	154	155	155	155	
 9 (97%)	sec	18								
 9	sec	16								
Temperatur / temperature	°C ¹⁾	-1								
Probenerwärmung/Sample temp. rise	K ²⁾	7								

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
 4) not closeable with lid 5053
 6) with decanting aid

4753	5051 + 5053								
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$									
	5264	5266		5267		5268			
									
									
Kapazität / capacity	ml	4 - 7	30	30	3	1,1 - 1,4	2,6 - 2,9	4,9	1,6 - 5
Maße / dimensions $\varnothing \times L$	mm	16 x 75	25 x 110	25 x 110	10 x 60	8 x 66	13 x 65	13 x 90	13 x 75
Anzahl p. Rotor / number p. rotor		48	20	20	80	80	48	48	48
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		2773	2755	2755	2737	2737	2808	2808	2808
Radius / radius	mm	155	154	154	153	153	157	157	157
 9 (97%)	sec	18							
 9	sec	16							
Temperatur / temperature	°C ¹⁾	-1							
Probenerwärmung/Sample temp. rise	K ²⁾	7							

4753	5051 + 5053							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$	 							
	5268	6306	6306	SK 113/88				
		 4)						
								
Kapazität / capacity	ml	4 - 7	15	12	25			
Maße / dimensions $\varnothing \times L$	mm	16 x 75	17 x 120	17 x 100	25 x 90			
Anzahl p. Rotor / number p. rotor		48	28	28	20			
Drehzahl / speed	RPM	4000	4000	4000	4000			
RZB / RCF		2808	2898	2898	2826			
Radius / radius	mm	157	162	162	158			
 9 (97%)	sec	18						
 9	sec	16						
Temperatur / temperature	$^\circ\text{C}^{1)}$	-1						
Probenerwärmung/Sample temp. rise	$\text{K}^{2)}$	7						

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








- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
- 3) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.


4753	5092 + 5093								
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°									
	mit Bioabdichtung / with bio-containment ¹⁰⁾								
	5123	5124	5125	5126	5128		5129	5134	
									
	0513	0521	0526	0523	0501		0509		
									
Kapazität / capacity	ml	50	50	100	100	6	5	15	25
Maße / dimensions Ø x L	mm	29 x 115	34 x 100	44 x 100	40 x 115	12 x 82	13 x 75	17 x 120	25 x 90
Anzahl p. Rotor / number p. rotor		8	4	4	4	48	48	28	12
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		3095	2952	2952	2952	3005	3005	3095	2826
Radius / radius	mm	173	165	165	165	168	168	173	158
 9 (97%)	sec	18							
 9	sec	16							
Temperatur / temperature	°C ¹⁾	-1							
Probenerwärmung/Sample temp. rise	K ²⁾	7							

4753	5092 + 5093								
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°									
	mit Bioabdichtung / with bio-containment 10)								
	5135	5136							
									
	---	2079	0507						
									
Kapazität / capacity	ml	50	10	15	10	4 – 4,5	7,5 – 8,2	9 - 10	10
Maße / dimensions Ø x L	mm	29 x 115	17 x 70	17 x 100	16 x 80	15 x 75	15 x 92	16 x 92	15 x 102
Anzahl p. Rotor / number p. rotor		8	32	32	32	32	32	32	32
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		3023	2952	2952	2952	2952	2952	2952	2952
Radius / radius	mm	169	165	165	165	165	165	165	165
 9 (97%)	sec	18							
 9	sec	16							
Temperatur / temperature	°C ¹⁾	-1							
Probenerwärmung/Sample temp. rise	K ²⁾	7							

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

- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C room temperature (only with cooling centrifuges)
 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".


4753	5092 + 5093								
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$	 								
	mit Bioabdichtung / with bio-containment ¹⁰⁾								
	5136				5137				
									
									
Kapazität / capacity	ml	8,5 – 10	4 - 7	15	5	6	1,1 – 1,4	2,7 - 3	2,6 – 2,9
Maße / dimensions $\varnothing \times L$	mm	16 x 100	16 x 75	17 x 100	12 x 75	12 x 82	8 x 66	11 x 66	13 x 65
Anzahl p. Rotor / number p. rotor		32	32	32	32	32	32	32	32
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		2952	2952	2952	2952	2952	2952	2952	2952
Radius / radius	mm	165	165	165	165	165	165	165	165
 9 (97%)	sec	18							
 9	sec	16							
Temperatur / temperature	°C ¹⁾	-1							
Probenerwärmung/Sample temp. rise	K ²⁾	7							

4753	5092 + 5093								
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$	 								
	mit Bioabdichtung / with bio-containment ¹⁰⁾								
	5137						SK 21.00		
									
									
Kapazität / capacity	ml	4,9	4,5 - 5	1,6 - 5	4 - 7	5	3	4	1,5 2,0
Maße / dimensions $\varnothing \times L$	mm	13 x 90	11 x 92	13 x 75	13 x 100	13 x 75	10 x 60	10 x 88	11 x 38
Anzahl p. Rotor / number p. rotor		32	32	32	32	32	40	40	40
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF		2952	2952	2952	2952	2952	2916	2916	2916
Radius / radius	mm	165	165	165	165	165	163	163	163
 9 (97%)	sec	18							
 9	sec	16							
Temperatur / temperature	°C ¹⁾	-1							
Probenerwärmung/Sample temp. rise	K ²⁾	7							

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4753	5092 + 5093									
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°										
	mit Bioabdichtung / with bio-containment 10)									
	SK 21.00		SK 94.89 4)	SK 94.89	5138 					
										
	Kapazität / capacity	ml	1,1 – 1,4	2,7 - 3	225	175	1,1 – 1,4	2,7-3	2,6 - 2,9	1,6 - 5
	Maße / dimensions $\varnothing \times L$	mm	8 x 66	11 x 66	61 x 137	61 x 118	8 x 66	11 x 66	13 x 65	13 x 75
Anzahl p. Rotor / number p. rotor		40	40	4	4	48	48	48	48	
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF		2916	2916	3095	3095	2540	2540	2540	2540	
Radius / radius	mm	163	163	173	173	142	142	142	142	
 9 (97%)	sec	18								
 9	sec	16								
Temperatur / temperature	°C ¹⁾	-1								
Probenerwärmung/Sample temp. rise	K ²⁾	7								

4754	1425								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	1431				1432		1433	1434	
									
	0500	0507		0518	Rhesus	---		0501	
									
Kapazität / capacity	ml	9	15	7,5 - 8,2	15	1	0,4	3	6
Maße / dimensions $\varnothing \times L$	mm	14 x 100	17 x 100	15 x 92	17 x 100	6 x 45	6 x 45	10 x 60	12 x 82
Anzahl p. Rotor / number p. rotor		28	28	28	28	144	144	56	48
Drehzahl / speed	RPM	4600	4600	4600	4600	4600	4600	4600	4600
RZB / RCF		3596	3596	3596	3596	3643	3643	3596	3596
Radius / radius	mm	152	152	152	152	154	154	152	152
 9 (97%)	sec	14							
 9	sec	14							
Temperatur / temperature	°C ¹⁾	0							
Probenerwärmung/Sample temp. rise	K ²⁾	9							

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









- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)


4754		1425						
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
		1441	1443	1444	1458	1736	1737	
								
								
Kapazität / capacity	ml	8,5 - 10	50	1,5	2,0	1,1 - 1,4	15	50
Maße / dimensions $\varnothing \times L$	mm	16 x 100	29 x 115	11 x 38	11 x 38	8 x 66	17 x 120	29 x 115
Anzahl p. Rotor / number p. rotor		28	4	36	36	36	16	4
Drehzahl / speed	RPM	4600	4600	4600	4600	4600	4600	4600
RZB / RCF		3596	3738	3572	3572	3690	3738	3738
Radius / radius	mm	152	158	151	151	156	158	158
 9 (97%)	sec	14						
 9	sec	14						
Temperatur / temperature	$^{\circ}\text{C}^{1)}$	0						
Probenerwärmung/Sample temp. rise	$\text{K}^{2)}$	9						

4754		1427 + 1421							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$		 							
		1326	1327	1357	1732				
									
				Rhesus					
									
Kapazität / capacity	ml	4	3	1	0,4	5	2,6 - 2,9	1,6 - 5	5
Maße / dimensions $\varnothing \times L$	mm	12 x 60	10 x 60	6 x 45	6 x 45	12 x 75	13 x 65	13 x 75	13 x 75
Anzahl p. Rotor / number p. rotor		48	48	120	120	32	32	32	32
Drehzahl / speed	RPM	4600	4600	4600	4600	4600	4600	4600	4600
RZB / RCF		3714	3714	3738	3738	3690	3690	3690	3690
Radius / radius	mm	157	157	158	158	156	156	156	156
 9 (97%)	sec	14							
 9	sec	14							
Temperatur / temperature	$^{\circ}\text{C}^{1)}$	0							
Probenerwärmung/Sample temp. rise	$\text{K}^{2)}$	9							

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




- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
5) Take the inserts out of the frame / adapter

4754	1427 + 1421						
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°	 + 						
	5278	5279					
							
							
Kapazität / capacity ml	1,1 – 1,4	4,5 - 5	7,5 – 8,2				
Maße / dimensions Ø x L mm	8 x 66	15 x 75	15 x 92				
Anzahl p. Rotor / number p. rotor	48	20	20				
Drehzahl / speed RPM	4600	4600	4600				
RZB / RCF	3761	3761	3761				
Radius / radius mm	159	159	159				
 9 (97%) sec	14						
 9 sec	14						
Temperatur / temperature °C ¹⁾	0						
Probenerwärmung/Sample temp. rise K ²⁾	9						

4754	1427						
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°							
	1732	1733	5230	5231	5233 ⁵⁾	5276	
					 2 x 6316		
							
Kapazität / capacity ml	4,9	4 - 7	15	4,5 - 5	8,5 - 10	50	50
Maße / dimensions Ø x L mm	13 x 90	13 x 100	17 x 120	11 x 92	16 x 100	29 x 115	29 x 115
Anzahl p. Rotor / number p. rotor	32	32	12	48	24	4	4
Drehzahl / speed RPM	4600	4600	4600	4600	4600	4600	4600
RZB / RCF	3690	3690	3809	3620	3620	3714	3714
Radius / radius mm	156	156	161	153	153	157	157
 9 (97%) sec	14						
 9 sec	14						
Temperatur / temperature °C ¹⁾	0						
Probenerwärmung/Sample temp. rise K ²⁾	9						

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








1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)


2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)

5) Take the inserts out of the frame / adapter

7) can not be centrifugated when plug is attached, Scal. 10µl-300µl, 15ml, 30ml



10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".


4754	1430 + 1382								
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°									
	mit Bioabdichtung / with bio-containment 10)								
	1348				1351		1363	1365	
									
									
Kapazität / capacity	ml	10	8	4 – 7	8,5 - 10	1,5	2,0	25	30
Maße / dimensions Ø x L	mm	16 x 80	16 x 81	16 x 75	16 x 100	11 x 38	11 x 38	25 x 90	25 x 110
Anzahl p. Rotor / number p. rotor		16	16	16	16	20	20	4	4
Drehzahl / speed	RPM	4600	4600	4600	4600	4600	4600	4600	4600
RZB / RCF		3690	3690	3690	3690	3549	3549	3407	3832
Radius / radius	mm	156	156	156	156	150	150	144	162
 9 (97%)	sec	14							
 9	sec	14							
Temperatur / temperature	°C ¹⁾	0							
Probenerwärmung/Sample temp. rise	K ²⁾	9							

4754	1430 + 1382								
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°									
	mit Bioabdichtung / with bio-containment 10)								
	1383								
									
									
Kapazität / capacity	ml	5	6	7	2,7 - 3	2,6 – 2,9	4,9	4,5 - 5	1,6 – 5
Maße / dimensions Ø x L	mm	12 x 75	12 x 82	12 x 100	11 x 66	13 x 65	13 x 90	11 x 92	13 x 75
Anzahl p. Rotor / number p. rotor		20	20	20	20	20	20	20	20
Drehzahl / speed	RPM	4600	4600	4600	4600	4600	4600	4600	4600
RZB / RCF		3690	3690	3690	3690	3690	3690	3690	3690
Radius / radius	mm	156	156	156	156	156	156	156	156
 9 (97%)	sec	14							
 9	sec	14							
Temperatur / temperature	°C ¹⁾	0							
Probenerwärmung/Sample temp. rise	K ²⁾	9							

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
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
 10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".


4754	1430 + 1382								
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$									
	mit Bioabdichtung / with bio-containment ¹⁰⁾								
	1383		1384	1396		1457	1459		
									
									
Kapazität / capacity	ml	4 - 7	5	50	85	85	1,1 – 1,4	4 – 5,5	7,5 – 8,2
Maße / dimensions $\varnothing \times L$	mm	16 x 75	13 x 75	29 x 115	38 x 106	38 x 106	8 x 66	15 x 75	15 x 92
Anzahl p. Rotor / number p. rotor		20	20	4	4	4	28	16	16
Drehzahl / speed	RPM	4600	4600	4600	4600	4600	4600	4600	4600
RZB / RCF		3690	3690	3832	3761	3761	3714	3714	3714
Radius / radius	mm	156	156	162	159	159	157	157	157
 9 (97%)	sec	14							
 9	sec	14							
Temperatur / temperature	°C ¹⁾	0							
Probenerwärmung/Sample temp. rise	K ²⁾	9							

4754	1430 + 1382							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
	mit Bioabdichtung / with bio-containment ¹⁰⁾							
	4416	4417	6311	6318				
								
	0546	0545						
								
Kapazität / capacity	ml	50	30	12	50			
Maße / dimensions $\varnothing \times L$	mm	29 x 107	26 x 95	17 x 100	29 x 115			
Anzahl p. Rotor / number p. rotor		4	4	4	4			
Drehzahl / speed	RPM	4600	4600	4600	4600			
RZB / RCF		3785	3549	3832	3832			
Radius / radius	mm	160	150	162	162			
 9 (97%)	sec	14						
 9	sec	14						
Temperatur / temperature	°C ¹⁾	0						
Probenerwärmung/Sample temp. rise	K ²⁾	9						

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






1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)

2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)

6) with decanting aid

11) Plastic suspension can only be used in temperatures up to a maximum of 40°C / 104°F.

4758	1741				1742							
Ausschwingrotor 8-fach / Swing out rotor 8-times   90°	 11)				 11)							
	0701							0716				
												
		0500			0507	 *)	0509		0518			
Kapazität / capacity	ml	4,9	4,5 - 5	9	1,1 – 1,4	15	15	15	2,6 - 2,9	4 - 4,5	1,6 - 5	4-7
Maße / dimensions Ø x L	mm	13 x 90	11 x 92	14 x 100	8 x 66	17 x 100	17 x 120	17 x 100	13 x 65	15 x 75	13 x 75	16x 75
Anzahl p. Rotor / number p. rotor		80		80	80	56	24	56	56		56	
Drehzahl / speed	RPM	3800		3800	3800	3800	3800	3800	3800		3800	
RZB / RCF		2906		2874	2874	2906	2906	2906	2793		2793	
Radius / radius	mm	180		178	178	180	180	180	173		173	
 9 (97%)	sec	19										
 9	sec	16										
Temperatur / temperature	°C ¹⁾	0										
Probenerwärmung/Sample temp. rise	K ²⁾	11										

4758	1742				1745	1746			SK 13.06 11)	
Ausschwingrotor 8-fach / Swing out rotor 8-times   90°	 11)				 11)	 11)				
							1641	1641		
								1462		
	   *)     				0545	0521	0513	0509	0519	
Kapazität / capacity	ml	7,5-8,2	9-10	10	8,5 - 10	30	50	50	15	25
Maße / dimensions Ø x L	mm	15 x 92	16 x 92	15 x 102	16 x 100	26 x 95	34 x 100	29 x 115	17 x 120	24 x 100
Anzahl p. Rotor / number p. rotor		56		24	56	16	8	8	8	8
Drehzahl / speed	RPM	3800		3800	3800	3800	3800	3800	3800	3800
RZB / RCF		2906		2906	2906	2906	2906	2906	2906	2793
Radius / radius	mm	180		180	180	180	180	180	180	173
 9 (97%)	sec	19								
 9	sec	16								
Temperatur / temperature	°C ¹⁾	0								
Proben erwärmung / Sample temp. rise	K ²⁾	11								

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













- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
 2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
 3) Centrifuge containers of glass will not stand RCF values exceeding 4000
 12) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.




4760 + 1162								
Winkelrotor 6-fach / Angle rotor 6-times  $\angle 45^\circ$								
		PCR-Strips	---					
								
Kapazität / capacity	ml	0,2	0,2					
Maße / dimensions $\varnothing \times L$	mm	----	6 x 18					
Anzahl p. Rotor / number p. rotor		6	48					
Drehzahl / speed	RPM	14000	14000					
RZB / RCF	³⁾	18845	18845					
Radius / radius	mm	86	86					
 9 (97%)	sec	15						
 9	sec	15						
 1	sec	285						
Temperatur / temperature	°C ¹⁾	1						
Probenerwärmung/sample temp. rise	K ²⁾	19						

4790-A		2031 ⁹⁾	2023	2024				
Winkelrotor 30-fach / Angle rotor 30-times  $\angle 45^\circ$ mit Bioabdichtung / with bio-containment 10)								
		0536	2078	---	---	---	---	
								
Kapazität / capacity	ml	2,0	1,5	0,8	0,5	0,4	0,2	
Maße / dimensions $\varnothing \times L$	mm	11 x 38	11 x 38	8 x 45	8 x 30	6 x 45	6 x 18	
Anzahl p. Rotor / number p. rotor		30	30	30	30	30	30	
Drehzahl / speed	RPM	15000	15000	15000	15000	15000	15000	
RZB / RCF	³⁾	24400	24400	24400	24400	24400	24400	
Radius / radius	mm	97	97	97	97	97	97	
 9 (97%)	sec	19						
 9	sec	24						
Temperatur / temperature	°C ¹⁾	4						
Probenerwärmung/sample temp. rise	K ²⁾	19						

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
2) Sample temp. rise during maximum speed and 1 h running time (only with centrifuges without cooling)
3) Centrifuge containers of glass will not stand RCF values exceeding 4000
9) recommended for high-speed centrifugation
10) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4794				1446	1447	1448	1449	
mit Bioabdichtung / with bio-containment 10)								
	0547	0538	0549	0546	0519	0545	---	---
								
Kapazität / capacity	85	94	85	50	25	30	10	3
Maße / dimensions Ø x L	38 x 106	38 x 106	38 x 106	29 x 107	24 x 100	26 x 95	16 x 80	10 x 60
Anzahl p. Rotor / number p. rotor	6	6	6	6	6	6	12	24
Drehzahl / speed	11000	11000	11000	11000	11000	11000	11000	11000
RZB / RCF	16504	16504	16504	15828	15151	15151	15557	15828
Radius / radius	122	122	122	117	112	112	115	117
 9 (97%)	40							
 9	63							
Temperatur / temperature	4							
Probenerwärmung/Sample temp. rise	19							

4794	1449		1451					1454		
 45° mit Bioabdichtung / with bio-containment 10)										
	2078	0536	0507					0518	0513	
										
Kapazität / capacity	ml	1,5	2	6	7,5 – 8,2	9 - 10	10	8,5 - 10	15	20
Maße / dimensions Ø x L	mm	11 x 38	11 x 38	17 x 100	15 x 92	16 x 92	15 x 102	16 x 100	17 x 100	29 x 115
Anzahl p. Rotor / number p. rotor		24	24	6	6	6		6	6	6
Drehzahl / speed	RPM	11000	11000	11000	11000	11000		11000	11000	11000
RZB / RCF	³⁾	15828	15828	15557	15557	15557		15828	15557	16098
Radius / radius	mm	117	117	115	115	115		117	115	119
 9 (97%)	sec	40								
 9	sec	63								
Temperatur / temperature	°C ¹⁾	4								
Probenerwärmung/Sample temp. rise	K ²⁾	19								

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All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.