

Glasses Reader Program Documentation and Help Summary page.

Click on one of the following links to go directly to the area of documentation which interests you.
Click on the title page of page 1 of each section to return to this menu.

[Build Inventory](#)

[Manual Entry](#)

[Barcode Entry](#)

[Lensmeter Entry](#)

[File Options](#)

[Reports Options](#)

[Options Window](#)

[Printer](#)

Building an inventory can be an extensive task. Typical inventories are from 2000-8000 pairs of glasses. The Glasses Reader program makes this a considerably less complex task.

The Glasses Reader program allows input from three different sources:

1. Lensmeters (currently there are 8 different models),
2. Glasses which are labeled with Barcodes (produced by this program), and
3. Manual data entry.

In order to build such a large inventory, it is necessary for several people to be working on it simultaneously. To understand how this might work, here is an example.

Suppose that 4 people are building an inventory of 4000 pairs of glasses. Prior to their starting their work, a team of people at their church have gotten together on a Saturday and placed labels on the temples and the bags of all the glasses. They printed these labels using the Excel spreadsheet named **Boxes and Glasses Labels.xls** on the Glasses Reader CD. The 4 people now have glasses which have already been assigned inventory numbers.

Here is the list:

Name of person	Range of numbers assigned	How he/she is building the inventory.
James	1-1000	Using a Reichert LensChek lensmeter.
Mary	1001-2000	All her glasses are barcoded.
George	2001-3000	George is entering the prescriptions manually.
Fred	3001-4000	Fred has a TopCon CL-100 Lensmeter and some of his are entered manually.

James has taken the responsibility of maintaining the master copy of the inventory. Mary, George and Fred will E-Mail their work to James to enter into his system. **Browse and load / merge other Inventory Files** is the Glasses Reader options that he uses:

James and the other people have set up their options Window as shown below.

Path to collection files.

C:\My Documents\Glasses_Data\ Total_Inventory.txt

OR→ No File Output or input. Check to read file into memory at startup.

File

- Load Inventory File from Default Location.
- Browse and load / merge other Inventory Files. ←
- Load Deleted & Rejected inventory from Glasses Inventory program.
- Load Update inventory into Vacant Locations.
- Transfer inventory to Glasses Inventory Program

Their inventory will be placed into the:

C:\My Documents\Glasses_Data
directory and is called: **Total_Inventory.txt**

[Go to Next Page →](#)

James receives a message from Mary. She has completed the entry of all her glasses. Barcoded glasses can be entered in 1/10 the time of any other form of data entry. James saves the attachment to C:\TEMP on his PC. Here are the steps that he follows to merge her inventory with his inventory:

1. Click: **File | Browse and load / merge other inventory files.**

Enter the Inventory Number Load Range:


Lowest Allowable Number. 1 Enter the lowest inventory number that you expect will be within this file that you are loading. Lower numbers will be discarded.	Highest Allowable Number. 32000 Enter the highest inventory number that you expect will be within this file that you are loading. Higher numbers will be discarded.
Exit and do nothing.	Continue & Select File.

2. Change the range of acceptable number to match what Mary is providing.


Enter the Inventory Number Load Range:

Lowest Allowable Number. 1001 Enter the lowest inventory number that you expect will be within this file that you are loading. Lower numbers will be discarded.	Highest Allowable Number. 2000 Enter the highest inventory number that you expect will be within this file that you are loading. Higher numbers will be discarded.
Exit and do nothing.	Continue & Select File.

3. James browses to where he saved Mary's file and clicks Open.



4. After the load is completed, the status window is below.




and the highest inventory number displays below.





5. James can then click the following button, both options are already selected as shown to the right and then the new combined inventory will be written and error checked 100 records at a time.

Write Sorted Inventory and Check for Errors.




[Go to Next Page](#) →

When the new inventory has been written, there is an indication as to how many errors were detected as shown below.



The report will list the numbers of the locations that have an error. Enter the location number of each location in error and the system will describe the error. You can then correct this error and save the record again. After all locations have been corrected, you should write the sorted inventory file again.

Should James try to exit without saving, he will see.



Click YES and the save will take place. Otherwise, the change is not saved.

As different team members complete their inventory, James using the same procedure to merge all their files together.

This help covers the part of the screen which is shown below. Click on blue underlined areas for the part of the screen for which you need help or information.

USING: MANUAL ENTRY PROCEDURE

Status: Complete.

Inventory Number: 15

PRINT Label to: DYMO LabelWriter 320

of Rejections: []

Manual Entry Help (click here)

Left (OS) Clear Left

Sphere: [] Both

Cylinder: [] Both

Axis: [] Both

Add: [] Both

Transpose Reading

Both the same (OU)

C [] Undo

1	2	3	.00
4	5	6	.25
7	8	9	.50
+	0	-	.75

Right (OD) Clear Right

Sphere: [] Both

Cylinder: [] Both

Axis: [] Both

Add: [] Both

Transpose Reading

Clear Record

FLAGS?

Sunglasses? Progressive? Scratched? For Male? Either Male or Female For Female? For Child?

Reading Glasses? Bi-Trifocal?

Frame Size? Small? Medium? Large? PD [] mm

Save & Clear

Exit

Comment on Glasses (Use NO commas)

Enter your Initials: []

Save your entries after completing manual entry of data.

PAGE 1

PAGE 2

PAGE 3

PAGE 4

Be sure you select your inventory number.

Inventory Number

15

Left (OS)	Both the same (OU)	Right (OD)
<input type="button" value="Clear Left"/>	<input type="button" value="C"/> <input type="text" value=""/> <input type="button" value="Undo"/>	<input type="button" value="Clear Right"/>
Sphere: <input type="text" value=""/> <input type="button" value="Both"/>	1 2 3 .00	Sphere: <input type="text" value=""/> <input type="button" value="Both"/>
Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>	4 5 6 .25	Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>
Axis: <input type="text" value=""/> <input type="button" value="Both"/>	7 8 9 .50	Axis: <input type="text" value=""/> <input type="button" value="Both"/>
Add: <input type="text" value=""/> <input type="button" value="Both"/>	+ 0 - .75	Add: <input type="text" value=""/> <input type="button" value="Both"/>
<input type="button" value="Transpose Reading"/>		<input type="button" value="Transpose Reading"/>

The above keypad is used for manual data entry. Below are the steps for entering the following prescription

<-----Left Eye ----->

Sphere Cylinder Axis Add

2.50 -1.25 87 2.25

<-----Right Eye ----->

Sphere Cylinder Axis Add

2.50 -.75 97 2.25

Click 2 .50

Left (OS)	Both the same (OU)	Right (OD)
<input type="button" value="Clear Left"/>	<input type="button" value="C"/> 2.50 <input type="button" value="Undo"/>	<input type="button" value="Clear Right"/>
Sphere: <input type="text" value=""/> <input type="button" value="Both"/>	1 2 3 .00	Sphere: <input type="text" value=""/> <input type="button" value="Both"/>
Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>	4 5 6 .25	Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>
Axis: <input type="text" value=""/> <input type="button" value="Both"/>	7 8 9 .50	Axis: <input type="text" value=""/> <input type="button" value="Both"/>
Add: <input type="text" value=""/> <input type="button" value="Both"/>	+ 0 - .75	Add: <input type="text" value=""/> <input type="button" value="Both"/>
<input type="button" value="Transpose Reading"/>		<input type="button" value="Transpose Reading"/>

Then click Either BOTH to copy it to both sides.

Left (OS)	Both the same (OU)	Right (OD)
<input type="button" value="Clear Left"/>	<input type="button" value="C"/> 2.50 <input type="button" value="Undo"/>	<input type="button" value="Clear Right"/>
Sphere: 2.50 <input type="button" value="Both"/>	1 2 3 .00	Sphere: 2.50 <input type="button" value="Both"/>
Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>	4 5 6 .25	Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>
Axis: <input type="text" value=""/> <input type="button" value="Both"/>	7 8 9 .50	Axis: <input type="text" value=""/> <input type="button" value="Both"/>
Add: <input type="text" value=""/> <input type="button" value="Both"/>	+ 0 - .75	Add: <input type="text" value=""/> <input type="button" value="Both"/>
<input type="button" value="Transpose Reading"/>		<input type="button" value="Transpose Reading"/>

Click 1 .25

Left (OS)	Both the same (OU)	Right (OD)
<input type="button" value="Clear Left"/>	<input type="button" value="C"/> 1.25 <input type="button" value="Undo"/>	<input type="button" value="Clear Right"/>
Sphere: 2.50 <input type="button" value="Both"/>	1 2 3 .00	Sphere: 2.50 <input type="button" value="Both"/>
Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>	4 5 6 .25	Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>
Axis: <input type="text" value=""/> <input type="button" value="Both"/>	7 8 9 .50	Axis: <input type="text" value=""/> <input type="button" value="Both"/>
Add: <input type="text" value=""/> <input type="button" value="Both"/>	+ 0 - .75	Add: <input type="text" value=""/> <input type="button" value="Both"/>
<input type="button" value="Transpose Reading"/>		<input type="button" value="Transpose Reading"/>

Then click into the Left Cylinder Field.

Left (OS)	Both the same (OU)	Right (OD)
<input type="button" value="Clear Left"/>	<input type="button" value="C"/> 1.25 <input type="button" value="Undo"/>	<input type="button" value="Clear Right"/>
Sphere: 2.50 <input type="button" value="Both"/>	1 2 3 .00	Sphere: 2.50 <input type="button" value="Both"/>
Cylinder: -1.25 <input type="button" value="Both"/>	4 5 6 .25	Cylinder: <input type="text" value=""/> <input type="button" value="Both"/>
Axis: <input type="text" value=""/> <input type="button" value="Both"/>	7 8 9 .50	Axis: <input type="text" value=""/> <input type="button" value="Both"/>
Add: <input type="text" value=""/> <input type="button" value="Both"/>	+ 0 - .75	Add: <input type="text" value=""/> <input type="button" value="Both"/>
<input type="button" value="Transpose Reading"/>		<input type="button" value="Transpose Reading"/>

Notice that minus numbers always show up as **yellow**.

Notice that when you enter a positive cylinder the system automatically makes it minus. (Option setting)

Click **8 7** and then click the **Left Axis Field**.

Left (OS)		Both the same (OU)		Right (OD)	
Clear Left		C 87 Undo		Clear Right	
Sphere:	2.50 Both	1	2	3	.00
Cylinder:	-1.25 Both	4	5	6	.25
Axis:	87 Both	7	8	9	.50
Add:	Both	+ 0 -	.75		
Transpose Reading		Transpose Reading			

Click **2 .25** and then click the **Left Add Field**.

Left (OS)		Both the same (OU)		Right (OD)	
Clear Left		C 2.25 Undo		Clear Right	
Sphere:	2.50 Both	1	2	3	.00
Cylinder:	-1.25 Both	4	5	6	.25
Axis:	87 Both	7	8	9	.50
Add:	2.25 Both	+ 0 -	.75		
Transpose Reading		Transpose Reading			

Notice how the **Left Add** is duplicated on both the left and right. If that is wrong, just overwrite the other **Add**.

Now click **.75** and then click the **right Cylinder field**. Click **97** then click the **right axis field**.

Left (OS)		Both the same (OU)		Right (OD)	
Clear Left		C .75 Undo		Clear Right	
Sphere:	2.50 Both	1	2	3	.00
Cylinder:	-1.25 Both	4	5	6	.25
Axis:	87 Both	7	8	9	.50
Add:	2.25 Both	+ 0 -	.75		
Transpose Reading		Transpose Reading			

Left (OS)		Both the same (OU)		Right (OD)	
Clear Left		C 97 Undo		Clear Right	
Sphere:	2.50 Both	1	2	3	.00
Cylinder:	-1.25 Both	4	5	6	.25
Axis:	87 Both	7	8	9	.50
Add:	2.25 Both	+ 0 -	.75		
Transpose Reading		Transpose Reading			

Now select the flags which apply for this pair of glasses (Prescription sunglasses for a Male, Medium size)

<input checked="" type="checkbox"/> Sunglasses? <input type="checkbox"/> Reading Glasses?		<input type="checkbox"/> Progressive? <input checked="" type="checkbox"/> Bi-Trifocal?		<input checked="" type="checkbox"/> For Male? <input type="checkbox"/> For Female?		<input type="checkbox"/> Either Male or Female <input type="checkbox"/> For Child?		Frame Size? <input type="checkbox"/> Small? <input checked="" type="checkbox"/> Medium? <input type="checkbox"/> Large?		PD <input type="text" value=""/> mm
---	--	---	--	---	--	---	--	--	--	-------------------------------------

Enter your comment (if any) on this pair of glasses and then click **Save & Clear**.

You can enter the mm size of the PD here.

Save & Clear	Comment on pair of glasses
	Comment on Glasses (Use NO commas)

You will see:

Save & Clear	15 saved successfully...
--------------	--------------------------

and you are ready to enter the next pair of glasses.

Notice that you can double click any field to clear it. Double click an **Add** field and both **Adds** are cleared.

Click the appropriate UP or DOWN arrows to change the inventory number up or down the indicated amount.

Inventory Number
15

10000 1000 100 10 1

<-Last Empty Next Empty ->

Click Last Empty to go backwards to the last empty location.

Click Next Empty to go forwards to the next empty location.

This displays the number of times this pair of glasses has been rejected by the Eyeglasses Inventory Program.

of Rejections
0

Click this button to empty a location.

Clear Record

Enter the initials of the operator and they will appear on the barcoded label in the comments field.


Enter your Initials:

PRINT Label to:
DYMO LabelWriter 320


Click to print the label (as selected in Options) to the indicated label printer.

This help covers the part of the screen which is shown below. Click on blue underlined areas for the part of the screen for which you need help or information.

Located at the bottom / left portion of the screen.




Located at the bottom / right portion of the screen



1. When you click the above checkmark,
2. The yellow button on the right bottom appears.
3. When you click the above yellow button,
4. The below window appears.



5. When you click the above Expand Window button.
6. The below window appears.



Below is a repeat of the larger barcode window and the very important inventory number on the main window.

Barcode Label

OD: -1.75 -0.75 89 1.25 [10453]
OS: -1.25 -1.25 103 1.25 PD: MD
[BiFocal, Either M/F (8/8/2004)]



1. The bottom barcode is scanned into the window.
(You may scan either barcode first)

Inventory Number
3976

5. And the Inventory number increments automatically to the next number.

Bar Code Reader Input.

BarCode Data: S1./-2./123/2.25/12

Shrink Window Save HIDE Clear Record To Keyboard

FLAGS? (Change to modify barcode data)

Sunglasses? Progressive? For Male? Either Male or Female

Reading Glasses? Bi-Trifocal? For Female? For Child?

Frame Size? (Change to modify barcode data)

Small? Medium? Large? [] mm

Inventory Number

10000 1000 100 10 1

<--Last Empty Next Empty-->

2. The data is extracted and placed in the left side of the main window.

4. When both barcodes are scanned, the SAVE barcode is scanned also.

Left (OS)		Right (OD)
Sphere: 1.00 Both		Sphere: [] Both
Cylinder: -2.00 Both		Cylinder: [] Both
Axis: 123 Both		Axis: [] Both
Add: 2.25 Both		Add: 2.25 Both

3. The Add is copied from left to right.



Under Options you will see the following section at the top of this window. This is where you choose the lensmeter you are using to measure your glasses.

You can also switch to lensmeter data entry at any time by clicking the checkmark as shown below on the bottom / left of the main window at an time.

Identify your Lensmeter

<input type="radio"/> Reichert LensChek	<input type="radio"/> Humphrey LA360	<input type="radio"/> Tomey TL-3000
<input checked="" type="radio"/> Reichert AL200	<input type="radio"/> Topcon CL-100	<input type="radio"/> Nikon NL-2
<input type="radio"/> Humphrey LA340	<input type="radio"/> Tomey TL-900	
<input type="radio"/> Barcode Reader	<input type="radio"/> Keyboard input	

Data Input Selection:

<input type="checkbox"/> Manual Data Entry?
<input type="checkbox"/> Barcode Data Entry?
<input checked="" type="checkbox"/> AL200 Lensmeter

The below shows where you select the Comm Port number on your computer to which your lensmeter is attached.

Select COMM Port number

None 1 2 3 4 5 6 7 8 9 10 11 12

Click on the picture of your lensmeter below for help on this unit.

- 1) Reichert LensCheck, 2) Reichert AL-200, 3) Humphrey LA360, 4) Tomey TL-900, 5) Tomey TL-3000, 6) Humphrey LA340, 7) Nikon NL2, and 8) Topcon CL-100.

-----1-----



[Click to Go to Page 2](#)

-----2-----



-----3-----



-----4-----



-----5-----



-----6-----



-----7-----



-----8-----



Reichert LensChek Help

Connect Kendall Optometry Ministries, Inc part number 30005 cable between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.


To Measure Single Vision and BiFocal Glasses:

1. Use top blue button and select the setting shown to the right.
2. Be sure that the correct inventory number is showing in the Glasses Reader program. Verify that this number matches the number you have on the glasses bag and frame.
3. Place Glasses on platform facing you with temples away from you. Bottom of frame is down. Right lens will be first.
4. Press the third blue button on the right until you see the word "Right" showing under the starburst pattern.
5. Make sure you have pressed the top right blue button to point to the center symbol (Single / BiFocal).
6. Adjust platform up and down and glasses left and right until you see the starburst pattern (to the right).
7. If these are BiFocals press the 2nd blue button until the pointer points to "Add". If these are single vision, skip to step 8.
8. For BiFocals move the pair of glasses until the Add segment is position over the hold in the rubber boot. The Add will be automatically read.
9. Press the Hold Button.
10. Press the third blue button on the right until you see the word "Left" showing under the starburst pattern.
11. Adjust platform up and down and glasses left and right until you see the starburst pattern (to the right).
12. If these are BiFocals press the 2nd blue button until the pointer points to "Add". If these are single vision, skip to step 8.
13. For BiFocals move the pair of glasses until the Add segment is position over the hold in the rubber boot. The Add will be automatically read.
14. Press the Hold Button.
15. Press the PRINT blue button (above the green button).
16. Measurement will be transferred to the glasses reader program.
17. Click SAVE on the Glasses Reader program.
18. Start at step 1 for the next pair of glasses.

Select Single/BiFocal



Starburst



Hold Button



To Measure Progressive Lens Glasses, select the top symbol in the above and refer to the user's manual to understand what you will see on the display. As before press PRINT to transfer the reading to the P.C.

[Go to page 1](#)

Lensmeter Data Entry Help

PAGE 3 12/28/2006

Reichert AL200 Lensmeter Help

Connect Kendall Optometry Ministries, Inc part number 30005 cable between the lensmeter and the P,C. The default settings from the factory are sufficient to operate this unit with this program.

Humphrey LA360 Lensmeter Help

Connect Kendall Optometry Ministries, Inc part number 30005 cable between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.

Detailed step-by-step operating instructions are not available for this unit. There are setup instructions available and they will be provided upon request.

Tomey TL-900 Lensmeter Help

Connect Kendall Optometry Ministries, Inc part number 30005 cable between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.

Detailed step-by-step operating instructions are not available for this unit. There are setup instructions available and they will be provided upon request.

Tomey TL-3000 Lensmeter Help

Connect Kendall Optometry Ministries, Inc part number 30005 cable between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.

Detailed step-by-step operating instructions are not available for this unit. There are setup instructions available and they will be provided upon request.

Humphrey LA-340 Lensmeter Help

Connect Kendall Optometry Ministries, Inc part number 30001 cable between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.

Proper operation of the lensmeter.

When you measure a pair of glasses you must make sure that the instrument displays RIGHT or LEFT on the panel when you do a STORE using the either the foot pedal or STORE button. You will see something like the below on the screen when (for example) you measure the right lens. Realize that you can cause the RIGHT light to come on by pressing the edge of the left lens while you are doing your measurement. Make sure you also see ADD when measured (see below/right). Press the STORE button or foot pedal after each operation.



Store Button

Press the PRINT button when you have measured both the RIGHT and the LEFT lens doing a STORE after each operation. 8 seconds later the computer will receive and process the data. The 8 second pause is due to the fact that the lensmeter is delaying the amount of time that it previously took to print data to the printer.

PRINT button




Nikon NL-2 Lensmeter Help

Connect Kendall Optometry Ministries, Inc part number 30004 cable between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.

For the following usage instructions refer the drawing to the right for the location of the HOLD and ADD buttons on the front of the instrument.

For bifocals and progressive lens glasses.

1. Put glasses in holder right lens first and with bottom of glasses against the holding plate.
2. Adjust until you get a centered position.
3. Press the ADD button to measure the power of the ADD segment.
4. Move glasses until you are viewing through the ADD part (BiFocal segment) of the lens.
5. Press HOLD button.
6. Press Right button until green right light comes on.
7. Press HOLD button again to free up instrument.
8. Put glasses in holder left lens first and with bottom of glasses against the holding plate.
9. Adjust until you get a centered position.
10. Press the ADD button to measure the power of the ADD segment.
11. Move glasses until you are viewing through the ADD part (BiFocal segment) of the lens.
12. Press HOLD button.
13. Press left button until red left light comes on
14. Press PRINT button and information is transmitted to the computer.



For single vision glasses.

1. Put glasses in holder right lens first and with bottom of glasses against the holding plate.
2. Adjust until you get a centered position.
3. Press HOLD button.
4. Press Right button until green right light comes on.
5. Press HOLD button again to free up instrument.
6. Put glasses in holder left lens first and with bottom of glasses against the holding plate.
7. Adjust until you get a centered position.
8. Press HOLD button.
9. Press left button until red left light comes on
10. Press PRINT button and information is transmitted to the computer.



The picture below on the right shows a completed reading where the PRINT key is being depressed to send the information to the computer.

Topcon CL-100 Lensmeter Help

Connect Topcon cable part number ZZ-1-T013308 between the lensmeter and the P.C. The default settings from the factory are sufficient to operate this unit with this program.

Detailed step-by-step operating instructions are not available for this unit. There are setup instructions available and they will be provided upon request.

Below is what you see when you click the **File** option in the Glasses Reader program. The numbers have been added as reference points for explanation. Click the Jump links to go to selections on different pages.

File

1. Load Inventory File from Default Location.
2. Browse and load / merge other Inventory Files.
3. Load Deleted & Rejected inventory from Glasses Inventory program. [JUMP PAGE 2](#)
4. Load Update inventory into Vacant Locations. [JUMP PAGE 3](#)
5. Transfer inventory to Glasses Inventory Program [JUMP PAGE 4](#)
6. Create Vacancy List File for Label Generation. [JUMP PAGE 5](#)

Selection 1: (Load Inventory File from Default Location)

When you choose selection number one you will load the inventory into memory from the default location (C:\My Documents\Glasses_Data) as defined in the options screen entry as shown below.

Path to collection files.


C:\My Documents\Glasses_Data\ Total_Inventory.txt

OR→ No File Output or input. Check to read file into memory at startup.


Notice that the checkmark in the above indicates that this load has already happened when Glasses Reader started.

Selection 2: (Browse and load / merge other Inventory Files.)

Selection number two allows you to combine the inventory built by several different people. For example one person may create records 1-1000, another 1001-2000 and the third 2001-3000. When you choose this option, you will get the following windows:



If you are loading inventory records 1001-2000 you will change the window to read:




You will then browse to the inventory file that you are loading, click OPEN and it will be loaded into memory.

Notice that this option is described in a separate help file which describes how to build your inventory.


Selection 3: (Load Deleted & Rejected inventory from Glasses Inventory Program).

The Glasses Reader program is used to build the inventory while the Glasses Inventory Program is used to dispense the inventory. This option provides an interface between the two programs. If on the Glasses Inventory program a person has used 199 pairs of glasses, these eyeglasses need to be removed from the inventory maintained by the Glasses Reader program. Here are the steps involved:


On Glasses Inventory select the first option:




and the following popup will appear. This indicates that the deleted inventory file is now in C:\



On Glasses Reader if you select the vacancy report from the below.




you will then see:




Notice that zero locations are empty.


If you use Glasses Reader selection 3 and load the Deleted & Rejected inventory file from the Glasses Inventory program you will see:



Click Open and you will in the status window:



and the vacancy report will now show:



Note that if any of the locations being deleted have previously already been filled with some more inventory, this new inventory will not be deleted by this procedure.

Selection 4: (Load Update inventory into Vacant Locations).

After the inventory removed by the Glasses Inventory program has been removed from the file maintained by the Glasses Reader program this selection can be used to introduce new inventory into the system and in the process label with inventory with appropriate labels.

While the mission team is on their trip, people at the church are busily inputting new inventory into a separate file. This file should start numbering at 1 and represents the "Update Inventory File". When the team returns and the deleted inventory is removed from the Glasses Reader file, this option can be used to easily and quickly remove the inventory. New labels can also be produced for the update inventory eyeglasses bags and temples. Here are the steps involved.

Select the option and you will see:

These same two numbers (circled) are retained in the options window.

Click Continue & Select File and you will see:

Click OPEN and you will see the status window to the left.

The first number is set to be +1 more than the last number used by this option. The second number should correctly indicate the maximum number of locations in your inventory system.


and you will see in the status window:

Note that in the Default location you will also find a file called "Vacancy_list.txt" which shows exactly how the update inventory has been used. To the far right is an excerpt of that file. This says that update number 1 went into location 22, 2 into location 39, etc.

Vacancy file



File	Edit	For
1,	22	
2,	39	
3,	40	
4,	41	
5,	76	
6,	80	
7,	86	
8,	90	
9,	105	
10,	110	
11,	113	
12,	129	
13,	136	
14,	138	
15,	145	
16,	149	
17,	157	
18,	167	
19,	170	

A spreadsheet called "Optical Labels.xls" provided with this program allows you to print labels to go over the existing update labels. To the right is a sample of one of those labels. In **Red** is the previous update number (10101) and in **black** is the new location number.




Selection 5: (Transfer inventory to Glasses Inventory Program).

When you have completed building your inventory, you now need to transfer this inventory over to the Eyeglasses Inventory Program. If you select this option you will see the following window. Read carefully what it says. In most cases both options are selected and the transfer takes place. This will appear in the status window upon completion.



You may instead receive this error message:



Note that any previous inventory and journal file kept by the Glasses Inventory program has been backed up.

if you have changed your inventory but not saved it.

NOTICE


In some cases the Glasses Inventory files are kept on an entirely different computer. To handle this condition the Glasses Reader program will execute a command file (if exists) called: MoveInventory.bat located in C:. The likely contents of this file might be:

```
@echo off
Echo Place a blank floppy in drive A:
pause
Copy /Y c:\total_inventory.txt a:\
Copy /Y c:\journalfile.txt a:\
```

This will place the files on a floppy. When you get to the target PC both files should be copied to C:\.

Selection 6: (Create Vacancy List for Label Generation).

After having used option #3 (Load Deleted and Rejected inventory from Glasses Inventory Program), there are a block of vacancies within the inventory managed by the Glasses Reader program. The purpose of this option is to create a file called “vacancies.txt” which can be used by the Excel spreadsheet to print labels for glasses which will be placed in these locations. After selecting this function, you will see the following graphic:




If the maximum number of locations in your inventory is correct, click continue. Otherwise, correct the number and click continue. The vacancies file will then be created. Click CANCEL if you have changed your mind.

When you have completed this function launch the Excel Spreadsheet named “Print Update Labels with Excel” from:

START | Programs | Glasses Reader

And you will see a button which looks like the below.



Click this button to print the labels from the vacancies.txt file created by running this option.

Below is what you see when you click the **Reports** option in the Glasses Reader program. The numbers have been added as reference points for explanation. Click the Jump links to go to selections on different pages.

Reports	Options	Help
1	Error Report	
2	Rejected Item Report	← JUMP PAGE 2
3	Vacancy Report	← JUMP PAGE 3
4	Gender Report	← JUMP PAGE 4
5	Inventory Analysis Report	← JUMP PAGE 5
6	Clear glass report (Run after above)	← JUMP PAGE 5 (Bottom)

Selection 1: (Display the error list report from a previous sorted inventory save.)

Here is a sample of a report:

```
----- Numbers of Erroneous Glasses in Inventory. -----  
-----  
As of: 10/7/2005 2:35:47 PM  
10092 14973  
  
Number of errors = 2
```

Location 14973 Error

Status:

Right: Add too high!

Date/Time: 5/7/2005 10:34:29 AM

If you change the inventory number on the screen to display the two numbers, you will see the error below for location 10092. The error to the right is for location number 14973.

Location 10092 Error

Left (OS)

Sphere:	4.25	Both
Cylinder:	-10.75 *	Both
Axis:	37	Both
Add:	3.50	Both

Right (OD)

Sphere:	11.50	Both
Cylinder:	-7.75	Both
Axis:	3	Both
Add:	7.75 *	Both

Selection 2: (Display the rejected item report from the Eyeglasses Inventory Program)

-Inventory Numbers and Occurrences of Rejected Glasses in Inventory. --

As of: 10/7/2005 2:55:36 PM

0076 01	0121 01	0129 01	0163 01	0164 01	0268 01	0271 01	0288 01
0290 01	0298 02	0304 01	0305 01	0310 01	0324 02	0458 01	0490 01
0546 01	0567 01	0568 01	0585 01	0586 01	0591 01	0784 01	0844 01
0856 01	0885 01	1092 01	1106 01	1141 01	1210 02	1254 01	1262 01
1293 01	1298 01	1456 01	1477 01	1480 01	1484 01	1566 01	1584 01
1586 01	1604 01	1673 01	1683 01	1692 01	1695 01	1757 01	1763 01
1853 01	1864 01	1928 01	1985 01	2142 02	2211 01	2284 01	2305 01
2318 01	2561 01	2568 01	2582 01	2640 01	2693 01	2722 01	2739 01
2787 01	2803 02	2819 01	2872 01	2881 01	2914 01	2932 01	2968 01
3037 01	3062 01	3225 01	3327 01	3331 01	3332 01	3526 01	3566 02
3576 01	3599 01	3616 01	3684 01	3743 01	3747 01	3779 01	3867 03
3884 01	3908 01	4014 01	4059 01	4083 01	4111 01	4247 01	4261 01
4275 01	4284 01	4292 01	4324 01	4371 01	4379 01	4380 02	4408 02
4421 01	4427 01	4430 01	4470 02	4484 01	4536 01	4567 02	4596 02
4609 01	4615 01	4626 01	4650 01	4664 01	4751 01	4891 01	4896 01
4899 01	4927 01	4967 01	5017 01	5022 01			

Number of rejections = 125

Please note that a “rejection” is defined as a pair of glasses which were the right prescription for the patient but the patient could not see properly through the glasses. It is not a pair of glasses which simply were not liked by the patient. This is a pair which you suspect may be improperly measured. Generally when occurrence count for a pair of glasses goes as high as 4 or 5, you remove it from the inventory as likely it has been improperly measured.

Realize a patient can reject a pair of glasses that are properly measured because they do not physically fit or the focal center of each lens is not positioned properly on the pupil of the patient (PD is not right). Other patients may find this pair to be totally acceptable.

Selection 3: (Display the list of vacant locations in inventory.)

Here is the list of vacant locations in inventory.

```
----- Open Location Numbers in Inventory. -----  
As of: 10/7/2005 2:58:29 PM  
0022 0039 0040 0041 0076 0080 0086 0090 0105 0110  
0113 0129 0136 0138 0145 0149 0157 0167 0170 0177  
0197 0202 0210 0219 0262 0268 0277 0282 0290 0295  
0297 0312 0314 0337 0375 0386 0391 0399 0406 0438  
0454 0469 0474 0489 0490 0511 0512 0520 0542 0547  
0581 0584 0593 0596 0601 0606 0608 0624 0625 0627  
0635 0642 0658 0660 0666 0669 0714 0721 0738 0748  
0790 0800 0822 0836 0843 0896 0902 0908 0917 0948  
0988 0992 1029 1050 1055 1124 1172 1187 1188 1211  
1214 1233 1239 1242 1245 1259 1271 1275 1283 1285  
1315 1326 1327 1476 1480 1520 1525 1547 1599 1608  
1609 1617 1647 1651 1665 1680 1686 1687 1735 1736  
1778 1788 1856 1860 1868 1874 1895 1912 1956 1963  
1999 2004 2028 2029 2032 2053 2068 2117 2126 2165  
2192 2219 2270 2282 2329 2336 2340 2414 2415 2420  
2426 2429 2438 2443 2468 2477 2497 2546 2562 2577  
2602 2621 2643 2735 2755 2780 2794 2795 2798 2828  
2829 2847 2854 2871 2889 2899 2904 3013 3020 3026  
3027 3033 3059 3078 3081 3161 3187 3205 3252 3388  
3422 3522 3536 3604 3632 3846 3849 3939 3959
```

Number of Openings = 199

Selection 4: (Display the number of eyeglasses with a specific gender.)

This is a typical display of the gender of the eyeglasses in inventory.

<p>GENDER OF EYEGLASSES IN INVENTORY As of: 10/7/2005 2:49:07 PM Eyeglasses for: Male = 2344 Female = 4833 Either = 1938 Child = 97 With a Total Count of 9212</p>
--

Selection 5: (Display the Inventory Analysis Report.)

This report provides a detailed analysis of the inventory. This report is very computer intensive as it will scan the entire inventory and indicate to you how many pairs of glasses you have with a certain range of sphere, cylinder or axis. Below is a sample of this report:

Circled in **Red** in the Total column you see an indication that the highest number of glasses are in the spherical range from **+1.25** to **+1.25**.

Circled in **green** you can see that most of your glasses have NO cylinder (<.25).

Looking at the **black** circled items you can see that most of your glasses have axis which are WTR (With The Rule) which means from 0-30 and 150-180 or on a horizontal axis.

Looking at the **orange** circled items you can see that you have nearly 80 pairs of NON BiFocal glasses which have clear lens on one side or the other or both. The **yellow** alert indicates you should look at the clear glasses report to see which pairs are clear on both lens. There are some glasses in your inventory which are clear glass and are only useful for eye protection.

Below is a copy of the clear glass report for this inventory. You should probably remove this list of glasses.

```

----- Clear Glasses in Inventory. -----
      As of: 10/22/2005 10:08:58 PM
0399 0843 1265 1285 1756 1970 2133 2330 2363 2613
2634 2877 3178 3542

Number of clear glasses = 14
    
```

[Go to Page 1](#)

Click **Display to PRINT** to produce a copy of this report in notepad that you can print.

Inventory Analysis Report As of: 10/22/2005 9:54:30

Left Eye		Cylinder Counts					Axis Counts		
Sphere values from A to B.		NO Cyl (<=.25)	Low Cyl .26->.75	Med Cyl .76->2.75	High Cyl (> 2.75)	Total	ATR	WTR	Oblique
+5.26	+1000	2	7	16	2	27	15	10	2
+4.26	+5.25	9	19	17	2	47	25	20	2
+3.26	+4.25	49	51	23	6	129	60	60	9
+2.26	+3.25	140	90	43	13	286	109	156	21
+1.26	+2.25	290	115	76	12	493	164	289	40
+0.26	+1.25	378	214	181	22	795	307	416	72
-0.25	+0.25	261	211	202	20	701	282	345	74
-1.25	-0.26	266	224	156	11	657	269	313	75
-2.25	-1.26	144	83	60	8	295	89	180	26
-3.25	-2.26	77	79	43	1	200	51	133	16
-4.25	-3.26	60	50	27	3	140	33	92	15
-5.25	-4.26	31	37	16	2	86	13	61	12
-1000	-5.26	45	26	33	6	110	14	84	12
Totals: ==>>>		1752	1206	900	108	3966	1431	2159	376

=> Clear glass has been found! Click to run Clear Glass report. <=

Display to PRINT
CLOSE

Right Eye		Cylinder Counts					Axis Counts		
Sphere values from A to B.		NO Cyl (<=.25)	Low Cyl .26->.75	Med Cyl .76->2.75	High Cyl (> 2.75)	Total	ATR	WTR	Oblique
+5.26	+1000	4	3	14	3	24	12	10	2
+4.26	+5.25	10	14	11	3	38	22	15	1
+3.26	+4.25	47	70	27	10	154	86	55	13
+2.26	+3.25	162	84	34	9	289	102	166	21
+1.26	+2.25	255	125	67	11	458	145	277	36
+0.26	+1.25	372	225	203	21	821	319	428	74
-0.25	+0.25	264	228	210	21	723	307	350	66
-1.25	-0.26	261	205	182	15	663	275	334	54
-2.25	-1.26	142	74	48	7	271	75	171	25
-3.25	-2.26	92	68	29	0	189	48	124	17
-4.25	-3.26	56	50	27	5	138	30	90	18
-5.25	-4.26	37	29	25	7	98	22	68	8
-1000	-5.26	39	28	28	5	100	18	63	19
Total: ==>>>		1741	1203	905	117	3966	1461	2151	354

Options Help

PAGE 1

5/22/2007

This help covers the part of the screen which is shown below. Click on blue boxed or circled areas or underlined text for the part of the screen for which you need help or information.

PAGE 2

PAGE 2

PAGE 2C

PAGE 2C

PAGE 2B

PAGE 2B

PAGE 3

PAGE 3

PAGE 4

PAGE 4

PAGE 4B

PAGE 4B

PAGE 5

PAGE 5

Set Data Collection Options. V3.0-x+

Help with options. (Click for Help)

Identify your Lensmeter

- Reichert LensChek
- Humphrey LA350/360
- Tomey TL-3000
- Reichert AL200
- Topcon CL-100
- Nikon NL-2
- Humphrey LA340
- Tomey TL-900
- Barcode Reader
- Keyboard input

Select COMM Port number

None 1 2 3 4 5 6 7 8 9 10 11 12

Path to collection files.

C:\My Documents\Glasses_Data\ Total_Inventory.txt

OR-> No File Output or input. Check to read file into memory at startup.

Printing Control

Choose Your Printer: Adobe PDF

Printing

ON OFF

Print After Measurement.

ON OFF

Show Inventory # on printout?

YES NO

Choose Printout and Printer Form

- Detailed text plus bar codes.
- Brief large text plus bar codes.
- Print labels for 2 uncut lens
- Brief Large text only - No barcode

30252 Address Label

Other Options

- Make ADDs Equal
- Add Minus "-" to cylinder
- Force ADDs to be PLUS.
- Right first then left on screen.
- Check enables automatic reading glasses detection. (NOT saved)

Default Axis for single Lens entry. 90 NOT saved.

Allowed Maximums and Minimums

MAX Sphere	MAX Cylinder	MAX Add	MIN Add	Maximum Number of locations.
18	8	6	.75	32000

Historical Information

Previous Next Number 00000 Last Used Update Number 00000

UnProtect UnProtect

Save to File and continue Continue

Select the model of your lensmeter or if you have no lensmeter select either Barcode Reader or Keyboard input.

Identify your Lensmeter

- Reichert LensChek
- Humphrey LA360
- Tomey TL-3000
- Reichert AL200
- Topcon CL-100
- Nikon NL-2
- Humphrey LA340
- Tomey TL-900
- Barcode Reader
- Keyboard input

Select the COMM Port number to which your Lensmeter is attached. If you have no lensmeter select None.

Select COMM Port number

- None
- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8
- 9
- 10
- 11
- 12

The default path can be changed but not the file name. If you click on the circled item, no file will be produced. If you click the item in a square then the file will be loaded into memory upon startup of the program.

Path to collection files.

C:\My Documents\Glasses_Data\ Total_Inventory.txt

OR→ No File Output or input. Check to read file into memory at startup.

Click here to turn printing on or off. PRINT button will appear or disappear to/from the main window.

Click this dropdown to choose your label printer from printers configured on your system.

Click this dropdown to select the type of label material you are using (if any).

Select your style of label here.

The numbers of the sample labels below correspond to the numbers of the selections above.

Click YES here so that the inventory number will appear on the top/right part of every label printed.

Turn ON this to automatically print a label after completion of a lensmeter measurement. If you do this, you should select your flags before entering data from the lensmeter.

1

OD Sph	Cyl	Axis	Add	OS Sph	Cyl	Axis	Add	PD
-0.75	-1.25	126		-0.50	-1.25	102		MD
[Either M/F (5/22/2007)]								
OD: -0.75 -1.25 126 FLAGS 5/22/2007								
OS: -0.50 -1.25 102 MD								
(Manually)								

2

OD:	-0.75	-1.25	126	PD:	MD
OS:	-0.50	-1.25	102	[Either M/F (5/22/2007)]	
(Manually)					

3

SPH:	-0.75	CYL:	-1.25	[1]	- R
(Manually)					
SPH:	-0.50	CYL:	-1.25	[1]	- L
(Manually)					

Above are two labels.

4

OD:	-0.75	-1.25	126	PD:	MD
OS:	-0.50	-1.25	102	[Either M/F (5/22/2007)]	
(Manually)					

Click this option to automatically put the same Add in both the Left and Right positions. If this is wrong the duplicated Add can be corrected.

Select this option so that every time you enter a number into a cylinder field a "-" sign is automatically added at the first. This save you the time required to enter a minus sign when all your cylinder values are minus.

Other Options

<input checked="" type="checkbox"/> Make ADDs Equal	<input checked="" type="checkbox"/> Add Minus "-" to cylinder
<input checked="" type="checkbox"/> Force ADDs to be PLUS.	<input type="checkbox"/> Right first then left on screen.
<input type="checkbox"/> Check enables automatic reading glasses detection. (NOT saved)	

Default Axis for single Lens entry.

90 **NOT saved.**

Click to reverse the position of the Right and Left areas of the main window.

This is used when you are measuring glasses and they are zero cylinder and the sphere is plus and the same for both sides. (Rarely Used)

When scanning in barcoded labels from lens which have been measured by Glasses Reader, enter here the default AXIS assumed for any lens entered which have some amount of cylinder.

Do not allow the entry of minus (-) Adds.

Enter here the maximum allows Sphere, Cylinder and Add.

Enter here the smallest add which can be entered without flagging an error.

Allowed Maximums and Minimums

MAX Sphere	MAX Cylinder	MAX Add	MIN Add	Maximum Number of locations.
18	8	6	.75	32000

Should the above be exceeded, you will see an error message in the status window (below) and a large red asterisk (*) to the right of the suspect field.

Enter the maximum number of locations you ever expect your inventory file to contain. 32000 is the architectural maximum of this program.

Click this button to save the options which you have chosen to file and then continue the operation of the Glasses Reader program.

Click this number to continue with the currently chosen options without saving them to file for use next time the program starts.

Historical Information	
Previous Next Number	Last Used Update Number
15	00000
<input type="checkbox"/> UnProtect	<input type="checkbox"/> UnProtect


Save to File and continue

Continue

Click UnProtect and the number will no longer be grey. You can then enter a new "Next Number" which will be put in the main window "Inventory Number" when it next starts up.

Click UnProtect and then the number will no longer be grey. You can then enter a new "Last Used Update Number". This number +1 will be the next update inventory number used by the system.

Currently it is recommended that the Dymo labelwriter series be used for printing the labels generated by this program. Here is the section of the options windows which is used to setup the Dymo Labelwriter.



The option which is circled above is used to select the printer on your system to which labels will be printed. Click the dropdown to select the right printer. This prevents you from having to always print to just the default printer. Be sure that the printer that you select is not greyed out (offline) in your printer folder.

The second areas (in a square box) allows you to select the type of labels. To the right if each selection you are allowed to select the type of label. There are three different types of Dymo paper involved.

Part number	Label Description	Program Setting
30256	White Shipping Labels	Brief Large Text Plus Barcodes and Detailed Text plus Barcodes
30252	White Address Labels	Brief Large Text only - No Barcode
30270	White Continuous Wide	All three types of labels.

NOTICE




Please understand that as of version 3.0-x of the Glasses Reader program, the selection of the **Type** of label stock does not consistently work. It is therefore necessary to manually set this label type on the Dymo Labelwriter printer icon. The following pages describe the steps to perform to select the paper type used by the Dymo Labelwriter printer.

[Go to first page of printer setup.](#)

Printer Paper Setup – Phase 1

1. Right click the printer icon and select **Properties** and you will see the window #1 to the right.
2. Click Printing Preferences (circled) and you will see window #2.
3. Click Advanced and you will see window #3 below.
4. Click the dropdown window (circled) and select the correct type of printer paper.
5. Click OK twice to return to the main window.




[Click to go to the next page.](#)



Printer Paper Setup – Phase 2

6. Now click the Advanced tab (circled) in Window # 4.
7. Next click the Printing Defaults (circled below) and you will see Windows #5.
8. Click the Advanced button (circled) and you will see Window #6.
9. Click the dropdown (circled) and select the type of printer paper which you wish to use.
10. Click OK twice to return to the original window.

[Click to go to the next page.](#)



Printer Paper Setup – Phase 3

11. Click the Device Settings tab (circled) and you will see window #7
12. Click the paper selection (circled) and then use the dropdown to select the correct type of paper.
13. Click OK and the printer is ready for use.

