

PicoLog

Data Logger Software

User Manual

Version 2.0 rev 3

by

M K Green

Pico Technology Limited

149-151 St Neots Road

Hardwick

Cambs CB3 7QJ

GB

Tel: 44-1954-211716

Fax: 44-1954-211880

E-mail: 100073.2365@compuserve.com

Copyright 1994-1995 Pico Technology Limited

Contents

1 Introduction	1
2 Basic principles	2
3 Quick Guide	5
Installing the PicoLog program	5
Setting up for your computer	5
Collecting data for the first time	5
Repeating a run	6
Reports from an old run	7
A simple example	7
Setting up the system	8
Checking the ADC	9
Collecting data	9
Monitoring during the run	10
Stopping the run early	11
Additional reports on a completed run	11
A more complex example	13
Checking the ADC	14
Sampling details	14
Specifying the channels	15
Thermistor Scaling	16
Flow meter scaling	18
Power- a derived parameter	19
Tabulation report	20
X vs Y graph	21
Notes	22
Collecting the data	22
Spreadsheet files	23
Repeating the run	23
4 Running the program	25
Menus	25
Forms	26
Graphic Reports	27

5 Main menu	28
6 Setup menu	30
Graphics screen	30
Interface options	31
Graphics Printer details	32
Text Printer Details	33
Data and backup paths	33
7 Display voltages	35
8 Data collection menu	36
Define new run	37
Repeat existing run	37
Start collecting	38
Display report	39
Stop collecting	39
Edit sampling rate	40
Add or Edit channels	41
Add or Edit parameters	44
Linear scaling	46
Table lookup scaling	48
Equation scaling	50
Add or Edit reports	51
Monitor report	52
Tabulation report	53
Y vs Time Report	55
Edit X vs Y report	58
Edit Notes	60
9 Generate reports	61
Add or Edit parameters	61
Add or Edit reports	61
Edit notes	61
Display report	62
Print report	62
Write text report to disk	63
Display run setup	63

- 10 File management 64
 - Load old run from disk 64
 - Display or print index 64
 - Delete one run 64
 - Copy one run 65
 - Backup all runs 65
 - Extract records 66

1 Introduction

The PicoLog data logger is a powerful and flexible program which enables you to carry out data logging and analysis. It is designed to make full use of the facilities offered by the Pico range of data acquisition products.

To use PicoLog, you will need an IBM compatible PC and a Pico analog to digital converter. If you wish to print reports, you will also need an HP Deskjet or Laserjet printer, or an Epson FX or LQ printer.

This manual starts off with an overview of the functions provided by PicoLog, then it explains how to set up a couple of typical data collection runs. Finally, there is a detailed description of each of the menus and forms in PicoLog.

The content of some forms depends upon the type of ADC: see the manual for your ADC for a detailed description of the special features of PicoLog for your ADC.

We have tried to use the following typographical conventions throughout the manual:

Bold text is used for terms which have some special meaning within the manual, when they first appear, and for menu and field names:


- Ⓒ ... has one or more **channels**
- Ⓒ the **Setup menu**

Courier text is used to indicate user input and also menu options, for example:

- Ⓒ ... enter a filename of TEST1
- Ⓒ select the *Display Voltages* option from the...

2 Basic principles

The PicoLog data logger collects sets of measurements from the channels of an analog to digital converter (ADC) and stores them on disk. It also provides tools to convert the ADC values to give parameters scaled in engineering units and to produce text and graphical reports.



An analog to digital converter has one or more **channels**: each channel is used to measure a voltage which is a representation of a real-world parameter like pressure or temperature.

The program takes a number of **readings** from a channel and combines them to form a **sample**. During each **sampling interval**, the program can do this on several channels. The program collects samples at regular time intervals (the sampling interval). It collects a specified maximum number of samples, then stops recording: you can stop it recording at any time

before the maximum number of samples. A sequence of samples is referred to as a **run**: each run is stored in a separate file on disk.


The program can scale the values for a channel to produce a **parameter**: a parameter is a real-world measurement derived from the channel values for a sample. You can use **linear** scaling, **table lookup** or an **equation** to convert the channel values into real-world units.

- C Linear conversion is suitable for use with sensors which have a linear response (for example strain gauges and thermocouples)
- C table lookup or equation scaling must be used when the sensor output is non-linear, for example a thermistor.
- C Equation scaling can be used both for scaling channel data and for combining parameters- for example, you could measure two temperatures, then use an equation to evaluate the temperature difference.

Note: you can also add the values from a channel directly onto a report, if the values are already scaled correctly. There are also a number of special parameters (sample number, time since start of run, etc) which are generated automatically by the program.

A **report** is a text or graphical representation of the data for a run. Text reports can be displayed, printed or written to disk: graphic reports can be displayed or printed.

Graphical report



Text report

Indoor Temp degC	Outdoor Temp degC	Wind Speed km/h	Humidity %RH
16.2	18.0	1492	3890
16.3	18.5	1400	3193
16.3	18.1	1431	3113
16.3	17.8	1883	3080
16.2	17.6	1439	3044
16.2	17.6	1254	2988

Before, during or after collection of data, you can define a number of parameters and reports. During recording, you can display a report to

monitor what is happening: once the run is complete, you can print out reports and output text reports to a disk file.

Each repeat run will be stored in a separate file. The run file contains all details about the run: this includes channels, sampling details, data values, parameters and reports. You can read in a run file at a later date, then display or print existing reports, and new parameters, reports or notes. You can also repeat the run to obtain another set of data, then generate the same reports as for the first run.

3 Quick Guide

The rest of the manual explains in detail the operation of each menu and form in the system. First, here is a quick guide to the things that you are most likely to want to do.

3.1 Installing the PicoLog program

Insert the distribution diskette into drive A and type in `a:install`. The computer will first ask you what directory to install into, then it will offer you a menu of modules to install (drivers, PicoScope, PicoLog et cetera). If you just want to install PicoLog, make sure that there is `Yes` next to PicoLog and `No` next to the other modules.

3.2 Setting up for your computer

Use each of the options on the **Setup menu** to define the type of screen, display, printer and ADC, and the default locations for data and program files.

3.3 Collecting data for the first time

All of the options that you require to start a new run are on the **Collect menu**.

First, select the `Define new run` option: the computer will first ask for a filename for the new run, then the computer will ask for the sampling details (duration, number of samples), then the channels to collect data from. Once you have specified a filename, the computer will automatically update the file each time you change any of the run details and at regular intervals whilst collecting data.

You may wish to define the parameters to be derived from the data you collect (for example to convert ADC readings to temperatures) and then to define text or graphical reports to display this data while you are collecting data. Use the `Add or edit Parameter` and `Add or edit Report` options.

Next, select `Start collecting`. If you selected a sample interval less than 50ms, the computer will just display a message saying `collecting...` until all samples have been taken. For lower sampling rates, the computer will display the collecting menu, so that you can display reports or modify parameters and reports whilst you are collecting data. For sampling rates faster than one per second, changing the reports and parameters whilst collecting may interfere with data collection.

If you wish to terminate the run before the computer has collected the maximum number of samples, select the `Stop collecting` option. The computer will not let you out of this menu until you have stopped collecting.

3.4 Repeating a run

Once you have successfully carried out a run, it is very easy to repeat the run with the same parameters and reports.

If you exited from the program since you collected the data, you must first read in the run that you wish to repeat, using the `Load old run from disk` option on the **File management menu**.

Once you have loaded the run that you wish to repeat, all of the functions that you require are on the **Collect menu**.

Next, select `Repeat existing run` and specify a new filename for the run. Once you have done this, you can either start collecting with exactly the same samples, parameters, and reports, or you can first edit the details using the `Add or edit...` options.

You start the repeat run exactly the same as the first run, and all subsequent steps are the same.

Use the `Edit Notes` option to annotate reports with information about each particular run. You can attach the notes to reports: this will also appear if you request a file index for a directory.

3.5 Reports from an old run

If you wish to examine or print reports for an old run, you should first read in the run using the `Load existing run` option on the **File Management menu**. If you are not sure which run you require, use the `File index` option to check the notes, date and time and number of records in each of the files in a directory.

Once you have loaded a run, all of the functions that you require are on the **Generate reports menu**.

You can define new parameters and reports, then to print or display the reports. You can also write text reports to disk in a format suitable for entry into a spreadsheet or a database.

3.6 A simple example

The following is an example of a simple experiment, in which you wish to record the voltage across a battery as it discharges in a circuit using an ADC-12.

After you have installed PicoLog on your computer (see section 3.1), you need to set up the details about your system.

1. Type `picolog` and press `Enter`. the **main menu** appears.

3.6.1 Setting up the system

1. Select `Setup` (display, ADC, printer, paths). The **setup menu** appears.

2. To use a colour display, select the Graphics Screen (mono/colour) option.
3. Press the space bar change the screen type from Black and White to Colour.
4. Press F10 to accept this change.

Note: this will not take effect until you next run the program. If you want to start using colour immediately, Press ESC twice to exit from the program, type in picolog again and then repeat steps 1 and 2.

5. Select Interface details (ADC type, port) the **SELECT ADC** form appears.
6. Press the space bar until ADC-12 appears.
7. Press F10 to accept this change. The **EDIT SETUP** form appears.
8. Press the space bar until the required printer port name appears.
9. Press F10 to accept this printer port: the **Setup menu** appears.

Note: There are three more options on the setup menu. If you do not need printed reports and all of your data files will be in the current directory, it is not necessary to use these options. We will ignore them for the purposes of this example.

3.6.2 Checking the ADC

1. Press `Esc` to return to the main menu
2. Select `Display Voltages`. the **DISPLAY VOLTAGES** form appears.

if you have set up the interface correctly, the computer will display the ADC value and voltage. Connect your battery and check that the voltage reading is correct.

3. press `Esc` to return to the main menu.

3.6.3 Collecting data

1. Select `Collect Data (new or repeat run)`: the **Collect menu** appears.
2. Select `Define new run`: the **DEFINE NEW RUN** form appears.
3. Type in a filename of `TEST1`
4. Press `F10` to accept the filename. The **EDIT SAMPLING RATE** form appears.
5. Set **sampling interval** to `1` and **units** to `minutes`: this will take one sample every minute. If the maximum time to discharge the battery will be 8 hours, set the maximum number of samples to `480` (8 hours * 60 minutes).
6. Press `F10` to accept the sampling details. The **EDIT CHANNEL** form appears.
7. Set the **Channel name** to `Battery Voltage` and change the **Measurement mode** to `DC Volts`: leave other fields as they are.
8. Press **F10**: the computer returns to the **Collect menu**.

9. Select **Start collecting**: the computer will say **press any key to start**.

10. Press a key: again, the computer goes back to the data collection menu. The number of samples in the top right window will increase once per minute.

3.6.4 Monitoring during the run

If you wish to see what data you are recording, you must first define a report. You can do this whilst the computer is still collecting.

1. Select `Add or edit report`: the **EDIT REPORTS menu** appears. This shows any reports that you have already defined and the option `<NEW>`.

2. Select `<NEW>` because we wish to create a new report: the **EDIT REPORT** form appears.

3. Set the **report type** to `monitor`. The report name is optional.

4. Press `F10`. The **EDIT MONITOR REPORT** form appears.

5. Move the cursor to the first parameter, then press the space bar until `sample no` is selected.

6. Move the cursor down, then press the space bar till `Battery Voltage` is selected.

7. Press `F10` to save these details. The **Collect menu** appears again.

8. Select `Display report`: the **Display report** menu appears.

9. Select the `monitor` report you have created: the computer will now display the current sample number and the voltage.

3.6.5 Stopping the run early

When the battery is discharged (the voltage has reached a suitably low value) you can stop the run early:

1. Press `Escape` to return to the **Collect menu**
2. Select `Stop collecting` to stop the run.

3.6.6 Additional reports on a completed run

You can define reports before, during or after you have collected the data. To add a graphical report once the run is completed:

1. Select `Add or edit report`. the **Edit Reports menu** appears.
2. Select `<NEW>` again.
3. Set the name to `Battery voltage vs time`
3. Change the **Report type** to `Y vs time`.
4. Press `F10`: the **EDIT Y VS TIME** form appears.

Note: the **title** for the graph has been set automatically to `Battery voltage vs time` (the name you entered for the report). You can change it to something different if you wish.

6. Move down to the left axis **units** field and enter `Volts`.
7. Set **min Y** to `0` and **max Y** to the nearest whole number above the starting battery voltage (say `4`).
8. Move down to the first `parameter` field and press the space bar until the parameter name is `Battery voltage`.
9. Press `F10` to save the details: the **Collect menu** appears.
10. Select `Display report`: the computer displays a graph of battery voltage against time.


3.7 A more complex example

We shall now work through a complete example which demonstrates most of the functions of PicoLog. Suppose that you wish to carry out a number of experiments to determine the efficiency of a boiler, using various different fuels and burners.

The boiler has two thermistors attached- one to measure the inlet temperature and one to measure the outlet temperature. It also has a flow rate sensor which reports the amount of water flowing through the boiler. All three sensors are connected to an ADC-11.

Each experiment will last five minutes, in the course of which the pump will be manually adjusted to run at five different rates for about a minute each.

Here is a diagram showing the setup:



3.7.1 Checking the ADC

First, we wish to check that the sensors are working correctly.

1. From the **Main menu**, select `Display voltages` option: the computer displays the ADC counts and voltages for each of the channels on the ADC-11.

Channel 1 changes when we increase the inlet temperature, Channel 2 changes when the outlet temperature changes and the frequency of Channel 3 changes when the pump is turned on. See section 3.7.4 for how to check that the values are accurate.

2. Press `Escape` to return to the main menu.

3.7.2 Sampling details

1. Select `Collect data (new or repeat run)`: the computer displays the **Collect menu**.

2. Select `Define new run`. The computer displays the **DEFINE NEW RUN** form.

3. Type in a **Filename** of `BOILER1`

Note: when we repeat this run for subsequent tests, the computer will automatically increment the filename- `BOILER2`, `BOILER3` et cetera.

4. Press `F10` to accept the filename: the computer displays the **EDIT SAMPLING RATE** form.

5. Set the **sampling interval** to `1`

6. Set the **units** to `seconds`

7. Set the **number of samples** to `300` (300 seconds is 5 minutes).

8. Press F10 to accept: the computer now displays the **Edit channel menu**.

3.7.3 Specifying the channels

1. Select **Channel 1**: the computer displays the **Edit channel form**.

2. Set the **Channel name** to `Inlet Thermistor`

3. Set the **measurement mode** to `ADC counts`.

4. Press **F10** to accept these details: the computer displays the **Edit channel menu**.

Note: the entry for channel 1 in the channel menu now says `Inlet Thermistor`. You can select this option again at any time to change the details for this channel.

5. Select **Channel 2**

6. Set the **channel name** to `Outlet thermistor`

7. Set the **measurement mode** to `ADC counts`.

8. Set the **number of readings** to `100`: the outlet temperature could vary a lot as the heated water will be at various different temperatures, so it makes sense to take lots of readings and then average them.

9. Press F10 to accept these details and return to the **Edit channel menu**.

10. Select **Channel 3**

11. Set the **Channel name** to `Flow meter`

12. Set the **measurement mode** to `Frequency`: the flow meter generates pulses whose speed depends on the flow rate. The frequency measurement is effectively an average, so only one reading will be enough.

13. Press `F10` to return to the **Edit channel menu**

13. Press `Escape` to return to the **Collect menu**.

3.7.4 Thermistor Scaling

We now need to define some parameters to convert the ADC readings into real-world values (degrees C, litres per second et cetera).

1. Select `Add or edit parameters` the **Edit parameters menu** appears.

2. Select `<NEW>`: the **EDIT PARAMETER form** appears.

3. Set the **Name** to `Inlet temperature`,

4. Set **Heading lines 1 and 2** to `Inlet` and `Temp`

5. Set **Units** to `degC`

Note: the degrees symbol `E` appears on the screen, but may not appear on all printers. 'deg' will be printed on all printers.

6. Set the **Scaling method** to `Table lookup`

7. Press **F10**: the **EDIT TABLE LOOKUP SCALING** form appears.

8. Set **Channel** to `Inlet thermistor`

9. We now need to do some sums. The thermistor input is a 30k ohm Betacurve device connected in series with a 30k resistor to 2.5V. The voltage input to the ADC-11 is $(2.5 * R_t)/(R_t+30k)$. The following table shows the resistance of the Betacurve resistor at temperatures from 0 to C, and the resulting ADC values. and enter the values from columns 4 and 1 into the form as **raw** and **scaled** values

Temperature (degrees C) SCALED	Resistance form datasheet	Voltage $2.5 \cdot R_t /$ $(R_t + 30k)$	ADC counts $V \cdot 1023 / 2.5$ RAW
0	105.3k	1.946	796
10	62.1k	1.686	690
20	38.1k	1.399	572
30	23.8k	1.106	453
40	15.3k	0.844	345
50	10.1k	0.630	258
60	6.78k	0.461	189
70	4.65k	0.335	137
80	3.24k	0.244	100

10. Press **F10** to save the table lookup values and return to the **Collect data** menu.
11. Select the `Add or edit parameters` menu option
12. Select the `<NEW>` menu option.
13. Set the **name** to `Outlet temperature`
14. Set **Heading lines 1 and 2** to `Outlet` and `Temp`
15. Set the **units** to `deg C`
16. Set the **Scaling method** to `Table lookup`
17. Press **F10** to save the details and go to the **EDIT TABLE LOOKUP SCALING** form.

18. Set the **Channel** to `Outlet thermistor`
19. Set the **raw** and **scaled values** to the same values as for the inlet temperature.
20. Press `F10` to save the scaling information and return to the **Collect data** menu.

3.7.5 Flow meter scaling

1. Select the `Add or edit parameters` menu option
2. Select the `<NEW>` menu option.
3. Set the **name** to `Flow Rate`
4. Set **Heading lines 1 and 2** to `Flow and Rate`
5. Set the **units** to `l/min`
6. Set the **Scaling method** to `Linear`
7. Press **F10** to save the parameter details: the computer will display the `EDIT LINEAR SCALING` form.
8. Set the **raw parameter** to `Flow Meter`.
9. Time for some more calculations. According to the data sheet, the flow rate sensor outputs a frequency which is 25Hz per litre per minute. Set the **offset** to 0 (litres per minute) and the **slope** to $1 / 25$, which is 0.04 (litres per Hz).
10. Press `F10` to return to the **Collect menu**.

3.7.6 Power- a derived parameter

The last parameter to specify is the power, which we calculate from the temperature difference and the flow rate.

1. Select the `Add or edit parameters` menu option
2. Select the `<NEW>` menu option.
3. Set the **name** to `Power`
4. Set **Heading line 1** to `Power`
5. Set the **units** to `kW`
6. Set the **Scaling method** to `Equation`
7. Press `F10` to go to the **EDIT EQUATION SCALING** form.
8. Move the highlight to **parameter A** and press the space bar until **Inlet temperature** appears.
9. Likewise, set **B** to `Outlet temperature`
10. Set **C** to **Flow rate**.
11. Type in an **Equation** which is $(B-A) * C * 4.192 / 60$: (B-A) is the temperature difference produced by the boiler, C/60 is the number of litres per second and 4.192 is the number of kiloWatt seconds required to heat a litre of water by one degree C.
12. Press `F10` to save the scaling information and return to the **Collect menu**.

3.7.7 Tabulation report

We are now ready to define some reports. We require a tabulation report to output the data to a spreadsheet, and a graph of power versus flow rate

to show how well the boiler works at each flow rate. Either report can be used to monitor progress whilst the experiment is running.

1. Select `Add or edit report` form the **Collect menu**: the **EDIT REPORT menu** appears.
2. Select `<NEW>`: the **EDIT REPORT** form appears.
3. Type in a **name** of Spreadsheet - just to remind your self what the report is for when you see it on menus.
4. Select the **Format** `Tabulation report`
5. Press **F10** to save these details: the **EDIT TABULATION REPORT** form appears.
6. There are two columns of fields below the line that says **Parameter**. Move the cursor to the first **Parameter** field, then select **Flow rate**.
7. Move to the next **Parameter** field and select **Power**.
8. Press **F10** to save the details and return to the **Collect Menu**.

3.7.8 X vs Y graph

1. Select `Add or edit report` form the **Collect menu**: the **EDIT REPORT menu** appears.
2. Select `<NEW>`: the **EDIT REPORT** form appears.
3. Type in a **name** of Spreadsheet - just to remind your self what the report is for when you see it on menus.
4. Select the **Format** `X vs Y graph`
5. Press **F10** to save: the **EDIT X VS Y REPORT** form appears.

6. Type in a **Title** of Power vs Flow Rate.
7. Move to **Notes** and select `Notes 1`: we can use notes 1 for details of each run on the graph.
8. Move to **Mark points with shapes** and select `Yes`
9. Move to **Join points** and select `No`: for this graph, the connection between points is less important than their position on the XY graph.
10. Move to the **X axis Param** field and select `Flow rate`.
11. Move to the **X axis Min** and **Max** and type in 0 and 50 respectively: this gives a range of 0 to 50 degrees C on the X axis.
12. Move to the **Y axis Param** field and select `Power`.
13. Move to the **Y axis Min** and `Max` and type in 0 and 5000. This gives a range of 0 to 5000 watts on the axis.
14. Press `F10` to save these details and return to the **Collect menu**.

3.7.9 Notes

To enter some notes for the first run:

1. Select `Edit notes` on the collect menu: this takes you to the **EDIT NOTES** form.
2. Move to the first line of **notes 1** and type in
burner type: Egremont AZ35
3. Move to the second line of notes 1 and type in
Fuel: MAXOIL 230
4. Press `F10` to save these details.

3.7.10 Collecting the data

1. Select `Start recording` on the **Collect menu**: the computer says

Press any key to start

2. press a key. The computer starts recording and returns to the **Collect menu**. You can see the sample number at the top of the screen increase once per second.

3. Select **Display report**

4. Select either **Spreadsheet** or **Power graph**: the computer will display the report you requested, and will update the report every time a new sample is taken.

5. Press `Esc` to exit from the report and return to the **collect menu**.

Once the full five minutes has elapsed, the computer will stop collecting data.

3.7.11 Spreadsheet files

Once the run is completed, you can write the spreadsheet to disk:

1. Press `Esc` to return to the **main menu**.

2. Select `Generate reports (old run)`.

3. Select `Write text report to disk`: the computer displays the **DISK REPORT** menu.

4. There is only one report that can be written to disk: select `Spreadsheet`.

3.7.12 Repeating the run

Once you have entered the details for a run, it is very easy to repeat the run.

1. Press `Esc` to return to the **Main menu**
2. Select `Collect data (new or repeat run)`
3. Select `Repeat current run`: the computer asks for a filename for the new run: it offers the old filename plus one (`BOILER1 + 1 = BOILER2`) as a default.

Note: if you have just completed a run, PicoLog will assume that you wish to repeat the same run. If, however, you exit from PicoLog between runs, the computer will ask you which run you wish to repeat at this point.

5. Press `F10` to save this filename.
6. Select **Edit notes**: the computer displays the **EDIT NOTES form**.
7. Move to the first line of **notes 1** and type in
`burner type: Gloucester AP20`
8. leave the second line of **notes 1** still saying
`Fuel: MAXOIL 230`
9. Press `F10` to save these details and return to the **Collect menu**.
10. Select `Start Collecting`
11. Proceed as for the first run.

6 Setup menu

The setup menu is usually used only when you are installing the system. It enables you to define the hardware and software configuration for the logger program.

```
?4444444444U SETUP MENU W4444444444@
* Graphics screen (mono/colour) *
* Interface details (ADC type, port)*
* Graphics Printer details *
* Text Printer details *
* Data and backup paths *
B4444444444444444444444444444444444A
```

6.1 Graphics screen

```
?4444U EDIT SCREEN SETUP W444@
* *
* Graphics display: Colour *
* *
B4444444444444444444444444444444444A
```

This option enables you to specify how graphs will appear on your screen. If you have a colour display, select `colour`: if you have a monochrome display, try both of the `black on white` and `white on black` options, then decide which looks best on your computer.

Note: that PicoLog uses CGA screens in two colour mode (black and white) to get maximum resolution.

This also affects the use of colour on text displays. No harm will come from using the wrong type, but some menus or graphic symbols may be difficult to see or invisible. If you run into trouble (for example the menus become invisible), delete the file **PICOLOG.SET**.

6.3 Graphics Printer details

```
?4444U EDIT PRINTER SETUP W4444@
*                                     *
* Printer:           Deskjet          *
* Port:              LPT1             *
*                                     *
* Resolution:        High             *
* Orientation:       Landscape        *
*                                     *
B44444444444444444444444444444444A
```

This menu option is used to specify the type of printer to be used for graphics output. You can also specify that you wish to write graphics output to a disk file (Warning! graphic files can be very big). it also enables you to select the quality of graphic reports.

If you have one of the printers listed, just select that option. If not, check your printer manual to see whether it is compatible with one of these types. It will do no harm to experiment, though you may end up with either no output or 30 pages of garbage. If you get really stuck, call Pico: we may be able to help.

In addition to the standard graphics drivers supplied with PicoLog, there is a disk containing additional graphics drivers for other printer types, plotters and word processor graphics formats.

The printer **port** can be one of LPT1, LPT2 and File.

The **resolution** can be high, medium or low. High resolution looks good in reports, but it can take a very long time to print. Low resolution is quicker to print, and its chunky appearance works well on overhead projector slides.

The **orientation** can be Landscape, Portrait or Half portrait. Landscape usually gives the best results.

directory, but it also offers the option to change to other directories or drives.

The **backup path** is used when you select the `Copy all runs to backup path` option on the **File management** menu: this option copies all run files in the data path to the backup path.

When you select this option, the computer offers the current run filename with a number at the end incremented (for example, TEST1 becomes TEST2). You can change to a different name if you wish. If you enter a filename that already exists, the program will ask you to confirm that you wish to overwrite the existing file.

This option clears the samples for the existing run, but retains all channel, parameter and report information.

8.3 Start collecting

This option is used to collect data. Before using this option, you must have entered the following details:

```
C    run filename
C    sampling rate
C    channels
```

If all of the details have been entered, the computer will display the following window:

```
?44444U COLLECT DATA W44444@
*                                     *
*Press any key to start             *
*                                     *
B444444444444444444444444444444A
```

When you are ready to start recording, press a key. If the sampling rate is too fast to update reports at the same time, it will say **too fast for reports** and the **collect data** window will remain in place until all data has been collected. If the sampling rate is slow, it will return to the **collect data** menu: you can then edit and display reports whilst data is collected. The number of samples will appear in the window at the top right of the screen, and will be updated unless a report is active.

To stop the run before all samples have been taken, select the **Stop collecting** option on the **Collect data** menu.

8.4 Display report

```
?444444U DISPLAY REPORT W44444444@
*Temp diff          Y vs Time Graph *
*All temps          Tabulation      *
*Alarms             Monitor         *
B444444444444444444444444444444A
```

When you select this option, the program gives you a menu containing the name and type of all of the reports suitable for display whilst data is being collected. Move the highlight to the report that you wish to display, and press Enter.

Note: Monitor reports do not appear on this menu once the run is completed.

8.5 Stop collecting

This option is used to terminate sampling before the maximum number of samples has been collected. Once you have started a run, the program will not let you out of the **Collect menu** until you have stopped the run.

8.6 Edit sampling rate

```
?444444U EDIT SAMPLING RATE 4444444444@
*                                     *
* Sampling interval:      1          *
* Units:                  Second     *
* Number of samples:     500        *
*                                     *
* Action at end of run:  Stop        *
*                                     *
B444444444444444444444444444444A
```

This option is used to set the interval between samples and the maximum number of samples to be taken. To select as sample every 5 minutes, set the **sampling interval** to 5 and the **units** to **Minutes**. The **sampling interval** must be a whole number, so if you wish to sample at half-minute intervals you must specify 30 seconds rather than 0.5 minutes.

The **units** can be milliseconds (ms), seconds, minutes, hours and days.

The **number of samples** defines the maximum number of samples to be taken: during data collection, you can stop sampling at any time before the maximum.

You must ensure that the **sampling interval** is greater than the total sampling time which appears at the bottom of the edit adc-.. channel forms.

Action at end of run defines what to do when all of the samples have been collected. You can select either `Stop`, which stops collecting, or `Repeat`, which start a new run. It automatically increments the filename, so that each run file is successively numbered.


If, for example, you use a thermocouple amplifier which outputs the temperature in degrees C at 10mV per degree C, with an ADC-11, the calculation is as follows.

The ADC-11 outputs an ADC value of 1023 at 2.5 volts (2500mV). so the slope is $2500/(10*1023)= 0.2444$. If, for some reason, the thermocouple registers 3 degrees C when immersed in ice, it needs an offset of -3 degrees. The offset is therefore -3 and the slope is 0.2444.

The maximum and minimum ADC values are as follows.

Note: bear in mind that Pico ADCs are typically 1% accuracy- see spec for more details.

ADC	Min voltage	Max voltage	ADC Range
ADC-10	0	5	0..255
ADC-11	0	2.5	0..1023
ADC-12	0	5	0..4095
ADC-16	-2.5	2.5	8 bit: -255..255 9 bit: -511..511 10 bit: -1023..1023 11bit: -2047..2047 12 bit: -4095..4095 13 bit: -8191..8191 14 bit: -16383..16383 15bit: -32767..32767 16bit: -65535..65535
ADC-22	0	2.5	0..1023
ADC-100	-range	+range	0..4095
TC-08	n/a	n/a	n/a



The raw and scaled values are used to convert the raw channel value from ADC counts into a real-world parameter value. The computer selects the raw values immediately above and below the raw value which must be scaled, then uses linear interpolation between the two corresponding scaled values to estimate the parameter value. You can specify up to 10 raw/scaled value pairs.

Note: next time you edit the values, the entries will have been sorted into increasing raw value order.

The raw and scaled values can either be calculated from details supplied with the sensor, or can be determined empirically by measuring a range of parameter values and noting the corresponding raw values.

Note: the raw value will change if you change the resolution for a channel.

The **title** will appear on each page of the tabulation report, as will the heading details for each of the parameters that you select.

The **notes** option can be either `None`, `Notes 1` or `Notes 2`. This specifies whether any of the notes for this run are to appear in the heading. This is useful if you carry out a number of similar runs: the report title can be the same for every run, but the notes contain information specific to each run.

First sample and **Last sample** specify the range of samples that you wish to appear in the report.

One sample every is useful if you have a lot of samples: you can generate a report which contains only one sample every 10.

Max columns specifies how wide the report is to be: it only has any effect if you specify `right` or `centre` justify.

Justify affects the positions of the columns for each parameter across the page. It can be one of:

- C `left` - display the parameters at the left hand side of the page
- C `centre` - put the parameters in the middle of the page
- C `right` - spread the columns evenly across the page.


The **timestamp** provides a facility to display the time and/or date at which a sample was taken, at the start of each line. It can be `None`, `Date`, `Time` or `Date and time`.

Set one or more entries in the **parameter** column to the parameters that you wish to appear on the report.

If the run is in progress or completed, press `F3` to display the report.

8.9.3 Y vs Time Report

A y vs time report is a graph of up to ten parameters against time. Here is a typical example:



Min and **max x** specify the range to be displayed on the graph, in the units of the sampling interval (eg for a sampling interval of 500ms, the units are milliseconds). Once the run is completed, you can alter the x scale whilst displaying the graph.

The **format** defines the format for the scale on the X (time) axis. It can be `linear` (for example 0..5000 minutes) or `date/time` (for example 16:20 to 21:45 or 03Jul92 to 12Jul92).

Note: the min and max X must still be specified in linear time, even if you select `date/time` format for the x axis scale.

The graph can have an axis scale on each side, so that parameters measured in different units can be displayed on the same graph. The parameters can then be assigned to either the left or the right axis. If all of the parameters have the same scale, only one axis is necessary.

The two **enabled** fields say whether an axis is to be displayed on the corresponding side. If, for example, all of your readings are temperatures, you can use the same scale for all parameters. If you want to display a voltage and a temperature use both scales.

Label and **units** are added as notes to the corresponding axis.


Min y appears at the bottom of the axis and **max y** appears at the top of the axis.

The four **parameter** options specify which parameters are to appear on the graph: unused parameters should be left at `<NONE>`. You can specify that the parameters are scaled using either the `left` axis or the `right` axis.

If data collection is completed, you can press `F3` to display the report. Press `F10` to save the details.

8.9.4 Edit X vs Y report

An X vs Y graph is a plot of one parameter against another. A typical graph looks like this:



9.4 Display report

```
?444444U DISPLAY REPORT W4444444@
*Temp diff      Y vs Time Graph *
*All temps      Tabulation      *
B444444444444444444444444444444A
```

When you select this option, the program gives you a menu containing the name and type of all of the reports suitable for display once a run is completed. Move the highlight to the report that you wish to display, and press **Enter**.

Note: Monitor reports do not appear on this menu once the run is completed.

If you display a graphical report, you can alter the scale of the graph and select which portion is displayed: press **F1** for help while the graph is displayed. See section 4.3 for a description of this function.

9.5 Print report

```
?444444U PRINT REPORT W4444444@
*Temp diff      Y vs Time Graph *
*All temps      Tabulation      *
B444444444444444444444444444444A
```

When you select this option, the program gives you a menu containing the name and type of all of the reports suitable for printing. Move the highlight to the report that you wish to display, and press **Enter**.

9.6 Write text report to disk

```
?444444U DISK REPORT W4444444@  
*All temps      Tabulation      *  
B44444444444444444444444444444444A
```

This option is used to write text to a disk file, so that you can import it into a word processor, spreadsheet or database.

When you select this option, the program gives you a menu containing the name and type of all text format reports. Move highlight to the report that you wish to display, and press `Enter`.

9.7 Display run setup

When you select this option, the program displays details of all channels, parameters and reports for this run. You can use it to check on the setup. If you press `F5`, the program will print out a hardcopy version of the setup.

10 File management

```
?44444U FILE MANAGEMENT W4444444@
* Load old run from disk          *
* Display or print index          *
* Delete one run                  *
* Copy one run                    *
* Copy all runs to backup path    *
* Extract records                 *
B444444444444444444444444444444A
```

10.1 Load old run from disk

When you select this option, the program displays a menu of all of the run files in the data directory defined in the setup. It also displays all of the directories in this directory and all of the drives. You can either select one of the files in the current directory, or select a different directory or drive to see what run files there are in that directory.

10.2 Display or print index

When you select this option, the program produces a report showing details for each of the files in the data directory. The details displayed are the filename, date, time, number of samples and the first line of notes 1.

Index

AC Volts	42	tabulation	51
tabulation	53	X vs Y graph	52
ADC	31	y vs time	56
alarm	52	Y vs time graph	52
channel	2	resolution	32
channel name	42	run	3
colour	30	copy	65
combination method	43	delete	65
data path	34	index	64
dB	42	load	64
decimal places	45	save	37
display values	35	setup	63
extract records	66	sample	2
field width	45	sampling interval	2, 40
Flow Meter	18	scaling	3, 45
form	26	equation	3, 45, 50
graphics printer	32	linear	3, 45, 46
graphics screen	30	lookup	48
interface	31	table lookup	3
measurement mode	42	setup	30
menu	25	spreadsheet	53
new run	37	text printer	33
notes	54, 56, 59, 60	thermistor	16
number of readings	43	timestamp	54, 57
number of samples	40	units	43, 44
orientation	32	Volts	42
parameter	3	word processor	32
parameter name	44		
port	31, 32		
raw parameter	18		
reading	2		
repeat	37		
report	3, 39, 61		
display	52, 62		
file	63		
graphical	27		
Monitor	51, 52		
print	62		