

USER MANUAL

LCD Array Kit

SPECIES 2.0

DNA-based identification of
Animal species

Code: A-750-04

Code: A-750-12

1	Product information	3
	1.1 Intended use	3
	1.2 Test principle	4
2	Storage conditions & kit content	5
3	Equipment and reagents - not supplied with the kit	6
4	Preparation of reagents	6
5	Warnings and precautions	7
6	General safety information	7
7	Method description	8
	7.1 PCR	8
	7.2 Hybridisation / Detection	8
8	Array description	9
9	Protocol	10
	9.1 PCR amplification	10
	9.1.1 PCR set-up	10
	9.1.2 Template	10
	9.1.3 Cyclor settings	11
	9.1.4 Agarose gel examination	11
	9.2 LCD Array hybridisation and detection	11
	9.2.1 Hybridisation	12
	9.2.2 Labelling	13
	9.2.3 Staining	13
	9.3 Short protocol	14
10	Analysis and interpretation of the results	15
	10.1 General remarks	15
	10.2 Methods of analysis	15
	10.2.1 Analysis by simple optical examination with naked eye	15
	10.2.2 Analysis by transmission light scanning and SlideReader software	15
11	Device performance	16
	11.1 Sensitivity	16
	11.2 Determination of device performance	17
12	Device limits	18
13	Trouble shooting	18
14	Literature	19
15	Order information	20
16	Contact manufacturer	20
17	Symbol key	20

Instructions for use

Instructions version: SPECIES 2.0 V-1.0-2012-ENG Revised : 26-July-2012

1 Product information

Product name: LCD Array Kit SPECIES 2.0

Order code	Kit size
A-750-04	32 Tests / 04 Chips
A-750-12	96 Tests / 12 Chips

1.1 Intended use

This array has been developed for rapid, easy and reliable identification of animal DNA in fresh meat preparations and products manufactured thereof. DNA of animal species most commonly used in food production will be detected in parallel using extracted DNA as starting material.

Beef	<i>Bos taurus</i>	Mallard Duck	<i>Anas platyrhynchos</i>
Buffalo	<i>Bubalus bubalis</i>	Muscovy Duck	<i>Cairina moschata</i>
Pork	<i>Sus scrofa</i>	Dog	<i>Canis lupus (familiaris)</i>
Sheep	<i>Ovis aries</i>	Cat	<i>Felis catus</i>
Goat	<i>Capra hircus</i>	Kangaroo	<i>M. giganteus / M. rufus</i>
Horse	<i>Equus caballus</i>	Ostrich	<i>Struthio camelus</i>
Donkey	<i>Equus asinus</i>	Bison	<i>Bos bison</i>
Rabbit	<i>Oryctolagus cuniculus</i>	Red Deer	<i>Cervus elaphus</i>
Hare	<i>Lepus europaeus</i>	Fallow Deer	<i>Dama dama</i>
Chicken	<i>Gallus gallus</i>	Roe Deer	<i>Capreolus capreolus</i>
Turkey	<i>Meleagris gallopavo</i>	Camel	<i>Camelus sp.</i>
Goose	<i>Ansa albifrons</i>	Seagull	<i>Larus sp.</i>

List of species which can be detected by hybridisation of specific capture probes.

1.2 Test principle

Using extracted DNA as starting material (1), small fragments of the mitochondrial 16S rRNA genes of the target organisms will be amplified by the initial PCR. During this amplification the generated PCR fragments are labelled with Biotin (2).

PCR amplicons are hybridised to specific capture probes on the surface of the array (3). Non-specifically bound PCR amplicons will be removed by short washes under high stringency. The remaining, specifically bound amplicons can be visualized by incubation with a Streptavidin-Peroxidase conjugate and the subsequent formation of a dark precipitate after incubation with peroxidase substrate (TMB) provided as staining solution (4).




Figure 1: Test principle

2 Storage conditions & kit content

When stored at the indicated temperature and not handled otherwise as specifically described, all components are stable until their expiry date. Expiry dates are printed on the components label or the labels of the surrounding packages. More than 6 freezing and thawing cycles should be avoided for components which require frozen storage. Prepare aliquots if necessary.

Component	A-300-04 (32 tests)	A-300-12 (96 tests)	Storage
LCD-Arrays	1 box, 04 chips	3 boxes, 12 chips	+4°C to +28°
Wash Powder	1 bottle	2 bottles	+4°C to +28°
Primer			
MEAT Primer Mix	1 x 100 µl	2 x 100 µl	-10°C to -25°C
Detection Kit I			
Modulator	1 x 300 µl	2 x 300 µl	+4°C to +28°
Dilution Buffer	1 x 2000 µl	2 x 2000 µl	+4°C to +28°
Stain	1 x 2000 µl	2 x 2000 µl	(*) +4°C to +28°
Detection Kit II			
Hybridisation Buffer B	1 x 1750 µl	2 x 1750 µl	-10°C to -25°C
Label	1 x 20 µl	2 x 20 µl	-10°C to -25°C
Additional supplies			
Analysis Matrix			---
Manual, Pattern File, MSDS	1 CD	1 CD	---
Chip-Box Connectors	2 connectors	6 connectors	---



(*) = Keep the ‚Stain‘ solution always dark, protected from direct light.

3 Equipment and reagents - not supplied with the kit

Instrumentation

- Thermal cycler
- Water bath (adjustable to 35°C)
- Centrifuge for LCD-Arrays (cat. no HS-500-01, Chipron GmbH)
- Micro pipettes (range from 2 µl to 1000 µl)

optional:

- SlideScanner PF3650u (cat. no HS-300-01, Chipron GmbH)
- SlideReader Analysis Software (cat. no HS-200-01, Chipron GmbH)

Reagents & Materials

- Reagents for DNA extraction
- PCR chemicals (Taq Polymerase, Buffer, dNTPs)
- PCR grade water
- Deionised water
- Disposable gloves
- Sterile filter tips
- PCR reaction vessels
- 3 wash containers (150 to 400 ml each)
- 1l-bottles

4 Preparation of reagents

Wash Buffer

To prepare the 20 x wash buffer concentrate dissolve the content of one wash powder bottle in 1 liter of deionised water. To prepare the 1 x wash buffer add 50 ml of 20 x wash buffer concentrate to 950 ml of deionised water. Prepare the 1 x wash buffer always fresh.

Stability of Wash-Buffer solutions:

- | | |
|----------------------|--|
| 20 x Concentrate | 2 month at room temperature, 12 month at +4°C |
| | If formation of precipitate occurs, warm the solution to 42°C and let equilibrate to room temperature prior to use |
| 1 x Working solution | min. 10 days at room temperature |

5 Warnings and precautions

- All steps described in the protocol should only be carried out by well trained lab personnel.
- Read the manual carefully and completely before starting.
- Avoid any exposure to light of the 'Stain' solution. Always keep it dark.
- All reagents in the kit are optimised for this particular test. Substitution of kit reagents may effect the performance.
- The same general safety guidelines which apply to the sample material should be followed during the whole protocol.
- The PCR products generated in the first protocol step have to be considered as contamination sources for further PCRs. Therefore, all hybridisation, washing, staining, drying and analysis steps should be carried out in the 'post-PCR' area.
- Observe the standard guidelines for working in a PCR molecular diagnostic laboratory to prevent contaminations.

6 General safety information

When working with chemicals, make sure that you always wear a suitable lab coat, disposable gloves and protective goggles. For more information about our products, please refer to the appropriate material safety data sheets (MSDS) which can be found on the CD provided with the kit.

Hybridisation Buffer

contains formamide (>50%) and N-Dodecanoyl-N-methylglycin-Sodium salt.

Hybridisation buffer should therefore be handled as formamide.

Toxic, harmful, irritant

Risk and safety phrases: R61, R36, R38, S24, S26

Wash Powder

contains N-Dodecanoyl-N-methylglycin-Sodium salt.

Harmful, irritant

Risk and safety phrases: R36, S24, S26

Stain

contains 3,3',5,5'-Tetramethylbenzidin (> 0.5 %).

Harmful, irritant

Risk and safety phrases: R20/21/22, R36/37/38, R40 (Self Assessment)

Modulator

Harmful, irritant

Risk and safety phrases: R36/37/38 S26/36 (Self Assessment)

7 Method description

The test consists of two main procedure steps.

- PCR amplification of DNA fragments with biotin incorporation
- Hybridisation of PCR fragments to LCD-Arrays and detection

7.1 PCR

One primer mix is provided with the kit. The multiple primer pairs in this mix target a portion of the 16S RNA gene of a broad range of Vertebrate species. When using the primer mix 'MEAT' distinct amplification products of 115 – 125 bp can be expected (species dependent):

Examples of fragment sizes:	Species	Size
	<i>Bos taurus</i>	116 bp
	<i>Sus scrofa</i>	119 bp
	<i>Equus caballus</i>	124 bp
	<i>Gallus gallus</i>	121 bp
	<i>Struthio camelus</i>	122 bp


7.2 Hybridisation / Detection

The labelled PCR fragments are combined with the hybridisation buffer (provided) and hybridised to the individual array fields of one chip. During hybridisation the labelled PCR fragments will bind to the specific capture probes immobilized as spots on the bottom of each field. Following a short washing procedure, each field is incubated with a secondary label solution (enzyme-conjugate). After a second washing step, the positions (spots) where PCR fragments and secondary label are bound can be visualized as blue precipitate formed by the enzyme substrate provided as "STAIN". The data read-out can either be done by simple 'naked-eyed' examination, using the 'Analysis Matrix' provided with the kit or, alternatively, with the scanner and software from the "Analysis-Package" which can be obtained from Chipron GmbH.

For experienced users the whole procedure will take 3-4 hrs (depending on the duration of the PCR amplification).

8 Array description

Each LCD-Chip contains eight identical arrays in rectangular reaction chambers which can be addressed individually. Functional controls to monitor hybridisation, secondary labeling and staining are located in three corners. The arrays of the SPECIES 2.0 kit consist of an 8 x 8 pattern with average spot diameters of 350 µm. The species specific capture probes are positioned as vertical duplicates..



No.	Name	Species	No.	Name	Species
01	Hyb-Ctrl	Hybridization Control			
02	Beef ¹⁾	<i>Bos taurus</i>	14	Mall. Duck	<i>Anas platyrhynchos</i>
03	Buffalo	<i>Bubalus bubalis</i>	15	Musc. Duck	<i>Cairina moschata</i>
04	Pork	<i>Sus scrofa</i>	16	Dog ³⁾	<i>Canis lupus (familiaris)</i>
05	Sheep	<i>Ovis aries</i>	17	Cat ⁴⁾	<i>Felis catus</i>
06	Goat	<i>Capra hircus</i>	18	Kangaroo	<i>M. giganteus / M. rufus</i>
07	Horse	<i>Equus caballus</i>	19	Ostrich	<i>Struthio camelus</i>
08	Donkey ²⁾	<i>Equus asinus</i>	20	Bison	<i>Bos bison</i>
09	Rabbit	<i>Oryctolagus cuniculus</i>	21	Red Deer	<i>Cervus elaphus</i>
10	Hare	<i>Lepus europaeus</i>	22	Fallow Deer	<i>Dama dama</i>
11	Chicken	<i>Gallus gallus</i>	23	Roe Deer	<i>Capreolus capreolus</i>
12	Turkey	<i>Meleagris gallopavo</i>	24	Camel ⁵⁾	<i>Camelus sp.</i>
13	Goose	<i>Ansa albifrons</i>	25	Seagull ⁶⁾	<i>Larus sp.</i>

- 1) Cross reactivity of the capture probe for Beef with 100% Bison may occur
- 2) Cross reactivity of the capture probe for Donkey with 100% Horse may occur
- 3) Capture probe detects *Canis lupus*, *C. lupus familiaris*, *C. lupus chanco* and *Canis indica*
- 4) Capture probe detects several species from the subfamily Felinae (*Felis*, *Puma*, *Leopardus*, *Lynx*)
- 5) Capture probe detects both species of the Genus *Camelus* (*C. bactrianus* & *C. dromedarius*).
- 6) Capture probe detects several species from the Genus *Larus* (*L. argentatus*, *L. cachinnans*, *L. occidentalis*, *L. dominicanus*, *L. glaucoides*).

9 Protocol

9.1 PCR amplification

9.1.1 PCR set-up

We recommend the HotStarTaq Plus Master Mix Kit [Qiagen, Code: 203645]

Using a 2 x PCR Master Mix (incl. Taq-Polymerase, dNTPs, Buffer and MgCl₂)

Number of reactions (25 µl each)	1	4 (set up 5)	8 (set up 9)
2x Master Mix (incl. MgCl ₂ final 1.5-2.0 mM)	12.5 µl	62.5 µl	112.5 µl
Primer Mix 'MEAT'	1.5 µl	7.5 µl	13.5 µl
PCR grade water	6.0 µl	30.0 µl	54.0 µl
Total volume	20.0 µl	100.0 µl	180.0 µl
Aliquot into		4 x 20.0 µl	8 x 20.0 µl

Using Taq-Polymerase, dNTPs, and Buffer as separate components

Number of reactions (25 µl each)	1	4 (set up 5)	8 (set up 9)
10x PCR Buffer (incl. 15-20 mM MgCl ₂)	2.5 µl	12.5 µl	22.5 µl
dNTP mix (10mM each)	1.0 µl	5.0 µl	9.0 µl
Primer Mix 'MEAT'	1.5 µl	7.5 µl	13.5 µl
Taq Polymerase (5 U/µl)	0.3 µl	1.5 µl	2.7 µl
PCR grade water	14.7 µl	73.5 µl	132.3 µl
Total volume	20.0 µl	100.0 µl	180.0 µl
Aliquot into		4 x 20.0 µl	8 x 20.0 µl

9.1.2 Template

Add 5.0 µl of extracted DNA as template to each reaction tube.

Note: If larger or smaller volumes of template need to be added, reduce the amount of added water accordingly

9.1.3 Cyclers settings

The cycle regime given below has been optimised for the use of

- HotStarTaq Plus Master Mix Qiagen GmbH [Code 203645]
- TProfessional Standard Cycler Analytik Jena AG [Code 070-951]

When other combinations of PCR cycler and enzyme will be used, slight modifications of the given protocol could become necessary. Please contact your local distributor or the Chipron support team if you wish assistance for any kind of assay optimization.

Step	Duration	Temperature	Ramp	Note
1) Initial denaturation	5:00 min	95°C		longer with some "Hot-Start" enzymes
2) * 35 repetitions of:				
Denaturation	0:30 min	94°C	3°C / sec	
Annealing	0:45 min	57°C	3°C / sec	
Elongation	0:45 min	72°C	3°C / sec	
3) Strand completion	2:00 min	72°C		

* Increase cycle number for higher sensitivity, or decrease cycle number for lower sensitivity.

9.1.4 Agarose gel examination

When the arrays are used for the first time or optimisation becomes necessary, it is recommended to analyse the PCR amplification by agarose gel electrophoreses (2.0 % agarose gel).

9.2 LCD Array hybridisation and detection

Detailed protocol for first time user – advanced users refer to short protocol (chapter 9.3.)

General Remarks

Since the working principle of LCD-Arrays is DNA/DNA hybridisation, the specificity and sensitivity of the assays are mainly controlled by the hybridisation stringency during the 30 minute hybridisation step. Apart from the concentration of the interacting molecules, the two factors with the biggest influence on hybridisation stringency are temperature and buffer composition (e.g. concentration of salts, formamide, urea etc.). It is crucial that the temperature during hybridisation and the pipetting volumes are precisely controlled. The use of calibrated thermometers and micro pipettes is strongly recommended.

Deviations of more than 1 °C or 1 µl can have severe effects.

9.2.1 Hybridisation

- 1) Set water bath temperature to 35 °C and add one droplet of water to each corner of the humidity chamber. Do NOT preheat the LCD-Chip. The chip should be equilibrated to room temperature.
- 2) Prepare the Hybridisation Mix.
Make sure that all components are equilibrated to room temperature.

Number of reactions	1	4 (set up 5)	8 (set up 9)
Hybridisation buffer B	22 µl	110 µl	198 µl
Modulator	2 µl	10 µl	18 µl
Aliquot per*	-	24 µl	24 µl

* 0.2 ml reaction vials or 8-well PCR strips are well suited for setting up the reactions

- 3) Combine the PCR product with the Hybridisation Mix.
Add 10 µl of the PCR product to the respective aliquot of Hybridisation Mix. Mix well by pipetting up and down for several times – do not vortex!
- 4) Initiate the hybridisation.
Place the slide in the humidity chamber and pipette 28 µl of the PCR-Hybridisation mixes to the respective fields. Make sure that the PCR-Hybridisation mixes come into contact with the entire reaction zone of the respective array field – avoid any contact of the pipetting tip with the chip surface. Transfer the closed humidity chamber as quickly as possible to the preset water bath (the chamber will float on the surface).
- 5) Incubate the slide at 35°C for 30 minutes
- 6) Prepare three wash containers with 1 x wash buffer (~ 150 ml each).
Prepare the 1 x wash buffer from the 20 x wash buffer concentrate, by adding 50 ml of concentrate to 950 ml of deionized water. Mix thoroughly.
- 7) Submerge the slide completely in the first container with the wash buffer and move it 3 times slowly backward and forward. Submerge the slide quickly, to avoid field to field cross talk. Repeat the procedure in wash container 2. Subsequently, transfer the slide to container 3 and incubate for 1 minute
- 8) Dry the slide by spinning for 15 seconds in the CHIP Spin FVL2400N (Cat.No. HS-500-01).

9.2.2 Labelling

- 8) Prepare the labelling mix

Number of reactions	1	4 (setup 5)	8 (setup 10)
Dilution Buffer	27.0 μ l	135 μ l	270 μ l
Modulator	3.0 μ l	15 μ l	30 μ l
Label	0.2 μ l	1 μ l	2 μ l
Total	30.2 μ l	151 μ l	302 μ l

Mix well by vortex or intensive pipetting

- 9) Apply 28 μ l of the label mix to each field of the slide.
Make sure that the label mix has contact with the entire reaction zone of the respective array field – avoid any contact of the pipette tip with the chip surface.
- 10) Incubate the slide at room temperature for 5 min
- 11) Replace the wash buffer in all three containers and repeat the wash procedure from step 7).
- 12) Dry the slide as in step 8).

9.2.3 Staining

- 13) Apply 28 μ l of 'Stain' solution to each field and incubate for 5 minutes at room temperature. Make sure that the fields do not contain any traces of washing solution – wait until fields are completely dry.



Avoid any contamination of the 'Stain' stock solution! Prepare aliquots if necessary.

- 14) Stop the staining process by rinsing the slide in the last wash container from step 11) for 15 seconds.
- 15) Dry the slide as in step 8).

The slides are now ready for read out and can be kept in the dark for several years without a significant loss of signal intensity.

9.3 Short protocol

- 1) Set water bath temperature to 35°C
- 2) For each reaction mix 22 µl Hybridisation buffer **B**, 2 µl of Modulator and 10 µl of PCR product. Apply 28 µl of this solution to the respective array field.
- 3) Incubate the slide at 35°C for 30 minutes (in humidity chamber)
- 4) Prepare 3 wash containers with ~ 150 ml each of 1x wash solution.
(1 x wash buffer freshly prepared from 20x concentrate)
- 5) Rinse slide in wash container 1 and 2 for 10 seconds each and incubate in wash container 3 for 1 minute.
- 6) Dry the slide by spinning it for 15 seconds in the CHIP-Spin Centrifuge
- 7) Prepare the labelling mix by combining 27 µl of Dilution buffer, 3 µl of Modulator and 0.2 µl of Label per array field.
- 8) Apply 28 µl of label mix to each field of the slide and incubate for 5 minutes at room temperature.
- 9) Replace the wash solution in all containers and repeat the wash procedure from step 5).
- 10) Dry the slide as in step 6).
- 11) Apply 28 µl of Stain solution to each field of the slide and incubate for 5 minutes at room temperature.



Avoid contamination of the STAIN solution with residues of Label solution

- 12) Stop the staining after 5 minutes by rinsing the slide in the last wash container from step 9) for 10 seconds.
- 13) Dry the slide as in step 6) – the slide is now ready for analysis.

10 Analysis and interpretation of the results

10.1 General remarks

LCD-Arrays generate qualitative results indicating the 'presence' or 'absence' of the respective parameter (species) within the sample material used for DNA extraction and PCR amplification. Since the assay principle is based on DNA detection, vital and dead organisms will be identified likewise. 'Presence' and 'absence' in this respect is defined by the assay specific detection range and limits for each parameter (see chapter 11 Device performance, 11.1 Sensitivity).

Although different signal intensities can be observed during the analysis of LCD-Arrays and these intensities are generally correlated with the amount of target copies in the starting material, it should be noted that LCD-Arrays have not been validated as tools for absolute quantification.

10.2 Methods of analysis

The formation of dark visible precipitates at positions (spots) where DNA/DNA hybridisation between the PCR amplicons and the immobilized capture probes took place, combined with the crystal clear polymer support (LCD chip), offers the opportunity to use two different analysis methods. Simple optical examination (naked eye, with or without magnifying lenses) or transmission light scanning followed by software-assisted image analysis (details are given below). The experiments of all validations and performance tests of the product have been analysed with both methods in parallel with identical results.

Regardless of the method chosen for analysis, negative controls should always be included into the experimental scheme to ensure that no artificial background 'signals' are detected due to cross contamination or 'over staining' (see chapter 13, Trouble shooting).

10.2.1 Analysis by simple optical examination with naked eye

Each LCD-Array Kit contains a pattern graphic with the imprint of an array specific 'Analysis Matrix' for simple data analysis. The 'Analysis Matrix' displays the array pattern for the specific array type, a table with the spot number code and the respective spot identities (Genera or species).

A simple comparison of all visible spots on the array, with the corresponding positions given in the 'Analysis Matrix', will lead to the correct data interpretation.

10.2.2 Analysis by transmission light scanning and SlideReader software

An alternative method for data analysis including the generation of comprehensive data reports in PDF format is the use of the SlideScanner PF3650u in combination with the SlideReader Software.

Product	Order No.	Company
SlideScanner PF3650u	HS-300-01	Chipron GmbH, Berlin, Germany
SlideReader Software	HS-200-01	Chipron GmbH, Berlin, Germany

Detailed instructions for the use of SlideScanner and SlideReader software are given in the SlideReader Manual.

The pattern file for automatic analysis with the SlideReader software is provided as data file on the CD which is part of every LCD-Array kit. Choose the directory \ Pattern File and dependent on the SlideReader software version in use one of the two following files.

SlideReader Vers 6 to 9:	use the file	Gal-00377 SPECIES 2.0.txt
SlideReader Version 11 (and higher)	use the file	0377 SPECIES 2.0.gal
<u>Required software settings:</u>	Cut off:	set to 2000

Please contact your local distributor or the Chipron support team to learn more about the automatic analysis options with the SlideReader software.

11 Device performance

All DNA sequences underlying the design of primer and probes have been carefully chosen and thoroughly confirmed by broad comparisons in publicly available databases (BLAST analysis). No indication for cross reactivity apart from those as described below the specificity table on page 9 could be revealed by such *in silico* analysis.

It should be noted that, as for all PCR and hybridisation-based molecular tests, it can not be excluded that rare genotypes, comprising variations within the sequence regions used for the assay design, may remain undetected. Additionally, other, yet unknown / yet unpublished species (especially among the *Bovidae*) might contain highly similar sequence regions in their genomes and could therefore interfere with the assay performance. The assay in its present form reflects the current state of knowledge of Chipron GmbH.

11.1 Sensitivity

The sensitivity of the assay is mainly determined by the number of PCR cycles and the binding affinity of the capture probes during hybridisation. Under the given protocol conditions (hybridisation buffer, temperature and duration), the binding affinity of the capture probes is well standardised and allows for the detection of ~ 10 fmol of single-stranded, biotinylated target sequences. The number of PCR cycles, on the other hand, can easily be adjusted to modulate the sensitivity range according to the experimental demands. The maximal sensitivity of the SPECIES 2.0 LCD-Array kit for 16S rRNA sequences of the target organisms represented by specific capture probes is approx. 0.1 Genome Equivalent (GE) / µl and will be reached with 42 PCR cycles.

Note: Detection limits below 1 GE / reaction are not unusual for mitochondrial target genes like the 16S rRNA gene which are present in high copy numbers / cell (tissue dependent).

12 Device limits

The following factors might limit the assay performance:

- total concentration of extracted DNA is too low or DNA is degraded
- PCR amplification is hampered by PCR inhibitors which have been co-extracted
- use of PCR additives which can interfere with the hybridisation (avoid additives like Urea, DMSO, Betaine etc.)
- use of non calibrated or functionally impaired instruments (pipettes, water bath, PCR cycler)
- any variation from the protocol given in this user instruction

13 Trouble shooting

General remarks

In cases where unexpected results are obtained or when signal intensities or background levels appear to be different from usual observations, check the following list first

- kit components and PCR chemicals were not expired
- agarose gel examination of PCR products is positive
- correct results with positive and negative controls
- correct combination of primer mixes, chip types, hybridisation buffer and protocol

Detailed trouble shooting

Observation	No signals detected (including guide dots)
Probable causes	Wrong reaction set-up, wrong temperature, use of wrong hybridisation buffer .
Solution	Compare the experimental set-up with the given protocol, check temperature of the water bath, compare the letters printed on the hybridisation buffer label and the chip label (should be B) .
Observation	Guide dots detected, but no other signals visible
Probable causes	PCR failure, DNA extraction inefficient, template DNA degraded, DNA content too low or test result for target organisms is negative.
Solution	Use fresh PCR chemicals, make sure to run the controls provided with kit, check the PCR cycler program, try an alternative DNA extraction procedure.
Observation	Signals detected , but background is too high
Probable causes	Wrong reaction set up or incubation times too long

Solution	Make sure the Modulator has been added at all stages according to the protocol, monitor incubation times and temperatures for hybridisation, labelling and staining carefully.
Observation	'Bleeding' guide dots and very strong control signals
Probable causes	Staining or secondary labeling time is too long
Solution	Reduce the incubation time gradually in 1 minute steps
Observation	Lots of small dark particles on chip surface
Probable causes	1 x wash buffer is too old (formation of precipitates)
Solution	Prepare a fresh solution of 1 x wash buffer from 20 x concentrate
Observation	Negative control shows species specific signals
Probable causes	Contamination of PCR chemicals or primer mixes or sample carryover during the hybridisation protocol
Solution	Use fresh PCR chemicals, make sure to follow the general guide lines for PCR set up to avoid laboratory born contaminations. If necessary, use new primer mixes. Make sure that no well to well cross talk occurs during the hybridisation steps.
Observation	Only single spots are seen for some features (not duplicates)
Probable causes	The missing spots have been mechanically removed during the pipetting or drying steps or by touching the surface of the reaction chambers.
Solution	Use solely the Chip-Spin FVL2400N centrifuge for drying steps. Higher g-forces can result in spot dislocations. Avoid any contact with the reaction zone surface of the array fields

If none of the above given solutions solves your problem, please get in contact with your local distributor or contact our support team at support@chipron.com .

14 Literature

- [1] Iwobi A., Huber I., Hauner G., Miller A. and Busch U., (2011):
Biochip Technology for the Detection of Animal Species in Meat Products;
Food Analytical Methods, *Vol. 4, No. 3, pp. 389-398*

15 Order information

Product:	Quantity	Order No.
SPECIES 2.0	32 Tests / 04 Chips	A-750-04
SPECIES 2.0	96 Tests / 12 Chips	A-750-12

Related Products:

Product	Description	Order No.
Slide Reader Software	1 License	HS-200-01
Scanner PF 3650	Slide Scanner w. 10µm resolution	HS-300-01
CHIP Spin FVL 2400N	Mini-Centrifuge for LCD-Arrays	HS-500-01

16 Contact manufacturer

Chipron GmbH	Tel.: + 49 (0) 30 787994 70
Eresburgstrasse 22-23	Fax: + 49 (0) 30 787994 99
D-12103 Berlin	Email: info@chipron.com
Germany	Support: support@chipron.com

17 Symbol key

The following symbols are used on labels of the kit box and components therein



Manufacturer



Number of tests



Lot number



Date of expiry [month – year]



Order number



Storage temperature [range °C]