

SybGREEN qPCR Primer Pairs and Panels

Tools for RT-PCR

Application Guide

Table of Contents

Package Contents and Storage Conditions	2
Introduction	3
SybGREEN qPCR Primer Pairs and Panels	3
Gene Specific qPCR Template Standards	3
Figure 1: Preparation of qPCR Standard	4
Figure 2: Preparation of qPCR Standard (dilution)	4
Experimental Procedures	5
qPCR protocol: Preparation of the qPCR template standard	5
Preparation of PCR mix:	5
Figure 3: Suggested Plate Layout	6
qPCR setup	6
Figure 4: Sample data (SYBR Green I Real -Time PCR on ABI9700HT)	7
Arraying 200 reactions (in matrix box) into PCR plates	7
Trouble-shooting Guide	7
Frequently Asked Questions	8

Package Contents and Storage Conditions

For SybGREEN primer only:

- 1 Vial SybGREEN qPCR primer mix: 2 nmol lyophilized (sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 μ l dH₂O to make a final concentration of 10 nM.
- Application Guide

For SybGREEN primer and qPCR standard kit:

- 1 Vial SybGREEN qPCR primer mix: 2 nmol lyophilized (sufficient for 200 reactions)
- 1 Vial Gene Specific qPCR template standard: 50 X 10⁷ Copies** (double-stranded DNA) lyophilized. Add 50 μ l of d H₂O to the tube and vortex briefly. The resuspended template is 10,000,000 copy/ μ l.
- 1 Vial containing 1 mL 10 x PCR template standard dilution buffer.
- Application guide

For SybGREEN primer panels (200 reactions):

- 1 matrix box containing 96 gene-specific qPCR primer pairs: 2nmol lyophilized (sufficient for 200 reactions)
- Flash Drive containing plate layout, gene symbol, locus ID, and primer sequences
- Application Guide

For SybGREEN primer panels (3 pre-arrayed panels):

- 3 PCR plates containing 96 pre-arrayed qPCR primers mix in each plate lyophilized
- Flash Drive containing plate layout, gene symbol, locus ID, and primer sequences
- Application Guide

For SybGREEN primer panels (10 pre-arrayed panels):

- 10 PCR plates containing 96 pre-arrayed qPCR primers mix in each plate lyophilized
- Flash Drive containing plate layout, gene symbol, locus ID, and primer sequences
- 1 Application Guide

For SybGREEN primer panels (20 pre-arrayed panels):

- 20 PCR plates containing 96 pre-arrayed qPCR primers mix in each plate lyophilized
- Flash Drive containing plate layout, gene symbol, locus ID, and primer sequences
- Application Guide

The above components are shipped at room temperature, but should be kept at -20°C for long-term storage. If properly stored, they have a 12-month shelf life.

*The polymerase chain reaction (PCR) is protected by patents held by Hoffmann-La Roche. Purchase of any of OriGene's PCR-related products does not convey a license to use the PCR process, covered by these patents. Purchasers of these products must obtain a license before performing PCR.

** see page 6 on how to calculate the copy number of an unknown sample.

NOTE: FOR RESEARCH PURPOSES ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USAGE.

Introduction

Real-time PCR is one of the most important approaches in detecting the transcript level of a gene in a particular cell or tissue sample. Two types of real time PCR protocols are currently used. One is probe-based and the other is SYBR Green I based. The two protocols utilize fluorescent probes or dyes that are proportional to the amount of the PCR products generated at the end of each cycle. Differences between the two methods involve the types of fluorescent signals used and the ways in which the signals are generated.

In the TaqMan protocol, a pair of gene-specific amplification primers and a sequence-specific, fluorogenic probe is present in the PCR mixture. During amplification, the probes anneal proportionally to the single strand amplicons and are subsequently removed base-by-base by the 5' exo-nuclease activity of Taq polymerase. Consequently, the released fluorescence moiety generates fluorescent signals that are proportional to the amount of PCR product generated. This method is very sensitive and reliable.

The SYBR Green I protocol, on the other hand, does not require a sequence-specific fluorescent probe. SYBR Green I binds to double stranded DNA to generate detectable fluorescence, and the amount of signal is proportional to the amount of double-stranded DNA present. Since SYBR Green I binds indiscriminately to double-stranded DNAs, it will generate false signals if non-specific DNA contamination exists. Nevertheless, it is frequently used in many laboratories due to its simple protocol and low cost.

SybGREEN qPCR Primer Pairs and Panels


The key components of a successful qPCR are the primers used to amplify the template. Only well-designed primers can give trustworthy expression data using qPCR. Poor primer design usually results in non-specific PCR products and/or primer dimers. These undesired DNA byproducts would interact equally well with SYBR Green reagents, and therefore generate a false positive signal. To design primers with a high specificity and a low incidence of primer dimer formation, OriGene conducted a large scale study of qPCR involving thousands of primers. Based on the vast PCR data acquired through this study, OriGene scientists were able to develop a superb primer design algorithm. With this algorithm, OriGene is able to provide customers with genome-wide pre-designed, ready-to-use qPCR primers. Besides high specificity and high efficiency, each primer pair is also selected to span exon junctions whenever possible, thus avoiding the amplification of genomic DNA. The same effective primer design algorithm is utilized to design oncology-focused qPCR primer panels covering a variety of pathways for cancer biomarker profiling.

Each primer pair is supplied as 2 nmol lyophilized powder. Before use, reconstitute the primer mix in 200µl dH₂O to make a final concentration of 10 µM. The primer panels are supplied either as 10pmol lyophilized pre-arrayed primers mix in PCR plates or 2nmol lyophilized pairs as 200 reactions in matrix box.

Gene Specific qPCR Template Standards


A known qPCR template standard serves two main purposes. It functions as a positive control and as a reference for measuring the exact copy number of a transcript in an unknown sample. OriGene has a large collection of full-length human and mouse cDNA constructs which are ideal PCR templates. To increase the accuracy of DNA quantification, a full-length cDNA is purified using high purity ion-exchange columns, linearized by a restriction enzyme and quantified by the Picogreen method. The Picogreen method relies upon an ultra-sensitive fluorescent nucleic acid stain and has a much higher accuracy than a UV spectrophotometer for determining DNA concentration. The process for the preparation of a qPCR standard is illustrated in Figure 1.

Figure 1: Preparation of qPCR Standard



A. An EcoRI linearized qPCR template standard is run on a 1% agarose gel. **B.** Plot of a template DNA and PicoGreen standard from a PicoGreen assay. The concentration of the template DNA is converted from ng per μl to copy number per μl using the following formula: $(C * 10^{-9} / \text{MW}) * N_A$ (C: template concentration ng/ μl , MW: template molecular weight in Daltons, N_A : Avogadro's constant, 6.022×10^{23})

Figure 2: Preparation of qPCR Standard (dilution)



Experimental Procedures

qPCR protocol: Preparation of the qPCR template standard

1. The qPCR template standard is provided as a dried pellet. Add 50 μl of d H_2O to the tube and vortex briefly. The resuspended template is 10,000,000 copy/ μl . Thaw the 10 x dilution buffer tube from the kit at room temperature.
2. If the SYBR Green I method is used, reconstitute SybGREEN primer by adding 200 μl d H_2O to the primer tube and vortex it briefly. The resuspended primers are 10 pmol/ μl each. (Omit this step if Taqman probe method is used)
3. Prepare 1.0 mL of 1 x dilution buffer by adding 100 μl of the 10 X buffer to 900 μl H_2O .
4. Label six (6) Eppendorf tubes from 1 to 6. Transfer 90 μl of 1x buffer to each pre-labeled eppendorf tube (1 to 6).
5. Transfer 10 μl of template standard from the stock tube to eppendorf tube #1 (Fig. 2). Mix the solution by briefly vortexing.
6. Transfer 10 μl solution from tube #1 to tube #2, then mix the solution by vortexing..
7. Perform remaining dilutions by repeating steps 5 and 6.
8. Transfer 5 μl of diluted template solution from each tube to a 96-well PCR plate (use 2 μl when a 384-well plate is used).

Preparation of PCR mix:

The qPCR template standard can be used with either the Taqman probe based method or the SYBR Green I method. Step one, calculate the total volume of the PCR mixture based on the number of the total samples. The following recipes are for setting up a cocktail for 60 samples (6 template standard dilutions in duplicate and 48 unknown samples). The qPCR reaction volume is 25 μl volume in 96-well format for both the SYBR Green method and the Taqman probe method: When the 384-well format is adopted, reduce all reagents proportionally. The recommended reaction volume for 384 is 10 μl .

SYBR Green I method

Prepare the cocktail as follows:

	64X solution
2 x Master mix	800 μl
Primer mix	64 μl
H_2O	416 μl
<hr/>	
	1,300 μl

Taqman probe method


Probe method

2 x Master mix	800 μl
20 x Taqman Probe	80 μl

* with the gene specific primers included

H_2O	420 μl
<hr/>	
	1,300 μl

Figure 3: Suggested Plate Layout




qPCR setup

1. Transfer 5 μ l of each of the six qPCR Template Standard dilutions to Row A in duplicate
2. Add 5 μ l of each unknown sample to row C, D, E and F. The amount of cDNA in each reaction is ideally around 1 to 10 ng, diluted with the dilution buffer included.
3. Add 20 μ l of PCR cocktail to each of well of row A, C, D, E and F.
4. Seal the PCR plate with a real time PCR sealer.
5. Mount the plate to a real time PCR machine, such as ABI 7900HT.
6. Set up the PCR run according to the machine's user manual. Be sure to input the copy number when setting up the PCR standards (using the standard curve method).
7. Start the PCR reaction.
8. Analyze the data according to the machine's user manual.

Calculation of the copy number of your gene of interest

1. If the cDNA templates in your samples are single-stranded such as cDNA from RT reactions, the actual copy number of your gene of interest is two times the number you got by comparing to the qPCR standards. For example, if the copy number of your gene is 5 copy/ μ l from the standard curve of a qPCR program, the actual number is 10 copy/ μ l. The reason is that the first cycle for a single-stranded sample is to make the complementary strand; therefore, there is one cycle delay in the PCR reaction compared to the double-stranded cDNA template standard's reaction.
2. If the templates in your samples are double-stranded, the copy number of your gene of interest is the same as that calculated according to the qPCR standards.

Figure 4: Sample data (SYBR Green I Real -Time PCR on ABI9700HT)
 qPCR Template Standard and primer pair for NM_004808, Human N-myristoyltransferase 2 (NMT2)



Real time PCR data were obtained by using serial dilutions of a PCR standard for NM_004808 as the template with specific primers. The PCR was performed on the ABI 7900HT utilizing the SYBR method. The plot of ΔCt vs. cycle number for qPCR Template Standard is shown on the left (with copy numbers from 10,000,000 to 10 per μl) and the concentration dilution plot is shown on the right. An unknown sample (breast cancer tissue cDNA) is detected to contain approximately 37 copies/ μl of NMT2 based on the standard curve results.

Arraying 200 reactions (in matrix box) into PCR plates

1. Re-suspend the lyophilized primer mix in matrix tubes with 200 μl dH₂O to make a final concentration of 10 μM .
2. Prepare a master mix of PCR reagents (SYBR Green I method) using the proportions in the table below.

Item	Vol/well
dH ₂ O	5 μl
SybGreen	12.5 μl
cDNA template	5 μl

3. Add the PCR mix prepared above into each well of PCR plate.
4. Add 1 μl of re-suspended primer pairs from matrix box into the PCR plate using a multi-channel. While adding the primer pairs, make certain that the plate layout of matrix box corresponds to the qPCR plate.
5. Set up the PCR run according to the machine's user manual. Be sure to input the copy number when setting up the PCR standards (using the standard curve method).
6. Start the PCR reaction.
7. Analyze the data according to the machine's user manual

Trouble-shooting Guide

No qPCR product detected

1. A qPCR component may have been omitted in the reaction. Be sure to use a checklist when assembling the reaction mix.
2. Make sure the master mix is working by using a positive template control.
3. The sample used may have very limited qPCR template. Add more sample or prepare efficient templates by using a better reverse transcriptase.

The standard curve generated is not linear

The dilution of the Template Standard may have been done incorrectly. Use an accurate pipette and mix each dilution sample thoroughly.

A Ct value is out of the range of the standard curve

If the Ct reading for an unknown sample is too small (high abundance of a template), dilute the sample further. If the Ct is too high (low abundance of a template), then dilute the template standard further or use more unknown sample.

Frequently Asked Questions

Q: What are SybGREEN qPCR Primer Pairs and panels?

A: SybGREEN qPCR Primer Pairs are pre-designed, qPCR tested and ready-to-use primer sets for SYBR Green based qPCR experiments. The primer pairs are also available in a 96-well panel format covering broad range of pathways for cancer biomarker profiling.

Q: Can SybGREEN qPCR Primer Pairs be used in probe-based qPCR?

A: No. This kit has not been designed to accommodate any commercial probe based qPCR.

Q: What is the T_m of the primer and what is a typical size of the amplicon?

A: The T_m of a SybGREEN qPCR Primer is around 60, and the amplicon is around 95 to 140 bp.

Q: Does a SybGREEN qPCR Primer Pair amplify an exon junction in a cDNA target?

A: Whenever possible, an amplicon by a SybGREEN qPCR Primer Pair covers exon junction or junctions.

Q: My qPCR machine is not on the list of compatible machines on the website. Do you provide qPCR primer panels in plates compatible with other machines?

A: Yes, we can provide our primer panels in other plates within approximately one week of receiving the custom order. Please email the catalog numbers of the qPCR primer panels needed plus the catalog number and manufacturer of the plates to techsupport@origene.com. Our technical support scientists will prepare a custom quote that can be used to order custom panels.

Q: If I provide a list of genes I want to assay, can OriGene prepare a custom qPCR primer panel?

A: Yes, OriGene can design a custom qPCR primer panel. Custom panels are offered for human and mouse genes in 200 reactions format for purchase, delivered in 2-D bar code matrix tubes.

Q: Where can I find the sequence of the primers in each panel?

A: Each panel is delivered with a USB drive that contains a spreadsheet of primer sequences and locations.

Q: How much primer pair is in each well?

A: Each well in 200 reactions format contains 2 nmol of lyophilized primer pair. Please re-dissolve in 200ul of water so the final concentration is 10uM. Each well in the PCR plate contains 10pmol of primer pairs. The primers will be re-suspended once the PCR reagents are added to the plate.

Q: What is a qPCR template standard?

A: A gene specific qPCR template standard is a tube of linear DNA made from a full-length cDNA plasmid. The copy number of a template standard solution is determined by the PicoGreen method and calculation based on MW.

Q: Can a template standard be used in a probe-based qPCR?

A: Yes, as long as a corresponding probe is used.

Q: In your protocol, it is recommended to make six serial dilutions with one log of difference, can I make more dilutions and with different dilution scheme?

A: Yes, as long as it is diluted in the qPCR detection range and each dilution is mixed thoroughly.

Q: Can I use my own buffer to dilute the template standard instead of the buffer in the kit?

A: Yes, as long as the buffer has no negative effects on qPCR.

Q: Does OriGene offer custom-made template standards?

A: Yes. Please contact our Tech Support at techsupport@origene.com

Q: What is the OriGene guarantee on qPCR primers and template standards?

A: OriGene qPCR primers and templates are warranted for SYBR Green qPCR experiments. If your experimental results are not satisfactory, our scientists will work with you to pinpoint the problem and, if it is determined that our products are at fault, OriGene will refund your money in the form of a credit.

Q: How should I calculate the copy number of my gene of interest?

A:

- If the cDNA templates in your samples are single-stranded such as cDNA from RT reactions, the actual copy number of your gene of interest is two times the number you got by comparing to the qPCR standards. For example, if the copy number of your gene is 5 copy/ μ l from the standard curve of a qPCR program, the actual number is 10 copy/ μ l. The reason is that the first cycle for a single-stranded sample is to make the complementary strand; therefore, there is one cycle delay in the PCR reaction compared to the double-stranded cDNA template standard's reaction.
- If the templates in your samples are double-stranded, the copy number of your gene of interest is the same as that calculated according to the qPCR standards.

Q: I am writing a paper for publication and need to describe this product. How should I cite?

A: We recommend that you refer to the product by its specific catalog number and refer to us as OriGene Technologies (Rockville, MD). Furthermore, we'd love to hear from you when your paper is published. Inform us and we will send a gift.