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
and service paper of your Siro II

Siro II
Typ 801

We are a company-member of rehaKIND e.V.

International Support Group for the
Rehabilitation of Children and Young People

Editor: Bernhard Wendel
Technical collaborator: Bernd Rothardt
As at March 13
Technical alterations and misprints excepted.

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© All texts and pictures are subject to copyright and may not be published - not even as an excerpt! - without the explicit authorization of SORG Rollstuhltechnik GmbH + Co.KG

© translated by: www.my-translator.biz

1. Contents

1. Contents	3
2. Preface	7
3. General advice	8
3.1 Functional purpose	8
3.2 Mark of SiroII	8
4. Safety precautions	9
5. Handling	13
5.1. Completion of use condition	13
5.2. Unfolding the chair	13
5.3. Folding up the chair	14
5.4. Loading and transport	15
5.5. Getting in and out the chair	15
5.6. Transporting the wheelchair	16
5.7. Rolling off	16
5.8. Carrying	16
6. Brakes	17
6.1. Standard knee lever brakes	17
6.2. Knee lever brakes - braked lockable	17
6.3. Knee lever brakes - braked and unbraked lockable	17
6.4. Brake lever lengthening	17
Safety precautions and advice for your orthopaedic mechanic master	18
6.5. Drum brakes	19
Advice for your orthopaedic mechanic master	19
7. Adjustment of wheel postions	21
7.1. Punch plate and camber	21
7.2. Alteration of the centre of gravity's position	22
7.3. Adjustment of seat-tilting	22
7.4. Critical height of obstacle	22
7.5. Static stability	22
Safety precautions and advice for your orthopaedic mechanic master	22
7.6. Adjustment of front wheel	23
Advice for your orthopaedic mechanic master	23
7.7. Supporting wheel (overturning protection)	24
Advice for your orthopaedic mechanic master	24

8. Wheels	25
8.1. Driving wheels	25
8.2. Grip rings	25
8.3. Spoke cover	26
8.4. Wheel cover	26
9. Seatback	27
9.1. Backrest	27
9.2. Seatback cover cloth standard	27
9.3. Seatback cover cloth restretchable	27
9.4. Fixed seatback	27
9.5. Lateral truss pads	28
9.6. Seatback inclination with notch system	28
9.7. Headrests	28
9.8. Stabilizing rod	28
10. Sliding handles	29
10.1. Extendable sliding handles / One-handed sliding handle	29
10.2. Sliding handles with sliding pole	29
Safety precautions	29
11. Seat	31
11.1. Fixed seatFix	31
11.2. Seat extension	31
11.3. Bolting device	31
11.4. Seat cover cloth	32
12. Side parts	33
12.1. Side parts	33
12.2. Arm pads	33

13. Leg supports	36
Advice for your orthopaedic mechanic master	36
14. Further options	37
14.1. Therapy table	37
14.2. Safety belt	37
14.3. Tilt clip	37
14.4. Extension parts of other concerns	39
15. Cleaning and maintenance	41
16. Warranty statement	43
16.1. Quality assurance system	43
16.2. Design life	43
16.3. Working life	43
16.4. Complaints	44
16.5. Data plate	44
17. Warranty clauses	45
18. Compliance statement	47
19. Warranty of delivery	48
19.1. Assemblage and inspection	48
19.2. User-chronology	48
19.3. Handing over	48
Please make sure that the concerning items are filled out	
20. List of maintenance	49
21. List of repair works	51

2. Preface

Dear Customer,

today, your child has met its new “daily companion” for the first time. A companion that would never let your child down and that will accompany him/her through thick and thin. Day after day, month after month, year after year. The most important precondition for such a cheerful and safe companionship is a sensitive reaction against all individual needs.

We want to be such a reliable companion to your child. Therefore, we use our years of experience and we take a lot of care over development and assemblage. This is the only way for us to guarantee the permanent best quality, individuality and safety. After having been certified according to DIN EN ISO 9001:2000, we are under the obligation to continue this standard of quality perpetually. With a measurable result!

Special people need special solutions. This is why we are developing mobilization solutions, all being as unique as their user. Therefore it may happen that some descriptions slightly differ from your child's wheelchair because we're on the road of permanent development and this individual variant may be not yet included in this user manual. Our specialist staff or your intermediary medical service centre will be pleased to answer all your questions about the handling.

Please keep this service paper in a safe place and have filled out the required data on page 44. It will be proof of regular inspections, which might be decisive in case of a warranty claim. Moreover, it will be an important source of information for the cost-bearer in the event of the wheelchair's further use.

For a reliable use of the new wheelchair, please read the following sections carefully and pay attention to our advice and instructions for your child's safety! More questions? Please consult our proficient staff, they'll be pleased to help you!
Service phone +49 07254 / 9279-0

Hugo Sorg, Managing director

3. General advice

This user manual has been drawn up not only for you but also for your orthopaedic mechanic master of the medical service centre – specialists who are able to adjust the wheelchair exactly to your individual needs and who will complete the wheelchair's use-state.

Therefore, this user manual includes instructions which you may in no circumstances carry out yourself – improper action will always be a risk to you, the user and your environment. Please take some time to read this user manual carefully and keep it in a safe place.

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4. Safety precautions

When developing and manufacturing our wheelchairs, we set great store by safety. In order to maintain this high quality of safety in everyday life, it is in your own and your child's interest to carefully pay attention to the following rules. Improper use could cause a hazardous situation to you and your child and our assurance of warranty would no longer be valid.

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Familiarizing with Siro II

- **Take your time to familiarize yourself and your child with the new situation. Be very patient and use caution to make sure that your child will feel safe in its new wheelchair. Because Muby is an activity-wheelchair whose limits you'll have to get to know. You should really let your child undergo these limits quite literally.**

Make sure that during the following steps both of you are safeguarded by an experienced and strong assistant!

- Lock the parking brakes and the drum brake and put your child into Siro II.
- Carefully practise the wheelchair's behaviour on shifting the centre of gravity in all 4 directions! Make your child reach for objects, holding them further and further.
- Sensitize your child for these changed situations and ask if something changes, and what exactly. And do not forget to express clearly that you will now intensify the situation.
- In any case, this must all be practised without fear. Do not force the learning process!
- Only after your child has been familiarized with the wheelchair's static stability in its standstill position, you should risk trying out the wheeling behaviour.
- Children soon understand how to drive a wheelchair, it nearly is an automaticity. The process of braking down however, needs thorough practice because it is dangerous! Tell your child in clear terms that the knee lever brake is **NOT** designed for braking down the wheelchair's speed and that it may fall out of the chair.

- Your child must learn that it should only use the grip rings for braking down the wheelchair. Set Siro II in a **slow** motion and make the child brake down the wheelchair to the standstill position with its palm on the grip rings. In doing so, the child should keep in mind what will happen as it reaches between tyre and grip ring (as far as the wheelchair is not equipped with our Revolution-wheels). If your child has got delicate hands, please give it some gloves, eg. for cycling, meaning with a leather inside – in no case mitts nor woolen gloves!

(Never forget that strong friction will heat up the grip rings. Therefore, either wear some gloves or use the grip ring coating. Your medical assistant will be pleased to give you some advice.)

- As your child has thoroughly been familiarized with the items „centre of gravity“ and „brakes“, you may carefully begin with the next items.
- At first, repeat the above exercises in the open country on even territory and only after that on gradients and slopes and always with another assistant.
- It is most important to create a situation your child is not afraid of. But: do never “play” dangerous situations, eg. by tilting the chair and then catching it meaning it to be a joke... Only real-life dangers may be mastered by your child.

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- plug-in axles have been thoroughly tightened.
- Your child may **never** use the **foot rest(s) or leg supports** for getting in or out the wheelchair. The chair may turn over forwards.
 - Lock the parking brake each time your child gets in or out the wheelchair.
 - If possible, always use the supporting wheel. Especially when wheeling on gradients or slopes, on uneven territory or thresholds, the danger of turning over is really high.
 - When wheeling in public traffic, keep in mind that you must act according to Road Traffic Act. In order to avoid accidents, you should always control if the reflectors are tidy and both of you should wear bright colours to be well seen.
 - The wheelchair may not be carried together with its occupant but may only be wheeled.
 - When transporting the wheelchair in a vehicle, your child may only sit in the chair as the chair is fixed onto the required devices. (see handling – transporting the wheelchair p.16).
 - Have the wheelchair serviced regularly by a qualified medical service centre. It is necessary to examine if all screws are properly tightened, if welded joints show any damage and if signs of wear can be detected.
 - Your child should get over steps and stairs **only** with an assistant. Take care that the supporting wheel is swung in and that the assistant(s) always uses the fixedly screwed **metal parts** of the wheelchair for carrying, because the sliding handles' rubber caps may unexpectedly come off and cause serious accidents!
 - In order to avoid injuries, your child should not take its hand between the driving wheels and knee lever brake or wheel cover.
 - Do not wheel on gradients over 4°; meaning ca. 7%. The danger of overturning increases on gradients over 4°.
 - An occupied wheelchair should never be left unsupervised on a slope, not even with locked knee lever brake.
 - With an active (swung out) supporting wheel, the max. critical obstacle height is ca. 3 cm.
 - The max. useful vehicle load capacity is 100 kg. The standard frame is not equipped for more weight.

5. Handling

5.1 Completion of use condition

The wheelchair will be in working order for everyday life after having examined the following basic items:

- do the parking brakes show functional performance?
- do the drum brake show functional performance?
- are the plug-in axles thoroughly snapped in the adapters?
- are the steering wheels in working order and may they be easily rotated?
- do the driving wheels show adequate air pressure?
- is the supporting wheel in working order?
- is the fixed seat shell thoroughly locked into the securing devices at the back?
- is the seat shell in perfect condition?
- **is the seat shell fest im Sitzschalen-Tragerahmen arretiert?**
- is the seat-filling unobjectionable working?
- is the gas pressure spring in working order?
- **Sitzt der Bowdenzug am Auslösehebel der Kantelung fest in seiner Aufhängung?**

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5.6 Transporting the wheelchair

First of all we have to point out that transporting a person sitting in a wheelchair signifies a considerable security risk!

In principle, a wheelchair may only be transported unoccupied!

Therefore, transporting and securing a wheelchair in a vehicle, with or without occupant, is always at your own risk. SORG Rollstuhltechnik GmbH + Co.KG will not be liable for injured people or damaged materials, resulting from such a transport. At its immovable points the unoccupied wheelchair must be secured in the vehicle with suitable securing devices, so it will not put passengers in danger in case of a sudden brake manoeuvre!

In order to transport a wheelchair as instructed, it must absolutely be equipped with an **AMF-system** or a **similar retention system with MOT-certificate (Technical Authority Inspection)**. In addition to this, the vehicle must be equipped with a **headrest, firmly connected with the vehicle**, as well as with an **integrated person retention system on the basis of three point seatbelt** if the wheelchair is to be transported together with its occupant.

More questions? Do not hesitate to contact us, just give us a call: service phone +49 7254-92790. We will be pleased to give you advice and to refit your Siro II if necessary.

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5.8 Carrying

Your child should, together with its wheelchair, be carried by an assistant **only** in extremely seldom exceptions. To the most part it will be easier rolling off the wheelchair over the obstacle with its driving wheels! And if it really should be unavoidable, only frame and back tubes are suited for carrying.

Safety precautions:

The extendable sliding handles, as well, may in no circumstances be used for carrying the wheelchair. When wheeling with your chair, never use escalators! Not even with an assistant and especially not being carried! Department stores and public institutions with escalators always have suitable elevators.

6. Brakes

6.1 The knee lever brake is a parking brake for blocking the wheels and may in no circumstances be used for braking down the wheelchair's speed!

Each wheelchair is equipped with two knee lever brakes. These only have the purpose of blocking the wheels in their standstill position. They are not designed for braking down the wheeling chair. For this, your child should always use the grip rings and it is also possible for you to deal with this, using the drum brake. Please, keep in mind that strong friction will considerably heat up the aluminium grip rings. For your child's and your own safety, always take care that the brakes are working perfectly. Low air pressure, rainy and wet days, bad treads and a distance too large between brakes and tyres, all may have an effect on the functional performance of the brakes.

6.2 Knee lever brake - braked lockable (as option)

The same function as above, yet, in addition lockable by a safety bolt. For securing the knee lever brake in braked-state, lock the brake. The bolt will be interlocked as far as it is not pulled out and turned by 90°. For releasing the brake again, pull out the safety bolt and at the same time open the brake.

6.3 Knee lever brake – braked and unbraked lockable (as option)

The same function as above, yet, in addition lockable by a red safety bolt in every position. This bolt locks automatically in braked- or unbraked-state. To actuate the brake, you only have to unburden the brake lever and pull the bolt outwards to get the brake in the needed position. After that the bolt again locks automatically and so it prevents a change of the brake position.

6.4 Brake lever lengthening (as option)

In case the brake is correctly adjusted but the operation lever may no longer be reliably operated by manual power, in no circumstances you may change the distance between brake dowel and tyre! In this case, your child will need our brake lever lengthening and your medical centre will be pleased to equip the wheelchair with it. Please take care that your child is not sitting on the lengthening – danger of breakage


Safety precautions:

- Take care that the tyres always show correct air pressure. Low tyre pressure is not only a needless obstruction for easy wheeling, it above all reduces the parking brakes' functional performance considerably and their operation function may even be eliminated completely.
- Always take care that the brakes have to operate reliably. Especially when wheeling on slopes, it must be ensured at any time that the wheelchair will stand still when the brakes are locked.
- Keep in mind that on rainy and wet days and as the chair is dirty, the braking operation is considerably changed and reduced.
- After each adjustment of the wheel positions as well as after each tyre-change, the brakes must be adjusted anew. Have the knee lever brakes adjusted by a qualified medical service centre.

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For adjusting the brake as the tyres are properly inflated, both filister head screws (A+B), fixing the brake on the brake support, must be loosened (max. air pressure on the outer tyre). The brake must be adjusted in such a way that with a released brake on 1"-tyres a distance (s) of ca. 15-18mm between brake dowel and tyre surface exists and a distance of ca. 12-15mm on 1 3/8"-tyres. Re-fix the brake in this position with the filister head screws.

Now do a test on it: the brakes should block the wheel thus far that the occupied chair can no longer be wheeled. The wheelchair will instead skid across the floor with blocked wheels.

6.5 Drum brakes


In contrast to the parking brake, the drum brake may be used for **braking down** your speed. In addition, our drum brakes have a locking-mechanism with whom they can be braked in their maximum braking-state. (as an option, a childproof safety catch is available as well)

- Please keep in mind that the drum brake on its own is **not suited** for a reliable **securing** of the wheelchair, eg. on slopes. The brake pull rod may unexpectedly be torn, resulting in bad accidents.
- Hence, for your child's safety, take special care that when locking the drum brake, the **parking brake** (knee lever brake) must always be activated **as well!**
- Drum brakes have delicate adjustment alternatives and they require some service and maintenance, however, they are not dependent on air pressure.
- In the second lock-in position of the brake lever, the brake should be so firmly fixed that the occupied wheelchair may no longer be wheeled (the wheelchair will instead skid across the floor with blocked wheels).
- Please have the drum brakes regularly serviced by your intermediary medical service centre.

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The adjustment of the drum brake is effected by a set screw on the lower end of the brake pull rod. The brake is adjusted by turning the set screw counter-clockwise.

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9.2 Alteration of the centre of gravity's position

Fundamentally, Siro's centre of gravity is situated in front of the wheelchair's rear axles. The distance between centre of gravity and rear axles is decisive on the handling and the safety of the chair. If the axles are right behind the centre of gravity, the child will wheel and manoeuvre much easier, yet the danger of overturning increases, **above all if the chair is tilted!**

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9.3 Critical height of obstacle

Siro's critical height of obstacle, with an activated supporting wheel is at the best ca. 2.5 cm.

9.4 Static stability

Siro's static stability is 4°, i.e. ca. 7%. On steep slopes, the wheelchair together with its occupant may tilt over and your child may be injured badly!

Safety precautions:

Fundamentally, your child may turn over backwards in its wheelchair as it is getting over an obstacle. This depends on the wheelchair's centre of gravity and on the height of obstacle. Please always make sure that an assistant is helping and safeguarding the child as it wants to get over high obstacles, e.g. stairs and thresholds.

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On trying out different positions of the rear wheels, use the supporting wheels on your first trials and test the behaviour of the centre of gravity at best together with a strong helper who will safeguard your tests. **Beachten Sie bitte auch die veränderte Situation bei Kantelung.**

9.5 Adjustment of front wheel


The following alterations should only be carried out by a qualified medical centre, in order to ensure the wheelchair's functional performance, its therapeutical usefulness and especially your safety.

On changing the extension of the steering wheel onto the frame, the seatangle may be changed as well. Take care on all alterations that exactly the same adjustments are carried out on both sides, that the centre of gravity is not shifted too far backwards (danger of overturning) and that the rotation performance is not obstructed by the leg support.

Please remember to adjust the correct position of the front wheel adapter after each modification of the wheel size, of the wheel base etc., otherwise the perfect wheeling function of the front wheels will no longer be ensured.

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For adjusting the steering wheel adapter, loosen both screws (A) which are fixing the front wheel adapter onto the frame. Now the eccentric disc (B) can be turned until the front wheel adapter is in an exact **perpendicular position** - with a protractor the angle may be rectified. **Do not forget to re-tighten all screws after having finished the adjustment!**

9.6 Supporting wheel (overturning protection)

The purpose of the supporting wheel is to avoid turning over backwards. As depicted in the illustration, the wheel should be adjusted with a distance of 2-5cm towards the floor.


For getting over a threshold, the supporting wheel must be swung in by the assistant and afterwards be swung out again; otherwise the wheelchair may not be wheeled properly.

The swing jib of the supporting wheel is pivoted in the fixed fitment. In order to swing out the supporting wheel, push it down with your foot or hand and turn it by 180°. Then release it, so it may be engaged into the position slot.

If you push down the jib in the fitment, it will come out of its position and may be swung back. Make sure that the jib is snapped in its end position because a loose wheel may result in an ankle injury of the assistant.

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For adjusting the supporting wheel **higher**, loosen screw A and pull out the jib from the spring and the fitment. Re-set screw B in the jib by the needed punch hole distance. Then assemble spring and jib into the fitment and re-tighten all screws properly.

10. Wheels

10.1 Driving wheels

Fundamentally, the driving wheels are nearly the most delicate parts of the wheelchair because many factors concerning your safety and the wheeling performance depend on them!

In this connection, once again we have to point out politely, yet insistent, that the safety precautions are most important!

- In case the wheelchair is equipped with pneumatic tyres, please regularly examine the air pressure (at least every 2 weeks) because the functional performance of the brakes depends on them! You'll find details concerning the required air pressure on the outer tyre.
- Examine the tyres' tread regularly because a bad tread will not only affect the functional performance of the knee lever brake, it may even be completely eliminated in unfavourable circumstances.
- This is also valid on puncture resistant tyres!
- Examine regularly (at least every 3 months) if the spokes are tightened properly. If necessary, have the spokes re-tightened by your medical service centre, otherwise the wheel will be a danger.
- The grip ring's screws may loosen themselves as they are strained too much. Then the grip rings as well will be a considerable danger because they may come off and the wheelchair may not be braked down. Hence, also examine if the grip rings are tightened properly at least every 3 months and if necessary, have them re-adjusted by your orthopaedic mechanic master.
- Please also examine if the spoke covers, if equipped with, are firmly fixed! Loose ends may result in injuries of your child's hand. If you need the required clips or cable fixers, ask your medical service centre or just give us a call, we'll send them right away; service phone +49 7254-92790.

10.2 Grip rings

Strong friction will easily heat up aluminium grip rings. Some gloves are surely a good aid yet your child's wheelchair may also be equipped with a grip ring coating, even if the chair is equipped with rehavolution wheels with integrated grip rings.

On standard wheels, the grip ring extension and its distance towards the wheel rim is a variable feature. For this, the tyre must be removed from the wheel and the grip ring screws on the rim must be loosened. Then re-assemble the screws in the needed position and put the tyre back onto the rim. Please take special care that the rubber band covers the screw heads otherwise the air tube may be damaged.

10.3 Spoke cover (as option)

Depending on the design, the spoke cover is attached to the spokes by clips or cable fixers. For the most part, these parts may not be used twice if the spoke cover has been removed.

10.4 Wheel cover (as option)

The most important purpose of the wheel cover is to protect your child's hands from getting injured, but it will also keep off the dirt caused by the rotating driving wheels. At our works, we have adjusted the wheel cover in the optimum position.

So, if you want to modify the wheel position, you will also have to re-adjust the wheel cover. For this, the screws on the wheel cover holder must be loosened, then get the wheel cover into the needed position by means of the slotted holes. Re-tighten the screws - and off it goes!

In case the wheel cover holder is bent out of shape, please ask your medical service centre, so that the aluminium parts will not be broken off or the paint layer will not peel off as the part is bent improperly.

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

12. Sliding handles

The sliding handles of all our wheelchairs have black handle caps at their end made out of PVC. Those are connected firmly with the sliding rod. In unfavourable circumstances however, the caps may loosen themselves from the tube – eg. in cases of heat, humidity or extreme strain.

With insistence we have to warn you:

Never use the sliding handle caps to roll off an occupied wheelchair over an obstacle, or even worse - on stairs! For this, always use parts which are firmly fixed onto the frame and make sure that another person is safeguarding the action.

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13. Leg supports


15.1

To your child's seating convenience, it is most important that its feet are resting on the leg supports in a relaxed manner. On leg supports that are not properly adjusted, your child's knees may fall outwards (due to insufficient support underneath the thighs). And unpleasant congestions may occur in the thighs as the leg supports are too far below.

Do not take an unnecessary risk and have the following adjustments carried out only by your experienced orthopaedist.

Advice for your orthopaedic mechanic master:

Siro's leg supports may be adjusted in height in several steps. Hence, the wheelchair may be adjusted exactly to the needed length not only of the thigh but also of the lower leg.

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Please also keep in mind that the foot rest/s may not bump into the front wheels! When changing of direction, the front wheels must always rotate and swing free underneath, i.e. behind the foot rest, i.e. behind the leg support fitment, otherwise you may hurt somebody or even have an accident!

Finally, do not forget to re-tighten all screws and keep in mind that self-locking nuts with a nylon-inset may be used only a few times and must then be replaced or else they will no longer serve their purpose in a reliable way!

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16 Further options

16.1 Headrests

Our headrests may be adjusted in height, in depth and in angle. For this, loosen the clamping levers of the connection you want to adjust by turning them slightly. Then re-tighten them after you have reached the needed position. The headrest's angle may be adjusted by slightly loosening the screwed connection and tightening them as you have reached the needed position.

16.2 Therapy tables

Our therapy tables are adjustable according to your child's individual needs. For this, loosen the locking screws of the table fitment at the back tube by slightly turning them and then adjust the table to the needed height or depth. Furthermore, you may modify the table-angle by adjusting the locking screw at the table's underside (depending on design).

16.3 Safety belt

The safety belt protects your child from falling out of the chair. It is fixed onto the frame, i.e. onto the seat pan. Always take care that the safety belt is fastened in a firm and secure manner.

16.4 Tilt clip

The tilt clip is an aid for tilting the wheelchair much easier as the child wants to get over small obstacles. Your assistant shall carefully operate the tilt clip using his/her foot and at the same time keep hold of the wheelchair at its sliding clamp.

16.5 Stick holder

The stick holder must only be used for transporting walking aids. Excessive use may turn over the wheelchair backwards.

16.6 Extension parts of other concerns

For using extension parts of other concerns, such as:

**AMF-system, Scalamobil, s-max,
viamobil, E-Fix E20 or e motion M12,**

it is necessary to carry out some modifications on the wheelchair frame. If you are interested, please consult your intermediary medical service centre or us, we'll be pleased to give you advice! Service-phone +49 7254-92790.

17. Cleaning and maintenance

Please pay attention to the following advice for cleaning and maintenance so that the new wheelchair will remain the reliable everyday therapeutical partner for a long time.

- Examine the tyres' air pressure at least every 2 weeks or according to demand at regular intervals, as well as the connected function of the brakes. It is absolutely necessary for your safety. The spokes of the driving wheels should be tightened once a year. For this, please consult your specialist of the medical service centre.
- Every quarter year you will have to control if all relevant screwed connections are still tight and if necessary, you will have to tighten them. In doing so, keep in mind that a safety nut can be used only once without reduction of its safety efficiency. Your specialist of the medical service centre will answer further questions.
- For cleaning the wheelchair's frame, only use a gentle cleaner and a soft cloth.
- If the wheelchair has become wet, please dry it as soon as possible with a clean, soft cloth. Sand, sea water and salt for icy roads will lead to damage and corrosion. In order to avoid this, clean the wheelchair carefully and thoroughly after contact with one of these or similar substances.
- All ball bearings are maintenance-free. For preserving the running properties of the plug-in axle permanently, clean the plug-in axles and the fitment bushes of the plug-in axles every quarter year with a clean and soft cloth and than use some lubricating oil as grease.
- Gradually dirt particles will gather between the steering wheel and the fork of the steering wheel, which will render the wheel heavy to move. For cleaning, remove the steering wheel out of the fork and free it from dirt.
- In case of a puncture, the tyre may be changed in the same way as on the usual bicycle. If you're on your way, it is recommendable to take a puncture spray or an emergency repair set, or else to have the wheelchair equipped with puncture resistant tyres

18. Warranty statement

18.1 Quality assurance system

When manufacturing our wheelchairs, we pay attention to the greatest possible quality and comply with the relevant German Industrial Standards (DIN) and EC-directives. Therefore we have been certified according to DIN EN ISO 9001:2000 in 2003 and we take the stipulated subsequent audits regularly. All our wheelchairs, mobile standing chairs and seat pan undercarriages meet the requirements of law concerning medical devices and are correspondingly CE-labelled (special designs excepted).

18.2 Design life

According to the criteria of  , the International Support Group for the Rehabilitation of Children and Young People, concerning quality and safety on the further use of rehabilitation devices, we can as manufacturer of Muby reliably guarantee a design life of 7 years, provided the wheelchair is used in a normal way and is regularly serviced by a qualified specialist. Until over this period of time we can guarantee supplying you with the according spare parts for another 2 years (special designs excepted).

18.3 Working life

Especially devices in further use demand high standards. Therefore we freely recommend you, in case of further use of Muby, not to exceed a further working life of 3 years within the period of design life.

„Working life“ shall be understood by the period of time during which a device has been permanently in use. Hence, storage periods with the customer and/or the specialist store can not be called working life. “Design life“ however, means the period of time which we have defined on the basis of our experience and technical facts and which should not be exceeded in further use. Yet, we emphasize that this life span may directly be influenced as it is clearly dependent on service and maintenance. Under the most favourable conditions the wheelchair’s functionality may be extended over the defined design life in a reliable way.

18.4 Complaints

If you really should have cause for complaint, it goes without saying that we will solve your problems as soon as possible and replace or repair the defective parts free of charge. We give you this guarantee for a period of 5 years, however wear and tear parts as well as all sewed parts are unfortunately excluded (cover cloth for seat and back, cushions, truss pads, headrest pads).

In order to be served promptly, please enclose with the return the invoice-no., the order-no. and/or the consignment-no., as well as a short description what has happened. So we may be able to avoid the failure in the future. We are grateful for your help. You'll find the required details on the data plate (example below).

Damage in transit must immediately after receipt and examination of the goods be communicated to the freight company, otherwise it shall not be liable for the defective consignment, which is in accordance with the transportation's general conditions.

The five-year guarantee period starts with the date of delivery. Please read the warranty clauses for an exact definition of our guarantee conditions.

18.5 Data plate

The data plate can be found on the lower frame tube.

It contains the following information:

manufacturer, order-no. (Auftr.-Nr.)

design name, (Modell)

type-no. (Typ)

aid index no. / licensing no. (Zul.-Nr.)

max. nominal weight in kg (max. G.)

seat width (SB), seat depth (ST) and back height (RH) in cm.

In order to be served promptly, please indicate the **order-no.** in case of a spare part order or complaint. Both last figures refer to the year of construction.

 Rollstuhltechnik GmbH + Co. KG	
Auftr.-Nr.	1744 04
Modell	SIRO II
Typ	801
Zul.-Nr.	26.99.01.3 xxx
max. G.	100 kg
SB	48
ST	56
RH	50
	

19. Warranty clauses

I. Subject of Warranty

These warranty conditions are valid for wheelchairs from the limited company Sorg Rollstuhltechnik GmbH & Co. KG. Legal claims out of deficiencies warranty of the end-user with regard to the medical service centre or other service companies, responsible for the patient's medical care, will not be affected.

II. Extent of Warranty

Sorg Rollstuhltechnik GmbH & Co. KG grants a warranty period of three years after first use on frame parts and cross braces in accordance with the conditions from section III and to the exclusion of section IV. With the warranty conditions, the end-user is granted the right to claim substitute delivery and improvement(s) according to paragraph 2.

In case production and/or material faults occur two years after first use, Sorg Rollstuhltechnik GmbH & Co. KG will either replace the wheelchair partly or as a whole or will repair the defective frame parts and the cross braces. After this period, Sorg Rollstuhltechnik GmbH & Co. KG shall only be committed to replacing the defective parts. Further claims for reduction, contract cancellation or compensation will not be granted.

III. Claim Conditions

In order to put in a claim, the following conditions have to be met:

- the wheelchair must be in first use
- the wheelchair must be adapted by a specialist centre from the National Health Service (NHS)
- the transmission of the full care must be carried out by the experts from the NHS with enclosed list of faults
- wheelchair adaptations may not exceed the designed completion of Sorg Rollstuhltechnik GmbH & Co. KG, yet, in exceptional cases they may, as far as an agreement with Sorg Rollstuhltechnik GmbH & Co. KG has been made
- Sorg Rollstuhltechnik GmbH & Co. KG's products may only be repaired and extended by specifically qualified employees with the exclusive use of original spare parts.

IV. Exclusion Clause

Sorg Rollstuhltechnik GmbH & Co. KG has the right to reject a claim on the following conditions:

- as the claim conditions according to section III are not met
- as the reduction of the wheelchair's use feature is based upon unprofessional use, especially after improper alterations and modifications
- as the reduction of the wheelchair's use feature is based upon the usual process of consumption
- as the reduction of the wheelchair's use feature is based upon faulty or improper storage, transportation or use
- as the reduction of the wheelchair's use feature is based upon a change of the patient's physical constitution, eg. in case of a considerable increase in weight
- as the reduction of the wheelchair's use feature is based upon Force Majeure or as the wheelchair was wilfully damaged
- as Sorg Rollstuhltechnik GmbH & Co. KG determines that the maintenance regulations and instructions have not been complied with.

V. Liability

Sorg Rollstuhltechnik GmbH & Co. KG shall not be liable in case of negligence and accepts no responsibility for intentional action of accomplices. Furthermore, entitlements to compensation are restricted to damages, foreseeable at the beginning of the obligations' performance.

VI. Additional Provisions

- Warranty cases must be communicated within appropriate period.
- Replaced parts must be returned to Sorg Rollstuhltechnik GmbH & Co. KG.
- Warranty obligations shall be fulfilled in Oberhausen-Rheinhausen, Germany.
- These warranty conditions are subject to the civil and commercial laws of Federal Government of Germany.

VII. Warranty Period for Repairs

The warranty period for repairs will end together with the wheelchair's warranty, however, will be at least one year.

20. Compliance Statement

EC – Compliance Statement
in accordance with Appendix VII of directive
93/42/EEC concerning medical devices

We, Company Sorg Rollstuhltechnik
GmbH & Co.KG, Benzstraße 3,
D 68794 Oberhausen-Rheinhausen,
in sole responsibility as a manufacturer,
declare that our product

Siro II

meets the basic requirements
of the directive 93/42/EEC, Appendix I.

Oberhausen-Rheinhausen, March 13
(signature)

21. Warranty of delivery

This wheelchair has been assembled and inspected with great care!
Our staff stands surety for that with its signature.

Mechanic:

Date _____ signature _____

Inspector:

Date _____ signature _____

User-chronology

name of user or its legal representative	date of birth	street	code	Place
1.				
2.				
3.				

Wheelchair handed over by _____ on _____

_____ stamp medical service centre

_____ date and signature

22. List of maintenance

maintenance due:	parking brakes and drum brakes rota of 3 months					
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						

maintenance due:	all frame parts and welded joints rota of 6 months					
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						

22. List of maintenance

maintenance due:	parking brakes and drum brakes rota of 3 months					
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						

maintenance due:	all frame parts and welded joints rota of 6 months					
date, stamp medical service centre, signature						
date, stamp medical service centre, signature						