

VPAP[®] III & III ST

Reorder number: 24807/1 03 06

USER'S MANUAL

English (USA)

RESMED

VPAP[®] III & III ST

Reorder number: 24807/1 03 06

USER'S MANUAL

English (USA)

Manufactured by:

ResMed Ltd 97 Waterloo Road North Ryde NSW 2113 Australia

Tel: +61 (2) 9886 5000 or 1 800 658 189 (toll free) Fax: +61 (2) 9878 0120 Email: reception@resmed.com.au

Distributed by:

ResMed Corp 14040 Danielson Street Poway CA 92064-6857 USA

Tel: +1 (858) 746-2400 or 1-800-424-0737 (toll free) Fax: +1 (858) 746-2900 Email: reception@resmed.com

ResMed (UK) Limited 65 Milton Park Abingdon Oxfordshire OX14 4RX UK

Tel: +44 (1235) 862 997 Fax: +44 (1235) 831 336 Email: reception@resmed.co.uk

ResMed Asia Pacific Ltd 97 Waterloo Road North Ryde NSW 2113 Australia

Tel: +61 (2) 9886 5000 or 1 800 991 900 (toll free) Fax: +61 (2) 9889 1471 Email: reception@resmed.com.au

ResMed Finland Niittykatu 6 FIN 02200 ESPOO Suomi

Puh: +358 9 8676820 Faksi: +358 9 86768222 Sähköposti: reception@resmed.fi

ResMed GmbH & Co. KG Rudolfstraße 10 D-41068 Mönchengladbach Deutschland

Tel: +49 (0) 2161-3521-0 (Reception), +49 (0) 180 22 22 668 (Service-Telefon; 0.06 €/Anruf),

+49 (0) 180 22 66 888 (Wartungstelefon; 0.06 €/Anruf) Fax: +49 (0) 2161-3521-1499 Email: reception@resmed.de

ResMed Japan Nihonbashi Hisamatsu Bldg. 4F, 2-28-1 Nihonbashi-Hamacho, Chuo-Ku, Tokyo 103-0007, Japan

Tel: +81 (3) 3662 5056 Fax: +81 (3) 3662 5040

ResMed Malaysia Sdn Bhd Suite E-10-20, Plaza Mon't Kiara No. 2 Jalan 1/70C Mon't Kiara 50480 Kuala Lumpur

Malaysia Tel: +60 3 6201 7177 Fax: +60 3 6201 2177 Email: reception@resmed.com.my

ResMed NZ Ltd PO Box 51-048 Pakuranga Auckland New Zealand

Tel: +64 274 737 633 Fax: +64 9 239 0193 Email: reception@resmed.co.nz

ResMed SA Parc de la Bandonnière 2, rue Maurice Audibert 69800 Saint-Priest France

Tél: +33 (0) 4 37 251 251 Fax: +33 (0) 4 37 251 260 Email: reception@resmed.fr

ResMed Singapore Pte Ltd 57 Ubi Ave 1 #07-09 Ubi Centre Singapore 408936

Tel: +65 284 7177 Fax: +65 284 7787 Email: reception@resmed.com.sg

ResMed Spain SL C/Arturo Soria, 245 28033 Madrid España

Tel: +34 (93) 5908154 Fax: +34 (93) 5908153 Email: angelo@resmed.es

ResMed Sweden AB Industrigatan 2 S-461 37 Trollhättan Sverige

Tel: +46 520 420 110 Fax: +46 520 397 15 Email: reception@resmed.se

Labhardt AG Thannerstrasse 57 CH-4054 Basel Schweiz

Tel: +41 (061) 307 9711 Fax: +41 (061) 307 9722 Email: info@labhardt.ch

US DESIGNATED AGENT: ResMed Corp

EU AUTHORIZED REPRESENTATIVE: ResMed (UK) Ltd

Internet: www.resmed.com

Protected by patents: AU 697652, AU 699726, AU 713679, EP 0661071, US 4944310, US 5199424, US 5522382, US 6213119, US 6240921. Other patents pending.

Protected by design registrations: AU 147283, AU 147335, AU 147336, CH 128.709, CH 128.710, CH 128.711, CH 128.712, GB 3001791, GB 3001819, GB 3001820, GB 3001821, JP 1164087, JP 1164265, JP 1164266, JP 1164267, US D467335, US D468011. Other designs pending.

VPAP, HumidAire, HumidAire 2i, IPAP Max, IPAP Min, Mirage, SmartStart, Smart Data, *Ultra* Mirage, and Vista are trademarks of ResMed Ltd.

CONTENTS

INTRODUCTION	1
DEFINITIONS	1
USER/OWNER RESPONSIBILITY	1
MEDICAL INFORMATION	1
WHAT THE VPAP® III AND VPAP III ST ARE FOR	1
TELL YOUR DOCTOR IF...	1
WARNINGS	2
CAUTIONS	3
THE VPAP SYSTEM	5
VPAP COMPONENTS	5
MASKS	6
HUMIDIFIER	6
ACCESSORIES	7
PREPARING FOR USE	9
SETTING UP THE VPAP	9
FEATURES OF THE VPAP	14
OPERATING INSTRUCTIONS	21
STARTING TREATMENT	21
STOPPING TREATMENT	23
USING THE MASK-FITTING FEATURE	23
HELPFUL HINTS	24
CLEANING AND MAINTENANCE	27
DAILY	27
WEEKLY	27
PERIODICALLY	27
REPLACING THE AIR FILTER	28
SERVICING	28
TROUBLESHOOTING	29
SYSTEM SPECIFICATIONS	31
LIMITED WARRANTY	39
INDEX	41

INTRODUCTION

DEFINITIONS

This manual contains special terms and icons that appear in the margins to draw your attention to specific and important information.



WARNING

Alerts you to possible injury.



CAUTION

Explains special measures for the safe and effective use of the device.

Note: *Is an informative or helpful note.*

USER/OWNER RESPONSIBILITY

The user or owner of this system shall have sole responsibility and liability for any injury to persons or damage to property resulting from:

- operation which is not in accordance with the operating instructions supplied
- maintenance or modifications carried out unless in accordance with authorized instructions and by authorized persons.

Please read this manual carefully before use.

MEDICAL INFORMATION

WHAT THE VPAP® III AND VPAP III ST ARE FOR

The VPAP® III and VPAP III ST systems are intended for the treatment of adult patients with obstructive sleep apnea (OSA).



CAUTION (USA ONLY)

Federal law restricts this device to sale by or on the order of a physician.

TELL YOUR DOCTOR IF...

You must tell your doctor, and CPAP (continuous positive airway pressure) or bilevel therapy must not be used, if you have any of the following conditions:

- pneumothorax or pneumomediastinum
- severe heart failure, low blood pressure, or dehydration
- surgery to the brain, middle or inner ear, pituitary gland, or sinuses
- respiratory distress syndrome

- middle ear infection or perforated ear drum
- severe nosebleed.

WARNINGS

CPAP or bilevel therapy should be used with caution if you have any of the following conditions:

- respiratory failure
- cavities or cysts in the lung, or previous history of pneumothorax
- previous history of severe nosebleed
- sinus infection.

Tell your doctor if you have any of these conditions. Your doctor will advise you whether the likely benefits of CPAP or bilevel therapy outweigh the expected risks.

Special care should be exercised if you are dehydrated, or may become dehydrated, for example as a result of fluid restriction or diuretic therapy (including changes in therapy).

Discontinue therapy and seek medical advice if, during therapy or when you start therapy each night, you feel faint or light-headed.

The following are general warnings that pertain to your use of a **VPAP** unit. Specific warnings appear next to the relevant instructions in the manual.



WARNING

- **This is NOT a life support device.** It may stop operating with power failure or if a fault occurs in the unit.
- The air flow for breathing produced by this device can be as much as 11°F (6°C) higher than the temperature of the room. Caution should be exercised if the room temperature is warmer than 90°F (32°C).
- The **VPAP** unit should only be used with masks (and connectors)* recommended by ResMed, or by your physician or respiratory therapist. A mask should not be used unless the **VPAP** unit is turned on and operating properly. The vent hole or holes associated with the mask should never be blocked.

Explanation: The **VPAP** unit is intended to be used with special masks (or connectors)* which have vent holes to allow continuous flow of air out of the mask. When the device is turned on and functioning properly, new air from the device flushes the exhaled air out through the mask vent holes. However, when the device is not operating, insufficient fresh air will be provided through the mask, and the exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation. This applies to most models of CPAP or bilevel devices.

* Ports may be incorporated in the mask or in connectors that are near the mask.

- The **VPAP** unit should only be connected to the components, humidifiers, or accessories specified in this manual. Connection of other items may result in injury, or damage to the **VPAP** unit.

- At low pressures, the flow through the exhalation ports of your mask may not clear all exhaled gas from the tubing. Some rebreathing may occur.
- Explosion hazard—do not use in the vicinity of flammable anesthetics.
- If oxygen is used with this device, the oxygen flow must be turned off when the device is not operating.

Explanation of the warning: When the device is not in operation and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire. This warning applies to most types of flow generators.

- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- Always ensure air flow is being generated by the device before the oxygen supply is turned on.
- Always turn the oxygen supply off before stopping the air flow from the device.

Note: At a fixed rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on where the oxygen is introduced, pressure settings, patient breathing pattern, mask selection, and leak rate. This applies to most types of flow generators.

WARNINGS RELATED TO TREATMENT



WARNING

- If you stop your CPAP or bilevel treatment, your sleep apnea will return immediately.
- Always consult your clinician if you expect to be in a situation where you cannot use your **VPAP** unit.
- If you are admitted to a hospital or prescribed any other form of medical treatment, always inform the medical staff that you are being treated with CPAP or bilevel. It is also important to contact the clinician who is treating you for sleep apnea.
- If you experience an infection of the upper respiratory tract, middle ear, or sinuses, contact your clinician before continuing your CPAP or bilevel treatment. You may be advised to stop treatment until the infection has cleared. If you continue with treatment during an infection, be sure to clean your mask and tubing after every use.

CAUTIONS

The following are general cautions. Specific cautions appear next to the relevant instructions in the manual.



CAUTION

- Do not open the **VPAP** case. There are no user serviceable parts inside. Repairs and internal servicing should only be performed by an authorized service agent.

POSSIBLE SIDE EFFECTS

The **VPAP** flow generator is designed to help you have a good night's sleep. However, you need to be aware of possible problems that may arise during CPAP or bilevel treatment.



WARNING

Consult your clinician immediately if you experience any of the following symptoms during your CPAP or bilevel treatment:

- headache
- middle ear or sinus discomfort
- chest pain
- dryness of the nose, mouth, or throat
- feeling bloated due to air swallowing
- air continually leaking out of the mouth while sleeping
- recurrence of any sleep apnea symptoms while on CPAP or bilevel.



Skin irritation may occur from sensitivity to the mask materials or from excessively tight headgear straps. A correctly fitted mask and appropriately adjusted straps will often prevent skin irritation. If problems persist, contact your clinician for advice.

THE VPAP SYSTEM

VPAP COMPONENTS

Please identify and familiarize yourself with the following components of the VPAP unit:


VPAP III or VPAP III ST Unit



Power cord



Carry bag



Air tubing
6ft6in (2m)




WARNING

- Do not connect any device to the auxiliary port. Although your health care provider may connect specially designed devices to the auxiliary port of the VPAP unit, connection of other devices could result in injury, or damage of the unit.
- In the home environment the only device that may be connected to the communications port is a modem that is locally approved. Locally approved modems may also be connected in the clinical environment.
- In the clinical environment any PC that is used with the VPAP system must be at least 1.5m (5ft) away from, or at least 2.5m (8ft) above the patient. It must also comply with IEC 60950 or equivalent.

MASKS

You will also need a **ResMed mask system** (supplied separately).


The following ResMed mask systems are recommended for use with the **VPAP**:




MIRAGE VISTA™
MASK




ULTRA MIRAGE™
MASK



MIRAGE® MASK



MIRAGE FULL FACE
MASK SERIES 2




MODULAR MASK


Note: ResMed VPAP devices have been designed and manufactured to provide optimum performance using ResMed vented mask systems. Other mask systems may be used, however performance and data outputs may be affected. To select an appropriate setting for another mask system, find the closest match to a ResMed mask in the “Mask Flow/Pressure Characteristics” on page 34.

HUMIDIFIER


A humidifier may be required if you are experiencing dryness of the nose, throat, or mouth. The **VPAP** is compatible for use with the following humidifiers.




HUMIDAIRE 2i™ heated humidifier



HUMIDAIRE 2iC passover humidifier



HUMIDAIRE® heated humidifier




ResMed PASSOVER humidifier

**WARNING**


Only the HumidAire 2i, HumidAire 2iC, HumidAire heated humidifier, and the ResMed Passover are compatible for use with the VPAP. Please refer to Warnings on page 2.

ACCESSORIES

The following accessories are available for use with the VPAP:

VPAP ACCESSORY

Air tubing,
9ft10in (3m)

HUMIDIFIER ACCESSORY (HUMIDAIRE AND RESMED PASSOVER ONLY)

Medium air tubing
21in (52cm)

PREPARING FOR USE

SETTING UP THE VPAP


- 1 Place the **VPAP** unit on a flat surface near the head of your bed. If the unit is placed on the floor, ensure that the area is free from dust and clear of bedding, clothes, or any other objects that could block the air inlet.



CAUTION

Be careful not to place the device where it can be bumped or where someone is likely to trip over the power cord.

- 2 Connect the power cord to the socket at the rear of the flow generator. Plug the other end of the power cord into a power outlet.




WARNING

Make sure the power cord and plug are in good condition and the equipment is not damaged.

The air filter cover protects the device in the event of accidental liquid spillage onto the device. Ensure that the air filter and air filter cover are fitted at all times.

- 3 Connect one end of the air tubing firmly onto the air outlet of the unit.






WARNING

Only ResMed air tubing should be used with your flow generator. A different type of air tubing may alter the pressure you actually receive reducing the effectiveness of your treatment.

-
- 4 Assemble your mask system according to the mask user instructions.
-
- 5 Connect your mask system to the free end of the air tubing.



The **VPAP** is now ready for use. To start treatment, see “Operating Instructions” on page 21.

HUMIDIFIER USE




WARNING

When using a humidifier, position it lower than you, and at the same level or lower than the VPAP.


HUMIDAIRE 2i™

The HUMIDAIRE 2i™ attaches to the front of the VPAP to provide heated humidification. No other accessories are required for its use. The VPAP automatically detects the presence of the HUMIDAIRE 2i. No menu changes are required. Please refer to the *HumidAire 2i User's Manual* for details.



HUMIDAIRE 2iC

The HUMIDAIRE 2iC attaches to the front of a VPAP unit to provide passover humidification. No other accessories are required for its use. Please refer to the *HumidAire 2iC User's Manual* for details.




HUMIDAIRE® AND RESMED PASSOVER HUMIDIFIERS


Medium size (21 in) air tubing is a necessary accessory for connecting the **VPAP** unit to the **HUMIDAIRE®** and ResMed **PASSOVER** humidifiers.

To set up the **VPAP** with the **HUMIDAIRE** or ResMed **PASSOVER**:

- 1 Fill the **HUMIDAIRE** or **PASSOVER** with water as described in the humidifier manual.
- 2 **HumidAire Users** Place the filled water chamber inside the **HUMIDAIRE**. Connect the medium (21 in) air tubing to the right connector port, and the long air tubing (6ft6in or 9ft10in) to the left connector port on the humidifier. Close the **HUMIDAIRE** lid.
ResMed Passover Users Connect the medium (21 in) air tubing to the right connector port, and the long air tubing (6ft6in or 9ft10in) to the left connector port on the humidifier.




HumidAire




ResMed Passover

- 3 Place the **VPAP** on top of the **HUMIDAIRE** or **PASSOVER**. Do not place the **VPAP** unit underneath the humidifier. (This is to avoid water spilling into the unit.)




HumidAire




ResMed Passover

- 4 Connect the free end of the medium air tubing to the air outlet of the **VPAP**.




HumidAire




ResMed Passover

- 5 Connect the mask system to the free end of the long air tubing. The final assembly should look like this:



HumidAire



ResMed Passover

- 6 **HumidAire Users** Plug the HUMIDAIRE power cord into a power outlet.

- 7 Connect the power cord to the socket at the rear of the **VPAP**. Plug the other end of the power cord into a power outlet.



WARNING

Make sure that the power cord and plug are in good condition and the equipment is not damaged.


- 8 Navigate to the humidifier setting (if available) in the **VPAP** menu. See “Detailed Menu” on page 16.

The **VPAP** is now ready for use with the **HUMIDAIRE** or the ResMed **PASSOVER**. To start treatment, see “Operating Instructions” on page 21.





FEATURES OF THE VPAP

LCD SCREEN AND KEYPAD

The control panel of the **VPAP** includes an LCD screen and keypad.



The **VPAP** keypad has the following keys:

Key	Function
Front 	<ul style="list-style-type: none"> Starts or stops treatment. Extended hold for at least 3 seconds starts the Mask-Fitting feature.
Up/Down 	<ul style="list-style-type: none"> Allows you to scroll through the VPAP menus, submenus, and setting options.
Left 	<ul style="list-style-type: none"> Performs the function indicated by the guiding text displayed above it on the LCD screen. Guiding text includes menu, enter, change, and apply.
Right 	<ul style="list-style-type: none"> Performs the function indicated by the guiding text displayed above it on the LCD screen. Guiding text includes exit and cancel.

To assist you in adjusting the **VPAP**, the keypad and LCD are equipped with a backlight. The LCD backlight comes on when the unit is turned on or when you press a key, and turns off after 2 minutes. The keypad backlight is on at all times when the **VPAP** is on.

USING THE MENUS

The **VPAP** unit provides a set of functions which are arranged in menus and submenus. Via the LCD screen, the menus and submenus allow you to view and change the settings for a particular function. You can access the menus at any time, regardless of whether the **VPAP** is in stand-by mode, delivering therapy, or downloading data. After the Welcome screen appears and device self-checks are complete, the **VPAP** (or Ramp) screen appears.

RAMP SCREEN

If your clinician has set a maximum ramp time, the Ramp screen is displayed after the Welcome screen. On the Ramp screen, you can immediately set a ramp time. Ramp time is the period during which the pressure increases from a low pressure to the prescribed treatment pressure. See “Ramp time” on page 22.

Ramp time can be altered in 5 minute increments (from OFF to a maximum ramp time set by your clinician) by the using the **Up/Down** key.

MENU TYPE

Depending on the type of menu that your clinician has set for your machine, either a *standard* menu or a *detailed* menu appears.

STANDARD MENU

The standard menu allows you to view details about the run hours and the current software version of your **VPAP**. The Run Hours screen displays the total number of hours the **VPAP** has been used.

Figure 1 summarizes the **VPAP** standard menu series.

- To access the **VPAP** menus:
Press the **Left** key (menu) while the **VPAP** (or Ramp) screen is displayed.
- To scroll through items within the menu:
Press the **Up/Down** key
- To exit out of the menu:
Press the **Right** key (exit)




Figure 1: VPAP Standard Menu Series

DETAILED MENU

The detailed menu allows you to view and change settings such as mask type, tube length, and the humidifier used. You can also view the run hours, serial number, and current software version of your **VPAP**.

Figure 2 summarizes the **VPAP** detailed menu series.

- To access the **VPAP** menus:
Press the **Left** key (menu) while the **VPAP** (or Ramp) screen is displayed.
- To scroll through items within a menu or submenu:
Press the **Up/Down** key
- To enter a submenu:
Press the **Left** key (enter).
- To change a setting option for a function:
 1. Press the **Left** key (change)
 2. Press the **Up/Down** key until the desired setting option appears.
 3. Press the **Left key** (apply) to select the setting option.
- To exit without changing options:
Press the **Right** key (cancel)
- To exit out of a menu or submenu:
Press the **Right** key (exit)

Note: You can return to the **VPAP** (or Ramp) screen at any time by holding the **Right** key for at least 3 seconds.




Figure 2: VPAP Detailed Menu Series

MENU FUNCTIONS (DETAILED MENU ONLY)

The **VPAP** menu functions are summarized in Tables 1–4 below with a brief description of what each function does and the available setting options. To access these functions see “Detailed Menu” on page 16.

SETTINGS MENU

The Settings Menu allows you to view and change certain operating features of the **VPAP** unit.

Table 1: Settings Menu Functions

Function	Function Description	Setting Options
Mask	Selects your mask type.	MIRAGE, MIR FULL (Full Face), STANDARD (Vista, Modular), ULTRA (Ultra Mirage)
Tube Length	Selects the length of air tubing connecting your mask to the VPAP.	2m / 3m (6ft6in / 9ft10in)
Humidifier	Selects the type of humidifier to be used with the VPAP.	NONE, PASSOVER, HUMIDAIRE If the HumidAire 2i is used, the above setting options do not appear and H2i is displayed.
SmartStart®	Turns the SmartStart function on or off. See “SmartStart®” on page 22.	ON/OFF
Mask Alarm	Turns the Mask Alarm signal on or off. If enabled, the Mask Alarm feature will alert you when a high mask leak is detected. An audible tone will sound and a high leak message will appear on the LCD screen.	ON/OFF

Note: If you select “Mir Full” as the mask option, SmartStop is automatically disabled. SmartStart may not work with a Mirage Full Face Mask due to safety features of the mask.

When Mask Alarm is set to ON, SmartStart/Stop automatically reverts to OFF. SmartStop cannot be used with Mask Alarm because if a high leak occurs, SmartStop will stop treatment before the Mask Alarm signal is activated.

RESULTS MENU

Note: These menus appear only if at least one Smart Data option has been enabled by the clinician. See the VPAP III Smart Data Diary for further details.

Table 2: Results Menu

Function	Function Description	Setting Options
Mask Fit (Smart Data)	Displays a star rating corresponding to the mask leak from the previous session.	View only
Usage (Smart Data)	Displays usage hours from the previous session.	View only

OPTIONS MENU

Table 3: Options Menu

Function	Function Description	Setting Options
Smart Data – Auto Appear	The Smart Data menu is displayed only if one or more of the options have been set to ON by the clinician. If Auto Appear is set to ON, the Smart Data screens are displayed in the morning if you reset the device. If Auto Appear is set to OFF, Smart Data is displayed in the Results menu only.	ON/OFF
Language	Selects the language the VPAP uses for all its display text. English is the default language.	English, German, French, Italian, Spanish, Portuguese, Swedish, Dutch.

SERVICING MENU

Table 4: Servicing Menu

Function	Function Description	Setting Options
Run hours*	Displays the total number of machine hours.	View only
Serial Number (SN)	Displays the serial number for the VPAP.	View only
Printed Circuit Board (PCB)	Displays the printed circuit board number.	View only
Software*	Displays the current software version installed in the VPAP.	View only


* These items also appear in the Standard Menu.

OPERATING INSTRUCTIONS

STARTING TREATMENT

The **VPAP** unit should be assembled beside your bed with the air tubing and mask system connected. See “Setting Up The VPAP” on page 9.

- 1 Turn the main power switch at the back of the unit to on (I).



When the **VPAP** is turned on, a welcome message is displayed on the LCD screen. The **VPAP** (or Ramp) screen then appears.

Note: If you have the HumidAire 2i attached, please refer to the HumidAire 2i User's Manual for operating instructions.


- 2 Fit your mask as described in the mask user instructions.
- 3 Lie down and arrange the air tubing so that it is free to move if you turn in your sleep.



CAUTION

- Do not leave long lengths of air tubing around the top of your bed. It could twist around your head or neck while you are sleeping.
- Make sure the area around the flow generator is dry and clean. It should also be clear of bedding, clothes, and other potential blockages.

- 4 To start treatment, press the **Front** key
or
if your clinician has enabled the SmartStart function, simply breathe into the mask and treatment will begin.
After starting treatment, you can display one of the two treatment screens below. Press the **Up/Down** key to switch between views.



Following are the descriptions of the treatment screens.

Treatment mode: Mode of treatment set by your clinician. Options include: CPAP, Spontaneous, Spontaneous/Timed (VPAP III ST), and Timed (VPAP III ST).

Ramping indicator: Appears if the VPAP is in ramp mode. This disappears once the ramp time has elapsed.

Set pressure(s): In CPAP mode, it is the set treatment pressure (centimetres of water). In other modes, it is exhalation and inhalation pressures (centimetres of water).

Trigger indicator: How the VPAP changes the pressure when you are inhaling. "S" (Spontaneous) indicates a patient triggered change and "T" (Timed) indicates a device triggered change.

Pressure bar graph: Graphical display of the changing pressure.

Leak: Current mask leak (litres per minute).

Respiratory rate: Number of breaths per minute.

Minute ventilation: Volume of air inhaled per minute (litres per minute). It is the product of respiratory rate and tidal volume.

Tidal volume: Volume of air inhaled per breath (millilitres per breath).

RAMP TIME

Ramp time is a feature which can be enabled by your clinician. If you have difficulty falling asleep with full pressure, select a ramp time. The airflow will start very gently while you fall asleep. The pressure will slowly increase to full operating pressure over the selected ramp time. The clinician has set a maximum ramp time; you may select any value up to the maximum.

SMARTSTART[®]

The VPAP has a function called SmartStart which can be enabled by your clinician. If SmartStart is enabled, VPAP will start automatically when you breathe into the mask and will stop automatically when you take your mask off. This means you do not have to press the **Front** key to begin or end treatment.

STOPPING TREATMENT


To stop treatment at any time, remove your mask and press the **Front** key
or

if your clinician has enabled the SmartStart function, simply remove your mask and treatment will end.

USING THE MASK-FITTING FEATURE

The **VPAP** Mask-Fitting feature can be used to help you fit your mask properly. The mask-fitting feature delivers air pressure for a three-minute period, prior to starting treatment, for checking and adjusting your mask fit to minimize leaks. If a Ramp time is selected, the mask can be adjusted at a pressure closer to the prescribed pressure. To use the mask-fitting feature:

- 1 Fit your mask as described in the user instructions.
-
- 2 Hold down the Front key for at least 3 seconds until air pressure delivery starts. The following display will appear on the LCD screen indicating that the Mask-Fitting feature is in operation. The flow generator will ramp to the Mask-Fit pressure and will remain at this pressure for 3 minutes. A Mask-Fit star rating is also displayed. See “Definitions of Mask-Fit Star Rating” on page 24.



Notes

- The Mask-Fit star rating display disappears after 3 minutes.
- The Mask-Fitting feature can only be started from the VPAP (or Ramp) screen.
- The Mask-Fit pressure is the set treatment pressure or 10 cm H₂O, whichever is greater.

- 3 Adjust your mask, mask cushion, and headgear until you have a secure and comfortable fit.

Once you have a secure and comfortable fit, check your Mask-Fit star rating on the LCD screen. Definitions of the Mask-Fit star ratings are presented in Table 4.

Note: If there is another person nearby to check your Mask-Fit star rating, you can adjust your mask, mask cushion, and headgear while lying down.

- 4 After 3 minutes, treatment will begin. Definitions of the Mask-Fit star ratings are presented in Table 4.
 - If you do not wish to wait 3 minutes, hold down the **Front** key for at least 3 seconds and treatment will begin immediately.

- If you press the **Front** key for less than 3 seconds, the unit will return to standby mode (the **VPAP** or Ramp screen is displayed).

Table 5: Definitions of Mask-Fit Star Rating

Star rating	Definition
*****	Excellent
****—	Very good
***—	Good
**——	Adjust mask
*——	Adjust mask
HIGH LEAK	Adjust mask

HELPFUL HINTS

STARTING OUT

MOUTH LEAKS

If using a nasal mask, try to keep your mouth closed during treatment. Air leaks from your mouth can decrease the effectiveness of your treatment. If mouth leaks are a problem, a full face mask or chin strap may help. Contact your clinician or equipment supplier for further details.

MASK FITTING

The flow generator delivers the most effective treatment when the mask is well fitted and comfortable. Treatment can be affected by leaks, so it is important to eliminate any leaks that may arise.

If you have problems trying to get a comfortable mask fit, contact your sleep clinic or equipment supplier. You may benefit from a different size or style of mask.

You can also use the Mask-Fitting feature to help you fit your mask properly. See “Using the Mask-Fitting Feature” on page 23.

Before wearing your mask, wash your face to remove excess facial oils. This will allow a better fit and prolong the life of the mask cushion.

NASAL IRRITATION

DRYNESS

You may experience dryness of the nose, mouth, and/or throat during the course of treatment, especially during winter. In many cases, a humidifier may resolve this discomfort. Contact your clinician for advice.

RUNNY OR BLOCKED NOSE

You may experience sneezing and/or a runny or blocked nose during the first few weeks of treatment. In many cases, nasal irritation can be resolved with a humidifier. Consult your clinician for advice.

TRAVELLING WITH THE VPAP

INTERNATIONAL USE

Your **VPAP** flow generator has an internal power adapter that enables it to operate in other countries. It will operate on power supplies of 100–240V and 50–60Hz. No special adjustment is necessary, but you may need a plug adapter for the power outlet.

CLEANING AND MAINTENANCE

You should regularly carry out the cleaning and maintenance described in this section.

DAILY

1. Disconnect the air tubing and hang it in a clean, dry place until next use. Do not hang the air tubing in direct sunlight as it may harden and crack over time.
2. Clean the mask according to the mask user instructions.
3. If you are using a humidifier, clean it according to the instructions in the manual.

WEEKLY

1. Remove the air tubing from the **VPAP** unit and the mask.
2. Wash the mask system according to the instructions supplied with it.
3. Wash the air tubing in warm water using mild detergent. Rinse thoroughly, hang and allow to dry.
4. Before next use, assemble the mask and headgear according to the mask user instructions.
5. Reconnect the air tubing to the air outlet and mask.



CAUTION

- Do not use bleach, chlorine-, alcohol- or aromatic-based solutions (including all scented oils), moisturizing, or antibacterial soaps to clean the cushion, mask, air tubing, or the VPAP. These solutions may cause hardening and reduce the life of the product.
- Do not wash or dry the mask frame at a temperature above 176°F (80°C). Exposure to higher temperatures may reduce the life of the product.
- Do not hang the air tubing in direct sunlight as the tubing may harden over time and eventually crack.

PERIODICALLY

1. The mask and air tubing are subject to normal wear and tear. Inspect them regularly for damage.
2. Clean the exterior of the flow generator with a damp cloth and mild detergent.
3. Inspect the air filter to check if it is blocked by dirt or contains holes. See “Replacing the Air Filter” on page 28.



WARNING

Beware of electric shock. Do not immerse the flow generator or power cord in water. Always unplug the flow generator before cleaning and be sure that it is dry before reconnecting.




CAUTION

Do not attempt to open the VPAP. There are no user serviceable parts inside. Repairs and internal servicing should only be performed by an authorized service agent.

REPLACING THE AIR FILTER

Inspect the air filter every month to check if it is blocked by dirt or contains holes. With normal use of a VPAP unit, the air filter needs to be replaced every six months (or more often if your unit is in a dusty environment). To replace the air filter:

1. Remove the air filter cover at the back of the VPAP.



2. Remove and discard the old air filter.
3. Insert a new filter with the blue tinted side facing out.
4. Replace the air filter cover.



WARNING

Do not wash the air filter. The air filter is not washable or reusable.

Note: The air filter should be inspected once a month.

SERVICING

Your VPAP flow generator is designed to give you years of trouble-free operation.

The flow generator should not require regular servicing if it is maintained according to the instructions in this manual. If you feel that your unit is not performing properly, see "Troubleshooting" on page 29.



CAUTION

Inspection and repair should only be performed by an authorized agent. Under no circumstances should you attempt to service or repair the flow generator yourself.

TROUBLESHOOTING

If there is a problem, try the following suggestions. If the problem cannot be solved, contact your equipment supplier or ResMed. Do not attempt to open the unit.

Problem	Possible Cause	Solution
<i>No display.</i>	Power not connected or switch at back is not on.	Ensure the power cable is connected and that the switch at the back of the unit is in the ON position.
<i>Insufficient air delivered from the VPAP.</i>	Ramp Time is in use.	Wait for air pressure to build up.
	Air filter is dirty.	Replace air filter.
	Air tubing is kinked or punctured.	Straighten or replace tubing.
	Air tubing not connected properly.	Check air tubing.
	Mask and headgear not positioned correctly.	Adjust position of mask and headgear.
	Plug(s) missing from access port(s) on mask.	Replace plug(s).
	Pressure required for treatment may have changed.	See your clinician to adjust the pressure.
<i>The VPAP does not start when you breathe into the mask.</i>	Power cord not connected properly.	Connect power cord firmly at both ends.
	Power outlet may be faulty.	Try another power outlet.
	The VPAP unit not switched on.	Switch power switch at rear of the VPAP to ON.
	SmartStart not on.	Enable SmartStart.
	Mask Alarm has been enabled; SmartStart has automatically been disabled.	Disable Mask Alarm to enable SmartStart.
	Use of a Mirage Full Face Mask.	SmartStart does not work with a full face mask as the anti asphyxia valve will not allow sufficiently high pressure on exhalation.

Problem	Possible Cause	Solution
	Breath is not deep enough to trigger SmartStart.	Take a deep breath in and out through the mask.
	There is excessive leak.	Adjust position of mask and headgear. Plugs may be missing from ports on mask. Replace them. Air tubing not connected properly. Connect firmly at both ends. Air tubing kinked or punctured. Straighten or replace.
VPAP unit does not stop when you remove your mask.	SmartStart/Stop is disabled.	Enable SmartStart/Stop.
SmartStart is enabled but the flow generator does not stop automatically when you remove your mask.	Incompatible humidifier or mask system being used.	Use only equipment as recommended and supplied by ResMed.
Display error message: Check tube!! Key if done	The air tubing is loose.	Check that the air tubing is connected securely to your mask and the air outlet on the front of the VPAP. To clear the error message, press any key on the VPAP keypad.
Displays error message: SYSTEM ERROR Call service!	Component failure.	Return your VPAP for servicing.
Excessive motor noise.	Component failure.	Return your unit for servicing.
Display error message: High leak in last session.	You have experienced excessively high leak levels during the night.	Check that your air tubing is connected properly and that your mask does not leak excessively. Use the mask-fitting feature to help you to fit your mask properly. If this message appears again, contact your clinician.

SYSTEM SPECIFICATIONS

Dynamic pressure characteristics

IPAP: 4 cm H₂O to 25 cm H₂O (measured at the end of standard 2m air tubing)

EPAP: 4 cm H₂O to 25 cm H₂O (measured at the end of standard 2m air tubing)

Maximum single fault pressure: 40 cmH₂O

Dynamic flow characteristics

130 L/min at 4 to 20 cm H₂O

Sound pressure level: <30 dB (tested in accordance with the requirements of ISO 17510-1:2002)

Dimensions (H x W x D): 5.6in x 9.1in x 10.6in

Weight: 5.1lb

Air outlet: 22mm taper, compatible with EN 1281-1:1997 Anaesthetic & Respiratory Equipment - Conical Connectors

Pressure measurement: Internally mounted pressure transducer

Flow measurement: Internally mounted flow transducer

Power Supply: Input range 100–240V, 50–60Hz, 40VA (typical power consumption), < 100VA (maximum power consumption)

Housing Construction: Flame retardant engineering thermoplastic

Environmental Conditions

Operating Temperature: +41°F to +104°F

Operating Humidity: 10%–95% non-condensing

Storage and Transport Temperature: -4°F to +140°F

Storage and Transport Humidity: 10%–95% non-condensing

Electromagnetic Compatibility

Product complies with all applicable electromagnetic compatibility requirements (EMC) according to IEC60601-1-2, for residential, commercial, and light industry environments. For further details, see “Guidance and Manufacturer’s Declaration - Electromagnetic Emissions and Immunity” on page 35.

Air Filter: Two-layered, powder-bonded, polyester non-woven fiber

Air Tubing: Flexible plastic, 6ft6in or 9ft10in length

IEC 60601-1 Classifications

Class II (double insulation)

Type CF

Table 6: Displayed values

Value	Range	Accuracy	Display Resolution
Pressure sensor at air outlet			
Pressure	-5 to 30 cm H ₂ O	±0.5 cm H ₂ O	0.1 cm H ₂ O
Flow sensor in flow generator*			
Leak	0–120 L/min	**	1 L/min
Tidal volume	100–3000 mL	**	1 mL
Respiratory rate	6–60 BPM	±0.5 BPM	0.1 BPM
Minute ventilation	0.6–60 L/min	**	0.25 L/min


* Results may be inaccurate in the presence of leaks.

** The displayed values are estimates. They are provided for trending purposes only.


Table 7: Maximum low impedance flow at stated pressures

Pressure (cm H ₂ O)	Flow (L/min)
6.6	177
13.2	202
20.0	202

PRESSURE VARIATION



PRESSURE VOLUME CURVE



Note: The manufacturer reserves the right to change these specifications without notice.

SYMBOLS WHICH APPEAR ON THE PRODUCT



Attention, consult accompanying documents



Class II equipment



Type CF equipment




IPX1

Drip Proof



Start/Stop or Mask-Fit

MASK FLOW/PRESSURE CHARACTERISTICS



GUIDANCE AND MANUFACTURER'S DECLARATION - ELECTROMAGNETIC EMISSIONS AND IMMUNITY

Guidance and manufacturer's declaration – electromagnetic emissions

The VPAP is intended for use in the electromagnetic environment specified below. The customer or the user of the VPAP should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11	Group 1	The VPAP uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The VPAP is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to EMC information provided in this document.

Warnings: The VPAP should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the VPAP should be observed to verify normal operation in the configuration in which it will be used.

The use of accessories (eg humidifiers) other than those specified in this manual is not recommended. They may result in increased emissions or immunity of the VPAP.

Guidance and manufacturer's declaration – electromagnetic immunity

The VPAP is intended for use in the electromagnetic environment specified below. The customer or the user of the VPAP should assure that it is used in such an environment.


Immunity test	IEC60601-1-2 test level	Compliance level	Electromagnetic environment –guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV Not Applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines. IEC 61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	< 12V (>95% dip in 240V) for 0.5 cycle 96V (60% dip in 240V) for 5 cycles 168V (30% dip in 240V) for 25 cycles <12V (>95% dip in 240V) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VPAP requires continued operation during power mains interruptions, it is recommended that the VPAP be powered from an uninterruptible power source
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

NOTE: Ut is the a.c. mains voltage prior to application of the test level.

(Continued next page)

Guidance and manufacturer’s declaration – electromagnetic immunity (Continued)

The VPAP is intended for use in the electromagnetic environment specified below. The customer or the user of the VPAP should assure that it is used in such an environment.

Immunity test	IEC60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>10 V/m 80 MHz to 2.5 GHz</p>	<p>3 Vrms</p> <p>10 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the VPAP, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.17 \sqrt{P}$</p> <p>$d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1: At 80 MHz and 800MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VPAP is used exceeds the applicable RF compliance level above, the VPAP should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VPAP.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the VPAP

The VPAP is intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the VPAP can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the VPAP as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800MHz to 2.5 GHz $d = 0.35 \sqrt{P}$
0.01	0.17	0.04	0.04
0.1	0.37	0.11	0.11
1	1.17	0.35	0.35
10	3.69	1.11	1.11
100	11.70	3.50	3.50

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

LIMITED WARRANTY

ResMed warrants that your ResMed product shall be free from defects in material and workmanship for the period specified below from the date of purchase by the initial consumer.

Product	Warranty Period
AutoSet CS™ flow generator, ResMed humidifiers, ResControl™, ResLink™.	1 Year
VPAP™ flow generator, CPAP flow generator, AutoSet T™ flow generator, AutoSet Spirit™ flow generator.	2 Years
Accessories, mask systems (including mask frame, cushion, headgear and tubing). Excludes single-use devices.	90 Days

If the product fails under conditions of normal use, ResMed will repair or replace, at its option, the defective product or any of its components. This Limited Warranty does not cover: a) any damage caused as a result of improper use, abuse, modification or alteration of the product; b) repairs carried out by any service organization that has not been expressly authorized by ResMed to perform such repairs; c) any damage or contamination due to cigarette, pipe, cigar or other smoke; d) any damage caused by water being spilled on or into the flow generator. Any product repaired or replaced under warranty will be returned, freight prepaid, to the dealer designated by the consumer. The cost of transporting the product to an authorized service organization will be borne by the consumer.

This warranty is in lieu of all other express or implied warranties, including any implied warranty of merchantability or fitness for a particular purpose. Some regions or states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you.

ResMed shall not be responsible for any incidental or consequential damages claimed to have occurred as a result of the sale, installation or use of any ResMed product. Some regions or states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from region to region.

For further information on your warranty rights, contact your local ResMed dealer or ResMed office.

INDEX

A

Accessories 7
 Air Filter 31
 Air filter 28
 Air Tubing 31

B

Back-light, LCD and Keypad 14

C

Cleaning and Maintenance 27
 Cleaning Periodically 27
 Components 5

D

Daily Cleaning 27
 Definitions 1
 Detailed Menu 15, 16, 17
 Dimensions 31
 Dryness 24
 Dynamic flow characteristics 31
 Dynamic pressure characteristics 31

E

Electromagnetic Compatibility 31
 Environmental Conditions 31
 Error Messages 30

F

Features of the VPAP 14
 Front Key 14

G

Glossary of Symbols 33

H

Helpful Hints 24
 High Leak Message 18
 Housing Construction 31
 HumidAire 2i 11
 HumidAire 2iC 11
 Humidifier 6, 11
 Humidifier Use 11

I

IEC 60601-1 Classifications 31
 International Use 25

K

Keys, functions 14

L

LCD Screen and Keypad 14
 Left Key 14

M

Mask Alarm 18
 Mask Fitting 24
 Mask Flow/Pressure Characteristics 34
 Mask-Fit Star Rating, definitions 24
 Mask-Fitting Feature 23
 Masks 6
 Medical Information 1
 Menu Functions 18
 Menu Type 15
 Menus, using the VPAP 15
 Mouth Leaks 24

N

Nasal Irritation 24

O

Operating Instructions 21
 Options Menu 19

P

Power Supply 31
 Preparing for Use 9

R

Ramp Screen 15
 Ramp time 15
 Replacing the Air Filter 28
 Responsibility, user/owner 1
 Results Menu 19
 Right Key 14
 Run Hours Screen 15
 Runny or Blocked Nose 24

S

Servicing 28
 Servicing Menu 19
 Setting Up 9
 Settings Menu 18
 Side Effects 4
 SmartStart 22

Standard Menu 15
Starting Treatment 21
Stopping Treatment 23
System Specifications 31

T

Travelling 25
Troubleshooting 29

U

Up/Down Key 14
User/Owner Responsibility 1
Using the Mask-Fitting feature 23
Using the Menus 15

W

Warnings Related to Treatment 3
Weekly Cleaning 27
Weight 31