



VivoVISION®
SYSTEMS

IVIS® Lumina System Manual

February 2006

XENOGEN®
Life. Changing.™

Xenogen Corporation
860 Atlantic Avenue
Alameda, California 94501
USA
Phone (510) 291-6100
Fax (510) 291-6196
www.xenogen.com



IVIS[®] Lumina System Manual

February 2006

For Research Use Only



© 2006 Xenogen Corporation. All rights reserved.

PN40150RevB

Xenogen Corporation

860 Atlantic Avenue
Alameda, California 94501, USA
Main Phone: 1.510.291.6100
Fax: 1.510.291.6196
E-mail: ivistechsupport@xenogen.com
www.xenogen.com

Xenogen Technical Support
1.888.810.8055 (US only)
1.510.291.6275 (International)

Discovery in the Living Organism, IVIS Imaging System and Living Image are either registered trademarks or trademarks of Xenogen Corporation. The names of companies and products mentioned herein may be the trademarks of their respective owners. Apple, Macintosh and QuickTime are registered trademarks of Apple Computer, Inc. Microsoft, PowerPoint and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Adobe and Illustrator are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries.

Contents

1 Welcome	1
1.1 Introduction	1
1.2 Xenogen Technical Support	2
2 Important Safety Instructions	3
2.1 Definitions	3
2.2 Instructions	4
2.3 Environmental Considerations for the System Components	4
2.4 Cleaning or Moving the System Components	5
2.5 Power Considerations	5
2.6 Servicing	6
2.7 Other Equipment	6
3 Warnings	7
3.1 Electrical Safety	7
3.2 Eye Safety and Burn Hazard	8
3.3 Mechanical Safety	8
3.4 Chemical & Biological Safety	9
3.5 Panels, Cover, and Modules	9
4 Legal Notices	11
4.1 Limited Warranty	11
4.2 Patents	12
4.3 Trademarks	13
5 IVIS® Lumina Components & Specifications	15
5.1 CCD Camera	16
5.2 Imaging Chamber	16
5.3 Optics	18
5.4 Optical Filter Wheel	18
5.5 Acquisition Computer	18
5.6 Environmental Requirements	19
6 Operating the IVIS® Lumina	21
6.1 System Start Up	21
6.2 Restarting the System After a Power Outage	22
6.3 Gas Plumbing	22
6.4 Door Operation	24
6.5 Imaging Basics	25
6.6 System Shut Down Procedure	26
7 Care & Maintenance	27
7.1 Cleaning the IVIS Lumina	27
8 Troubleshooting	29
8.1 Measured Temperature Does Not Equal the Demand Temperature	29

8.2 Photographic Image Is Unacceptable	30
8.3 Luminescent Image Is Unacceptable	31
8.4 No Image Is Produced	32
9 XFO-12 Fluorescence Equipment	33
9.1 Installation Requirements	34
9.2 Specifications	34
9.3 Description and Theory of Operation	35
9.4 Acquiring Fluorescent Images	40
9.5 Troubleshooting	42
9.6 Care and Maintenance of the XFO-12 Fluorescence Equipment	46
9.7 XFO-12 Hardware Replacement Parts List	47
Appendix A: Moving IVIS® Lumina With the XWS-260 Workstation	49
Index	51

1 Welcome

Introduction	1
Xenogen Technical Support	2

1.1 Introduction

The Xenogen IVIS® Lumina is a high-sensitivity, low noise, *in vivo* imaging technology platform that enables noninvasive visualization and tracking of cellular and genetic activity within a living organism, in real time (Figure 1.1).



Figure 1.1 IVIS Lumina hardware

Specific genes, cells, or organisms are “tagged” with a gene encoding one of the luciferase enzymes that enable some bacteria, insects, and animals to glow. When the tagged entity is active, it glows. The emitted light corresponds to the number and location of the tagged entities. This information enables scientists to observe the spread of disease or the effects of a drug throughout the system, non-invasively.

The IVIS Lumina takes very low light level images that can be saved and displayed for subsequent analysis. In addition to bioluminescent imaging, the IVIS Lumina incorporates the XFO-12 fluorescence equipment that enables fluorescent imaging capabilities for both *in vitro* and *in vivo* applications. The sensitive range of the CCD camera sets the wavelength range for fluorescence applications from 400 to 950nm.

The IVIS Lumina is an integrated imaging system that includes:

- A scientific grade thermoelectrically cooled CCD camera mounted on a light-tight imaging chamber
- A camera power supply
- A Windows®-based computer system for data acquisition and analysis
- The XFO-12 fluorescence equipment

This manual explains how to operate and maintain the equipment, and provides guidelines for obtaining the best bioluminescent or fluorescent images. Before using the IVIS® Lumina, please read this manual carefully to obtain safe, optimum performance and a maximum service life from the unit.

For instructions on how to acquire fluorescent images using the XFO-12 fluorescence equipment, please see [Chapter 9, page 33](#). For instructions on using the system software, please see the *Living Image® Software User's Manual*.

If you have questions regarding this manual or the IVIS Lumina, please call Xenogen technical support.

1.2 Xenogen Technical Support

Technical Support:	1.888.810.8055 (US) 1.510.291.6275 (International)
Main Telephone:	1.510.291.6100
E-mail:	IVISTechSupport@xenogen.com
Fax:	1.510.291.6196
Address:	Xenogen Corporation 860 Atlantic Avenue Alameda, CA 94501

2 Important Safety Instructions

Definitions	3
Instructions	4
Environmental Considerations for the System Components	4
Cleaning or Moving the System Components	5
Power Considerations	5
Servicing	6
Other Equipment	6

2.1 Definitions



CAUTION

CAUTION! This symbol indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury and/or mechanical damage. It is also used to alert against unsafe practices. It reminds you that all safety instructions should be read and understood before installation, operation, maintenance, or repair of this instrument. When you see this symbol, pay particular attention to the safety information presented. Observance of safety precautions will help avoid actions that could damage or adversely affect the performance of the IVIS® Lumina. If the equipment is used in a manner not specified in this manual, the protection provided by the equipment may be impaired.



WARNING

WARNING! Used when an action or condition may potentially cause serious personal injury or loss of life. Mechanical damage may also result.



VOLTAGE

VOLTAGE! Indicates high voltage or risk of electric shock.

2.2 Instructions

Read Instructions	Read and understand all the safety and operating instructions before you install, operate, or perform maintenance on this product. Make sure that you fully understand the following safety instructions, warnings, and disclaimers before proceeding to the rest of the manual.
Retain Instructions	Retain the safety and operating instructions for future reference.
Follow Instructions	Follow all operating and handling instructions. Failure to follow operating or handling instructions may void any warranty covering this product.
Heed Warnings	Abide by all warnings on the product and in the operating instructions. Failure to adhere to warnings or safety precautions may void any warranty covering this product.

2.3 Environmental Considerations for the System Components

Location for the IVIS® Lumina	<p>Before the IVIS Lumina is installed, consider the proper environment for the components.</p> <p>Install the equipment in an environment where:</p> <ul style="list-style-type: none">• The temperature does not fluctuate widely and is maintained between 15-25° C.• The humidity does not exceed 80%.• No strong electric or magnetic fields exist.• No vibrations are present.• No corrosive gases are present.• High amounts of dust are not present.• No open flame is present.• There is sufficient space behind the IVIS Lumina equipment. A minimum space of four inches from the flat surface of the rear panel should be provided behind the IVIS Lumina to provide unobstructed air flow and access to the main power on/off switch.• The work space is level.
Heat	The system should be situated away from heat sources such as open flames, radiators, heat registers, stoves, and other heat-generating electrical equipment.

Water & Moisture



VOLTAGE

VOLTAGE! Do not use this product near water (for example, near a sink or wet room) due to risk of electric shock, electrical damage, and/or failure.

2.4 Cleaning or Moving the System Components

Cleaning/Liquid Entry



VOLTAGE

VOLTAGE! Do not use liquid or aerosol cleaners and never spill liquid of any kind on any of the IVIS® Lumina components. Sprays and liquids that come into contact with the IVIS Lumina hardware may result in damage to the system or electrocution. For more details on proper care of the system, see [Cleaning the IVIS Lumina, page 27](#).

Moving the IVIS Lumina

You can move the light source module on the laboratory bench within the extent of the fiber optic cable. Be careful not to bend the fiber optic cable (minimum radius of bending: 3 inches/7.6 cm).



CAUTION

CAUTION! The IVIS Lumina is sensitive, scientific equipment and should not be moved by any user unless the system is located on the XWS-260 workstation. Due to the risk of potential damage, it is critical that only a trained technician moves the IVIS Lumina. If it is necessary to relocate the IVIS Lumina, contact Xenogen technical support. For more details on moving the IVIS Lumina on the XWS-260 workstation, see [page 49](#).

2.5 Power Considerations

Power Sources

The IVIS Lumina is configured for the voltage requirements of the installation locality that was specified at the time of order. If the IVIS Lumina is moved to another area, check to make sure that the same voltage requirements exist.



VOLTAGE

VOLTAGE! The IVIS Lumina can operate at multiple voltages (100-240 VAC); however, you are **not permitted** to change the input voltage to any of the system components. Several internal modifications are required for voltage change when carried out by Xenogen technical support. If the operating voltage must be changed, contact Xenogen technical support.

Power Cord Protection

Power supply cords should be routed so that they are unlikely to be walked on or pinched by items placed upon or against them. Pay close attention to receptacles and to points of connection between cords and equipment.

Lightning & Power Line Surges

The IVIS® Lumina is supplied with a surge protector. All components should be connected to this device to protect against electrical transient events. Failure to isolate the camera from electrical transients may result in damage to the CCD camera.

Power Outages

If the IVIS Lumina experiences a loss of supply power, turn off the power switch for all components and do not restart the system until reliable power has been restored.

Overloading



WARNING

WARNING! Do not overload wall outlets, extension cords, or integral convenience receptacles as this can result in a risk of fire or electric shock. For more details on the power requirements of the IVIS Lumina equipment, see [Chapter 5, page 15](#).

Facilities should be adequately wired according to local building codes.

2.6 Servicing

Refer all servicing to Xenogen technical support. If the product is damaged and requires service, unplug the product from the outlet and contact Xenogen technical support. Servicing by anyone other than those authorized by Xenogen voids any warranty covering this product.

2.7 Other Equipment

Use of any equipment other than that recommended by this manual has not been evaluated for safety and, therefore, is the sole responsibility of the user.

3 Warnings

Electrical Safety	7
Eye Safety and Burn Hazard	8
Mechanical Safety	8
Chemical & Biological Safety	9
Panels, Cover, and Modules	9

3.1 Electrical Safety



VOLTAGE

VOLTAGE! DO NOT attempt to service the IVIS Lumina yourself. Although there are no voltages in excess of 24V inside the imaging chamber, local line voltages (110VAC or 230VAC) are present in the lower electronics tray. The light source module may be user serviced for line fuse and lamp replacement only. There are no other user serviceable electrical parts in the light source module with the exception of the line fuse and tungsten halogen lamp. See [Troubleshooting, page 42](#), for instructions on these procedures. Contact Xenogen Technical Service for other electrical service needs.



WARNING

WARNING! If necessary, wipe exterior surfaces of the light source module with a soft cloth. DO NOT use fluids to clean the exterior of the module. Do not allow fluids of any kind to enter the light source module under any circumstances.



WARNING

WARNING! When the power is on, DO NOT disconnect or reattach the electrical control cable that connects the XFO-12 fluorescence equipment (excitation filter assembly) to the light source module and the IVIS® Lumina (electronics tray). See [System Components, page 35](#) for photographs of these components. Disconnecting or reconnecting the control cable when the system has electrical power will damage the system. Always turn off the switch on the front panel of the light source module and the rear-mounted ON/OFF switch on the IVIS Lumina before you connect or disconnect any of these cable connections.

3.2 Eye Safety and Burn Hazard

The light source module and the connecting fiber optic cables produce intense light that can cause eye damage. The module uses a tungsten halogen lamp bulb that operates at a high temperature, which if exposed to a user's skin, could cause a burn.



WARNING

WARNING! DO NOT operate the light source module or the XFO-12 fluorescence equipment without all of the fiber optic cables connected at both of their end connections.



WARNING

WARNING! If you need to replace the tungsten halogen lamp or lamp assembly in the light source module, wait until the lamp and socket are cool. Confirm that the components are cool before proceeding. See [Lamp Replacement, page 44](#).

3.3 Mechanical Safety

The imaging chamber of the IVIS Lumina is heavy and weighs 90 lbs. (41 kg).

IMPORTANT



Do not move the IVIS Lumina after installation. If you need to move the IVIS Lumina, contact Xenogen technical support for assistance.

The IVIS Lumina has many internal motorized components that can move at any time. The imaging stage can move when the door is open. Care should be taken to keep hands and equipment away from the sides of the platform when it is moving. Never place anything underneath the platform.

Do not attempt to put anything into the lens opening of the camera as there are optical components that can be compromised or damaged.

If the imaging chamber makes an unusual noise or appears to be jammed, turn off the power switch located on the back of the chamber.



WARNING

WARNING! If the gas hoses become caught, kinked, or disconnected, do not operate the instrument. Overexposure to anesthesia gas may occur.

**CAUTION**

CAUTION! DO NOT touch or expose the four diffusing reflectors and the exposed emission filter to contaminants (described in [System Components, page 35](#) and shown in [Figure 9.2](#)), otherwise imaging performance may be impaired. The reflectors' surfaces have been surface treated for optimum light diffusion.

**CAUTION**

CAUTION! The Living Image® software controls excitation filter selection. Do not manually turn the numbered knob on the Excitation Filter Wheel Assembly (see [Figure 9.5, page 37](#)).

3.4 Chemical & Biological Safety

Normal operation may involve the use of test samples that are pathogenic, toxic, or radioactive. It is your responsibility to ensure that all necessary safety precautions are taken before such materials are used.

Dispose of all waste materials according to appropriate environmental health and safety guidelines.

It is your responsibility to decontaminate the IVIS® Lumina before requesting service by Xenogen technical support. Ask your laboratory safety officer to advise you about the level of containment required for your application and about the proper decontamination or sterilization procedures to follow.

Handle all infectious samples according to good laboratory procedures and methods to prevent the spread of disease.

3.5 Panels, Cover, and Modules

There are no user serviceable components in the lower electronics tray of the IVIS® Lumina. Do not remove the electronics tray from the IVIS Lumina or the cover from the light source module unless you are instructed by and under the supervision of a Xenogen technical service representative.

The only exception to this is removing the light source module front panel mounted lamp assembly for tungsten halogen lamp replacement (see [page 44](#)).

[This page intentionally left blank.]

4 Legal Notices

Limited Warranty	11
Patents	12
Trademarks	13

4.1 Limited Warranty

Xenogen Corporation ("Xenogen") provides the following limited warranty for each new IVIS® Lumina ("System") purchased from it as follows ("Limited Warranty"):

- i. This Limited Warranty for the System extends for a period of one (1) year following delivery to, and installation of, the System to the original customer, purchaser, or user ("Customer") and is not assignable or transferable to any successor.
- ii. During the Limited Warranty period, Xenogen will repair or replace, at Xenogen's sole option, any defective parts if such repair or replacement is needed because of System malfunction or failure during normal usage in accordance with the instructions in this manual. Repairs and replacements under the Limited Warranty will be made at Xenogen's expense. Xenogen's limit of liability under the Limited Warranty shall be the purchase price of the Imaging System. Xenogen shall not be liable for any other losses or damages. These remedies are the Customer's exclusive remedies for breach of warranty.
- iii. No coverage or benefits shall be provided under this Limited Warranty if any of the following conditions apply:
 - a. The System has been subjected to abnormal use, abnormal conditions, unauthorized modifications (e.g., unauthorized installation of hardware or software), unauthorized repair or servicing, misuse, neglect, abuse, accident, alteration, any use inconsistent with or in contradiction to the instructions in this manual, or other acts which are not the fault of Xenogen.
 - b. Xenogen was not advised in writing by the Customer of the alleged defect or malfunction of the System within ten (10) days after the expiration of the applicable limited warranty period.
- iv. If a problem develops during the Limited Warranty period, the Customer shall contact Xenogen technical support for assistance.
- v. THE FOREGOING LIMITED WARRANTY IS THE CUSTOMER'S SOLE AND EXCLUSIVE REMEDY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. XENOGEN SHALL NOT BE LIABLE FOR SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOSS OF ANTICIPATED BENEFITS OR PROFITS, LOSS OF SAVINGS

OR REVENUE, PUNITIVE DAMAGES, LOSS OF USE OF THE SYSTEM OR ANY ASSOCIATED EQUIPMENT, COST OF CAPITAL, COST OF ANY SUBSTITUTE EQUIPMENT OR FACILITIES, DOWNTIME, THE CLAIMS OF ANY THIRD PARTIES, INCLUDING CUSTOMERS, AND INJURY TO PROPERTY, RESULTING FROM THE PURCHASE OR USE OF THE SYSTEM OR ARISING FROM BREACH OF THE WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, STRICT TORT, OR ANY OTHER LEGAL OR EQUITABLE THEORY, EVEN IF XENOVEN KNEW OF THE LIKELIHOOD OF SUCH DAMAGES. XENOVEN SHALL NOT BE LIABLE FOR DELAY IN RENDERING SERVICE UNDER THE LIMITED WARRANTY, OR LOSS OF USE DURING THE PERIOD THAT THE SYSTEM IS BEING REPAIRED.

- vi. Some countries, states or provinces do not allow the exclusion or limitation of implied warranties or the limitation of incidental or consequential damages for certain products or the limitation of liability for personal injury, so the above limitations and exclusions may be limited in their application to you. When any implied warranties are not allowed to be excluded in their entirety, they will be limited to the duration of the applicable written warranty. This Limited Warranty gives you specific legal rights which may vary depending on local law.
- vii. This Limited Warranty shall be governed by the laws of the State of California, U.S.A., excluding its conflicts of laws principles and excluding the United Nations Convention on Contracts for the International Sale of Goods.

4.2 Patents

The detection and imaging of light originating within mammals is the subject of several issued patents and pending patent applications in the United States and around the world, including U.S. Patent Numbers 5,650,135, 6,217,847, 6,649,143, 6,890,515, and European Patent Commission Number EP0861093, for which Xenogen Corporation is the exclusive licensor. The use of an IVIS® Imaging System for such applications requires a sublicense from Xenogen Corporation.

In addition, many of the hardware and software components of the Imaging System are the subject of various issued patents and pending patent applications owned by Xenogen, including: United States Patent Number 6,614,452 (Graphical User Interface for In Vivo Imaging) and 6,775,567 (Improved Imaging Apparatus); and United States Patent Applications 09/905668 (Multi-view Imaging Systems), 10/606976 (Method and Apparatus for 3-D Reconstruction of Light Emitting Sources), 10/151463 (Method and Apparatus for Determining Target Depth, Brightness, and Size Within a Body Region), 10/189886 (Fluorescence illumination assembly for an imaging apparatus), and 10/068573 (Light calibration device for use in low level light imaging systems).

4.3 Trademarks

IVIS and Living Image are registered trademarks of Xenoven Corporation. The names of companies and products mentioned herein may be the trademarks of their respective owners. Microsoft and Windows are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries. Pentium III is a registered trademark of Intel Corporation.

[This page intentionally left blank.]

5 IVIS® Lumina Components & Specifications

CCD Camera	16
Imaging Chamber	16
Optics	18
Optics	18
Acquisition Computer	18
Environmental Requirements	19

The IVIS® Lumina consists of a charge coupled device (CCD) camera, an imaging chamber, the XFO-12 fluorescence equipment, and a preconfigured computer. There are no user-serviceable parts in the IVIS Lumina. For more details on the XFO-12 fluorescence equipment, see [Chapter 9, page 33](#).



Figure 5.1 IVIS Lumina, computer, and monitor



VOLTAGE

VOLTAGE! While the IVIS® Lumina can operate at multiple voltages (100 - 240 VAC), users are **not permitted** to change the input voltage to any of the system components. Several internal modifications are required for voltage change when carried out by XENOGEN technical support. If you want to change the operating voltage, contact Xenogen technical support.

IMPORTANT



If you modify the IVIS® Lumina in any way, without prior approval from Xenogen, all warranties that cover this product are void. In addition, the computer included with the IVIS Lumina is specifically configured to run all system-related applications. Any modification of existing software or hardware voids all warranties.

If you have any questions, please contact Xenogen technical support.

5.1 CCD Camera

The camera is a scientific grade, back-thinned, back-illuminated, large format CCD manufactured by Andor Technologies ([Figure 5.1](#)).

CCD Camera Features

- Low dark current
- Thermoelectrically cooled
- 16 bit CCD digitization
- Low-noise electronic readout for extremely low-background images

CCD Camera Specifications

CCD Camera	Specification
Sensor Type	Back illuminated
CCD Format	1024 x 1024
Pixel Dimensions	13 x 13 μm
Quantum Efficiency	>85% 400-700nm >50% 350-900nm 550-650nm
Readout Noise - bin 1	≤ 2 e ⁻ RMS
Dark Charge	<0.0015 e ⁻ /pixel/sec

5.2 Imaging Chamber

The imaging chamber ([Figure 5.2](#)) is a highly specialized device consisting of the imaging chamber housing, a heated, moveable platform, an auto focusing lens system with f/stop control, a filter wheel, and sample illumination LEDs. All adjustable components are motorized and computer-controlled, including the illumination system.

The imaging chamber is *light tight*, so that no light penetrates from the outside. The interior of the imaging chamber is constructed from materials that are non-phosphorescent and non-fluorescent to prevent internal light contamination that could compromise sample measurements.



WARNING

WARNING! Under no circumstances should you attempt to make any mechanical modifications to the imaging chamber.



Figure 5.2 IVIS® Lumina imaging chamber

Imaging Chamber Features

- Custom zero-background imaging chamber
- Eight position optical filter emission wheel and a 12 position excitation filter wheel
- High-efficiency lens assembly
- Sample illumination system
- f/stop control
- Heated and regulated sample shelf temperature to reduce stress on an animal under anesthesia
- Gas anesthesia manifold, including gas delivery and exhaust plumbing
- Software-controlled field of view, f-stop, focus, and optical filter wheel

Imaging Chamber Specifications

IVIS® Lumina Imaging Chamber	Description		
Power requirements	4.0 A max at 120 V	2.0 A max at 240 V	50-60 Hz
Dimensions	19" x 28" x 39"	48 cm x 71 cm x 100 cm	
Door opening dimensions	15" x 20.25"	38 cm x 51 cm	
Weight	90 lbs	41 Kg	

5.3 Optics

Optics	Specification
Lens f/stop	f/.95- f/16
Field of View	5-12 cm square

5.4 Optical Filter Wheel

The imaging chamber has a 8-position computer-controlled optical filter wheel that is located at the top of the imaging chamber, in front of the imaging lens. The system is supplied with four optical emission filters for GFP, DsRed, Cy5.5, and ICG emission. Contact Xenogen technical support for information on installing additional optical filters.

A 12-position excitation filter wheel is attached to the back of the imaging chamber. A set of four excitation and four background subtraction filters that correspond to the standard emission filter set is included.

The filter wheel settings are selected in the Living Image® software. For more details, see the *Living Image Software Manual*.

5.5 Acquisition Computer

The computer contains a Intel family processor and Windows® operating system (Figure 5.1). Microsoft® Office is installed as well as the Xenogen Corporation Living Image® software that controls the IVIS® Lumina.

The computer controls the IVIS Lumina hardware, including the CCD camera. A printer can be connected to the computer.

Computer Features

- Modern, high speed Windows®-based PC
- Microsoft Windows family operating system
- Living Image® software installed. This software controls the IVIS® Lumina, and displays and analyzes image data
- CD-burner installed for data storage and transport
- Network ready
- 20" high-resolution flat screen monitor for image viewing
- Microsoft® Office installed

Computer Specifications

Computer	Description
Power requirements	1.0 A at 120 V 0.5 A at 240 V 50-60 Hz
Dimensions	15.75" x 17" x 4.75" 40 cm x 44 cm x 12 cm

Computer	Description
Weight	22 lbs 10 Kg

Computer Monitor Specifications

Computer Monitor (Flat screen)	Description
Power requirements	0.6 A at 120 V 0.35 A at 240 V 50-60 Hz
Dimensions with stand	17.5" x 17.5" x 9" 45 cm x 45 cm x 23 cm
Weight with stand	33 lbs 15 Kg
Dimensions without stand	17.5" x 17.5" x 2.5" 45 cm x 37 cm x 6.5 cm
Weight without stand	20 lbs 9 Kg

5.6 Environmental Requirements

Environmental Requirements	Specification
Temperature	15° C to 25° C (50° F to 78° F)
Humidity	0% to 80% non-condensing
Type of use	Indoor
Imaging chamber shelf temperature	Ambient to 37° C
Altitude rating	<2000 meters (6560 ft.)
Pollution degree	2
Installation category	II

[This page intentionally left blank.]

6 Operating the IVIS® Lumina

System Start Up	21
Restarting the System After a Power Outage	22
Gas Plumbing	22
Door Operation	24
Imaging Basics	25
System Shut Down Procedure	26

The IVIS® Lumina is intended for use in biophotonic imaging procedures. The system is designed to detect extremely low-level light emissions that are orders of magnitude dimmer than can be seen with the naked eye. The IVIS Lumina allows you to monitor and record cellular and genetic activity within a living organism, in real time. The imaging system captures, quantifies, and images the light emitted by a sample.


6.1 System Start Up

NOTE



All components of the IVIS Lumina should be left on at all times. Periodically rebooting the computer is permissible and does not affect the camera operation.

1. Turn on the computer and monitor.
2. Turn on the IVIS Lumina imaging chamber (the On/Off switch is located on the back of the unit).
3. After the desktop screen is displayed, start the Living Image® software.
4. Enter a User ID (up to three letters) when prompted, then click **Done**.
5. To initialize the system, click **Initialize IVIS System** in the System Control panel (Figure 6.1).



The temperature square is green when the operating temperature is reached. Click the temperature square to display the current

Figure 6.1 IVIS System Control panel in the Living Image software

6. Allow the system to initialize.
 - You will hear the motors move. The System Status box displays the current changes.
- The Temperature square in the IVIS System Control panel is red at startup and turns green when the operating temperature is reached. The IVIS System Control panel displays the current temperature (Figure 6.1).
7. When the temperature is locked at -90°C, as indicated by the green light in the System Control panel, the instrument is ready for operation. (For operating instructions, see the *Living Image® Software Manual*.)

6.2 Restarting the System After a Power Outage

If the IVIS® Lumina experiences a loss of supply power, turn off the power switch on each component. Do not restart the system until reliable line power has been restored.

Turn off all components and wait until stable line power has been restored.

1. Turn on the computer
2. Turn on the imaging chamber.
3. Start the Living Image® software and click **Initialize IVIS System** in the IVIS System Control panel.

6.3 Gas Plumbing

Anesthesia gas tubing is built into the IVIS® Lumina imaging chamber. On the back of the imaging chamber are 0.25" hose barbs that are marked "GAS IN" and "GAS OUT" (Figure 6.2).



Figure 6.2 Imaging chamber, external gas ports

"GAS IN" means the direction of flow is into that port. Similarly, the port labeled "GAS OUT" means that flow can be exhausted out of this port. The valve at the right side of the imaging chamber is used to turn on or off the flow of gas to the animal manifold.



WARNING

WARNING! Use only isoflurane with the IVIS® Lumina. **DO NOT USE FLAMMABLE ANESTHESIA GAS.**



CAUTION

CAUTION! Xenogen recommends using the XGI-8 Gas Anesthesia System when imaging small animals ([Figure 6.3](#)). This system supplies a controlled amount of isoflurane to the imaging chamber and continuously reduces the build-up of isoflurane in the chamber. If you want to use a gas other than the recommended isoflurane/oxygen gas mixture or pure air, contact Xenogen technical support.



Figure 6.3 XGI-8 Gas Anesthesia Delivery System



Be careful to use only tubing and other plumbing fixtures that do not fluoresce or phosphoresce (glow) in the imaging chamber. Contact Xenogen technical support for a list of acceptable materials.

6.4 Door Operation

The IVIS® Lumina imaging chamber door has custom designed hinges, seals, and a four-point closure mechanism. The door is designed to provide a light-tight seal over numerous opening/closing cycles and should close easily without excessive handle turning resistance.

WARNING

WARNING! Do not open the door when the status light is red. If the door is opened when the status light is red, the acquisition is immediately discontinued. This protects the CCD from over exposure.

Stage Curtain

A *stage curtain* is attached to the IVIS imaging chamber platform and covers the empty space beneath the platform (Figure 6.4). The stage curtain serves as a reminder to not place anything below the platform. The curtain attaches to the platform by a rod that is held in place by clips at three locations. If it is necessary to access this area, the curtain can be easily removed and replaced.

To release the rod and curtain from the clips, pull the stage rod slightly forward. Do not allow the curtain to retract around the roller. To reattach the rod, push the rod back into the clips.




Figure 6.4 Stage curtain

6.5 Imaging Basics

Black Paper

Although the platform is black anodized, it is recommended that you image items on a high quality black paper, especially biological specimens. Xenogen Corporation has surveyed many types of paper and recommends Swarthmore, Artagain, Black, part no. 445-109, size 8.5 inch x 11 inch. This paper prevents illumination reflections and helps keep the stage clean.

Centering a Subject

It is recommended that you confirm that the subject is centered on the stage before acquiring an image.

1. Place the subject on the center of the imaging stage and close the imaging chamber door.
2. In the IVIS System Control panel, place a check mark in the **Live** check box.
— The monitor displays an image of the subject.
3. If the subject is centered properly, uncheck the **Live** check box and proceed with imaging.

If the subject is not centered, uncheck the **Live** check box and repeat the step 1 and step 2. For more details, refer to the *Living Image® Software User's Manual*.

Glowing Materials

Always keep in mind that nearly EVERYTHING glows (that is, has the potential to phosphoresce and contaminate the image). Most plastics, almost all tape, plants, paint, rodent food (mostly plants), mouse urine, and animal bedding have been found to glow.

Use caution when introducing materials into the IVIS® Lumina. It is advisable to pre-screen all items by imaging them alone, before imaging them with samples under study. Xenogen recommends using non-powdered gloves when working with IVIS® Lumina equipment.

6.6 System Shut Down Procedure

Xenogen Corporation does not recommend power cycling the IVIS Lumina (turning the system components on and off). If it is necessary to shut down the imaging system for any reason, it is important to follow the procedure below.

1. Close the Living Image® software and save any information of interest at the prompt.
2. Turn off the IVIS® Lumina imaging chamber.
3. Turn off the computer using the standard Windows® shut down procedure.

If you have any problems during the shut down or start up procedure, please contact Xenogen technical support for assistance.

7 Care & Maintenance

7.1 Cleaning the IVIS Lumina

Approved Cleaning Solutions

The compounds shown in [Table 7.1](#) do not damage the internal finish of the IVIS® Lumina imaging chamber and are suitable as cleaners, if required. Do not use any solution not included in this list. In particular, avoid strong bases, bleach, or acids that may potentially damage the unit and compromise its operation.

IMPORTANT



Do not spray cleaning solutions in the imaging chamber.

Table 7.1 Acceptable cleaning solutions for the IVIS Lumina imaging chamber

Cleaning Solution	Manufacturer
Cidexplus® Solution (3.4% glutaraldehyde)	Johnson & Johnson Medical
70% methyl alcohol/30% deionized water solution	
70% ethyl alcohol/30% deionized water solution	
Sporicidin® Sterilizing Solution (1.56% phenol)	Sporicidin International
Clidox-s® Disinfectant	Pharmaca Research Laboratories, Inc.

NOTE



Xenogen makes no claims as to the sterility of the IVIS Lumina imaging chamber after using the solutions in [Table 7.1](#). Please refer to the manufacturer's literature for information as to the applicability of the compound for the organism of interest.

It is recommended that you use a lint-free wipe, such as Scott Pure® wipe or a Kaydry EX-L® wipe to minimize the presence of particulate matter in the imaging chamber.

After saturating a lint-free wipe, clean the internal surfaces using a gentle circular motion. Do not pour or spray the solution directly onto internal surfaces. Rinse surfaces using a wipe saturated with sterile deionized water. Do not allow puddles of water to remain on the surfaces. To avoid any phosphorescence from the cleaner, be sure that

the surfaces are dry before using the imaging chamber. Be careful not to smudge the camera lens and optical filters.

Consider dedicating an IVIS® Lumina for immunodeficient animals to remove the risk of cross contamination.

8 Troubleshooting

Measured Temperature Does Not Equal the Demand Temperature . . .	29
Photographic Image Is Unacceptable	30
Luminescent Image Is Unacceptable	31
No Image Is Produced	32

8.1 Measured Temperature Does Not Equal the Demand Temperature

At start up, the Living Image® software programs the CCD camera to maintain the CCD at -90° C. If the camera power supply remains on (IVIS® Lumina box), the system maintains this temperature regardless of whether the Living Image software is open or the computer is turned on.

- To check the temperature of the CCD, click the Temperature square (red or green) in the IVIS System Control panel of the Living Image software (Figure 8-1).




Figure 8-1 IVIS System Control panel in the Living Image software

Problem	Possible Cause	Corrective Action
Measured temperature is warmer than the demand temperature.	Ambient temperature may be too high, camera air vents may be blocked, or system needs service. A problem may exist with the camera.	Verify that the room temperature is within operational limits. Check air vents in the camera head by removing protective cover and inspecting. Contact Xenogen technical support for assistance.


Problem	Corrective Action
Something is obstructing the stage.	<ol style="list-style-type: none"> 1. Open the door to the imaging chamber and visually inspect the stage. 2. Remove anything that is physically obstructing the stage. If there is no obstruction and/or the stage still does not move, contact Xenogen technical support for assistance.

8.2 Photographic Image Is Unacceptable

Photograph imaging parameters are automatically controlled and generally produce a good quality photo. If shiny objects are imaged, creating specular reflections, the automatic algorithm may get confused and produce an underexposed image.

Also, refer to the *Living Image Software Manual* for further details on acquiring images.

Problem	Possible Cause	Corrective Action
Image is streaked.	Subject moved during the exposure.	Check to see if the subject may have moved. If the subject is not on the sample stage, it is probably on the floor of the imaging chamber. If the sample has moved, locate and reanesthetize it. If gas anesthesia is being used, confirm that the anesthesia is turned on and the flow rate is appropriate.
Image is blurry.	Subject height is significantly less than or greater than 1.5 cm.	The focus is set for a sample height of 1.5 cm. Significant deviation from this height results in an out-of-focus photograph.
	Incorrect f/stop setting.	The f/stop for photographs should be set to f/8 or f/16. An f/stop smaller than 8 reduces the depth of field in the photograph.
A white spot appears in the center of the field of view.	An excessively moist environment in the imaging chamber can result in condensation on the CCD window (Figure 8-2).	Turn off the entire system and remove excess moisture in the imaging chamber. Allow the chamber to thoroughly dry. If the problem persists, contact Xenogen technical support for assistance.



Click #: DW20010705122916
Acq Date: Thu, Jul 05, 2001
Acq Time: 12:29:29
Bin: M (2), FOV: 15, 6
Camera: IVIS 1356, LN1300EB
Series: Plate Test
Experiment:
Label:
Comment:
Analysis Comment:

Figure 8-2 Condensation example

8.3 Luminescent Image Is Unacceptable

Binning, f/stop, and exposure time affect the appearance of a luminescent image. Please refer to the *Living Image® Software Manual* for instructions on setting binning, exposure time, and f/stop values.

In order to function properly and reduce camera noise, the CCD camera must be cooled to the demand temperature before acquiring an image. If the camera is not cooled to the demand temperature, imaging may result in false positive signals.

Problem	Corrective Action
Light contamination- Internal	Check to see that there are no extraneous light sources in the imaging chamber. Many substances phosphoresce when exposed to light. Be especially cautious of plastics and substances that contain pigment. Be sure to pre-screen any substance or material before performing actual experiments.
Light contamination- External	<p>A 2" diameter high reflectance hemisphere (Xenogen part no. 30098) is used to help check for light leaks (Figure 8-3). Contact Xenogen technical support to purchase this accessory.</p> <p>To check for light leaks:</p> <ol style="list-style-type: none"> 1. Place the high reflectance hemisphere in the imaging chamber on the stage using a subject height of 3.5 cm at field of view D. 2. Take a luminescent image of the hemisphere using the luminescent settings: f/stop = 1, Binning = Large (high sensitivity), and exposure time = 5 minutes. <p>If the hemisphere can be easily seen, there is a light leak. Contact Xenogen technical support for assistance.</p>
Camera noise	<p>Verify that the camera is cooled to the demand temperature.</p> <ol style="list-style-type: none"> 1. Check the measured temperature in the IVIS System Control panel to ensure that it is locked. If the camera temperature is locked, the camera temperature box is green. 2. If the camera temperature box is red, click the red box to display the actual temperature. See Measured temperature is warmer than the demand temperature., page 29.



Figure 8-3 High reflectance hemisphere

Wear powderless gloves when handling the hemisphere to prevent surface contamination that may glow when using the IVIS[®] Lumina.

8.4 No Image Is Produced

If no image is produced, there may be an error in the Living Image software, a problem with the physical connections to the camera, or a hardware failure.

1. Close the Living Image[®] software and restart the computer.
2. Restart the Living Image software and try to acquire an image.
3. If after restarting the computer, you are still unable to produce an image, contact Xenogen technical support for assistance.

9 XFO-12 Fluorescence Equipment

Installation Requirements	34
Specifications	34
Description and Theory of Operation	35
Acquiring Fluorescent Images	40
Troubleshooting	42
Care and Maintenance of the XFO-12 Fluorescence Equipment	46
XFO-12 Hardware Replacement Parts List	47

The XFO-12 fluorescence equipment provides the IVIS® Lumina with fluorescent imaging capability. The XFO-12 equipment can be used for *in vitro* or *in vivo* applications. The sensitive range of the IVIS Lumina CCD camera sets the wavelength range for fluorescence applications, which is approximately 400-950 nm (30-85% quantum efficiency). As with bioluminescent imaging, wavelengths greater than 600 nm are preferred for *in vivo* fluorescent applications due to lower absorption in tissue. The Living Image software controls fluorescent image acquisition, including lamp power, level, and filter selection.

This chapter explains how to operate the XFO-12 fluorescence equipment with the IVIS Lumina. It also provides important safety and maintenance information.

IMPORTANT



To ensure optimum and safe performance of the XFO-12 fluorescence equipment with maximum service life, read this chapter carefully before you use IVIS Lumina with the fluorescence equipment.

You should also be familiar with and refer to the other chapters of this manual for the IVIS Lumina.

For instructions on how to use the Living Image software that controls fluorescent image acquisition, see the *Living Image® Software Manual* that is provided with the software.

The *Schott Fostec DCR® III Direct Current Regulated Light Source User's Manual and Technical Reference* is a separate manual for the light source module. It provides additional useful information on the light source module and its safe operation.

9.1 Installation Requirements

The XFO-12 fluorescence equipment requires **100 - 240 VAC** electrical power. The system automatically accepts the required voltage.

IMPORTANT



The XFO-12 fluorescence equipment operates at the same voltage as the IVIS® Lumina chamber and must not be used at other than its labeled voltage.

For more details on the XFO-12 operating requirements, see [Specifications, page 34](#).

9.2 Specifications

Electrical Power and Fuses

Voltage Available	100VAC, 50/60 Hz - 240VAC, 50/60Hz
Power Consumption	200 W

NOTE



For more on additional electrical power requirements of the IVIS Lumina, see [Chapter 5, page 15](#).

Environmental - Temperature and Humidity

Temperature	15-25° C
Humidity	0% to 80% Non-Condensing
Type of Use	Indoor

Lamp and Fuse

Lamp	150 W Quartz Tungsten Halogen 21 V EKE lamp life: 400 hrs ("High" setting), 6,000 hrs ("Low" setting) Color temperature 3250 K°
Fuse	3A(10) 250VAC

Ventilation Requirements

Provide sufficient space (minimum 6 inches) behind the fan of the light source module so that airflow is unobstructed. Provide a similar minimum distance behind the IVIS Lumina to enable fan cooling as well as adequate room for fiber optic cable routing.

Chemicals Required for Operation

No chemicals are required for the operation of the IVIS Lumina or the XFO-12 fluorescence equipment. Other user supplied chemicals or materials may be required for your specific biological testing procedures.

Weight and Dimensions of the XFO-12 Fluorescence Light Source Module

Weight	4.94 lbs (2.24 kg)
Depth	8.61 inches (21.9 cm)
Width	7.27 inches (18.5 cm)
Height	4.6 inches (17.7 cm)

9.3 Description and Theory of Operation

System Components

The Xenogen XFO-12 fluorescence equipment provides fluorescent imaging capability. You can conveniently switch between bioluminescent and fluorescent imaging applications.

Figure 9.1 shows the XFO-12 equipment and Figure 9.2 depicts the interior of the IVIS Lumina imaging chamber (green indicates the excitation light).




Figure 9.1 IVIS Lumina and XFO-12 fluorescence equipment




Figure 9.2 IVIS Lumina and XFO-12 fluorescence equipment.

The XFO-12 fluorescence light source module (Figure 9.3) provides the fluorescence excitation light. This light source consists of a 150-Watt quartz tungsten halogen lamp with a dichroic reflector. Figure 9.4 shows the relative spectral radiance output of the lamp/reflector combination and indicates emission throughout the IVIS Lumina wavelength range of 400-950 nm.




Figure 9.3 XFO-12 fluorescence light source module, front view (left) and rear view (right).




Figure 9.4 Relative spectral radiance for the quartz halogen lamp

The dichroic reflector reduces infrared transmission (>700 nm) enough to prevent overheating of fiber optic bundles, while still allowing sufficient infrared light throughput to permit imaging at these wavelengths. The lamp intensity level is computer controlled by the Living Image® software. The user can adjust the lamp output intensity by means of software to settings of low or high. The XFO-12 fluorescence light source module operates under software control, therefore manual adjustment of the front panel lamp potentiometer is disabled.

The lamp output is delivered to the excitation filter wheel assembly (Figure 9.5) on the back of the IVIS® Lumina imaging chamber using a 1/4 inch diameter fiber optic bundle. Figure 9.6 shows a cross section of this excitation filter wheel assembly. Light from the input fiber optic bundle passes through a collimating lens followed by a 25 mm diameter excitation filter. Twelve filter wheel locations allow you to choose up to eleven excitation filters. A light block is provided in one filter slot and is used during bioluminescent imaging to prevent external light from entering the IVIS Lumina imaging chamber. The twelve-position excitation filter wheel is motor-controlled through the Living Image software.




Figure 9.5 Excitation filter wheel assembly




Figure 9.6 The cross-sectional area of the excitation filter wheel.

Following the excitation filter, a second lens then focuses light into a one quarter inch fused silica fiber optic bundle inside the IVIS® Lumina imaging chamber. Fused silica (core and clad) fibers are used in this bundle to avoid the generation of autofluorescence in the fiber, as is the case with ordinary glass fibers.

The fused silica fiber bundle splits into four separate legs that deliver filtered light to four reflectors located on the ceiling of the IVIS Lumina imaging chamber (Figure 9.2, page 35). Typical illumination profiles for stage locations A-D (fields of view 5-12 cm respectively) are shown in Figure 9.7. Note that the profiles for all the stage locations are peaked

near their center. The non-uniformity of the illumination pattern is compensated for in Living Image when units of efficiency are selected (for more details, see the *Living Image® Software User's Manual*). When imaging 96-well plates, the lower stage positions (C and D) are preferred to minimize shadowing effects due to the off-axis illumination.

Fluorescent emission from the target fluorophore is collected through an emission filter wheel located at the top of the IVIS® Lumina imaging chamber and then focused into the CCD Camera. The emission filter wheel contains eight openings. Users have the ability to choose up to seven emission filters (60 mm diameter), leaving one position open for bioluminescent imaging.




Figure 9.7 Typical illumination profiles for stage locations (FOV of 5, 7, 10 and 12 cm, respectively) measured from the FOV center.

Understanding Filter Spectra

The use of high quality filters is essential for obtaining good signal-to-background levels (contrast) in fluorescence measurements, particularly in a high sensitivity instrument such as the IVIS® Lumina. Typical excitation and emission fluorophore spectra are shown in [Figure 9.8](#), along with idealized excitation and emission filter transmission curves shown as rectangles. The excitation and emission filters are called "bandpass" filters; ideally they transmit all the wavelengths within the bandpass region and block (absorb or reflect) all wavelengths outside the bandpass. This spectral band is like a window, characterized by its central wavelength and its width at 50% peak transmission (full width half maximum, FWHM).




Figure 9.8 Typical excitation and emission spectra for a fluorescent compound. Included in the picture are two idealized bandpass filters that would be used with this fluorescent compound.

Real filters are not perfect and have some leakage outside of the bandpass region and can also exhibit autofluorescence, depending on the materials used in the filter construction. More realistic filter transmission curves are shown in [Figure 9.9](#). Referring to this figure, we can define several terms: "bandgap" is the spacing between the transmission regions of the excitation and emission filters; "transmission" is the fraction of light that passes through the filter bandpass region; "blockage" is the light rejection in the non-transmitting region (or the region outside the band) of the filter spectrum; "leakage" is undesirable light that is not blocked properly by the filter and is detected by the camera.




Figure 9.9 Illustration of typical excitation and emission curves.

The vertical axis in [Figure 9.9](#) is *optical density*, defined as $OD = -\log(T)$ where T is the transmission. An optical density of 0 indicates 100% transmission, whereas OD7 indicates a reduction of the transmission to 10^{-7} . Typical transmission of an XFO-12 filter in the bandpass region is about 0.7 (OD0.15) and typical blocking outside of the bandpass region is about OD7. The band gap between the two filters is usually defined as the gap at 50% transmission (OD0.3). There is a slope in the

transition region from bandpass to blocking, as indicated in [Figure 9.9](#). A steep slope is required to avoid overlap between the two filters.

The XFO-12 fluorescence filters are high quality interference filters, constructed from alternating layers of dielectric films on a substrate of low autofluorescent glass. Care has been taken to minimize autofluorescence of the filter so that its level is below OD7. Filter passbands for the standard set are listed in [Table 9.1](#).

Table 9.1 Standard filter sets, fluorescent dyes, and proteins used with the IVIS® Lumina.

Filter Set	Label	Background Passband	Excitation Passband	Emission Passband	Fluorescent Dyes & Proteins
1	GFP	410-440	445-490	515-575	GFP, EGFP, FITC
2	DsRed	460-490	500-550	575-650	DsRed2-1, PKH26, CellTracker™ Orange, 605 Qdot® Bioconjugate
3	Cy5.5	580-610	615-665	695-770	Cy5.5, Alexa Fluor® 660, Alexa Fluor® 680, 705 Qdot® Bioconjugate
4	ICG	665-695	710-760	810-875	Indocyanine green (ICG), 805 Qdot Bioconjugate

9.4 Acquiring Fluorescent Images

NOTE



This chapter provides a quick start guide to acquiring fluorescent images. See the *Living Image® Software User Manual* for complete details.

IVIS® System Control Panel Overview

Acquiring fluorescent images using Living Image software is controlled through the IVIS System Control Panel ([Figure 9.10](#)). To acquire a fluorescent image, the user must check the **Fluorescent** box on the left side of the panel. Once selected, controls for the illumination lamp – **Fluor Lamp Level**, and **Filter Lock** – will appear in the top half of the panel. Checking the **Filter Lock** box ensures that the excitation and emission filters are properly paired. During image acquisition, the lamp is computer-controlled through Living Image software. The **Fluorescent level** pull-down menu controls the illumination intensity level of the lamp with options – **Off**, **Low**, **High**, and **Inspect**. The **Low** setting is approximately 18% of the **High** setting. **Inspect** turns on the illumination lamp, allowing the user to manually inspect the excitation lamp.

NOTE



Make sure the desired filters are currently selected in the filter popup menus before selecting Inspect. The Inspect operation will automatically position the filters currently selected in the filter popup menus before turning on the lamp. Subsequent changes to the filter popup menus will have no effect until another Inspect operation is done.




Figure 9.10 IVIS System Control panel.

Procedure for Acquiring Fluorescent Images

1. If it is not already on, start up the acquisition computer and start the Living Image® software.
— The System Control panel appears ((Figure 9.10)).
2. Click **Initialize IVIS system**.
— After initialization, the **Temperature** box in the center of the panel should be green, indicating that the CCD camera is adequately cooled. (Allow 10 to 15 minutes for the camera to reach the proper temperature.) The **Temperature** box changes from red to green when the CCD camera has reached the proper operating temperature.
3. Place the sample to be imaged in the center of the stage in the imaging chamber and close the door.
4. Select the **Fluorescent** check box.
5. Make a selection from the **Field of View** drop-down menu on the left side of the Control Panel.
6. Enter the approximate (0.5 cm) **Subject Height** (height) in the lower left entry box (or focus manually).
7. Select the **Emission Filter** and **Excitation Filter**. If the **Filter Lock** box is checked, select the excitation or emission filter of interest. Select only one, as the other filter will be selected automatically.
8. Select **High** or **Low** from the Fluor Lamp Level drop-down menu. **High** is the recommended setting
9. Set the **Exposure Time**, **Binning**, and **f/stop**.

Fluorescence is generally brighter than bioluminescence, so the exposure time is shorter and f/stop higher (smaller lens opening). Typical fluorescent image camera settings might be 10 sec exposure time, Binning=small, and f/2.

10. Click **Acquire**.

- After the exposure is completed, the overlay image is displayed.

NOTE



During the fluorescent image acquisition, the **Acquire** button becomes a **Stop** button, which can be used to terminate the exposure, if necessary.

11. To save the data, select **Living Image → Save Living Image Data** from the main menu bar.

This completes the data acquisition. To obtain additional images, repeat the process, beginning with step 3.

The Image image window that displays the fluorescent image includes annotations specific to fluorescence (including emission filter, excitation filter, and fluorescence level) as well as standard annotations such as exposure time, f/stop, FOV, and date/time of exposure.

9.5 Troubleshooting

Hardware Problems

If you have difficulty during fluorescent imaging, the trouble may be due to the lack of excitation light. Loss of excitation light can result from either a burned out quartz tungsten halogen lamp or a blown line fuse. The following procedure describes a troubleshooting process for determining the problem. [Figure 9.3, page 36](#) shows the XFO-12 fluorescence light source module.)

1. Verify that the fiber optic cables are not loosened or disconnected from their proper connectors.
2. Adjust "Fluor Light Level" (excitation lamp) to a value of high in the IVIS System Control Panel.
3. Check the "Change Lamp When Red" indicator on the front of the light source module. If indicator is red, proceed to [Lamp Replacement, page 44](#).
4. Take a fluorescent image and check for light and fan operation by looking through the fan guard on the rear panel of the light source module. If there is neither, a blown fuse is the probable cause for light loss. See [Fuse Replacement](#) below.
5. Observe the operation of the system by selecting "Inspect" in the Fluorescent Lamp Level drop down list. This causes the selected filter to move into place and the lamp to turn on. Open the chamber and visually inspect to see if the excitation light is incident on the sample stage.

**WARNING**

WARNING! DO NOT disconnect or reattach the electrical control cable that connects the XFO-12 fluorescence equipment (excitation filter assembly) to the light source module and IVIS® Lumina (electronics tray) when the power is on. See [System Components, page 35, Figure 9.1](#) and [Figure 9.2](#) for photographs of these components. Disconnecting or reconnecting the control cable when the system has electrical power will damage the system. Always turn off the front panel light source module switch and the rear-mounted ON/OFF switch on the IVIS Lumina before making or breaking any of these cable connections.

**WARNING****Fuse Replacement**

WARNING! The following procedure can expose the user to hazardous voltages unless the electrical power to the XFO-12 fluorescence light source module is completely eliminated. As instructed in the procedure, turn off the lamp from the module front panel ON/OFF switch and remove the electrical power to the module by disconnecting the plug from the surge protector and the back of the module.

1. See the *Fuse Replacement* section in the separate manual for the light source module (*SCHOTT Fostec DCR® III Direct Current Regulated Light Source User's Manual and Technical Reference*).
2. Press the ON(1)/OFF(0) switch to the OFF(0) position on the light source module.
3. Remove AC line cord from the AC power receptacle on the surge protector.
4. Remove the AC line cord from the power entry module on the rear panel of the light source module.
5. Pull out the fuse drawer. The fuse drawer is part of the power entry module, and is located directly beneath where the AC line cord plugs in. Use a thin bladed screwdriver or penknife blade if necessary. Use care not to damage the drawer. Remove the blown fuse that is positioned closest to the light source module. The second fuse in front of the holder is a spare. Discard the blown fuse.
6. Examine the spare fuse and verify it is the correct rating. See [Electrical Power and Fuses, page 34](#).
7. Place the replacement fuse into the fuse drawer. The fuse will work in either direction. Fill the spares slot in the drawer holder if another good fuse of correct rating is available. Fuse value: 3.0A, 5 x 20, 250V (Xenogen part no. 20395).
8. Push the fuse drawer until it clicks into position.
9. Reattach the AC line cord to the power entry module, then to the AC receptacle on the surge protector.
10. Return the ON/OFF switch to the ON(1) position (the red indicator will be visible on the switch).

NOTE



For the software to control lamp functions, the front panel switch on the light source module must be in the ON position.

11. Resume normal operation.

Lamp Replacement



WARNING

WARNING! The following procedure can expose the user to hazardous voltages unless the electrical power to the light source module of the XFO-12 fluorescence equipment is completely eliminated. As instructed in the procedure, turn off the lamp from the module front panel ON/OFF switch, and remove the electrical power to the module by disconnecting the plug from the surge protector and the back of the module.



WARNING

WARNING! Before you replace the tungsten halogen lamp or lamp assembly in the XFO-12 fluorescence light source module, verify that the lamp and socket are cool before proceeding.

1. Review the sections *Lamp Module Replacement* and *Lamp Replacement* in the separate manual for the XFO-12 fluorescence light source module (*SCHOTT Fostec DCR® III Direct Current Regulated Light Source User's Manual and Technical Reference*). It contains additional useful information and warnings regarding the light source.
2. In the IVIS System Control panel of the Living Image software, set "Fluorescence On" and adjust "Fluor lamp Level" to OFF ([Figure 9.11](#)). Allow the fan on the light source module to run for several minutes until the fiber optic connector (nosepiece with thumb screw) is cool to touch.




Figure 9.11 IVIS System Control panel

3. Press the ON (1)/OFF(0) switch to the OFF(0) position on the light source module.
4. Remove AC line cord from the AC power receptacle on the surge protector and from the back of the module.
5. Release the lamp assembly from the front panel of the light source module by turning the retaining screw counterclockwise using a straight blade screwdriver. The retaining screw will be disengaged from the front panel, but will remain in place in the lamp assembly.
6. Remove the lamp assembly from the light source module by slowly pulling outward.
7. **Check the lamp assembly to verify that the lamp and socket are cool before proceeding.**
8. Push back on the lampholder release lever to raise the lamp (L) (Figure 9.12).
9. Gently unplug the ceramic lamp socket (K).




Figure 9.12 Light source module

10. Lift the release lever. Insert a new lamp into the lamp holder. Push the lamp all the way down. Be careful not to touch the bulb or the inside of the reflector. Fingerprints may affect the light output.
11. Push the ceramic socket onto the two lamp pins until flush with the lamp base.
12. Return the “Modulamp” assembly into the recess and tighten the captive thumbscrew.
13. For complete “Modulamp” assembly removal and installation, see sections 5.1 and 5.2 in the *SCHOTT Fostec DCR® III Direct Current Regulated Light Source User’s Manual and Technical Reference*.
14. Reattach the AC line cord to the back of the module and then to the AC receptacle on the surge protector.

15. Return the ON/OFF switch to the ON(1) position (the red indicator will be visible on the switch).

NOTE



For the software to control lamp functions, the front panel switch located on the light source module must be in the “On” position.

16. Resume normal operation.

9.6 Care and Maintenance of the XFO-12 Fluorescence Equipment

Cleaning the XFO-12 Fluorescence Light Source Module

If necessary, wipe the exterior surfaces of the XFO-12 light source module with a soft cloth.



WARNING

WARNING! DO NOT use fluids to clean the exterior (or interior) of the module. Do not allow fluids of any kind to enter the light source module under any circumstances. Sprays and liquids that come into contact with the light source module or the IVIS® Lumina may result in damage to the system or electrocution.

If the XFO-12 fluorescence light source module requires more aggressive cleaning or sterilization, contact Xenogen technical support.

Cleaning the IVIS Lumina and Other Components of the XFO-12 Fluorescence Equipment



WARNING

WARNING! DO NOT use fluids or moistened towels to clean the any part of the IVIS Lumina where electrical or fiber optic cables make connections. Do not use fluids of any kind in the vicinity of the XFO-12 Excitation Filter Wheel Assembly (mounted on the rear of the imaging chamber). Turn off electrical power to the IVIS Lumina before engaging in cleaning operations using fluids. The Imaging Chamber power switch is located in the rear on the electronics tray.

For more details on how to clean the imaging chamber, see [Cleaning the IVIS Lumina, page 27](#).

Cleaning the Optical Components and Filter Replacement

Contact Xenogen technical support for information about cleaning or sterilizing any of the optical components or the optical filter replacement.

9.7 XFO-12 Hardware Replacement Parts List

Item	Xenogen Part No.
Quartz tungsten halogen lamp for XFO-12 light source module	20396
Fuse, 3.0A, 5x20, 250V	20395
Low fluorescent imaging mat	60104

[This page intentionally left blank.]

Appendix A: Moving IVIS® Lumina With the XWS-260 Workstation

This appendix explains how to move the IVIS® Lumina that is configured with the XWS-260 workstation (Figure A.1).

CAUTION

CAUTION! Xenogen recommends that you do not move an IVIS Lumina that is not located on the XWS-260 workstation. If you need to move an imaging system that is not on the workstation or cart, contact Xenogen technical support for assistance.

The procedure involves shutting down the system components and unplugging them from the line power, moving the workstation with the components, and restarting the system components after the power is restored.

CAUTION

CAUTION! The IVIS Lumina includes many cables and lines. To avoid damaging the system components, it is very important to closely follow all directions.



Figure A.1 IVIS Lumina configured with the XWS-260 Workstation

Shutting Down the Imaging System

1. Close the Living Image® software after you save any important information.
2. Turn off the computer using the standard shut down procedure.
3. Turn off the imaging chamber.

Moving the Imaging System on the XWS-260 Workstation

1. Unlock the wheels on the XWS-260 workstation and carefully roll the workstation to the new location.

CAUTION

CAUTION! When moving the imaging system on the XWS-260 workstation, be sure to grasp and push the workstation from as low as comfortably possible. Xenogen advises having two people present during any major move to minimize the risk of damage to the imaging system. If the IVIS® Lumina is damaged during movement, Xenogen cannot be held responsible.

2. When you are finished moving the workstation, lock the wheels.
If you moved the imaging system to a new location, verify that the power outlets meet the power requirements for the system. For more details on equipment power requirements, see [Chapter 5, page 15](#).

Starting the Imaging System

1. Turn on the computer and monitor.
2. Turn on the imaging chamber (the power switch is on the back of the unit).
3. After the desktop screen appears on the monitor, start the Living Image® software.
4. When prompted, enter a user ID (up to three letters) and click **Done**.
5. To initialize the system, click **Initialize IVIS System** in the IVIS System Control panel ([Figure A.2](#)).
— As the system initializes, you will hear the motors moving. The System Status box displays the current changes.




Figure A.2 IVIS System Control panel in the Living Image software

6. When the measured temperature reaches the demand temperature, the imaging system is ready to use. This requires about 10 minutes.

Index

A

- acquiring fluorescent images 41–42
- anesthesia delivery system 23

B

- bandgap 39
- bandpass filter 38
- black paper 25
- blockage 39

C

- camera noise 31
- caution definition 3
- CCD camera 16
- centering a subject 25
- chemical & biological safety 9
- cleaning
 - imaging system 5
- cleaning solutions 27
- cleaning the system 27–28
- computer 18
 - specifications 18
- computer monitor
 - specifications 19
- cover 9

D

- definitions
 - caution 3
 - voltage 3
 - warning 3
- door operation 24

E

- electrical power 34
- electrical safety 7
- emission bandpass 40
- emission filter transmission curve 38
- emission filter wheel 38
- emission fluorophore spectra 38
- environmental considerations
 - heat 4

- water and moisture 5
- environmental requirements 19
- excitation filter 40
- excitation filter transmission curve 38
- excitation filter wheel 37
- excitation fluorophore spectra 38
- eye and burn hazard 8

F

- filter lock 40
- filter replacement 46
- filter spectra 38
- fluorescent dyes 40
- fluorescent images
 - acquiring 41–42
- fluorescent proteins 40
- fuse 34, 47
- fuse replacement 43–44

G

- gas plumbing 22
- glowing materials 25

H

- hardware problems 42
- heat safety 4
- humidity 34

I

- imaging chamber 16–17
 - door 24
 - gas plumbing 22
 - specifications 17
 - stage curtain 24
- imaging system
 - cleaning 5
 - servicing 6
- installation requirements 34
- isoflurane gas anesthesia 23
- IVIS System Control panel 40

L

- lamp 34
- lamp replacement 44–46
- leakage 39
- light block 37
- Light Source Module 36
 - dimensions 35
 - weight 35
- lightning 6
- limited warranty 11–12
- Living Image software 37
- low fluorescent imaging mat 47
- luminescent image unacceptable 31

M

- manual information 4
- mechanical safety 8
- modules 9
- moving the system 49–50

N

- no image produced 32

O

- optical density 39
- optical filter wheel 18
- other equipment 6

P

- panels 9
- patents 12
- photographic image unacceptable 30
- power cord 6
- power line surges 6
- power outage
 - restarting the system 22
- power outages 6
- power overloading 6
- power sources 5

R

- replace
 - filter 46
 - fuse 43–44

- lamp 44–46
- replacement parts
 - fuse 47
 - low fluorescent imaging mat 47
 - quartz tungsten halogen lamp 47
- requirements
 - installation 34
 - voltage 34
- restarting after power outage 22

S

- safety
 - chemical & biological 9
 - electrical 7
 - electrical safety 7
 - eye and burn hazard 8
 - mechanical 8
 - panels, cover, and modules 9
- safety instructions 3–6
- servicing 6
- servicing the equipment 7
- shut down procedure 26
- specifications
 - computer 18
 - computer monitor 19
 - electrical power 34
 - environmental 34
 - fuse 34
 - humidity 34
 - imaging chamber 17
 - lamp 34
 - temperature 34
- stage curtain 24
- start up 21–22
- startup
 - restarting after power outage 22
- subject placement 25
- system
 - cleaning 27–28
 - environmental requirements 19
 - moving with XWS-260 workstation 49–50
 - restarting after power outage 22
 - shut down 26
 - start up 21–22
- system components
 - XFO-12 35

T

- technical support 2
- temperature 34
 - measured unequal to demand 29
- trademarks 13
- transmission 39
- troubleshooting
 - camera noise 31
 - measured and demand unequal 29
 - no image produced 32
 - unacceptable luminescent image 31
 - unacceptable photographic image 30

V

- ventilation requirements 34
- voltage definition 3
- voltage requirements 34

W

- warning definition 3
- water and moisture safety 5

X

- Xenogen technical support 2
- XFO-12 fluorescence equipment
 - installation requirements 34
 - system components 35
- XWS-260 workstation 49–50

[This page intentionally left blank.]