



50 (cat. no. 762164)

PAXgene®

Blood RNA Kit Handbook

Version 2

IVD

For in vitro diagnostic use

The PAXgene Blood RNA System consists of a blood collection tube (PAXgene Blood RNA Tube) and nucleic acid purification kit (PAXgene Blood RNA Kit). It is intended for the collection, storage, and transport of blood and stabilization of intracellular RNA in a closed tube and subsequent isolation and purification of intracellular RNA from whole blood for RT-PCR used in molecular diagnostic testing.

Performance characteristics for the PAXgene Blood RNA System have only been established with FOS and IL1B gene transcripts. The user is responsible for establishing appropriate PAXgene Blood RNA System performance characteristics for other target transcripts.

REF

762164



PreAnalytiX GmbH,
Feldbachstrasse, CH-8634 Hombrechtikon
Produced by QIAGEN GmbH for PreAnalytiX

March 2009

 PreAnalytiX
A QIAGEN / BD Company

Trademarks: PAXgene®, PreAnalytiX® (PreAnalytiX GmbH); QIAGEN®, QIAcube® (QIAGEN Group); BD Vacutainer®, BD Hemogard™, Safety-Lok™ (Becton, Dickinson and Company, Franklin Lakes, NJ, USA); Eppendorf® (Eppendorf-Netheler-Hinz GmbH).

Limited License Agreement

Use of this product signifies the agreement of any purchaser or user of the PAXgene Blood RNA Kit to the following terms:

1. The PAXgene Blood RNA Kit may be used solely in accordance with the *PAXgene Blood RNA Kit Handbook* and for use with components contained in the Kit only. PreAnalytiX grants no license under any of its intellectual property to use or incorporate the enclosed components of this Kit with any components not included within this Kit except as described in the *PAXgene Blood RNA Kit Handbook* and additional protocols available at www.preanalytix.com.
2. Other than expressly stated licenses, PreAnalytiX makes no warranty that this Kit and/or its use(s) do not infringe the rights of third-parties.
3. This Kit and its components are licensed for one-time use and may not be reused, refurbished, or resold.
4. PreAnalytiX specifically disclaims any other licenses, expressed or implied other than those expressly stated.
5. The purchaser and user of the Kit agree not to take or permit anyone else to take any steps that could lead to or facilitate any acts prohibited above. PreAnalytiX may enforce the prohibitions of this Limited License Agreement in any Court, and shall recover all its investigative and Court costs, including attorney fees, in any action to enforce this Limited License Agreement or any of its intellectual property rights relating to the Kit and/or its components.

For updated license terms, see www.preanalytix.com.

© 2005–2009 PreAnalytiX GmbH, all rights reserved.

Contents

Explanation of Symbols	4
Kit Contents	5
Storage Conditions	5
Intended Use	6
Product Use Limitations	6
Quality Control	6
Technical Assistance	6
Safety Information	7
Introduction	9
Principle and procedure	9
Sample collection and stabilization	9
RNA concentration and purification	12
Manual RNA purification	14
Automated RNA purification	22
Equipment and Reagents to Be Supplied by User	26
Important Notes	27
Using the QIAcube	27
Starting the QIAcube	27
Installing protocols on the QIAcube	27
Loading the QIAcube	29
Protocol	
■ Manual Purification of Total RNA from Human Whole Blood Collected into PAXgene Blood RNA Tubes	35
■ Automated Purification of Total RNA from Human Whole Blood Collected into PAXgene Blood RNA Tubes	40
Troubleshooting Guide	45
Appendix A: General Remarks on Handling RNA	47
Appendix B: Quantification and Determination of Quality of Total RNA	48
Appendix C: Handling PAXgene Blood RNA Tubes	50
Ordering Information	51

Explanation of Symbols



Contains sufficient for <N> tests



Consult instructions for use



Use by



In vitro diagnostic medical device



Catalog number



Batch code



Temperature limitation



Upper limit of temperature



Manufacturer

Kit Contents

PAXgene Blood RNA Kit	(50)
Catalog no.	762164
Number of preps	50
Buffer BR1 (Resuspension Buffer)	20 ml
Buffer BR2 (Binding Buffer) *	18 ml
Buffer BR3 (Wash Buffer) *	45 ml
Buffer BR4 (Wash Buffer) †	11 ml
Buffer BR5 (Elution Buffer)	6 ml
RNase-Free Water (bottle)	2 x 125 ml
Proteinase K (green lid)	2 x 1.4 ml
PAXgene RNA Spin Columns (red)	5 x 10
Processing Tubes (2 ml)	6 x 50
Secondary BD Hemogard™ Closures	50
Microcentrifuge Tubes (1.5 ml)	3 x 50, 1 x 10
DNase I, RNase-Free (lyophilized)	1500 Kunitz units‡
Buffer RDD (white lid)	2 x 2 ml
RNase-Free Water (tube, lilac lid)	2 ml
PAXgene Shredder Spin Columns (lilac)	5 x 10
<i>PAXgene Blood RNA Kit Handbook (Version 2)</i>	1

* Not compatible with disinfecting reagents containing bleach. Contains a guanidine salt. See page 7 for safety information.

† Buffer BR4 is supplied as a concentrate. Before using for the first time, add 4 volumes of ethanol (96–100%, purity grade p.a.) as indicated on the bottle to obtain a working solution.

‡ Kunitz units are the commonly used units for measuring DNase I; see page 36 or 41 for definition.

Storage Conditions

PAXgene RNA spin columns, PAXgene Shredder spin columns, proteinase K, and Buffers BR1 to BR5 can be stored dry at room temperature (15–25°C).

The RNase-Free DNase Set, which contains DNase I, Buffer RDD, and RNase-free water (tube), is shipped at ambient temperature. Store all components of the RNase-Free DNase Set immediately upon receipt at 2–8°C.

Intended Use

The PAXgene Blood RNA Kit is for the purification of intracellular RNA from whole blood collected in the PAXgene Blood RNA Tube. When the kit is used in conjunction with the PAXgene Blood RNA Tube, the system provides purified intracellular RNA from whole blood for RT-PCR used in molecular diagnostic testing. See the *PAXgene Blood RNA Tube Product Circular* for information about the use of PAXgene Blood RNA Tubes.

Performance characteristics for the PAXgene Blood RNA System have only been established with FOS and IL1B gene transcripts. The user is responsible for establishing appropriate PAXgene Blood RNA System performance characteristics for other target transcripts.

Product Use Limitations

The PAXgene Blood RNA Kit is intended for purification of intracellular RNA from human whole blood ($4.8 \times 10^6 - 1.1 \times 10^7$ leukocytes/ml) for in vitro diagnostics applications. It is not for the purification of genomic DNA or viral nucleic acids from human whole blood. Due to the limited number of transcripts validated for stabilization specifications (FOS and IL1B gene transcripts), the performance characteristics have not been established for all transcripts. Laboratory personnel should review the manufacturer's data and their own data to determine whether validation is necessary for other transcripts.

Quality Control

In accordance with QIAGEN's ISO-certified Total Quality Management System, each lot of PAXgene Blood RNA Kit is tested against predetermined specifications to ensure consistent product quality.

Technical Assistance

At QIAGEN we pride ourselves on the quality and availability of our technical support. Our Technical Service Departments are staffed by experienced scientists with extensive practical and theoretical expertise in molecular biology and the use of PreAnalytiX products. If you have any questions regarding the PAXgene Blood RNA Kit, please do not hesitate to contact us.

For technical assistance and more information please call QIAGEN Technical Services (see page 55).

Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.

To avoid the risk of infection (e.g., from HIV or hepatitis B viruses) or injury when working with biological and chemical materials, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDSs). These are available online in convenient and compact PDF format at www.preanalytix.com/rna_msdss.asp where you can find, view, and print the MSDSs for this kit.



CAUTION: DO NOT add bleach or acidic solutions directly to the sample-preparation waste.

Buffers BR2 and BR3 contain guanidine thiocyanate, which can form highly reactive compounds when combined with bleach. If Buffer BR2 or Buffer BR3 is spilt, clean with suitable laboratory detergent and water. If liquid containing potentially infectious agents is spilt, clean the affected area first with laboratory detergent and water, and then with 1% (v/v) sodium hypochlorite.

The RNA stabilizing solution and blood mixture from the PAXgene Blood RNA Tube can be disinfected using 1 volume of commercial bleach solution (5% sodium hypochlorite) per 9 volumes of the RNA stabilizing solution and blood mixture.

Sample-preparation waste, such as supernatants from centrifugation steps in the RNA purification procedure, is to be considered potentially infectious. Before disposal, the waste must be autoclaved or incinerated to destroy any infectious material.

The following risk and safety phrases apply to components of the PAXgene Blood RNA Kit. See the *PAXgene Blood RNA Tube Product Circular* for safety information about PAXgene Blood RNA Tubes.

Buffer BR2



X_n

Contains guanidine thiocyanate: harmful (X_n). Risk and safety phrases: *
R20/21/22-32, S13-26-36-46

* R10: Flammable; R20/21/22: Harmful by inhalation, in contact with skin and if swallowed;
R32: Contact with acids liberates very toxic gas; R36/37/38: Irritating to eyes, respiratory system and skin;
R42/43: May cause sensitization by inhalation and skin contact; S13: Keep away from food, drink and animal feedingstuffs; S22: Do not breathe dust; S23: Do not breathe spray; S24: Avoid contact with skin;
S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice;
S36: Wear suitable protective clothing; S36/37: Wear suitable protective clothing and gloves;
S46: If swallowed, seek medical advice immediately and show this container or label.

Buffer BR3

Contains ethanol: flammable. Risk phrase:* R10

Proteinase K



Xn

Contains proteinase K (*Tritirachium album*): sensitizer, irritant. Risk and safety phrases:* R36/37/38-42/43, S23-24-26-36/37

DNase I



Xn

Contains deoxyribonuclease (bovine): sensitizer. Risk and safety phrases:* R42/43, S22-24-26-36/37

24-hour emergency information

Emergency medical information in English, French, and German can be obtained 24 hours a day from:

Poison Information Center Mainz, Germany

Tel: +49-6131-19240

* R10: Flammable; R20/21/22: Harmful by inhalation, in contact with skin and if swallowed;
R32: Contact with acids liberates very toxic gas; R36/37/38: Irritating to eyes, respiratory system and skin;
R42/43: May cause sensitization by inhalation and skin contact; S13: Keep away from food, drink and animal feedingstuffs; S22: Do not breathe dust; S23: Do not breathe spray; S24: Avoid contact with skin;
S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice;
S36: Wear suitable protective clothing; S36/37: Wear suitable protective clothing and gloves;
S46: If swallowed, seek medical advice immediately and show this container or label.

Introduction

Collection of whole blood is the first step in many molecular assays used to study cellular RNA. However, a major problem in such experiments is the instability of the cellular RNA profile *in vitro*. Studies at PreAnalytiX have shown that the copy numbers of individual mRNA species in whole blood can change more than 1000-fold during storage or transport at room temperature.* This is caused both by rapid RNA degradation and by induced expression of certain genes after the blood is drawn. Such changes in the RNA expression profile make reliable studies of gene expression impossible. A method that preserves the RNA expression profile during and after phlebotomy is therefore essential for accurate analysis of gene expression in human whole blood.

Principle and procedure

PreAnalytiX has developed a system that enables the collection, stabilization, storage, and transportation of whole blood specimens, together with a rapid and efficient protocol for purification of intracellular RNA. The system requires the use of PAXgene Blood RNA Tubes (US Patents 6,602,718 and 6,617,170) for blood collection and RNA stabilization, followed by manual or automated RNA purification using the PAXgene Blood RNA Kit. Both manual and automated protocols provide substantially equivalent performance with regards to RNA quality and yield. Performance data for the manual protocol (pages 14–21) and the automated protocol (pages 22–25) are included in this handbook.

Sample collection and stabilization

PAXgene Blood RNA Tubes contain a proprietary reagent composition based on a patented RNA stabilization technology. This reagent composition protects RNA molecules from degradation by RNases and minimizes induction of gene expression. PAXgene Blood RNA Tubes are intended for the collection of whole blood and stabilization of cellular RNA for up to 3 days at 18–25°C (Figure 1, page 10) or up to 5 days at 2–8°C (Figure 2, page 11). Currently available data shows stabilization of cellular RNA for at least 50 months at –20°C or –70°C. For more information from ongoing studies evaluating stability for longer time periods, please contact QIAGEN Technical Services.

The actual duration of RNA stabilization may vary depending upon the species of cellular RNA and the downstream application used. Due to the limited number of transcripts validated for stabilization specifications (FOS and IL1B gene transcripts), the performance characteristics have not been established for all transcripts. Laboratory personnel should review the manufacturer's data and their own data to determine whether validation is necessary for other transcripts.

* Rainen, L. et al. (2002) Stabilization of mRNA expression in whole blood samples. *Clin. Chem.* **48**, 1883.

RNA Stability in Blood Samples at 18–25°C




Figure 1 Blood was drawn from 10 donors, with duplicate samples, and stored at 18–25°C for the indicated number of days, followed by total RNA purification. Blood was collected and stored in PAXgene Blood RNA Tubes, and total RNA was purified using the PAXgene Blood RNA Kit. Relative transcript levels of **A** FOS and **B** IL1B were determined by real-time, duplex RT-PCR, using 18S rRNA as an internal standard. The values for all samples are plotted (20 data sets for each gene, blue lines), with means (red lines) and standard deviations (black bars) of all samples shown. The dashed lines indicate the $\pm 3x$ total precision of the assays (FOS: 2.34 C_t ; IL1B: 1.93 C_t).

RNA Stability in Blood Samples at 2–8°C





Figure 2 Blood was drawn and total RNA purified, after storage at 2–8°C, as described in Figure 1. Relative transcript levels of **A** FOS and **B** IL1B were determined by real-time, duplex RT-PCR, using 18S rRNA as an internal standard. The values for all samples are plotted (20 data sets for each gene, blue lines), with means (red lines) and standard deviations (black bars) of all samples shown. The dashed lines indicate the $\pm 3x$ total precision of the assays (FOS: 2.34 C_t ; IL1B: 1.93 C_t).


RNA concentration and purification

The PAXgene Blood RNA Kit is for the purification of total RNA from 2.5 ml human whole blood collected in a PAXgene Blood RNA Tube. The procedure is simple and can be performed using manual or automated procedures (see flowcharts). In both protocols, purification begins with a centrifugation step to pellet nucleic acids in the PAXgene Blood RNA Tube. The pellet is washed and resuspended, followed by manual or automated RNA purification. In principle, both protocols follow the same protocol steps with the same kit components.

The Manual PAXgene Blood RNA Procedure



The Automated PAXgene Blood RNA Procedure



Manual RNA purification

In detail, the resuspended pellet is incubated in optimized buffers together with proteinase K to bring about protein digestion. An additional centrifugation through the PAXgene Shredder spin column is carried out to homogenize the cell lysate and remove residual cell debris, and the supernatant of the flow-through fraction is transferred to a fresh microcentrifuge tube. Ethanol is added to adjust binding conditions, and the lysate is applied to a PAXgene RNA spin column. During a brief centrifugation, RNA is selectively bound to the PAXgene silica membrane as contaminants pass through. Remaining contaminants are removed in several efficient wash steps. Between the first and second wash steps, the membrane is treated with DNase I to remove trace amounts of bound DNA. After the wash steps, RNA is eluted in elution buffer (Buffer BR5) and heat-denatured.

Total RNA purified using the PAXgene Blood RNA System is highly pure. Using the manual protocol, A_{260}/A_{280} values are between 1.8 and 2.2, and $\leq 1\%$ (w/w) genomic DNA is present in $\geq 95\%$ of all samples, as measured by quantitative, real-time PCR of a sequence of the beta-actin gene. At least 95% of samples show no inhibition in RT-PCR, when using up to 30% of the eluate.

Using the manual protocol, average sample preparation time (based on data from 12 sample preps) is approximately 90 minutes, with only 40 minutes of hands-on time. RNA yields from 2.5 ml healthy human whole blood are $\geq 3 \mu\text{g}$ for $\geq 95\%$ of the samples processed. Since yields are highly donor-dependent, individual yields may vary. For individual donors, the PAXgene Blood RNA system provides highly reproducible and repeatable yields (Figures 3 and 4, pages 15 and 16) and reproducible and repeatable RT-PCR (Figures 5 and 6, pages 19 and 20), making it highly robust for clinical diagnostic tests.

Figure 3 indicates the overall repeatability and reproducibility of the PAXgene Blood RNA System. Additional studies were conducted to show the influence of different PAXgene Blood RNA Kit lots and different operators on the reproducibility of RNA yield and real time RT-PCR performance. As pooled blood samples instead of individual PAXgene tubes were used for these studies, the results do not reflect the system repeatability, including fluctuation between individual blood draws, but only the repeatability of the sample preparation (see Figure 4).

Reproducible and Repeatable RNA Purification




Figure 3 Quadruplicate blood samples from 14 donors were manually processed by each of 3 technicians (A, B, C). Three sets of equipment were used, and all samples prepared by a single technician were processed using the same equipment. **A** Means and standard deviations of RNA yield per replicate samples from the same donors and different technicians are shown. **B** Twelve replicate blood samples from each of 14 donors were processed by the 3 different technicians. Means and standard deviations of RNA yield per samples from the same donors and all technicians are presented. For all RNA samples, A_{260}/A_{280} ratios ranged from 1.8 to 2.2.

Repeatability and Reproducibility of RNA Yield for Different Operators and PAXgene Blood RNA Kit Lots Using Pooled Blood Samples




Figure 4 Blood samples from 30 different donors were collected in PAXgene Blood RNA tubes (12 tubes per donor, 360 tubes in total). The contents of the tubes from 3 donors were pooled and subsequently re-aliquoted into 36 samples. These 36 samples per 3-donor-pool were manually processed by 3 different operators. Each operator used 3 different PAXgene Blood RNA Kit lots for the extraction and processed quadruplicate samples from each of the 10 donor pools. **A** RNA yield and standard deviation for every operator–lot combination. Quadruplicate blood samples from 10 donor pools were processed by 3 different operators (A, B, C) with each of 3 kit lots (1, 2, 3). The mean yields (columns) and standard deviations (error bars) per quadruplicate sample from the same donor pool for different operator and different kit lot are presented. **B** CV of RNA yield per donor pool for all operator–lot combinations (A, B, C; 1, 2, 3) as calculated from the mean yield and standard deviation of the yield shown in Figure 4A.

Table 1A. Reproducibility within Each Lot and within Each User for Selected Donor Pools (1, 6, 9, 10)

Combination of data	Donor pool 1 5.1 x 10 ⁶ cells/ml			Donor pool 6 6.5 x 10 ⁶ cells/ml			Donor pool 9 8.4 x 10 ⁶ cells/ml			Donor pool 10 10.2 x 10 ⁶ cells/ml		
	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)
Lot 1, user A	8.03	0.42	5	9.55	0.99	10	7.52	0.41	6	7.96	0.49	6
Lot 1, user B	7.98	1.17	15	9.38	1.94	21	8.82	1.72	19	8.90	0.76	9
Lot 1, user C	7.87	0.45	6	10.71	0.65	6	10.14	1.46	14	10.22	1.29	13
Lot 2, user A	7.32	0.98	13	9.78	1.89	19	6.92	0.27	4	7.63	1.23	16
Lot 2, user B	6.09	1.04	17	9.82	2.83	29	7.20	0.71	10	7.00	0.56	8
Lot 2, user C	6.87	0.31	4	10.37	0.74	7	9.14	1.52	17	11.56	1.21	10
Lot 3, user A	7.04	0.90	13	8.96	0.68	8	8.18	0.76	9	7.85	0.82	10
Lot 3, user B	6.98	1.22	17	7.73	0.97	13	6.41	0.88	14	8.88	2.17	24
Lot 3, user C	8.78	0.89	10	10.59	1.94	18	10.78	0.56	5	10.88	0.37	3

Table 1B. Reproducibility within Each User and between All Lots for Selected Donor Pools (1, 6, 9, 10)

Combination of data	Donor pool 1 5.1 x 10 ⁶ cells/ml			Donor pool 6 6.5 x 10 ⁶ cells/ml			Donor pool 9 8.4 x 10 ⁶ cells/ml			Donor pool 10 10.2 x 10 ⁶ cells/ml		
	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)
User A, all lots	7.46	0.85	11	9.43	1.22	13	7.54	0.72	10	7.81	0.82	11
User B, all lots	7.02	1.31	19	8.98	2.09	23	7.48	1.50	20	8.26	1.54	19
User C, all lots	7.84	0.98	13	10.56	1.15	11	10.02	1.34	13	10.89	1.10	10

Table 1C. Reproducibility within Each Lot and between All Users for Selected Donor Pools (1, 6, 9, 10)

Combination of data	Donor pool 1			Donor pool 6			Donor pool 9			Donor pool 10		
	5.1 x 10 ⁶ cells/ml			6.5 x 10 ⁶ cells/ml			8.4 x 10 ⁶ cells/ml			10.2 x 10 ⁶ cells/ml		
	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)
Lot 1, all users	7.96	0.69	9	9.88	1.34	14	8.83	1.63	19	9.02	1.27	14
Lot 2, all users	6.76	0.93	14	9.99	1.84	18	7.75	1.36	18	8.73	2.31	26
Lot 3, all users	7.60	1.27	17	9.09	1.71	19	8.46	1.99	24	9.20	1.80	20

Table 1D. Reproducibility between All Lots and All Users for Selected Donor Pools (1, 6, 9, 10)

Combination of data	Donor pool 1			Donor pool 6			Donor pool 9			Donor pool 10		
	5.1 x 10 ⁶ cell/ml			6.5 x 10 ⁶ cell/ml			8.4 x 10 ⁶ cell/ml			10.2 x 10 ⁶ cell/ml		
	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)	Mean yield (µg)	SD (µg)	CV (%)
All lots and all users	7.44	1.09	15	9.66	1.65	17	8.35	1.70	20	8.99	1.80	20

Detailed analysis of 4 representative donor pools. The pools were selected according to the white blood cell count and reflect the upper, medium, and the lower values of the normal range of white blood cell counts (4.8 x 10⁶ – 1.1 x 10⁷ leukocytes/ml). The white blood cell count represents the mean value of the 3 white blood cell counts from the 3 donors per donor pool.

Reproducibility of RT-PCR — between Users




Figure 5 RNA purified in the experiment described in Figure 4 was used for real-time RT-PCR. Relative transcript levels of **A** FOS and **B** IL1B were determined by real-time, duplex RT-PCR using 18S rRNA as an internal standard. The values for all samples are plotted, relative to the values for user 1 (10 donor pools \times 3 kit lots \times 4 replicates = 120 data sets for each gene), with means (red lines) and standard deviations (black bars) for all samples shown. The dashed lines indicate the $\pm 3x$ total precision of the assays (FOS: 2.34 C_t ; IL1B, 1.93 C_t).

Reproducibility of RT-PCR — between Kit Lots




Figure 6 RNA purified in the experiment described in Figure 4 was used for real-time RT-PCR. Relative transcript levels of **A** FOS and **B** IL1B were determined by real-time, duplex RT-PCR using 18S rRNA as an internal standard. The values for all samples are plotted, relative to the values for kit lot 1 (10 donor pools \times 3 users \times 4 replicates = 120 data sets for each gene), with means (red lines) and standard deviations (black bars) for all samples shown. The dashed lines indicate the $\pm 3x$ total precision of the assays (FOS: $2.34 C_t$; IL1B, $1.93 C_t$).

Table 2. Summary of RT-PCR Data from Figures 5 and 6

Test system	FOS/18S rRNA assay		IL1B/18S rRNA assay	
	Mean ($\Delta\Delta C_T$)	\pm SD ($\Delta\Delta C_T$)	Mean ($\Delta\Delta C_T$)	\pm SD ($\Delta\Delta C_T$)
Comparison of data				
Reproducibility within each user and between all lots				
All users, lot 1 – lot 1	0.00	0.00	0.00	0.00
All users, lot 1 – lot 2	-0.03	0.48	-0.07	0.66
All users, lot 1 – lot 3	-0.21	0.52	0.11	0.71
Reproducibility within each lot and between all users				
All lots, user A – user A	0.00	0.00	0.00	0.00
All lots, user A – user B	-0.46	0.44	-0.06	0.69
All lots, user A – user C	-0.31	0.60	-0.15	0.71

User: Technician, performed the study.

Lot: Number of kit lot used in this study.

SD: Standard deviation.

Mean $\Delta\Delta C_T$ values (N = 120) and standard deviations are shown for the data presented in Figures 5 and 6.

Automated RNA purification

Sample preparation using the QIAcube® follows the same steps as the manual procedure, enabling you to continue using the PAXgene Blood RNA Kit for purification of high-quality RNA. See the *QIAcube User Manual* and www.qiagen.com/MyQIAcube for more information about the QIAcube.

The automated RNA purification protocol consists of 2 parts (or protocols), “PAXgene Blood RNA Part A” and “PAXgene Blood RNA Part B”, with a brief manual intervention between the 2 parts.

The centrifuged, washed, and resuspended nucleic acid pellet (see “RNA concentration and purification”, page 12) is transferred from the PAXgene Blood RNA Tube into processing tubes, which are placed into the thermoshaker unit on the QIAcube worktable. The operator selects and starts the “PAXgene Blood RNA Part A” protocol from the menu. The QIAcube performs the steps of the protocol through to elution of RNA in elution buffer. The operator transfers the microcentrifuge tubes, containing the purified RNA, into the thermoshaker unit of the QIAcube. The operator selects and starts the “PAXgene Blood RNA Part B” protocol from the menu, and heat denaturation is performed by the QIAcube.

Average sample preparation time (based on data from 12 sample preps) is 125 minutes, with only approximately 20 minutes of hands-on time.

RNA yields from 2.5 ml healthy human whole blood are $\geq 3 \mu\text{g}$ for $\geq 95\%$ of the samples processed. Figure 7 indicates the RNA yields from a total of 288 samples prepared using the automated protocol with 3 kit lots by 3 operators. As pooled blood samples instead of individual PAXgene Blood RNA Tubes were used for these studies, the results do not reflect the RNA yield expected from single samples of individual blood draws. Since yields are highly donor-dependent, individual yields may vary (Figure 7).

RNA Yield — Automated Processing




Figure 7 Blood samples from 48 different donors were collected in PAXgene Blood RNA Tubes (6 tubes per donor, 288 tubes in total). The contents of the tubes from 6 donors were pooled and subsequently realiquoted into 36 samples. These 36 samples per 6-donor-pool were processed by 3 different operators (A, B, C). Each operator used 3 different lots (1, 2, 3) of the PAXgene Blood RNA Kit for automated extraction and processed quadruplicate samples from each of the 8 donor pools. RNA yields of all individual samples are shown for every operator–lot combination.

At least 95% of samples show no inhibition in RT-PCR, when using up to 30% of the eluate.

RNA purified with the PAXgene Blood RNA System and the automated protocol is highly pure, as shown by lack of RT-PCR inhibition (see above) and A_{260}/A_{280} values between 1.8 and 2.2. Genomic DNA is present at $\leq 1\%$ (w/w) in $\geq 95\%$ of all samples, as measured by quantitative, real-time PCR of a sequence of the beta-actin gene. Figures 8 and 9 show the A_{260}/A_{280} values and relative genomic DNA of a total of 288 samples prepared using the automated protocol with 3 kit lots by 3 operators.

RNA Purity (A_{260}/A_{280} Values) — Automated Processing




Figure 8 RNA was purified by 3 different operators (A, B, C) using 3 different lots (1, 2, 3) of the PAXgene Blood RNA Kit in the experiment described in Figure 7. A_{260}/A_{280} values of all individual samples are shown for every operator–lot combination.

RNA Purity (% Genomic DNA Contamination) — Automated Processing




Figure 9 RNA was purified by 3 different operators (A, B, C) using 3 different lots (1, 2, 3) of the PAXgene Blood RNA Kit in the experiment described in Figure 7. Genomic DNA amounts (w/w) in all individual samples are shown for every operator–lot combination.

The automated protocol of RNA purification using the PAXgene Blood RNA System provides highly reproducible and repeatable RT-PCR results, as shown in Figure 10, making it highly robust for clinical diagnostic tests.

Reproducibility of RT-PCR — between Automated and Manual Protocols




Figure 10 RNA was purified by 3 different operators (A, B, C) using 3 different lots (1, 2, 3) of the PAXgene Blood RNA Kit using the automated protocol in the experiment described in Figure 7. In parallel, RNA was purified from the corresponding replicate tubes using the manual protocol. Relative transcript levels of **A** FOS and **B** IL1B were determined by real-time, duplex RT-PCR using 18S rRNA as an internal standard. Possible differences of transcript levels between RNA prepared from paired blood samples using both extraction protocols (automated and manual protocol) were calculated by the $\Delta\Delta C_t$ method. Individual $\Delta\Delta C_t$ values for all sample pairs (4 replicates \times 8 donor pools \times 3 kit lots \times 3 operators = 288 pairs for each gene) are plotted as single dots with means (larger dots) and standard deviations (black bars) for all samples shown. The dashed lines indicate the $\pm 3x$ total precision of the assays (FOS: 2.34 C_t ; IL1B, 1.93 C_t).

Equipment and Reagents to Be Supplied by User

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate material safety data sheets (MSDSs), available from the product supplier.

For all protocols

- PAXgene Blood RNA Tubes (cat. no. 762165)
- Ethanol (96–100%)
- Pipets* (10 μ l – 4 ml)
- Sterile, aerosol-barrier, RNase-free pipet tips[†]
- Graduated cylinder[†]
- Centrifuge* capable of attaining 3000–5000 $\times g$, and equipped with a swing-out rotor and buckets to hold PAXgene Blood RNA Tubes
- Vortex mixer*
- Crushed ice
- Permanent pen for labeling

For the manual protocol

- Variable-speed microcentrifuge* capable of attaining 1000–8000 $\times g$, and equipped with a rotor for 2 ml microcentrifuge tubes
- Shaker-incubator* capable of incubating at 55°C and 65°C and shaking at \geq 400 rpm, not exceeding 1400 rpm (e.g., Eppendorf[®] Thermomixer Compact, or equivalent)

For the automated protocol

- QIAcube* (QIAGEN, cat. no. 9001292)
- QIAcube consumables
 - Filter-Tips, 1000 μ l (1024) (QIAGEN, cat. no. 990352)[§]
 - Reagent Bottles, 30 ml (6) (QIAGEN, cat. no. 990393)[§]
 - Rotor Adapters (10 \times 24) (QIAGEN, cat. no. 990394)[§]
- QIAcube accessories
 - Reagent Bottle Rack (QIAGEN, cat. no. 990390)[§]
 - Rotor Adapter Holder (QIAGEN, cat. no. 990392)[§]
- Scissors

* Ensure that instruments have been checked, maintained, and calibrated regularly according to the manufacturer's recommendations.

[†] Ensure that you are familiar with the guidelines on handling RNA (Appendix A, page 47).

[‡] For the addition of ethanol to Buffer BR4 concentrate.

[§] Also included in the Starter Pack, QIAcube (QIAGEN, cat. no. 990395).

Important Notes

Using the QIAcube

Ensure that you are familiar with operating the QIAcube. Please read the *QIAcube User Manual* and any additional information supplied with the QIAcube, paying careful attention to the safety information, before beginning the automated PAXgene Blood RNA protocols.

Starting the QIAcube

Close the QIAcube door, and switch on the QIAcube with the power switch (see Figure 11).

A beeper sounds and the startup screen appears. The instrument automatically performs initialization tests.

Installing protocols on the QIAcube

An initial protocol installation is required before the first RNA preparation run on the QIAcube can be performed. Install both "PAXgene Blood RNA Part A" and "PAXgene Blood RNA Part B" protocols.

Protocols are provided at www.qiagen.com/MyQIAcube and need to be downloaded to the USB stick supplied with the QIAcube and transferred to the QIAcube via the USB port.

The USB port, located behind the protective panel (see Figure 11), allows connection of the QIAcube to a USB stick (supplied with the QIAcube). Data files, such as log files or report files can also be transferred via the USB port from the QIAcube to the USB stick.

Note: The USB port is only for use with the USB stick provided by QIAGEN. Do not connect other devices to this port.

Note: Do not remove the USB stick while downloading protocols or transferring data files or during a protocol run.

Front View of the QIAcube



Figure 11

1 Touchscreen

2 Door

3 RS232 serial port behind protective panel (for use by QIAGEN Instrument Service Specialists only)

4 USB port behind protective panel

5 Power switch

6 Waste drawer

Loading the QIAcube

To save time, loading can be performed during one or both of the 10-minute centrifugation steps (steps 3 and 5) in “Protocol: Automated Purification of Total RNA from Human Whole Blood Collected into PAXgene Blood RNA Tubes”, page 42.

Reagent bottles

Carefully fill 4 QIAcube reagent bottles with the reagents listed in Table 3 (fill the bottles up to the indicator level on the reagent bottles). Label the bottles and lids clearly with buffer names and place into the appropriate position in the reagent bottle rack (see Figures 12 and 13).

Before every run on the QIAcube, make sure that the reagent bottles are filled up to the indicator levels (remaining volumes in the original kit buffer bottles should be used to fill the reagent bottles). Position the reagent rack with filled reagent bottles onto the QIAcube worktable as shown (Figures 12 and 13).

Note: Be sure to remove lids from the bottles before placing onto the worktable.

Note: Buffer volumes provided in the PAXgene Blood RNA Kit (50) are sufficient for a maximum of 7 RNA preparation runs on the QIAcube. Multiple runs with few samples should be avoided in order allow sufficient buffer volumes for processing the full 50 samples.

Table 3. Positions in the reagent bottle rack

Position	Reagent
1	Binding buffer (BR2)
2	96–100% ethanol
3	Wash buffer 1 (BR3)
4	Wash buffer 2 (BR4) *
5	– (leave empty)
6	– (leave empty)

* Wash buffer 2 (BR4) is supplied as a concentrate. Before using for the first time, add 4 volumes of ethanol (96–100%) as indicated on the bottle to obtain a working solution.

Loading the Reagent Bottle Rack




Figure 12 **A** Schematic of positions and contents of bottles in the reagent bottle rack. **B** Loading the rack onto the QIAcube.

Internal View of the QIAcube




Figure 13

- 1** Centrifuge lid
- 2** Centrifuge
- 3** Shaker
- 4** Reagent bottle rack
- 5** Tip sensor
- 6** Microcentrifuge tube slots
- 7** Tip racks
- 8** Disposal slots for tips and columns
- 9** Robotic arm

Spin columns, microcentrifuge tubes, and QIAcube plasticware

Place 2 tip racks filled with Filter-Tips 1000 µl onto the QIAcube (see Figure 13). Refill racks with tips when necessary.

Note: Only use 1000 µl filter-tips designed for use with the QIAcube.

Label rotor adapters and microcentrifuge tubes for each sample using a permanent pen. Open the PAXgene Shredder spin columns to be used, and cut the lids off completely using scissors (see Figure 14).

Note: For proper operation of the QIAcube robotic gripper, completely remove (cut off) the lids and all plastic parts connecting the lid to the PAXgene Shredder spin columns (see Figure 14). Otherwise, the robotic gripper cannot grip the spin columns properly. Load the PAXgene RNA spin column, PAXgene Shredder spin column (without lid), and labeled microcentrifuge tube into the appropriate positions in each labeled rotor adapter as shown in Table 4 and Figure 15 (page 32).

Note: Make sure that the spin column and microcentrifuge tube lids are pushed all the way down to the bottom of the slots at the edge of the rotor adapter otherwise the lids will break off during centrifugation.

Loading a PAXgene Shredder Spin Column



Figure 14 The PAXgene Shredder spin column is loaded into the middle position of the rotor adapter. Cut off the lid before loading the column.

Table 4. Labware in the rotor adapter

Position	Labware	Lid position
1	PAXgene RNA spin column (red)	L1
2	PAXgene Shredder spin column (lilac) (cut off lid before placing in rotor adapter)	–
3	Microcentrifuge tube*	L3

* Use the microcentrifuge tubes (1.5 ml) included in the PAXgene Blood RNA Kit.

Positions in the Rotor Adapter




Figure 15 The rotor adapter has three tube positions (1–3) and three lid positions (L1–L3).

Loading the centrifuge

Load the assembled rotor adapters into the centrifuge buckets as shown in Figure 16.

Note: If processing fewer than 12 samples, make sure to load the centrifuge rotor balanced radially (see Figure 17). All centrifuge buckets must be mounted before starting a protocol run, even if fewer than 12 samples are to be processed. A single (one) sample or 11 samples cannot be processed.

Loading the Centrifuge



Figure 16 Load the assembled rotor adapters into the centrifuge buckets.

Loading the Centrifuge and Shaker




Figure 17 Centrifuge and shaker positions are shown for processing from two to ten samples. One or 11 samples cannot be processed.

Processing tubes

Remove any processing tubes left in the microcentrifuge tube slots from previous runs (see Figure 13, page 30). Fill 3 processing tubes with the amount of reagents given in Table 5, according to the number of samples in the run. Label the tubes clearly with reagent names, and place them into the appropriate position in the microcentrifuge tube slots, as indicated in Table 6.

For DNase I incubation mix, pipet the indicated volume of DNA digestion buffer into a processing tube, and add the indicated volume of DNase I stock solution. Mix by gently pipeting the complete mixture up and down 3 times using a 1000 µl pipet tip.

Note: DNase I is especially sensitive to physical denaturation. Mix only by pipeting, using wide-bore pipet tips to reduce shearing. Do not vortex.

Use the 2 ml processing tubes included in the PAXgene Blood RNA Kit.

Note: Be sure to only pipet the required volume as indicated in Table 5.

Table 5. Volume of reagents required in processing tubes for the microcentrifuge tube slots

Number of samples	Volume of reagent required for the indicated number of samples (µl)		
	Proteinase K	DNase I incubation mix	Elution buffer (Buffer BR5)
2	126	187 (23 DNase I + 164 Buffer RDD)	313
3	170	261 (33 DNase I + 228 Buffer RDD)	399
4	213	334 (42 DNase I + 292 Buffer RDD)	486
5	256	407 (51 DNase I + 356 Buffer RDD)	572
6	299	481 (60 DNase I + 421 Buffer RDD)	658
7	342	554 (69 DNase I + 485 Buffer RDD)	745
8	386	627 (78 DNase I + 549 Buffer RDD)	831
9	429	701 (88 DNase I + 613 Buffer RDD)	918
10	472	775 (97 DNase I + 678 Buffer RDD)	1004
12	558	921 (115 DNase I + 806 Buffer RDD)	1177

Table 6. Microcentrifuge tube slots

	Position		
	A	B	C
Reagent	Proteinase K	DNase I incubation mix	Elution buffer (Buffer BR5)
Container	Processing tube*	Processing tube*	Processing tube*

* Use the 2 ml processing tubes included in the PAXgene Blood RNA Kit.

Protocol: Manual Purification of Total RNA from Human Whole Blood Collected into PAXgene Blood RNA Tubes

Important points before starting

- Make sure that the kit box is intact and undamaged, and that buffers have not leaked. Do not use a kit that has been damaged.
- When using a pipet, ensure that it is set to the correct volume, and that liquid is carefully and completely aspirated and dispensed.
- To avoid transferring samples to the wrong tube or spin column, ensure that all tubes and spin columns are properly labeled using a permanent pen. Label the lid and the body of each tube. For spin columns, label the body of its processing tube. Close each tube or spin column after liquid is transferred to it.
- Spillages of samples and buffers during the procedure may reduce the yield and purity of RNA.
- Unless otherwise indicated, all steps of this protocol, including centrifugation steps, should be carried out at room temperature (15–25°C).
- Because of the sensitivity of nucleic acid amplification technologies, the following precautions are necessary when handling samples to avoid cross-contamination:
 - Carefully pipet the sample into the spin column without moistening the rim of the column.
 - Always change pipet tips between liquid transfers. Use aerosol-barrier pipet tips.
 - Avoid touching the spin column membrane with the pipet tip.
 - After vortexing or heating a microcentrifuge tube, briefly centrifuge it to remove drops from the inside of the lid.
 - Wear gloves throughout the entire procedure. In case of contact between gloves and sample, change gloves immediately.
 - Close the spin column before placing it in the microcentrifuge. Centrifuge as described in the procedure.
 - Open only one spin column at a time, and take care to avoid generating aerosols.
 - For efficient parallel processing of multiple samples, fill a rack with processing tubes to which the spin columns can be transferred after centrifugation. Discard the used processing tubes containing flow-through, and place the new processing tubes containing spin columns directly in the microcentrifuge.

Things to do before starting

- Blood must be collected in PAXgene Blood RNA Tubes according to the instructions in the *PAXgene Blood RNA Tube Product Circular*. If necessary, see Appendix C (page 50) for recommendations on handling PAXgene Blood RNA Tubes.
- Ensure that the PAXgene Blood RNA Tubes are incubated for at least 2 hours at room temperature after blood collection to ensure complete lysis of blood cells. Incubation of the PAXgene Blood RNA Tube overnight may increase yields. If the PAXgene Blood RNA Tube was stored at 2–8°C or –20°C or –70°C after blood collection, first equilibrate it to room temperature, and then store it at room temperature for 2 hours before starting the procedure.
- Read the safety information on page 7.
- Read the guidelines on handling RNA (Appendix A, page 47).
- Ensure that instruments, such as pipets and the shaker–incubator, have been checked and calibrated regularly according to the manufacturer’s recommendations.
- A shaker–incubator is required in steps 5 and 20. Set the temperature of the shaker–incubator to 55°C.
- Buffer BR2 may form a precipitate upon storage. If necessary, warm to 37°C to dissolve.
- Buffer BR4 is supplied as a concentrate. Before using for the first time, add 4 volumes of ethanol (96–100%) as indicated on the bottle to obtain a working solution.
- If using the RNase-Free DNase Set for the first time, prepare DNase I stock solution. Dissolve the solid DNase I (1500 Kunitz units)* in 550 µl of the RNase-free water provided with the set. Take care that no DNase I is lost when opening the vial. Do not vortex the reconstituted DNase I. DNase I is especially sensitive to physical denaturation. Mixing should only be carried out by gently inverting the tube.
- Current data shows that reconstituted DNase I can be stored at 2–8°C for up to 6 weeks. For long-term storage of DNase I, remove the stock solution from the glass vial, divide it into single-use aliquots (use the 1.5 ml microcentrifuge tubes supplied with the kit; there are enough for 5 aliquots), and store at –20°C for up to 9 months. Thawed aliquots can be stored at 2–8°C for up to 6 weeks. Do not refreeze the aliquots after thawing.
- When reconstituting and aliquoting DNase I, ensure that you follow the guidelines for handling RNA (Appendix A, page 47).

* Kunitz units are the commonly used units for measuring DNase I, defined as the amount of DNase I that causes an increase in A_{260} of 0.001 per minute per milliliter at 25°C, pH 5.0, with highly polymerized DNA as the substrate (Kunitz, M. (1950) *J. Gen. Physiol.* **33**, 349 and 363).

Procedure

1. **Centrifuge the PAXgene Blood RNA Tube for 10 minutes at 3000–5000 x g using a swing-out rotor.**

Note: Ensure that the blood sample has been incubated in the PAXgene Blood RNA Tube for a minimum of 2 hours at room temperature (15–25°C), in order to achieve complete lysis of blood cells.

Note: The rotor must contain tube adapters for round-bottom tubes. If other types of tube adapter are used, the tubes may break during centrifugation.

2. **Remove the supernatant by decanting or pipetting. Add 4 ml RNase-free water to the pellet, and close the tube using a fresh secondary BD Hemogard closure (supplied with the kit).**

If the supernatant is decanted, take care not to disturb the pellet, and dry the rim of the tube with a clean paper towel.

3. **Vortex until the pellet is visibly dissolved, and centrifuge for 10 minutes at 3000–5000 x g using a swing-out rotor. Remove and discard the entire supernatant.**

Small debris remaining in the supernatant after vortexing but before centrifugation will not affect the procedure.

Note: Incomplete removal of the supernatant will inhibit lysis and dilute the lysate, and therefore affect the conditions for binding RNA to the PAXgene membrane.

4. **Add 350 µl Buffer BR1, and vortex until the pellet is visibly dissolved.**
5. **Pipet the sample into a 1.5 ml microcentrifuge tube. Add 300 µl Buffer BR2 and 40 µl proteinase K. Mix by vortexing for 5 seconds, and incubate for 10 minutes at 55°C using a shaker–incubator at 400–1400 rpm. After incubation, set the temperature of the shaker–incubator to 65°C (for step 20).**

Note: Do not mix Buffer BR2 and proteinase K together before adding them to the sample.

6. **Pipet the lysate directly into a PAXgene Shredder spin column (lilac) placed in a 2 ml processing tube, and centrifuge for 3 minutes at maximum speed (but not to exceed 20,000 x g).**

Note: Carefully pipet the lysate into the spin column and visually check that the lysate is completely transferred to the spin column.

To prevent damage to columns and tubes, do not exceed 20,000 x g.

Note: Some samples may flow through the PAXgene Shredder spin column without centrifugation. This is due to low viscosity of some samples and should not be taken as an indication of product failure.

7. **Carefully transfer the entire supernatant of the flow-through fraction to a fresh 1.5 ml microcentrifuge tube without disturbing the pellet in the processing tube.**

8. Add 350 μ l ethanol (96–100%). Mix by vortexing, and centrifuge briefly (1–2 seconds at 500–1000 \times g) to remove drops from the inside of the tube lid.

Note: The length of the centrifugation **must not** exceed 1–2 seconds, as this may result in pelleting of nucleic acids and reduced yields of total RNA.

9. Pipet 700 μ l sample into the PAXgene RNA spin column (red) placed in a 2 ml processing tube, and centrifuge for 1 minute at 8000–20,000 \times g. Place the spin column in a new 2 ml processing tube, and discard the old processing tube containing flow-through.
10. Pipet the remaining sample into the PAXgene RNA spin column, and centrifuge for 1 minute at 8000–20,000 \times g. Place the spin column in a new 2 ml processing tube, and discard the old processing tube containing flow-through.

Note: Carefully pipet the sample into the spin column and visually check that the sample is completely transferred to the spin column.

11. Pipet 350 μ l Buffer BR3 into the PAXgene RNA spin column. Centrifuge for 1 minute at 8000–20,000 \times g. Place the spin column in a new 2 ml processing tube, and discard the old processing tube containing flow-through.
12. Add 10 μ l DNase I stock solution to 70 μ l Buffer RDD in a 1.5 ml microcentrifuge tube. Mix by gently flicking the tube, and centrifuge briefly to collect residual liquid from the sides of the tube.

If processing, for example, 10 samples, add 100 μ l DNase I stock solution to 700 μ l Buffer RDD. Use the 1.5 ml microcentrifuge tubes supplied with the kit.

Note: DNase I is especially sensitive to physical denaturation. Mixing should only be carried out by gently flicking the tube. Do not vortex.

13. Pipet the DNase I incubation mix (80 μ l) directly onto the PAXgene RNA spin column membrane, and place on the benchtop (20–30°C) for 15 minutes.

Note: Ensure that the DNase I incubation mix is placed directly onto the membrane. DNase digestion will be incomplete if part of the mix is applied to and remains on the walls or the O-ring of the spin column.

14. Pipet 350 μ l Buffer BR3 into the PAXgene RNA spin column, and centrifuge for 1 minute at 8000–20,000 \times g. Place the spin column in a new 2 ml processing tube, and discard the old processing tube containing flow-through.
15. Pipet 500 μ l Buffer BR4 to the PAXgene RNA spin column, and centrifuge for 1 minute at 8000–20,000 \times g. Place the spin column in a new 2 ml processing tube, and discard the old processing tube containing flow-through.

Note: Buffer BR4 is supplied as a concentrate. Ensure that ethanol is added to Buffer BR4 before use (see “Things to do before starting”, page 36).

16. Add another 500 μ l Buffer BR4 to the PAXgene RNA spin column. Centrifuge for 3 minutes at 8000–20,000 \times g.

17. Discard the tube containing the flow-through, and place the PAXgene RNA spin column in a new 2 ml processing tube. Centrifuge for 1 minute at 8000–20,000 x g.
18. Discard the tube containing the flow-through. Place the PAXgene RNA spin column in a 1.5 ml microcentrifuge tube, and pipet 40 µl Buffer BR5 directly onto the PAXgene RNA spin column membrane. Centrifuge for 1 minute at 8000–20,000 x g to elute the RNA.

It is important to wet the entire membrane with Buffer BR5 in order to achieve maximum elution efficiency.

19. Repeat the elution step (step 18) as described, using 40 µl Buffer BR5 and the same microcentrifuge tube.
20. Incubate the eluate for 5 minutes at 65°C in the shaker-incubator (from step 5) without shaking. After incubation, chill immediately on ice.

Note: This incubation at 65°C denatures the RNA for downstream applications. Do not exceed the incubation time or temperature.

21. If the RNA samples will not be used immediately, store at –20°C or –70°C. Since the RNA remains denatured after repeated freezing and thawing, it is not necessary to repeat the incubation at 65°C. If using the RNA samples in a diagnostic assay, follow the instructions supplied by the manufacturer.

Note: For quantification in Tris buffer, use the relationship $A_{260} = 1 \Rightarrow 44 \mu\text{g/ml}$. See Appendix B, page 48.

Protocol: Automated Purification of Total RNA from Human Whole Blood Collected into PAXgene Blood RNA Tubes

Important points before starting

- Make sure that the kit box is intact and undamaged, and that buffers have not leaked. Do not use a kit that has been damaged.
- When using a pipet, ensure that it is set to the correct volume, and that liquid is carefully and completely aspirated and dispensed.
- To avoid transferring samples to the wrong tubes and plastic consumables, ensure that all processing tubes, microcentrifuge tubes, and rotor adapters are properly labeled using a permanent pen. Label the lid and the body of each microcentrifuge tube, the body of each processing tube, and the outer wall of each rotor adapter.
- Spillages of samples and buffers during the procedure may reduce the yield and purity of RNA.
- Unless otherwise indicated, all steps of this protocol, including centrifugation steps, should be carried out at room temperature (15–25°C).
- Because of the sensitivity of nucleic acid amplification technologies, the following precautions are necessary when handling samples to avoid cross-contamination:
 - Carefully pipet the sample into the processing tube, on the bottom of the tube without moistening the rim of the tube.
 - Always change pipet tips between liquid transfers. Use aerosol-barrier pipet tips.
 - Avoid touching the spin column membrane with the pipet tip.
 - After vortexing or heating a microcentrifuge tube, briefly centrifuge it to remove drops from the inside of the lid.
 - Wear gloves throughout the entire procedure. In case of contact between gloves and sample, change gloves immediately.

Things to do before starting

- Blood must be collected in PAXgene Blood RNA Tubes according to the instructions in the *PAXgene Blood RNA Tube Product Circular*. If necessary, see Appendix C (page 50) for recommendations on handling PAXgene Blood RNA Tubes.
- Ensure that the PAXgene Blood RNA Tubes are incubated for at least 2 hours at room temperature after blood collection to ensure complete lysis of blood cells. Incubation of the PAXgene Blood RNA Tube overnight may increase yields. If the PAXgene Blood RNA Tube was stored at 2–8°C or –20°C or –70°C after blood collection, first equilibrate it to room temperature, and then store it at room temperature for 2 hours before starting the procedure.

- Read the safety information on page 7.
- Read “Important Notes”, pages 27–34.
- Read the guidelines on handling RNA (Appendix A, page 47).
- Read the *QIAcube User Manual* and any additional information supplied with the QIAcube, paying careful attention to the safety information.
- Ensure that instruments, such as pipets and the QIAcube, have been checked and calibrated regularly according to the manufacturer’s recommendations.
- Binding buffer (Buffer BR2) may form a precipitate upon storage. If necessary, warm to 37°C to dissolve.
- Wash buffer 2 (Buffer BR4) is supplied as a concentrate. Before using for the first time, add 4 volumes of ethanol (96–100%) as indicated on the bottle to obtain a working solution.
- If using the RNase-Free DNase Set for the first time, prepare DNase I stock solution. Dissolve the solid DNase I (1500 Kunitz units)* in 550 µl of the RNase-free water provided with the set. Take care that no DNase I is lost when opening the vial. Do not vortex the reconstituted DNase I. DNase I is especially sensitive to physical denaturation. Mixing should only be carried out by gently inverting the tube.
- Current data shows that reconstituted DNase I can be stored at 2–8°C for up to 6 weeks. For long-term storage of DNase I, remove the stock solution from the glass vial, divide it into single-use aliquots (use the 1.5 ml microcentrifuge tubes supplied with the kit; there are enough for 5 aliquots), and store at –20°C for up to 9 months. Thawed aliquots can be stored at 2–8°C for up to 6 weeks. Do not refreeze the aliquots after thawing.
- When reconstituting and aliquoting DNase I, ensure that you follow the guidelines for handling RNA (Appendix A, page 47).
- Install the correct shaker adapter (included with the QIAcube; use the adapter for 2 ml safe-lock tubes, marked with a “2”), and place the shaker rack on top of the adapter.
- Check the waste drawer and discard any contents so that it is empty.
- Install the protocols if not already done for previous runs. Install both “PAXgene Blood RNA Part A” and “PAXgene Blood RNA Part B” protocols. See “Installing protocols on the QIAcube”, page 27.

* Kunitz units are the commonly used units for measuring DNase I, defined as the amount of DNase I that causes an increase in A_{260} of 0.001 per minute per milliliter at 25°C, pH 5.0, with highly polymerized DNA as the substrate (Kunitz, M. (1950) *J. Gen. Physiol.* **33**, 349 and 363).

Procedure

1. **Close the QIAcube door, and switch on the QIAcube with the power switch (see Figure 11, page 28).**

A beeper sounds and the startup screen appears. The instrument automatically performs initialization tests.

2. **Open the QIAcube door, and load the necessary reagents and plasticware into the QIAcube. See "Loading the QIAcube", pages 29–34.**

To save time, loading can be performed during one or both of the following 10-minute centrifugation steps (steps 3 and 5).

3. **Centrifuge the PAXgene Blood RNA Tube for 10 minutes at 3000–5000 x g using a swing-out rotor.**

Note: Ensure that the blood sample has been incubated in the PAXgene Blood RNA Tube for a minimum of 2 hours at room temperature (15–25°C), in order to achieve complete lysis of blood cells.

Note: The rotor must contain tube adapters for round-bottom tubes. If other types of tube adapter are used, the tubes may break during centrifugation.

4. **Remove the supernatant by decanting or pipetting. Add 4 ml RNase-free water to the pellet, and close the tube using a fresh secondary BD Hemogard closure (supplied with the kit).**

If the supernatant is decanted, take care not to disturb the pellet, and dry the rim of the tube with a clean paper towel.

5. **Vortex until the pellet is visibly dissolved, and centrifuge for 10 minutes at 3000–5000 x g using a swing-out rotor. Remove and discard the entire supernatant.**

Small debris remaining in the supernatant after vortexing but before centrifugation will not affect the procedure.

Note: Incomplete removal of the supernatant will inhibit lysis and dilute the lysate, and therefore affect the conditions for binding RNA to the PAXgene membrane.

6. **Add 350 µl resuspension buffer (Buffer BR1), and vortex until the pellet is visibly dissolved.**

7. **Pipet the sample into a 2 ml processing tube.**

Note: Use the 2 ml processing tubes included in the PAXgene Blood RNA Kit.

8. **Load the open 2 ml processing tubes containing sample into the QIAcube shaker (see Figure 13, page 30, and Figure 17, page 33). The sample positions are numbered for ease of loading. Insert shaker rack plugs (included with the QIAcube) into the slots at the edge of the shaker rack next to each processing tube. This enables detection of samples during the load check.**

Note: Make sure that the correct shaker adapter (Shaker Adapter, 2 ml, safe-lock tubes, marked with a "2", included with the QIAcube) is installed.

Note: If processing fewer than 12 samples, make sure to load the shaker rack as shown in Figure 17, page 33. One or 11 samples cannot be processed. The position numbers in the shaker rack correspond to the position numbers in the centrifuge.

9. **Close the QIAcube instrument door (see Figure 11, page 28).**
10. **Select the "PAXgene Blood RNA Part A" protocol, and start the protocol.**

Follow the instructions given on the QIAcube touchscreen.

Note: Make sure that both program parts (part A and part B) are installed on the QIAcube instrument (see "Installing protocols on the QIAcube", page 27).

Note: The QIAcube will perform load checks for samples, tips, rotor adapters, and reagent bottles.

11. **After the "PAXgene Blood RNA Part A" protocol is finished, as indicated by a display message, open the QIAcube instrument door (see Figure 11, page 28). Remove and discard the PAXgene RNA spin columns from rotor adapters and the empty processing tubes from the shaker.**

Note: During the run, spin columns are transferred from the rotor adapter position 1 (lid position L1) to rotor adapter position 3 (lid position L2) by the instrument (see Figure 15, page 32).

12. **Close the lids of all 1.5 ml microcentrifuge tubes containing the purified RNA in the rotor adapters (position 3, lid position L3, see Figure 15, page 32). Transfer the 1.5 ml microcentrifuge tubes into the corresponding positions in the QIAcube shaker adapter (see Figure 13, page 30, and Figure 17, page 33).**
13. **Close the QIAcube instrument door (see Figure 11, page 28).**
14. **Select the "PAXgene Blood RNA Part B" protocol, and start the protocol.**

Follow the instructions given on the QIAcube touchscreen.

Note: This program incubates the samples at 65°C and denatures the RNA for downstream applications. Even if the downstream application includes a heat denaturation step, do not omit this step. Sufficient RNA denaturation is essential for maximum efficiency in downstream applications.

15. After the “PAXgene Blood RNA Part B” program is finished, as indicated by a display message, open the QIAcube instrument door (see Figure 11, page 28). Immediately place the microcentrifuge tubes containing the purified RNA on ice.

WARNING**Hot surface**

The shaker can reach temperatures of up to 70°C (158°F).
Avoid touching it when it is hot.

Note: Do not let the purified RNA remain in the QIAcube. Since the samples are not cooled, the purified RNA can be degraded. Unattended overnight sample preparation runs are therefore not recommended.

16. If the RNA samples will not be used immediately, store at -20°C or -70°C .

Since the RNA remains denatured after repeated freezing and thawing, it is not necessary to repeat the heat incubation protocol (“PAXgene Blood RNA Part B”).

Note: For quantification in Tris buffer, use the relationship $A_{260} = 1 \rightarrow 44 \mu\text{g}/\text{ml}$. See Appendix B, page 48.

17. Remove the reagent bottle rack from QIAcube worktable (see Figure 13, page 30), and close all bottles with the appropriately labeled lids. Buffer in bottles can be stored at room temperature (15–25°C) for up to 3 months. Remove and discard remaining reagents in the processing tubes in the QIAcube microcentrifuge tube slots (see Figure 13, page 30). Remove and discard rotor adapters from the centrifuge (see Figure 13, page 30). Empty the QIAcube waste drawer (see Figure 11, page 28). Close the QIAcube instrument door, and switch off the instrument with the power switch (see Figure 11, page 28).

Troubleshooting Guide

This troubleshooting guide may be helpful in explaining any questions that may arise. The scientists in QIAGEN Technical Services are always happy to answer any questions you may have about the information and the manual and automated protocols in this handbook (for contact information, see page 55 or visit www.qiagen.com).

Comments and suggestions

RNA degraded

RNase contamination Be careful not to introduce any RNases into the reagents during the procedure or later handling (see Appendix A, page 47).

Low RNA yield

- a) Less than 2.5 ml blood collected in the PAXgene Blood RNA Tube Ensure that 2.5 ml blood is collected in the PAXgene Blood RNA Tube (see *PAXgene Blood RNA Tube Product Circular*).
- b) RNA concentration measured in water RNA concentration must be measured in 10 mM Tris-Cl, pH 7.5 for accurate quantification (Appendix B, page 48).
- c) Cell debris transferred to the PAXgene RNA spin column in steps 9 and 10 of the manual protocol Avoid transferring large particles when pipetting the supernatant in step 7 of the manual protocol. (Transfer of small debris will not affect the procedure.)
- d) Supernatant not completely removed in step 3 Ensure the entire supernatant is removed. If the supernatant is decanted, remove drops from the rim of the tube by dabbing onto a paper towel. Take appropriate precautions to prevent cross-contamination.
- e) After collection in the PAXgene Blood RNA Tube, blood is incubated for less than 2 hours Incubate blood in the PAXgene Blood RNA Tube for at least 2 hours after collection.

Low A_{260}/A_{280} ratio

- | | |
|---|---|
| a) RNA diluted in water before purity is measured | Use 10 mM Tris-Cl, pH 7.5 to dilute RNA before measuring purity* (Appendix B, page 48). |
| b) Spectrophotometer not properly zeroed | To zero the spectrophotometer, use a blank containing the same proportion of elution and dilution buffers as in the samples to be measured. Components in buffers may affect the A_{260}/A_{280} ratio. |

Instrument malfunction

- | | |
|-------------------------------|--|
| QIAcube not properly operated | Read the <i>QIAcube User Manual</i> , paying careful attention to the Troubleshooting section. Make sure that the QIAcube is properly maintained, as described in the <i>QIAcube User Manual</i> . |
|-------------------------------|--|

* Wilfinger, W.W., Mackey, M., and Chomczynski, P. (1997) Effect of pH and ionic strength on the spectrophotometric assessment of nucleic acid purity. *BioTechniques* **22**, 474.

Appendix A: General Remarks on Handling RNA

Handling RNA

Ribonucleases (RNases) are very stable and active enzymes that generally do not require cofactors to function. Since RNases are difficult to inactivate and even minute amounts are sufficient to destroy RNA, do not use any plasticware or glassware without first eliminating possible RNase contamination. Great care should be taken to avoid inadvertently introducing RNases into the RNA sample during or after the purification procedure. In order to create and maintain an RNase-free environment, precautions must be taken during pretreatment and use of disposable and non-disposable vessels and solutions while working with RNA.

General handling

Proper microbiological, aseptic technique should always be used when working with RNA. Hands and dust particles carry bacteria and molds and are the most common sources of RNase contamination. Always wear latex or vinyl gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin or from dusty laboratory equipment. Change gloves frequently and keep tubes closed whenever possible. Keep purified RNA on ice when aliquots are pipetted for downstream applications.

Protocols for removing RNase-contamination from glassware and solutions can be found in general molecular biology guides, such as Sambrook, J. and Russell, D. W. (2001) *Molecular Cloning: A Laboratory Manual*, 3rd ed. Cold Spring Harbor, NY: Cold Spring Harbor Laboratory Press.

Appendix B: Quantification and Determination of Quality of Total RNA

Quantification of RNA

The concentration of RNA should be determined by measuring the absorbance at 260 nm (A_{260}) in a spectrophotometer. To ensure significance, readings should be in the linear range of the spectrophotometer. An absorbance of 1 unit at 260 nm corresponds to 44 μg of RNA per ml ($A_{260} = 1 \Rightarrow 44 \mu\text{g/ml}$). This relation is valid only for measurements in 10 mM Tris-Cl, * pH 7.5. Therefore, if it is necessary to dilute the RNA sample, this should be done in 10 mM Tris-Cl. As discussed below (see "Purity of RNA"), the ratio between the absorbance values at 260 and 280 nm gives an estimate of RNA purity.

When measuring RNA samples, be certain that cuvettes are RNase-free. Use the buffer in which the RNA is diluted to zero the spectrophotometer, and make sure to add the same volume of Buffer BR5 as the volume of eluted RNA to be diluted.

An example of the calculation involved in RNA quantification is shown below:

Volume of RNA sample = 80 μl

Dilution = 10 μl of RNA sample + 140 μl 10 mM Tris-Cl, pH 7.5 (1/15 dilution)

Measure absorbance of diluted sample in a cuvette (RNase-free).

$A_{260} = 0.3$

Concentration of RNA sample = $44 \times A_{260} \times \text{dilution factor}$

= $44 \times 0.3 \times 15$

= 198 $\mu\text{g/ml}$

Total yield = concentration x volume of sample in milliliters

= $198 \mu\text{g/ml} \times 0.08 \text{ ml}$

= 15.8 μg RNA

* When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, consult the appropriate material safety data sheets (MSDSs), available from the product supplier.

Purity of RNA

The ratio of the readings at 260 nm and 280 nm (A_{260}/A_{280}) provides an estimate of the purity of RNA with respect to contaminants that absorb in the UV, such as protein. However, the A_{260}/A_{280} ratio is influenced considerably by pH. Lower pH results in a lower A_{260}/A_{280} ratio and reduced sensitivity to protein contamination.* For accurate values, we recommend measuring absorbance in 10 mM Tris-Cl, pH 7.5. Pure RNA has an A_{260}/A_{280} ratio of 1.8–2.2 in 10 mM Tris-Cl, pH 7.5. Always calibrate the spectrophotometer with the same solution.

* Wilfinger, W.W., Mackey, M., and Chomczynski, P. (1997) Effect of pH and ionic strength on the spectrophotometric assessment of nucleic acid purity. *BioTechniques* **22**, 474.

Appendix C: Handling PAXgene Blood RNA Tubes

The following recommendations from BD may be helpful when handling PAXgene Blood RNA Tubes. See the *PAXgene Blood RNA Tube Product Circular* for more information about PAXgene Blood RNA Tubes.

Instructions for removal of BD Hemogard Closure

1. Grasp the PAXgene Blood RNA Tube with one hand, placing the thumb under the BD Hemogard closure. (For added stability, place arm on solid surface.) With the other hand, twist the BD Hemogard closure while simultaneously pushing up with the thumb of the other hand **ONLY UNTIL THE TUBE STOPPER IS LOOSENED**.
2. Move thumb away before lifting closure. **DO NOT** use thumb to push closure off tube. Caution: If the tube contains blood, an exposure hazard exists. To help prevent injury during closure removal, it is important that the thumb used to push upward on the closure be removed from contact with the tube as soon as the BD Hemogard closure is loosened.
3. Lift closure off tube. In the unlikely event of the plastic shield separating from the rubber stopper, **DO NOT REASSEMBLE CLOSURE**. Carefully remove rubber stopper from tube.

Instructions for insertion of Secondary BD Hemogard Closure

1. Replace closure over tube.
2. Twist and push down firmly until stopper is fully reseated. Complete reinsertion of the stopper is necessary for the closure to remain securely on the tube during handling.

Ordering Information

Product	Contents	Cat. no.
PAXgene Blood RNA System		
Products that can be ordered from QIAGEN		
PAXgene Blood RNA Kit (50)	50 PAXgene RNA Spin Columns, 50 PAXgene Shredder Spin Columns, Processing Tubes, Microcentrifuge Tubes, Secondary BD Hemogard Closures, RNase-Free DNase Set, Proteinase K, RNase-Free Reagents and Buffers. To be used in conjunction with the PAXgene Blood RNA Tubes	762164*
QIAcube (110 V)	Robotic workstation for automated purification of DNA, RNA, or proteins using QIAGEN spin-column kits, 1-year warranty on parts and labor	9001292
Warranty PLUS 2 Full, QIAcube	3-year warranty, 48-hour (2 working days) priority response, all labor, travel, and repair parts	9240834
Starter Pack, QIAcube	Pack includes: reagent bottle racks (3); rack labeling strips (8); 200 µl filter-tips (1024); 1000 µl filter-tips (1024); 1000 µl filter-tips, wide-bore (1024); 30 ml reagent bottles (18); rotor adapters (240); rotor adapter holder	990395
Filter-Tips, 1000 µl (1024)	Sterile, Disposable Filter-Tips, racked;	990352
Reagent Bottles, 30 ml (6)	Reagent Bottles (30 ml) with lids; pack of 6; for use with the QIAcube reagent bottle rack	990393
Rotor Adapters (10 x 24)	For 240 preps: 240 Disposable Rotor Adapters; for use with the QIAcube	990394

* USA only.

Ordering Information

Product	Contents	Cat. no.
Reagent Bottle Rack	Rack for accommodating 6 x 30 ml reagent bottles on the QIAcube worktable	990390
Rotor Adapter Holder	Holder for 12 disposable rotor adapters; for use with the QIAcube	990392
Products that can be ordered from BD and BD authorized distributors*		
PAXgene Blood RNA Tubes (100)	100 Blood Collection Tubes. To be used in conjunction with the PAXgene Blood RNA Kit (50)	762165
Blood Collection Set	BD Vacutainer® Safety-Lok™ Blood Collection Set: 21G, 0.75 inch needle, 12 inch tubing with luer adapter; 50 per box, 200 per case	367281
BD Vacutainer One-Use Holder	Case only for 13 mm and 16 mmdiameter; 1000/case	364815
BD Vacutainer Plus Serum Tubes	13 x 75 mm 4.0 ml draw with Red BD Hemogard closure and paper label; 100/box, 1000/case	367812

* These blood collection accessories represent typical products that can be used with PAXgene Blood RNA Tubes. To find out more about these accessories, including how to order, visit www.bd.com/vacutainer/products/venous.

This page intentionally left blank

This page intentionally left blank

PreAnalytiX products are distributed in the USA by

QIAGEN Inc. ■ 27220 Turnberry Lane ■ Suite 200 ■ Valencia ■ CA 91355

Orders 888-387-2381 or 800426-8157 ■ Fax 800-718-2056 ■ Technical 888-387-2381 or 800-DNA-PREP (800-362-7737)

www.qiagen.com

BD Diagnostics — Preanalytical Systems ■ 1 Becton Drive ■ Franklin Lakes ■ NJ 07417

Orders 888-237-2762 ■ Fax 800-847-2220 ■ Technical support 800-631-0174

www.bd.com

www.PreAnalytiX.com

1051084 127153018



A QIAGEN / BD Company

