

Vacuum Oven

VAC·LCV



Vacuum ovens and vacuum driers that serve a wide array of uses

ESPEC vacuum ovens are

designed with easy operation on production lines in mind.

They provide a variety of features such as six application-specific operation modes including an Expert Mode that effectively adjusts the depressurization rate during high-volume processing of identical specimens, and a Program Mode that can remember up to nine steps.

Direct-heated vacuum driers speedily perform tasks that use nitrogen such as anaerobic baking, defoaming, hardening and deaerating.

VAC - 100PR



VAC - 200PR



VAC - 300PR



LCV - 233P



Can control depressurization rate adjusting to specimens



Jog dial for fine depressurization control

- **Expert Mode demonstrates its capabilities in repeated high-volume processing**

A jog dial is provided for fine control of the depressurization rate.

The depressurization schedule used is stored and can be called up for subsequent operations to ensure accurate processing. Expert Mode eliminates the fussing with valve controls for each process, and is ideally suited for repeated high-volume processing of identical specimens.

- **Torr-Pa automatic selection function**

Display readings are given in Pa in accordance with the shift toward SI pressure units. Torr-Pa automatic selection function is equipped as a standard; however, the display panel can be switched to the previously used Torr units by touching the button on the panel.



Pa readout



Torr readout

- **Viewing window for full view of specimen**

The viewing window is curved very slightly to eliminate exterior reflections. In Expert Mode, the jog dial can be used to adjust the depressurization rate accordingly while checking the defoaming of the resin inside the oven, preventing problems in the specimen caused by rapid pressure changes.



Curved viewing window

Vacuum control modes to suit a wide range of applications

● Six operation modes to be chosen

The ovens feature a selection of six operating modes, including constant operation, allowing timer-controlled start and stop (on/off); programs 1 and 2, allowing programmed operation of up to nine steps; gas exchange mode; vacuum gradient control mode; and expert mode, allowing repeated high-volume processing of identical specimens. These modes eliminate troublesome operations and use program controls to meet operations to particular applications.

TEMP. & PRESSURE INDICATOR CONTEOLLER

Operation mode	Constant value, Programs 1, 2, Gas exchange, Vacuum gradient control, Expert
Setting range	Temperature: 0 to +200 Pressure: 0 to 1013 × 10 ² Pa Time: 1 min. to 99 hours 59 min.
Setting and indication resolution	Temperature: 1 Pressure: 1 × 10 ² Pa Time: 1 minute
Communications function	E-BUS
Alarm functions	Upper and lower temp. and pressure limit alarm Temp. and pressure sensor disconnection alarm Overheat protector operation alarm Thermal fuse disconnection alarm CPU memory error alarm Motor valve malfunction alarm Vacuum pump thermal relay operation alarm

Six operation modes

Program and Typical pattern	Details	Main applications
<p>Constant value</p>	<ul style="list-style-type: none"> Sets the constant operation temp. and pressure. Timer-controlled setting "on/off" is also possible. 	Vacuum drying
<p>Program 1, 2</p>	<ul style="list-style-type: none"> Two temp. and pressure program patterns can be set in up to nine steps. Pressure decrease and increase time cannot be controlled (ramp operation). 	Vacuum hardening
<p>Gas exchange</p>	<ul style="list-style-type: none"> Gas exchange is performed three times in step one. Temp. cannot be controlled, however. Temp. and pressure can be programmed and controlled from step 2 to step 9. 	Drying in N ₂ gas
<p>Vacuum gradient control</p>	<ul style="list-style-type: none"> The pressure decrease and increase times can be controlled (ramp operation). Temp. can be programmed and controlled in up to step 9. 	Defoaming vacuum drying
<p>Expert</p>	<ul style="list-style-type: none"> The jog dial can be used to control, record, and reproduce depressurization. The depressurization process can be recorded in step 1 only. (Pressure increase cannot be controlled.) Temp. and pressure can be programmed and controlled in up to step 9. 	Defoaming

Excellent temperature uniformity and ease of operation




Inside chamber

- **Double-layered interior construction for great temperature uniformity**

The vacuum chamber interior features a second chamber, forming a double-layered construction. A heater is installed on the exterior of the inner chamber to minimize heat loss and to improve temperature uniformity (± 4 for VAC-300PR). This allows even more uniform heat treatment, and also improves machine efficiency by dramatically reducing the time required to heat up.


Double-skin interior construction



- **Design emphasizing ease of use**

A space at least 110 mm high is provided underneath the stand to allow easy loading and unloading of specimens using a hand lift during high-volume processing. The design also includes ease-of-use features, such as door handles with a recoil-free locking mechanism for smooth opening and closing.

Interior chamber temperature uniformity measurement data



Front

Model: VAC-300R
 Temperature setting: + 200°C
 Pressure: 400 Pa vacuum
 Ambient temperature: + 27°C

*Temperature uniformity for interior center (point 5) and eight interior points with no specimen.

Point	Temp.
1	+ 181.7
2	+ 189.1
3	+ 188.6
4	+ 185.6
5	+ 186.6
6	+ 186.2
7	+ 190.5
8	+ 188.6
9	+ 189.7
Uniformity	± 4

The measurements above are typical values, and are presented for reference only.

- **Equipment range to suit production scale**

Ovens are available in three types, with internal dimensions of 45cm x 45cm, 60cm x 60cm, and 80cm x 80cm. Rotary vacuum pumps are used.

- **Suitable for a wide range of usages**

The ovens are ideal for wide range of applications, particularly in electronic component production processes, including defoaming when mixing silicone rubber or resins in LED production, deaerating during various types of resin forming, hardening when injecting epoxy for hybrid ICs, and drying electronic components after washing.

SPECIFICATIONS

Model	VAC-100PR	VAC-200PR	VAC-300PR
Power supply *1	200V AC 1 2W 50/60Hz, 200V AC 3 3W 50/60Hz 220V AC 1 2W 60Hz, 230V AC 1 2W 50/60Hz 240V AC 1 2W 50/60Hz		200V AC 3 3W 50/60Hz
Maximum power consumption	2.75kVA	3.65kVA	2.75kVA
Pressure control system	Fuzzy control		
Operating temperature	+ 5 to + 35 (+ 41 to + 95°F)		
Temperature control range	+ 40 to + 200 (+ 104 to + 392°F)		
Temperature fluctuations	±0.5 (vacuum), ±1 (atmospheric)		
Temperature heat-up rate *2	Within 50 min.	Within 70 min.	Within 80 min.
Pressure control range	933 to 1 [$\times 10^2$ Pa]		
Attainment pressure *3	Below 133 Pa		
Pressure pull-down rate *3	Within 7 min.	From atmospheric pressure to 133 Pa Within 15 min.	Within 30 min.
Pressure recovery time	Within 4 min.	Inlet open to atmosphere Within 8 min.	Within 15 min.
Construction	Exterior material	Enameled cold-rolled steel plate	
	Vacuum chamber	Stainless steel plate (SUS304)	
	Interior material	Stainless steel plate (NSS430M3)	
	Viewing window	W324 × H336 mm (W13.0 × H13.4 inch)	
	Insulation	Glass wool	
	Heater	Mica heater	
	Inlet	R 1/4 inch, max. operating pressure 0.05 MPa (0.5 kg/cm ² G)	
	Outlet	28 mm external dia. rubber hose connection	
Oil rotary vacuum pump	Motor	200V AC 1 50/60Hz 550W	200V AC 3 50/60Hz 550W
	Design exhaust speed	200L/min. (50Hz), 240L/min. (60Hz)	
	Attainment pressure	6.7×10^{-2} Pa	
	Auxiliary functions	Gas ballast valve, oil mist trap	
Fittings	Adjuster feet and casters (free moving) (× 4 each)		
Effective inside capacity (L)	91	216	512
Effective inside dimensions	W450 × H450 × D450 mm (W18.0 × H18.0 × D18.0 inch)	W600 × H600 × D600 mm (W24.0 × H24.0 × D24.0 inch)	W800 × H800 × D800 mm (W32.0 × H32.0 × D32.0 inch)
Outside dimensions *4	W870 × H1450 × D662 mm (W34.8 × H58.0 × D26.5 inch)	W1020 × H1600 × D812 mm (W40.8 × H64.0 × D32.5 inch)	W1220 × H1800 × D1012 mm (W48.8 × H72.0 × D40.5 inch)
Weight (kg)	320 (328 for 220, 230, 240V)	400 (408 for 220, 230, 240V)	610

*1 Voltage fluctuations within ± 10% of rated voltage.

*2 Time to attain stable temperature the center of chamber with no specimen, under vacuum with ambient temperature of +23°C and temperature setting of +200°C.

*3 Constant temperature inside chamber with no gas generation from specimen.

*4 Excluding protrusions.



DANGER

Do not use specimens which are explosive or flammable, or which contain such substances. To do so could be hazardous, as this may lead to fire or explosion.



CAUTION

Read the User's manual thoroughly prior to use to ensure correct operation of the vacuum pump.

SAFETY DEVICES

Leakage breaker
Overheat protector
Thermal fuse
Thermal relay

ACCESSORIES

Shelves

Stainless steel punched trays 2

VAC-100 PR:

W435 x H13.5 x D435 (mm)

Maximum load: 30kg per shelf
(evenly-distributed load)

VAC-200 PR:

W585 x H13.5 x D585 (mm)

Maximum load: 30kg per shelf
(evenly-distributed load)

VAC-300 PR:

W785 x H13.5 x D785 (mm)

Maximum load: 20kg per shelf
(evenly-distributed load)

* Two shelves are included as standard, and up to five can be fitted.

Total specimen weight must not exceed 100kg.

User's manual 1

Vacuum pump instruction manual 1

OPTIONS

Hermetic terminals for voltage application

Used when applying to specimens voltage.

- Specifications: Hermetic terminal (four-core)
- Max. current: 6 A
- Max. voltage: 200V AC, 250V DC
- Mounted location: Oven rear

* Up to four hermetic connectors can be connected to terminals for voltage application and thermocouples.

Hermetic terminals for thermocouples

Used for connecting to thermocouples from specimens or interior chamber.

Specifications: Hermetic terminal (eight-core, four pairs)

Mounted location: Oven rear

* Up to four hermetic connectors can be connected to terminals for thermo-couples and voltage application.



for voltage application

for thermocouples

Terminal for recorder

Output interior temperature and pressure via 1 to 5V DC linear output.

Temperature: +20 to +220, 1V to 5V

- Pressure: 0 to 106.7 kPa, 1V to 5V
- Mounted location: Oven rear (above inlet)

Pirani vacuum gauge

Pressure is displayed digitally, while this gauge is used to measure pressure accurately below 2,700 Pa.

- Measuring range: 0.4 to 2,700 Pa
- Measuring accuracy: within $\pm 3\%$ of full-scale 100% equivalent on linear scale

* The temperature and pressure recorder or the paperless recorder cannot be fitted to the VAC-100PR if the Pirani vacuum gauge is installed.



Paperless recorder

Records temperature and pressure inside the chamber. Additional inputs may also be recorded.

Temperature range: +20 to +220

Pressure range: 0 to 106.8kPa

Number of inputs (Initial setting):

Temperature 1

Pressure 1

(4 more channels can be turned ON)

Data saving cycle: 5 sec

External recording media:

CF memory card (32MB)

Language Support: ENG, JPN

* The Pirani vacuum gauge cannot be fitted to the VAC-100PR if the paperless recorder is installed.

Temperature and pressure recorder

Records the interior oven temperature and pressure.

- Temperature range: +20 to +220
- Pressure range: 0 to 106.7 kPa
- Inputs: Temperature ($\times 1$), Pressure ($\times 1$)
- Recording method: Intermittent recording

* The Pirani vacuum gauge cannot be fitted to the VAC-100PR if the temperature and pressure recorder is installed.

OPTIONS

Integrating hour meter

Indicates the total integrated operating time.

This is used as a guide for time recording during continuous operation, as well as for maintenance and inspection timing.

- Mounting location: Bottom of operating panel

External alarm terminal

If the safety device of the chamber activates, an error is notified to a distance via the external alarm terminal.

- Power capacity: 250V AC, 3A
- Operation: Connection output when error occurs (closed)
- Mounted location: Oven rear (above inlet)

Signal tower

Illuminates to indicate errors when the safety device activates.

- Color: Red
- Mounting location: Top panel

Inlet filter

Filters the air drawn into the depressurized interior.

- Pore size: 0.2 μm
- Max. pressure: 411.9 kPa (4.2 kg/cm²)
- Connector: NPT 1/8, male screw
- Mounting location: Inlet

Cold trap

Cools and removes moisture and organic solvents contained in the outlet air before being drawn into the vacuum pump.

(Separate from oven)

- Outside dimensions: W306 x H700 x D355mm

Vacuum pump outlet port

Vents gas from the vacuum pump externally.

- Outside connector: NW25 (ISO standard)
- Connection: Quick coupling Center ring with O-ring (not provided)
- Mounting location: Shelf rear

Vacuum pump oil (one-liter can)

Used when maintaining the vacuum pump.

Shelves

Stainless steel punched trays

* Up to five can be fitted inside the oven.



E-BUS cable

- 5, 10m

Power cord

Length from oven: 5 and 10m (two extra cords provided)

* The standard cord provided is 2.5 m from the oven.



LCV - 233P



Direct heating system for fast vacuum-dry.

In addition to the gas exchange function, it can treat specimens in oxygen-free atmospheres using nitrogen or other gases, and supports baking, degassing, hardening, deaeration and numerous other applications.


Easy operation

Temperature setting and upper/lower temperature limit alarm setting can be done with simple key operation.

TEMP. PROGRAM INDICATOR CONTROLLER

Operation mode	Program operation, Constant operation
Program capacity	9 steps / 1 pattern (Number of repetition: 1 to 99)
Setting and indication ranges	Temperature : 0 to + 215 Time : 0 to 99hours 59min., 100 to 999hours
Setting and indication resolution	Temperature : 1 Time : 1min.
Input	Thermocouple type K (Nickel-Chromium/ Nickel-Aluminum)
Control	PID control
Auxiliary functions	Input burn-out detection function Upper and lower temp. limit alarm function Self-diagnostic function (Watchdog timer) Alarm indication function Power failure protection function Timer function (automatic start/ stop)

Temperature heat-up rate (Example)



* Measurement results above are shown as an example.

SPECIFICATIONS

Model	LCV-233	LCV-243	
System	Direct PID control		
Vacuum control	Manual LEAK-VACUUM balance system		
Power supply	AC200V 3 50/60Hz		
Maximum current	8A	9A	
Performance *1	Temperature range	(Ambient + 20) to + 200 (± 392°F)	
	Pressure range	0 to - 101kPa (Gauge)	
	Temperature fluctuation	± 1.0 (± 1.8°F)	
	Temperature heat-up rate	Ambient temperature to + 200 (+ 392°F) 70 min.	110 min.
Construction	External material	Painted steel (melamine coating)	
	Internal material	18-8 Cr-Ni stainless steel plate (2B polish)	
	Viewing window	Tempered glass	
	Vacuum gauge	Bourdon tube vacuum gauge	
Heater	Mica heater		
Capacity	90L	165L	
Inside dimensions *2	W450 × H450 × D450 mm (W17.7 × H17.7 × D17.7 in.)	W550 × H550 × D550 mm (W21.7 × H21.7 × D21.7 in.)	
Outside dimensions *2	W670 × H890 × D700 mm (W26.4 × H35.0 × D27.6 in.)	W770 × H990 × D800 mm (W30.3 × H39.0 × D31.5 in.)	
Weight	170kg	250kg	

*1 Figures for an ambient temperature of + 23 with no specimen in the chamber.

The performance is according to JTM K 05-2000 of Japan Testing Machinery Association.

*2 Excluding protrusions

A separate type transformer for voltage modification is available upon request.

Vacuum Oven with vacuum pump

(Specification for Vacuum Oven is the same as stated above.)

Model	LCV-233P	LCV-243P
Vacuum pump performance *	Direct coupled oil - rotating vacuum pump 6.7 × 10 ⁻² Pa (abs) with gas ballast valve closed 0.67Pa (abs) with gas ballast valve open	
Power supply	AC200V 3 50/60Hz	
Discharge speed *	253/ 309L/ min.	
Outside dimensions *2	W670 × H1540 × D700 mm (W26.4 × H60.6 × D27.6 in.)	W770 × H1640 × D800 mm (W30.3 × H64.6 × D31.5 in.)
Weight	240kg	320kg

* Individual performance rate of vacuum pump.

SAFETY DEVICES

Leakage breaker for power supply	Overheat protector (independent type)
Thermal fuse	Upper and lower temperature limit alarms
Watchdog timer	Sensor burn-out detection circuit

ACCESSORIES

Shel/ Shelf bracket (Stainless steel plate)	5 sets
User's manual	1 set

OPTIONS

Hermetic terminal

The terminals are used to apply voltage to specimen inside chamber and to measure in-chamber temperatures.

- for thermocouple 8P (× 4 pairs)
- for voltage impression 4P


Reverse flow prevention valve

The valve prevents lubricating oil inside vacuum pump from reverse flow when chamber is vacuum state.

* LCV-233P, 243P models only.

Shelf, Shelf bracket

Standard specification shelves and shelf brackets are added as required.



Chamber stand

The stand is equipped with casters enabling chamber to move.

* LCV-233, 243 models only.

* Standard equipment in LCV-233P, 243P models.

Communication function

Connects chamber to a PC, enabling operation control of the chamber.

- RS-485
- GPIB
- RS-232C
- E-BUS

Communication cable

- RS-485 cable (5, 10 m)
- GPIB cable (2, 4 m)
- RS-232C cable (1.5, 3, 5 m)
- E-BUS cable (5, 10 m)

Special spare parts

In addition to standard parts, control parts are installed on request.

ESPEC CORP. <http://www.espec.co.jp/english>

Head Office

3-5-6, Tenjinbashi, Kita-ku, Osaka 530-8550, Japan
Tel: 81-6-6358-4741 Fax: 81-6-6358-5500

ESPEC NORTH AMERICA, INC.

Tel: 1-616-896-6100 Fax: 1-616-896-6150

ESPEC EVALUATION & TEST SYSTEMS, INC.

Tel: 1-408-592-4059 Fax: 1-408-778-4353

ESPEC EUROPE GmbH

Tel: 49-89-18939630 Fax: 49-89-189396379

ESPEC ENVIRONMENTAL EQUIPMENT (SHANGHAI) CO., LTD.

Head Office

Tel: 86-21-51036677 Fax: 86-21-63372237

BEIJING Rep. Office

Tel: 86-10-64627025 Fax: 86-10-64627036

GUANGZHOU Rep. Office

Tel: 86-20-83317826 Fax: 86-20-83317825

SHENZHEN Rep. Office

Tel: 86-755-83674422 Fax: 86-755-83674228

SUZHOU Rep. Office

Tel: 86-512-68664007 Fax: 86-512-68601994

WUXI Rep. Office

Tel: 86-510-82735036 Fax: 86-510-82735039

ESPEC TEST TECHNOLOGY (SHANGHAI) CO., LTD.

Tel: 86-21-68798008 Fax: 86-21-68798088

ESPEC (MALAYSIA) SDN. BHD.

Tel: 60-3-89451377 Fax: 60-3-89451287



QMS
JIS Q 9001:2000
JSAQ 004



JAB
QMS Accreditation
R001



JAB
EMS Accreditation
RE 009



ISO 9001/JIS Q 9001

Quality Management System Assessed and Registered

ESPEC CORP. has been assessed by and registered in the Quality Management System based on the International Standard ISO 9001:2000 (JIS Q 9001:2000) through the Japanese Standards Association (JSA).

ISO 14001 (JIS Q 14001)

Environmental Management System Assessed and Registered

ESPEC GROUP ESPEC CORP.
ESPEC ENGINEERING CORP.
ESPEC TEST CENTER CORP.

Specifications are subject to change without notice due to design improvements.

Corporate names and trade names mentioned in this catalog are trademarks or registered trademarks.

Printed on recycled paper.