


User's Manual

Chart Projector GCP-7000



IMPORTANT NOTICE

This product may malfunction due to electromagnetic waves caused by portable personal telephones, transceivers, radio-controlled toys, etc. Be sure to avoid having objects such as, which affect this product, brought near the product.

The information in this publication has been carefully checked and is believed to be entirely accurate at the time of publication. GILRAS assumes no responsibility, however, for possible errors or omissions, or for any consequences resulting from the use of the information contained herein.

GILRAS reserves the right to make changes in its products or product specifications at any time and without prior notice, and is not required to update this documentation to reflect such changes.

Ver. 1.0

Gilras, 9990 NW 14th St. Unit 105 Miami, FL 33172, USA

Phone: +1 305 722-7321, Fax: +1 786 664-3347

info@gilras.com, www.gilras.com

All rights are reserved. Under copyright laws, this manual may not be copied, in whole or in part, without the prior written consent of GILRAS Co., Ltd.

CONTENTS

1. Introduction	4
1.1. Outline of product	4
1.2. Classifications.....	4
2. Safety Information	5
2.1. Introduction.....	5
2.2. Safety Symbols.....	6
2.3. Environmental factors	7
2.4. Safety Precautions.....	8
3. Features	10
4. Notes for Using the Instrument	11
5. Configurations	12
5.1. Main Body	12
5.2. Main Body cont.....	13
5.3. Accessories	14
6. Installation	15
6.1. Position of Chart Projector	15
6.2. How to obtain correct Letter Size and focus (Zooming Type only)	16
7. Input and Use of Program	17
7.1. How to input user's program	17
7.2. How to use your stored programs.	18

8. Types of Chart (Type C is Standard)	19
9. Maintenance and Repair.....	20
9.1. Lamp Replacement	20
9.2. Fuse Replacement	22
10. Specifications	23
11. Service Information	24

1. Introduction

1.1. Outline of product

GILRAS GCP-7000 is an automatic chart projector which provides 41 different charts including red/green and polarized filters for the eye optometry. Users can select the desired charts instantaneously by cordless remote controller. Charts move quickly and quietly. The customized two programs will help you to perform more convenient and effective examination. You can also execute the binocular balance tests, stereo tests, fixation disparity tests, fusion tests and aniseikonia tests with GCP-7000. This automatic chart project is compatible with GILRAS digital refractor, CDR-7000.

1.2. Classifications

GCP-7000 is classified into Class I system

2. Safety Information

2.1. Introduction

Safety is everyone's responsibility. The safe use of this equipment is largely dependent upon the installer, user, operator, and maintainer. It is imperative that personnel study and become familiar with this entire manual before attempting to install, use, clean, service or adjust this equipment and any associated accessories. It is paramount that the instructions contained in this manual are fully understood and followed to enhance safety to the patient and the user/operator. It is for this reason that the following safety notices have been placed appropriately within the text of this manual to highlight safety related information or information requiring special emphasis. All users, operators, and maintainers must be familiar with and pay particular attention to all Warnings and Cautions incorporated herein.

WARNING

“Warning ” indicates the presence of a hazard that could result in severe personal injury, death or substantial property damage if ignored.

NOTE

“Note” describes information for the installation, operation, or maintenance of which is important but hazard related if ignored.







CAUTION

“Caution” indicates the presence of a hazard that could result in minor injury, or property damaged if ignored.

2.2. Safety Symbols









The International Electrotechnical Commission (IEC) has established a set of symbols for medical electronic equipment, which classify a connection or warn of any potential hazards. The classifications and symbols are shown below.

Save these instructions

	<p>I and O on the power switch represent ON and OFF respectively.</p>
	<p>Signal Input/Output Connection.</p>
	<p>This symbol identifies a safety note. Ensure you understand the function of this control before using it. Control function is described in the appropriate User's or Service Manual.</p>
	<p>Manufactured by xxxx (year).</p>
	<p>Hot surface.</p>
	<p>Identifies the point where the system safety ground is fastened to the chassis. Protective earth connected to conductive parts of Class I equipment for safety purposes.</p>

2.3. Environmental factors

Avoid the following environments for operation or storage:

	<p>Where the equipment is exposed to water vapor. Don't operate an equipment with a wet hand.</p>
	<p>Where the equipment is exposed to direct sunlight.</p>
	<p>Where the temperature changes extremely. Normal operating temperature range is from 10 °C to 40 °C, Humidity is from 30% to 75%.</p>
	<p>Where it is near the heat equipment.</p>
	<p>Where the humidity is extremely high or there is a ventilation problem.</p>
	<p>Where the equipment is subject to excessive shocks or vibrations.</p>
	<p>Where equipment is exposed to chemical material or explosive gas.</p>
	<p>Be careful not to be inserted dust, especially, metal</p>

8 Chart Projector GCP-7000 -----

Avoid places where the ambient temperature falls below 10°C or exceeds 40°C for normal operation, or below -10°C or exceeds 40°C (14°F-104°F) for transportation and storage. Humidity should be maintained between 30 and 75% for normal operation, transportation and storage. Avoid environments where the equipment is exposed to excessive shocks or vibrations.

2.4. Safety Precautions

This equipment has been developed and tested according to national as well as international safety standards. This guarantees a very high degree of safety for this device. The correct handling of this equipment is imperative for its safe operation. Therefore, please read carefully all instructions before switching on this device. For more detailed information, please contact our Customer Service Department or one of our authorized representatives.

1. This equipment must not be used (a) in an area that is in danger of explosions and (b) in the presence of flammable, explosive, or volatile solvents such as alcohol, benzene or similar chemicals.
2. Do not put or use this device in humid rooms. Humidity should be maintained between 30 and 75% for normal operation. Do not expose the device to water splashes, dripping water, or sprayed water. Do not place containers containing fluids, liquids, or gases on top of any electrical equipment or devices.
3. The equipment must be operated only by, or under direct supervision of a properly trained and qualified person.
4. Customer maintenance of this equipment may only be performed as stated in the User's Manual. Any additional maintenance may only be performed by GILRAS service technicians or other authorized persons.

5. The manufacturer is only responsible for effects on safety, reliability, and performance of this equipment when the following requirements are fulfilled: (1) The electrical installation in the respective room corresponds to the specifications stated in this manual and (2) This equipment is used, operated, and maintained according to this manual.
6. The manufacturer is not liable for damage caused by unauthorized tampering with the device(s). Such tampering will forfeit any rights to claim under warranty.
7. Only persons who have undergone proper training and instructions are authorized to install, use, operate, and maintain this equipment.
8. Keep the User's Manual in a place easily accessible at all times for persons operating and maintaining the equipment.
9. Please do not pull on any cable. Always hold on to the plug when disconnecting cables.
10. Before every operation, visually check the equipment for exterior mechanical damage(s) and for proper function.
11. Do not cover any ventilation grids or slits.

10 Chart Projector GCP-7000 -----

-

3. Features

Bright (halogen), high quality images

Quiet, instantaneous chart selection

41 different charts including: red/green and polarizing filters, binocular balance tests, stereo tests, fixation disparity tests, fusion tests and Aniseikonia tests.

34 Masks (vertical line, horizontal line, and single character)

2.5 ~ 8 m (8 ~ 26 feet) projection range

Cordless Remote Control

User Programmable (2 programs up to 30 charts each)

Easily networked with GILRAS's Digital Refractor

4. Notes for Using the Instrument

Do not hit or drop the instrument. The instrument may be damaged if it receives a strong impact. The impact can damage the function of this instrument. Handle it with care.

Exposure to direct sunlight or very bright indoor lights can influence the results of measurements.

If you want to connect this with other equipment, consult the dealer.

Sudden heating of the room in cold areas will cause condensation of vapor on the protective glass in the measurement window and on optical parts inside the instrument. In this case, wait until the condensation disappears before performing measurements.

Keep the objective glass of the examinee side clean. If it is stained with other substances or a fingerprint, it may cause an error or inaccurate measurements.

Disconnect the power supply and consult the dealer when there is smoke, strange odors, or noise while working.

Don't use organic solution such as alcohol, thinner, benzene, etc. to clean the surface of this instrument. It may damage the instrument.


When moving the GCP-7000, always make sure the power supply is off, and then lift the bottom of the unit with both hands.

If you leave the GCP-7000 without using it for long periods, disconnect the power supply and protect the unit with dust cover.

5. Configurations

5.1. Main Body

The GCP-7000 Chart Projector is mainly composed of the body, stand, screen and remote control.




[Fig. 1 Front View with Remote Control]

1. IR receiving window
3. Stand (table mount)

2. Projection lens
4. Wireless remote control


5.2. Main Body cont.



[Fig. 3 Rear View]

-
- 5. Tilting Adjustment
 - 6. Fuse Holder & 110/220 Voltage selector
 - 7. AC Power Inlet
 - 8. Power On/Off Switch
 - 9. CAN (Controller Area Network)

5.3. Accessories




[Fig.4 Accessories]

-
- | | | | |
|--------------------------|------------------|---------------|---------|
| 1. Screen | 2. T wrench(5mm) | 3. lamp | 4. fuse |
| 5. Batteries (2 size AA) | 6. Power cable | 7. Dust Cover | |
| 8. Screen Screw | | | |

6. Installation

6.1. Position of Chart Projector

The single tube projector must be installed at the same distance from the screen as the patient.



[Fig. 6 Installation without mirrors]

1. Install chart projector at desired position within 2.5~8m (8~26 feet).
2. Set the projection angle of chart processor based on Screen.
3. Project 0.05(20/400, 6/120) "C" chart onto the screen.
4. Adjust it for a bright clear image using the focusing wheel.

6.2. How to obtain correct Letter Size and focus (Zooming Type only)

1. Install chart projector at desired position of 2.5~8m (8~26 feet)
2. Measure the distance from patient to screen
3. Refer to the test scale shown in the Appendix A
4. Have someone hold the test scale at the screen
5. Project the 20/400(6/120, 0.05) alphabet "C"
6. Adjust the position of the projector forward or backward and turn the focusing wheel until "C" is brought in sharp focus within the respective distance scale

7. Input and Use of Program

7.1. How to input user's program

- Step 1. Press the "LIGHT" button **once**.
Results: The light goes out.
- Step 2. Press the "PGM 1" button **three** times.
Results: The light will come on and the chart projector will display the first chart of the default program.
- Step 3. Select the first chart for your new program.
Press the "Program ►" button **once**.
- Step 4. Select the second chart for your new program.
Press the "Program ►" button **once**.
- Step 5. Repeat steps 3 & 4 as desired for a total of thirty charts per program.
Note: If you want to check or modify the previously programmed chart, press the "Program ◀" button.
- Step 6. Press the "PGM 1" button **one** time to save your new program.
Caution : If you do not press the PGM1 button, your program will not be saved.

Program No. 2

Repeat Steps 1 through 6 substituting "PGM 2" for "PGM 1".


-

7.2. How to use your stored programs.


- Step 1. Press “ PGM 1 ” or “PGM 2” to select the program of your choice.
- Step 2. Use the “Program ►” or “Program ◀” to scroll through your stored program.

8. Types of Chart (Type C is Standard)


A TYPE



C TYPE



D TYPE



9. Maintenance and Repair


9.1. Lamp Replacement

1. Turn off power and disconnect the power cord.
2. Allow the projector to cool down for about 5 minutes before replacing the lamp.
3. Using the cross screwdriver remove the five screws from the bottom of the cover.
4. Remove the top cover.
5. Loosen the screw-A and remove steel cover
6. Remove and Replace the lamp.
7. Connect the power cord and check the lamp.
8. Reinstall the steel cover



Do not touch the lamp while it is hot. (Risk of burns)

Do not touch the projector with wet hands. (Risk of electric shock)




[Fig. 8 Lamp assembly]

9.2. Fuse Replacement

1. Turn off power, and disconnect the power cord.
2. Pull out the plastic cover at the lower part of power socket, to replace the fuse.



Do not touch the projector with wet hands. (Risk of electric shock)



[Fig. 9 Fuse assembly]

10. Specifications

Charts	41charts -34 mask, Red/green& Polarization Filters
Distance of projection	2.5~8m(Optional: 8m screen)
Chart Rotation Speed	Average 0.15 sec
Chart magnification	30X (When 5m distance)
Power saving	Automatic Lamp Off (10 min)
Program	2 programs with a maximum of 30 charts each
Tilt angle	15 degrees
Power Supply	100-120V~, 200-240~, 50/60 Hz, 0.6A
Lamp	6V/30W(Halogen)
Fuse	5x20mm, 250V, 630mA
Dimension	362x210x265mm
Weight	6.7kg (Body: 5.9kg + Stand: 0.8kg)



11. Service Information

Repair: If the problem is not solved in spite of the settlement according to the contents of chapter 9, please contact to GILRAS's agent with the information on the following items.

- Name of Equipment Type : GCP-7000
- Typical No. of Equipment : Typical number consisted of 8 digits and characters written on its name plate
- Explanation on its symptom : Description in details

Supply of parts required for repair:

- The preservation period of parts required for repair of this machine is by seven(7) years after stopping to produce the product.

Parts to be repaired by qualified service manpower :

- Parts below are consumable in their characteristics, or the quality of them shall be degraded after the long time use. User should not replace them by him or herself. Please contact to GILRAS's agent for the replacement if these parts are consumed enough or degraded by the long time use.
- Back-up battery for clerk and data



CAUTION

As this machine use battery(AA), the reckless abandon of the machine itself or the battery(AA) can cause the environmental pollution. Please contact to the professional waste disposal company.