

**510(k) Summary
for the TeleEMG, LLC
Focus**

MAR - 4 2011

(per 21 CFR 807.92 and <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/default.htm>)

1. SUBMITTER/510(K) HOLDER

TeleEMG, LLC
65 Arlington Road
Woburn, MA 01801, USA

Contact Person: Joe F. Jabre, M.D.
Telephone: 617-840-3253
Date Prepared: August 10, 2010

2. DEVICE NAME

Proprietary Name: Focus EMG Device
Common/Usual Name: Diagnostic Electromyograph
Classification Name: Diagnostic Electromyograph
Product Code: IKN, JXE

3. PREDICATE DEVICES

- Nemus System, K073415
- Synergy LT, K981405

4. DEVICE DESCRIPTION

Physical Description

The Focus is a 2-channel neurodiagnostic testing system designed to enable reliable recording, display and documentation of electrophysiological information from the human nervous and muscular systems in a clinical environment. The system enables the healthcare provider to perform evoked potentials, nerve conductions and needle EMG studies as an aid in the evaluation of patients with central and peripheral nervous system symptoms.

The device provides functionality for the Electromyography (EMG), Nerve Conduction Studies (NCS), and Evoked Potential (EP) testing for a range of clinical applications and has a fast, intuitive and flexible graphical user interface (GUI) that conforms to a Windows philosophy.

The Focus provides the healthcare professional with:

- A means of adjusting or activating all of the unit's settings and controls.
- Multiple ways in which a user can carry a task, such as from a menu, using the mouse or the keyboard or by double-clicking on an icon.
- Configurable function keys and fields that enable the user to customize the interface to meet their requirements.
- Single-stroke keyboard function keys for key tasks carried out during a patient examination.

In addition to using the GUI, key tasks can be carried out via control switches on the Control Panel, Patient Interface Unit, Amplifiers, and Footswitch.

Intended Use

The Focus is intended for use by a healthcare provider to perform nerve conductions and EMG studies as an aid in the evaluation of patients with diseases of muscle and nerves. The machine can also use electrical stimulus or sound stimulus for evoked potentials (EP) studies.

How the Device Functions

The device's principle of operation is based on the recording and input of electrical physiological signals (biopotentials) to a computer with the purpose of its analysis. The functional scheme of the device is represented in the figure below:




Figure 5-1. Device Functional Scheme

EMG biopotentials are amplified by the channels amplifier (CA) and quantized by time and level with the use of the analog-digital converter (ADC) under the control of microcontroller MC1. The commands of the microcontroller MC1 control are transferred from the microcontroller MC4 via the galvanic isolation device (GID), the counts of the input signal are delivered in the reverse direction.

The microcontroller MC1 controls the gain and bandpass of EMG channels amplifiers, is used in the calibration procedures and electrode impedance measurement, and performs the internal diagnostic operations.

The power supply of the amplifier unit is done via the galvanically isolated direct voltage transducer DVT1. The microcontroller MC4 controls the DVT1.

The galvanically isolated electrical stimulator unit includes the electrical stimulus generator (ESG). The microcontroller MC1 specifies the pulse amplitude and duration for the electrical stimulus generator. ESG unit measures the value of the real stimulus current via a patient and transfers it to the microcontroller MC1.

The control of the auditory stimulator, photic stimulator and checkerboard stimulator

is performed via the corresponding auditory, photic, checkerboard pattern generators with the use of the microcontroller MC2 which receives the commands from the computer via MC4.

The microcontroller MC3 performs the displaying of the information to the indicator and entering of information from the buttons, encoders, footswitch and patient button. The given block is located constructively on a separate board.

The microcontroller MC4 also generates the synchro-signal input/output.

The supply unit (SU) transforms the supplying voltage of USB bus to 3.3 V voltage which is required for MC4 operation.

The device operates under control of a PC (IBM PC type) with the mouse, keyboard, laser or jet printer and an installed licensed Windows operational system.

The electronic unit is attached to the USB connector of a computer via the interface cable.

Significant Physical and Performance Characteristics of the Device, such as Device Design, Material Used, and Physical Properties

The Focus EMG Device consists of a handheld electronic unit, a stimulator unit with steel electrodes and felt pads, and a footswitch.























Figure 5-2. The Electronic Unit

Table 5-1. Controls for Handheld Unit

Button/Encoder	Menu Command/Action
	Acquisition Impedance
	Key [Tab] (intended for switching between input and traces review windows)
	Test Next, Test Previous
	Trace Next
	Acquisition Accumulation/averaging

	Acquisition Stimulus protocol (autoincrement)
 <p>Esc</p>	Acquisition Cancel, same as [Esc] key
	Acquisition Monitoring
	Acquisition Repetitive stimulation
	Acquisition Pause
	Acquisition Start acquisition/stimulus
	Acquisition Stop
	Trace Marker Move left and Trace Marker Move right

**Table 5-1. Controls for Handheld Unit
(Continued)**

Button/Encoder	Menu Command/Action
 (pressing)	Selection of next marker
 (rotation)	View Sweep Increase and View Sweep Decrease
 (rotation)	View Sensitivity Increase and View Sensitivity Decrease
 (rotation)	Acquisition Stimulus Increase stimulus duration and Acquisition Stimulus Decrease Stimulus duration
 (rotation)	For tests with stimulation: Acquisition Stimulus Stimulus increase and Acquisition Stimulus Stimulus decrease. For tests without stimulation: Acquisition Sound Volume up and Acquisition Sound Volume down
 (pressing)	For tests without stimulation: Acquisition Sound Sound on/off

The connectors that attach electrodes for measuring potentials are shown in Figure 5-3.




Figure 5-3. Panel View showing Electrode Connectors

View of the panel with connectors for the stimulators is given in Figure 5-4.




Figure 5-4. Panel View showing Stimulator Connectors

The electronic unit measures 92x167x46 mm and weighs .35kg.




Figure 5-5. The Stimulator Unit

The footswitch allows the operator to control the electronic unit. The footswitch measures 103x273x43 mm and weighs 1 kg.



Figure 5-6. Foot Switch


The auditory stimulator is used for auditory Evoked Potentials testing. It is connected to the device via the auditory stimulator connection.



Multiple Transducer Calibration	Supports up to 6 different transducers with saved calibration values
Intensity	0-126 dB SPL, 0-96 dB HL
Polarity	-ve, +ve or alternating
Stimulation Frequency	0.05-30 Hz
Tone Frequency	100-8000 Hz
Click Duration	0.1-6 ms
Pips/Tone Bursts Intensity	0-126 dB SPL, 0-96 dB HL
Presentation	Left / Right / Binaural
Loudspeakers	External (computer)

Figure 5-6. Foot Switch

The LED goggles are used in the “Flash Visual EP Test”. They are connected via the connector for photic or checkerboard stimulation pattern.



Maximum brightness of LED goggles	$(1100 \pm 110) \text{ kd/m}^2$
Flash Duration	1-500 ms
Stimulation Frequency	0.05-50 Hz
Presentation	Left / Right / Both

Figure 5-6. Technical Characteristics of the LED Goggles

5. INTENDED USE

The Focus is intended for use by a healthcare provider to perform nerve conductions and EMG studies as an aid in the evaluation of patients with diseases of muscle and nerves. The machine can also use electrical stimulus or sound stimulus for evoked potentials (EP) studies.

6. SUMMARY OF TECHNOLOGICAL CHARACTERISTICS COMPARED TO THE PREDICATE DEVICE

Table 5-2. Side-by-Side Comparison of the Proposed Device with Cited Predicate Devices

PRODUCT FEATURES	FOCUS EMG DEVICE NEUROSOFT, LTD. PROPOSED	NEMUS SYSTEM EBNEURO S.P.A. K073415	SYNERGY LT OXFORD INSTRUMENTS K981405
Intended Use	The Focus is intended for use by a healthcare provider to perform nerve conductions and EMG studies as an aid in the evaluation of patients with diseases of muscle and nerves. The machine can also use electrical stimulus or sound stimulus for evoked potentials (EP) studies.	The NEMUS system is intended to monitor, record and display the bioelectric signal produced by the muscles, to stimulate peripheral nerves, and to monitor, record and display the electrical activity produced by nerves to aid the clinician in the diagnosis and prognosis of neuro-muscular diseases (EMG). The device may use electrical stimulus or sound stimulus for use in evoked response measurements (EP).	2 channel electromyograph which provides facilities for EMG and Evoked Potentials testing for a range of clinical application. Synergy LT is designed to enable reliable recording display and documentation of electrophysiological information from the human nervous and muscular system in a clinical environment.
Warnings	Items related to off-label use.	Items related to off-label use.	Items related to off-label use.
Contraindications	Items related to design and indicated use limitations, such as not for use in the presence of flammable anesthetics or in conjunction with defibrillation equipment.	Items related to design and indicated use limitations, such as not for use in the presence of flammable anesthetics or in conjunction with defibrillation equipment.	Items related to design and indicated use limitations, such as not for use in the presence of flammable anesthetics
Design - General			
General system approach	Computer based equipment with dedicated hardware peripherals/components	Computer based equipment with dedicated hardware peripherals/components	Computer based equipment with dedicated hardware peripherals/components
User input device	Microsoft Windows mouse/keyboard driven graphic interface. Built-in keyboard	Microsoft Windows mouse keyboard driven graphic interface. Dedicated external keyboard (optional)	Microsoft Windows mouse/keyboard driven graphic interface.
User output device	Digital color display Commercial printers	Digital color display Commercial printers	Digital color display Commercial printers
Patient inputs	2 channels preamplifiers, isolated	2 channels preamplifiers, isolated	2 channels preamplifiers, isolated
Signal acquisition	Analog-to-digital conversion at variable sampling rate	Analog-to-digital conversion at variable sampling rate	Analog-to-digital conversion at variable sampling rate
Trigger input (synchronization to external events)	Yes	Yes	Yes
Trigger output (synchronization for external devices)	Yes	Yes	Yes
Footswitch for hands-free operation	Yes	Yes (Optional)	Yes
Use of standard software platform (Operating System)	YES — Microsoft Windows	YES — Microsoft Windows	YES — Microsoft Windows

Table 5-2. Side-by-Side Comparison of the Proposed Device with Cited Predicate Devices (Continued)

PRODUCT FEATURES	FOCUS EMG DEVICE NEUROSOFT, LTD. PROPOSED	NEMUS SYSTEM EBNEURO S.P.A. K073415	SYNERGY LT OXFORD INSTRUMENTS K981405
Customization of clinical protocols	Via storage / retrieval of user-defined settings	Via storage / retrieval of user-defined settings	Via storage / retrieval of user-defined settings
Application flexibility /expandability	Via software update	Via software update	Via software update
Safety Standards	60601-1-1-2000 60601-1-2-2001 60601-1-2-40-1998 62471-2006	IEC 60601-1 IEC 60601-1-1 IEC 60601-2-26 IEC 60601-2-40 IEC 60601-1-2	IEC 60601-1
Patient circuitry isolation	Optic/transformer	Optic	Optic
System Components	Focus EMG Device USB Cable Footswitch Electrical Stimulator Notebook PC based system	Nemus I Base Unit (Amplifier) Host Computer (PC) Printer (optional) Cart (optional) Dedicated Keyboard (optional)	Headbox (Amplifier) Computer (PC) Cart Printer (optional)
Amplifier- Computer interface	USB	LAN Ethernet 100 Mbit	Proprietary
System Power Supply	Electronic unit: 5V DC System with Personal Computer: 100-120/200-240 V AC (50-60 Hz) System with Notebook: 100-120/200-240 V AC (50-60 Hz)	From mains (110-240 VAC) through isolation transformer	From mains (110-240 VAC)
Amplifier Power Supply	5V DC from USB	15 VDC from medical approved AC/DC converter	Internal power supply
Size (H/W/D) mm	Electronic Unit: 90x270x50 mm	Nemus System: 1250/750/850 (complete system, with cart, monitor, arm)	500/209/321 (cart : 720/640/780)
Weight (complete system-kg)	Electronic Unit: 0.35 System: 3.5	Nemus System: 90 (complete system, with cart, monitor, arm)	55 (with cart)
Design - Acquisition			
Number of channels	2	2	2
CMRR	>100 dB	>100 dB	>110dB
Noise	< 0. 6 μ Vrms.	< 0.3 μ Vrms (0.1 — 00 Hz) < 20 nV/ \sqrt Hz (10 Hz — 10 kHz)	<0.7 μ Vrms (0.1 Hz — 10 kHz)
Input impedance	> 100 MOhm < 25pF	> 1000 MOhm / 8 pF	> 1000 MOhm / 30 pF
Low pass filter	(-12dB/octave) 10, 20, 35, 50, 75, 100, 150, 200, 300, 500, 1000, 2000, 3000, 5000, 10000 Hz	20 Hz — 16 kHz (15 step)	50 Hz — 16 kHz (10 step)
High pass filter	(-6dB/octave) 0.05, 0.1, 0.2, 0.3, 0.5, 1, 2, 3, 5, 10, 20, 30, 50, 100, 200, 300, 500, 1000, 2000, 3000 Hz	0.01 — 500 Hz (11 step)	DC, 30 — 2000 Hz
Notch filter	50/60 Hz selectable	50/60 Hz selectable	50/60 Hz selectable
A/D conversion	16 Bit ADC	24 bit Sigma-Delta	16 bit
Sampling rate	200-80000 Hz	4.194 MHz	50 kHz
Analysis time	2 ms - 5 min	5 ms — 10 s	5 ms — 10s
Time base	Single	Single	Single

Table 5-2. Side-by-Side Comparison of the Proposed Device with Cited Predicate Devices (Continued)

PRODUCT FEATURES	FOCUS EMG DEVICE NEUROSOFT, LTD. PROPOSED	NEMUS SYSTEM EBNEURO S.P.A. K073415	SYNERGY LT OXFORD INSTRUMENTS K981405
Trigger mode	Free, Auto, Internal, External	Free, Auto, Internal, External	Auto, Internal, External
Signal delay (pre/post)	0-10 s	0-5 div	0 - 9 div
Ohmmeter	0-100 kOhm	0-100 kOhm (auto full scale)	0-32 kOhm
Design – Stimulators			
Somatosensory (electrical) Stimulator	Type: constant current N. output 1 Max output: 100 mA Pulse width: 0.1 - 5 ms Mode: single, train (50 Hz)	Type: constant current N. output 1 Max output: 100 mA Pulse width: 0.05 - 1 ms Mode: single, train	Type: constant current/voltage N. output : 1 Max output: 100 mA Pulse width: 0.05-1 ms Mode: single, double, train (200Hz)
Audio Stimulator	Output mode: click, tone Sound pressure: 0-126 dB SPL (TDH 39) Phase: condens., raref, alternate Signal frequency: 1008000 Hz Click width: 100-5000 μ s Stimulus presen. Left, right, binaural Headset: TDH 39	Output mode: click, tone Sound pressure: 0-132 dB SPL Phase: condens., raref, alternate Signal frequency: 1258000 Hz Plateau time: 1-200 ms Rise/fall time: 1-100 ms Mask level: -40 + +10 dB (relative) Click width: 1-100 μ s Stimulus presen. Left, right, binaural Headset: TDH 39	Output mode: click, pip, tone Sound pressure: 0-122 dB SPL Phase: condens., raref, alternate Signal frequency: 1258000 Hz Plateau time: 1-999 ms Rise/fall time: 1-255 ms Mask level: -40 + +10 dB (relative) Click width: 0.05 - 1 ms Stimulus presen. Left, right, binaural Headset: TDH 39
Basic EMG Application Modules			
Free run acquisition	Yes	Yes	Yes
Nerve conduction study	Yes	Yes	Yes
Self triggered acquisition	Yes	Yes	Yes
Spontaneous activity	Yes	Yes	Yes
Single fiber EMG	Yes	Yes	Yes
Motor Unit Analysis	Yes	Yes	Yes
F wave analysis	Yes	Yes	Yes
H-Reflex analysis	Yes	No	Yes
Galvanic Skin Response Also known as Sympathetic Skin Response (SSR)	Yes	No	Yes (Referred to as Sympathetic Skin Response or SSR)
Basic EP Application Modules			
Somatosensory EP	Yes	Yes	Yes
Auditory EP	Yes	Yes	Yes
Visual EP (flash)	Yes	No	Yes
P300 EP	Yes	No	Yes
Visual EP (pattern)	Yes	No	Yes

7. SUMMARY OF NON-CLINICAL PERFORMANCE TESTING AS BASIS FOR SUBSTANTIAL EQUIVALENCE

Performance Testing

Performance evaluation of the features described in the Focus user's manual has been successfully completed utilizing hardware and software tests and validations. Hardware qualification is performed using the following industry standards:

- IEC 60601-1-1:2000 Medical electric equipment Part 1-1: Safety requirements for medical electrical systems
- IEC 60601-1-2:2001 Medical electric equipment Part 1-2: General requirements for safety – Collateral standard: Electromagnetic compatibility –Requirements and tests
- IEC 60601-2-40:1998 Medical electric equipment Part 2-40: Particular requirements for the safety of electromyographs and evoked response equipment
- IEC 62471:2006 The TeleEMG Focus EMG Device has met the IEC 62471:2006: Photobiological safety of lamps and lamp systems. According to IEC 62471, our LED goggles luminance is 1,203 cd/m² that is much less than the admissible level 10,000 cd/m² specified in clause 4.1 of this standard. That is why it is not required to make additional spectral calculations according to IEC 62471.

The TeleEMG Focus is of a "moderate" level of concern according to the FDA draft software guidance document ("Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices"), and meets the requirements set forth in the FDA draft software guidance document for a "moderate" in the level of concern device.

Biocompatibility Testing

The only patient contacting surfaces of the Focus EMG device are the stainless steel electrodes and felt pads applied during testing. The electrodes and pads were cleared under the submission of the NC Stat, K060584.

8. SUMMARY OF CLINICAL TESTING AS BASIS FOR SUBSTANTIAL EQUIVALENCE

No clinical testing was conducted to support this submission.

9. SUMMARY OF OTHER INFORMATION

No other information is available.

10. CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL TESTS

Based on the information and supporting documentation provided in the premarket notification, the Focus EMG Device is substantially equivalent to the cited predicate device. Testing demonstrates that the Focus EMG Device fulfills prospectively defined design and performance specifications.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

TeleEMG, LLC
c/o Joe F. Jabre, M.D.
65 Arlington Road
Woburn, MA 01801

MAR - 4 2011

Re: K102610

Trade/Device Name: Focus EMG
Regulation Number: 21 CFR 890.1375
Regulation Name: Diagnostic Electromyograph
Regulatory Class: Class II
Product Code: IKN, JXE
Dated: January 22, 2011
Received: January 24, 2011

Dear Dr. Jabre:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in cursive script that reads "Malvina B. Eydelman". The signature is written in dark ink and is positioned to the left of the typed name and title.

Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic, Neurological,
and Ear, Nose and Throat Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K102610

Device Name: Focus EMG Device

Indications for Use:

The Focus is intended for use by a healthcare provider to perform nerve conductions and EMG studies as an aid in the evaluation of patients with diseases of muscle and nerves. The machine can also use electrical stimulus or sound stimulus for evoked potentials (EP) studies.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

KRISTEN BOWSHER
(Division Sign-Off)
Division of Ophthalmic, Neurological and Ear,
Nose and Throat Devices

510(k) Number K102610