

Factura™

Feature Identification Software

User's Manual

© Copyright 2001, Applied Biosystems. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

Notice to Purchaser: License Disclaimer

Purchase of this product alone does not imply any license under any process, instrument or other apparatus, system, composition, reagent or kit rights under patent claims owned or otherwise controlled by PE Corporation (NY), either expressly, or by estoppel.

ABI PRISM and its design and Applied Biosystems are registered trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries. ABI, AutoAssembler, GeneAssist and Sequence Navigator are trademarks of Applied Biosystems Corporation or its subsidiaries in the U.S. and certain other countries.

All other trademarks are the sole property of their respective owners.

Contents

1 Introduction

About the Factura Software User's Manual	1-3
What is Factura Software?	1-4
Special Text Usage	1-5
Technical Support	1-6

2 Getting Started

Factura Macintosh Installation	2-3
Factura Program Files	2-3
Factura Input and Output Files	2-4
Factura Windows and Views	2-7

3 Working With the Batch Worksheet

Specifying Settings for the Batch Worksheet	3-3
Opening the Batch Worksheet	3-10
Adding and Removing Sequences	3-14
Changing Batch Worksheet Parameters	3-18
Submitting the Batch Worksheet	3-21
Closing the Batch Worksheet	3-25

4 Working With Sequences in Factura

Opening the Sequence Window	4-3
Description of the Sequence Window	4-8
Editing in the Sequence Window	4-14
Viewing Options in the Sequence Window	4-21
Closing the Sequence Window	4-42

5	Printing and Saving in Factura	
	Printing Batch Reports & Electropherograms	5-3
	Saving the Batch Worksheet	5-5
	Saving Results in a Batch Report	5-7
	Saving Results to Sequence Data Source Files	5-9
	Saving Changes Made in the Sequence Window	5-12
	Exporting Results to Other Formats	5-13
6	Setting up Factura Libraries	
	Library Files	6-3
	Setting Up the Libraries	6-5
7	Overview—Factura Menu Commands	
	Factura Main Menu	7-3
	Apple Menu	7-3
	File Menu	7-4
	Edit Menu	7-6
	Library Menu	7-7
	Sequence Menu	7-8
	Worksheet Menu	7-10
	Window Menu	7-12

1 Introduction

This section, *Introduction*, provides a general introduction to Fatura Software. It also provides information about the organization of this manual, about special text usage in the manual, and instructions on how get help from Applied Biosystems.

Contents

About the Factura Software User's Manual	1-3
What is Factura Software?	1-4
Factura (Feature Identification)	1-4
Factura Dataflow	1-4
Factura Setup	1-4
Special Text Usage	1-5
Attention Users	1-5
Technical Support	1-6
Contacting Technical Support	1-6
To Contact Technical Support by E-Mail	1-6
Hours for Telephone Technical Support	1-6
To Contact Technical Support by Telephone or Fax	1-7
To Reach Technical Support Through the Internet	1-11
To Obtain Documents on Demand	1-12

About the Fatura Software User's Manual

This manual provides you with detailed information about the Fatura application program. For information on how to use Fatura, consult Tutorial 6, *Using the Fatura Program* in the *GeneAssist Applications Tutorials* manual.

- This section, *Introduction*, provides a general introduction to Fatura. It also provides information about the organization of this manual, about special text usage in the manual, and instructions on how get help from Applied Biosystems.
- Section 2, *Getting Started*, lists the files provided with the Fatura application. It also describes Fatura input and output files, and Fatura windows and views.
- Section 3, *Working With the Batch Worksheet*, explains the use of the batch worksheet. This section provides information on how to open, close, and set up the batch worksheet, as well as how to add and remove sequences, change the processing parameters, and submit the batch of sequences for processing.
- Section 4, *Working With Sequences in Fatura*, provides descriptions of several ways to open a sequence window. It includes a brief description of the sequence window and its various views and provides procedures for editing and changing the appearance of a sequence in the sequence window.
- Section 5, *Printing and Saving in Fatura*, explains how to save the batch worksheet, view and save results in a batch report, save Fatura-identified features to the individual sequence files, and export the contents of a sequence file into other formats.
- Section 6, *Setting up Fatura Libraries*, explains the difference between Main and Custom libraries and describes how to set up custom libraries of vectors, enzymes, and primers to be used with your Fatura processing.
- Section 7, *Overview—Fatura Menu Commands*, provides tables with brief descriptions of the Fatura main menus along with cross references to other sections that provide more detail.
- The *Fatura Software User's Manual* concludes with an index.

What is Factura Software?

Factura (Feature Identification)

The Factura application enables you to clean up ABI 373 or ABI PRISM 377 or 310 Sequencer sample files or other data prior to assembly or prior to analysis. Factura does this by automatically identifying designated sequence features, such as vector sequences on both ends of the sequence fragment and ambiguous regions at the ends of the fragment. After identification of features, Factura takes advantage of the 373, 377, or 310 sample files (Applied Biosystems automated sequencer integrated data files—ABI 373 or ABI PRISM 377 or 310) to add identified features to sequence files for use by the AutoAssembler and GeneAssist applications.

Factura Dataflow

Factura is used to clean up sample files or other data. Input data often contains vector sequences and might have ambiguously called bases at both ends of the sequence that the user has to remove prior to assembly or analysis. Factura functions by identifying designated features, such as vectors and ambiguous regions, and flagging the features in the sequence file. This allows the AutoAssembler and GeneAssist applications to ignore designated sequence features, effectively cleaning up the data.

Factura Setup

To enable you to clean up data input from ABI 373 or ABI PRISM 377 or 310 sequencer files, Factura allows you to specify information needed to quickly access the location of the vector sequence as well as the information needed to remove ambiguous regions. Information required to locate the vector sequence includes the names of the vector, primer, and the cloning site. Information needed to remove ambiguities includes the specification of the amount of ambiguity acceptable at both the 5' and 3' ends and the number of Ns acceptable in the sequence. You can also specify the valid length of the sequence (possible maximum range of base positions).

Special Text Usage

Attention Users

Four User Attention words appear in the text of this manual. They are designed to draw your attention to safety issues or to issues relevant to proper operation of the instrument or software programs. Each one requires a certain level of observation or action as follows:

<i>Note</i>	<i>Used to call attention to information.</i>
-------------	---

IMPORTANT	<i>Indicates information that is necessary for proper instrument or program operation.</i>
------------------	--

Caution	<i>Damage to the instrument or data could result if you do not comply with this information.</i>
----------------	--

WARNING	<i>Physical injury to you or other people could result if these required precautions are not taken.</i>
----------------	---

Technical Support

Contacting Technical Support

You can contact Applied Biosystems for technical support by telephone or fax, by e-mail, or through the Internet. You can order Applied Biosystems user documents, MSDSs, certificates of analysis, and other related documents 24 hours a day. In addition, you can download documents in PDF format from the Applied Biosystems Web site (please see the section “To Obtain Documents on Demand” following the telephone information below).

To Contact Technical Support by E-Mail

Contact technical support by e-mail for help in the following product areas:

Product Area	E-mail address
Genetic Analysis (DNA Sequencing)	galab@appliedbiosystems.com
Sequence Detection Systems and PCR	pclab@appliedbiosystems.com
Protein Sequencing, Peptide and DNA Synthesis	corelab@appliedbiosystems.com
Biochromatography, PerSeptive DNA, PNA and Peptide Synthesis systems, CytoFluor [®] , FMAT [™] , Voyager [™] , and Mariner [™] Mass Spectrometers	tsupport@appliedbiosystems.com
LC/MS (Applied Biosystems/MDS Sciex)	apisupport@sciex.com or api3-support@sciex.com
Chemiluminescence (Tropix)	tropix@appliedbiosystems.com

Hours for Telephone Technical Support

In the United States and Canada, technical support is available at the following times:

Product	Hours
Chemiluminescence	8:30 a.m. to 5:30 p.m. Eastern Time
Framingham support	8:00 a.m. to 6:00 p.m. Eastern Time
All Other Products	5:30 a.m. to 5:00 p.m. Pacific Time

To Contact Technical Support by Telephone or Fax

In North America

To contact Applied Biosystems Technical Support, use the telephone or fax numbers given below. (To open a service call for other support needs, or in case of an emergency, dial **1-800-831-6844** and press 1.)

Product or Product Area	Telephone Dial...	Fax Dial...
ABI PRISM® 3700 DNA Analyzer	1-800-831-6844 , then press 8	1-650-638-5981
DNA Synthesis	1-800-831-6844 , then press 21	1-650-638-5981
Fluorescent DNA Sequencing	1-800-831-6844 , then press 22	1-650-638-5981
Fluorescent Fragment Analysis (includes GeneScan® applications)	1-800-831-6844 , then press 23	1-650-638-5981
Integrated Thermal Cyclers (ABI PRISM® 877 and Catalyst 800 instruments)	1-800-831-6844 , then press 24	1-650-638-5981
ABI PRISM® 3100 Genetic Analyzer	1-800-831-6844 , then press 26	1-650-638-5981
Bioinformatics (includes BioLIMS®, BioMerge™, and SQL GT™ applications)	1-800-831-6844 , then press 25	1-505-982-7690
Peptide Synthesis (433 and 43X Systems)	1-800-831-6844 , then press 31	1-650-638-5981
Protein Sequencing (Procise® Protein Sequencing Systems)	1-800-831-6844 , then press 32	1-650-638-5981
PCR and Sequence Detection	1-800-762-4001 , then press 1 for PCR, 2 for the 7700 or 5700, 6 for the 6700 or dial 1-800-831-6844 , then press 5	1-240-453-4613

Product or Product Area	Telephone Dial...	Fax Dial...
Voyager™ MALDI-TOF Biospectrometry and Mariner™ ESI-TOF Mass Spectrometry Workstations	1-800-899-5858 , then press 13	1-508-383-7855
Biochromatography (BioCAD® Workstations and Poros® Perfusion Chromatography Products)	1-800-899-5858 , then press 14	1-508-383-7855
Expedite™ Nucleic acid Synthesis Systems	1-800-899-5858 , then press 15	1-508-383-7855
Peptide Synthesis (Pioneer™ and 9050 Plus Peptide Synthesizers)	1-800-899-5858 , then press 15	1-508-383-7855
PNA Custom and Synthesis	1-800-899-5858 , then press 15	1-508-383-7855
FMAT™ 8100 HTS System and Cytofluor® 4000 Fluorescence Plate Reader	1-800-899-5858 , then press 16	1-508-383-7855
Chemiluminescence (Tropix)	1-800-542-2369 (U.S. only), or 1-781-271-0045	1-781-275-8581
Applied Biosystems/MDS Sciex	1-800-952-4716	1-650-638-6223

Outside North America

Region	Telephone Dial...	Fax Dial...
Africa and the Middle East		
Africa (English Speaking) and West Asia (Fairlands, South Africa)	27 11 478 0411	27 11 478 0349
South Africa (Johannesburg)	27 11 478 0411	27 11 478 0349
Middle Eastern Countries and North Africa (Monza, Italia)	39 (0)39 8389 481	39 (0)39 8389 493

Region	Telephone Dial...	Fax Dial...
Eastern Asia, China, Oceania		
Australia (Scoresby, Victoria)	61 3 9730 8600	61 3 9730 8799
China (Beijing)	86 10 64106608	86 10 64106617
Hong Kong	852 2756 6928	852 2756 6968
Korea (Seoul)	82 2 593 6470/6471	82 2 593 6472
Malaysia (Petaling Jaya)	60 3 758 8268	60 3 754 9043
Singapore	65 896 2168	65 896 2147
Taiwan (Taipei Hsien)	886 2 2358 2838	886 2 2358 2839
Thailand (Bangkok)	66 2 719 6405	66 2 319 9788
Europe		
Austria (Wien)	43 (0)1 867 35 75 0	43 (0)1 867 35 75 11
Belgium	32 (0)2 712 5555	32 (0)2 712 5516
Czech Republic and Slovakia (Praha)	420 2 61 222 164	420 2 61 222 168
Denmark (Naerum)	45 45 58 60 00	45 45 58 60 01
Finland (Espoo)	358 (0)9 251 24 250	358 (0)9 251 24 243
France (Paris)	33 (0)1 69 59 85 85	33 (0)1 69 59 85 00
Germany (Weiterstadt)	49 (0) 6150 101 0	49 (0) 6150 101 101
Hungary (Budapest)	36 (0)1 270 8398	36 (0)1 270 8288
Italy (Milano)	39 (0)39 83891	39 (0)39 838 9492
Norway (Oslo)	47 23 12 06 05	47 23 12 05 75
Poland, Lithuania, Latvia, and Estonia (Warszawa)	48 (22) 866 40 10	48 (22) 866 40 20
Portugal (Lisboa)	351 (0)22 605 33 14	351 (0)22 605 33 15
Russia (Moskva)	7 095 935 8888	7 095 564 8787
South East Europe (Zagreb, Croatia)	385 1 34 91 927	385 1 34 91 840
Spain (Tres Cantos)	34 (0)91 806 1210	34 (0)91 806 1206
Sweden (Stockholm)	46 (0)8 619 4400	46 (0)8 619 4401
Switzerland (Rotkreuz)	41 (0)41 799 7777	41 (0)41 790 0676

Region	Telephone Dial...	Fax Dial...
The Netherlands (Nieuwerkerk a/d IJssel)	31 (0)180 331400	31 (0)180 331409
United Kingdom (Warrington, Cheshire)	44 (0)1925 825650	44 (0)1925 282502
All other countries not listed (Warrington, UK)	44 (0)1925 282481	44 (0)1925 282509
Japan		
Japan (Hacchobori, Chuo-Ku, Tokyo)	81 3 5566 6230	81 3 5566 6507
Latin America		
Del.A. Obregon, Mexico	305-670-4350	305-670-4349

To Reach Technical Support Through the Internet

We strongly encourage you to visit our Web site for answers to frequently asked questions and for more information about our products. You can also order technical documents or an index of available documents and have them faxed or e-mailed to you through our site. The Applied Biosystems Web site address is

<http://www.appliedbiosystems.com/techsupp>

To submit technical questions from North America or Europe:

Step	Action
1	Access the Applied Biosystems Technical Support Web site.
2	Under the Troubleshooting heading, click Support Request Forms , then select the relevant support region for the product area of interest.
3	Enter the requested information and your question in the displayed form, then click Ask Us RIGHT NOW (blue button with yellow text).
4	Enter the required information in the next form (if you have not already done so), then click Ask Us RIGHT NOW . You will receive an e-mail reply to your question from one of our technical experts within 24 to 48 hours.

To Obtain Documents on Demand

Free, 24-hour access to Applied Biosystems technical documents, including MSDSs, is available by fax or e-mail or by download from our Web site.

To order documents...	Then...
by index number	a. Access the Applied Biosystems Technical Support Web site at http://www.appliedbiosystems.com/techsupp b. Click the Index link for the document type you want, then find the document you want and record the index number. c. Use the index number when requesting documents following the procedures below.
by phone for fax delivery	a. From the U.S. or Canada, call 1-800-487-6809 , or from outside the U.S. and Canada, call 1-858-712-0317 . b. Follow the voice instructions to order the documents you want. Note There is a limit of five documents per request.
through the Internet for fax or e-mail delivery	a. Access the Applied Biosystems Technical Support Web site at http://www.appliedbiosystems.com/techsupp b. Under Resource Libraries , click the type of document you want. c. Enter or select the requested information in the displayed form, then click Search . d. In the displayed search results, select a check box for the method of delivery for each document that matches your criteria, then click Deliver Selected Documents Now (or click the PDF icon for the document to download it immediately). e. Fill in the information form (if you have not previously done so), then click Deliver Selected Documents Now to submit your order. Note There is a limit of five documents per request for fax delivery but no limit on the number of documents you can order for e-mail delivery.

2 Getting Started

This section provides:

- A list of the files provided with the Factura application
- Descriptions of Factura input and output files
- Descriptions of Factura windows and views

These descriptions do not walk you through the use of the Factura application. See Section 3, *Working With the Batch Worksheet* through Section 6, *Setting up Factura Libraries* for more detail.

Contents

Factura Macintosh Installation	2-3
Factura Program Files	2-3
Factura Input and Output Files	2-4
The Sequence File	2-4
The Batch Worksheet File	2-5
The Batch Report File	2-5
Diagram of File Input And Output	2-5
Factura Windows and Views	2-7
The Batch Worksheet Window	2-7
The Batch Report Window	2-7
The Sequence Window	2-7
Sequence View	2-8
Annotation View	2-9
Feature View	2-10
Electropherogram View	2-11

Factura Macintosh Installation

To install the Factura application and related tutorial files, refer to the *User's Manual* that accompanies this software.

Factura Program Files

Factura is bundled with either GeneAssist, AutoAssembler, or Sequence Navigator. The two applications are installed simultaneously from one set of disks. Figure 2-1 illustrates an example of the folder and file structure following a typical Factura and GeneAssist installation.




Figure 2-1. Folders and files with typical Factura and GeneAssist installation

The primary application (in this case, GeneAssist) folder is installed at the top level on the hard disk. The ABI Folder is installed inside the System folder. For a more information on other Applied Biosystems, Applied Biosystems software applications, see the *User's Manual* that accompanies this software.

Factura Input and Output Files

The Sequence File

A sequence file contains the nucleic acid characters of a sequence and any annotation associated with the sequence.

Factura uses sequence files as input and creates them as output. Sequences processed in the batch worksheet can be created externally or created within the Factura application. Factura accepts several different types of sequence files as input and can create the same types as output:

- *Sample* files are data files created by the ABI 373 Sequencing Analysis application. They contain base calls, peak locations, electropherograms, and other information. After a sequence is processed in Factura, you can save the identified features (vector, ambiguity, confidence range, heterozygote positions, and so on) to a feature table in the sample file.

Note *The original ABI 373/377/310-produced sequence data is maintained in its unmodified state in a sample file. A copy of the data is stored in the file as editable data. When you save a sequence, only the editable data is changed. The Factura sequence windows display the editable data. The Settings dialog box in Factura allows you to revert sequences to the original data for feature identification. If you choose to revert them, the editable data is overwritten with original data, and any editing you have performed is lost.*

- *Text* files contain a string of characters. An input sequence can be a standard word processing text file or any of the standard text formats for sequences, such as Staden or GCG. You can also type a new text sequence manually in Factura application. Sequences output by Factura as text files are created in this format for easy export into other applications.

Individual sequence files can be viewed on the screen in the sequence window (see page 2-7).

The Batch Worksheet File

The batch worksheet is the main file created by Factura. It summarizes a batch of sequences and the feature identification parameters you have applied to them using the application. You can save this information in a batch worksheet file, and import the batch of sequences into another application, such as the Sequence Navigator application.

When you open a batch worksheet file, the contents are displayed on the screen in a batch worksheet window (described on page 2-7). Information about working with a batch worksheet is provided in Section 3.

The Batch Report File

The batch report is created by Factura to summarize the results of processing. For each sequence the report describes the identified features (vector range, ambiguity range, confidence range), and details the clear range of data, the number of bases in the clear range, the percentage of ambiguities in the clear data, and the original length of the sequence. This report is a text file that can be opened in any word processor or in Factura.

The batch report appears on the screen in the batch report window.

For a description of the batch report, refer to Section 3. To print or save a batch report, refer to Section 5.

Diagram of File Input And Output

Figure 2-2 shows the relationship of the Factura input and output files.




Figure 2-2. Factura input and output files

You can import and export sequence files that are in Sample or text file format. They are displayed and processed using the batch worksheet. After feature identification, you can save the results in a batch report file, and into the individual sequence files. You can also save the entire batch worksheet in a batch worksheet file.

Factura Windows and Views

The Batch Worksheet Window

The batch worksheet window displays a worksheet that summarizes the sequences and feature identification parameters you have chosen to work with in Factura. It allows you quick access to the sequence files for editing and display. In it you can easily edit the identification parameters either individually for each sequence or for the entire batch of sequences prior to processing.

Information about working with a batch worksheet is provided in Section 3.

The Batch Report Window

After Factura processing, you can choose to display the batch report, detailing the results of processing. It appears on the screen in a batch report window.

For a description of the batch report, refer to Section 3. To print or save a batch report, refer to Section 5.

The Sequence Window

The sequence window displays the contents of a single sequence file. Sequences from Applied Biosystems genetic analyzer instruments produce information in four views, described under the headings below. You can change to the different views by clicking the appropriate buttons in the bottom left corner of the window.

When you open a sequence window sample file, a sequence created using the New Sequence command, or a text sequence entered on a word processor, information in the sequence window is only available in three views (Sequence, Annotation, and Feature). Electropherogram View is available only for sequences that have electropherograms.

For a full description of the Factura sequence window, refer to Section 4.



Sequence View

Sequence View displays the listing of individual bases in the sequence (see Figure 2-3). To change to Sequence View from any of the other views, click the button shown here.




Figure 2-3. Sequence View of the sequence window

In sequence view you can search for specified patterns, modify features, and edit the sequence. To learn more about finding patterns, refer to page 4-21. For information about modifying features and editing, refer to page 4-14.



Annotation View

Annotation View shows information stored in the file about the run that produced the sequence data, as well as annotation from a database entry. If the file is not an Applied Biosystems genetic analyzer instrument sample file, some basic information is presented. Click the button shown here to display Annotation View. Figure 2-4 shows an example of Annotation View.




Figure 2-4. Annotation View of the sequence window



Feature View

Feature View allows you to view identified features in sequences produced by an Applied Biosystems genetic analyzer and database sequences. To display Feature View, click the button shown here.

Figure 2-5 shows an example of Feature View.

Feature key:	Range(s):	Description:
Δ *ABI_Vector	1 17	This range of bases indicates the vector
Δ *ABI_Ambiguity	1 3	This range of bases indicates ambiguity
+ *ABI_Limits	1 450	This is the confidence range
*ABI_Multibase	2	0.89 N
*ABI_Multibase	3	1.00 V
*ABI_Multibase	14	0.51 Y
*ABI_Multibase	420	0.79 B
*ABI_Multibase	423	0.74 B
*ABI_Multibase	435	0.66 B
*ABI_Multibase	440	0.90 B
*ABI_Multibase	465	0.72 N
*ABI_Multibase	472	0.54 N
*ABI_Multibase	473	0.57 B

Figure 2-5. Feature View of the sequence window

In Feature View you can add, modify, and remove features.

Features identified by Factura might include vector segments, ambiguous regions, and positions that represent multiple bases (IUB code calls).



Electropherogram View

Electropherogram View shows a four-color picture of a sequence, with peaks that represent the bases or amino acids. This view is available only with Applied Biosystems genetic analyzer instrument data files. Click the button shown here to display Electropherogram View.

Figure 2-6 shows an example:




Figure 2-6. Electropherogram View of the sequence window

In Electropherogram View you can edit individual bases in the sequence, and display the original data for comparison.

3 Working With the Batch Worksheet

The batch worksheet is the main window in the Factura application. It summarizes the sequences and feature-identification parameters you have chosen to work with and allows you quick access to the sequences for editing and display. On the batch worksheet you can easily edit the identification parameters prior to processing, either individually or for the entire batch of sequences.

This section provides the following information about the batch worksheet.

- How to open it, close it, and set it up
- What it looks like
- How to add and remove sequences
- How to change the processing parameters displayed on it
- How to submit the batch of sequences for processing

Contents

Specifying Settings for the Batch Worksheet	3-3
Settings Dialog Box Options	3-4
Identify Vector Sequence	3-4
Identify Ambiguity	3-6
Identify Confidence Range	3-7
Identify IUB/Heterozygous Bases	3-7
Automatic Options	3-9
Opening the Batch Worksheet	3-10
Opening a New Batch Worksheet	3-10
Opening an Existing Batch Worksheet	3-10
From the Factura Application	3-11
From the Finder	3-11
Description of the Batch Worksheet	3-11
Adding and Removing Sequences	3-14
Adding Sequences to the Batch Worksheet	3-14
Moving in the Batch Worksheet	3-16
Removing Sequences from the Batch Worksheet	3-17
Changing Batch Worksheet Parameters	3-18
Submitting the Batch Worksheet	3-21
The Batch Report	3-22
Closing the Batch Worksheet	3-25

Specifying Settings for the Batch Worksheet

The Settings command allows you to specify the vector sequence, ambiguity, confidence range, heterozygote/IUB code threshold, and other parameters to use for the submission of a worksheet with sequences from the ABI 373, ABI PRISM 377 or 310 for processing in Factura.

IMPORTANT Set up your custom libraries before you designate Settings. The choices in the Settings dialog box (Figure 3-1) for vector, enzyme, and primer draw from the lists you generate when you set up the libraries. To set up the libraries (which you need to do only once), refer to Section 6.

To designate Settings:

1. Choose Settings from the Worksheet menu. The Settings dialog box appears.

Settings

Identify Vector Sequence
Vector: M13MP18
Primer: M13-21
Cloning Site: SmaI Match parameters...

Identify Ambiguity
Remove bases from the beginning and end until no more than
1 ambiguities remain out of 20 bases.
 Reject sequences with > 10% ambiguities.

Identify Confidence Range
Keep the range from 1 to 450

Identify IUB/Heterozygous Bases
Threshold: 50% Update edited bases

Automatically save to sequence file
 Revert sequences to original base calls
 Use these settings as default value

OK
Cancel

Figure 3-1. Settings dialog box

Note *The names listed for vector and cloning site initially reflect the top names in the Vectors Used and Cloning Sites Used list in each library (see Setting Up the Libraries on page 6-5). The name listed for primer is the default or top name in the Primer library list (see Setting up the Primer Library on page 6-10).*

2. Enter the proper parameters (descriptions are provided below).
3. Click OK to accept your settings.

Settings Dialog Box Options

The identification parameters you choose in the Settings dialog box (except Vector, Primer and Cloning Site) depend on the quality of the data you plan to analyze. Familiarize yourself with the descriptions in this section, and try various settings to determine what values give the best results.

Identify Vector Sequence

Identify Vector Sequence

Vector: M13MP18

Primer: M13-21

Cloning Site: SmaI Match parameters...

This checkbox is selected by default. It allows you to specify the vector, cloning site, primer, and match parameters. Deselecting the box disables vector identification.

- Use the Vector, Primer, and Cloning Site pop-up menus to choose parameter settings from the respective libraries.
- Click the Match Parameters button to display the dialog box shown in Figure 3-2.




Figure 3-2. Match Parameters dialog box

- Use the Maximum permitted gaps field to specify the maximum number of gaps allowable in the overlapping sequences when the comparison is made between the vector data in the sample file and actual vector sequence data. If you reduce this number, the identification is performed more quickly, but is not as sensitive for a noisy vector.
- Use the Minimum match between vector & sequence field to specify how closely the vector data included in the sample file must match actual vector sequence data for the vector to be identified. Choices range from 50% to 100%.

For example, if a sample file contains a vector that is 20 bases long before the insert, specifying a percentage of 90% would mean that a match of 18 bases would identify the vector, but a match of only 17 bases would fail.

- Select the checkbox labeled Search for short inserts to perform a more rigorous search for the vector, both at the 5' and 3' ends of sequences from the ABI 373, ABI PRISM 377 or 310 sequences. You should always select this checkbox if your sample data might contain short inserts (shorter than the ABI 373, ABI PRISM 377 or 310 read, for example, <400 bases) in addition to long inserts (>800 bases).

If you don't select this parameter, a faster algorithm is used. Selecting the parameter causes the use of a more detailed algorithm and slightly increases processing time.

- Select the checkbox labeled Ignore leading N's to instruct the application to ignore Ns at the beginning of a sequence when identifying the vector. If this checkbox is selected, vector identification begins at the first non-N base. **We recommend that you always check this checkbox.**

Identify Ambiguity

Identify Ambiguity
 Remove bases from the beginning and end until no more than
 N's remains out of bases.
 Reject sequences with > % N's.

This checkbox is selected by default. It enables identification of regions of ambiguous sequence. Deselecting the box turns off ambiguity identification.

- When you enter values in the first two fields, the application removes one base at a time from each end of the data (beginning and end) until less than the specified number of ambiguities remain in the specified beginning and ending ranges.

With the default setting, 1 out of 20 bases, the application examines the first and last 20 bases of the sequence. If two Ns are located in either end, the application moves in one base and checks again. This continues until each end of the sequence consists of 20 bases with no more than one ambiguity.

The default setting is equivalent to specifying that the beginning and ending ranges of the data must contain no more than 5% ambiguity.

IMPORTANT *This method looks at a sliding window of data; it does not check all possible windows. It stops after it finds a series of unambiguous bases that matches the parameter specified. It might therefore miss a region in the middle that has more ambiguity. If you know that substantial ambiguity exists in the middle of the sequence, you might specify a larger window size as shown in the example.*

Example:



In this example, the window stops at each end when 20 unambiguous bases are located, and does not identify the ambiguous region (all Ns) in the middle of the sequence. Specifying a larger window size such as 1 N out of 50 bases might identify the ambiguity.

- Select the checkbox and enter a value in the third entry field to specify the total percentage of ambiguity to be allowed in the data *after* the beginning and ending range specification (immediately above) is met.

Identify Confidence Range

Identify Confidence Range
 Keep the range from to

This parameter allows you to adjust the length of the sequence fragment that you will accept as containing good data. Sequence data tends to deteriorate after a certain number of bases are sequenced. The first item under *Identify Ambiguity* can limit the ends of the data when many ambiguities exist. The Identify Confidence Range parameter limits the range of the data to what is likely to be accurate even if few ambiguities are identified.

IMPORTANT *This value depends on what sequence you are using. If you use “stretched version” of the ABI 373A or ABI PRISM 377, which generate longer reads, increase the confidence range to 600–900 bases.*

Identify IUB/Heterozygous Bases

Identify IUB/Heterozygous Bases
 Threshold
 Update edited bases

- Use the pop-up menu to specify a threshold used to assign IUB codes to mixed base positions. The application uses this ratio to compare the highest peak to each of the other three peaks in the same location. If the ratio between any of the three lower peaks and the higher one is above the threshold percentage, an IUB code is assigned. If not, the existing base assignment is retained.

For example, assume the highest peak at an ambiguous location is C and the application calculates ratios as follows:

- A to C is 85%
- T to C is 60%
- G to C is 30%

If the threshold parameter is set at 80%, the IUB code assigned is M (representing A or C). If the threshold is set at 60%, the IUB code assigned is H (indicating A, C, or T).

- Select the checkbox labeled Update edited bases to store the IUB codes in the sequence files when you use the Save to Sequence Files command in the Worksheet menu (refer to page 5-9). If you deselect this option, you must manually edit the data in the sequence files to clarify ambiguities.

Note *IUB/heterozygote identification operates with sequences from ABI PRISM 310, ABI 373, and ABI PRISM 377 sample files, since it requires four-color electropherogram data.*

Automatic Options

<input type="checkbox"/>	Automatically save to sequence file
<input type="checkbox"/>	Revert sequences to original base calls
<input checked="" type="checkbox"/>	Use these settings as default value

- Select the checkbox labeled Automatically save to sequence file to save the results of batch processing to the sample file feature tables after you submit the batch worksheet.
- Select the checkbox labeled Revert sequences to original base calls if you want to restore original data to ABI 373 and ABI PRISM 377 or 310 sample files that you have edited. This might be useful if your data was substantially edited and you want to return it to the original form before identifying features. The reversion occurs when you submit the batch worksheet for processing.

IMPORTANT *If you select the Revert parameter, the edited sequence is overwritten with the original data. The edited information is lost. If the sequence has a feature table saved from a previous submission, the feature table remains with the file. It is not deleted.*


The Revert parameter is normally disabled so that you can edit sample files, either in Factura or in another GeneAssist application, and submit the edited files to Factura for processing.

- Select the checkbox labeled Use these settings as default value to specify that the settings you choose in the Settings dialog box take precedence over any information stored in the sequence files when sequences are imported into the batch worksheet.

When you finish specifying parameters on the Settings dialog box, click OK.

Opening the Batch Worksheet

When you start Fatura, a blank batch worksheet (Figure 3-3) appears on the screen:



#	File name	Vector	Site	Primer	Match	Ambiguity	Limits	IUB/Hetero
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								

Figure 3-3. The new batch worksheet

The application opens with Untitled-1 as the worksheet name and the number increments with each new worksheet opened.

If you have not yet set up libraries or specified settings, close this worksheet and perform those tasks. For a brief description, refer to Section 6, *Setting up Fatura Libraries* and to page 3-3 to specify settings. When you open a new batch worksheet the libraries and settings you specify are automatically invoked.

Opening a New Batch Worksheet

If you are starting a new application with Fatura and have set up the libraries and worksheet identification parameters, you can open a new batch worksheet by choosing New Batch Worksheet from the File menu.

Opening an Existing Batch Worksheet

You can open a previously saved batch worksheet either before or after you start Fatura.

From the Factura Application

To open an existing batch worksheet, choose Open Batch Worksheet from the File menu. A standard file dialog box appears. Select the file you wish to open. The file dialog box shows folders at different levels but shows only batch worksheet file names.

From the Finder



Batch Worksheet File icon

Batch worksheet files are identified by the icon shown here. When you double-click a batch worksheet icon, Factura starts automatically and displays the batch worksheet just as it was saved to the file after you were last using it.

Description of the Batch Worksheet

The gray area of the batch worksheet beneath the standard Macintosh title line displays pop-up menus and entry fields when you want to change the parameters in the columns of the worksheet. It is more fully described in *Changing Batch Worksheet Parameters* on page 3-18.

The nine columns of the batch worksheet reflect information about the sequences and identifications you are working with. Use the Tab key to move from field to field or click on a field to select it. The columns are described as follows:

- The first column of the worksheet lists the numerical order (1–10) in which sequences were added. Once you have added sequences (detailed on page 3-14), this column displays a symbol that indicates the present state of the project. A diamond symbol (◆) next to a sequence number indicates that the sequence is not processed. When the sequences are processed, a triangle symbol (▲) appears next to each of the sequence numbers.

Figure 3-4 shows an example of a batch worksheet with unprocessed sequences on it.






Figure 3-4. Batch worksheet containing unprocessed sequences.

- The File name column lists the file name of each added sequence and other information about it. Since the default name field contains only nine characters, a longer name than this is only partially shown unless you widen the column.

Note Widen a column by placing the cursor over the vertical line to the right of the column. When the cursor symbol changes from  to , hold down the mouse key and drag the line to the right.

- The Vector column lists the vector name in each sequence file. If none is assigned, or if you chose to use the parameters in the Settings dialog box as defaults, this column shows the name of the vector you specified.
- The Site column lists the cloning site in each sequence file. If none is assigned, or if you chose to use the parameters in the Settings dialog box as defaults, this column shows the cloning site you designated.
- The Primer column lists the primer in each sequence file. If none is assigned, or if you chose to use the parameters in the Settings dialog box as defaults, this column shows the primer you designated.

-
- The Match column shows the current match parameters value you specified in the Settings dialog box.
 - The Ambiguity column lists the three current ambiguity values you specified in the Settings dialog box.
 - The Limits column shows the current limits of the confidence range you specified in the Settings dialog box.
 - The IUB/Hetero column shows the threshold percentage you chose for heterozygote identification in the Settings dialog box. This identification operates only with sequences from ABI 373 and ABI PRISM 377 or 310 sample files, since it requires electropherogram data.

Note *If you save results to a sample file after processing, the settings saved with it are those used for processing.*

Adding and Removing Sequences

Note To see the contents of a sequence file before adding it to the batch worksheet, choose one of the Open Sequence commands or the New Sequence command from the File menu (refer to page 4-4).

Adding Sequences to the Batch Worksheet

Adding sequences to a batch worksheet fills out the worksheet with sequence and Settings information (see Figure 3-4 on page 3-12). The parameters you specified when you set up the libraries and chose settings for the batch worksheet are reflected.

To add sequences to a new batch worksheet:

1. Choose Add Sequences from the Worksheet menu.

A file dialog box appears.




Figure 3-5. Add Sequences file dialog box


2. Open the folder containing your sequences.

To limit the file list, select and deselect the checkboxes under the heading File Type. The file list only shows files of the type you select.


3. Add files in one of the following ways:

- Add any single file by selecting it and clicking Add.
- With one file selected, click Add All to add all the files in the folder.

If you add multiple sequences using Add All, the following dialog box appears while the sequences are loading into the batch worksheet:



If you have deselected the Exclude Files Already in List checkbox (see Figure 3-5) and attempt to add a sequence already in your batch worksheet, a dialog box appears that allows you to reconsider.



If you like the default parameters on the worksheet after you have added sequences, the worksheet is ready to submit (described on page 3-21). If you want to change the parameters, see *Changing Batch Worksheet Parameters* on page 3-18.

Moving in the Batch Worksheet

Click the number to select a sequence

Click any field to edit it, or click the title to select an entire column

When the cursor is this symbol, drag it to widen a column

Click the number column of an empty row to de-select sequences

Use the scroll bar to scroll and the size box to stretch the window

#	File name	Vector	Site	Primer	Match	Ambiguity	Limits	IUB/Hetero
1	Sample 1...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
2	Sample 2...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
3	Sample 3...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
4	Sample 4...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
5	Sample 5...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
6	Sample 6...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
7								
8								
9								
10								

You can click any field in the batch worksheet to select it for editing. Click the number column on the left to select an entire sequence row.

To move from column to column within one row:

Press the Tab key or the Right Arrow key (→) to move to the right one field. Press the Left Arrow key (←) to move to the left one field.

To move from row to row within one column:

Press the Return key or the Down Arrow key (↓) to move down one field. Press the Up Arrow key (↑) to move up one field.

To select sequences:

When you select a sequence in the batch worksheet, the row containing the sequence is shown against a reversed background (highlighted). Select sequences as follows:



- To select a specific sequence and to deselect all other sequences, click the number of the sequence you wish to select.

-
- To select a consecutive range of sequences, click the number of the first sequence in the group, then hold down the shift key and click the number of the last sequence in the group.
 - To select a sequence and leave other (discontinuous) sequences selected, hold down the command (⌘) key and click the sequence.
 - To deselect sequences, click the number column of an empty row.

To select entire columns:


You might want to quickly select an entire column to use the Fill Down command, which is described on page 3-19. To do so, click the title of the column.

To change column width:

Place the cursor over the vertical line to the right of the column. When the cursor symbol changes from  to , hold down the mouse button and drag the line to the right to widen the column or to the left to decrease the column width.

You can widen the File name column to see more information about the sequence.

To see more of the sequences in the worksheet:

If some of the rows are not visible on the worksheet, you can use the size box () in the bottom right corner of the worksheet to stretch the window, or click the scroll bar on the right of the worksheet to move down the worksheet.

To scroll one window at a time, click in the gray region of the scroll bar. The window shifts so the top sequence becomes the bottom sequence or vice versa. If you click on an arrow at the end of the scroll bar, the window shifts one sequence at a time.

Removing Sequences from the Batch Worksheet

To remove a sequence from the batch worksheet, select the sequence on the worksheet by clicking in the left column between the diamond and the assigned sequence number. With the sequence highlighted, choose Remove Sequence from the Worksheet menu.


Changing Batch Worksheet Parameters

You can apply changes in the batch worksheet to either a single sequence or the entire batch of sequences you plan to process. The procedure is described first in this section, followed by more detail about the individual parameters.

To change identification parameters:

1. Click the field you wish to change. To change the same parameter for all sequences in the batch, click the parameter field for the first sequence in the column.

The editing field at the top of the worksheet displays either a pop-up menu or entry fields, depending on the parameter you wish to change.



After you have made the change, click the title row to select the column

2. Make the desired changes, using the pop-up menu or entry fields.
The change is automatically applied to the selected field.
3. To change the parameter for all sequences, hold down the mouse button and drag to select all table cells in the column, or select the column by clicking the title row.

Note *To make the change for selected sequences, hold down the command (⌘) key and click the cells to change.*


4. Choose Fill Down from the Worksheet menu.

Fill Down applies the change to all the selected cells in the column.

Note *Skip Step 3 and Step 4 to apply the change to only a single table cell.*

The following describes each of the pop-up menus and fields that can appear in the editing field at the top of the batch worksheet.

- When you click a cell in the Vector column, the following pop-up menu appears in the editing field:



Change the vector in the selected cell by clicking the menu and choosing a vector from the pop-up list, which shows all the vectors in your custom library.

- When you click a cell in the Site column, the following pop-up menu appears in the editing field:



Change the site in the selected cell by clicking the menu and choosing a cloning site from the pop-up list, which shows all the cloning sites in your custom library.

- When you click a cell in the Primer column, the following pop-up menu appears in the editing field:



Change the primer in the selected cell by clicking the menu and choosing a primer from the pop-up list, which shows all primers in the primer library.

- When you click a cell in the Match column, the following appears in the

editing field:

Select % match: %

You can choose a new percentage value in two ways. Either click the arrow to display a pop-up menu from which you can choose the new value, or select the existing value in the entry field and type a new value.

- When you click a cell in the Ambiguity column, the following editable entry fields appear in the editing field:

Specify %N Match: No more than N's out of bases.
Sequence < % N's.

To change the selected ambiguity setting, highlight an existing entry and type a new value.

- When you click a cell in the Limits column, the following editable entry fields appear in the editing field:

Specify the confidence range: to

To change a setting in the selected cell, highlight an existing entry and type a new value.

- When you click a cell in the IUB/Hetero column, the following appears in the editing field:

Select % match: %

You can choose a new threshold percentage value in two ways. Either click the arrow to display a pop-up menu from which you can choose the new value, or select the existing value in the entry field and type a new value.

Note *IUB/heterozygote identification operates only with sequences from ABI 373, ABI PRISM 377 or 310 sample files because it requires electropherogram data.*

Submitting the Batch Worksheet

When you submit a batch worksheet in Fatura, all sequences on the worksheet are processed according to the parameters shown. For each sequence, Fatura identifies the vector sequence, ambiguity, confidence range, and heterozygote features, and calculates the clear range, clear length and percentage of ambiguities. To store this information in the feature tables of the sequences, you must save the information to the sample files (see Section 5).

Note *IUB/heterozygote identification operates only with sequences from ABI genetic analyzer instrument sample files because it requires electropherogram data.*

To process a batch worksheet in Fatura:

1. Choose Submit from the Worksheet menu (see *Worksheet Menu* on page 7-10).

A progress indicator appears while the application is processing the sequences. A dialog box like that in Figure 3-6 appears after processing is complete.




Figure 3-6. Batch report dialog box

If you click Yes, a batch report like that shown in Figure 3-8 appears. The batch report is described on page 3-22.

After processing, the left column of the batch worksheet shows a triangle symbol (▲) for each sequence, indicating that the worksheet was processed (Figure 3-7).

#	File name	Vector	Site	Primer	Match	Ambiguity	Limits	IUB/Hetero
1	Sample 1 ...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
2	Sample 2 ...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
3	Sample 3 ...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
4	Sample 4 ...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
5	Sample 5 ...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
6	Sample 6 ...	M13MP19	SmaI	M13-21	80%	1/20/10%	1/450	50%
7								
8								
9								
10								

Figure 3-7. Batch worksheet after processing

The Batch Report

The batch report (Figure 3-8) repeats the Vector, Site, and Primer information from the batch worksheet and provides additional information described in Table 3-1.

#	Sequence	Vector	Site	Primer	Vector range	Ambiguity range	Confidence range	Clear range	Clear length	%N	Original length
1	Sample 1.short Sample 1.short	M13MP19	SmaI	M13-21	1-15	*	1-193	16-193	178	0%	193
2	Sample 2 Sample 2	M13MP19	SmaI	M13-21	*	*	1-402	1-402	402	1%	402
3	Sample 3 Sample 3	M13MP19	SmaI	M13-21	1-17	1-3 465-472	1-450	18-450	433	1%	472
4	Sample 4 Sample 4	M13MP19	SmaI	M13-21	*	415-434	1-434	1-414	414	1%	434
5	Sample 5 Sample 5	M13MP19	SmaI	M13-21	1-15	533-560	1-450	16-450	435	0%	560
6	Sample 6 Sample 6	M13MP19	SmaI	M13-21	1-12	503-527	1-450	13-450	438	1%	527

Figure 3-8. Batch report example

Table 3-1. Batch Report Columns

Column Heading	Description
Vector range	The portion or portions of the sequence data that correspond to the vector. As shown in the example report, one or two vector ranges might exist for any one sequence.
Ambiguity range	The portion or portions of the sequence data containing ambiguously called bases.
Confidence range	The largest sequence fragment that you will accept as containing good data. This parameter addresses the fact that sequence data tends to deteriorate after a certain number of bases are sequenced.
Clear range	A range of bases in the sequence data that represents clear data, that is, data not containing portions identified as vector or ambiguous and which is within the confidence range.
Clear length	The number of clear bases in the Clear range.
%N	The percentage of Ns (ambiguities) remaining in the clear data.
Original length	The length of the sequence data.

You can also observe the results of Fatura processing by opening individual sequence windows and examining the Sequence and Feature Views. For more information about viewing results in sequence windows, refer to Section 4.

To save the batch report:


You can save the batch report by making the batch report window active and choosing any of the three Save commands from the File menu (see *Saving the Batch Worksheet* on page 5-5).

To close the batch report window:

You can close a batch report window by clicking the close box in the upper left corner or by choosing Close from the File menu.

Note *It is a good idea to save the batch report in a different folder from the sequences. This allows you to use the Add All button when adding sequences to a new batch worksheet without accidentally adding the batch report.*

If you close the batch report before saving it, the following dialog box appears.



- To save the report in a batch report file, click Save. The application asks you for a file name and location.
- To continue to close the report without saving it, click Don't Save.
- To cancel closing the report, click Cancel.

To open a batch report file:

The batch report is a delimited text file, and can be opened by many different applications. In Factura, you can open a batch report file in two ways.

IMPORTANT *Do NOT use Open Sequence in the File menu. Doing so causes the application to try to open the batch report as a DNA sequence.*

- If Factura is active, choose Open Batch Report from the Worksheet menu.
- If the Finder is active, double-click the icon of the batch report you wish to open. This starts Factura automatically and displays the batch report just as it was saved to the file.




Batch report icon

Closing the Batch Worksheet


Note *It is a good idea to save your work before you close the batch worksheet. See Section 5 for instructions on saving.*

You can close a batch worksheet in two ways. You can either close the window by clicking the close box, or choose Close from the File menu. If you processed the batch worksheet and did not save it, you are given the opportunity to save the identified features on a sequence-by-sequence basis. The following dialog box appears for *each* sequence in the worksheet.



- To save identified features to the sample file, click Yes. The application asks you for a file name and location. The features are stored in the feature table of the sample file.
- To continue to close the worksheet without saving features for the named sequence, click No.
- To cancel closing the worksheet, click Cancel.
- To close the worksheet without saving features and to bypass additional dialog boxes, click the button labeled "No and no more questions!"

If you have modified the batch worksheet since it was last saved and have not saved it, the following message appears, allowing you to save the changes if you wish.



4 Working With Sequences in Factura

The Factura sequence window allows you to display information from your sequence files in four different formats. The window can be opened in several different ways, with or without a batch worksheet displayed.

This section provides:

- Descriptions of several ways to open a sequence window
- A brief description of the sequence window and its various views
- Procedures for editing and changing the appearance of a sequence in the sequence window

Contents

Opening the Sequence Window	4-3
Opening a Sequence File Independent of the Batch Worksheet	4-4
Opening an Existing Sequence File	4-4
Creating a New Sequence File	4-5
Opening a Sequence File from the Batch Worksheet	4-7
Description of the Sequence Window	4-8
Sequence View	4-9
Annotation View	4-10
Feature View	4-11
Electropherogram View	4-13
Editing in the Sequence Window	4-14
Adding, Deleting or Changing Bases	4-14
Editing Bases In Sequence View	4-14
Editing Bases In Electropherogram View	4-15
Adding, Removing, or Changing Features	4-16
Joining Features	4-20
Viewing Options in the Sequence Window	4-21
Finding Patterns in a Sequence	4-21
Displaying the Complement of a Sequence	4-25
Changing the Font and Grouping of a Sequence	4-26
Changing the Appearance of Features	4-27
Applying a Marking Style to a Feature	4-27
Changing the Definition of A Marking Style	4-28
Marking Features Automatically by Renaming Mark Styles	4-31
Changing the Appearance of Electropherogram View	4-32
Closing the Sequence Window	4-42

Opening the Sequence Window

Open a sequence window to view or edit the contents of an individual sequence file, or to see the results of batch processing on an individual sequence. You can open the window in several ways, as described below.

Two menus become available whenever a sequence window is active. The Edit menu allows you to undo your last action, manipulate text and graphics, and select the entire contents of a window view. The Sequence menu allows you to change the visual display of the information in the sequence window.

Note *The Edit menu and the Sequence menu are only available when a sequence window is active. You cannot use either menu when the batch worksheet is the active window.*

Sequences open in the Sequence View (Figure 4-1). Other possible views are discussed beginning on page 4-9.




Figure 4-1. Sequence View

You can create a new sequence file or open and work with sequence files that already exist, whether or not you have a batch worksheet open.

Opening a Sequence File Independent of the Batch Worksheet

You might wish to open sequence files independent of the batch worksheet to examine the contents of the files or to copy information from text files into new sequence files.

Because Fatura does not save feature table information into text files, your text files must be copied to a different format so that Fatura-identified features can be stored in them (for a procedure, see page 4-7).

Opening an Existing Sequence File

Two Open Sequence commands are provided in Fatura. The first command, which is followed by an ellipsis (...), displays the standard dialog box allowing you to open sequence files of all three types (Sample, Inherit, and Text).

If you have a large number of different types of sequence files, the second Open Sequence command can simplify file selection by filtering for the selected file type. It displays a submenu, as shown in Figure 4-2.




Figure 4-2. Open Sequence submenu

You must choose the type of sequence file (Sample, Inherit, Text) while selecting the command. Only files of the chosen type appear in the standard file dialog box.

The file types are as follows:

- *Sample* files contain sequence data in four different formats: Sequence View presents a listing of data; Annotation View may contain annotation from an automated sequencer; Feature View conveys the Factura-generated results to the Sequence Navigator application in the form of defined vector, ambiguity, confidence range, and IUB/heterozygote features; and Electropherogram View displays 373 sequence data in electropherogram format.
- *Inherit* files are files exported from any of the applications of the INHERIT sequence analysis system. When opened within Factura they contain only Sequence, Annotation and Feature Views.
- *Text* files contain the same views as the INHERIT type files when you open them in Factura. When you save from Factura to a text file and open it in a word processor, only the sequence listing view is available.

Creating a New Sequence File

Creating a new sequence file is a useful tool that enables Factura to process information you have drawn from a text file and store the results in a file with a feature table. Since Factura does not save feature table information into text files, your text files must be copied to a different format so that Factura-identified features can be stored in them. The new sequence file created by Factura is like an ABI 373 sample file except that it contains no electropherogram data.

To create a new sequence file:

1. Select New Sequence from the File menu.

An empty sequence window appears (Figure 4-3).




Figure 4-3. Empty sequence window

2. Create a new sequence in one of the following ways:
 - Place the insertion point in the central box and type a sequence.
 - Use Finder and the Macintosh Copy (or Cut) and Paste commands to move sequence data from another application.
3. Choose Save As from the File menu to save the new sequence under a new file name.

To copy a sequence from a text file into a new sequence file:

1. Open the existing text file using Open Sequence, as described on page 4-4.
This opens a sequence window containing the sequence data stored in the text file.
2. Choose New Sequence from the File menu.
A new, blank sequence window appears.
3. Click the sequence window containing the text file to make it active.
4. Choose Select All from the Edit menu.
The entire sequence is highlighted.
5. Choose Copy (⌘-C) from the Edit menu.
6. Click the new (empty) sequence window to make it active.
7. Choose Paste (⌘-V) from the Edit menu.
8. Choose Save As from the File menu.
9. Enter a name for the new sequence.
10. Click Save (⌘-S).

Opening a Sequence File from the Batch Worksheet

You can open sequence windows for any number of selected sequences in the batch worksheet.

To open a sequence window from the batch worksheet:

1. Select a sequence.
2. Choose Show Sequence from the Worksheet menu.
A sequence window appears for the selected sequence.

Note *If you select more than one sequence before choosing the command, a sequence window appears for each selected sequence.*

Description of the Sequence Window

Following is a brief description of the sequence window displays and buttons (Figure 4-4).




Figure 4-4. The sequence window displays and buttons

Immediately below the standard Macintosh title line and close box is a display window to the right of a lock image. The horizontal line in this summary graphic represents the length of the sequence and reflects the cursor position as you move it to different places in the sequence. The arrowhead on the summary graphic indicates the orientation of the sequence.

If you click the lock image, the sequence is protected from edits. You cannot Cut from or Paste to the sequence (using the Edit menu) as long as the lock is closed. Click the image a second time to unlock it.

The main portion of the sequence window contains the information from the sequence file. You can display up to four different views using the buttons located in the bottom left corner of the window. The function of each button is described under a heading below.

Note *If you open a database sequence saved in the GeneAssist application, a sequence created with the New Sequence command, or a Text sequence entered on a word processor, the sequence window shows only three of the views. Electropherogram View is available only for sequences produced by Applied Biosystems genetic analysis instruments.*



Sequence View

Sequence View (Figure 4-4) is the default view when you open a sequence window in Factura. The center column contains the sequence list. The left and right columns show the base positions at the beginning and end of each row. To change to Sequence View from any of the other views, click the button shown here.

In sequence view you can search for specified patterns, show the complement of the sequence, add features, and change the way the sequence shows on the screen. More information about these options is available as follows:

- To find a specified pattern, see page 4-21
- To complement the sequence, see page 4-25.
- To add or remove features, see page 4-16.
- To edit the display of the sequence, see page 4-26.



Annotation View

Annotation View (Figure 4-5) shows information stored in the file about the run that produced the sequence data, as well as annotation from a database entry. If the file is not from an Applied Biosystems genetic analysis instrument sample file, some basic information is presented. Click the button shown here to display Annotation View. Figure 4-5 shows an example of Annotation View.




Figure 4-5. Sequence Annotation View



Feature View

Feature View (Figure 4-6) displays Factura-identified features only after you have submitted the batch worksheet for processing (refer to page 3-21) and updated the sequence files (refer to page 5-9). It also shows the features for a database entry. To display Feature View, click the button shown here.

The triangle indicates a hidden feature

The plus sign indicates an active feature

Feature key:	Range(s):	Description:
△ *ABI_Vector	1 17	This range of bases indicates the vector
△ *ABI_Ambiguity	1 3	This range of bases indicates ambiguity
+ *ABI_Limits	465 472	This is the confidence range
*ABI_Multibase	2	0.89 M
*ABI_Multibase	3	1.00 V
*ABI_Multibase	14	0.51 Y
*ABI_Multibase	420	0.79 B
*ABI_Multibase	423	0.74 B
*ABI_Multibase	435	0.66 B
*ABI_Multibase	440	0.90 B
*ABI_Multibase	465	0.72 M
*ABI_Multibase	472	0.54 M
*ABI_Multibase	473	0.57 B

Ratio of peak heights (refer to Identify IUB/Heterozygous Bases on page 3-7)

Figure 4-6. Sequence Feature View

After you update the sequence files with the results of batch worksheet processing, Factura adds feature ranges to the view, identifying portions of the data that represent vector, ambiguity, and confidence range. When you use these sequences with other ABI PRISM data products, data in the vector and ambiguity ranges can be ignored or hidden, effectively eliminating poor quality data.

Note All of the information is maintained in the original data. The data that does not use these feature ranges is an editable copy of the original data.

In Feature View you can add, change, or remove features, and change the colors or borders that mark the features. For more information about these procedures, see the following references:

- To add, remove, or edit features, see page 4-16.
- To change colors or borders that mark the features, see page 4-27.

Note *Factura does not save feature table information to text files. Feature View for text files does not display Factura-identified features. To change a text file so that it can display Factura-identified features, see the copying procedure on page 4-7.*



Electropherogram View

Electropherogram View (Figure 4-7) is available only with ABI genetic analysis instrument data files. It is helpful to compare several files of the same data in Electropherogram View if you want to edit ambiguously called bases. Click the right-most button to display Electropherogram View.




Figure 4-7. Sequence Electropherogram View

The default view shows a range of bases in the region of the sequence where you have placed the insertion point. You can zoom in or out, change the horizontal and vertical ruler display, and selectively turn off the drawing of particular bases. For more information about working with Electropherogram View, see the following references:

- To edit bases in Electropherogram View, see page 4-15.
- To change the display of information in Electropherogram View see page 4-32.

Editing in the Sequence Window

You can perform editing functions in a sequence window, including:

- Adding, deleting, or changing bases in the sequence
- Adding, removing, or changing the features of the sequence.

You can edit individual bases in Sequence View or Electropherogram View. You must be in Feature View to remove or change features.

For information about complementing the sequence, finding patterns, and changing how the sequence displays on the screen, see *Viewing Options in the Sequence Window* on page 4-21.

Adding, Deleting or Changing Bases

You can add, delete, or change bases in both Electropherogram View and Sequence View. For information about changing how the sequence appears on the screen, see page 4-26.

Editing Bases In Sequence View

In Sequence View you can use the standard Edit menu commands to cut, copy, paste, and clear bases or ranges of the sequence in the active window. You can select the entire sequence (including marked features) with Select All in the Edit menu. The Edit menu commands operate as described in the *Apple System Software User's Guide*.

Following are brief descriptions of how to add, delete, and replace bases.

Note *Editing bases in Electropherogram View (see below) is more precise than editing in Sequence View. If you edit ABI PRISM sample files in Sequence View, you should check Electropherogram View when you are finished.*

To add a base or range of bases to the sequence:

1. Place the insertion point at the position in the sequence where you wish to add a base.
2. Type the characters you wish to insert.

To delete a base or range of bases from the sequence:

1. Select the base or range of bases.
2. Press the Delete key or choose Clear from the Edit menu.

To change a base in the sequence:

1. Select the base you wish to change.
2. Type the new character you wish to be in that position.

Note *You can also place the insertion point to the right of the character you wish to replace, press the Delete key, then type the character you want in that position.*

Editing Bases In Electropherogram View

In Electropherogram View the Edit menu commands are not available and you can edit only one base at a time. You can add, delete, or change bases in much the same way as described above for Sequence View; however, the spacing of the characters is much more precise.

Multiple base positions (approximately ten) are available between the displayed bases in Electropherogram View. If you place the insertion point between two characters and click, a position is selected. Following are some hints about moving the selection from one position to another:

- To move from base to base, use the Left Arrow and Right Arrow keys.
- To move from position to position (often pixel-by-pixel) hold down the Option key while you use the Left and Right Arrow keys.

IMPORTANT *Because the available base positions are so close together, it is possible to select a position very close to one of the bases when you are actually trying to select the base itself. If you do so, you might insert a character when you intend to change an existing character. To be sure you are selecting the base, use the following procedure.*

To select a base in Electropherogram View to delete or change it:

1. Place the insertion point to the left or the right of the character you wish to select and click.
2. Press the Right Arrow key or Left Arrow key to move to the base you wish to select.

This procedure ensures that you have selected the base, not a position only one pixel away from it. Once you have selected the base, you can delete it by pressing Delete or replace it by typing a new character.

To add a base in Electropherogram View:

1. Place the insertion point between bases where you wish to insert the character and click.
2. To move the insertion point closer to one of the flanking bases, hold down the Option key and use the Left or Right Arrow key to move to the position where you wish to insert the base.
3. Type the new character.

Note *If you add bases in Sequence View then change to Electropherogram View, the new bases are spaced as evenly as possible between the two previously existing bases.*

Adding, Removing, or Changing Features

After you have processed a batch worksheet and saved the features identified by Factura into the feature tables of the individual sequences, you can add other specialized features, and change or delete the features identified by Factura. You can add a feature in Sequence View, but you must be in Feature View to remove or change a feature, or to specify how it is marked for display.

Note *For information about how to mark features with color and borders, see page 4-27.*

To add a feature to the feature table of a sequence:

1. Select a feature, if you are in Feature View. If you are in Sequence View, select a range of characters.
2. Choose Feature from the Sequence menu, and choose Add from the submenu that appears.




Figure 4-8. Feature submenu

A dialog box appears:




Figure 4-9. Add/Edit Feature dialog box

Note *If you chose a range of characters in Sequence view, it is entered as a default in the Feature range table.*

3. Type a name for the new feature (up to 15 characters long) in the Feature Key entry field.

Note *To comply with the GenBank feature table definition the name should start with an asterisk (*).*

4. Enter text describing the new feature in the Description entry field.
5. Choose a marking style from the Style pull-down menu.

Using the Style menu, you can choose one of eight marking styles to help you identify your feature. The color of the style you choose identifies your feature in Sequence View and Feature View, and in the summary graphic in all views. For information about defining the marking styles, see page 4-28.

6. Check the range of your new feature.

If you selected a range prior to opening the dialog box, it shows as the default range. You can change the range using the two entry fields under Feature Range(s). Select either the beginning or ending value in the appropriate entry field, type a new value, and click Replace.

7. Note the following, and perform any of the options you wish:
 - You can add multiple ranges to your feature. For each new range, enter the beginning and ending values in the entry fields and click Add. If you do not change the Feature Key and Description entry fields, the new range is added to the same feature.
 - If you wish to remove a range, select it in the feature window on the right side of the dialog box, and click Delete.
 - If you wish to replace a range with a different one, select the range you wish to replace, type a new range in the Feature Range(s) entry fields, and click Replace.
8. When you are finished defining the feature, click OK.

IMPORTANT *If you reprocess the batch worksheet and save the reprocessed features to the sequence files, the new features overwrite the feature tables. Any special features you added disappear. You must add them again after reprocessing. Refer also to Saving Results to Sequence Data Source Files on page 5-9.*

To remove a feature from the feature table of a sequence:

1. Be sure you are in Feature View.
2. Select the feature you wish to remove.
3. Choose Feature from the Sequence menu, and choose Remove from the submenu that appears (see Figure 4-8 on page 4-17).

To change a feature in the feature table of a sequence:

1. Be sure you are in Feature View.
2. Select the feature you wish to change.
3. Choose Feature from the Sequence menu, and choose Modify from the submenu that appears (Figure 4-8 on page 4-17).

The Add/Edit Feature dialog box (shown in Figure 4-9 on page 4-17) appears.

4. Make changes as described in Step 5 through Step 7 under *To add a feature to the feature table of a sequence:* (the steps are on page 4-18).

You can remove one or more of multiple ranges, add additional ranges, change the values of existing ranges, or change the style that marks the feature. You can also change the name or description of the feature by typing new information in the Feature Key or Description entry fields.

5. Click OK to implement the changes.

Joining Features

You can join two or more individual features in the feature table into a single feature with multiple ranges. The following procedure creates a new feature that contains all the ranges included in the individual features you select. The original features remain in the feature table. If you wish to delete any of them, see the procedure under *To remove a feature from the feature table of a sequence*: on page 4-19.

To join features in the feature table of a sequence:

1. In Feature View select the features you wish to join by clicking each while holding down the command key (⌘) key.

<i>Note</i>	<i>To select contiguous features, click the first feature in the group, then hold down the shift key and click the last feature in the group.</i>
-------------	---

2. Choose Feature from the Sequence menu, and choose Join from the submenu that appears.

The Add/Edit dialog box appears, with blank Feature Key and Description entry fields. The Feature Ranges box at the right shows all ranges from the features you selected.

3. Type a new name for the joined feature in the Feature Key entry field.

<i>Note</i>	<i>To comply with the GenBank feature table definition the name should start with an asterisk (*).</i>
-------------	--

4. Type text describing the joined feature in the Description entry field.
5. Choose a new marking style from the Style pull-down menu.
6. Click OK.

Viewing Options in the Sequence Window

Among the viewing options available when the sequence window is active are:

- Finding patterns in the sequence
- Displaying the complement of the sequence
- Changing how the sequence appears on the screen

You must be in Sequence View to find patterns and display the complement of the sequence. You can change the appearance of the screen in Electropherogram View and Feature View as well as in Sequence View.

Finding Patterns in a Sequence

In addition to the standard Apple commands, the Edit menu contains commands for the following tasks.

- Finding a particular pattern of bases in the Sequence view of an active sequence window.
- Finding the same pattern again as many times as necessary.

Note *You cannot use the Find and Find Again commands in Electropherogram View. Use the Find command in Sequence view, and when the pattern is highlighted, switch to Electropherogram View. The electropherogram shows a range of bases in the highlighted region of the sequence.*


To find a pattern in a sequence:

1. Make sure you are in Sequence View in the sequence window.

Note *The search begins at the cursor position. If the pattern is before the cursor, it is only found if Wrap around is turned on (see page 4-25).*

2. Choose Find from the Edit menu.

The Find dialog box appears:



3. In the Find what? entry field, enter the sequence of characters you wish to locate.
4. Click the appropriate radio buttons and checkboxes. They are described here:
 - Select Literal to search for patterns that match exactly what you have typed.
 - Select IUPAC/IUB if you have entered an IUB character as part of the pattern. The Find command locates all possible matches. For instance, if the pattern you enter is TAR, the Find command locates either TAG or TAA.
 - Select Grep to set an expression for the search. Table 4-1 on page 4-23 describes some of the expressions you can use and how they function.

Table 4-1. Selection Expressions

Expression	Match Performed	Example
[] (brackets)	Any character inside the brackets	AA [AC] [GT] matches AAAG, AAAT, AACG, or AACT. [AGC] matches A, G or C.
[^] (brackets with ^ as first character inside)	Any character EXCEPT the character(s) inside the brackets	A [^AG] C matches ACC or ATC.
* after character	Zero or more such characters	AT [CG] * T matches ATT or ATCT or ATGGT, and so on.
. (period)	Any character	AA . A matches AAAA, AACAA, AAGA, AATA, AANA, and so on.
– (dash) enclosed by brackets	A range of characters	AA [A–z] matches AAA, AAC, AAG, AAz, and so on.

- Select Offset to move the cursor to the position or range of positions you specify. If you simply enter a number in the Find what? entry field, the insertion point is moved to that base position. If you enter a range of numbers, the whole range is highlighted.

Example:

If you enter 123 in the Find what? entry field, Figure 4-10 shows where the insertion point is positioned.




Figure 4-10. Example of single position found by Find/Offset

If you enter the range 123...250 (using Option-semicolon to create the ellipsis), Figure 4-11 shows how the range is highlighted after the Find is completed.

The summary graphic shows the relative position of the highlighted range




Figure 4-11. Example of range found by Find/Offset

-
- Select Case sensitive if you want upper- and lower-case variants of a letter to be recognized as different symbols.
 - Select Wrap around if you want the search to start again at the beginning of the sequence after it has reached the end. If the Wrap around checkbox is not selected, the search stops at the end of the sequence.

5. Click Find to perform the search.

The first instance of the specified pattern is highlighted and its position is marked in the summary graphic at the top of the sequence window (note examples on page 4-24).

Note *If you only want to find a pattern in the valid range, place the insertion point just before this range in the sequence.*

To find other occurrences of the same pattern:

Choose Find Again (⌘-G) from the Edit menu to bypass the dialog box and use the pattern defined in the previous Find command. Each time you use it, the next occurrence of the specified pattern is located.

Displaying the Complement of a Sequence

You can display the complement of a sequence in all sequence window views. When vector and ambiguity features have been identified, the characters in these features as well as feature locations in the Feature table are changed to reflect the complement. The display of original data in the electropherogram is complemented as well (see *Changing the Appearance of Electropherogram View* on page 4-32).

To display the complement of a sequence:

1. Make sure the relevant sequence window is active.
2. Choose Complement from the Sequence menu.

If you wish to return to the original data, simply choose Complement a second time.

Changing the Font and Grouping of a Sequence

In Sequence View you have the option to change the font size and style in which the sequence is displayed, as well as the grouping of the characters. If you want to change the way features are marked in Sequence View, see *Changing the Appearance of Features* on page 4-27.

To change the font and grouping of a sequence:

1. Make the sequence window active.
2. Choose Options from the Sequence menu.

The dialog box shown in Figure 4-12 appears.




Figure 4-12. Display Options dialog box

3. Indicate your choices as follows:
 - Use the Font pull-down menu to choose either Monaco (the default) or Courier. Both fonts have fixed-width letters. As you change the parameter, the changes show in the Sample box.
 - Use the Size pull-down menu to select one of seven values in a range from 6 to 24 points. As you change the size, you'll see the change in the Sample box.

-
- You do not need to select Prefer Outline Fonts unless you plan to increase the font size substantially (to the point it would otherwise look jagged). It specifies that dynamically scaled fonts be used.
 - The last entry field allows you to specify the number of characters displayed in each group in Sequence view.
4. Click OK to accept your choices.

Changing the Appearance of Features

To easily locate the features identified by Factura and the special features you add to sequence feature tables, you can mark them with defined styles. The marking appears in Sequence View, in Feature View, and in the summary graphic at the top of all views. The styles specify colors and borders that mark the features to identify them (the borders appear only in Sequence View). You can change the way features are marked and the definitions of the marking styles when you are in Feature View of a sequence window.

Applying a Marking Style to a Feature

You can apply any of eight styles to a previously created feature. Each style has a default color and border, but you can assign new parameters using the Edit Mark Styles command (see page 4-28).


Table 4-2 shows the defaults for the eight marking styles.

Table 4-2. Default Marking Styles

Style Name	Color	Border
Blue	Blue	No Underline
*ABI_Multibase	Red	Light Underline
Green Bold	Bright Green	Heavy Underline
Gray Double	Gray	Double Underline
Brown	Brown	No Underline
D Green Single	Dark Green	Light Underline
D Blue Bold	Dark Blue	Heavy Underline
Purple Double	Purple	Double Underline

To change the marking style that is applied to a feature:

1. Make sure you are in Feature View.
2. Select the feature.
3. Choose Mark Style from the Sequence menu.
4. From the submenu that appears, choose the style you wish to apply to the selected feature.



Changing the Definition of A Marking Style

You can change the color or border assigned to any of the eight marking styles that you can apply to identify a feature. You can do this any time a sequence window is active (in any of the four views).

If you mark a feature with one of the styles then edit the style, the change occurs immediately in the marked feature.

To edit the marking styles:

1. Make sure a sequence window is active.
2. Choose Edit Mark Styles from the Sequence menu.

A dialog box appears (Figure 4-13).




Figure 4-13. Edit Mark Styles dialog box

3. In the Available Styles list on the left, select the style you want to edit.

Note When the dialog box appears, the first style in the list is selected. Before making any change, be sure you select the style you want to edit. If you make any entries or choices in the right side of the dialog box, they are applied to the style selected in the list on the left.

4. Make changes or choices as follows:
 - To change the name of the style selected in the list on the left of the dialog box, type a new name in the Name entry field.
 - To assign a new color to a style, click the color box.

The color picker (Figure 4-14) appears:




Figure 4-14. Color picker

You can assign a new color to a marking style using either the scrolling windows to the left or the color wheel to the right.

To use the scrolling fields, place the cursor on one of the up or down scrolling arrows and hold down the mouse button. This scrolls through available colors and color parameters. While scrolling, observe the changes in the color box above. The upper half of the box shows the changes and the lower half retains the original color until you accept a new color by clicking OK.

To use the color wheel, place the insertion point on the desired color in the wheel and click. The new color appears in the upper half of the color box. You can continue moving the insertion point around and clicking until you obtain the desired color.

In either case, click OK when you are satisfied with your new color.

- Click the radio button identifying the border of your choice. Four choices are available to differentiate the display of a feature in Sequence View. The default places no underline under the marked feature. The second option, shown in the upper right corner of the Border box, places a line under the feature. The third option, directly below the default in the

Border box, places a bold line under the feature. The fourth option (in the lower right corner of the Border box) places a double line under the feature.

5. When you have finished editing marking styles, click OK in the Edit Mark Styles dialog box (shown in Figure 4-13 on page 4-29).

Marking Features Automatically by Renaming Mark Styles

You can have Factura automatically mark features as they are identified by renaming the marking styles with feature names. To see the results, you must reprocess the batch worksheet.

Example:

If you rename the Blue marking style to *ABI_Ambiguity, the ambiguity range identified by Factura is automatically marked blue when you process the batch worksheet.

To automatically mark a feature:

1. Choose Edit Mark Styles from the Sequence menu.

The Edit Mark Styles dialog box appears (see Figure 4-13 on page 4-29).

2. Select the marking style that you wish to use for the feature.
3. Type the name of the feature in the Name entry field.

IMPORTANT *The name you type must be exactly the same as the feature name.*

4. Click OK.
5. Choose Submit from the Worksheet menu.

After the application finishes processing, the features are marked with the color and border of the marking style you renamed.

Note *If you wish to stop the automatic marking of the *ABI_Multibase features, change the name of the *ABI_Multibase marking style.*

Changing the Appearance of Electropherogram View

You can change the appearance of Electropherogram View in several ways. You can choose to display the original form of the sequence data to compare it with editable data that automatically shows in the window. You can also change the scaling of the electropherogram, hide the display of certain bases or change the vertical scale displayed by the window.

To show the original sequence data on the electropherogram:

1. Make sure you are in Electropherogram View in the sequence window.
2. Choose Electropherogram in the Sequence menu, and choose Show Original from the submenu that appears (Figure 4-15).




Figure 4-15. Electropherogram submenu

A second line of data that represents the original form of the data appears at the top of the window. The line below it represents the data you can edit. Figure 4-16 is an example of Electropherogram View showing the original data.

Note *If you complement the sequence, the original form of the data is also displayed as a complement.*




Figure 4-16. Electropherogram View with original data

By comparing the two lines of data, you can keep track of any edits you make to the sequence.

To change the scaling of the electropherogram:

1. Choose Electropherogram from the Sequence menu.

The submenu choices other than Show Original are related to scaling of the electropherogram (see Figure 4-15).

Compare Figure 4-17, which shows the actual size electropherogram, to Figure 4-18, Figure 4-19, and Figure 4-20 to see the results of the scaling options.




Figure 4-17. Normal Electropherogram View (actual size)

Figure 4-17 shows the normal size of the electropherogram in Electropherogram View, when it is first displayed.

- To see views with ever greater detail, choose **Zoom In**. Click the region you want to zoom before choosing the command.




Figure 4-18. Electropherogram View after Zoom In command

After you have used the command three times in succession, it becomes grayed out and you can use either *Zoom Out* or *Actual Size* to reduce the size (both are described below). Compare Figure 4-17 with Figure 4-18 to see the effect of using *Zoom In* once. Note the difference in the horizontal scale at the top of the window.

- To see successively smaller scale views of the electropherogram, choose *Zoom Out*. Click the area you want to zoom out before choosing the command.




Figure 4-19. Electropherogram View after Zoom Out command

After you have used the command three times in succession, it becomes grayed out. Use either Zoom In or Actual Size to return the Electropherogram to its original size. Compare Figure 4-17 with Figure 4-19 to see the effect of using Zoom Out once.

- To scale the electropherogram so that the entire electropherogram fits within the standard size Electropherogram View, choose Full View. Figure 4-20 shows an example of the result. Use the Actual Size command to return to the original scale.




Figure 4-20. Electropherogram View after Full View command

- To return the Electropherogram to its original size after using Zoom In, Zoom Out, or Full View, choose Actual Size. Click the sequence at the area you wish to view in its actual size.

To change how data is displayed in Electropherogram View:

1. Choose Display Option from the Sequence menu.

The dialog box shown in Figure 4-21 allows you to change a number of settings for display of electropherograms.




Figure 4-21. Display Options Dialog Box

Set your preferred options as follows:

- Select and deselect the Show Data checkboxes to selectively turn off the display of one or more of the electropherograms for the four bases. The data shown on the line at the top of the Electropherogram View window is not affected.

Example:

If you deselect “A”, the electropherogram might appear like that shown in Figure 4-22. On your color monitor the green portion of the electropherogram corresponding to “A” is no longer visible.




Figure 4-22. Electropherogram with “A” deselected in Show Data

- The checkbox labeled “Show data points” is selected by default, so that the window shows vertical and horizontal rulers with values and tick points. If you are using relative values for the vertical scale (Show Relative Values parameter described below), the maximum vertical value is about 1200 full scale, so you may want to disable the display of data points.
- Use the Counts Per Tick entry fields to specify the horizontal and vertical indexing of the rulers on the Electropherogram display. The unit of measure (count) on the horizontal axis is the number of scans performed in the ABI Sequencers to obtain the data. The vertical axis indicates signal intensity.
- Select a radio button in the Vertical Display box to change the scaling of the vertical display. The default setting is Show Relative Values.

The default, when selected, compresses the scale of the electropherogram display vertically so that the electropherogram fits within a standard size sequence window.

As you can see by examining Figure 4-22, the Counts Per Tick value is not initially applied when the vertical display is set to Show relative values. In this mode the default Counts Per Tick value (20) would show too many

tick marks. If you desire, you can change the Counts Per Tick value to recalibrate the vertical scale.

Example:

If you change the vertical Counts Per Tick value from 20 to 125, the vertical scale shown in Figure 4-23 results:




Figure 4-23. Electropherogram (Recalibrated Vertical Scale)

Selecting the radio button labeled Show Real Values displays the real scale of the data, as shown in Figure 4-24 (the window was stretched vertically to show a larger portion of the electropherogram). It shows the original value (20) for vertical scale.

Note

The real value electropherogram shown in Figure 4-24 would have to be stretched quite a bit further vertically to show the tops of the highest peaks, which are at a scale value of approximately 1200.





Figure 4-24. Real Value electropherogram

Closing the Sequence Window

Note *It is a good idea to save any editing before you close the sequence window. See Section 5 for instructions on saving. To print any of the four views of the sequence window, refer to page 5-3 for instructions.*

You can close a sequence window in two ways. You can either click the close box, or choose Close from the File menu while the window is active.

If you have modified the sequence and have not saved it, the following dialog box appears, allowing you to save the changes.



- To save changes, click Save. Changes are stored in the sequence file.
- To continue to close the window without saving, click Don't Save.
- To cancel closing the window, click Cancel.

Note *When you save changes, the original sequence remains unmodified. Feature changes are stored in the feature table, and edits are stored in the editable copy of the sequence.*

5 Printing and Saving in Factura

This section describes:

- How to save the batch worksheet
- How to view and save results in a batch report
- How to save Factura-identified features to the individual sequence files
- How to export the contents of a sequence file into other formats

Contents

Printing Batch Reports & Electropherograms	5-3
Saving the Batch Worksheet	5-5
Saving a New Worksheet	5-5
Saving As a Different File Name	5-5
Saving a Copy	5-5
Saving Results in a Batch Report	5-7
Saving Results to Sequence Data Source Files	5-9
Saving Changes Made in the Sequence Window	5-12
Exporting Results to Other Formats	5-13

Printing Batch Reports & Electropherograms

You can print copies of the batch report created after Fectura processing, or any of the four sequence window views (Sequence, Electropherogram, Feature, or Annotation). The batch report has a header showing the date it was printed. All of the sequence window views are printed without a header and without any of the borders you see on the screen.

To print the batch report:

1. Click the batch report window to make it active.
2. Choose Print from the File menu.

The standard Macintosh print dialog box appears.

3. Click OK to start printing.

To print a view in the sequence window:

1. Click the sequence window to make it active.
2. Choose Page Setup from the File menu.

The dialog box allows additional options to the standard page setup.

Select Landscape orientation to print an electropherogram

These settings apply when you print Electropherogram View

Select this radio button to print the entire electropherogram on one page

LaserWriter Page Setup 7.0

Paper: US Letter R4 Letter
 US Legal B5 Letter Tabloid

Reduce or Enlarge: 100%

Orientation

Printer Effects:
 Font Substitution?
 Text Smoothing?
 Graphics Smoothing?
 Faster Bitmap Printing?

Chromatogram Settings

Single Page Variable Size

Number of Panes per Page 5

Number of Points per Pane 350

OK Cancel Options

Change these settings to specify custom print parameters. The options are described below.


3. Click Landscape Orientation, if you are printing Electropherogram View.
4. Change the Chromatogram Settings options as desired.

Select the radio button labeled “Single Page” to print the entire electropherogram on one page. Selecting the radio button labeled “Variable Size” causes the electropherogram to print on several pages, depending on the settings in the two entry fields.

In most cases the default settings are sufficient, although you can fine tune the print by changing the entry fields. They specify the number of times the electropherogram wraps down the page and the number of data points displayed within each wrap.

5. Choose Print from the File menu.

A dialog box appears.



6. Select the checkboxes indicating the views you wish to print.
 - If you select all four checkboxes, the first three views are printed in landscape orientation on a single page. The electropherogram is longer and prints on several pages.
 - If you select an individual view or only the first three, they print in portrait orientation.
 - Select the checkbox labeled “Allow for 3-hole punch” if you want to punch holes in the page and file it in a three-ring binder.
7. Click OK to start printing.

Saving the Batch Worksheet

When you save a batch worksheet, only the worksheet information (sequence names and identification parameters) is stored. Saving the worksheet does not save the Factura-identified features into the sample file feature tables.

IMPORTANT *If you want to preserve the results of processing a batch worksheet, you must either save the batch report (see page 5-7), or save the features to the sequence data files (see page 5-9), or both.*

If you submit a batch worksheet for processing, save it, then open it again, the indicator in the left column shows the sequences as unprocessed, indicating that you have not submitted the worksheet during the current use. If you previously saved Factura-identified features to the sample files, those features are still in the feature tables of the files.

The File menu contains three options for saving a batch worksheet. All are briefly described here.

Saving a New Worksheet

You can save a newly created batch worksheet using the standard Macintosh Save command in the File menu. A standard dialog box allows you to name the file when you save it.


Saving As a Different File Name

The Save As command is another standard Macintosh command that allows you to save the worksheet under a file name different from the original. When you use Save As, the name of the currently open worksheet is changed to the new file name.

Saving a Copy

The Save a Copy In command also saves results under a file name other than that of the original batch worksheet. The difference from the Save As command is that the name of your current worksheet is not changed. A copy is saved to the file you name, but the original remains on the screen.

Just as with the other two Save commands, when you choose this option to initiate saving, the dialog box shown below allows you to assign the file name and location.



After you click Save, a status indicator appears while the file is being saved.

Saving Results in a Batch Report

After you have processed a batch worksheet in Fatura, you can generate a batch report that shows the following:

- Vector, Site, Primer, and Confidence range (Limits)
- Vector range, Ambiguity range, Confidence range, Clear range.
- Clear length, Original length
- %Ns (ambiguities)


The batch report is fully described on page 3-22.

To save the batch report:

Save the results of batch processing by making the batch report window active and choosing any of the three Save commands from the File menu (see page 5-5).

Note *It is a good idea to save the batch report in a different folder from the sequences. This allows you to use the Add All button to add sequences to a file, such as a batch worksheet, without accidentally adding the batch report at the same time.*

If you close the batch report before you save it, the following dialog box appears:



- To save the report in a batch report file, click Save. The program asks you for a file name and location.
- To continue to close the report without saving it, click Don't Save.
- To cancel closing the report, click Cancel.

To open a saved batch report file:

You can open a batch report and review the results at a later time in two ways.

- If Factura is active, choose Open Batch Report from the Worksheet menu.
- If the Finder is active, double-click the icon of the batch report you wish to open. This starts Factura automatically and displays the batch report just as it was saved to the file.



Batch Report Icon

Note *Saved batch reports are tab-delimited text files and can be opened by many text editor programs.*


Saving Results to Sequence Data Source Files

If you are satisfied with the results from processing your data files, the Save to Sequence File command in the Worksheet menu puts the results into Sequence feature tables. You must do this to use the results in other software programs or to view and edit the newly identified features in Feature View (refer to page 4-11).

To save results to sequence files:

1. Choose Save to Sequence File from the Worksheet menu.

The following dialog box appears.



Note *This dialog box automatically appears after the batch worksheet is processed if you select the checkbox labeled "Automatically save to sequence file" when you specify settings for the batch worksheet (refer to page 3-3).*

The default stores batch processing results into the feature table of each sequence on the worksheet. You can then see them in Feature View of the sequence window (illustrated in Figure 5-1 on page 5-10).


Feature key:	Range(s):	Description:
Δ *ABI_Vector	1 17	This range of bases indicates the vector
Δ *ABI_Ambiguity	1 3	This range of bases indicates ambiguity
+ *ABI_Limits	465 472	This is the confidence range
*ABI_Multibase	1 450	
*ABI_Multibase	2	0.89 M
*ABI_Multibase	3	1.00 V
*ABI_Multibase	14	0.51 V
*ABI_Multibase	420	0.79 B
*ABI_Multibase	423	0.74 B
*ABI_Multibase	435	0.66 B
*ABI_Multibase	440	0.90 B
*ABI_Multibase	465	0.72 M
*ABI_Multibase	472	0.54 M
*ABI_Multibase	473	0.57 B

Figure 5-1. Feature table in Feature View of sequence window

2. Click the checkbox labeled “Overwrite existing feature table” to specify that the features identified by Factura be written over the existing feature table in the sequence file.

IMPORTANT *Any special features you might have previously added to the sequence are overwritten if you select this option.*

If you do not select the checkbox, the following dialog box is displayed for each sequence in the worksheet.



This dialog box allows you to overwrite the feature tables of specific sequences, while maintaining the existing information in others. The option is useful if you have added special features to a sequence, then reprocessed the batch worksheet using the same parameters as for a previous submission (for instance, if you have added more sequences to the worksheet). You can choose not to overwrite the sequence containing the special feature, because the Factura-identified features already exist in it from a previous submission.

Note *Overwriting the feature table does not destroy your original data. The original data and electropherogram are stored separately from edited versions and can be easily retrieved.*

3. Select the radio button labeled “Save to text file(s)” to save the data in each sequence to a text file.

If you save to text files when your original data is in sample files, text files are created that start with the original name and have “.seq” appended. These text files contain only the text of the sequence; they do not contain feature tables, electropherograms, or annotation data.

IMPORTANT *When you save data to a text file, the portions of the sequence outside of the confidence range are deleted. Only confidence range data is saved to a text file.*

The options allow you to (in order):

- Store the information in the original text files, overwriting the previous information.
- Rename the existing text files by adding “.old” to the file name, and save the data to a file with the same name as the existing file.
- Save the file in the current folder, where the existing text files are stored.
- Select a different folder, if the text files are not stored in the current one.

Saving Changes Made in the Sequence Window


If you have edited any bases or marked features in the sequence window, you should save the changes before closing the sequence window.

To save changes made in the sequence window:

1. If the window in which you made changes is not the active window, click it to make it active.
2. Choose Save (⌘-S) from the File menu.

The changes you made in that sequence window are saved to the sequence file.

If you click the close box in the sequence window without saving the changes, the following dialog box allows you to save them if you wish.



- To save the changes made in the sequence window, click Save.
- To close the window without saving changes, click Don't Save.
- To cancel closing the window, click Cancel.

Exporting Results to Other Formats

You should use the Export command (⌘-E) only after you have saved the results of Factura processing into the source files (see *Saving Results to Sequence Data Source Files* on page 5-9). The command is available when a sequence window is the active window.

The Export command exports the contents of the sequence file displayed in the active window into any of three formats:

- *Sample* files contain sequence data in four different formats: Sequence View presents a listing of data; Annotation View may contain annotation from an ABI PRISM sequencer; Feature View conveys the Factura-generated results to downstream programs in the form of defined vector, ambiguity, confidence range, and heterozygote/IUB code features; and Electropherogram View displays ABI PRISM sequence data in electropherogram format.
- *Inherit* files are compatible with the INHERIT software. When opened within Factura they contain only Sequence, Annotation and Feature Views.
- *Text* files contain the same views as the INHERIT type files when you open them in Factura. When you export from Factura to a text file, the file contains only the text of the sequence; it does not contain a feature table, electropherogram information, or annotations.

Note When you export to a text file using *Export* in the *File* menu, the entire sequence is exported. If you save to a text file using *Save to Sequence Files* in the *Worksheet* menu, the portions of the sequence outside of the confidence range are deleted.

6 Setting up Factura Libraries

This section explains the difference between Main and Custom libraries and describes how to set up custom libraries of vectors, enzymes, and primers to be used with your Factura processing.

Note *If you have previously set up libraries using an earlier version of Factura, they are still valid. You need not repeat the setup.*

Contents

Library Files	6-3
Main Libraries	6-4
Customized Libraries	6-4
Setting Up the Libraries	6-5
Setting up the Vector Library	6-5
Setting up the Enzyme Library	6-7
Setting up the Primer Library	6-10

Library Files

Factura provides three libraries (Vector, Enzyme, and Primer) to be used as resources with a batch worksheet. The following pages describe the library files and how to set up the libraries so they are customized for your purposes.

When you install Factura, the library files are automatically installed in the ABI folder inside the System folder on your hard disk, as shown in Figure 6-1.




Figure 6-1. Location of Library files

The SeqEdRanks file shown in Figure 6-1 is used during vector searching. The other files are library files. Three of the library files are standard libraries provided

with Factura. The other two are files that will contain only the libraries you choose to have immediately available to you.

IMPORTANT *If you remove, rename or replace any of the library files, the sequences in those files are no longer available for use with Factura.*

Main Libraries

The following files contain the main libraries provided with Factura:

- *VecBase* contains the Factura-provided library of vectors, which consists of more than 140 vectors. From it you can choose the subset that becomes your customized library.
- *PrimerBank* contains the library of primers provided with Factura. The Primer library is small, so it requires no set up unless you wish to add your own special primers.
- *All Commercial* contains the Factura-provided list of over 300 enzymes. From it you can choose the subset that becomes your customized enzyme library.

Using Factura with the large vector and enzyme libraries can be cumbersome, so you are allowed to set up customized vector and enzyme libraries, choosing small subsets that represent the options you use most frequently.

Customized Libraries

By setting up libraries in Factura, you create customized libraries that contain only the vectors and enzymes you use most frequently. The subsets you select provide the alternatives in the Settings dialog box when you specify Factura settings (see Section 3). Because the primer library is small, you need only set it up if you wish to add primers to it.

When you set up the Vector and Enzyme libraries, the following files hold your specified choices:

- *VecSelect* in the Libraries folder stores the sequences you choose to be a part of your customized vector library.
- *All Commercial Select* in the Restriction Enzymes folder stores the enzymes you choose for your customized enzyme library.

Setting Up the Libraries

Setting up the Vector Library


The vector library you set up should contain any vectors that you might want to identify using Factura. You can copy vectors from the Factura VecBase library or from your own files. Once you have added them, you can also rename them, as described below.

To set up the Vector library:

1. Choose Vector Library Setup from the Library Menu.

The following dialog box appears.

This list contains your customized list of vectors that are stored in the VecSelect file described on page 6-4



-
2. Perform any of the following you wish:

To copy vectors from the VecBase library:

- a. In the left (All Vectors in VecBase) list, select the vector you wish to copy.

As soon as you select a vector, the vector sequence appears in the table in the bottom of the dialog box. This enables you to determine if the selected vector is the one you want.

- b. Click the Copy button in the center of the dialog box.

To remove a vector you don't want:

- a. Select the vector in the Vectors Used list (on the right side of the dialog box).

- b. Click Remove.

- c. Click Delete in the dialog box that appears.

To copy vectors from your own files:

- a. Click the Add button.

A standard file dialog box appears so you can choose your vector file.




Figure 6-2. Standard file dialog box

Note that this dialog box allows you to filter for a specific file type. Only the selected file types show in the dialog box.

- b. Open the folder that contains the vector you wish to add.
- c. Click the vector in the dialog box to select it.
- d. Click Add to add it to the Vectors Used list.

To rename a vector in the Vectors Used list:

- a. Select the vector you wish to rename.

The current name appears in the entry field under the Vectors Used list.

- b. Select the existing name in the entry field and type a new name.

You may enter a name with up to 20 characters.

- c. Click Rename.

3. When you have finished, click OK to accept your customized library.


Setting up the Enzyme Library

Use Enzyme Library Setup to set up a list of the enzymes you plan to use with the batch worksheet. You can add enzymes from the All Commercial library listed in Factura. Your customized library is stored in the All Commercial Select file described on page 6-4.

To set up the Enzyme Library:

1. Choose Enzyme Library Setup from the Library menu.

The setup dialog box appears, with entry fields similar to those in the Vector library setup dialog box.



2. Perform any of the following you wish.

To copy enzymes from the All Commercial library:

- a. Select the enzyme you wish to copy from the All Commercial list.

The cloning sites appear in the table in the bottom of the dialog box, so you can determine if the selected enzyme is the one you want.

- b. Click the Copy button in the center of the dialog box.


To remove an enzyme:

- a. Select the enzyme you wish to remove from the Cloning Sites Used list.
- b. Click Remove.
- c. Click Delete in the dialog box that appears.

To create new enzymes:

- a. Click the New button.

The table at the bottom of the dialog box is cleared and a new file name (Untitled 1) appears in the Cloning Sites Used list.



- b. Type the appropriate sequence in the upper entry field of the table in the lower part of the dialog box. The complementary sequence appears beneath it as you type.
- c. Click above and below the table to place the arrows in the correct positions to indicate where the site cuts on the target sequence.
- d. Rename the new enzyme as described below.

To rename an enzyme in the Cloning Sites Used list:

- a. Select the enzyme you wish to rename.

The name appears in the entry field under the Cloning Sites Used list.

- b. Select the existing name in the entry field and type a name of your choice.

You may enter a name with up to 20 characters.

- c. Click Rename.

3. When you have finished, click OK to accept your customized library.


Setting up the Primer Library

The Primer Library setup dialog box differs from the other two setup dialog boxes in that it contains only a single list you can use to add, remove, or rename primers.

To set up the Primer Library:

1. Choose Primer Library Setup from the Library menu.

The following dialog box appears.



2. Review the Primer list.
3. Add or rename primers as follows:

To add a primer:

- a. Click the Add button.

A standard file dialog box like the one shown in Figure 6-2 on page 6-6 appears. Note that the dialog box allows you to filter for the type of file you add. Only the selected file types show in the dialog box.

- b. Open the folder that contains the desired primer file.
- c. Click the primer you wish to add so it is highlighted.
- d. Click Add.

Note *If you add a primer you don't want, select the primer in the list, click Remove, then click Delete in the dialog box that appears.*

To rename a primer:

- a. Choose Primer Library Setup from the Library menu.
- b. Select the primer you wish to rename.

The current name appears in the entry field under the Primer list.

- c. Select the existing name in the entry field and type a new name of your choice.

You may enter a name with up to 20 characters.

- d. Click Rename.
4. When you have finished, click OK to accept your customized library.

7 Overview—Factura Menu Commands

This section provides:

- Tables with brief descriptions of the Factura main menus
- Cross references to other sections with more detail

These descriptions do not walk you through the use of the Factura program. See Section 3, *Working With the Batch Worksheet* through Section 6, *Setting up Factura Libraries* for more detail.

Contents

Factura Main Menu	7-3
Apple Menu	7-3
File Menu	7-4
Edit Menu	7-6
Library Menu	7-7
Sequence Menu	7-8
Worksheet Menu	7-10
Window Menu	7-12

Factura Main Menu

The Factura main menu bar is arranged into six menus. To perform any task, pull down the menu that is most related to the task and choose the appropriate command from the list under the menu.

Note *Instead of choosing a command from a menu using the mouse, you can choose a command using the Command-key equivalent listed on the menu next to the command. For example, you can choose the New Sequence command from the File menu by holding down the ⌘ key and pressing N.*

This section briefly describes each of the menus on the Factura Main Menu bar (Figure 7-1).





Figure 7-1. Factura Main Menu Bar

For greater detail on specific commands, please see the expanded descriptions cross-referenced in the following tables.

Apple Menu

The Apple menu (🍏) allows you to use desk accessories installed in the system file on your hard disk. It is a standard feature of Macintosh applications. Refer to the *Macintosh System Software User's Guide* for details.

File Menu



The File menu provides access to tasks that you do routinely with active windows and files. This includes opening a new file, opening and closing an existing file or window, saving your work, getting information about a file, printing and setting up the page, and quitting the program.

Table 7-1. File Menu Commands

Command	Description	See for detail
<i>New Batch Worksheet and Open Batch Worksheet...</i>	Open a batch worksheet for working with sequences. Use New Batch Worksheet to open a new blank batch worksheet or Open Batch Worksheet to open a previously stored one.	Section 3, page 3-10
<i>New Sequence</i>	Allows you to type a sequence or to move sequence data from another application using the clipboard and the Macintosh Copy (or Cut) and Paste commands.	Section 4, page 4-5
<i>Open Sequence...</i>	Displays a dialog box that allows you to open any of your sequence files, regardless of the format (Sample, GeneAssist, and/or text files).	Section 4, page 4-4
<i>Open Sequence</i>	Asks you to specify the type of sequence file you wish to open before presenting a list from which to choose.	Section 4, page 4-4
<i>Save a Copy in...</i>	Allows you to save a backup copy of your results in a file other than the original batch worksheet. The name of the batch worksheet that you are working with does not change.	Section 5, page 5-5
<i>Export</i>	Exports the contents of sequence files into ABI 373 sequence (sample) files, GeneAssist sequence files, or plain text files.	Section 5, page 5-13
<i>Close Save, Save As... Revert, Page Setup... Print..., Quit</i>	Standard Macintosh commands. See the Apple <i>System Software User's Guide</i> .	

Edit Menu

Edit	
Redo Cut	⌘Z
Cut	⌘H
Copy	⌘C
Paste	⌘U
Clear	
Select All	⌘A
Find...	⌘F
Find Again	⌘G

The Edit menu is available when a sequence window is active. It provides a means of changing sequences by cutting and pasting, and includes a quick way to select the entire contents of a window.

Table 7-2. Edit Menu Commands

Command	Description	See for detail
<i>Select All</i>	Selects the entire contents of an active sequence window.	Section 4, page 4-7
<i>Find...</i>	Allows you to locate a particular pattern (a specified sequence of bases or amino acids) in a sequence.	Section 4, page 4-21
<i>Find Again</i>	Allows you to continue the search begun with the Find command without having to re-type the information.	Section 4, page 4-25
<i>Undo (Redo), Cut, Copy, Paste, and Clear</i>	Standard Macintosh commands. See <i>Apple System Software User's Guide</i> .	

Library Menu




The Library menu contains commands for setting up customized libraries of the vectors, enzymes, and primers you plan to use with the Factura program.

Table 7-3. Library Menu Commands

Command	Description	See for detail
<i>Vector Library Setup...</i>	Allows you to set up a list of vectors, drawn from either the VecBase library or your own files, for later use with a batch worksheet.	Section 6, page 6-5
<i>Enzyme Library Setup...</i>	Allows you to create a list of enzymes, selected from a library of all commercially available enzymes, for later use with a batch worksheet.	Section 6, page 6-7
<i>Primer Library Setup...</i>	Allows you to create a list of primers for later use with a batch worksheet. You can add to an existing list.	Section 6, page 6-10

Sequence Menu



The Sequence menu is available when a sequence window is active. It includes options for changing how the sequence is displayed and editing the features of the sequence.

Table 7-4. Sequence Menu Commands

Command	Description	See for detail
<i>Complement</i>	Replaces the sequence in the selected sequence window with its complement.	Section 4, page 4-25
<i>Options...</i>	Provides pull-down menus for changing the character display (font, size, grouping) in the sequence window.	Section 4, page 4-26
<i>Feature</i>	Allows you to add, modify, join, or remove features in the feature table of the displayed sequence.	Section 4, page 4-16
<i>Mark Style</i>	Enables you to change the style (color and border) applied to a feature or to assign a style for the first time.	Section 4, page 4-27
<i>Edit Mark Styles...</i>	Lets you change the name, color, or border assigned to any of the eight styles that can be applied to features.	Section 4, page 4-28
<i>Electropherogram</i>	Is only available when a sequence window is active and the electropherogram is displayed. The command allows you to zoom in or out, display the original sequence as a reference for edits, or display a full view of the electropherogram.	Section 4, page 4-32
<i>Display Option...</i>	Lets you change the display parameters for the electropherogram view. You can change the ruler indexing, show relative or real values for the vertical display, and selectively turn off the display of individual base electropherograms.	Section 4, page 4-38

Worksheet Menu


Worksheet	
Settings...	⌘H
Add Sequences...	⌘T
Remove Sequence	⌘R
Show Sequence	⌘D
Fill Down	
Submit	⌘Y
Save to Sequence Files...	
Open Batch Report...	

The Worksheet menu provides the commands that control the batch worksheet.

Table 7-5. Worksheet Menu Commands

Command	Description	See for detail
<i>Settings...</i>	Allows you to define the parameters most commonly used for processing the sequences with Factura (for example, vector sequence, ambiguity, confidence range). You can selectively enable or disable identification of vector, ambiguity, confidence range, or IUB/heterozygote calling.	Section 3, page 3-3
<i>Add Sequences...</i>	Adds the sequences you select (ABI PRISM sample files, GeneAssist files, or text files) to the batch worksheet for processing.	Section 3, page 3-14
<i>Remove Sequence</i>	Allows you to remove one or more selected sequences from the batch worksheet.	Section 3, page 3-17
<i>Show Sequence</i>	Opens a separate sequence window for each sequence you have selected on the batch worksheet. The sequence window allows up to four views of the data: Sequence, Annotation, Electropherogram, and Feature views.	Section 4, page 4-7
<i>Fill Down</i>	Allows you to quickly apply a parameter change in the batch worksheet to all sequences listed on the worksheet.	Section 3, page 3-18
<i>Submit</i>	Processes the sequences in the batch worksheet according to the parameters you have specified.	Section 3, page 3-21
<i>Save to Sequence Files...</i>	Allows you to save batch results back to sequence files. You can save to sample files (feature tables are updated) or to text files.	Section 5, page 5-9
<i>Open Batch Report...</i>	Allows you to review a saved batch report. The report displays vector range, ambiguity range, confidence range, clear range, clear length, % N, and original length, as well as the parameters used for processing.	Section 5, page 5-8

Window Menu



This menu lists all open windows associated with the current Factura program. Choose a window name from the menu to make the window that displays that data the active window. The active window is preceded by a check mark.

Index

A

ABI folder 2-3, 6-3
Add Sequence(s) command 7-11
adding
 bases in Electropherogram View 4-16
 features 4-16
 primers to primer library 6-11
 ranges to features 4-18
 sequences to batch
 worksheet 3-14, 7-11
All Commercial (Primers) library 6-7
All Commercial file 6-4
All Commercial Select file 6-4
ambiguity
 designating in settings 3-3, 3-6
 editing on batch worksheet 3-20
ambiguous bases, editing using
 Electropherogram View 4-13
Annotation View
 button to display 2-9
 described 2-9, 4-10
 example 2-9
 in text, Inherit, new sequences 4-9
Apple menu 7-3
automatic options
 Revert the sequences 3-9
 Save to Sequence File 3-9

B

bases
 adding in Sequence View 4-14
 changing in Electropherogram
 View 4-16
 changing in Sequence View 4-14
 deleting in Electropherogram View 4-16
 deleting in Sequence View 4-14

 editing in sequence window 4-14 to 4-16
 multiple positions in
 electropherograms 4-15
 selecting in electropherograms 4-15 to
 4-16
batch report
 closing 3-23
 described 2-5, 3-23
 example 3-22
 icon 3-24, 5-8
 opening 3-24, 5-8, 7-11
 printing 5-3
 saving 3-23, 5-7
 viewing 3-21
batch worksheet
 adding sequences 3-14, 7-11
 changing column width 3-17
 closing 3-25
 column description 3-11
 column selection 3-16
 description of 2-5, 3-1, 3-11
 editing parameters 7-11
 icon 3-11
 moving around 3-16
 opening 3-10, 7-5
 opening sequence windows from 4-7
 removing sequences 3-14, 3-17, 7-11
 reprocessing 4-19
 saving 3-26, 5-5, 7-5
 selecting sequences 3-16
 setting parameters 3-3, 3-10, 7-11
 submitting for processing 3-21, 7-11
 Worksheet menu commands 7-10
borders, marking features 4-27

buttons

- to display Annotation View 2-9
- to display Electropherogram View 2-11
- to display Feature View 2-10
- to display Sequence View 2-8

C

Caution, purpose of 1-5

changing

- bases in Electropherogram View 4-15
- bases in Sequence View 4-14
- feature appearance 4-27
- feature range 4-16
- font and grouping of a sequence 4-26
- marking style 4-28
- scaling of electropherograms in sequence window 4-33

characters

see also bases; symbols

Clear command 7-6

cloning sites. *see* enzymes

Close command 7-5

closing

- batch report 3-23
- batch worksheet 3-25
- sequence windows 4-42

codes, IUB. *see* IUB codes

color

- changing marking style of features 4-28
- marking features 4-27

columns

- described (batch worksheet) 3-11
- moving around batch worksheet 3-16
- widening 3-12, 3-17

command key 4-20

command-key equivalents (hot keys) 7-3

commands

- cross-reference tables 7-3
- Edit menu 7-6
- File menu 7-4
- keyboard 7-3
- Library menu 7-7
- Sequence menu 7-8
- tables of menu commands 7-3 to 7-11
- Worksheet menu 7-10

complement

- displaying for sequence 4-25

Complement command 7-9

complementing a sequence 4-25, 7-9

confidence range

- designating in settings 3-3, 3-7

Copy command 7-6

copying

- a sequence from a text file 4-7
- enzymes 6-8
- Save a Copy In command 5-5, 7-5
- vectors 6-6

creating

- new enzyme for library 6-9
- new sequence 4-5, 7-5

cross-reference tables 7-3 to 7-11

Cut and Paste commands 7-6

D

data

- electropherogram in sample file 2-4
- original data 2-4
- reverting to original sequences 2-4, 3-9
- saving to sequence files 5-9
- showing original on
 - electropherogram 4-32
 - unmodified in sample file 2-4

defaults

- marking styles 4-27

deleting
 bases in Electropherogram View 4-16
 bases in Sequence View 4-15

diagram
 Factura input and output files 2-5
 location of library files 6-3

Display Option... command 7-9

displaying
 complement of a sequence 4-25
 sequence windows 4-3

DNA Technical Support 1-6

E

Edit Mark Styles command 7-9

Edit menu

 Factura 4-3, 7-6
 not available in Electropherogram
 View 4-15

editable data 2-4

editing

 ambiguity column on batch
 worksheet 3-20
 bases in Electropherogram View 4-15
 bases in sequence window 4-14 to 4-16
 the batch worksheet 7-11
 features 4-16
 marking style 4-28
 protecting sequence from edits 4-8
 saving changes in the sequence
 window 4-42

Electropherogram command 7-9

electropherograms

 adding bases 4-16
 button to display 2-11
 changing appearance 4-32, 7-9
 changing scaling in Electropherogram
 View 4-33
 deleting and changing bases 4-16

 editing 4-15
 Electropherogram View 2-11, 4-13
 finding patterns 4-21
 in files 2-4
 in Sample files only 4-9
 multiple base positions 4-15
 selecting bases 4-15, 4-16
 showing original data 4-32
 turning off display 4-38
 zooming in sequence window 4-35, 7-9

ellipsis (...)

 creating for Find command 4-24

Enzyme Library, setting up 6-7

enzymes

 copying 6-8
 creating new 6-9
 renaming 6-9

Export command 5-13, 7-5

exporting

 Factura results 5-13

F

Factura

 description 1-4
 input files 2-4
 output files 2-4
 saving results 5-7
 settings 3-3
 windows and views 2-7

Feature command 7-9

feature key 4-18

feature tables

 adding features to 4-17
 Factura results in 5-5
 marking features 4-27
 saving to 3-21, 5-9
 see also features; Feature View

Feature View

- button to display 2-10
- described 2-10, 4-11
- example 2-10, 4-11
- in text, Inherit, new sequences 4-9
- marking features 4-27

features

- adding ranges 4-18
- adding, removing, or changing 4-16
- changing appearance 4-27
- editing 4-16, 7-9
- identified by Factura 3-21
- joining 4-20
- marking 4-18, 4-27, 7-9
- overwriting 5-10
- removing ranges 4-18
- special 4-16, 5-10
- see also* feature tables; Feature View

File menu

- Factura 7-4

files

- All Commercial file 6-4
- All Commercial Select file 6-4
- batch report file 2-5, 3-24, 5-8
- batch worksheet file 2-5, 3-11
- copying text file to new sequence 4-7
- diagram of Factura input and output files 2-5
- editable data in 2-4
- electropherogram data in 2-4
- Factura input and output files 2-4
- INHERIT 4-5
- opening. *see type of file* (batch worksheet; sequences)
- PrimerBank file 6-4
- sample files. *see* sample files
- saving a copy 5-5
- saving as text 5-11
- saving automatically 3-9

- text files. *see* text files
- unmodified data in 2-4
- VecBase file 6-4
- VecSelect file 6-4

Fill Down command 3-19, 7-11

Find command 4-22

finding

- a range of bases 4-23
- Find Again command 4-25, 7-6
- Find command 7-6
- IUB codes 4-22
- patterns 4-21
- selection expressions 4-22

folders

- ABI 2-3, 6-3
- GeneAssist 2-3
- system 6-3

font, changing 4-26

G

GeneAssist folder 2-3

grep 4-22

H

help. *see* technical support

heterozygote

- identification threshold 3-7

hot keys 7-3

How to get help 1-6 to 1-12

I

icons

- batch report file 3-24, 5-8
- batch worksheet file 3-11
- diamond shapes by sequence names 3-12

identification parameters

- designating in settings 3-3

IMPORTANT, purpose of 1-5

INHERIT files
 description 4-5, 5-13
 opening 4-4
input files
 Factura 2-4
 Factura diagram 2-5
installation, top level folder 2-3
IUB codes
 designating threshold 3-7
 finding 4-22
 identified only for sample
 files 3-8, 3-20, 3-21
IUPAC 4-22

J
joining features 4-20

K
key, feature 4-18
keyboard keys
 command key 4-20
 Shift key 4-20
 Tab key 3-11

L
Left-Arrow key
 moving in electropherograms 4-15
libraries
 All Commercial (Primers) 6-7
 Enzyme 6-7
 Primer 6-4, 6-10
 setting up 3-10, 6-5, 7-7
 VecBase 6-5
 Vector 6-5
Library menu commands 7-7
lock image 4-8

M
main menu
 bar 7-3
 tables of commands 7-3
mark style
 choosing 4-18
 defaults 4-27
 editing 4-28, 7-9
Mark Style command 7-9
marking
 features 4-18, 4-27, 7-9
menus
 Apple menu 7-3
 cross-reference tables 7-3
 Edit menu commands 7-6
 File menu commands 7-4
 Library menu commands 7-7
 main menu bar 7-3
 Sequence menu commands 7-8
 tables of commands 7-3
 Window menu 7-12
 Worksheet menu commands 7-10
multiple base positions 4-15

N
New Batch Worksheet command 7-5
New Sequence command 4-6, 7-5
Note, purpose of 1-5

O
Open Batch Report command 7-11
Open Batch Worksheet command 3-11, 7-5
Open Sequence commands 4-4, 7-5
opening
 batch report 3-24, 5-8, 7-11
 batch worksheet 3-10, 7-5
 sequence from the batch worksheet 4-7
 sequence in Factura 3-14
 sequence window 4-3, 7-5

Option key
 creating ellipsis 4-24
 moving cursor in
 electropherograms 4-15

Options... command 7-9

output files
 Factura 2-4
 Factura diagram 2-5

overwriting features 5-10

P

Page Setup... command 7-5

parameters
 identification 3-4

Paste command 7-6

patterns
 finding in electropherograms 4-21
 finding in sequences 4-21

primer
 adding to primer library 6-11
 editing on batch worksheet 3-19
 renaming 6-11

Primer Library 6-4, 6-10

PrimerBank file 6-4

Print command 7-5

printing
 batch report 5-3
 menu commands for 7-5
 sequence window views 5-3

processing
 batch worksheet 3-21
 exporting Factura results 5-13
 saving results 5-9

program files
 Factura 2-3

Q

Quit command 7-5

R

range
 adding to features 4-18
 changing range of a feature 4-19
 finding a range of bases 4-23
 removing from features 4-18

Redo command
 Factura 7-6

Remove Sequences command 7-11

removing
 bases in Factura 4-14
 features 4-16
 ranges from features 4-18
 sequences from batch
 worksheet 3-14, 3-17, 7-11

renaming
 enzyme 6-9
 primer 6-11
 vector 6-7

reports
 batch. *see* batch report

results
 the batch report 3-22
 exporting 5-13
 in feature tables 5-5
 saving batch report 5-7
 saving to sequence files 5-9
 viewing Factura results 3-23

Revert command 7-5

reverting sequences to original data 3-9

Right-Arrow key
 moving in electropherograms 4-15

S

sample files

- described 2-4, 4-5, 5-13
- IUB/heterozygote
 - assignment 3-8, 3-20, 3-21
- opening 4-4
- restoring original data 3-9
- saving changes 5-9
- viewing information from run 4-10

saving

- automatically to sequences 3-9
- batch report 3-23, 5-7
- batch worksheet 3-26, 5-5, 7-5
- Factura results to sequence
 - files 5-9, 7-11
- modifications made in sequence
 - window 4-42
- Save a Copy in... command 7-5
- Save As... command 7-5
- Save command 7-5
- Save to Sequence File(s) command 7-11
 - to feature tables 5-9
 - to text files 5-11

scaling

- electropherograms in the sequence
 - window 4-33

Select All command 4-7, 7-6

selecting

- sequences 3-16

selection expressions, tables 4-23

Sequence menu

- available with sequence window 4-3
- commands 7-8

Sequence View

- button to display 2-8
- changing font and grouping 4-26
- described 2-8, 4-9
- example 2-8, 4-3
- in text, Inherit, new sequences 4-9
- marking features 4-27
- See also* sequence windows

sequence windows

- changing bases 4-14
- closing 4-42
- deleting bases 4-14
- described 2-7, 4-1, 4-8
- editing 4-14, 7-9
- opening 4-3, 7-5, 7-11
- printing 5-3
- saving modifications 4-42
- viewing options 4-21
- Views 4-9
- zooming electropherograms 4-33 to 4-37

sequences

- adding bases in Sequence View 4-14
- adding to a batch worksheet 3-14, 7-11
- changing bases in Sequence View 4-14
- changing font and grouping 4-26
- closing the sequence window 4-42
- complementing 4-25, 7-9
- copying from text files 4-7
- creating new 4-5, 7-5
- deleting bases in Sequence View 4-15
- displaying complement 4-25, 7-9
- editable data in sample file 2-4
- editing in Electropherogram View 4-15
- editing in Sequence View 4-14
- exporting to other formats 5-13
- Factura results in feature tables 5-5
- finding patterns 4-21

sequences (*continued*)

- locking 4-8
 - opening in Factura 4-7
 - opening sequence windows 3-14, 4-3
 - processing in Factura 3-21
 - protecting from edits 4-8
 - removing from batch
 - worksheet 3-14, 3-17
 - reverting to original data 3-9
 - saving to sample file 5-9, 7-11
 - selecting in Factura 3-16
 - showing original 4-32
- service specialists 1-6
- setting up
 - batch worksheet 3-3, 3-10, 7-11
 - Enzyme library 6-7
 - libraries 6-5, 7-7
 - Primer library 6-10
 - Vector library 6-5
- settings
 - designating ambiguity 3-6
 - designating confidence range 3-7
 - designating IUB code threshold 3-7
 - designating vector sequence 3-4
 - dialog box options 3-4
 - set up libraries before 3-3
- Settings command 7-11
- Shift key 4-20
- Show Sequence command 4-7
- Show Sequence(s) command 7-11
- special text
 - Caution 1-5
 - IMPORTANT 1-5
 - Note 1-5
 - WARNING 1-5
- specialized features 4-16, 5-10
- Style menu 4-18

styles

- changing definition 4-28
 - defaults for marking features 4-27
 - marking features 4-18, 4-27
- Submit command 7-11
- submitting a batch worksheet 3-21, 7-11
- summary graphic 4-8, 4-18, 4-25, 4-27
- support. *see* technical support
- symbols
 - see also* characters; bases
- system folder 6-3

T

- Tab key 3-11
- tables
 - batch report columns 3-23
 - command cross-reference 7-3
 - default marking styles 4-27
 - Edit menu commands 7-6
 - File menu commands 7-4
 - Library menu commands 7-7
 - main menu commands 7-3
 - selection expressions for Find
 - command 4-23
 - Sequence menu commands 7-8
 - Worksheet menu commands 7-10
- technical support 1-6 to 1-12
 - contacting 1-6
 - service and applications support 1-6
 - telephone numbers 1-7 to 1-10
- text files
 - and feature table information 4-4
 - copying to sequence files 4-7
 - described 2-4, 4-5, 5-13
 - opening 4-4
 - saving to 5-11
- threshold
 - IUB/heterozygote identification 3-7

U

Undo command
 Factura 7-6

V

VecBase file 6-4
VecBase library 6-5
VecSelect file 6-4
vector
 copying 6-6
 designating in settings 3-4
 editing on batch worksheet 3-19
 renaming 6-7
Vector Library, setting up 6-5
viewing
 Factura results 3-23
 options in sequence window 4-21
views
 associated with file types 4-5
 Factura windows and views 2-7
 marking features 4-27
 sequence window views 4-9

W

WARNING, purpose of 1-5
Window menu 7-12
windows
 batch report 2-7
 batch worksheet 2-7
 closing the sequence window 4-42
 Factura windows and views 2-7
 sequence window. *see* sequence windows
word processing
 opening batch reports in 5-8
Worksheet menu commands 7-10
Wrap checkbox 4-25

Z

zooming
 electropherograms in the sequence
 window 4-33 to 4-37

Headquarters

850 Lincoln Centre Drive
Foster City, CA 94404 USA
Phone: +1 650.638.5800
Toll Free: +1 800.345.5224
Fax: +1 650.638.5884

Worldwide Sales Offices

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches into 150 countries on six continents. For international office locations, please call our local office or refer to our web site at www.appliedbiosystems.com.

www.appliedbiosystems.com



Applera Corporation is committed to providing the world's leading technology and information for life scientists. Applera Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Printed in the USA, 07/2001
Part Number 903696 Rev. B

an **Applera** business