

Stabilizer® T1 Quick Guide


- 1 Turn on power and insert a USB flash drive



- 2 System ready after 3-5 min



- 3 Enter ▶ Sample name
▶ Maintainer® card number
Select ▶ Method



- 4 ▶ Instrument ready
▶ Start ex vivo time and immediately collect the sample




- 5 Place sample centered in Maintainer card



- 6 Close Maintainer card



- 7 Place Maintainer card in Stabilizer T1



- 8 Press start



- 9 Green illumination = stabilization finished



- 10 Remove Maintainer card



- 11 The sample is ready for further preparation/analysis (see back page)

or

long time storage in freezer

- 12 Remove the USB flash drive with the log file.
View the log-file on any computer with Microsoft Excel.



Please refer to Stabilizer T1 User Manual for safety instructions and more detailed information.

Stabilizor® T1 Quick Guide

Sample stabilization methods

- ▶ **Structure Preserving**
For preserving the anatomical structures of the sample. Measures sample and adjusts the heating time.
- ▶ **Fresh Auto**
For most fresh samples. Measures sample and adjusts the heating time.
- ▶ **Frozen Auto**
For frozen samples. Measures sample and adjusts the heating time (includes extra time due to low starting temperature).
- ▶ **Custom**
Manual mode where all parameters can be set by the user.

Advice for multiple samples

- ▶ **Similar *ex vivo* time** between samples in an experimental series minimizes result variation due to e.g. degradation.
- ▶ **Plan your experiment** to combine short *ex vivo* time and careful sample collection. Perform sample collection efficiently, but without stress.
- ▶ **Make all preparations before starting** sample collection, i.e. prepare Maintainor cards, enter data into Stabilizor T1 instrument.
- ▶ **Perform a mock run before starting** to ensure that everything is prepared before risking any samples.

Maintainor Tissue

— for tissues, up to 7 mm thick

Ex vivo time

The time the sample has been outside the living organism.

Further sample preparations

After the stabilization, the sample can, for example, be:

- ▶ Dissected
- ▶ Homogenized for further analysis
- ▶ Microtome sectioned for imaging
- ▶ Stored frozen for later preparation/analysis

Denator AB

Arvid Wallgrens Backe 20
SE 413 46 Gothenburg, Sweden
Phone: +46 (0)31 412 841
Fax: +46 (0)31 412 840
order@denator.com
support@denator.com

www.denator.com