



# **MPR 010**

# **Mobile Phase Recycler**

Version 1.00

## **User's Manual**

Revision A, May 2000

Copyright © 2000 LabAlliance.

All rights reserved. This manual or any of its parts must not be electronically, mechanically or photographically copied, reproduced or stored otherwise without permission.

## TABLE OF CONTENTS

<b>TABLE OF CONTENTS</b> .....	<b>1</b>
<b>LIST OF FIGURES</b> .....	<b>2</b>
<b>INTRODUCTION</b> .....	<b>3</b>
SPECIFICATION .....	4
<b>INSTALLATION</b> .....	<b>6</b>
<i>Liquid connection</i> .....	6
<i>Electrical connection</i> .....	6
<i>Connection with computer (optional)</i> .....	8
<b>HOW MPR 010 WORKS?</b> .....	<b>9</b>
MOBILE PHASE RECYCLING .....	9
INPUT RANGE .....	9
AUTO ZERO VERSUS ACTIVEZERO™ .....	9
DELAY AND TRANSPORT DELAY .....	9
DIRECT CONTROL MODE .....	10
FAIL SAFE DESIGN .....	10
MPR 010 AND GLP .....	10
MPR 010 FACTORY SETTING .....	10
<b>MPR 010 START UP</b> .....	<b>11</b>
MPR 010 CONTROLS .....	11
FIRST STEPS WITH THE MPR 010: .....	12
<b>CONTACTS FOR TECHNICAL SUPPORT</b> .....	<b>13</b>

## LIST OF FIGURES

<u>Figure 1: MPR 010 Mobile Phase Recycler</u> .....	3
<u>Figure 2: MPR 010 Position in the HPLC System</u> .....	4
<u>Figure 3: MPR 010 Operational Principle</u> .....	4
<u>Figure 4: Flange-Free Fitting</u> .....	6
<u>Figure 5: MPR 010 Rear Panel</u> .....	7
<u>Figure 6: MPR 010 Front Panel</u> .....	11

## INTRODUCTION

MPR 010 Mobile Phase Recycler (Figure 1) is a useful cost-saving device for isocratic HPLC. With efficient recycling of mobile phase, one can save up to 90% of solvent. Usually, most of the mobile phase is free of contaminants and need not be wasted. MPR 010 allows the user to redirect the uncontaminated mobile phase to the solvent reservoir, saving solvent and waste-disposal costs.



**Figure 1: MPR 010 Mobile Phase Recycler**

MPR 010 continuously monitors the detector signal. In response to signal level, the 3-way valve is switched to waste or recycle position. The switching valve is located just after the detector output (Figure 2). If the signal level detected is greater than the set signal threshold, the valve is switched to waste position. The switching may be delayed slightly to account for transport time from the detector cell to 3-way valve (Figure 3).




Figure 2: MPR 010 Position in the HPLC System




Figure 3: MPR 010 Operational Principle

MPR 010 can be connected to any HPLC detector with an integrator output (1 V/Full Scale). There are two major benefits of using the integrator output: 1) better resistance against radio-frequency interference, and 2) constant signal to voltage ratio independent of the actual range setting.

The analytical fitting kit (1/16") is a standard part of delivery. Optional semi-prep (1/8") fitting kit is also available.

### Specification

Input Range	$\pm 1V$
Threshold	0.1 - 99 % of input range (1 ÷ 990 mV)
Delay	1 - 99 s with 1 s increment
Transport Delay	0 - 99 s with 1 s increment
Autozero	100% of input range
Valve	3-way electrically actuated valve ¼" – 28 flat bottom fitting
Working Pressure	0.2 MPa (30 psi) maximum
Wetted materials	PTFE


Display	3 digit LED
Valve position indication	Red LED - WASTE Green LED – RECYCLE
Inputs (TTL compatible)	Analog signal, autozero, waste external control
Outputs	ActiveZERO™ (contact closure)
PC communication	RS 232, control program included
Power requirements	115 V/5 VA (power adapter, MPR 010 itself uses the safe voltage)
Environmental requirements	10 - 35 °C, up to 60% relative non-condensing humidity
Dimensions	7.9 × 5.5 × 2.2 inch (200 × 140 × 55 mm)
Weight	600 g

## INSTALLATION

It is easy to install MPR-010 for the experienced HPLC user. No special tools are required. Once the tubings are connected to the valve, and detector signal output is connected with MPR-010, the recycler is ready for use. An optional computer connection can be made.

### ***Liquid connections***

The recycler valve should be placed in the HPLC system according to Figure 2 (page 4). The Detector liquid output should be connected with the valve port "D", the Waste bottle with port "W", and mobile phase Reservoir with "R" port. Use the flange-free fittings (Figure 4) and capillaries from the start-up kit for the connections.



**Figure 4: Flange-Free Fitting**

### ***Electrical connections***

They are different ways to make electrical connections from MPR 010 to the chromatographic system. The schematic depends on your specific requirements (see following examples).



**Figure 5: MPR 010 Rear Panel**

Signal (input, yellow)	Input from the integrator output of HPLC detector (1Volt/Full Scale)
Zero (input, black)	External autozero command (e.g. injection valve, autosampler etc...)
Waste (input, red)	External signal controlled waste position
ActiveZERO™ (output, green)	Allows the synchronized zeroing of detector and MPR 010 prior to sample injection. When the ZERO button is pressed (or external zero input is shorted), the HPLC detector is zeroed immediately, after ActiveZERO™ Delay MPR 010 is zeroed itself.
PC	PC connector allows user to connect the MPR 010 to a computer with Control Panel Software. Additional functions are available with PC control.
Power	MPR 010 is powered by safe voltage from the power adapter.

### Connection Examples:

- 1) MPR 010 will be used independent of the other HPLC system components. Using the "Universal Cable" (P/N DY2K012507) connect the detector 1V/FS output with the yellow connector on the MPR 010 rear panel.
- 2) MPR 010 Auto-zero function should be activated at the moment of injection. Connect the analog signal as described in example 1. The black MPR 010 connector should be



connected to the contact closure of the manual injection valve or autosampler. MPR 010 will be zeroed at the injection.

- 3) Both detector and MPR 010 should be zeroed manually before an injection. Connect the detector signal output to the yellow MPR 010 connector. Green ActiveZERO™ output should be connected to the detector external zero input. ZERO button of MPR 010 cause the auto zero function on both instruments.

### **Connection with computer (optional)**

If advanced functions of MPR-010 are required, connect the recycler to a computer.

Advanced functions:

- Transport Delay setting
- Using of Delay Calculator
- Setting of external ZERO mode
- Detect negative peak setting
- GLP functions
- PQ functions (for authorized technicians only)

RS232 cable should be connected to the free communication port of your computer. COM1 – COM 4 can be used. Then install the *MPR 010 Control Panel Software*. The program runs under Windows 95/98/NT.

## HOW DOES THE MPR-010 RECYCLER WORK?

### Mobile phase recycling

MPR 010 continuously watches the detector signal. If the signal level exceeds the *Threshold*, the 3-way valve switches the mobile phase flow to the waste bottle. After the signal decreases to a value below the *Threshold* level, the valve still holds in the WASTE position for *Delay* time. Then, it switches to the RECYCLE position (and the mobile phase is recycled).

Depending on the chromatographic analysis, the optimal *Threshold* setting usually is between 0.2 - 1%. The *Delay* parameter should match the internal volume of tubing between detector and MPR 010. A useful tool for *Delay* determination is the *Delay Calculator* in the *MPR 010 Control Panel* software.

### Input range

MPR 010 Input Range is equal to the standard integrator output range of HPLC detectors,  $\pm 1V$ . MPR-010 analog input accepts both positive and negative peaks. However, the negative peak detection can be suppressed using MPR 010 Control Panel Options Setting.

### Auto Zero versus ActiveZERO™

The Auto-zero function is activated by pressing the ZERO button (or external ZERO signal). The current voltage monitored by the analog IN is then used as the MPR-010 baseline signal. ZERO also activates the ActiveZERO™ function: At the moment of ZERO the ActiveZERO™ contact closure is closed for the pre-defined period (3 seconds by default). This (ActiveZERO) contact closure when connected to the auto zero input of detector sets the detector output to zero. After this period when the detector signal is stabilized, MPR 010 sets itself to zero. The ActiveZERO™ function allows the user to synchronize the zero level of both detector and MPR 010 with a single button command.

### Delay and Transport Delay

The functions of *Delay* and *Transport Delay* are very similar. The *Delay* is applied after the peak, whereas *Transport Delay* prior to peak. The *Transport Delay* factory setting is 0 and can be changed from MPR 010 Control Panel only.

## Direct Control Mode

The basic method for switching valve position is signal dependent. In special cases MPR 010 can be switched directly to WASTE/RECYCLE positions. They are two ways: using front panel buttons WASTE and RECYCLE, or external signal WASTE. When the MPR 010 valve is under direct control, the position LEDs are blinking. The direct control mode can be escaped by pressing MODE button.

## Fail Safe Design

The MPR 010 was designed to recycle the mobile phase only if active and properly initialized by user.

- The normal position of the MPR-010 3-way valve is WASTE. In the case of power fail the valve remains in the WASTE position and the mobile phase remains uncontaminated.
- **After MPR-010 is Reset (Power On), the detector signal is ignored and the ZERO button must be pressed for initialization.**

## MPR 010 and GLP

MPR-010 GLP functions are available from *MPR 010 Control Panel* software. The Installation Qualification, Operational Qualification and Performance Qualification procedures and protocols are available upon request.

## MPR 010 Factory Setting

Threshold	0.5 %
Delay	3 s
Transport Delay	0 s
ActiveZERO™ Delay	3 s
External ZERO	negative edge (or contact closure)
Detect Peak	Both positive and negative peaks

## MPR 010 START UP

### MPR 010 Controls




Figure 6: MPR 010 Front Panel

Display	Display the current signal level (in % of full scale), Threshold and Delay. 000 indicates the zero state
Display mode LED	SIGNAL – signal level [%] THRESHOLD – switching level [%] DELAY – WASTE to RECYCLE switching delay [s]
WASTE/RECYCLE LED	Continuous light – input signal level dependent position blinking – direct control mode
MODE button	Switches the display mode
ZERO button	Activates the ActiveZERO™ output and performs the MPR 010 auto zero function
UP and DOWN	Allows the parameter setting

arrows	
WASTE button	Direct WASTE position
RECYCLE button	Direct RECYCLE position

## First steps with the MPR 010:

- Connect the power adapter
- **Press the ZERO button**
- Press MODE button to access the Threshold setting (THRESHOLD LED is lit)
- Using UP and DOWN arrows set the Threshold value. If you are new with using MPR 010, the default value of 0.5% is good choice.
- Press MODE button again to access the Delay setting. It can be calculated by formula:

$$\text{Delay}[s] = \text{Tubing volume [mL]} / \text{Flow Rate [mL/min]}$$

Delay Calculator in the MPR 010 Control Panel Software can be also used for the calculation.

- Press MODE button. Display now shows the input signal.

## CONTACT FOR TECHNICAL SUPPORT

Should you have further questions, comments or suggestions concerning the MPR 010, contact your distributor below.

LabAlliance  
349 N.Science Park Road  
State College  
PA 16803  
USA

Phone: 800-345-5557 (toll-free within U.S.)  
814-234-7317  
Fax: 814-236-6072  
Email: sales@laballiance.com

When contacting technical support, please provide the following information:

Product name	MPR 010 Mobile Phase Recycler
Serial number	Accessible on the instrument bottom