# Waters® Heat Sealer



**User Manual** 

Waters

# Waters® Heat Sealer

# **Contents Page**

Introduction
Safety
Declaration of Conformity
Caution and Safety Signs in Accordance with IEC 417
Technical Specifications
Environmental Requirements
Installation and Set-up
Operation Procedure (Sealing & Removing Seal)
Cleaning Your Thermo-Sealer 10
Guarantee 10
Contact Information
Technical Information
Circuit Diagram for 110v unit
Circuit Diagram for 240v unit
Fault Finding
Accessories

CAUTION: Avoid contact with the heater plate, which can cause serious injury.

#### Introduction



The Waters® Heat Sealer offers a convenient method for heat sealing collection plates made from various types of plastic.

The heat sealer provides heat to form a seal between the seal and plate, preventing evaporative loss during exposure to the atmosphere or for long-term storage.

## Safety

The Waters® Heat Sealer must only be used in accordance with proper safety standards and procedures, together with the instructions contained in this manual.

The unit has been designed and manufactured to conform to international safety specifications. Under normal conditions, once the aluminum block is up to temperature, the amber heater lamp will flicker on and off to maintain temperature. In the unlikely event of an over-temperature problem with the heat sealer, the unit is fitted with a thermal cut out which removes power to the heater should the block temperature exceed 220°C. If this is triggered the amber heater lamp will go out and will not come back on. The block will gradually cool to room temperature. The thermal cut out cannot be reset (for safety reasons) and the unit should be returned for repair.

If there is any doubt relating to the proper use of this equipment, refer to this manual or call your Waters® representative. Only properly trained personnel should service this equipment.

## **Declaration of Conformity**

The Waters® Heat Sealer, Model Numbers SP-0669/110 & SP-0669/220 have been designed in accordance with, and satisfies the requirements of, article 11 of the Low Voltage Directive 73/23/EEC as realigned by 93/68/EEC on the harmonization of the laws of the Member States relating to electrical equipment designed for use within certain voltage limits, to the essential requirements of BS EN 61010-1: 1993 & BS EN 61010-2-010:1995.

The Waters® Heat Sealer, Model Numbers SP-0669/110 & SP-0669/220 have been type tested by EMC Projects NAWAS approved test house, and issued a certificate of testing No 3418/98 to the following E.M.C. Standards:

EN 50081-1:	1992
EN 55022(B):	1994
EN 50082-1:	1997
EN 61000-4-3:	1997
EN 61000-4-2:	1995
EN 61000-4-4:	1995
EN 61000-4-5:	1995
EN 61000-4-6:	1996
EN 61000-4-11:	1994

Covering: Conducted Emissions, Radiated Emissions, Electrostatic Discharge, Radiated Immunity, Surge, Fast Transient Bursts and Voltage Dips satisfying the EMC Directive(s) 89/336/EEC and 92/31/EEC as realigned by 93/68/EEC.

# Caution and Safety Signs in Accordance with IEC 417

Caution



Refer to manual

Caution



Risk of electrical shock

Warning



Hot surface (this instrument is fitted with an internal heater with a maximum temperature of 220°C when switched on)

Only refit correct type of fuse

Must be IEC 127 approved for use in EU countries.

Must be CSA or UL listed or recognized for use in USA & Canada.

Mains Switch Symbols

I = ON

O = OFF

## **Technical Specifications**

Input Voltage 230 (nominal) Volts

110 (nominal) Volts

Input Frequency 47Hz - 63Hz

Power Consumption 110V Version 500 Watts

240V Version 400 Watts

Fuse Rating 5 Amp(T) 110 Volts

3.15 Amp(T) 240 Volts

Internal Heater 350 Watts (Ambient to 220°C max)

Instrument Housing Sheet Metal

Dimensions  $140 \times 330 \times 150$ mm

Instrument Weight 3.4kg

## **Environmental Requirements**

**Temperature Range** - Operating 10 to 35°C

- Storage -20 to  $40^{\circ}\text{C}$ 

**Relative Humidity** - Operating 10% to 80% non-condensing

- Storage 10% to 80% non-condensing

# **Installation and Set-up**



**WARNING -** Do not operate this instrument in an atmosphere containing explosive gases.

**WARNING -** Only approved power cord set must be used with this instrument.

**WARNING -** If it is required to use an extension lead, the lead must be earthed.


Before attempting to connect the unit to the electricity supply, or to operate the unit, please read the operating information and pay particular attention to the following:

1. Check the voltage marked on the unit's rating plate against that of the mains supply. The unit is supplied with a pre-wired plug and socket.

The fused plug supplied with the mains cable for use in the UK is fitted with the following value fuse to protect the cable.

Supply	Fuse Fitted in Plug	Fuse Fitted in Sealer
240/220V	5 Amp	3.15 Amp
120/110V	No Fuse	5 Amp

The fuses in the unit protect the unit and the operator. **Note:** The second fuseholder fitted to the rear of the sealer is a spare in the case of 110V units, but both fuseholders are operational in the 240V units. Note that units marked 240V on the rating plate work at 220V; units marked 120V on the rating plate work at 110V. In both these cases, however, the heating rate will degrade by approximately 16%.



- 2. Connect the unit to the mains electricity supply and switch on the power switch on the side of the unit. If electricity is being supplied to the unit, the green power light will illuminate. If the heater within the unit is receiving power then the amber light will also illuminate.
- 3. When switching on, it is necessary to leave the unit on for 20–25 minutes to allow the unit to come to the proper temperature before carrying out the first sealing operation.
  - On initial power-up, both green and amber lights will illuminate. After approximately 10–15 minutes the amber light will go out. **DO NOT USE** the Thermo-Sealer until 20-25 minutes has elapsed from power-up. By this time the amber light will blink on and off constantly. The amber light, blinking on and off, indicates the unit is at the correct working condition.
- 4. The heater plate is set to approximately 170°C. Take care to avoid skin contact with the heater plate and to prevent contact between the heater plate and any material which may be damaged or ignited.
- 5. In the event of spillage or other contamination of the unit, switch off and disconnect the unit from the mains power supply, and do not attempt to clean the unit until it has thoroughly cooled.
  If the equipment is used in a manner NOT specified by the manufacturer, the protection provided by the equipment may be impaired.

## **Operation Procedure**

#### Sealing

After installation, follow the procedures set out below:

Check that the unit is connected to an appropriate power supply and leave for 20–25 minutes to settle down to correct operating temperature. The green power lamp should be illuminated and the amber heater lamp flicker on and off.

Slide the plate to be sealed onto the platform of the sealer, making sure that it is flat between the two centralizing springs. You may also choose to slide the plate holder (black anodized metal plate) between the springs.

Place the sealing film over the plate, ensuring that the sealing surface is face down.

Place one hand above the other on the handle on the top of the unit and push the heater head down onto the seal and plate. Once the heater head touches, continue to push down to compress the springs that are behind the heated plate. Once these are fully compressed (a distance of 5mm) the force applied to the plate is 10kg. Each plate should be sealed in both orientations to achieve uniform sealing across the plate. Hold the heater head down for 3 seconds. Allow the heater head to return to the top of the track.



Rotate the plate (180 degrees – turn it end for end) in the unit and seal again for 3 seconds.

Allow the heater head to return to the top of the track.

Remove the plate.

**NOTE:** Foil backed film may remain hot for a number of seconds after sealing. Care should be taken in handling these plates.

#### **Removing Seal**

Depending on the seal type, samples can be accessed by either piercing the seal or peeling it off.

#### Seals

After piercing the seal, if storage and reuse are desired, peal off the seal and apply a new seal for best storage results.

## **Cleaning Your Thermo-Sealer**



Before using any cleaning or decontamination method other than that recommended by the manufacturer, users should check with the manufacturer that the proposed method will not damage the equipment.



**Do not** use acetone or abrasive cleaners.

The Thermo-Sealer cover can be cleaned with a cloth dipped in water or ethanol (methanol or formaldehyde can also be used). The unit should **not** be immersed in solvents.

#### **Guarantee**

This instrument is guaranteed against any defect in material or workmanship for a period as specified on the enclosed guarantee. This period is from the date of purchase, and within this period all defective parts will be replaced free of charge provided that the defect is not the result of misuse, accident or negligence. Servicing under this guarantee should be obtained from the supplier.

This User's Manual has been prepared solely for the convenience of customers and nothing in this User's Manual shall be taken as a warranty, condition or representation concerning the description, the merchantability, fitness for purpose or otherwise of the units or their components.


#### **Contact Information**


For technical, sales or servicing information, contact Waters Corporation.

#### **IMPORTANT NOTE:**

THIS EQUIPMENT SHOULD ONLY BE DISMANTLED BY PROPERLY TRAINED PERSONNEL. REMOVING THE TOP CASE EXPOSES POTENTIALLY LETHAL MAINS VOLTAGES.

THERE ARE NO USER MAINTENANCE PARTS WITHIN THE EQUIPMENT.





# **Fault Finding**

#### **General**

In the unlikely event that you experience any problems with your Heat Sealer which cannot be easily remedied, you should contact Waters Corporation and return the unit in its original packing.

#### **Accessories**

For Waters selection of collection plates and seals, visit us at www.waters.com

720001330EN

©2005 Waters Corporation.

Waters is a registered trademark of Waters Corporation.