# Afinion™ HbA1c

REF 1115017

Hemoglobin A1c test	GB	3...10
Hemoglobin A1c test	NO	11...18
Hemoglobin A1c test	SE	19...26
Hæmoglobin A1c test	DK	27...34
Hemoglobiini A1c testi	FI	35...42

**afinion**™  
POINT OF CARE CENTER

For use with the Afinion™ AS100 Analyzer:

## PRODUCT DESCRIPTION

### Intended use

Afinion™ HbA1c is an *in vitro* diagnostic medical device for quantitative determination of glycated hemoglobin (% hemoglobin A1c, HbA1c) in human whole blood. The measure of % HbA1c is recommended as a marker of long-term metabolic control in persons with diabetes mellitus<sup>1</sup>.

### Summary and explanation of the test

The human erythrocyte is freely permeable to glucose. Within each erythrocyte a slow, continuous, non-enzymatic process between hemoglobin A and various sugars takes place. The product formed is known as glycated hemoglobin, or glycohemoglobin<sup>2</sup>.

The chronic elevated blood sugar level of persons with diabetes mellitus will over time cause damage to the small vessels of the body. This damage develops slowly over years and is known to cause late complications<sup>3</sup>.

Good metabolic control, i.e. lowering the % HbA1c, has proven to delay the onset and slowing the progression of diabetes late complications<sup>3,4</sup>.

### Principle of the assay

Afinion™ HbA1c is a fully automated boronate affinity assay for the determination of the percentage of hemoglobin A1c in human whole blood.

The Afinion™ HbA1c Test Cartridge contains all of the reagents necessary for the determination of % HbA1c. The sample material is collected using the sampling device integrated in the Test Cartridge. The blood sample is then diluted and mixed with a liquid that releases hemoglobin from the erythrocytes. The hemoglobin precipitates. This sample mixture is transferred to a blue boronic acid conjugate, which binds to the cis-diols of glycated hemoglobin. This reaction mixture is soaked through a filter membrane and all precipitated hemoglobin, conjugate-bound and unbound (i.e. glycated and non-glycated hemoglobin) remains on the membrane. Any excess of conjugate is removed with a washing reagent.

The Afinion™ instrument evaluates the precipitate on the membrane. By measuring the reflectance, the blue (glycated hemoglobin) and the red (total hemoglobin) colour intensities are evaluated, the ratio between them being proportional to the percentage of HbA1c in the sample. The % HbA1c is displayed on the Afinion™ AS100 Analyzer.

### Standardization

Afinion™ HbA1c is traceable to the IFCC Reference Method for Measurement of HbA1c<sup>6</sup>. HbA1c values are reported according to national recommendations at DCCT level. Please contact your local supplier if the national recommendation is unknown.

The Afinion™ HbA1c meets the performance standards established by NGSP (National Glycohemoglobin Standardization Program, USA).

### Kit contents (per 15 tests unit)


- 15 Test Cartridges packaged separately in foil pouches with a desiccant bag
- 1 Package Insert

### Materials required but not provided with the kit

- Afinion™ AS100 Analyzer
- Afinion™ HbA1c Control
- Standard blood collection equipment

### Description of the Afinion™ HbA1c Test Cartridge

The main components of the Test Cartridge are the sampling device and the reagent container. The Test Cartridge has a handle, a barcode label with lot specific information and an ID area for sample ID.

Component	Function/composition
1 Sampling device a. Closed position b. Lifted position	For collection of patient sample or control.
2 Capillary	1.5 µL glass capillary to be filled with sample material.
3 Reaction wells a. Conjugate b. Membrane tube c. Washing solution  d. Reconstitution reagent e. Empty	Contains reagents necessary for one test: Patented blue boronic acid conjugate. Tube with a polyethersulfone membrane. Morpholine buffered sodium chloride with detergents and preservative. HEPES buffered sodium chloride with lysis and precipitation agents. N/A
4 Handle	For correct finger grip.
5 Barcode label	Contains assay specific information for the Analyzer.
6 Optical reading area	Area for transmission measurement.
7 ID area	Space for written or labelled sample identification.

### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Do not use Test Cartridges after the expiration date or if the Test Cartridges have not been stored in accordance with recommendations.
- Do not use the Test Cartridge if the foil pouch or the Test Cartridge itself has been damaged.
- Do not use the Test Cartridge if the desiccant bag is damaged and desiccant particles are found on the Test Cartridge. Do not wipe off.
- Each foil pouch contains a desiccant bag with 1 g silica gel. This material shall not be used in the assay. Discard the desiccant in a suitable container. Do not swallow.
- The Afinion™ HbA1c Test Cartridge contains sodium azide as a preservative. The concentration is <0.1%, which is below what is considered hazardous in normal use<sup>5</sup>. In case of leakage, avoid contact with eyes and skin. Wash with plenty of water.
- Do not re-use any part of the Test Cartridge.
- Handle and dispose the Test Cartridges and sampling equipment as potentially bio hazardous materials. Use gloves.

### STORAGE INSTRUCTIONS

#### Refrigerated storage (2-8°C, 36-46°F)

- The Afinion™ HbA1c Test Cartridges are stable until the expiry date only when stored refrigerated in sealed foil pouches. The expiry date is the last day of the month stated on the foil pouches and outer container.
- The Afinion™ HbA1c Test Cartridge must reach room temperature before use. Upon removal from refrigerated storage, leave the Test Cartridge in unopened foil pouches for at least 15 minutes.
- Do not freeze.

#### Room temperature storage (15-25°C, 59-77°F)

- The Afinion™ HbA1c Test Cartridges can be stored in unopened foil pouches at room temperature for 90 days. Note the date of removal from the refrigerator on the kit container.
- Avoid exposure to direct sunlight.

#### Opened foil pouch

- When the foil pouch is opened, the Test Cartridge must be used within 1 hour.
- Avoid exposure to direct sunlight.
- Avoid relative humidity above 90% (non-condensing).

### SAMPLE MATERIAL

The following sample materials can be used with the Afinion™ HbA1c test:

- Capillary blood sample (from finger stick)
- Venous whole blood without anticoagulants
- Venous whole blood with anticoagulants (EDTA, heparin, citrate or NaF)
- Afinion™ HbA1c Control

## Sample storage

- Capillary blood samples and venous blood samples without anticoagulants cannot be stored.
- Venous whole blood with anticoagulants (EDTA, heparin or citrate) may be stored refrigerated for 10 days. Do not freeze.
- Consult the Afinion™ HbA1c Control Package Insert for storage of control materials.


## SAMPLE COLLECTION

 Consult the Afinion™ AS100 Analyzer User Manual for general Test Cartridge handling instructions.

When collecting a sample for HbA1c analysis, take note of the following.


Blood sampling from finger (a):

- Make sure that the finger is clean, warm and dry.
- Use a suitable lancet to prick the finger.
- Hold the sampling device at an acute angle and bring the tip of the capillary just beneath the surface of the blood drop.



Sampling from tube (b):

- Patient samples stored refrigerated can be used without equilibration to room temperature.
- Do not use the sample if the blood is hemolysed.
- Mix the sample material thoroughly, but gently, by inverting the tube 8-10 times before collecting a sample.
- Tilt the tube, hold the sampling device at an acute angle and bring the tip of the capillary just beneath the surface of the sample material. The sample can also be extracted from the tube cap.

Sampling from control vial (c):

- Equilibrate the sample material to room temperature (15-30°C, 59-86°F) before use. Equilibration to room temperature takes approximately 30 minutes.
- Mix the control material thoroughly by shaking the vial for 30 sec. before extracting a sample.
- Tilt the vial, hold the sampling device at an acute angle and bring the tip of the capillary just beneath the surface of the control material. The sample can also be extracted from the vial cap.
- Wipe off any material on the outside of the vial neck and replace the cap.

### Important!

- Be sure that the capillary is completely filled, see arrow (d).
- Avoid air bubbles in the capillary.
- Do not wipe off the capillary.
- Do not use the Test Cartridge if dropped on the bench or floor after sample collection.

- Once the capillary is filled with sample, analysis of the Test Cartridge in the Afinion™ AS100 Analyzer must start within 1 minute.

## ANALYSING A SAMPLE

 Consult the Afinion™ AS100 Analyzer User Manual for detailed instructions on how to analyse a patient or control sample. An illustrated step-by-step procedure is also given in the Afinion™ HbA1c Quick Guide.

## INTERPRETATION OF RESULTS

Interpret the Afinion™ HbA1c test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the Afinion™ AS100 Analyzer System.

### Measuring range

The Afinion™ HbA1c measuring range is 4.0-18.0% HbA1c. The hemoglobin measuring range is 6-20 g/dL. The HbA1c results are displayed in 0.1% intervals.

 If the patient's HbA1c or hemoglobin value is outside the measuring range, no test result will be reported and an information code will be displayed. Consult the Afinion™ AS100 Analyzer User Manual.

### Analytical specificity

Afinion™ HbA1c measures the total glycated hemoglobin but reports the % HbA1c value.

The following hemoglobin (Hb) variants have been analysed and found not to affect the Afinion™ HbA1c test result: HbAS, HbAE, HbF and HbA<sub>J</sub><sup>10</sup>.

Carbamylated hemoglobin and pre-HbA1c do not affect the Afinion™ HbA1c test result<sup>10</sup>.

### Interference

No significant interference was observed up to the following concentrations:

- |  |            |
|--|------------|
| • Bilirubin  | 0.2 mg/mL  |
| • Triglycerides  | 5.6 mmol/L |
| • Cholesterol  | 5.8 mmol/L |
| • Glucose  | 5.0 mg/mL  |
| • Fructosamine   | 680 µmol/L |
| • Hemolysis  | <5%        |
| • Anticoagulants: EDTA, heparin, citrate and NaF in normal concentrations in blood collection tubes. |            |

### Limitations of the test

Diluted samples *cannot* be used in Afinion™ HbA1c. For accurate results outside the measuring range the sample should be analysed using another method.

### Reference range

For methods reporting DCCT traceable values, the upper limit of non-diabetic reference range is approximately 6% HbA1c<sup>7,8,9</sup>.

## QUALITY CONTROL

Quality control testing should be done to confirm that your Afinion™ AS100 Analyzer System is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results for patient samples be assured.

### Choosing control material

 Afinion™ HbA1c Control from Axis-Shield is recommended for routine quality control testing. Consult the Afinion™ HbA1c Control Package Insert.

If controls from another supplier are used, the precision must be determined and acceptable ranges for the Afinion™ AS100 Analyzer System must be established.

Lyophilised cell free lysate or hemolysed whole blood control materials cannot be used with Afinion™ HbA1c.

### Frequency of control testing

It is recommended analysing controls:

- With each shipment of Afinion™ HbA1c test kits.
- With each new lot number of Afinion™ HbA1c test kits.
- If the Afinion™ HbA1c Test Cartridge has not been stored or handled in accordance with specifications (see "Storage Instructions").
- If an unexpected test result is obtained.
- When training new personnel in correct use of Afinion™ HbA1c and the Afinion™ AS100 Analyzer.
- In compliance with national or local regulations for your facility.

### Verifying the control results

 The measured value should be within the acceptable limits stated for the control material. Consult the Afinion™ HbA1c Control package insert.

Re-test the control material if the control does not perform as expected. Re-test a sample from a new control vial if the measured value is still outside acceptable limits. Patients' test results must be declared invalid if the control does not perform as expected. If the problem persists contact your local Afinion™ supplier for advice before analysing patient samples.

## PERFORMANCE CHARACTERISTICS

### Method comparison

Fresh whole blood samples with % HbA1c within the range 4.5-11.4% were analysed using the affinity HPLC system Primus CLC385 (x) and Afinion™ HbA1c. The correlation data are summarised in table 1.

Table 1: Correlation to Primus CLC385 (x)

n	Bias	95% CI*
39	0.00% HbA1c	-0.31 - 0.30% HbA1c

\* CI = Confidence Interval

### Precision

The between-day and total precision values were determined according to NCCLS Protocol EP5-A. Two levels of HbA1c samples were assayed for 20 days. Two replicates were performed per day at each level.

Within-run precision was evaluated by running two levels of HbA1c samples, 20 determinations of each sample in sequence.

Precision data are summarised in table 2 and table 3.

Table 2: Between-day and total precision (% CV\*)

Sample	n	Average % HbA1c	Between day	Total
Level 1	17	5.6	0.2%	1.1%
Level 2	20	10.0	0.0%	1.1%


Table 3: Within-run precision (% CV\*)

Sample	n	Average % HbA1c	Within run
Level 1	17	5.6	0.9%
Level 2	20	10.0	0.7%

\* CV = Coefficient of Variation

## SYMBOLS

The following symbols are used in the packaging material for Afinion™ HbA1c.

	Conformity to the European directive 98/79/EC on <i>in vitro</i> diagnostic medical devices
	In Vitro Diagnostic Medical Device
	Catalogue number
	Lot number
	Test Cartridge
	Contents sufficient for "Σ" number of tests
	Consult the Afinion™ user instructions
	Caution, consult instructions for use
	Expiration date (year-month)
	Storage temperature (store at 2-8°C, 36-46°F)
	Manufacturer

## BIBLIOGRAPHY

- 1 Lenzi S *et al.*, The Clinical Usefulness of Glycated Hemoglobin in Monitoring Diabetes Mellitus: A long-Term Study. *Clin Chem* 1987; 33:55-56.
- 2 Frantzen F, Studies on glycated hemoglobin and new analytical methods for its determination. Dr. philos. thesis. NTNU, Trondheim 1998. ISBN 82-7861-116-5.
- 3 The Diabetes Control and Complications Trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. *N Engl J Med* 1993; 329:977-986.
- 4 Sacks DB *et al.*, Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. *Clin Chem* 2002; 48:436-472.
- 5 Council Directive 67/548/EEC, Council Directive 76/769/EEC, Council Directive 1999/45/EC and Commission Directive 2001/60/EC.
- 6 Jeppsson, JO *et al.*, Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, *Clin Chem Lab Med* 2002; 40(1):78-89.
- 7 Little, RR *et al.*, Relationship of glycosylated hemoglobin to oral glucose tolerance. Implications for diabetes screening. *Diabetes* 1988;37(1):60-64.
- 8 Little, RR *et al.*, Glycated haemoglobin predicts progression to diabetes mellitus in Pima Indians with impaired glucose tolerance. *Diabetologia* 1994; 37:252-256.
- 9 Rohlfing, CL *et al.*, Use of GHb (HbA1c) in Screening for Undiagnosed Diabetes in the U.S. Population. *Diabetes Care* 2000; 23(2):187-191.
- 10 European Reference Laboratory, Manufacturer Check-up Certification of Afinion™ HbA1c, 2004.

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## REFERANSER

- Lenzi S *et al.*, The Clinical Usefulness of Glycated Hemoglobin in Monitoring Diabetes Mellitus: A long-Term Study. *Clin Chem* 1987; 33:55-56.
- Frantzen F, Studies on glycated hemoglobin and new analytical methods for its determination. Dr. philos. thesis. NTNU, Trondheim 1998. ISBN 82-7861-116-5.
- The Diabetes Control and Complications Trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. *N Engl J Med* 1993; 329:977-986.
- Sacks DB *et al.*, Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. *Clin Chem* 2002; 48:436-472.
- Council Directive 67/548/EEC, Council Directive 76/769/EEC, Council Directive 1999/45/EC and Commission Directive 2001/60/EC.
- Jeppsson, JO *et al.*, Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, *Clin Chem Lab Med* 2002; 40(1):78-89.
- Little, RR *et al.*, Relationship of glycosylated hemoglobin to oral glucose tolerance. Implications for diabetes screening. *Diabetes* 1988;37(1):60-64.
- Little, RR *et al.*, Glycated haemoglobin predicts progression to diabetes mellitus in Pima Indians with impaired glucose tolerance. *Diabetologia* 1994; 37:252-256.
- Rohlfing, CL *et al.*, Use of GHB (HbA1c) in Screening for Undiagnosed Diabetes in the U.S. Population. *Diabetes Care* 2000; 23(2):187-191.
- European Reference Laboratory, Manufacturer Check-up Certification of Afinion™ HbA1c, 2004.

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## BIBLIOGRAFI

- 1 Lenzi S *et al.*, The Clinical Usefulness of Glycated Hemoglobin in Monitoring Diabetes Mellitus: A long-Term Study. *Clin Chem* 1987; 33:55-56.
- 2 Frantzen F, Studies on glycated hemoglobin and new analytical methods for its determination. Dr. philos. thesis. NTNU, Trondheim 1998. ISBN 82-7861-116-5.
- 3 The Diabetes Control and Complications Trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. *N Engl J Med* 1993; 329:977-986.
- 4 Sacks DB *et al.*, Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. *Clin Chem* 2002; 48:436-472.
- 5 Council Directive 67/548/EEC, Council Directive 76/769/EEC, Council Directive 1999/45/EC and Commission Directive 2001/60/EC.
- 6 Jeppsson, JO *et al.*, Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, *Clin Chem Lab Med* 2002; 40(1):78-89.
- 7 Little, RR *et al.*, Relationship of glycosylated hemoglobin to oral glucose tolerance. Implications for diabetes screening. *Diabetes* 1988;37(1):60-64.
- 8 Little, RR *et al.*, Glycated haemoglobin predicts progression to diabetes mellitus in Pima Indians with impaired glucose tolerance. *Diabetologia* 1994; 37:252-256.
- 9 Rohlfing, CL *et al.*, Use of GHB (HbA1c) in Screening for Undiagnosed Diabetes in the U.S. Population. *Diabetes Care* 2000; 23(2):187-191.
- 10 European Reference Laboratory, Manufacturer Check-up Certification of Afinion™ HbA1c, 2004.

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