GD2000 Slab Gel Dryer

Vacuum Gel Dryer System

Contents

Important Information	ii
Waste Electrical and Electronic Equipment (WEEE)	vii
1. Gel Dryer function and description	1
Design Features	2
2. Unpacking the GD2000	3
Setting Up the Vacuum System	4
Gel dryer controls	5
3. Specifications	7
4. Operating instructions.....	8
Step 1: Prepare the dryer.....	8
Step 2: Prepare the gel drying stack	9
Step 3: Set the temperature	12
Step 4: Set the timer.....	13
Step 5: Create a vacuum seal	15
Step 6: Disassembly	15
Options for manual operations	16
5. Care and maintenance	17
Cleaning	17
Replacing fuses.....	18
Replacing the rubber sealing sheet	19
6. Troubleshooting	20
7. Ordering information.....	21

Important Information – English

- If this equipment is used in a manner not specified by Hoefel, Inc. the protection provided by the equipment may be impaired.
- This instrument is designed for indoor laboratory use only.
- Only accessories and parts approved or supplied by Hoefel, Inc. may be used for operating, maintaining, and servicing this product.
- **Warning!** Because this instrument can develop sufficient voltage and current to produce a lethal shock, care must be exercised in its operation.
- This instrument is designed in accordance with the EN61010-1:2001 electrical safety standard. Nevertheless, it should be used only by properly trained operators. Read this entire manual before using the instrument and use only according to the instructions.
- The instrument must always be used with the earth lead of the power cord correctly grounded to earth at the mains outlet.
- Use only undamaged electrical wire and equipment specific for the voltages you will use. All equipment connected to high voltage should be in accordance with EN61010-1:2001.
- Keep the instrument as dry and clean as possible. Wipe regularly with a soft, damp cloth. Let the instrument dry completely before use.
- Do not operate the instrument in extreme humidity (above 80%). Avoid condensation by letting the unit equilibrate to ambient temperature when taking the instrument from a colder to a warmer environment.
- To permit sufficient cooling, ensure that the vents of the instrument are not covered.

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Waste Electrical and Electronic Equipment (WEEE)

English

This symbol indicates that the waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Please contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.

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1. Gel Dryer function and description

The Hofer® GD2000 Slab Gel dryer rapidly dries acrylamide and agarose gels and permanently bonds them to filter paper or transparent porous cellophane. This is accomplished by heating the gel slab while simultaneously drawing away released moisture with an external vacuum pump. The dryer accommodates one large (34 × 44 cm) gel, up to four standard (14 × 16 cm) gels, or twelve mini (8 × 10 cm) gels.


Fig 1. Features of the GD2000 Slab Gel Dryer.

- For a description of the dryer controls, see page 5.
- For diagrams of recommended drying stacks, see Fig 3 on page 9, Fig 4 on page 10, and Fig 5 on page 11.

Design Features

Drying surface

The PTFE coated, cast aluminum platen is resistant to acidic fumes that may be released while drying certain types of gels. A grid of vacuum conduits on the surface allows released moisture to be pulled away.

Vacuum port

Accepts vacuum tubing of 9-mm (3/8") ± 1 mm i.d. The port for the external vacuum source is at the back of the instrument, at the center of the platen.

Mains power

The mains power module houses the power cord receptacle and one or two input fuses.

115 V~ One F 12 A, 250 V, 3 AG fuse

230 V~ Two T 6.3 A, 250 V, 5 \times 20 mm

See "Care and maintenance" on page 17 for an illustration of the mains power module.

Vacuum power cord receptacle

Connects the vacuum pump to the vacuum timer.

115 V~, 50/60 Hz. Accommodates pumps that draw up to 5 A.

230 V~, 50/60 Hz. Accommodates pumps that draw up to 2.5 A.

2. Unpacking the GD2000

Carefully unwrap all packages and compare the items received with the packing list, making sure all items arrived. If any part is missing, contact your local Hoefer, Inc. sales office. Inspect all components for damage that may have occurred while the unit was in transit. If any part appears damaged, contact the carrier immediately. Be sure to keep all packing material for damage claims or to use in case you need to return the unit.

Setting Up the Vacuum System

The vacuum system must include a vacuum pump capable of moving an air volume of at least 1.5 m³/h, but not more than 6 m³/h. Rotary-vane-type pumps require both a chemical vapor trap and a cold trap (set to maximum cooling) to remove vapors that may damage the pump.

Note: We recommend a diaphragm vacuum pump such as the VP200 because the pump is chemically resistant to the liquids and vapors removed from the gel during drying. A water aspirator or house vacuum is insufficient for gel drying.

If you are using the VP200 Chemistry Diaphragm Pump, no traps are required because the pump is constructed of chemically resistant materials and is equipped with two vapor trap flasks. A cold trap is optional but can be installed to control the amount of vapor released to the atmosphere.

1

Attach the vacuum tubing (9 mm ± 1 mm i.d.) from the vacuum pump to the vacuum port on the GD2000.

2


Plug the vacuum pump either into the dryer control cabinet (the power cord receptacle is on the left side of the control panel) or into a lab power receptacle.

When the pump is plugged into the dryer control cabinet, the vacuum timer automatically turns the pump on and off. The 115 V~ model pump receptacle can accommodate any pump that draws less than 5 A. The maximum rating for the 230 V~ model is 2.5 A. If your pump exceeds this rating, plug it into a lab receptacle and manually turn the power on and off.

Gel dryer controls

The GD2000 control panel is located on the top of the gel dryer. (See Fig 1). A detailed view of the features of the control panel is shown in Fig 2.

Fig 2. Gel dryer controls.

number	control feature	function
1	LED display	<p>Shows Set Temperature, Set Time or Elapsed Time. On the left of the LED display, a light indicates whether the displayed value is °C (temperature) or Hrs (time). When the displayed value is time and the LED colon blinks, the value is Elapsed Time. When the colon is not blinking, the value is Set Time.</p> <p>Press either Adjust key once to go from Elapsed Time to Set Time.</p> <p>On the right side of the display, two red lights, labeled "Heater On" and "Vacuum On," indicate the status of the heater and vacuum.</p>
2	Time/Temperature key	Toggles the LED display between showing time or temperature.
3	Vacuum key	Toggles the vacuum outlet On or Off. A red light appears on the right side of the LED display when the vacuum is On.
4	Adjust keys	<p>Adjust Set Temperature and Set Time; reset Elapsed Time.</p> <p>Press a key once briefly to move one interval. Press and hold a key down to move in larger intervals. Press both keys simultaneously to reset Temperature or Time to the Reset values.</p>
5	Reset values	Lists the Reset values for Elapsed Time, Set Time and Set Temperature: Elapsed Time 0:00; Set Time OFF; Temperature 40 °C.

This declaration of conformity is only valid for the instrument when it is:

- used in laboratory locations,
- used as delivered from Hoefer, Inc., except for alterations described in the user manual, and
- connected to other CE-labeled instruments or products recommended or approved by Hoefer, Inc.

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Maximum Operating Ratings

Wattage	Heater: 800 W, thermostatically controlled
Pump Outlet	575 W
Power requirements	Model: 80-6428-84: 115 V~, 50/60 Hz 80-6429-03: 230 V~, 50/60 Hz

Environment

Operating environment:	Indoor use, 15–40 °C ambient Relative humidity ≤ 80% for 15–31 °C, decreasing linearly to 50% for 31–40 °C Altitude ≤ 2000 m Installation category II Pollution Degree 2
------------------------	--

Dimensions (w × d × h)	55.0 × 43.5 × 8.5 cm (21.6 × 17 × 3.4 in.)
----------------------------------	---

Weight	8 kg (17.6 lbs)
---------------	-----------------

Product certifications	CE, UL61010A-1, CSA
-------------------------------	---------------------

Important! The plate surface of the GD2000 Gel Dryer reaches high temperatures during operation. Do not touch the plate surface during operation!

4. Operating instructions

After you have attached the vacuum pump to the GD2000, follow these instructions to prepare the dryer and gel drying stack. Once you set the temperature and the timer, the GD2000 automatically starts the vacuum pump and turns on the heat after 10 seconds. At the end of a timed run, the heat is turned off first and the vacuum ten minutes later.

Step 1: Prepare the dryer

Wipe away all contaminants with a soft damp cloth. See “Care and maintenance” on page 17 for recommendations on how to remove accumulations of radioactive materials. Fit the stainless steel screen into the recess on the platen and then place a sheet of filter paper on the screen slightly larger than the surface area required by the gel(s). The paper should not extend over the ridge that surrounds the platen.

Step 2: Prepare the gel drying stack

The slab dryer accommodates both agarose and polyacrylamide gels. The configuration of the gel stack layers depends on the next processing step and the thickness of the gel(s).

Drying between cellophane sheets

1

Immerse two sheets of porous cellophane in water.

2

Lay one cellophane sheet smoothly on top of the filter paper liner. Carefully center the gel on the cellophane. Cover the gel with the second sheet of wet cellophane.

3

The stiff cover sheet produces a smooth gel surface that reduces scan irregularities.

4

Check that the edges of all the sheets fall within the recess of the platen. If necessary, trim the corners of the sheets to fit within the recess.

5

Cover the stack with the silicone rubber sealing sheet.


Note: Cellophane provides a transparent support for densitometric scanning. Remove all air pockets at every step of this procedure to avoid scanning distortions.

Drying between cellophane sheets

A filter paper liner and a stainless steel screen lie under each stack. A translucent rubber sealing sheet lies on top.

The gel type and thickness determine the cover sheet and layers used around the gel.

Fig 3. Gel drying stack for drying to cellophane. Use this configuration for gels to be scanned, imaged, and stored.

Drying thin or low-concentration polyacrylamide gels (≤ 1.5 mm) and agarose gels on paper

1

Lay a sheet of filter paper on top of the filter paper liner and position the gel on this sheet, taking care to avoid trapping air beneath the gel.

2

Cover the gel with thin plastic wrap. Do not leave wrinkles in the plastic wrap.

3

Use the stiff plastic cover sheet with polyacrylamide gels but not with agarose gels.

4


Check that the edges of all the sheets fall within the recess of the platen.

If necessary, trim the sheets to fit within the recess.

5

Cover the stack with the silicone rubber sealing sheet.

Fig 4. Stacks for drying thin or low-concentration polyacrylamide gels and agarose gels to filter paper for autoradiography.

*not used with agarose gels

Drying thick gels (> 1.5 mm), high concentration gels, gradient gels on paper

1

Lay a sheet of filter paper on top of the filter paper liner and then position the gel on this sheet, taking care to avoid trapping air beneath the gel.

2

Cover the gel with the porous polyethylene cover sheet, with the smooth side toward the gel.

3


Check that the edges of all the sheets fall within the recess of the platen.

If necessary, trim the sheets to fit within the recess.

4

Cover the stack with the clear silicone rubber sealing sheet.

Fig 5. Stacks for drying thick gels to filter paper for autoradiography.

Note: The melting temperature of an agarose gel is dependent on its concentration and properties. The drying temperature should not exceed the melting temperature. We recommend a drying temperature of 50 °C for most agarose gels.

Step 3: Set the temperature

The highest setting, 80 °C, can be used for many types of gels for reliable, rapid drying. Use a temperature of 50 °C for agarose gels. Polyacrylamide gels prepared for fluorography may require a setting of 60 °C to protect the floors. Follow the manufacturer's instructions. If your gels tend to crack, slower drying at a lower temperature setting may be indicated.

1

To set the temperature, make sure the red light labeled "°C" is lit. If necessary, press the Time/Temp key to toggle between the time and temperature display.

2

Use the Adjust keys to change the temperature. You can set the heater to any temperature from 40 to 80 °C or to OFF (room temperature).

Press the Up or Down key once briefly to change by one degree. Press and hold the key down to count five 1-degree intervals, followed by 5-degree intervals.

When the display reads 40 °C, press the Down key once to go to OFF. When the temperature is 80 °C, press the Up key once to go to OFF.

Press the Up and Down arrows simultaneously to reset the temperature to 40 °C. To use the vacuum at room temperature, turn the heater off.

Note: When the heater attains the temperature setting, the red heater light goes off. The heater light blinks when the heater comes on to maintain the set temperature.

The platen begins to heat at the same time that the vacuum pump starts— 10 seconds after you have finished setting the time. A red light next to the words "Heater On" indicates when the heating element is on.

Step 4: Set the timer

The amount of time required for a gel to dry depends on such factors as gel thickness, gel concentration, drying temperature, and vacuum applied. A typical sequencing or 1.5 mm 10% T gel can be expected to dry in approximately 45 minutes at 80 °C. Larger gels may take 2 to 3 hours. When dry, the thickness of an agarose gel seen through the silicone flap decreases to about 1 mm.

1

Press the Time/Temp key to go to the Time mode.

The LED display reads 0:FF and the LED colon blinks to indicate Elapsed Time.

2

Press either the Up or Down key once to go from Elapsed Time to Set Time.

3

Press either the Up or Down key to change the Set Time.

Press Up once to set time for a continuous run. On a continuous run, the LED display reads r:un and the heater and vacuum run continuously until you manually turn them both off.

Press Up again to count up in 15-minute intervals.

Press the Up and Down keys simultaneously to reset the time to 0. When the Set Time is set to 0, the display reads 0:FF.

When you have finished setting the time, you can start the heater and vacuum.

Note: Agarose gels become brittle when over-dried.

To start the heater and the vacuum pump

- You can press the Time/Temp key to immediately start the heater and vacuum pump.
- If you don't press the Time/Temp key, ten seconds after you stop adjusting the time, the heater and vacuum pump start automatically.

As the heater and vacuum pump start, the timer begins to count Elapsed Time (Hours:Minutes). The LED colon between the hours and minutes blinks each second when the timer is in Elapsed Time.

You may change the Set Time at any point after the LED begins to count Elapsed Time.

To change the set time while the heater is on

1

Make sure the LED displays time.

The red light labeled Hrs is lit when the LED displays time. Press the Time/Temp key to toggle between temperature and time display.

The LED colon blinks when the display shows Elapsed Time.

2

Press the Up or Down arrow key to change the display from Elapsed Time to Set Time.

3

Press the Up or Down arrow key to change the Set Time.

Note: If using a cold trap with inline valves, close the valve between the trap and gel dryer and then open the valve between the trap and the pump. After the trap is pumped down, open the valve to the gel dryer. The extra vacuum should quickly pull the sealing sheet down and accelerate the sealing process.

Important! Once the gel has begun to dry, do not break the vacuum seal until the gel is completely dry. Gels may crack if you turn off the vacuum before the gel is dry.

Note: Some gels may curl as they dry. To minimize curling, the vacuum continues for 10 minutes after the heating timer shuts off.

Note: If you manually turn off the vacuum during a continuous run, the heat remains on until you also turn it off manually.

Note: If the gel contains radioactive materials and was covered with plastic wrap, dispose of the wrap according to local regulations pertaining to radioactive waste.

Note: Fluorescent compounds, such as ethidium bromide, cannot be visualized after drying.

Step 5: Create a vacuum seal

Ten seconds after you set the timer, the vacuum starts automatically if it is connected through the vacuum receptacle on the dryer. Watch for a seal to form between the rubber overlay and the platen.

If the seal does not form almost immediately, check stack for misalignment. No materials should extend beyond the edge of the recess. Assist seal formation by pressing gently at each corner to ensure that the rubber sheet is pulled into the recess.

Without removing the rubber sealing sheet, periodically inspect the gel as it dries. When the gel appears dry, check the temperature of the gel by briefly touching the sealing sheet over the gel. Wet gels feel cold compared to the platen. Typically, the gel is dry when the drying surface has become evenly hot. A markedly flattened gel contour also indicates that the gel has dried.

In automatic mode, when the set time is reached, the dryer beeps once and the heater turns off. The vacuum power remains on for 10 minutes. During this cooling down period, the Elapsed Time display counts from “C:00” to “C:10”. After ten minutes, the vacuum power automatically turns off and the dryer beeps once.

Step 6: Disassembly

Remove each layer of the drying stack and clean the dryer according to instructions in “Care and maintenance” on page 17.

Note: When you manually preheat the dryer, be sure the set time includes the time needed to prepare the gel stack as well as the time needed to dry the gel.

Options for manual operations

To apply a vacuum without heat. Set the Temperature to 0:FF, then set the timer. Ten seconds after you set the time, the vacuum starts without heat. The vacuum stops when the Elapsed Time equals the Set Time.

To preheat the dryer. Set the Temperature and set the Time. After ten seconds, the vacuum starts. Press the Vacuum key to turn off the vacuum and leave the heat on. To restart the vacuum, press the Vacuum key again.

5. Care and maintenance

Cleaning

1

Turn mains power switch off and unplug the power cord.

2

Remove the stainless steel screen and cover sheets and wash separately with a mild laboratory detergent. Do not use abrasives or solvents on any part of the dryer.

Periodically remove accumulations left by autoradiography reagents from the platen and stainless steel screen. Apply a strong detergent, such as Contrad™ 70 or Decon™ 90, for no longer than 5 minutes and rinse thoroughly.

3

Dry with a soft towel.

Important! Fuses protect equipment by disconnecting loads too large for the circuit design. Always replace fuses with those that conform to the specified fuse rating.

Important! Detach the power cord before replacing fuses.

Replacing fuses

115 V~ Model. The fuse drawer holds one F 12A 250V 3AG fuse and one shorting coil.

230 V~ Model. The fuse drawer holds two F 6.3A 250V 5 × 20 mm fuses.

The fuse drawer is in the power entry module, found on the left side of the control panel (See Fig 6).

1

Insert a small, flat-blade screwdriver into the slot below the fuse drawer. (See Fig 6.) Push in the direction of the arrow to release the drawer. Grasp the fuse drawer with your fingers and pull it out.

2

Pull the fuse out of the drawer to inspect it. If the fuse element is burned or broken, replace it.

If the fuse appears to be intact, check it with a multimeter. A reading of 1Ω or less indicates the fuse is still usable.

3


Push the fuse drawer back into the power entry module until it snaps into place.

4


Plug the power cord in and turn the mains power switch on.

Fig 6. Mains power module.

Insert a small screwdriver blade in the slot. Push in the direction of the arrow to release the fuse drawer.

Use your fingers to pull out the fuse drawer

Replacing the rubber sealing sheet

Occasionally, a tear or nick in the rubber sealing sheet (SE1143) may inhibit the formation of the vacuum seal. The rubber sealing sheet is attached to the dryer by four screws in a retaining bar along the edge of the platen. To change the sealing sheet, you also need a small tube of clear silicone adhesive, available at hardware stores.

1

Use a Phillips-head screwdriver to remove the four screws on top of the retaining bar that holds the rubber sealing sheet in place.

2

Lift off the retaining bar and the damaged rubber sealing sheet.

If necessary, use a sharp edge to clean out any old silicone adhesive on the inside of the retaining bar.

3

Position the replacement sealing sheet on the platen, aligning the four holes in the overlay with the four holes for the screws.

4

Lay a bead of clear silicone adhesive along the inside edge of the retaining bar. Replace the retaining bar, aligning the four holes of the retaining bar over the holes in the platen and sealing sheet.

5

Screw the retaining bar and sealing sheet in place.

6. Troubleshooting

problem	solution
No power or LED Display	Check that the power switch is turned on. Check that the dryer is plugged into a working receptacle. Check the fuse(s).
No heat	Make sure time is set and the LED colon is blinking. Make sure the temperature is not set to OFF. If still no heat, contact your local Hoefer, Inc. distributor for service.
No vacuum	Make sure the tubing connects the vacuum port to the vacuum pump. Make sure the vacuum pump is plugged into the vacuum receptacle on the dryer. Check that time is set. Vacuum automatically starts 10 seconds after time is set. Check switches on the vacuum source.
Can't create vacuum seal	Make sure the rubber sealing sheet seals around the entire inside edge of the recess. Omit the mylar sheet and use plastic wrap on top of gel. Check for tears or punctures in the rubber sealing sheet. Replace it, if necessary.
Gels crack	Use thinner gels (≤ 0.75 mm), if possible. Thin gels rarely crack. Reduce % T. Equilibrate gels with 30% ethanol, 2% glycerol for one hour before drying. Make sure gel is completely dry before turning off vacuum.
Gels do not dry	Do not use >5% glycerol during pre-drying treatment. Empty the liquid or cold trap. Replenish dry ice in the cold trap. Make sure to place only porous cellophane or filter paper under the gel. Do not use plastic film or the stiff plastic cover sheet in the stack below the gel.
Fluors become degraded	Follow manufacturer's handling instructions, paying close attention to recommended temperature exposure.

7. Ordering information

product	quantity	code number
GD2000 Vacuum Gel Dryer System Includes: stainless steel screen, VP200 Vacuum pump, Vacuum tubing, 10 sheets of filter paper, 50 sheets of porous cellophane, one mylar sheet and one porous polyethylene sheet		
115 V~	1	GD2001
230 V~	1	GD2002

Replacement Parts

Filter paper, 35 × 44 cm	25	SE1141
Porous cellophane, 35 × 44 cm	50	SE1142
Clear silicone rubber sealing sheet	1	SE1143
Stiff plastic cover sheet	1	SE1144
Porous polyethylene cover sheet	1	SE1145
Stainless steel screen	1	SE1146
Vacuum tubing, 8 mm i.d., 3 m	1	VT3

115 V~ model

Detachable power cord, 115 V~, 15 A	1	PSCORD15A-115V
Fuses, F 12 A, 250 V, 3AG	5	PSF12A-FB-3AG

230 V~ model

Detachable power cord, 230 V~	1	PSCORD-2230V
Fuses, T 6.3 A, 250 V, 5 × 20 mm	5	PSF6.3A-SB-SX20

Hoefler, Inc.

84 October Hill Road
Holliston, MA 01746

Toll Free: 1-800-227-4750

Phone: 1-508-893-8999

Fax: 1-508-893-0176

E-mail: support@hoeflerinc.com

Web: www.hoeflerinc.com

Hoefler is a registered trademark
of Hoefler, Inc.

Coomassie is a trademark of
ICI plc.

Conrad 70 and Decon 90 are
trademarks of Decon Lab.

© 2012 Hoefler, Inc. —
All rights reserved.

Printed in the USA.