

QuantorDent Vet User Manual

Image Acquisition and Diagnosis Software

Doc No.: TM-806-EN

Rev 0.2.2 February 2014

Part No.: CR-FPM-11-006-EN

3DISC, FireCR, Quantor and the **3D Cube** are trademarks of **3D Imaging & Simulations Corp.**, South Korea, and its affiliates. All other trademarks are held by their respective owners and are used in an editorial fashion with no intention of infringement. The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications, which must be met by **3D Imaging & Simulations Corp.** All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available in your local area. Please contact your local sales representative for availability information. **3D Imaging & Simulations Corp.** strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© Copyright 2014 **3D Imaging & Simulations Corp.**, all rights reserved, printed, and published in South Korea by **3D Imaging & Simulations Corp.**

3DISC
I M A G I N G

Contact



815, Tamnip-Dong, Yuseong-Gu,
Daejeon, Korea

Tel : 82-42-931-2100
Fax : 82-42-931-2299

Website : www.3DISCimaging.com
E-mail : info@3DISCimaging.com

3DISC Americas
22560 Glenn Dr, Suite 116
Sterling, VA 20164 USA
Tel : 1-703-430-6080
E-mail : info@3DISCimaging.com



3DISC Europe
Gydevang, 39-41, 3450 Alleroed, Denmark
Tel : 45-88-276-650
E-mail : info@3DISCimaging.com

Warnings and Symbols used

To ensure the safety of patients, staff and other persons, any changes to software and hardware delivered by **3D Imaging & Simulations Corp** may only be made with prior written permission from **3D Imaging & Simulations Corp**.

Please read the respective manuals of the connected devices, such as of the X-ray generator, sensor, or reader, before starting to use the **QuantorDentVet** software.

Following symbols will be used throughout this manual.



DANGER

The functionality of the software can be destroyed in the case of incorrect use.

If unauthorized changes have been made to delivered software and hardware components, the warranty by **3D Imaging & Simulations Corp.** becomes void. **3D Imaging & Simulations Corp.** will not accept any responsibility or liability for the proper functioning of the product in such a case.



CAUTION

The functionality of the software can be limited in the case of incorrect use. Hints that require special attention.



NOTE

Notes represent information that is important to know but which do not affect the functionality of the software.

Medical Device Security

Users must take steps to secure their networks and protect their Medical Information Systems which includes a risk assessment strategy, network defense in depth strategy, business continuity planning, etc.

- ✓ **User Authentication**
Only authorized users should log on to computers on which medical information systems are installed.
- ✓ **Password Security**
In today's world, passwords can be compromised in literally seconds by using a wide variety of tools and techniques. To lower the possibility of a compromised password, it is vital to adhere to a set of protocols.
 - Choose a password between 7 ~ 10 characters using both alpha and numeric characters.
 - Do not share the password.
 - Do not base the password on a pet's name, a relative's name or any dictionary word.
 - Do not write down the password.
 - Do not leave the account logged on.
- ✓ **User Access Control**
Configure the workstation to prompt for logon after coming out of stand-by mode.
- ✓ **Internet Usage**
Accessing to the Internet exposes the computer to a plethora of vulnerabilities such as:
 - Viruses
 - Spyware
 - Trojans
 - Hostile CodesIt is not recommended to install any unauthorized software on the computer. Peer-to-peer software can expose your entire hard drive to any individual running the same type of software.
- ✓ **Antivirus Products**
Use of antivirus software can increase CPU and memory usage, which can cause a slight degradation in the performance of the system. However, functionality should not be affected.
- ✓ **Physical Security**
It is recommended that the user employs some method of physical security when dealing with the system to ensure that only authorized personnel have access to the product.

There are several vulnerabilities a malicious user could exploit locally. Some examples are:

- Theft of equipment
- Local password cracking
- Installation of hardware key loggers
-

Table of Contents

1. <i>Introduction</i>	10
1.1. <i>Major Features</i>	11
2. <i>Installation</i>	12
2.1. <i>Installation Environment</i>	12
2.1.1. Recommended Computer Environment	12
2.1.2. Minimum Computer Requirements	12
2.2. <i>Installation</i>	13
2.2.1. Software Installation	13
2.3. <i>Start and Exit</i>	17
2.3.1. Start system	17
2.3.2. Start program	17
2.3.3. Exit program	18
2.4. <i>FireCR Dental Calibration</i>	19
2.4.1. Step 1: Erase	20
2.4.2. Step 2: Scan Blank	20
2.4.3. Step 3: Scan Low Dose	20
2.4.4. Step 4: Scan Mid Dose	20
2.4.5. Step 2: Scan High Dose	20
2.4.6. Calibration	20
2.4.7. Upload	20
2.4.8. Download	20
3. <i>Software Overview</i>	21
3.1. <i>Supported Resolutions</i>	21
3.2. <i>Home</i>	21
3.3. <i>Navigation Bar</i>	22
3.4. <i>Study list</i>	22
3.5. <i>Create New Study</i>	23
3.6. <i>Patient Information</i>	23
3.7. <i>Scan</i>	24
3.8. <i>Review</i>	24

4. Home Screen	25
4.1. Screen Layout	25
4.2. System menu.....	26
5. Create New Study	27
5.1. Screen Layout	27
5.2. Edit Column Dialog Box.....	29
6. Scan.....	30
6.1. Screen Layout	30
6.2. Tool Box.....	31
6.2.1. Image.....	31
6.2.2. Image+ Dentition	32
6.2.3. Image + Chart.....	32
6.2.4. Chart.....	33
6.2.5. Template	33
6.3. Image Scan	34
6.3.1. For Standalone version	34
6.3.2. For Client version	35
7. Image.....	36
7.1. Screen Layout	36
7.2. ROI.....	37
7.3. Marking.....	38
7.4. Rotation/Flip	39
7.5. 180-degree rotation	39
7.6. Auto Window	39
7.7. Reject	40
7.8. QA – Image Adjustment.....	41
7.9. QA – Image Processing.....	42
7.10. Image Transport	43
7.11. Image Comment.....	43
8. Task List.....	44
8.1. Screen Layout	44
8.2. Print.....	45

9. Review	46
9.1. Screen Layout	46
9.2. View Caption Bar	48
9.3. Full Screen	49
9.4. Image Processing	50
9.5. Measure	51
9.6. Marking	54
9.7. File	54
9.8. Tool	54
9.9. Apply	55
10. Study list	56
10.1. Screen Layout	56
10.2. List	57
10.3. Edit Column	58
10.4. Search	59
10.5. Export	59
11. Confirm Transport Results	61
11.1. Screen Layout	61
11.2. List	62
12. Settings	63
12.1. Screen Layout	63
12.2. General	64
12.2.1. Options	64
12.2.2. User Information	64
12.2.3. Password for Processing Parameter	64
12.2.4. Display Language	64
12.3. Fonts	65
12.3.1. General	65
12.3.2. Marking	65
12.3.3. DICOM Printer	65
12.4. Network / Export	66
12.4.1. Image Server	66

12.4.2. Options	67
12.5. Overlay	67
12.5.1. Marking	67
12.5.2. Auto Marking Margin	67
12.5.3. Overlay Contents	68
12.6. System	68
12.6.1. File Management	68
12.6.2. System Options	68
12.6.3. Date	69
12.6.4. Dental	69

1. Introduction

QuantorDentVet is an IP reader for intraoral X-ray imaging. **QuantorDentVet** is used to scan and review intraoral X-ray images. When scanning an IP radiated with X-ray, **FireCR Dental** shows the preview image on the LCD screen and the image is saved in the built-in memory storage. Shortly after the end of scanning, the IP is deleted and it can be used for the next imaging.

Images stored in the scanner are transferred to **QuantorDentVet** in the image workstation via network or USB. The images are post-processed through **QuantorDentVet** and saved in the database that allows searching in the future.

QuantorDentVet can transport the acquired images to the hospital image storage system using DICOM Send.

The scanner can be connected to the in-hospital network or directly to the image workstation.

The system has following components:

- FireCR Dental
- Fire ID
- QuantorDentVet

Safety Instruction

To ensure the safety of patients, staff and other people, any changes to software and hardware delivered by **3D Imaging & Simulations Corp** may only be made with prior written permission from **3D Imaging & Simulations Corp**.

Responsibilities

If unauthorized changes have been made to software and hardware components delivered by **3D Imaging & Simulations Corp**, the warranty by **3D Imaging & Simulations Corp** becomes void. **3D Imaging & Simulations Corp** will not accept any responsibility or liability for the improper functioning of the product in such a case.

1.1. *Major Features*

Image Acquisition and Management

QuantorDent supports various image acquisition methods that meet the needs of dental hospital. It can be used in a simple structure where one workstation is connected to a scanner for basic use as well as in a multi-structure where multiple image workstations are connected to multiple scanners. When multiple systems are used, it supports **Fire ID**, which enables automatic image classification to designated patient and test information regardless of the scanning location.

Image Processing

To create an optimal image, you can apply different functions, including cropping, enhancing edges, and increasing contrast. You can also add or delete images and modify the image sequence at your discretion.

Built-in Viewer

QuantorDentVet not just simple image acquisition software, but can also be used as a diagnostic viewer. Like other viewers, it can manage patient and test information and contains all basic viewer functions.

2. *Installation*

2.1. *Installation Environment*

2.1.1. *Recommended Computer Environment*

Operating system	Microsoft Window XP, 7 and 8
CPU	Intel CORE i5
Memory	4G or higher
Hard disk	500GB or higher
Network	1Gbps Ethernet
Video output	32-bit color display
Screen resolution	1920 x 1080, 1440 x 900

2.1.2. *Minimum Computer Requirements*

Operating system	Microsoft Window XP
CPU	Intel Processor
Memory	2GB
Hard disk	80GB
Network	1Gbps Ethernet
Video output	32 bit-color display
Screen resolution	1366 x 768, 1280 x 800, 1600 x 900, 1280 x 1024

2.2. Installation

2.2.1. Software Installation

- Log on to a PC with an administrator account.
- Connect the USB dongle to a USB port.
- Insert the installation CD. The installation program runs automatically.
- If it does not automatically run, manually find a file named “QuantorDentVet 2.x.x.Setup.exe” in the CD and run the file.
- Click the “Next” button as shown in the Figure.




Figure 2.1 Welcome Dialog Box

- Select components and click “Next”.




Figure 2.2 Components Dialog Box

- As shown in the Figure, select language and click the “Next” button.




Figure 2.3 Select Language Dialog Box

- As shown in the Figure, select the application type. If you do not have QuantorDentVet Server, select Standard. Otherwise, select Client of QuantorDentVet Server.




Figure 2.4 Select Application Dialog Box

- As shown in the Figure, select the installation path and click the “Next” button.




Figure 2.5 Select Installation Path

- The installation status is shown like the Figure.




Figure 2.6 Display Installation Status

- The dongle driver installation status is shown like the Figure.




Figure 2.7 Installation of dongle Driver

- The **FireCR** USB driver
The installation status is shown
like the Figure.
like the Figure.




Figure 2.8 Driver Installation Status

- When the installation is
successfully completed,
shown like the Figure.
Click "Finish" to finish the
installation.




Figure 2.9 Finish Installation

2.3. Start and Exit

2.3.1. Start system

1. Turn on the Fire CR device.
2. Turn on the PC and monitor.



CAUTION

Do not run any application software other than **QuantorDentVet** during operation of the scanner. This may slow the scanner response.

2.3.2. Start program

After successfully booting Windows, double-click the **QuantorDent** icon on the desktop screen as shown in Figure 2.10 to start the **QuantorDent** program.




Figure 2.10 Home Screen



QuantorDentVet

Figure 2.11 Desktop Icon

2.3.3. *Exit program*

Select the system button that looks like Figure 2.12 on the Home screen.




Figure 2.12 System button

You can select "Exit" in the system menu to close the program.





Figure 2.13 System menu

2.4. FireCR Dental Calibration

In the system menu, select IP Calibration to open a window as shown in Figure 2.14. Process calibration as instructed in the window.



CAUTION

Calibrate the scanner before use. The “Please calibrate system” message pops up if the scanner is not calibrated or calibration data does not exist.

It is recommended to calibrate the reader after moving it and as a part of regular maintenance..


Device Calibration

The X-ray beam should cover the entire IP area.

needs SIZE4C IP

ScanBlank	ScanLowDose	ScanMidDose	ScanHighDose
224	3529	9055	34136
(50~500) 0.00kVp 0.00mAs SID 280 mm	(1500~5000) 70.00kVp 0.49mAs SID 280 mm	(8000~16000) 70.00kVp 1.40mAs SID 280 mm	(30000~45000) 70.00kVp 3.50mAs SID 280 mm

Figure 2.14 Calibration



CAUTION

SIZE 4C IP must be used for calibration.

2.4.1. *Step 1: Erase*

At first, you can see the Erase button flickering. Place the IP on the tray and insert it. When the erase is completed, the Scan Blank part starts flickering.

2.4.2. *Step 2: Scan Blank*

Place an IP that has not been radiated with X-Rays on the tray and insert it. At the scan completion, the Scan Low Dose part starts flickering.

2.4.3. *Step 3: Scan Low Dose*

Adjust the X-Ray imaging settings to 70kVp, 0.49mAs and scan the IP from the recommended X-Ray imaging distance.

Place the IP on the tray and insert it. Once the scan is completed, the measurement is displayed below the button. It should be between 1500 and 5000. If so, click Accept to process the next step. If not, click Reject to repeat Step 3.

2.4.4. *Step 4: Scan Mid Dose*

Adjust the X-Ray imaging settings to 70kVp, 1.4mAs and scan the IP from the recommended X-Ray imaging distance. Place the IP on the tray and insert it. Once the scan is completed, the measurement is displayed below the button. It should be between 8000 and 16000. If so, click Accept to process the next step. If not, click Reject to repeat Step 3.

2.4.5. *Step 2: Scan High Dose*

Adjust the X-Ray imaging settings to 70kVp, 3.5mAs and scan the IP from the recommended X-Ray imaging distance. Place the IP on the tray and insert it. Once the scan is completed, the measurement is displayed below the button. It should be between 30000 and 45000. If so, click Accept to process the next step. If not, click Reject to repeat Step 3.

2.4.6. *Calibration*

This creates a Calibration data file.

2.4.7. *Upload*

If you click the Upload button, the calibration data is transferred to the scanner for storage.

2.4.8. *Download*

If you click the Download button, stored calibration data in the scanner is downloaded to the PC.

3. Software Overview

3.1. Supported Resolutions

Landscape resolutions with 16:9, 4:3 and 5:4 ratios are supported. Portrait ratios are not supported. The horizontal resolution of the monitor must be at least 1280 pixels, and the vertical resolution of the monitor must be at least 768 pixels. Table 3.1 provides the recommended resolution.

Table 3.1. Supported Resolutions

Screen ratio	Screen resolution
16:9	1366 x 768 1280 x 800 1400 x 900 1600 x 900 1920 x 1080
4:3, 5:4	1280 x 1024

3.2. Home

The Home screen is the first screen that appears when the program is run. From the Home screen, you can create a new study or view existing studies.




Figure 3.1 Home Screen

3.3. *Navigation Bar*

You can see the Navigation Bar in the bottom of screen. The Navigation Bar is always displayed while the program is running, allowing the user to easily switch screens during a task.




Figure 3.2 Navigation Bar

- ① Home
- ② Studies
- ③ Scan
- ④ Review

3.4. *Study list*

This displays a collection of windows that are related to the study. You can go to this screen by clicking on the Study button in the Navigation Bar. On the left side, the study tool bar is located. In the study tool bar, you can select study-related windows for the study list, work list and transport.




Figure 3.3 Studies Screen

- ① You can perform various tasks by searching existing studies and selecting a study.
- ② You can confirm the transmission results of the acquired images.

3.5. Create New Study

The New Study screen is the first page to start a scan. After entering the necessary patient information and study information, and adding a new study, you can begin a scan.




Figure 3.4 Create New Study

3.6. Patient Information

On the top of all screens except the Home and New Study screens, there is a patient Information Bar. The Patient Information Bar displays patient information and other related information.



Figure 3.5 Patient Information

- ① Patient Name
- ② Patient ID
- ③ Patient's Date of Birth
- ④ Patient's Age
- ⑤ Patient Sex (M: Male, F: Female, O: Other, SM: Sterile Male, SF: Sterile Female, S: Sterile)
- ⑥ Current Time
- ⑦ Model Name
- ⑧ System

3.7. Scan

In the Scan screen, you can insert an IP into the **FireCR Dental** and process image acquisition.




Figure 3.6 Scan Screen

3.8. Review

You can search existing images in the Review screen. It is useful for diagnosis after acquiring the images.




Figure 3.7 Review Screen

4. Home Screen

4.1. Screen Layout

The Home screen is the first screen that appears when the program is run. From the Home screen, you can create a new study or view existing studies. You can go to the Home screen by:

- Clicking the Home button in the Navigation Bar;
- Clicking the Study Exit button in the Scan screen; or
- Clicking the Study Exit button in the Review screen.




Figure 4.1 Home Screen

- ① Model name: Displays the model name of **QuantorDent**.
- ② Current time: Displays current time in hours and minutes.
- ③ System: Displays the System menu.
- ④ Create New Study: Displays the New Study screen.
- ⑤ Study List: Displays the Study List screen.

4.2. *System menu*

The System button is always displayed in the upper-right corner of the program. When you click the button, the following menu appears on the screen.




Figure 4.2 System menu

- Settings: Displays a setting dialog box.
- IP calibration: Runs IP calibration program.
- Scanner control: Runs scanner control program.
- Check update: Checks new version. Checks for a new version. If a new version is available, it is automatically downloaded and installed.
- About this software: Checks the current software version.
- Exit: Exits the program.

5. Create New Study

5.1. Screen Layout

The New Study screen is the first page to start a scan. Enter necessary patient and study information. You can easily go to the next input field by using the Tab key. Fields marked with * are mandatory fields that must be filled in.

Figure 5.1 Create New Study

- ① Emergency scan: Automatically fill in all the mandatory fields. Used in a case of emergency where patient information has not been identified or when scanning needs to be performed first and the patient information will be entered later.
- ② Edit Input Field: It displays the Edit Input Field dialog box. In the Edit Input Field dialog box, you can change the order of input fields or add or delete input fields.
- ③ System Menu: Displays the system menu.
- ④ Input Field: Table 5.1 provides the input information.

Table 5.1. Input Field

Title	Description	Mandatory
Name	Enter the patient's name.	O
ID	Enter the patient ID.	O
Birth	Enter the patient's date of birth. You can select a date from the drop-down box or manually enter numbers. When you enter a date of birth only without an age, the age is calculated automatically.	O
Sex	Select the patient's sex. You can choose one of Male, Female, NA, Male Sterile, Female Sterile or Sterile.	O
Age	Enter the patient's age. Leave the Year field blank if the patient is less than 1 year old. When you enter age only without a date of birth, the date of birth are calculated automatically. In that case, month and day will be designated as January 1.	X
Ref. Phy.	Enter the ref. physician's name.	X
Description	Enter study information.	X
Accession No.	Enter the accession number. If it is not entered, although it is a mandatory field, a random number is automatically assigned.	O
Species	Enter the patient species. You can choose either Canine or Feline.	
Breed	Enter the patient breed.	
Guardian	Enter the guardian information.	

- ⑤ Next/Previous Page: Move between pages.
- ⑥ Next: Go to the next screen of site selection.

5.2. *Edit Column Dialog Box*

In the Input Field Edit dialog box, you can change the order of the input fields, and add or remove input fields. They can be adjusted depending on your environment. However, mandatory input fields cannot be removed.




Figure 5.2 Edit Input Fields

- ① Lists input fields that will appear on the screen.
- ② Lists input fields that will not appear on the screen.
- ③ Moves the selected fields to the left or right. Fields moved to the left will appear on the screen and those moved to the right will not appear.
- ④ The buttons are used to move between pages.
- ⑤ Moves the selected field upward or down to change the order.
- ⑥ Saves changes and closes the dialog box.
- ⑦ Cancels changes and closes the dialog box.
- ⑧ Resets the program to default settings.

6. Scan

In the Scan screen, you can insert the IP into **FireCR Dental** and acquire the image.

6.1. Screen Layout




Figure 6.1 Scan Screen

- ① Image window: Displays an image selected in the Task List.
- ② Support window: Displays various information related to the teeth.
- ③ Tool box: You can change contents shown in the image view window and support window.

6.2. *Tool Box*



Figure 6.2 Tool Box

- ① Image window: Displays images only.
- ② Image window + Dentition: Displays images and dentition.
- ③ Image window + chart: Displays images and chart.
- ④ Chart: Displays chart only.
- ⑤ Template: Displays template instead of chart.

6.2.1. *Image*

When you click the Image button in the tool box, the image window only displays an image.




Figure 6.3 Image only

6.2.2. *Image+ Dentition*

When you click the "Image + Dentition" button in the tool box, it displays an image on the top and dentition on the bottom. Using dentition, you can record which teeth have been scanned after scanning completion.




Figure 6.4 Image+ Dentition

6.2.3. *Image + Chart*

When you click the "Image + Chart" button in the tool box, it displays an image on the top and the chart on the bottom. The chart provides the international standard teeth number in the center. You can drag an image from the image window or task bar and align them on the chart.




Figure 6.5 Image + Chart

6.2.4. *Chart*

When you click the "Chart" button in the tool box, the chart is provided in full screen.




Figure 6.6 Chart

6.2.5. *Template*

When you click the Template button in the tool box, templates are provided that can scan the Full Mouth Series. Drag images from the image field or task bar and align them on the template.




Figure 6.7 Template

6.3. *Image Scan*

6.3.1. *For Standalone version*

1. Enter study information in the New Study screen and click "Next" . Or select a study in

- the Study List screen and click “Scan”.
2. Place an IP on the scanner tray and push it into the scanner.
 3. The screen displays the scanning progression.




Figure 6.8 Scan in Standalone version

4. When the scan is completed, the image is provided on the screen.

Note: When you insert the tray into a scanner without opening the study, the scanner conducts scanning but does not transport the image to **QuantorDent**.

6.3.2. *For Client version*

1. Enter study information in the New Study screen and click “Next” . Or select a study in the Study List screen and click “Scan”.
2. Tag an IP to get the image on the Fire ID. Confirm that a task is added to the task bar.




Figure 6.9 Scan in Client version

3. Radiate the X-ray on an IP.
4. Place the IP on the scanner tray and push it into the scanner.
5. The navigation bar displays the scanning progression.




Figure 6.10 Indication of scanning progression

6. When the scan is completed, the image is provided on the screen.

Note: In the Client version, you must open a study for tagging but you can use other studies during scanning.

7. Image

7.1. Screen Layout

When an IP scan is completed, the Scan window displays the acquired image as shown in Figure 7.1.




Figure 7.1 Image output

- ① ROI box: A tool for cropping the valid region of the scanned image.
- ② ROI: Displays or hides the ROI box.
- ③ Marking: Adds marking on the image.
- ④ Rotation/flip: It rotates the image or flips it.
- ⑤ 180 degree rotation: Rotates the image 180 degrees.
- ⑥ Auto window: Automatically adjusts the image brightness and contrast.
- ⑦ Reject: It adds Rejects marking the image.
- ⑧ QA: Calls tools for adjusting image details.
- ⑨ Transport: Transports the image.
- ⑩ Memo: Used to enter image-related contents.

7.2. ROI

After completion of scanning, the ROI box is automatically provided as shown in Figure 7.1 along with the image. When you click and drag the ROI Box to the desired location and then double-click the ROI Box or click on the ROI button, the ROI is cropped and fitted to the screen as shown in Figure 7.2. To adjust the size of the ROI Box, drag the small rectangles at each corner of the ROI. When transmitting or printing the image, only the ROI are used. To bring up the ROI Box again, double-click on the image or click on the ROI button.




Figure 7.2 Image with cropped ROI region


7.3. Marking

This adds markings on the image. You can add markings to the image by selecting a predefined marking or manually typing text.

Add predefined marking


- ① Click on the Marking button.
- ② Select a marking.
- ③ Click a point on the image.

You can edit predefined markings in the [Overlay - Marking] under Settings.







Direct input marking

- ① Click on the Marking button.
- ② Select Direct Input.
- ③ Click on a position where you want to add the marking and a text box appears. Enter your text. Click outside the text box when you are done.



7.4. *Rotation/Flip*

You can rotate or flip the image.

Icon	Action
	Rotate counter-clockwise 90 degrees
	Rotate clockwise 90 degrees
	Horizontal flip
	Vertical flip

7.5. *180-degree rotation*

Rotate the image 180 degrees.

7.6. *Auto Window*

Automatically adjusts the image brightness and contrast.

7.7. Reject

In spite of poor quality of a scanned image, Reject preserves the current image with a “Reject” marking without rescanning. When you click the Reject button, a “Reject” marking is added on the upper-left corner of the image as shown in the Figure, and a “Reason” text box appears. Type in the reason for rejection, and click outside the text box when you are done. If you click the Reject button again, Reject is inactivated.




Figure 7.3 Reject

7.8. QA – Image Adjustment

This can adjust the LUT curve of the image. Adjusting the image requires expert knowledge about images. Adjusting the image without proper knowledge may result in lower image quality.





Figure 7.4 LUT Curve

- ① LUT curve: The LUT curve is applied to the image. A default curve for each region is provided. The user can adjust the curve to their preference. To adjust the curve, locate a point with the mouse and move it. Clicking a point on the line adds a new point. When the a point is dragged out of the box, it is removed.
- ② Histogram: Shows the image histogram. The bright region at the center is the region currently displayed on the screen. The left part of the histogram indicates the frequency where the X-ray penetration level is low and image is bright. The right part of the histogram indicates the frequency where X-ray penetration level is high and the image is dark. Images have low pixel values in bright regions and high pixel values in dark regions.
- ③ Window Center/Width: Adjusts the window center and width of the image. The window center is adjusted when you drag the image up or down, and window width is adjusted when you drag the image left or right. Dragging the sliders has the same effect as dragging the image.



- ④ Set: Sets the curve and window value adjusted by the user as the default value for images of the current region.
- ⑤ Load: Calls and applies default curve and window values of the current region.
- ⑥ Reset: Resets the curve.
- ⑦ Processing: Displays the Image Processing Filters screen.

7.9. QA – Image Processing

This is used to adjust image processing filters. Adjusting the image requires expert knowledge about images. Adjusting the image without proper knowledge may result in lower image quality.




Figure 7.5 Image processing filter

- ① Contrast: Increases the image contrast. The level can be adjusted from 0 to 10. Higher levels increase the intensity. Filter is not applied when 0 is selected.
- ② Detail contrast: Enhances contrast of regions with a slight change in brightness. The level can be adjusted from 0 to 10. Higher levels increase the intensity. Filter is not applied when 0 is selected.
- ③ Edge Enhancement: Emphasizes the border lines. It is used to increase visibility of blood vessels or small calcifications. The level can be adjusted from 0 to 5. Higher levels increase the intensity. Filter is not applied when 0 is selected.
- ④ Edge Frequency: Designates the range of the edges to emphasize. The level can be adjusted from 0 to 5. Lower level emphasizes the smaller edges. Higher levels emphasize the bigger edges.
- ⑤ Latitude Reduction: Compresses the low frequency regions to improve the overall visibility of the image. The level can be adjusted from 0 to 10. Higher levels increase the intensity. Filter is not applied when 0 is selected.
- ⑥ Noise Reduction: Reduces the image noise. The level can be adjusted from 0 to 10. Higher levels increase the intensity. Filter is not applied when 0 is selected.
- ⑦ Set: Sets the filter values adjusted by the user as the default value for images of the current region.
- ⑧ Load: Loads and applies filter values of the current region.
- ⑨ Reset: Resets the filter.
- ⑩ LUT: Displays the LUT Curve adjustment screen.

7.10. *Image Transport*

Immediately transmits the current image to the image server designated in Settings.

7.11. *Image Comment*

You can add a memo about image-related findings and extra information.

8. Task List

8.1. Screen Layout

This is a list of 'to-be-scanned' tasks. Tasks are group by unit of study, and arranged in chronological order. The user can select the desired region and perform a scan. The task list consists of lists of studies and each study is +made up of lists of tasks. The check boxes located to the right of the task list are used to select tasks.




Figure 8.1 Task List

- ① Study bar: Displays date, time and description of the study. Studies on the day of scanning display time only, while studies before the scanning day display date only.
- ② Task bar: Displays image preview and scanning region.
- ③ Select all: Selects all tasks.
- ④ Delete: Deletes the selected tasks.
- ⑤ Transport: Transports the selected studies. The transport destination is added in Settings.
- ⑥ Print: Prints the image using the Window Printer.
- ⑦ Previous/Next page: Moves the page.

8.2. *Print*

You can print out images with the Window Printer by selecting images in the task bar and clicking the Print button. The printed page can be previewed on the left part of the screen.





Figure 8.2 Print

- Layout: Selects the layout of the image to be printed on film or paper.



- Print: Starts printing.
- Setting: Displays the Print Setting screen.
- Close: Closes the dialog box.

9. Review

9.1. Screen Layout

This supports various image comparisons, processing and measurement tools necessary for diagnosis. The image view screen is composed with multiple views and each view consists of multiple images. You can customize the view layout and image layout according to your needs.




Figure 9.1 View – Image layout

To display an image on the view, you can simply drag an image from the task list. When an image is dragged from the task list, the view displays the selected image only. Dragging a study displays all images of the study on the view.




Figure 9.2 Review screen

- ① View caption bar: The view caption bar.
- ② Full screen: Provides the review window in full screen.
- ③ View layout: Adjusts the view layout.
- ④ Zoom: Used to zoom in or out of the image. Drag the image left, right, up, or down using the mouse.
- ⑤ Pan: Pans the image. Drag the image left, right, up, or down using the mouse.
- ⑥ Partial zoom in: Zooms in a part of the image.
- ⑦ Image: Opens image manipulation tools for rotating, flipping and inverting images.
- ⑧ Measure: Provides various measurement tools for measuring distance and angle.
- ⑨ Marking: Provides marking tools.
- ⑩ File: Displays DICOM file on the screen.
- ⑪ Tools: Provides other tools.
- ⑫ Apply: Designates the range in which image adjustments apply.

9.2. *View Caption Bar*



Figure 9.3 View Caption Bar

- ① Image layout: Modifies the image layout. When you click on the button, a 3x3 grid appears on the screen as shown in Figure 9.4. Select your desired layout size. You can select up to a 3x3 layout.




Figure 9.4 Image Layout

- ② Previous Page
③ Next Page
④ Close









9.3. Full Screen

This extends the review window to full screen as shown in Figure 9.5.











Figure 9.5 Full Screen


9.4. *Image Processing*

Icon	Action
	Rotate counter-clockwise 90 degrees
	Rotate clockwise 90 degrees
	Bilateral (right and left) flip
	Up and down flip
	Reset the image to the initial status.
	Negative image
	Fit to the screen size
	100% size


9.5. Measure

Icon	Action
	Distance
	Angle
	Rectangle
	Circle
	Polygon
	Free curve
	Pixel value
	Profile

- ① Distance: Measures distance. Click the start and end point of the measurement. A straight line is displayed between the two points and the measurement value is displayed at the end of the straight line. You can move the measuring line by dragging it. You can move either one of the two control points by dragging them with your mouse. You can also move the position of the measurement values by dragging them with your mouse.




- ② Angle: Measures an angle. Click on the center point of the angle you wish to measure, and then click on two control points. You can move both lines at once by dragging. You can also move either one of the two control points. You can also move the position of the measurement values by dragging them with your mouse.




- ③ Rectangle: Draws a rectangle. Click the mouse at the desired position of the upper left corner of the rectangle, and then drag the mouse to the desired position of the lower right corner of the rectangle and release the mouse button. Click on an area within the rectangle to move it. Drag the control points to resize the rectangle.




- ④ Ellipse : Draws an ellipse. Click the mouse at the desired position where the top and left-most point of the ellipse would intersect, drag the mouse toward the bottom right and release the mouse after reaching a desired size and shape. Click on an area within the ellipse to move it. Drag the control points to resize the ellipse.



- ⑤ Polygon: Draws a polygon. A control point is added every time you click the mouse. Double-click when you are done drawing the polygon. Click on an area within the polygon to move it. Drag the control points to change the corners of the polygon.




- ⑥ Free Draw: Creates a free draw selection. Click and drag the mouse to the desired shape. Release the mouse when you are done drawing. You can click inside the selection to move the free draw selection.



- ⑦ Pixel Value: Measures pixel value. Click and drag the mouse and the coordinates (x, y) and pixel values (p) of the pixel at the cursor position will be displayed on the screen.





- ⑧ Profile: Draw a straight line in the same manner as measuring a distance, and a profile dialog box will be displayed automatically.






1. Distance
2. Start point
3. End point
4. Minimum value
5. Maximum value
6. Average value
7. Pixel value: Pixel value of the current position when the mouse was clicked over the graph.





9.6. Marking

Icon	Action
	Arrow
	Text

9.7. File

Icon	Action
	Open DICOM file
	Open a general image file.
	Export

9.8. Tool

Icon	Action
	Copy to the clip board
	<p>Set a pixel size.</p> <ol style="list-style-type: none"> 1. Scan a measuring tool with a determined length. 2. Select a scanned image and click the Pixel Gap Calibration button. 3. The following dialog box will be displayed on the screen. Enter the actual length of the measuring tool. <div data-bbox="750 1464 1152 1677" data-label="Image">  </div> <ol style="list-style-type: none"> 4. Click on the start point of the measuring tool over the image and start dragging. 5. Release the mouse at the end point of the measuring tool to complete the calibration.
	Box Auto Window

9.9. *Apply*

Set the range in which image processing will be applied. When comparing two images, you can apply same filters to both windows for easier comparison.





Icon	Action
	Active: The last-clicked image is the active image. Active image is indicated with orange edges.
	Select: Multiple images can be selected by clicking images while holding the Ctrl key. The selected image is shown with blue edges.
	Study: Image processing is applied to all images in the study that the active image belongs to.
	All: Image processing is applied to all images.

Image filters that can be used with Apply Range are as follows.

- Zoom, move, rotate, flip, window

10. Study list

10.1. Screen Layout

In the study list, you can manage scanned images.





Figure 10.1 Study List


- ① Edit column: It displays a dialogue box that can edit the list column.
- ② Today: Searches studies from today only.
- ③ Search: Searches studies with various conditions.
- ④ Filter: Searches patients with the same ID.
- ⑤ Edit: Edits the selected studies.
- ⑥ Delete: Deletes the selected studies.
- ⑦ Add a study: Adds a new study.
- ⑧ Scan: Opens the selected study and goes to the Scan screen.
- ⑨ Review: Opens the selected study and goes to the Review screen.
- ⑩ Transport: Transports the selected studies.
- ⑪ Export: Exports the selected studies to a CD, DVD or a local drive.
- ⑫ Previous/Next Page: Moves between pages.
- ⑬ Preview: Previews images in the selected studies.

10.2. List

- Sorting: When clicking the each column header, the list is sorted by the selected column. Each time the column header is clicked, the sorting order changes.



Sort in ascending order



Sort in descending order

- Multiple selections: You can use the check boxes for multiple selections. To select all lists, use the check box located in the column header.
- Move between pages: The mouse wheel or page button can be used to move between pages.
- Column width change: To change the column width, hold and move the mouse from the edge of the column header.

Table 10.1 provides the column of study list.

Table 10.1. Study List Column

Column title	Description
Date	Study date
Time	Study time
Name	Patient name
ID	Patient ID
Sex	Patient sex
Birth	Patient's date of birth
Age	Patient's age
Description	Study description
Ref. Phy.	Ref. physician's name
Acc. No.	Accession No.
Image	Number of acquired images
Transport status	Results of image transport
Guardian	Name of patient's guardian
Species	Patient species
Breed	Patient breed
W./S Name	Workstation name (when run in Client mode)

10.3. *Edit Column*




Figure 10.2 Edit Column

- ① The columns to use in the Study List.
- ② The columns not used in the Study List.
- ③ Moves the selected fields to the left or right.
- ④ The buttons are used to move between pages.
- ⑤ Moves the selected field upward or down to change the order.
- ⑥ Saves changes and closes the dialog box.
- ⑦ Resets the program to default settings.

10.4. Search



- ① Filter: It selects the column where searching is applied. When a filter has not been set and All is selected, all columns are searched.
- ② Keyword: Enters keywords to search.
- ③ Search: Starts searching.

10.5. Export




Figure 10.3 Export


- ① Device: Selects a storage device to export the images to. Optical disk attached to the PC are listed. When selecting Local, you can use other storage devices attached to the PC.
- ② Label: The label is recorded on the disk.
- ③ Media Type: Displays the type of media inserted in the drive.
- ④ Local Folder: This is activated when you select Local in the Device. Select a destination folder to save images.
- ⑤ Format: Select a format to save images. Supported formats are as follows.
 - A. DICOM
 - B. Bitmap: Windows bit map file

- C. Raw: Original file
- D. Jpg: Jpeg
- E. Tiff: Tiff
- ⑥ Log message: Displays the optical disk recording progression status.
- ⑦ Progression bar: Indicates the image export process.
- ⑧ Study List: Displays the list of selected studies.
- ⑨ Image List: Displays the selected fields from the Study List. If you click and deactivate an image, the image will not be exported.
- ⑩ Select all: Activates all images.
- ⑪ Deselect all: Deactivates all images.
- ⑫ Anonymize: It deletes patient name and ID when saving as DICOM images.
- ⑬ Close: Closes the dialog box.
- ⑭ Abort: Stops recording.
- ⑮ Burn: Starts recording.

11. Confirm Transport Results

In the Transport screen, you can see the transport results of images sent to the image server.

11.1. Screen Layout



Date	Time	Lap	AE Title	Host	Port	Name	ID	Status	
2014/02/07	16:45:16	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:14	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:13	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:12	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:10	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:09	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:07	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:06	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:05	0	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:03	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:01	1	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>
2014/02/07	16:45:00	0	bunny	127.0.0.1	3000	Emergency	ID000001	Sent	<input type="checkbox"/>

Figure 11.1 Transport Information

- ① Today: Searches studies from today only.
- ② Study: Searches studies with various conditions. The same method as using the Study List is used.
- ③ Delete: Deletes the selected lines.
- ④ Previous/Next Page: Moves between pages.

11.2. *List*

Column title	Description
Date	Date to start
Time	Time to start
Lap	Time taken
AE Title	Destination AE title
IP address	Destination IP address
Port	Destination port
Name	Patient name
ID	Patient ID
Status	Transport status

12. Settings

12.1. Screen Layout

The Settings screen consists of a main category, sub-category, and specific item. Clicking on the main category item changes the sub-category items. Clicking on the sub-category item changes the specific item in the selected sub-category. Setting values can be changed in the specific item.




Figure 12.1 Settings

- ① Main category item list
- ② Sub-category item list
- ③ Specific item list

12.2. General

12.2.1. Options

Item	Description	Default value
Selected ROI Color	Select the ROI box color. You can change the display color by clicking the item.	Orange
Apply Auto Window After Cropping	Determines whether to apply Auto Window after cropping the ROI box.	Enabled

12.2.2. User Information

Item	Description
Institution name	Designate the hospital name. The entered values are included in the DICOM information.
Dept. name	The department name where the system is installed. The entered values are included in the DICOM information.

12.2.3. Password for Processing Parameter

Item	Description	Default value
New password	Set a password.	
Confirm password	Enter the password again.	
Require Password for Modifying Image Processing	Set the system to request a password when attempting to modify LUT and processing parameter.	Enabled
Require Password for Saving Image Process Values	It activates password screen when saving default values in the LUT and processing windows.	Enabled

12.2.4. Display Language

Item	Description	Default value
Language	Select language.	Language selected upon product installation

12.3. *Fonts*

12.3.1. *General*

Item	Description	Default value
Font name	Name of font used on the screen.	Malgun Gothic
Font size	Size of screen font. Adjustable from 8 to 12 points.	10

12.3.2. *Marking*

Item	Description	Default value
Font name	Name font used for marking.	ARIAL
Font size	Size of marking font. Adjustable from 16 to 30 points.	22

12.3.3. *DICOM Printer*

Item	Description	Default value
Font name	Name of font to use for printing.	ARIAL
Font size	Size of printing font. Adjustable from 8 to 12 points.	10

12.4. Network / Export

12.4.1. Image Server

This manages lists of DICOM servers to transport images.




Figure 12.2 Network / Export

- ① Add: Adds an image server.

Item	Description
Local AE Title	Set AE title of Quantor.
Server AE Title	Set the AE Title of image server.
IP address	IP address of image server
Port	Port of image server
Description	Description of the image server

- ② Edit: Edits the selected studies. This displays the same screen as adding an image server.
- ③ Delete: Deletes the selected servers.
- ④ Network test: Checks the connection by sending DICOM Echo to the selected server.
- ⑤ Move up: Moves the selected server up.

- ⑥ Move down: Moves the selected server down.
- ⑦ Timeout: The connection is dropped if a response is not received within the set amount of time.

Once a server is added, the check boxes for each item are automatically activated. Quantor sends images to all checked image servers in the list order.

12.4.2. Options

Item	Description	Default value
Automatic transfer of studies	Study images are transferred automatically after completion of the study.	Disabled
Series Encoding	Set DICOM Encoding method. <ul style="list-style-type: none"> ● A series has one image. ● A series has multiple images. 	A series has one image.

12.5. Overlay

12.5.1. Marking

Here you can edit predefined markings that appear when adding marking on the image in the Scan screen. After entering text in the input field, you can perform the following actions:

- Add: Add a new marking.
- Modify: Modify the selected marking.
- Delete: Deletes the selected marking.
- Move up: Moves the selected marking to one line up.
- Move down: Moves the selected marking to one line down.

12.5.2. Auto Marking Margin

This sets the margins when markings are automatically added to an image. Each setting value shows margins by % compared to the width and height of the image.

Item	Description	Default value
Left	The left margin.	10
Top	The top margin.	10
Right	The right margin.	10
Bottom	The bottom margin.	10

12.5.3. *Overlay Contents*

This selects items to show in the image when printing.

12.6. *System*

12.6.1. *File Management*

Item	Description	Default value
Data path	This sets the destination folder where images are stored.	Installation Folder/Database
Data space	This allocates disk space for saving images. Auto deletion is activated when the allocated disk space becomes full.	30GB
Auto Deletion	This sets whether to use Auto Deletion and Auto Deletion intervals. When the user does not take any action, and the auto deletion interval is exceeded, auto deletion starts.	Enabled 60 sec
Delete studies that are not sent	This defines the action to take when deleting an image that has not been transmitted. <ul style="list-style-type: none"> ● Never delete ● Check before deletion ● Always delete 	Never delete

12.6.2. *System Options*

They are other system options.

Item	Description	Default value
Temporary File Path	This sets a folder to save temporary files used in the program.	Installation Folder/Temp
Virtual Keyboard Enabled	This defines whether to use the Windows virtual keyboard when entering data in the screen.	Disabled

12.6.3. *Date*

Item	Description	Default value
Date format	This sets the date display format. You can choose one of following three formats. <ul style="list-style-type: none"> ● YYYY/MM/DD ● DD/MM/YYYY ● MM/DD/YYYY 	YYYY/MM/DD
Separator	This selects a separator for year, month and day. Choose one from following three. <ul style="list-style-type: none"> ● / ● - ● . 	/

12.6.4. *Dental*

Item	Description	Default value
Workstation name	Enter the workstation name.	

Appendix A: Canine Dental Chart

G	110	109	108	107	106	105	104	103	102	101	201	202	203	204	205	206	207	208	209	210	G
411	410	409	408	407	406	405	404	403	402	401	301	302	303	304	305	306	307	308	309	310	311

Appendix B: Feline Dental Chart

109	108	107	106		104	103	102	101	201	202	203	204		206	207	208	209
409	408	407			404	403	402	401	301	302	303	304			307	308	309

Appendix C: Template

1	2	3	4	5
6	7	8	9	10